

Memo To:

Jonathan Young

7/6/88

From:

Glenn W. Calvert

Subject:

Leak Test of Transfer Cask

On July 5, 1988 a Model L-3 Leak Test Kit was used to perform a leak test on a Transfer Cask located at the Parsippany, NJ facility for detection of any removable radioactive contamination. Samples taken by you were brought to Dosimetry Services for analysis.

The samples were taken at seven (7) locations: three on the lid of the cask, the vent hole top and bottom and two on the sides of the cask.

Measurements for removable contamination were made using an Eberline Mini-Scaler (Model MS-2) Detector with a Model HP-21 Probe. An efficiency of 4.54% was determined with the high voltage set at 721. Counts per minute (cpm) of the Cobalt 60 Standard were 388 and counter background cpm was 37. Counts for each of the seven locations are listed below:

Location	Counts per Minute
 Lid Lid Lid Vent Hole (top) Vent Hole (bottom) Side of Cask Side of Cask 	39 29 40 34 33 38 38

Note: Remonable contamination was < 005 plie problem mphi. Malt 1/1/26 cc:M. Saylor

ITEM #

ISOMEDIX INC.

CORPORATE OFFICES • 11 APOLLO DRIVE, WHIPPANY, NEW JERSEY 07981 • (201) 887-4700 • FAX 887-1476 • TELEX 317361