

METAL MANAGEMENT, INC.

Metal Management – Northeast

Connecticut Operations Schiavone Operations 234 Universal Drive North Haven, CT 06473 Fax: 203/782-4272 Fax: 203/782-4210 Voice: 203/777-2591

Fax Cover Page



Metal Management Michael Schiavone & Sons Tony Altieri 234 Universal Dr North Haven, CT 06473



Most regulatory agencies require that you leak-test your sealed radioactive source before you use it & at six month intervals thereafter. You may be asked to show proof of compliance during an inspection by the regulatory agency.

***Leak Test Certificate ***

Leak test samples are obtained using a" KVXL-TERT" or "JH-7" leak test kit as authorized by California Radioactive Materials License, # 1655-31.

The sealed source(s) noted below was/were checked for leakage by wiping the source(s) and the holder assembly with a moistened swab which was then analyzed for contamination using a window-less gas-flow proportional counter. The test showed less than 0.001 microCuries of leakage (contamination). The maximum allowable leakage is 0.005 microCuries.

----- This radioactive sealed source passes the leak-test

Source	Activity, Isotope,	Activity, Isotope,	Leak test	Leak test	Leak test
serial #	and Assay Date	and Assay Date	Collected by:	Analyzed by:	Due Date
3713A	21 mCi of Cd109 as of 15 Sep 98	21 mCi of Am241 as of 15 Nov 82	Ed Ammerer on 28 Dec 99	John Harrison on 07 Jan 00	28 Jun 2000

rison, Radiation Sylety Office

0202 Source



Michael Schiavone & Sons John Sartori 234 Universal Dt North Haven, CT 06473



Most regulatory agencies require that you leak-test your sealed radioactive source before you use it & at six month intervals thereafter. You may be asked to show proof of compliance during an inspection by the regulatory agency.

***Leak Test Certificate ***

Losk test samples are obtained using a" KVXL-TEST" or "HLT" leak test kit as authorized by California Radioactive Materials License, # 1655-31.

The sealed source(s) noted below was/were checked for leakage by wiping the source(s) and the holder assembly with a moistened swab which was then analyzed for contamination using a window-less gas-flow proportional counter. The test showed less than 0.001 microCuries of leakage (contamination). The maximum allowable leakage is 0.005 microCuries.

----- This radioactive sealed source passes the leak-test

Source	Activity, Isotope,	Activity, Isotope,	Leak test	Leak test	Leak test
serial #	and Assay Date	and Assay Date	Collected by:	Analyzed by:	Due Date
3713A	21 mCi of Cd109 as of 15 Sep 98	21 mCi of Am241 as of 15 Nov 82	Ed Ammerer on 14 Jul 99	John Harrison on 04 Aug 99	14 Jan 2000

John Harrison, Radiation Safety Offices



Metal Management Tony Altieri 234 Universal Dr North Haven, CT 06473



Most regulatory agencies require that you leak-test your sealed radioactive source before you use it & at six month intervals thereafter. You may be asked to show proof of compliance during an inspection by the regulatory agency.

***Leak Test Certificate ***

Leak test samples are obtained using a" KVRL-IEST" or "H-7" leak test kit as authorized by California Radioactive Materials License, # 1655-31.

The sealed source(s) noted below was/were checked for leakage by wiping the source(s) and the holder assembly with a moistened swab which was then analyzed for contamination using a window-less gas-flow proportional counter. The test showed less than 0.001 microCuries of leakage (contamination). The maximum allowable leakage is 0.005 microCuries.

----- This radioactive sealed source **passes** the leak-test ------

Source	Activity, Isotope,	Activity, Isotope,	Leak test	Leak test	Leak test
serial #	and Assay Date	and Assay Date	Collected by:	Analyzed by:	Due Date
3713A	21 mCi of Cd109 as of 15 Sep 98	21 mCi of Am241 as of 15 Nov 82	Ed Ammerer on 01 Aug 00	John Harrison on 09 Aug 00	01 Feb 2001

John Harrison, Radiation S

RSA LABORATORIES, INC W-100 Leak Test Kit

SAMPLE DOCUMENTATION AND REPORT OF RESULTS

Company Name: MTLM Acrospace

Radioactive Materials License #: State of CT Ionizing Reg # 917

Company Address: 500 Flatbush Ave. Hartford, CT 06106

Person to be contacted: Eric R. Hopkins Phone: (860) 777-2591

For the source to be leak tested:

Source Manufacturer: Kevux

Source Model No.: 6700

Source Serial No.: 18917

Source Isotope: Cd-109

Source Quantity: 8 mCi

Date of Wipe Test: 03/05/2001

Person performing wipe test: Eric R. Hopkins Return these samples, the empty vial and this data sheet to:

> RSA Laboratories, Inc. 21 Pendleton Drive, P.O. Box 61 Hebron, CT 06248

If you have any questions about how to perform this sample procedure, call RSA Laboratories, Inc. at (860) 228-0721

SEALED SOURCE LEAK TEST CERTIFICATION

This is to certify that a wipe test performed on the above-named source indicates the presence of $<.005 \ \mu\text{Ci}$ of total removable radioactivity.

	SAMP	LE ANALYSIS RE	SULTS	
Total cpin	Background cpm	Efficiency	dpm	Total Activity
Gamma Spec	Gamma Spec	Gamina Spec	Gaiuma Spec	<4.001/-05

Person Performing the Analysis: Jay R. Dockendorff

Date 9hAnalysis: 03/26/2001

12101

Original to customer Copy to ESA file

Rev. 1/8/93

VHJr 8899 | e.mail: jvhje@neworld.net JohnHarrison 140 Seriedand Drive Auburn, California, 95603

916, 782, 8899 or: 530, 906.

٢,

Eric Hopkins, Metals Management

203 782 4272

31 Oct, 2001

Re Your question regarding the leak testing of radioactive sealed sources.

Eric,

months. An x-ray tube based analyzer contains no "radioactive sealed source" and does not require a "leak test". Radioactive sealed sources (in most cases) need to be leak tested every 6

Sincerely, John Harrison

 ${}^{\prime\prime}$

Jun. 26 2001 01:49PM

ñ.

. . FRX ND. Valley Safety Services Associates 330 Old Enfield Road, Belchertown, MA 01007 (413) 323-9571

Source/Device Description

Source tested prior to installation

Installed in XL Serial #: V1296 Source #: 9980LG Activity: 14 mCi

Device: XRF Analyzer Make/Model : XL 300, 700 or 800 Manuf: Niton

Radioisotope: Am-241 Source Model # : AMC.P4 Manufacturer: Amersham/AEA

The leak testing of this source/device was performed as indicated in manufacturer recommendations for LT-952 kit instructions using kit # 2.5 8 2/1 by: Jim Failla _____ Date: 2.5 6 1

Company Name:Niton Corporation (Ken Martin, RSO)Address:1130 Ten Rod Road Suite C-207North Kingstown, RI 02852

Telephone: (401) 294-1234

CERTIFICATE #: 958211

Do Not Write Below This Line

Fax: (401) 295-2090

Leak Test Analysis Result for above source/device

Analysis of the above test on $2\frac{1}{2600} ol$ yielded the following:

Statistical analysis of the radioactive count data of this leak test specimen indicated any activity present is less than 0.005uCi. The source may be used as authorized.

Statistical analysis of the radioactive count data of this leak test specimen indicated any activity present is <u>greater than 0.005 microcuries of activity present</u>. This source should be considered leaking. Consult your device operations manual; place this unit in storage and make the required notification to your regulatory agency.

YOUR NEXT REQUIRED LEAK TEST FOR THIS DEVICE/SOURCE IS DUE ON OR BEFORE _______ AUGUST 01

THIS CERTIFICATE IS AN ESSENTIAL RECORD AND SHOULD BE MAINTAINED FOR INSPECTION BY THE REGULATORY AGENCY.

BY

MA MAL DATE DEBOI

Valley Safety Services Associates 330 Old Enfield Road, Belchertown, MA 01007 (413) 323-9571

Source/Device Description

Source tested prior to installation

Installed in XL Serial #: 1/296 Source #: 9980 LQ Activity: 14 mCl

Manuf: Niton Make/Model : XL 300, 700 or 800 Device: XRF Analyzer

Radioisotope: Am-241 Source Model # ; AMC.P4 Manufacturer; Amersham/AEA

The leak testing of this source/device was performed as indicated in manufacturer recommendations for LT-952 kit instructions using kit # 759320 9/7/01 by: Jim Failla Date:

Company Na	me: Niton Corporat	tion (Ken Martin, RSO)
Address:	1130 Ten Rod I	Road Suite C-207
	North Kingstov	vn, RI 02852
Telephone:	(401) 294-1234	Fax: (401) 295-2090

Telephone: (401) 294-1234

CERTIFICATE #: <u>459.320</u> BY:

Do Not Write Below This Line

Leak Test Analysis Result for above source/device

Analysis of the above test on 12367-01 yielded the following:

Statistical analysis of the radioactive count data of this leak test specimen indicated any activity present is less than 0.005uCi. The source may be used as authorized.

Statistical analysis of the radioactive count data of this leak test specimen indicated any activity present is greater than 0,005 microcuries of activity present. This source should be considered leaking. Consult your device operations manual; place this unit in storage and make the required notification to your regulatory agency.

YOUR NEXT REQUIRED LEAK TEST FOR THIS DEVICE/SOURCE IS DUE ON 12 MARC 2002 OR BEFORE

THIS CERTIFICATE IS AN ESSENTIAL RECORD AND SHOULD BE MAINTAINED FOR INSPECTION BY THE REGULATORY AGENCY.

emitai DATE: 12 Sept 01



Niton XL# : U12964951LY

Page: 1 of 2

	SERV	ICE RE	PORT	Ship Date:
		PO	Number:	
Ship To:	JOHN SARTORI METAL MANAGEMENT- NORTHEAST 234 UNIVERSAL DRIVE		Customer# : 2193	
	NORTH HAVEN USA	СТ	06473-3105	•

RA:

Niton XL# :U12964951LYModel: 821S

Detector Type: T-08 Features : Alloy Analyzer

Speed: Seconds

Hours Used:

Leak Test Certificate: Yes

Leak Test Performed on: 6/23/00

Other :See SO# 3042 for upgrade info

Description	QTY	Serial#
T-08 DETECTOR	1	577-T
AMERSHAM-10 Cadmium-109 Pin Source	1	4951LY
AMERSHAM-14 Americium-241 Pin Source	1	9980LQ
821S METAL ANALYZER (W/cd-109 & AM-241	1	U1296
BATTERY PACK - STANDARD #1	1	
LEATHER WRIST STRAP	1	
GRAY CARRYING CASE #1500	1	
CASE LABEL 1500	1	
CASE FOAM 1500	1	
CASE FOAM SOURCE GUARD	1	
BAR CODE WAND (light pen) W/ LEMO END	1	
SPARE SAPHIRE TIP FOR BAR CODE WAND	1	
12 VOLT AUXILLARY POWER CABLE	1	
RS232 DOWNLOADING CABLE	1	

900 Middlesex Turnpike, Building 8 • Billerica, MA 01821-3926 USA Tel: 978-670-7460 • 800-875-1578 • Fax: 978-670-7430 • www.niton.com • email: xrf@niton.com



Niton XL# : U12964951LY

Page: 2 of 2

Description	QTY	Serial#
RS232 CABLE ADAPTOR FOR PC	1	
ALLOY USER MANUAL	1	
XL CONVERSION FROM 801S TO 821S	1	

Problems & Work Performed Detail

	Problem	Work performed
1	Upgrade to an 821S	Upgraded
2	Windows Upgrade from earlier version to 5.5	Supplied new diskette in customers case.
3	XL Software upgrade to current version.	Upgraded XL software to current version.
4	Calibration & Quality Check.	Calibrated XL, performed a quality check on standards, downloaded to XTRAS Data System.

EVALUATION OF POSSESSION AND USE OF BYPRODUCT MATERIAL

Jame: <u>Michael Schiavone de Sons</u> , Tur Address: <u>234 Univ</u> ersal Drive <u>North Haven</u> , CT 08473 FAX Number: () -	1 1 1
. Name and Title of person responsible for radiation safety program: Anthony Altieri, Radiation Safety Officer	
. Describe how you safeguard the byproduct material from:	
(a) use by unauthorized personnel: The portable Niton unit	
(b) loss or that T) here when not in use	
(b) iss of men. The two Kever units are kept in a	
Representation file bada and capiest any	
. Describe controls that prevent individuals who work in the area around the watthe	<
the material from becoming exposed to radiation: The sector	2.
are kept in self-shieldert devices with turntables for	-
the samples so the sources are not exposed.	-
	-
I. Do you have a personal monitoring program for your employees, such as	
film badges, dosimeters, etc.? Yes ☑ No □	
It yes, what was the maximum dose received since <u>1985</u> ? (year of last	
sephone contact or inspection) <u>20-50 milliven</u> in a calendar year.	
b. Do you perform surveys to detect external radiation in the area around	
If ves how often are the surveys performed?	
What instrument is used to perform the surveys? Whe	n l
/as this instrument last calibrated?	
5. On what date was the last physical inventory of all byproduct material	
in your possession performed? <u>/ol_l2c</u> o/	
Were all sources accounted for? $(1,2)$ Yes \square No \square N/A \square	
I DO YOU PERFORM IEAK TESTS ON THE SEALED SOURCE? YES NO DO N/A D	
Who evaluates the leak test resulte? Manufacture for the	
0. describe the provisions you have made to have the lask tasts	
done: Note: Niton leakstested 9/1999 6/7mm + /2001	
ind 9/2001	
3. Describe your provisions for repair and maintenance of your device or	11
3. Describe your provisions for repair and maintenance of your device or source holder: <u>Send to manufacturer</u>	
3. Describe your provisions for repair and maintenance of your device or source holder: <u>Send to manufacturer</u>	
3. Describe your provisions for repair and maintenance of your device or source holder: <u>Send to manufacturer</u>	
3. Describe your provisions for repair and maintenance of your device or source holder: <u>Send to manufacturer</u> 3. Describe any unusual events involving the byproduct material or device(a) in which it is used.	
B. Describe your provisions for repair and maintenance of your device or source holder: <u>Send to manufacturer</u> D. Describe any unusual events involving the byproduct material or device(s) in which it is used: <u>None</u>	
B. Describe your provisions for repair and maintenance of your device or source holder: <u>Send to manufacturer</u> . Describe any unusual events involving the byproduct material or device(s) in which it is used: <u>None</u>	
B. Describe your provisions for repair and maintenance of your device or source holder: <u>Send to manufacturer</u> . Describe any unusual events involving the byproduct material or device(s) in which it is used: <u>None</u>	
B. Describe your provisions for repair and maintenance of your device or source holder: <u>Send to manufacturer</u> . Describe any unusual events involving the byproduct material or device(s) in which it is used: <u>None</u>	

Issue Date: 12/01/00

· :- :-

in the second .