CENTER FOR NUCLEAR WASTE REGULATORY ANALYSES

Proc. <u>TOP-18</u>	
Revision Change	
Page of	

TECHNICAL OPERATING PROCEDURE

SOFTWARE RELEASE NOTICE PA-SRN-287 1. SRN Number: **Project** 2. Project Title: Diablo Canyon Spent Fuel Storage Facility No.20.06003.03.31 3. SRN Title: MCNP 4C2 4. Originator/Requestor: Oleg Povetko Date: 01/27/2003 5. Summary of Actions ☐ Change of access software Release of new software □ Software Retirement Release of modified software: □ Enhancements made □ Corrections made 6. Validation Status Validated Limited Validation Not Validated Explain: To be determined 7. Persons Authorized Access Addition/Change/ Read Only/Read-Write Name Delete **Read Only** N/A Oleg Povetko 8. Element Manager Approval: 9. Remarks:

SOFTWARE SUMMARY FORM

01. Summary Date: January 27, 2003	02. Summary prepared by (Name and phone) Oleg Povetko (210)522-5258		03. Summary Action:	
04. Software Date:	05. Short Title: MCNP 4C2			
06. Software Title: MCNP, version 4C2		07. Internal Software ID: None		
08. Software Type:	09. Processing Mode:	10. Application Area		
☐ Automated Data System	√ Interactive	a. General: √ Scientific/Engineering □ Auxiliary Analyses □ Total System PA □ Subsystem PA □ Other b. Specific:		
√ Computer Program	□ Batch			
□ Subroutine/Module	☐ Combination			
		12. Technical Contact(s) and Phone:		
CNWRA/SwRI 6220 Culebra Road San Antonio, TX 78228		Oleg Povetko (CNWRA) (210) 522-5258		
13. Software Application:MCNP is a general-purpose, continuous-energy, generalized geometry, time-dependent, coupled neutron-photon-electron Monte Carlo transport code system. Energy ranges are 0-60 MeV for neutrons (data generally only available up to 20 MeV) and 1 keV - 1 GeV for photons and electrons. MCNP4C2 is a version upgrade for MCNP4A which is currently controlled under TOP-018. The changes from version 4A to 4C2 include: photonuclear physics model addition, interactive plotting, implementing macrobody surfaces and others. The code utilizes the ASCII DLC-200/MCNPDATA data library which is included in the software package and installed along with the code.				
14. Computer Platform Windows-based personal computer	15. Computer Operating System: Windows-95/NT or later	16. Programming Language(s): N/A	17. Number of Source Program Statements: N/A	
18. Computer Memory Requirements: Not available	19. Tape Drives: N/A	20. Disk Units: 940 MB of hard disk space	21. Graphics: Windows-based Minimum VGA graphics card	
22. Other Operational Requirements None				
23. Software Availability: √ Available □ Limited				
25. Software User: Oleg Povetko	My 6 1	Inn//h	Date: 01/29/2003	

CENTER FOR NUCLEAR WASTE REGULATORY ANALYSES QA VERIFICATION REPORT

FOR

→ ACQUIRED SOFTWARE <u>NOT</u> TO BE MODIFIED ←

Software Title/Name: Version: Demonstration workstation: Operating System: User:	MCNP 4CZ 4CZ Palandin Windows NT Oleg Povetko	
	ot meet all requirements and will be evaluated	l on a case-by-case basis.
Installation Testing [TOP-018, Sec	tion 5.6]	
Has installation testing been conduc	operating Systems: Windows See en elosed memo Mabrillo, January 21	Yes: B No: L N/A: L
Software Output [TOP-018, Sectio	n 5.5.4]	
	dual runs are uniquely identified by d	_
Date and Time Displayed:		Yes: 🗖 No: 🗖 N/A: 🗖
NOTE: Output identification content and	format is typically taken as is.	
Medium Documentation [TOP-013		
The physical labeling of software m Module Revision, File type (ASCII,	nedium (tapes, disks, etc.) contains: Pro , OBJ, EXE), Recording Date, and Ope	gram Name, Module/Name/Title, rating System(s)? Yes: No: □ N/A: □

CENTER FOR NUCLEAR WASTE REGULATORY ANALYSES QA VERIFICATION REPORT

FOR

→ ACOUIRED SOFTWARE NOT TO BE MODIFIED ←

User Documentation [TOP-018, Section 5.5.7]				
Is there a Users' Manual for the software and is it up-to-date?	Yes:	2	No: □	N/A: □
User's Manual Version and Date: Comments:				
Are there basic instructions for the installation and use of the software?			No: □	
Location of Instructions: <u>See enclosed</u> meno Povetko Comments: January 21, 2003	+» ±	Mads.	vito da NP4CS	tel 2 Nutes
Configuration Control [TOP-018, Section 5.7, 5.9.3]				
Is the Software Summary Form (Form TOP-4-1) completed and signed?	Yes:	1 /	No: □	N/A: □
Date of Approval: 01/29(03				
Is the list of files attached to the Software Summary Form complete and acc	urate? Yes:	?	No: □	N/A: □
Comments:				
Is the source code available or, is the executable code available in the case of	of (acc Yes:	quired/o	commercial No: 🏻	codes)? N/A: □
Location of Source Code: <u>See enclosed</u> CD - Comments:				
•				
Have all the script/make files and executable files been submitted to the Sof	tware	Custo	dian?	
Only the executable files are being submitted.		_/	·· -	 -
Location of executable files: See enclosed CD. Comments:	Yes:	: (1)	No: □	N/A: □

CENTER FOR NUCLEAR WASTE REGULATORY ANALYSES QA VERIFICATION REPORT

FOR

→ ACQUIRED SOFTWARE <u>NOT</u> TO BE MODIFIED ←				
Software Release [TOP-018, Section 5.9]				
Upon acceptance of the software as verified above, has a Software Release Notice (SRN), Form TOP-6 been issued and does the version number of the software match the documentation? Yes: No: No: NA: SRN Number: 287 Comments:				
Software Validation [TOP-018, Section 5.10]				
Has a Software Validation Test Plan (SVTP) been prepared for the <i>range of application</i> of the software? Yes: □ No: ► N/A: □				
Version and Date of SVTP:				
Date Reviewed and Approved via QAP-002:				
Comments: NCMP, 4A classified as cot. II with validation due Sept. 2008.				
Has a Software Validation Test Report (SVTR) been prepared that documents the results of the validation cases, interpretation of the results, and determination if the software has been validated? Yes: □ No: ► N/A: □				
Version and Date of SVTR:				
Date Reviewed and Approved via QAP-002:				
Comments.:				
Additional Comments: Market Market				

Volume in drive D is MCNP4C2 Volume Serial Number is ECA0-7969

Directory of D:\

```
701ALLCP 00 4,122 06-06-01 6:20a 701ALLCP.00 C701DOS0 EXE 16,903,243 06-04-01 7:08a C701DOS0.EXE C701TARO GZ 13,330,060 06-04-01 8:11a C701TARO.GZ D200DOS2 EXE 268,238,353 03-01-01 8:03a D200DOS2.EXE D200TAR2 GZ 274,994,967 02-12-01 7:09a D200TAR2.GZ DOC CDIR> 06-11-01 1:22a DOC READ4C2 TXT 18,592 06-01-01 5:48a READ4C2.TXT README 1ST 8,180 06-06-01 6:11a README.1ST README RSI 4,096 12-31-99 6:00a README.RSI 8 file(s) 573,501,613 bytes
```

Directory of D:\DOC

Total files listed:

10 file(s) 589,698,632 bytes 3 dir(s) 0 bytes free TO:

Bruce Mabrito

FROM:

Oleg Povetko

SUBJECT:

TOP-018 for MCNP-4C2

DATE:

January 21, 2003

Software Application: MCNP is a general-purpose, continuous-energy, generalized geometry, time-dependent, coupled neutron-photon-electron Monte Carlo transport code system. Energy ranges are 0-60 MeV for neutrons (data generally only available up to 20 MeV) and 1 keV - 1 GeV for photons and electrons. MCNP4C3 is a version upgrade for MCNP4A which is currently controlled under TOP-018. The changes from version 4A to 4C2 include: photonuclear physics model addition, interactive plotting, implementing macrobody surfaces and others. The code utilizes the ASCII DLC-200/MCNPDATA data library which is included in the software package and installed along with the code.

It was acquired from ORNL (Oak Ridge National Laboratory) RSICC (Radiation Safety Information Computational Center), P.O. Box 2008, Oak Ridge, TN 37831-6362. Phone: (865)574-6182. Fax: (865)574-6182. Email: pdc@ornl.gov. URL: http://www-rsicc.ornl.gov/rsic.html.

Installation of MCNP4C2.

MCNP4C2/MCNPDATA software package was successfully installed on Paladin machine following included instructions:

Expanding files from the distribution CD will create a subdirectory called "mcnp4c2" that includes the code system plus the new la150u photonuclear library and associated documentation in the "xs" subdirectory under "mcnp4c2". The mcnp4c2\xsdir file includes lines for the new la150u as well as the other MCNP libraries distributed in the DLC-200 package.

DLC-200 cross sections will expand to a "xs" subdirectory under "mcnp".

QUICK START on PC

The MCNP4C2 code system is transmitted in a compressed file for PC users running Windows 9x or NT or 2000. PKware 2.60 was used to create the self-extracting compressed files. Note that the mcnp executables are not Windows applications; they must be run from a DOS prompt.

The PC distribution contains executables for MCNP4C2 and MAKXSF, which were built with the Digital Visual Fortran 6.0 Professional Edition and Lahey/Fujitsu Fortran 95 Compiler Release 5.50h Professional edition. Both X11 graphics and DVF QuickWin and Lahey Winteracter graphics are supported. See the Readme4c2.txt file for details details on compiling and on software

requirements.

Insert distribution CD-ROM and access CD drive.

Double-click on the file c701dos0.exe to expand mcnp4c2 + la150u to mcnp4c2. Double-click on the file d200dos2.exe to expand cross sections to mcnp\xs. Ensure that "Recreate Subdirectories" option is checked and change drives if desired. Then click OK to create an MCNP4C2 subdirectory and extract files.

Open a DOS window and compile as instucted in README.TXT or copy an MCNP executable from the EXE subdirectory and run test cases by invoking the DOS commands below:

cd \mcnp4c2 copy exe\mcnp4.dvq mcnp.exe (copy DVF-compiled MCNP with dynamically allcoated storage & QuickWin graphics0)

One must move the DLC-200 mcnp\xs cross sections to mcnp4c2\xs. You may edit the datapath statement in mcnp4c2\xsdir to tell MCNP where to find DLC-200 cross sections.

Installation test results.

MCNP4C2 comes with 34 test files (contrary to a number of 32 indicated in README.1ST file) that implement a wide range of the code features when executed.

Test cases use TESTDIR and TESTLIB1. Test files are stored in /samples directory.

In order to run test cases batch file runsamps was executed.

Two output files were generated for each test case.

Results of test runs are programmatically compared with results of the same cases executed by the developer and provided on the installation CD. Results of file comparison are documented in dif* files. After execution 68 dif* files were examined to verify that the test cases were run successfully. Additional visual inspection of various portions of generated and provided output files verified that corresponding output files contained the same results.

Additional test case from Reused Soil Scenario Analysis project was run using both MCNP4A and MCNP4C2 codes.

<u>Input file:</u>

tr137i*

Description of the case:

Objective:

To calculate gamma component of the external dose to the truck operator caused by irradiation of gamma rays from the contaminated soil in cargo box.

Model description:

Truck cargo box is modeled as stainless steel box, truck cab is modeled as a stainless steel sphere. Truck operator is modeled as a 70-kg sphere positioned in the center of the cab for truck types 2-4. Truck dimensions are presented in Table 1. Truck type 4 was examined in the test case (file tr137i).

Based on visual inspection of various portions of the output files, the two versions of the software produced the same results. Results are not identical due to Monte Carlo sampling of the input parameters.

Table 1. Truck Parameters for the Truck Operator Scenario				
	P2 Truck Operator Truck Types			
Parameter	Type 1 (Pickup truck)	Type 2 (Light/Medium Truck)	Type 3 (Standard dump truck)	Type 4 (Heavy dump truck)
Cargo box length	6.5'	8'	10'	17'
Cargo box width	5.025'	6'	7'	7'
Height of soil in cargo box front "shining plane"	1.2'	2.5'	3'	4'
Cargo box wall thickness	0.089 mm	0.476 mm	0.476 mm	0.476 mm
Volume of soil in cargo box	1306997 cm3	3754568.594 cm3	6225079.679 cm3	15357142 cm3
Mass of soil in cargo box	2,091 kg	6,007 kg	9,960 kg	24,571 kg
Sphere cab diameter	80 cm	80 cm	100 cm	100 cm
Sphere cab wall thickness	0.089 mm	0.089 mm	0.089 mm	0.089 mm
Position of the operator in the cab	Shifted 30 cm down and 5.02 cm back from the cab center towards cargo box	Center of the cab	Center of the cab	Center of the cab
Distance from the center of the modeled operator to the "shining plane" ¹	30.4 cm	81.476 cm	101.476 cm	101.484 cm

¹Front plain surface of the hauled soil closest to the receptor

Attached:

1 cdrom, MCNP4C2/MCNPDATA software

1 cdrom, installation test results