

**SOUTHWEST RESEARCH INSTITUTE
NUCLEAR PROJECT**

CLIENT: Div. 20

TASK ORDER: 031106-12

SRR: 25200

SDG: 237578

CASE: CNWRA

VTSR: November 06, 2003

PROJECT#: 06002.01.141

FINAL REPORT

SOUTHWEST RESEARCH INSTITUTE

SAMPLE ANALYSIS DATA SHEET

010001

Sample ID

16P-1003-FA-5

Lab Name: Southwest Research Institute

Client: Division 20

Lab Code: SwRI

Date Received: 11/06/03

Matrix: Water

Project No.: 20.06002.01.141

Lab System ID: 237578

SRR: 25200

TO: 031106-12

Analysis	Sample Result (mg/L)	Reporting Limit (mg/L)
Aluminum	0.620	0.05
Antimony	<0.01	0.01
Arsenic	0.014	0.005
Barium	0.022	0.005
Beryllium	<0.005	0.005
Bismuth	<0.015	0.015
Boron	0.204	0.025
Cadmium	<0.005	0.005
Calcium	2.17	0.05
Chromium	<0.005	0.005
Cobalt	<0.005	0.005
Copper	<0.005	0.005
Iron	1.14	0.05
Lanthanum	<0.005	0.005
Lead	<0.005	0.005
Lithium	0.081	0.005
Magnesium	0.388	0.05
Manganese	0.024	0.005
Molybdenum	0.063	0.01
Nickel	<0.005	0.005
Palladium	<0.01	0.01
Phosphorus	0.043	0.04
Potassium	1.92	0.2
Selenium	<0.015	0.015
Silicon	21.1	0.05
Silver	<0.005	0.005
Sodium	110	0.2
Strontium	0.028	0.005
Sulfur	20.2	0.025
Thallium	<0.025	0.025
Thorium	<0.01	0.01
Tin	<0.01	0.01
Titanium	<0.005	0.005
Tungsten	<0.015	0.015
Uranium	<0.1	0.1
Vanadium	<0.005	0.005
Yttrium	<0.005	0.005
Zinc	0.015	0.005
Zirconium	<0.005	0.005

SOUTHWEST RESEARCH INSTITUTE

DUPLICATE SUMMARY

010002

Sample ID
16P-1003-FA-5

Lab Name: Southwest Research Institute

Client: Division 20

Lab Code: SwRI

Date Received: 11/06/03

Matrix: Water

Project No.: 20.06002.01.141

Lab System ID: 237578

SRR: 25200

TO: 031106-12

Analysis	Sample Result (mg/L)	Duplicate Result (mg/L)	RPD
Aluminum	0.620	0.659	6.19%
Antimony	<0.01	<0.01	0.00%
Arsenic	0.014	0.008	54.9%
Barium	0.022	0.020	10.3%
Beryllium	<0.005	<0.005	0.00%
Bismuth	<0.015	<0.015	0.00%
Boron	0.204	0.206	0.95%
Cadmium	<0.005	<0.005	0.00%
Calcium	2.17	2.18	0.63%
Chromium	<0.005	<0.005	0.00%
Cobalt	<0.005	<0.005	0.00%
Copper	<0.005	<0.005	0.00%
Iron	1.14	1.23	7.02%
Lanthanum	<0.005	0.005	200%
Lead	<0.005	<0.005	0.00%
Lithium	0.081	0.081	0.43%
Magnesium	0.388	0.373	3.89%
Manganese	0.024	0.025	0.12%
Molybdenum	0.063	0.064	2.48%
Nickel	<0.005	<0.005	0.00%
Palladium	<0.01	<0.01	0.00%
Phosphorus	0.043	0.040	6.90%
Potassium	-----	-----	-----
Selenium	<0.015	<0.015	0.00%
Silicon	21.1	21.4	1.55%
Silver	<0.005	<0.005	0.00%
Sodium	-----	-----	-----
Strontium	0.028	0.028	0.40%
Sulfur	20.2	19.6	2.61%
Thallium	<0.025	<0.025	0.00%
Thorium	<0.01	<0.01	0.00%
Tin	<0.01	<0.01	0.00%
Titanium	<0.005	<0.005	0.00%
Tungsten	<0.015	<0.015	0.00%
Uranium	<0.1	<0.1	0.00%
Vanadium	<0.005	<0.005	0.00%
Yttrium	<0.005	<0.005	0.00%
Zinc	0.015	0.015	0.34%
Zirconium	<0.005	<0.005	0.00%

SOUTHWEST RESEARCH INSTITUTE

SAMPLE ANALYSIS DATA SHEET

010003

Sample ID

16P-1003-FUA-5

Lab Name: Southwest Research Institute

Client: Division 20

Lab Code: SwRI

Date Received: 11/06/03

Matrix: Water

Project No.: 20.06002.01.141

Lab System ID: 237579

SRR: 25200

TO: 031106-12

Analysis	Sample Result (mg/L)	Reporting Limit (mg/L)
Aluminum	0.340	0.05
Antimony	<0.01	0.01
Arsenic	0.011	0.005
Barium	0.008	0.005
Beryllium	<0.005	0.005
Bismuth	<0.015	0.015
Boron	0.210	0.025
Cadmium	<0.005	0.005
Calcium	1.85	0.05
Chromium	<0.005	0.005
Cobalt	<0.005	0.005
Copper	<0.005	0.005
Iron	0.262	0.05
Lanthanum	<0.005	0.005
Lead	<0.005	0.005
Lithium	0.082	0.005
Magnesium	0.269	0.05
Manganese	0.014	0.005
Molybdenum	0.066	0.01
Nickel	<0.005	0.005
Palladium	<0.01	0.01
Phosphorus	<0.04	0.04
Potassium	1.79	0.2
Selenium	<0.015	0.015
Silicon	20.7	0.05
Silver	<0.005	0.005
Sodium	105	0.2
Strontium	0.022	0.005
Sulfur	19.8	0.025
Thallium	<0.025	0.025
Thorium	<0.01	0.01
Tin	<0.01	0.01
Titanium	<0.005	0.005
Tungsten	<0.015	0.015
Uranium	<0.1	0.1
Vanadium	<0.005	0.005
Yttrium	<0.005	0.005
Zinc	0.007	0.005
Zirconium	<0.005	0.005

SOUTHWEST RESEARCH INSTITUTE

MATRIX SPIKE SUMMARY

010004

Sample ID

16P-1003-FUA-5

Lab Name: Southwest Research Institute

Client: Division 20

Lab Code: SwRI

Date Received: 11/06/03

Matrix: Water

Project No.: 20.06002.01.141

Lab System ID: 237579

SRR: 25200

TO: 031106-12

Analysis	Sample Result (mg/L)	Spike Result (mg/L)	Spike Added (mg/L)	Recovery
Aluminum	0.340	2.30	2.00	98.2%
Antimony	<0.01	0.370	0.500	74.0%
Arsenic	0.011	2.12	2.00	105.4%
Barium	0.008	2.05	2.00	102.0%
Beryllium	<0.005	0.053	0.050	105.9%
Bismuth	NA	NA	NA	NA
Boron	NA	NA	NA	NA
Cadmium	<0.005	0.052	0.050	103.5%
Calcium	1.85	22.6	20.0	103.9%
Chromium	<0.005	0.184	0.200	91.8%
Cobalt	<0.005	0.520	0.500	103.9%
Copper	<0.005	0.259	0.250	103.6%
Iron	0.262	1.23	1.00	97.2%
Lanthanum	NA	NA	NA	NA
Lead	<0.005	0.523	0.500	104.5%
Lithium	NA	NA	NA	NA
Magnesium	0.269	21.2	20.0	104.6%
Manganese	0.014	0.523	0.500	102.0%
Molybdenum	NA	NA	NA	NA
Nickel	<0.005	0.509	0.500	101.9%
Palladium	NA	NA	NA	NA
Phosphorus	NA	NA	NA	NA
Potassium	----	----	----	----
Selenium	<0.015	2.32	2.00	116.0%
Silicon	NA	NA	NA	NA
Silver	<0.005	0.042	0.050	84.5%
Sodium	----	----	----	----
Strontium	NA	NA	NA	NA
Sulfur	NA	NA	NA	NA
Thallium	<0.025	2.28	2.00	114.2%
Thorium	NA	NA	NA	NA
Tin	NA	NA	NA	NA
Titanium	NA	NA	NA	NA
Tungsten	NA	NA	NA	NA
Uranium	NA	NA	NA	NA
Vanadium	<0.005	0.514	0.500	102.9%
Yttrium	NA	NA	NA	NA
Zinc	0.007	0.536	0.500	105.8%
Zirconium	NA	NA	NA	NA

NA- Not Applicable.

SOUTHWEST RESEARCH INSTITUTE

SAMPLE ANALYSIS DATA SHEET

010005

Sample ID

16P-1003-UFUA-5

Lab Name: Southwest Research Institute

Client: Division 20

Lab Code: SwRI

Date Received: 11/06/03

Matrix: Water

Project No.: 20.06002.01.141

Lab System ID: 237580

SRR: 25200

TO: 031106-12

Analysis	Sample Result (mg/L)	Reporting Limit (mg/L)
Aluminum	0.354	0.05
Antimony	<0.01	0.01
Arsenic	0.019	0.005
Barium	0.038	0.005
Beryllium	<0.005	0.005
Bismuth	<0.015	0.015
Boron	0.273	0.025
Cadmium	<0.005	0.005
Calcium	6.53	0.05
Chromium	<0.005	0.005
Cobalt	0.010	0.005
Copper	0.009	0.005
Iron	1.07	0.05
Lanthanum	<0.005	0.005
Lead	0.005	0.005
Lithium	0.112	0.005
Magnesium	0.901	0.05
Manganese	0.054	0.005
Molybdenum	0.063	0.01
Nickel	<0.005	0.005
Palladium	<0.01	0.01
Phosphorus	0.061	0.04
Potassium	1.98	0.2
Selenium	<0.015	0.015
Silicon	27.1	0.05
Silver	<0.005	0.005
Sodium	117	0.2
Strontium	0.078	0.005
Sulfur	27.7	0.025
Thallium	<0.025	0.025
Thorium	<0.01	0.01
Tin	0.011	0.01
Titanium	<0.005	0.005
Tungsten	<0.015	0.015
Uranium	<0.1	0.1
Vanadium	<0.005	0.005
Yttrium	<0.005	0.005
Zinc	0.045	0.005
Zirconium	<0.005	0.005

SOUTHWEST RESEARCH INSTITUTE

DUPLICATE SUMMARY

010006

Sample ID

16P-1003-UFUA-5

Lab Name: Southwest Research Institute

Client: Division 20

Lab Code: SwRI

Date Received: 11/06/03

Matrix: Water

Project No.: 20.06002.01.141

Lab System ID: 237580

SRR: 25200

TO: 031106-12

Analysis	Sample Result (mg/L)	Duplicate Result (mg/L)	RPD
Aluminum	-----	-----	-----
Antimony	-----	-----	-----
Arsenic	-----	-----	-----
Barium	-----	-----	-----
Beryllium	-----	-----	-----
Bismuth	-----	-----	-----
Boron	-----	-----	-----
Cadmium	-----	-----	-----
Calcium	-----	-----	-----
Chromium	-----	-----	-----
Cobalt	-----	-----	-----
Copper	-----	-----	-----
Iron	-----	-----	-----
Lanthanum	-----	-----	-----
Lead	-----	-----	-----
Lithium	-----	-----	-----
Magnesium	-----	-----	-----
Manganese	-----	-----	-----
Molybdenum	-----	-----	-----
Nickel	-----	-----	-----
Palladium	-----	-----	-----
Phosphorus	-----	-----	-----
Potassium	1.98	1.89	4.72%
Selenium	-----	-----	-----
Silicon	-----	-----	-----
Silver	-----	-----	-----
Sodium	117	117	0.31%
Strontium	-----	-----	-----
Sulfur	-----	-----	-----
Thallium	-----	-----	-----
Thorium	-----	-----	-----
Tin	-----	-----	-----
Titanium	-----	-----	-----
Tungsten	-----	-----	-----
Uranium	-----	-----	-----
Vanadium	-----	-----	-----
Yttrium	-----	-----	-----
Zinc	-----	-----	-----
Zirconium	-----	-----	-----

SOUTHWEST RESEARCH INSTITUTE

SAMPLE ANALYSIS DATA SHEET

010007

Sample ID

19D-1003-FA-5

Lab Name: Southwest Research Institute

Client: Division 20

Lab Code: SwRI

Date Received: 11/06/03

Matrix: Water

Project No.: 20.06002.01.141

Lab System ID: 237581

SRR: 25200

TO: 031106-12

Analysis	Sample Result (mg/L)	Reporting Limit (mg/L)
Aluminum	<0.05	0.05
Antimony	<0.01	0.01
Arsenic	0.007	0.005
Barium	<0.005	0.005
Beryllium	<0.005	0.005
Bismuth	<0.015	0.015
Boron	0.232	0.025
Cadmium	<0.005	0.005
Calcium	2.23	0.05
Chromium	<0.005	0.005
Cobalt	<0.005	0.005
Copper	<0.005	0.005
Iron	0.275	0.05
Lanthanum	<0.005	0.005
Lead	<0.005	0.005
Lithium	0.249	0.005
Magnesium	0.103	0.05
Manganese	0.040	0.005
Molybdenum	0.039	0.01
Nickel	<0.005	0.005
Palladium	<0.01	0.01
Phosphorus	<0.04	0.04
Potassium	2.99	0.2
Selenium	<0.015	0.015
Silicon	7.10	0.05
Silver	<0.005	0.005
Sodium	173	0.2
Strontium	0.006	0.005
Sulfur	9.97	0.025
Thallium	<0.025	0.025
Thorium	<0.01	0.01
Tin	0.010	0.01
Titanium	<0.005	0.005
Tungsten	<0.015	0.015
Uranium	<0.1	0.1
Vanadium	<0.005	0.005
Yttrium	<0.005	0.005
Zinc	<0.005	0.005
Zirconium	<0.005	0.005

SOUTHWEST RESEARCH INSTITUTE

MATRIX SPIKE SUMMARY

010008

Sample ID

19D-1003-FA-5

Lab Name: Southwest Research Institute

Client: Division 20

Lab Code: SwRI

Date Received: 11/06/03

Matrix: Water

Project No.: 20.06002.01.141

Lab System ID: 237581

SRR: 25200

TO: 031106-12

Analysis	Sample Result (mg/L)	Spike Result (mg/L)	Spike Added (mg/L)	Recovery
Aluminum	-----	-----	-----	-----
Antimony	-----	-----	-----	-----
Arsenic	-----	-----	-----	-----
Barium	-----	-----	-----	-----
Beryllium	-----	-----	-----	-----
Bismuth	-----	-----	-----	-----
Boron	-----	-----	-----	-----
Cadmium	-----	-----	-----	-----
Calcium	-----	-----	-----	-----
Chromium	-----	-----	-----	-----
Cobalt	-----	-----	-----	-----
Copper	-----	-----	-----	-----
Iron	-----	-----	-----	-----
Lanthanum	-----	-----	-----	-----
Lead	-----	-----	-----	-----
Lithium	-----	-----	-----	-----
Magnesium	-----	-----	-----	-----
Manganese	-----	-----	-----	-----
Molybdenum	-----	-----	-----	-----
Nickel	-----	-----	-----	-----
Palladium	-----	-----	-----	-----
Phosphorus	-----	-----	-----	-----
Potassium	2.99	23.0	20.0	100.1%
Selenium	-----	-----	-----	-----
Silicon	-----	-----	-----	-----
Silver	-----	-----	-----	-----
Sodium	173	189	20.0	84.5%
Strontium	-----	-----	-----	-----
Sulfur	-----	-----	-----	-----
Thallium	-----	-----	-----	-----
Thorium	-----	-----	-----	-----
Tin	-----	-----	-----	-----
Titanium	-----	-----	-----	-----
Tungsten	-----	-----	-----	-----
Uranium	-----	-----	-----	-----
Vanadium	-----	-----	-----	-----
Yttrium	-----	-----	-----	-----
Zinc	-----	-----	-----	-----
Zirconium	-----	-----	-----	-----

NA- Not Applicable.

SOUTHWEST RESEARCH INSTITUTE

SAMPLE ANALYSIS DATA SHEET

010009

Sample ID

19D-1003-FUA-5

Lab Name: Southwest Research Institute

Client: Division 20

Lab Code: SwRI

Date Received: 11/06/03

Matrix: Water

Project No.: 20.06002.01.141

Lab System ID: 237582

SRR: 25200

TO: 031106-12

Analysis	Sample Result (mg/L)	Reporting Limit (mg/L)
Aluminum	<0.05	0.05
Antimony	<0.01	0.01
Arsenic	0.008	0.005
Barium	<0.005	0.005
Beryllium	<0.005	0.005
Bismuth	<0.015	0.015
Boron	0.233	0.025
Cadmium	<0.005	0.005
Calcium	2.24	0.05
Chromium	<0.005	0.005
Cobalt	<0.005	0.005
Copper	<0.005	0.005
Iron	0.202	0.05
Lanthanum	<0.005	0.005
Lead	<0.005	0.005
Lithium	0.247	0.005
Magnesium	0.106	0.05
Manganese	0.039	0.005
Molybdenum	0.038	0.01
Nickel	<0.005	0.005
Palladium	<0.01	0.01
Phosphorus	<0.04	0.04
Potassium	3.08	0.2
Selenium	<0.015	0.015
Silicon	7.13	0.05
Silver	<0.005	0.005
Sodium	173	0.2
Strontium	0.006	0.005
Sulfur	9.89	0.025
Thallium	<0.025	0.025
Thorium	<0.01	0.01
Tin	0.025	0.01
Titanium	<0.005	0.005
Tungsten	<0.015	0.015
Uranium	<0.1	0.1
Vanadium	<0.005	0.005
Yttrium	<0.005	0.005
Zinc	<0.005	0.005
Zirconium	<0.005	0.005

SOUTHWEST RESEARCH INSTITUTE

SAMPLE ANALYSIS DATA SHEET

010010

Sample ID

19D-1003-UFUA-5

Lab Name: Southwest Research Institute

Client: Division 20

Lab Code: SwRI

Date Received: 11/06/03

Matrix: Water

Project No.: 20.06002.01.141

Lab System ID: 237583

SRR: 25200

TO: 031106-12

Analysis	Sample Result (mg/L)	Reporting Limit (mg/L)
Aluminum	0.096	0.05
Antimony	<0.01	0.01
Arsenic	0.007	0.005
Barium	<0.005	0.005
Beryllium	<0.005	0.005
Bismuth	<0.015	0.015
Boron	0.234	0.025
Cadmium	<0.005	0.005
Calcium	2.18	0.05
Chromium	<0.005	0.005
Cobalt	<0.005	0.005
Copper	<0.005	0.005
Iron	1.33	0.05
Lanthanum	<0.005	0.005
Lead	<0.005	0.005
Lithium	0.256	0.005
Magnesium	0.111	0.05
Manganese	0.049	0.005
Molybdenum	0.035	0.01
Nickel	<0.005	0.005
Palladium	<0.01	0.01
Phosphorus	0.046	0.04
Potassium	2.93	0.2
Selenium	<0.015	0.015
Silicon	7.80	0.05
Silver	<0.005	0.005
Sodium	173	0.2
Strontium	0.007	0.005
Sulfur	9.85	0.025
Thallium	<0.025	0.025
Thorium	<0.01	0.01
Tin	0.024	0.01
Titanium	<0.005	0.005
Tungsten	<0.015	0.015
Uranium	<0.1	0.1
Vanadium	<0.005	0.005
Yttrium	<0.005	0.005
Zinc	0.014	0.005
Zirconium	<0.005	0.005

SOUTHWEST RESEARCH INSTITUTE

SAMPLE ANALYSIS DATA SHEET

010011

Sample ID

19IM1-25-1003-FA-5

Lab Name: Southwest Research Institute

Client: Division 20

Lab Code: SwRI

Date Received: 11/06/03

Matrix: Water

Project No.: 20.06002.01.141

Lab System ID: 237584

SRR: 25200

TO: 031106-12

Analysis	Sample Result (mg/L)	Reporting Limit (mg/L)
Aluminum	0.167	0.05
Antimony	<0.01	0.01
Arsenic	0.035	0.005
Barium	<0.005	0.005
Beryllium	<0.005	0.005
Bismuth	<0.015	0.015
Boron	0.179	0.025
Cadmium	<0.005	0.005
Calcium	0.543	0.05
Chromium	<0.005	0.005
Cobalt	<0.005	0.005
Copper	<0.005	0.005
Iron	0.074	0.05
Lanthanum	<0.005	0.005
Lead	<0.005	0.005
Lithium	0.125	0.005
Magnesium	0.071	0.05
Manganese	0.015	0.005
Molybdenum	<0.01	0.01
Nickel	<0.005	0.005
Palladium	<0.01	0.01
Phosphorus	<0.04	0.04
Potassium	3.04	0.2
Selenium	<0.015	0.015
Silicon	27.9	0.05
Silver	<0.005	0.005
Sodium	94.7	0.2
Strontium	<0.005	0.005
Sulfur	5.56	0.025
Thallium	<0.025	0.025
Thorium	<0.01	0.01
Tin	<0.01	0.01
Titanium	<0.005	0.005
Tungsten	<0.015	0.015
Uranium	<0.1	0.1
Vanadium	0.005	0.005
Yttrium	<0.005	0.005
Zinc	<0.005	0.005
Zirconium	<0.005	0.005

SOUTHWEST RESEARCH INSTITUTE

SAMPLE ANALYSIS DATA SHEET

010012

Sample ID

19IM1-25-1003-FUA-5

Lab Name: Southwest Research Institute

Client: Division 20

Lab Code: SwRI

Date Received: 11/06/03

Matrix: Water

Project No.: 20.06002.01.141

Lab System ID: 237585

SRR: 25200

TO: 031106-12

Analysis	Sample Result (mg/L)	Reporting Limit (mg/L)
Aluminum	0.132	0.05
Antimony	<0.01	0.01
Arsenic	0.034	0.005
Barium	<0.005	0.005
Beryllium	<0.005	0.005
Bismuth	<0.015	0.015
Boron	0.162	0.025
Cadmium	<0.005	0.005
Calcium	0.538	0.05
Chromium	<0.005	0.005
Cobalt	<0.005	0.005
Copper	<0.005	0.005
Iron	<0.05	0.05
Lanthanum	<0.005	0.005
Lead	<0.005	0.005
Lithium	0.125	0.005
Magnesium	0.068	0.05
Manganese	0.014	0.005
Molybdenum	<0.01	0.01
Nickel	<0.005	0.005
Palladium	<0.01	0.01
Phosphorus	0.040	0.04
Potassium	2.92	0.2
Selenium	<0.015	0.015
Silicon	28.2	0.05
Silver	<0.005	0.005
Sodium	94.4	0.2
Strontium	<0.005	0.005
Sulfur	5.48	0.025
Thallium	<0.025	0.025
Thorium	<0.01	0.01
Tin	<0.01	0.01
Titanium	<0.005	0.005
Tungsten	<0.015	0.015
Uranium	<0.1	0.1
Vanadium	<0.005	0.005
Yttrium	<0.005	0.005
Zinc	<0.005	0.005
Zirconium	<0.005	0.005

SOUTHWEST RESEARCH INSTITUTE.

SAMPLE ANALYSIS DATA SHEET

010013

Sample ID

19IM1-25-1003-UFUA-5

Lab Name: Southwest Research Institute

Client: Division 20

Lab Code: SwRI

Date Received: 11/06/03

Matrix: Water

Project No.: 20.06002.01.141

Lab System ID: 237586

SRR: 25200

TO: 031106-12

Analysis	Sample Result (mg/L)	Reporting Limit (mg/L)
Aluminum	0.204	0.05
Antimony	<0.01	0.01
Arsenic	0.033	0.005
Barium	<0.005	0.005
Beryllium	<0.005	0.005
Bismuth	<0.015	0.015
Boron	0.163	0.025
Cadmium	<0.005	0.005
Calcium	0.761	0.05
Chromium	<0.005	0.005
Cobalt	<0.005	0.005
Copper	<0.005	0.005
Iron	0.195	0.05
Lanthanum	<0.005	0.005
Lead	<0.005	0.005
Lithium	0.128	0.005
Magnesium	0.094	0.05
Manganese	0.025	0.005
Molybdenum	<0.01	0.01
Nickel	<0.005	0.005
Palladium	<0.01	0.01
Phosphorus	0.045	0.04
Potassium	3.04	0.2
Selenium	<0.015	0.015
Silicon	28.7	0.05
Silver	<0.005	0.005
Sodium	94.8	0.2
Strontium	<0.005	0.005
Sulfur	5.63	0.025
Thallium	<0.025	0.025
Thorium	<0.01	0.01
Tin	<0.01	0.01
Titanium	<0.005	0.005
Tungsten	<0.015	0.015
Uranium	<0.1	0.1
Vanadium	<0.005	0.005
Yttrium	<0.005	0.005
Zinc	0.019	0.005
Zirconium	<0.005	0.005

SOUTHWEST RESEARCH INSTITUTE

SAMPLE ANALYSIS DATA SHEET

010014

Sample ID

19IM2-25-1003-FA-5

Lab Name: Southwest Research Institute

Client: Division 20

Lab Code: SwRI

Date Received: 11/06/03

Matrix: Water

Project No.: 20.06002.01.141

Lab System ID: 237587

SRR: 25200

TO: 031106-12

Analysis	Sample Result (mg/L)	Reporting Limit (mg/L)
Aluminum	0.158	0.05
Antimony	<0.01	0.01
Arsenic	0.037	0.005
Barium	<0.005	0.005
Beryllium	<0.005	0.005
Bismuth	<0.015	0.015
Boron	0.157	0.025
Cadmium	<0.005	0.005
Calcium	0.360	0.05
Chromium	<0.005	0.005
Cobalt	<0.005	0.005
Copper	<0.005	0.005
Iron	0.074	0.05
Lanthanum	<0.005	0.005
Lead	<0.005	0.005
Lithium	0.126	0.005
Magnesium	<0.05	0.05
Manganese	<0.005	0.005
Molybdenum	<0.01	0.01
Nickel	<0.005	0.005
Palladium	<0.01	0.01
Phosphorus	0.040	0.04
Potassium	2.96	0.2
Selenium	<0.015	0.015
Silicon	28.5	0.05
Silver	<0.005	0.005
Sodium	95.6	0.2
Strontium	<0.005	0.005
Sulfur	5.69	0.025
Thallium	<0.025	0.025
Thorium	<0.01	0.01
Tin	<0.01	0.01
Titanium	<0.005	0.005
Tungsten	<0.015	0.015
Uranium	<0.1	0.1
Vanadium	<0.005	0.005
Yttrium	<0.005	0.005
Zinc	<0.005	0.005
Zirconium	<0.005	0.005

SOUTHWEST RESEARCH INSTITUTE

SAMPLE ANALYSIS DATA SHEET

010015

Sample ID

19IM2-25-1003-FUA-5

Lab Name: Southwest Research Institute

Client: Division 20

Lab Code: SwRI

Date Received: 11/06/03

Matrix: Water

Project No.: 20.06002.01.141

Lab System ID: 237588

SRR: 25200

TO: 031106-12

Analysis	Sample Result (mg/L)	Reporting Limit (mg/L)
Aluminum	0.120	0.05
Antimony	<0.01	0.01
Arsenic	0.035	0.005
Barium	<0.005	0.005
Beryllium	<0.005	0.005
Bismuth	<0.015	0.015
Boron	0.157	0.025
Cadmium	<0.005	0.005
Calcium	0.359	0.05
Chromium	<0.005	0.005
Cobalt	<0.005	0.005
Copper	<0.005	0.005
Iron	<0.05	0.05
Lanthanum	<0.005	0.005
Lead	<0.005	0.005
Lithium	0.127	0.005
Magnesium	<0.05	0.05
Manganese	<0.005	0.005
Molybdenum	<0.01	0.01
Nickel	<0.005	0.005
Palladium	<0.01	0.01
Phosphorus	<0.04	0.04
Potassium	3.09	0.2
Selenium	<0.015	0.015
Silicon	28.5	0.05
Silver	<0.005	0.005
Sodium	97.1	0.2
Strontium	<0.005	0.005
Sulfur	5.81	0.025
Thallium	<0.025	0.025
Thorium	<0.01	0.01
Tin	<0.01	0.01
Titanium	<0.005	0.005
Tungsten	<0.015	0.015
Uranium	<0.1	0.1
Vanadium	0.005	0.005
Yttrium	<0.005	0.005
Zinc	<0.005	0.005
Zirconium	<0.005	0.005

SOUTHWEST RESEARCH INSTITUTE

SAMPLE ANALYSIS DATA SHEET

010016

Sample ID

19IM2-25-1003-UFUA-5

Lab Name: Southwest Research Institute

Client: Division 20

Lab Code: SwRI

Date Received: 11/06/03

Matrix: Water

Project No.: 20.06002.01.141

Lab System ID: 237589

SRR: 25200

TO: 031106-12

Analysis	Sample Result (mg/L)	Reporting Limit (mg/L)
Aluminum	0.174	0.05
Antimony	<0.01	0.01
Arsenic	0.032	0.005
Barium	<0.005	0.005
Beryllium	<0.005	0.005
Bismuth	<0.015	0.015
Boron	0.157	0.025
Cadmium	<0.005	0.005
Calcium	0.409	0.05
Chromium	<0.005	0.005
Cobalt	<0.005	0.005
Copper	<0.005	0.005
Iron	0.190	0.05
Lanthanum	<0.005	0.005
Lead	<0.005	0.005
Lithium	0.127	0.005
Magnesium	<0.05	0.05
Manganese	0.034	0.005
Molybdenum	<0.01	0.01
Nickel	<0.005	0.005
Palladium	<0.01	0.01
Phosphorus	<0.04	0.04
Potassium	3.11	0.2
Selenium	<0.015	0.015
Silicon	28.7	0.05
Silver	<0.005	0.005
Sodium	97.3	0.2
Strontium	<0.005	0.005
Sulfur	5.69	0.025
Thallium	<0.025	0.025
Thorium	<0.01	0.01
Tin	<0.01	0.01
Titanium	<0.005	0.005
Tungsten	<0.015	0.015
Uranium	<0.1	0.1
Vanadium	0.005	0.005
Yttrium	<0.005	0.005
Zinc	0.006	0.005
Zirconium	<0.005	0.005

SOUTHWEST RESEARCH INSTITUTE

SAMPLE ANALYSIS DATA SHEET

010017

Sample ID

7SC-24-1003-FA-5

Lab Name: Southwest Research Institute

Client: Division 20

Lab Code: SwRI

Date Received: 11/06/03

Matrix: Water

Project No.: 20.06002.01.141

Lab System ID: 237590

SRR: 25200

TO: 031106-12

Analysis	Sample Result (mg/L)	Reporting Limit (mg/L)
Aluminum	<0.05	0.05
Antimony	<0.01	0.01
Arsenic	<0.005	0.005
Barium	0.056	0.005
Beryllium	<0.005	0.005
Bismuth	<0.015	0.015
Boron	0.259	0.025
Cadmium	<0.005	0.005
Calcium	28.1	0.05
Chromium	<0.005	0.005
Cobalt	<0.005	0.005
Copper	<0.005	0.005
Iron	<0.05	0.05
Lanthanum	<0.005	0.005
Lead	<0.005	0.005
Lithium	0.105	0.005
Magnesium	28.9	0.05
Manganese	0.075	0.005
Molybdenum	<0.01	0.01
Nickel	<0.005	0.005
Palladium	<0.01	0.01
Phosphorus	<0.04	0.04
Potassium	8.08	0.2
Selenium	<0.015	0.015
Silicon	13.7	0.05
Silver	<0.005	0.005
Sodium	82.2	0.2
Strontium	0.413	0.005
Sulfur	45.0	0.025
Thallium	<0.025	0.025
Thorium	<0.01	0.01
Tin	<0.01	0.01
Titanium	<0.005	0.005
Tungsten	<0.015	0.015
Uranium	<0.1	0.1
Vanadium	<0.005	0.005
Yttrium	<0.005	0.005
Zinc	<0.005	0.005
Zirconium	<0.005	0.005

SOUTHWEST RESEARCH INSTITUTE

SAMPLE ANALYSIS DATA SHEET

010018

Sample ID

7SC-24-1003-FUA-5

Lab Name: Southwest Research Institute

Client: Division 20

Lab Code: SwRI

Date Received: 11/06/03

Matrix: Water

Project No.: 20.06002.01.141

Lab System ID: 237591

SRR: 25200

TO: 031106-12

Analysis	Sample Result (mg/L)	Reporting Limit (mg/L)
Aluminum	<0.05	0.05
Antimony	<0.01	0.01
Arsenic	0.005	0.005
Barium	0.056	0.005
Beryllium	<0.005	0.005
Bismuth	<0.015	0.015
Boron	0.260	0.025
Cadmium	<0.005	0.005
Calcium	28.3	0.05
Chromium	<0.005	0.005
Cobalt	0.005	0.005
Copper	<0.005	0.005
Iron	<0.05	0.05
Lanthanum	<0.005	0.005
Lead	<0.005	0.005
Lithium	0.104	0.005
Magnesium	28.9	0.05
Manganese	0.071	0.005
Molybdenum	<0.01	0.01
Nickel	<0.005	0.005
Palladium	<0.01	0.01
Phosphorus	0.042	0.04
Potassium	8.16	0.2
Selenium	<0.015	0.015
Silicon	13.7	0.05
Silver	<0.005	0.005
Sodium	82.0	0.2
Strontium	0.413	0.005
Sulfur	43.7	0.025
Thallium	<0.025	0.025
Thorium	<0.01	0.01
Tin	0.016	0.01
Titanium	<0.005	0.005
Tungsten	<0.015	0.015
Uranium	<0.1	0.1
Vanadium	<0.005	0.005
Yttrium	<0.005	0.005
Zinc	<0.005	0.005
Zirconium	<0.005	0.005

SOUTHWEST RESEARCH INSTITUTE

SAMPLE ANALYSIS DATA SHEET

010019

Sample ID

7SC-24-1003-UFUA-5

Lab Name: Southwest Research Institute

Client: Division 20

Lab Code: SwRI

Date Received: 11/06/03

Matrix: Water

Project No.: 20.06002.01.141

Lab System ID: 237592

SRR: 25200

TO: 031106-12

Analysis	Sample Result (mg/L)	Reporting Limit (mg/L)
Aluminum	<0.05	0.05
Antimony	<0.01	0.01
Arsenic	<0.005	0.005
Barium	0.057	0.005
Beryllium	<0.005	0.005
Bismuth	<0.015	0.015
Boron	0.260	0.025
Cadmium	<0.005	0.005
Calcium	28.2	0.05
Chromium	<0.005	0.005
Cobalt	<0.005	0.005
Copper	<0.005	0.005
Iron	0.081	0.05
Lanthanum	<0.005	0.005
Lead	<0.005	0.005
Lithium	0.105	0.005
Magnesium	28.9	0.05
Manganese	0.071	0.005
Molybdenum	<0.01	0.01
Nickel	<0.005	0.005
Palladium	<0.01	0.01
Phosphorus	<0.04	0.04
Potassium	8.12	0.2
Selenium	<0.015	0.015
Silicon	13.7	0.05
Silver	<0.005	0.005
Sodium	81.6	0.2
Strontium	0.416	0.005
Sulfur	44.3	0.025
Thallium	<0.025	0.025
Thorium	<0.01	0.01
Tin	<0.01	0.01
Titanium	<0.005	0.005
Tungsten	<0.015	0.015
Uranium	<0.1	0.1
Vanadium	<0.005	0.005
Yttrium	<0.005	0.005
Zinc	0.016	0.005
Zirconium	<0.005	0.005

SOUTHWEST RESEARCH INSTITUTE

SAMPLE ANALYSIS DATA SHEET

010020

Sample ID

98Q-1003-FA-5

Lab Name: Southwest Research Institute

Client: Division 20

Lab Code: SwRI

Date Received: 11/06/03

Matrix: Water

Project No.: 20.06002.01.141

Lab System ID: 237593

SRR: 25200

TO: 031106-12

Analysis	Sample Result (mg/L)	Reporting Limit (mg/L)
Aluminum	<0.05	0.05
Antimony	<0.01	0.01
Arsenic	<0.005	0.005
Barium	<0.005	0.005
Beryllium	<0.005	0.005
Bismuth	<0.015	0.015
Boron	<0.025	0.025
Cadmium	<0.005	0.005
Calcium	10.1	0.05
Chromium	<0.005	0.005
Cobalt	<0.005	0.005
Copper	<0.005	0.005
Iron	<0.05	0.05
Lanthanum	<0.005	0.005
Lead	<0.005	0.005
Lithium	<0.005	0.005
Magnesium	0.065	0.05
Manganese	<0.005	0.005
Molybdenum	<0.01	0.01
Nickel	<0.005	0.005
Palladium	<0.01	0.01
Phosphorus	<0.04	0.04
Potassium	5.25	0.2
Selenium	<0.015	0.015
Silicon	19.4	0.05
Silver	<0.005	0.005
Sodium	27.0	0.2
Strontium	<0.005	0.005
Sulfur	0.114	0.025
Thallium	<0.025	0.025
Thorium	<0.01	0.01
Tin	<0.01	0.01
Titanium	<0.005	0.005
Tungsten	<0.015	0.015
Uranium	<0.1	0.1
Vanadium	<0.005	0.005
Yttrium	<0.005	0.005
Zinc	<0.005	0.005
Zirconium	0.008	0.005

SOUTHWEST RESEARCH INSTITUTE

SAMPLE ANALYSIS DATA SHEET

010021

Sample ID

98Q-1003-FA-6

Lab Name: Southwest Research Institute

Client: Division 20

Lab Code: SwRI

Date Received: 11/06/03

Matrix: Water

Project No.: 20.06002.01.141

Lab System ID: 237594

SRR: 25200

TO: 031106-12

Analysis	Sample Result (mg/L)	Reporting Limit (mg/L)
Aluminum	<0.05	0.05
Antimony	<0.01	0.01
Arsenic	<0.005	0.005
Barium	<0.005	0.005
Beryllium	<0.005	0.005
Bismuth	<0.015	0.015
Boron	<0.025	0.025
Cadmium	<0.005	0.005
Calcium	9.91	0.05
Chromium	<0.005	0.005
Cobalt	<0.005	0.005
Copper	<0.005	0.005
Iron	<0.05	0.05
Lanthanum	<0.005	0.005
Lead	<0.005	0.005
Lithium	<0.005	0.005
Magnesium	<0.05	0.05
Manganese	<0.005	0.005
Molybdenum	<0.01	0.01
Nickel	<0.005	0.005
Palladium	<0.01	0.01
Phosphorus	<0.04	0.04
Potassium	5.32	0.2
Selenium	<0.015	0.015
Silicon	19.4	0.05
Silver	<0.005	0.005
Sodium	26.6	0.2
Strontium	<0.005	0.005
Sulfur	<0.025	0.025
Thallium	<0.025	0.025
Thorium	<0.01	0.01
Tin	<0.01	0.01
Titanium	<0.005	0.005
Tungsten	<0.015	0.015
Uranium	<0.1	0.1
Vanadium	<0.005	0.005
Yttrium	<0.005	0.005
Zinc	<0.005	0.005
Zirconium	0.005	0.005

SOUTHWEST RESEARCH INSTITUTE

LABORATORY CONTROL SAMPLE

010022

Sample ID

LCSW - L11W1 / L24W1

Lab Name: Southwest Research Institute

Client: Division 20

Lab Code: SwRI

Date Received: NA

Matrix: Water

Project No.: 20.06002.01.141

Lab System ID: NA

SRR: 25200

TO: 031106-12

Analysis	Sample Result (mg/L)	True Value (mg/L)	Recovery
Aluminum	1.92	2.00	95.9%
Antimony	0.515	0.500	102.9%
Arsenic	2.08	2.00	104.1%
Barium	2.02	2.00	101.2%
Beryllium	0.054	0.050	108.4%
Bismuth	NA	NA	NA
Boron	NA	NA	NA
Cadmium	0.052	0.050	104.9%
Calcium	20.8	20.0	104.0%
Chromium	0.203	0.200	101.7%
Cobalt	0.515	0.500	103.0%
Copper	0.251	0.250	100.4%
Iron	1.03	1.00	103.1%
Lanthanum	NA	NA	NA
Lead	0.527	0.500	105.4%
Lithium	NA	NA	NA
Magnesium	20.4	20.0	101.9%
Manganese	0.516	0.500	103.1%
Molybdenum	NA	NA	NA
Nickel	0.511	0.500	102.2%
Palladium	NA	NA	NA
Phosphorus	NA	NA	NA
Potassium	18.3	20.0	91.4%
Selenium	2.15	2.00	107.3%
Silicon	NA	NA	NA
Silver	0.052	0.050	104.4%
Sodium	18.6	20.0	93.2%
Strontium	NA	NA	NA
Sulfur	NA	NA	NA
Thallium	2.22	2.00	110.9%
Thorium	NA	NA	NA
Tin	NA	NA	NA
Titanium	NA	NA	NA
Tungsten	NA	NA	NA
Uranium	NA	NA	NA
Vanadium	0.510	0.500	102.0%
Yttrium	NA	NA	NA
Zinc	0.531	0.500	106.2%
Zirconium	NA	NA	NA

NA- Not Applicable.

SOUTHWEST RESEARCH INSTITUTE

BLANK SUMMARY

010023

Sample ID

PBW - L11W1 / L24W2

Lab Name: Southwest Research Institute

Client: Division 20

Lab Code: SwRI

Date Received: NA

Matrix: Water

Project No.: 20.06002.01.141

Lab System ID: NA

SRR: 25200

TO: 031106-12

Analysis	Sample Result (mg/L)	Reporting Limit (mg/L)
Aluminum	<0.05	0.05
Antimony	<0.01	0.01
Arsenic	<0.005	0.005
Barium	<0.005	0.005
Beryllium	<0.005	0.005
Bismuth	<0.015	0.015
Boron	<0.025	0.025
Cadmium	<0.005	0.005
Calcium	<0.05	0.05
Chromium	<0.005	0.005
Cobalt	<0.005	0.005
Copper	<0.005	0.005
Iron	<0.05	0.05
Lanthanum	<0.005	0.005
Lead	<0.005	0.005
Lithium	<0.005	0.005
Magnesium	<0.05	0.05
Manganese	<0.005	0.005
Molybdenum	<0.01	0.01
Nickel	<0.005	0.005
Palladium	<0.01	0.01
Phosphorus	<0.04	0.04
Potassium	<0.2	0.2
Selenium	<0.015	0.015
Silicon	<0.05	0.05
Silver	<0.005	0.005
Sodium	<0.2	0.2
Strontium	<0.005	0.005
Sulfur	<0.025	0.025
Thallium	<0.025	0.025
Thorium	<0.01	0.01
Tin	<0.01	0.01
Titanium	<0.005	0.005
Tungsten	<0.015	0.015
Uranium	<0.1	0.1
Vanadium	<0.005	0.005
Yttrium	<0.005	0.005
Zinc	<0.005	0.005
Zirconium	<0.005	0.005

NA- Not Applicable.

PAGE 23 OF 23

This report may not be reproduced except in its entirety without the written approval of SwRI.

**SOUTHWEST RESEARCH INSTITUTE
NUCLEAR PROJECT
CLIENT: Div. 20
TASK ORDER: 031106-12
SRR: 25200
SDG: 237578
CASE: CNWRA
VTSR: November 06, 2003
PROJECT#: 06002.01.141**

Task Orders/01-QPP-015

Laboratory Task Order

TO #: 031106-12 Revision: 1

010024SDG: 237578
VTSR: 11/06/03
CASE: DIV 20 - CNWRASRR #s: 25200
Client(s): DIV 20Project(s): 06002.01.141
Manager(s): DAMMANN, MIKE
To PM: 12/02/03
To QA: 12/31/03
To Client: 01/02/04**Instructions**

DIVISION 20 - CNWRA. 28-day TAT listed on COC (subject to change). Point of Contact is Bradley Werling (ext. 6565). Work is 10 CFR 50, Part 21, Appendix B. Contact Charlie Butcher (ext. 5928, pager 271-5172) before starting ANY WORK on this task order. ICP analysis is for MAJOR and MINOR elements. CONTACT PM WITH ANY QUESTIONS.
revision 1: updating task order. (dri23003)

Documents Related to this task order: 7451

Test: DIL-DILUTION
Section: METALPREP

Holding: 28 days from CED

Prep, Dilution

Cnt: 17

System ID	Type	Cont	Matrix	Customer ID	CED	Method Date
237578		1	Water	16P-1003-FA-5	05 Nov 03	03 Dec 03
237579		1	Water	16P-1003-FUA-5	05 Nov 03	03 Dec 03
237580		1	Water	16P-1003-UFUA-5	05 Nov 03	03 Dec 03
237581		1	Water	19D-1003-FA-5	05 Nov 03	03 Dec 03
237582		1	Water	19D-1003-FUA-5	05 Nov 03	03 Dec 03
237583		1	Water	19D-1003-UFUA-5	05 Nov 03	03 Dec 03
237584		1	Water	19IM1-25-1003-FA-5	05 Nov 03	03 Dec 03
237585		1	Water	19IM1-25-1003-FUA-5	05 Nov 03	03 Dec 03
237586		1	Water	19IM1-25-1003-UFUA-5	05 Nov 03	03 Dec 03
237587		1	Water	19IM2-25-1003-FA-5	05 Nov 03	03 Dec 03
237588		1	Water	19IM2-25-1003-FUA-5	05 Nov 03	03 Dec 03
237589		1	Water	19IM2-25-1003-UFUA-5	05 Nov 03	03 Dec 03
237590		1	Water	7SC-24-1003-FA-5	05 Nov 03	03 Dec 03
237591		1	Water	7SC-24-1003-FUA-5	05 Nov 03	03 Dec 03
237592		1	Water	7SC-24-1003-UFUA-5	05 Nov 03	03 Dec 03
237593		1	Water	98Q-1003-FA-5	05 Nov 03	03 Dec 03
237594		1	Water	98Q-1003-FA-6	05 Nov 03	03 Dec 03

Test: ICP-SWRI
Section: METALS

Holding: 180 days from CED

ICP Analysis by SwRI Method

Cnt: 17

System ID	Type	Cont	Matrix	Customer ID	CED	Method Date
237578		1	Water	16P-1003-FA-5	05 Nov 03	03 May 04
237579		1	Water	16P-1003-FUA-5	05 Nov 03	03 May 04
237580		1	Water	16P-1003-UFUA-5	05 Nov 03	03 May 04
237581		1	Water	19D-1003-FA-5	05 Nov 03	03 May 04
237582		1	Water	19D-1003-FUA-5	05 Nov 03	03 May 04
237583		1	Water	19D-1003-UFUA-5	05 Nov 03	03 May 04
237584		1	Water	19IM1-25-1003-FA-5	05 Nov 03	03 May 04
237585		1	Water	19IM1-25-1003-FUA-5	05 Nov 03	03 May 04
237586		1	Water	19IM1-25-1003-UFUA-5	05 Nov 03	03 May 04
237587		1	Water	19IM2-25-1003-FA-5	05 Nov 03	03 May 04
237588		1	Water	19IM2-25-1003-FUA-5	05 Nov 03	03 May 04
237589		1	Water	19IM2-25-1003-UFUA-5	05 Nov 03	03 May 04
237590		1	Water	7SC-24-1003-FA-5	05 Nov 03	03 May 04
237591		1	Water	7SC-24-1003-FUA-5	05 Nov 03	03 May 04
237592		1	Water	7SC-24-1003-UFUA-5	05 Nov 03	03 May 04
237593		1	Water	98Q-1003-FA-5	05 Nov 03	03 May 04
237594		1	Water	98Q-1003-FA-6	05 Nov 03	03 May 04

Southwest Research Institute

Laboratory Task Order

TO #: 031106-12 Revision: 1

SDG: 237578
VTSR: 11/06/03
CASE: DIV 20 - CNWRA

SRR #'s: 25200
Client(s): DIV 20

010028
Project(s): 06002.01.141
Manager(s): DAMMANN, MIKE
To PM: 12/02/03
To QA: 12/31/03
To Client: 01/02/04

010028

01-QPP-015
Division 01
Revision 4
November 2002

CONTROLLED COPY
IF STAMP IS NOT RED, THIS DOCUMENT IS UNCONTROLLED

Document No. 3



Chemistry and Chemical
Engineering Division

QUALITY PROJECT PLAN FOR

**PERFORMANCE OF CHEMICAL ANALYSES
FOR COMMERCIAL NUCLEAR POWER PLANTS
WITHIN THE DEPARTMENT OF ANALYTICAL
AND ENVIRONMENTAL CHEMISTRY**


SOUTHWEST RESEARCH INSTITUTE
Chemistry and Chemical Engineering Division
6220 CULEBRA ROAD, SAN ANTONIO, TEXAS 78238

810027


**QUALITY PROJECT PLAN FOR PERFORMANCE OF CHEMICAL ANALYSES
FOR COMMERCIAL NUCLEAR POWER PLANTS
WITHIN THE DEPARTMENT OF ANALYTICAL AND ENVIRONMENTAL CHEMISTRY**

SwRI AUTHORIZATION SIGNATORIES


This is to certify that this Quality Project Plan of Southwest Research Institute (SwRI) has been reviewed and approved by the following personnel:


JOANN BOYD (210) 522-2169
Quality Assurance Manager


10/30/02
DATE


REZA KARIMI (210) 522-2412
Director, Department of Analytical and Environmental Chemistry

10/30/02
DATE


for MICHAEL G. MACNAUGHTON (210) 522-5162
Vice President, Chemistry and Chemical Engineering Division

10/30/02
DATE


CHRISTOPHER HOBSON (210) 522-5838
Quality Assurance Engineer

10/30/02
DATE

TABLE OF CONTENTS

	Page
1.0 INTRODUCTION	1
2.0 SCOPE.....	1
3.0 REFERENCES.....	1
4.0 APPLICABLE SECTIONS OF SwRI PROGRAM QUALITY PLAN (PQP-NUCLEAR)	1
4.1 Indoctrination and Training	1
4.2 Qualification of Personnel.....	2
4.3 Design Control.....	2
4.4 Right of Access	2
4.5 Control of Supplier-Generated Documents	2
4.6 Acceptance of Services Only.....	2
4.7 Commercial Grade Items.....	3
4.8 Inspection.....	4
4.9 Inspection and Testing	5
4.10 Handling, Storage, Packaging, Preservation, and Delivery	5
4.11 Quality Assurance Records	5
4.12 10 CFR, Part 21	6
4.13 Certified Test Report.....	6
4.14 Valid Documents List.....	6
5.0 HISTORY OF REVISIONS	6

**PERFORMANCE OF CHEMICAL ANALYSES
FOR COMMERCIAL NUCLEAR POWER PLANTS WITHIN THE
DEPARTMENT OF ANALYTICAL AND ENVIRONMENTAL CHEMISTRY**

1.0 INTRODUCTION

This Quality Project Plan (QPP) defines the Quality Assurance (QA) program requirements for personnel providing the chemical analyses for commercial nuclear power plants. Southwest Research Institute (SwRI) **Program Quality Plan (PQP-Nuclear)**, *Nuclear Services* shall implement the QA requirements. Project activities controlled by the PQP-Nuclear shall be accomplished as specified by the appropriate sections of **01-QAP-004**, *Quality Assurance Plan for Analytical and Environmental Services* and/or nationally recognized testing methods as specified on individual purchase orders. This QPP shall be applied to all projects initiated for nuclear utilities in the Department of Analytical and Environmental Chemistry. If, as a result of complexity, duration, or other factors, it is determined that a unique, project-specific quality plan is required, the project QAE shall notify the Project Manager and a project-specific quality plan shall be generated in accordance with **SOP-01-5.2**, *Preparation and Revision of Plans*.

2.0 SCOPE

This Quality Project Plan shall be applied to the chemical analyses performed for commercial nuclear power plants by the Department of Analytical and Environmental Chemistry within the Chemistry and Chemical Engineering Division. Although the majority of the work performed for nuclear facilities resides within the Department of Analytical and Environmental Chemistry, other departments within the division may utilize this Quality Project Plan as deemed necessary when nuclear projects are conducted.

3.0 REFERENCES

- 3.1 *SwRI Quality System Manual – 2000*
- 3.2 *10 CFR 50, Appendix B, ASME NQA-1*
- 3.3 *SwRI Program Quality Plan (PQP-Nuclear), Nuclear Services*
- 3.4 *01-QAP-004, Quality Assurance Plan for Analytical and Environmental Services*

4.0 APPLICABLE SECTIONS OF SwRI PROGRAM QUALITY PLAN (PQP-NUCLEAR)

4.1 Indoctrination and Training

- 4.1.1 Personnel performing duties affecting quality shall receive quality training to the *SwRI Program Quality Plan (PQP-Nuclear), Nuclear Services* prior to performing any work on projects for nuclear utilities. Institute Quality Systems (IQS) personnel shall perform this training and documentation shall be evident in the personnel training files maintained in Division Quality Assurance.
- 4.1.2 Indoctrination and training of personnel shall be conducted in accordance with **SOP-01-18.1**, *Qualification and Training*.

4.2 Qualification of Personnel

- 4.2.1 Testing personnel shall be designated as qualified to perform applicable project activities as specified in **SOP-01-18.1, *Qualification and Training***.
- 4.2.2 During the performance of each testing process, testing personnel shall have access to the necessary documented procedures, i.e., QPP, QAP, Work Order, Division Quality System Standard Operating Procedures (SOPs), and applicable test/analytical procedures (TAPs) available for ready reference.
- 4.2.3 Any person who has not performed testing activities associated with any particular method being used for nuclear utilities projects for a period of one year shall be reevaluated prior to the conduct of the test.
- 4.2.4 Quality Assurance personnel witnessing the testing process for nuclear utilities shall have documented evidence of qualifications maintained by Institute Quality Systems.

4.3 Design Control

Not applicable to activities conducted within the Department of Analytical and Environmental Chemistry.

4.4 Right of Access

- 4.4.1 Procurement documents shall provide for access to the suppliers' facilities and records for surveillance, inspection, or audit by SwRI and clients.
- 4.4.2 Where appropriate, quality clause **Q32** shall be noted on the procurement documents to indicate that right of access for inspection and surveillance of activities associated with the order shall be afforded to SwRI and clients.

4.5 Control of Supplier-Generated Documents

- 4.5.1 Client documents shall be controlled in accordance with procedures under **SOP-01-5.0, *Document and Data Control***, depending on the type of document supplied. These procedures provide the requirements for the preparation, review, approval, issue, distribution, and revision of documents controlled by the Chemistry and Chemical Engineering Division.
- 4.5.2 Documents may be controlled as Plans or Work Instructions and shall be accessible through the Division Intranet link, ***Contract Requirements*** as PDF files.
- 4.5.3 Nationally recognized test methods shall be of the most current issue or as specified in the purchase order. Work orders shall identify the applicable test methods to be used on the nuclear project.

4.6 Acceptance of Services Only

Not applicable to activities conducted within the Department of Analytical and Environmental Chemistry.

4.7 Commercial Grade Items

- 4.7.1 Where an item is to be incorporated into a test or deliverable to a client, and that item is not subject to design or specification requirements that are unique to nuclear facilities, used in applications other than nuclear facilities, and procured from the supplier on the specifications set forth in the manufacturers' published product and description, the item shall be considered "commercial grade".
- 4.7.2 Chemical reagents and standards used for testing purposes shall be ordered to specific chemical grades and certificates of analysis shall be required with each lot.
- 4.7.3 Controls for procurement planning, supplier selection, supplier performance evaluation, and acceptance of procured items and services other than chemical reagents and standards shall be as identified in **SOP-01-6.1, Purchasing**, and any referenced document within that procedure.
- 4.7.4 Receipt inspection of chemical reagents, standards, and test items for use on nuclear safety-related projects shall be performed by department personnel and documented on the *SwRI Receipt Traveler* or **FRM-109, Item Receipt Report**, as specified in **SOP-01-10.1, Inspection and Test Conduct**. Any discrepancy such as a damaged container or container label shall be documented on the form and the client shall be contacted for disposition.
- 4.7.5 Prior to acceptance of a commercial grade item, the receipt inspection shall determine the following:
- (a) Damage was not sustained during shipment;
 - (b) The item has satisfied the specified acceptance criteria; and
 - (c) Specified documentation, as applicable to the item, was received and is acceptable.
- 4.7.6 Receipt inspection of chemical reagents and standards shall also consist of verification of chemical type, grade, container integrity, certificate of analysis, and shelf life, where applicable. Upon acceptance of chemical reagents and standards, the containers shall be labeled with the following:
- (a) Chemical name;
 - (b) Chemical grade;
 - (c) Lot code;
 - (d) Date received; and
 - (e) Shelf life, when applicable.

010032

- 4.7.7 Expired shelf life items shall not be used for testing purposes.
- 4.7.8 Lot codes of chemical reagents and standards used during equipment standardization and testing shall be recorded on the individual testing data sheets to provide traceability.
- 4.7.9 Samples supplied to SwRI for testing shall be received by the Sample Custodian and logged into the laboratory logbook. Sample documentation and sample custody shall be maintained in accordance with **TAP-01-0407-001**, *Sample Receipt and Login*, and **TAP-01-0407-035**, *Sample/Extract Storage and Custody*.
- 4.7.10 Samples supplied to SwRI for testing shall be labeled with the following:
 - (a) Sample control number;
 - (b) Purchase order number;
 - (c) Purchase order line item number, as applicable;
 - (d) Work order number;
 - (e) Nuclear QA label; and
 - (f) Sample retention date, when applicable.
- 4.7.11 In the event that samples are damaged upon receipt, a **Sample Discrepancy Record** shall be generated from the Division Intranet.
- 4.7.12 The testing work order shall list the project number, tests required, test methods required, and shall be labeled *Nuclear Quality*.
- 4.7.13 Identification and traceability shall be maintained in accordance with **SOP-01-8.1**, *Item Identification and Traceability*.

4.8 Inspection

- 4.8.1 Inspection for acceptance shall be performed by qualified persons other than those who conduct or directly supervise the work being inspected.
- 4.8.2 Institute Quality System (IQS) personnel shall perform surveillance activities as required to ensure compliance with the contract and this Quality Project Plan. Specific areas in which IQS may perform surveillance activities include, but are not limited to, the following:
 - (a) Receiving inspection and labeling of chemical reagents, standards, and testing samples;
 - (b) Testing processes;
 - (c) Calibration and major equipment;
 - (d) Sample and record retention; and

(e) Test records.

4.9 Inspection and Testing

- 4.9.1 Required tests for acceptance shall be conducted under appropriate environmental conditions using the tools and equipment necessary to conduct the test in a manner to fulfill test requirements and acceptance criteria.
- 4.9.2 Tests shall be conducted, controlled, and verified in accordance with **SOP-01-10.1, *Inspection and Test Conduct***.
- 4.9.3 Controls for measuring and test equipment shall be as specified in **SOP-01-11.1, *Control of Measuring and Test Equipment***.
- 4.9.4 Controls for identification, segregation, reporting, and resolution of nonconforming items and conditions shall be as specified in **SOP-01-13.1, *Nonconformance Reporting***.

4.10 Handling, Storage, Packaging, Preservation, and Delivery

- 4.10.1 Controls for handling, storage, packaging, preservation, and delivery of items are identified in **SOP-01-15.1, *Handling, Storage, Packaging, Preservation, and Delivery of Items***.
- 4.10.2 Samples specified on the purchase order to be returned to the client shall be prepared and packaged as specified on the purchase order. Each package shall be marked legibly and indelibly with the purchase order/release number and line item number(s) relevant to the package.

4.11 Quality Assurance Records

- 4.11.1 Quality assurance records shall furnish documentary evidence that items or activities meet specified quality requirements. Documents that ensure this evidence include **TAP-01-0407-014, *Inventory of Case File Purges***, and **SOP-01-16.1, *Storage and Maintenance of Quality Records***. These documents and this QPP ensure that QA records shall be legible, identifiable, retrievable, and maintained in dual storage.
- 4.11.2 Records shall be traceable to associated items and activities and shall accurately reflect the work accomplished or information required.
- 4.11.3 Documents shall be considered valid records only if stamped, initialed or signed and dated by authorized personnel or otherwise authenticated.
- ~~4.11.4~~ Records of test analyses performed by the Department of Analytical and Environmental Chemistry are classified as *nonpermanent* and shall be retained for a minimum of five years. Nonpermanent records are those required to show evidence that an activity was performed in accordance with the applicable requirements, but need not be retained for the life of the item. Based on the use of the final data, the client shall be responsible for determining and implementing permanent storage requirements.

4.11.5 In order to satisfy duplicate storage requirements, one copy of the QA record shall be maintained by the Project Manager in Building 70 and a separate copy shall be maintained in the Division Quality Assurance Archives in Building 201. Storage requirements shall be as stated in **SOP-01-16.1**, *Storage and Maintenance of Quality Records*, to ensure protection against the risk of damage or destruction.

4.12 10 CFR, Part 21

4.12.1 SwRI procurement documents shall include requirements for reporting and approving disposition of supplier nonconformances and, when required, compliance to 10 CFR, Part 21.

4.12.2 The Manager of Institute Quality Assurance or Director of Institute Quality Systems shall determine if a nonconforming condition is reportable under 10 CFR, Part 21, and initiate reporting and condition in accordance with the SwRI Operating Policies and Procedures (OPP). Safety hazards or defects that could create a substantial safety hazard shall be reported. Substantial safety hazard means a loss of safety function to the extent that there is a major reduction in the degree of protection provided to public health and safety.

4.13 Certified Test Report

The Project Manager and Institute Quality Assurance Manager as complying with all contractual requirements shall certify test reports. The certified test report shall reference the purchase order/release number, the test methods performed, and the purchase order/release line item number.

4.14 Valid Documents List

The Department of Analytical and Environmental Chemistry work order shall specify all applicable documents and appropriate document revision level for each document. The work order shall then serve as the Valid Documents List (VDL) for each individual project.

5.0 HISTORY OF REVISIONS

Revision 4

Title of document changed from the Standard Project Quality Plan *SPQP-CH/AN* to Quality Project Plan, *QPP-015*

Extensive revision to comply with Project Quality Plan PQP-Nuclear, *Nuclear Services*, which replaces SwRI NQAPM, *Nuclear Quality Assurance Program Manual*.

**SOUTHWEST RESEARCH INSTITUTE
NUCLEAR PROJECT
CLIENT: Div. 20
TASK ORDER: 031106-12
SRR: 25200
SDG: 237578
CASE: CNWRA
VTSR: November 06, 2003
PROJECT#: 06002.01.141**

Chain of Custody/Login Paperwork

Shipper Name/Address		SAMPLE LIST/CHAIN OF CUSTODY													Requested Turnaround:		
Client		Southwest Research Institute Chemistry and Chemical Engineering Division 6220 Culebra Road San Antonio, Texas 78238-5166													☐ 2 Weeks ☐ 3 Weeks ☒ Other: <u>4 wks</u>		
		Client Purchase Order/Other ID						Site/Zone ID							SwRI Contact		
		Analyses Requested													Mike Dammann		
Sample ID	Sample Collection Date (mm/dd/yy)	Sample Collection Time	Matrix Type	Sample Type	# of Containers	Analysis of major & minor elements by ICP											REMARKS
19IM1-25-1003-UFUA-5	11-5-03		W		1	X									intact	✓	Nuclear Safety Related-use
19IM2-25-1003-UFUA-5	11-5-03		W		1	X									✓		appropriate QA procedures
7SC-24-1003-UFUA-5	11-5-03		W		1	X									✓		FOC - Brad Werling
16P-1003-UFUA-5	11-5-03		W		1	X									✓		phone - 6565 fax - 5184
19D-1003-UFUA-5	11-5-03		W		1	X									✓		
Matrix Types: A - Air B - Biota D - Dust E - Emission/Stack L - Liquid P - Product Sd - Solid S - Soil SED - Sediment T - Tissue W - Water WP - Wipe		Sample Types: D - Duplicate ER - Equipment Rinsate ES - Environmental Sample FB - Field Blank FD - Field Duplicate MS - Matrix Spike MSD - Matrix Spike Dup TB - Trip Blank		Relinquished by (Print/Signature) <u>Brad Werling/Brad Werling</u>		Date 11-6-03		Time 0916		SwRI Project#: 20.06002.01.141							
Temp: <u>22°C</u>		Therm #: <u>026</u>		Received by (Print/Signature)		Date		Time		Received by SwRI Lab. (Signature) <u>[Signature]</u>							
Comments: <u>610/81-85 EWDP 27-Oct</u>				Relinquished by (Print/Signature)		Date		Time		Date 11/6/03		Time 0916					
				Received by (Print/Signature)		Date		Time		Samples Disposed: Date		Time					
				Relinquished by (Print/Signature)		Date		Time		Samples Disposed by:							

Shipper Name/Address BRAD WERLING CWRA - DIV 20 BLD - 57		SAMPLE LIST/CHAIN OF CUSTODY Southwest Research Institute Chemistry and Chemical Engineering Division 6220 Culebra Road San Antonio, Texas 78238-5166										Requested Turnaround: ☐ 2 Weeks ☐ 3 Weeks ☒ Other: <u>4 wks</u>					
Client		Client Purchase Order/Other ID					Site/Zone ID					SwRI Contact Mike Damann					
		Analyses Requested										REMARKS Preservation a = HCl to pH <2 <u>b = HNO₃ to pH <2</u> c = H ₂ SO ₄ to pH <2 d = NaOH to pH >12 e = Cool (4°C±2°C) f = Other (specify)					
Sample ID	Sample Collection Date (mm/dd/yy)	Sample Collection Time	Matrix Type	Sample Type	# of Containers	Analysis of Major + minor elements by ICP											
98Q-1003 - FA-6	11-5-03		W		1	X									intact	✓	Nuclear Safety Related-use appropriate QA procedures
98Q-1003- FA-5	11-5-03		W		1	X									✓		
19IM1-25-1003- FUA-5	11-5-03		W		1	X									✓		
19IM2-25-1003- FUA-5	11-5-03		W		1	X									✓		POL - Brad Werling phone - 6565 fax - 5184
19D-1003- FUA-5	11-5-03		W		1	X									✓		
Matrix Types: A - Air B - Biota D - Dust E - Emission/Stack L - Liquid P - Product Sd - Solid S - Soil SED - Sediment T - Tissue W - Water WP - Wipe		Sample Types: D - Duplicate ER - Equipment Rinsate ES - Environmental Sample FB - Field Blank FD - Field Duplicate MS - Matrix Spike MSD - Matrix Spike Dup TB - Trip Blank		Relinquished by (Print/Signature) Brad Werling / Brad Werling					Date 11/6/03		Time 0916		SwRI Project#: 20.06002.01.141				
		Received by (Print/Signature)					Date		Time		Received by SwRI Lab: (Signature) re p [Signature]						
		Relinquished by (Print/Signature)					Date		Time		Date 11/6/03		Time 0916				
		Received by (Print/Signature)					Date		Time		Samples Disposed: Date Time						
Temp: 22C Therm #: 026		Relinquished by (Print/Signature)					Date		Time		Samples Disposed by:						
Comments: 610/81-85 EWDP 27-30-Oct																	

Shipper Name/Address Brad Werling CNWRA-DIV 20 BLD-57		SAMPLE LIST/CHAIN OF CUSTODY Southwest Research Institute Chemistry and Chemical Engineering Division 6220 Culebra Road San Antonio, Texas 78238-5166										Requested Turnaround: ☐ 2 Weeks ☐ 3 Weeks ☒ Other: <u>4 wks</u>					
Client		Client Purchase Order/Other ID					Site/Zone ID					SwRI Contact Mike Dumann					
		Analyses Requested												REMARKS Preservation a = HCl to pH <2 b = HNO₃ to pH <2 c = H ₂ SO ₄ to pH <2 d = NaOH to pH >12 e = Cool (4°C±2°C) f = Other (specify)			
Sample ID	Sample Collection Date (mm/dd/yyyy)	Sample Collection Time	Matrix Type	Sample Type	# of Containers	Analysis of major + minor elements by ICP											
16P-1003-FUA-5	11-5-03		W		1	X									intact	✓	Nuclear Safety Related-use appropriate QA procedures POC-Brad Werling phone- 6565 Fax - 5184
75C-24-1003-FUA-5	11-5-03		W		1	X								/	/		
19IM1-25-1003-FA-5	11-5-03		W		1	X								/	/		
19IM2-25-1003-FA-5	11-5-03		W		1	X								/	/		
16P-1003-FA-5	11-5-03		W		1	X								/	/		
Matrix Types: A - Air B - Biota D - Dust E - Emission/Stack L - Liquid P - Product Sd - Solid S - Soil SED - Sediment T - Tissue W - Water WP - Wipe		Sample Types: D - Duplicate ER - Equipment Rinsate ES - Environmental Sample FB - Field Blank FD - Field Duplicate MS - Matrix Spike MSD - Matrix Spike Dup TB - Trip Blank		Relinquished by (Print/Signature) Brad Werling					Date 11-6-03		Time 0916		SwRI Project#: <u>20,06002.01,141</u>				
Temp: <u>22°C</u>		Therm #: <u>026</u>		Received by (Print/Signature)					Date		Time		Received by SwRI Lab: (Signature) Ken E				
Comments: 610/81-85 EWDP 27-30-05		Relinquished by (Print/Signature)					Date		Time		Date 11/6/03		Time 0916				
Samples Disposed:		Date		Time		Samples Disposed by:		Date		Time		Samples Disposed by:					

Shipper Name/Address BRAD WERLING CNRA - DIV 20 BLD - 57		SAMPLE LIST/CHAIN OF CUSTODY Southwest Research Institute Chemistry and Chemical Engineering Division 6220 Culebra Road San Antonio, Texas 78238-5166												Requested Turnaround: ☐ 2 Weeks ☐ 3 Weeks ☒ Other: 4 WKS			
		Client Purchase Order/Other ID						Site/Zone ID						SwRI Contact			
Client		Analyses Requested												Mike Dummam			
Sample ID	Sample Collection Date (mm/dd/yyyy)	Sample Collection Time	Matrix Type	Sample Type	# of Containers	Analysis of major + minor elements by ICP										REMARKS Preservation a = HCl to pH <2 <u>b = HNO₃ to pH <2</u> c = H ₂ SO ₄ to pH <2 d = NaOH to pH >12 e = Cool (4°C ± 2°C) f = Other (specify)	
19D-1003-FA-5	11-5-03		W		1	X									intact	Nuclear Safety Related - use appropriate QA procedures	
75C-24-1003-FA-5	11-5-03		W		1	X										appropriate QA procedures	
Matrix Types: A - Air B - Biota D - Dust E - Emission/Stack L - Liquid P - Product Sd - Solid S - Soil SED - Sediment T - Tissue W - Water WP - Wipe Temp:		Sample Types: D - Duplicate ER - Equipment Rinsate ES - Environmental Sample FB - Field Blank FD - Field Duplicate MS - Matrix Spike MSD - Matrix Spike Dup TB - Trip Blank Therm #:		Relinquished by (Print/Signature) <i>Brad Werling / Brad Werling</i>						Date		Time		SwRI Project#:			
				Received by (Print/Signature)						Date		Time		Received by SwRI Lab: (Signature)			
				Relinquished by (Print/Signature)						Date		Time		<i>Kenel</i>			
				Received by (Print/Signature)						Date		Time		Date: 11/6/03 Time: 0916			
				Relinquished by (Print/Signature)						Date		Time		Samples Disposed: Date Time			
Comments: 610/81-85 FWDP 27-30-02														Samples Disposed by:			



Southwest Research Institute
LOG-IN CHECKLIST
 Sample Receipt Report No.: 25200

010039

Project Information			
All product information shall be derived from a contract, statement of work (SOW), or project manager (PM).			
Previous PM information	Yes ☐	No ☒	
Current PM information	Yes ☒	No ☐	
Sponsor's Documentation			
All original documentation must be attached to the work order and sent to Div 01 Quality Assurance (DQA).			
Delivery documentation (*) = (Tape on 8 1/2 x 11 white paper)	Yes ☐	No ☒	
Hand carried	Yes ☒	No ☐	
Correspondence	Yes ☐	No ☒	
Custody seals*	Yes ☐	No ☒	
Chain-of-custody/Traffic report	Yes ☒	No ☐	
Sample tags* (Confirm with container information before removal)	Yes ☐	No ☒	
Hazard notifications	Yes ☐	No ☒	
Product Condition			
Sample(s) shall not be left unattended once an ice chest has been opened.			
Sample(s) received intact	Yes ☒	No ☐	
Container(s) received leaking	Yes ☐	No ☒	
Container(s) received damaged	Yes ☐	No ☒	
Lid(s) received damaged	Yes ☐	No ☒	
Product Preparation			
FRM-115 (pH) Preservation Report / Metals & Wetchem H2O	Yes ☒	No ☐	
Product Storage			
Sample(s) shall be stored in a designated location with their corresponding sample control records (SCR).			
Sample(s) in designated location	Yes ☒	No ☐	
Sample Receipt Report			
All copies must be legible.			
Computerized sample receipt report with original documents	Yes ☒	No ☐	
Sample custodian signature and date on all original documents.	Yes ☒	No ☐	
Copy of sample receipt report to PM	Yes ☒	No ☐	
PM sample receipt report approval with initial and date	Yes ☒	No ☐	
Case file documents to DQA	Yes ☒	No ☐	
Print sample log-in page and place in log-in notebook	Yes ☒	No ☐	
Chain-of-Custody/Traffic Report			
Sample Custodian's signature, date, and time	Yes ☒	No ☐	
Document thermometer no. / temperature (Confirm calibration is current)	Yes ☒	No ☐	
Document condition of sample(s)	Yes ☒	No ☐	
FRM-109			
Discrepancies	Yes ☒	No ☐	
FRM-109 to PM	Yes ☒	No ☐	
FRM-109 to Client (Only for required contracts)	Yes ☐	No ☒	
FRM-109 to Inorganics for discrepancies of FRM-115	Yes ☒	No ☐	
Ice Chest Shipment			
Clean ice chest/Discard used bags and ice/Return Sponsor's packing material	Yes ☐	No ☒	
SwRI shipping ticket with project identification	Yes ☐	No ☒	
Log-In Closure			
Sample Custodian's Signature: <i>[Signature]</i>	Date: 11/06/03	Time: 11:57 AM	

SAMPLE LOG-IN SHEET

000040

Lab Name

Southwest Research Institute

Page 1 of 1

Received By (Print Name)

Log-in Date

11/06/2003

KHALED EDRISI

Received By (Signature)

Case Number

CNWRA

Sample Delivery Group No.

SAS Number

N/A

Remarks: 06002.01.141

Remarks:
Condition of Sample
Shipment, etc

		Corresponding			
		EPA Sample #	Sample Tag #	Assigned Lab #	
1.	Custody Seal(s) Present/ ~~Absent~~ Intact/ ~~Broken~~	16P-1003-FA-5	None	237578	Intact
2.	Custody Seal Nos. <u>N/A</u>	16P-1003-FUA-5	None	237579	Intact
		16P-1003-UFUA-5	None	237580	Intact
3.	Chain-of Custody Records ~~Present~~ /Absent*	19D-1003-FA-5	None	237581	Intact
4.	Traffic Reports or Packing Lists Present/ ~~Absent~~	19D-1003-FUA-5	None	237582	Intact
5.	Airbill Airbill/Sticker ~~Present~~ /Absent*	19D-1003-UFUA-5	None	237583	Intact
6.	Airbill No. HAND DELIVERED	19IM1-25-1003-FA-5	None	237584	Intact
		19IM1-25-1003-FUA-5	None	237585	Intact
7.	Sample Tags Present/ ~~Absent~~	19IM1-25-1003-UFUA-5	None	237586	Intact
	Sample Tag Numbers Listed/ ~~Not~~ listed on Chain of Custody	19IM2-25-1003-FA-5	None	237587	Intact
		19IM2-25-1003-FUA-5	None	237588	Intact
8.	Sample Condition ~~Intact~~ /Broken*/Leaking	19IM2-25-1003-UFUA-5	None	237589	Intact
9.	Cooler Temperature 22.0C	7SC-24-1003-FA-5	None	237590	Intact
10.	Does Information on custody records, traffic reports, and sample tags agree? ~~Yes~~ /No*	7SC-24-1003-FUA-5	None	237591	Intact
		7SC-24-1003-UFUA-5	None	237592	Intact
11.	Date Received at Lab 11/06/2003	98Q-1003-FA-5	None	237593	Intact
12.	Time Received 09:16:00	98Q-1003-FA-6	None	237594	Intact

Sample Transfer

Fraction

Inorganics

Area #

RF 13

Fraction

Area #

By

KHALED EDRISI

By

On

11/06/2003

On

* Contact SMO and attach record of resolution

Reviewed By

CYNTHIA A. SAUCEON

Logbook No.

Sample Receipt (25200)

Date

11/07/2003

Logbook Page No.

4882 (SECTION 2 OF 4)

**SOUTHWEST RESEARCH INSTITUTE
NUCLEAR PROJECT
CLIENT: Div. 20
TASK ORDER: 031106-12
SRR: 25200
SDG: 237578
CASE: CNWRA
VTSR: November 06, 2003
PROJECT#: 06002.01.141**

Copies of Login Book

Sample Login Book

Nov 06, 2003

060041

SwRI Login Area
Division 1

Sample Receipt: 25199		Project: 06002.01.141	Client: DIV 20
VTSR Date: Nov 06, 2003		VTSR Time: 09:16:00	Manager: DAMMANN, MIKE
System ID	Customer Sample ID	Matrix	
237571	19IM1-25-1003-UFUA-1	Water	
237572	19IM2-25-1003-FUA-1	Water	
237573	19IM2-25-1003-UFUA-1	Water	
237574	7SC-24-1003-FUA-1	Water	
237575	7SC-24-1003-UFUA-1	Water	
237576	7SC-24-1003-FUA-2	Water	
237577	99Q-1003-FUA-1	Water	

Sample Receipt: 25200		Project: 06002.01.141	Client: DIV 20
VTSR Date: Nov 06, 2003		VTSR Time: 09:16:00	Manager: DAMMANN, MIKE
System ID	Customer Sample ID	Matrix	
237578	16P-1003-FA-5	Water	
237579	16P-1003-FUA-5	Water	
237580	16P-1003-UFUA-5	Water	
237581	19D-1003-FA-5	Water	
237582	19D-1003-FUA-5	Water	
237583	19D-1003-UFUA-5	Water	
237584	19IM1-25-1003-FA-5	Water	
237585	19IM1-25-1003-FUA-5	Water	
237586	19IM1-25-1003-UFUA-5	Water	
237587	19IM2-25-1003-FA-5	Water	
237588	19IM2-25-1003-FUA-5	Water	
237589	19IM2-25-1003-UFUA-5	Water	
237590	7SC-24-1003-FA-5	Water	
237591	7SC-24-1003-FUA-5	Water	
237592	7SC-24-1003-UFUA-5	Water	
237593	98Q-1003-FA-5	Water	
237594	98Q-1003-FA-6	Water	

Sample Receipt: 25201		Project: 05748.01.002	Client: EGLIN A.F.B.
VTSR Date: Nov 06, 2003		VTSR Time: 09:35:00	Manager: SCOTT, JIM
System ID	Customer Sample ID	Matrix	
237596	1008 A	Water	

**SOUTHWEST RESEARCH INSTITUTE
NUCLEAR PROJECT**

CLIENT: Div. 20

TASK ORDER: 031106-12

SRR: 25200

SDG: 237578

CASE: CNWRA

VTSR: November 06, 2003

PROJECT#: 06002.01.141

RAW DATA

TRACE METALS PREPARATORY LABORATORY DIGESTION LOG

SOUTHWEST RESEARCH INSTITUTE
SAN ANTONIO, TEXAS 78228BOOK / PAGE: **54 294**

CLIENT(S): SW 11-11-03 TRC DIV 20 **010042**

TASK ORDER(S): 031106-12 SDG(S): 237578

PROJECT NO(S): 06002.01.141

METHOD: 3005A 3050B 3050B-7.5 3010A 3020A 7760A 7740A HClO₄ H₂SO₄ Sb

Microwave Fusion Teflon Rock OTHER AL Quots

MATRIX: Water ✓ Soil Biota Solid Liquid TCLP Ext OTHER

INSTRUMENT: GFAA ICP ICP-MS IC FLAA HYDRIDE OTHER

ACID INORG #: HNO₃# XX HCl# XX H₂SO₄# HF HClO₄

INTERNAL STD: Sc @ 10 PPM ✓ Be @ 10 PPM SOURCE: 10 INORG# 4261 EXP: 10/04 AMT: 50NL

Oven/Hotplate ID: Temperature (°C):

SAMPLE IDENTIFICATION	pH	WT (g)	I.V. (mL)	F.V. (mL)	
PBW - L11W1			5	5	
LCSW - L11W1 *					
237578					
75781					
7579					
75796 *					
7580					
7581					
7582					
7583					
7584					
7585					
7586					
7587					
7588					
7589					
7590					
7591					
7592					
7593					
7594					
* Spike 50NL Sp: KCl	Spent	4306	Exp. 10/04		
↓ 20NL B-CAL	Spent	4245	Exp. 9/04		
XX - RLB-02-032-04	190	HNO ₃ / 590	HCL	Exp. 1/04	
PBW and LCSW are prepared as 5mLs - XX					
					LOCATION: N/A

PREPARED BY: John WilksDATE: 11-11-03REVIEWED BY: DDATE: 11-11-03

DISPOSAL INT/DATE/LOC: _____

TRACE METALS PREPARATORY LABORATORY DIGESTION LOG

SOUTHWEST RESEARCH INSTITUTE
SAN ANTONIO, TEXAS 78228

BOOK / PAGE: **55 124**

CLIENT(S): DIV 20 **010043**
 TASK ORDER(S): 031106-12 SDG(S): 237578
 PROJECT NO(S): 06002.01.141
 METHOD: 3005A ☐ 3050B ☐ 3050B-7.5 ☐ 3010A ☐ 3020A ☐ 7760A ☐ 7740A ☐ HClO₄ ☐ H₂SO₄ ☐ Sb ☐
 Microwave ☐ Fusion ☐ Teflon ☐ Rock ☐ OTHER ALiquots
 MATRIX: Water ☒ Soil ☐ Biota ☐ Solid ☐ Liquid ☐ TCLP Ext ☐ OTHER ☐
 INSTRUMENT: GFAA ☐ ICP ☒ ICP-MS ☐ IC ☐ FLAA ☐ HYDRIDE ☐ OTHER ☐
 ACID INORG #: HNO₃# ** HCl# ** H₂SO₄# HClO₄# HF#
 INTERNAL STD: Sc @ 10 PPM ☒ Be @ 10 PPM ☐ SOURCE: 10 INORG# 4261 EXP: 10/04 AMT: 50NL
 Oven/Hotplate ID: Temperature (°C):

SAMPLE IDENTIFICATION	pH	WT (g)	I.V. (mL)	F.V. (mL)	
PBW-L24W2			5	50	11-24-03
LCSW-L24W1 *					
237578					
7579					
7580					
7580d					
7581					
7581⑤					
~~7582~~ 7582 8 ⁰⁰ 11-24-03					
7583					
7584					
7585					
7586					
7587					
7588					
7589					
7590					
7591					
7592					
7593					
✓ 7594			↓	↓	
* Spike 50NL Spike1	Spect#	4306	Exp	10/04	
↓ 20NL I-CAL 1	Spect#	4245	Exp	9/04	
** - RLB-02-032-04	1% HNO ₃ / 5% HCL				
PBW and LCSW are prepared as 5 mLs-*					
					LOCATION: N/A

PREPARED BY: John Wills

REVIEWED BY:

DISPOSAL INT/DATE/LOC:

DATE: 11-24-03

DATE: 11-24-03

Trace Metals Reagent Logbook

SOUTHWEST RESEARCH INSTITUTE
SAN ANTONIO, TEXAS 78228

BOOK / PAGE: 02 032

010044

Reagent I.D.:	Preparation Description:	Prep Date:	Exp Date:	Initials:
RLB-02-032-01	10% HCL	10-28-03	1-28-04	JW
↓	add 180 mLs di water to a 200 mL Nalgene bottle, add 20 mLs HCL invy # 4291	↓	↓	↓
↓	F.V. is 200 mLs	↓	↓	↓
RLB02-032-02	10 x sol. As			
	9.50 L Ascorbic #4272 + 18.75 G KI	10-28-03	10-29-03	BL
	# 5758 + 18.75 mL HCL # 4292 to	↓	↓	↓
	100 mL with Reag. water	↓	↓	↓
RLB-02-032-03	1% HNO ₃ / 5% HCL	10-29-03	1-29-04	JW
↓	add 100 mLs di water to a 250 mL Nalgene bottle, add 2.5 mLs HNO ₃ invy # 4295,	↓	↓	↓
↓	add 12.5 mLs HCL invy # 4292,	↓	↓	↓
↓	F.V. to 250 mLs with di water	↓	↓	↓
RLB-02-032-04	1% HNO ₃ / 5% HCL	11-4-03	1-4-04	JW
↓	add 100 mLs di water to 250 mL Nalgene bottle, add 2.5 mLs HNO ₃ invy # 4296,	↓	↓	↓
↓	add 12.5 mLs HCL invy # 4292, F.V. to 250 mLs with di water	↓	↓	↓
RLB2-032-05	5% KMNO ₃			
	150 G to 3 L reag. water	↓	↓	↓
	# 3663	11-4-03	11-4-04	BL

SOUTHWEST RESEARCH INSTITUTE®

6220 Culebra Rd
San Antonio, Texas 78228

010043

TJA TRACE ICP DAILY LOG


ANALYST JLDATE 11-12-03As 189.042 Profile Line As Intensity: 234.469PEAK POSITION -0.26293VERNIER POSITION 243AUTOSAMPLER Yes X No

QC PREP DATE:

CCV/ICV	03K12
CRI	03K04
ICSA	03K11
ICSAB	03K04

STDs PREP DATE:

CLP_STD1_SC	03K05
CLP_STD2_SC	03K05
CLP_STD3_SC	03K05
CLP_STD4_SC	03K05
CLP_STD5_SC	03301
CLP_STD6_SC	03K05
BLK_SC	03329

COMMENTS	FILE	CLIENT	TO#	PROJECT NO.	METHOD
A3110612	DA11Y32	Div 20	031106-12	06007.01.141	DA11Y154294
A310104	DA11Y32	INEEL	031010-4 031013-11	06355.04.006	54256
I	I	I	031024-5	I	54256
					

COMMENTS: _____

MAINTENANCE:

REVIEWED BY: T. OhniDATE: 12/3/03

SOUTHWEST RESEARCH INSTITUTE

6220 Culebra Rd
San Antonio, Texas 78228

SPECTRO ICP DAILY LOG

ANALYST: DP

DATE 11-24-03

POWER: 1700

FLOWS:

Aux 39
Coolant 60
Mass Flow Controller 1054

CURRENT	PROPOSED	
5006/5073	5073 / X 5063	Na
4761/4800	4800 / 4795	Fe
4913/4991	4991 / 4975	Sr

QC PREP DATE:	
CCV/ICV	03X12
CRI	~~011210~~
ICSA	
ICSAB	

CLP STD1 SC	03205
CLP STD2 SC	11-203
CLP STD3 SC	
CLP STD4 SC	
CLP STD5 SC	
BLK SC	03221

[illegible]

COMMENTS: _____

MAINTENANCE:

OTHER: _____

Cleaned Torch: _____ YES
 Changed Pump Tubing: X YES
 Cleaned Optics: _____ YES
 Polished Optics: _____ YES

REVIEWED BY:

DATE: 11/24/03

ICP Calibration Blank/ICB/CCB Solution

ID: BLK-03302

Date Prepared: 9-2-03

Prepared By: DL

Make up as needed in 1000ml volumetric flask.

Added ☒ 10 ml HNO₃ INORG #: 4197

Added ☒ 50 ml HCL INORG #: 4190

Added ☒ 1000ul of 10000ppm Sc (INORG. VENT.) EXP. Date: 10-1-03 INORG #: 3683

020047

ICP Calibration Blank/ICB/CCB Solution

ID: BLK-03325

Date Prepared: 9-25-03

Prepared By: DL

Make up as needed in 1000ml volumetric flask.

Added ☒ 10 ml HNO₃ INORG #: 4237

Added ☒ 50 ml HCL INORG #: 4247

Added ☒ 1000ul of 10000ppm Sc (INORG. VENT.) EXP. Date: 10-1-04 INORG #: 4262

ICP Calibration Blank/ICB/CCB Solution

ID: BLK- ^{04 10-1-03} ~~03331~~ 03301

Date Prepared: ~~9-3-03~~ 10-1-03

Prepared By: DL

Make up as needed in 1000ml volumetric flask.

Added ☒ 10 ml HNO₃ INORG #: 4237

Added ☒ 50 ml HCL INORG #: 4249

Added ☒ 1000ul of 10000ppm Sc (INORG. VENT.) EXP. Date: 10-1-04 INORG #: 4262

ICP Calibration Blank/ICB/CCB Solution

ID: BLK-03329

Date Prepared: 10-29-03

Prepared By: DL

Make up as needed in 1000ml volumetric flask.

Added ☒ 10 ml HNO₃ INORG #: 4295

Added ☒ 50 ml HCL INORG #: 4292

Added ☒ 1000ul of 10000ppm Sc (INORG. VENT.) EXP. Date: 10-1-04 INORG #: 4262

ICP Calibration Blank/ICB/CCB Solution

ID: BLK- 03221

Date Prepared: 11-21-03

Prepared By: DL

Make up as needed in 1000ml volumetric flask.

Added ☒ 10 ml HNO3 INORG #: 4299

Added ☒ 50 ml HCL INORG #: 4289

Added ☒ 1000ul of 10000ppm Sc (INORG. VENT.) EXP. Date: 10-1-04 INORG #: 4262

032048

ICP Calibration Blank/ICB/CCB Solution

ID: BLK-03M04

Date Prepared: 12/4/03

Prepared By: T. OH

Make up as needed in 1000ml volumetric flask.

Added ☒ 10 ml HNO3 INORG #: 4299

Added ☒ 50 ml HCL INORG #: 4335

Added ☒ 1000ul of 10000ppm Sc (INORG. VENT.) EXP. Date: 10-1-04 INORG #: 4262

ICP Calibration Blank/ICB/CCB Solution

ID: BLK-

Date Prepared: _____

Prepared By: _____

Make up as needed in 1000ml volumetric flask.

Added _____ 10 ml HNO3 INORG #: _____

Added _____ 50 ml HCL INORG #: _____

Added _____ 1000ul of 10000ppm Sc (INORG. VENT.) EXP. Date: _____ INORG #: _____

ICP Calibration Blank/ICB/CCB Solution

ID: BLK-

Date Prepared: _____

Prepared By: _____

Make up as needed in 1000ml volumetric flask.

Added _____ 10 ml HNO3 INORG #: _____

Added _____ 50 ml HCL INORG #: _____

Added _____ 1000ul of 10000ppm Sc (INORG. VENT.) EXP. Date: _____ INORG #: _____

ICP ICV/CCV SOLUTION

CCV-03212

010049

Date Prepared: 11-12-03

Prepared By: JD

HNO3 INORG #: 4292

HCl INORG #: 4294

Make up as needed in 1000ml volumetric flask in 1% HNO3 AND 5% HCl.

Element	Std Conc (ppm)	Amt added	Check	Source	Inorg #	Stock Conc (ppm)	Exp Date
Sc	10	1ml	✓	INORGVENT	43064262	10000	10-30-04
B	5	5ml	✓	SPEX	4139	1000	6-30-04
Li	5	5ml	✓	SPEX	4065	1000	5-15-04
Mo	5	5ml	✓	SPEX	4066	1000	5-15-04
P	5	5ml	✓	SPEX	4307	1000	10-30-04
Si	5	5ml	✓	SPEX	4232	1000	8-30-04
Ti	5	5ml	✓	SPEX	4234	1000	8-30-04
Sr	5	5ml	✓	SPEX	4308	1000	10-30-04
Sn	5	5ml	✓	SPEX	4141	1000	6-30-04
Bi	5	5ml	✓	SPEX	4138	1000	6-30-04
La	5	5ml	✓	SPEX	4064	1000	5-15-04
Y	5	5ml	✓	SPEX	4067	1000	5-15-04
Pd	1	1ml	✓	SPEX	3928	1000	1-15-04
S	1	1ml	✓	SPEX	4140	1000	6-30-04
Th	1	1ml	✓	SPEX	4233	1000	8-30-04
U	1	1ml	✓	SPEX	4142	1000	6-30-04
W	1	1ml	✓	SPEX	4212	1000	8-15-04
Zr	5	5ml	✓	SPEX	4213	1000	8-15-04
Na	10	1ml	✓	SPEX	3951	10000	2-15-04
ICV-2A	vary	10ml	✓	SPEX	3773	mix	11-30-03
ICV-2B	vary	1ml	✓	SPEX	3774	mix	11-30-03
ICV-2C	vary	10ml	✓	SPEX	3775	mix	11-30-03

Expiration Date: 11-30-03

ICP Calibration Standards

Date Prepared: 11-5-03

Prepared By: [Signature]

HNO3 INORG #: 4296

HCl INORG #: 4293

040030

Make up as needed in 500 ml volumetric flasks in 1% HNO3 and 5% HCl.

Prepared	Standard Name	Element	Std Conc (ppm)	Added ml	Check	Source	INORG #	Stock Conc (ppm)	Exp Date
11-5-03	STD1- 03405 m	Al	50	2.50	✓	INORVENT	4220	10000	9-1-04
		Ca	50	2.50	✓	INORVENT	3900	10000	2-1-04
		Fe	50	2.50	✓	INORVENT	3963	10000	3-1-04
		K	50	2.50	✓	INORVENT	4320	10000	12-1-04
		Mg	25	1.25	✓	INORVENT	4204	10000	8-1-04
		Na	50	2.50	✓	INORVENT	4205	10000	9-1-04
		Li	10	5.00	✓	INORVENT	4149	1000	7-1-04
		Sc	10	0.500	✓	INORVENT	4262	10000	10-1-04
11-5-03	STD2- 03405	Ba	10	5.00	✓	INORVENT	3971	1000	3-1-04
		Be	5	2.50	✓	INORVENT	4062	1000	6-1-04
		Cr	10	5.00	✓	INORVENT	4314	1000	12-1-04
		Cu	10	5.00	✓	INORVENT	3955	1000	3-1-04
		Ni	10	5.00	✓	INORVENT	3957	1000	3-1-04
		Sc	10	0.500	✓	INORVENT	4262	10000	10-1-04
11-5-03	STD3- 03405	Cd	10	5.00	✓	INORVENT	3972	1000	3-1-04
		Co	10	5.00	✓	INORVENT	3933	1000	3-1-04
		Mn	10	5.00	✓	INORVENT	3921	1000	2-1-04
		V	10	5.00	✓	INORVENT	4321	1000	12-1-04
		Zn	10	5.00	✓	INORVENT	4319	1000	12-1-04
		Sc	10	0.500	✓	INORVENT	4262	10000	10-1-04
11-5-03	STD4- 03405	Ag	2	1.00	✓	INORVENT	4222	1000	9-1-04
		As	10	5.00	✓	INORVENT	3999	1000	2-1-04
		Pb	10	5.00	✓	INORVENT	4313	1000	11-1-04
		Sb	10	5.00	✓	INORVENT	3953	1000	3-1-04
		Se	10	5.00	✓	INORVENT	4152	1000	7-1-04
		TL	10	5.00	✓	INORVENT	3922	1000	2-1-04
		Sc	10	0.500	✓	INORVENT	4262	10000	10-1-04
	STD5-	B	10	5.00		INORVENT		1000	
		Mo	10	5.00		INORVENT		1000	
		P	10	5.00		INORVENT		1000	
		Si	10	5.00		INORVENT		1000	
		Ti	10	5.00		INORVENT		1000	
		Sr	10	5.00		INORVENT		1000	
		Sn	10	5.00		INORVENT		1000	
		Bi	5	2.50		INORVENT		1000	
		Sc	10	0.500		INORVENT		10000	
11-5-03	STD6- 03405	La	10	5.00	✓	INORVENT	4221	1000	9-1-04
		Na	1	0.05	✓	INORVENT	4205	10000	8-1-04
		Pd	10	5.00	✓	INORVENT	3964	1000	3-1-04
		S	10	5.00	✓	INORVENT	4317	1000	12-1-04
		Th	10	5.00	✓	INORVENT	4283	1000	11-1-04
		U	10	5.00	✓	INORVENT	3974	1000	3-1-04
		W	5	2.50	✓	INORVENT	4203	1000	9-1-04
		Y	10	5.00	✓	INORVENT	4063	1000	6-1-04
		Zr	10	5.00	✓	INORVENT	3954	1000	3-1-04
		SC	10	0.500	✓	INORVENT	4262	10000	10-1-04

Expiration Dates:

STD1: 2-1-04

STD4: 2-1-04

STD2: 3-1-04

STD5:

STD3: 2-1-04

STD6: 3-1-04

ICP Calibration Standards

Date Prepared: 9-31-03 ^{at} 10-1-03

Prepared By: [Signature]

HNO3 INORG #: 4237

HCl INORG #: 4249

#10051

Make up as needed in 500 ml volumetric flasks in 1% HNO3 and 5% HCl.

Prepared	Standard Name	Element	Std Conc (ppm)	Added ml	Check	Source	INORG #	Stock Conc (ppm)	Exp Date
<u>9-31-03</u> <u>10-1-03</u>	STD1- <u>03301</u>	Al	50	2.50	✓	INORVENT	4220	10000	9-1-04
		Ca	50	2.50	✓	INORVENT	3900	10000	2-1-04
		Fe	50	2.50	✓	INORVENT	3963	10000	3-1-04
		K	50	2.50	✓	INORVENT	3234	10000	11-1-03
		Mg	25	1.25	✓	INORVENT	4204	10000	9-1-04
		Na	50	2.50	✓	INORVENT	4205	10000	9-1-04
		Li	10	5.00	✓	INORVENT	4149	1000	7-1-04
		Sc	10	0.500	✓	INORVENT	4262	10000	10-1-04
<u>9-31-03</u> <u>10-1-03</u>	STD2- <u>03301</u>	Ba	10	5.00	✓	INORVENT	3931	1000	3-1-04
		Be	5	2.50	✓	INORVENT	4062	1000	6-1-04
		Cr	10	5.00	✓	INORVENT	3732	1000	11-1-03
		Cu	10	5.00	✓	INORVENT	3955	1000	3-1-04
		Ni	10	5.00	✓	INORVENT	3957	1000	3-1-04
		Sc	10	0.500	✓	INORVENT	4262	10000	10-1-04
<u>9-31-03</u> <u>10-1-03</u>	STD3- <u>03301</u>	Cd	10	5.00	✓	INORVENT	3922	1000	3-1-04
		Co	10	5.00	✓	INORVENT	3933	1000	3-1-04
		Mn	10	5.00	✓	INORVENT	3921	1000	2-1-04
		V	10	5.00	✓	INORVENT	3740	1000	11-1-03
		Zn	10	5.00	✓	INORVENT	3741	1000	11-1-03
		Sc	10	0.500	✓	INORVENT	4262	10000	10-1-04
<u>9-31-03</u> <u>10-1-03</u>	STD4- <u>03301</u>	Ag	2	1.00	✓	INORVENT	4222	1000	9-1-04
		As	10	5.00	✓	INORVENT	3499	1000	2-1-04
		Pb	10	5.00	✓	INORVENT	3735	1000	11-1-03
		Sb	10	5.00	✓	INORVENT	3953	1000	3-1-04
		Se	10	5.00	✓	INORVENT	4152	1000	7-1-04
		TL	10	5.00	✓	INORVENT	3922	1000	2-1-04
		Sc	10	0.500	✓	INORVENT	4262	10000	10-1-04
<u>9-31-03</u> <u>10-1-03</u>	STD5- <u>03301</u>	B	10	5.00	✓	INORVENT	3954	1000	3-1-04
		Mo	10	5.00	✓	INORVENT	3956	1000	3-1-04
		P	10	5.00	✓	INORVENT	4049	1000	6-1-04
		Si	10	5.00	✓	INORVENT	4153	1000	7-1-04
		Ti	10	5.00	✓	INORVENT	3262	1000	12-1-03
		Sr	10	5.00	✓	INORVENT	4154	1000	7-1-04
		Sn	10	5.00	✓	INORVENT	4026	1000	4-1-04
		Bi	5	2.50	✓	INORVENT	4200	1000	9-1-04
		Sc	10	0.500	✓	INORVENT	4262	10000	10-1-04
<u>10-1-03</u> <u>03301</u>	STD6- <u>03301</u>	La	10	5.00	✓	INORVENT	4221	1000	9-1-04
		Na	1	0.05	✓	INORVENT	4205	10000	9-1-04
		Pd	10	5.00	✓	INORVENT	3964	1000	3-1-04
		S	10	5.00	✓	INORVENT	3739	1000	11-1-03
		Th	10	5.00	✓	INORVENT	3681	1000	10-1-03
		U	10	5.00	✓	INORVENT	3924	1000	3-1-04
		W	5	2.50	✓	INORVENT	4203	1000	8-1-04
		Y	10	5.00	✓	INORVENT	4063	1000	6-1-04
		Zr	10	5.00	✓	INORVENT	3958	1000	3-1-04
		SC	10	0.500	✓	INORVENT	4262	10000	10-1-04

Expiration Dates:

STD1: 11-1-03

STD4: 11-1-03

STD2: 11-1-03

STD5: 12-1-03

STD3: 11-1-03

STD6: 11-1-03 10-1-03

Exp. order

Div 20
to#031106-12
06002.01.141

Handwritten: 12/18/03

Handwritten: 12/22/03



010052

system id	elem	result	qual	units	rl	tv	rec	ug/ml	date	time
237578	Ag3280	0.005	U	mg/L	0.005			0.00051	11/12/03	11:51
237578	Al3082	0.620		mg/L	0.05			0.61976	11/12/03	11:51
237578	As1890	0.014		mg/L	0.005			0.01445	11/12/03	11:51
237578	B 2496	0.204		mg/L	0.025			0.20366	11/12/03	11:51
237578	Ba4934	0.022		mg/L	0.005			0.02162	11/12/03	11:51
237578	Be3130	0.005	U	mg/L	0.005			0.00037	11/12/03	11:51
237578	Bi2230	0.015	U	mg/L	0.015			0.00519	11/12/03	11:51
237578	Ca3179	2.17		mg/L	0.05			2.16732	11/12/03	11:51
237578	Cd2265	0.005	U	mg/L	0.005			0.00024	11/12/03	11:51
237578	Co2286	0.005	U	mg/L	0.005			0.00499	11/12/03	11:51
237578	Cr2677	0.005	U	mg/L	0.005			0.00416	11/12/03	11:51
237578	Cu3247	0.005	U	mg/L	0.005			0.00272	11/12/03	11:51
237578	Fe2714	1.14		mg/L	0.05			1.14329	11/12/03	11:51
237578	K 7664	1.92		mg/L	0.2			1.9201	11/24/03	
237578	La4086	0.005	U	mg/L	0.005			0.00486	11/12/03	11:51
237578	Li6707	0.081		mg/L	0.005			0.08147	11/12/03	11:51
237578	Mg2790	0.388		mg/L	0.05			0.38805	11/12/03	11:51
237578	Mn2576	0.024		mg/L	0.005			0.02448	11/12/03	11:51
237578	Mo2020	0.063		mg/L	0.01			0.06260	11/12/03	11:51
237578	Na3302	110		mg/L	0.2			110.1565	11/24/03	
237578	Ni2316	0.005	U	mg/L	0.005			0.00148	11/12/03	11:51
237578	P 1782	0.043		mg/L	0.04			0.04318	11/12/03	11:51
237578	Pb220	0.005	U	mg/L	0.005			0.00177	11/12/03	11:51
237578	Pd3404	0.01	U	mg/L	0.01			0.00411	11/12/03	11:51
237578	S 1820	20.2		mg/L	0.025			20.16623	11/12/03	11:51
237578	Sb2068	0.01	U	mg/L	0.01			0.00324	11/12/03	11:51
237578	Se196	0.015	U	mg/L	0.015			0.00737	11/12/03	11:51
237578	Si2881	21.1		mg/L	0.05			21.06237	11/12/03	11:51
237578	Sn1899	0.01	U	mg/L	0.01			0.00231	11/12/03	11:51
237578	Sr4215	0.028		mg/L	0.005			0.02773	11/12/03	11:51
237578	Th2837	0.01	U	mg/L	0.01			0.00066	11/12/03	11:51
237578	Ti3372	0.005	U	mg/L	0.005			0.00405	11/12/03	11:51
237578	Ti1908	0.025	U	mg/L	0.025			0.00128	11/12/03	11:51
237578	U 3859	0.1	U	mg/L	0.1			0.07833	11/12/03	11:51
237578	V 2924	0.005	U	mg/L	0.005			0.00275	11/12/03	11:51
237578	W 2079	0.015	U	mg/L	0.015			-0.00026	11/12/03	11:51
237578	Y 3710	0.005	U	mg/L	0.005			0.00069	11/12/03	11:51
237578	Zn2062	0.015		mg/L	0.005			0.01472	11/12/03	11:51
237578	Zr3496	0.005	U	mg/L	0.005			0.00073	11/12/03	11:51
237578d	Ag3280	0.005	U	mg/L	0.005			-0.00017	11/12/03	11:56
237578d	Al3082	0.659		mg/L	0.05			0.65934	11/12/03	11:56
237578d	As1890	0.008		mg/L	0.005			0.00823	11/12/03	11:56
237578d	B 2496	0.206		mg/L	0.025			0.20561	11/12/03	11:56
237578d	Ba4934	0.020		mg/L	0.005			0.01951	11/12/03	11:56
237578d	Be3130	0.005	U	mg/L	0.005			0.00021	11/12/03	11:56
237578d	Bi2230	0.015	U	mg/L	0.015			0.00382	11/12/03	11:56
237578d	Ca3179	2.18		mg/L	0.05			2.18110	11/12/03	11:56
237578d	Cd2265	0.005	U	mg/L	0.005			0.00000	11/12/03	11:56
237578d	Co2286	0.005	U	mg/L	0.005			0.00450	11/12/03	11:56
237578d	Cr2677	0.005	U	mg/L	0.005			0.00455	11/12/03	11:56
237578d	Cu3247	0.005	U	mg/L	0.005			0.00274	11/12/03	11:56
237578d	Fe2714	1.23		mg/L	0.05			1.22652	11/12/03	11:56
237578d	La4086	0.005		mg/L	0.005			0.00506	11/12/03	11:56
237578d	Li6707	0.081		mg/L	0.005			0.08112	11/12/03	11:56

Div 20
to#031106-12
06002.01.141

010053

system id	elem	result	qual	units	rl	tv	rec	ug/ml	date	time
237578d	Mg2790	0.373		mg/L	0.05			0.37323	11/12/03	11:56
237578d	Mn2576	0.025		mg/L	0.005			0.02451	11/12/03	11:56
237578d	Mo2020	0.064		mg/L	0.01			0.06417	11/12/03	11:56
237578d	Ni2316	0.005	U	mg/L	0.005			0.00246	11/12/03	11:56
237578d	P 1782	0.040		mg/L	0.04			0.04030	11/12/03	11:56
237578d	Pb220	0.005	U	mg/L	0.005			0.00336	11/12/03	11:56
237578d	Pd3404	0.01	U	mg/L	0.01			0.00365	11/12/03	11:56
237578d	S 1820	19.6		mg/L	0.025			19.64684	11/12/03	11:56
237578d	Sb2068	0.01	U	mg/L	0.01			0.00094	11/12/03	11:56
237578d	Se196	0.015	U	mg/L	0.015			0.00428	11/12/03	11:56
237578d	Si2881	21.4		mg/L	0.05			21.39072	11/12/03	11:56
237578d	Sn1899	0.01	U	mg/L	0.01			0.00362	11/12/03	11:56
237578d	Sr4215	0.028		mg/L	0.005			0.02784	11/12/03	11:56
237578d	Th2837	0.01	U	mg/L	0.01			-0.00157	11/12/03	11:56
237578d	Ti3372	0.005	U	mg/L	0.005			0.00438	11/12/03	11:56
237578d	Tl1908	0.025	U	mg/L	0.025			0.00641	11/12/03	11:56
237578d	U 3859	0.1	U	mg/L	0.1			0.08619	11/12/03	11:56
237578d	V 2924	0.005	U	mg/L	0.005			0.00285	11/12/03	11:56
237578d	W 2079	0.015	U	mg/L	0.015			0.00089	11/12/03	11:56
237578d	Y 3710	0.005	U	mg/L	0.005			0.00096	11/12/03	11:56
237578d	Zn2062	0.015		mg/L	0.005			0.01467	11/12/03	11:56
237578d	Zr3496	0.005	U	mg/L	0.005			0.00161	11/12/03	11:56
237579	Ag3280	0.005	U	mg/L	0.005			0.00058	11/12/03	12:01
237579	Al3082	0.340		mg/L	0.05			0.34002	11/12/03	12:01
237579	As1890	0.011		mg/L	0.005			0.01112	11/12/03	12:01
237579	B 2496	0.210		mg/L	0.025			0.20976	11/12/03	12:01
237579	Ba4934	0.008		mg/L	0.005			0.00804	11/12/03	12:01
237579	Be3130	0.005	U	mg/L	0.005			0.00011	11/12/03	12:01
237579	Bi2230	0.015	U	mg/L	0.015			0.00672	11/12/03	12:01
237579	Ca3179	1.85		mg/L	0.05			1.84598	11/12/03	12:01
237579	Cd2265	0.005	U	mg/L	0.005			0.00014	11/12/03	12:01
237579	Co2286	0.005	U	mg/L	0.005			0.00390	11/12/03	12:01
237579	Cr2677	0.005	U	mg/L	0.005			-0.00067	11/12/03	12:01
237579	Cu3247	0.005	U	mg/L	0.005			0.00132	11/12/03	12:01
237579	Fe2714	0.262		mg/L	0.05			0.26168	11/12/03	12:01
237579	K 7664	1.79		mg/L	0.2			1.7867	11/24/03	
237579	La4086	0.005	U	mg/L	0.005			0.00225	11/12/03	12:01
237579	Li6707	0.082		mg/L	0.005			0.08247	11/12/03	12:01
237579	Mg2790	0.269		mg/L	0.05			0.26938	11/12/03	12:01
237579	Mn2576	0.014		mg/L	0.005			0.01352	11/12/03	12:01
237579	Mo2020	0.066		mg/L	0.01			0.06640	11/12/03	12:01
237579	Na3302	105		mg/L	0.2			105.1779	11/24/03	
237579	Ni2316	0.005	U	mg/L	0.005			0.00021	11/12/03	12:01
237579	P 1782	0.04	U	mg/L	0.04			0.03875	11/12/03	12:01
237579	Pb220	0.005	U	mg/L	0.005			-0.00006	11/12/03	12:01
237579	Pd3404	0.01	U	mg/L	0.01			0.00430	11/12/03	12:01
237579	S 1820	19.8		mg/L	0.025			19.77896	11/12/03	12:01
237579	Sb2068	0.01	U	mg/L	0.01			0.00279	11/12/03	12:01
237579	Se196	0.015	U	mg/L	0.015			0.00559	11/12/03	12:01
237579	Si2881	20.7		mg/L	0.05			20.72975	11/12/03	12:01
237579	Sn1899	0.01	U	mg/L	0.01			0.00986	11/12/03	12:01
237579	Sr4215	0.022		mg/L	0.005			0.02184	11/12/03	12:01
237579	Th2837	0.01	U	mg/L	0.01			-0.00301	11/12/03	12:01
237579	Ti3372	0.005	U	mg/L	0.005			0.00251	11/12/03	12:01

system id	elem	result	qual	units	rl	tv	rec	ug/ml	date	time
237579	Tl1908	0.025	U	mg/L	0.025			0.00126	11/12/03	12:01
237579	U 3859	0.1	U	mg/L	0.1			0.08519	11/12/03	12:01
237579	V 2924	0.005	U	mg/L	0.005			0.00186	11/12/03	12:01
237579	W 2079	0.015	U	mg/L	0.015			0.00260	11/12/03	12:01
237579	Y 3710	0.005	U	mg/L	0.005			0.00009	11/12/03	12:01
237579	Zn2062	0.007		mg/L	0.005			0.00739	11/12/03	12:01
237579	Zr3496	0.005	U	mg/L	0.005			0.00059	11/12/03	12:01
237579s	Ag3280	0.042		mg/L	0.005	0.05	84.5%	0.04223	11/12/03	12:06
237579s	Al3082	2.30		mg/L	0.05	2	98.2%	2.30374	11/12/03	12:06
237579s	As1890	2.12		mg/L	0.005	2	105.4%	2.11930	11/12/03	12:06
237579s	B 2496	0.197		mg/L	0.025			0.19743	11/12/03	12:06
237579s	Ba4934	2.05		mg/L	0.005	2	102.0%	2.04871	11/12/03	12:06
237579s	Be3130	0.053		mg/L	0.005	0.05	105.9%	0.05294	11/12/03	12:06
237579s	Bi2230	0.015	U	mg/L	0.015			0.00713	11/12/03	12:06
237579s	Ca3179	22.6		mg/L	0.05	20	103.9%	22.63233	11/12/03	12:06
237579s	Cd2265	0.052		mg/L	0.005	0.05	103.5%	0.05174	11/12/03	12:06
237579s	Co2286	0.520		mg/L	0.005	0.5	103.9%	0.51951	11/12/03	12:06
237579s	Cr2677	0.184		mg/L	0.005	0.2	91.8%	0.18356	11/12/03	12:06
237579s	Cu3247	0.259		mg/L	0.005	0.25	103.6%	0.25908	11/12/03	12:06
237579s	Fe2714	1.23		mg/L	0.05	1	97.2%	1.23375	11/12/03	12:06
237579s	La4086	0.005	U	mg/L	0.005			0.00204	11/12/03	12:06
237579s	Li6707	0.083		mg/L	0.005			0.08274	11/12/03	12:06
237579s	Mg2790	21.2		mg/L	0.05	20	104.6%	21.18984	11/12/03	12:06
237579s	Mn2576	0.523		mg/L	0.005	0.5	101.9%	0.52340	11/12/03	12:06
237579s	Mo2020	0.068		mg/L	0.01			0.06783	11/12/03	12:06
237579s	Ni2316	0.509		mg/L	0.005	0.5	101.9%	0.50942	11/12/03	12:06
237579s	P 1782	0.053		mg/L	0.04			0.05273	11/12/03	12:06
237579s	Pb220	0.523		mg/L	0.005	0.5	104.5%	0.52253	11/12/03	12:06
237579s	Pd3404	0.01	U	mg/L	0.01			0.00482	11/12/03	12:06
237579s	S 1820	18.8		mg/L	0.025			18.80460	11/12/03	12:06
237579s	Sb2068	0.370		mg/L	0.01	0.5	74.0%	0.36982	11/12/03	12:06
237579s	Se196	2.32		mg/L	0.015	2	116.0%	2.31989	11/12/03	12:06
237579s	Si2881	19.7		mg/L	0.05			19.72542	11/12/03	12:06
237579s	Sn1899	0.01	U	mg/L	0.01			0.00353	11/12/03	12:06
237579s	Sr4215	0.021		mg/L	0.005			0.02092	11/12/03	12:06
237579s	Th2837	0.01	U	mg/L	0.01			-0.00588	11/12/03	12:06
237579s	Ti3372	0.005	U	mg/L	0.005			0.00201	11/12/03	12:06
237579s	Tl1908	2.28		mg/L	0.025	2	114.2%	2.28312	11/12/03	12:06
237579s	U 3859	0.1	U	mg/L	0.1			0.05602	11/12/03	12:06
237579s	V 2924	0.514		mg/L	0.005	0.5	102.9%	0.51431	11/12/03	12:06
237579s	W 2079	0.015	U	mg/L	0.015			-0.00183	11/12/03	12:06
237579s	Y 3710	0.005	U	mg/L	0.005			0.00015	11/12/03	12:06
237579s	Zn2062	0.536		mg/L	0.005	0.5	105.9%	0.53649	11/12/03	12:06
237579s	Zr3496	0.005	U	mg/L	0.005			0.00114	11/12/03	12:06
237580	Ag3280	0.005	U	mg/L	0.005			0.00031	11/12/03	12:10
237580	Al3082	0.354		mg/L	0.05			0.35365	11/12/03	12:10
237580	As1890	0.019		mg/L	0.005			0.01909	11/12/03	12:10
237580	B 2496	0.273		mg/L	0.025			0.27345	11/12/03	12:10
237580	Ba4934	0.038		mg/L	0.005			0.03849	11/12/03	12:10
237580	Be3130	0.005	U	mg/L	0.005			0.00001	11/12/03	12:10
237580	Bi2230	0.015	U	mg/L	0.015			0.00637	11/12/03	12:10
237580	Ca3179	6.53		mg/L	0.05			6.52550	11/12/03	12:10
237580	Cd2265	0.005	U	mg/L	0.005			0.00017	11/12/03	12:10
237580	Co2286	0.010		mg/L	0.005			0.00978	11/12/03	12:10

Div 20
to#031106-12
06002.01.141

010055

system id	elem	result	qual	units	rl	tv	rec	ug/ml	date	time
237580	Cr2677	0.005	U	mg/L	0.005			0.00148	11/12/03	12:10
237580	Cu3247	0.009		mg/L	0.005			0.00943	11/12/03	12:10
237580	Fe2714	1.07		mg/L	0.05			1.06622	11/12/03	12:10
237580	K 7664	1.98		mg/L	0.2			1.9776	11/24/03	
237580	La4086	0.005	U	mg/L	0.005			0.00225	11/12/03	12:10
237580	Li6707	0.112		mg/L	0.005			0.11168	11/12/03	12:10
237580	Mg2790	0.901		mg/L	0.05			0.90135	11/12/03	12:10
237580	Mn2576	0.054		mg/L	0.005			0.05414	11/12/03	12:10
237580	Mo2020	0.063		mg/L	0.01			0.06266	11/12/03	12:10
237580	Na3302	117		mg/L	0.2			117.4242	11/24/03	
237580	Ni2316	0.005	U	mg/L	0.005			0.00457	11/12/03	12:10
237580	P 1782	0.061		mg/L	0.04			0.06107	11/12/03	12:10
237580	Pb220	0.005		mg/L	0.005			0.00545	11/12/03	12:10
237580	Pd3404	0.01	U	mg/L	0.01			0.00517	11/12/03	12:10
237580	S 1820	27.7		mg/L	0.025			27.65908	11/12/03	12:10
237580	Sb2068	0.01	U	mg/L	0.01			0.00376	11/12/03	12:10
237580	Se196	0.015	U	mg/L	0.015			0.01308	11/12/03	12:10
237580	Si2881	27.1		mg/L	0.05			27.09569	11/12/03	12:10
237580	Sn1899	0.011		mg/L	0.01			0.01054	11/12/03	12:10
237580	Sr4215	0.078		mg/L	0.005			0.07765	11/12/03	12:10
237580	Th2837	0.01	U	mg/L	0.01			-0.00404	11/12/03	12:10
237580	Ti3372	0.005	U	mg/L	0.005			0.00114	11/12/03	12:10
237580	Ti1908	0.025	U	mg/L	0.025			-0.00774	11/12/03	12:10
237580	U 3859	0.1	U	mg/L	0.1			0.07550	11/12/03	12:10
237580	V 2924	0.005	U	mg/L	0.005			0.00391	11/12/03	12:10
237580	W 2079	0.015	U	mg/L	0.015			0.00130	11/12/03	12:10
237580	Y 3710	0.005	U	mg/L	0.005			0.00013	11/12/03	12:10
237580	Zn2062	0.045		mg/L	0.005			0.04497	11/12/03	12:10
237580	Zr3496	0.005	U	mg/L	0.005			0.00064	11/12/03	12:10
237580d	K 7664	1.89		mg/L	0.2			1.8865	11/24/03	
237580d	Na3302	117		mg/L	0.2			117.0635	11/24/03	
237581	Ag3280	0.005	U	mg/L	0.005			-0.00060	11/12/03	12:15
237581	Al3082	0.05	U	mg/L	0.05			0.03098	11/12/03	12:15
237581	As1890	0.007		mg/L	0.005			0.00684	11/12/03	12:15
237581	B 2496	0.232		mg/L	0.025			0.23215	11/12/03	12:15
237581	Ba4934	0.005	U	mg/L	0.005			0.00350	11/12/03	12:15
237581	Be3130	0.005	U	mg/L	0.005			-0.00002	11/12/03	12:15
237581	Bi2230	0.015	U	mg/L	0.015			0.00396	11/12/03	12:15
237581	Ca3179	2.23		mg/L	0.05			2.22925	11/12/03	12:15
237581	Cd2265	0.005	U	mg/L	0.005			0.00005	11/12/03	12:15
237581	Co2286	0.005	U	mg/L	0.005			0.00348	11/12/03	12:15
237581	Cr2677	0.005	U	mg/L	0.005			-0.00101	11/12/03	12:15
237581	Cu3247	0.005	U	mg/L	0.005			0.00086	11/12/03	12:15
237581	Fe2714	0.275		mg/L	0.05			0.27546	11/12/03	12:15
237581	K 7664	2.99		mg/L	0.2			2.9881	11/24/03	
237581	La4086	0.005	U	mg/L	0.005			0.00140	11/12/03	12:15
237581	Li6707	0.249		mg/L	0.005			0.24912	11/12/03	12:15
237581	Mg2790	0.103		mg/L	0.05			0.10256	11/12/03	12:15
237581	Mn2576	0.040		mg/L	0.005			0.04014	11/12/03	12:15
237581	Mo2020	0.039		mg/L	0.01			0.03947	11/12/03	12:15
237581	Na3302	173		mg/L	0.2			172.5599	11/24/03	
237581	Ni2316	0.005	U	mg/L	0.005			-0.00007	11/12/03	12:15
237581	P 1782	0.04	U	mg/L	0.04			0.02545	11/12/03	12:15
237581	Pb220	0.005	U	mg/L	0.005			0.00000	11/12/03	12:15

Div 20
to#031106-12
06002.01.141

010038

system id	elem	result	qual	units	rl	tv	rec	ug/ml	date	time
237581	Pd3404	0.01	U	mg/L	0.01			0.00002	11/12/03	12:15
237581	S 1820	9.97		mg/L	0.025			9.96728	11/12/03	12:15
237581	Sb2068	0.01	U	mg/L	0.01			-0.00130	11/12/03	12:15
237581	Se196	0.015	U	mg/L	0.015			0.00785	11/12/03	12:15
237581	Si2881	7.10		mg/L	0.05			7.09798	11/12/03	12:15
237581	Sn1899	0.010		mg/L	0.01			0.01041	11/12/03	12:15
237581	Sr4215	0.006		mg/L	0.005			0.00643	11/12/03	12:15
237581	Th2837	0.01	U	mg/L	0.01			-0.00192	11/12/03	12:15
237581	Ti3372	0.005	U	mg/L	0.005			-0.00029	11/12/03	12:15
237581	Tl1908	0.025	U	mg/L	0.025			-0.02088	11/12/03	12:15
237581	U 3859	0.1	U	mg/L	0.1			0.04897	11/12/03	12:15
237581	V 2924	0.005	U	mg/L	0.005			0.00011	11/12/03	12:15
237581	W 2079	0.015	U	mg/L	0.015			-0.00309	11/12/03	12:15
237581	Y 3710	0.005	U	mg/L	0.005			0.00006	11/12/03	12:15
237581	Zn2062	0.005	U	mg/L	0.005			0.00229	11/12/03	12:15
237581	Zr3496	0.005	U	mg/L	0.005			0.00092	11/12/03	12:15
237581s	K 7664	23.0		mg/L	0.2	20	100.1%	23.0114	11/24/03	
237581s	Na3302	189		mg/L	0.2	20	82.3%	189.4536	11/24/03	
237582	Ag3280	0.005	U	mg/L	0.005			-0.00001	11/12/03	12:20
237582	Al3082	0.05	U	mg/L	0.05			0.02816	11/12/03	12:20
237582	As1890	0.008		mg/L	0.005			0.00770	11/12/03	12:20
237582	B 2496	0.233		mg/L	0.025			0.23320	11/12/03	12:20
237582	Ba4934	0.005	U	mg/L	0.005			0.00312	11/12/03	12:20
237582	Be3130	0.005	U	mg/L	0.005			-0.00001	11/12/03	12:20
237582	Bi2230	0.015	U	mg/L	0.015			0.00610	11/12/03	12:20
237582	Ca3179	2.24		mg/L	0.05			2.24230	11/12/03	12:20
237582	Cd2265	0.005	U	mg/L	0.005			0.00023	11/12/03	12:20
237582	Co2286	0.005	U	mg/L	0.005			0.00370	11/12/03	12:20
237582	Cr2677	0.005	U	mg/L	0.005			-0.00123	11/12/03	12:20
237582	Cu3247	0.005	U	mg/L	0.005			-0.00065	11/12/03	12:20
237582	Fe2714	0.202		mg/L	0.05			0.20248	11/12/03	12:20
237582	K 7664	3.08		mg/L	0.2			3.0764	11/24/03	
237582	La4086	0.005	U	mg/L	0.005			0.00219	11/12/03	12:20
237582	Li6707	0.247		mg/L	0.005			0.24661	11/12/03	12:20
237582	Mg2790	0.106		mg/L	0.05			0.10601	11/12/03	12:20
237582	Mn2576	0.039		mg/L	0.005			0.03876	11/12/03	12:20
237582	Mo2020	0.038		mg/L	0.01			0.03837	11/12/03	12:20
237582	Na3302	173		mg/L	0.2			172.8243	11/24/03	
237582	Ni2316	0.005	U	mg/L	0.005			-0.00012	11/12/03	12:20
237582	P 1782	0.04	U	mg/L	0.04			0.03902	11/12/03	12:20
237582	Pb220	0.005	U	mg/L	0.005			-0.00150	11/12/03	12:20
237582	Pd3404	0.01	U	mg/L	0.01			0.00464	11/12/03	12:20
237582	S 1820	9.89		mg/L	0.025			9.88691	11/12/03	12:20
237582	Sb2068	0.01	U	mg/L	0.01			-0.00021	11/12/03	12:20
237582	Se196	0.015	U	mg/L	0.015			0.00788	11/12/03	12:20
237582	Si2881	7.13		mg/L	0.05			7.13240	11/12/03	12:20
237582	Sn1899	0.025		mg/L	0.01			0.02542	11/12/03	12:20
237582	Sr4215	0.006		mg/L	0.005			0.00641	11/12/03	12:20
237582	Th2837	0.01	U	mg/L	0.01			-0.00171	11/12/03	12:20
237582	Ti3372	0.005	U	mg/L	0.005			0.00021	11/12/03	12:20
237582	Tl1908	0.025	U	mg/L	0.025			-0.02226	11/12/03	12:20
237582	U 3859	0.1	U	mg/L	0.1			0.06348	11/12/03	12:20
237582	V 2924	0.005	U	mg/L	0.005			0.00000	11/12/03	12:20
237582	W 2079	0.015	U	mg/L	0.015			-0.00088	11/12/03	12:20

system id	elem	result	qual	units	rl	tv	rec	ug/ml	date	time
237582	Y 3710	0.005	U	mg/L	0.005			0.00005	11/12/03	12:20
237582	Zn2062	0.005	U	mg/L	0.005			0.00118	11/12/03	12:20
237582	Zr3496	0.005	U	mg/L	0.005			0.00074	11/12/03	12:20
237583	Ag3280	0.005	U	mg/L	0.005			0.00005	11/12/03	12:25
237583	Al3082	0.096		mg/L	0.05			0.09599	11/12/03	12:25
237583	As1890	0.007		mg/L	0.005			0.00715	11/12/03	12:25
237583	B 2496	0.234		mg/L	0.025			0.23412	11/12/03	12:25
237583	Ba4934	0.005	U	mg/L	0.005			0.00423	11/12/03	12:25
237583	Be3130	0.005	U	mg/L	0.005			-0.00001	11/12/03	12:25
237583	Bi2230	0.015	U	mg/L	0.015			0.00111	11/12/03	12:25
237583	Ca3179	2.18		mg/L	0.05			2.17952	11/12/03	12:25
237583	Cd2265	0.005	U	mg/L	0.005			-0.00004	11/12/03	12:25
237583	Co2286	0.005	U	mg/L	0.005			0.00363	11/12/03	12:25
237583	Cr2677	0.005	U	mg/L	0.005			-0.00174	11/12/03	12:25
237583	Cu3247	0.005	U	mg/L	0.005			0.00005	11/12/03	12:25
237583	Fe2714	1.33		mg/L	0.05			1.32630	11/12/03	12:25
237583	K 7664	2.93		mg/L	0.2			2.9317	11/24/03	
237583	La4086	0.005	U	mg/L	0.005			0.00198	11/12/03	12:25
237583	Li6707	0.256		mg/L	0.005			0.25611	11/12/03	12:25
237583	Mg2790	0.111		mg/L	0.05			0.11050	11/12/03	12:25
237583	Mn2576	0.049		mg/L	0.005			0.04894	11/12/03	12:25
237583	Mo2020	0.035		mg/L	0.01			0.03477	11/12/03	12:25
237583	Na3302	173		mg/L	0.2			173.4907	11/24/03	
237583	Ni2316	0.005	U	mg/L	0.005			-0.00058	11/12/03	12:25
237583	P 1782	0.046		mg/L	0.04			0.04621	11/12/03	12:25
237583	Pb220	0.005	U	mg/L	0.005			-0.00124	11/12/03	12:25
237583	Pd3404	0.01	U	mg/L	0.01			0.00293	11/12/03	12:25
237583	S 1820	9.85		mg/L	0.025			9.85413	11/12/03	12:25
237583	Sb2068	0.01	U	mg/L	0.01			0.00353	11/12/03	12:25
237583	Se196	0.015	U	mg/L	0.015			0.01123	11/12/03	12:25
237583	Si2881	7.80		mg/L	0.05			7.80068	11/12/03	12:25
237583	Sn1899	0.024		mg/L	0.01			0.02423	11/12/03	12:25
237583	Sr4215	0.007		mg/L	0.005			0.00653	11/12/03	12:25
237583	Th2837	0.01	U	mg/L	0.01			-0.00086	11/12/03	12:25
237583	Ti3372	0.005	U	mg/L	0.005			0.00095	11/12/03	12:25
237583	Tl1908	0.025	U	mg/L	0.025			-0.02756	11/12/03	12:25
237583	U 3859	0.1	U	mg/L	0.1			0.02347	11/12/03	12:25
237583	V 2924	0.005	U	mg/L	0.005			0.00003	11/12/03	12:25
237583	W 2079	0.015	U	mg/L	0.015			-0.00304	11/12/03	12:25
237583	Y 3710	0.005	U	mg/L	0.005			-0.00005	11/12/03	12:25
237583	Zn2062	0.014		mg/L	0.005			0.01397	11/12/03	12:25
237583	Zr3496	0.005	U	mg/L	0.005			0.00010	11/12/03	12:25
237584	Ag3280	0.005	U	mg/L	0.005			0.00040	11/12/03	12:46
237584	Al3082	0.167		mg/L	0.05			0.16689	11/12/03	12:46
237584	As1890	0.035		mg/L	0.005			0.03486	11/12/03	12:46
237584	B 2496	0.179		mg/L	0.025			0.17861	11/12/03	12:46
237584	Ba4934	0.005	U	mg/L	0.005			0.00264	11/12/03	12:46
237584	Be3130	0.005	U	mg/L	0.005			-0.00016	11/12/03	12:46
237584	Bi2230	0.015	U	mg/L	0.015			0.00826	11/12/03	12:46
237584	Ca3179	0.543		mg/L	0.05			0.54330	11/12/03	12:46
237584	Cd2265	0.005	U	mg/L	0.005			0.00013	11/12/03	12:46
237584	Co2286	0.005	U	mg/L	0.005			0.00395	11/12/03	12:46
237584	Cr2677	0.005	U	mg/L	0.005			-0.00070	11/12/03	12:46
237584	Cu3247	0.005	U	mg/L	0.005			0.00018	11/12/03	12:46

Div 20
to#031106-12
06002.01.141

010058

system id	elem	result	qual	units	rl	tv	rec	ug/ml	date	time
237584	Fe2714	0.074		mg/L	0.05			0.07445	11/12/03	12:46
237584	K 7664	3.04		mg/L	0.2			3.0363	11/24/03	
237584	La4086	0.005	U	mg/L	0.005			0.00259	11/12/03	12:46
237584	Li6707	0.125		mg/L	0.005			0.12491	11/12/03	12:46
237584	Mg2790	0.071		mg/L	0.05			0.07063	11/12/03	12:46
237584	Mn2576	0.015		mg/L	0.005			0.01529	11/12/03	12:46
237584	Mo2020	0.01	U	mg/L	0.01			0.00422	11/12/03	12:46
237584	Na3302	94.7		mg/L	0.2			94.7452	11/24/03	
237584	Ni2316	0.005	U	mg/L	0.005			0.00015	11/12/03	12:46
237584	P 1782	0.04	U	mg/L	0.04			0.02897	11/12/03	12:46
237584	Pb220	0.005	U	mg/L	0.005			-0.00260	11/12/03	12:46
237584	Pd3404	0.01	U	mg/L	0.01			0.00373	11/12/03	12:46
237584	S 1820	5.56		mg/L	0.025			5.56027	11/12/03	12:46
237584	Sb2068	0.01	U	mg/L	0.01			0.00289	11/12/03	12:46
237584	Se196	0.015	U	mg/L	0.015			0.00640	11/12/03	12:46
237584	Si2881	27.9		mg/L	0.05			27.90196	11/12/03	12:46
237584	Sn1899	0.01	U	mg/L	0.01			0.00219	11/12/03	12:46
237584	Sr4215	0.005	U	mg/L	0.005			0.00192	11/12/03	12:46
237584	Th2837	0.01	U	mg/L	0.01			-0.00133	11/12/03	12:46
237584	Ti3372	0.005	U	mg/L	0.005			0.00116	11/12/03	12:46
237584	Ti1908	0.025	U	mg/L	0.025			0.01142	11/12/03	12:46
237584	U 3859	0.1	U	mg/L	0.1			0.07460	11/12/03	12:46
237584	V 2924	0.005		mg/L	0.005			0.00506	11/12/03	12:46
237584	W 2079	0.015	U	mg/L	0.015			0.00119	11/12/03	12:46
237584	Y 3710	0.005	U	mg/L	0.005			0.00013	11/12/03	12:46
237584	Zn2062	0.005	U	mg/L	0.005			0.00344	11/12/03	12:46
237584	Zr3496	0.005	U	mg/L	0.005			0.00123	11/12/03	12:46
237585	Ag3280	0.005	U	mg/L	0.005			0.00045	11/12/03	12:51
237585	Al3082	0.132		mg/L	0.05			0.13187	11/12/03	12:51
237585	As1890	0.034		mg/L	0.005			0.03385	11/12/03	12:51
237585	B 2496	0.162		mg/L	0.025			0.16160	11/12/03	12:51
237585	Ba4934	0.005	U	mg/L	0.005			0.00044	11/12/03	12:51
237585	Be3130	0.005	U	mg/L	0.005			-0.00018	11/12/03	12:51
237585	Bi2230	0.015	U	mg/L	0.015			0.00592	11/12/03	12:51
237585	Ca3179	0.538		mg/L	0.05			0.53801	11/12/03	12:51
237585	Cd2265	0.005	U	mg/L	0.005			0.00026	11/12/03	12:51
237585	Co2286	0.005	U	mg/L	0.005			0.00433	11/12/03	12:51
237585	Cr2677	0.005	U	mg/L	0.005			-0.00125	11/12/03	12:51
237585	Cu3247	0.005	U	mg/L	0.005			-0.00095	11/12/03	12:51
237585	Fe2714	0.05	U	mg/L	0.05			0.04260	11/12/03	12:51
237585	K 7664	2.92		mg/L	0.2			2.9233	11/24/03	
237585	La4086	0.005	U	mg/L	0.005			0.00239	11/12/03	12:51
237585	Li6707	0.125		mg/L	0.005			0.12451	11/12/03	12:51
237585	Mg2790	0.068		mg/L	0.05			0.06789	11/12/03	12:51
237585	Mn2576	0.014		mg/L	0.005			0.01392	11/12/03	12:51
237585	Mo2020	0.01	U	mg/L	0.01			0.00292	11/12/03	12:51
237585	Na3302	94.4		mg/L	0.2			94.418	11/24/03	
237585	Ni2316	0.005	U	mg/L	0.005			-0.00037	11/12/03	12:51
237585	P 1782	0.040		mg/L	0.04			0.04036	11/12/03	12:51
237585	Pb220	0.005	U	mg/L	0.005			-0.00285	11/12/03	12:51
237585	Pd3404	0.01	U	mg/L	0.01			0.00433	11/12/03	12:51
237585	S 1820	5.48		mg/L	0.025			5.48117	11/12/03	12:51
237585	Sb2068	0.01	U	mg/L	0.01			0.00226	11/12/03	12:51
237585	Se196	0.015	U	mg/L	0.015			0.00556	11/12/03	12:51

Div 20
to#031106-12
06002.01.141

010059

system id	elem	result	qual	units	rl	tv	rec	ug/ml	date	time
237585	Si2881	28.2		mg/L	0.05			28.16215	11/12/03	12:51
237585	Sn1899	0.01	U	mg/L	0.01			0.00780	11/12/03	12:51
237585	Sr4215	0.005	U	mg/L	0.005			0.00182	11/12/03	12:51
237585	Th2837	0.01	U	mg/L	0.01			-0.00327	11/12/03	12:51
237585	Ti3372	0.005	U	mg/L	0.005			0.00130	11/12/03	12:51
237585	Ti1908	0.025	U	mg/L	0.025			0.00850	11/12/03	12:51
237585	U 3859	0.1	U	mg/L	0.1			0.09038	11/12/03	12:51
237585	V 2924	0.005	U	mg/L	0.005			0.00410	11/12/03	12:51
237585	W 2079	0.015	U	mg/L	0.015			0.00274	11/12/03	12:51
237585	Y 3710	0.005	U	mg/L	0.005			-0.00007	11/12/03	12:51
237585	Zn2062	0.005	U	mg/L	0.005			0.00220	11/12/03	12:51
237585	Zr3496	0.005	U	mg/L	0.005			0.00040	11/12/03	12:51
237586	Ag3280	0.005	U	mg/L	0.005			0.00014	11/12/03	12:56
237586	Al3082	0.204		mg/L	0.05			0.20392	11/12/03	12:56
237586	As1890	0.033		mg/L	0.005			0.03324	11/12/03	12:56
237586	B 2496	0.163		mg/L	0.025			0.16308	11/12/03	12:56
237586	Ba4934	0.005	U	mg/L	0.005			0.00108	11/12/03	12:56
237586	Be3130	0.005	U	mg/L	0.005			-0.00013	11/12/03	12:56
237586	Bi2230	0.015	U	mg/L	0.015			0.00608	11/12/03	12:56
237586	Ca3179	0.761		mg/L	0.05			0.76082	11/12/03	12:56
237586	Cd2265	0.005	U	mg/L	0.005			-0.00003	11/12/03	12:56
237586	Co2286	0.005	U	mg/L	0.005			0.00328	11/12/03	12:56
237586	Cr2677	0.005	U	mg/L	0.005			-0.00093	11/12/03	12:56
237586	Cu3247	0.005	U	mg/L	0.005			-0.00027	11/12/03	12:56
237586	Fe2714	0.195		mg/L	0.05			0.19455	11/12/03	12:56
237586	K 7664	3.04		mg/L	0.2			3.0399	11/24/03	
237586	La4086	0.005	U	mg/L	0.005			0.00217	11/12/03	12:56
237586	Li6707	0.128		mg/L	0.005			0.12799	11/12/03	12:56
237586	Mg2790	0.094		mg/L	0.05			0.09406	11/12/03	12:56
237586	Mn2576	0.025		mg/L	0.005			0.02485	11/12/03	12:56
237586	Mo2020	0.01	U	mg/L	0.01			0.00261	11/12/03	12:56
237586	Na3302	94.8		mg/L	0.2			94.7542	11/24/03	
237586	Ni2316	0.005	U	mg/L	0.005			0.00009	11/12/03	12:56
237586	P 1782	0.045		mg/L	0.04			0.04487	11/12/03	12:56
237586	Pb220	0.005	U	mg/L	0.005			-0.00102	11/12/03	12:56
237586	Pd3404	0.01	U	mg/L	0.01			0.00333	11/12/03	12:56
237586	S 1820	5.63		mg/L	0.025			5.63160	11/12/03	12:56
237586	Sb2068	0.01	U	mg/L	0.01			-0.00049	11/12/03	12:56
237586	Se196	0.015	U	mg/L	0.015			0.00396	11/12/03	12:56
237586	Si2881	28.7		mg/L	0.05			28.74402	11/12/03	12:56
237586	Sn1899	0.01	U	mg/L	0.01			0.00617	11/12/03	12:56
237586	Sr4215	0.005	U	mg/L	0.005			0.00294	11/12/03	12:56
237586	Th2837	0.01	U	mg/L	0.01			-0.00192	11/12/03	12:56
237586	Ti3372	0.005	U	mg/L	0.005			0.00131	11/12/03	12:56
237586	Ti1908	0.025	U	mg/L	0.025			0.01264	11/12/03	12:56
237586	U 3859	0.1	U	mg/L	0.1			0.07838	11/12/03	12:56
237586	V 2924	0.005	U	mg/L	0.005			0.00491	11/12/03	12:56
237586	W 2079	0.015	U	mg/L	0.015			0.00053	11/12/03	12:56
237586	Y 3710	0.005	U	mg/L	0.005			0.00011	11/12/03	12:56
237586	Zn2062	0.019		mg/L	0.005			0.01937	11/12/03	12:56
237586	Zr3496	0.005	U	mg/L	0.005			0.00112	11/12/03	12:56
237587	Ag3280	0.005	U	mg/L	0.005			-0.00031	11/12/03	13:00
237587	Al3082	0.158		mg/L	0.05			0.15800	11/12/03	13:00
237587	As1890	0.037		mg/L	0.005			0.03730	11/12/03	13:00

Div 20
to#031106-12
06002.01.141

918820

system id	elem	result	qual	units	rl	tv	rec	ug/ml	date	time
237587	B 2496	0.157		mg/L	0.025			0.15683	11/12/03	13:00
237587	Ba4934	0.005	U	mg/L	0.005			0.00025	11/12/03	13:00
237587	Be3130	0.005	U	mg/L	0.005			-0.00017	11/12/03	13:00
237587	Bi2230	0.015	U	mg/L	0.015			0.00726	11/12/03	13:00
237587	Ca3179	0.360		mg/L	0.05			0.35965	11/12/03	13:00
237587	Cd2265	0.005	U	mg/L	0.005			0.00000	11/12/03	13:00
237587	Co2286	0.005	U	mg/L	0.005			0.00387	11/12/03	13:00
237587	Cr2677	0.005	U	mg/L	0.005			-0.00129	11/12/03	13:00
237587	Cu3247	0.005	U	mg/L	0.005			-0.00073	11/12/03	13:00
237587	Fe2714	0.074		mg/L	0.05			0.07428	11/12/03	13:00
237587	K 7664	2.96		mg/L	0.2			2.9562	11/24/03	
237587	La4086	0.005	U	mg/L	0.005			0.00145	11/12/03	13:00
237587	Li6707	0.126		mg/L	0.005			0.12569	11/12/03	13:00
237587	Mg2790	0.05	U	mg/L	0.05			0.04042	11/12/03	13:00
237587	Mn2576	0.005	U	mg/L	0.005			0.00383	11/12/03	13:00
237587	Mo2020	0.01	U	mg/L	0.01			0.00303	11/12/03	13:00
237587	Na3302	95.6		mg/L	0.2			95.6233	11/24/03	
237587	Ni2316	0.005	U	mg/L	0.005			-0.00029	11/12/03	13:00
237587	P 1782	0.040		mg/L	0.04			0.04010	11/12/03	13:00
237587	Pb220	0.005	U	mg/L	0.005			-0.00119	11/12/03	13:00
237587	Pd3404	0.01	U	mg/L	0.01			0.00419	11/12/03	13:00
237587	S 1820	5.69		mg/L	0.025			5.69496	11/12/03	13:00
237587	Sb2068	0.01	U	mg/L	0.01			-0.00043	11/12/03	13:00
237587	Se196	0.015	U	mg/L	0.015			0.00007	11/12/03	13:00
237587	Si2881	28.5		mg/L	0.05			28.50282	11/12/03	13:00
237587	Sn1899	0.01	U	mg/L	0.01			0.00083	11/12/03	13:00
237587	Sr4215	0.005	U	mg/L	0.005			0.00085	11/12/03	13:00
237587	Th2837	0.01	U	mg/L	0.01			-0.00099	11/12/03	13:00
237587	Ti3372	0.005	U	mg/L	0.005			0.00110	11/12/03	13:00
237587	Tl1908	0.025	U	mg/L	0.025			0.00640	11/12/03	13:00
237587	U 3859	0.1	U	mg/L	0.1			0.07039	11/12/03	13:00
237587	V 2924	0.005	U	mg/L	0.005			0.00486	11/12/03	13:00
237587	W 2079	0.015	U	mg/L	0.015			0.00234	11/12/03	13:00
237587	Y 3710	0.005	U	mg/L	0.005			-0.00002	11/12/03	13:00
237587	Zn2062	0.005	U	mg/L	0.005			0.00339	11/12/03	13:00
237587	Zr3496	0.005	U	mg/L	0.005			0.00138	11/12/03	13:00
237588	Ag3280	0.005	U	mg/L	0.005			0.00001	11/12/03	13:05
237588	Al3082	0.120		mg/L	0.05			0.12017	11/12/03	13:05
237588	As1890	0.035		mg/L	0.005			0.03467	11/12/03	13:05
237588	B 2496	0.157		mg/L	0.025			0.15656	11/12/03	13:05
237588	Ba4934	0.005	U	mg/L	0.005			0.00029	11/12/03	13:05
237588	Be3130	0.005	U	mg/L	0.005			-0.00004	11/12/03	13:05
237588	Bi2230	0.015	U	mg/L	0.015			0.00851	11/12/03	13:05
237588	Ca3179	0.359		mg/L	0.05			0.35877	11/12/03	13:05
237588	Cd2265	0.005	U	mg/L	0.005			0.00012	11/12/03	13:05
237588	Co2286	0.005	U	mg/L	0.005			0.00382	11/12/03	13:05
237588	Cr2677	0.005	U	mg/L	0.005			-0.00112	11/12/03	13:05
237588	Cu3247	0.005	U	mg/L	0.005			-0.00096	11/12/03	13:05
237588	Fe2714	0.05	U	mg/L	0.05			0.04585	11/12/03	13:05
237588	K 7664	3.09		mg/L	0.2			3.0895	11/24/03	
237588	La4086	0.005	U	mg/L	0.005			0.00193	11/12/03	13:05
237588	Li6707	0.127		mg/L	0.005			0.12680	11/12/03	13:05
237588	Mg2790	0.05	U	mg/L	0.05			0.03637	11/12/03	13:05
237588	Mn2576	0.005	U	mg/L	0.005			0.00406	11/12/03	13:05

Div 20
to#031106-12
06002.01.141

010081

system id	elem	result	qual	units	rl	tv	rec	ug/ml	date	time
237588	Mo2020	0.01	U	mg/L	0.01			0.00272	11/12/03	13:05
237588	Na3302	97.1		mg/L	0.2			97.0979	11/24/03	
237588	Ni2316	0.005	U	mg/L	0.005			-0.00051	11/12/03	13:05
237588	P 1782	0.04	U	mg/L	0.04			0.03274	11/12/03	13:05
237588	Pb220	0.005	U	mg/L	0.005			-0.00244	11/12/03	13:05
237588	Pd3404	0.01	U	mg/L	0.01			0.00332	11/12/03	13:05
237588	S 1820	5.81		mg/L	0.025			5.80537	11/12/03	13:05
237588	Sb2068	0.01	U	mg/L	0.01			0.00023	11/12/03	13:05
237588	Se196	0.015	U	mg/L	0.015			0.00136	11/12/03	13:05
237588	Si2881	28.5		mg/L	0.05			28.50225	11/12/03	13:05
237588	Sn1899	0.01	U	mg/L	0.01			0.00611	11/12/03	13:05
237588	Sr4215	0.005	U	mg/L	0.005			0.00087	11/12/03	13:05
237588	Th2837	0.01	U	mg/L	0.01			-0.00190	11/12/03	13:05
237588	Ti3372	0.005	U	mg/L	0.005			0.00090	11/12/03	13:05
237588	Ti1908	0.025	U	mg/L	0.025			0.01041	11/12/03	13:05
237588	U 3859	0.1	U	mg/L	0.1			0.07138	11/12/03	13:05
237588	V 2924	0.005		mg/L	0.005			0.00528	11/12/03	13:05
237588	W 2079	0.015	U	mg/L	0.015			0.00175	11/12/03	13:05
237588	Y 3710	0.005	U	mg/L	0.005			0.00001	11/12/03	13:05
237588	Zn2062	0.005	U	mg/L	0.005			0.00262	11/12/03	13:05
237588	Zr3496	0.005	U	mg/L	0.005			0.00204	11/12/03	13:05
237589	Ag3280	0.005	U	mg/L	0.005			0.00029	11/12/03	13:10
237589	Al3082	0.174		mg/L	0.05			0.17385	11/12/03	13:10
237589	As1890	0.032		mg/L	0.005			0.03159	11/12/03	13:10
237589	B 2496	0.157		mg/L	0.025			0.15731	11/12/03	13:10
237589	Ba4934	0.005	U	mg/L	0.005			0.00049	11/12/03	13:10
237589	Be3130	0.005	U	mg/L	0.005			-0.00008	11/12/03	13:10
237589	Bi2230	0.015	U	mg/L	0.015			0.00836	11/12/03	13:10
237589	Ca3179	0.409		mg/L	0.05			0.40912	11/12/03	13:10
237589	Cd2265	0.005	U	mg/L	0.005			0.00004	11/12/03	13:10
237589	Co2286	0.005	U	mg/L	0.005			0.00418	11/12/03	13:10
237589	Cr2677	0.005	U	mg/L	0.005			-0.00053	11/12/03	13:10
237589	Cu3247	0.005	U	mg/L	0.005			-0.00038	11/12/03	13:10
237589	Fe2714	0.190		mg/L	0.05			0.18990	11/12/03	13:10
237589	K 7664	3.11		mg/L	0.2			3.1053	11/24/03	
237589	La4086	0.005	U	mg/L	0.005			0.00193	11/12/03	13:10
237589	Li6707	0.127		mg/L	0.005			0.12693	11/12/03	13:10
237589	Mg2790	0.05	U	mg/L	0.05			0.04487	11/12/03	13:10
237589	Mn2576	0.034		mg/L	0.005			0.03400	11/12/03	13:10
237589	Mo2020	0.01	U	mg/L	0.01			0.00247	11/12/03	13:10
237589	Na3302	97.3		mg/L	0.2			97.3213	11/24/03	
237589	Ni2316	0.005	U	mg/L	0.005			0.00003	11/12/03	13:10
237589	P 1782	0.04	U	mg/L	0.04			0.02992	11/12/03	13:10
237589	Pb220	0.005	U	mg/L	0.005			-0.00105	11/12/03	13:10
237589	Pd3404	0.01	U	mg/L	0.01			0.00234	11/12/03	13:10
237589	S 1820	5.69		mg/L	0.025			5.68979	11/12/03	13:10
237589	Sb2068	0.01	U	mg/L	0.01			0.00081	11/12/03	13:10
237589	Se196	0.015	U	mg/L	0.015			0.00498	11/12/03	13:10
237589	Si2881	28.7		mg/L	0.05			28.72847	11/12/03	13:10
237589	Sn1899	0.01	U	mg/L	0.01			0.00963	11/12/03	13:10
237589	Sr4215	0.005	U	mg/L	0.005			0.00107	11/12/03	13:10
237589	Th2837	0.01	U	mg/L	0.01			-0.00255	11/12/03	13:10
237589	Ti3372	0.005	U	mg/L	0.005			0.00105	11/12/03	13:10
237589	Ti1908	0.025	U	mg/L	0.025			-0.01542	11/12/03	13:10

Div 20
to#031106-12
06002.01.141

010062

system id	elem	result	qual	units	rl	tv	rec	ug/ml	date	time
237589	U 3859	0.1	U	mg/L	0.1			0.08163	11/12/03	13:10
237589	V 2924	0.005		mg/L	0.005			0.00545	11/12/03	13:10
237589	W 2079	0.015	U	mg/L	0.015			-0.00159	11/12/03	13:10
237589	Y 3710	0.005	U	mg/L	0.005			0.00011	11/12/03	13:10
237589	Zn2062	0.006		mg/L	0.005			0.00616	11/12/03	13:10
237589	Zr3496	0.005	U	mg/L	0.005			0.00137	11/12/03	13:10
237590	Ag3280	0.005	U	mg/L	0.005			0.00012	11/12/03	13:15
237590	Al3082	0.05	U	mg/L	0.05			0.04353	11/12/03	13:15
237590	As1890	0.005	U	mg/L	0.005			0.00035	11/12/03	13:15
237590	B 2496	0.259		mg/L	0.025			0.25874	11/12/03	13:15
237590	Ba4934	0.056		mg/L	0.005			0.05606	11/12/03	13:15
237590	Be3130	0.005	U	mg/L	0.005			-0.00035	11/12/03	13:15
237590	Bi2230	0.015	U	mg/L	0.015			0.00460	11/12/03	13:15
237590	Ca3179	28.1		mg/L	0.05			28.14830	11/12/03	13:15
237590	Cd2265	0.005	U	mg/L	0.005			0.00018	11/12/03	13:15
237590	Co2286	0.005	U	mg/L	0.005			0.00462	11/12/03	13:15
237590	Cr2677	0.005	U	mg/L	0.005			-0.00132	11/12/03	13:15
237590	Cu3247	0.005	U	mg/L	0.005			-0.00044	11/12/03	13:15
237590	Fe2714	0.05	U	mg/L	0.05			0.02468	11/12/03	13:15
237590	K 7664	8.08		mg/L	0.2			8.079	11/24/03	
237590	La4086	0.005	U	mg/L	0.005			0.00227	11/12/03	13:15
237590	Li6707	0.105		mg/L	0.005			0.10455	11/12/03	13:15
237590	Mg2790	28.9		mg/L	0.05			28.89275	11/12/03	13:15
237590	Mn2576	0.075		mg/L	0.005			0.07487	11/12/03	13:15
237590	Mo2020	0.01	U	mg/L	0.01			0.00552	11/12/03	13:15
237590	Na3302	82.2		mg/L	0.2			82.1895	11/24/03	
237590	Ni2316	0.005	U	mg/L	0.005			0.00118	11/12/03	13:15
237590	P 1782	0.04	U	mg/L	0.04			0.02175	11/12/03	13:15
237590	Pb220	0.005	U	mg/L	0.005			-0.00189	11/12/03	13:15
237590	Pd3404	0.01	U	mg/L	0.01			0.00464	11/12/03	13:15
237590	S 1820	45.0		mg/L	0.025			45.01411	11/12/03	13:15
237590	Sb2068	0.01	U	mg/L	0.01			0.00131	11/12/03	13:15
237590	Se196	0.015	U	mg/L	0.015			0.00631	11/12/03	13:15
237590	Si2881	13.7		mg/L	0.05			13.68723	11/12/03	13:15
237590	Sn1899	0.01	U	mg/L	0.01			0.00742	11/12/03	13:15
237590	Sr4215	0.413		mg/L	0.005			0.41273	11/12/03	13:15
237590	Th2837	0.01	U	mg/L	0.01			-0.00739	11/12/03	13:15
237590	Ti3372	0.005	U	mg/L	0.005			0.00064	11/12/03	13:15
237590	Ti1908	0.025	U	mg/L	0.025			0.00545	11/12/03	13:15
237590	U 3859	0.1	U	mg/L	0.1			0.05746	11/12/03	13:15
237590	V 2924	0.005	U	mg/L	0.005			-0.00103	11/12/03	13:15
237590	W 2079	0.015	U	mg/L	0.015			0.00481	11/12/03	13:15
237590	Y 3710	0.005	U	mg/L	0.005			0.00003	11/12/03	13:15
237590	Zn2062	0.005	U	mg/L	0.005			0.00189	11/12/03	13:15
237590	Zr3496	0.005	U	mg/L	0.005			0.00188	11/12/03	13:15
237591	Ag3280	0.005	U	mg/L	0.005			0.00202	11/12/03	13:20
237591	Al3082	0.05	U	mg/L	0.05			0.04931	11/12/03	13:20
237591	As1890	0.005		mg/L	0.005			0.00506	11/12/03	13:20
237591	B 2496	0.260		mg/L	0.025			0.26003	11/12/03	13:20
237591	Ba4934	0.056		mg/L	0.005			0.05610	11/12/03	13:20
237591	Be3130	0.005	U	mg/L	0.005			-0.00024	11/12/03	13:20
237591	Bi2230	0.015	U	mg/L	0.015			0.01187	11/12/03	13:20
237591	Ca3179	28.3		mg/L	0.05			28.28334	11/12/03	13:20
237591	Cd2265	0.005	U	mg/L	0.005			0.00059	11/12/03	13:20

Div 20
to#031106-12
06002.01.141

010063

system id	elem	result	qual	units	rl	tv	rec	ug/ml	date	time
237591	Co2286	0.005		mg/L	0.005			0.00509	11/12/03	13:20
237591	Cr2677	0.005	U	mg/L	0.005			-0.00203	11/12/03	13:20
237591	Cu3247	0.005	U	mg/L	0.005			0.00192	11/12/03	13:20
237591	Fe2714	0.05	U	mg/L	0.05			0.01057	11/12/03	13:20
237591	K 7664	8.16		mg/L	0.2			8.1617	11/24/03	
237591	La4086	0.005	U	mg/L	0.005			0.00375	11/12/03	13:20
237591	Li6707	0.104		mg/L	0.005			0.10407	11/12/03	13:20
237591	Mg2790	28.9		mg/L	0.05			28.90833	11/12/03	13:20
237591	Mn2576	0.071		mg/L	0.005			0.07132	11/12/03	13:20
237591	Mo2020	0.01	U	mg/L	0.01			0.00484	11/12/03	13:20
237591	Na3302	82.0		mg/L	0.2			82.0032	11/24/03	
237591	Ni2316	0.005	U	mg/L	0.005			-0.00047	11/12/03	13:20
237591	P 1782	0.042		mg/L	0.04			0.04198	11/12/03	13:20
237591	Pb220	0.005	U	mg/L	0.005			-0.00234	11/12/03	13:20
237591	Pd3404	0.01	U	mg/L	0.01			0.00656	11/12/03	13:20
237591	S 1820	43.7		mg/L	0.025			43.74258	11/12/03	13:20
237591	Sb2068	0.01	U	mg/L	0.01			0.00859	11/12/03	13:20
237591	Se196	0.015	U	mg/L	0.015			0.00491	11/12/03	13:20
237591	Si2881	13.7		mg/L	0.05			13.65035	11/12/03	13:20
237591	Sn1899	0.016		mg/L	0.01			0.01634	11/12/03	13:20
237591	Sr4215	0.413		mg/L	0.005			0.41301	11/12/03	13:20
237591	Th2837	0.01	U	mg/L	0.01			-0.00822	11/12/03	13:20
237591	Ti3372	0.005	U	mg/L	0.005			0.00147	11/12/03	13:20
237591	Tl1908	0.025	U	mg/L	0.025			0.00383	11/12/03	13:20
237591	U 3859	0.1	U	mg/L	0.1			0.09174	11/12/03	13:20
237591	V 2924	0.005	U	mg/L	0.005			-0.00123	11/12/03	13:20
237591	W 2079	0.015	U	mg/L	0.015			0.01042	11/12/03	13:20
237591	Y 3710	0.005	U	mg/L	0.005			-0.00005	11/12/03	13:20
237591	Zn2062	0.005	U	mg/L	0.005			0.00132	11/12/03	13:20
237591	Zr3496	0.005	U	mg/L	0.005			0.00091	11/12/03	13:20
237592	Ag3280	0.005	U	mg/L	0.005			-0.00164	11/12/03	13:25
237592	Al3082	0.05	U	mg/L	0.05			0.03167	11/12/03	13:25
237592	As1890	0.005	U	mg/L	0.005			0.00053	11/12/03	13:25
237592	B 2496	0.260		mg/L	0.025			0.25969	11/12/03	13:25
237592	Ba4934	0.057		mg/L	0.005			0.05670	11/12/03	13:25
237592	Be3130	0.005	U	mg/L	0.005			-0.00018	11/12/03	13:25
237592	Bi2230	0.015	U	mg/L	0.015			0.00404	11/12/03	13:25
237592	Ca3179	28.2		mg/L	0.05			28.23118	11/12/03	13:25
237592	Cd2265	0.005	U	mg/L	0.005			-0.00013	11/12/03	13:25
237592	Co2286	0.005	U	mg/L	0.005			0.00361	11/12/03	13:25
237592	Cr2677	0.005	U	mg/L	0.005			-0.00040	11/12/03	13:25
237592	Cu3247	0.005	U	mg/L	0.005			0.00001	11/12/03	13:25
237592	Fe2714	0.081		mg/L	0.05			0.08110	11/12/03	13:25
237592	K 7664	8.12		mg/L	0.2			8.1205	11/24/03	
237592	La4086	0.005	U	mg/L	0.005			0.00072	11/12/03	13:25
237592	Li6707	0.105		mg/L	0.005			0.10467	11/12/03	13:25
237592	Mg2790	28.9		mg/L	0.05			28.88541	11/12/03	13:25
237592	Mn2576	0.071		mg/L	0.005			0.07123	11/12/03	13:25
237592	Mo2020	0.01	U	mg/L	0.01			0.00626	11/12/03	13:25
237592	Na3302	81.6		mg/L	0.2			81.5987	11/24/03	
237592	Ni2316	0.005	U	mg/L	0.005			0.00073	11/12/03	13:25
237592	P 1782	0.04	U	mg/L	0.04			0.03380	11/12/03	13:25
237592	Pb220	0.005	U	mg/L	0.005			-0.00061	11/12/03	13:25
237592	Pd3404	0.01	U	mg/L	0.01			0.00278	11/12/03	13:25

Div 20
to#031106-12
06002.01.141

system id	elem	result	qual	units	rl	tv	rec	ug/ml	date	time
237592	S 1820	44.3		mg/L	0.025			44.27149	11/12/03	13:25
237592	Sb2068	0.01	U	mg/L	0.01			-0.00229	11/12/03	13:25
237592	Se196	0.015	U	mg/L	0.015			0.00239	11/12/03	13:25
237592	Si2881	13.7		mg/L	0.05			13.73487	11/12/03	13:25
237592	Sn1899	0.01	U	mg/L	0.01			0.00984	11/12/03	13:25
237592	Sr4215	0.416		mg/L	0.005			0.41608	11/12/03	13:25
237592	Th2837	0.01	U	mg/L	0.01			-0.00623	11/12/03	13:25
237592	Ti3372	0.005	U	mg/L	0.005			-0.00045	11/12/03	13:25
237592	Tl1908	0.025	U	mg/L	0.025			-0.01200	11/12/03	13:25
237592	U 3859	0.1	U	mg/L	0.1			0.04218	11/12/03	13:25
237592	V 2924	0.005	U	mg/L	0.005			0.00053	11/12/03	13:25
237592	W 2079	0.015	U	mg/L	0.015			0.00616	11/12/03	13:25
237592	Y 3710	0.005	U	mg/L	0.005			0.00025	11/12/03	13:25
237592	Zn2062	0.016		mg/L	0.005			0.01636	11/12/03	13:25
237592	Zr3496	0.005	U	mg/L	0.005			0.00223	11/12/03	13:25
237593	Ag3280	0.005	U	mg/L	0.005			0.00061	11/12/03	13:30
237593	Al3082	0.05	U	mg/L	0.05			0.04240	11/12/03	13:30
237593	As1890	0.005	U	mg/L	0.005			0.00245	11/12/03	13:30
237593	B 2496	0.025	U	mg/L	0.025			0.00276	11/12/03	13:30
237593	Ba4934	0.005	U	mg/L	0.005			0.00035	11/12/03	13:30
237593	Be3130	0.005	U	mg/L	0.005			-0.00024	11/12/03	13:30
237593	Bi2230	0.015	U	mg/L	0.015			0.00456	11/12/03	13:30
237593	Ca3179	10.1		mg/L	0.05			10.06773	11/12/03	13:30
237593	Cd2265	0.005	U	mg/L	0.005			0.00008	11/12/03	13:30
237593	Co2286	0.005	U	mg/L	0.005			0.00376	11/12/03	13:30
237593	Cr2677	0.005	U	mg/L	0.005			-0.00117	11/12/03	13:30
237593	Cu3247	0.005	U	mg/L	0.005			0.00354	11/12/03	13:30
237593	Fe2714	0.05	U	mg/L	0.05			-0.00052	11/12/03	13:30
237593	K 7664	5.25		mg/L	0.2			5.251	11/24/03	
237593	La4086	0.005	U	mg/L	0.005			0.00243	11/12/03	13:30
237593	Li6707	0.005	U	mg/L	0.005			0.00028	11/12/03	13:30
237593	Mg2790	0.065		mg/L	0.05			0.06452	11/12/03	13:30
237593	Mn2576	0.005	U	mg/L	0.005			0.00007	11/12/03	13:30
237593	Mo2020	0.01	U	mg/L	0.01			-0.00126	11/12/03	13:30
237593	Na3302	27.0		mg/L	0.2			27.0296	11/24/03	
237593	Ni2316	0.005	U	mg/L	0.005			0.00058	11/12/03	13:30
237593	P 1782	0.04	U	mg/L	0.04			0.01365	11/12/03	13:30
237593	Pb220	0.005	U	mg/L	0.005			-0.00062	11/12/03	13:30
237593	Pd3404	0.01	U	mg/L	0.01			0.00295	11/12/03	13:30
237593	S 1820	0.114		mg/L	0.025			0.11411	11/12/03	13:30
237593	Sb2068	0.01	U	mg/L	0.01			0.00315	11/12/03	13:30
237593	Se196	0.015	U	mg/L	0.015			0.00461	11/12/03	13:30
237593	Si2881	19.4		mg/L	0.05			19.37984	11/12/03	13:30
237593	Sn1899	0.01	U	mg/L	0.01			0.00298	11/12/03	13:30
237593	Sr4215	0.005	U	mg/L	0.005			0.00127	11/12/03	13:30
237593	Th2837	0.01	U	mg/L	0.01			-0.00160	11/12/03	13:30
237593	Ti3372	0.005	U	mg/L	0.005			0.00018	11/12/03	13:30
237593	Tl1908	0.025	U	mg/L	0.025			0.00549	11/12/03	13:30
237593	U 3859	0.1	U	mg/L	0.1			0.08671	11/12/03	13:30
237593	V 2924	0.005	U	mg/L	0.005			-0.00008	11/12/03	13:30
237593	W 2079	0.015	U	mg/L	0.015			0.00127	11/12/03	13:30
237593	Y 3710	0.005	U	mg/L	0.005			0.00004	11/12/03	13:30
237593	Zn2062	0.005	U	mg/L	0.005			0.00202	11/12/03	13:30
237593	Zr3496	0.008		mg/L	0.005			0.00762	11/12/03	13:30

010054

Div 20
to#031106-12
06002.01.141

010063

system id	elem	result	qual	units	rl	tv	rec	ug/ml	date	time
237594	Ag3280	0.005	U	mg/L	0.005			0.00008	11/12/03	13:51
237594	Al3082	0.05	U	mg/L	0.05			0.03466	11/12/03	13:51
237594	As1890	0.005	U	mg/L	0.005			0.00297	11/12/03	13:51
237594	B 2496	0.025	U	mg/L	0.025			0.00672	11/12/03	13:51
237594	Ba4934	0.005	U	mg/L	0.005			0.00033	11/12/03	13:51
237594	Be3130	0.005	U	mg/L	0.005			-0.00026	11/12/03	13:51
237594	Bi2230	0.015	U	mg/L	0.015			0.00643	11/12/03	13:51
237594	Ca3179	9.91		mg/L	0.05			9.91204	11/12/03	13:51
237594	Cd2265	0.005	U	mg/L	0.005			-0.00004	11/12/03	13:51
237594	Co2286	0.005	U	mg/L	0.005			0.00355	11/12/03	13:51
237594	Cr2677	0.005	U	mg/L	0.005			-0.00149	11/12/03	13:51
237594	Cu3247	0.005	U	mg/L	0.005			-0.00047	11/12/03	13:51
237594	Fe2714	0.05	U	mg/L	0.05			-0.01388	11/12/03	13:51
237594	K 7664	5.32		mg/L	0.2			5.3166	11/24/03	
237594	La4086	0.005	U	mg/L	0.005			0.00178	11/12/03	13:51
237594	Li6707	0.005	U	mg/L	0.005			0.00000	11/12/03	13:51
237594	Mg2790	0.05	U	mg/L	0.05			0.01025	11/12/03	13:51
237594	Mn2576	0.005	U	mg/L	0.005			-0.00010	11/12/03	13:51
237594	Mo2020	0.01	U	mg/L	0.01			0.00021	11/12/03	13:51
237594	Na3302	26.6		mg/L	0.2			26.6394	11/24/03	
237594	Ni2316	0.005	U	mg/L	0.005			-0.00015	11/12/03	13:51
237594	P 1782	0.04	U	mg/L	0.04			0.01818	11/12/03	13:51
237594	Pb220	0.005	U	mg/L	0.005			-0.00106	11/12/03	13:51
237594	Pd3404	0.01	U	mg/L	0.01			0.00366	11/12/03	13:51
237594	S 1820	0.025	U	mg/L	0.025			0.00859	11/12/03	13:51
237594	Sb2068	0.01	U	mg/L	0.01			-0.00077	11/12/03	13:51
237594	Se196	0.015	U	mg/L	0.015			0.00639	11/12/03	13:51
237594	Si2881	19.4		mg/L	0.05			19.40569	11/12/03	13:51
237594	Sn1899	0.01	U	mg/L	0.01			0.00162	11/12/03	13:51
237594	Sr4215	0.005	U	mg/L	0.005			0.00056	11/12/03	13:51
237594	Th2837	0.01	U	mg/L	0.01			-0.00252	11/12/03	13:51
237594	Ti3372	0.005	U	mg/L	0.005			0.00005	11/12/03	13:51
237594	Ti1908	0.025	U	mg/L	0.025			0.02433	11/12/03	13:51
237594	U 3859	0.1	U	mg/L	0.1			0.06365	11/12/03	13:51
237594	V 2924	0.005	U	mg/L	0.005			-0.00018	11/12/03	13:51
237594	W 2079	0.015	U	mg/L	0.015			0.00096	11/12/03	13:51
237594	Y 3710	0.005	U	mg/L	0.005			0.00004	11/12/03	13:51
237594	Zn2062	0.005	U	mg/L	0.005			0.00217	11/12/03	13:51
237594	Zr3496	0.005		mg/L	0.005			0.00520	11/12/03	13:51
lcsW-L11W1	Ag3280	0.052		mg/L	0.005	0.05	104.4%	0.05221	11/12/03	11:46
lcsW-L11W1	Al3082	1.92		mg/L	0.05	2	95.9%	1.91868	11/12/03	11:46
lcsW-L11W1	As1890	2.08		mg/L	0.005	2	104.1%	2.08211	11/12/03	11:46
lcsW-L11W1	B 2496	0.025	U	mg/L	0.025			0.00598	11/12/03	11:46
lcsW-L11W1	Ba4934	2.02		mg/L	0.005	2	101.2%	2.02382	11/12/03	11:46
lcsW-L11W1	Be3130	0.054		mg/L	0.005	0.05	108.4%	0.05419	11/12/03	11:46
lcsW-L11W1	Bi2230	0.015	U	mg/L	0.015			0.00685	11/12/03	11:46
lcsW-L11W1	Ca3179	20.8		mg/L	0.05	20	104.0%	20.79470	11/12/03	11:46
lcsW-L11W1	Cd2265	0.052		mg/L	0.005	0.05	104.9%	0.05243	11/12/03	11:46
lcsW-L11W1	Co2286	0.515		mg/L	0.005	0.5	103.0%	0.51509	11/12/03	11:46
lcsW-L11W1	Cr2677	0.203		mg/L	0.005	0.2	101.7%	0.20337	11/12/03	11:46
lcsW-L11W1	Cu3247	0.251		mg/L	0.005	0.25	100.4%	0.25093	11/12/03	11:46
lcsW-L11W1	Fe2714	1.03		mg/L	0.05	1	103.1%	1.03091	11/12/03	11:46
lcsW-L24W1	K 7664	18.3		mg/L	0.2	20	91.4%	18.282	11/24/03	
lcsW-L11W1	La4086	0.005	U	mg/L	0.005			0.00144	11/12/03	11:46

Div 20
to#031106-12
06002.01.141


010060

system id	elem	result	qual	units	rl	tv	rec	ug/ml	date	time
lcsW-L11W1	Li6707	0.005	U	mg/L	0.005			0.00030	11/12/03	11:46
lcsW-L11W1	Mg2790	20.4		mg/L	0.05	20	101.9%	20.38276	11/12/03	11:46
lcsW-L11W1	Mn2576	0.516		mg/L	0.005	0.5	103.1%	0.51556	11/12/03	11:46
lcsW-L11W1	Mo2020	0.01	U	mg/L	0.01			0.00140	11/12/03	11:46
lcsW-L24W1	Na3302	18.6		mg/L	0.2	20	93.2%	18.6344	11/24/03	
lcsW-L11W1	Ni2316	0.511		mg/L	0.005	0.5	102.2%	0.51096	11/12/03	11:46
lcsW-L11W1	P 1782	0.04	U	mg/L	0.04			0.01782	11/12/03	11:46
lcsW-L11W1	Pb220	0.527		mg/L	0.005	0.5	105.4%	0.52704	11/12/03	11:46
lcsW-L11W1	Pd3404	0.01	U	mg/L	0.01			0.00705	11/12/03	11:46
lcsW-L11W1	S 1820	0.025	U	mg/L	0.025			0.02038	11/12/03	11:46
lcsW-L11W1	Sb2068	0.515		mg/L	0.01	0.5	102.9%	0.51453	11/12/03	11:46
lcsW-L11W1	Se196	2.15		mg/L	0.015	2	107.3%	2.14626	11/12/03	11:46
lcsW-L11W1	Si2881	0.05	U	mg/L	0.05			0.01599	11/12/03	11:46
lcsW-L11W1	Sn1899	0.01	U	mg/L	0.01			0.00294	11/12/03	11:46
lcsW-L11W1	Sr4215	0.005	U	mg/L	0.005			0.00030	11/12/03	11:46
lcsW-L11W1	Th2837	0.01	U	mg/L	0.01			-0.00654	11/12/03	11:46
lcsW-L11W1	Ti3372	0.005	U	mg/L	0.005			0.00071	11/12/03	11:46
lcsW-L11W1	Ti1908	2.22		mg/L	0.025	2	110.9%	2.21897	11/12/03	11:46
lcsW-L11W1	U 3859	0.1	U	mg/L	0.1			0.07891	11/12/03	11:46
lcsW-L11W1	V 2924	0.510		mg/L	0.005	0.5	102.0%	0.51012	11/12/03	11:46
lcsW-L11W1	W 2079	0.015	U	mg/L	0.015			0.00184	11/12/03	11:46
lcsW-L11W1	Y 3710	0.005	U	mg/L	0.005			0.00015	11/12/03	11:46
lcsW-L11W1	Zn2062	0.531		mg/L	0.005	0.5	106.2%	0.53080	11/12/03	11:46
lcsW-L11W1	Zr3496	0.005	U	mg/L	0.005			0.00140	11/12/03	11:46
pbw-L11W1	Ag3280	0.005	U	mg/L	0.005			0.00070	11/12/03	11:42
pbw-L11W1	Al3082	0.05	U	mg/L	0.05			0.01881	11/12/03	11:42
pbw-L11W1	As1890	0.005	U	mg/L	0.005			0.00151	11/12/03	11:42
pbw-L11W1	B 2496	0.025	U	mg/L	0.025			0.00927	11/12/03	11:42
pbw-L11W1	Ba4934	0.005	U	mg/L	0.005			0.00012	11/12/03	11:42
pbw-L11W1	Be3130	0.005	U	mg/L	0.005			-0.00005	11/12/03	11:42
pbw-L11W1	Bi2230	0.015	U	mg/L	0.015			0.00563	11/12/03	11:42
pbw-L11W1	Ca3179	0.05	U	mg/L	0.05			0.03755	11/12/03	11:42
pbw-L11W1	Cd2265	0.005	U	mg/L	0.005			0.00018	11/12/03	11:42
pbw-L11W1	Co2286	0.005	U	mg/L	0.005			0.00306	11/12/03	11:42
pbw-L11W1	Cr2677	0.005	U	mg/L	0.005			-0.00164	11/12/03	11:42
pbw-L11W1	Cu3247	0.005	U	mg/L	0.005			0.00019	11/12/03	11:42
pbw-L11W1	Fe2714	0.05	U	mg/L	0.05			-0.00499	11/12/03	11:42
pbw-L24W2	K 7664	0.2	U	mg/L	0.2			-0.1632	11/24/03	
pbw-L11W1	La4086	0.005	U	mg/L	0.005			0.00104	11/12/03	11:42
pbw-L11W1	Li6707	0.005	U	mg/L	0.005			0.00016	11/12/03	11:42
pbw-L11W1	Mg2790	0.05	U	mg/L	0.05			0.01320	11/12/03	11:42
pbw-L11W1	Mn2576	0.005	U	mg/L	0.005			0.00012	11/12/03	11:42
pbw-L11W1	Mo2020	0.01	U	mg/L	0.01			0.00290	11/12/03	11:42
pbw-L24W2	Na3302	0.2	U	mg/L	0.2			-0.0131	11/24/03	
pbw-L11W1	Ni2316	0.005	U	mg/L	0.005			-0.00034	11/12/03	11:42
pbw-L11W1	P 1782	0.04	U	mg/L	0.04			0.01860	11/12/03	11:42
pbw-L11W1	Pb220	0.005	U	mg/L	0.005			-0.00158	11/12/03	11:42
pbw-L11W1	Pd3404	0.01	U	mg/L	0.01			0.00226	11/12/03	11:42
pbw-L11W1	S 1820	0.025	U	mg/L	0.025			0.00864	11/12/03	11:42
pbw-L11W1	Sb2068	0.01	U	mg/L	0.01			0.00358	11/12/03	11:42
pbw-L11W1	Se196	0.015	U	mg/L	0.015			0.00139	11/12/03	11:42
pbw-L11W1	Si2881	0.05	U	mg/L	0.05			0.01859	11/12/03	11:42
pbw-L11W1	Sn1899	0.01	U	mg/L	0.01			0.00268	11/12/03	11:42
pbw-L11W1	Sr4215	0.005	U	mg/L	0.005			0.00021	11/12/03	11:42

Div 20
to#031106-12
06002.01.141

410067

system id	elem	result	qual	units	rl	tv	rec	ug/ml	date	time
pbw-L11W1	Th2837	0.01	U	mg/L	0.01			-0.00142	11/12/03	11:42
pbw-L11W1	Ti3372	0.005	U	mg/L	0.005			0.00086	11/12/03	11:42
pbw-L11W1	Ti1908	0.025	U	mg/L	0.025			0.00399	11/12/03	11:42
pbw-L11W1	U 3859	0.1	U	mg/L	0.1			0.04330	11/12/03	11:42
pbw-L11W1	V 2924	0.005	U	mg/L	0.005			-0.00032	11/12/03	11:42
pbw-L11W1	W 2079	0.015	U	mg/L	0.015			0.00468	11/12/03	11:42
pbw-L11W1	Y 3710	0.005	U	mg/L	0.005			-0.00002	11/12/03	11:42
pbw-L11W1	Zn2062	0.005	U	mg/L	0.005			0.00193	11/12/03	11:42
pbw-L11W1	Zr3496	0.005	U	mg/L	0.005			0.00021	11/12/03	11:42



I C P ANALYSIS

PROJ. NO.	PROJECT	TO#	DATE	MATRIX	LOGBK PG
06002.01.141	Div 20	031106-12	11-12-03	Water	S4294
_____	_____	_____	_____	_____	_____
_____	_____	_____	_____	_____	_____
_____	_____	_____	_____	_____	_____
_____	_____	_____	_____	_____	_____
_____	_____	_____	_____	_____	_____
_____	_____	_____	_____	_____	_____

INSTRUMENT: Trace 1 FILENAME: A3110612

INSTRUMENT DL: _____

010008

Analysis Report

Averages

11/12/03 02:24:14 PM

page 1

Method: DAILY1 Sample Name: blk Operator:
 Run Time: 11/12/03 11:03 Filename: DAILY32
 Mode: IR Type: X Corr. Factor: 1.00000
 Lab ID.: Cust. Smpl. ID.: Cust. ID.:

Elems	Ag3280	Al3082	As1890	B_2496	Ba4934	Be3130
Units	count	count	count	count	count	count
Avg	-.00001	.00054	-.00005	.00007	-.00003	-.00008
SDev	.00003	.00002	.00002	.00003	0	0
%RSD	246.223	5.10779	33.5174	53.1128	14.7535	1.99075

Elems	Bi2230	Ca3179	Cd2265	Co2286	Cr2677	Cu3247
Units	count	count	count	count	count	count
Avg	-.00005	.00001	0	-.00002	.00003	.00042
SDev	.00001	0	.00002	0	.00001	0
%RSD	17.3184	33.0883	497.925	30.5879	34.2539	1.29706

Elems	Fe2714	K_7664	La4086	Li6707	Mg2790	Mn2576
Units	count	count	count	count	count	count
Avg	-.00003	.00012	.00043	.00002	-.00002	0
SDev	0	.00222	.00005	.00007	.00002	0
%RSD	4.06668	1722.11	13.8519	294.31	107.74	64.7716

Elems	Mo2020	Na5889	Na3302	Ni2316	P_1782	2203/1
Units	count	count	count	count	count	count
Avg	-.00001	-.00834	-.00028	0	.00003	-.00004
SDev	.00001	.00164	.00012	.00001	.00001	.00015
%RSD	112.914	19.6707	44.1359	488.129	48.9545	340.458

Elems	2203/2	Pd3404	S_1820	Sb2068	Sc3613	1960/1
Units	count	count	count	count	count	count
Avg	.00011	-.00005	0	0	61.0989	-.00016
SDev	.00012	.00002	.00001	.00006	.37731	.00009
%RSD	107.093	43.4153	602.047	812.529	.61754	58.1811

Elems	1960/2	Si2881	Sn1899	Sr4215	Th2837	Ti3372
Units	count	count	count	count	count	count
Avg	.00017	.00152	-.00001	0	.00001	-.00142
SDev	.00017	.00006	.00002	.00001	.00001	.00006
%RSD	101.696	4.1703	129.537	1324.24	124.646	4.51407

Elems	Tl1908	U_3859	V_2924	W_2079	Y_3710	Zn2062
Units	count	count	count	count	count	count
Avg	-.00005	.00383	0	.00005	0	0
SDev	.00001	.00007	0	.00007	.00001	0
%RSD	32.5196	1.8867	85.2472	132.906	308.777	85.2472

Elems	Zr3496	*Sc	*TIME
Units	count		
Avg	.00012	610950	10000
SDev	.00001	3774.54	0
%RSD	11.5619	.61781	0



11-12-03

12/15/03

Analysis Report

Averages

11/12/03 02:24:14 PM

page 2

Method: DAILY1 Sample Name: clp_std4 Operator:
Run Time: 11/12/03 11:08 Filename: DAILY32
Mode: IR Type: X Corr. Factor: 1.00000
Lab ID.: Cust. Smpl. ID.: Cust. ID.:

Elems	Ag3280	As1890	2203/1	2203/2	Sb2068	1960/1
Units	count	count	count	count	count	count
Avge	.05223	.09466	.37759	.38988	.10509	.20881
SDev	.00028	.00005	.00176	.00262	.00064	.00068
%RSD	.54034	.05799	.46751	.67257	.61124	.32904

Elems	1960/2	Tl1908	*Sc
Units	count	count	
Avge	.35814	.02603	598604
SDev	.00456	.00014	5673.12
%RSD	1.27563	.54955	.94772

Analysis Report

Averages

11/12/03 02:24:14 PM

page 3

Method: DAILY1 Sample Name: clp_std1 Operator:
Run Time: 11/12/03 11:12 Filename: DAILY32
Mode: IR Type: X Corr. Factor: 1.00000
Lab ID.: Cust. Smpl. ID.: Cust. ID.:

Elms	Al3082	Ca3179	Fe2714	K_7664	Li6707	Mg2790
Units	count	count	count	count	count	count
Avge	.12804	.13188	.09492	1.71313	3.37792	.0516
SDev	.00107	.003	.00172	.01911	.02184	.00047
%RSD	.836	2.27879	1.81885	1.11561	.64675	.92183

Elms	Na3302	*Sc
Units	count	
Avge	.01212	591915
SDev	.00016	4258.2
%RSD	1.36953	.71939

Analysis Report

Averages

11/12/03 02:24:14 PM

Method: DAILY1 Sample Name: clp_std5 Operator:
Run Time: 11/12/03 11:15 Filename: DAILY32
Mode: IR Type: X Corr. Factor: 1.00000
Lab ID.: Cust. Smpl. ID.: Cust. ID.:

Elems	B_2496	Bi2230	Mo2020	P_1782	Si2881	Sn1899
Units	count	count	count	count	count	count
Avge	.36537	.05776	.13862	.02937	.12354	.16023
SDev	.00299	.00043	.00173	.00024	.00062	.00031
%RSD	.82074	.75345	1.24979	.82635	.50543	.19755

Elems	Sr4215	Ti3372	*Sc
Units	count	count	
Avge	1.58791	1.0464	608022
SDev	.00334	.00253	1980.61
%RSD	.21072	.24183	.32574

Method: DAILY1 Sample Name: clp_std2 Operator:
Run Time: 11/12/03 11:19 Filename: DAILY32
Mode: IR Type: X Corr. Factor: 1.00000
Lab ID.: Cust. Smpl. ID.: Cust. ID.:

Elms	Ba4934	Be3130	Cr2677	Cu3247	Ni2316	*Sc
Units	count	count	count	count	count	
Avge	.66296	.18354	.26688	.32022	.67386	580391
SDev	.00348	.00332	.00278	.00197	.00648	23689.5
%RSD	.5255	1.81218	1.04246	.61674	.96189	4.08164

Method: DAILY1 Sample Name: clp_std3 Operator:
 Run Time: 11/12/03 11:22 Filename: DAILY32
 Mode: IR Type: X Corr. Factor: 1.00000
 Lab ID.: Cust. Smpl. ID.: Cust. ID.:

Elements	Cd2265	Co2286	Mn2576	V_2924	Zn2062	*Sc
Units	count	count	count	count	count	
Avge	1.00697	.24642	.37331	.07038	.07214	594958
SDev	.00259	.00017	.001	.00002	.00017	6829.94
%RSD	.25792	.06953	.26864	.03238	.23738	1.14797

Analysis Report

Averages

11/12/03 02:24:14 PM

page 7

Method: DAILY1 Sample Name: clp_std6 Operator:
Run Time: 11/12/03 11:26 Filename: DAILY32
Mode: IR Type: X Corr. Factor: 1.00000
Lab ID.: Cust. Smpl. ID.: Cust. ID.:

Elms	La4086	Na5889	Pd3404	S_1820	Th2837	U_3859
Units	count	count	count	count	count	count
Avge	.39382	.12868	.15273	.03741	.05715	.03296
SDev	.0001	.00008	.00004	.0001	.00033	.00005
%RSD	.02664	.06931	.02785	.27022	.59116	.17054

Elms	W_2079	Y_3710	Zr3496	*Sc
Units	count	count	count	
Avge	.10741	1.01665	1.17921	611792
SDev	.00112	.00008	.00042	6528.72
%RSD	1.04502	.00793	.03584	1.06714

Method: DAILY1 Sample Name: icv/ccv Operator:
 Run Time: 11/12/03 11:31 Filename: DAILY32
 Mode: CONC Type: Q Corr. Factor: 1.00000
 Lab ID.: Cust. Smpl. ID.: Cust. ID.:

Elms	Ag3280	Al3082	As1890	B_2496	Ba4934	Be3130
Units	ppm	ppm	ppm	ppm	ppm	ppm
Avg	1.0182	9.887	5.132	4.859	10.07	1.086
SDev	.00445	.0151	.0208	.02	.0207	.0025
%RSD	.43771	.1528	.4047	.4114	.2056	.2322

Elms	Bi2230	Ca3179	Cd2265	Co2286	Cr2677	Cu3247
Units	ppm	ppm	ppm	ppm	ppm	ppm
Avg	5.111	20.35	1.034	5.081	2.013	2.053
SDev	.0166	.0781	.0011	.005	.0064	.0007
%RSD	.3253	.3839	.1067	.0992	.3204	.0336

Elms	Fe2714	K_7664	La4086	Li6707	Mg2790	Mn2576
Units	ppm	ppm	ppm	ppm	ppm	ppm
Avg	10.24	18.51	5.014	4.790	20.13	1.036
SDev	.011	.1362	.0028	.037	.0593	.0039
%RSD	.1071	.7357	.0561	.7726	.2944	.3768

Elms	Mo2020	Na5889	Na3302	Ni2316	P_1782	2203/1
Units	ppm	ppm	ppm	ppm	ppm	ppm
Avg	5.135	046.32	27.19	5.015	4.946	5.242
SDev	.0563	.3808	.0577	.0049	.1098	.0451
%RSD	1.097	.8222	.2122	.0986	2.22	.8599

Elms	2203/2	Pd3404	S_1820	Sb2068	Sc3613	1960/1
Units	ppm	ppm	ppm	ppm	%R	ppm
Avg	5.272	1.045	1.048	1.039	99.64	5.513
SDev	.0341	.0111	.0237	.0018	.8645	.0589
%RSD	.6475	1.064	2.262	.1728	.8676	1.068

Elms	1960/2	Si2881	Pb220	Se196	Sn1899	Sr4215
Units	ppm	ppm	ppm	ppm	ppm	ppm
Avg	5.427	5.003	5.257	5.450	5.066	5.223
SDev	.046	.0202	.0077	.011	.0055	.0038
%RSD	.8469	.4032	.147	.2018	.108	.0724

Elms	Th2837	Ti3372	Tl1908	U_3859	V_2924	W_2079
Units	ppm	ppm	ppm	ppm	ppm	ppm
Avg	.9637	4.906	05.510	1.048	5.092	01.109
SDev	.0164	.0167	.055	.0141	.0001	.0053
%RSD	1.697	.341	.9986	1.344	.0014	.4739

Elms	Y_3710	Zn2062	Zr3496	*Sc	*TIME
Units	ppm	ppm	ppm		
Avg	5.057	1.055	4.865	608723	10000
SDev	.0091	.0016	.0006	5235.42	0
%RSD	.1808	.1497	.0132	.86006	0

Method: DAILY1 Sample Name: icb/ccb Operator:
 Run Time: 11/12/03 11:37 Filename: DAILY32
 Mode: CONC Type: B Corr. Factor: 1.00000
 Lab ID.: Cust. Smpl. ID.: Cust. ID.:

Elems	Ag3280	Al3082	As1890	B_2496	Ba4934	Be3130
Units	ppm	ppm	ppm	ppm	ppm	ppm
Avg	.00010	.0022	.0030	.0178	.0033	.0004
SDev	.00026	.0086	.0012	.0033	.0025	.0002
%RSD	251.39	397.7	38.22	18.57	75.99	58.3

Elems	Bi2230	Ca3179	Cd2265	Co2286	Cr2677	Cu3247
Units	ppm	ppm	ppm	ppm	ppm	ppm
Avg	H.0052	.0093	.0004	.0023	-.0006	.0005
SDev	.0015	.0053	.0003	.001	.0012	.0005
%RSD	28.04	57.47	57.98	45.08	180.6	91.89

Elems	Fe2714	K_7664	La4086	Li6707	Mg2790	Mn2576
Units	ppm	ppm	ppm	ppm	ppm	ppm
Avg	-.0060	.0437	.0021	.0011	.0097	.0001
SDev	.01	.0088	.0011	.0007	.005	.0003
%RSD	167.9	20.11	54.01	64.21	51.2	220.8

Elems	Mo2020	Na5889	Na3302	Ni2316	P_1782	2203/1
Units	ppm	ppm	ppm	ppm	ppm	ppm
Avg	H.0092	.0133	H.3280	.0012	H.0186	.0078
SDev	.0027	.0064	.1192	.0016	.0165	.0008
%RSD	29.77	47.79	36.35	136.3	88.4	9.853

Elems	2203/2	Pd3404	S_1820	Sb2068	Sc3613	1960/1
Units	ppm	ppm	ppm	ppm	%R	ppm
Avg	-.0020	.0029	.0077	.0026	102.0	.0068
SDev	.0002	.0001	.0002	.0007	.5923	.0021
%RSD	11.25	5.055	3.228	28.19	.5805	30.23

Elems	1960/2	Si2881	Pb220	Se196	Sn1899	Sr4215
Units	ppm	ppm	ppm	ppm	ppm	ppm
Avg	.0028	.0008	.0013	.0042	.0018	.0018
SDev	.0049	.0035	.0004	.0026	.0042	.0013
%RSD	171.8	460.8	31.91	61.21	236.1	72

Elems	Th2837	Ti3372	Tl1908	U_3859	V_2924	W_2079
Units	ppm	ppm	ppm	ppm	ppm	ppm
Avg	-.0013	.0019	.0044	.0471	.0015	.0041
SDev	.0022	.0009	.0014	.0031	.0015	.0037
%RSD	160.4	48.63	31.87	6.496	100.5	89.33

Elems	Y_3710	Zn2062	Zr3496	*Sc	*TIME
Units	ppm	ppm	ppm		
Avg	.0016	.0005	.0019	623254	10000
SDev	.0013	.0001	.0015	3585.74	0
%RSD	83.63	33.73	78.44	.57532	0

Analysis Report

Averages

11/12/03 02:24:14 PM

page 10

Method: DAILY1 Sample Name: pbw-L11W1 pg54-294 Operator:
 Run Time: 11/12/03 11:42 Filename: DAILY32
 Mode: CONC Type: S Corr. Factor: 1.00000
 Lab ID.: Cust. Smpl. ID.: Cust. ID.:

Elems	Ag3280	Al3082	As1890	B_2496	Ba4934	Be3130
Units	ppm	ppm	ppm	ppm	ppm	ppm
Avg	.00070	.0188	.0015	.0093	.0001	-.0001
SDev	.00011	.0071	.006	.0007	.0001	.0001
%RSD	16.304	37.67	393.9	7.358	53.33	211.4

Elems	Bi2230	Ca3179	Cd2265	Co2286	Cr2677	Cu3247
Units	ppm	ppm	ppm	ppm	ppm	ppm
Avg	H.0056	.0376	.0002	.0031	-.0016	.0002
SDev	.0021	.0002	0	.0001	.0001	.0003
%RSD	37.26	.442	23.56	2.515	6.07	138.9

Elems	Fe2714	K_7664	La4086	Li6707	Mg2790	Mn2576
Units	ppm	ppm	ppm	ppm	ppm	ppm
Avg	-.0050	.0432	.0010	.0002	.0132	.0001
SDev	.0014	.0035	.0008	0	.0011	0
%RSD	28.38	8.101	77.94	8.044	8.427	30.76

Elems	Mo2020	Na5889	Na3302	Ni2316	P_1782	2203/1
Units	ppm	ppm	ppm	ppm	ppm	ppm
Avg	.0029	.0273	H.3231	-.0003	H.0186	.0031
SDev	.0006	.0005	.1124	.0003	.0036	.0021
%RSD	21.1	2.001	34.79	91.11	19.17	68.81

Elems	2203/2	Pd3404	S_1820	Sb2068	Sc3613	1960/1
Units	ppm	ppm	ppm	ppm	%R	ppm
Avg	-.0039	.0023	.0086	.0036	99.31	.0052
SDev	.0001	.0002	.0076	.0002	1.349	.0049
%RSD	3.304	7.945	87.87	4.971	1.358	93.63

Elems	1960/2	Si2881	Pb220	Se196	Sn1899	Sr4215
Units	ppm	ppm	ppm	ppm	ppm	ppm
Avg	-.0005	H.0186	-.0016	.0014	.0027	.0002
SDev	.0031	.0009	.0006	.0037	.0025	0
%RSD	584.9	4.608	38.67	266.4	93.91	1.812

Elems	Th2837	Ti3372	Tl1908	U_3859	V_2924	W_2079
Units	ppm	ppm	ppm	ppm	ppm	ppm
Avg	-.0014	.0009	.0040	.0433	-.0003	.0047
SDev	.0018	.0005	.0007	.0167	0	.0017
%RSD	127.4	52.69	16.63	38.52	3.635	36.26

Elems	Y_3710	Zn2062	Zr3496	*Sc	*TIME
Units	ppm	ppm	ppm		
Avg	-.0000	.0019	.0002	606746	10000
SDev	0	.0003	.0002	8264.66	0
%RSD	188	17.65	115.1	1.36212	0

Analysis Report

Averages

11/12/03 02:24:14 PM

page 11

Method: DAILY1 Sample Name: lcsw-L11W1 Operator:
 Run Time: 11/12/03 11:46 Filename: DAILY32
 Mode: CONC Type: S Corr. Factor: 1.00000
 Lab ID.: Cust. Smpl. ID.: Cust. ID.:

Elms	Ag3280	Al3082	As1890	B_2496	Ba4934	Be3130
Units	ppm	ppm	ppm	ppm	ppm	ppm
Avg	.05221	1.919	2.082	.0060	2.024	.0542
SDev	.00049	.0009	.0055	.0011	.0056	.0007
%RSD	.94608	.0473	.2638	17.99	.2776	1.274

Elms	Bi2230	Ca3179	Cd2265	Co2286	Cr2677	Cu3247
Units	ppm	ppm	ppm	ppm	ppm	ppm
Avg	.0069	20.79	.0524	.5151	.2034	.2509
SDev	.0045	.0752	0	.0012	.0004	.0003
%RSD	65.96	.3615	.0058	.2344	.1805	.1154

Elms	Fe2714	K_7664	La4086	Li6707	Mg2790	Mn2576
Units	ppm	ppm	ppm	ppm	ppm	ppm
Avg	1.031	16.37	.0014	.0003	20.38	.5156
SDev	.0032	.1185	.0004	0	.0158	.0013
%RSD	.3067	.7234	26.67	.3516	.0774	.2478

Elms	Mo2020	Na5889	Na3302	Ni2316	P_1782	2203/1
Units	ppm	ppm	ppm	ppm	ppm	ppm
Avg	.0014	28.79	17.38	.5110	.0178	.5318
SDev	0	.1332	.1794	.0003	.0149	.0126
%RSD	.0554	.4628	1.032	.0531	83.45	2.366

Elms	2203/2	Pd3404	S_1820	Sb2068	Sc3613	1960/1
Units	ppm	ppm	ppm	ppm	%R	ppm
Avg	.5255	.0071	.0204	.5145	98.67	2.132
SDev	.0049	.0013	.0029	.0037	.2754	.0208
%RSD	.9383	18.12	14.17	.7139	.2792	.976

Elms	1960/2	Si2881	Pb220	Se196	Sn1899	Sr4215
Units	ppm	ppm	ppm	ppm	ppm	ppm
Avg	2.156	.0160	.5270	2.146	.0029	.0003
SDev	.0106	.0003	.0075	.014	.0007	0
%RSD	.491	1.852	1.418	.6516	22.23	7.538

Elms	Th2837	Ti3372	Tl1908	U_3859	V_2924	W_2079
Units	ppm	ppm	ppm	ppm	ppm	ppm
Avg	-.0065	.0007	2.219	.0789	.5101	.0018
SDev	.0045	0	.0302	.0077	.0001	.0023
%RSD	69.37	1.02	1.36	9.789	.0209	125.7

Elms	Y_3710	Zn2062	Zr3496	*Sc	*TIME
Units	ppm	ppm	ppm		
Avg	.0001	.5308	.0014	602832	10000
SDev	0	.0004	.0002	1710.49	0
%RSD	29.39	.0683	17.15	.28374	0

Method: DAILY1 Sample Name: 237578 Operator:
 Run Time: 11/12/03 11:51 Filename: DAILY32
 Mode: CONC Type: S Corr. Factor: 1.00000
 Lab ID.: Cust. Smpl. ID.: Cust. ID.:

Elems	Ag3280	Al3082	As1890	B_2496	Ba4934	Be3130
Units	ppm	ppm	ppm	ppm	ppm	ppm
Avg	.00051	.6198	.0145	.2037	.0216	.0004
SDev	.00012	.0024	.0027	.0003	.0017	.0001
%RSD	23.876	.3792	18.73	.1581	7.836	18.82

Elems	Bi2230	Ca3179	Cd2265	Co2286	Cr2677	Cu3247
Units	ppm	ppm	ppm	ppm	ppm	ppm
Avg	.0052	2.167	.0002	.0050	.0042	.0027
SDev	.0002	.0148	0	.0003	.0003	.0004
%RSD	4.107	.6838	10.27	5.464	7.248	16.35

Elems	Fe2714	K_7664	La4086	Li6707	Mg2790	Mn2576
Units	ppm	ppm	ppm	ppm	ppm	ppm
Avg	1.143	2.341	.0049	.0815	.3881	.0245
SDev	.0098	.0125	.0005	.0001	.0115	.0006
%RSD	.8614	.5333	10.95	.1787	2.957	2.402

Elems	Mo2020	Na5889	Na3302	Ni2316	P_1782	2203/1
Units	ppm	ppm	ppm	ppm	ppm	ppm
Avg	.0626	168.2	107.8	.0015	.0432	.0073
SDev	0	.6254	.297	.0005	.0006	.0015
%RSD	.0346	.3717	.2756	33.62	1.486	20.3

Elems	2203/2	Pd3404	S_1820	Sb2068	Sc3613	1960/1
Units	ppm	ppm	ppm	ppm	%R	ppm
Avg	-.0010	.0041	20.17	.0032	97.49	.0117
SDev	.0015	.001	.188	.0016	.2953	.0014
%RSD	156.3	23.29	.9322	48.63	.3029	12.2

Elems	1960/2	Si2881	Pb220	Se196	Sn1899	Sr4215
Units	ppm	ppm	ppm	ppm	ppm	ppm
Avg	.0052	21.06	.0018	.0074	.0023	.0277
SDev	.0004	.0324	.0005	.0007	.0014	.0001
%RSD	7.67	.1537	29.43	10.05	61.22	.392

Elems	Th2837	Ti3372	Tl1908	U_3859	V_2924	W_2079
Units	ppm	ppm	ppm	ppm	ppm	ppm
Avg	.0007	.0041	.0013	.0783	.0028	-.0003
SDev	.0011	.0003	.0182	.0193	.0003	.0004
%RSD	159.2	7.542	1418	24.68	12.33	164.7

Elems	Y_3710	Zn2062	Zr3496	*Sc	*TIME
Units	ppm	ppm	ppm		
Avg	.0007	.0147	.0007	595628	10000
SDev	.0001	.0007	.0004	1833.53	0
%RSD	12.09	4.773	54.81	.30783	0

Method: DAILY1 Sample Name: 237578d Operator:
 Run Time: 11/12/03 11:56 Filename: DAILY32
 Mode: CONC Type: S Corr. Factor: 1.00000
 Lab ID.: Cust. Smpl. ID.: Cust. ID.:

Elms	Ag3280	Al3082	As1890	B_2496	Ba4934	Be3130
Units	ppm	ppm	ppm	ppm	ppm	ppm
Avg	-.00017	.6593	.0082	.2056	.0195	.0002
SDev	.00029	.0317	.001	.0009	.0003	.0001
%RSD	165.31	4.813	11.59	.446	1.601	37.02

Elms	Bi2230	Ca3179	Cd2265	Co2286	Cr2677	Cu3247
Units	ppm	ppm	ppm	ppm	ppm	ppm
Avg	.0038	2.181	.0000	.0045	.0046	.0027
SDev	.0029	.0057	0	.0004	.0021	.0018
%RSD	76.73	.2624	607.1	8.931	46.34	65.93

Elms	Fe2714	K_7664	La4086	Li6707	Mg2790	Mn2576
Units	ppm	ppm	ppm	ppm	ppm	ppm
Avg	1.227	2.278	.0051	.0811	.3732	.0245
SDev	.031	.0083	.0008	.0006	.0039	.0001
%RSD	2.528	.366	16.62	.7205	1.054	.5058

Elms	Mo2020	Na5889	Na3302	Ni2316	P_1782	2203/1
Units	ppm	ppm	ppm	ppm	ppm	ppm
Avg	.0642	167.7	108.2	.0025	.0403	.0107
SDev	.003	1.351	.7028	.0011	.0004	.0004
%RSD	4.647	.8056	.6497	46.49	.9893	3.399

Elms	2203/2	Pd3404	S_1820	Sb2068	Sc3613	1960/1
Units	ppm	ppm	ppm	ppm	%R	ppm
Avg	-.0003	.0037	19.65	.0009	98.61	.0126
SDev	.0011	.0002	.4706	.0015	.7603	.0006
%RSD	386.4	6.796	2.395	156.4	.771	5.063

Elms	1960/2	Si2881	Pb220	Se196	Sn1899	Sr4215
Units	ppm	ppm	ppm	ppm	ppm	ppm
Avg	.0001	21.39	.0034	.0043	.0036	.0278
SDev	.0018	.1564	.0009	.001	.0023	.0001
%RSD	1321	.7312	26.09	22.28	64.2	.2985

Elms	Th2837	Ti3372	Tl1908	U_3859	V_2924	W_2079
Units	ppm	ppm	ppm	ppm	ppm	ppm
Avg	-.0016	.0044	.0064	.0862	.0029	.0009
SDev	.0039	.0003	.0019	.0024	.0018	.0019
%RSD	244.6	7.88	28.82	2.784	64.47	210.8

Elms	Y_3710	Zn2062	Zr3496	*Sc	*TIME
Units	ppm	ppm	ppm		
Avg	.0010	.0147	.0016	602458	10000
SDev	.0005	.0003	.0019	4671.85	0
%RSD	50.95	1.831	119.7	.77546	0

Method: DAILY1 Sample Name: 237579 Operator:
 Run Time: 11/12/03 12:01 Filename: DAILY32
 Mode: CONC Type: S Corr. Factor: 1.00000
 Lab ID.: Cust. Smpl. ID.: Cust. ID.:

Elems	Ag3280	Al3082	As1890	B_2496	Ba4934	Be3130
Units	ppm	ppm	ppm	ppm	ppm	ppm
Avg	.00058	.3400	.0111	.2098	.0080	.0001
SDev	.00045	.0302	.0003	.0021	.0001	.0001
%RSD	77.164	8.889	3.056	.9863	1.594	88.71

Elems	Bi2230	Ca3179	Cd2265	Co2286	Cr2677	Cu3247
Units	ppm	ppm	ppm	ppm	ppm	ppm
Avg	.0067	1.846	.0001	.0039	-.0007	.0013
SDev	.0003	.006	.0001	.0005	.0002	.0002
%RSD	4.971	.323	50.07	12.29	25.81	15.19

Elems	Fe2714	K_7664	La4086	Li6707	Mg2790	Mn2576
Units	ppm	ppm	ppm	ppm	ppm	ppm
Avg	.2617	2.291	.0023	.0825	.2694	.0135
SDev	.0017	.0046	.0001	.0002	.0037	.0001
%RSD	.6545	.2003	3.295	.2619	1.386	.739

Elems	Mo2020	Na5889	Na3302	Ni2316	P_1782	2203/1
Units	ppm	ppm	ppm	ppm	ppm	ppm
Avg	.0664	171.7	110.2	.0002	.0388	.0085
SDev	.0019	.629	.1004	.0005	.0108	.0065
%RSD	2.903	.3664	.0911	231.4	27.85	76.52

Elems	2203/2	Pd3404	S_1820	Sb2068	Sc3613	1960/1
Units	ppm	ppm	ppm	ppm	%R	ppm
Avg	-.0044	.0043	19.78	.0028	95.96	.0109
SDev	.0049	.0019	.2997	.0013	1.378	.0008
%RSD	113.3	43.11	1.515	45.14	1.436	7.621

Elems	1960/2	Si2881	Pb220	Se196	Sn1899	Sr4215
Units	ppm	ppm	ppm	ppm	ppm	ppm
Avg	.0029	20.73	-.0001	.0056	.0099	.0218
SDev	.0054	.2227	.0011	.0033	.0035	0
%RSD	181.8	1.074	1644	58.88	35.57	.1282

Elems	Th2837	Ti3372	Tl1908	U_3859	V_2924	W_2079
Units	ppm	ppm	ppm	ppm	ppm	ppm
Avg	-.0030	.0025	.0013	.0852	.0019	.0026
SDev	.0016	.0001	.015	.0023	.0003	.0023
%RSD	52.46	4.901	1182	2.706	18.61	86.66

Elems	Y_3710	Zn2062	Zr3496	*Sc	*TIME
Units	ppm	ppm	ppm		
Avg	.0001	.0074	.0006	586314	10000
SDev	.0001	.0008	.0002	8406.79	0
%RSD	78.27	10.35	30.56	1.43383	0

Analysis Report

Averages

11/12/03 02:24:14 PM

Method: DAILY1 Sample Name: 237579s Operator:
 Run Time: 11/12/03 12:06 Filename: DAILY32
 Mode: CONC Type: S Corr. Factor: 1.00000
 Lab ID.: Cust. Smpl. ID.: Cust. ID.:

Elms	Ag3280	Al3082	As1890	B_2496	Ba4934	Be3130
Units	ppm	ppm	ppm	ppm	ppm	ppm
Avg	.04223	2.304	2.119	.1974	2.049	.0529
SDev	.00001	.0141	.0045	.001	.0159	.0005
%RSD	.04598	.6135	.2109	.5038	.7753	1.007

Elms	Bi2230	Ca3179	Cd2265	Co2286	Cr2677	Cu3247
Units	ppm	ppm	ppm	ppm	ppm	ppm
Avg	.0071	22.63	.0517	.5195	.1836	.2591
SDev	.0037	.1079	.0002	.0031	.0015	.0016
%RSD	51.3	.4769	.3557	.6	.8403	.6359

Elms	Fe2714	K_7664	La4086	Li6707	Mg2790	Mn2576
Units	ppm	ppm	ppm	ppm	ppm	ppm
Avg	1.234	26.88	.0020	.0827	21.19	.5234
SDev	.0117	.1006	.0003	0	.1199	.003
%RSD	.9493	.3744	17.06	.0489	.5657	.5772

Elms	Mo2020	Na5889	Na3302	Ni2316	P_1782	2203/1
Units	ppm	ppm	ppm	ppm	ppm	ppm
Avg	.0678	202.6	132.5	.5094	.0527	.5257
SDev	.0018	.2619	.0545	.0028	.0159	.0047
%RSD	2.693	.1293	.0411	.5414	30.16	.8887

Elms	2203/2	Pd3404	S_1820	Sb2068	Sc3613	1960/1
Units	ppm	ppm	ppm	ppm	%R	ppm
Avg	.5217	.0048	18.80	.3698	100.2	2.331
SDev	.0069	.0003	.2966	.0089	.2412	.0283
%RSD	1.322	5.954	1.577	2.403	.2407	1.213

Elms	1960/2	Si2881	Pb220	Se196	Sn1899	Sr4215
Units	ppm	ppm	ppm	ppm	ppm	ppm
Avg	2.318	19.73	.5225	2.320	.0035	.0209
SDev	.0151	.0401	.0061	.0194	.002	.0001
%RSD	.6494	.2035	1.177	.8381	57.02	.4694

Elms	Th2837	Ti3372	Tl1908	U_3859	V_2924	W_2079
Units	ppm	ppm	ppm	ppm	ppm	ppm
Avg	-.0059	.0020	2.283	.0560	.5143	-.0018
SDev	.0018	.0005	.0166	.0025	.0019	.0011
%RSD	31.28	23.78	.7282	4.543	.3659	60.43

Elms	Y_3710	Zn2062	Zr3496	*Sc	*TIME
Units	ppm	ppm	ppm		
Avg	.0001	.5365	.0011	612252	10000
SDev	0	.0067	.0001	1495.53	0
%RSD	29.47	1.248	13.52	.24426	0

Method: DAILY1 Sample Name: 237580 Operator:
 Run Time: 11/12/03 12:10 Filename: DAILY32
 Mode: CONC Type: S Corr. Factor: 1.00000
 Lab ID.: Cust. Smpl. ID.: Cust. ID.:

Elems	Ag3280	Al3082	As1890	B_2496	Ba4934	Be3130
Units	ppm	ppm	ppm	ppm	ppm	ppm
Avg	.00031	.3537	.0191	.2735	.0385	.0000
SDev	.00019	.0129	.0005	.003	.0026	.0001
%RSD	60.799	3.642	2.699	1.109	6.742	768

Elems	Bi2230	Ca3179	Cd2265	Co2286	Cr2677	Cu3247
Units	ppm	ppm	ppm	ppm	ppm	ppm
Avg	.0064	6.526	.0002	.0098	.0015	.0094
SDev	.0016	.0322	.0001	.0012	.0001	.0005
%RSD	25.54	.4942	72.93	12.47	4.155	5.656

Elems	Fe2714	K_7664	La4086	Li6707	Mg2790	Mn2576
Units	ppm	ppm	ppm	ppm	ppm	ppm
Avg	1.066	3.190	.0023	.1117	.9014	.0541
SDev	.0072	.0133	.0004	.001	.026	.0004
%RSD	.671	.4169	17.78	.9363	2.887	.8199

Elems	Mo2020	Na5889	Na3302	Ni2316	P_1782	2203/1
Units	ppm	ppm	ppm	ppm	ppm	ppm
Avg	.0627	228.7	147.9	.0046	.0611	.0146
SDev	.002	2.307	.3329	.001	.0076	.0041
%RSD	3.197	1.009	.2252	22.56	12.43	28.02

Elems	2203/2	Pd3404	S_1820	Sb2068	Sc3613	1960/1
Units	ppm	ppm	ppm	ppm	%R	ppm
Avg	.0009	.0052	27.66	.0038	74.09	.0211
SDev	.0012	.0001	.4256	.0046	.1556	.0033
%RSD	132	1.169	1.539	121.7	.2101	15.63

Elems	1960/2	Si2881	Pb220	Se196	Sn1899	Sr4215
Units	ppm	ppm	ppm	ppm	ppm	ppm
Avg	.0091	27.10	.0055	.0131	.0105	.0777
SDev	.0119	.0631	.0006	.0068	.0029	.0006
%RSD	130.8	.2328	10.83	52.1	27.78	.7234

Elems	Th2837	Ti3372	Tl1908	U_3859	V_2924	W_2079
Units	ppm	ppm	ppm	ppm	ppm	ppm
Avg	-.0040	.0011	-.0077	.0755	.0039	.0013
SDev	.0032	.0008	.0041	.0087	.0013	.0035
%RSD	78.06	72.02	53.33	11.54	33.65	267.4

Elems	Y_3710	Zn2062	Zr3496	*Sc	*TIME
Units	ppm	ppm	ppm		
Avg	.0001	.0450	.0006	452996	100000
SDev	.0001	.0007	.0007	912.168	0
%RSD	100.4	1.638	112.1	.20136	0

Method: DAILY1 Sample Name: 237581 Operator:
 Run Time: 11/12/03 12:15 Filename: DAILY32
 Mode: CONC Type: S Corr. Factor: 1.00000
 Lab ID.: Cust. Smpl. ID.: Cust. ID.:

Elems	Ag3280	Al3082	As1890	B_2496	Ba4934	Be3130
Units	ppm	ppm	ppm	ppm	ppm	ppm
Avg	-.00060	.0310	.0068	.2322	.0035	-.0000
SDev	.00115	.0043	.0034	.0003	.0001	0
%RSD	191.81	13.99	49	.1145	2.801	66.76

Elems	Bi2230	Ca3179	Cd2265	Co2286	Cr2677	Cu3247
Units	ppm	ppm	ppm	ppm	ppm	ppm
Avg	.0040	2.229	.0001	.0035	-.0010	.0009
SDev	.0051	.0429	.0005	.0012	.0001	.0001
%RSD	129.6	1.923	1032	33.38	15.49	7.006

Elems	Fe2714	K_7664	La4086	Li6707	Mg2790	Mn2576
Units	ppm	ppm	ppm	ppm	ppm	ppm
Avg	.2755	4.071	.0014	.2491	.1026	.0401
SDev	.0038	.0046	.0011	.001	.0103	.0003
%RSD	1.373	.1131	77.91	.4032	10.03	.6428

Elems	Mo2020	Na5889	Na3302	Ni2316	P_1782	2203/1
Units	ppm	ppm	ppm	ppm	ppm	ppm
Avg	.0395	586.18	190.4	-.0001	.0255	-.0017
SDev	.0011	.3641	1.396	.0006	.0006	.0076
%RSD	2.791	.4225	.7329	839	2.508	446.2

Elems	2203/2	Pd3404	S_1820	Sb2068	Sc3613	1960/1
Units	ppm	ppm	ppm	ppm	%R	ppm
Avg	.0008	.0000	9.967	-.0013	102.9	.0086
SDev	.0052	.0042	.1976	.0058	.7184	.0127
%RSD	609.7	16140	1.983	445	.6981	148.6

Elems	1960/2	Si2881	Pb220	Se196	Sn1899	Sr4215
Units	ppm	ppm	ppm	ppm	ppm	ppm
Avg	.0075	7.098	.0000	.0079	.0104	.0064
SDev	.0012	.0591	.0009	.0034	.0046	.0001
%RSD	16.03	.8326	18690	43.9	44.56	1.477

Elems	Th2837	Ti3372	Tl1908	U_3859	V_2924	W_2079
Units	ppm	ppm	ppm	ppm	ppm	ppm
Avg	-.0019	-.0003	-.0209	.0490	.0001	-.0031
SDev	.0013	.0014	.0154	.0288	.0002	.004
%RSD	67.99	484.7	73.58	58.73	150.9	128.3

Elems	Y_3710	Zn2062	Zr3496	*Sc	*TIME
Units	ppm	ppm	ppm		
Avg	.0001	.0023	.0009	628706	10000
SDev	0	.0015	.0003	4357.19	0
%RSD	16.91	67.47	31.16	.69304	0

Method: DAILY1 Sample Name: 237582 Operator:
 Run Time: 11/12/03 12:20 Filename: DAILY32
 Mode: CONC Type: S Corr. Factor: 1.00000
 Lab ID.: Cust. Smpl. ID.: Cust. ID.:

Elems	Ag3280	Al3082	As1890	B_2496	Ba4934	Be3130
Units	ppm	ppm	ppm	ppm	ppm	ppm
Avg	-.00001	.0282	.0077	.2332	.0031	-.0000
SDev	.00044	.0002	.0028	.0001	.0001	0
%RSD	2747.4	.718	36.36	.0629	4.054	72.56

Elems	Bi2230	Ca3179	Cd2265	Co2286	Cr2677	Cu3247
Units	ppm	ppm	ppm	ppm	ppm	ppm
Avg	.0061	2.242	.0002	.0037	-.0012	-.0007
SDev	.0002	.0209	.0001	.0006	.0007	.0007
%RSD	3.556	.9306	24.33	16.57	53.76	109.5

Elems	Fe2714	K_7664	La4086	Li6707	Mg2790	Mn2576
Units	ppm	ppm	ppm	ppm	ppm	ppm
Avg	.2025	4.073	.0022	.2466	.1060	.0388
SDev	.0151	.0187	.0004	.0002	.001	.0003
%RSD	7.481	.46	17.68	.0949	.9072	.8092

Elems	Mo2020	Na5889	Na3302	Ni2316	P_1782	2203/1
Units	ppm	ppm	ppm	ppm	ppm	ppm
Avg	.0384	585.78	189.2	-.0001	.0390	.0061
SDev	.002	1.076	.3177	.0011	.0033	.0002
%RSD	5.163	1.255	.168	865.5	8.549	3.742

Elems	2203/2	Pd3404	S_1820	Sb2068	Sc3613	1960/1
Units	ppm	ppm	ppm	ppm	%R	ppm
Avg	-.0053	.0046	9.887	-.0002	103.4	.0248
SDev	.0019	.0029	.0524	.0029	.9925	.0032
%RSD	35.01	61.59	.5297	1374	.96	12.93

Elems	1960/2	Si2881	Pb220	Se196	Sn1899	Sr4215
Units	ppm	ppm	ppm	ppm	ppm	ppm
Avg	-.0006	7.132	-.0015	.0079	.0254	.0064
SDev	.0002	.024	.0013	.0012	.0009	0
%RSD	38.65	.3359	87.08	15.49	3.358	.1862

Elems	Th2837	Ti3372	Tl1908	U_3859	V_2924	W_2079
Units	ppm	ppm	ppm	ppm	ppm	ppm
Avg	-.0017	.0002	-.0223	.0635	.0000	-.0009
SDev	.0038	.0001	.0206	.0105	.0009	.0014
%RSD	221.5	75.14	92.56	16.6	9450	161.5

Elems	Y_3710	Zn2062	Zr3496	*Sc	*TIME
Units	ppm	ppm	ppm		
Avg	.0001	.0012	.0007	631618	10000
SDev	.0002	.0003	.0007	6053.54	0
%RSD	386.1	25.6	90.39	.95841	0

Analysis Report

Averages

11/12/03 02:24:14 PM

page 19

Method: DAILY1 Sample Name: 237583 Operator:
 Run Time: 11/12/03 12:25 Filename: DAILY32
 Mode: CONC Type: S Corr. Factor: 1.00000
 Lab ID.: Cust. Smpl. ID.: Cust. ID.:

Elms	Ag3280	Al3082	As1890	B_2496	Ba4934	Be3130
Units	ppm	ppm	ppm	ppm	ppm	ppm
Avg	.00005	.0960	.0072	.2341	.0042	-.0000
SDev	.00032	.0005	.0001	.0037	.0001	0
%RSD	595.15	.5399	1.301	1.601	2.903	61.07

Elms	Bi2230	Ca3179	Cd2265	Co2286	Cr2677	Cu3247
Units	ppm	ppm	ppm	ppm	ppm	ppm
Avg	.0011	2.180	-.0000	.0036	-.0017	.0001
SDev	.0008	.0019	.0001	.0003	.0002	.0001
%RSD	74.05	.0858	127.3	9.557	10.24	225.5

Elms	Fe2714	K_7664	La4086	Li6707	Mg2790	Mn2576
Units	ppm	ppm	ppm	ppm	ppm	ppm
Avg	1.326	4.214	.0020	.2561	.1105	.0489
SDev	.0058	.0437	.0004	.0005	.0002	.0001
%RSD	.4368	1.037	18.32	.1995	.2198	.1998

Elms	Mo2020	Na5889	Na3302	Ni2316	P_1782	2203/1
Units	ppm	ppm	ppm	ppm	ppm	ppm
Avg	.0348	587.57	193.3	-.0006	.0462	.0044
SDev	.001	1.071	.898	.0009	.0011	.0012
%RSD	2.961	1.223	.4645	146.8	2.301	27.49

Elms	2203/2	Pd3404	S_1820	Sb2068	Sc3613	1960/1
Units	ppm	ppm	ppm	ppm	%R	ppm
Avg	-.0041	.0029	9.854	.0035	101.1	.0205
SDev	.0012	.0009	.3171	.0004	.8749	.0008
%RSD	30.1	30.31	3.218	12.59	.8651	3.872

Elms	1960/2	Si2881	Pb220	Se196	Sn1899	Sr4215
Units	ppm	ppm	ppm	ppm	ppm	ppm
Avg	.0066	7.801	-.0012	.0112	.0242	.0065
SDev	.0019	.0597	.0012	.0015	.0001	0
%RSD	28.75	.7651	97.09	13.63	.3635	.2072

Elms	Th2837	Ti3372	Tl1908	U_3859	V_2924	W_2079
Units	ppm	ppm	ppm	ppm	ppm	ppm
Avg	-.0009	.0010	-.0276	.0235	.0000	-.0030
SDev	.0033	.0004	.0005	.0049	.0002	.0043
%RSD	376.3	42.53	1.728	20.85	537.8	141.7

Elms	Y_3710	Zn2062	Zr3496	*Sc	*TIME
Units	ppm	ppm	ppm		
Avg	-.0001	.0140	.0001	617850	10000
SDev	.0001	.0002	.0004	5342.19	0
%RSD	165.1	1.428	350.3	.86464	0

Analysis Report

QC Standard

11/12/03 12:39:08 PM

page 1

Method: DAILY1 Sample Name: icv/ccv

Operator:

Run Time: 11/12/03 12:34:23

Comment:

Mode: CONC Corr. Factor: 1

Elem	Ag3280	Al3082	As1890	B_2496	Ba4934	Be3130	Bi2230
Units	ppm	ppm	ppm	ppm	ppm	ppm	ppm
Avg	1.0232	10.03	5.187	4.861	10.16	1.055	5.102
SDev	.0043	.04	.013	.011	.02	.001	.005
%RSD	.42133	.3592	.2516	.2274	.2151	.0572	.1048
#1	1.0263	10.05	5.178	4.853	10.17	1.055	5.105
#2	1.0202	10.000	5.197	4.868	10.14	1.054	5.098
Errors	QC Pass	QC Pass	QC Pass	QC Pass	QC Pass	QC Pass	QC Pass
Value	1.0000	10.00	5.000	5.000	10.00	1.000	5.000
Range	10.000	10.00	10.00	10.00	10.00	10.00	10.00
Elem	Ca3179	Cd2265	Co2286	Cr2677	Cu3247	Fe2714	K_7664
Units	ppm	ppm	ppm	ppm	ppm	ppm	ppm
Avg	20.32	1.022	5.051	2.009	2.076	10.20	18.81
SDev	.03	.002	.016	.006	.007	.06	.08
%RSD	.1470	.2266	.3077	.3003	.3379	.5614	.4178
#1	20.34	1.024	5.062	2.013	2.081	10.24	18.75
#2	20.29	1.020	5.040	2.005	2.071	10.16	18.86
Errors	QC Pass	QC Pass	QC Pass	QC Pass	QC Pass	QC Pass	QC Pass
Value	20.00	1.000	5.000	2.000	2.000	10.00	20.00
Range	10.00	10.00	10.00	10.00	10.00	10.00	10.00
Elem	La4086	Li6707	Mg2790	Mn2576	Mo2020	Na5889	Na3302
Units	ppm	ppm	ppm	ppm	ppm	ppm	ppm
Avg	5.060	4.846	20.12	1.036	5.149	Q47.02	27.25
SDev	.001	.000	.04	.002	.044	.04	.33
%RSD	.0098	.0079	.2161	.2195	.8475	.0812	1.204
#1	5.060	4.846	20.15	1.038	5.118	Q46.99	27.02
#2	5.059	4.847	20.09	1.035	5.180	Q47.05	27.48
Errors	QC Pass	QC Pass	QC Pass	QC Pass	QC Pass	QC Fail	QC Pass
Value	5.000	5.000	20.00	1.000	5.000	30.00	30.00
Range	10.00	10.00	10.00	10.00	10.00	10.00	10.00
Elem	Ni2316	P_1782	2203/1	2203/2	Pd3404	S_1820	Sb2068
Units	ppm	ppm	ppm	ppm	ppm	ppm	ppm
Avg	4.970	5.207	5.201	5.242	1.036	1.024	1.040
SDev	.009	.034	.044	.011	.004	.021	.008
%RSD	.1750	.6568	.8491	.2086	.3628	2.073	.7389
#1	4.976	5.231	5.232	5.234	1.039	1.009	1.035
#2	4.964	5.182	5.170	5.249	1.034	1.039	1.046
Errors	QC Pass	QC Pass	NOCHECK	NOCHECK	QC Pass	QC Pass	QC Pass
Value	5.000	5.000			1.000	1.000	1.000
Range	10.00	10.00			10.00	10.00	10.00
Elem	Sc3613	1960/1	1960/2	Si2881	Pb220	Se196	Sn1899

Analysis Report

QC Standard

11/12/03 12:39:08 PM

Units	%R	ppm	ppm	ppm	ppm	ppm	ppm
Avge	100.5	5.543	5.502	4.999	5.223	Q5.510	5.009
SDev	.7	.041	.035	.018	.007	.010	.000
%RSD	.7369	.7449	.6354	.3582	.1421	.1730	.0032
#1	99.99	5.572	5.477	5.012	5.228	Q5.503	5.009
#2	101.0	5.514	5.527	4.987	5.218	Q5.517	5.010
Errors	NOCHECK	NOCHECK	NOCHECK	QC Pass	QC Pass	QC Fail	QC Pass
Value				5.000	5.000	5.000	5.000
Range				10.00	10.00	10.00	10.00
Elem	Sr4215	Th2837	Ti3372	Tl1908	U_3859	V_2924	W_2079
Units	ppm	ppm	ppm	ppm	ppm	ppm	ppm
Avge	5.256	.9572	4.956	Q5.556	1.029	5.110	1.076
SDev	.018	.0085	.010	.036	.046	.019	.015
%RSD	.3384	.8867	.2063	.6455	4.515	.3655	1.428
#1	5.269	.9632	4.963	Q5.530	.9960	5.124	1.066
#2	5.244	.9512	4.948	Q5.581	1.062	5.097	1.087
Errors	QC Pass	QC Pass	QC Pass	QC Fail	QC Pass	QC Pass	QC Pass
Value	5.000	1.000	5.000	5.000	1.000	5.000	1.000
Range	10.00	10.00	10.00	10.00	10.00	10.00	10.00
Elem	Y_3710	Zn2062	Zr3496				
Units	ppm	ppm	ppm				
Avge	5.082	1.032	4.828				
SDev	.012	.005	.018				
%RSD	.2273	.4463	.3744				
#1	5.074	1.029	4.815				
#2	5.090	1.036	4.841				
Errors	QC Pass	QC Pass	QC Pass				
Value	5.000	1.000	5.000				
Range	10.00	10.00	10.00				

010089

Analysis Report

QC Standard

11/12/03 12:39:08 PM

page 3

IntStd	1	2	3	4	5	6	7
Mode	*Counts	Time	NOTUSED	NOTUSED	NOTUSED	NOTUSED	NOTUSED
Elem	Sc	--	--	--	--	--	--
Wavlen	361.384	--	--	--	--	--	--
Avge	614067	10000	--	--	--	--	--
SDev	4485.885	.00000000	--	--	--	--	--
%RSD	.7305205	.00000000	--	--	--	--	--
#1	610895	10000	--	--	--	--	--
#2	617239	10000	--	--	--	--	--

Analysis Report

Blank Sample

11/12/03 12:46:03 PM

page 1

Method: DAILY1 Sample Name: icb/ccb

Operator:

Run Time: 11/12/03 12:41:18

Comment:

Mode: CONC Corr. Factor: 1

Elem	Ag3280	Al3082	As1890	B_2496	Ba4934	Be3130	Bi2230
Units	ppm	ppm	ppm	ppm	ppm	ppm	ppm
Avge	.00110	.0231	.0045	.0156	.0022	.0001	H.0074
SDev	.00027	.0066	.0003	.0037	.0016	.0002	.0022
%RSD	24.196	28.60	6.213	23.95	71.54	347.2	29.31

#1	.00091	.0278	.0043	.0183	.0034	.0002	H.0059
#2	.00129	.0184	.0047	.0130	.0011	-.0001	H.0089

Errors	LC Pass	LC Pass	LC Pass	LC Pass	LC Pass	LC Pass	LC High
High	.00500	.0500	.0050	.0500	.0050	.0050	.0050
Low	-.00500	-.0500	-.0050	-.0500	-.0050	-.0050	-.0050

Elem	Ca3179	Cd2265	Co2286	Cr2677	Cu3247	Fe2714	K_7664
Units	ppm	ppm	ppm	ppm	ppm	ppm	ppm
Avge	.0089	.0004	.0018	-.0009	-.0001	-.0041	.0527
SDev	.0026	.0001	.0009	.0007	.0001	.0009	.0028
%RSD	29.28	30.90	49.82	77.40	110.0	20.58	5.231

#1	.0107	.0003	.0024	-.0004	-.0000	-.0035	.0546
#2	.0070	.0004	.0012	-.0014	-.0001	-.0047	.0507

Errors	LC Pass	LC Pass	LC Pass	LC Pass	LC Pass	LC Pass	LC Pass
High	.0500	.0050	.0050	.0050	.0050	.0250	.1000
Low	-.0500	-.0050	-.0050	-.0050	-.0050	-.0250	-.1000

Elem	La4086	Li6707	Mg2790	Mn2576	Mo2020	Na5889	Na3302
Units	ppm	ppm	ppm	ppm	ppm	ppm	ppm
Avge	.0029	.0008	.0104	.0001	H.0076	.0181	H.4028
SDev	.0004	.0005	.0017	.0001	.0034	.0056	.0302
%RSD	12.35	57.15	16.15	218.7	45.20	30.80	7.503

#1	.0032	.0011	.0092	.0001	H.0100	.0220	H.4241
#2	.0027	.0005	.0116	-.0000	H.0052	.0141	H.3814

Errors	LC Pass	LC Pass	LC Pass	LC Pass	LC High	LC Pass	LC High
High	.0050	.0050	.0500	.0050	.0050	.0500	.0500
Low	-.0050	-.0050	-.0500	-.0050	-.0050	-.0500	-.0500

Elem	Ni2316	P_1782	2203/1	2203/2	Pd3404	S_1820	Sb2068
Units	ppm	ppm	ppm	ppm	ppm	ppm	ppm
Avge	.0009	H.0182	.0035	-.0023	.0040	.0015	.0032
SDev	.0003	.0056	.0019	.0007	.0012	.0092	.0039
%RSD	36.90	30.67	55.26	31.52	29.65	611.6	120.4

#1	.0012	H.0222	.0048	-.0028	.0048	-.0050	.0059
#2	.0007	H.0143	.0021	-.0018	.0031	.0080	.0005

Errors	LC Pass	LC High	NOCHECK	NOCHECK	LC Pass	LC Pass	LC Pass
High	.0050	.0100			.0050	.0100	.0100
Low	-.0050	-.0100			-.0050	-.0100	-.0100

Elem	Sc3613	1960/1	1960/2	Si2881	Pb220	Se196	Sn1899

Analysis Report

Blank Sample

11/12/03 12:46:03 PM

page 2

Units	%R	ppm	ppm	ppm	ppm	ppm	ppm
Avg	100.7	.0043	-.0001	.0055	-.0004	.0014	.0026
SDev	.4	.0047	.0007	.0003	.0002	.0011	.0014
%RSD	.3658	108.4	1170.	6.220	44.05	75.33	54.74
#1	100.5	.0076	-.0006	.0057	-.0003	.0022	.0037
#2	101.0	.0010	.0005	.0053	-.0005	.0007	.0016
Errors	NOCHECK	NOCHECK	NOCHECK	LC Pass	LC Pass	LC Pass	LC Pass
High				.0100	.0030	.0050	.0050
Low				-.0100	-.0030	-.0050	-.0050
Elem	Sr4215	Th2837	Ti3372	Ti1908	U_3859	V_2924	W_2079
Units	ppm	ppm	ppm	ppm	ppm	ppm	ppm
Avg	.0012	-.0004	.0014	.0092	.0423	.0012	.0024
SDev	.0008	.0008	.0009	.0083	.0084	.0005	.0026
%RSD	69.81	186.4	60.88	90.29	19.92	42.18	108.0
#1	.0017	-.0010	.0020	H.0151	.0364	.0016	.0006
#2	.0006	.0001	.0008	.0033	.0483	.0009	.0042
Errors	LC Pass	LC Pass	LC Pass	LC Pass	LC Pass	LC Pass	LC Pass
High	.0050	.0100	.0050	.0100	.1000	.0050	.0100
Low	-.0050	-.0100	-.0050	-.0100	-.1000	-.0050	-.0100
Elem	Y_3710	Zn2062	Zr3496				
Units	ppm	ppm	ppm				
Avg	.0011	.0011	.0021				
SDev	.0008	.0002	.0009				
%RSD	71.95	14.36	44.85				
#1	.0016	.0012	.0028				
#2	.0005	.0010	.0015				
Errors	LC Pass	LC Pass	LC Pass				
High	.0050	.0050	.0050				
Low	-.0050	-.0050	-.0050				

010092

Analysis Report

Blank Sample

11/12/03 12:46:03 PM

page 3

IntStd	1	2	3	4	5	6	7
Mode	*Counts	Time	NOTUSED	NOTUSED	NOTUSED	NOTUSED	NOTUSED
Elem	Sc	--	--	--	--	--	--
Wavlen	361.384	--	--	--	--	--	--
Avge	615314	10000	--	--	--	--	--
SDev	2253.549	.0000000	--	--	--	--	--
%RSD	.3662441	.0000000	--	--	--	--	--
#1	613720	10000	--	--	--	--	--
#2	616907	10000	--	--	--	--	--

Analysis Report

11/12/03 12:50:55 PM

page 1

Method: DAILY1 Sample Name: 237584

Operator:

Run Time: 11/12/03 12:46:13

Comment:

Mode: CONC Corr. Factor: 1

Elem	Ag3280	Al3082	As1890	Br_2496	Ba4934	Be3130	Bi2230
Units	ppm	ppm	ppm	ppm	ppm	ppm	ppm
Avge	.00040	.1669	.0349	.1786	.0026	-.0002	.0083
SDev	.00002	.0059	.0012	.0009	.0000	.0000	.0001
%RSD	3.9637	3.525	3.455	.5124	.3433	24.30	.8548

#1	.00042	.1627	.0340	.1780	.0027	-.0002	.0082
#2	.00039	.1711	.0357	.1793	.0026	-.0001	.0083

Elem	Ca3179	Cd2265	Co2286	Cr2677	Cu3247	Fe2714	K_7664
Units	ppm	ppm	ppm	ppm	ppm	ppm	ppm
Avge	.5433	.0001	.0040	-.0007	.0002	.0745	3.681
SDev	.0008	.0001	.0002	.0003	.0003	.0096	.032
%RSD	.1372	90.31	4.908	37.12	146.4	12.84	.8699

#1	.5438	.0000	.0041	-.0009	.0004	.0812	3.658
#2	.5428	.0002	.0038	-.0005	-.0000	.0677	3.704

Elem	La4086	Li6707	Mg2790	Mn2576	Mo2020	Na5889	Na3302
Units	ppm	ppm	ppm	ppm	ppm	ppm	ppm
Avge	.0026	.1249	.0706	.0153	.0042	149.5	96.10
SDev	.0000	.0005	.0031	.0000	.0004	.3	.25
%RSD	.7505	.4089	4.311	.1701	10.19	.1733	.2592

#1	.0026	.1246	.0728	.0153	.0045	149.3	95.92
#2	.0026	.1253	.0685	.0153	.0039	149.6	96.27

Elem	Ni2316	P_1782	2203/1	2203/2	Pd3404	S_1820	Sb2068
Units	ppm	ppm	ppm	ppm	ppm	ppm	ppm
Avge	.0002	.0290	.0009	-.0044	.0037	5.560	.0029
SDev	.0006	.0006	.0015	.0023	.0006	.121	.0004
%RSD	372.6	1.913	161.1	51.33	17.23	2.180	14.33

#1	-.0003	.0294	-.0001	-.0060	.0042	5.475	.0026
#2	.0006	.0286	.0020	-.0028	.0033	5.646	.0032

Elem	Sc3613	1960/1	1960/2	Si2881	Pb220	Se196	Sn1899
Units	%R	ppm	ppm	ppm	ppm	ppm	ppm
Avge	102.2	.0145	.0024	27.90	-.0026	.0064	.0022
SDev	.4	.0022	.0008	.03	.0020	.0013	.0020
%RSD	.4255	15.26	33.72	.1033	77.01	19.74	89.77

#1	101.9	.0129	.0018	27.92	-.0040	.0055	.0036
#2	102.5	.0160	.0029	27.88	-.0012	.0073	.0008

Elem	Sr4215	Th2837	Ti3372	Tl1908	U_3859	V_2924	W_2079
Units	ppm	ppm	ppm	ppm	ppm	ppm	ppm
Avge	.0019	-.0013	.0012	.0114	.0746	.0051	.0012
SDev	.0000	.0031	.0004	.0104	.0093	.0002	.0012
%RSD	.7829	228.1	36.26	90.82	12.52	3.588	102.5

#1	.0019	.0008	.0015	.0188	.0812	.0052	.0021
#2	.0019	-.0035	.0009	.0041	.0680	.0049	.0003

010094

Analysis Report

11/12/03 12:50:55 PM

page 2

Elem	Y_3710	Zn2062	Zn3496
Units	ppm	ppm	ppm
Avg	.0001	.0034	.0012
SDev	.0000	.0003	.0001
%RSD	16.89	9.423	8.315

#1	.0002	.0037	.0013
#2	.0001	.0032	.0012

IntStd	1	2	3	4	5	6	7
Mode	*Counts	Time	NOTUSED	NOTUSED	NOTUSED	NOTUSED	NOTUSED
Elem	Sc	--	--	--	--	--	--
Wavlen	361.384	--	--	--	--	--	--
Avg	624375	10000	--	--	--	--	--
SDev	2623.366	.0000000	--	--	--	--	--
%RSD	.4201587	.0000000	--	--	--	--	--
#1	622520	10000	--	--	--	--	--
#2	626230	10000	--	--	--	--	--

Analysis Report

11/12/03 12:55:50 PM

Method: DAILY1 Sample Name: 237585

Operator:

Run Time: 11/12/03 12:51:06

Comment:

Mode: CONC Corr. Factor: 1

Elem	Ag3280	Al3082	As1890	B_2496	Ba4934	Be3130	Bi2230
Units	ppm	ppm	ppm	ppm	ppm	ppm	ppm
Avg	.00045	.1319	.0339	.1616	.0004	-.0002	.0059
SDev	.00009	.0058	.0017	.0023	.0002	.0000	.0018
%RSD	20.613	4.368	5.102	1.418	37.72	14.99	30.75

#1	.00039	.1359	.0351	.1632	.0003	-.0002	.0072
#2	.00052	.1278	.0326	.1600	.0006	-.0002	.0046

Elem	Ca3179	Cd2265	Co2286	Cr2677	Cu3247	Fe2714	K_7664
Units	ppm	ppm	ppm	ppm	ppm	ppm	ppm
Avg	.5380	.0003	.0043	-.0013	-.0010	.0426	3.658
SDev	.0056	.0001	.0005	.0010	.0004	.0250	.036
%RSD	1.044	43.71	11.44	79.82	45.33	58.67	.9885

#1	.5420	.0003	.0047	-.0020	-.0013	.0249	3.633
#2	.5340	.0002	.0040	-.0005	-.0006	.0603	3.684

Elem	La4086	Li6707	Mg2790	Mn2576	Mo2020	Na5889	Na3302
Units	ppm	ppm	ppm	ppm	ppm	ppm	ppm
Avg	.0024	.1245	.0679	.0139	.0029	149.0	95.88
SDev	.0001	.0001	.0012	.0002	.0015	.1	.63
%RSD	5.574	.0895	1.777	1.191	50.41	.0811	.6575

#1	.0023	.1244	.0670	.0138	.0019	148.9	95.43
#2	.0025	.1246	.0687	.0140	.0040	149.1	96.33

Elem	Ni2316	P_1782	2203/1	2203/2	Pd3404	S_1820	Sb2068
Units	ppm	ppm	ppm	ppm	ppm	ppm	ppm
Avg	-.0004	.0404	.0063	-.0075	.0043	5.481	.0023
SDev	.0006	.0081	.0005	.0009	.0011	.061	.0032
%RSD	169.8	19.98	8.317	12.73	25.52	1.108	141.4

#1	-.0008	.0347	.0067	-.0068	.0051	5.438	.0000
#2	.0001	.0461	.0060	-.0081	.0035	5.524	.0045

Elem	Sc3613	1960/1	1960/2	Si2881	Pb220	Se196	Sn1899
Units	%R	ppm	ppm	ppm	ppm	ppm	ppm
Avg	103.4	.0171	-.0002	28.16	-.0029	.0056	.0078
SDev	.1	.0041	.0051	.17	.0008	.0048	.0030
%RSD	.1230	23.80	2346.	.5989	28.27	85.92	38.35

#1	103.5	.0142	-.0039	28.28	-.0023	.0022	.0057
#2	103.3	.0200	.0034	28.04	-.0034	.0089	.0099

Elem	Sr4215	Th2837	Ti3372	Tl1908	U_3859	V_2924	W_2079
Units	ppm	ppm	ppm	ppm	ppm	ppm	ppm
Avg	.0018	-.0033	.0013	.0085	.0904	.0041	.0027
SDev	.0000	.0019	.0001	.0103	.0039	.0011	.0023
%RSD	.2711	58.00	3.820	121.5	4.311	27.44	83.18

#1	.0018	-.0046	.0013	.0158	.0876	.0033	.0044
#2	.0018	-.0019	.0013	.0012	.0931	.0049	.0011

Analysis Report

11/12/03 12:55:50 PM

page 2

Elem	Y_3710	Zn2062	Zr3496
Units	ppm	ppm	ppm
Avg	-.0001	.0022	.0004
SDev	.0002	.0001	.0012
%RSD	245.6	6.902	301.7

#1	-.0002	.0023	-.0005
#2	.0001	.0021	.0013

IntStd	1	2	3	4	5	6	7
Mode	*Counts	Time	NOTUSED	NOTUSED	NOTUSED	NOTUSED	NOTUSED
Elem	Sc	--	--	--	--	--	--
Wavlen	361.384	--	--	--	--	--	--
Avg	631880	10000	--	--	--	--	--
SDev	716.2992	.00000000	--	--	--	--	--
%RSD	.1133601	.00000000	--	--	--	--	--
#1	632386	10000	--	--	--	--	--
#2	631373	10000	--	--	--	--	--

Analysis Report

11/12/03 01:00:46 PM

page 1

Method: DAILY1 Sample Name: 237586

Operator:

Run Time: 11/12/03 12:56:01

Comment:

Mode: CONC Corr. Factor: 1

Elem	Ag3280	Al3082	As1890	B_2496	Ba4934	Be3130	Bi2230
Units	ppm	ppm	ppm	ppm	ppm	ppm	ppm
Avge	.00015	.2039	.0332	.1631	.0011	-.0001	.0061
SDev	.00097	.0051	.0033	.0006	.0001	.0001	.0003
%RSD	649.91	2.518	9.960	.3860	6.551	67.26	5.178

#1	.00083	.2076	.0356	.1635	.0011	-.0001	.0063
#2	-.00053	.2003	.0309	.1626	.0010	-.0002	.0059

Elem	Ca3179	Cd2265	Co2286	Cr2677	Cu3247	Fe2714	K_7664
Units	ppm	ppm	ppm	ppm	ppm	ppm	ppm
Avge	.7608	-.0000	.0033	-.0009	-.0003	.1946	3.807
SDev	.0034	.0001	.0003	.0004	.0000	.0110	.014
%RSD	.4429	305.4	8.587	44.32	.1328	5.642	.3619

#1	.7632	.0000	.0035	-.0006	-.0003	.2023	3.798
#2	.7584	-.0001	.0031	-.0012	-.0003	.1868	3.817

Elem	La4086	Li6707	Mg2790	Mn2576	Mo2020	Na5889	Na3302
Units	ppm	ppm	ppm	ppm	ppm	ppm	ppm
Avge	.0022	.1280	.0941	.0249	.0026	153.4	98.51
SDev	.0011	.0001	.0086	.0001	.0008	.5	.54
%RSD	49.70	.1100	9.173	.1981	29.03	.3003	.5495

#1	.0029	.1281	.1002	.0249	.0032	153.8	98.13
#2	.0014	.1279	.0880	.0248	.0021	153.1	98.90

Elem	Ni2316	P_1782	2203/1	2203/2	Pd3404	S_1820	Sb2068
Units	ppm	ppm	ppm	ppm	ppm	ppm	ppm
Avge	.0001	.0449	.0040	-.0035	.0033	5.632	-.0005
SDev	.0003	.0011	.0026	.0001	.0015	.051	.0033
%RSD	349.7	2.380	65.58	1.400	43.75	.9097	669.8

#1	-.0001	.0441	.0059	-.0036	.0044	5.595	.0018
#2	.0003	.0456	.0022	-.0035	.0023	5.668	-.0028

Elem	Sc3613	1960/1	1960/2	Si2881	Pb220	Se196	Sn1899
Units	%R	ppm	ppm	ppm	ppm	ppm	ppm
Avge	100.1	.0124	-.0003	28.74	-.0010	.0040	.0062
SDev	.2	.0013	.0024	.01	.0008	.0020	.0038
%RSD	.2043	10.46	809.6	.0343	82.27	51.08	62.05

#1	99.93	.0134	.0014	28.75	-.0004	.0054	.0089
#2	100.2	.0115	-.0020	28.74	-.0016	.0025	.0035

Elem	Sr4215	Th2837	Ti3372	Tl1908	U_3859	V_2924	W_2079
Units	ppm	ppm	ppm	ppm	ppm	ppm	ppm
Avge	.0029	-.0019	.0013	.0126	.0784	.0049	.0005
SDev	.0000	.0000	.0008	.0013	.0192	.0002	.0039
%RSD	1.180	2.187	64.58	10.35	24.44	3.696	721.5

#1	.0030	-.0019	.0019	.0136	.0919	.0050	.0033
#2	.0029	-.0020	.0007	.0117	.0648	.0048	-.0022

Elem	Y_3710	Zn2062	Zr3496
Units	ppm	ppm	ppm
Avge	.0001	.0194	.0011
SDev	.0001	.0010	.0007
%RSD	70.65	5.156	59.12

#1	.0002	.0201	.0016
#2	.0001	.0187	.0007

IntStd	1	2	3	4	5	6	7
Mode	*Counts	Time	NOTUSED	NOTUSED	NOTUSED	NOTUSED	NOTUSED
Elem	Sc	--	--	--	--	--	--
Wavlen	361.384	--	--	--	--	--	--
Avge	611354	10000	--	--	--	--	--
SDev	1260.771	.00000000	--	--	--	--	--
%RSD	.2062259	.00000000	--	--	--	--	--
#1	610463	10000	--	--	--	--	--
#2	612246	10000	--	--	--	--	--

Analysis Report

11/12/03 01:05:41 PM

page 1

Method: DAILY1 Sample Name: 237587

Operator:

Run Time: 11/12/03 13:00:57

Comment:

Mode: CONC Corr. Factor: 1

Elem	Ag3280	Al3082	As1890	B_2496	Ba4934	Be3130	Bi2230
Units	ppm	ppm	ppm	ppm	ppm	ppm	ppm
Avge	-.00031	.1580	.0373	.1568	.0003	-.0002	.0073
SDev	.00059	.0040	.0004	.0016	.0000	.0001	.0064
%RSD	188.77	2.545	1.106	1.032	12.94	55.86	88.27

#1	.00010	.1552	.0370	.1557	.0003	-.0002	.0118
#2	-.00073	.1608	.0376	.1580	.0002	-.0001	.0027

Elem	Ca3179	Cd2265	Co2286	Cr2677	Cu3247	Fe2714	K_7664
Units	ppm	ppm	ppm	ppm	ppm	ppm	ppm
Avge	.3597	-.0000	.0039	-.0013	-.0007	.0743	3.781
SDev	.0014	.0001	.0002	.0002	.0000	.0004	.037
%RSD	.3858	1338.	5.768	12.26	5.675	.4963	.9919

#1	.3587	-.0001	.0037	-.0012	-.0008	.0740	3.754
#2	.3606	.0001	.0040	-.0014	-.0007	.0745	3.808

Elem	La4086	Li6707	Mg2790	Mn2576	Mo2020	Na5889	Na3302
Units	ppm	ppm	ppm	ppm	ppm	ppm	ppm
Avge	.0015	.1257	.0404	.0038	.0030	151.8	98.19
SDev	.0004	.0006	.0017	.0000	.0003	1.5	.91
%RSD	25.52	.5173	4.245	.1704	8.458	.9599	.9296

#1	.0017	.1252	.0416	.0038	.0032	150.7	97.54
#2	.0012	.1262	.0392	.0038	.0029	152.8	98.83

Elem	Ni2316	P_1782	2203/1	2203/2	Pd3404	S_1820	Sb2068
Units	ppm	ppm	ppm	ppm	ppm	ppm	ppm
Avge	-.0003	.0401	.0041	-.0039	.0042	5.695	-.0004
SDev	.0005	.0012	.0033	.0014	.0003	.013	.0025
%RSD	184.3	3.078	80.51	35.41	7.098	.2253	581.7

#1	.0001	.0410	.0018	-.0029	.0040	5.686	.0014
#2	-.0007	.0392	.0065	-.0048	.0044	5.704	-.0022

Elem	Sc3613	1960/1	1960/2	Si2881	Pb220	Se196	Sn1899
Units	%R	ppm	ppm	ppm	ppm	ppm	ppm
Avge	104.0	.0088	-.0043	28.50	-.0012	.0001	.0008
SDev	.2	.0044	.0028	.06	.0002	.0004	.0026
%RSD	.1849	49.58	65.67	.2235	16.52	582.9	306.4

#1	103.8	.0057	-.0023	28.46	-.0013	.0004	.0026
#2	104.1	.0119	-.0063	28.55	-.0011	-.0002	-.0010

Elem	Sr4215	Th2837	Ti3372	Tl1908	U_3859	V_2924	W_2079
Units	ppm	ppm	ppm	ppm	ppm	ppm	ppm
Avge	.0009	-.0010	.0011	.0064	.0704	.0049	.0023
SDev	.0000	.0002	.0008	.0004	.0080	.0000	.0014
%RSD	3.440	24.06	69.16	6.968	11.31	.2288	57.85

#1	.0009	-.0012	.0016	.0061	.0760	.0049	.0033
#2	.0008	-.0008	.0006	.0067	.0648	.0049	.0014

Elem	Y_3710	Zn2062	Zn3496
Units	ppm	ppm	ppm
Avge	-.0000	.0034	.0014
SDev	.0001	.0003	.0002
%RSD	308.4	9.326	17.23
#1	.0000	.0036	.0016
#2	-.0001	.0032	.0012

IntStd	1	2	3	4	5	6	7
Mode	*Counts	Time	NOTUSED	NOTUSED	NOTUSED	NOTUSED	NOTUSED
Elem	Sc	--	--	--	--	--	--
Wavlen	361.384	--	--	--	--	--	--
Avge	635058	10000	--	--	--	--	--
SDev	1194.303	.0000000	--	--	--	--	--
%RSD	.1880622	.0000000	--	--	--	--	--
#1	634213	10000	--	--	--	--	--
#2	635902	10000	--	--	--	--	--

Analysis Report

11/12/03 01:10:36 PM

page 1

Method: DAILY1 Sample Name: 237588

Operator:

Run Time: 11/12/03 13:05:52

Comment:

Mode: CONC Corr. Factor: 1

Elem	Ag3280	Al3082	As1890	B_2496	Ba4934	Be3130	Bi2230
Units	ppm	ppm	ppm	ppm	ppm	ppm	ppm
Avge	.00001	.1202	.0347	.1566	.0003	-.0000	.0085
SDev	.00030	.0037	.0056	.0017	.0000	.0002	.0004
%RSD	2554.8	3.071	16.29	1.071	5.052	343.4	4.204

#1	-.00020	.1228	.0307	.1578	.0003	.0001	.0083
#2	.00023	.1176	.0387	.1554	.0003	-.0002	.0088

Elem	Ca3179	Cd2265	Co2286	Cr2677	Cu3247	Fe2714	K_7664
Units	ppm	ppm	ppm	ppm	ppm	ppm	ppm
Avge	.3588	.0001	.0038	-.0011	-.0010	.0459	3.816
SDev	.0006	.0000	.0000	.0000	.0002	.0053	.057
%RSD	.1767	18.13	.6734	.9531	25.16	11.48	1.493

#1	.3583	.0001	.0038	-.0011	-.0008	.0496	3.776
#2	.3592	.0001	.0038	-.0011	-.0011	.0421	3.857

Elem	La4086	Li6707	Mg2790	Mn2576	Mo2020	Na5889	Na3302
Units	ppm	ppm	ppm	ppm	ppm	ppm	ppm
Avge	.0019	.1268	.0364	.0041	.0027	153.1	98.99
SDev	.0004	.0005	.0003	.0000	.0001	1.2	.67
%RSD	21.99	.4168	.8642	.5215	5.207	.7855	.6807

#1	.0016	.1264	.0362	.0040	.0026	152.3	98.51
#2	.0022	.1272	.0366	.0041	.0028	154.0	99.47

Elem	Ni2316	P_1782	2203/1	2203/2	Pd3404	S_1820	Sb2068
Units	ppm	ppm	ppm	ppm	ppm	ppm	ppm
Avge	-.0005	.0327	.0041	-.0058	.0033	5.805	.0002
SDev	.0001	.0173	.0009	.0006	.0031	.004	.0007
%RSD	13.07	52.92	21.75	10.09	91.95	.0680	307.5

#1	-.0006	.0205	.0048	-.0062	.0012	5.808	.0008
#2	-.0005	.0450	.0035	-.0053	.0055	5.803	-.0003

Elem	Sc3613	1960/1	1960/2	Si2881	Pb220	Se196	Sn1899
Units	%R	ppm	ppm	ppm	ppm	ppm	ppm
Avge	104.7	.0109	-.0034	28.50	-.0024	.0014	.0061
SDev	.9	.0015	.0024	.21	.0001	.0011	.0027
%RSD	.8654	13.55	69.53	.7295	3.528	80.10	43.35

#1	104.1	.0098	-.0017	28.65	-.0025	.0021	.0042
#2	105.3	.0119	-.0051	28.36	-.0024	.0006	.0080

Elem	Sr4215	Th2837	Ti3372	Tl1908	U_3859	V_2924	W_2079
Units	ppm	ppm	ppm	ppm	ppm	ppm	ppm
Avge	.0009	-.0019	.0009	.0104	.0714	.0053	.0018
SDev	.0000	.0013	.0003	.0043	.0077	.0000	.0024
%RSD	.6928	70.68	36.95	41.54	10.77	.7442	135.7

#1	.0009	-.0029	.0007	.0074	.0659	.0053	.0001
#2	.0009	-.0010	.0011	.0135	.0768	.0053	.0034

Analysis Report

11/12/03 01:10:36 PM

page 2

Elem	Y_3710	Zn2062	Zr3496
Units	ppm	ppm	ppm
Avge	.0000	.0026	.0020
SDev	.0000	.0002	.0012
%RSD	190.4	6.577	57.05

#1	-.0000	.0027	.0012
#2	.0000	.0025	.0029

IntStd	1	2	3	4	5	6	7
Mode	*Counts	Time	NOTUSED	NOTUSED	NOTUSED	NOTUSED	NOTUSED
Elem	Sc	--	--	--	--	--	--
Wavlen	361.384	--	--	--	--	--	--
Avge	639535	10000	--	--	--	--	--
SDev	5552.203	.0000000	--	--	--	--	--
%RSD	.8681624	.0000000	--	--	--	--	--
#1	635609	10000	--	--	--	--	--
#2	643461	10000	--	--	--	--	--

Analysis Report

11/12/03 01:15:28 PM

page 1

Method: DAILY1 Sample Name: 237589
 Run Time: 11/12/03 13:10:46
 Comment:
 Mode: CONC Corr. Factor: 1

Operator:

Elem	Ag3280	Al3082	As1890	B_2496	Ba4934	Be3130	Bi2230
Units	ppm	ppm	ppm	ppm	ppm	ppm	ppm
Avg	.00030	.1739	.0316	.1573	.0005	-.0001	.0084
SDev	.00090	.0057	.0018	.0020	.0000	.0001	.0010
%RSD	301.97	3.278	5.550	1.254	6.843	125.5	12.28

#1	-.00034	.1779	.0304	.1587	.0005	-.0000	.0076
#2	.00093	.1698	.0328	.1559	.0005	-.0002	.0091

Elem	Ca3179	Cd2265	Co2286	Cr2677	Cu3247	Fe2714	K_7664
Units	ppm	ppm	ppm	ppm	ppm	ppm	ppm
Avg	.4091	.0000	.0042	-.0005	-.0004	.1899	3.834
SDev	.0017	.0001	.0001	.0006	.0002	.0084	.046
%RSD	.4238	277.6	3.707	107.9	47.77	4.440	1.211

#1	.4079	-.0000	.0041	-.0010	-.0003	.1839	3.801
#2	.4104	.0001	.0043	-.0001	-.0005	.1959	3.867

Elem	La4086	Li6707	Mg2790	Mn2576	Mo2020	Na5889	Na3302
Units	ppm	ppm	ppm	ppm	ppm	ppm	ppm
Avg	.0019	.1269	.0449	.0340	.0025	153.4	98.94
SDev	.0006	.0002	.0029	.0000	.0005	.1	.99
%RSD	33.22	.1439	6.510	.1114	21.97	.0845	1.002

#1	.0015	.1271	.0428	.0340	.0021	153.5	98.24
#2	.0024	.1268	.0469	.0340	.0029	153.3	99.64

Elem	Ni2316	P_1782	2203/1	2203/2	Pd3404	S_1820	Sb2068
Units	ppm	ppm	ppm	ppm	ppm	ppm	ppm
Avg	.0000	.0299	.0021	-.0027	.0023	5.690	.0008
SDev	.0001	.0019	.0023	.0021	.0011	.025	.0018
%RSD	322.6	6.430	106.5	81.00	45.87	.4459	223.9

#1	-.0000	.0286	.0005	-.0042	.0016	5.672	-.0005
#2	.0001	.0313	.0037	-.0011	.0031	5.708	.0021

Elem	Sc3613	1960/1	1960/2	Si2881	Pb220	Se196	Sn1899
Units	%R	ppm	ppm	ppm	ppm	ppm	ppm
Avg	103.2	.0134	.0007	28.73	-.0011	.0050	.0096
SDev	.9	.0010	.0038	.16	.0022	.0029	.0031
%RSD	.9054	7.687	506.8	.5443	207.5	57.45	32.01

#1	102.6	.0127	-.0019	28.84	-.0026	.0030	.0075
#2	103.9	.0142	.0034	28.62	.0005	.0070	.0118

Elem	Sr4215	Th2837	Ti3372	Tl1908	U_3859	V_2924	W_2079
Units	ppm	ppm	ppm	ppm	ppm	ppm	ppm
Avg	.0011	-.0026	.0011	-.0154	.0816	.0055	-.0016
SDev	.0000	.0004	.0001	.0003	.0006	.0004	.0019
%RSD	2.361	16.60	4.485	1.947	.7879	8.046	121.8

#1	.0011	-.0023	.0011	-.0156	.0821	.0051	-.0002
#2	.0011	-.0028	.0010	-.0152	.0812	.0058	-.0030

Elem	Y_3710	Zn2062	Zr3496
Units	ppm	ppm	ppm
Avge	.0001	.0062	.0014
SDev	.0001	.0001	.0001
%RSD	89.48	1.740	8.199

#1	.0000	.0061	.0013
#2	.0002	.0062	.0015

IntStd	1	2	3	4	5	6	7
Mode	*Counts	Time	NOTUSED	NOTUSED	NOTUSED	NOTUSED	NOTUSED
Elem	Sc	--	--	--	--	--	--
Wavlen	361.384	--	--	--	--	--	--
Avge	630548	10000	--	--	--	--	--
SDev	5670.289	.0000000	--	--	--	--	--
%RSD	.8992630	.0000000	--	--	--	--	--
#1	626539	10000	--	--	--	--	--
#2	634558	10000	--	--	--	--	--

Analysis Report

11/12/03 01:20:24 PM

page 1

Method: DAILY1 Sample Name: 237590

Operator:

Run Time: 11/12/03 13:15:38

Comment:

Mode: CONC Corr. Factor: 1

Elem	Ag3280	Al3082	As1890	B_2496	Ba4934	Be3130	Bi2230
Units	ppm	ppm	ppm	ppm	ppm	ppm	ppm
Avge	.00012	.0435	.0004	.2587	.0561	-.0004	.0046
SDev	.00052	.0021	.0002	.0028	.0007	.0000	.0040
%RSD	428.16	4.885	60.11	1.067	1.187	9.771	87.92

#1	.00049	.0420	.0005	.2568	.0556	-.0003	.0075
#2	-.00025	.0450	.0002	.2607	.0565	-.0004	.0017

Elem	Ca3179	Cd2265	Co2286	Cr2677	Cu3247	Fe2714	K_7664
Units	ppm	ppm	ppm	ppm	ppm	ppm	ppm
Avge	28.15	.0002	.0046	-.0013	-.0004	.0247	9.669
SDev	.08	.0002	.0003	.0005	.0002	.0068	.027
%RSD	.2830	129.1	7.660	34.14	43.91	27.38	.2823

#1	28.20	.0004	.0049	-.0010	-.0006	.0295	9.650
#2	28.09	.0000	.0044	-.0016	-.0003	.0199	9.689

Elem	La4086	Li6707	Mg2790	Mn2576	Mo2020	Na5889	Na3302
Units	ppm	ppm	ppm	ppm	ppm	ppm	ppm
Avge	.0023	.1046	28.89	.0749	.0055	133.3	84.93
SDev	.0006	.0002	.10	.0001	.0003	.4	.20
%RSD	25.57	.1772	.3298	.0740	5.704	.3297	.2349

#1	.0027	.1047	28.83	.0749	.0057	133.6	85.07
#2	.0019	.1044	28.96	.0748	.0053	133.0	84.79

Elem	Ni2316	P_1782	2203/1	2203/2	Pd3404	S_1820	Sb2068
Units	ppm	ppm	ppm	ppm	ppm	ppm	ppm
Avge	.0012	.0218	.0041	-.0049	.0046	45.01	.0013
SDev	.0006	.0129	.0002	.0007	.0015	1.20	.0020
%RSD	51.47	59.36	5.814	13.92	31.37	2.657	154.0

#1	.0016	.0309	.0039	-.0044	.0057	44.17	.0027
#2	.0008	.0126	.0042	-.0054	.0036	45.86	-.0001

Elem	Sc3613	1960/1	1960/2	Si2881	Pb220	Se196	Sn1899
Units	%R	ppm	ppm	ppm	ppm	ppm	ppm
Avge	103.4	.0168	.0010	13.69	-.0019	.0063	.0074
SDev	.2	.0013	.0032	.02	.0004	.0025	.0022
%RSD	.1645	7.469	305.1	.1151	19.70	40.02	29.35

#1	103.5	.0177	.0033	13.68	-.0016	.0081	.0059
#2	103.3	.0160	-.0012	13.70	-.0022	.0045	.0090

Elem	Sr4215	Th2837	Ti3372	Tl1908	U_3859	V_2924	W_2079
Units	ppm	ppm	ppm	ppm	ppm	ppm	ppm
Avge	.4127	-.0074	.0006	.0055	.0575	-.0010	.0048
SDev	.0035	.0024	.0008	.0009	.0089	.0010	.0001
%RSD	.8419	32.33	118.6	16.05	15.46	92.51	2.469

#1	.4103	-.0057	.0012	.0061	.0637	-.0004	.0047
#2	.4152	-.0091	.0001	.0048	.0512	-.0017	.0049

Elem	Y_3710	Zn2062	Zr3496
Units	ppm	ppm	ppm
Avge	.0000	.0019	.0019
SDev	.0001	.0003	.0007
%RSD	166.3	16.27	38.85

#1	.0001	.0021	.0024
#2	-.0000	.0017	.0014

IntStd	1	2	3	4	5	6	7
Mode	*Counts	Time	NOTUSED	NOTUSED	NOTUSED	NOTUSED	NOTUSED
Elem	Sc	--	--	--	--	--	--
Wavlen	361.384	--	--	--	--	--	--
Avge	631833	10000	--	--	--	--	--
SDev	1042.275	.00000000	--	--	--	--	--
%RSD	.1649606	.00000000	--	--	--	--	--
#1	632570	10000	--	--	--	--	--
#2	631096	10000	--	--	--	--	--

Analysis Report

11/12/03 01:25:16 PM

page 1

Method: DAILY1 Sample Name: 237591

Operator:

Run Time: 11/12/03 13:20:34

Comment:

Mode: CONC Corr. Factor: 1

Elem	Ag3280	Al3082	As1890	B_2496	Ba4934	Be3130	Bi2230
Units	ppm	ppm	ppm	ppm	ppm	ppm	ppm
Avge	.00202	.0493	.0051	.2600	.0561	-.0002	.0119
SDev	.00310	.0069	.0011	.0007	.0006	.0000	.0072
%RSD	153.09	14.05	22.68	.2828	1.071	6.826	60.83

#1	-.00017	.0444	.0043	.2595	.0557	-.0002	.0068
#2	.00421	.0542	.0059	.2606	.0565	-.0003	.0170

Elem	Ca3179	Cd2265	Co2286	Cr2677	Cu3247	Fe2714	K_7664
Units	ppm	ppm	ppm	ppm	ppm	ppm	ppm
Avge	28.28	.0006	.0051	-.0020	.0019	.0106	9.673
SDev	.00	.0005	.0019	.0007	.0002	.0157	.086
%RSD	.0085	84.89	36.49	34.43	8.137	148.2	.8869

#1	28.29	.0002	.0038	-.0025	.0018	-.0005	9.612
#2	28.28	.0009	.0064	-.0015	.0020	.0216	9.733

Elem	La4086	Li6707	Mg2790	Mn2576	Mo2020	Na5889	Na3302
Units	ppm	ppm	ppm	ppm	ppm	ppm	ppm
Avge	.0038	.1041	28.91	.0713	.0048	132.9	83.61
SDev	.0029	.0000	.02	.0001	.0002	.4	.83
%RSD	77.66	.0303	.0711	.1517	4.436	.2764	.9922

#1	.0017	.1041	28.92	.0713	.0047	132.7	83.02
#2	.0058	.1041	28.89	.0714	.0050	133.2	84.20

Elem	Ni2316	P_1782	2203/1	2203/2	Pd3404	S_1820	Sb2068
Units	ppm	ppm	ppm	ppm	ppm	ppm	ppm
Avge	-.0005	.0420	.0113	-.0092	.0066	43.74	.0086
SDev	.0002	.0031	.0076	.0049	.0059	.22	.0079
%RSD	37.92	7.403	67.53	53.11	89.52	.4929	91.51

#1	-.0006	.0398	.0059	-.0057	.0024	43.59	.0030
#2	-.0004	.0442	.0167	-.0126	.0107	43.90	.0142

Elem	Sc3613	1960/1	1960/2	Si2881	Pb220	Se196	Sn1899
Units	%R	ppm	ppm	ppm	ppm	ppm	ppm
Avge	104.5	.0182	-.0018	13.65	-.0023	.0049	.0163
SDev	.6	.0069	.0064	.04	.0007	.0019	.0054
%RSD	.5462	37.72	362.7	.2624	30.03	39.59	32.86

#1	104.1	.0134	.0027	13.68	-.0018	.0063	.0126
#2	104.9	.0231	-.0062	13.63	-.0028	.0035	.0201

Elem	Sr4215	Th2837	Ti3372	Tl1908	U_3859	V_2924	W_2079
Units	ppm	ppm	ppm	ppm	ppm	ppm	ppm
Avge	.4130	-.0082	.0015	.0038	.0917	-.0012	.0104
SDev	.0003	.0020	.0021	.0008	.0509	.0007	.0048
%RSD	.0659	23.82	139.2	20.18	55.49	52.63	45.74

#1	.4128	-.0068	.0000	.0044	.0557	-.0017	.0071
#2	.4132	-.0096	.0029	.0033	.1277	-.0008	.0138

Analysis Report

11/12/03 01:25:16 PM

page 2

Elem	Y_3710	Zn2062	Zr3496
Units	ppm	ppm	ppm
Avge	-.0001	.0013	.0009
SDev	.0002	.0014	.0009
%RSD	288.9	103.5	103.5

#1	-.0002	.0004	.0002
#2	.0001	.0023	.0016

IntStd	1	2	3	4	5	6	7
Mode	*Counts	Time	NOTUSED	NOTUSED	NOTUSED	NOTUSED	NOTUSED
Elem	Sc	--	--	--	--	--	--
Wavlen	361.384	--	--	--	--	--	--
Avge	638193	10000	--	--	--	--	--
SDev	3464.823	.0000000	--	--	--	--	--
%RSD	.5429115	.0000000	--	--	--	--	--
#1	635743	10000	--	--	--	--	--
#2	640643	10000	--	--	--	--	--

Analysis Report

11/12/03 01:30:08 PM

page 1

Method: DAILY1 Sample Name: 237592

Operator:

Run Time: 11/12/03 13:25:27

Comment:

Mode: CONC Corr. Factor: 1

Elem	Ag3280	Al3082	As1890	B_2496	Ba4934	Be3130	Bi2230
Units	ppm	ppm	ppm	ppm	ppm	ppm	ppm
Avge	-.00164	.0317	.0005	.2597	.0567	-.0002	.0040
SDev	.00106	.0013	.0001	.0015	.0008	.0000	.0044
%RSD	64.311	4.094	21.14	.5813	1.432	7.928	109.0

#1	-.00089	.0326	.0005	.2608	.0561	-.0002	.0072
#2	-.00239	.0308	.0006	.2586	.0573	-.0002	.0009

Elem	Ca3179	Cd2265	Co2286	Cr2677	Cu3247	Fe2714	K_7664
Units	ppm	ppm	ppm	ppm	ppm	ppm	ppm
Avge	28.23	-.0001	.0036	-.0004	.0000	.0811	9.735
SDev	.47	.0002	.0000	.0016	.0007	.0431	.082
%RSD	1.681	170.2	.6565	384.2	3641.	53.09	.8402

#1	28.57	.0000	.0036	.0007	.0005	.1116	9.677
#2	27.90	-.0003	.0036	-.0015	-.0005	.0507	9.793

Elem	La4086	Li6707	Mg2790	Mn2576	Mo2020	Na5889	Na3302
Units	ppm	ppm	ppm	ppm	ppm	ppm	ppm
Avge	.0007	.1047	28.89	.0712	.0063	133.9	83.95
SDev	.0015	.0006	.12	.0008	.0017	.5	.21
%RSD	208.9	.5375	.4152	1.056	27.38	.3454	.2549

#1	.0018	.1043	28.97	.0718	.0075	133.6	84.11
#2	-.0003	.1051	28.80	.0707	.0051	134.2	83.80

Elem	Ni2316	P_1782	2203/1	2203/2	Pd3404	S_1820	Sb2068
Units	ppm	ppm	ppm	ppm	ppm	ppm	ppm
Avge	.0007	.0338	.0036	-.0027	.0028	44.27	-.0023
SDev	.0006	.0149	.0004	.0009	.0003	.70	.0061
%RSD	81.63	44.00	11.96	31.00	10.52	1.590	263.8

#1	.0012	.0443	.0033	-.0021	.0030	43.77	.0020
#2	.0003	.0233	.0039	-.0033	.0026	44.77	-.0066

Elem	Sc3613	1960/1	1960/2	Si2881	Pb220	Se196	Sn1899
Units	%R	ppm	ppm	ppm	ppm	ppm	ppm
Avge	104.1	.0123	-.0026	13.73	-.0006	.0024	.0098
SDev	.6	.0011	.0015	.08	.0004	.0014	.0009
%RSD	.6090	8.869	58.77	.5751	68.38	56.87	8.688

#1	104.6	.0115	-.0036	13.79	-.0003	.0014	.0105
#2	103.7	.0131	-.0015	13.68	-.0009	.0034	.0092

Elem	Sr4215	Th2837	Ti3372	Tl1908	U_3859	V_2924	W_2079
Units	ppm	ppm	ppm	ppm	ppm	ppm	ppm
Avge	.4161	-.0062	-.0005	-.0120	.0422	.0005	.0062
SDev	.0031	.0015	.0013	.0070	.0207	.0019	.0035
%RSD	.7557	23.87	286.8	57.98	49.06	355.9	56.93

#1	.4139	-.0073	.0005	-.0071	.0568	.0019	.0087
#2	.4183	-.0052	-.0014	-.0169	.0276	-.0008	.0037

Analysis Report

11/12/03 01:30:08 PM

page 2

Elem	Y_3710	Zn2062	Zr3496
Units	ppm	ppm	ppm
Avg	.0003	.0164	.0022
SDev	.0004	.0007	.0014
%RSD	162.2	4.069	60.79

#1	.0005	.0168	.0032
#2	-.0000	.0159	.0013

IntStd	1	2	3	4	5	6	7
Mode	*Counts	Time	NOTUSED	NOTUSED	NOTUSED	NOTUSED	NOTUSED
Elem	Sc	--	--	--	--	--	--
Wavlen	361.384	--	--	--	--	--	--
Avg	636258	10000	--	--	--	--	--
SDev	3877.066	.0000000	--	--	--	--	--
%RSD	.6093549	.0000000	--	--	--	--	--
#1	638999	10000	--	--	--	--	--
#2	633516	10000	--	--	--	--	--

Analysis Report

11/12/03 01:35:00 PM

page 1

Method: DAILY1 Sample Name: 237593

Operator:

Run Time: 11/12/03 13:30:18

Comment:

Mode: CONC Corr. Factor: 1

Elem	Ag3280	Al3082	As1890	B_2496	Ba4934	Be3130	Bi2230
Units	ppm	ppm	ppm	ppm	ppm	ppm	ppm
Avge	.00062	.0424	.0025	.0028	.0004	-.0002	.0046
SDev	.00057	.0014	.0000	.0013	.0004	.0000	.0009
%RSD	92.255	3.255	.0608	48.32	104.3	11.39	19.21

#1	.00102	.0434	.0025	.0037	.0006	-.0003	.0052
#2	.00021	.0414	.0025	.0018	.0001	-.0002	.0039

Elem	Ca3179	Cd2265	Co2286	Cr2677	Cu3247	Fe2714	K_7664
Units	ppm	ppm	ppm	ppm	ppm	ppm	ppm
Avge	10.07	.0001	.0038	-.0012	.0035	-.0005	5.243
SDev	.18	.0003	.0008	.0013	.0008	.0320	.112
%RSD	1.809	295.0	21.52	112.7	21.64	6061.	2.128

#1	10.20	.0003	.0043	-.0002	.0041	.0221	5.164
#2	9.939	-.0001	.0032	-.0021	.0030	-.0231	5.322

Elem	La4086	Li6707	Mg2790	Mn2576	Mo2020	Na5889	Na3302
Units	ppm	ppm	ppm	ppm	ppm	ppm	ppm
Avge	.0024	.0003	.0645	.0001	-.0013	39.00	24.59
SDev	.0006	.0002	.0472	.0003	.0009	.46	.17
%RSD	26.79	69.60	73.19	372.9	75.37	1.180	.7001

#1	.0029	.0004	.0979	.0003	-.0006	38.67	24.47
#2	.0020	.0001	.0311	-.0001	-.0019	39.32	24.72

Elem	Ni2316	P_1782	2203/1	2203/2	Pd3404	S_1820	Sb2068
Units	ppm	ppm	ppm	ppm	ppm	ppm	ppm
Avge	.0006	.0137	.0015	-.0017	.0030	.1141	.0032
SDev	.0012	.0169	.0001	.0004	.0024	.0768	.0005
%RSD	213.9	123.9	8.580	26.02	82.06	67.32	15.65

#1	.0015	.0256	.0016	-.0020	.0047	.1684	.0035
#2	-.0003	.0017	.0014	-.0014	.0012	.0598	.0028

Elem	Sc3613	1960/1	1960/2	Si2881	Pb220	Se196	Sn1899
Units	%R	ppm	ppm	ppm	ppm	ppm	ppm
Avge	105.3	.0120	.0009	19.38	-.0006	.0046	.0030
SDev	1.3	.0004	.0035	.05	.0003	.0024	.0001
%RSD	1.278	3.234	382.9	.2657	39.79	52.73	1.896

#1	106.2	.0123	.0034	19.42	-.0008	.0063	.0030
#2	104.3	.0117	-.0015	19.34	-.0005	.0029	.0029

Elem	Sr4215	Th2837	Ti3372	Tl1908	U_3859	V_2924	W_2079
Units	ppm	ppm	ppm	ppm	ppm	ppm	ppm
Avge	.0013	-.0016	.0002	.0055	.0867	-.0001	.0013
SDev	.0006	.0009	.0002	.0028	.0339	.0013	.0003
%RSD	48.08	53.37	122.8	51.06	39.05	1574.	19.83

#1	.0017	-.0010	.0003	.0035	.1107	.0008	.0015
#2	.0008	-.0022	.0000	.0075	.0628	-.0010	.0011

Elem	Y_3710	Zn2062	Zr3496
Units	ppm	ppm	ppm
Avg	.0000	.0020	.0076
SDev	.0003	.0007	.0038
%RSD	676.6	36.49	49.31

#1	.0003	.0025	.0103
#2	-.0002	.0015	.0050

IntStd	1	2	3	4	5	6	7
Mode	*Counts	Time	NOTUSED	NOTUSED	NOTUSED	NOTUSED	NOTUSED
Elem	Sc	--	--	--	--	--	--
Wavlen	361.384	--	--	--	--	--	--
Avg	643068	10000	--	--	--	--	--
SDev	8172.033	.0000000	--	--	--	--	--
%RSD	1.270787	.0000000	--	--	--	--	--
#1	648847	10000	--	--	--	--	--
#2	637290	10000	--	--	--	--	--

Analysis Report

QC Standard

11/12/03 01:43:55 PM

page 1

Method: DAILY1 Sample Name: icv/ccv

Operator:

Run Time: 11/12/03 13:39:10

Comment:

Mode: CONC Corr. Factor: 1

Elem	Ag3280	Al3082	As1890	B_2496	Ba4934	Be3130	Bi2230
Units	ppm	ppm	ppm	ppm	ppm	ppm	ppm
Avg	1.0296	10.08	5.221	4.925	10.22	01.101	5.156
SDev	.0003	.04	.006	.037	.08	.011	.024
%RSD	.02710	.3743	.1090	.7546	.7355	.9566	.4711

#1	1.0294	10.05	5.217	4.899	10.17	1.093	5.139
#2	1.0298	10.11	5.225	4.951	10.28	01.108	5.173

Errors	QC Pass	QC Pass	QC Pass	QC Pass	QC Pass	QC Fail	QC Pass
Value	1.0000	10.00	5.000	5.000	10.00	1.000	5.000
Range	10.000	10.00	10.00	10.00	10.00	10.00	10.00

Elem	Ca3179	Cd2265	Co2286	Cr2677	Cu3247	Fe2714	K_7664
Units	ppm	ppm	ppm	ppm	ppm	ppm	ppm
Avg	20.36	1.027	5.084	2.019	2.097	10.23	18.89
SDev	.11	.003	.013	.006	.011	.01	.17
%RSD	.5416	.2519	.2595	.3181	.5425	.0778	.8972

#1	20.28	1.025	5.075	2.015	2.089	10.22	18.77
#2	20.44	1.029	5.094	2.024	2.105	10.24	19.01

Errors	QC Pass	QC Pass	QC Pass	QC Pass	QC Pass	QC Pass	QC Pass
Value	20.00	1.000	5.000	2.000	2.000	10.00	20.00
Range	10.00	10.00	10.00	10.00	10.00	10.00	10.00

Elem	La4086	Li6707	Mg2790	Mn2576	Mo2020	Na5889	Na3302
Units	ppm	ppm	ppm	ppm	ppm	ppm	ppm
Avg	5.103	4.855	20.31	1.043	5.193	047.24	27.66
SDev	.047	.035	.05	.005	.056	.22	.18
%RSD	.9286	.7140	.2607	.4323	1.075	.4600	.6572

#1	5.070	4.831	20.27	1.040	5.153	047.09	27.54
#2	5.137	4.880	20.34	1.046	5.232	047.39	27.79

Errors	QC Pass	QC Pass	QC Pass	QC Pass	QC Pass	QC Fail	QC Pass
Value	5.000	5.000	20.00	1.000	5.000	30.00	30.00
Range	10.00	10.00	10.00	10.00	10.00	10.00	10.00

Elem	Ni2316	P_1782	2203/1	2203/2	Pd3404	S_1820	Sb2068
Units	ppm	ppm	ppm	ppm	ppm	ppm	ppm
Avg	4.999	5.120	5.173	5.226	1.045	1.042	1.054
SDev	.019	.079	.001	.018	.002	.006	.009
%RSD	.3884	1.545	.0243	.3389	.2318	.6159	.8298

#1	4.986	5.064	5.174	5.238	1.043	1.047	1.048
#2	5.013	5.176	5.172	5.213	1.047	1.037	1.060

Errors	QC Pass	QC Pass	NOCHECK	NOCHECK	QC Pass	QC Pass	QC Pass
Value	5.000	5.000			1.000	1.000	1.000
Range	10.00	10.00			10.00	10.00	10.00

Elem	Sc3613	1960/1	1960/2	Si2881	Pb220	Se196	Sn1899

Analysis Report

QC Standard

11/12/03 01:43:55 PM

Units	%R	ppm	ppm	ppm	ppm	ppm	ppm
Avge	101.2	5.530	5.494	5.015	5.203	Q5.501	5.061
SDev	.6	.014	.023	.013	.012	.011	.026
%RSD	.5833	.2470	.4123	.2619	.2348	.1916	.5201
#1	100.7	5.520	5.510	5.006	5.212	Q5.508	5.042
#2	101.6	5.539	5.478	5.025	5.194	5.493	5.079
Errors	NOCHECK	NOCHECK	NOCHECK	QC Pass	QC Pass	QC Fail	QC Pass
Value				5.000	5.000	5.000	5.000
Range				10.00	10.00	10.00	10.00
Elem	Sr4215	Th2837	Ti3372	Tl1908	U_3859	V_2924	W_2079
Units	ppm	ppm	ppm	ppm	ppm	ppm	ppm
Avge	5.284	.9618	4.983	Q5.572	1.066	5.137	1.099
SDev	.042	.0057	.015	.037	.014	.010	.007
%RSD	.8025	.5959	.3013	.6722	1.339	.1910	.5948
#1	5.254	.9577	4.972	Q5.546	1.056	5.130	1.094
#2	5.314	.9658	4.993	Q5.598	1.076	5.144	Q1.103
Errors	QC Pass	QC Pass	QC Pass	QC Fail	QC Pass	QC Pass	QC Pass
Value	5.000	1.000	5.000	5.000	1.000	5.000	1.000
Range	10.00	10.00	10.00	10.00	10.00	10.00	10.00
Elem	Y_3710	Zn2062	Zr3496				
Units	ppm	ppm	ppm				
Avge	5.108	1.039	4.848				
SDev	.017	.003	.034				
%RSD	.3245	.3111	.6970				
#1	5.097	1.037	4.824				
#2	5.120	1.042	4.872				
Errors	QC Pass	QC Pass	QC Pass				
Value	5.000	1.000	5.000				
Range	10.00	10.00	10.00				

Analysis Report

QC Standard

11/12/03 01:43:55 PM

page 3

IntStd	1	2	3	4	5	6	7
Mode	*Counts	Time	NOTUSED	NOTUSED	NOTUSED	NOTUSED	NOTUSED
Elem	Sc	--	--	--	--	--	--
Wavlen	361.384	--	--	--	--	--	--
Avge	618050	10000	--	--	--	--	--
SDev	3573.718	.0000000	--	--	--	--	--
%RSD	.5782247	.0000000	--	--	--	--	--
#1	615523	10000	--	--	--	--	--
#2	620577	10000	--	--	--	--	--

Analysis Report

Blank Sample

11/12/03 01:50:50 PM

page 1

Method: DAILY1

Sample Name: icb/ccb

Operator:

Run Time: 11/12/03 13:46:05

Comment:

Mode: CONC Corr. Factor: 1

Elem	Ag3280	Al3082	As1890	B_2496	Ba4934	Be3130	Bi2230
Units	ppm	ppm	ppm	ppm	ppm	ppm	ppm
Avge	.00015	.0239	.0013	.0140	.0026	-.0001	.0029
SDev	.00153	.0069	.0023	.0013	.0019	.0001	.0037
%RSD	1004.6	28.92	168.3	9.052	73.23	164.5	127.6

#1	-.00093	.0190	-.0003	.0149	.0040	.0000	.0003
#2	.00124	.0288	.0030	.0131	.0013	-.0001	H.0056

Errors	LC Pass	LC Pass	LC Pass	LC Pass	LC Pass	LC Pass	LC Pass
High	.00500	.0500	.0050	.0500	.0050	.0050	.0050
Low	-.00500	-.0500	-.0050	-.0500	-.0050	-.0050	-.0050

Elem	Ca3179	Cd2265	Co2286	Cr2677	Cu3247	Fe2714	K_7664
Units	ppm	ppm	ppm	ppm	ppm	ppm	ppm
Avge	.0095	.0003	.0018	-.0010	-.0002	-.0061	.0332
SDev	.0015	.0004	.0006	.0006	.0005	.0071	.0285
%RSD	15.36	135.9	30.98	55.60	303.8	116.2	85.67

#1	.0106	.0000	.0022	-.0006	.0002	-.0011	.0131
#2	.0085	.0005	.0014	-.0015	-.0006	-.0111	.0534

Errors	LC Pass	LC Pass	LC Pass	LC Pass	LC Pass	LC Pass	LC Pass
High	.0500	.0050	.0050	.0050	.0050	.0250	.1000
Low	-.0500	-.0050	-.0050	-.0050	-.0050	-.0250	-.1000

Elem	La4086	Li6707	Mg2790	Mn2576	Mo2020	Na5889	Na3302
Units	ppm	ppm	ppm	ppm	ppm	ppm	ppm
Avge	.0020	.0009	.0057	.0001	H.0081	.0186	H.2237
SDev	.0002	.0004	.0072	.0002	.0029	.0004	.2019
%RSD	11.59	53.20	126.5	286.5	35.76	2.004	90.24

#1	.0018	.0012	.0006	.0002	H.0101	.0184	H.0810
#2	.0021	.0005	.0108	-.0001	H.0061	.0189	H.3665

Errors	LC Pass	LC Pass	LC Pass	LC Pass	LC High	LC Pass	LC High
High	.0050	.0050	.0500	.0050	.0050	.0500	.0500
Low	-.0050	-.0050	-.0500	-.0050	-.0050	-.0500	-.0500

Elem	Ni2316	P_1782	2203/1	2203/2	Pd3404	S_1820	Sb2068
Units	ppm	ppm	ppm	ppm	ppm	ppm	ppm
Avge	.0008	.0037	.0085	-.0052	.0047	H.0156	-.0023
SDev	.0008	.0054	.0106	.0073	.0018	.0016	.0020
%RSD	99.88	146.5	124.8	141.7	37.83	10.11	87.52

#1	.0014	-.0001	.0010	.0000	.0034	H.0167	-.0038
#2	.0002	.0075	.0161	-.0104	H.0059	H.0144	-.0009

Errors	LC Pass	LC Pass	NOCHECK	NOCHECK	LC Pass	LC High	LC Pass
High	.0050	.0100			.0050	.0100	.0100
Low	-.0050	-.0100			-.0050	-.0100	-.0100

Elem	Sc3613	1960/1	1960/2	Si2881	Pb220	Se196	Sn1899

Analysis Report

Blank Sample

11/12/03 01:50:50 PM

page 2

Units	%R	ppm	ppm	ppm	ppm	ppm	ppm
Avg	101.5	.0097	-.0017	-.0008	-.0006	.0021	.0016
SDev	.3	.0086	.0081	.0048	.0013	.0026	.0026
%RSD	.2555	88.42	469.5	592.7	220.7	123.1	158.7
#1	101.4	.0036	.0040	-.0042	.0003	.0039	-.0002
#2	101.7	.0157	-.0075	.0026	-.0016	.0003	.0034
Errors	NOCHECK	NOCHECK	NOCHECK	LC Pass	LC Pass	LC Pass	LC Pass
High				.0100	.0030	.0050	.0050
Low				-.0100	-.0030	-.0050	-.0050
Elem	Sr4215	Th2837	Ti3372	Ti1908	U_3859	V_2924	W_2079
Units	ppm	ppm	ppm	ppm	ppm	ppm	ppm
Avg	.0013	-.0039	.0011	.0025	.0323	.0012	.0023
SDev	.0009	.0050	.0002	.0065	.0271	.0018	.0021
%RSD	66.96	129.1	19.28	260.0	83.77	145.9	93.66
#1	.0020	-.0003	.0009	.0071	.0132	.0025	.0008
#2	.0007	-.0074	.0012	-.0021	.0515	-.0000	.0037
Errors	LC Pass	LC Pass	LC Pass	LC Pass	LC Pass	LC Pass	LC Pass
High	.0050	.0100	.0050	.0100	.1000	.0050	.0100
Low	-.0050	-.0100	-.0050	-.0100	-.1000	-.0050	-.0100
Elem	Y_3710	Zn2062	Zr3496				
Units	ppm	ppm	ppm				
Avg	.0013	.0008	.0020				
SDev	.0010	.0003	.0013				
%RSD	77.87	39.88	66.09				
#1	.0020	.0006	.0029				
#2	.0006	.0010	.0011				
Errors	LC Pass	LC Pass	LC Pass				
High	.0050	.0050	.0050				
Low	-.0050	-.0050	-.0050				

010118

Analysis Report

Blank Sample

11/12/03 01:50:50 PM

page 3

IntStd	1	2	3	4	5	6	7
Mode	*Counts	Time	NOTUSED	NOTUSED	NOTUSED	NOTUSED	NOTUSED
Elem	Sc	--	--	--	--	--	--
Wavlen	361.384	--	--	--	--	--	--
Avge	620384	10000	--	--	--	--	--
SDev	1571.898	.0000000	--	--	--	--	--
%RSD	.2533753	.0000000	--	--	--	--	--
#1	619272	10000	--	--	--	--	--
#2	621495	10000	--	--	--	--	--

Analysis Report

11/12/03 01:55:42 PM

page 1

Method: DAILY1 Sample Name: 237594

Operator:

Run Time: 11/12/03 13:51:01

Comment:

Mode: CONC Corr. Factor: 1

Elem	Ag3280	Al3082	As1890	B_2496	Ba4934	Be3130	Bi2230
Units	ppm	ppm	ppm	ppm	ppm	ppm	ppm
Avg	.00009	.0347	.0030	.0067	.0003	-.0003	.0064
SDev	.00025	.0019	.0010	.0003	.0001	.0001	.0030
%RSD	279.58	5.559	35.22	4.073	24.67	36.95	46.09

#1	.00027	.0360	.0037	.0069	.0004	-.0003	.0043
#2	-.00009	.0333	.0022	.0065	.0003	-.0002	.0085

Elem	Ca3179	Cd2265	Co2286	Cr2677	Cu3247	Fe2714	K_7664
Units	ppm	ppm	ppm	ppm	ppm	ppm	ppm
Avg	9.912	-.0000	.0036	-.0015	-.0005	-.0139	5.238
SDev	.032	.0001	.0003	.0003	.0001	.0005	.109
%RSD	.3265	235.7	8.739	19.45	29.81	3.739	2.078

#1	9.889	.0000	.0033	-.0013	-.0006	-.0135	5.161
#2	9.935	-.0001	.0038	-.0017	-.0004	-.0143	5.315

Elem	La4086	Li6707	Mg2790	Mn2576	Mo2020	Na5889	Na3302
Units	ppm	ppm	ppm	ppm	ppm	ppm	ppm
Avg	.0018	-.0000	.0103	-.0001	.0002	39.20	24.37
SDev	.0002	.0002	.0006	.0000	.0005	.16	.72
%RSD	9.958	2436.	5.359	28.65	220.5	.3983	2.969

#1	.0017	-.0001	.0106	-.0001	-.0001	39.09	23.85
#2	.0019	.0001	.0099	-.0001	.0006	39.31	24.88

Elem	Ni2316	P_1782	2203/1	2203/2	Pd3404	S_1820	Sb2068
Units	ppm	ppm	ppm	ppm	ppm	ppm	ppm
Avg	-.0002	.0182	.0061	-.0046	.0037	.0086	-.0008
SDev	.0002	.0003	.0016	.0005	.0004	.0030	.0013
%RSD	125.3	1.708	25.90	11.66	12.26	34.99	166.7

#1	-.0003	.0180	.0072	-.0042	.0040	.0107	-.0017
#2	-.0000	.0184	.0050	-.0050	.0033	.0065	.0001

Elem	Sc3613	1960/1	1960/2	Si2881	Pb220	Se196	Sn1899
Units	%R	ppm	ppm	ppm	ppm	ppm	ppm
Avg	103.7	.0135	.0028	19.41	-.0011	.0064	.0016
SDev	.2	.0070	.0051	.01	.0009	.0011	.0006
%RSD	.2383	51.63	183.0	.0324	83.18	17.19	38.48

#1	103.5	.0086	.0064	19.41	-.0004	.0072	.0021
#2	103.8	.0185	-.0008	19.40	-.0017	.0056	.0012

Elem	Sr4215	Th2837	Ti3372	Tl1908	U_3859	V_2924	W_2079
Units	ppm	ppm	ppm	ppm	ppm	ppm	ppm
Avg	.0006	-.0025	.0001	.0243	.0637	-.0002	.0010
SDev	.0000	.0005	.0001	.0052	.0021	.0002	.0021
%RSD	5.453	18.87	287.1	21.15	3.321	82.41	215.2

#1	.0006	-.0029	.0002	.0280	.0652	-.0001	-.0005
#2	.0005	-.0022	-.0001	.0207	.0622	-.0003	.0024

Elem	Y_3710	Zn2062	Zr3496
Units	ppm	ppm	ppm
Avg	.0000	.0022	.0052
SDev	.0001	.0009	.0005
%RSD	162.1	42.98	8.902

#1	.0001	.0028	.0055
#2	-.0000	.0015	.0049

IntStd	1	2	3	4	5	6	7
Mode	*Counts	Time	NOTUSED	NOTUSED	NOTUSED	NOTUSED	NOTUSED
Elem	Sc	--	--	--	--	--	--
Wavlen	361.384	--	--	--	--	--	--
Avg	633282	10000	--	--	--	--	--
SDev	1501.895	.0000000	--	--	--	--	--
%RSD	.2371605	.0000000	--	--	--	--	--
#1	632220	10000	--	--	--	--	--
#2	634344	10000	--	--	--	--	--

Analysis Report

QC Standard

11/12/03 02:04:35 PM

page 1

Method: DAILY1 Sample Name: icv/ccv

Operator:

Run Time: 11/12/03 13:59:51

Comment:

Mode: CONC Corr. Factor: 1

Elem	Ag3280	Al3082	As1890	B_2496	Ba4934	Be3130	Bi2230
Units	ppm	ppm	ppm	ppm	ppm	ppm	ppm
Avg	1.0212	10.12	5.171	4.904	10.32	1.072	5.177
SDev	.0071	.02	.017	.027	.00	.005	.001
%RSD	.69636	.2253	.3325	.5490	.0066	.4490	.0131

#1	1.0162	10.11	5.159	4.885	10.32	1.068	5.177
#2	1.0262	10.14	5.184	4.923	10.32	1.075	5.178

Errors	QC Pass	QC Pass	QC Pass	QC Pass	QC Pass	QC Pass	QC Pass
Value	1.0000	10.00	5.000	5.000	10.00	1.000	5.000
Range	10.000	10.00	10.00	10.00	10.00	10.00	10.00

Elem	Ca3179	Cd2265	Co2286	Cr2677	Cu3247	Fe2714	K_7664
Units	ppm	ppm	ppm	ppm	ppm	ppm	ppm
Avg	20.08	1.016	5.041	2.005	2.117	10.16	19.20
SDev	.08	.005	.026	.011	.008	.05	.13
%RSD	.4069	.4691	.5107	.5517	.3961	.5154	.6906

#1	20.02	1.013	5.023	1.997	2.112	10.13	19.11
#2	20.14	1.020	5.059	2.013	2.123	10.20	19.30

Errors	QC Pass	QC Pass	QC Pass	QC Pass	QC Pass	QC Pass	QC Pass
Value	20.00	1.000	5.000	2.000	2.000	10.00	20.00
Range	10.00	10.00	10.00	10.00	10.00	10.00	10.00

Elem	La4086	Li6707	Mg2790	Mn2576	Mo2020	Na5889	Na3302
Units	ppm	ppm	ppm	ppm	ppm	ppm	ppm
Avg	5.136	4.937	20.31	1.040	5.155	048.09	27.63
SDev	.035	.015	.12	.005	.025	.19	.25
%RSD	.6807	.3019	.5882	.5289	.4773	.4047	.9227

#1	5.160	4.926	20.22	1.036	5.137	047.95	27.45
#2	5.111	4.947	20.39	1.044	5.172	048.22	27.81

Errors	QC Pass	QC Pass	QC Pass	QC Pass	QC Pass	QC Fail	QC Pass
Value	5.000	5.000	20.00	1.000	5.000	30.00	30.00
Range	10.00	10.00	10.00	10.00	10.00	10.00	10.00

Elem	Ni2316	P_1782	2203/1	2203/2	Pd3404	S_1820	Sb2068
Units	ppm	ppm	ppm	ppm	ppm	ppm	ppm
Avg	4.961	5.172	5.210	5.273	1.047	1.057	1.058
SDev	.013	.011	.059	.022	.004	.008	.002
%RSD	.2532	.2192	1.131	.4252	.4124	.7765	.1447

#1	4.952	5.180	5.168	5.257	1.044	1.051	1.059
#2	4.970	5.164	5.252	5.289	1.050	1.063	1.057

Errors	QC Pass	QC Pass	NOCHECK	NOCHECK	QC Pass	QC Pass	QC Pass
Value	5.000	5.000			1.000	1.000	1.000
Range	10.00	10.00			10.00	10.00	10.00

Elem	Sc3613	1960/1	1960/2	Si2881	Pb220	Se196	Sn1899

Units	%R	ppm	ppm	ppm	ppm	ppm	ppm
Avge	99.00	5.597	5.583	5.021	5.247	Q5.582	5.010
SDev	.10	.072	.054	.032	.035	.060	.029
%RSD	.0992	1.283	.9698	.6326	.6585	1.074	.5764
#1	98.93	5.546	5.545	4.998	5.223	Q5.540	4.990
#2	99.07	5.647	5.621	5.043	5.271	Q5.624	5.031
Errors	NOCHECK	NOCHECK	NOCHECK	QC Pass	QC Pass	QC Fail	QC Pass
Value				5.000	5.000	5.000	5.000
Range				10.00	10.00	10.00	10.00
Elem	Sr4215	Th2837	Ti3372	Ti1908	U_3859	V_2924	W_2079
Units	ppm	ppm	ppm	ppm	ppm	ppm	ppm
Avge	5.313	.9606	4.961	Q5.583	1.052	5.119	1.084
SDev	.034	.0028	.010	.092	.007	.025	.000
%RSD	.6349	.2936	.1966	1.652	.6815	.4982	.0139
#1	5.337	.9586	4.954	Q5.518	1.047	5.101	1.084
#2	5.289	.9626	4.968	Q5.648	1.057	5.137	1.084
Errors	QC Pass	QC Pass	QC Pass	QC Fail	QC Pass	QC Pass	QC Pass
Value	5.000	1.000	5.000	5.000	1.000	5.000	1.000
Range	10.00	10.00	10.00	10.00	10.00	10.00	10.00
Elem	Y_3710	Zn2062	Zr3496				
Units	ppm	ppm	ppm				
Avge	5.103	1.020	4.889				
SDev	.016	.007	.025				
%RSD	.3182	.6929	.5179				
#1	5.115	1.015	4.871				
#2	5.092	1.025	4.907				
Errors	QC Pass	QC Pass	QC Pass				
Value	5.000	1.000	5.000				
Range	10.00	10.00	10.00				

010123

Analysis Report

QC Standard

11/12/03 02:04:35 PM

page 3

IntStd	1	2	3	4	5	6	7
Mode	*Counts	Time	NOTUSED	NOTUSED	NOTUSED	NOTUSED	NOTUSED
Elem	Sc	--	--	--	--	--	--
Wavlen	361.384	--	--	--	--	--	--
Avge	604808	10000	--	--	--	--	--
SDev	625.7895	.0000000	--	--	--	--	--
%RSD	.1034690	.0000000	--	--	--	--	--
#1	604366	10000	--	--	--	--	--
#2	605251	10000	--	--	--	--	--

Analysis Report

Blank Sample

11/12/03 02:11:29 PM

page 1

Method: DAILY1 Sample Name: icb/ccb

Operator:

Run Time: 11/12/03 14:06:45

Comment:

Mode: CONC Corr. Factor: 1

Elem	Ag3280	Al3082	As1890	B_2496	Ba4934	Be3130	Bi2230
Units	ppm	ppm	ppm	ppm	ppm	ppm	ppm
Avge	.00078	.0276	.0010	.0158	.0029	.0001	H.0095
SDev	.00019	.0107	.0022	.0031	.0021	.0002	.0023
%RSD	24.181	38.72	224.0	19.79	70.71	420.0	24.31

#1	.00091	.0351	.0026	.0180	.0044	.0002	H.0111
#2	.00064	.0200	-.0006	.0136	.0015	-.0001	H.0079

Errors	LC Pass	LC Pass	LC Pass	LC Pass	LC Pass	LC Pass	LC High
High	.00500	.0500	.0050	.0500	.0050	.0050	.0050
Low	-.00500	-.0500	-.0050	-.0500	-.0050	-.0050	-.0050

Elem	Ca3179	Cd2265	Co2286	Cr2677	Cu3247	Fe2714	K_7664
Units	ppm	ppm	ppm	ppm	ppm	ppm	ppm
Avge	.0125	.0004	.0021	-.0009	-.0000	-.0050	.0477
SDev	.0050	.0001	.0014	.0003	.0006	.0002	.0061
%RSD	40.07	36.28	68.09	34.26	3659.	4.767	12.72

#1	.0161	.0004	.0031	-.0007	.0004	-.0049	.0434
#2	.0090	.0003	.0011	-.0011	-.0005	-.0052	.0520

Errors	LC Pass	LC Pass	LC Pass	LC Pass	LC Pass	LC Pass	LC Pass
High	.0500	.0050	.0050	.0050	.0050	.0250	.1000
Low	-.0500	-.0050	-.0050	-.0050	-.0050	-.0250	-.1000

Elem	La4086	Li6707	Mg2790	Mn2576	Mo2020	Na5889	Na3302
Units	ppm	ppm	ppm	ppm	ppm	ppm	ppm
Avge	.0028	.0010	.0109	.0001	H.0086	.0208	H.4180
SDev	.0014	.0006	.0055	.0003	.0037	.0030	.0104
%RSD	48.59	55.50	50.47	261.6	43.10	14.34	2.484

#1	.0038	.0014	.0147	.0003	H.0112	.0229	H.4106
#2	.0019	.0006	.0070	-.0001	H.0060	.0187	H.4253

Errors	LC Pass	LC Pass	LC Pass	LC Pass	LC High	LC Pass	LC High
High	.0050	.0050	.0500	.0050	.0050	.0500	.0500
Low	-.0050	-.0050	-.0500	-.0050	-.0050	-.0500	-.0500

Elem	Ni2316	P_1782	2203/1	2203/2	Pd3404	S_1820	Sb2068
Units	ppm	ppm	ppm	ppm	ppm	ppm	ppm
Avge	.0007	.0062	.0063	-.0025	.0030	.0014	.0019
SDev	.0013	.0199	.0002	.0009	.0017	.0097	.0001
%RSD	179.8	323.2	2.924	38.12	56.27	669.1	2.918

#1	.0016	H.0202	.0062	-.0018	.0042	-.0054	.0019
#2	-.0002	-.0079	.0064	-.0031	.0018	.0083	.0018

Errors	LC Pass	LC Pass	NOCHECK	NOCHECK	LC Pass	LC Pass	LC Pass
High	.0050	.0100			.0050	.0100	.0100
Low	-.0050	-.0100			-.0050	-.0100	-.0100

Elem	Sc3613	1960/1	1960/2	Si2881	Pb220	Se196	Sn1899

Analysis Report

Blank Sample

11/12/03 02:11:29 PM

Units	%R	ppm	ppm	ppm	ppm	ppm	ppm
Avge	102.5	.0072	-.0003	.0044	.0004	.0022	.0034
SDev	1.6	.0018	.0026	.0036	.0006	.0023	.0045
%RSD	1.525	24.24	784.9	81.62	126.4	105.7	132.3
#1	101.4	.0085	.0015	.0069	.0009	.0038	H.0066
#2	103.6	.0060	-.0022	.0019	.0000	.0006	.0002
Errors	NOCHECK	NOCHECK	NOCHECK	LC Pass	LC Pass	LC Pass	LC Pass
High				.0100	.0030	.0050	.0050
Low				-.0100	-.0030	-.0050	-.0050
Elem	Sr4215	Th2837	Ti3372	Ti1908	U_3859	V_2924	W_2079
Units	ppm	ppm	ppm	ppm	ppm	ppm	ppm
Avge	.0016	-.0019	.0016	H.0149	.0493	.0011	.0037
SDev	.0011	.0019	.0011	.0047	.0024	.0010	.0032
%RSD	69.20	104.6	69.10	31.40	4.906	88.74	87.18
#1	.0023	-.0005	.0024	H.0182	.0476	.0018	.0060
#2	.0008	-.0032	.0008	H.0116	.0510	.0004	.0014
Errors	LC Pass	LC Pass	LC Pass	LC High	LC Pass	LC Pass	LC Pass
High	.0050	.0100	.0050	.0100	.1000	.0050	.0100
Low	-.0050	-.0100	-.0050	-.0100	-.1000	-.0050	-.0100
Elem	Y_3710	Zn2062	Zr3496				
Units	ppm	ppm	ppm				
Avge	.0015	.0008	.0022				
SDev	.0010	.0003	.0010				
%RSD	69.86	39.90	47.51				
#1	.0022	.0010	.0029				
#2	.0007	.0006	.0014				
Errors	LC Pass	LC Pass	LC Pass				
High	.0050	.0050	.0050				
Low	-.0050	-.0050	-.0050				

010120

Analysis Report

Blank Sample

11/12/03 02:11:29 PM

page 3

IntStd	1	2	3	4	5	6	7
Mode	*Counts	Time	NOTUSED	NOTUSED	NOTUSED	NOTUSED	NOTUSED
Elem	Sc	--	--	--	--	--	--
Wavlen	361.384	--	--	--	--	--	--
Avge	626253	10000	--	--	--	--	--
SDev	9543.113	.00000000	--	--	--	--	--
%RSD	1.523843	.00000000	--	--	--	--	--
#1	619505	10000	--	--	--	--	--
#2	633001	10000	--	--	--	--	--

✓

ICP ANALYSIS

K + Na

PROJ. NO.	PROJECT	TO#	DATE	MATRIX	LOGBK PG
06002.01.141	Div 20	031106-12	11-24-03	Water	SS 124
18023.01.078	Div 18	031113-5	11-24-03	Water	NA

INSTRUMENT: Spectro FILENAME: 031124

INSTRUMENT DL: _____

Starting run ...

Creating high priority queue entries ...

BACKGROUND CORRECTED INTENSITIES

Identity 1 : BLK_SC Identity 2 : Type : STD

Weight : 1.0000 Volume : 1.00 Printed : 2:42:30 PM November 24, 2003

	K_766	Li670	Na589	Sc361
# 1	30.0	-14.5	35.5	4974.0
# 2	28.0	-7.5	38.5	4722.0
Mean	29.0	-11.0	37.0	4848.0
SD	1.4	4.9	2.1	178.2
%RSD	4.9	45.0	5.7	3.7

INTENSITIES

Identity 1 : BLK_SC Identity 2 : Type : STD

Weight : 1.0000 Volume : 1.00 Printed : 2:42:30 PM November 24, 2003

	K_766	Li670	Na589	Sc361
# 1	0.0	-0.0	0.0	4974.0
# 2	0.0	-0.0	0.0	4722.0
Mean	0.0	-0.0	0.0	4848.0
SD	0.0	0.0	0.0	178.2
%RSD	1.2	41.7	9.4	3.7

BACKGROUND CORRECTED INTENSITIES

Identity 1 : CLP_STD1_SC Identity 2 : Type : STD

Weight : 1.0000 Volume : 1.00 Printed : 2:45:36 PM November 24, 2003

	K_766	Li670	Na589
# 1	4846.5	39923.0	24769.0
# 2	4736.5	39034.0	24189.0
Mean	4791.5	39478.5	24479.0
SD	77.8	628.6	410.1
%RSD	1.6	1.6	1.7

INTENSITIES

Identity 1 : CLP_STD1_SC Identity 2 : Type : STD

Weight : 1.0000 Volume : 1.00 Printed : 2:45:36 PM November 24, 2003

	K_766	Li670	Na589
# 1	1.0	8.1	5.0
# 2	1.0	8.1	5.0
Mean	1.0	8.1	5.0
SD	0.0	0.0	0.0
%RSD	0.0	0.0	0.1

BACKGROUND CORRECTED INTENSITIES

Identity 1 : CLP_CCV_SC Identity 2 : Type : CV

Weight : 1.0000 Volume : 1.00 Printed : 2:48:44 PM November 24, 2003

	K_766	Li670	Na589	Sc	Sc361
# 1	1957.5	19923.5	14779.0	4981.0	4981.0
# 2	1906.5	19654.5	14581.0	4914.0	4914.0
Mean	1932.0	19789.0	14680.0	4947.5	4947.5

010127

11/15/03
300cc 03

11-24-03

%RSD 1.9 1.0 1.0 1.0 1.0

APPARENT CONCENTRATIONS

Identity 1 : CLP_CCV_SC Identity 2 : Type : CV

EVOLUTION by Micro-Active Australia Pty Ltd 2:48:56 PM November 24, 2003

010128

Weight : 1.0000 Volume : 1.00 Printed : 2:48:44 PM November 24, 2003

	K_766	Li670	Na589	Sc	Sc361
	ppm	ppm	ppm		ppm
# 1	19.8002	4.9387	29.5014	4981.000 H	102.7471
# 2	19.5433	4.9384	29.5030	4914.000 H	101.3632
Mean	19.6718	4.9385	29.5022	4947.500 H	102.0551
SD	0.1816	0.0002	0.0011	47.376	0.9785
%RSD	0.9232	0.0036	0.0039	0.958	0.9588

Checking calibration verification ...

Identity 1 : CLP_CCV_SC Identity 2 :

Report name	Low limit	Value	High limit
K_766	18.000	19.672	22.000
Li670	4.500	4.939	5.500
Na589	27.000	29.502	33.000

EVOLUTION by Micro-Active Australia Pty Ltd 2:58:58 PM November 24, 2003

BACKGROUND CORRECTED INTENSITIES

Identity 1 : Calibration blank Identity 2 : Type : CB

Weight : 1.0000 Volume : 1.00 Printed : 2:52:52 PM November 24, 2003

	K_766	Li670	Na589	Sc	Sc361
# 1	33.5	-16.0	5.5	5058.5	5058.5
# 2	16.5	-18.0	38.5	4935.5	4935.5
Mean	25.0	-17.0	22.0	4997.0	4997.0
SD	12.0	1.4	23.3	87.0	87.0
%RSD	48.1	8.3	106.1	1.7	1.7

APPARENT CONCENTRATIONS

Identity 1 : Calibration blank Identity 2 : Type : CB

Weight : 1.0000 Volume : 1.00 Printed : 2:52:52 PM November 24, 2003

	K_766	Li670	Na589	Sc	Sc361
	ppm	ppm	ppm		ppm
# 1	0.0328 L	-0.0011 L	-0.0652	5058.500 H	104.3478
# 2 L	-0.1350 L	-0.0017	0.0017	4935.500 H	101.8073
Mean L	-0.0511 L	-0.0014 L	-0.0318	4997.000 H	103.0776
SD	0.1186	0.0004	0.0473	86.974	1.7964
%RSD	232.1138	30.1289	148.9011	1.741	1.7428

Checking calibration blank ...

Identity 1 : Calibration blank Identity 2 :

Report name	CRDL	Value
K_766	0.100	-0.051
Li670	0.010	-0.001
Na589	0.050	-0.032
Sc361	0.000	103.078

BACKGROUND CORRECTED INTENSITIES

Identity 1 : pbw-124v2 Identity 2 : pg55-124 Type : SAMPLE

Weight : 1.0000 Volume : 1.00 Printed : 2:56:00 PM November 24, 2003

	K_766	Li670	Na589	Sc	Sc361
# 1	10.5	-13.0	44.0	5263.0	5263.0
# 2	18.5	-17.0	22.0	5154.0	5154.0
Mean	14.5	-15.0	33.0	5208.5	5208.5
SD	5.7	2.8	15.6	77.1	77.1
%RSD	39.0	18.9	47.1	1.5	1.5

APPARENT CONCENTRATIONS

Identity 1 : pbw-124v2 Identity 2 : pg55-124 Type : SAMPLE

Weight : 1.0000 Volume : 1.00 Printed : 2:56:00 PM November 24, 2003

	K_766	Li670	Na589	Sc	Sc361
	ppm	ppm	ppm		ppm
# 1 L	-0.2040 L	-0.0002	0.0073	5263.000 H	108.5717
# 2 L	-0.1224 L	-0.0013 L	-0.0335	5154.000 H	106.3204
Mean L	-0.1632 L	-0.0008 L	-0.0131	5208.500 H	107.4460
SD	0.0577	0.0007	0.0288	77.075	1.5920
%RSD	35.3462	95.1983	219.5836	1.480	1.4816

010130

EVOLUTION by Micro-Active Australia Pty Ltd 3:05:40 PM November 24, 2003

BACKGROUND CORRECTED INTENSITIES

Identity 1 : lcsw-L24w1 Identity 2 : pg55-124 Type : SAMPLE

Weight : 1.0000 Volume : 1.00 Printed : 2:59:08 PM November 24, 2003

	K_766	Li670	Na589	Sc	Sc361
# 1	1925.5	-24.5	9947.5	5305.0	5305.0
# 2	1891.5	-23.5	9771.5	5201.0	5201.0
Mean	1908.5	-24.0	9859.5	5253.0	5253.0
SD	24.0	0.7	124.5	73.5	73.5
%RSD	1.3	2.9	1.3	1.4	1.4

APPARENT CONCENTRATIONS

Identity 1 : lcsw-L24w1 Identity 2 : pg55-124 Type : SAMPLE

Weight : 1.0000 Volume : 1.00 Printed : 2:59:08 PM November 24, 2003

	K_766	Li670	Na589	Sc	Sc361
	ppm	ppm	ppm		ppm
# 1	18.2636 L	-0.0029	18.6162	5305.000 H	109.4392
# 2	18.3005 L	-0.0028	18.6526	5201.000 H	107.2911
Mean	18.2820 L	-0.0028	18.6344	5253.000 H	108.3652
SD	0.0261	0.0001	0.0258	73.539	1.5189
%RSD	0.1426	3.0729	0.1383	1.400	1.4017

BACKGROUND CORRECTED INTENSITIES

Identity 1 : 237578 Identity 2 : pg55-124 Type : SAMPLE

Weight : 1.0000 Volume : 1.00 Printed : 3:02:16 PM November 24, 2003

	K_766	Li670	Na589	Sc	Sc361
# 1	240.0	342.0	58705.0	5322.0	5322.0
# 2	220.0	332.0	58165.0	5247.0	5247.0
Mean	230.0	337.0	58435.0	5284.5	5284.5
SD	14.1	7.1	381.8	53.0	53.0
%RSD	6.1	2.1	0.7	1.0	1.0

APPARENT CONCENTRATIONS

Identity 1 : 237578 Identity 2 : pg55-124 Type : SAMPLE

Weight : 1.0000 Volume : 1.00 Printed : 3:02:16 PM November 24, 2003

	K_766	Li670	Na589	Sc	Sc361
	ppm	ppm	ppm		ppm
# 1	2.0011	0.0821 H	109.8836	5322.000 H	109.7904
# 2	1.8391	0.0809 H	110.4294	5247.000 H	108.2412
Mean	1.9201	0.0815 H	110.1565	5284.500 H	109.0158
SD	0.1146	0.0009	0.3860	53.033	1.0954
%RSD	5.9671	1.0572	0.3504	1.004	1.0048

BACKGROUND CORRECTED INTENSITIES

Identity 1 : 237579 Identity 2 : pg55-124 Type : SAMPLE

Weight : 1.0000 Volume : 1.00 Printed : 3:05:24 PM November 24, 2003

	K_766	Li670	Na589	Sc	Sc361
# 1	219.5	335.0	58304.0	5564.0	5564.0
# 2	232.5	337.0	58394.0	5489.0	5489.0
Mean	226.0	336.0	58349.0	5526.5	5526.5

%RSD 4.1 0.4 0.1 1.0 1.0

APPARENT CONCENTRATIONS

Identity 1 : 237579 Identity 2 : pg55-124 Type : SAMPLE

010131

EVOLUTION by Micro-Active Australia Pty Ltd 3:15:04 PM November 24, 2003

Weight : 1.0000 Volume : 1.00 Printed : 3:05:24 PM November 24, 2003

	K_766	Li670	Na589	Sc	Sc361
	ppm	ppm	ppm		ppm
# 1	1.7123	0.0771 H	104.3825	5564.000 H	114.7888
# 2	1.8610	0.0786 H	105.9733	5489.000 H	113.2397
Mean	1.7867	0.0778 H	105.1779	5526.500 H	114.0143
SD	0.1052	0.0010	1.1248	53.033	1.0954
%RSD	5.8870	1.3308	1.0694	0.960	0.9607

BACKGROUND CORRECTED INTENSITIES

Identity 1 : 237580 Identity 2 : pg55-124 Type : SAMPLE

Weight : 1.0000 Volume : 1.00 Printed : 3:08:30 PM November 24, 2003

	K_766	Li670	Na589	Sc	Sc361
# 1	218.5	340.0	58357.0	4993.0	4993.0
# 2	223.5	333.0	58374.0	4911.0	4911.0
Mean	221.0	336.5	58365.5	4952.0	4952.0
SD	3.5	4.9	12.0	58.0	58.0
%RSD	1.6	1.5	0.0	1.2	1.2

APPARENT CONCENTRATIONS

Identity 1 : 237580 Identity 2 : pg55-124 Type : SAMPLE

Weight : 1.0000 Volume : 1.00 Printed : 3:08:30 PM November 24, 2003

	K_766	Li670	Na589	Sc	Sc361
	ppm	ppm	ppm		ppm
# 1	1.9329	0.0868 H	116.4343	4993.000 H	102.9949
# 2	2.0223	0.0865 H	118.4142	4911.000 H	101.3012
Mean	1.9776	0.0867 H	117.4242	4952.000 H	102.1481
SD	0.0633	0.0003	1.4000	57.983	1.1976
%RSD	3.1992	0.2904	1.1923	1.171	1.1724

BACKGROUND CORRECTED INTENSITIES

Identity 1 : 237580d Identity 2 : pg55-124 Type : SAMPLE

Weight : 1.0000 Volume : 1.00 Printed : 3:11:38 PM November 24, 2003

	K_766	Li670	Na589	Sc	Sc361
# 1	207.5	337.5	58017.0	4984.0	4984.0
# 2	215.5	333.5	57989.0	4889.0	4889.0
Mean	211.5	335.5	58003.0	4936.5	4936.5
SD	5.7	2.8	19.8	67.2	67.2
%RSD	2.7	0.8	0.0	1.4	1.4

APPARENT CONCENTRATIONS

Identity 1 : 237580d Identity 2 : pg55-124 Type : SAMPLE

Weight : 1.0000 Volume : 1.00 Printed : 3:11:38 PM November 24, 2003

	K_766	Li670	Na589	Sc	Sc361
	ppm	ppm	ppm		ppm
# 1	1.8240	0.0864 H	115.9646	4984.000 H	102.8090
# 2	1.9491	0.0870 H	118.1623	4889.000 H	100.8468
Mean	1.8865	0.0867 H	117.0635	4936.500 H	101.8279
SD	0.0885	0.0004	1.5540	67.175	1.3875

BACKGROUND CORRECTED INTENSITIES

Identity 1 : 237581 Identity 2 : pg55-124 Type : SAMPLE
 Weight : 1.0000 Volume : 1.00 Printed : 3:14:46 PM November 24, 2003
 EVOLUTION by Micro-Active Australia Pty Ltd 3:21:28 PM November 24, 2003

010132

	K_766	Li670	Na589	Sc	Sc361
# 1	351.5	990.5	95345.5	5524.0	5524.0
# 2	357.5	999.5	95358.5	5488.0	5488.0
Mean	354.5	995.0	95352.0	5506.0	5506.0
SD	4.2	6.4	9.2	25.5	25.5
%RSD	1.2	0.6	0.0	0.5	0.5

APPARENT CONCENTRATIONS

Identity 1 : 237581 Identity 2 : pg55-124 Type : SAMPLE
 Weight : 1.0000 Volume : 1.00 Printed : 3:14:46 PM November 24, 2003

	K_766	Li670	Na589	Sc	Sc361
	ppm	ppm	ppm		ppm
# 1	2.9495	0.2241 H	171.9838	5524.000 H	113.9626
# 2	3.0267	0.2275 H	173.1361	5488.000 H	113.2190
Mean	2.9881	0.2258 H	172.5599	5506.000 H	113.5908
SD	0.0547	0.0025	0.8148	25.456	0.5258
%RSD	1.8290	1.0882	0.4722	0.462	0.4629

BACKGROUND CORRECTED INTENSITIES

Identity 1 : 237581s Identity 2 : pg55-124 Type : SAMPLE
 Weight : 1.0000 Volume : 1.00 Printed : 3:17:52 PM November 24, 2003

	K_766	Li670	Na589	Sc	Sc361
# 1	2517.5	994.0	104922.0	5546.5	5546.5
# 2	2508.5	988.0	104744.0	5481.5	5481.5
Mean	2513.0	991.0	104833.0	5514.0	5514.0
SD	6.4	4.2	125.9	46.0	46.0
%RSD	0.3	0.4	0.1	0.8	0.8

APPARENT CONCENTRATIONS

Identity 1 : 237581s Identity 2 : pg55-124 Type : SAMPLE
 Weight : 1.0000 Volume : 1.00 Printed : 3:17:54 PM November 24, 2003

	K_766	Li670	Na589	Sc	Sc361
	ppm	ppm	ppm		ppm
# 1	22.9157	0.2239 H	188.4974	5546.500 H	114.4273
# 2	23.1071	0.2252 H	190.4098	5481.500 H	113.0848
Mean	23.0114	0.2246 H	189.4536	5514.000 H	113.7561
SD	0.1353	0.0009	1.3523	45.962	0.9493
%RSD	0.5880	0.4004	0.7138	0.834	0.8345

BACKGROUND CORRECTED INTENSITIES

Identity 1 : 237582 Identity 2 : pg55-124 Type : SAMPLE
 Weight : 1.0000 Volume : 1.00 Printed : 3:21:02 PM November 24, 2003

	K_766	Li670	Na589	Sc	Sc361
# 1	366.0	1032.0	96850.5	5637.0	5637.0
# 2	376.0	1010.0	97816.5	5587.0	5587.0
Mean	371.0	1021.0	97333.5	5612.0	5612.0
SD	7.1	15.6	683.1	35.4	35.4
%RSD	1.9	1.5	0.7	0.6	0.6

Identity 1 : 237582 Identity 2 : pg55-124 Type : SAMPLE
 Weight : 1.0000 Volume : 1.00 Printed : 3:21:02 PM November 24, 2003

K_766 Li670 Na589 Sc Sc361
 EVOLUTION by Micro-Active Australia Pty Ltd 3:29:36 PM November 24, 2003

010133

	ppm	ppm	ppm		ppm
# 1	3.0158	0.2287 H	171.1961	5637.000 H	116.2966
# 2	3.1371	0.2259 H	174.4525	5587.000 H	115.2639
Mean	3.0764	0.2273 H	172.8243	5612.000 H	115.7802
SD	0.0858	0.0020	2.3026	35.355	0.7303
%RSD	2.7880	0.8827	1.3323	0.630	0.6307

BACKGROUND CORRECTED INTENSITIES

Identity 1 : 237583 Identity 2 : pg55-124 Type : SAMPLE
 Weight : 1.0000 Volume : 1.00 Printed : 3:24:08 PM November 24, 2003

	K_766	Li670	Na589	Sc	Sc361
# 1	351.0	1028.5	97329.0	5645.5	5645.5
# 2	360.0	1017.5	98296.0	5590.5	5590.5
Mean	355.5	1023.0	97812.5	5618.0	5618.0
SD	6.4	7.8	683.8	38.9	38.9
%RSD	1.8	0.8	0.7	0.7	0.7

APPARENT CONCENTRATIONS

Identity 1 : 237583 Identity 2 : pg55-124 Type : SAMPLE
 Weight : 1.0000 Volume : 1.00 Printed : 3:24:08 PM November 24, 2003

	K_766	Li670	Na589	Sc	Sc361
# 1	2.8749	0.2276 H	171.7832	5645.500 H	116.4722
# 2	2.9885	0.2274 H	175.1982	5590.500 H	115.3362
Mean	2.9317	0.2275 H	173.4907	5618.000 H	115.9042
SD	0.0804	0.0002	2.4148	38.891	0.8033
%RSD	2.7414	0.0672	1.3919	0.692	0.6931

BACKGROUND CORRECTED INTENSITIES

Identity 1 : CLP_CCV_SC Identity 2 : Type : CV
 Weight : 1.0000 Volume : 1.00 Printed : 3:29:26 PM November 24, 2003

	K_766	Li670	Na589	Sc	Sc361
# 1	1940.0	19925.5	14742.5	4999.5	4999.5
# 2	1893.0	19471.5	14387.5	4848.5	4848.5
Mean	1916.5	19698.5	14565.0	4924.0	4924.0
SD	33.2	321.0	251.0	106.8	106.8
%RSD	1.7	1.6	1.7	2.2	2.2

APPARENT CONCENTRATIONS

Identity 1 : CLP_CCV_SC Identity 2 : Type : CV
 Weight : 1.0000 Volume : 1.00 Printed : 3:29:26 PM November 24, 2003

	K_766	Li670	Na589	Sc	Sc361
# 1	19.5467	4.9209	29.3192	4999.500 H	103.1292
# 2	19.6690	4.9585	29.5048	4848.500 H	100.0103
Mean	19.6079	4.9397	29.4120	4924.000 H	101.5698
SD	0.0865	0.0266	0.1312	106.773	2.2054
%RSD	0.4412	0.5385	0.4462	2.168	2.1713

Identity 1 : CLP_CCV_SC Identity 2 :
 Report name Low limit Value High limit
 K_766 18.000 19.608 22.000
 Li670 4.500 4.940 5.500
 EVOLUTION by Micro-Active Australia Pty Ltd 3:40:12 PM November 24, 2003

010134

Na589 27.000 29.412 33.000

BACKGROUND CORRECTED INTENSITIES

Identity 1 : Calibration blank Identity 2 : Type : CB
 Weight : 1.0000 Volume : 1.00 Printed : 3:33:40 PM November 24, 2003

	K_766	Li670	Na589	Sc	Sc361
# 1	29.5	-13.0	42.5	5111.0	5111.0
# 2	25.5	-10.0	39.5	5071.0	5071.0
Mean	27.5	-11.5	41.0	5091.0	5091.0
SD	2.8	2.1	2.1	28.3	28.3
%RSD	10.3	18.4	5.2	0.6	0.6

APPARENT CONCENTRATIONS

Identity 1 : Calibration blank Identity 2 : Type : CB
 Weight : 1.0000 Volume : 1.00 Printed : 3:33:40 PM November 24, 2003

	K_766	Li670	Na589	Sc	Sc361
	ppm	ppm	ppm		ppm
# 1 L	-0.0107 L	-0.0003	0.0068	5111.000 H	105.4322
# 2 L	-0.0488	0.0004	0.0016	5071.000 H	104.6060
Mean L	-0.0298	0.0000	0.0042	5091.000 H	105.0191
SD	0.0269	0.0005	0.0037	28.284	0.5842
%RSD	90.3636	3604.7288	88.4777	0.556	0.5563

Checking calibration blank ...

Identity 1 : Calibration blank Identity 2 :

Report name	CRDL	Value
K_766	0.100	-0.030
Li670	0.010	0.000
Na589	0.050	0.004
Sc361	0.000	105.019

BACKGROUND CORRECTED INTENSITIES

Identity 1 : 237584 Identity 2 : pg55-124 Type : SAMPLE
 Weight : 1.0000 Volume : 1.00 Printed : 3:36:46 PM November 24, 2003

	K_766	Li670	Na589	Sc	Sc361
# 1	373.5	559.5	54836.0	5795.5	5795.5
# 2	377.5	547.5	54519.0	5701.5	5701.5
Mean	375.5	553.5	54677.5	5748.5	5748.5
SD	2.8	8.5	224.2	66.5	66.5
%RSD	0.8	1.5	0.4	1.2	1.2

APPARENT CONCENTRATIONS

Identity 1 : 237584 Identity 2 : pg55-124 Type : SAMPLE
 Weight : 1.0000 Volume : 1.00 Printed : 3:36:46 PM November 24, 2003

	K_766	Li670	Na589	Sc	Sc361
	ppm	ppm	ppm		ppm
# 1	2.9912	0.1219 H	94.2448	5795.500 H	119.5704
# 2	3.0814	0.1213 H	95.2456	5701.500 H	117.6288
Mean	3.0363	0.1216 H	94.7452	5748.500 H	118.5996
SD	0.0638	0.0004	0.7077	66.468	1.3729
%RSD	2.1019	0.3681	0.7469	1.156	1.1576

Identity 1 : 237585 Identity 2 : pg55-124 Type : SAMPLE
Weight : 1.0000 Volume : 1.00 Printed : 3:39:54 PM November 24, 2003

010135

	K_766	Li670	Na589	Sc	Sc361
EVOLUTION by Micro-Active Australia Pty Ltd 3:40:12 PM November 24, 2003					
# 1	358.5	555.0	54291.5	5753.5	5753.5
# 2	364.5	543.0	54290.5	5701.5	5701.5
Mean	361.5	549.0	54291.0	5727.5	5727.5
SD	4.2	8.5	0.7	36.8	36.8
%RSD	1.2	1.5	0.0	0.6	0.6

APPARENT CONCENTRATIONS

Identity 1 : 237585 Identity 2 : pg55-124 Type : SAMPLE
Weight : 1.0000 Volume : 1.00 Printed : 3:39:54 PM November 24, 2003

	K_766	Li670	Na589	Sc	Sc361
	ppm	ppm	ppm		ppm
# 1	2.8818	0.1218 H	93.9899	5753.500 H	118.7029
# 2	2.9648	0.1203 H	94.8461	5701.500 H	117.6288
Mean	2.9233	0.1211 H	94.4180	5727.500 H	118.1659
SD	0.0586	0.0011	0.6054	36.770	0.7595
%RSD	2.0056	0.8828	0.6412	0.642	0.6427

BACKGROUND CORRECTED INTENSITIES

Identity 1 : 237586 Identity 2 : pg55-124 Type : SAMPLE

Weight : 1.0000 Volume : 1.00 Printed : 3:43:02 PM November 24, 2003

	K_766	Li670	Na589	Sc	Sc361
# 1	368.0	554.0	53777.0	5689.5	5689.5
# 2	373.0	544.0	54011.0	5641.5	5641.5
Mean	370.5	549.0	53894.0	5665.5	5665.5
SD	3.5	7.1	165.5	33.9	33.9
%RSD	1.0	1.3	0.3	0.6	0.6

APPARENT CONCENTRATIONS

Identity 1 : 237586 Identity 2 : pg55-124 Type : SAMPLE

Weight : 1.0000 Volume : 1.00 Printed : 3:43:02 PM November 24, 2003

	K_766	Li670	Na589	Sc	Sc361
	ppm	ppm	ppm		ppm
# 1	3.0031	0.1230 H	94.1466	5689.500 H	117.3810
# 2	3.0766	0.1218 H	95.3618	5641.500 H	116.3895
Mean	3.0399	0.1224 H	94.7542	5665.500 H	116.8853
SD	0.0520	0.0008	0.8592	33.941	0.7010
%RSD	1.7097	0.6732	0.9068	0.599	0.5998

BACKGROUND CORRECTED INTENSITIES

Identity 1 : 237587 Identity 2 : pg55-124 Type : SAMPLE

Weight : 1.0000 Volume : 1.00 Printed : 3:46:10 PM November 24, 2003

	K_766	Li670	Na589	Sc	Sc361
# 1	360.5	541.5	53959.5	5645.5	5645.5
# 2	355.5	531.5	53837.5	5583.5	5583.5
Mean	358.0	536.5	53898.5	5614.5	5614.5
SD	3.5	7.1	86.3	43.8	43.8
%RSD	1.0	1.3	0.2	0.8	0.8

APPARENT CONCENTRATIONS

Identity 1 : 237587 Identity 2 : pg55-124 Type : SAMPLE

Weight : 1.0000 Volume : 1.00 Printed : 3:46:10 PM November 24, 2003

	K_766	Li670	Na589	Sc	Sc361
	ppm	ppm	ppm		ppm
# 1	2.9609	0.1212 H	95.2032	5645.500 H	116.4722
# 2	2.9514	0.1203 H	96.0434	5583.500 H	115.1916
Mean	2.9562	0.1207 H	95.6233	5614.500 H	115.8319
SD	0.0067	0.0006	0.5941	43.841	0.9055
%RSD	0.2281	0.5247	0.6213	0.781	0.7818

BACKGROUND CORRECTED INTENSITIES

Identity 1 : 237588 Identity 2 : pg55-124 Type : SAMPLE

Weight : 1.0000 Volume : 1.00 Printed : 3:49:18 PM November 24, 2003

	K_766	Li670	Na589	Sc	Sc361
# 1	368.5	561.0	55038.5	5689.0	5689.0
# 2	382.5	559.0	55275.5	5628.0	5628.0
Mean	375.5	560.0	55157.0	5658.5	5658.5

%RSD 2.6 0.3 0.3 0.8 0.8

APPARENT CONCENTRATIONS

Identity 1 : 237588 Identity 2 : pg55-124 Type : SAMPLE

010137

EVOLUTION by Micro-Active Australia Pty Ltd 3:59:00 PM November 24, 2003

Weight : 1.0000 Volume : 1.00 Printed : 3:49:18 PM November 24, 2003

	K_766	Li670	Na589	Sc	Sc361
	ppm	ppm	ppm		ppm
# 1	3.0079	0.1245 H	96.3654	5689.000 H	117.3706
# 2	3.1711	0.1254 H	97.8304	5628.000 H	116.1107
Mean	3.0895	0.1249 H	97.0979	5658.500 H	116.7407
SD	0.1154	0.0006	1.0360	43.134	0.8909
%RSD	3.7349	0.4983	1.0669	0.762	0.7632

BACKGROUND CORRECTED INTENSITIES

Identity 1 : 237589 Identity 2 : pg55-124 Type : SAMPLE

Weight : 1.0000 Volume : 1.00 Printed : 3:52:24 PM November 24, 2003

	K_766	Li670	Na589	Sc	Sc361
# 1	377.0	538.0	55187.0	5697.0	5697.0
# 2	376.0	535.0	55152.0	5597.0	5597.0
Mean	376.5	536.5	55169.5	5647.0	5647.0
SD	0.7	2.1	24.7	70.7	70.7
%RSD	0.2	0.4	0.0	1.3	1.3

APPARENT CONCENTRATIONS

Identity 1 : 237589 Identity 2 : pg55-124 Type : SAMPLE

Weight : 1.0000 Volume : 1.00 Printed : 3:52:26 PM November 24, 2003

	K_766	Li670	Na589	Sc	Sc361
	ppm	ppm	ppm		ppm
# 1	3.0796	0.1193 H	96.4898	5697.000 H	117.5359
# 2	3.1309	0.1208 H	98.1528	5597.000 H	115.4704
Mean	3.1053	0.1200 H	97.3213	5647.000 H	116.5031
SD	0.0363	0.0010	1.1759	70.711	1.4605
%RSD	1.1693	0.8368	1.2083	1.252	1.2536

BACKGROUND CORRECTED INTENSITIES

Identity 1 : 237590 Identity 2 : pg55-124 Type : SAMPLE

Weight : 1.0000 Volume : 1.00 Printed : 3:55:32 PM November 24, 2003

	K_766	Li670	Na589	Sc	Sc361
# 1	919.5	426.5	46156.0	5639.5	5639.5
# 2	905.5	434.5	45734.0	5496.5	5496.5
Mean	912.5	430.5	45945.0	5568.0	5568.0
SD	9.9	5.7	298.4	101.1	101.1
%RSD	1.1	1.3	0.6	1.8	1.8

APPARENT CONCENTRATIONS

Identity 1 : 237590 Identity 2 : pg55-124 Type : SAMPLE

Weight : 1.0000 Volume : 1.00 Printed : 3:55:32 PM November 24, 2003

	K_766	Li670	Na589	Sc	Sc361
	ppm	ppm	ppm		ppm
# 1	8.0357	0.0961 H	81.5109	5639.500 H	116.3482
# 2	8.1224	0.1003 H	82.8681	5496.500 H	113.3946
Mean	8.0790	0.0982 H	82.1895	5568.000 H	114.8714
SD	0.0613	0.0030	0.9597	101.116	2.0885

BACKGROUND CORRECTED INTENSITIES

Identity 1 : 237591 Identity 2 : pg55-124 Type : SAMPLE
 Weight : 1.0000 Volume : 1.00 Printed : 3:58:38 PM November 24, 2003

010138

EVOLUTION by Micro-Active Australia Pty Ltd 4:05:02 PM November 24, 2003

	K_766	Li670	Na589	Sc	Sc361
# 1	930.5	448.0	46316.0	5673.5	5673.5
# 2	930.5	445.0	46266.0	5571.5	5571.5
Mean	930.5	446.5	46291.0	5622.5	5622.5
SD	0.0	2.1	35.4	72.1	72.1
%RSD	0.0	0.5	0.1	1.3	1.3

APPARENT CONCENTRATIONS

Identity 1 : 237591 Identity 2 : pg55-124 Type : SAMPLE
 Weight : 1.0000 Volume : 1.00 Printed : 3:58:38 PM November 24, 2003

	K_766 ppm	Li670 ppm	Na589 ppm	Sc	Sc361 ppm
# 1	8.0849	0.1002 H	81.3031	5673.500 H	117.0505
# 2	8.2385	0.1014 H	82.7034	5571.500 H	114.9437
Mean	8.1617	0.1008 H	82.0032	5622.500 H	115.9971
SD	0.1086	0.0008	0.9902	72.125	1.4897
%RSD	1.3309	0.7853	1.2075	1.283	1.2843

BACKGROUND CORRECTED INTENSITIES

Identity 1 : 237592 Identity 2 : pg55-124 Type : SAMPLE
 Weight : 1.0000 Volume : 1.00 Printed : 4:01:46 PM November 24, 2003

	K_766	Li670	Na589	Sc	Sc361
# 1	929.5	452.0	46705.0	5769.0	5769.0
# 2	951.5	450.0	46875.0	5654.0	5654.0
Mean	940.5	451.0	46790.0	5711.5	5711.5
SD	15.6	1.4	120.2	81.3	81.3
%RSD	1.7	0.3	0.3	1.4	1.4

APPARENT CONCENTRATIONS

Identity 1 : 237592 Identity 2 : pg55-124 Type : SAMPLE
 Weight : 1.0000 Volume : 1.00 Printed : 4:01:46 PM November 24, 2003

	K_766 ppm	Li670 ppm	Na589 ppm	Sc	Sc361 ppm
# 1	7.9371	0.0995 H	80.6281	5769.000 H	119.0230
# 2	8.3039	0.1010 H	82.5693	5654.000 H	116.6477
Mean	8.1205	0.1002 H	81.5987	5711.500 H	117.8354
SD	0.2593	0.0011	1.3726	81.317	1.6796
%RSD	3.1934	1.0792	1.6822	1.424	1.4254

BACKGROUND CORRECTED INTENSITIES

Identity 1 : 237593 Identity 2 : pg55-124 Type : SAMPLE
 Weight : 1.0000 Volume : 1.00 Printed : 4:04:52 PM November 24, 2003

	K_766	Li670	Na589	Sc	Sc361
# 1	610.0	-20.0	15493.5	5731.0	5731.0
# 2	626.0	-19.0	15451.5	5650.0	5650.0
Mean	618.0	-19.5	15472.5	5690.5	5690.5
SD	11.3	0.7	29.7	57.3	57.3
%RSD	1.8	3.6	0.2	1.0	1.0

Identity 1 : 237593 Identity 2 : pg55-124 Type : SAMPLE
 Weight : 1.0000 Volume : 1.00 Printed : 4:04:52 PM November 24, 2003

K_766 Li670 Na589 Sc Sc361
 EVOLUTION by Micro-Active Australia Pty Ltd 4:14:40 PM November 24, 2003

010139

	ppm	ppm	ppm		ppm
# 1	5.1395 L	-0.0015	26.8735	5731.000 H	118.2381
# 2	5.3625 L	-0.0013	27.1858	5650.000 H	116.5651
Mean	5.2510 L	-0.0014	27.0296	5690.500 H	117.4016
SD	0.1577	0.0001	0.2208	57.276	1.1830
%RSD	3.0023	7.7570	0.8169	1.007	1.0077

BACKGROUND CORRECTED INTENSITIES

Identity 1 : CLP_CCV_SC Identity 2 : Type : CV
 Weight : 1.0000 Volume : 1.00 Printed : 4:10:10 PM November 24, 2003

	K_766	Li670	Na589	Sc	Sc361
# 1	1937.5	19830.5	14637.5	4978.5	4978.5
# 2	1899.5	19534.5	14493.5	4872.5	4872.5
Mean	1918.5	19682.5	14565.5	4925.5	4925.5
SD	26.9	209.3	101.8	75.0	75.0
%RSD	1.4	1.1	0.7	1.5	1.5

APPARENT CONCENTRATIONS

Identity 1 : CLP_CCV_SC Identity 2 : Type : CV
 Weight : 1.0000 Volume : 1.00 Printed : 4:10:10 PM November 24, 2003

	K_766	Li670	Na589	Sc	Sc361
	ppm	ppm	ppm		ppm
# 1	19.6047	4.9181	29.2330	4978.500 H	102.6954
# 2	19.6389	4.9501	29.5760	4872.500 H	100.5060
Mean	19.6218	4.9341	29.4045	4925.500 H	101.6007
SD	0.0241	0.0226	0.2425	74.953	1.5481
%RSD	0.1231	0.4581	0.8248	1.522	1.5238

Checking calibration verification ...

Identity 1 : CLP_CCV_SC Identity 2 :

Report name	Low limit	Value	High limit
K_766	18.000	19.622	22.000
Li670	4.500	4.934	5.500
Na589	27.000	29.404	33.000

BACKGROUND CORRECTED INTENSITIES

Identity 1 : Calibration blank Identity 2 : Type : CB
 Weight : 1.0000 Volume : 1.00 Printed : 4:14:22 PM November 24, 2003

	K_766	Li670	Na589	Sc	Sc361
# 1	19.5	-7.0	15.0	5110.0	5110.0
# 2	14.5	5.0	28.0	4995.0	4995.0
Mean	17.0	-1.0	21.5	5052.5	5052.5
SD	3.5	8.5	9.2	81.3	81.3
%RSD	20.8	848.5	42.8	1.6	1.6

APPARENT CONCENTRATIONS

Identity 1 : Calibration blank Identity 2 : Type : CB
 Weight : 1.0000 Volume : 1.00 Printed : 4:14:24 PM November 24, 2003

	K_766	Li670	Na589	Sc	Sc361
	ppm	ppm	ppm		ppm
# 1 L	-0.1108	0.0011 L	-0.0468	5110.000 H	105.4115

Mean	L	-0.1342	0.0026	L	-0.0335	5052.500	H	104.2239
SD		0.0330	0.0021		0.0188	81.317		1.6796
%RSD		24.6224	80.4230		56.1685	1.609		1.6115

EVOLUTION by Micro-Active Australia Pty Ltd 4:24:02 PM November 24, 2003

010140

Checking calibration blank ...

Identity 1 : Calibration blank Identity 2 :

Report name	CRDL	Value	
K_766	0.100	-0.134	Contaminated
Li670	0.010	0.003	
Na589	0.050	-0.034	
Sc361	0.000	104.224	

BACKGROUND CORRECTED INTENSITIES

Identity 1 : 237594 Identity 2 : pg55-124 Type : SAMPLE

Weight : 1.0000 Volume : 1.00 Printed : 4:17:30 PM November 24, 2003

	K_766	Li670	Na589	Sc	Sc361
# 1	652.0	-11.0	15874.5	5947.5	5947.5
# 2	645.0	-12.0	15753.5	5854.5	5854.5
Mean	648.5	-11.5	15814.0	5901.0	5901.0
SD	4.9	0.7	85.6	65.8	65.8
%RSD	0.8	6.1	0.5	1.1	1.1

APPARENT CONCENTRATIONS

Identity 1 : 237594 Identity 2 : pg55-124 Type : SAMPLE

Weight : 1.0000 Volume : 1.00 Printed : 4:17:30 PM November 24, 2003

	K_766	Li670	Na589	Sc	Sc361
	ppm	ppm	ppm		ppm
# 1	5.3026	0.0005	26.5311	5947.500 H	122.7099
# 2	5.3305	0.0003	26.7477	5854.500 H	120.7890
Mean	5.3166	0.0004	26.6394	5901.000 H	121.7495
SD	0.0197	0.0002	0.1532	65.761	1.3583
%RSD	0.3714	44.3237	0.5750	1.114	1.1156

BACKGROUND CORRECTED INTENSITIES

Identity 1 : 237907 Identity 2 : pg55-124 Type : SAMPLE

Weight : 1.0000 Volume : 1.00 Printed : 4:20:38 PM November 24, 2003

	K_766	Li670	Na589	Sc	Sc361
# 1	372.5	55.0	36126.5	5791.5	5791.5
# 2	381.5	55.0	35878.5	5691.5	5691.5
Mean	377.0	55.0	36002.5	5741.5	5741.5
SD	6.4	0.0	175.4	70.7	70.7
%RSD	1.7	0.0	0.5	1.2	1.2

APPARENT CONCENTRATIONS

Identity 1 : 237907 Identity 2 : pg55-124 Type : SAMPLE

Weight : 1.0000 Volume : 1.00 Printed : 4:20:38 PM November 24, 2003

	K_766	Li670	Na589	Sc	Sc361
	ppm	ppm	ppm		ppm
# 1	2.9846	0.0145 H	62.1064	5791.500 H	119.4878
# 2	3.1233	0.0147 H	62.7646	5691.500 H	117.4223
Mean	3.0540	0.0146 H	62.4355	5741.500 H	118.4550
SD	0.0981	0.0001	0.4654	70.711	1.4605
%RSD	3.2119	0.9957	0.7454	1.232	1.2330

Identity 1 : 237907d Identity 2 : pg55-124 Type : SAMPLE
Weight : 1.0000 Volume : 1.00 Printed : 4:23:46 PM November 24, 2003

	K_766	Li670	Na589	Sc	Sc361
EVOLUTION by Micro-Active Australia Pty Ltd 4:32:32 PM November 24, 2003					
# 1	365.5	51.5	36117.5	5844.5	5844.5
# 2	384.5	54.5	36024.5	5725.5	5725.5
Mean	375.0	53.0	36071.0	5785.0	5785.0
SD	13.4	2.1	65.8	84.1	84.1
%RSD	3.6	4.0	0.2	1.5	1.5

010141

APPARENT CONCENTRATIONS

Identity 1 : 237907d Identity 2 : pg55-124 Type : SAMPLE
Weight : 1.0000 Volume : 1.00 Printed : 4:23:46 PM November 24, 2003

	K_766	Li670	Na589	Sc	Sc361
	ppm	ppm	ppm		ppm
# 1	2.8935	0.0137 H	61.5271	5844.500 H	120.5825
# 2	3.1298	0.0145 H	62.6456	5725.500 H	118.1245
Mean	3.0116	0.0141 H	62.0864	5785.000 H	119.3535
SD	0.1671	0.0006	0.7909	84.146	1.7380
%RSD	5.5477	4.3729	1.2738	1.455	1.4562

BACKGROUND CORRECTED INTENSITIES

Identity 1 : 237907s Identity 2 : pg55-124 Type : SAMPLE
Weight : 1.0000 Volume : 1.00 Printed : 4:26:54 PM November 24, 2003

	K_766	Li670	Na589	Sc	Sc361
# 1	2547.0	67.5	45924.0	5704.0	5704.0
# 2	2569.0	53.5	45827.0	5641.0	5641.0
Mean	2558.0	60.5	45875.5	5672.5	5672.5
SD	15.6	9.9	68.6	44.5	44.5
%RSD	0.6	16.4	0.1	0.8	0.8

APPARENT CONCENTRATIONS

Identity 1 : 237907s Identity 2 : pg55-124 Type : SAMPLE
Weight : 1.0000 Volume : 1.00 Printed : 4:26:54 PM November 24, 2003

	K_766	Li670	Na589	Sc	Sc361
	ppm	ppm	ppm		ppm
# 1	22.5391	0.0174 H	80.1828	5704.000 H	117.6805
# 2	22.9938	0.0145 H	80.9078	5641.000 H	116.3792
Mean	22.7665	0.0160 H	80.5453	5672.500 H	117.0298
SD	0.3215	0.0021	0.5126	44.548	0.9201
%RSD	1.4122	12.8517	0.6364	0.785	0.7862

See sample spiked with
2001 SCAL-1 # 4245
DR 11-24-03

BACKGROUND CORRECTED INTENSITIES

Identity 1 : CLP_CCV_SC Identity 2 : Type : CV
Weight : 1.0000 Volume : 1.00 Printed : 4:32:12 PM November 24, 2003

	K_766	Li670	Na589	Sc	Sc361
# 1	1935.0	19797.0	14632.0	4959.5	4959.5
# 2	1923.0	19514.0	14478.0	4871.5	4871.5
Mean	1929.0	19655.5	14555.0	4915.5	4915.5
SD	8.5	200.1	108.9	62.2	62.2
%RSD	0.4	1.0	0.7	1.3	1.3

APPARENT CONCENTRATIONS

Identity 1 : CLP_CCV_SC Identity 2 : Type : CV

	K_766 ppm	Li670 ppm	Na589 ppm	Sc	Sc361 ppm
# 1	19.6552	4.9286	29.3342	4959.500 H	102.3030
EVOLUTION by Micro-Active Australia Pty Ltd 4:42:34 PM November 24, 2003					
# 2	19.8898	4.9459	29.5503	4871.500 H	100.4854
Mean	19.7725	4.9372	29.4423	4915.500 H	101.3942
SD	0.1659	0.0122	0.1528	62.225	1.2853
%RSD	0.8388	0.2477	0.5191	1.266	1.2676

010142

Checking calibration verification ...

Identity 1 : CLP_CCV_SC Identity 2 :

Report name	Low limit	Value	High limit
K_766	18.000	19.773	22.000
Li670	4.500	4.937	5.500
Na589	27.000	29.442	33.000

BACKGROUND CORRECTED INTENSITIES

Identity 1 : Calibration blank Identity 2 : Type : CB

Weight : 1.0000 Volume : 1.00 Printed : 4:36:24 PM November 24, 2003

	K_766	Li670	Na589	Sc	Sc361
# 1	20.5	-11.5	45.0	5087.0	5087.0
# 2	16.5	-12.5	62.0	5008.0	5008.0
Mean	18.5	-12.0	53.5	5047.5	5047.5
SD	2.8	0.7	12.0	55.9	55.9
%RSD	15.3	5.9	22.5	1.1	1.1

APPARENT CONCENTRATIONS

Identity 1 : Calibration blank Identity 2 : Type : CB

Weight : 1.0000 Volume : 1.00 Printed : 4:36:24 PM November 24, 2003

	K_766 ppm	Li670 ppm	Na589 ppm	Sc	Sc361 ppm
# 1 L	-0.0999	0.0000	0.0121	5087.000 H	104.9365
# 2 L	-0.1375 L	-0.0003	0.0473	5008.000 H	103.3048
Mean L	-0.1187 L	-0.0001	0.0297	5047.500 H	104.1206
SD	0.0266	0.0002	0.0249	55.861	1.1538
%RSD	22.4109	152.1709	83.8274	1.107	1.1081

Checking calibration blank ...

Identity 1 : Calibration blank Identity 2 :

Report name	CRDL	Value
K_766	0.100	-0.119 Contaminated
Li670	0.010	-0.000
Na589	0.050	0.030
Sc361	0.000	104.121

**SOUTHWEST RESEARCH INSTITUTE
NUCLEAR PROJECT
CLIENT: Div. 20
TASK ORDER: 031106-12
SRR: 25200
SDG: 237578
CASE: CNWRA
VTSR: November 06, 2003
PROJECT#: 06002.01.141**

Certificates of Analysis

010143

**FISHER SCIENTIFIC
TRACEMETAL GRADE NITRIC ACID
CERTIFICATE OF ANALYSIS**

Catalog No. A509


Lot No: 1103070

Release Date: July, 2003

Expiry Date: July, 2006

<u>Tests</u>	<u>Units</u>	<u>Value</u>
Assay	%	69%
Color	APHA	<10
Aluminum	ppb	<0.2
Antimony	ppb	<0.1
Arsenic	ppb	<0.1
Barium	ppb	<0.1
Beryllium	ppb	<0.1
Bismuth	ppb	<0.1
Boron	ppb	<1
Cadmium	ppb	<0.1
Calcium	ppb	<0.5
Chromium	ppb	<0.1
Cobalt	ppb	<0.1
Copper	ppb	<0.1
Iron	ppb	<0.5
Lead	ppb	<0.1
Lithium	ppb	<0.1
Magnesium	ppb	<0.2
Manganese	ppb	<0.1
Mercury	ppb	<0.2
Molybdenum	ppb	<0.1
Nickel	ppb	<0.1
Potassium	ppb	<0.2
Selenium	ppb	<0.1
Silver	ppb	<0.1
Sodium	ppb	<0.2
Strontium	ppb	<0.1
Thorium	ppb	<0.1
Tin	ppb	<0.1
Titanium	ppb	<0.1
Uranium	ppb	<0.1
Vanadium	ppb	<0.1
Zinc	ppb	<0.2
Zirconium	ppb	<0.1

Element concentrations are at the point of bottling.
Concentrations of some elements in particular, Ca, Si,
K, Na, B, Al, Mg & Mn will increase due to storage in
glass bottles.



Dr. B. McKelvey
QA/QC Manager

Fisher Scientific Chemical Division
Pittsburgh, PA., 15275 Phone (412) 490-8300



Fisher Chemical

A Fisher Scientific Company

INORGANIC LABS/RADCHEM LABS
DATE RECEIVED: 10/10/03
DATE EXPIRED: 10/10/2013
DATE OPENED: 10/13/03
INORG: 9995-4309 PO: F52242

010144

**FISHER SCIENTIFIC
TRACEMETAL GRADE NITRIC ACID
CERTIFICATE OF ANALYSIS**


Catalog No. A509

Lot No: 1103050

Date: May, 2003

<u>Tests</u>	<u>Units</u>	<u>Value</u>
Assay	%	69%
Color	APHA	<10
Aluminum	ppb	<0.2
Antimony	ppb	<0.1
Arsenic	ppb	<0.1
Barium	ppb	<0.1
Beryllium	ppb	<0.1
Bismuth	ppb	<0.1
Boron	ppb	>1
Cadmium	ppb	<0.1
Calcium	ppb	<0.5
Chromium	ppb	<0.1
Cobalt	ppb	<0.1
Copper	ppb	<0.1
Iron	ppb	<0.5
Lead	ppb	<0.1
Lithium	ppb	<0.1
Magnesium	ppb	<0.2
Manganese	ppb	<0.1
Mercury	ppb	<0.2
Molybdenum	ppb	<0.1
Nickel	ppb	<0.1
Potassium	ppb	<0.2
Selenium	ppb	<0.1
Silver	ppb	<0.1
Sodium	ppb	<0.2
Strontium	ppb	<0.1
Thorium	ppb	<0.1
Tin	ppb	<0.1
Titanium	ppb	<0.1
Uranium	ppb	<0.1
Vanadium	ppb	<0.1
Zinc	ppb	<0.2
Zirconium	ppb	<0.1

Element concentrations are at the point of bottling.
Concentrations of some elements in particular, Ca, Si,
K, Na, B, Al, Mg & Mn will increase due to storage in
glass bottles.



Dr. B. McKelvey

QA/QC Manager

Fisher Scientific Chemical Division
Pittsburgh, PA., 15275 Phone (412) 490-8300

**Fisher Chemical**

A Fisher Scientific Company

INORGANIC LABS/RADCHEM LABS
DATE RECEIVED: 09/19/03
DATE EXPIRED: 09/19/2013
DATE OPENED: 09/19/03
INORG: 4935-4040 PO: F500005

010143

FISHER SCIENTIFIC
TRACEMETAL GRADE HYDROCHLORIC ACID

CERTIFICATE OF ANALYSIS

Catalog No. A508

Lot No: 4103030

Date: March, 2003

<u>Tests</u>	<u>Units</u>	<u>Value</u>
Assay	%	36%
Color	APHA	<10
Aluminum	ppb	<0.5
Antimony	ppb	<0.1
Arsenic	ppb	<0.1
Barium	ppb	<0.1
Beryllium	ppb	<0.1
Bismuth	ppb	<0.1
Boron	ppb	<1
Cadmium	ppb	<0.1
Calcium	ppb	<0.5
Chromium	ppb	<0.1
Cobalt	ppb	<0.1
Copper	ppb	<0.1
Iron	ppb	<0.5
Lead	ppb	<0.1
Lithium	ppb	<0.1
Magnesium	ppb	<0.5
Manganese	ppb	<0.1
Mercury	ppb	<0.2
Molybdenum	ppb	<0.1
Nickel	ppb	<0.1
Potassium	ppb	<0.1
Selenium	ppb	<0.1
Silver	ppb	<0.1
Sodium	ppb	<0.5
Strontium	ppb	<0.1
Thorium	ppb	<0.1
Tin	ppb	<0.1
Titanium	ppb	<0.1
Uranium	ppb	<0.1
Vanadium	ppb	<0.1
Zinc	ppb	<0.5
Zirconium	ppb	<0.1

Element concentrations are at the point of bottling.
 Concentrations of some elements in particular, Ca, Si, K,
 Na, B, Al, Mg & Mn will increase due to storage in glass
 bottles.

B McKelvey

Dr. B. McKelvey
 QA/QC Manager

Fisher Scientific Chemical Division
 Pittsburgh, PA., 15275 Phone (412) 490-8300



Fisher Chemical

A Fisher Scientific Company

INORGANIC LABS/RADCHEM LABS
 DATE RECEIVED: 10/10/03
 DATE EXPIRED: 10/10/2013
 DATE OPENED: 10/13/03
 INORG: 4289-4994 PO: F502412

010148

FISHER SCIENTIFIC
TRACEMETAL GRADE HYDROCHLORIC ACID

CERTIFICATE OF ANALYSIS

Catalog No. A508

Lot No: 4102110

Date: December, 2002

<u>Tests</u>	<u>Units</u>	<u>Value</u>
Assay	%	36%
Color	APHA	<10
Aluminum	ppb	<0.5
Antimony	ppb	<0.1
Arsenic	ppb	<0.1
Barium	ppb	>0.1
Beryllium	ppb	>0.1
Bismuth	ppb	>0.1
Boron	ppb	>0.1
Cadmium	ppb	>0.1
Calcium	ppb	>0.5
Chromium	ppb	>0.1
Cobalt	ppb	>0.1
Copper	ppb	>0.1
Iron	ppb	>0.5
Lead	ppb	>0.1
Lithium	ppb	>0.1
Magnesium	ppb	>0.5
Manganese	ppb	>0.1
Mercury	ppb	>0.2
Molybdenum	ppb	>0.1
Nickel	ppb	>0.1
Potassium	ppb	>0.1
Selenium	ppb	>0.1
Silver	ppb	>0.1
Sodium	ppb	>0.5
Strontium	ppb	>0.1
Thorium	ppb	>0.1
Tin	ppb	>0.1
Titanium	ppb	>0.1
Uranium	ppb	>0.1
Vanadium	ppb	>0.1
Zinc	ppb	>0.5
Zirconium	ppb	>0.1

Element concentrations are at the point of bottling.
 Concentrations of some elements in particular,
 Ca, Si, K, Na, B, Al, Mg & Mn will increase due to
 storage in glass bottles.

B McKelvey

Dr. B. McKelvey
 QA/QC Manager

Fisher Scientific Chemical Division
 Pittsburgh, PA., 15275 Phone (412) 490-8300



INORGANIC LABS/RADCHEM LABS
 DATE RECEIVED: 09/10/03
 DATE EXPIRED: 09/10/2013
 DATE OPENED: 09/17/03
 INORG: 4947-4952 PO: F50008

INORGANIC LABS/RADCHEM LABS
 DATE RECEIVED: 12/31/03
 DATE EXP. 12/31/04
 DATE OPENED: 11/3/03
 INORG: 4306 PO: PS201
 010147
 DK

SPEXertificate™

Certificate of Reference Material

Catalog Number: SPIKE-1 **Lot No.:** 25-23AS
Description: Spike Sample Standard 1
Matrix: 5% Nitric Acid/tr Tartaric Acid - HF

This ASSURANCE® certified reference material, CRM, is intended primarily for use as a calibration standard or quality control standard for inorganic spectroscopic instrumentation such as ICPOES, DCP, AA, ICPMS, and XRF. It can be employed in USEPA, ASTM and other methods relevant to the certified properties listed below.

The CRM is prepared from high purity single element concentrates of individual elements using Class A laboratory ware to give precise concentration.

Refer to side 2 for details of measurement uncertainties.

Instrumental Analysis by ICP Spectrometer:

Element	Labeled (mg/L)	Measured (mg/L)	NIST SRM	Element	Labeled (mg/L)	Measured (mg/L)	NIST SRM
Al	200	199.51	3101a	Pb	50	49.98	3128
As	200	199.89	3103a	Sb	50	50.02	3102a
Ba	200	199.68	3104a	V	50	49.95	3165
Se	200	200.10	3149	Zn	50	50.02	3168a
TL	200	200.07	3158	Cu	25	25.34	3114
Fe	100	99.91	3126a	Cr	20	20.04	3112a
Co	50	50.25	3113	Ag	5	5.00	3151
Mn	50	49.98	3132	Be	5	5.00	3105a
Ni	50	50.11	3136	Cd	5	4.99	3108

Spex Reference Multi: Lot #2-61BD, 17-55AS, 19-85ASREF

Balances are calibrated regularly with weight sets traceable to NIST#s 32856, 32867 and others. This CRM is guaranteed stable and accurate to +/- 0.5% on the average of all the certified concentrations with no single component exceeding +/- 2%. This guarantee is valid for a period of one year from the date of certification only when the material is kept tightly capped and transported and stored under laboratory conditions.

Date of Certification: OCT - - 2003 Certifying Officer: N. Kocherakota

SPEXertificate™

Certificate of Reference Material

Catalog Number: ICAL-1 **Lot No.:** 25-176AS
Description: Instrument Calibration Standard 1
Matrix: 5% Nitric Acid

This ASSURANCE® certified reference material, CRM, is intended primarily for use as a calibration standard or quality control standard for inorganic spectroscopic instrumentation such as ICPOES, DCP, AA, ICPMS, and XRF. It can be employed in USEPA, ASTM and other methods relevant to the certified properties listed below.

The CRM is prepared from high purity single element concentrates of individual elements using Class A laboratory ware to give precise concentration.

Refer to side 2 for details of measurement uncertainties.

Instrumental Analysis by ICP Spectrometer:

Element	Labeled (mg/L)	Measured (mg/L)	NIST SRM
Ca	5,000	4,984.92	3109a
K	5,000	4,990.26	3141a
Mg	5,000	4,991.82	3131a
Na	5,000	4,998.07	3152a

Spex Reference Multi: Lot #10-100AS, 12-113AS, 5-198VY, 6-28VY-REF

Balances are calibrated regularly with weight sets traceable to NIST#s 32856, 32867 and others. This CRM is guaranteed stable and accurate to +/- 0.5% on the average of all the certified concentrations with no single element exceeding +/-2%. This includes uncertainty of measurements and other effects, such as transpiration losses. This guarantee is valid for a period of one year from the date of certification only when the material is kept tightly capped and transported and stored under laboratory conditions.

Date of Certification: SEP - 2003, Certifying Officer: N. Kocherakota

INORGANIC LABS/RADCHEM LABS
DATE RECEIVED: 09/15/03
DATE EXP. DATE: 09/15/2004
DATE OPENED: 09/16/03
INORG: 4945 PO: E50008



inorganic ventures / iv labs

195 lehigh avenue, suite 4, lakewood, nj 08701 usa
phone: 800-669-6799 • 732-901-1900 • fax: 732-901-1903
e-mail: ivsales@ivstandards.com • website: www.ivstandards.com

010149

certificate of analysis

CUSTOM-GRADE SOLUTION

1000 µg/mL Scandium IN 5% HNO₃ (abs)

Catalog Number: CGSC1-1and CGSC1-5

Lot Number: **W-SC02055**

Starting Material:

Starting Material Purity:

Starting Material Lot No:

Sc₂O₃
99.99%
632-5721

INORGANIC LABS/RADCHEM LABS

DATE RECEIVED: 09/04/03

DATE EXPIRED: 10/01/2004

DATE OPENED: 09/24/03

INORG: 4201 PO: F52232

CERTIFIED CONCENTRATION: 1000 ± 5 µg/mL

The Certified Value is based upon the most precise method used to analyze this CRM. The following equations are used in the calculation of the certified value and the uncertainty:

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

$$\text{Uncertainty } (\pm) = \frac{2(\sum s_i^2)^{1/2}}{(n)^{1/2}}$$

($\bar{x}$) = mean

x_i = individual results

n = number of measurements

$\sum s_i$ = The summation of all significant estimated errors.

Classical Wet Assay: 1000 ± 5 µg/mL

Method: EDTA Titration vs NIST SRM 928 Lead Nitrate.

Instrument Analysis: 995 ± 5 µg/mL

Method: Inductively Coupled Plasma Spectroscopy (ICP) vs NIST SRM 3148a.

The independent samples t-test was used to determine if there is agreement between the above assay methods at the 95% confidence interval. Both methods were compared and showed agreement within the stated uncertainties. This agreement is a confirmation of the accuracy of this CRM.

TRACE METALLIC IMPURITIES DETERMINED BY ICP-MS AND ICP-OES IN µg/mL:

Custom-Grade solutions tested for trace metallic impurities by ICP-MS were analyzed in an ULPA-Filtered Clean Room. An ULPA-Filter is 99.9985% efficient for the removal of particles down to 0.3 µm.

<u>Q</u>	Al	<0.070	<u>M</u>	Dy	<0.00060	<u>M</u>	Li	<0.0010	<u>M</u>	Pr	<0.000030	<u>M</u>	Te	<0.0030
<u>M</u>	Sb	<0.000050	<u>M</u>	Er	<0.00050	<u>M</u>	Lu	<0.000040	<u>M</u>	Re	<0.00010	<u>M</u>	Tb	<0.000030
<u>M</u>	As	<0.0010	<u>M</u>	Eu	<0.00030	<u>M</u>	Mg	<0.0030	<u>M</u>	Rh	<0.00010	<u>M</u>	Ti	<0.00010
<u>M</u>	Ba	<0.0010	<u>M</u>	Gd	<0.00010	<u>M</u>	Mn	<0.00040	<u>M</u>	Rb	<0.00010	<u>M</u>	Th	0.0028
<u>M</u>	Be	<0.000050	<u>M</u>	Ga	<0.00010	<u>Q</u>	Hg	i	<u>M</u>	Ru	<0.00020	<u>M</u>	Tm	<0.000040
<u>M</u>	Bi	0.0043	<u>M</u>	Ge	<0.00060	<u>M</u>	Mo	<0.00020	<u>M</u>	Sm	<0.00010	<u>M</u>	Sn	<0.00050
<u>Q</u>	B	<0.020	<u>M</u>	Au	<0.00030	<u>M</u>	Nd	<0.00020	<u>s</u>	Sc		<u>M</u>	Tl	<0.0050
<u>M</u>	Cd	<0.00030	<u>M</u>	Hf	0.0030	<u>Q</u>	Ni	<0.084	<u>Q</u>	Se	<0.67	<u>M</u>	W	<0.0010
<u>Q</u>	Ca	0.016	<u>M</u>	Ho	<0.000050	<u>Q</u>	Nb	<0.000050	<u>Q</u>	Si	<0.034	<u>M</u>	U	<0.00020
<u>M</u>	Ce	<0.00050	<u>M</u>	In	<0.00010	<u>n</u>	Os		<u>M</u>	Ag	0.0050	<u>M</u>	V	<0.00020
<u>M</u>	Cs	<0.000030	<u>M</u>	Ir	<0.00050	<u>M</u>	Pd	<0.00050	<u>Q</u>	Na	<0.16	<u>M</u>	Yb	<0.00010
<u>M</u>	Cr	<0.00050	<u>Q</u>	Fe	<0.16	<u>Q</u>	P	i	<u>M</u>	Sr	<0.000050	<u>M</u>	Y	<0.0040
<u>M</u>	Co	<0.00030	<u>M</u>	La	<0.000050	<u>M</u>	Pt	<0.00020	<u>n</u>	S		<u>M</u>	Zn	0.0075
<u>M</u>	Cu	<0.00060	<u>M</u>	Pb	0.00050	<u>Q</u>	K	<5.01	<u>M</u>	Ta	<0.00070	<u>M</u>	Zr	0.032

M - checked by ICP-MS

O - checked by ICP-OES

i - spectral interference

n - not checked for

s - solution standard element

ANALYZED DENSITY OF SOLUTION (measured at 22°C): 1.038 g/mL


QA:KL Rev.031303/JTS

Paul R. Gaines

Quality Assurance Manager

Expires:

EXPIRES
01/22/04



inorganic ventures / iv labs

195 lehigh avenue, suite 4, lakewood, nj 08701 usa
phone: 800-669-6799 • 732-901-1900 • fax: 732-901-1903
e-mail: ivsales@ivstandards.com • website: www.ivstandards.com

certificate of analysis

010150

CUSTOM-GRADE SOLUTION

Catalog Number: CGSC10-1 and CGSC10-5

10,000 µg/mL Scandium IN 5% HNO₃ (abs)

Lot Number: T-SC02053

Starting Material:
Starting Material Purity:
Starting Material Lot No:

Sc₂O₃
99.999%
632-5721

INORGANIC LABS/RADCHEM LABS

DATE RECEIVED: 09/24/03
DATE EXPIRED: 10/01/2004
DATE OPENED: 09/24/03
INORG: 4262 PO: F50030

CERTIFIED CONCENTRATION: 10,047 ± 29 µg/mL

The Certified Value is based upon the most precise method used to analyze this CRM. The following equations are used in the calculation of the certified value and the uncertainty:

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

$$\text{Uncertainty } (\pm) = \frac{2[(\sum s_i)^2]^{1/2}}{(n)^{1/2}}$$

($\bar{x}$) = mean

x_i = individual results

n = number of measurements

$\sum S_i$ = The summation of all significant estimated errors.

Classical Wet Assay: 10,047 ± 29 µg/mL

Method: EDTA Titration vs NIST SRM 928 Lead Nitrate.

Instrument Analysis: 9994 ± 41 µg/mL

Method: Inductively Coupled Plasma Spectroscopy (ICP) vs NIST SRM 3148a.

The independent samples t-test was used to determine if there is agreement between the above assay methods at the 95% confidence interval. Both methods were compared and showed agreement within the stated uncertainties. This agreement is a confirmation of the accuracy of this CRM.

TRACE METALLIC IMPURITIES DETERMINED BY ICP-MS AND ICP-OES IN µg/mL:

Custom-Grade solutions tested for trace metallic impurities by ICP-MS were analyzed in an ULPA-Filtered Clean Room.
An ULPA-Filter is 99.9985% efficient for the removal of particles down to 0.3 µm.

<u>Q</u>	Al	<0.070	<u>M</u>	Dy	<0.0060	<u>M</u>	Li	<0.010	<u>M</u>	Pr	<0.00030	<u>M</u>	Te	<0.030
<u>M</u>	Sb	<0.00050	<u>M</u>	Er	<0.0050	<u>M</u>	Lu	<0.00040	<u>M</u>	Re	<0.0010	<u>M</u>	Tb	<0.00030
<u>M</u>	As	<0.010	<u>M</u>	Eu	<0.0030	<u>M</u>	Mg	<0.030	<u>M</u>	Rh	<0.0010	<u>M</u>	Ti	<0.0010
<u>M</u>	Ba	<0.010	<u>M</u>	Gd	<0.0010	<u>M</u>	Mn	<0.0040	<u>M</u>	Rb	<0.0010	<u>M</u>	Th	0.028
<u>M</u>	Be	<0.00050	<u>M</u>	Ga	<0.0010	<u>i</u>	Hg		<u>M</u>	Ru	<0.0020	<u>M</u>	Tm	<0.00040
<u>M</u>	Bi	0.043	<u>M</u>	Ge	<0.0060	<u>M</u>	Mo	<0.0020	<u>M</u>	Sm	<0.0010	<u>M</u>	Sn	<0.0050
<u>Q</u>	B	<0.034	<u>M</u>	Au	<0.0030	<u>M</u>	Nd	<0.0020	<u>s</u>	Sc		<u>n</u>	Ti	
<u>M</u>	Cd	<0.0030	<u>M</u>	Hf	0.030	<u>Q</u>	Ni	<0.084	<u>Q</u>	Se	<0.67	<u>M</u>	W	<0.010
<u>Q</u>	Ca	0.17	<u>M</u>	Ho	<0.00050	<u>M</u>	Nb	<0.00050	<u>Q</u>	Si	<0.034	<u>M</u>	U	<0.0020
<u>M</u>	Ce	<0.0050	<u>M</u>	In	<0.0010	<u>n</u>	Os		<u>M</u>	Ag	0.0050	<u>M</u>	V	<0.0020
<u>M</u>	Cs	<0.00030	<u>M</u>	Ir	<0.0050	<u>n</u>	Pd	<0.0050	<u>Q</u>	Na	<0.16	<u>M</u>	Yb	<0.0010
<u>M</u>	Cr	<0.0050	<u>Q</u>	Fe	<0.16	<u>i</u>	P		<u>M</u>	Sr	<0.00050	<u>M</u>	Y	<0.040
<u>M</u>	Co	<0.0030	<u>M</u>	La	<0.00050	<u>M</u>	Pt	<0.0020	<u>n</u>	S		<u>M</u>	Zn	0.075
<u>M</u>	Cu	<0.0060	<u>M</u>	Pb	0.0050	<u>Q</u>	K	<5.01	<u>M</u>	Ta	<0.0070	<u>M</u>	Zr	0.32

M - checked by ICP-MS

O - checked by ICP-OES

i - spectral interference

n - not checked for

s - solution standard element

ANALYZED DENSITY OF SOLUTION (measured at 22°C): 1.073 g/mL

(over)

QA:KL Rev.0212030N

Paul R. Gaines

Quality Assurance Manager

Expires:

EXPIRES

01/2004

SPEXTM Certificate

Certificate of Reference Material

010151

INORGANIC LABS/RADCHEM LABS
DATE RECEIVED: 06/09/03
DATE EXP. 06/30/04
DATE OPENED: 06/03/03
INORG: 4139 PD: F52370

Catalog Number: PLB9-2X/2Y/2T**Lot No.** 9-143B**Description:** 1000 mg/L Boron**Matrix:** H₂O

This ASSURANCE[®] certified reference material, CRM, is intended primarily for use as a calibration standard or quality control standard for inorganic spectroscopic instrumentation such as ICPOES, DCP, AA, ICPMS, and XRF. It can be employed in USEPA, ASTM and other methods relevant to the certified properties listed below.

Certified Value: 999.5 mg/L**Uncertainty Associated with Measurement:** +/-3.0mg/L**Certified Value is Traceable to:** NIST SRM 3107

The CRM is prepared gravimetrically using high purity (NH₄)₂(B₄O₇)-4H₂O Lot# 08001E. The certified value listed is the average of values obtained by classical wet assay and ICP spectrometer analysis

Refer to side 2 for details of measurement uncertainties.

Classical Wet Assay: 1001 mg/L

Method: Titration with Sodium Hydroxide using Phenolphthalein as indicator. Sodium Hydroxide standardized against Potassium Biphthalate NIST SRM #84K.

Instrumentation Analysis By ICP spectrometer: 998 mg/L**Uncertified Properties:****Density:** 1.000 @ 23.7 Degrees Celsius**Trace Metallic Impurities in the Actual Solution via ICP / ICPMS Analysis:**

Element	mg/L	Element	mg/L	Element	mg/L
Al	0.06	Cu	<0.001	Pb	<0.001
As	<0.001	Fe	0.001	Re	<0.001
Ag	<0.002	Ga	<0.001	Rb	<0.001
Ba	<0.001	In	<0.001	Sr	<0.001
Be	<0.001	K	<0.01	Sb	<0.001
Bi	<0.001	Li	<0.001	Sn	<0.001
Cd	<0.001	Mn	<0.001	Ti	<0.001
Co	<0.001	Mo	<0.001	Tl	<0.001
Ca	0.003	Mg	<0.001	V	0.004
Cr	<0.003	Na	0.026	Zr	<0.001
		Ni	<0.001	Zn	<0.007

Balances are calibrated regularly with weight sets traceable to NIST #32856, #32857 and others. This CRM is guaranteed stable to +/-0.5% of the certified concentration inclusive of uncertainty of measurements and other effects, such as transpiration losses, for a period of one year from the date of certification. This guarantee is valid only when the material is kept tightly capped and transported and stored under laboratory conditions.

Date of Certification: JUN '03

Certifying Officer: N. Kocherakota

SPEXTM Certificate

Certificate of Reference Material

010152

Catalog Number: PLLI2-2X/2Y

Lot No. 9-102LI

Description: 1000 mg/L Lithium

Matrix: 2% HNO₃

This ASSURANCE[®] certified reference material, CRM, is intended primarily for use as a calibration standard or quality control standard for inorganic spectroscopic instrumentation such as ICPOES, DCP, AA, ICPMS, and XRF. It can be employed in USEPA, ASTM and other methods relevant to the certified properties listed below.

Certified Value: 1002 mg/L

Uncertainty Associated with Measurement: +/- 3 mg/L

Certified Value is Traceable to: NIST SRM 3129a

The CRM is prepared gravimetrically using high purity Lithium Carbonate Lot# 03021A. The certified value listed is the average of values obtained by classical wet assay and ICP spectrometer analysis

Refer to side 2 for details of measurement uncertainties.

Classical Wet Assay: 1002 mg/L

Method: Evaporate to dryness. Fume with Sulfuric Acid. Ignite and weigh as Li₂SO₄.

Instrumentation Analysis By ICP spectrometer: 1002 mg/L

Uncertified Properties:

Density: 1.014 @ 24.5 Degrees Celsius

Trace Metallic Impurities in the Actual Solution via ICP / ICPMS Analysis:

Element	mg/L	Element	mg/L	Element	mg/L
Al	0.009	Cu	0.01	Pb	0.001
As	0.004	Fe	0.045	Re	0.002
Ag	<0.001	Ga	<0.001	Rb	<0.001
B	<0.004	In	<0.001	Sr	0.001
Ba	0.004	K	0.05	Sb	0.002
Be	0.002	Mn	<0.001	Sn	<0.001
Bi	<0.001	Mo	0.002	Ti	0.005
Ca	0.03	Mg	0.003	Tl	<0.001
Cr	<0.001	Na	0.02	V	<0.001
Cd	<0.001	Ni	<0.001	Zr	<0.001
Co	0.002			Zn	0.20

Balances are calibrated regularly with weight sets traceable to NIST #32856, #32857 and others. This CRM is guaranteed stable to +/-0.5% of the certified concentration inclusive of uncertainty of measurements and other effects, such as transpiration losses, for a period of one year from the date of certification. This guarantee is valid only when the material is kept tightly capped and transported and stored under laboratory conditions.

Date of Certification: MAY '03

Certifying Officer: N. Kocherakota

INORGANIC LABS/RADCHEM LABS
 DATE RECEIVED: 05/12/03
 DATE EXPIRED: 05/15/2004
 DATE OPENED: 05/12/03
 INORG: 4065
 PO: F52181

SPExertificate™

Certificate of Reference Material

010153

INORGANIC LABS/RADCHEM LABS
DATE RECEIVED: 05/12/03
DATE EXPIRED: 05/12/04
DATE OPENED: 05/12/03
INORG: 4066 PO: F52181

Catalog Number: PLMO9-2X/2Y/2T **Lot No.** 9-49MO
Description: 1000 mg/L Molybdenum
Matrix: H₂O

This ASSURANCE® certified reference material, CRM, is intended primarily for use as a calibration standard or quality control standard for inorganic spectroscopic instrumentation such as ICPOES, DCP, AA, ICPMS, and XRF. It can be employed in USEPA, ASTM and other methods relevant to the certified properties listed below.

Certified Value: 1000 mg/L
Uncertainty Associated with Measurement: +/- 3 mg/L
Certified Value is Traceable to: NIST SRM 3134

The CRM is prepared gravimetrically using high purity Ammonium Molybdate Lot# 03011C. The certified value listed is the average of values obtained by classical wet assay and ICP spectrometer analysis. Refer to side 2 for details of measurement uncertainties.

Classical Wet Assay: 999 mg/L

Method: Precipitation using 8-Hydroxyquinoline. Filter, dry, and weigh as MoO₂(C₉H₆NO)₂

Instrumentation Analysis By ICP spectrometer: 1001 mg/L

Uncertified Properties:

Density: 0.9987 @ 24.6 Degrees Celsius

Trace Metallic Impurities in the Actual Solution via ICP / ICPMS Analysis:

Element	mg/L	Element	mg/L	Element	mg/L
Al	0.001	Cu	0.004	Pb	0.001
As	0.04	Fe	<0.01	Re	0.03
Ag	<0.001	Ga	<0.001	Rb	<0.001
B	<0.003	In	<0.001	Sr	<0.001
Ba	<0.001	K	0.02	Sb	<0.001
Be	<0.001	Li	<0.001	Sn	<0.001
Bi	<0.001	Mg	0.007	Ti	0.004
Ca	0.006	Mn	<0.001	Tl	<0.001
Cr	<0.005	Na	0.009	V	0.004
Cd	<0.05	Ni	<0.001	Zr	<0.001
Co	<0.001			Zn	<0.001

Balances are calibrated regularly with weight sets traceable to NIST #32856, #32857 and others. This CRM is guaranteed stable to +/-0.5% of the certified concentration inclusive of uncertainty of measurements and other effects, such as transpiration losses, for a period of one year from the date of certification. This guarantee is valid only when the material is kept tightly capped and transported and stored under laboratory conditions.

Date of Certification: MAY '03 **Certifying Officer:** N. Kocherlakota

SPEXTM Certificate

Certificate of Reference Material

010154**Catalog Number:** PLP9-2X/2Y/2T**Lot No.** 9-150P**Description:** 1000 mg/L Phosphorus**Matrix:** H₂O

This ASSURANCE[®] certified reference material, CRM, is intended primarily for use as a calibration standard or quality control standard for inorganic spectroscopic instrumentation such as ICPOES, DCP, AA, ICPMS, and XRF. It can be employed in USEPA, ASTM and other methods relevant to the certified properties listed below.

Certified Value: 1002.5 mg/L**Uncertainty Associated with Measurement:** +/- 3 mg/L**Certified Value is Traceable to:** NIST SRM 3139a

The CRM is prepared gravimetrically using high purity (NH₄)H₂(PO₄) Lot# W1002B. The certified value listed is the average of values obtained by classical wet assay and ICP spectrometer analysis.

Refer to side 2 for details of measurement uncertainties.

Classical Wet Assay: 1003 mg/L**Method:** Precipitation using Magnesia Mixture. Filter, ignite, and weigh as Mg₂P₂O₇.**Instrumentation Analysis By ICP spectrometer:** 1002 mg/L**Uncertified Properties:****Density:** 0.9996 @ 24.0 Degrees Celsius**Trace Metallic Impurities in the Actual Solution via ICP / ICPMS Analysis:**

Element	mg/L	Element	mg/L	Element	mg/L
Al	0.002	Cu	<0.001	Pb	<0.001
As	0.001	Fe	<0.001	Rb	<0.001
Ag	<0.002	Ga	<0.001	Re	<0.001
B	<0.002	In	<0.001	Sn	<0.001
Ba	<0.001	K	0.006	Sr	<0.001
Be	<0.001	Li	<0.001	Sb	0.004
Bi	<0.001	Mg	<0.001	Ti	0.004
Ca	0.004	Mn	<0.001	Tl	<0.001
Cr	<0.008	Mo	<0.001	V	<0.006
Cd	<0.001	Na	0.003	Zr	<0.001
Co	<0.001	Ni	<0.001	Zn	0.07

Balances are calibrated regularly with weight sets traceable to NIST #32856, #32857 and others. This CRM is guaranteed stable to +/-0.5% of the certified concentration inclusive of uncertainty of measurements and other effects, such as transpiration losses, for a period of one year from the date of certification. This guarantee is valid only when the material is kept tightly capped and transported and stored under laboratory conditions.

Date of Certification: OCT -- 2003**Certifying Officer:** N. Kocherakota

INORGANIC LABS/RADCHEM LABS
DATE RECEIVED: 10/31/03
DATE EXPIRED: 10/31/04
DATE OPENED: 11/30/03
INORG: 4305 PO: F53261

SPExertificate™

Certificate of Reference Material

010153**Catalog Number:** PLSI9-2X/2Y/2T**Lot No.** 10-07SI**Description:** 1000 mg/L Silicon**Matrix:** H₂O / 0.4% F-

This ASSURANCE® certified reference material, CRM, is intended primarily for use as a calibration standard or quality control standard for inorganic spectroscopic instrumentation such as ICPOES, DCP, AA, ICPMS, and XRF. It can be employed in USEPA, ASTM and other methods relevant to the certified properties listed below.

Certified Value: 998.5 mg/L**Uncertainty Associated with Measurement:** +/- 3 mg/L**Certified Value is Traceable to:** NIST SRM #3150

The CRM is prepared gravimetrically using high purity (NH₄)₂SiF₆ Lot# 02021D. The certified value listed is the average of values obtained by classical wet assay and ICP spectrometer analysis.

Refer to side 2 for details of measurement uncertainties.

Classical Wet Assay: 997 mg/L

Method: Precipitation using Ammonium Molybdate and 8-Hydroxy Quinoline. Filter, dry, and weigh as (C₉H₇ON)₄(H₄)[Si(Mo₁₂O₄₀)]

Instrumentation Analysis By ICP spectrometer: 1000 mg/L**Uncertified Properties:****Density:** 1.010 @ 26.5 Degrees Celsius**Trace Metallic Impurities in the Actual Solution via ICP / ICPMS Analysis:**

Element	mg/L	Element	mg/L	Element	mg/L
Al	0.002	Cu	<0.001	Pb	<0.001
As	<0.001	Fe	0.020	Rb	<0.001
Ag	<0.001	Ga	<0.001	Re	<0.001
B	<0.003	In	<0.001	Sr	<0.001
Ba	<0.001	K	<0.010	Sb	0.03
Be	<0.001	Li	<0.001	Ti	<0.001
Bi	<0.001	Mg	<0.001	Tl	<0.001
Ca	0.018	Mn	<0.001	V	<0.001
Cr	<0.002	Mo	<0.001	Zr	0.05
Cd	<0.001	Na	0.02	Zn	0.06
Co	<0.001	Ni	0.001		

Balances are calibrated regularly with weight sets traceable to NIST #32856, #32857 and others. This CRM is guaranteed stable to +/-0.5% of the certified concentration inclusive of uncertainty of measurements and other effects, such as transpiration losses, for a period of one year from the date of certification. This guarantee is valid only when the material is kept tightly capped and transported and stored under laboratory conditions.

Date of Certification: AUG '03 Certifying Officer: N. Kocherlakota

INORGANIC LABS/RADCHEM LABS
DATE RECEIVED: 09/05/03
DATE EXPIRED: 08/30/2004
DATE OPENED: 09/05/03
INORG: 4232 PO: 552225

*SPEx*ertificate™

Certificate of Reference Material

010156**Catalog Number:** PLTI9-2X/2Y/2T**Lot No.** 10-38TI**Description:** 1000 mg/L Titanium**Matrix:** H₂O/ 0.24% F-

This ASSURANCE® certified reference material, CRM, is intended primarily for use as a calibration standard or quality control standard for inorganic spectroscopic instrumentation such as ICPOES, DCP, AA, ICPMS, and XRF. It can be employed in USEPA, ASTM and other methods relevant to the certified properties listed below.

Certified Value: 1001 mg/L**Uncertainty Associated with Measurement:** 3.0mg/L**Certified Value is Traceable to:** NIST SRM #3162a

The CRM is prepared gravimetrically using high purity (NH₄)₂TiF₆ Lot# 02021E. The certified value listed is the average of values obtained by classical wet assay and ICP spectrometer analysis.

Refer to side 2 for details of measurement uncertainties.

Classical Wet Assay: 1003 mg/L**Method:** Precipitation using Ammonium Hydroxide. Filter, ignite, and weigh as TiO₂.**Instrumentation Analysis By ICP spectrometer:** 999 mg/L**Uncertified Properties:****Density:** 1.001 @ 22.5 Degrees Celsius**Trace Metallic Impurities in the Actual Solution via ICP / ICPMS Analysis:**

Element	mg/L	Element	mg/L	Element	mg/L
Al	0.006	Cu	<0.10	Pb	<0.001
As	<0.001	Fe	<0.01	Rb	<0.001
Ag	<0.001	Ga	<0.001	Re	<0.001
B	0.003	In	<0.001	Si	0.52
Ba	<0.001	K	<0.01	Sr	0.001
Be	<0.001	Li	<0.001	Sb	<0.001
Bi	<0.001	Mg	<0.001	Tl	<0.001
Ca	0.013	Mn	<0.001	V	<0.001
Cr	<0.003	Mo	<0.001	Zr	0.01
Cd	<0.001	Na	0.02	Zn	0.03
Co	0.002	Ni	<0.001		

Balances are calibrated regularly with weight sets traceable to NIST #32856, #32857 and others. This CRM is guaranteed stable to +/-0.5% of the certified concentration inclusive of uncertainty of measurements and other effects, such as transpiration losses, for a period of one year from the date of certification. This guarantee is valid only when the material is kept tightly capped and transported and stored under laboratory conditions.

Date of Certification: AUG '03**Certifying Officer:** N. Kocherakota

INORGANIC LABS/RADCHEM LABS
DATE RECEIVED: 09/05/03
DATE EXPIRED: 08/30/2004
DATE OPENED: 09/05/03
INORG: 4204 PO: 550005

SPEX CertificateTM

Certificate of Reference Material

010157

Catalog Number: PLSR2-2X/2Y/2T **Lot No.** 9-166SR
Description: 1000 mg/L Strontium in 2% HNO₃
Matrix: 2% HNO₃

This ASSURANCE[®] certified reference material, CRM, is intended primarily for use as a calibration standard or quality control standard for inorganic spectroscopic instrumentation such as ICPOES, DCP, AA, ICPMS, and XRF. It can be employed in USEPA, ASTM and other methods relevant to the certified properties listed below.

Certified Value: 1002.5 mg/L
Uncertainty Associated with Measurement: +/- 3 mg/L
Certified Value is Traceable to: NIST SRM 3153a

The CRM is prepared gravimetrically using high purity Strontium Carbonate Lot# 02001B. The certified value listed is the average of values obtained by classical wet assay and ICP spectrometer analysis. Refer to side 2 for details of measurement uncertainties.

Classical Wet Assay: 1002 mg/L

Method: EDTA titration using Methyl Thymol Blue as indicator. EDTA standardized against Pb(NO₃)₂ NIST SRM #928.

Instrumentation Analysis By ICP spectrometer: 1003 mg/L

Uncertified Properties:

Density: 1.010 @ 22.7 Degrees Celsius

Trace Metallic Impurities in the Actual Solution via ICP / ICPMS Analysis:

Element	mg/L	Element	mg/L	Element	mg/L
Al	0.02	Cu	<0.001	Pb	<0.001
As	<0.001	Fe	0.001	Rb	<0.001
Ag	<0.002	Ga	<0.001	Re	<0.001
B	<0.003	In	<0.001	Si	0.043
Ba	0.008	K	0.10	Sb	<0.001
Be	<0.001	Li	0.007	Ti	<0.002
Bi	<0.001	Mg	<0.003	Tl	<0.001
Ca	0.014	Mn	<0.001	V	<0.001
Cr	0.001	Mo	<0.001	Zr	<0.001
Cd	<0.001	Na	0.01	Zn	0.04
Co	<0.001	Ni	<0.001		

Balances are calibrated regularly with weight sets traceable to NIST #32856, #32857 and others. This CRM is guaranteed stable to +/-0.5% of the certified concentration inclusive of uncertainty of measurements and other effects, such as transpiration losses, for a period of one year from the date of certification. This guarantee is valid only when the material is kept tightly capped and transported and stored under laboratory conditions.

Date of Certification: OCT -- 2003

Certifying Officer: N. Kocherakota

INORGANIC LABS/RADCHEM LABS
DATE RECEIVED: 10/31/03
DATE EXPIRED: 10/31/04
DATE OPENED: 11/3/03
INORG: 4308 PO: F50061

SPEXTM Certificate

Certificate of Reference Material

010158

INORGANIC LABS/REACHCHEM LABS
 DATE RECEIVED: 06/30/03
 DATE EXPIRED: 06/30/2004
 DATE OPENED: 06/23/03
 INRG: 411 PD: F52370

Catalog Number: PLSN5-2X/2Y/2T

Lot No. 9-62SN

Description: 1000 mg/L Tin

Matrix: 20% HCL

This ASSURANCE[®] certified reference material, CRM, is intended primarily for use as a calibration standard or quality control standard for inorganic spectroscopic instrumentation such as ICPOES, DCP, AA, ICPMS, and XRF. It can be employed in USEPA, ASTM and other methods relevant to the certified properties listed below.

Certified Value: 1000 mg/L

Uncertainty Associated with Measurement: +/- 3 mg/L

Certified Value is Traceable to: NIST SRM 3161a

The CRM is prepared gravimetrically using high purity Tin Metal Lot# 02891N. The certified value listed is the average of values obtained by classical wet assay and ICP spectrometer analysis.

Refer to side 2 for details of measurement uncertainties.

Classical Wet Assay: 999 mg/L

Method: Precipitation using Ammonium Hydroxide. Filter, ignite, and weigh as SnO₂.

Instrumentation Analysis By ICP spectrometer: 1001 mg/L

Uncertified Properties:

Density: 1.034 @ 24.8 Degrees Celsius

Trace Metallic Impurities in the Actual Solution via ICP / ICPMS Analysis:

Element	mg/L	Element	mg/L	Element	mg/L
Al	0.007	Cu	<0.001	Pb	0.001
As	0.01	Fe	0.02	Rb	<0.001
Ag	0.002	Ga	<0.001	Re	<0.001
B	<0.03	In	<0.001	Sr	<0.001
Ba	<0.001	K	0.10	Sb	0.002
Be	<0.001	Li	<0.001	Ti	<0.001
Bi	<0.001	Mg	<0.001	Tl	<0.001
Ca	0.004	Mn	<0.001	V	<0.001
Cr	<0.005	Mo	<0.001	Zr	<0.001
Cd	<0.001	Na	0.02	Zn	0.03
Co	0.008	Ni	<0.01		

Balances are calibrated regularly with weight sets traceable to NIST #32856, #32857 and others. This CRM is guaranteed stable to +/-0.5% of the certified concentration inclusive of uncertainty of measurements and other effects, such as transpiration losses, for a period of one year from the date of certification. This guarantee is valid only when the material is kept tightly capped and transported and stored under laboratory conditions.

Date of Certification: JUN 03

Certifying Officer: N. Kocherakota

SPEXTM Certificate

Certificate of Reference Material

010159

INORGANIC LABS/RADCHEM LABS
DATE RECEIVED: 06/29/03
DATE EXPIRED: 06/30/04
DATE OPENED: 06/23/03
INORG: 4138 PO: F52370

Catalog Number: PLBI4-2X/2Y**Lot No.** 9-36BI**Description:** 1000 mg/L Bismuth**Matrix:** 10% HNO₃

This ASSURANCE[®] certified reference material, CRM, is intended primarily for use as a calibration standard or quality control standard for inorganic spectroscopic instrumentation such as ICPOES, DCP, AA, ICPMS, and XRF. It can be employed in USEPA, ASTM and other methods relevant to the certified properties listed below.

Certified Value: 1001 mg/L**Uncertainty Associated with Measurement:** +/- 3 mg/L**Certified Value is Traceable to:** NIST SRM 3106

The CRM is prepared gravimetrically using high purity Bismuth Metal Lot# 04941B. The certified value listed is the average of values obtained by classical wet assay and ICP spectrometer analysis.

Refer to side 2 for details of measurement uncertainties.

Classical Wet Assay: 1001 mg/L

Method: EDTA titration using Xylenol Orange as indicator. EDTA standardized against Pb(NO₃)₂ NIST SRM #928.

Instrumentation Analysis By ICP spectrometer: 1001 mg/L**Uncertified Properties:****Density:** 1.052 @ 23.1 Degrees Celsius**Trace Metallic Impurities in the Actual Solution via ICP / ICPMS Analysis:**

Element	mg/L	Element	mg/L	Element	mg/L
Al	0.002	Cu	0.034	Pb	0.006
As	<0.001	Fe	0.001	Re	<0.001
Ag	0.002	Ga	<0.001	Rb	<0.001
B	<0.04	In	<0.001	Sr	<0.001
Ba	<0.001	K	0.001	Sb	<0.001
Be	<0.001	Li	<0.001	Sn	<0.001
Cd	<0.001	Mn	<0.001	Ti	<0.001
Co	<0.001	Mo	<0.001	Tl	<0.001
Ca	0.006	Mg	<0.001	V	<0.002
Cr	<0.005	Na	0.005	Zr	<0.001
		Ni	<0.001	Zn	0.02

Balances are calibrated regularly with weight sets traceable to NIST #32856, #32857 and others. This CRM is guaranteed stable to +/-0.5% of the certified concentration inclusive of uncertainty of measurements and other effects, such as transpiration losses, for a period of one year from the date of certification. This guarantee is valid only when the material is kept tightly capped and transported and stored under laboratory conditions.

Date of Certification: JUN '03**Certifying Officer:** N. Kocherakota

SPEX Certificate TM

Certificate of Reference Material

010160**Catalog Number:** PLLA2-2X/2Y**Lot No.** 10-17LA**Description:** 1000 mg/L Lanthanum**Matrix:** 2% HNO₃

This **ASSURANCE**® certified reference material, CRM, is intended primarily for use as a calibration standard or quality control standard for inorganic spectroscopic instrumentation such as ICPOES, DCP, AA, ICPMS, and XRF. It can be employed in USEPA, ASTM and other methods relevant to the certified properties listed below.

Certified Value: 996.5 mg/L**Uncertainty Associated with Measurement:** +/- 3 mg/L**Certified Value is Traceable to:** NIST SRM 3127a

The CRM is prepared gravimetrically using high purity Lanthanum Oxide Lot# 06981D. The certified value listed is the average of values obtained by classical wet assay and ICP spectrometer analysis

Refer to side 2 for details of measurement uncertainties.

Classical Wet Assay: 998 mg/L

Method: EDTA titration using Methyl Thymol Blue as indicator. EDTA standardized against Pb(NO₃)₂ NIST SRM #928.

Instrumentation Analysis By ICP spectrometer: 995 mg/L**Uncertified Properties:****Density:** 1.010 @ 23.8 Degrees Celsius**Trace Metallic Impurities in the Actual Solution via ICP / ICPMS Analysis:**

Element	mg/L	Element	mg/L	Element	mg/L
Ce	0.009	Lu	<0.001	Tm	<0.001
Ca	0.05	Mn	<0.001	Ti	<0.001
Dy	<0.001	Mo	<0.001	Tb	<0.001
Er	<0.001	Nd	<0.001	Ta	<0.001
Eu	<0.001	Ni	<0.001	Tl	<0.001
Fe	0.006	Na	0.006	V	<0.001
Gd	<0.2	Pr	<0.001	W	<0.001
Ga	<0.001	Rb	<0.001	Y	<0.001
Hf	<0.001	Sc	<0.003	Yb	<0.001
Ho	<0.001	Sm	<0.001	Zr	<0.001
In	<0.001	Th	<0.001		

Balances are calibrated regularly with weight sets traceable to NIST #32856, #32857 and others. This CRM is guaranteed stable to +/-0.5% of the certified concentration inclusive of uncertainty of measurements and other effects, such as transpiration losses, for a period of one year from the date of certification. This guarantee is valid only when the material is kept tightly capped and transported and stored under laboratory conditions.

Date of Certification: MAY '03 Certifying Officer: N. Kocherakota

INORGANIC LABS/KADCHEM LABS
DATE RECEIVED: 05/12/03
DATE EXP. 05/15/2004
DATE OPENED: 05/12/03
INDRG: 4064 PO: F52181

SPEXTM Certificate

Certificate of Reference Material

010161

Catalog Number: PLY2-2X/2Y/2T

Lot No. 9-152Y

Description: 1,000 mg/L Yttrium

Matrix: 2% HNO₃

This ASSURANCE[®] certified reference material, CRM, is intended primarily for use as a calibration standard or quality control standard for inorganic spectroscopic instrumentation such as ICPOES, DCP, AA, ICPMS, and XRF. It can be employed in USEPA, ASTM and other methods relevant to the certified properties listed below.

Certified Value: 1001.5 mg/L

Uncertainty Associated with Measurement: +/- 3 mg/L

Certified Value is Traceable to: NIST SRM 3167a.

The CRM is prepared gravimetrically using high purity Yttrium Oxide Lot# 08001A. The certified value listed is the average of values obtained by classical wet assay and ICP spectrometer analysis.

Refer to side 2 for details of measurement uncertainties.

Classical Wet Assay: 1002 mg/L

Method: EDTA titration using Methyl Thymol Blue as indicator. EDTA standardized against Pb(NO₃)₂ NIST SRM #928.

Instrumentation Analysis By ICP spectrometer: 1001 mg/L

Uncertified Properties:

Density: 1.010 @ 24.8 Degrees Celsius

Trace Metallic Impurities in the Actual Solution via ICP / ICPMS Analysis:

Element	mg/L	Element	mg/L	Element	mg/L
Ce	<0.001	La	<0.001	Tb	<0.001
Ca	0.007	Lu	<0.001	Tm	<0.001
Dy	<0.001	Mn	<0.001	Tl	<0.001
Er	<0.001	Mo	<0.001	Th	<0.001
Eu	<0.001	Nd	<0.001	Ta	<0.001
Fe	0.003	Ni	<0.001	Ti	<0.001
Gd	<0.001	Na	0.005	V	<0.001
Ga	<0.001	Pr	<0.001	W	<0.001
Hf	<0.001	Rb	<0.001	Yb	<0.001
Ho	<0.001	Sc	<0.001	Zr	0.003
In	<0.001	Sm	<0.001		

Balances are calibrated regularly with weight sets traceable to NIST #32856, #32857 and others. This CRM is guaranteed stable to +/-0.5% of the certified concentration inclusive of uncertainty of measurements and other effects, such as transpiration losses, for a period of one year from the date of certification. This guarantee is valid only when the material is kept tightly capped and transported and stored under laboratory conditions.

Date of Certification: MAY '03

Certifying Officer: N. Kocherakota

INORGANIC LABS/RADCHEM LABS
 DATE RECEIVED: 05/12/03
 DATE EXPIRED: 05/15/2004
 DATE OPENED: 05/12/03
 INORG: 4067 PD: F52181

SPEXTM Certificate

Certificate of Reference Material

010162

INORGANIC LABS/RADCHEM LABS

DATE RECEIVED: 01/16/03

DATE EXPIRED: 01/15/2004

DATE OPENED: 01/23/03

INORG: 3928 PO: F52013

Catalog Number: PLPD3-2X/2Y

Lot No. 9-74PD

Description: 1000 mg/L Palladium

Matrix: 10% HCl

This ASSURANCE[®] certified reference material, CRM, is intended primarily for use as a calibration standard or quality control standard for inorganic spectroscopic instrumentation such as ICPOES, DCP, AA, ICPMS, and XRF. It can be employed in USEPA, ASTM and other methods relevant to the certified properties listed below.

Certified Value: 995 mg/L

Uncertainty Associated with Measurement: +/- 3 mg/L

Certified Value is Traceable to: NIST SRM 3138

The CRM is prepared gravimetrically using high purity Palladium Powder Lot# 07991E. The certified value listed is the average of values obtained by classical wet assay and ICP spectrometer analysis

Refer to side 2 for details of measurement uncertainties.

Classical Wet Assay: 990 mg/L

Method: Precipitation using dimethyl glyoxime, filter, dry and weigh as Pd(C₄H₇O₂N₂)₂.

Instrumentation Analysis By ICP spectrometer: 1000 mg/L

Uncertified Properties:

Density: 1.025 @ 26.3 Degrees Celsius

Trace Metallic Impurities in the Actual Solution via ICP / ICPMS Analysis:

Element	mg/L	Element	mg/L	Element	mg/L
Al	<0.004	Cr	<0.009	Pt	<0.004
Au	0.03	Fe	0.02	Re	<0.001
Ag	0.007	Ga	<0.001	Rh	0.002
B	<0.05	Ir	<0.001	Rb	<0.001
Be	<0.002	In	<0.001	Ru	<0.001
Bi	<0.001	Mg	<0.005	Sn	0.08
Ca	0.02	Mn	<0.001	Te	<0.001
Cd	<0.001	Na	<0.03	Ti	<0.001
Co	<0.001	Ni	0.004	W	<0.001
Cu	<0.006	Pb	<0.04	Zr	<0.001
				Zn	0.16

Balances are calibrated regularly with weight sets traceable to NIST #32856, #32857 and others. This CRM is guaranteed stable to +/-0.5% of the certified concentration inclusive of uncertainty of measurements and other effects, such as transpiration losses, for a period of one year from the date of certification. This guarantee is valid only when the material is kept tightly capped and transported and stored under laboratory conditions.

Date of Certification: JAN - - '03 **Certifying Officer:** N. Kocherakota

SPExertificate™

Certificate of Reference Material

020163

INORGANIC LABS/RADCHEM LABS
DATE RECEIVED: 06/20/03
DATE EXPIRED: 06/30/2004
DATE OPENED: 06/23/03
INORG: 4140 PO: F52370

Catalog Number: PLS9-2X/2Y/2T**Lot No.** 8-74S**Description:** 1000 mg/L Sulfur**Matrix:** H2O

This ASSURANCE® certified reference material, CRM, is intended primarily for use as a calibration standard or quality control standard for inorganic spectroscopic instrumentation such as ICPOES, DCP, AA, ICPMS, and XRF. It can be employed in USEPA, ASTM and other methods relevant to the certified properties listed below.

Certified Value: 1003 mg/L**Uncertainty Associated with Measurement:** +/- 3 mg/L**Certified Value is Traceable to:** NIST SRM 3154

The CRM is prepared gravimetrically using high purity Ammonium Sulfate Lot# 05891M. The certified value listed is the average of values obtained by classical wet assay and ICP spectrometer analysis.

Refer to side 2 for details of measurement uncertainties.

Classical Wet Assay: 1003 mg/L**Method:** Precipitation using barium chloride, filter, ignite and weigh as BaSO4.**Instrumentation Analysis By ICP spectrometer:** 1003 mg/L**Uncertified Properties:****Density:** 1.007 @ 23.6 Degrees Celsius**Trace Metallic Impurities in the Actual Solution via ICP / ICPMS Analysis:**

Element	mg/L	Element	mg/L	Element	mg/L
Al	<0.001	Cu	<0.001	Pb	0.002
As	<0.001	Fe	0.008	Rb	<0.001
Ag	<0.001	Ga	<0.001	Re	<0.001
B	<0.004	In	<0.001	Sn	<0.001
Ba	<0.001	K	<0.001	Sr	<0.001
Be	<0.001	Li	<0.001	Sb	<0.001
Bi	<0.001	Mg	0.005	Ti	<0.002
Ca	0.009	Mn	<0.001	Tl	<0.001
Cr	<0.004	Mo	<0.001	V	<0.001
Cd	<0.001	Na	0.02	Zr	<0.001
Co	<0.001	Ni	<0.001	Zn	0.0075

Balances are calibrated regularly with weight sets traceable to NIST #32856, #32857 and others. This CRM is guaranteed stable to +/-0.5% of the certified concentration inclusive of uncertainty of measurements and other effects, such as transpiration losses, for a period of one year from the date of certification. This guarantee is valid only when the material is kept tightly capped and transported and stored under laboratory conditions.

Date of Certification: JUN '03**Certifying Officer:** N. Kocherakota

SPEXTM Certificate

Certificate of Reference Material

040164

Catalog Number: PLTH2-2X/2Y

Lot No. 10-24TH

Description: 1000 mg/L Thorium

Matrix: 2% HNO₃

This ASSURANCE[®] certified reference material, CRM, is intended primarily for use as a calibration standard or quality control standard for inorganic spectroscopic instrumentation such as ICPOES, DCP, AA, ICPMS, and XRF. It can be employed in USEPA, ASTM and other methods relevant to the certified properties listed below.

Certified Value: 999 mg/L

Uncertainty Associated with Measurement: +/- 3.0 mg/L

Certified Value is Traceable to: NIST SRM #3159

The CRM is prepared gravimetrically using high purity Th(NO₃)₄·4H₂O Lot# 01851R. The certified value listed is the average of values obtained by classical wet assay and ICP spectrometer analysis

Refer to side 2 for details of measurement uncertainties.

Classical Wet Assay: 1000 mg/L

Method: EDTA titration using Xylenol Orange as indicator. EDTA syandardized against Pb(NO₃)₂ NIST SRM #928.

Instrumentation Analysis By ICP spectrometer: 998 mg/L

Uncertified Properties:

Density: 1.010 @ 22.0 Degrees Celsius

Trace Metallic Impurities in the Actual Solution via ICP / ICPMS Analysis:

Element	mg/L	Element	mg/L	Element	mg/L
Ce	0.01	La	0.003	Tb	<0.001
Ca	0.27	Lu	<0.001	Tm	<0.001
Dy	<0.001	Mn	<0.001	Ti	<0.002
Er	<0.001	Mo	<0.001	Ta	<0.001
Eu	<0.001	Nd	0.003	Tl	<0.001
Fe	<0.01	Ni	<0.001	V	<0.001
Gd	<0.001	Na	0.04	W	<0.001
Ga	<0.001	Pr	<0.001	Y	0.002
Hf	<0.001	Rb	<0.001	Yb	<0.001
Ho	<0.001	Sc	<0.03	Zr	<0.001
In	<0.001	Sm	<0.001		

Balances are calibrated regularly with weight sets traceable to NIST #32856, #32857 and others. This CRM is guaranteed stable to +/-0.5% of the certified concentration inclusive of uncertainty of measurements and other effects, such as transpiration losses, for a period of one year from the date of certification. This guarantee is valid only when the material is kept tightly capped and transported and stored under laboratory conditions.

Date of Certification: AUG '03

Certifying Officer: N. Kocherakota

INORGANIC LABS/RADCHEM LABS
 DATE RECEIVED: 09/05/03
 DATE EXPIRED: 08/30/2004
 DATE OPENED: 09/05/03
 INDRG: 4233
 FO: F52225

SPEXcertificate™

Certificate of Reference Material

010165**Catalog Number:** PLU2-2X/2Y**Lot No.** 9-179U**Description:** 1000 mg/L Uranium**Matrix:** 2% HNO₃

This ASSURANCE® certified reference material, CRM, is intended primarily for use as a calibration standard or quality control standard for inorganic spectroscopic instrumentation such as ICPOES, DCP, AA, ICPMS, and XRF. It can be employed in USEPA, ASTM and other methods relevant to the certified properties listed below.

Certified Value: 999.5 mg/L**Uncertainty Associated with Measurement:** +/- 3 mg/L**Certified Value is Traceable to:** NIST SRM 3164.

The CRM is prepared gravimetrically using high purity Uranium Oxide Lot# 04001D. The certified value listed is the average of values obtained by classical wet assay and ICP spectrometer analysis.

Refer to side 2 for details of measurement uncertainties.**Classical Wet Assay:** 999 mg/L**Method:** Evaporate to dryness. Ignite and weigh as U₃O₈.**Instrumentation Analysis By ICP spectrometer:** 1000 mg/L**Uncertified Properties:****Density:** 1.010 @ 23.6 Degrees Celsius**Trace Metallic Impurities in the Actual Solution via ICP / ICPMS Analysis:**

Element	mg/L	Element	mg/L	Element	mg/L
Al	0.005	Cu	0.02	Pb	0.004
As	0.06	Fe	0.011	Rb	<0.001
Ag	<0.001	Ga	<0.001	Re	<0.001
B	<0.005	In	<0.001	Si	<0.10
Ba	0.004	K	0.008	Sr	0.003
Be	<0.001	Li	<0.001	Sb	0.003
Bi	<0.001	Mg	0.003	Ti	<0.001
Ca	0.012	Mn	0.003	Tl	<0.001
Cr	<0.010	Mo	0.006	V	<0.003
Cd	<0.001	Na	0.10	Zr	<0.001
Co	<0.001	Ni	<0.001	Zn	0.008

Balances are calibrated regularly with weight sets traceable to NIST #32856, #32857 and others. This CRM is guaranteed stable to +/-0.5% of the certified concentration inclusive of uncertainty of measurements and other effects, such as transpiration losses, for a period of one year from the date of certification. This guarantee is valid only when the material is kept tightly capped and transported and stored under laboratory conditions.

Date of Certification: JUN '03**Certifying Officer:** N. Kocherakota

INORGANIC LABS/RADCHEM LABS
DATE RECEIVED: 06/20/03
DATE EXPIRED: 06/30/2004
DATE OPENED: 06/23/03
INORG: 4152 PO: F52370

SPEXTM Certificate

Certificate of Reference Material

010168

Catalog Number: PLW9-2X/2Y

Lot No. 9-177W

Description: 1000 mg/L Tungsten

Matrix: H₂O

This ASSURANCE[®] certified reference material, CRM, is intended primarily for use as a calibration standard or quality control standard for inorganic spectroscopic instrumentation such as ICPOES, DCP, AA, ICPMS, and XRF. It can be employed in USEPA, ASTM and other methods relevant to the certified properties listed below.

Certified Value: 1,000 mg/L

Uncertainty Associated with Measurement: +/- 3 mg/L

Certified Value is Traceable to: NIST SRM 3163

The CRM is prepared gravimetrically using high purity Ammonium Tungstate Lot# 02001H. The certified value listed is the average of values obtained by classical wet assay and ICP spectrometer analysis. Refer to side 2 for details of measurement uncertainties.

Classical Wet Assay: 1000 mg/L

Method: Fume with Sulfuric Acid to dryness. Ignite and weigh as WO₃.

Instrumentation Analysis By ICP spectrometer: 1000 mg/L

Uncertified Properties:

Density: 0.9979 @ 23.7 Degrees Celsius

Trace Metallic Impurities in the Actual Solution via ICP / ICPMS Analysis:

Element	mg/L	Element	mg/L	Element	mg/L
Al	0.002	Cu	<0.001	Pb	<0.001
As	0.01	Fe	<0.01	Rb	<0.001
Ag	<0.003	Ga	<0.001	Re	0.004
B	<0.005	In	<0.001	Si	.56
Ba	<0.001	K	0.05	Sr	<0.001
Be	<0.001	Li	<0.001	Sb	0.001
Bi	<0.001	Mg	<0.001	Ti	<0.001
Ca	0.009	Mn	<0.001	Tl	<0.001
Cr	<0.001	Mo	0.005	V	0.001
Cd	<0.001	Na	0.03	Zr	<0.001
Co	0.001	Ni	<0.001	Zn	0.01

Balances are calibrated regularly with weight sets traceable to NIST #32856, #32857 and others. This CRM is guaranteed stable to +/-0.5% of the certified concentration inclusive of uncertainty of measurements and other effects, such as transpiration losses, for a period of one year from the date of certification. This guarantee is valid only when the material is kept tightly capped and transported and stored under laboratory conditions.

Date of Certification: AUG '03

Certifying Officer: N. Kocherakota

INORGANIC LABS/RADCHEM LABS
DATE RECEIVED: 08/11/03
DATE EXPIRED: 08/15/2004
DATE OPENED: 08/13/03
INORG: 4212 PD: F52218

SPEXTM Certificate

Certificate of Reference Material

010167

Catalog Number: PLZR2-2X/2Y/2T

Lot No. 10-05ZR

Description: 1000 mg/L Zirconium

Matrix: 2% HNO₃

This ASSURANCE[®] certified reference material, CRM, is intended primarily for use as a calibration standard or quality control standard for inorganic spectroscopic instrumentation such as ICPOES, DCP, AA, ICPMS, and XRF. It can be employed in USEPA, ASTM and other methods relevant to the certified properties listed below.

Certified Value: 997 mg/L

Uncertainty Associated with Measurement: +/- 3.0 mg/L

Certified Value is Traceable to: NIST SRM 3169

The CRM is prepared gravimetrically using high purity Zirconyl Nitrate Lot# 11011C. The certified value listed is the average of values obtained by classical wet assay and ICP spectrometer analysis.

Refer to side 2 for details of measurement uncertainties.

Classical Wet Assay: 997 mg/L

Method: Evaporate to dryness. Fume with Sulfuric Acid. Ignite and weigh as ZrO₂.

Instrumentation Analysis By ICP spectrometer: 997 mg/L

Uncertified Properties:

Density: 1.010 @ 23.6 Degrees Celsius

Trace Metallic Impurities in the Actual Solution via ICP / ICPMS Analysis:

Element	mg/L	Element	mg/L	Element	mg/L
Al	0.03	Cu	0.002	Pb	0.002
As	<0.001	Fe	0.017	Rb	<0.001
Ag	<0.05	Ga	<0.001	Re	<0.001
B	<0.004	In	<0.001	Si	0.10
Ba	<0.001	K	0.10	Sr	<0.001
Be	<0.001	Li	0.002	Sb	<0.001
Bi	<0.001	Mg	0.003	Ti	<0.001
Ca	0.11	Mn	<0.001	Tl	<0.001
Cr	<0.009	Mo	<0.001	V	<0.001
Cd	0.004	Na	0.04	Zn	0.02
Co	<0.001	Ni	<0.001		

Balances are calibrated regularly with weight sets traceable to NIST #32856, #32857 and others. This CRM is guaranteed stable to +/-0.5% of the certified concentration inclusive of uncertainty of measurements and other effects, such as transpiration losses, for a period of one year from the date of certification. This guarantee is valid only when the material is kept tightly capped and transported and stored under laboratory conditions.

Date of Certification: AUG 03 **Certifying Officer:** N. Kocherakota

DATE RECEIVED: 08/11/03
DATE EXP. 08/15/2004
DATE OPENED: 08/13/03
INDRG: 4913 PO: F52218

*SPEX*ertificateTM

Certificate of Reference Material

010168**Catalog Number:** PLNA2-3X/3Y**Lot No.** T8-73NA**Description:** 10,000 mg/L Sodium**Matrix:** 5% HNO₃

This ASSURANCE[®] certified reference material, CRM, is intended primarily for use as a calibration standard or quality control standard for inorganic spectroscopic instrumentation such as ICPOES, DCP, AA, ICPMS, and XRF. It can be employed in USEPA, ASTM and other methods relevant to the certified properties listed below.

Certified Value: 10,036.5 mg/L**Uncertainty Associated with Measurement:** +/- 30 mg/L**Certified Value is Traceable to:** NIST SRM 3152a

The CRM is prepared gravimetrically using high purity Sodium Carbonate Lot# 02021A. The certified value listed is the average of values obtained by classical wet assay and ICP spectrometer analysis. Refer to side 2 for details of measurement uncertainties.

Classical Wet Assay: 10,039 mg/L**Method:** Evaporate to dryness. Fume with Sulfuric Acid. Ignite and weigh as Na₂SO₄.**Instrumentation Analysis By ICP spectrometer:** 10,034 mg/L**Uncertified Properties:****Density:** 1.048 @ 23.9 Degrees Celsius**Trace Metallic Impurities in the Actual Solution via ICP / ICPMS Analysis:**

Element	mg/L	Element	mg/L	Element	mg/L
Al	0.02	Cu	<0.002	Pb	<0.001
As	<0.08	Fe	0.03	Re	<0.001
Ag	<0.02	Ga	<0.001	Rb	<0.001
B	<0.1	In	<0.001	Sr	<0.001
Ba	0.008	K	1.36	Sb	<0.001
Be	<0.01	Li	<0.002	Sn	<0.001
Bi	<0.001	Mg	0.60	Ti	<0.03
Ca	0.60	Mn	<0.02	Tl	<0.001
Cr	0.002	Mo	<0.001	V	<0.001
Cd	<0.01	Ni	<0.003	Zr	<0.001
Co	<0.001			Zn	<0.05

Balances are calibrated regularly with weight sets traceable to NIST #32856, #32857 and others. This CRM is guaranteed stable to +/-0.5% of the certified concentration inclusive of uncertainty of measurements and other effects, such as transpiration losses, for a period of one year from the date of certification. This guarantee is valid only when the material is kept tightly capped and transported and stored under laboratory conditions.

Date of Certification: FEB -- '03 Certifying Officer: N. Kocherakota

INORGANIC LABS/RADCHEM LABS
DATE RECEIVED: 02/13/03
DATE EXPIRED: 02/15/2004
DATE OPENED: 02/13/03
INORG: 3951 PO: F52060

SPE Certificate™

Certificate of Reference Material

010169

Catalog Number: ICV-2A Lot No.: 22-12AS
Description: Initial Calibration Verification Standard II
Matrix: 5% Nitric Acid

This ASSURANCE® certified reference material, CRM, is intended primarily for use as a calibration standard or quality control standard for inorganic spectroscopic instrumentation such as ICPOES, DCP, AA, ICPMS, and XRF. It can be employed in USEPA, ASTM and other methods relevant to the certified properties listed below.

The CRM is prepared from high purity single element concentrates of individual elements using Class A laboratory ware to give precise concentration.

Refer to side 2 for details of measurement uncertainties.

Instrumental Analysis by ICP Spectrometer:

Element	Labeled (mg/L)	Measured (mg/L)	NIST SRM	Element	Labeled (mg/L)	Measured (mg/L)	NIST SRM
Ca	2,000	2,005.40	3109a	Ni	500	500.58	3136
K	2,000	1,997.89	3141a	V	500	504.23	3165
Mg	2,000	1,992.26	3131a	Cr	200	203.21	3112a
Na	2,000	1,992.99	3152a	Cu	200	199.75	3114
Al	1,000	1,005.90	3101a	Ag	100	100.46	3151
Ba	1,000	1,001.51	3104a	Be	100	100.04	3105a
Fe	1,000	1,003.17	3126a	Mn	100	100.64	3132
Co	500	505.10	3113	Zn	100	100.52	3168a

Spex Reference Multi: Lot #4-63BD, 14-125AS

Balances are calibrated regularly with weight sets traceable to NIST#s 32856, 32867 and others. This CRM is guaranteed stable and accurate to +/- 0.5% on the average of all the certified concentrations with no single component exceeding +/- 2%. This guarantee is valid for a period of one year from the date of certification only when the material is kept tightly capped and transported and stored under laboratory conditions.

NOV -- '02

Date of Certification: _____ Certifying Officer: N. Kocherakota

INORGANIC LABS/RADCHEM LABS
DATE RECEIVED: 11/19/02
DATE EXPIRED: 11/30/03
DATE OPENED: 11/19/02
INORG: 3923 PO: F58064
AR

SPE Certificate™

Certificate of Reference Material

010170

Catalog Number: PLSB7-2X/2Y/2T **Lot No.:** 8-175SB-X/Y/T
Description: 1000 mg/L Antimony
Matrix: Water/0.6% Tartaric Acid/tr HNO₃

This ASSURANCE® certified reference material, CRM, is intended primarily for use as a calibration standard or quality control standard for inorganic spectroscopic instrumentation such as ICPOES, DCP, AA, ICPMS, and XRF. It can be employed in USEPA, ASTM and other methods relevant to the certified properties listed below.

Certified Value: Antimony(Sb): 999 mg/L $\pm$ 3 mg/L
Traceable to: NIST SRM 3102a

The CRM is prepared gravimetrically using high purity Antimony Metal (Sb) Lot#R1198A. The certified value listed is the average of values obtained by classical wet assay and ICP spectrometer analysis.

Refer to side 2 for details of measurement uncertainties.

Classical Wet Assay: 998 mg/L

Method: Gravimetric analysis by evaporating to dryness, fuming with Nitric Acid, igniting and weighing as Sb₂O₄.

Instrumental Analysis by ICP spectrometer: 1000 mg/L

Uncertified Properties:

Density: 1.007 at 24.0°C

Trace Metallic Impurities in the Actual Solution via ICP Analysis:

Element	mg/L	Element	mg/L	Element	mg/L
Ag	<0.001	Cu	<0.001	Pb	<0.001
Al	0.030	Fe	0.013	Rb	<0.001
As	0.001	Ga	<0.001	Re	<0.001
B	<0.001	In	<0.001	Sn	<0.001
Ba	<0.001	K	0.030	Sr	<0.001
Be	<0.001	Li	<0.001	Ti	0.001
Bi	<0.001	Mg	0.003	Tl	<0.001
Ca	0.012	Mn	<0.001	V	<0.001
Cd	<0.001	Mo	<0.001	Zn	<0.010
Co	<0.001	Na	0.005	Zr	<0.001
Cr	<0.001	Ni	<0.001		

Balances are calibrated regularly with weight sets traceable to NIST#s 32856, 32867 and others. This CRM is guaranteed stable to $\pm$ 0.5% of the certified concentration inclusive of uncertainty of measurements and other effects, such as transpiration losses, for a period of one year from the date of certification. This guarantee is valid only when the material is kept tightly capped and transported and stored under laboratory conditions.

Date of Certification: NOV -- '02 Certifying Officer: N. Kocherakota

INORGANIC LABS/RADCHEM LABS
DATE RECEIVED: 11/19/02
DATE EXPIRED: 11/30/03
DATE OPENED: 11/19/02
INORG: 3224 PO: F53064

010171

SPE Certificate™**Certificate of Reference Material**

INORGANIC LABS/RADCHEM LABS
 DATE RECEIVED: 11/16/02
 DATE EXPIRED: 11/30/2003
 DATE OPENED: 11/19/02
 INORG: 3775 PO: F58664

Catalog Number: ICV-2C **Lot No.:** 22-13AS
Description: Initial Calibration Verification Standard II
Matrix: 5% Nitric Acid

This ASSURANCE® certified reference material, CRM, is intended primarily for use as a calibration standard or quality control standard for inorganic spectroscopic instrumentation such as ICPOES, DCP, AA, ICPMS, and XRF. It can be employed in USEPA, ASTM and other methods relevant to the certified properties listed below.

The CRM is prepared from high purity single element concentrates of individual elements using Class A laboratory ware to give precise concentration.

Refer to side 2 for details of measurement uncertainties.

Instrumental Analysis by ICP Spectrometer:

Element	Labeled (mg/L)	Measured (mg/L)	NIST SRM
As	500	497.85	3103a
Pb	500	495.41	3128
Se	500	501.98	3149
TL	500	501.89	3158
Cd	100	99.77	3108

Spex Reference Multi: Lot #4-51BDREF, 15-39AS, 11-173AS

Balances are calibrated regularly with weight sets traceable to NIST#s 32856, 32867 and others. This CRM is guaranteed stable and accurate to +/- 0.5% on the average of all the certified concentrations with no single element exceeding +/- 2%. This includes uncertainty of measurements and other effects, such as transpiration losses. This guarantee is valid for a period of one year from the date of certification only when the material is kept tightly capped and transported and stored under laboratory conditions.

Date of Certification: NOV -- '02 Certifying Officer: N. Kocherakota



inorganic ventures / iv labs

195 lehigh avenue, suite 4, lakewood, nj 08701 usa
phone: 800-669-6799 • 732-901-1900 • fax: 732-901-1903
e-mail: ivsales@ivstandards.com • website: www.ivstandards.com

certificate of analysis

010172

- 1.0 Inorganic Ventures / IV Labs is an ISO Guide 34-2000 Certified Reference Material (CRM) Manufacturer: Certificate #883-02.** The certificate is designed and the certified value(s) and uncertainty(ies) are determined in accordance with ISO Guide 31-2000 (Reference Materials - Contents of certificates and label(s)), ISO Guide 34-2000 "Quality System Guidelines for the Production of Reference Materials," and ISO Guide 35-1989 "Certification of Reference Materials - General and Statistical Principles."

2.0 DESCRIPTION OF CRM Custom-Grade 10000 µg/mL Aluminum in 5% (abs) HNO₃

Catalog Number: CGAL10-1 and CGAL10-5

Lot Number: W-AL04008

Starting Material: Al metal

Starting Material Purity (%): 99.998460

Starting Material Lot No 607116

Matrix: 5% (abs) HNO₃

INORGANIC LABS/RADCHEM LABS Pg. 1 of 2

DATE RECEIVED: 08/26/03

DATE EXPIRED: 09/01/2004

DATE OPENED: 08/26/03

INORG: 4220 PO: F52224

3.0 CERTIFIED VALUES AND UNCERTAINTIES

Certified Concentration: 10070 ± 31 µg/mL

Certified Density: 1.059 g/mL (measured at 22° C)

The Certified Value is based upon the most precise method used to analyze this CRM. The following equations are used in the calculation of the certified value and the uncertainty:

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

$$\text{Uncertainty } (\pm) = \frac{2[(\sum s_i)^2]^{1/2}}{(n)^{1/2}}$$

($\bar{x}$) = mean

x_i = individual results

n = number of measurements

$\sum s_i$ = The summation of all significant estimated errors.

(Most common are the errors from instrumental measurement, weighing, dilution to volume, and the fixed error reported on the NIST SRM certificate of analysis.)

The independent samples t-test was used to determine if there is agreement between the above assay methods at the 95% confidence interval. Both methods were compared and showed agreement within the stated uncertainties. This agreement is a confirmation of the accuracy of this CRM.

4.0 TRACEABILITY TO NIST AND VALUES OBTAINED BY INDEPENDENT METHODS

□ "Property of the result of a measurement or the value of a standard whereby it can be related to stated references, usually national or international standards, through an unbroken chain of comparisons all having stated uncertainties." (ISO VIM, 2nd ed., 1993, definition 6.10)

□ This IV product is Traceable to NIST via direct comparison to NIST SRMs. The uncertainties for each certified value are reported, taking into account the SRM uncertainty error and the measurement, weighing and volume dilution errors.

4.1 Assay Method #1 10006 ± 55 µg/mL

ICP Assay NIST SRM 3101a Lot Number: 992003

Assay Method #2 10070 ± 31 µg/mL

EDTA NIST SRM 928 Lot Number: 880710

4.2 BALANCE CALIBRATION - All balances are checked daily using in-house procedure number 6-IMM-001. The weights used for testing are annually compared to Gerhart Scale Corporation's master weights and are traceable to the National Institute of Standards and Technology (NIST). The NIST Traceability numbers are 692476 - Class 1 and 692476A - Class 2. The NIST test number is 822/260017-98. All analytical balances are calibrated every 4 months by Gerhart Scale Corp. of South Amboy. The balances are calibrated with a class 1 and/or class 2 analytical weight set. These weights are tested annually by a NIST / NVLAP accredited calibration lab. The NIST test number is 822/260017-98.

4.3 THERMOMETER CALIBRATION - The thermometers used in the determination of the final densities are calibrated vs standard thermometer No. 903-2680 which was certified in accordance with the procedures outlined by ASTM E77-87 and NIST Monograph 150 using NIST Test Nos. and Std Nos.: 769543, 217368/769543, 217368/P14452, 176240/P14452, 176240. The in-house procedure No. is 2-QC-001. Thermometers which are not calibrated vs standard thermometer No. 903-2680 are traceable to NIST Identification Nos. 92564, 119016, 471047 and NIST test report Nos. 811/258522, 811/2557078, and 236090.

4.4 GLASSWARE CALIBRATION - In-house procedure 3-QC-002 is used to calibrate all Class A Glassware used in the manufacture and quality control of Custom Grade Standards.

5.0 TRACE METALLIC IMPURITIES (TMI) DETERMINED BY ICP/MS AND ICP-OES IN µg/mL

Custom-Grade solutions are tested for trace metallic impurities by Axial ICP-OES and ICP-MS. The result from the most sensitive method for each element, is reported below. Solutions tested by ICP-MS were analyzed in an ULPA-Filtered Clean Room. An ULPA-Filter is 99.9985% efficient for the removal of particles down to 0.3 µm.

<u>S</u> Al	<u>M</u> Dy < 0.02695	<u>O</u> Li 0.00011	<u>M</u> Pr < 0.00135	<u>M</u> Te < 0.13473
<u>M</u> Sb < 0.00225	<u>M</u> Er < 0.02245	<u>M</u> Lu < 0.00180	<u>M</u> Re < 0.00449	<u>M</u> Tb < 0.00135
<u>M</u> As < 0.04491	<u>M</u> Eu < 0.01347	<u>O</u> Mg 0.00470	<u>M</u> Rh < 0.00449	<u>M</u> Tl < 0.00449
<u>M</u> Ba < 0.04491	<u>M</u> Gd < 0.00449	<u>M</u> Mn < 0.01796	<u>M</u> Rb < 0.00449	<u>M</u> Th < 0.00449
<u>O</u> Be < 0.00017	<u>M</u> Ga < 0.00449	<u>O</u> Hg < 0.00700	<u>M</u> Ru < 0.00898	<u>M</u> Tm < 0.00180
<u>M</u> Bi < 0.00180	<u>M</u> Ge < 0.02695	<u>M</u> Mo < 0.00898	<u>M</u> Sm < 0.00449	<u>M</u> Sn < 0.02245
<u>O</u> B 0.01164	<u>M</u> Au < 0.01347	<u>M</u> Nd < 0.00898	<u>M</u> Sc < 0.04491	<u>M</u> Ti < 0.22454
<u>M</u> Cd < 0.01347	<u>M</u> Hf < 0.00898	<u>O</u> Ni < 0.00600	<u>M</u> Se < 0.03593	<u>M</u> W < 0.04491
<u>O</u> Ca 0.01903	<u>M</u> Ho < 0.00225	<u>M</u> Nb < 0.00225	<u>O</u> Si 0.07389	<u>M</u> U < 0.00898
<u>M</u> Ce < 0.02245	<u>O</u> In < 0.03000	<u>n</u> Os	<u>M</u> Ag < 0.00898	<u>M</u> V < 0.00898
<u>M</u> Cs < 0.00135	<u>M</u> Ir < 0.02245	<u>M</u> Pd < 0.02245	<u>O</u> Na 0.03359	<u>M</u> Yb < 0.00449
<u>O</u> Cr 0.00336	<u>O</u> Fe 0.00493	<u>O</u> P < 0.03000	<u>M</u> Sr < 0.00225	<u>M</u> Y < 0.17963
<u>M</u> Co < 0.01347	<u>M</u> La < 0.00225	<u>M</u> Pt < 0.00898	<u>O</u> S < 0.10000	<u>M</u> Zn < 0.08982
<u>M</u> Cu < 0.02695	<u>M</u> Pb < 0.01347	<u>O</u> K 0.02911	<u>M</u> Ta < 0.03144	<u>M</u> Zr < 0.02245

M - Checked by ICP-MS O - Checked by ICP-OES i - Spectral Interference n - Not Checked For s - Solution Standard Element

6.0 INTENDED USE

For the calibration of analytical instruments including but not limited to the following:
 ICP-MS, ICP-OES, FAAS, GFAA, XRF, and DCP
 For the validation of analytical methods
 For the preparation of "working reference samples"
 For interference studies and the determination of correction coefficients
 For detection limit and linearity studies
 For additional intended uses, contact IV Technical Staff

7.0 INSTRUCTIONS FOR THE CORRECT USE OF THIS REFERENCE MATERIAL

Storage & Handling - Keep tightly sealed when not in use. Store and use at $20 \pm 4^\circ\text{C}$. Do not pipet from container. Do not return portions removed for pipetting to container.

Atomic Weight; Valence; Coordination Number; Chemical Form in Solution - 26.98154; +3, 6; $\text{Al}(\text{H}_2\text{O})_6^{3+}$

Chemical Compatibility - Soluble in HCl , HNO_3 , HF and H_2SO_4 . Avoid neutral media. Soluble in strongly basic NaOH forming the $\text{Al}(\text{OH})_4^-$ species. Stable with most metals and inorganic anions. The phosphate is insoluble in water and only slightly soluble in acid.

Stability - 2-100 ppb levels stable for months in 1% HNO_3 / LDPE container. 1-10,000 ppm solutions chemically stable for years in 2-5% HNO_3 / LDPE container.

Al Containing Samples (Preparation and Solution) - Metal (Best dissolved in HCl / HNO_3); α - Al_2O_3 (Na_2CO_3 fusion in Pt°); γ - Al_2O_3 (Soluble in acids such as HCl); Ores (Carbonate fusion in Pt° followed by HCl dissolution); Organic Matrices (sulfuric/peroxide digestion or nitric / sulfuric / perchloric acid decomposition, or dry ash and dissolution in dilute HCl).

Atomic Spectroscopic Information (ICP-OES D.L.s are given as radial/axial view):

Technique/Line	Estimated D.L.	Order	Type	Interferences (underlined indicates severe at 400nm)
ICP-OES 394.401 nm	0.05 / 0.006 $\mu\text{g/mL}$	1	atom	U, Ce
ICP-OES 396.152 nm	0.03 / 0.006 $\mu\text{g/mL}$	1	atom	Mo, Zr, Ce
ICP-OES 167.078 nm	0.1 / 0.009 $\mu\text{g/mL}$	1	ion	Fe
ICP-MS 27 amu	30 ppt	n/a	M'	$^{12}\text{C}^{14}\text{N}$, $^{13}\text{C}^{14}\text{N}$, $^{14}\text{C}^{14}\text{N}$, ^{18}O , ^{52}Cr , ^{56}Fe

8.0 **HAZARDOUS INFORMATION** - Please refer to the enclosed Material Safety Data sheet for information regarding this CRM.

9.0 **HOMOGENEITY** - This solution was mixed according to procedure IV-MPM-004 and is guaranteed to be homogeneous.

10.0 QUALITY STANDARD DOCUMENTATION

10.1 **ISO 9001:2000 Quality Management System Registration - QMI Certificate Number 010105**

Recognized by:

Registrar Accreditation Board (ANSI-RAB)

Standards Council of Canada (SCC)

Dutch Council for Accreditation (RVA)

Entidad Mexicana de Acreditación, a.c.(EMA)

Members of IQ Net International Certification Network:

Argentina (IRAM), Australia (QAS), Austria (ÖQS), Belgium (Avinter), Brazil (FCAV), Canada (QMI), Hong Kong (HKQAA), Columbia (ICONTEC), Czech Republic (CQS), Denmark (DS), Finland (SFS), France (AFAQ), Germany (DQS), Greece (ELOT), Hungary (MSZT), Ireland (NSAI), Israel (SII), Italy (CISQ), Japan (JQA), Korea (KSA-QA), Netherlands (KEMA), Norway (NCS), Poland (PCBC), Portugal (APCER), Singapore (PSB), Slovenia (SIQ), Spain (AENOR), Switzerland (SQS)

10.2 **ISO/IEC 17025 - 1999 "General Requirements for the Competence of Testing and Calibration"**

- Chemical Testing - Accredited A2LA Certificate Number 883.01



10.3 **ISO/IEC Guide 34 - 2000 "General Requirements for the Competence of Reference Material Producers"**

- Reference Materials Production - Accredited A2LA Certificate Number 883.02

A2LA Mutual Recognition Agreement Partners:

Australia (NATA), Austria (BmWA), Belgium (BELTEST) (BKO-OBE), Canada (SCC), Chinese Taipei (CNLA), Czech Republic (NAO), Denmark (DANAK), Finland (FINAS), France (COFRAC), Germany (DAR), Hong Kong (HKAS), Ireland (NAB), Italy (SIT) (SINAL), Japan (JAB) (JNLA), Republic of Korea (KOLAS), The Netherlands (RvA), New Zealand (IANZ), Norway (NA), Portugal (IPQ), Singapore (SAC-SINGLAS), Spain (ENAC), Sweden (SWEDAC), Switzerland (SAS), United Kingdom (UKAS) and United States (NVLAP) (ICBO ES)

10.4 **10CFR50 Appendix B - Nuclear Regulatory Commission**

- Domestic Licensing of Production and Utilization Facilities

10.5 **10CFR21 - Nuclear Regulatory Commission - Reporting Defects and Non-Compliance**

10.6 **MIL-STD-45662A (Obsolete/Observed)**

INORGANIC LABS/RADCHEM LABS Pg. 2 of 2
 DATE RECEIVED: 08/26/03
 DATE EXPIRED: 09/01/2004
 DATE OPENED: 08/26/03
 INORG: 4020 PO: F52224

11.0 DATE OF CERTIFICATION AND PERIOD OF VALIDITY



11.1 IV Shelf Life - The period of time during which the concentration of the analyte(s) in a properly packaged, unopened, and unused standard stored under environmentally controlled and monitored conditions will remain within the specified uncertainty range. Shelf life is limited primarily by transpiration (loss of water from the solution) and infrequently, by chemical instability. Transpiration studies (P-SP01020) of chemically-stable solutions performed at Inorganic Ventures / IV Labs indicate a CRM shelf-life of four years for solutions packaged in 500-mL low density polyethylene bottles. When stored under special conditions that minimize transpiration and instability, the shelf life can be extended past this limit.

11.2 Expiration Date - The date after which a CRM should not be used. Routine laboratory use of a CRM increases transpiration losses and the chance of contamination which affect the integrity of the CRM and limit its useful life. Inorganic Ventures / IV Labs concurs with state and federal regulatory agencies' recommendations that solution standards be assigned a one-year expiration date.

Certification Date: February 13, 2003

Expiration Date:

EXPIRES
1/1/2004

12.0 NAMES AND SIGNATURES OF CERTIFYING OFFICERS

Certificate Prepared By: JoAnn Struthers, QA Administrative Assistant

JoAnn Struthers

Certificate Approved By: Katalin Le, QC Supervisor

Katalin Le

Certifying Officer: Paul Gaines, Chemist, Senior Technical Director

Paul Gaines



inorganic ventures / iv labs

195 lehigh avenue, suite 4, lakewood, nj 08701 usa
phone: 800-669-6799 • 732-901-1900 • fax: 732-901-1903
e-mail: ivsales@ivstandards.com • website: www.ivstandards.com

010176

certificate of analysis

1.0 Inorganic Ventures / IV Labs is an ISO Guide 34-2000 Certified Reference Material (CRM) Manufacturer: Certificate #883-02. The certificate is designed and the certified value(s) and uncertainty(ies) are determined in accordance with ISO Guide 31-2000 (Reference Materials - Contents of certificates and label(s)), ISO Guide 34-2000 "Quality System Guidelines for the Production of Reference Materials," and ISO Guide 35-1989 "Certification of Reference Materials - General and Statistical Principles."

2.0 DESCRIPTION OF CRM Custom-Grade 10000 µg/mL Calcium in 1.4% (abs) HNO₃

Catalog Number: CGCA10-1 AND CGCA10-5

Lot Number: T-CA03010

Starting Material: CaO

Starting Material Purity (%): 99.9981

Starting Material Lot No C27L01

Matrix: 1.4% (abs) HNO₃

INORGANIC LABS/RADCHEM LABS Pg. 1 of 2

DATE RECEIVED: 01/09/03

DATE EXPIRED: 02/01/2004

DATE OPENED: 01/10/03

INORG: 3900 PO: F50071

3.0 CERTIFIED VALUES AND UNCERTAINTIES

Certified Concentration: 10,007 ± 22 µg/mL

Certified Density: 1.037 g/mL (measured at 22° C)

The Certified Value is based upon the most precise method used to analyze this CRM. The following equations are used in the calculation of the certified value and the uncertainty:

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

($\bar{x}$) = mean

x_i = individual results

n = number of measurements

$$\text{Uncertainty } (\pm) = \frac{2(\sum s_i^2)^{1/2}}{(n)^{1/2}}$$

$\sum s_i^2$ = The summation of all significant estimated errors.

(Most common are the errors from instrumental measurement, weighing, dilution to volume, and the fixed error reported on the NIST SRM certificate of analysis.)

The independent samples t-test was used to determine if there is agreement between the above assay methods at the 95% confidence interval. Both methods were compared and showed agreement within the stated uncertainties. This agreement is a confirmation of the accuracy of this CRM.

4.0 TRACEABILITY TO NIST AND VALUES OBTAINED BY INDEPENDENT METHODS

"Property of the result of a measurement or the value of a standard whereby it can be related to stated references, usually national or international standards, through an unbroken chain of comparisons all having stated uncertainties." (ISO VIM, 2nd ed., 1993, definition 6.10)

This IV product is Traceable to NIST via direct comparison to NIST SRMs. The uncertainties for each certified value are reported, taking into account the SRM uncertainty error and the measurement, weighing and volume dilution errors.

4.1 Assay Method #1 10,034 ± 25 µg/mL

ICP Assay NIST SRM 3109a Lot Number: 000622

Assay Method #2 10,007 ± 22 µg/mL

EDTA NIST SRM 928 Lot Number: 880710

4.2 BALANCE CALIBRATION - All balances are checked daily using in-house procedure number 6-IMM-001. The weights used for testing are annually compared to Gerhart Scale Corporation's master weights and are traceable to the National Institute of Standards and Technology (NIST). The NIST Traceability numbers are 692476 - Class 1 and 692476A - Class 2. The NIST test number is 822/260017-98. All analytical balances are calibrated every 4 months by Gerhart Scale Corp. of South Amboy. The balances are calibrated with a class 1 and/or class 2 analytical weight set. These weights are tested annually by a NIST / NVLAP accredited calibration lab. The NIST test number is 822/260017-98.

010177

4.3 THERMOMETER CALIBRATION - The thermometers used in the determination of the final densities are calibrated vs standard thermometer No. 903-2680 which was certified in accordance with the procedures outlined by ASTM E77-87 and NIST Monograph 150 using NIST Test Nos. and Std Nos.: 769543, 217368/769543, 217368/P14452, 176240/P14452, 176240. The in-house procedure No. is 2-QC-001. Thermometers which are not calibrated vs standard thermometer No. 903-2680 are traceable to NIST Identification Nos. 92564, 119016, 471047 and NIST test report Nos. 811/258522, 811/2557078, and 236090.

4.4 GLASSWARE CALIBRATION - In-house procedure 3-QC-002 is used to calibrate all Class A Glassware used in the manufacture and quality control of Custom Grade Standards.

5.0 TRACE METALLIC IMPURITIES (TMI) DETERMINED BY ICP/MS AND ICP-OES IN µg/mL

Custom-Grade solutions are tested for trace metallic impurities by Axial ICP-OES and ICP-MS. The result from the most sensitive method for each element, is reported below. Solutions tested by ICP-MS were analyzed in an ULPA-Filtered Clean Room. An ULPA-Filter is 99.9985% efficient for the removal of particles down to 0.3 µm.

Q Al	0.00064	M Dy	< 0.02560	Q Li	< 0.00002	M Pr	< 0.00128	M Te	< 0.12798
M Sb	< 0.00213	M Er	< 0.02133	M Lu	< 0.00171	M Re	< 0.00427	M Tb	< 0.00128
M As	< 0.04266	M Eu	< 0.01280	Q Mg	0.07143	M Rh	< 0.00427	M Tl	< 0.00427
Q Ba	0.00071	M Gd	< 0.00427	Q Mn	0.00041	M Rb	< 0.00427	M Th	< 0.00427
Q Be	< 0.00009	M Ga	< 0.00427	Q Hg	< 0.01100	M Ru	< 0.00853	M Tm	< 0.00171
M Bi	< 0.00171	M Ge	< 0.02560	M Mo	< 0.00853	M Sm	< 0.00427	M Sn	< 0.02133
Q B	< 0.00054	M Au	< 0.01280	M Nd	< 0.00853	Q Sc	< 0.00002	M Ti	< 0.21329
Q Cd	< 0.00450	M Hf	< 0.00853	Q Ni	< 0.00230	Q Se	< 0.00620	M W	< 0.04266
S Ca		M Ho	< 0.00213	M Nb	< 0.00213	Q Si	0.00214	M U	< 0.00853
M Ce	< 0.02133	Q In	< 0.00200	n Os		Q Ag	< 0.04000	Q V	< 0.00090
M Cs	< 0.00128	M Ir	< 0.02133	M Pd	< 0.02133	Q Na	0.00571	M Yb	< 0.00427
Q Cr	0.00238	Q Fe	< 0.00110	Q P	< 0.00480	Q Sr	0.08095	M Y	< 0.17064
Q Co	< 0.00120	M La	< 0.00213	M Pt	< 0.00853	Q S	0.04048	Q Zn	0.07381
Q Cu	0.00405	M Pb	< 0.01280	Q K	< 0.00170	M Ta	< 0.02986	M Zr	< 0.02133

M - Checked by ICP-MS O - Checked by ICP-OES i - Spectral Interference n - Not Checked For s - Solution Standard Element

6.0 INTENDED USE

For the calibration of analytical instruments including but not limited to the following:

ICP-MS, ICP-OES, FAAS, GFAA, XRF, and DCP

For the validation of analytical methods

For the preparation of "working reference samples"

For interference studies and the determination of correction coefficients

For detection limit and linearity studies

For additional intended uses, contact IV Technical Staff

7.0 INSTRUCTIONS FOR THE CORRECT USE OF THIS REFERENCE MATERIAL

010178

Storage & Handling - Keep tightly sealed when not in use. Store and use at $20 \pm 4^\circ\text{C}$. Do not pipet from container. Do not return portions removed for pipetting to container.

Atomic Weight; Valence; Coordination Number; Chemical Form in Solution - 40.078; +2; 6; $\text{Ca}(\text{H}_2\text{O})_6^{++}$

Chemical Compatibility - Soluble in HCl and HNO_3 . Avoid H_2SO_4 , HF, H_3PO_4 , and neutral to basic media. Stable with most metals and inorganic anions forming insoluble silicate, carbonate, hydroxide, oxide, fluoride, sulfate, oxalate, chromate, arsenate and tungstate in neutral aqueous media.

Stability - 2-100 ppb levels stable for months in 1% HNO_3 / LDPE container. 1-10,000 ppm solutions chemically stable for years in 1-10% HNO_3 / LDPE container.

Ca Containing Samples (Preparation and Solution) - Metal (best dissolved in diluted HNO_3), Ores (Carbonate fusion in P^{1+} followed by HCl dissolution); Organic Matrices (dry ash and dissolution in dilute HCl. Do not heat when dissolving to avoid precipitation of SiO_2). The oxide, hydroxide, carbonate, phosphate, and fluoride of calcium are soluble in % levels of HCl or HNO_3 . The sulfates (gypsum, anhydrite, etc.), certain silicates and complex compounds require fusion with Na_2CO_3 , followed by HCl / water dissolution. Contamination is a very real problem when analyzing for trace levels.

Atomic Spectroscopic Information (ICP-OES D.L.s are given as radial/axial view):

Technique/Line	Estimated D.L.	Order	Type	Interferences (underlined indicates severe at concs.)
ICP-OES 393.366 nm	0.0002 / 0.00004 $\mu\text{g/mL}$	1	ion	U, Ce
ICP-OES 396.847 nm	0.0005 / 0.00006 $\mu\text{g/mL}$	1	ion	Th
ICP-OES 422.673 nm	0.01 / 0.001 $\mu\text{g/mL}$	1	atom	Ge
ICP-MS 44 amu	1200 ppt	n/a	M'	^{40}O , ^{12}C , ^{28}Si , ^{40}Ar , ^{88}Sr

8.0 HAZARDOUS INFORMATION - Please refer to the enclosed Material Safety Data sheet for information regarding this CRM.

9.0 HOMOGENEITY - This solution was mixed according to procedure IV-MPM-004 and is guaranteed to be homogeneous.

10.0 QUALITY STANDARD DOCUMENTATION

10.1 ISO 9001:2000 Quality Management System Registration - QMI Certificate Number 010105

Recognized by:

Registrar Accreditation Board (ANSI-RAB)

Standards Council of Canada (SCC)

Dutch Council for Accreditation (RVA)

Entidad Mexicana de Acreditacion, a.c.(EMA)

Members of IQ Net International Certification Network:

Argentina (IRAM), Australia (QAS), Austria (ÖQS), Belgium (Avinter), Brazil (FCAV), Canada (QMI), Hong Kong (HKQAA), Columbia (ICONTEC), Czech Republic (CQS), Denmark (DS), Finland (SFS), France (AFAQ), Germany (DQS), Greece (ELOT), Hungary (MSZT), Ireland (NSAI), Israel (SII), Italy (CISQ), Japan (JQA), Korea (KSA-QA), Netherlands (KEMA), Norway (NCS), Poland (PCBC), Portugal (APCER), Singapore (PSB), Slovenia (SIQ), Spain (AENOR), Switzerland (SQS)

10.2 ISO/IEC 17025 - 1999 "General Requirements for the Competence of Testing and Calibration"

- Chemical Testing - Accredited A2LA Certificate Number 883.01



10.3 ISO/IEC Guide 34 - 2000 "General Requirements for the Competence of Reference Material Producers"

- Reference Materials Production - Accredited A2LA Certificate Number 883.02

A2LA Mutual Recognition Agreement Partners:

Australia (NATA), Austria (BmWA), Belgium (BELTEST) (BKO-OBE), Canada (SCC), Chinese Taipei (CNLA), Czech Republic (NAO), Denmark (DANAK), Finland (FINAS), France (COFRAC), Germany (DAR), Hong Kong (HKAS), Ireland (NAB), Italy (SIT) (SINAL), Japan (JAB) (JNLA), Republic of Korea (KOLAS), The Netherlands (RvA), New Zealand (IANZ), Norway (NA), Portugal (IPQ), Singapore (SAC-SINGLAS), Spain (ENAC), Sweden (SWEDAC), Switzerland (SAS), United Kingdom (UKAS) and United States (NVLAP) (ICBO ES)

10.4 10CFR50 Appendix B - Nuclear Regulatory Commission

- Domestic Licensing of Production and Utilization Facilities

10.5 10CFR21 - Nuclear Regulatory Commission - Reporting Defects and Non-Compliance

10.6 MIL-STD-45662A (Obsolete/Observed)

INORGANIC LABS/RADCHEM LABS Pg. 2 of 2


DATE RECEIVED: 01/09/03

DATE EXPIRED: 02/01/2004

DATE OPENED: 01/10/03

INORG: 3900 PO: F52071

11.0 DATE OF CERTIFICATION AND PERIOD OF VALIDITY



11.1 IV Shelf Life - The period of time during which the concentration of the analyte(s) in a properly packaged, unopened, and unused standard stored under environmentally controlled and monitored conditions will remain within the specified uncertainty range. Shelf life is limited primarily by transpiration (loss of water from the solution) and infrequently, by chemical instability. Transpiration studies (P-SP01020) of chemically-stable solutions performed at Inorganic Ventures / IV Labs indicate a CRM shelf-life of four years for solutions packaged in 500-mL low density polyethylene bottles. When stored under special conditions that minimize transpiration and instability, the shelf life can be extended past this limit.

11.2 Expiration Date - The date after which a CRM should not be used. Routine laboratory use of a CRM increases transpiration losses and the chance of contamination which affect the integrity of the CRM and limit its useful life. Inorganic Ventures / IV Labs concurs with state and federal regulatory agencies' recommendations that solution standards be assigned a one-year expiration date.

Certification Date: August 26, 2002

Expiration Date: **EXPIRES**
1st 2004

12.0 NAMES AND SIGNATURES OF CERTIFYING OFFICERS

Certificate Prepared By: Debbie Newman, QA Administrator

Certificate Approved By: Katalin Le, QC Supervisor

Certifying Officer: Paul Gaines, Chemist, Senior Technical Director

Debbie Newman
Katalin Le
Paul Gaines

certificate of analysis

1.0 Inorganic Ventures / IV Labs is an ISO Guide 34-2000 Certified Reference Material (CRM) Manufacturer: Certificate #883-02. The certificate is designed and the certified value(s) and uncertainty(ies) are determined in accordance with ISO Guide 31-2000 (Reference Materials - Contents of certificates and label(s)), ISO Guide 34-2000 "Quality System Guidelines for the Production of Reference Materials," and ISO Guide 35-1989 "Certification of Reference Materials - General and Statistical Principles."

2.0 DESCRIPTION OF CRM Custom-Grade 10000 µg/mL Iron in 3.5% (abs) HNO₃

Catalog Number: CGFE10-1 and CGFE10-5
Lot Number: T-FE03028
Starting Material: Fe metal
Starting Material Purity (%): 99.9992
Starting Material Lot No: 23024
Matrix: 3.5% (abs) HNO₃

INORGANIC LABS/RADCHEM LABS Pg. 1 of 2
DATE RECEIVED: 2/28/03
DATE EXPIRED: 3/01/2004
DATE OPENED: 2/28/03
INORG: 3963 PO: E52064

3.0 CERTIFIED VALUES AND UNCERTAINTIES

Certified Concentration: 9920 ± 22 µg/mL

Certified Density: 1.037 g/mL (measured at 22° C)

The Certified Value is based upon the most precise method used to analyze this CRM. The following equations are used in the calculation of the certified value and the uncertainty:

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

($\bar{x}$) = mean

x_i = individual results

n = number of measurements

$$\text{Uncertainty } (\pm) = \frac{2[(\sum s_i^2)^{1/2}]}{(n)^{1/2}}$$

$\sum S$ = The summation of all significant estimated errors.

(Most common are the errors from instrumental measurement, weighing, dilution to volume, and the fixed error reported on the NIST SRM certificate of analysis.)

The independent samples t-test was used to determine if there is agreement between the above assay methods at the 95% confidence interval. Both methods were compared and showed agreement within the stated uncertainties. This agreement is a confirmation of the accuracy of this CRM.

4.0 TRACEABILITY TO NIST AND VALUES OBTAINED BY INDEPENDENT METHODS

• "Property of the result of a measurement or the value of a standard whereby it can be related to stated references, usually national or international standards, through an unbroken chain of comparisons all having stated uncertainties." (ISO VIM, 2nd ed., 1993, definition 6.10)

• This IV product is Traceable to NIST via direct comparison to NIST SRMs. The uncertainties for each certified value are reported, taking into account the SRM uncertainty error and the measurement, weighing and volume dilution errors.

4.1 Assay Method #1 9920 ± 22 µg/mL

ICP Assay NIST SRM 3126a Lot Number: 000606

Assay Method #2 9962 ± 41 µg/mL

EDTA NIST SRM 928 Lot Number: 880710

4.2 BALANCE CALIBRATION - All balances are checked daily using in-house procedure number 6-IMM-001. The weights used for testing are annually compared to Gerhart Scale Corporation's master weights and are traceable to the National Institute of Standards and Technology (NIST). The NIST Traceability numbers are 692476 - Class 1 and 692476A - Class 2. The NIST test number is 822/260017-98. All analytical balances are calibrated every 4 months by Gerhart Scale Corp. of South Amboy. The balances are calibrated with a class 1 and/or class 2 analytical weight set. These weights are tested annually by a NIST / NVLAP accredited calibration lab. The NIST test number is 822/260017-98.

010181

4.3 THERMOMETER CALIBRATION - The thermometers used in the determination of the final densities are calibrated vs standard thermometer No. 903-2680 which was certified in accordance with the procedures outlined by ASTM E77-87 and NIST Monograph 150 using NIST Test Nos. and Std Nos.: 769543, 217368/769543, 217368/P14452, 176240/P14452, 176240. The in-house procedure No. is 2-QC-001. Thermometers which are not calibrated vs standard thermometer No. 903-2680 are traceable to NIST Identification Nos. 92564, 119016, 471047 and NIST test report Nos. 811/258522, 811/2557078, and 236090.

4.4 GLASSWARE CALIBRATION - In-house procedure 3-QC-002 is used to calibrate all Class A Glassware used in the manufacture and quality control of Custom Grade Standards.

5.0 TRACE METALLIC IMPURITIES (TMI) DETERMINED BY ICP/MS AND ICP-OES IN µg/mL

Custom-Grade solutions are tested for trace metallic impurities by Axial ICP-OES and ICP-MS. The result from the most sensitive method for each element, is reported below. Solutions tested by ICP-MS were analyzed in an ULPA-Filtered Clean Room. An ULPA-Filter is 99.9985% efficient for the removal of particles down to 0.3 µm.

<u>Q</u> Al < 0.00270	<u>M</u> Dy < 0.02421	<u>Q</u> Li < 0.00003	<u>M</u> Pr < 0.00121	<u>M</u> Te < 0.12103
<u>M</u> Sb < 0.00202	<u>M</u> Er < 0.02017	<u>M</u> Lu < 0.00161	<u>M</u> Re < 0.00403	<u>M</u> Tb < 0.00121
<u>M</u> As < 0.04034	<u>M</u> Eu < 0.01210	<u>Q</u> Mg < 0.00006	<u>M</u> Rh < 0.00403	<u>M</u> Tl < 0.00403
<u>M</u> Ba < 0.04034	<u>M</u> Gd < 0.00403	<u>Q</u> Mn < 0.05000	<u>M</u> Rb < 0.00403	<u>M</u> Th < 0.00403
<u>Q</u> Be < 0.00005	<u>M</u> Ga 0.00394	<u>Q</u> Hg < 0.01100	<u>M</u> Ru < 0.00807	<u>M</u> Tm < 0.00161
<u>M</u> Bi < 0.00161	<u>i</u> Ge	<u>M</u> Mo < 0.00807	<u>M</u> Sm < 0.00403	<u>M</u> Sn 0.04920
<u>Q</u> B < 0.00090	<u>M</u> Au < 0.01210	<u>M</u> Nd < 0.00807	<u>M</u> Sc < 0.04034	<u>M</u> Ti < 0.20172
<u>M</u> Cd < 0.01210	<u>M</u> Hf < 0.00807	<u>Q</u> Ni < 0.00230	<u>M</u> Se < 0.03228	<u>M</u> W < 0.04034
<u>Q</u> Ca 0.00707	<u>M</u> Ho < 0.00202	<u>M</u> Nb < 0.00202	<u>Q</u> Si 0.00781	<u>M</u> U < 0.00807
<u>M</u> Ce < 0.02017	<u>M</u> In < 0.04034	<u>n</u> Os	<u>M</u> Ag < 0.00807	<u>M</u> V < 0.00807
<u>M</u> Cs < 0.00121	<u>M</u> Ir < 0.02017	<u>M</u> Pd < 0.02017	<u>Q</u> Na 0.00756	<u>M</u> Yb < 0.00403
<u>M</u> Cr 0.00541	<u>s</u> Fe	<u>i</u> P	<u>M</u> Sr < 0.00202	<u>M</u> Y < 0.16138
<u>Q</u> Co < 0.00110	<u>M</u> La < 0.00202	<u>M</u> Pt < 0.00807	<u>Q</u> S < 0.07200	<u>M</u> Zn 0.03739
<u>M</u> Cu < 0.02421	<u>M</u> Pb < 0.01210	<u>Q</u> K < 0.00170	<u>M</u> Ta < 0.02824	<u>M</u> Zr < 0.02017

M - Checked by ICP-MS O - Checked by ICP-OES i - Spectral Interference n - Not Checked For s - Solution Standard Element

6.0 INTENDED USE

For the calibration of analytical instruments including but not limited to the following:
 ICP-MS, ICP-OES, FAAS, GFAA, XRF, and DCP
 For the validation of analytical methods
 For the preparation of "working reference samples"
 For interference studies and the determination of correction coefficients
 For detection limit and linearity studies
 For additional intended uses, contact IV Technical Staff

7.0 INSTRUCTIONS FOR THE CORRECT USE OF THIS REFERENCE MATERIAL

Storage & Handling - Keep tightly sealed when not in use. Store and use at $20 \pm 4^\circ\text{C}$. Do not pipet from container. Do not return portions removed for pipetting to container.

Atomic Weight; Valence; Coordination Number; Chemical Form in Solution - 55.847; +3; 6; $\text{Fe}(\text{H}_2\text{O})_6^{3+}$

Chemical Compatibility - Stable in HCl , HNO_3 , H_2SO_4 , HF and H_3PO_4 . Avoid basic media. Stable with most metals and inorganic anions in acidic media.

Stability - 2-100 ppb levels stable for months in 1% HNO_3 / LDPE container. 1-10,000 ppm solutions chemically stable for years in 1-5% HNO_3 / LDPE container.

Fe Containing Samples (Preparation and Solution) - Metal (Soluble in HCl); Oxides (If the oxide has been at a high temperature then Na_2CO_3 fusion in Pt° followed by HCl dissolution otherwise dissolve in dilute HCl); Ores (See Oxides above using only the fusion approach).

Atomic Spectroscopic Information (ICP-OES D.L.s are given as radial/axial view):

Technique/Line	Estimated D.L.	Order	Type	Interferences (underlined indicates severe at = concs.)
ICP-OES 238.204 nm	0.005 / 0.001 $\mu\text{g/mL}$	1	ion	Ru, Co
ICP-OES 239.562 nm	0.005 / 0.001 $\mu\text{g/mL}$	1	ion	Co, W, Cr
ICP-OES 259.940 nm	0.006 / 0.001 $\mu\text{g/mL}$	1	ion	Hf, Nb
ICP-MS 56 amu	970 ppt	n/a	M ⁺	$^{40}\text{Ar}^{15}\text{N}^+\text{H}$, $^{40}\text{Ar}^{16}\text{O}$, $^{36}\text{Ar}^{17}\text{O}^+\text{H}$, $^{38}\text{Ar}^{18}\text{O}$, $^{37}\text{Cl}^{18}\text{O}^+\text{H}$, $^{40}\text{Ca}^{16}\text{O}$

8.0 **HAZARDOUS INFORMATION** - Please refer to the enclosed Material Safety Data sheet for information regarding this CRM.

9.0 **HOMOGENEITY** - This solution was mixed according to procedure IV-MPM-004 and is guaranteed to be homogeneous.

10.0 QUALITY STANDARD DOCUMENTATION

10.1 ISO 9001:2000 Quality Management System Registration - QMI Certificate Number 010105

Recognized by:

Registrar Accreditation Board (ANSI-RAB)

Standards Council of Canada (SCC)

Dutch Council for Accreditation (RVA)

Entidad Mexicana de Acreditacion, a.c.(EMA)

Members of IQ Net International Certification Network:

Argentina (IRAM), Australia (QAS), Austria (ÖQS), Belgium (Avinter), Brazil (FCAV), Canada (QMI), Hong Kong (HKQAA), Columbia (ICONTEC), Czech Republic (CQS), Denmark (DS), Finland (SFS), France (AFAQ), Germany (DQS), Greece (ELOT), Hungary (MSZT), Ireland (NSAI), Israel (SII), Italy (CISQ), Japan (JQA), Korea (KSA-QA), Netherlands (KEMA), Norway (NCS), Poland(PCBC), Portugal (APCER), Singapore (PSB), Slovenia (SIQ), Spain (AENOR), Switzerland (SQS)



10.2 ISO/IEC 17025 - 1999 "General Requirements for the Competence of Testing and Calibration"

- Chemical Testing - Accredited A2LA Certificate Number 883.01



10.3 ISO/IEC Guide 34 - 2000 "General Requirements for the Competence of Reference Material Producers"

- Reference Materials Production - Accredited A2LA Certificate Number 883.02

A2LA Mutual Recognition Agreement Partners:

Australia (NATA), Austria (BmWA), Belgium (BELTEST) (BKO-OBE), Canada (SCC), Chinese Taipei (CNLA), Czech Republic (NAO), Denmark (DANAK), Finland (FINAS), France (COFRAC), Germany (DAR), Hong Kong (HKAS), Ireland (NAB), Italy (SIT) (SINAL), Japan (JAB) (JNLA), Republic of Korea (KOLAS), The Netherlands (RvA), New Zealand (IANZ), Norway (NA), Portugal (IPQ), Singapore (SAC-SINGLAS), Spain (ENAC), Sweden (SWEDAC), Switzerland (SAS), United Kingdom (UKAS) and United States (NVLAP) (ICBO ES)

10.4 10CFR50 Appendix B - Nuclear Regulatory Commission

- Domestic Licensing of Production and Utilization Facilities

10.5 10CFR21 - Nuclear Regulatory Commission - Reporting Defects and Non-Compliance

10.6 MIL-STD-45662A (Obsolete/Observed)

INORGANIC LABS/RADCHEM LABS Pg. 2 of 2

DATE RECEIVED: 2/28/03


DATE EXPIRED: 3/01/2004

DATE OPENED: 2/28/03

INORG: 3963 PO: F52064

11.0 DATE OF CERTIFICATION AND PERIOD OF VALIDITY

010183



11.1 IV Shelf Life - The period of time during which the concentration of the analyte(s) in a properly packaged, unopened, and unused standard stored under environmentally controlled and monitored conditions will remain within the specified uncertainty range. Shelf life is limited primarily by transpiration (loss of water from the solution) and infrequently, by chemical instability. Transpiration studies (P-SP01020) of chemically-stable solutions performed at Inorganic Ventures / IV Labs indicate a CRM shelf-life of four years for solutions packaged in 500-mL low density polyethylene bottles. When stored under special conditions that minimize transpiration and instability, the shelf life can be extended past this limit.

11.2 Expiration Date - The date after which a CRM should not be used. Routine laboratory use of a CRM increases transpiration losses and the chance of contamination which affect the integrity of the CRM and limit its useful life. Inorganic Ventures / IV Labs concurs with state and federal regulatory agencies' recommendations that solution standards be assigned a one-year expiration date.

Certification Date: October 11, 2002

Expiration Date:

EXPIRES
12/2004

12.0 NAMES AND SIGNATURES OF CERTIFYING OFFICERS

Certificate Prepared By: Debbie Newman, QA Administrator

Debbie Newman

Certificate Approved By: Katalin Le, QC Supervisor

Katalin Le

Certifying Officer: Paul Gaines, Chemist, Senior Technical Director

Paul Gaines

certificate of analysis

1.0 Inorganic Ventures / IV Labs is an ISO Guide 34-2000 Certified Reference Material (CRM) Manufacturer: Certificate #883-02. The certificate is designed and the certified value(s) and uncertainty(ies) are determined in accordance with ISO Guide 31-2000 (Reference Materials - Contents of certificates and label(s)), ISO Guide 34-2000 "Quality System Guidelines for the Production of Reference Materials," and ISO Guide 35-1989 "Certification of Reference Materials - General and Statistical Principles."

2.0 DESCRIPTION OF CRM Custom-Grade 10000 µg/mL Potassium in 1.4% (abs) HNO₃

Catalog Number: CGK10-1, CGK10-2, and CGK10-5
Lot Number: W-K02111
Starting Material: KNO₃
Starting Material Purity (%): 99.997230
Starting Material Lot No: K18J19
Matrix: 1.4% (abs) HNO₃

INORGANIC LABS/RADCHEM LABS
DATE RECEIVED: 11/5/03 1/2
DATE EXPIRED: 12/1/2004 OF
DATE OPENED: 11/5/03
INORG: 4320 PO: FS2256

3.0 CERTIFIED VALUES AND UNCERTAINTIES

Certified Concentration: 9930 ± 9 µg/mL

Certified Density: 1.024 g/mL (measured at 22° C)

The Certified Value is based upon the most precise method used to analyze this CRM. The following equations are used in the calculation of the certified value and the uncertainty:

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

($\bar{x}$) = mean

x_i = individual results

n = number of measurements

$$\text{Uncertainty } (\pm) = \frac{2[(\sum s_i^2)^{1/2}]}{(n)^{1/2}}$$

$\sum S$ = The summation of all significant estimated errors.

(Most common are the errors from instrumental measurement, weighing, dilution to volume, and the fixed error reported on the NIST SRM certificate of analysis.)

The independent samples t-test was used to determine if there is agreement between the above assay methods at the 95% confidence interval. Both methods were compared and showed agreement within the stated uncertainties. This agreement is a confirmation of the accuracy of this CRM.

4.0 TRACEABILITY TO NIST AND VALUES OBTAINED BY INDEPENDENT METHODS

□ "Property of the result of a measurement or the value of a standard whereby it can be related to stated references, usually national or international standards, through an unbroken chain of comparisons all having stated uncertainties." (ISO VIM, 2nd ed., 1993, definition 6.10)

□ This IV product is Traceable to NIST via direct comparison to NIST SRMs. The uncertainties for each certified value are reported, taking into account the SRM uncertainty error and the measurement, weighing and volume dilution errors.

4.1 Assay Method #1 9926 ± 62 µg/mL

ICP Assay NIST SRM 3141a Lot Number: 891312

Assay Method #2 9930 ± 9 µg/mL

Gravimetric NIST SRM Lot Number: See Sec. 4.2

4.2 BALANCE CALIBRATION - All balances are checked daily using in-house procedure number 6-IMM-001. The weights used for testing are annually compared to Gerhart Scale Corporation's master weights and are traceable to the National Institute of Standards and Technology (NIST). The NIST Traceability numbers are 692476 - Class 1 and 692476A - Class 2. The NIST test number is 822/260017-98. All analytical balances are calibrated every 4 months by Gerhart Scale Corp. of South Amboy. The balances are calibrated with a class 1 and/or class 2 analytical weight set. These weights are tested annually by a NIST / NVLAP accredited calibration lab. The NIST test number is 822/260017-98.

4.3 THERMOMETER CALIBRATION - The thermometers used in the determination of the final densities are calibrated vs standard thermometer No. 903-2680 which was certified in accordance with the procedures outlined by ASTM E77-87 and NIST Monograph 150 using NIST Test Nos. and Std Nos.: 769543, 217368/769543, 217368/P14452, 176240/P14452, 176240. The in-house procedure No. is 2-QC-001. Thermometers which are not calibrated vs standard thermometer No. 903-2680 are traceable to NIST Identification Nos. 92564, 119016, 471047 and NIST test report Nos. 811/258522, 811/2557078, and 236090.

4.4 GLASSWARE CALIBRATION - In-house procedure 3-QC-002 is used to calibrate all Class A Glassware used in the manufacture and quality control of Custom Grade Standards.

5.0 TRACE METALLIC IMPURITIES (TMI) DETERMINED BY ICP/MS AND ICP-OES IN µg/mL

Custom-Grade solutions are tested for trace metallic impurities by Axial ICP-OES and ICP-MS. The result from the most sensitive method for each element, is reported below. Solutions tested by ICP-MS were analyzed in an ULPA-Filtered Clean Room. An ULPA-Filter is 99.9985% efficient for the removal of particles down to 0.3 µm.

<u>O</u> Al < 0.00090	<u>M</u> Dy < 0.02400	<u>O</u> Li < 0.00003	<u>M</u> Pr < 0.00120	<u>M</u> Te < 0.11998
<u>M</u> Sb < 0.00200	<u>M</u> Er < 0.02000	<u>M</u> Lu < 0.00160	<u>M</u> Re < 0.00400	<u>M</u> Tb < 0.00120
<u>M</u> As < 0.03999	<u>M</u> Eu < 0.01200	<u>O</u> Mg 0.00100	<u>M</u> Rh < 0.00400	<u>M</u> Tl < 0.00400
<u>M</u> Ba < 0.03999	<u>M</u> Gd < 0.00400	<u>O</u> Mn < 0.00003	<u>M</u> Rb 0.49948	<u>M</u> Th < 0.00400
<u>O</u> Be < 0.00020	<u>M</u> Ga < 0.00400	<u>O</u> Hg < 0.01500	<u>M</u> Ru < 0.00800	<u>M</u> Tm < 0.00160
<u>M</u> Bi < 0.00160	<u>O</u> Ge < 0.00150	<u>M</u> Mo < 0.00800	<u>M</u> Sm < 0.00400	<u>M</u> Sn < 0.02000
<u>O</u> B < 0.00060	<u>O</u> Au < 0.00300	<u>M</u> Nd < 0.00800	<u>O</u> Sc < 0.00002	<u>O</u> Ti < 0.00070
<u>M</u> Cd < 0.01200	<u>M</u> Hf < 0.00800	<u>O</u> Ni < 0.00230	<u>O</u> Se < 0.05000	<u>M</u> W < 0.03999
<u>O</u> Ca 0.00075	<u>M</u> Ho < 0.00200	<u>M</u> Nb < 0.00200	<u>O</u> Si < 0.00340	<u>M</u> U < 0.00800
<u>M</u> Ce < 0.02000	<u>M</u> In < 0.03999	<u>n</u> Os	<u>M</u> Ag < 0.00800	<u>O</u> V < 0.00090
<u>M</u> Cs < 0.00120	<u>M</u> Ir < 0.02000	<u>M</u> Pd < 0.02000	<u>O</u> Na 0.21730	<u>M</u> Yb < 0.00400
<u>M</u> Cr < 0.02000	<u>O</u> Fe 0.00212	<u>O</u> P < 0.00250	<u>M</u> Sr < 0.00200	<u>M</u> Y < 0.15998
<u>M</u> Co < 0.01200	<u>M</u> La < 0.00200	<u>M</u> Pt < 0.00800	<u>O</u> S < 0.07200	<u>O</u> Zn 0.00050
<u>M</u> Cu < 0.02400	<u>M</u> Pb < 0.01200	<u>S</u> K	<u>M</u> Ta < 0.02800	<u>M</u> Zr < 0.02000

M - Checked by ICP-MS O - Checked by ICP-OES I - Spectral Interference n - Not Checked For s - Solution Standard Element

6.0 INTENDED USE

For the calibration of analytical instruments including but not limited to the following:
 ICP-MS, ICP-OES, FAAS, GFAA, XRF, and DCP
 For the validation of analytical methods
 For the preparation of "working reference samples"
 For interference studies and the determination of correction coefficients
 For detection limit and linearity studies
 For additional intended uses, contact IV Technical Staff

7.0 INSTRUCTIONS FOR THE CORRECT USE OF THIS REFERENCE MATERIAL

Storage & Handling - Keep tightly sealed when not in use. Store and use at 20 ± 4°C. Do not pipet from container. Do not return portions removed for pipetting to container.

Atomic Weight; Valence; Coordination Number; Chemical Form in Solution - 29.0983; +1; (6); K'(aq)
 (Coordination Number in parentheses is assumed, not certain.)

Chemical Compatibility - Soluble in HCl, HNO₃, H₂SO₄, and HF aqueous matrices. Avoid use of HClO₄ due to insolubility of the perchlorate. Stable with all metals and inorganic anions except ClO₄⁻.

Stability - 2-100 ppb levels stable for months in 1% HNO₃ / LDPE container. 1-10,000 ppm solutions chemically stable for years in 1-5% HNO₃ / LDPE container.

K Containing Samples (Preparation and Solution) - Metal (Dissolves very rapidly in water). Ores (Sodium carbonate fusion in Pt⁺ followed by HCl dissolution-blank levels of K in sodium carbonate critical). Organic Matrices (Sulfuric/peroxide digestion)

Atomic Spectroscopic Information (ICP-OES D.L.s are given as radial/axial view):

Technique/Line	Estimated D.L.	Order	Type	Interferences (underlined indicates severe at concs.)
ICP-OES 766.490 nm	0.4 / 0.001 µg/mL	1	atom	2 nd order radiation from R.E.s on some optical designs
ICP-OES 771.531 nm	1.0 / 0.03 µg/mL	1	atom	2 nd order radiation from R.E.s on some optical designs
ICP-OES 404.721 nm	1.1 / 0.05 µg/mL	1	atom	<u>U</u> , <u>Ce</u> , <u>³³ArH</u> , <u>²³Na</u> ¹⁰ O, <u>⁷⁸Se</u> ¹²
ICP-MS 39 amu	10 ppt	na	M'	

- 8.0 **HAZARDOUS INFORMATION** - Please refer to the enclosed Material Safety Data sheet for information regarding this CRM.
- 9.0 **HOMOGENEITY** - This solution was mixed according to procedure IV-MPM-004 and is guaranteed to be homogeneous.
- 10.0 **QUALITY STANDARD DOCUMENTATION**

010186

10.1 **ISO 9001:2000 Quality Management System Registration - QMI Certificate Number 010105**

Recognized by:

Registrar Accreditation Board (ANSI-RAB)
Standards Council of Canada (SCC)
Dutch Council for Accreditation (RVA)
Entidad Mexicana de Acreditacion, a.c.(EMA)

Members of IQ Net International Certification Network:

Argentina (IRAM), Australia (QAS), Austria (ÖQS), Belgium (Avinter), Brazil (FCAV), Canada (QMI), Hong Kong (HKQAA), Columbia (ICONTEC), Czech Republic (CQS), Denmark (DS), Finland (SFS), France (AFAQ), Germany (DQS), Greece (ELOT), Hungary (MSZT), Ireland (NSAI), Israel (SII), Italy (CISQ), Japan (JQA), Korea (KSA-QA), Netherlands (KEMA), Norway (NCS), Poland(PCBC), Portugal (APCER), Singapore (PSB), Slovenia (SIQ), Spain (AENOR), Switzerland (SQS)



10.2 **ISO/IEC 17025 - 1999 "General Requirements for the Competence of Testing and Calibration"**

- Chemical Testing - Accredited A2LA Certificate Number 883.01



10.3 **ISO/IEC Guide 34 - 2000 "General Requirements for the Competence of Reference Material Producers"**

- Reference Materials Production - Accredited A2LA Certificate Number 883.02

A2LA Mutual Recognition Agreement Partners:

Australia (NATA), Austria (BmWA), Belgium (BELTEST) (BKO-OBE), Canada (SCC), Chinese Taipei (CNLA), Czech Republic (NAO), Denmark (DANAK), Finland (FINAS), France (COFRAC), Germany (DAR), Hong Kong (HKAS), Ireland (NAB), Italy (SIT) (SINAL), Japan (JAB) (JNLA), Republic of Korea (KOLAS), The Netherlands (RvA), New Zealand (IANZ), Norway (NA), Portugal (IPQ), Singapore (SAC-SINGLAS), Spain (ENAC), Sweden (SWEDAC), Switzerland (SAS), United Kingdom (UKAS) and United States (NVLAP) (ICBO ES)

10.4 **10CFR50 Appendix B - Nuclear Regulatory Commission**

- Domestic Licensing of Production and Utilization Facilities

10.5 **10CFR21 - Nuclear Regulatory Commission - Reporting Defects and Non-Compliance**

10.6 **MIL-STD-45662A (Obsolete/Observed)**

11.0 **DATE OF CERTIFICATION AND PERIOD OF VALIDITY**



11.1 **IV Shelf Life** - The period of time during which the concentration of the analyte(s) in a properly packaged, unopened, and unused standard stored under environmentally controlled and monitored conditions will remain within the specified uncertainty range. Shelf life is limited primarily by transpiration (loss of water from the solution) and infrequently, by chemical instability. Transpiration studies (P-SP01020) of chemically-stable solutions performed at Inorganic Ventures / IV Labs indicate a CRM shelf-life of four years for solutions packaged in 500-mL low density polyethylene bottles. When stored under special conditions that minimize transpiration and instability, the shelf life can be extended past this limit.

11.2 **Expiration Date** - The date after which a CRM should not be used. Routine laboratory use of a CRM increases transpiration losses and the chance of contamination which affect the integrity of the CRM and limit its useful life. Inorganic Ventures / IV Labs concurs with state and federal regulatory agencies' recommendations that solution standards be assigned a one-year expiration date.

INORGANIC LABS/RADCHEM LABS
DATE RECEIVED: 11/5/03
DATE EXPIRED: 12/1/04
DATE OPENED: 11/5/03
INORG: 4320 PO: R52258

Certification Date: January 30, 2003

Expiration Date:

EXPIRES
122004

12.0 NAMES AND SIGNATURES OF CERTIFYING OFFICERS

Certificate Prepared By: JoAnn Struthers, QA Administrative Assistant

Certificate Approved By: Katalin Le, QC Supervisor

Certifying Officer: Paul Gaines, Chemist, Senior Technical Director

010187
JoAnn Struthers
Katalin Le
Paul Gaines

inorganic ventures / iv labs

195 lehigh avenue, suite 4, lakewood, nj 08701 usa
phone: 800-669-6799 • 732-901-1900 • fax: 732-901-1903
e-mail: ivsales@ivstandards.com • website: www.ivstandards.com

certificate of analysis

010188

1.0 Inorganic Ventures / IV Labs is an ISO Guide 34-2000 Certified Reference Material (CRM) Manufacturer: Certificate #883-02. The certificate is designed and the certified value(s) and uncertainty(ies) are determined in accordance with ISO Guide 31-2000 (Reference Materials - Contents of certificates and label(s)), ISO Guide 34-2000 "Quality System Guidelines for the Production of Reference Materials," and ISO Guide 35-1989 "Certification of Reference Materials - General and Statistical Principles."

2.0 DESCRIPTION OF CRM Custom-Grade 10000 µg/mL Magnesium In 1.4% (abs) HNO₃

Catalog Number: CGMG10-1 and CGMG10-5

Lot Number: T-MG03006

Starting Material: Mg metal

Starting Material Purity (%): 99.9968

Starting Material Lot No RML91191

Matrix: 1.4% (abs) HNO₃

INORGANIC LABS/RADCHEM LABS Pg. 1 of 2
DATE RECEIVED: 07/31/03
DATE EXPIRED: 08/01/2004
DATE OPENED: 08/01/03
INORG: 4204 PO: F52391

3.0 CERTIFIED VALUES AND UNCERTAINTIES

Certified Concentration: 9921 ± 20 µg/mL

Certified Density: 1.050 g/mL (measured at 22° C)

The Certified Value is the instrument analysis value. The following equations are used in the calculation of the certified value and the uncertainty:

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

$$\text{Uncertainty } (\pm) = \frac{2(\sum s_i)^{1/2}}{(n)^{1/2}}$$

($\bar{x}$) = mean

x_i = individual results

n = number of measurements

$\sum S$ = The summation of all significant estimated errors

(Most common are the errors from instrumental measurement, weighing, dilution to volume, and the fixed error reported on the NIST SRM certificate of analysis.)

4.0 TRACEABILITY TO NIST AND VALUES OBTAINED BY INDEPENDENT METHODS

"Property of the result of a measurement or the value of a standard whereby it can be related to stated references, usually national or international standards, through an unbroken chain of comparisons all having stated uncertainties." (ISO VIM, 2nd ed., 1993, definition 6.10)

This IV product is Traceable to NIST via direct comparison to NIST SRMs. The uncertainties for each certified value are reported, taking into account the SRM uncertainty error and the measurement, weighing and volume dilution errors.

4.1 Assay Method #1 9998 ± 20 µg/mL

EDTA NIST SRM 928 Lot Number: 880710

Assay Method #2 9921 ± 20 µg/mL

ICP Assay NIST SRM 3131a Lot Number: 991107

4.2 BALANCE CALIBRATION - All balances are checked daily using in-house procedure number 6-IMM-001. The weights used for testing are annually compared to Gerhart Scale Corporation's master weights and are traceable to the National Institute of Standards and Technology (NIST). The NIST Traceability numbers are 692476 - Class 1 and 692476A - Class 2. The NIST test number is 822/260017-98. All analytical balances are calibrated every 4 months by Gerhart Scale Corp. of South Amboy. The balances are calibrated with a class 1 and/or class 2 analytical weight set. These weights are tested annually by a NIST / NVLAP accredited calibration lab. The NIST test number is 822/260017-98.

010189

4.3 THERMOMETER CALIBRATION - The thermometers used in the determination of the final densities are calibrated vs standard thermometer No. 903-2680 which was certified in accordance with the procedures outlined by ASTM E77-87 and NIST Monograph 150 using NIST Test Nos. and Std Nos.: 769543, 217368/769543, 217368/P14452, 176240/P14452, 176240. The in-house procedure No. is 2-QC-001. Thermometers which are not calibrated vs standard thermometer No. 903-2680 are traceable to NIST Identification Nos. 92564, 119016, 471047 and NIST test report Nos. 811/258522, 811/2557078, and 236090.

4.4 GLASSWARE CALIBRATION - In-house procedure 3-QC-002 is used to calibrate all Class A Glassware used in the manufacture and quality control of Custom Grade Standards.

5.0 TRACE METALLIC IMPURITIES (TMI) DETERMINED BY ICP/MS AND ICP-OES IN µg/mL

Custom-Grade solutions are tested for trace metallic impurities by Axial ICP-OES and ICP-MS. The result from the most sensitive method for each element, is reported below. Solutions tested by ICP-MS were analyzed in an ULPA-Filtered Clean Room. An ULPA-Filter is 99.9985% efficient for the removal of particles down to 0.3 µm.

Q Al 0.02454	M Dy < 0.02455	Q Li 0.00797	M Pr < 0.00123	M Te < 0.12275
M Sb 0.00306	M Er < 0.02046	M Lu < 0.00164	M Re < 0.00409	M Tb < 0.00123
M As < 0.04092	M Eu < 0.01228	S Mg	M Rh < 0.00409	M Tl < 0.00409
M Ba < 0.04092	M Gd < 0.00409	M Mn < 0.01637	M Rb < 0.00409	M Th < 0.00409
Q Be < 0.00017	M Ga < 0.00409	Q Hg < 0.00900	M Ru < 0.00818	M Tm < 0.00164
M Bi < 0.00164	M Ge < 0.02455	M Mo < 0.00818	M Sm < 0.00409	M Sn < 0.02046
Q B 0.00871	M Au < 0.01228	M Nd < 0.00818	M Sc < 0.04092	Q Ti 0.10206
M Cd < 0.01228	M Hf < 0.00818	Q Ni 0.01404	M Se < 0.03273	M W < 0.04092
Q Ca 0.01070	M Ho < 0.00205	M Nb < 0.00205	Q Si 0.03186	M U < 0.00818
M Ce < 0.02046	M In < 0.04092	n Os	M Ag < 0.00818	M V < 0.00818
M Cs < 0.00123	M Ir < 0.02046	M Pd < 0.02046	Q Na 0.01817	M Yb < 0.00409
Q Cr 0.02315	Q Fe 0.02467	Q P < 0.01600	M Sr < 0.00205	M Y < 0.16367
M Co < 0.01228	M La < 0.00205	M Pt < 0.00818	n S	Q Zn 0.01892
Q Cu 0.00672	Q Pb 0.03236	Q K < 0.05000	M Ta < 0.02864	M Zr < 0.02046

M - Checked by ICP-MS O - Checked by ICP-OES I - Spectral Interference n - Not Checked For s - Solution Standard Element

6.0 INTENDED USE

For the calibration of analytical instruments including but not limited to the following:

ICP-MS, ICP-OES, FAAS, GFAA, XRF, and DCP

For the validation of analytical methods

For the preparation of "working reference samples"

For interference studies and the determination of correction coefficients

For detection limit and linearity studies

For additional intended uses, contact IV Technical Staff

7.0 INSTRUCTIONS FOR THE CORRECT USE OF THIS REFERENCE MATERIAL

Storage & Handling - Keep tightly sealed when not in use. Store and use at $20 \pm 4^\circ\text{C}$. Do not pipet from container. Do not return portions removed for pipetting to container.

Atomic Weight; Valence; Coordination Number; Chemical Form in Solution - 24.305; +2; 6; $\text{Mg}(\text{H}_2\text{O})_6^{+2}$

Chemical Compatibility - Soluble in HCl , HNO_3 , and H_2SO_4 ; avoid HF , H_3PO_4 , and neutral to basic media. Stable with most metals and inorganic anions forming insoluble silicates, carbonates, hydroxides, oxides, and tungstates in neutral and slightly acidic media.

Stability - 2-100 ppb levels stable for months in 1% HNO_3 / LDPE container. 1-10,000 ppm solutions chemically stable for years in 1-10% HNO_3 / LDPE container.

Mg Containing Samples (Preparation and Solution) - Metal (Best dissolved in diluted HNO_3); Oxide (Readily soluble in above compatible aqueous acidic solutions); Ores (Carbonate fusion in Pt° followed by HCl dissolution); Organic Matrices (Sulfuric / peroxide digestion or nitric / sulfuric / perchloric acid decomposition, or dry ash and dissolution in dilute HCl).

Atomic Spectroscopic Information (ICP-OES D.L.s are given as radial/axial view):

Technique/Line	Estimated D.L.	Order	Type	Interferences
ICP-OES 279.553 nm	0.0002 / 0.00003 $\mu\text{g/mL}$	1	ion	Th
ICP-OES 280.270 nm	0.0003 / 0.00005 $\mu\text{g/mL}$	1	ion	U, V
ICP-OES 285.213 nm	0.002 / 0.00003 $\mu\text{g/mL}$	1	atom	U, Hf, Cr, Zr
ICP-MS 24 amu	42 ppt	n/a	M	$^6\text{Li}^{+10}$, $^{47}\text{Ti}^{+2}$, $^{40}\text{Ca}^{+2}$

010190

8.0 HAZARDOUS INFORMATION - Please refer to the enclosed Material Safety Data sheet for information regarding this CRM.

9.0 HOMOGENEITY - This solution was mixed according to procedure IV-MPM-004 and is guaranteed to be homogeneous.

10.0 QUALITY STANDARD DOCUMENTATION

10.1 ISO 9001:2000 Quality Management System Registration - QMI Certificate Number 010105

Recognized by:

Registrar Accreditation Board (ANSI-RAB)

Standards Council of Canada (SCC)

Dutch Council for Accreditation (RVA)

Entidad Mexicana de Acreditacion, a.c.(EMA)

Members of **IQ Net International Certification Network**:

Argentina (IRAM), Australia (QAS), Austria (ÖQS), Belgium (Avinter), Brazil (FCAV), Canada (QMI), Hong Kong (HKQAA), Columbia (ICONTEC), Czech Republic (CQS), Denmark (DS), Finland (SFS), France (AFAQ), Germany (DQS), Greece (ELOT), Hungary (MSZT), Ireland (NSAI), Israel (SII), Italy (CISQ), Japan (JQA), Korea (KSA-QA), Netherlands (KEMA), Norway (NCS), Poland (PCBC), Portugal (APCER), Singapore (PSB), Slovenia (SIQ), Spain (AENOR), Switzerland (SQS)

10.2 ISO/IEC 17025 - 1999 "General Requirements for the Competence of Testing and Calibration"

- Chemical Testing - Accredited A2LA Certificate Number 883.01



10.3 ISO/IEC Guide 34 - 2000 "General Requirements for the Competence of Reference Material Producers"

- Reference Materials Production - Accredited A2LA Certificate Number 883.02

A2LA Mutual Recognition Agreement Partners:

Australia (NATA), Austria (BmWA), Belgium (BELTEST) (BKO-OBE), Canada (SCC), Chinese Taipei (CNLA), Czech Republic (NAO), Denmark (DANAK), Finland (FINAS), France (COFRAC), Germany (DAR), Hong Kong (HKAS), Ireland (NAB), Italy (SIT) (SINAL), Japan (JAB) (JNLA), Republic of Korea (KOLAS), The Netherlands (RvA), New Zealand (IANZ), Norway (NA), Portugal (IPQ), Singapore (SAC-SINGLAS), Spain (ENAC), Sweden (SWEDAC), Switzerland (SAS), United Kingdom (UKAS) and United States (NVLAP) (ICBO ES)

10.4 10CFR50 Appendix B - Nuclear Regulatory Commission


- Domestic Licensing of Production and Utilization Facilities

10.5 10CFR21 - Nuclear Regulatory Commission - Reporting Defects and Non-Compliance

10.6 MIL-STD-45662A (Obsolete/Observed)

INORGANIC LABS/RADCHEM LABS Pg. 2 of 2
DATE RECEIVED: 07/31/03
DATE EXPIRED: 08/01/2004
DATE OPENED: 08/01/03
INORG: 4204 PO: E52391

11.0 DATE OF CERTIFICATION AND PERIOD OF VALIDITY



11.1 IV Shelf Life - The period of time during which the concentration of the analyte(s) in a properly packaged, unopened, and unused standard stored under environmentally controlled and monitored conditions will remain within the specified uncertainty range. Shelf life is limited primarily by transpiration (loss of water from the solution) and infrequently, by chemical instability. Transpiration studies (P-SP01020) of chemically-stable solutions performed at Inorganic Ventures / IV Labs indicate a CRM shelf-life of four years for solutions packaged in 500-mL low density polyethylene bottles. When stored under special conditions that minimize transpiration and instability, the shelf life can be extended past this limit.

11.2 Expiration Date - The date after which a CRM should not be used. Routine laboratory use of a CRM increases transpiration losses and the chance of contamination which affect the integrity of the CRM and limit its useful life. Inorganic Ventures / IV Labs concurs with state and federal regulatory agencies' recommendations that solution standards be assigned a one-year expiration date.

Certification Date: August 28, 2002

Expiration Date: **EXPIRES**
01 02 00 4

12.0 NAMES AND SIGNATURES OF CERTIFYING OFFICERS

Certificate Prepared By: Debbie Newman, QA Administrator

Certificate Approved By: Katalin Le, QC Supervisor

Certifying Officer: Paul Gaines, Chemist, Senior Technical Director

Debbie Newman
Katalin Le
Paul Gaines



inorganic ventures / iv labs

195 lehigh avenue, suite 4, lakewood, nj 08701 usa
phone: 800-669-6799 • 732-901-1900 • fax: 732-901-1903
e-mail: ivsales@ivstandards.com • website: www.ivstandards.com

010192 certificate of analysis

1.0 Inorganic Ventures / IV Labs is an ISO Guide 34-2000 Certified Reference Material (CRM) Manufacturer: Certificate #883-02. The certificate is designed and the certified value(s) and uncertainty(ies) are determined in accordance with ISO Guide 31-2000 (Reference Materials - Contents of certificates and label(s)), ISO Guide 34-2000 "Quality System Guidelines for the Production of Reference Materials," and ISO Guide 35-1989 "Certification of Reference Materials - General and Statistical Principles."

2.0 DESCRIPTION OF CRM Custom-Grade 10000 µg/mL Sodium In 1.4% (abs) HNO₃

Catalog Number: CGNA10-1, CGNA10-2, and CGNA10-5
Lot Number: T-NA03006
Starting Material: Na₂CO₃
Starting Material Purity (%): 99.999936
Starting Material Lot No 42095
Matrix: 1.4% (abs) HNO₃

INORGANIC LABS/RADCHEM LABS Pg 1 of 2
DATE RECEIVED: 07/31/03
DATE EXPIRED: 08/01/2004
DATE OPENED: 08/01/03
INORG: 4205 PO: F52391

3.0 CERTIFIED VALUES AND UNCERTAINTIES

Certified Concentration: 10,005 ± 7 µg/mL

Certified Density: 1.032 g/mL (measured at 22° C)

The Certified Value is based upon the most precise method used to analyze this CRM. The following equations are used in the calculation of the certified value and the uncertainty:

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

($\bar{x}$) = mean

x_i = individual results

n = number of measurements

$$\text{Uncertainty } (\pm) = \frac{2[(\sum s_i)^2]^{1/2}}{(n)^{1/2}}$$

$\sum s_i$ = The summation of all significant estimated errors.

(Most common are the errors from instrumental measurement, weighing, dilution to volume, and the fixed error reported on the NIST SRM certificate of analysis.)

The independent samples t-test was used to determine if there is agreement between the above assay methods at the 95% confidence interval. Both methods were compared and showed agreement within the stated uncertainties. This agreement is a confirmation of the accuracy of this CRM.

4.0 TRACEABILITY TO NIST AND VALUES OBTAINED BY INDEPENDENT METHODS

□ "Property of the result of a measurement or the value of a standard whereby it can be related to stated references, usually national or international standards, through an unbroken chain of comparisons all having stated uncertainties." (ISO VIM, 2nd ed., 1993, definition 6.10)

□ This IV product is Traceable to NIST via direct comparison to NIST SRMs. The uncertainties for each certified value are reported, taking into account the SRM uncertainty error and the measurement, weighing and volume dilution errors.

4.1 Assay Method #1 10,067 ± 75 µg/mL

ICP Assay NIST SRM 3152a Lot Number: 990907

Assay Method #2 10,005 ± 7 µg/mL

Gravimetric NIST SRM Lot Number: See Sec. 4.2

4.2 BALANCE CALIBRATION - All balances are checked daily using in-house procedure number 6-IMM-001. The weights used for testing are annually compared to Gerhart Scale Corporation's master weights and are traceable to the National Institute of Standards and Technology (NIST). The NIST Traceability numbers are 692476 - Class 1 and 692476A - Class 2. The NIST test number is 822/260017-98. All analytical balances are calibrated every 4 months by Gerhart Scale Corp. of South Amboy. The balances are calibrated with a class 1 and/or class 2 analytical weight set. These weights are tested annually by a NIST / NVLAP accredited calibration lab. The NIST test number is 822/260017-98.

010193

4.3 THERMOMETER CALIBRATION - The thermometers used in the determination of the final densities are calibrated vs standard thermometer No. 903-2680 which was certified in accordance with the procedures outlined by ASTM E77-87 and NIST Monograph 150 using NIST Test Nos. and Std Nos.: 769543, 217368/769543, 217368/P14452, 176240/P14452, 176240. The in-house procedure No. is 2-QC-001. Thermometers which are not calibrated vs standard thermometer No. 903-2680 are traceable to NIST Identification Nos. 92564, 119016, 471047 and NIST test report Nos. 811/258522, 811/2557078, and 236090.

4.4 GLASSWARE CALIBRATION - In-house procedure 3-QC-002 is used to calibrate all Class A Glassware used in the manufacture and quality control of Custom Grade Standards.

5.0 TRACE METALLIC IMPURITIES (TMI) DETERMINED BY ICP/MS AND ICP-OES IN µg/mL

Custom-Grade solutions are tested for trace metallic impurities by Axial ICP-OES and ICP-MS. The result from the most sensitive method for each element, is reported below. Solutions tested by ICP-MS were analyzed in an ULPA-Filtered Clean Room. An ULPA-Filter is 99.9985% efficient for the removal of particles down to 0.3 µm.

<u>Q</u> Al < 0.00090	<u>M</u> Dy < 0.02499	<u>Q</u> Li < 0.00003	<u>M</u> Pr < 0.00125	<u>M</u> Te < 0.12494
<u>M</u> Sb < 0.00208	<u>M</u> Er < 0.02082	<u>M</u> Lu < 0.00167	<u>M</u> Re < 0.00417	<u>M</u> Tb < 0.00125
<u>M</u> As < 0.04165	<u>M</u> Eu < 0.01249	<u>Q</u> Mg 0.00015	<u>M</u> Rh < 0.00417	<u>M</u> Tl < 0.00417
<u>M</u> Ba < 0.04165	<u>M</u> Gd < 0.00417	<u>Q</u> Mn < 0.00003	<u>M</u> Rb < 0.00417	<u>M</u> Th < 0.00417
<u>Q</u> Be < 0.00020	<u>M</u> Ga < 0.00417	<u>Q</u> Hg < 0.01500	<u>M</u> Ru < 0.00833	<u>M</u> Tm < 0.00167
<u>M</u> Bi < 0.00167	<u>Q</u> Ge < 0.00150	<u>M</u> Mo < 0.00833	<u>M</u> Sm < 0.00417	<u>M</u> Sn < 0.02082
<u>Q</u> B < 0.00060	<u>Q</u> Au < 0.00300	<u>M</u> Nd < 0.00833	<u>Q</u> Sc < 0.00002	<u>Q</u> Ti < 0.00070
<u>M</u> Cd < 0.01249	<u>M</u> Hf < 0.00833	<u>Q</u> Ni < 0.00230	<u>Q</u> Se < 0.05000	<u>M</u> W < 0.04165
<u>Q</u> Ca 0.00160	<u>M</u> Ho < 0.00208	<u>M</u> Nb < 0.00208	<u>Q</u> Si < 0.00340	<u>M</u> U < 0.00833
<u>M</u> Ce < 0.02082	<u>M</u> In < 0.04165	<u>n</u> Os	<u>M</u> Ag < 0.00833	<u>Q</u> V < 0.00090
<u>M</u> Cs 0.00104	<u>M</u> Ir < 0.02082	<u>M</u> Pd < 0.02082	<u>S</u> Na	<u>M</u> Yb < 0.00417
<u>M</u> Cr < 0.02082	<u>Q</u> Fe < 0.00110	<u>Q</u> P < 0.04000	<u>M</u> Sr < 0.00208	<u>M</u> Y < 0.16658
<u>M</u> Co < 0.01249	<u>M</u> La < 0.00208	<u>M</u> Pt < 0.00833	<u>Q</u> S < 0.07200	<u>Q</u> Zn 0.00130
<u>Q</u> Cu < 0.00140	<u>M</u> Pb < 0.01249	<u>Q</u> K 0.00873	<u>M</u> Ta < 0.02915	<u>M</u> Zr < 0.02082

M - Checked by ICP-MS

Q - Checked by ICP-OES

i - Spectral Interference

n - Not Checked For

s - Solution Standard Element

6.0 INTENDED USE

For the calibration of analytical instruments including but not limited to the following:

ICP-MS, ICP-OES, FAAS, GFAA, XRF, and DCP

For the validation of analytical methods

For the preparation of "working reference samples"

For interference studies and the determination of correction coefficients

For detection limit and linearity studies

For additional intended uses, contact IV Technical Staff

7.0 INSTRUCTIONS FOR THE CORRECT USE OF THIS REFERENCE MATERIAL

Storage & Handling - Keep tightly sealed when not in use. Store and use at 20 ± 4°C. Do not pipet from container. Do not return portions removed for pipetting to container.

Atomic Weight; Valence; Coordination Number; Chemical Form in Solution - 22.98977; +1; (6); Na⁺(aq) largely ionic in nature (Coordination Number in parentheses is assumed, not certain.)

Chemical Compatibility - Soluble in HCl, HNO₃, H₂SO₄, and HF aqueous matrices. Stable with all metals and inorganic anions.

Stability - 2-100 ppb levels stable for months in 1% HNO₃ / LDPE container. 1-10,000 ppm solutions chemically stable for years in 1-5% HNO₃ / LDPE container.

Na Containing Samples (Preparation and Solution) - Metal (Dissolves very rapidly in water). Ores (Lithium carbonate fusion in graphite crucible followed by HCl dissolution - blank levels of Na in lithium carbonate critical). Organic Matrices (Sulfuric / peroxide digestion or nitric/sulfuric/perchloric acid decomposition).

Atomic Spectroscopic Information (ICP-OES D.L.s are given as radial/axial view):

Technique/Line	Estimated D.L.	Order	Type	Interferences (underlined indicates severe at 1000°C)
ICP-OES 589.595 nm	0.07 / 0.00009 µg/mL	1	atom	2 nd order radiation from R.E.s on some optical designs
ICP-OES 588.995 nm	0.03 / 0.006 µg/mL	1	atom	2 nd order radiation from R.E.s on some optical designs
ICP-OES 330.237 nm	2.0 / 0.09 µg/mL	1	atom	<u>Pd, Zn</u>
ICP-MS 23 amu	310 ppt	n/a	M ⁺	⁴⁵ Ti ¹² , ⁴⁴ Ca ¹²

- 8.0 **HAZARDOUS INFORMATION** - Please refer to the enclosed Material Safety Data sheet for information regarding this CRM.
- 9.0 **HOMOGENEITY** - This solution was mixed according to procedure IV-MPM-004 and is guaranteed to be homogeneous.
- 10.0 **QUALITY STANDARD DOCUMENTATION**

010194

10.1 **ISO 9001:2000 Quality Management System Registration - QMI Certificate Number 010105**

Recognized by:

Registrar Accreditation Board (ANSI-RAB)
Standards Council of Canada (SCC)
Dutch Council for Accreditation (RVA)
Entidad Mexicana de Acreditacion, a.c.(EMA)

Members of IQ Net International Certification Network:

Argentina (IRAM), Australia (QAS), Austria (ÖQS), Belgium (Avinter), Brazil (FCAV), Canada (QMI), Hong Kong (HKQAA), Columbia (ICONTEC), Czech Republic (CQS), Denmark (DS), Finland (SFS), France (AFAQ), Germany (DQS), Greece (ELOT), Hungary (MSZT), Ireland (NSAI), Israel (SII), Italy (CISQ), Japan (JQA), Korea (KSA-QA), Netherlands (KEMA), Norway (NCS), Poland(PCBC), Portugal (APCER), Singapore (PSB), Slovenia (SIQ), Spain (AENOR), Switzerland (SQS)

10.2 **ISO/IEC 17025 - 1999 "General Requirements for the Competence of Testing and Calibration"**

- Chemical Testing - Accredited A2LA Certificate Number 883.01



10.3 **ISO/IEC Guide 34 - 2000 "General Requirements for the Competence of Reference Material Producers"**

- Reference Materials Production - Accredited A2LA Certificate Number 883.02

A2LA Mutual Recognition Agreement Partners:

Australia (NATA), Austria (BmWA), Belgium (BELTEST) (BKO-OBE), Canada (SCC), Chinese Taipei (CNLA), Czech Republic (NAO), Denmark (DANAK), Finland (FINAS), France (COFRAC), Germany (DAR), Hong Kong (HKAS), Ireland (NAB), Italy (SIT) (SINAL), Japan (JAB) (JNLA), Republic of Korea (KOLAS), The Netherlands (RvA), New Zealand (IANZ), Norway (NA), Portugal (IPQ), Singapore (SAC-SINGLAS), Spain (ENAC), Sweden (SWEDAC), Switzerland (SAS), United Kingdom (UKAS) and United States (NVLAP) (ICBO ES)

10.4 **10CFR50 Appendix B - Nuclear Regulatory Commission**

- Domestic Licensing of Production and Utilization Facilities

10.5 **10CFR21 - Nuclear Regulatory Commission - Reporting Defects and Non-Compliance**

10.6 **MIL-STD-45662A (Obsolete/Observed)**

11.0 **DATE OF CERTIFICATION AND PERIOD OF VALIDITY**



11.1 **IV Shelf Life** - The period of time during which the concentration of the analyte(s) in a properly packaged, unopened, and unused standard stored under environmentally controlled and monitored conditions will remain within the specified uncertainty range. Shelf life is limited primarily by transpiration (loss of water from the solution) and infrequently, by chemical instability. Transpiration studies (P-SP01020) of chemically-stable solutions performed at Inorganic Ventures / IV Labs indicate a CRM shelf-life of four years for solutions packaged in 500-mL low density polyethylene bottles. When stored under special conditions that minimize transpiration and instability, the shelf life can be extended past this limit.

11.2 **Expiration Date** - The date after which a CRM should not be used. Routine laboratory use of a CRM increases transpiration losses and the chance of contamination which affect the integrity of the CRM and limit its useful life. Inorganic Ventures / IV Labs concurs with state and federal regulatory agencies' recommendations that solution standards be assigned a one-year expiration date.

Certification Date: January 24, 2003

Expiration Date:

EXPIRES
01/24/04

INORGANIC LABS/RADCHEM LABS P. 2042
DATE RECEIVED: 07/31/03
DATE EXPIRED: 08/01/2004
DATE OPENED: 08/01/03
INORG: 4205 PO: F52391

12.0 NAMES AND SIGNATURES OF CERTIFYING OFFICERS

010195

Certificate Prepared By: Debbie Newman, LIMS Administrator

Debbie Newman
Katalin Le

Certificate Approved By: Katalin Le, QC Supervisor

Certifying Officer: Paul Gaines, Chemist, Senior Technical Director

Paul Gaines



Certificate of Analysis

010196

CUSTOM-GRADE SOLUTION

1000 µg/mL Lithium in 0.1% HNO₃ (abs)

Catalog Number: CGLI1-1, CGLI1-2 and CGLI1-5

Lot Number: **W-LI02066**

Starting Material:
Starting Material Purity:
Starting Material Lot No:

Li₂CO₃
99.999%
1053

INORGANIC LABS/RADCHEM LABS

DATE RECEIVED: 06/20/03

DATE EXPIRED: 07/01/2004

DATE OPENED: 06/23/03

INORG: 4149 PO: F52370

CERTIFIED CONCENTRATION: 998 ± 2 µg/mL

The Certified Value is based upon the most precise method used to analyze this CRM. The following equations are used in the calculation of the certified value and the uncertainty:

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

$$\text{Uncertainty } (\pm) = \frac{2(\sum s_i)^{2/3}}{(n)^{1/2}}$$

($\bar{x}$) = mean

x_i = individual results

n = number of measurements

$\sum S_i$ = The summation of all significant estimated errors.

Classical Wet Assay: 998 ± 2 µg/mL

Method: Gravimetric as the Sulfate vs NIST weights #822/254143-94.

Instrument Analysis: 1000 ± 4 µg/mL

Method: Inductively Coupled Plasma Spectroscopy (ICP) vs NIST SRM 3129a.

The independent samples t-test was used to determine if there is agreement between the above assay methods at the 95% confidence interval. Both methods were compared and showed agreement within the stated uncertainties. This agreement is a confirmation of the accuracy of this CRM.

TRACE METALLIC IMPURITIES DETERMINED BY ICP-MS AND ICP-OES IN µg/mL:

Custom Grade solutions tested for trace metallic impurities by ICP-MS were analyzed in an ULPA-Filtered Clean Room.

An ULPA-Filter is 99.9985% efficient for the removal of particles down to 0.3 µm.

<u>O</u> Al <0.010	<u>M</u> Dy <0.00060	<u>s</u> Li	<u>M</u> Pr <0.000030	<u>O</u> Te <0.0090
<u>M</u> Sb <0.000050	<u>M</u> Er <0.00050	<u>M</u> Lu <0.000040	<u>M</u> Re <0.00010	<u>M</u> Tb <0.000030
<u>O</u> As <0.044	<u>M</u> Eu <0.00030	<u>O</u> Mg <0.00010	<u>M</u> Rh <0.00010	<u>M</u> Tl <0.00010
<u>M</u> Ba <0.0010	<u>M</u> Gd <0.00010	<u>O</u> Mn <0.00020	<u>M</u> Rb <0.00010	<u>M</u> Th <0.00010
<u>O</u> Be <0.000050	<u>M</u> Ga <0.00010	<u>O</u> Hg <0.0070	<u>M</u> Ru <0.00020	<u>M</u> Tm <0.000040
<u>M</u> Bi <0.000040	<u>M</u> Ge <0.00060	<u>M</u> Mo <0.00020	<u>M</u> Sm <0.00010	<u>M</u> Sn <0.00050
<u>O</u> B <0.0060	<u>O</u> Au <0.010	<u>M</u> Nd <0.00020	<u>M</u> Sc <0.0010	<u>O</u> Ti <0.00030
<u>O</u> Cd <0.0018	<u>M</u> Hf <0.00020	<u>O</u> Ni <0.0040	<u>O</u> Se <0.020	<u>M</u> W <0.0010
<u>O</u> Ca 0.051	<u>M</u> Ho <0.000050	<u>M</u> Nb <0.000050	<u>O</u> Si 0.023	<u>M</u> U <0.00020
<u>M</u> Ce <0.00050	<u>O</u> In <0.030	<u>n</u> Os	<u>O</u> Ag <0.0040	<u>O</u> V <0.0010
<u>M</u> Cs 0.0018	<u>M</u> Ir <0.00050	<u>M</u> Pd <0.00050	<u>O</u> Na <0.10	<u>M</u> Yb <0.00010
<u>O</u> Cr <0.0020	<u>O</u> Fe <0.0020	<u>O</u> P <0.030	<u>O</u> Sr <0.0010	<u>M</u> Y <0.0040
<u>M</u> Co <0.00030	<u>M</u> La <0.000050	<u>M</u> Pt <0.00020	<u>O</u> S <0.050	<u>O</u> Zn <0.030
<u>M</u> Cu <0.00060	<u>M</u> Pb <0.00030	<u>O</u> K 0.0070	<u>M</u> Ta <0.00070	<u>M</u> Zr <0.00050

M - checked by ICP-MS

O - checked by ICP-OES

i - spectral interference

n - not checked for

s - solution standard element

ANALYZED DENSITY OF SOLUTION (measured at 22°C): 1.004 g/mL

(over)

QA:KL Rev.0224090W

Inorganic Ventures, Inc.

195 Lehigh Avenue • Suite 4 • Lakewood, NJ 08701

Orders: 800-669-6799 • FAX (732) 901-1903

Technical Support: 800-569-6799

Quality Assurance Manager

EXPIRES

01/22/04



Certificate of Analysis

010197**CUSTOM-GRADE SOLUTION****1000 µg/mL Barium in 0.1% HNO₃ (abs)**

Catalog Number: CGBA1-1, CGBA1-2, and CGBA1-5

Lot Number: **T-BA02021**Starting Material:
Starting Material Purity:
Starting Material Lot No:Ba(NO₃)₂
99.999%
21879

INORGANIC LABS/RADCHEM LABS

DATE RECEIVED: 3/03/03
DATE EXPIRED: 3/01/2004
DATE OPENED: 3/04/03
INORG: 3971 PO: F520064**CERTIFIED CONCENTRATION: 1003 ± 2 µg/mL**

The Certified Value is based upon the most precise method used to analyze this CRM. The following equations are used in the calculation of the certified value and the uncertainty:

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

$$\text{Uncertainty } (\pm) = \frac{2(\sum s_i^2)^{1/2}}{(n)^{1/2}}$$

 $(\bar{x})$ = mean x_i = individual results n = number of measurements $\sum S_i$ = The summation of all significant estimated errors.**Classical Wet Assay: 1003 ± 2 µg/mL**

Method: Gravimetric as the Sulfate vs NIST weights #822/254143-94.

Instrument Analysis: 1002 ± 3 µg/mL

Method: Inductively Coupled Plasma Spectroscopy (ICP) vs NIST SRM 3104a.

The independent samples t-test was used to determine if there is agreement between the above assay methods at the 95% confidence interval. Both methods were compared and showed agreement within the stated uncertainties. This agreement is a confirmation of the accuracy of this CRM.

TRACE METALLIC IMPURITIES DETERMINED BY ICP-MS AND ICP-OES IN µg/mL:

Custom-Grade solutions tested for trace metallic impurities by ICP-MS were analyzed in an ULPA-Filtered Clean Room.

An ULPA-Filter is 99.9985% efficient for the removal of particles down to 0.3 µm.

<u>Q</u>	Al	<0.040	<u>M</u>	Dy	<0.00060	<u>M</u>	Li	<0.0010	<u>n</u>	Pr		<u>i</u>	Te	
<u>M</u>	Sb	<0.000050	<u>M</u>	Er	<0.00050	<u>M</u>	Lu	<0.000040	<u>M</u>	Re	<0.00010	<u>M</u>	Tb	<0.000030
<u>M</u>	As	<0.0010	<u>M</u>	Eu	<0.00030	<u>M</u>	Mg	<0.0030	<u>M</u>	Rh	<0.00010	<u>M</u>	Tl	<0.00010
<u>s</u>	Ba		<u>M</u>	Gd	0.0020	<u>M</u>	Mn	<0.00040	<u>M</u>	Rb	<0.00010	<u>M</u>	Th	<0.00010
<u>M</u>	Be	<0.000050	<u>M</u>	Ga	<0.00010	<u>i</u>	Hg		<u>M</u>	Ru	<0.00020	<u>M</u>	Tm	<0.000040
<u>M</u>	Bi	<0.000040	<u>M</u>	Ge	0.0010	<u>M</u>	Mo	<0.00020	<u>M</u>	Sm	<0.00010	<u>M</u>	Sn	<0.00050
<u>M</u>	B	<0.0070	<u>M</u>	Au	<0.00030	<u>M</u>	Nd	0.0020	<u>M</u>	Sc	<0.0010	<u>M</u>	Ti	<0.0050
<u>M</u>	Cd	<0.00030	<u>M</u>	Hf	<0.00020	<u>i</u>	Ni		<u>Q</u>	Se	<0.40	<u>M</u>	W	<0.0010
<u>Q</u>	Ca	<0.010	<u>M</u>	Ho	<0.000050	<u>M</u>	Nb	<0.000050	<u>Q</u>	Si	<0.020	<u>M</u>	U	<0.00020
<u>n</u>	Ce		<u>i</u>	In		<u>n</u>	Os		<u>M</u>	Ag	<0.00020	<u>M</u>	V	<0.00020
<u>n</u>	Cs		<u>M</u>	Ir	<0.00050	<u>M</u>	Pd	<0.00050	<u>Q</u>	Na	<0.090	<u>M</u>	Yb	<0.00010
<u>M</u>	Cr	<0.00050	<u>Q</u>	Fe	<0.050	<u>Q</u>	P	<0.050	<u>M</u>	Sr	0.0040	<u>M</u>	Y	<0.0040
<u>M</u>	Co	<0.00030	<u>n</u>	La		<u>M</u>	Pt	<0.00020	<u>n</u>	S		<u>Q</u>	Zn	<0.030
<u>M</u>	Cu	<0.00060	<u>M</u>	Pb	<0.00030	<u>n</u>	K		<u>M</u>	Ta	<0.00070	<u>M</u>	Zr	<0.00050

M - checked by ICP-MS

O - checked by ICP-OES

i - spectral interference

n - not checked for

s - solution standard element

ANALYZED DENSITY OF SOLUTION (measured at 22°C): 1.000 g/mL**(over)**

QA:KL Rev. 112702DN

Inorganic Ventures, Inc.

195 Lehigh Avenue • Suite 4 • Lakewood, NJ 08701

Orders: 800-669-6799 • FAX (732) 901-1903

Technical Support: 800-569-6799

Quality Assurance Manager

EXPIRES
1/2004



inorganic ventures / iv labs

195 lehigh avenue, suite 4, lakewood, nj 08701 usa
phone: 800-669-6799 • 732-901-1900 • fax: 732-901-1903
e-mail: ivsales@ivstandards.com • website: www.ivstandards.com

010198

certificate of analysis

1.0 Inorganic Ventures / IV Labs is an ISO Guide 34-2000 Certified Reference Material (CRM) Manufacturer: Certificate #883-02. The certificate is designed and the certified value(s) and uncertainty(ies) are determined in accordance with ISO Guide 31-2000 (Reference Materials - Contents of certificates and label(s)), ISO Guide 34-2000 "Quality System Guidelines for the Production of Reference Materials," and ISO Guide 35-1989 "Certification of Reference Materials - General and Statistical Principles."

2.0 DESCRIPTION OF CRM Custom-Grade 1000 µg/mL Beryllium In 2% (abs) HNO₃

Catalog Number: CGBE1-1, CGBE1-2, and CGBE1-5
Lot Number: W-BE01100
Starting Material: Be(OOCCH₃)₂
Starting Material Purity (%): 99.999897
Starting Material Lot No 01-10-01
Matrix: 2% (abs) HNO₃

3.0 CERTIFIED VALUES AND UNCERTAINTIES

Certified Concentration: 991 ± 3 µg/mL

Certified Density: 1.022 g/mL (measured at 22° C)

The Certified Value is the instrument analysis value. The following equations are used in the calculation of the certified value and the uncertainty:

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

$$\text{Uncertainty } (\pm) = \frac{2[(\sum s_i)^2]^{1/2}}{(n)^{1/2}}$$

($\bar{x}$) = mean

x_i = individual results

n = number of measurements

$\sum s_i$ = The summation of all significant estimated errors.

(Most common are the errors from instrumental measurement, weighing, dilution to volume, and the fixed error reported on the NIST SRM certificate of analysis.)

4.0 TRACEABILITY TO NIST AND VALUES OBTAINED BY INDEPENDENT METHODS

☐ "Property of the result of a measurement or the value of a standard whereby it can be related to stated references, usually national or international standards, through an unbroken chain of comparisons all having stated uncertainties." (ISO VIM, 2nd ed., 1993, definition 6.10)

☐ This IV product is Traceable to NIST via direct comparison to NIST SRMs. The uncertainties for each certified value are reported, taking into account the SRM uncertainty error and the measurement, weighing and volume dilution errors.

4.1 Assay Method #1 991 ± 3 µg/mL (Avg 2 runs)

ICP Assay NIST SRM 3105a Lot Number: 892707

INORGANIC LABS/RADCHEM LABS Pg. 1 of 2

DATE RECEIVED: 05/09/03

DATE EXPIRED: 06/01/2004

DATE OPENED: 05/12/03

INORG: 4062 PO: F52181

4.2 BALANCE CALIBRATION - All balances are checked daily using in-house procedure number 6-IMM-001. The weights used for testing are annually compared to Gerhart Scale Corporation's master weights and are traceable to the National Institute of Standards and Technology (NIST). The NIST Traceability numbers are 692476 - Class 1 and 692476A - Class 2. The NIST test number is 822/260017-98. All analytical balances are calibrated every 4 months by Gerhart Scale Corp. of South Amboy. The balances are calibrated with a class 1 and/or class 2 analytical weight set. These weights are tested annually by a NIST / NVLAP accredited calibration lab. The NIST test number is 822/260017-98.

010199

4.3 THERMOMETER CALIBRATION - The thermometers used in the determination of the final densities are calibrated vs standard thermometer No. 903-2680 which was certified in accordance with the procedures outlined by ASTM E77-87 and NIST Monograph 150 using NIST Test Nos. and Std Nos.: 769543, 217368/769543, 217368/P14452, 176240/P14452, 176240. The in-house procedure No. is 2-QC-001. Thermometers which are not calibrated vs standard thermometer No. 903-2680 are traceable to NIST Identification Nos. 92564, 119016, 471047 and NIST test report Nos. 811/258522, 811/2557078, and 236090.

4.4 GLASSWARE CALIBRATION - In-house procedure 3-QC-002 is used to calibrate all Class A Glassware used in the manufacture and quality control of Custom Grade Standards.

5.0 TRACE METALLIC IMPURITIES (TMI) DETERMINED BY ICP/MS AND ICP-OES IN µg/mL

Custom-Grade solutions are tested for trace metallic impurities by Axial ICP-OES and ICP-MS. The result from the most sensitive method for each element, is reported below. Solutions tested by ICP-MS were analyzed in an ULPA-Filtered Clean Room. An ULPA-Filter is 99.9985% efficient for the removal of particles down to 0.3 µm.

<u>O</u> Al < 0.00800	<u>M</u> Dy < 0.01305	<u>O</u> Li < 0.00002	<u>M</u> Pr < 0.00065	<u>M</u> Te < 0.06525
<u>M</u> Sb < 0.00109	<u>M</u> Er < 0.01087	<u>M</u> Lu < 0.00087	<u>M</u> Re < 0.00218	<u>M</u> Tb < 0.00065
<u>M</u> As < 0.02175	<u>M</u> Eu < 0.00652	<u>O</u> Mg < 0.00003	<u>M</u> Rh < 0.00218	<u>M</u> Tl < 0.00218
<u>M</u> Ba < 0.02175	<u>M</u> Gd < 0.00218	<u>O</u> Mn < 0.00002	<u>M</u> Rb < 0.00218	<u>M</u> Th < 0.00218
<u>s</u> Be	<u>M</u> Ga < 0.00218	<u>O</u> Hg < 0.01500	<u>M</u> Ru < 0.00435	<u>M</u> Tm < 0.00087
<u>M</u> Bi < 0.00087	<u>M</u> Ge < 0.01305	<u>M</u> Mo < 0.00435	<u>M</u> Sm < 0.00218	<u>M</u> Sn < 0.01087
<u>O</u> B < 0.01200	<u>M</u> Au < 0.00652	<u>M</u> Nd < 0.00435	<u>O</u> Sc < 0.00009	<u>M</u> Ti < 0.10874
<u>M</u> Cd < 0.00652	<u>M</u> Hf < 0.00435	<u>M</u> Ni < 0.65245	<u>M</u> Se < 0.01740	<u>M</u> W < 0.02175
<u>O</u> Ca 0.00164	<u>M</u> Ho < 0.00109	<u>M</u> Nb < 0.00109	<u>O</u> Si 0.00649	<u>M</u> U < 0.00435
<u>M</u> Ce < 0.01087	<u>M</u> In < 0.02175	<u>n</u> Os	<u>M</u> Ag < 0.00435	<u>M</u> V < 0.00435
<u>M</u> Cs < 0.00065	<u>M</u> Ir < 0.01087	<u>M</u> Pd < 0.01087	<u>O</u> Na 0.00368	<u>M</u> Yb < 0.00218
<u>O</u> Cr < 0.00900	<u>O</u> Fe 0.00268	<u>n</u> P	<u>M</u> Sr < 0.00109	<u>M</u> Y < 0.08699
<u>M</u> Co < 0.00652	<u>M</u> La < 0.00109	<u>M</u> Pt < 0.00435	<u>!</u> S	<u>M</u> Zn < 0.04350
<u>M</u> Cu < 0.01305	<u>M</u> Pb < 0.00652	<u>O</u> K < 0.10000	<u>M</u> Ta < 0.01522	<u>M</u> Zr < 0.01087

M - Checked by ICP-MS O - Checked by ICP-OES ! - Spectral Interference n - Not Checked For s - Solution Standard Element

6.0 INTENDED USE

For the calibration of analytical instruments including but not limited to the following:
 ICP-MS, ICP-OES, FAAS, GFAA, XRF, and DCP
 For the validation of analytical methods
 For the preparation of "working reference samples"
 For interference studies and the determination of correction coefficients
 For detection limit and linearity studies
 For additional intended uses, contact IV Technical Staff

7.0 INSTRUCTIONS FOR THE CORRECT USE OF THIS REFERENCE MATERIAL

Storage & Handling - Keep tightly sealed when not in use. Store and use at 20 ± 4°C. Do not pipet from container. Do not return portions removed for pipetting to container.

Atomic Weight; Valence; Coordination Number; Chemical Form in Solution - 9.01218; +2; 4; Be(H₂O)₄.¹²

Chemical Compatibility - Soluble in HCl, HNO₃, H₂SO₄, and HF aqueous matrices. Stable with all metals and inorganic anions.

Stability - 2-100 ppb levels stable for months in 1% HNO₃ / LDPE container. 1-10,000 ppm solutions chemically stable for years in 5-10% HNO₃ / LDPE container.

Be Containing Samples (Preparation and Solution) - Metal (is best dissolved in diluted H₂SO₄). BeO (boiling nitric, hydrochloric, or sulfuric acids or KHSO₄ fusion). Ores (H₂SO₄/HF digestion or carbonate fusion in Pt⁶). Organic Matrices (sulfuric/peroxide digestion or nitric/sulfuric/perchloric acid decomposition, or dry ash and dissolution according to the BeO procedure above).

Atomic Spectroscopic Information (ICP-OES D.L.s are given as radial/axial view):

Technique/Line	Estimated D.L.	Order	Type	Interferences (underlined indicates severe at 1000x)
ICP-OES 313.042 nm	0.0003 / 0.00009 µg/mL	1	ion	V, Ce, U
ICP-OES 234.861 nm	0.0003 / 0.00016 µg/mL	1	atom	Fe, Ta, Mo
ICP-OES 313.107 nm	0.0007 / 0.0005 µg/mL	1	ion	Ce, Th, Tm
ICP-MS 9 amu	4 ppt	n/a	M'	

8.0 **HAZARDOUS INFORMATION** - Please refer to the enclosed Material Safety Data sheet for information regarding this CRM.

9.0 **HOMOGENEITY** - This solution was mixed according to procedure IV-MPM-004 and is guaranteed to be homogeneous.

10.0 QUALITY STANDARD DOCUMENTATION

10.1 ISO 9001:2000 Quality Management System Registration - QMI Certificate Number 010105

Recognized by:

Registrar Accreditation Board (ANSI-RAB)

Standards Council of Canada (SCC)

Dutch Council for Accreditation (RVA)

Entidad Mexicana de Acreditacion, a.c.(EMA)

Members of IQ Net International Certification Network:

Argentina (IRAM), Australia (QAS), Austria (ÖQS), Belgium (Avinter), Brazil (FCAV), Canada (QMI), Hong Kong (HKQAA), Columbia (ICONTEC), Czech Republic (CQS), Denmark (DS), Finland (SFS), France (AFAQ), Germany (DQS), Greece (ELOT), Hungary (MSZT), Ireland (NSAI), Israel (SII), Italy (CISQ), Japan (JQA), Korea (KSA-QA), Netherlands (KEMA), Norway (NCS), Poland (PCBC), Portugal (APCER), Singapore (PSB), Slovenia (SIQ), Spain (AENOR), Switzerland (SQS)

10.2 ISO/IEC 17025 - 1999 "General Requirements for the Competence of Testing and Calibration"

- Chemical Testing - Accredited A2LA Certificate Number 883.01

10.3 ISO/IEC Guide 34 - 2000 "General Requirements for the Competence of Reference Material Producers"

- Reference Materials Production - Accredited A2LA Certificate Number 883.02

A2LA Mutual Recognition Agreement Partners:

Australia (NATA), Austria (BmWA), Belgium (BELTEST) (BKO-OBE), Canada (SCC), Chinese Taipei (CNLA), Czech Republic (NAO), Denmark (DANAK), Finland (FINAS), France (COFRAC), Germany (DAR), Hong Kong (HKAS), Ireland (NAB), Italy (SIT) (SINAL), Japan (JAB) (JNLA), Republic of Korea (KOLAS), The Netherlands (RvA), New Zealand (IANZ), Norway (NA), Portugal (IPQ), Singapore (SAC-SINGLAS), Spain (ENAC), Sweden (SWEDAC), Switzerland (SAS), United Kingdom (UKAS) and United States (NVLAP) (ICBO ES)

10.4 10CFR50 Appendix B - Nuclear Regulatory Commission

- Domestic Licensing of Production and Utilization Facilities

10.5 10CFR21 - Nuclear Regulatory Commission - Reporting Defects and Non-Compliance

10.6 MIL-STD-45662A (Obsolete/Observed)

11.0 DATE OF CERTIFICATION AND PERIOD OF VALIDITY



11.1 **IV Shelf Life** - The period of time during which the concentration of the analyte(s) in a properly packaged, unopened, and unused standard stored under environmentally controlled and monitored conditions will remain within the specified uncertainty range. Shelf life is limited primarily by transpiration (loss of water from the solution) and infrequently, by chemical instability. Transpiration studies (P-SP01020) of chemically-stable solutions performed at Inorganic Ventures / IV Labs indicate a CRM shelf-life of four years for solutions packaged in 500-mL low density polyethylene bottles. When stored under special conditions that minimize transpiration and instability, the shelf life can be extended past this limit.

11.2 **Expiration Date** - The date after which a CRM should not be used. Routine laboratory use of a CRM increases transpiration losses and the chance of contamination which affect the integrity of the CRM and limit its useful life. Inorganic Ventures / IV Labs concurs with state and federal regulatory agencies' recommendations that solution standards be assigned a one-year expiration date.

Certification Date: April 14, 2003

Expiration Date:

EXPIRES
12/2004

INORGANIC LABS/RADCHEM LABS Pg. 2 of 2
DATE RECEIVED: 05/09/03
DATE EXPIRED: 06/01/2004
DATE OPENED: 05/12/03
INORG: 4062 PO: F52181

12.0 NAMES AND SIGNATURES OF CERTIFYING OFFICERS

Certificate Prepared By: JoAnn Struthers, QA Administrative Assistant

Certificate Approved By: Katalin Le, QC Supervisor

Certifying Officer: Paul Gaines, Chemist, Senior Technical Director

JoAnn Struthers

Katalin Le

Paul Gaines

010201



inorganic ventures / iv labs

195 lehigh avenue, suite 4, lakewood, nj 08701 usa
phone: 800-669-6799 • 732-901-1900 • fax: 732-901-1903
e-mail: ivsales@ivstandards.com • website: www.ivstandards.com

certificate of analysis

CUSTOM-GRADE SOLUTION 1000 µg/mL Chromium⁺³ in 1.4% HNO₃ (abs)

Catalog Number: CGCR(3)1-1, CGCR(3)1-2 and CGCR(3)1-5

Lot Number: **W-QCR02033**

Starting Material:
Starting Material Purity:
Starting Material Lot No:

Chromium Metal
99.995%
F16122

INORGANIC LABS/RADCHEM LABS
DATE RECEIVED: 11/5/03
DATE EXPIRED: 12/1/04
DATE OPENED: 11/5/03
INORG: 4318 PO: F52258

CERTIFIED CONCENTRATION: 995 ± 3 µg/mL

The Certified Value is the instrument analysis value. The following equations are used in the calculation of the certified value and the uncertainty:

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

$$\text{Uncertainty } (\pm) = \frac{2[(\sum s_i^2)^{1/2}]}{(n)^{1/2}}$$

($\bar{x}$) = mean

x_i = individual results

n = number of measurements

$\sum s_i$ = The summation of all significant estimated errors.

Instrument Analysis: 995 ± 3 µg/mL (Avg of 3 runs)

Method: Inductively Coupled Plasma Spectroscopy (ICP) vs NIST SRM 3112a.

Calculated Value: 1002 µg/mL

Method: Calculated, based on starting material.

TRACE METALLIC IMPURITIES DETERMINED BY ICP-MS AND ICP-OES IN µg/mL:

Custom-Grade solutions tested for trace metallic impurities by ICP-MS were analyzed in an ULPA-Filtered Clean Room. An ULPA-Filter is 99.9985% efficient for the removal of particles down to 0.3 µm.

<u>Q</u>	Al	0.0028	<u>M</u>	Dy	<0.00060	<u>M</u>	Li	<0.0010	<u>M</u>	Pr	<0.000030	<u>M</u>	Te	<0.0030
<u>M</u>	Sb	<0.000050	<u>M</u>	Er	<0.00050	<u>M</u>	Lu	<0.000040	<u>M</u>	Re	<0.00010	<u>M</u>	Tb	<0.000030
<u>M</u>	As	<0.0010	<u>M</u>	Eu	<0.00030	<u>Q</u>	Mg	<0.010	<u>M</u>	Rh	<0.00010	<u>M</u>	Ti	<0.00010
<u>M</u>	Ba	<0.0010	<u>M</u>	Gd	<0.00010	<u>Q</u>	Mn	<0.050	<u>M</u>	Rb	0.0066	<u>M</u>	Th	<0.00010
<u>M</u>	Be	<0.000050	<u>M</u>	Ga	0.00070	<u>Q</u>	Hg	<0.10	<u>M</u>	Ru	0.017	<u>M</u>	Tm	<0.000040
<u>M</u>	Bi	<0.000040	<u>M</u>	Ge	<0.00060	<u>M</u>	Mo	<0.00020	<u>M</u>	Sm	<0.00010	<u>M</u>	Sn	<0.00050
<u>M</u>	B	<0.0070	<u>M</u>	Au	<0.00030	<u>M</u>	Nd	<0.00020	<u>M</u>	Sc	<0.0010	<u>M</u>	Ti	<0.0050
<u>M</u>	Cd	<0.00030	<u>M</u>	Hf	<0.00020	<u>Q</u>	Ni	<0.10	<u>i</u>	Se		<u>M</u>	W	<0.0010
<u>Q</u>	Ca	0.0011	<u>M</u>	Ho	<0.000050	<u>M</u>	Nb	<0.000050	<u>Q</u>	Si	<0.10	<u>M</u>	U	<0.00020
<u>M</u>	Ce	<0.00050	<u>M</u>	In	<0.10	<u>n</u>	Os		<u>M</u>	Ag	0.00070	<u>i</u>	V	
<u>M</u>	Cs	<0.000030	<u>M</u>	Ir	<0.00050	<u>M</u>	Pd	<0.00050	<u>Q</u>	Na	0.016	<u>M</u>	Yb	<0.00010
<u>s</u>	Cr		<u>Q</u>	Fe	<0.10	<u>i</u>	P		<u>M</u>	Sr	<0.000050	<u>M</u>	Y	<0.0040
<u>Q</u>	Co	<0.10	<u>M</u>	La	<0.000050	<u>M</u>	Pt	<0.00020	<u>n</u>	S		<u>Q</u>	Zn	<0.10
<u>M</u>	Cu	<0.00060	<u>M</u>	Pb	0.00039	<u>n</u>	K		<u>M</u>	Ta	<0.00070	<u>M</u>	Zr	<0.00050

M - checked by ICP-MS

Q - checked by ICP-OES

i - spectral interference

n - not checked for

s - solution standard element

ANALYZED DENSITY OF SOLUTION (measured at 22°C): 1.010 g/mL

(over)

QA:KSL Rev.090403.JTS

Paul R. Gaines

Quality Assurance Manager

Expires:

EXPIRES
12/2004



Certificate of Analysis

010203

CUSTOM-GRADE SOLUTION 1000 µg/mL Copper in 2% HNO₃ (abs)

Catalog Number: CGCU1-1, CGCU1-2 and CGCU1-5

Lot Number: **T-CU02060**

Starting Material:
Starting Material Purity:
Starting Material Lot No:

Copper Metal
99.999%
K09C13

INORGANIC LABS/RADCHEM LABS

DATE RECEIVED: 02/13/03

DATE EXPIRED: 03/01/2004

DATE OPENED: 02/13/03

INORG: 3955 PO: F52060

CERTIFIED CONCENTRATION: 1005 ± 2 µg/mL

The Certified Value is based upon the most precise method used to analyze this CRM. The following equations are used in the calculation of the certified value and the uncertainty:

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

$$\text{Uncertainty } (\pm) = \frac{2[(\sum s_i)^2]^{1/2}}{(n)^{1/2}}$$

($\bar{x}$) = mean

x_i = individual results

n = number of measurements

$\sum S_i$ = The summation of all significant estimated errors.

Classical Wet Assay: 1005 ± 2 µg/mL

Method: EDTA Titration vs NIST SRM 928 Lead Nitrate.

Instrument Analysis: 1003 ± 5 µg/mL

Method: Inductively Coupled Plasma Spectroscopy (ICP) vs NIST SRM 3114.

The independent samples t-test was used to determine if there is agreement between the above assay methods at the 95% confidence interval. Both methods were compared and showed agreement within the stated uncertainties. This agreement is a confirmation of the accuracy of this CRM.

TRACE METALLIC IMPURITIES DETERMINED BY ICP-MS AND ICP-OES IN µg/mL:

Custom-Grade solutions tested for trace metallic impurities by ICP-MS were analyzed in an ULPA-Filtered Clean Room.

An ULPA-Filter is 99.9985% efficient for the removal of particles down to 0.3 µm.

<u>Q</u> Al <0.010	<u>M</u> Dy <0.00060	<u>Q</u> Li <0.0050	<u>M</u> Pr <0.000030	<u>Q</u> Te <0.031
<u>M</u> Sb <0.000050	<u>M</u> Er <0.00050	<u>M</u> Lu <0.00004	<u>M</u> Re <0.00010	<u>M</u> Tb <0.000030
<u>Q</u> As <0.044	<u>M</u> Eu <0.00030	<u>Q</u> Mg 0.00080	<u>M</u> Rh 0.0010	<u>M</u> Tl <0.00010
<u>M</u> Ba <0.0010	<u>M</u> Gd <0.00010	<u>Q</u> Mn <0.00020	<u>M</u> Rb <0.00010	<u>M</u> Th <0.00010
<u>Q</u> Be <0.00020	<u>Q</u> Ga <0.0070	<u>Q</u> Hg <0.0070	<u>M</u> Ru <0.00020	<u>M</u> Tm <0.000040
<u>M</u> Bi 0.000090	<u>Q</u> Ge <0.010	<u>M</u> Mo <0.00020	<u>M</u> Sm <0.00010	<u>M</u> Sn 0.0013
<u>Q</u> B <0.0060	<u>Q</u> Au <0.010	<u>M</u> Nd <0.00020	<u>M</u> Sc <0.0010	<u>Q</u> Ti <0.00030
<u>Q</u> Cd <0.0018	<u>M</u> Hf <0.00020	<u>Q</u> Ni <0.0060	<u>Q</u> Se <0.020	<u>M</u> W <0.0010
<u>Q</u> Ca 0.00058	<u>M</u> Ho <0.000050	<u>M</u> Nb <0.00005	<u>Q</u> Si <0.0033	<u>M</u> U <0.00020
<u>M</u> Ce <0.00050	<u>Q</u> In <0.030	<u>n</u> Os	<u>Q</u> Ag <0.0040	<u>M</u> V <0.00020
<u>M</u> Cs <0.000030	<u>M</u> Ir <0.00050	<u>M</u> Pd 0.0024	<u>Q</u> Na <0.10	<u>M</u> Yb <0.00010
<u>Q</u> Cr <0.0020	<u>Q</u> Fe <0.0020	<u>Q</u> P <0.030	<u>M</u> Sr <0.000050	<u>M</u> Y <0.0040
<u>M</u> Co <0.00030	<u>M</u> La <0.000050	<u>M</u> Pt <0.00020	<u>Q</u> S <0.050	<u>Q</u> Zn i
<u>s</u> Cu	<u>M</u> Pb 0.0029	<u>Q</u> K <0.0060	<u>M</u> Ta <0.00070	<u>M</u> Zr <0.00050

M - checked by ICP-MS

O - checked by ICP-OES

i - spectral interference

n - not checked for

s - solution standard element

ANALYZED DENSITY OF SOLUTION (measured at 22°C): 1.015 g/mL

QA:KSL Rev.070302DN

(over)

Inorganic Ventures, Inc.

195 Lehigh Avenue • Suite 4 • Lakewood, NJ 08701

Orders: 800-669-6799 • FAX (732) 901-1903

Technical Support: 800-569-6799

Quality Assurance Manager

EXPIRES

12/2004

010204

certificate of analysis

1.0 Inorganic Ventures / IV Labs is an ISO Guide 34-2000 Certified Reference Material (CRM) Manufacturer: Certificate #883-02. The certificate is designed and the certified value(s) and uncertainty(ies) are determined in accordance with ISO Guide 31-2000 (Reference Materials - Contents of certificates and label(s)), ISO Guide 34-2000 "Quality System Guidelines for the Production of Reference Materials," and ISO Guide 35-1989 "Certification of Reference Materials - General and Statistical Principles."

2.0 DESCRIPTION OF CRM Custom-Grade 1000 µg/mL Nickel in 1.4% (abs) HNO₃

Catalog Number: CGNI1-1, CGNI1-2, and CGNI1-5

Lot Number: T-NI02028

Starting Material: Ni pieces

Starting Material Purity (%): 99.9994

Starting Material Lot No L06L02

Matrix: 1.4% (abs) HNO₃

INORGANIC LABS/RADCHEM LABS Pg. 1 of 2
DATE RECEIVED: 02/13/03
DATE EXPIRED: 03/01/2004
DATE OPENED: 02/13/03
INORG: 3957 PO: ES2060

3.0 CERTIFIED VALUES AND UNCERTAINTIES

Certified Concentration: 998 ± 2 µg/mL

Certified Density: 1.011 g/mL (measured at 22° C)

The Certified Value is based upon the most precise method used to analyze this CRM. The following equations are used in the calculation of the certified value and the uncertainty:

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

($\bar{x}$) = mean

x_i = individual results

n = number of measurements

$$\text{Uncertainty } (\pm) = \frac{2[(\sum s_i^2)^{1/2}]}{(n)^{1/2}}$$

$\sum S$ = The summation of all significant estimated errors.

(Most common are the errors from instrumental measurement, weighing, dilution to volume, and the fixed error reported on the NIST SRM certificate of analysis.)

The independent samples t-test was used to determine if there is agreement between the above assay methods at the 95% confidence interval. Both methods were compared and showed agreement within the stated uncertainties. This agreement is a confirmation of the accuracy of this CRM.

4.0 TRACEABILITY TO NIST AND VALUES OBTAINED BY INDEPENDENT METHODS

• "Property of the result of a measurement or the value of a standard whereby it can be related to stated references, usually national or international standards, through an unbroken chain of comparisons all having stated uncertainties." (ISO VIM, 2nd ed., 1993, definition 6.10)

• This IV product is Traceable to NIST via direct comparison to NIST SRMs. The uncertainties for each certified value are reported, taking into account the SRM uncertainty error and the measurement, weighing and volume dilution errors.

4.1 Assay Method #1 1002 ± 3 µg/mL

ICP Assay NIST SRM 3136 Lot Number: 000612

Assay Method #2 998 ± 2 µg/mL

EDTA NIST SRM 928 Lot Number: 880710

4.2 BALANCE CALIBRATION - All balances are checked daily using in-house procedure number 6-IMM-001. The weights used for testing are annually compared to Gerhart Scale Corporation's master weights and are traceable to the National Institute of Standards and Technology (NIST). The NIST Traceability numbers are 692476 - Class 1 and 692476A - Class 2. The NIST test number is 822/260017-98. All analytical balances are calibrated every 4 months by Gerhart Scale Corp. of South Amboy. The balances are calibrated with a class 1 and/or class 2 analytical weight set. These weights are tested annually by a NIST / NVLAP accredited calibration lab. The NIST test number is 822/260017-98.

010205

4.3 THERMOMETER CALIBRATION - The thermometers used in the determination of the final densities are calibrated vs standard thermometer No. 903-2680 which was certified in accordance with the procedures outlined by ASTM E77-87 and NIST Monograph 150 using NIST Test Nos. and Std Nos.: 769543, 217368/769543, 217368/P14452, 176240/P14452, 176240. The in-house procedure No. is 2-QC-001. Thermometers which are not calibrated vs standard thermometer No. 903-2680 are traceable to NIST Identification Nos. 92564, 119016, 471047 and NIST test report Nos. 811/258522, 811/2557078, and 236090.

4.4 GLASSWARE CALIBRATION - In-house procedure 3-QC-002 is used to calibrate all Class A Glassware used in the manufacture and quality control of Custom Grade Standards.

5.0 TRACE METALLIC IMPURITIES (TMI) DETERMINED BY ICP/MS AND ICP-OES IN µg/mL

Custom-Grade solutions are tested for trace metallic impurities by Axial ICP-OES and ICP-MS. The result from the most sensitive method for each element, is reported below. Solutions tested by ICP-MS were analyzed in an ULPA-Filtered Clean Room. An ULPA-Filter is 99.9985% efficient for the removal of particles down to 0.3 µm.

<u>Q</u> Al < 0.00938	<u>M</u> Dy < 0.06577	<u>Q</u> Li < 0.00006	<u>M</u> Pr < 0.00329	<u>M</u> Te < 0.32886
<u>M</u> Sb < 0.00548	<u>M</u> Er < 0.05481	<u>M</u> Lu < 0.00439	<u>M</u> Re < 0.01096	<u>M</u> Tb < 0.00329
<u>Q</u> As < 0.01689	<u>M</u> Eu < 0.03289	<u>Q</u> Mg 0.00002	<u>M</u> Rh < 0.01096	<u>M</u> Tl < 0.01096
<u>M</u> Ba < 0.10962	<u>M</u> Gd < 0.01096	<u>M</u> Mn < 0.04385	<u>M</u> Rb < 0.01096	<u>M</u> Th < 0.01096
<u>Q</u> Be < 0.00626	<u>M</u> Ga < 0.01096	<u>Q</u> Hg < 0.03441	<u>M</u> Ru < 0.02192	<u>M</u> Tm < 0.00439
<u>M</u> Bi < 0.00439	<u>M</u> Ge < 0.06577	<u>M</u> Mo < 0.02192	<u>M</u> Sm < 0.01096	<u>M</u> Sn < 0.05481
<u>Q</u> B < 0.03097	<u>M</u> Au < 0.03289	<u>M</u> Nd < 0.02192	<u>M</u> Sc < 0.10962	<u>M</u> Ti < 0.54811
<u>M</u> Cd < 0.03289	<u>M</u> Hf < 0.02192	<u>S</u> Ni	<u>Q</u> Se < 0.01877	<u>M</u> W < 0.10962
<u>Q</u> Ca < 0.01157	<u>M</u> Ho < 0.00548	<u>M</u> Nb < 0.00548	<u>Q</u> Si 0.00188	<u>M</u> U < 0.02192
<u>M</u> Ce < 0.05481	<u>M</u> In < 0.10962	<u>n</u> Os	<u>M</u> Ag < 0.02192	<u>M</u> V < 0.02192
<u>M</u> Cs < 0.00329	<u>M</u> Ir < 0.05481	<u>M</u> Pd < 0.05481	<u>Q</u> Na 0.00102	<u>M</u> Yb < 0.01096
<u>M</u> Cr < 0.05481	<u>Q</u> Fe 0.00156	<u>Q</u> P < 0.31280	<u>M</u> Sr < 0.00548	<u>M</u> Y < 0.43849
<u>Q</u> Co 0.00182	<u>M</u> La < 0.00548	<u>M</u> Pt < 0.02192	<u>Q</u> S < 0.07820	<u>M</u> Zn 0.00189
<u>M</u> Cu < 0.06577	<u>M</u> Pb < 0.03289	<u>Q</u> K 0.00043	<u>M</u> Ta < 0.07674	<u>M</u> Zr < 0.05481

M - Checked by ICP-MS O - Checked by ICP-OES i - Spectral Interference n - Not Checked For s - Solution Standard Element

6.0 INTENDED USE

For the calibration of analytical instruments including but not limited to the following:

ICP-MS, ICP-OES, FAAS, GFAA, XRF, and DCP

For the validation of analytical methods

For the preparation of "working reference samples"

For interference studies and the determination of correction coefficients

For detection limit and linearity studies

For additional intended uses, contact IV Technical Staff

7.0 INSTRUCTIONS FOR THE CORRECT USE OF THIS REFERENCE MATERIAL

Storage & Handling - Keep tightly sealed when not in use. Store and use at 20 ± 4°C. Do not pipet from container. Do not return portions removed for pipetting to container.

Atomic Weight; Valence; Coordination Number; Chemical Form in Solution - 58.69; +2; 6; Ni(H₂O)₆²⁺

Chemical Compatibility - Stable in HCl, HNO₃, H₂SO₄, HF, H₃PO₄. Avoid basic media. Stable with most metals and inorganic anions in acidic media.

Stability - 2-100 ppb levels stable for months in 1% HNO₃ / LDPE container. 1-10,000 ppm solutions chemically stable for years in 1-5% HNO₃ / LDPE container.

Ni Containing Samples (Preparation and Solution) - Metal (Soluble in HNO₃); Oxides (Soluble in HCl); Ores (Dissolve in HCl / HNO₃).

Atomic Spectroscopic Information (ICP-OES D.L.s are given as radial/axial view):

Technique/Line	Estimated D.L.	Order	Type	Interferences (underlined indicates severe at ≈ concs.)
ICP-OES 221.647 nm	0.01 / 0.0009 µg/mL	1	ion	Si
ICP-OES 232.003 nm	0.02 / 0.006 µg/mL	1 atom		Cr, Re, Os, Nb, Ag, Pt, Fe
ICP-OES 231.604 nm	0.02 / 0.002 µg/mL	1 ion		Sb, Ta, Co
ICP-MS 60 amu	100 ppt	n/a	M ⁺	⁴³ Ca ¹⁶ O ¹ H, ⁴⁴ Ca ¹⁶ O, ²³ Na ³⁷ Cl

8.0 **HAZARDOUS INFORMATION** - Please refer to the enclosed Material Safety Data sheet for information regarding this CRM.

9.0 **HOMOGENEITY** - This solution was mixed according to procedure IV-MPM-004 and is guaranteed to be homogeneous.

10.0 QUALITY STANDARD DOCUMENTATION

010206

10.1 ISO 9001:2000 Quality Management System Registration - QMI Certificate Number 010105

Recognized by:

Registrar Accreditation Board (ANSI-RAB)

Standards Council of Canada (SCC)

Dutch Council for Accreditation (RVA)

Entidad Mexicana de Acreditacion, a.c.(EMA)

Members of IQ Net International Certification Network:

Argentina (IRAM), Australia (QAS), Austria (OQS), Belgium (Avinter), Brazil (FCAV), Canada (QMI), Hong Kong (HKQAA), Columbia (ICONTEC), Czech Republic (CQS), Denmark (DS), Finland (SFS), France (AFAQ), Germany (DQS), Greece (ELOT), Hungary (MSZT), Ireland (NSAI), Israel (SII), Italy (CISQ), Japan (JQA), Korea (KSA-QA), Netherlands (KEMA), Norway (NCS), Poland(PCBC), Portugal (APCER), Singapore (PSB), Slovenia (SIQ), Spain (AENOR), Switzerland (SQS)



10.2 ISO/IEC 17025 - 1999 "General Requirements for the Competence of Testing and Calibration"

- Chemical Testing - Accredited A2LA Certificate Number 883.01



10.3 ISO/IEC Guide 34 - 2000 "General Requirements for the Competence of Reference Material Producers"

- Reference Materials Production - Accredited A2LA Certificate Number 883.02

A2LA Mutual Recognition Agreement Partners:

Australia (NATA), Austria (BmWA), Belgium (BELTEST) (BKO-OBE), Canada (SCC), Chinese Taipei (CNLA), Czech Republic (NAO), Denmark (DANAK), Finland (FINAS), France (COFRAC), Germany (DAR), Hong Kong (HKAS), Ireland (NAB), Italy (SIT) (SINAL), Japan (JAB) (JNLA), Republic of Korea (KOLAS), The Netherlands (RvA), New Zealand (IANZ), Norway (NA), Portugal (IPQ), Singapore (SAC-SINGLAS), Spain (ENAC), Sweden (SWEDAC), Switzerland (SAS), United Kingdom (UKAS) and United States (NVLAP) (ICBO ES)


10.4 10CFR50 Appendix B - Nuclear Regulatory Commission

- Domestic Licensing of Production and Utilization Facilities

10.5 10CFR21 - Nuclear Regulatory Commission - Reporting Defects and Non-Compliance

10.6 MIL-STD-45662A (Obsolete/Observed)

11.0 DATE OF CERTIFICATION AND PERIOD OF VALIDITY



11.1 **IV Shelf Life** - The period of time during which the concentration of the analyte(s) in a properly packaged, unopened, and unused standard stored under environmentally controlled and monitored conditions will remain within the specified uncertainty range. Shelf life is limited primarily by transpiration (loss of water from the solution) and infrequently, by chemical instability. Transpiration studies (P-SP01020) of chemically-stable solutions performed at Inorganic Ventures / IV Labs indicate a CRM shelf-life of four years for solutions packaged in 500-mL low density polyethylene bottles. When stored under special conditions that minimize transpiration and instability, the shelf life can be extended past this limit.

11.2 **Expiration Date** - The date after which a CRM should not be used. Routine laboratory use of a CRM increases transpiration losses and the chance of contamination which affect the integrity of the CRM and limit its useful life. Inorganic Ventures / IV Labs concurs with state and federal regulatory agencies' recommendations that solution standards be assigned a one-year expiration date.

Certification Date: October 11, 2002

Expiration Date:

EXPIRES
11/2004

INORGANIC LABS/RADCHEM LABS Pg. 2 of 2

DATE RECEIVED: 02/13/03

DATE EXPIRED: 03/01/2004

DATE OPENED: 02/13/03

INORG: 3957 PO: F52060

12.0 NAMES AND SIGNATURES OF CERTIFYING OFFICERS

010207

Certificate Prepared By: Debbie Newman, QA Administrator

Debbie Newman
Katalin Le

Certificate Approved By: Katalin Le, QC Supervisor

Certifying Officer: Paul Gaines, Chemist, Senior Technical Director

Paul Gaines



Certificate of Analysis



CUSTOM-GRADE SOLUTION 1000 µg/mL Cadmium in 2% HNO₃ (abs)

Catalog Number: CGCD1-1, CGCD1-2 and CGCD1-5

010208

Lot Number: **T-CD01125**

Starting Material:
Starting Material Purity:
Starting Material Lot No:

Cadmium Metal
99.999%
H02K40A

INORGANIC LABS/RADCHEM LABS
DATE RECEIVED: 3/03/03
DATE EXPIRED: 3/01/2004
DATE OPENED: 3/04/03
INORG: 3972 PO: F52004

CERTIFIED CONCENTRATION: 1003 ± 3 µg/mL

The Certified Value is based upon the most precise method used to analyze this CRM. The following equations are used in the calculation of the certified value and the uncertainty:

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

$$\text{Uncertainty } (\pm) = \frac{2[(\sum s_i)^2]^{1/2}}{(n)^{1/2}}$$

($\bar{x}$) = mean

x_i = individual results

n = number of measurements

$\sum S_i$ = The summation of all significant estimated errors.

Classical Wet Assay: 1003 ± 3 µg/mL

Method: EDTA Titration vs NIST SRM 928 Lead Nitrate.

Instrument Analysis: 1003 ± 3 µg/mL

Method: Inductively Coupled Plasma Spectroscopy (ICP) vs NIST SRM 3108.

The independent samples t-test was used to determine if there is agreement between the above assay methods at the 95% confidence interval. Both methods were compared and showed agreement within the stated uncertainties. This agreement is a confirmation of the accuracy of this CRM.

TRACE METALLIC IMPURITIES DETERMINED BY ICP-MS AND ICP-OES IN µg/mL:

Custom-Grade solutions tested for trace metallic impurities by ICP-MS were analyzed in an ULPA-Filtered Clean Room.

An ULPA-Filter is 99.9985% efficient for the removal of particles down to 0.3 µm.

<u>Q</u> Al <0.00030	<u>M</u> Dy <0.00060	<u>Q</u> Li <0.0050	<u>M</u> Pr <0.000030	<u>Q</u> Te <0.031
<u>M</u> Sb <0.000050	<u>M</u> Er <0.00050	<u>M</u> Lu <0.000040	<u>M</u> Re <0.00010	<u>M</u> Tb <0.000030
<u>Q</u> As <0.12	<u>M</u> Eu <0.00030	<u>Q</u> Mg <0.00010	<u>M</u> Rh <0.00010	<u>M</u> Tl <0.00010
<u>M</u> Ba <0.0010	<u>M</u> Gd <0.00010	<u>Q</u> Mn <0.00020	<u>M</u> Rb <0.00010	<u>M</u> Th <0.00010
<u>Q</u> Be <0.000050	<u>M</u> Ga <0.00010	<u>Q</u> Hg <0.0070	<u>M</u> Ru <0.00020	<u>M</u> Tm <0.000040
<u>M</u> Bi <0.000040	<u>M</u> Ge <0.00060	<u>M</u> Mo <0.00020	<u>M</u> Sm <0.00010	<u>M</u> Sn <0.00050
<u>Q</u> B <0.0060	<u>Q</u> Au <0.010	<u>M</u> Nd <0.00020	<u>M</u> Sc <0.0010	<u>Q</u> Ti <0.00030
<u>s</u> Cd	<u>M</u> Hf <0.00020	<u>Q</u> Ni <0.0060	<u>Q</u> Se <0.020	<u>M</u> W <0.0010
<u>Q</u> Ca 0.0018	<u>M</u> Ho <0.000050	<u>M</u> Nb <0.000050	<u>Q</u> Si <0.010	<u>M</u> U <0.00020
<u>M</u> Ce <0.00050	<u>Q</u> In <0.030	<u>n</u> Os	<u>Q</u> Ag <0.27	<u>M</u> V <0.00020
<u>M</u> Cs <0.000030	<u>M</u> Ir <0.00050	<u>Q</u> Pd <0.053	<u>Q</u> Na <0.10	<u>M</u> Yb <0.00010
<u>Q</u> Cr <0.0040	<u>Q</u> Fe <0.0020	<u>Q</u> P <0.030	<u>M</u> Sr <0.000050	<u>M</u> Y <0.0040
<u>M</u> Co <0.00030	<u>M</u> La <0.000050	<u>M</u> Pt <0.00020	<u>Q</u> S <0.020	<u>Q</u> Zn <0.030
<u>M</u> Cu <0.00060	<u>M</u> Pb <0.00030	<u>Q</u> K <0.0060	<u>M</u> Ta <0.00070	<u>M</u> Zr <0.00050

M - checked by ICP-MS

Q - checked by ICP-OES

i - spectral interference

n - not checked for

s - solution standard element

ANALYZED DENSITY OF SOLUTION (measured at 22°C): 1.015 g/mL

(over)

QA:KL Rev.080102DN

Inorganic Ventures, Inc.

195 Lehigh Avenue • Suite 4 • Lakewood, NJ 08701

Orders: 800-669-6799 • FAX (732) 901-1903

Technical Support: 800-569-6799

Quality Assurance Manager

EXPIRES
1/2004



inorganic ventures / iv labs

195 lehigh avenue, suite 4, lakewood, nj 08701 usa
phone: 800-669-6799 • 732-901-1900 • fax: 732-901-1903
e-mail: ivsales@ivstandards.com • website: www.ivstandards.com

010209

certificate of analysis

1.0 Inorganic Ventures / IV Labs is an ISO Guide 34-2000 Certified Reference Material (CRM) Manufacturer: Certificate #883-02. The certificate is designed and the certified value(s) and uncertainty(ies) are determined in accordance with ISO Guide 31-2000 (Reference Materials - Contents of certificates and label(s)), ISO Guide 34-2000 "Quality System Guidelines for the Production of Reference Materials," and ISO Guide 35-1989 "Certification of Reference Materials - General and Statistical Principles."

2.0 DESCRIPTION OF CRM Custom-Grade 1000 µg/mL Cobalt in 2% (abs) HNO₃

Catalog Number: CGCO1-1, CGCO1-2, and CGCO1-5
Lot Number: T-CO01120
Starting Material: Co powder
Starting Material Purity (%): 99.9957
Starting Material Lot No 22897
Matrix: 2% (abs) HNO₃

INORGANIC LABS/RADCHEM LABS Pg. 1 of 2
DATE RECEIVED: 3/03/03
DATE EXPIRED: 3/01/2004
DATE OPENED: 3/04/03
INORG: 3973 PO: E52064

3.0 CERTIFIED VALUES AND UNCERTAINTIES

Certified Concentration: 1002 ± 3 µg/mL

Certified Density: 1.016 g/mL (measured at 22° C)

The Certified Value is based upon the most precise method used to analyze this CRM. The following equations are used in the calculation of the certified value and the uncertainty:

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

($\bar{x}$) = mean

x_i = individual results

n = number of measurements

$$\text{Uncertainty } (\pm) = \frac{2[(\sum s_i^2)^{1/2}]}{(n)^{1/2}}$$

$\sum S$ = The summation of all significant estimated errors.

(Most common are the errors from instrumental measurement, weighing, dilution to volume, and the fixed error reported on the NIST SRM certificate of analysis.)

The independent samples t-test was used to determine if there is agreement between the above assay methods at the 95% confidence interval. Both methods were compared and showed agreement within the stated uncertainties. This agreement is a confirmation of the accuracy of this CRM.

4.0 TRACEABILITY TO NIST AND VALUES OBTAINED BY INDEPENDENT METHODS

• "Property of the result of a measurement or the value of a standard whereby it can be related to stated references, usually national or international standards, through an unbroken chain of comparisons all having stated uncertainties." (ISO VIM, 2nd ed., 1993, definition 6.10)

• This IV product is Traceable to NIST via direct comparison to NIST SRMs. The uncertainties for each certified value are reported, taking into account the SRM uncertainty error and the measurement, weighing and volume dilution errors.

4.1 Assay Method #1 998 ± 4 µg/mL

ICP Assay NIST SRM 3181 Lot Number: 000630

Assay Method #2 1002 ± 3 µg/mL

EDTA NIST SRM 928 Lot Number: 880710

4.2 BALANCE CALIBRATION - All balances are checked daily using in-house procedure number 6-IMM-001. The weights used for testing are annually compared to Gerhart Scale Corporation's master weights and are traceable to the National Institute of Standards and Technology (NIST). The NIST Traceability numbers are 692476 - Class 1 and 692476A - Class 2. The NIST test number is 822/260017-98. All analytical balances are calibrated every 4 months by Gerhart Scale Corp. of South Amboy. The balances are calibrated with a class 1 and/or class 2 analytical weight set. These weights are tested annually by a NIST / NVLAP accredited calibration lab. The NIST test number is 822/260017-98.

010210

4.3 THERMOMETER CALIBRATION - The thermometers used in the determination of the final densities are calibrated vs standard thermometer No. 903-2680 which was certified in accordance with the procedures outlined by ASTM E77-87 and NIST Monograph 150 using NIST Test Nos. and Std Nos.: 769543, 217368/769543, 217368/P14452, 176240/P14452, 176240. The in-house procedure No. is 2-QC-001. Thermometers which are not calibrated vs standard thermometer No. 903-2680 are traceable to NIST Identification Nos. 92564, 119016, 471047 and NIST test report Nos. 811/258522, 811/2557078, and 236090.

4.4 GLASSWARE CALIBRATION - In-house procedure 3-QC-002 is used to calibrate all Class A Glassware used in the manufacture and quality control of Custom Grade Standards.

5.0 TRACE METALLIC IMPURITIES (TMI) DETERMINED BY ICP/MS AND ICP-OES IN µg/mL

Custom-Grade solutions are tested for trace metallic impurities by Axial ICP-OES and ICP-MS. The result from the most sensitive method for each element, is reported below. Solutions tested by ICP-MS were analyzed in an ULPA-Filtered Clean Room. An ULPA-Filter is 99.9985% efficient for the removal of particles down to 0.3 µm.

<u>Q</u> Al 0.00025	<u>M</u> Dy < 0.02419	<u>Q</u> Li 0.00001	<u>M</u> Pr < 0.00121	<u>M</u> Te < 0.12097
<u>M</u> Sb < 0.00202	<u>M</u> Er < 0.02016	<u>M</u> Lu < 0.00161	<u>M</u> Re < 0.00403	<u>M</u> Tb < 0.00121
<u>Q</u> As < 0.10000	<u>M</u> Eu < 0.01210	<u>Q</u> Mg 0.00045	<u>M</u> Rh < 0.00403	<u>M</u> Tl < 0.00403
<u>M</u> Ba < 0.04032	<u>M</u> Gd < 0.00403	<u>Q</u> Mn 0.00003	<u>M</u> Rb < 0.00403	<u>M</u> Th < 0.00403
<u>M</u> Be < 0.00202	<u>M</u> Ga < 0.00403	<u>Q</u> Hg < 0.05000	<u>M</u> Ru < 0.00807	<u>M</u> Tm < 0.00161
<u>M</u> Bi < 0.00161	<u>M</u> Ge < 0.02419	<u>M</u> Mo < 0.00807	<u>M</u> Sm < 0.00403	<u>M</u> Sn < 0.02016
<u>Q</u> B < 0.04000	<u>M</u> Au < 0.01210	<u>M</u> Nd < 0.00807	<u>M</u> Sc < 0.04032	<u>M</u> Ti < 0.20162
<u>M</u> Cd < 0.01210	<u>M</u> Hf < 0.00807	<u>Q</u> Ni < 0.02000	<u>M</u> Se < 0.03226	<u>M</u> W < 0.04032
<u>Q</u> Ca 0.00325	<u>M</u> Ho < 0.00202	<u>M</u> Nb < 0.00202	<u>Q</u> Si < 0.00400	<u>M</u> U < 0.00807
<u>M</u> Ce < 0.02016	<u>M</u> In < 0.04032	<u>n</u> Os	<u>M</u> Ag < 0.00807	<u>M</u> V < 0.00807
<u>M</u> Cs < 0.00121	<u>M</u> Ir < 0.02016	<u>M</u> Pd < 0.02016	<u>Q</u> Na 0.00138	<u>M</u> Yb < 0.00403
<u>M</u> Cr < 0.02016	<u>Q</u> Fe 0.00875	<u>n</u> P	<u>M</u> Sr < 0.00202	<u>M</u> Y < 0.16129
<u>s</u> Co	<u>M</u> La < 0.00202	<u>M</u> Pt < 0.00807	<u>n</u> S	<u>M</u> Zn < 0.08065
<u>M</u> Cu < 0.02419	<u>M</u> Pb < 0.01210	<u>Q</u> K 0.03000	<u>M</u> Ta < 0.02823	<u>M</u> Zr < 0.02016

M - Checked by ICP-MS O - Checked by ICP-OES i - Spectral Interference n - Not Checked For s - Solution Standard Element

6.0 INTENDED USE

For the calibration of analytical instruments including but not limited to the following:

ICP-MS, ICP-OES, FAAS, GFAA, XRF, and DCP

For the validation of analytical methods

For the preparation of "working reference samples"

For interference studies and the determination of correction coefficients

For detection limit and linearity studies

For additional intended uses, contact IV Technical Staff

7.0 INSTRUCTIONS FOR THE CORRECT USE OF THIS REFERENCE MATERIAL

Storage & Handling - Keep tightly sealed when not in use. Store and use at 20 ± 4°C. Do not pipet from container. Do not return portions removed for pipetting to container.

Atomic Weight; Valence; Coordination Number; Chemical Form in Solution - 58.9332; +2; 6; Co(H₂O)₆²⁺

Chemical Compatibility - Stable in HCl, HNO₃, H₂SO₄, HF, H₃PO₄. Avoid basic media. Stable with most metals and inorganic anions in acidic media.

Stability - 2-100 ppb levels stable for months in 1% HNO₃ / LDPE container. 1-10,000 ppm solutions chemically stable for years in 1-5% HNO₃ / LDPE container.

Co Containing Samples (Preparation and Solution) - Metal (soluble in HNO₃); Oxides (Soluble in HCl); Ores (Dissolve in HCl / HNO₃).

Atomic Spectroscopic Information (ICP-OES D.L.s are given as radial/axial view):

Technique/Line	Estimated D.L.	Order	Type	Interferences (underlined indicates severe at ≈ concs.)
ICP-OES 238.892 nm	0.01/.002 µg/mL	1	ion	Fe, W, Ta
ICP-OES 228.616 nm	0.01/.001 µg/mL	1	ion	
ICP-OES 237.862 nm	0.01/.002 µg/mL	1	ion	W, Re, Al, Ta
ICP-MS 59 amu	2 ppt	n/a	M ⁺	⁴² Ca ¹⁸ O ¹ H, ⁴⁰ Ar ¹⁸ O ¹ H, ³⁸ Ar ²³ Na, ⁴⁰ Ca ¹⁶ O, ²⁴ Mg ³⁵ Cl

8.0 HAZARDOUS INFORMATION - Please refer to the enclosed Material Safety Data sheet for information regarding this CRM.

9.0 HOMOGENEITY - This solution was mixed according to procedure IV-MPM-004 and is guaranteed to be homogeneous.

10.0 QUALITY STANDARD DOCUMENTATION

010211

10.1 ISO 9001:2000 Quality Management System Registration - QMI Certificate Number 010105

Recognized by:

Registrar Accreditation Board (ANSI-RAB)

Standards Council of Canada (SCC)

Dutch Council for Accreditation (RVA)

Entidad Mexicana de Acreditacion, a.c.(EMA)

Members of IQ Net International Certification Network:

Argentina (IRAM), Australia (QAS), Austria (OQS), Belgium (Avinter), Brazil (FCAV), Canada (QMI), Hong Kong (HKQAA), Columbia (ICONTEC), Czech Republic (CQS), Denmark (DS), Finland (SFS), France (AFAQ), Germany (DQS), Greece (ELOT), Hungary (MSZT), Ireland (NSAI), Israel (SII), Italy (CISQ), Japan (JQA), Korea (KSA-QA), Netherlands (KEMA), Norway (NCS), Poland(PCBC), Portugal (APCER), Singapore (PSB), Slovenia (SIQ), Spain (AENOR), Switzerland (SQS)



10.2 ISO/IEC 17025 - 1999 "General Requirements for the Competence of Testing and Calibration"

- Chemical Testing - Accredited A2LA Certificate Number 883.01



10.3 ISO/IEC Guide 34 - 2000 "General Requirements for the Competence of Reference Material Producers"

- Reference Materials Production - Accredited A2LA Certificate Number 883.02

A2LA Mutual Recognition Agreement Partners:

Australia (NATA), Austria (BmWA), Belgium (BELTEST) (BKO-OBE), Canada (SCC), Chinese Taipei (CNLA), Czech Republic (NAO), Denmark (DANAK), Finland (FINAS), France (COFRAC), Germany (DAR), Hong Kong (HKAS), Ireland (NAB), Italy (SIT) (SINAL), Japan (JAB) (JNLA), Republic of Korea (KOLAS), The Netherlands (RvA), New Zealand (IANZ), Norway (NA), Portugal (IPQ), Singapore (SAC-SINGLAS), Spain (ENAC), Sweden (SWEDAC), Switzerland (SAS), United Kingdom (UKAS) and United States (NVLAP) (ICBO ES)


10.4 10CFR50 Appendix B - Nuclear Regulatory Commission

- Domestic Licensing of Production and Utilization Facilities

10.5 10CFR21 - Nuclear Regulatory Commission - Reporting Defects and Non-Compliance

10.6 MIL-STD-45662A (Obsolete/Observed)

11.0 DATE OF CERTIFICATION AND PERIOD OF VALIDITY



11.1 IV Shelf Life - The period of time during which the concentration of the analyte(s) in a properly packaged, unopened, and unused standard stored under environmentally controlled and monitored conditions will remain within the specified uncertainty range. Shelf life is limited primarily by transpiration (loss of water from the solution) and infrequently, by chemical instability. Transpiration studies (P-SP01020) of chemically-stable solutions performed at Inorganic Ventures / IV Labs indicate a CRM shelf-life of four years for solutions packaged in 500-mL low density polyethylene bottles. When stored under special conditions that minimize transpiration and instability, the shelf life can be extended past this limit.

11.2 Expiration Date - The date after which a CRM should not be used. Routine laboratory use of a CRM increases transpiration losses and the chance of contamination which affect the integrity of the CRM and limit its useful life. Inorganic Ventures / IV Labs concurs with state and federal regulatory agencies' recommendations that solution standards be assigned a one-year expiration date.

Certification Date: September 05, 2002

Expiration Date:

EXPIRES

12/2004

INORGANIC LABS/RADCHEM LABS Pg. 2 of 2

DATE RECEIVED: 3/03/03

DATE EXPIRED: 3/01/2004

DATE OPENED: 3/04/03

INORG: 3973 PO: F52064

12.0 NAMES AND SIGNATURES OF CERTIFYING OFFICERS

Certificate Prepared By: Debbie Newman, QA Administrator

Certificate Approved By: Katalin Le, QC Supervisor

Certifying Officer: Paul Gaines, Chemist, Senior Technical Director

010212

Debbie Newman
Katalin Le

Paul Gaines



Certificate of Analysis

CUSTOM-GRADE SOLUTION

1000 µg/mL Manganese in 2% HNO₃ (abs)

Catalog Number: CGMN1-1, CGMN1-2, and CGMN1-5

010213

Lot Number: **T-MN02033**

Starting Material:
Starting Material Purity:
Starting Material Lot No:

Manganese Metal
99.998%
21563

INORGANIC LABS/RADCHEM LABS

DATE RECEIVED: 01/16/03

DATE EXPIRED: 02/01/2004

DATE OPENED: 01/16/03

INORG: 3921 PO: F52073

CERTIFIED CONCENTRATION: 1004 ± 3 µg/mL

The Certified Value is the wet assay value. The following equations are used in the calculation of the certified value and the uncertainty:

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

$$\text{Uncertainty } (\pm) = \frac{2((\sum s_i)^2)^{1/2}}{(n)^{1/2}}$$

($\bar{x}$) = mean

x_i = individual results

n = number of measurements

$\sum S_i$ = The summation of all significant estimated errors.

Classical Wet Assay: 1004 ± 3 µg/mL

Method: EDTA Titration vs NIST SRM 928 Lead Nitrate.

Instrument Analysis: 1001 ± 3 µg/mL

Method: Inductively Coupled Plasma Spectroscopy (ICP) vs NIST SRM 3132.

The independent samples t-test was used to determine if there is agreement between the above assay methods at the 95% confidence interval. Both methods were compared and showed agreement within the stated uncertainties. This agreement is a confirmation of the accuracy of this CRM.

TRACE METALLIC IMPURITIES DETERMINED BY ICP-MS AND ICP-OES IN µg/mL:

Custom-Grade solutions tested for trace metallic impurities by ICP-MS were analyzed in an ULPA-Filtered Clean Room. An ULPA-Filter is 99.9985% efficient for the removal of particles down to 0.3 µm.

<u>M</u> Al <0.00090	<u>M</u> Dy <0.00060	<u>M</u> Li <0.0010	<u>M</u> Pr <0.000030	<u>M</u> Te <0.0030
<u>M</u> Sb <0.000050	<u>M</u> Er <0.00050	<u>M</u> Lu <0.000040	<u>M</u> Re <0.00010	<u>M</u> Tb <0.000030
<u>M</u> As <0.0010	<u>M</u> Eu <0.00030	<u>M</u> Mg 0.0040	<u>M</u> Rh <0.00010	<u>M</u> Tl <0.00010
<u>M</u> Ba <0.0010	<u>M</u> Gd <0.00010	<u>s</u> Mn	<u>M</u> Rb <0.00010	<u>M</u> Th <0.00010
<u>M</u> Be <0.000050	<u>i</u> Ga	<u>O</u> Hg <0.030	<u>M</u> Ru <0.00020	<u>M</u> Tm <0.000040
<u>M</u> Bi <0.000040	<u>M</u> Ge <0.00060	<u>M</u> Mo 0.00024	<u>M</u> Sm <0.00010	<u>M</u> Sn <0.00050
<u>M</u> B <0.0070	<u>M</u> Au <0.00030	<u>M</u> Nd <0.00020	<u>M</u> Sc <0.0010	<u>M</u> Ti <0.0050
<u>M</u> Cd <0.00030	<u>M</u> Hf <0.00020	<u>O</u> Ni <0.050	<u>O</u> Se <0.40	<u>M</u> W <0.0010
<u>O</u> Ca <0.010	<u>M</u> Ho <0.000050	<u>M</u> Nb <0.000050	<u>O</u> Si <0.020	<u>M</u> U <0.00020
<u>M</u> Ce <0.00050	<u>i</u> In	<u>n</u> Os	<u>M</u> Ag <0.00020	<u>i</u> V
<u>M</u> Cs <0.000030	<u>M</u> Ir <0.00050	<u>M</u> Pd <0.00050	<u>O</u> Na <0.090	<u>M</u> Yb <0.00010
<u>M</u> Cr <0.00050	<u>i</u> Fe	<u>i</u> P	<u>M</u> Sr <0.000050	<u>M</u> Y <0.0040
<u>O</u> Co <0.050	<u>M</u> La <0.000050	<u>M</u> Pt <0.00020	<u>n</u> S	<u>M</u> Zn <0.0020
<u>M</u> Cu <0.00060	<u>M</u> Pb <0.00030	<u>n</u> K	<u>M</u> Ta <0.00070	<u>M</u> Zr <0.00050

M - checked by ICP-MS

O - checked by ICP-OES

i - spectral interference

n - not checked for

s - solution standard element

ANALYZED DENSITY OF SOLUTION (measured at 22°C): 1.015 g/mL

(over)

QA:KL Rev.082202DN

Inorganic Ventures, Inc.

195 Lehigh Avenue • Suite 4 • Lakewood, NJ 08701

Orders: 800-669-6799 • FAX (732) 901-1903

Technical Support: 800-569-6799

Quality Assurance Manager

EXPIRES

1/2004

certificate of analysis

CUSTOM-GRADE SOLUTION

Catalog Number: CGV1-1, CGV1-2 and CGV1-5

1000 µg/mL Vanadium in 1.4% HNO₃ (abs)

Lot Number: **T-V02032**

Starting Material:
Starting Material Purity:
Starting Material Lot No:

Vanadium Pentoxide
99.999%
46

INORGANIC LABS/RADCHEM LABS

DATE RECEIVED: 11/5/13
DATE EXPIRED: 12/1/2014 *OR*
DATE OPENED: 11/5/13
INORG: 4321 PO: F52258

CERTIFIED CONCENTRATION: 990 ± 2 µg/mL

The Certified Value is based upon the most precise method used to analyze this CRM. The following equations are used in the calculation of the certified value and the uncertainty:

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

$$\text{Uncertainty } (\pm) = \frac{2[(\sum s_i)^2]^{1/2}}{(n)^{1/2}}$$

($\bar{x}$) = mean

x_i = individual results

n = number of measurements

$\sum S_i$ = The summation of all significant estimated errors.

Classical Wet Assay: 993 ± 4 µg/mL

Method: EDTA Titration vs NIST SRM 928 Lead Nitrate.

Instrument Analysis: 990 ± 2 µg/mL

Method: Inductively Coupled Plasma Spectroscopy (ICP) vs NIST SRM 3165.

The independent samples t-test was used to determine if there is agreement between the above assay methods at the 95% confidence interval. Both methods were compared and showed agreement within the stated uncertainties. This agreement is a confirmation of the accuracy of this CRM.

TRACE METALLIC IMPURITIES DETERMINED BY ICP-MS AND ICP-OES IN µg/mL:

Custom-Grade solutions tested for trace metallic impurities by ICP-MS were analyzed in an ULPA-Filtered Clean Room.
An ULPA-Filter is 99.9985% efficient for the removal of particles down to 0.3 µm.

<u>M</u>	Al	0.0095	<u>M</u>	Dy	<0.00060	<u>M</u>	Li	<0.0010	<u>M</u>	Pr	<0.000030	<u>M</u>	Te	<0.0030
<u>M</u>	Sb	0.042	<u>M</u>	Er	<0.00050	<u>M</u>	Lu	<0.000040	<u>M</u>	Re	<0.00010	<u>M</u>	Tb	<0.000030
<u>M</u>	As	<0.0010	<u>M</u>	Eu	<0.00030	<u>M</u>	Mg	0.0089	<u>M</u>	Rh	<0.00010	<u>M</u>	Tl	<0.00010
<u>M</u>	Ba	<0.0010	<u>M</u>	Gd	<0.00010	<u>i</u>	Mn		<u>M</u>	Rb	<0.00010	<u>M</u>	Th	<0.00010
<u>M</u>	Be	<0.000050	<u>M</u>	Ga	<0.00010	<u>i</u>	Hg		<u>M</u>	Ru	<0.00020	<u>M</u>	Tm	<0.000040
<u>M</u>	Bi	<0.000040	<u>M</u>	Ge	<0.00060	<u>M</u>	Mo	0.016	<u>M</u>	Sm	<0.00010	<u>M</u>	Sn	<0.00050
<u>M</u>	B	<0.0070	<u>M</u>	Au	<0.00030	<u>M</u>	Nd	<0.00020	<u>M</u>	Sc	<0.0010	<u>M</u>	Ti	<0.0050
<u>M</u>	Cd	<0.00030	<u>M</u>	Hf	<0.00020	<u>O</u>	Ni	<0.050	<u>O</u>	Se	<0.40	<u>M</u>	W	0.00055
<u>O</u>	Ca	<0.010	<u>M</u>	Ho	<0.000050	<u>M</u>	Nb	0.00024	<u>O</u>	Si	<0.030	<u>M</u>	U	0.0011
<u>M</u>	Ce	<0.00050	<u>O</u>	In	<0.070	<u>n</u>	Os		<u>M</u>	Ag	0.00044	<u>s</u>	V	
<u>M</u>	Cs	<0.000030	<u>M</u>	Ir	<0.00050	<u>M</u>	Pd	<0.00050	<u>O</u>	Na	<0.090	<u>M</u>	Yb	<0.00010
<u>O</u>	Cr	<0.020	<u>O</u>	Fe	<0.050	<u>i</u>	P		<u>M</u>	Sr	<0.000050	<u>M</u>	Y	<0.0040
<u>O</u>	Co	<0.050	<u>M</u>	La	<0.000050	<u>M</u>	Pt	<0.00020	<u>n</u>	S		<u>M</u>	Zn	0.0041
<u>M</u>	Cu	<0.00060	<u>M</u>	Pb	<0.00030	<u>n</u>	K		<u>M</u>	Ta	<0.00070	<u>M</u>	Zr	<0.00050

M - checked by ICP-MS

O - checked by ICP-OES

i - spectral interference

n - not checked for

s - solution standard element

ANALYZED DENSITY OF SOLUTION (measured at 22°C): 1.015 g/mL

(over)


QA:KLRev.082202DH

Paul R. Gaines

Quality Assurance Manager

EXPIRES

Expires: 12/2004



inorganic ventures / iv labs

195 lehigh avenue, suite 4, lakewood, nj 08701 usa
phone: 800-669-6799 • 732-901-1900 • fax: 732-901-1903
e-mail: ivsales@ivstandards.com • website: www.ivstandards.com

certificate of analysis

010213

1.0 Inorganic Ventures / IV Labs is an ISO Guide 34-2000 Certified Reference Material (CRM) Manufacturer: Certificate #883-02. The certificate is designed and the certified value(s) and uncertainty(ies) are determined in accordance with ISO Guide 31-2000 (Reference Materials - Contents of certificates and label(s), ISO Guide 34-2000 "Quality System Guidelines for the Production of Reference Materials," and ISO Guide 35-1989 "Certification of Reference Materials - General and Statistical Principles."

2.0 DESCRIPTION OF CRM Custom-Grade 1000 µg/mL Zinc in 1.4% (abs) HNO₃

Catalog Number: CGZN1-1, CGZN1-2, and CGZN1-5

Lot Number: **W-ZN02018**

Starting Material: Zn shot

Starting Material Purity (%): 99.999889

Starting Material Lot No J17L26

Matrix: 1.4% (abs) HNO₃

INORGANIC LABS/RADCHEM LABS
DATE RECEIVED: 11/5/03
DATE EXPIRED: 12/1/2004
DATE OPENED: 11/5/03
INORG: 4319 PO: F52256

3.0 CERTIFIED VALUES AND UNCERTAINTIES

Certified Concentration: 1006 ± 3 µg/mL

Certified Density: 1.011 g/mL (measured at 22° C)

The Certified Value is based upon the most precise method used to analyze this CRM. The following equations are used in the calculation of the certified value and the uncertainty:

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

$$\text{Uncertainty } (\pm) = \frac{2(\sum s_i^2)^{1/2}}{(n)^{1/2}}$$

$\bar{x}$ = mean

x_i = individual results

n = number of measurements

$\sum s_i^2$ = The summation of all significant estimated errors

(Most common are the errors from instrumental measurement, weighing, dilution to volume, and the fixed error reported on the NIST SRM certificate of analysis.)

The independent samples t-test was used to determine if there is agreement between the above assay methods at the 95% confidence interval. Both methods were compared and showed agreement within the stated uncertainties. This agreement is a confirmation of the accuracy of this CRM.

4.0 TRACEABILITY TO NIST AND VALUES OBTAINED BY INDEPENDENT METHODS

[] "Property of the result of a measurement or the value of a standard whereby it can be related to stated references, usually national or international standards, through an unbroken chain of comparisons all having stated uncertainties." (ISO VIM, 2nd ed., 1993, definition 6.10)

[] This IV product is Traceable to NIST via direct comparison to NIST SRMs. The uncertainties for each certified value are reported, taking into account the SRM uncertainty error and the measurement, weighing and volume dilution errors.

4.1 Assay Method #1 1002 ± 6 µg/mL

ICP Assay NIST SRM 3168a Lot Number: 001402

Assay Method #2 1006 ± 3 µg/mL

EDTA NIST SRM 928 Lot Number: 880710

4.2 BALANCE CALIBRATION - All balances are checked daily using in-house procedure number 6-IMM-001. The weights used for testing are annually compared to Gerhart Scale Corporation's master weights and are traceable to the National Institute of Standards and Technology (NIST). The NIST Traceability numbers are 692476 - Class 1 and 692476A - Class 2. The NIST test number is 822/260017-98. All analytical balances are calibrated every 4 months by Gerhart Scale Corp. of South Amboy. The balances are calibrated with a class 1 and/or class 2 analytical weight set. These weights are tested annually by a NIST / NVLAP accredited calibration lab. The NIST test number is 822/260017-98.

010216

4.3 THERMOMETER CALIBRATION - The thermometers used in the determination of the final densities are calibrated vs standard thermometer No. 903-2680 which was certified in accordance with the procedures outlined by ASTM E77-87 and NIST Monograph 150 using NIST Test Nos. and Std Nos.: 769543, 217368/769543, 217368/P14452, 176240/P14452, 176240. The in-house procedure No. is 2-QC-001. Thermometers which are not calibrated vs standard thermometer No. 903-2680 are traceable to NIST Identification Nos. 92564, 119016, 471047 and NIST test report Nos. 811/258522, 811/2557078, and 236090.

4.4 GLASSWARE CALIBRATION - In-house procedure 3-QC-002 is used to calibrate all Class A Glassware used in the manufacture and quality control of Custom Grade Standards.

5.0 TRACE METALLIC IMPURITIES (TMI) DETERMINED BY ICP/MS AND ICP-OES IN µg/mL

Custom-Grade solutions are tested for trace metallic impurities by Axial ICP-OES and ICP-MS. The result from the most sensitive method for each element, is reported below. Solutions tested by ICP-MS were analyzed in an ULPA-Filtered Clean Room. An ULPA-Filter is 99.9985% efficient for the removal of particles down to 0.3 µm.

<u>O</u> Al < 0.00200	<u>M</u> Dy < 0.02440	<u>O</u> Li 0.00001	<u>M</u> Pr < 0.00122	<u>M</u> Te < 0.12198
<u>M</u> Sb < 0.00203	<u>M</u> Er < 0.02033	<u>M</u> Lu < 0.00163	<u>M</u> Re < 0.00407	<u>M</u> Tb < 0.00122
<u>M</u> As < 0.04066	<u>M</u> Eu < 0.01220	<u>O</u> Mg 0.00011	<u>M</u> Rh < 0.00407	<u>M</u> Tl < 0.00407
<u>M</u> Ba < 0.04066	<u>M</u> Gd < 0.00407	<u>M</u> Mn < 0.01626	<u>M</u> Rb < 0.00407	<u>M</u> Th < 0.00407
<u>M</u> Be < 0.00203	<u>M</u> Ga < 0.00407	<u>O</u> Hg < 0.01000	<u>M</u> Ru < 0.00813	<u>M</u> Tm < 0.00163
<u>M</u> Bi < 0.00163	<u>M</u> Ge < 0.02440	<u>M</u> Mo < 0.00813	<u>M</u> Sm < 0.00407	<u>M</u> Sn < 0.02033
<u>O</u> B 0.00015	<u>M</u> Au < 0.01220	<u>M</u> Nd < 0.00813	<u>M</u> Sc < 0.04066	<u>M</u> Ti < 0.20331
<u>M</u> Cd < 0.01220	<u>M</u> Hf < 0.00813	<u>O</u> Ni 0.00009	<u>M</u> Se < 0.03253	<u>M</u> W < 0.04066
<u>O</u> Ca 0.00022	<u>M</u> Ho < 0.00203	<u>M</u> Nb < 0.00203	<u>O</u> Si < 0.00400	<u>M</u> U < 0.00813
<u>M</u> Ce < 0.02033	<u>M</u> In < 0.04066	<u>n</u> Os	<u>M</u> Ag < 0.00813	<u>M</u> V < 0.00813
<u>M</u> Cs < 0.00122	<u>M</u> Ir < 0.02033	<u>M</u> Pd < 0.02033	<u>O</u> Na 0.00055	<u>M</u> Yb < 0.00407
<u>O</u> Cr < 0.00100	<u>O</u> Fe 0.00005	<u>O</u> P < 0.00300	<u>M</u> Sr < 0.00203	<u>M</u> Y < 0.16264
<u>M</u> Co < 0.01220	<u>M</u> La < 0.00203	<u>M</u> Pt < 0.00813	<u>O</u> S < 0.02000	<u>s</u> Zn
<u>O</u> Cu < 0.00050	<u>M</u> Pb < 0.01220	<u>O</u> K 0.00018	<u>M</u> Ta < 0.02846	<u>M</u> Zr < 0.02033

M - Checked by ICP-MS O - Checked by ICP-OES I - Spectral Interference n - Not Checked For s - Solution Standard Element

6.0 INTENDED USE

For the calibration of analytical instruments including but not limited to the following:
 ICP-MS, ICP-OES, FAAS, GFAA, XRF, and DCP
 For the validation of analytical methods
 For the preparation of "working reference samples"
 For interference studies and the determination of correction coefficients
 For detection limit and linearity studies
 For additional intended uses, contact IV Technical Staff

7.0 INSTRUCTIONS FOR THE CORRECT USE OF THIS REFERENCE MATERIAL

Storage & Handling - Keep tightly sealed when not in use. Store and use at $20 \pm 4^\circ\text{C}$. Do not pipet from container. Do not return portions removed for pipetting to container.

Atomic Weight; Valence; Coordination Number; Chemical Form in Solution - 65.39; +2; 4; $\text{Zn}(\text{OH})(\text{aq})^{++}$

Chemical Compatibility - Stable in HCl , HNO_3 , H_2SO_4 , HF , H_3PO_4 . Avoid basic media that promotes the formation of insoluble carbonate and hydroxide. Stable with most metals and inorganic anions in acidic media.

Stability - 2-100 ppb levels stable for months in 1% HNO_3 / LDPE container. 1-10,000 ppm solutions chemically stable for years in 1-5% HNO_3 / LDPE container.

Zn Containing Samples (Preparation and Solution) - Metal (Soluble in HNO_3); Oxides (Soluble in HCl); Ores (Dissolve in HCl / HNO_3); Organic based (Dry ash at 450°C and dissolve ash in HCl) (Sulfuric/peroxide acid digestion)

Atomic Spectroscopic Information (ICP-OES D.L.s are given as radial/axial view):

Technique/Line	Estimated D.L.	Order	Type	Interferences (underlined indicates severe at 100ppb)
ICP-OES 213.856 nm	0.002 / 0.0004 $\mu\text{g/mL}$	1 atom	Ni, Cu, V	
ICP-OES 202.548 nm	0.004 / 0.0002 $\mu\text{g/mL}$	1 ion	Nb, Cu, Co, Hf	
ICP-OES 206.200 nm	0.006 / 0.0006 $\mu\text{g/mL}$	1 ion	Sb, Ta, Bi, Os	
ICP-MS 66 amu	7 ppt	n/a	M'	$^{44}\text{Ti}^{+}\text{O}$, $^{52}\text{Cr}^{+}\text{O}$, $^{50}\text{V}^{+}\text{O}$, $^{34}\text{S}^{+}\text{O}_2$, $^{32}\text{S}^{+}\text{O}^{+}\text{O}$, $^{33}\text{S}^{+}\text{O}_2$, $^{34}\text{S}^{+}\text{O}^{+}\text{O}$, $^{32}\text{S}^{+}\text{S}$, $^{34}\text{S}_2$

8.0 HAZARDOUS INFORMATION - Please refer to the enclosed Material Safety Data sheet for information regarding this CRM.

9.0 HOMOGENEITY - This solution was mixed according to procedure IV-MPM-004 and is guaranteed to be homogeneous.

10.0 QUALITY STANDARD DOCUMENTATION

10.1 ISO 9001:2000 Quality Management System Registration - QMI Certificate Number 010105

Recognized by:

Registrar Accreditation Board (ANSI-RAB)

Standards Council of Canada (SCC)

Dutch Council for Accreditation (RVA)

Entidad Mexicana de Acreditacion, a.c.(EMA)

Members of **IQ Net International Certification Network**:

Argentina (IRAM), Australia (QAS), Austria (OQS), Belgium (Avinter), Brazil (FCAV), Canada (QMI), Hong Kong (HKQAA), Columbia (ICONTEC), Czech Republic (CQS), Denmark (DS), Finland (SFS), France (AFAQ), Germany (DQS), Greece (ELOT), Hungary (MSZT), Ireland (NSAI), Israel (SII), Italy (CISQ), Japan (JQA), Korea (KSA-QA), Netherlands (KEMA), Norway (NCS), Poland (PCBC), Portugal (APCER), Singapore (PSB), Slovenia (SIQ), Spain (AENOR), Switzerland (SQS)

10.2 ISO/IEC 17025 - 1999 "General Requirements for the Competence of Testing and Calibration"

- Chemical Testing - Accredited A2LA Certificate Number 883.01



10.3 ISO/IEC Guide 34 - 2000 "General Requirements for the Competence of Reference Material Producers"

- Reference Materials Production - Accredited A2LA Certificate Number 883.02

A2LA Mutual Recognition Agreement Partners:

Australia (NATA), Austria (BmWA), Belgium (BELTEST) (BKO-OBE), Canada (SCC), Chinese Taipei (CNLA), Czech Republic (NAO), Denmark (DANAK), Finland (FINAS), France (COFRAC), Germany (DAR), Hong Kong (HKAS), Ireland (NAB), Italy (SIT) (SINAL), Japan (JAB) (JNLA), Republic of Korea (KOLAS), The Netherlands (RvA), New Zealand (IANZ), Norway (NA), Portugal (IPQ), Singapore (SAC-SINGLAS), Spain (ENAC), Sweden (SWEDAC), Switzerland (SAS), United Kingdom (UKAS) and United States (NVLAP) (ICBO ES)

10.4 10CFR50 Appendix B - Nuclear Regulatory Commission

- Domestic Licensing of Production and Utilization Facilities

10.5 10CFR21 - Nuclear Regulatory Commission - Reporting Defects and Non-Compliance

10.6 MIL-STD-45662A (Obsolete/Observed)

11.0 DATE OF CERTIFICATION AND PERIOD OF VALIDITY

010218



11.1 IV Shelf Life - The period of time during which the concentration of the analyte(s) in a properly packaged, unopened, and unused standard stored under environmentally controlled and monitored conditions will remain within the specified uncertainty range. Shelf life is limited primarily by transpiration (loss of water from the solution) and infrequently, by chemical instability. Transpiration studies (P-SP01020) of chemically-stable solutions performed at Inorganic Ventures / IV Labs indicate a CRM shelf-life of four years for solutions packaged in 500-mL low density polyethylene bottles. When stored under special conditions that minimize transpiration and instability, the shelf life can be extended past this limit.

11.2 Expiration Date - The date after which a CRM should not be used. Routine laboratory use of a CRM increases transpiration losses and the chance of contamination which affect the integrity of the CRM and limit its useful life. Inorganic Ventures / IV Labs concurs with state and federal regulatory agencies' recommendations that solution standards be assigned a one-year expiration date.

Certification Date: May 02, 2003

Expiration Date:

EXPIRES
1st 2004

12.0 NAMES AND SIGNATURES OF CERTIFYING OFFICERS

Certificate Prepared By: JoAnn Struthers, QA Administrative Assistant

JoAnn Struthers

Certificate Approved By: Katalin Le, QC Supervisor

Katalin Le

Certifying Officer: Paul Gaines, Chemist, Senior Technical Director

Paul Gaines

INORGANIC LABS/RADCHEM LABS
DATE RECEIVED: 11/5/03
DATE EXPIRED: 12/1/2004
DATE OPENED: 11/5/03
INORG: 4319 PO: F52258

inorganic ventures / iv labs

195 lehigh avenue, suite 4, lakewood, nj 08701 usa
phone: 800-669-6799 • 732-901-1900 • fax: 732-901-1903
e-mail: ivsales@ivstandards.com • website: www.ivstandards.com

certificate of analysis 010219

CUSTOM-GRADE SOLUTION 1000 µg/mL Silver in 3.5% HNO₃ (abs)

Catalog Number: CGAG1-1, CGAG1-2 and CGAG1-5

Lot Number: **T-AG02015**

Starting Material:
Starting Material Purity:
Starting Material Lot No:

Silver Metal
99.999%
F15I02

INORGANIC LABS/RADCHEM LABS
DATE RECEIVED: 08/26/03
DATE EXPIRED: 09/01/2004
DATE OPENED: 08/26/03
INORG: 4222 PO: F52224

CERTIFIED CONCENTRATION: 1001 ± 2 µg/mL

The Certified Value is based upon the most precise method used to analyze this CRM. The following equations are used in the calculation of the certified value and the uncertainty:

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

$$\text{Uncertainty } (\pm) = \frac{2((\sum s_i)^2)^{1/2}}{(n)^{1/2}}$$

($\bar{x}$) = mean

x_i = individual results

n = number of measurements

$\sum S_i$ = The summation of all significant estimated errors.

Classical Wet Assay: 1004 ± 3 µg/mL

Method: Volhard Titration vs NIST SRM 999a Potassium Chloride

Instrument Analysis: 1001 ± 2 µg/mL

Method: Inductively Coupled Plasma Spectroscopy (ICP) vs NIST SRM 3151.

The independent samples t-test was used to determine if there is agreement between the above assay methods at the 95% confidence interval. Both methods were compared and showed agreement within the stated uncertainties. This agreement is a confirmation of the accuracy of this CRM.

TRACE METALLIC IMPURITIES DETERMINED BY ICP-MS AND ICP-OES IN µg/mL:

Custom-Grade solutions tested for trace metallic impurities by ICP-MS were analyzed in an ULPA-Filtered Clean Room.

An ULPA-Filter is 99.9985% efficient for the removal of particles down to 0.3 µm.

<u>Q</u> Al <0.00010	<u>M</u> Dy <0.00060	<u>Q</u> Li <0.000030	<u>M</u> Pr <0.000030	<u>Q</u> Te <0.030
<u>M</u> Sb <0.000050	<u>M</u> Er <0.00050	<u>M</u> Lu <0.000040	<u>M</u> Re <0.00010	<u>M</u> Tb <0.000030
<u>Q</u> As <0.0050	<u>M</u> Eu <0.00030	<u>Q</u> Mg <0.000040	<u>M</u> Rh <0.00010	<u>M</u> Tl <0.00010
<u>M</u> Ba <0.0010	<u>M</u> Gd <0.00010	<u>Q</u> Mn <0.00030	<u>M</u> Rb <0.00010	<u>M</u> Th <0.00010
<u>Q</u> Be <0.00050	<u>M</u> Ga <0.00010	<u>Q</u> Hg 0.00090	<u>M</u> Ru <0.00020	<u>M</u> Tm <0.000040
<u>M</u> Bi <0.000040	<u>M</u> Ge <0.00060	<u>M</u> Mo <0.00020	<u>M</u> Sm <0.00010	<u>M</u> Sn <0.00050
<u>Q</u> B <0.0020	<u>Q</u> Au <0.012	<u>M</u> Nd <0.00020	<u>M</u> Sc <0.0010	<u>Q</u> Ti <0.00070
<u>Q</u> Cd <0.0020	<u>M</u> Hf <0.00020	<u>Q</u> Ni <0.0070	<u>Q</u> Se <0.036	<u>M</u> W <0.0010
<u>Q</u> Ca <0.000050	<u>M</u> Ho <0.000050	<u>Q</u> Nb <0.000050	<u>Q</u> Si <0.0030	<u>M</u> U <0.00020
<u>M</u> Ce <0.00050	<u>Q</u> In <0.020	<u>n</u> Os	<u>s</u> Ag	<u>M</u> V <0.00020
<u>M</u> Cs <0.000030	<u>M</u> Ir <0.00050	<u>M</u> Pd <0.00050	<u>Q</u> Na <0.090	<u>M</u> Yb <0.00010
<u>Q</u> Cr <0.0020	<u>Q</u> Fe <0.00070	<u>Q</u> P <0.030	<u>M</u> Sr <0.000050	<u>M</u> Y <0.0040
<u>M</u> Co <0.00030	<u>M</u> La <0.000050	<u>M</u> Pt <0.00020	<u>Q</u> S <0.020	<u>Q</u> Zn <0.0010
<u>M</u> Cu <0.00060	<u>M</u> Pb <0.00030	<u>Q</u> K <0.0060	<u>M</u> Ta <0.00070	<u>M</u> Zr <0.00050

M - checked by ICP-MS

Q - checked by ICP-OES

i - spectral interference

n - not checked for

s - solution standard element

ANALYZED DENSITY OF SOLUTION (measured at 22°C): 1.024 g/mL

QA:KL Rev.082102DN

Paul R. Laines

Quality Assurance Manager

Expires:

EXPIRES

122004



Certificate of Analysis

010220

CUSTOM-GRADE SOLUTION 1000 µg/mL Arsenic in 1.4% HNO₃ (abs)

Catalog Number: CGAS1-1, CGAS1-2 and CGAS1-5

Lot Number: **T-AS02019**

Starting Material:
Starting Material Purity:
Starting Material Lot No:

Arsenic Metal
99.999%
23014

INORGANIC LABS/RADCHEM LABS

DATE RECEIVED: 01/09/03

DATE EXPIRED: 02/01/2004

DATE OPENED: 01/10/03

INORG: 3899 PO: F52071

CERTIFIED CONCENTRATION: 1006 ± 3 µg/mL

The Certified Value is based upon the most precise method used to analyze this CRM. The following equations are used in the calculation of the certified value and the uncertainty:

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

$$\text{Uncertainty } (\pm) = \frac{2[(\sum s_i)^2]^{1/2}}{(n)^{1/2}}$$

($\bar{x}$) = mean

x_i = individual results

n = number of measurements

$\sum S_i$ = The summation of all significant estimated errors.

Calculated Value: 1003 µg/mL

Method: Calculated, based on starting material.

Instrument Analysis: 1006 ± 3 µg/mL (Average of 3 runs)

Method: Inductively Coupled Plasma Spectroscopy (ICP) vs NIST SRM 3103a.

TRACE METALLIC IMPURITIES DETERMINED BY ICP-MS AND ICP-OES IN µg/mL:

Custom-Grade solutions tested for trace metallic impurities by ICP-MS were analyzed in an ULPA-Filtered Clean Room.

An ULPA-Filter is 99.9985% efficient for the removal of particles down to 0.3 µm.

<u>Q</u>	Al	0.00030	<u>M</u>	Dy	<0.012	<u>Q</u>	Li	<0.000020	<u>M</u>	Pr	<0.00060	<u>M</u>	Te	<0.060
<u>Q</u>	Sb	<0.010	<u>M</u>	Er	<0.010	<u>M</u>	Lu	<0.00080	<u>Q</u>	Re	<0.010	<u>M</u>	Tb	<0.00060
<u>s</u>	As		<u>M</u>	Eu	<0.0060	<u>Q</u>	Mg	<0.000030	<u>M</u>	Rh	<0.0020	<u>M</u>	Tl	<0.0020
<u>M</u>	Ba	<0.020	<u>M</u>	Gd	<0.0020	<u>Q</u>	Mn	<0.000030	<u>M</u>	Rb	<0.0020	<u>M</u>	Th	<0.0020
<u>M</u>	Be	<0.0010	<u>M</u>	Ga	<0.0020	<u>Q</u>	Hg	<0.012	<u>M</u>	Ru	<0.0040	<u>M</u>	Tm	<0.00080
<u>M</u>	Bi	<0.00080	<u>M</u>	Ge	<0.012	<u>M</u>	Mo	<0.0040	<u>M</u>	Sm	<0.0020	<u>Q</u>	Sn	0.00053
<u>Q</u>	B	<0.012	<u>M</u>	Au	<0.0060	<u>M</u>	Nd	<0.0040	<u>M</u>	Sc	<0.020	<u>M</u>	Ti	<0.10
<u>M</u>	Cd	<0.0060	<u>M</u>	Hf	<0.0040	<u>M</u>	Ni	<0.016	<u>M</u>	Se	<0.016	<u>M</u>	W	<0.020
<u>Q</u>	Ca	<0.0015	<u>M</u>	Ho	<0.0010	<u>Q</u>	Nb	<0.0020	<u>Q</u>	Si	0.0077	<u>M</u>	U	<0.0040
<u>M</u>	Ce	<0.010	<u>M</u>	In	<0.020	<u>n</u>	Os		<u>M</u>	Ag	<0.0040	<u>M</u>	V	<0.0040
<u>M</u>	Cs	<0.00060	<u>M</u>	Ir	<0.010	<u>M</u>	Pd	<0.010	<u>Q</u>	Na	0.00096	<u>M</u>	Yb	<0.0020
<u>M</u>	Cr	<0.010	<u>Q</u>	Fe	<0.0011	<u>Q</u>	P	<0.0026	<u>M</u>	Sr	<0.0010	<u>M</u>	Y	<0.080
<u>M</u>	Co	<0.0060	<u>M</u>	La	<0.0010	<u>M</u>	Pt	<0.0040	<u>Q</u>	S	<0.025	<u>Q</u>	Zn	0.00019
<u>M</u>	Cu	<0.012	<u>M</u>	Pb	<0.0060	<u>Q</u>	K	0.00096	<u>M</u>	Ta	<0.014	<u>M</u>	Zr	<0.010

M - checked by ICP-MS

O - checked by ICP-OES

i - spectral interference

n - not checked for

s - solution standard element

ANALYZED DENSITY OF SOLUTION (measured at 22°C): 1.010 g/mL

QA:KSL Rev.061202DN

(over)

Inorganic Ventures, Inc.

195 Lehigh Avenue • Suite 4 • Lakewood, NJ 08701

Orders: 800-669-6799 • FAX (732) 901-1903

Technical Support: 800-569-6799

Quality Assurance Manager

EXPIRES

1/1/2004



inorganic ventures / iv labs

195 lehigh avenue, suite 4, lakewood, nj 08701 usa
phone: 800-669-6799 • 732-901-1900 • fax: 732-901-1903
e-mail: ivsales@ivstandards.com • website: www.ivstandards.com

010221 certificate of analysis

- 1.0 Inorganic Ventures / IV Labs is an ISO Guide 34-2000 Certified Reference Material (CRM) Manufacturer: Certificate #883-02.** The certificate is designed and the certified value(s) and uncertainty(ies) are determined in accordance with ISO Guide 31-2000 (Reference Materials - Contents of certificates and label(s)), ISO Guide 34-2000 "Quality System Guidelines for the Production of Reference Materials," and ISO Guide 35-1989 "Certification of Reference Materials - General and Statistical Principles."

2.0 DESCRIPTION OF CRM Custom-Grade 1000 µg/mL Lead in 0.35% (abs) HNO₃

Catalog Number: CGPB1-1, CGPB1-2, and CGPB1-5
Lot Number: **W-PB02114**
Starting Material: Pb(NO₃)₂
Starting Material Purity (%): 99.999974
Starting Material Lot No 22150
Matrix: 0.35% (abs) HNO₃

3.0 CERTIFIED VALUES AND UNCERTAINTIES

Certified Concentration: 1006 ± 2 µg/mL
Certified Density: 1.002 g/mL (measured at 22° C)

The Certified Value is the wet assay value. The following equations are used in the calculations of the certified value and the uncertainty:

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

$$\text{Uncertainty } (\pm) = \frac{2[(\sum s_i^2)]^{1/2}}{(n)^{1/2}}$$

($\bar{x}$) = mean

x_i = individual results

n = number of measurements

$\sum s_i^2$ = The summation of all significant estimated errors.

(Most common are the errors from instrumental measurement, weighing, dilution to volume, and the fixed error reported on the NIST SRM certificate of analysis.)

The independent samples t-test was used to determine if there is agreement between the above assay methods at the 95% confidence interval. Both methods were compared and showed agreement within the stated uncertainties. This agreement is a confirmation of the accuracy of this CRM.

4.0 TRACEABILITY TO NIST AND VALUES OBTAINED BY INDEPENDENT METHODS

□ "Property of the result of a measurement or the value of a standard whereby it can be related to stated references, usually national or international standards, through an unbroken chain of comparisons all having stated uncertainties." (ISO VIM, 2nd ed., 1993, definition 6.10)

□ This IV product is Traceable to NIST via direct comparison to NIST SRMs. The uncertainties for each certified value are reported, taking into account the SRM uncertainty error and the measurement, weighing and volume dilution errors.

4.1 Assay Method #1 1005 ± 2 µg/mL
ICP Assay NIST SRM 3128 Lot Number: 991504

Assay Method #2 1006 ± 2 µg/mL
EDTA NIST SRM 928 Lot Number: 880710

INORGANIC LABS/RADCHEM LABS 142
DATE RECEIVED: 11/3/03
DATE EXPIRED: 11/1/04
DATE OPENED: 11/3/03
INORG: 4313
PO: F52258

4.2 BALANCE CALIBRATION - All balances are checked daily using in-house procedure number 6-IMM-001. The weights used for testing are annually compared to Gerhart Scale Corporation's master weights and are traceable to the National Institute of Standards and Technology (NIST). The NIST Traceability numbers are 692476 - Class 1 and 692476A - Class 2. The NIST test number is 822/260017-98. All analytical balances are calibrated every 4 months by Gerhart Scale Corp. of South Amboy. The balances are calibrated with a class 1 and/or class 2 analytical weight set. These weights are tested annually by a NIST / NVLAP accredited calibration lab. The NIST test number is 822/260017-98.

010222

4.3 THERMOMETER CALIBRATION - The thermometers used in the determination of the final densities are calibrated vs standard thermometer No. 903-2680 which was certified in accordance with the procedures outlined by ASTM E77-87 and NIST Monograph 150 using NIST Test Nos. and Std Nos.: 769543, 217368/769543, 217368/P14452, 176240/P14452, 176240. The in-house procedure No. is 2-QC-001. Thermometers which are not calibrated vs standard thermometer No. 903-2680 are traceable to NIST Identification Nos. 92564, 119016, 471047 and NIST test report Nos. 811/258522, 811/2557078, and 236090.

4.4 GLASSWARE CALIBRATION - In-house procedure 3-QC-002 is used to calibrate all Class A Glassware used in the manufacture and quality control of Custom Grade Standards.

5.0 TRACE METALLIC IMPURITIES (TMI) DETERMINED BY ICP/MS AND ICP-OES IN µg/mL

Custom-Grade solutions are tested for trace metallic impurities by Axial ICP-OES and ICP-MS. The result from the most sensitive method for each element, is reported below. Solutions tested by ICP-MS were analyzed in an ULPA-Filtered Clean Room. An ULPA-Filter is 99.9985% efficient for the removal of particles down to 0.3 µm.

<u>Q</u> Al < 0.00270	<u>M</u> Dy < 0.01193	<u>Q</u> Li < 0.00002	<u>M</u> Pr < 0.00060	<u>M</u> Te < 0.05965
<u>M</u> Sb < 0.00099	<u>M</u> Er < 0.00994	<u>M</u> Lu < 0.00080	<u>M</u> Re < 0.00199	<u>M</u> Tb < 0.00060
<u>M</u> As < 0.01989	<u>M</u> Eu < 0.00597	<u>Q</u> Mg 0.00008	<u>Q</u> Rh < 0.00900	<u>Q</u> Tl 0.00130
<u>M</u> Ba < 0.01989	<u>M</u> Gd < 0.00199	<u>M</u> Mn < 0.00795	<u>M</u> Rb < 0.00199	<u>M</u> Th < 0.00199
<u>M</u> Be < 0.00099	<u>M</u> Ga < 0.00199	<u>Q</u> Hg < 0.01500	<u>M</u> Ru < 0.00398	<u>M</u> Tm < 0.00080
<u>Q</u> Bi < 0.02000	<u>M</u> Ge < 0.01193	<u>M</u> Mo < 0.00398	<u>M</u> Sm < 0.00199	<u>M</u> Sn < 0.00994
<u>Q</u> B < 0.04000	<u>M</u> Au < 0.00597	<u>M</u> Nd < 0.00398	<u>M</u> Sc < 0.01989	<u>M</u> Ti < 0.09942
<u>M</u> Cd < 0.00597	<u>M</u> Hf < 0.00398	<u>M</u> Ni < 0.01591	<u>M</u> Se < 0.01591	<u>M</u> W < 0.01989
<u>Q</u> Ca 0.00009	<u>M</u> Ho < 0.00099	<u>M</u> Nb < 0.00099	<u>Q</u> Si < 0.00340	<u>M</u> U < 0.00398
<u>M</u> Ce < 0.00994	<u>M</u> In < 0.01989	<u>n</u> Os	<u>M</u> Ag < 0.00398	<u>M</u> V < 0.00398
<u>M</u> Cs < 0.00060	<u>M</u> Ir < 0.00994	<u>M</u> Pd < 0.00994	<u>Q</u> Na < 0.00600	<u>M</u> Yb < 0.00199
<u>M</u> Cr < 0.00994	<u>Q</u> Fe 0.00011	<u>Q</u> P < 0.00500	<u>M</u> Sr < 0.00099	<u>M</u> Y < 0.07954
<u>M</u> Co < 0.00597	<u>M</u> La < 0.00099	<u>M</u> Pt < 0.00398	<u>Q</u> S < 0.10000	<u>M</u> Zn < 0.03977
<u>M</u> Cu < 0.01193	<u>s</u> Pb	<u>Q</u> K < 0.00180	<u>M</u> Ta < 0.01392	<u>M</u> Zr < 0.00994

M - Checked by ICP-MS O - Checked by ICP-OES i - Spectral Interference n - Not Checked For s - Solution Standard Element

6.0 INTENDED USE

For the calibration of analytical instruments including but not limited to the following:

ICP-MS, ICP-OES, FAAS, GFAA, XRF, and DCP

For the validation of analytical methods

For the preparation of "working reference samples"

For interference studies and the determination of correction coefficients

For detection limit and linearity studies

For additional intended uses, contact IV Technical Staff

7.0 INSTRUCTIONS FOR THE CORRECT USE OF THIS REFERENCE MATERIAL

Storage & Handling - Keep tightly sealed when not in use. Store and use at $20 \pm 4^\circ\text{C}$. Do not pipet from container. Do not return portions removed for pipetting to container.

Atomic Weight; Valence; Coordination Number; Chemical Form in Solution - 207.2; +2; 6; $\text{Pb}(\text{H}_2\text{O})_6^{2+}$

010223

Chemical Compatibility - Soluble in HCl, HF and HNO_3 . Avoid H_2SO_4 . Stable with most metals and inorganic anions forming insoluble carbonate, borate, sulfate, sulfite, sulfide, phosphate, oxalate, chromate, tennate, iodate, and cyanide in neutral aqueous media.

Stability - 2-100 ppb levels stable for months in 1% HNO_3 / LDPE container. 1-10,000 ppm solutions chemically stable for years in 2-5% HNO_3 / LDPE container.

Pb Containing Samples (Preparation and Solution) - Metal (Best dissolved in 1:1 H_2O / HNO_3). Oxides (The many different Pb oxides are soluble in HNO_3 with the exception of PbO_2 which is soluble in HCl or HF); Ores and Alloys (Best attacked using 1:1 H_2O / HNO_3); Organic Matrices (Dry ash and dissolve in dilute HCl).

Atomic Spectroscopic Information (ICP-OES D.L.s are given as radial/axial view):

Technique/Line	Estimated D.L.	Order	Type	Interferences (underlined indicates severe at λ concs.)
ICP-OES 168.215 nm	0.03 / 0.003 $\mu\text{g/mL}$	1	ion	Co
ICP-OES 220.353 nm	0.04 / 0.006 $\mu\text{g/mL}$	1	ion	Bi, Nb
ICP-OES 217.000 nm	0.09 / 0.03 $\mu\text{g/mL}$	1	atom	W, Ir, Hf, Sb, Th
ICP-MS 208 amu	5 ppt	n/a	M'	^{202}Pb , ^{208}Pb

8.0 **HAZARDOUS INFORMATION** - Please refer to the enclosed Material Safety Data sheet for information regarding this CRM.

9.0 **HOMOGENEITY** - This solution was mixed according to procedure IV-MPM-004 and is guaranteed to be homogeneous.

10.0 QUALITY STANDARD DOCUMENTATION

10.1 ISO 9001:2000 Quality Management System Registration - QMI Certificate Number 010105

Recognized by:

Registrar Accreditation Board (ANSI-RAB)

Standards Council of Canada (SCC)

Dutch Council for Accreditation (RVA)

Entidad Mexicana de Acreditacion, a.c.(EMA)

Members of IQ Net International Certification Network:

Argentina (IRAM), Australia (QAS), Austria (ÖQS), Belgium (Avinter), Brazil (FCAV), Canada (QMI), Hong Kong (HKQAA), Columbia (ICONTEC), Czech Republic (CQS), Denmark (DS), Finland (SFS), France (AFAQ), Germany (DQS), Greece (ELOT), Hungary (MSZT), Ireland (NSAI), Israel (SII), Italy (CISQ), Japan (JQA), Korea (KSA-QA), Netherlands (KEMA), Norway (NCS), Poland (PCBC), Portugal (APCER), Singapore (PSB), Slovenia (SIQ), Spain (AENOR), Switzerland (SQS)

10.2 ISO/IEC 17025 - 1999 "General Requirements for the Competence of Testing and Calibration"

- Chemical Testing - Accredited A2LA Certificate Number 883.01



10.3 ISO/IEC Guide 34 - 2000 "General Requirements for the Competence of Reference Material Producers"

- Reference Materials Production - Accredited A2LA Certificate Number 883.02

A2LA Mutual Recognition Agreement Partners:

Australia (NATA), Austria (BmWA), Belgium (BELTEST) (BKO-OBE), Canada (SCC), Chinese Taipei (CNLA), Czech Republic (NAO), Denmark (DANAK), Finland (FINAS), France (COFRAC), Germany (DAR), Hong Kong (HKAS), Ireland (NAB), Italy (SIT) (SINAL), Japan (JAB) (JNLA), Republic of Korea (KOLAS), The Netherlands (RvA), New Zealand (IANZ), Norway (NA), Portugal (IPQ), Singapore (SAC-SINGLAS), Spain (ENAC), Sweden (SWEDAC), Switzerland (SAS), United Kingdom (UKAS) and United States (NVLAP) (ICBO ES)

10.4 10CFR50 Appendix B - Nuclear Regulatory Commission

- Domestic Licensing of Production and Utilization Facilities

10.5 10CFR21 - Nuclear Regulatory Commission - Reporting Defects and Non-Compliance

10.6 MIL-STD-45662A (Obsolete/Observed)

11.0 DATE OF CERTIFICATION AND PERIOD OF VALIDITY



11.1 IV Shelf Life - The period of time during which the concentration of the analyte(s) in a properly packaged, unopened, and unused standard stored under environmentally controlled and monitored conditions will remain within the specified uncertainty range. Shelf life is limited primarily by transpiration (loss of water from the solution) and infrequently, by chemical instability. Transpiration studies (P-SP01020) of chemically-stable solutions performed at Inorganic Ventures / IV Labs indicate a CRM shelf-life of four years for solutions packaged in 500-mL low density polyethylene bottles. When stored under special conditions that minimize transpiration and instability, the shelf life can be extended past this limit.

11.2 Expiration Date - The date after which a CRM should not be used. Routine laboratory use of a CRM increases transpiration losses and the chance of contamination which affect the integrity of the CRM and limit its useful life. Inorganic Ventures / IV Labs concurs with state and federal regulatory agencies' recommendations that solution standards be assigned a one-year expiration date.

INORGANIC LABS/RADCHEM LABS 2 of 2
 DATE RECEIVED: 11/13/03
 DATE EXPIRED: 11/11/04 DL
 DATE OPENED: 11/13/03
 INORG: 4313 PO: F52258

Certification Date: January 23, 2003

Expiration Date: **EXPIRES**
 11/2004

12.0 NAMES AND SIGNATURES OF CERTIFYING OFFICERS

Certificate Prepared By: JoAnn Struthers, QA Administrative Assistant

Certificate Approved By: Katalin Le, QC Supervisor

Certifying Officer: Paul Gaines, Chemist, Senior Technical Director

JoAnn Struthers
Katalin Le
Paul Gaines



010225

Certificate of Analysis

CUSTOM-GRADE SOLUTION 1000 µg/mL Antimony in 0.7% HNO₃ (abs) / 3% Tartaric Acid
Catalog Number: CGSB1-1, CGSB1-2 and CGSB1-5

Lot Number: **T-SB02074**

Starting Material:
Starting Material Purity:
Starting Material Lot No:

Antimony Metal
99.995%
E14H25

INORGANIC LABS/RADCHEM LABS

DATE RECEIVED: 02/13/03

DATE EXPIRED: 03/01/2004

DATE OPENED: 02/13/03

INORG: 3953 PO: F52060

CERTIFIED CONCENTRATION: 1003 ± 2 µg/mL

The Certified Value is the instrument analysis value. The following equations are used in the calculation of the certified value and the uncertainty:

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

$$\text{Uncertainty } (\pm) = \frac{2[(\sum s_i)^2]^{1/2}}{(n)^{1/2}}$$

($\bar{x}$) = mean

x_i = individual results

n = number of measurements

$\sum S_i$ = The summation of all significant estimated errors.

Calculated Value: 999 µg/mL

Method: Calculated, based on starting material.

Instrument Analysis: 1003 ± 2 µg/mL (Average of 2 runs)

Method: Inductively Coupled Plasma Spectroscopy (ICP) vs NIST SRM 3102a.

TRACE METALLIC IMPURITIES DETERMINED BY ICP-MS AND ICP-OES IN µg/mL:

Custom-Grade solutions tested for trace metallic impurities by ICP-MS were analyzed in an ULPA-Filtered Clean Room.

An ULPA-Filter is 99.9985% efficient for the removal of particles down to 0.3 µm.

<u>M</u> Al <0.0050	<u>M</u> Dy <0.00060	<u>M</u> Li <0.0010	<u>M</u> Pr <0.000030	<u>n</u> Te
<u>s</u> Sb	<u>M</u> Er <0.00050	<u>M</u> Lu <0.000040	<u>M</u> Re <0.00010	<u>M</u> Tb <0.000030
<u>M</u> As <0.0010	<u>M</u> Eu <0.00030	<u>M</u> Mg <0.0030	<u>M</u> Rh <0.00010	<u>M</u> Tl 0.00030
<u>Q</u> Ba <0.020	<u>M</u> Gd <0.00010	<u>M</u> Mn 0.0080	<u>M</u> Rb <0.00010	<u>M</u> Th <0.00010
<u>M</u> Be <0.000050	<u>M</u> Ga <0.00010	<u>Q</u> Hg <0.030	<u>M</u> Ru <0.00020	<u>M</u> Tm <0.000040
<u>M</u> Bi 0.00070	<u>M</u> Ge <0.00060	<u>M</u> Mo <0.00020	<u>M</u> Sm <0.00010	<u>n</u> Sn
<u>M</u> B <0.0070	<u>M</u> Au <0.00030	<u>M</u> Nd <0.00020	<u>M</u> Sc <0.0010	<u>M</u> Ti <0.0050
<u>M</u> Cd <0.00030	<u>M</u> Hf <0.00020	<u>Q</u> Ni <0.050	<u>Q</u> Se <0.40	<u>M</u> W <0.0010
<u>Q</u> Ca <0.020	<u>M</u> Ho <0.000050	<u>M</u> Nb <0.000050	<u>Q</u> Si <0.020	<u>M</u> U <0.00020
<u>M</u> Ce 0.0018	<u>Q</u> In <0.030	<u>n</u> Os	<u>M</u> Ag <0.00020	<u>M</u> V <0.00020
<u>M</u> Cs <0.000030	<u>M</u> Ir <0.00050	<u>M</u> Pd <0.00050	<u>Q</u> Na <0.090	<u>M</u> Yb <0.00010
<u>M</u> Cr 0.031	<u>Q</u> Fe <0.050	<u>Q</u> P <0.050	<u>M</u> Sr <0.000050	<u>M</u> Y <0.0040
<u>M</u> Co <0.00030	<u>M</u> La <0.040	<u>M</u> Pt <0.00020	<u>n</u> S	<u>M</u> Zn <0.0020
<u>M</u> Cu 0.0016	<u>M</u> Pb 0.0011	<u>n</u> K	<u>M</u> Ta <0.00070	<u>M</u> Zr <0.00050

M - checked by ICP-MS

O - checked by ICP-OES

i - spectral interference

n - not checked for

s - solution standard element

ANALYZED DENSITY OF SOLUTION (measured at 22°C): 1.019 g/mL

(over)

QA:KL Rev.0822020H

Inorganic Ventures, Inc.

195 Lehigh Avenue • Suite 4 • Lakewood, NJ 08701

Orders: 800-669-6799 • FAX (732) 901-1903

Technical Support: 800-569-6799

Quality Assurance Manager

Expires:

EXPIRES

12/2004



inorganic ventures / iv labs

195 lehigh avenue, suite 4, lakewood, nj 08701 usa
phone: 800-669-6799 • 732-901-1900 • fax: 732-901-1903
e-mail: ivsales@ivstandards.com • website: www.ivstandards.com

certificate of analysis 010228

1.0 Inorganic Ventures / IV Labs is an ISO Guide 34-2000 Certified Reference Material (CRM) Manufacturer: Certificate #883-02. The certificate is designed and the certified value(s) and uncertainty(ies) are determined in accordance with ISO Guide 31-2000 (Reference Materials - Contents of certificates and label(s)), ISO Guide 34-2000 "Quality System Guidelines for the Production of Reference Materials," and ISO Guide 35-1989 "Certification of Reference Materials - General and Statistical Principles."

2.0 DESCRIPTION OF CRM Custom-Grade 1000 µg/mL Selenium in 1.4% (abs) HNO₃

Catalog Number: CGSE1-1, CGSE1-2, and CGSE1-5

Lot Number: T-SE01102

Starting Material: Se shot

Starting Material Purity (%): 99.9971

Starting Material Lot No C09L08

Matrix: 1.4% (abs) HNO₃

INORGANIC LABS/RADCHEM LABS Pg. 1 of 2

DATE RECEIVED: 06/20/03

DATE EXPIRED: 07/01/2004

DATE OPENED: 06/23/03

INORG: 4152 PO: T52370

3.0 CERTIFIED VALUES AND UNCERTAINTIES

Certified Concentration: 995 ± 3 µg/mL

Certified Density: 1.010 g/mL (measured at 22° C)

The Certified Value is based upon the most precise method used to analyze this CRM. The following equations are used in the calculation of the certified value and the uncertainty:

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

$\bar{x}$ = mean

x_i = individual results

n = number of measurements

$$\text{Uncertainty } (\pm) = \frac{2[(\sum s_i)^2]^{1/2}}{(n)^{1/2}}$$

$\sum s_i$ = The summation of all significant estimated errors

(Most common are the errors from instrumental measurement, weighing, dilution to volume, and the fixed error reported on the NIST SRM certificate of analysis.)

4.0 TRACEABILITY TO NIST AND VALUES OBTAINED BY INDEPENDENT METHODS

□ "Property of the result of a measurement or the value of a standard whereby it can be related to stated references, usually national or international standards, through an unbroken chain of comparisons all having stated uncertainties." (ISO VIM, 2nd ed., 1993, definition 6.10)

□ This IV product is Traceable to NIST via direct comparison to NIST SRMs. The uncertainties for each certified value are reported, taking into account the SRM uncertainty error and the measurement, weighing and volume dilution errors.

4.1 Assay Method #1 995 ± 3 µg/mL (Avg. 2 runs)

ICP Assay NIST SRM 3149 Lot Number: 992106

Assay Method #2 1002 µg/mL

Gravimetric NIST SRM Lot Number: See Sec. 4.2

4.2 BALANCE CALIBRATION - All balances are checked daily using in-house procedure number 6-IMM-001. The weights used for testing are annually compared to Gerhart Scale Corporation's master weights and are traceable to the National Institute of Standards and Technology (NIST). The NIST Traceability numbers are 692476 - Class 1 and 692476A - Class 2. The NIST test number is 822/260017-98. All analytical balances are calibrated every 4 months by Gerhart Scale Corp. of South Amboy. The balances are calibrated with a class 1 and/or class 2 analytical weight set. These weights are tested annually by a NIST / NVLAP accredited calibration lab. The NIST test number is 822/260017-98.

4.3 THERMOMETER CALIBRATION - The thermometers used in the determination of the final densities are calibrated vs standard thermometer No. 903-2680 which was certified in accordance with the procedures outlined by ASTM E77-87 and NIST Monograph 150 using NIST Test Nos. and Std Nos.: 769543, 217368/769543, 217368/P14452, 176240/P14452, 176240. The in-house procedure No. is 2-QC-001. Thermometers which are not calibrated vs standard thermometer No. 903-2680 are traceable to NIST Identification Nos. 92564, 119016, 471047 and NIST test report Nos. 811/258522, 811/2557078, and 236090.

4.4 GLASSWARE CALIBRATION - In-house procedure 3-QC-002 is used to calibrate all Class A Glassware used in the manufacture and quality control of Custom Grade Standards.

5.0 TRACE METALLIC IMPURITIES (TMI) DETERMINED BY ICP/MS AND ICP-OES IN µg/mL

Custom-Grade solutions are tested for trace metallic impurities by Axial ICP-OES and ICP-MS. The result from the most sensitive method for each element, is reported below. Solutions tested by ICP-MS were analyzed in an ULPA-Filtered Clean Room. An ULPA-Filter is 99.9985% efficient for the removal of particles down to 0.3 µm.

<u>Q</u> Al 0.00017	<u>M</u> Dy < 0.01196	<u>Q</u> Li < 0.00003	<u>M</u> Pr < 0.00060	<u>M</u> Te < 0.05981
<u>M</u> Sb 0.00160	<u>M</u> Er < 0.00997	<u>M</u> Lu < 0.00080	<u>Q</u> Re < 0.00900	<u>M</u> Tb < 0.00060
<u>Q</u> As < 0.00500	<u>M</u> Eu < 0.00598	<u>Q</u> Mg < 0.00003	<u>M</u> Rh < 0.00199	<u>M</u> Tl < 0.00199
<u>M</u> Ba < 0.01994	<u>M</u> Gd < 0.00199	<u>M</u> Mn < 0.00798	<u>M</u> Rb < 0.00199	<u>M</u> Th < 0.00199
<u>Q</u> Be < 0.00009	<u>M</u> Ga < 0.00199	<u>Q</u> Hg 0.01950	<u>Q</u> Ru 0.00220	<u>M</u> Tm < 0.00080
<u>M</u> Bi < 0.00080	<u>M</u> Ge < 0.01196	<u>Q</u> Mo < 0.00400	<u>M</u> Sm < 0.00199	<u>M</u> Sn < 0.00997
<u>Q</u> B < 0.00006	<u>M</u> Au < 0.00598	<u>M</u> Nd < 0.00399	<u>M</u> Sc < 0.01994	<u>M</u> Ti < 0.09969
<u>M</u> Cd < 0.00598	<u>M</u> Hf < 0.00399	<u>Q</u> Ni < 0.00090	<u>S</u> Se	<u>M</u> W < 0.01994
<u>Q</u> Ca 0.00200	<u>M</u> Ho < 0.00100	<u>Q</u> Nb < 0.00400	<u>Q</u> Si 0.00055	<u>M</u> U < 0.00399
<u>M</u> Ce < 0.00997	<u>M</u> In < 0.01994	<u>n</u> Os	<u>M</u> Ag 0.00070	<u>M</u> V < 0.00399
<u>M</u> Cs < 0.00060	<u>M</u> Ir < 0.00997	<u>M</u> Pd < 0.00997	<u>Q</u> Na 0.00355	<u>M</u> Yb < 0.00199
<u>M</u> Cr < 0.00997	<u>Q</u> Fe 0.00060	<u>Q</u> P < 0.00300	<u>M</u> Sr < 0.00100	<u>M</u> Y < 0.07975
<u>M</u> Co < 0.00598	<u>M</u> La < 0.00100	<u>M</u> Pt < 0.00399	<u>Q</u> S 0.00500	<u>M</u> Zn < 0.03988
<u>M</u> Cu < 0.01196	<u>M</u> Pb < 0.00598	<u>Q</u> K 0.00070	<u>M</u> Ta < 0.01396	<u>Q</u> Zr < 0.00040

M - Checked by ICP-MS O - Checked by ICP-OES i - Spectral Interference n - Not Checked For s - Solution Standard Element

6.0 INTENDED USE

For the calibration of analytical instruments including but not limited to the following:
 ICP-MS, ICP-OES, FAAS, GFAA, XRF, and DCP
 For the validation of analytical methods
 For the preparation of "working reference samples"
 For interference studies and the determination of correction coefficients
 For detection limit and linearity studies
 For additional intended uses, contact IV Technical Staff

7.0 INSTRUCTIONS FOR THE CORRECT USE OF THIS REFERENCE MATERIAL

Storage & Handling - Keep tightly sealed when not in use. Store and use at $20 \pm 4^\circ\text{C}$. Do not pipet from container. Do not return portions removed for pipetting to container.

Atomic Weight; Valence; Coordination Number; Chemical Form in Solution - 78.96; +4; 6; H_2SeO_4

Chemical Compatibility - Soluble in HCl , HNO_3 , H_2PO_4 , H_2SO_4 and HF aqueous matrices and water. It is stable with most inorganic anions but many cationic metals form the insoluble selenites under pH neutral conditions. When fluorinated and/or under acidic conditions precipitation is typically not a problem at moderate to low concentrations.

Stability - 2-100 ppb levels - stable for months alone or mixed with other elements at equivalent levels - in 1 % HNO_3 / LDPE container. 1-10,000 ppm solutions chemically stable for years in 1-5% HNO_3 / LDPE container.

Se Containing Samples (Preparation and Solution) - Metal (Soluble in HNO_3); Oxides (Readily soluble in water); Minerals and alloys (Acid digestion with HNO_3 or HNO_3 / HF); Organic Matrices (Acid digestion with hot concentrated H_2SO_4 accompanied by the careful dropwise addition of H_2O_2 until clear)

Atomic Spectroscopic Information (ICP-OES D.L.s are given as radial/axial view):

Technique/Line	Estimated D.L.	Order	Type	Interferences (underlined indicates severe at λ boncs.)
ICP-OES 196.026 nm	0.08 / 0.006 $\mu\text{g/mL}$	1	atom	Fe
ICP-OES 203.985 nm	0.2 / 0.05 $\mu\text{g/mL}$	1	atom	<u>Sb</u> , <u>Ir</u> , <u>Cr</u> , <u>La</u>
ICP-OES 206.279 nm	0.3 / 0.16 $\mu\text{g/mL}$	1	atom	<u>Cr</u> , <u>Pt</u>
ICP-MS 82 amu	200 ppt	n/a	M'	^{13}C , ^{35}Cl

8.0 **HAZARDOUS INFORMATION** - Please refer to the enclosed Material Safety Data sheet for information regarding this CRM.

9.0 **HOMOGENEITY** - This solution was mixed according to procedure IV-MPM-004 and is guaranteed to be homogeneous.

10.0 QUALITY STANDARD DOCUMENTATION

10.1 ISO 9001:2000 Quality Management System Registration - QMI Certificate Number 010105

Recognized by:

Registrar Accreditation Board (ANSI-RAB)

Standards Council of Canada (SCC)

Dutch Council for Accreditation (RVA)

Entidad Mexicana de Acreditacion, a.c.(EMA)

Members of IQ Net International Certification Network:

Argentina (IRAM), Australia (QAS), Austria (ÖQS), Belgium (Avinter), Brazil (FCAV), Canada (QMI), Hong Kong (HKQAA), Columbia (ICONTEC), Czech Republic (CQS), Denmark (DS), Finland (SFS), France (AFAQ), Germany (DQS), Greece (ELOT), Hungary (MSZT), Ireland (NSAI), Israel (SII), Italy (CISQ), Japan (JQA), Korea (KSA-QA), Netherlands (KEMA), Norway (NCS), Poland (PCBC), Portugal (APCER), Singapore (PSB), Slovenia (SIQ), Spain (AENOR), Switzerland (SQS)

10.2 ISO/IEC 17025 - 1999 "General Requirements for the Competence of Testing and Calibration"

- Chemical Testing - Accredited A2LA Certificate Number 883.01

10.3 ISO/IEC Guide 34 - 2000 "General Requirements for the Competence of Reference Material Producers"

- Reference Materials Production - Accredited A2LA Certificate Number 883.02

A2LA Mutual Recognition Agreement Partners:

Australia (NATA), Austria (BmWA), Belgium (BELTEST) (BKO-OBE), Canada (SCC), Chinese Taipei (CNLA), Czech Republic (NAO), Denmark (DANAK), Finland (FINAS), France (COFRAC), Germany (DAR), Hong Kong (HKAS, Ireland (NAB), Italy (SIT) (SINAL), Japan (JAB) (JNLA), Republic of Korea (KOLAS), The Netherlands (RvA), New Zealand (IANZ), Norway (NA), Portugal (IPQ), Singapore (SAC-SINGLAS), Spain (ENAC), Sweden (SWEDAC), Switzerland (SAS), United Kingdom (UKAS) and United States (NVLAP) (ICBO ES)

10.4 10CFR50 Appendix B - Nuclear Regulatory Commission

- Domestic Licensing of Production and Utilization Facilities

10.5 10CFR21 - Nuclear Regulatory Commission - Reporting Defects and Non-Compliance

10.6 MIL-STD-45662A (Obsolete/Observed)



INORGANIC LABS/RADCHEM LABS 2 of 2

DATE RECEIVED: 06/20/03

DATE EXPIRED: 07/01/2004

DATE OPENED: 06/23/03

INORG: 4152 PO: F52370

11.0 DATE OF CERTIFICATION AND PERIOD OF VALIDITY

010229



11.1 IV Shelf Life - The period of time during which the concentration of the analyte(s) in a properly packaged, unopened, and unused standard stored under environmentally controlled and monitored conditions will remain within the specified uncertainty range. Shelf life is limited primarily by transpiration (loss of water from the solution) and infrequently, by chemical instability. Transpiration studies (P-SP01020) of chemically-stable solutions performed at Inorganic Ventures / IV Labs indicate a CRM shelf-life of four years for solutions packaged in 500-mL low density polyethylene bottles. When stored under special conditions that minimize transpiration and instability, the shelf life can be extended past this limit.

11.2 Expiration Date - The date after which a CRM should not be used. Routine laboratory use of a CRM increases transpiration losses and the chance of contamination which affect the integrity of the CRM and limit its useful life. Inorganic Ventures / IV Labs concurs with state and federal regulatory agencies' recommendations that solution standards be assigned a one-year expiration date.

Certification Date: November 27, 2002

Expiration Date: **EXPIRES**

01 22 00 4


12.0 NAMES AND SIGNATURES OF CERTIFYING OFFICERS

Certificate Prepared By: Debbie Newman, QA Administrator

Certificate Approved By: Katalin Le, QC Supervisor

Certifying Officer: Paul Gaines, Chemist, Senior Technical Director

Debbie Newman
Katalin Le
Paul Gaines



Certificate of Analysis



CUSTOM-GRADE SOLUTION **1000 µg/mL Thallium in 0.5% HNO₃ (abs)**

010230

Catalog Number: CGTL1-1, CGTL1-2 and CGTL1-5

Lot Number: **T-TL01075**

Starting Material:
Starting Material Purity:
Starting Material Lot No:

Thallium Nitrate
99.999%
59001

INORGANIC LABS/RADCHEM LABS

DATE RECEIVED: 01/16/03

DATE EXPIRED: 02/01/2004

DATE OPENED: 01/16/03

INORG: 3922 PO: F52073

CERTIFIED CONCENTRATION: 997 ± 1 µg/mL

The Certified Value is based upon the most precise method used to analyze this CRM. The following equations are used in the calculation of the certified value and the uncertainty:

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

$$\text{Uncertainty } (\pm) = \frac{2[(\sum s_i)^2]^{1/2}}{(n)^{1/2}}$$

($\bar{x}$) = mean

x_i = individual results

n = number of measurements

$\sum S_i$ = The summation of all significant estimated errors.

Calculated Value: 1000 µg/mL

Method: Calculated, based on starting material.

Instrument Analysis: 997 ± 1 µg/mL (Average of 2 runs)

Method: Inductively Coupled Plasma Spectroscopy (ICP) vs NIST SRM 3158.

TRACE METALLIC IMPURITIES DETERMINED BY ICP-MS AND ICP-OES IN µg/mL:

Custom-Grade solutions tested for trace metallic impurities by ICP-MS were analyzed in an ULPA-Filtered Clean Room.

An ULPA-Filter is 99.9985% efficient for the removal of particles down to 0.3 µm.

<u>Q</u> Al <0.040	<u>M</u> Dy <0.00060	<u>M</u> Li <0.0010	<u>M</u> Pr <0.000030	<u>M</u> Te <0.0030
<u>M</u> Sb 0.000080	<u>M</u> Er <0.00050	<u>M</u> Lu <0.000040	<u>M</u> Re <0.00010	<u>M</u> Tb <0.000030
<u>M</u> As <0.0010	<u>M</u> Eu <0.00030	<u>M</u> Mg <0.0030	<u>M</u> Rh <0.00010	<u>s</u> Tl
<u>M</u> Ba <0.0010	<u>M</u> Gd <0.00010	<u>M</u> Mn <0.00040	<u>M</u> Rb <0.00010	<u>M</u> Th <0.00010
<u>M</u> Be <0.000050	<u>M</u> Ga <0.00010	<u>Q</u> Hg i	<u>M</u> Ru <0.00020	<u>M</u> Tm <0.000040
<u>M</u> Bi <0.000040	<u>M</u> Ge <0.00060	<u>M</u> Mo <0.00020	<u>M</u> Sm <0.00010	<u>M</u> Sn 0.00060
<u>Q</u> B <0.020	<u>M</u> Au <0.00030	<u>M</u> Nd <0.00020	<u>M</u> Sc <0.0010	<u>M</u> Ti <0.0050
<u>M</u> Cd 0.00090	<u>M</u> Hf 0.0010	<u>Q</u> Ni <0.050	<u>Q</u> Se <0.40	<u>M</u> W <0.0010
<u>Q</u> Ca <0.010	<u>M</u> Ho <0.000050	<u>M</u> Nb <0.000050	<u>Q</u> Si <0.020	<u>M</u> U <0.00020
<u>M</u> Ce <0.00050	<u>Q</u> In i	<u>n</u> Os	<u>M</u> Ag <0.00020	<u>M</u> V <0.00020
<u>M</u> Cs 0.00010	<u>M</u> Ir <0.00050	<u>M</u> Pd <0.00050	<u>Q</u> Na <0.090	<u>M</u> Yb <0.00010
<u>M</u> Cr <0.00050	<u>Q</u> Fe <0.050	<u>Q</u> P i	<u>M</u> Sr <0.000050	<u>M</u> Y <0.0040
<u>M</u> Co <0.00030	<u>M</u> La <0.000050	<u>M</u> Pt <0.00020	<u>n</u> S	<u>M</u> Zn 0.015
<u>M</u> Cu <0.00060	<u>M</u> Pb <0.00030	<u>n</u> K	<u>M</u> Ta <0.00070	<u>M</u> Zr <0.00050

M - checked by ICP-MS

O - checked by ICP-OES

i - spectral interference

n - not checked for

s - solution standard element

ANALYZED DENSITY OF SOLUTION (measured at 22°C): 1.002 g/mL

QA:KL Rev.031402DN

(over)

Inorganic Ventures, Inc.

195 Lehigh Avenue • Suite 4 • Lakewood, NJ 08701

Orders: 800-669-6799 • FAX (732) 901-1903

Technical Support: 800-569-6799

Quality Assurance Manager





Certificate of Analysis

CUSTOM-GRADE SOLUTION

1000 µg/mL Lanthanum in 1.4% HNO₃ (abs)

010231

Catalog Number: CGLA1-1 and CGLA1-5

Lot Number: **T-QLA01057**

Starting Material:
Starting Material Purity:
Starting Material Lot No:

Lanthanum Oxide
99.999%
LA-0-5-017

INORGANIC LABS/RADCHEM LABS

DATE RECEIVED: 08/26/03

DATE EXPIRED: 09/01/2004

DATE OPENED: 08/26/03

INORG: 4221 PD: E52224

CERTIFIED CONCENTRATION: 1002 ± 3 µg/mL

The Certified Value is based upon the most precise method used to analyze this CRM. The following equations are used in the calculation of the certified value and the uncertainty:

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

$$\text{Uncertainty } (\pm) = \frac{2[(\sum s_i)^2]^{1/2}}{(n)^{1/2}}$$

($\bar{x}$) = mean

x_i = individual results

n = number of measurements

$\sum s_i$ = The summation of all significant estimated errors.

Classical Wet Assay: 1002 ± 3 µg/mL

Method: EDTA Titration vs NIST SRM 928 Lead Nitrate.

Instrument Analysis: 1007 ± 3 µg/mL

Method: Inductively Coupled Plasma Spectroscopy (ICP) vs NIST SRM 3127a.

The independent samples t-test was used to determine if there is agreement between the above assay methods at the 95% confidence interval. Both methods were compared and showed agreement within the stated uncertainties. This agreement is a confirmation of the accuracy of this CRM.

TRACE METALLIC IMPURITIES DETERMINED BY ICP-MS AND ICP-OES IN µg/mL:

Custom-Grade solutions tested for trace metallic impurities by ICP-MS were analyzed in an ULPA-Filtered Clean Room.

An ULPA-Filter is 99.9985% efficient for the removal of particles down to 0.3 µm.

<u>Q</u> Al <0.040	<u>M</u> Dy <0.00060	<u>M</u> Li <0.0010	<u>Q</u> Pr <0.020	<u>M</u> Te <0.0030
<u>M</u> Sb <0.000050	<u>M</u> Er 0.0010	<u>M</u> Lu 0.000040	<u>M</u> Re <0.00010	<u>M</u> Tb <0.000030
<u>M</u> As <0.0010	<u>M</u> Eu <0.00030	<u>M</u> Mg <0.0030	<u>M</u> Rh <0.00010	<u>M</u> Tl <0.00010
<u>Q</u> Ba <0.020	<u>M</u> Gd 0.039	<u>M</u> Mn <0.00040	<u>M</u> Rb <0.00010	<u>M</u> Th <0.00010
<u>M</u> Be <0.000050	<u>M</u> Ga <0.00010	<u>Q</u> Hg <0.030	<u>M</u> Ru <0.00020	<u>M</u> Tm <0.000040
<u>M</u> Bi <0.000040	<u>M</u> Ge <0.00060	<u>M</u> Mo <0.00020	<u>M</u> Sm 0.00040	<u>M</u> Sn <0.00050
<u>Q</u> B <0.020	<u>M</u> Au <0.00030	<u>M</u> Nd 0.00020	<u>M</u> Sc <0.0010	<u>M</u> Ti <0.0050
<u>M</u> Cd <0.00030	<u>M</u> Hf <0.00020	<u>Q</u> Ni <0.050	<u>Q</u> Se <0.40	<u>M</u> W <0.0010
<u>Q</u> Ca <0.010	<u>M</u> Ho 0.00010	<u>M</u> Nb <0.000050	<u>Q</u> Si <0.020	<u>M</u> U <0.00020
<u>i</u> Ce	<u>Q</u> In <0.030	<u>n</u> Os	<u>M</u> Ag <0.00020	<u>M</u> V <0.00020
<u>n</u> Cs	<u>M</u> Ir <0.00050	<u>M</u> Pd <0.00050	<u>Q</u> Na <0.090	<u>M</u> Yb <0.00010
<u>M</u> Cr <0.00050	<u>Q</u> Fe <0.050	<u>Q</u> P <0.050	<u>M</u> Sr <0.000050	<u>M</u> Y <0.0040
<u>M</u> Co <0.00030	<u>s</u> La	<u>M</u> Pt <0.00020	<u>n</u> S	<u>M</u> Zn <0.0020
<u>M</u> Cu <0.00060	<u>M</u> Pb <0.00030	<u>n</u> K	<u>M</u> Ta <0.00070	<u>M</u> Zr <0.00050

M - checked by ICP-MS

O - checked by ICP-OES

i - spectral interference

n - not checked for

s - solution standard element

ANALYZED DENSITY OF SOLUTION (measured at 22°C): 1.009 g/mL

(over)

QA:KSL Rev.121702DM

Inorganic Ventures, Inc.

195 Lehigh Avenue • Suite 4 • Lakewood, NJ 08701

Orders: 800-669-6799 • FAX (732) 901-1903

Technical Support: 800-569-6799

Quality Assurance Manager

EXPIRES

1/2004



Certificate of Analysis

010232

CUSTOM-GRADE SOLUTION

1000 µg/mL Palladium in 3.3% HCl (abs)

Catalog Number: CGPD1-1 and CGPD1-5

Lot Number: T-PD02014

Starting Material:
Starting Material Purity:
Starting Material Lot No:

Palladium Nitrate
99.999%
00614

INORGANIC LABS/RADCHEM LABS

DATE RECEIVED: 2/28/03

DATE EXPIRED: 3/01/2004

DATE OPENED: 2/28/03

INORG: 3964 PO: F52064

CERTIFIED CONCENTRATION: 1001 ± 3 µg/mL

The Certified Value is based upon the most precise method used to analyze this CRM. The following equations are used in the calculation of the certified value and the uncertainty:

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

$$\text{Uncertainty } (\pm) = \frac{2(\sum s_i^2)^{1/2}}{(n)^{1/2}}$$

($\bar{x}$) = mean

x_i = individual results

n = number of measurements

$\sum s_i$ = The summation of all significant estimated errors.

Calculated Value: 1002 µg/mL

Method: Calculated, based on starting material.

Instrument Analysis: 1001 ± 3 µg/mL

Method: Inductively Coupled Plasma Spectroscopy (ICP) vs NIST SRM 3138.

TRACE METALLIC IMPURITIES DETERMINED BY ICP-MS AND ICP-OES IN µg/mL:

Custom-Grade solutions tested for trace metallic impurities by ICP-MS were analyzed in an ULPA-Filtered Clean Room.
An ULPA-Filter is 99.9985% efficient for the removal of particles down to 0.3 µm.

M	Al	<0.00090	M	Dy	<0.00060	M	Li	<0.0010	M	Pr	<0.000030	M	Te	<0.0030
M	Sb	0.0017	M	Er	<0.00050	M	Lu	<0.000040	M	Re	<0.00010	M	Tb	<0.000030
M	As	<0.060	M	Eu	<0.00030	M	Mg	0.0090	O	Rh	<0.020	M	Ti	0.00010
M	Ba	<0.0010	M	Gd	<0.00010	M	Mn	<0.00040	M	Rb	<0.00010	M	Th	<0.00010
M	Be	<0.000050	M	Ga	0.00026	O	Hg	<0.030	O	Ru	<0.010	M	Tm	<0.000040
M	Bi	<0.000040	M	Ge	0.0058	M	Mo	<0.00020	M	Sm	<0.00010	M	Sn	<0.00050
M	B	<0.0070	M	Au	0.00090	M	Nd	0.00090	M	Sc	<0.0010	M	Tl	<0.0050
O	Cd	<0.010	M	Hf	<0.00020	O	Ni	<0.050	O	Se	<0.40	M	W	<0.0010
O	Ca	0.040	M	Ho	<0.000050	M	Nb	<0.000050	O	Si	<0.020	M	U	<0.00020
M	Ce	<0.00050	O	In	<0.030	n	Os		O	Ag	0.26	O	V	<0.010
i	Cs		M	Ir	<0.00050	s	Pd		O	Na	<0.090	M	Yb	<0.00010
O	Cr	<0.020	O	Fe	<0.050	O	P	<0.050	M	Sr	<0.000050	M	Y	<0.0040
M	Co	<0.00030	M	La	<0.000050	M	Pt	0.016	n	S		M	Zn	0.0030
M	Cu	0.0023	M	Pb	0.0022	n	K		M	Ta	<0.00070	M	Zr	<0.00050

M - checked by ICP-MS

O - checked by ICP-OES

i - spectral interference

n - not checked for

s - solution standard element

ANALYZED DENSITY OF SOLUTION (measured at 22°C): 1.022 g/mL

QA:KSL Rev.050202DM

(over)

Inorganic Ventures, Inc.

195 Lehigh Avenue • Suite 4 • Lakewood, NJ 08701

Orders: 800-669-6799 • FAX (732) 901-1903

Technical Support: 800-569-6799

Paul R. Gaines

Quality Assurance Manager

EXPIRES

1/2004



inorganic ventures / iv labs

195 lehigh avenue, suite 4, lakewood, nj 08701 usa
phone: 800-669-6799 • 732-901-1900 • fax: 732-901-1903
e-mail: ivsales@ivstandards.com • website: www.ivstandards.com

certificate of analysis ⁰¹⁰²³³

1.0 Inorganic Ventures / IV Labs is an ISO Guide 34-2000 Certified Reference Material (CRM) Manufacturer: Certificate #883-02. The certificate is designed and the certified value(s) and uncertainty(ies) are determined in accordance with ISO Guide 31-2000 (Reference Materials - Contents of certificates and label(s), ISO Guide 34-2000 "Quality System Guidelines for the Production of Reference Materials," and ISO Guide 35-1989 "Certification of Reference Materials - General and Statistical Principles."

2.0 DESCRIPTION OF CRM Custom-Grade 1000 µg/mL Sulfur in H₂O

Catalog Number: CGS1-1 and CGS1-5
Lot Number: **W-QS01098**
Starting Material: H₂SO₄
Starting Material Purity (%): 99.999965
Starting Material Lot No: N38818
Matrix: H₂O

INORGANIC LABS/RADCHEM LABS
DATE RECEIVED: 11/5/03 1062
DATE EXPIRED: 12/1/2004 DR
DATE OPENED: 11/5/03
INORG: 4317 PO: F52258

3.0 CERTIFIED VALUES AND UNCERTAINTIES

Certified Concentration: 1010 ± 2 µg/mL

Certified Density: 1.000 g/mL (measured at 22° C)

The Certified Value is based upon the most precise method used to analyze this CRM. The following equations are used in the calculation of the certified value and the uncertainty:

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

($\bar{x}$) = mean

x_i = individual results

n = number of measurements

$$\text{Uncertainty } (\pm) = \frac{2[(\sum s_i^2)]^{1/2}}{(n)^{1/2}}$$

$\sum s_i^2$ = The summation of all significant estimated errors

(Most common are the errors from instrumental measurement, weighing, dilution to volume, and the fixed error reported on the

NIST SRM certificate of analysis.)

The independent samples t-test was used to determine if there is agreement between the above assay methods at the 95% confidence interval. Both methods were compared and showed agreement within the stated uncertainties. This agreement is a confirmation of the accuracy of this CRM.

4.0 TRACEABILITY TO NIST AND VALUES OBTAINED BY INDEPENDENT METHODS

□ "Property of the result of a measurement or the value of a standard whereby it can be related to stated references, usually national or international standards, through an unbroken chain of comparisons all having stated uncertainties." (ISO VIM, 2nd ed., 1993, definition 6.10)

□ This IV product is Traceable to NIST via direct comparison to NIST SRMs. The uncertainties for each certified value are reported, taking into account the SRM uncertainty error and the measurement, weighing and volume dilution errors.

4.1 Assay Method #1 998 ± 8 µg/mL

ICP Assay NIST SRM 3154 Lot Number: 892205

Assay Method #2 1010 ± 2 µg/mL

Acidimetric NIST SRM 84k Lot Number: 84k

4.2 BALANCE CALIBRATION - All balances are checked daily using in-house procedure number 6-IMM-001. The weights used for testing are annually compared to Gerhart Scale Corporation's master weights and are traceable to the National Institute of Standards and Technology (NIST). The NIST Traceability numbers are 692476 - Class 1 and 692476A - Class 2. The NIST test number is 822/260017-98. All analytical balances are calibrated every 4 months by Gerhart Scale Corp. of South Amboy. The balances are calibrated with a class 1 and/or class 2 analytical weight set. These weights are tested annually by a NIST / NVLAP accredited calibration lab. The NIST test number is 822/260017-98.

010234

4.3 THERMOMETER CALIBRATION - The thermometers used in the determination of the final densities are calibrated vs standard thermometer No. 903-2680 which was certified in accordance with the procedures outlined by ASTM E77-87 and NIST Monograph 150 using NIST Test Nos. and Std Nos.: 769543, 217368/769543, 217368/P14452, 176240/P14452, 176240. The in-house procedure No. is 2-QC-001. Thermometers which are not calibrated vs standard thermometer No. 903-2680 are traceable to NIST Identification Nos. 92564, 119016, 471047 and NIST test report Nos. 811/258522, 811/2557078, and 236090.

4.4 GLASSWARE CALIBRATION - In-house procedure 3-QC-002 is used to calibrate all Class A Glassware used in the manufacture and quality control of Custom Grade Standards.

5.0 TRACE METALLIC IMPURITIES (TMI) DETERMINED BY ICP/MS AND ICP-OES IN µg/mL

Custom-Grade solutions are tested for trace metallic impurities by Axial ICP-OES and ICP-MS. The result from the most sensitive method for each element, is reported below. Solutions tested by ICP-MS were analyzed in an ULPA-Filtered Clean Room. An ULPA-Filter is 99.9985% efficient for the removal of particles down to 0.3 µm.

<u>Q</u> Al 0.00025	<u>M</u> Dy < 0.01197	<u>Q</u> Li < 0.00016	<u>M</u> Pr < 0.00060	<u>M</u> Te < 0.05984
<u>M</u> Sb < 0.00100	<u>M</u> Er < 0.00997	<u>M</u> Lu < 0.00080	<u>M</u> Re < 0.00200	<u>M</u> Tb < 0.00060
<u>M</u> As < 0.01995	<u>M</u> Eu < 0.00598	<u>Q</u> Mg < 0.00004	<u>M</u> Rh < 0.00200	<u>M</u> Tl < 0.00200
<u>M</u> Ba < 0.01995	<u>M</u> Gd < 0.00200	<u>M</u> Mn < 0.00798	<u>M</u> Rb < 0.00200	<u>M</u> Th < 0.00200
<u>Q</u> Be < 0.00200	<u>M</u> Ga < 0.00200	<u>Q</u> Hg < 0.01100	<u>M</u> Ru < 0.00399	<u>M</u> Tm < 0.00080
<u>M</u> Bi < 0.00080	<u>M</u> Ge < 0.01197	<u>M</u> Mo < 0.00399	<u>M</u> Sm < 0.00200	<u>M</u> Sn < 0.00997
<u>Q</u> B < 0.00990	<u>M</u> Au < 0.00598	<u>M</u> Nd < 0.00399	<u>M</u> Sc < 0.01995	<u>M</u> Ti < 0.09974
<u>M</u> Cd < 0.00598	<u>M</u> Hf < 0.00399	<u>Q</u> Ni < 0.00230	<u>Q</u> Se < 0.00620	<u>M</u> W < 0.01995
<u>Q</u> Ca 0.00020	<u>M</u> Ho < 0.00100	<u>M</u> Nb < 0.00100	<u>Q</u> Si < 0.00410	<u>M</u> U < 0.00399
<u>M</u> Ce < 0.00997	<u>M</u> In < 0.01995	<u>n</u> Os	<u>M</u> Ag < 0.00399	<u>M</u> V < 0.00399
<u>M</u> Cs < 0.00060	<u>M</u> Ir < 0.00997	<u>M</u> Pd < 0.00997	<u>Q</u> Na < 0.00010	<u>M</u> Yb < 0.00200
<u>M</u> Cr < 0.00997	<u>Q</u> Fe 0.00015	<u>Q</u> P < 0.00480	<u>M</u> Sr < 0.00100	<u>M</u> Y < 0.07979
<u>M</u> Co < 0.00598	<u>M</u> La < 0.00100	<u>M</u> Pt < 0.00399	<u>s</u> S	<u>Q</u> Zn 0.00125
<u>M</u> Cu < 0.01197	<u>M</u> Pb < 0.00598	<u>Q</u> K < 0.00170	<u>M</u> Ta < 0.01396	<u>M</u> Zr < 0.00997

M - Checked by ICP-MS O - Checked by ICP-OES i - Spectral Interference n - Not Checked For s - Solution Standard Element

6.0 INTENDED USE

For the calibration of analytical instruments including but not limited to the following:

ICP-MS, ICP-OES, FAAS, GFAA, XRF, and DCP

For the validation of analytical methods

For the preparation of "working reference samples"

For interference studies and the determination of correction coefficients

For detection limit and linearity studies

For additional intended uses, contact IV Technical Staff

7.0 INSTRUCTIONS FOR THE CORRECT USE OF THIS REFERENCE MATERIAL

Storage & Handling - Keep tightly sealed when not in use. Store and use at $20 \pm 4^\circ\text{C}$. Do not pipet from container. Do not return portions removed for pipetting to container.

Atomic Weight; Valence; Coordination Number; Chemical Form in Solution - 32.066; +6; 6 (O=), S(OH),

010235

Chemical Compatibility - Soluble in HCl, HNO₃, H₃PO₄, and HF aqueous matrices water and NH₄OH. Stable with all metals and inorganic anions at low to moderate ppm levels under acidic conditions except Ba and Pb and to a lesser extent Sr, and Ca.

Stability - 2-100 ppb levels - stability unknown - in 1% HNO₃ / LDPE container. 1-10,000 ppm solutions chemically stable for years in LDPE container.

S Containing Samples (Preparation and Solution) - We most often get questions about the determination of S in Rocks, Silicates and insoluble sulfates (the finely powdered sample is fused in a Pt* crucible with 20 times its weight of Na₂CO₃ + 0.5 grams KNO₃. The fuseate is extracted with water. Any BaSO₄ present in the sample is transposed by the carbonate fusion to the BaCO₃ which is left behind in the water-insoluble residue. If PbSO₄ is present the fuseate should be boiled with a sodium carbonate saturated with CO₂ solution for 1 hour or more where the PbSO₄ will be transposed to the water insoluble carbonate which can be filtered off. Boiling the fuseate with a saturated carbonate solution is good insurance for samples containing Ba, Sr, and Ca. The Ba, Pb, Sr, Ca, free filtrate can be acidified and measured by ICP.)

Atomic Spectroscopic Information (ICP-OES D.L.s are given as radial/axial view):

Technique/Line	Estimated D.L.	Order	Type	Interferences (underlined indicates severe at 4000s.)
ICP-OES 166.669nm	0.2 / 0.19 µg/mL	1	atom	Si, B
ICP-OES 182.034 nm	0.3 / 0.024 µg/mL	1	atom	
ICP-OES 143.328 nm	0.4 / 0.035 µg/mL	1	atom	
ICP-MS 32 amu	30,000 ppt	n/a	M	¹⁸ O ₂ , ¹⁴ N ¹⁶ O, ¹⁴ N ¹⁷ O, ¹⁴ N ¹⁸ O ² H, ¹⁴ N ¹⁶ O ² H

8.0 HAZARDOUS INFORMATION - Please refer to the enclosed Material Safety Data sheet for information regarding this CRM.

9.0 HOMOGENEITY - This solution was mixed according to procedure IV-MPM-004 and is guaranteed to be homogeneous.

10.0 QUALITY STANDARD DOCUMENTATION

10.1 ISO 9001:2000 Quality Management System Registration - QMI Certificate Number 010105

Recognized by:

Registrar Accreditation Board (ANSI-RAB)

Standards Council of Canada (SCC)

Dutch Council for Accreditation (RVA)

Entidad Mexicana de Acreditacion, a.c.(EMA)

Members of IQ Net International Certification Network:

Argentina (IRAM), Australia (QAS), Austria (ÖQS), Belgium (Avinter), Brazil (FCAV), Canada (QMI), Hong Kong (HKQAA), Columbia (ICONTEC), Czech Republic (CQS), Denmark (DS), Finland (SFS), France (AFAQ), Germany (DQS), Greece (ELOT), Hungary (MSZT), Ireland (NSAI), Israel (SII), Italy (CISQ), Japan (JQA), Korea (KSA-QA), Netherlands (KEMA), Norway (NCS), Poland(PCBC), Portugal (APCER), Singapore (PSB), Slovenia (SIQ), Spain (AENOR), Switzerland (SQS)

10.2 ISO/IEC 17025 - 1999 "General Requirements for the Competence of Testing and Calibration"

- Chemical Testing - Accredited A2LA Certificate Number 883.01



10.3 ISO/IEC Guide 34 - 2000 "General Requirements for the Competence of Reference Material Producers"

- Reference Materials Production - Accredited A2LA Certificate Number 883.02

A2LA Mutual Recognition Agreement Partners:

Australia (NATA), Austria (BmWA), Belgium (BELTEST) (BKO-OBE), Canada (SCC), Chinese Taipei (CNLA), Czech Republic (NAO), Denmark (DANAK), Finland (FINAS), France (COFRAC), Germany (DAR), Hong Kong (HKAS), Ireland (NAB), Italy (SIT) (SINAL), Japan (JAB) (JNLA), Republic of Korea (KOLAS), The Netherlands (RvA), New Zealand (IANZ), Norway (NA), Portugal (IPQ), Singapore (SAC-SINGLAS), Spain (ENAC), Sweden (SWEDAC), Switzerland (SAS), United Kingdom (UKAS) and United States (NVLAP) (ICBO ES)

10.4 10CFR50 Appendix B - Nuclear Regulatory Commission

- Domestic Licensing of Production and Utilization Facilities

10.5 10CFR21 - Nuclear Regulatory Commission - Reporting Defects and Non-Compliance

10.6 MIL-STD-45662A (Obsolete/Observed)

11.0 DATE OF CERTIFICATION AND PERIOD OF VALIDITY



11.1 IV Shelf Life - The period of time during which the concentration of the analyte(s) in a properly packaged, unused standard stored under environmentally controlled and monitored conditions will remain within the specified uncertainty range. Shelf life is limited primarily by transpiration (loss of water from the solution) and infrequently, by chemical instability. Transpiration studies (P-SP01020) of chemically-stable solutions performed at Inorganic Ventures / IV Labs indicate a CRM shelf-life of four years for solutions packaged in 500-mL low density polyethylene bottles. When stored under special conditions that minimize transpiration and instability, the shelf life can be extended past this limit.

11.2 Expiration Date - The date after which a CRM should not be used. Routine laboratory use of a CRM increases transpiration losses and the chance of contamination which affect the integrity of the CRM and limit its useful life. Inorganic Ventures / IV Labs concurs with state and federal regulatory agencies' recommendations that solution standards be assigned a one-year expiration date.

INORGANIC LABS/RADCHEM LABS
DATE RECEIVED: 11/5/03
DATE EXPIRED: 12/1/2004
DATE OPENED: 11/5/03
INORG: 4317 PO: F50058

Certification Date: August 27, 2003

Expiration Date: **EXPIRES**
12/2004

12.0 NAMES AND SIGNATURES OF CERTIFYING OFFICERS

Certificate Prepared By: JoAnn Struthers, QA Administrative Assistant

Certificate Approved By: Katalin Le, QC Supervisor

Certifying Officer: Paul Gaines, Chemist, Senior Technical Director

John Struthers
Katalin Le
Paul Gaines



Certificate of Analysis

010237

CUSTOM-GRADE SOLUTION

Catalog Number: CGTH1-1 and CGTH1-5

1000 µg/mL Thorium in 3% HNO₃ (abs)

Lot Number: T-TH01059

Starting Material:
Starting Material Purity:
Starting Material Lot No:

Thorium Nitrate
99.999%
C01L32

INORGANIC LABS/RADCHEM LABS

DATE RECEIVED: 10/08/03
DATE EXPIRED: 11/01/2004
DATE OPENED: 10/08/03
INORG: 4283 PO: F53240

CERTIFIED CONCENTRATION: 1001 ± 3 µg/mL

The Certified Value is based upon the most precise method used to analyze this CRM. The following equations are used in the calculation of the certified value and the uncertainty:

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

$$\text{Uncertainty } (\pm) = \frac{2(\sum s_i)^{2/3}}{(n)^{1/2}}$$

($\bar{x}$) = mean

x_i = individual results

n = number of measurements

$\sum s_i$ = The summation of all significant estimated errors.

Classical Wet Assay: 1001 ± 3 µg/mL

Method: EDTA Titration vs NIST SRM Lead Nitrate.

Instrument Analysis: 1002 ± 4 µg/mL

Method: Inductively Coupled Plasma Spectroscopy (ICP) vs NIST SRM 3159.

The independent samples t-test was used to determine if there is agreement between the above assay methods at the 95% confidence interval. Both methods were compared and showed agreement within the stated uncertainties. This agreement is a confirmation of the accuracy of this CRM.

TRACE METALLIC IMPURITIES DETERMINED BY ICP-MS AND ICP-OES IN µg/mL:

Custom-Grade solutions tested for trace metallic impurities by ICP-MS were analyzed in an ULPA-Filtered Clean Room.

An ULPA-Filter is 99.9985% efficient for the removal of particles down to 0.3 µm.

<u>Q</u> Al <0.00090	<u>M</u> Dy 0.0062	<u>Q</u> Li <0.000030	<u>M</u> Pr 0.00037	<u>Q</u> Te <0.031
<u>M</u> Sb <0.000050	<u>M</u> Er <0.00050	<u>M</u> Lu <0.000040	<u>M</u> Re <0.00010	<u>M</u> Tb <0.000030
<u>Q</u> As <0.014	<u>M</u> Eu <0.00030	<u>Q</u> Mg <0.000060	<u>M</u> Rh <0.00010	<u>M</u> Tl <0.00010
<u>M</u> Ba 0.0050	<u>M</u> Gd 0.0054	<u>Q</u> Mn <0.0000030	<u>M</u> Rb <0.00010	<u>s</u> Th
<u>Q</u> Be <0.00020	<u>M</u> Ga <0.00010	<u>i</u> Hg	<u>M</u> Ru <0.00020	<u>M</u> Tm <0.000040
<u>M</u> Bi <0.000040	<u>M</u> Ge <0.00060	<u>M</u> Mo <0.00020	<u>M</u> Sm 0.0095	<u>M</u> Sn <0.00050
<u>Q</u> B <0.00060	<u>M</u> Au <0.00030	<u>M</u> Nd 0.0026	<u>M</u> Sc <0.0010	<u>Q</u> Ti <0.00092
<u>Q</u> Cd <0.0045	<u>M</u> Hf <0.00020	<u>Q</u> Ni <0.0023	<u>M</u> Se <0.010	<u>M</u> W <0.0010
<u>Q</u> Ca <0.030	<u>M</u> Ho 0.00022	<u>M</u> Nb <0.000050	<u>Q</u> Si <0.0034	<u>M</u> U 0.074
<u>M</u> Ce <0.00050	<u>Q</u> In <0.0020	<u>n</u> Os	<u>M</u> Ag <0.00020	<u>M</u> V <0.00020
<u>M</u> Cs <0.000030	<u>M</u> Ir <0.00050	<u>M</u> Pd <0.00050	<u>Q</u> Na <0.00010	<u>M</u> Yb <0.00010
<u>Q</u> Cr <0.00080	<u>Q</u> Fe <0.0011	<u>i</u> P	<u>M</u> Sr <0.000050	<u>M</u> Y <0.0040
<u>M</u> Co <0.00030	<u>M</u> La <0.000050	<u>M</u> Pt <0.00020	<u>Q</u> S <0.072	<u>Q</u> Zn <0.00058
<u>M</u> Cu <0.00060	<u>M</u> Pb <0.00030	<u>Q</u> K <0.0017	<u>M</u> Ta <0.00070	<u>M</u> Zr 0.0085

M - checked by ICP-MS

O - checked by ICP-OES

i - spectral interference

n - not checked for

s - solution standard element

ANALYZED DENSITY OF SOLUTION (measured at 22°C): 1.022 g/mL

QA:KL Rev.050802DN

(over)

Inorganic Ventures, Inc.

195 Lehigh Avenue • Suite 4 • Lakewood, NJ 08701

Orders: 800-669-6799 • FAX (732) 901-1903

Technical Support: 800-569-6799

Paul R. Gaines

Quality Assurance Manager

EXPIRES

01/01/2004



Certificate of Analysis

CUSTOM-GRADE SOLUTION

1000 µg/mL Uranium in 1% HNO₃ (abs)

Catalog Number: CGU1-1 and CGU1-5

010238

Lot Number: T-U01056

Starting Material:
Starting Material Purity:
Starting Material Lot No:

Uranium Oxynitrate*
99.95%
E14L40

INORGANIC LABS/RADCHEM LABS
DATE RECEIVED: 3/03/03
DATE EXPIRED: 3/01/2004
DATE OPENED: 3/04/03
INORG: 3974 PO: F52064

CERTIFIED CONCENTRATION: 1001 ± 1 µg/mL

The Certified Value is based upon the most precise method used to analyze this CRM. The following equations are used in the calculation of the certified value and the uncertainty:

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

$$\text{Uncertainty } (\pm) = \frac{2(\sum S_i)^2)^{1/2}}{(n)^{1/2}}$$

($\bar{x}$) = mean

x_i = individual results

n = number of measurements

$\sum S_i$ = The summation of all significant estimated errors.

Instrument Analysis: 1001 ± 1 µg/mL (Avg. 2 runs)

Method: Inductively Coupled Plasma Spectroscopy (ICP) vs NIST SRM 3164.

Calculated Value: 1000 µg/mL

Method: Calculated, based on starting material.

TRACE METALLIC IMPURITIES DETERMINED BY ICP-MS AND ICP-OES IN µg/mL:

Custom-Grade solutions tested for trace metallic impurities by ICP-MS were analyzed in an ULPA-Filtered Clean Room.
An ULPA-Filter is 99.9985% efficient for the removal of particles down to 0.3 µm.

<u>O</u> Al <0.20	<u>M</u> Dy <0.0060	<u>O</u> Li <0.010	<u>M</u> Pr <0.00030	<u>M</u> Te <0.030
<u>M</u> Sb 0.00020	<u>M</u> Er <0.0050	<u>M</u> Lu <0.00040	<u>M</u> Re <0.0010	<u>M</u> Tb <0.00030
<u>M</u> As <0.010	<u>M</u> Eu <0.0030	<u>O</u> Mg 0.012	<u>M</u> Rh <0.0010	<u>M</u> Tl <0.0010
<u>M</u> Ba 0.0022	<u>M</u> Gd <0.0010	<u>M</u> Mn <0.0040	<u>M</u> Rb <0.0010	<u>M</u> Th <0.0010
<u>O</u> Be <0.10	<u>M</u> Ga <0.0010	<u>i</u> Hg	<u>M</u> Ru <0.0020	<u>M</u> Tm <0.00040
<u>M</u> Bi 0.00030	<u>M</u> Ge <0.0060	<u>M</u> Mo <0.0020	<u>M</u> Sm <0.0010	<u>M</u> Sn 0.043
<u>O</u> B 0.0020	<u>M</u> Au <0.0030	<u>M</u> Nd <0.0020	<u>M</u> Sc 0.00080	<u>M</u> Ti 0.029
<u>M</u> Cd <0.0030	<u>M</u> Hf <0.0020	<u>i</u> Ni	<u>M</u> Se <0.0080	<u>M</u> W <0.010
<u>O</u> Ca 0.11	<u>M</u> Ho <0.00050	<u>M</u> Nb <0.00050	<u>i</u> Si	<u>s</u> U
<u>M</u> Ce <0.0050	<u>M</u> In <0.010	<u>n</u> Os	<u>M</u> Ag <0.0020	<u>M</u> V 0.00040
<u>M</u> Cs <0.00030	<u>M</u> Ir <0.0050	<u>M</u> Pd <0.0050	<u>O</u> Na 0.19	<u>M</u> Yb <0.0010
<u>M</u> Cr 0.0047	<u>i</u> Fe	<u>i</u> P	<u>M</u> Sr 0.0011	<u>M</u> Y 0.56
<u>M</u> Co 0.014	<u>M</u> La 0.0014	<u>M</u> Pt <0.0020	<u>i</u> S	<u>M</u> Zn 0.020
<u>M</u> Cu 0.022	<u>M</u> Pb 0.0025	<u>i</u> K	<u>M</u> Ta <0.0070	<u>M</u> Zr 0.00060

M - checked by ICP-MS

O - checked by ICP-OES

i - spectral interference

n - not checked for

s - solution standard element

****NOTICE TO ICP-MS USERS:** The ²³⁵U in this standard is depleted. The certified abundances in Atom % are as follows:

Isotope	Natural Abundance		IV's Certified Abundance	
	Atom %		Atom %	
Uranium	²³⁸ U	99.3	99.8 ± 0.1	
	²³⁵ U	0.70	0.195 ± 0.005	

ANALYZED DENSITY OF SOLUTION (measured at 22°C): 1.008 g/mL

(over) QA:KSL Rev.092502DN

Inorganic Ventures, Inc.

195 Lehigh Avenue • Suite 4 • Lakewood, NJ 08701

Orders: 800-669-6799 • FAX (732) 901-1903

Technical Support: 800-569-6799

Quality Assurance Manager

EXPIRES
1/2004



inorganic ventures / iv labs

195 lehigh avenue, suite 4, lakewood, nj 08701 usa
phone: 800-669-6799 • 732-901-1900 • fax: 732-901-1903
e-mail: ivsales@ivstandards.com • website: www.ivstandards.com

010239

certificate of analysis

1.0 Inorganic Ventures / IV Labs is an ISO Guide 34-2000 Certified Reference Material (CRM) Manufacturer: Certificate #883-02. The certificate is designed and the certified value(s) and uncertainty(ies) are determined in accordance with ISO Guide 31-2000 (Reference Materials - Contents of certificates and label(s)), ISO Guide 34-2000 "Quality System Guidelines for the Production of Reference Materials," and ISO Guide 35-1989 "Certification of Reference Materials - General and Statistical Principles."

2.0 DESCRIPTION OF CRM Custom-Grade 1000 µg/mL Tungsten in 1% (abs) HNO₃/1% (abs) HF

Catalog Number: CGW1-1 and CGW1-5

Lot Number: W-W01080

Starting Material: W Powder

Starting Material Purity (%): 99.990703

Starting Material Lot No 21418,C31H46,D02J21,E03K06,D11F29

Matrix: 1% (abs) HNO₃/1% (abs) HF

INORGANIC LABS/RADCHEM LABS Pg. 1 of 2

DATE RECEIVED: 07/31/03

DATE EXPIRED: 08/01/2004

DATE OPENED: 08/01/03

INORG: 4203 PO: E52383

3.0 CERTIFIED VALUES AND UNCERTAINTIES

Certified Concentration: 1001 ± 2 µg/mL

Certified Density: 1.006 g/mL (measured at 22° C)

The Certified Value is based upon the most precise method used to analyze this CRM. The following equations are used in the calculation of the certified value and the uncertainty:

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

($\bar{x}$) = mean

x_i = individual results

n = number of measurements

$$\text{Uncertainty } (\pm) = \frac{2(\sum s_i^2)^{1/2}}{(n)^{1/2}}$$

$\sum S$ = The summation of all significant estimated errors.

(Most common are the errors from instrumental measurement, weighing, dilution to volume, and the fixed error reported on the NIST SRM certificate of analysis.)

4.0 TRACEABILITY TO NIST AND VALUES OBTAINED BY INDEPENDENT METHODS

□ "Property of the result of a measurement or the value of a standard whereby it can be related to stated references, usually national or international standards, through an unbroken chain of comparisons all having stated uncertainties." (ISO VIM, 2nd ed., 1993, definition 6.10)

□ This IV product is Traceable to NIST via direct comparison to NIST SRMs. The uncertainties for each certified value are reported, taking into account the SRM uncertainty error and the measurement, weighing and volume dilution errors.

4.1 Assay Method #1 1001 ± 2 µg/mL (Avg 2 runs)

ICP Assay NIST SRM 3163 Lot Number: 990209

Assay Method #2 1000 µg/mL

Gravimetric NIST SRM Lot Number: See Sec. 4.2

4.2 BALANCE CALIBRATION - All balances are checked daily using in-house procedure number 6-IMM-001. The weights used for testing are annually compared to Gerhart Scale Corporation's master weights and are traceable to the National Institute of Standards and Technology (NIST). The NIST Traceability numbers are 692476 - Class 1 and 692476A - Class 2. The NIST test number is 822/260017-98. All analytical balances are calibrated every 4 months by Gerhart Scale Corp. of South Amboy. The balances are calibrated with a class 1 and/or class 2 analytical weight set. These weights are tested annually by a NIST / NVLAP accredited calibration lab. The NIST test number is 822/260017-98.

4.3 THERMOMETER CALIBRATION - The thermometers used in the determination of the final densities are calibrated vs standard thermometer No. 903-2680 which was certified in accordance with the procedures outlined by ASTM E77-87 and NIST Monograph 150 using NIST Test Nos. and Std Nos.: 769543, 217368/769543, 217368/P14452, 176240/P14452, 176240. The in-house procedure No. is 2-QC-001. Thermometers which are not calibrated vs standard thermometer No. 903-2680 are traceable to NIST Identification Nos. 92564, 119016, 471047 and NIST test report Nos. 811/258522, 811/2557078, and 236090.

4.4 GLASSWARE CALIBRATION - In-house procedure 3-QC-002 is used to calibrate all Class A Glassware used in the manufacture and quality control of Custom Grade Standards.

5.0 TRACE METALLIC IMPURITIES (TMI) DETERMINED BY ICP/MS AND ICP-OES IN µg/mL

Custom-Grade solutions are tested for trace metallic impurities by Axial ICP-OES and ICP-MS. The result from the most sensitive method for each element, is reported below. Solutions tested by ICP-MS were analyzed in an ULPA-Filtered Clean Room. An ULPA-Filter is 99.9985% efficient for the removal of particles down to 0.3 µm.

<u>Q</u> Al 0.01792	<u>M</u> Dy < 0.00595	<u>Q</u> Li < 0.00008	<u>M</u> Pr < 0.00030	<u>M</u> Te < 0.02974
<u>M</u> Sb < 0.00050	<u>M</u> Er < 0.00496	<u>M</u> Lu < 0.00040	<u>I</u> Re	<u>M</u> Tb < 0.00030
<u>M</u> As < 0.00991	<u>M</u> Eu < 0.00297	<u>Q</u> Mg 0.00120	<u>M</u> Rh < 0.00099	<u>M</u> Tl < 0.00099
<u>M</u> Ba < 0.00991	<u>M</u> Gd < 0.00099	<u>M</u> Mn < 0.00397	<u>M</u> Rb < 0.00099	<u>M</u> Th < 0.00099
<u>M</u> Be < 0.00050	<u>M</u> Ga < 0.00099	<u>Q</u> Hg < 0.04778	<u>M</u> Ru < 0.00198	<u>M</u> Tm < 0.00040
<u>M</u> Bi < 0.00040	<u>M</u> Ge < 0.00595	<u>M</u> Mo 0.00050	<u>M</u> Sm < 0.00099	<u>M</u> Sn < 0.00496
<u>Q</u> B < 1.19460	<u>M</u> Au < 0.00297	<u>M</u> Nd < 0.00198	<u>Q</u> Sc < 0.00036	<u>M</u> Ti 0.00198
<u>M</u> Cd < 0.00297	<u>M</u> Hf < 0.00198	<u>M</u> Ni < 0.00793	<u>M</u> Se < 0.00793	<u>S</u> W
<u>Q</u> Ca 0.00080	<u>M</u> Ho < 0.00050	<u>Q</u> Nb < 0.06371	<u>Q</u> Si < 0.01354	<u>M</u> U < 0.00198
<u>M</u> Ce < 0.00496	<u>M</u> In < 0.00991	<u>n</u> Os	<u>M</u> Ag < 0.00198	<u>M</u> V < 0.00198
<u>M</u> Cs < 0.00030	<u>M</u> Ir < 0.00496	<u>M</u> Pd < 0.00496	<u>Q</u> Na 0.04778	<u>M</u> Yb < 0.00099
<u>M</u> Cr < 0.00496	<u>Q</u> Fe < 0.03982	<u>n</u> P	<u>M</u> Sr < 0.00050	<u>M</u> Y < 0.03965
<u>M</u> Co < 0.00297	<u>M</u> La < 0.00050	<u>M</u> Pt < 0.00198	<u>n</u> S	<u>M</u> Zn < 0.01983
<u>M</u> Cu < 0.00595	<u>M</u> Pb 0.00060	<u>Q</u> K 0.03146	<u>Q</u> Ta < 0.39820	<u>M</u> Zr 0.00079

M - Checked by ICP-MS

O - Checked by ICP-OES

I - Spectral Interference

n - Not Checked For

s - Solution Standard Element

6.0 INTENDED USE

For the calibration of analytical instruments including but not limited to the following:
ICP-MS, ICP-OES, FAAS, GFAA, XRF, and DCP

For the validation of analytical methods

For the preparation of "working reference samples"

For interference studies and the determination of correction coefficients

For detection limit and linearity studies

For additional intended uses, contact IV Technical Staff

7.0 INSTRUCTIONS FOR THE CORRECT USE OF THIS REFERENCE MATERIAL

Storage & Handling - Keep tightly sealed when not in use. Store and use at $20 \pm 4^\circ\text{C}$. Do not pipet from container. Do not return portions removed for pipetting to container.

Atomic Weight; Valence; Coordination Number; Chemical Form in Solution - 183.85; +6; 6, 7, 8, 9 WOF₆ (chemical form as received)

Chemical Compatibility - W is very readily hydrolyzed requiring 0.1 to 1% HF solutions for stable acidic solutions. The WOF₆ is soluble in % levels of HCl and HNO₃ provided it is in the WOF₆ form. Stable at ppm levels with some metals provided it is fluorinated. Do not mix with Alkaline or Rare Earths. Is best to be mixed only with other fluorinated metals (Ti, Zr, Hf, Nb, Ta, Mo, Si, Sn, Ge). Look for yellow WVO₃ precipitate if mixed with other transition elements at higher levels indicating instability. The yellow WVO₃ will form over a period of weeks even in trace HF, therefore, HF levels of W multi-element blends should be ~ 1%.

Stability - 2-100 ppb levels stable (Alone or mixed with all other metals that are at comparable levels) as the WOF₆ for months in 1% HNO₃ / LDPE container. 1-10,000 ppm single element solutions as the WOF₆ chemically stable for years in 1% HF in an LDPE container.

W Containing Samples (Preparation and Solution) - Metal (Soluble in HF / HNO₃); Oxide (Soluble in HF or NH₄OH); Organic Matrices (Dry ash at 450 °C in Pt⁺ and dissolve oxide with HF).

Atomic Spectroscopic Information (ICP-OES D.L.s are given as radial/axial view):

Technique/Line	Estimated D.L.	Order	Type	Interferences (underlined indicates severe at λ Doncs.)
ICP-OES 207.911 nm	0.03 / 0.001 $\mu\text{g/mL}$	1	ion	Ru, In
ICP-OES 224.875 nm	0.05 / 0.005 $\mu\text{g/mL}$	1	ion	Co, Rh, Ag
ICP-OES 209.475 nm	0.05 / 0.005 $\mu\text{g/mL}$	1	ion	Mo
ICP-MS 182 amu	5 ppt	n/a	M ⁺	¹⁸² Er ⁺ O

HF Note: This standard should not be prepared or stored in glass.

8.0 HAZARDOUS INFORMATION - Please refer to the enclosed Material Safety Data sheet for information regarding this CRM.

9.0 HOMOGENEITY - This solution was mixed according to procedure IV-MPM-004 and is guaranteed to be homogeneous.

10.0 QUALITY STANDARD DOCUMENTATION

10.1 ISO 9001:2000 Quality Management System Registration - QMI Certificate Number 010105

Recognized by:

Registrar Accreditation Board (ANSI-RAB)

Standards Council of Canada (SCC)

Dutch Council for Accreditation (RVA)

Entidad Mexicana de Acreditacion, a.c.(EMA)

Members of IQ Net International Certification Network:

Argentina (IRAM), Australia (QAS), Austria (ÖQS), Belgium (Avinter), Brazil (FCAV), Canada (QMI), Hong Kong (HKQAA), Columbia (ICONTEC), Czech Republic (CQS), Denmark (DS), Finland (SFS), France (AFAQ), Germany (DQS), Greece (ELOT), Hungary (MSZT), Ireland (NSAI), Israel (SII), Italy (CISQ), Japan (JQA), Korea (KSA-QA), Netherlands (KEMA), Norway (NCS), Poland (PCBC), Portugal (APCER), Singapore (PSB), Slovenia (SIQ), Spain (AENOR), Switzerland (SQS)

10.2 ISO/IEC 17025 - 1999 "General Requirements for the Competence of Testing and Calibration"

- Chemical Testing - Accredited A2LA Certificate Number 883.01



10.3 ISO/IEC Guide 34 - 2000 "General Requirements for the Competence of Reference Material Producers"

- Reference Materials Production - Accredited A2LA Certificate Number 883.02

A2LA Mutual Recognition Agreement Partners:

Australia (NATA), Austria (BmWA), Belgium (BELTEST) (BKO-OBE), Canada (SCC), Chinese Taipei (CNLA), Czech Republic (NAO), Denmark (DANAK), Finland (FINAS), France (COFRAC), Germany (DAR), Hong Kong (HKAS), Ireland (NAB), Italy (SIT) (SINAL), Japan (JAB) (JNLA), Republic of Korea (KOLAS), The Netherlands (RvA), New Zealand (IANZ), Norway (NA), Portugal (IPQ), Singapore (SAC-SINGLAS), Spain (ENAC), Sweden (SWEDAC), Switzerland (SAS), United Kingdom (UKAS) and United States (NVLAP) (ICBO ES)

10.4 10CFR50 Appendix B - Nuclear Regulatory Commission

- Domestic Licensing of Production and Utilization Facilities

10.5 10CFR21 - Nuclear Regulatory Commission - Reporting Defects and Non-Compliance

10.6 MIL-STD-45662A (Obsolete/Observed)

INORGANIC LABS/RADCHEM LABS Pg. 2 of 2

DATE RECEIVED: 07/31/03

DATE EXPIRED: 08/01/2004

DATE OPENED: 08/01/03

INORG: 4203 PO: F52383

11.0 DATE OF CERTIFICATION AND PERIOD OF VALIDITY

010242



11.1 IV Shelf Life - The period of time during which the concentration of the analyte(s) in a properly packaged, unopened, and unused standard stored under environmentally controlled and monitored conditions will remain within the specified uncertainty range. Shelf life is limited primarily by transpiration (loss of water from the solution) and infrequently, by chemical instability. Transpiration studies (P-SP01020) of chemically-stable solutions performed at Inorganic Ventures / IV Labs indicate a CRM shelf-life of four years for solutions packaged in 500-mL low density polyethylene bottles. When stored under special conditions that minimize transpiration and instability, the shelf life can be extended past this limit.

11.2 Expiration Date - The date after which a CRM should not be used. Routine laboratory use of a CRM increases transpiration losses and the chance of contamination which affect the integrity of the CRM and limit its useful life. Inorganic Ventures / IV Labs concurs with state and federal regulatory agencies' recommendations that solution standards be assigned a one-year expiration date.

Certification Date: February 10, 2003

Expiration Date: **EXPIRES**
12/2004

12.0 NAMES AND SIGNATURES OF CERTIFYING OFFICERS

Certificate Prepared By: JoAnn Struthers, QA Administrative Assistant

JoAnn Struthers

Certificate Approved By: Katalin Le, QC Supervisor

Katalin Le

Certifying Officer: Paul Gaines, Chemist, Senior Technical Director

Paul Gaines



Certificate of Analysis

CUSTOM-GRADE SOLUTION

1000 µg/mL Yttrium in 1.4% HNO₃ (abs)

Catalog Number: CGY1-1 and CGY1-5

010243

Lot Number: **T-Y01086**

INORGANIC LABS/RADCHEM LABS

Starting Material:
Starting Material Purity:
Starting Material Lot No:

Y₂O₃
99.999%
Y-O-5-9503201

DATE RECEIVED: 05/09/03
DATE EXPIRED: 06/01/2004
DATE OPENED: 05/12/03
INORG: 4063 PD: F52181

CERTIFIED CONCENTRATION: 1004 ± 3 µg/mL

The Certified Value is based upon the most precise method used to analyze this CRM. The following equations are used in the calculation of the certified value and the uncertainty:

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

$$\text{Uncertainty } (\pm) = \frac{2[(\sum s_i)^2]^{1/2}}{(n)^{1/2}}$$

($\bar{x}$) = mean

x_i = individual results

n = number of measurements

$\sum s_i$ = The summation of all significant estimated errors.

Classical Wet Assay: 1004 ± 3 µg/mL

Method: EDTA Titration vs NIST SRM 928 Lead Nitrate.

Instrument Analysis: 1004 ± 4 µg/mL

Method: Inductively Coupled Plasma Spectroscopy (ICP) vs NIST SRM 3167a.

The independent samples t-test was used to determine if there is agreement between the above assay methods at the 95% confidence interval. Both methods were compared and showed agreement within the stated uncertainties. This agreement is a confirmation of the accuracy of this CRM.

TRACE METALLIC IMPURITIES DETERMINED BY ICP-MS AND ICP-OES IN µg/mL:

Custom-Grade solutions tested for trace metallic impurities by ICP-MS were analyzed in an ULPA-Filtered Clean Room.

An ULPA-Filter is 99.9985% efficient for the removal of particles down to 0.3 µm.

<u>Q</u> Al 0.00021	<u>M</u> Dy <0.0060	<u>M</u> Li <0.010	<u>M</u> Pr <0.00030	<u>M</u> Te <0.030
<u>M</u> Sb <0.00050	<u>M</u> Er <0.0050	<u>Q</u> Lu <0.000020	<u>M</u> Re <0.0010	<u>M</u> Tb <0.00030
<u>M</u> As <0.010	<u>M</u> Eu <0.0030	<u>Q</u> Mg <0.000030	<u>M</u> Rh <0.0010	<u>M</u> Tl <0.0010
<u>M</u> Ba <0.010	<u>M</u> Gd <0.0010	<u>Q</u> Mn <0.000020	<u>M</u> Rb <0.0010	<u>M</u> Th <0.0010
<u>M</u> Be <0.00050	<u>M</u> Ga <0.0010	<u>Q</u> Hg <0.020	<u>M</u> Ru <0.0020	<u>M</u> Tm <0.00040
<u>M</u> Bi <0.00040	<u>M</u> Ge <0.0060	<u>M</u> Mo <0.0020	<u>M</u> Sm <0.0010	<u>M</u> Sn <0.0050
<u>Q</u> B <0.00014	<u>M</u> Au <0.0030	<u>M</u> Nd <0.0020	<u>Q</u> Sc 0.00024	<u>M</u> Ti <0.050
<u>M</u> Cd <0.0030	<u>M</u> Hf <0.0020	<u>M</u> Ni <0.080	<u>M</u> Se <0.0080	<u>M</u> W <0.010
<u>Q</u> Ca 0.00014	<u>M</u> Ho <0.00050	<u>M</u> Nb <0.00050	<u>Q</u> Si 0.000077	<u>M</u> U <0.0020
<u>M</u> Ce <0.0050	<u>M</u> In <0.010	<u>n</u> Os	<u>M</u> Ag <0.0020	<u>Q</u> V <0.00080
<u>M</u> Cs <0.00030	<u>M</u> Ir <0.0050	<u>Q</u> Pd <0.10	<u>Q</u> Na <0.050	<u>M</u> Yb <0.0010
<u>M</u> Cr <0.0050	<u>Q</u> Fe 0.0015	<u>Q</u> P <0.070	<u>M</u> Sr <0.00050	<u>s</u> Y
<u>M</u> Co <0.0030	<u>M</u> La <0.00050	<u>M</u> Pt <0.0020	<u>Q</u> S <0.043	<u>M</u> Zn <0.020
<u>M</u> Cu <0.0060	<u>M</u> Pb <0.0030	<u>Q</u> K <0.10	<u>M</u> Ta <0.0070	<u>Q</u> Zr <0.00070

M - checked by ICP-MS

O - checked by ICP-OES

i - spectral interference

n - not checked for

s - solution standard element

ANALYZED DENSITY OF SOLUTION (measured at 22°C): 1.010 g/mL

(over)

QA:KL Rev.061802DN

Inorganic Ventures, Inc.

195 Lehigh Avenue • Suite 4 • Lakewood, NJ 08701

Orders: 800-669-6799 • FAX (732) 901-1903

Technical Support: 800-569-6799


Quality Assurance Manager

EXPIRES

12/2004



Certificate of Analysis



CUSTOM-GRADE SOLUTION

1000 µg/mL Zirconium in H₂O tr. HF tr. HNO₃ **010244**

Catalog Number: CGZR1-1 and CGZR1-5

This standard should not be prepared or stored in glass.

Lot Number: **T-QZR01062**

Starting Material:
Starting Material Purity:
Starting Material Lot No:

Zirconium Metal
99.999%
SP4811A

INORGANIC LABS/RADCHEM LABS

DATE RECEIVED: 02/13/03

DATE EXPIRED: 03/01/2004

DATE OPENED: 02/13/03

INORG: 3958 PO: F52060

CERTIFIED CONCENTRATION: 1005 ± 6 µg/mL

The Certified Value is based upon the most precise method used to analyze this CRM. The following equations are used in the calculation of the certified value and the uncertainty:

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

$$\text{Uncertainty } (\pm) = \frac{2[(\sum s_i)^2]^{1/2}}{(n)^{1/2}}$$

($\bar{x}$) = mean

x_i = individual results

n = number of measurements

$\sum S_i$ = The summation of all significant estimated errors.

Instrument Analysis: 1005 ± 6 µg/mL (Average of 2 runs)

Method: Inductively Coupled Plasma Spectroscopy (ICP) vs NIST SRM 3169

Calculated Value: 1000 µg/mL (For informational purposes)

Method: Calculated based on starting material.

TRACE METALLIC IMPURITIES DETERMINED BY ICP-MS AND ICP-OES IN µg/mL:

Custom-Grade solutions tested for trace metallic impurities by ICP-MS were analyzed in an ULPA-Filtered Clean Room.
An ULPA-Filter is 99.9985% efficient for the removal of particles down to 0.3 µm.

<u>O</u>	Al	0.055	<u>M</u>	Dy	<0.012	<u>O</u>	Li	<0.000030	<u>M</u>	Pr	<0.00059	<u>M</u>	Te	<0.059
<u>M</u>	Sb	0.00017	<u>M</u>	Er	<0.0098	<u>M</u>	Lu	<0.00079	<u>M</u>	Re	<0.0020	<u>M</u>	Tb	<0.00059
<u>M</u>	As	<0.020	<u>M</u>	Eu	<0.0059	<u>O</u>	Mg	<0.000030	<u>M</u>	Rh	<0.0020	<u>M</u>	Tl	<0.0020
<u>M</u>	Ba	0.013	<u>M</u>	Gd	<0.0020	<u>O</u>	Mn	0.0020	<u>M</u>	Rb	<0.0020	<u>M</u>	Th	<0.0020
<u>O</u>	Be	<0.10	<u>M</u>	Ga	<0.0020	<u>O</u>	Hg	<0.011	<u>M</u>	Ru	<0.0040	<u>M</u>	Tm	<0.00079
<u>M</u>	Bi	<0.00079	<u>M</u>	Ge	<0.012	<u>O</u>	Mo	<0.10	<u>M</u>	Sm	<0.0020	<u>M</u>	Sn	0.0011
<u>i</u>	B		<u>M</u>	Au	<0.0059	<u>M</u>	Nd	<0.0040	<u>O</u>	Sc	<0.00016	<u>O</u>	Ti	0.0031
<u>O</u>	Cd	<0.0053	<u>M</u>	Hf	0.087	<u>O</u>	Ni	0.014	<u>M</u>	Se	<0.016	<u>M</u>	W	0.0016
<u>O</u>	Ca	0.0020	<u>M</u>	Ho	<0.00098	<u>i</u>	Nb		<u>i</u>	Si		<u>M</u>	U	0.00040
<u>M</u>	Ce	<0.0098	<u>M</u>	In	<0.020	<u>n</u>	Os		<u>i</u>	Ag		<u>M</u>	V	0.0016
<u>M</u>	Cs	<0.00059	<u>M</u>	Ir	<0.0098	<u>i</u>	Pd		<u>O</u>	Na	0.041	<u>M</u>	Yb	<0.0020
<u>O</u>	Cr	0.034	<u>O</u>	Fe	0.45	<u>O</u>	P	<0.0048	<u>M</u>	Sr	<0.00098	<u>O</u>	Y	<0.0010
<u>M</u>	Co	<0.0059	<u>M</u>	La	<0.00098	<u>M</u>	Pt	<0.0040	<u>i</u>	S		<u>O</u>	Zn	0.0057
<u>M</u>	Cu	0.013	<u>M</u>	Pb	<0.0059	<u>O</u>	K	0.013	<u>M</u>	Ta	0.0011	<u>s</u>	Zr	

M - checked by ICP-MS

O - checked by ICP-OES

i - spectral interference

n - not checked for

s - solution standard element

ANALYZED DENSITY OF SOLUTION (measured at 22°C): **1.000 g/mL**

QA:KL Rev.100302DN

(over)

Inorganic Ventures, Inc.

195 Lehigh Avenue • Suite 4 • Lakewood, NJ 08701

Orders: 800-669-6799 • FAX (732) 901-1903

Technical Support: 800-569-6799

Quality Assurance Manager

EXPIRES

1A2004



Certificate of Analysis

CUSTOM GRADE SOLUTION 1000 µg/mL Boron in H₂O

Catalog Number: CGB1-1, CGB1-2 and CGB1-5

010245

Lot Number: **T-B02038**

Starting Material:
Starting Material Purity: #
Starting Material Lot No:

H₃BO₃
99.999%
S-B02032

INORGANIC LABS/RADCHEM LABS

DATE RECEIVED: 02/13/03

DATE EXPIRED: 03/01/2004

DATE OPENED: 02/13/03

INORG: 3954 PO: F52060

CERTIFIED CONCENTRATION: 1001 ± 3 µg/mL

The Certified Value is the instrument analysis value. The following equations are used in the calculation of the certified value and the uncertainty:

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

$$\text{Uncertainty } (\pm) = \frac{2(\sum s_i^2)^{1/2}}{(n)^{1/2}}$$

($\bar{x}$) = mean

x_i = individual results

n = number of measurements

$\sum S_i$ = The summation of all significant estimated errors.

Instrument Analysis: 1001 ± 3 µg/mL

Method: Inductively Coupled Plasma Spectroscopy (ICP) vs NIST SRM 3107.

TRACE METALLIC IMPURITIES DETERMINED BY ICP-MS AND ICP-OES IN µg/mL:

Custom Grade solutions tested for trace metallic impurities by ICP-MS were analyzed in an ULPA-Filtered Clean Room.
An ULPA-Filter is 99.9985% efficient for the removal of particles down to 0.3 µm.

<u>Q</u> Al <0.00090	<u>M</u> Dy <0.0060	<u>Q</u> Li <0.000024	<u>M</u> Pr <0.00030	<u>M</u> Te <0.030
<u>M</u> Sb <0.00050	<u>M</u> Er <0.0050	<u>M</u> Lu <0.00040	<u>M</u> Re <0.0010	<u>M</u> Tb <0.00030
<u>M</u> As <0.010	<u>M</u> Eu <0.0030	<u>Q</u> Mg <0.000060	<u>M</u> Rh <0.0010	<u>M</u> Tl <0.0010
<u>Q</u> Ba <0.00010	<u>M</u> Gd <0.0010	<u>Q</u> Mn <0.000023	<u>M</u> Rb <0.0010	<u>M</u> Th <0.0010
<u>Q</u> Be <0.00017	<u>Q</u> Ga <0.0016	<u>Q</u> Hg <0.015	<u>M</u> Ru <0.0020	<u>M</u> Tm <0.00040
<u>M</u> Bi <0.00040	<u>M</u> Ge <0.0060	<u>M</u> Mo <0.0020	<u>M</u> Sm <0.0010	<u>M</u> Sn <0.0050
<u>s</u> B	<u>M</u> Au <0.0030	<u>M</u> Nd <0.0020	<u>Q</u> Sc <0.000023	<u>M</u> Ti <0.050
<u>M</u> Cd <0.0030	<u>M</u> Hf <0.0020	<u>Q</u> Ni <0.0023	<u>Q</u> Se <0.0062	<u>M</u> W <0.010
<u>Q</u> Ca <0.000070	<u>M</u> Ho <0.00050	<u>M</u> Nb <0.00050	<u>Q</u> Si 0.00067	<u>M</u> U <0.0020
<u>Q</u> Ce <0.0030	<u>M</u> In <0.010	<u>n</u> Os	<u>M</u> Ag <0.0020	<u>Q</u> V <0.00083
<u>M</u> Cs <0.00030	<u>M</u> Ir <0.0050	<u>M</u> Pd <0.0050	<u>Q</u> Na <0.00010	<u>M</u> Yb <0.0010
<u>M</u> Cr <0.0050	<u>Q</u> Fe <0.0011	<u>Q</u> P <0.0025	<u>M</u> Sr <0.00050	<u>M</u> Y <0.040
<u>Q</u> Co <0.0011	<u>M</u> La <0.00050	<u>M</u> Pt <0.0020	<u>Q</u> S <0.10	<u>Q</u> Zn <0.00019
<u>M</u> Cu <0.0060	<u>M</u> Pb <0.0030	<u>Q</u> K <0.0030	<u>M</u> Ta <0.0070	<u>M</u> Zr <0.0050

M - checked by ICP-MS

O - checked by ICP-OES

i - spectral interference

n - not checked for

s - solution standard element

ANALYZED DENSITY OF SOLUTION (measured at 22°C): 0.999 g/mL

QA:KL Rev.073102DN

(over)

Inorganic Ventures, Inc.

195 Lehigh Avenue • Suite 4 • Lakewood, NJ 08701


Orders: 800-669-6799 • FAX (732) 901-1903

Technical Support: 800-569-6799

Quality Assurance Manager

EXPIRES

1/2004



Certificate of Analysis

CUSTOM-GRADE SOLUTION 1000 µg/mL Molybdenum in H₂O tr. NH₄OH

Catalog Number: CGMO1-1, CGMO1-2 and CGMO1-5

010248

Lot Number: **T-MO01130**

Starting Material: Ammonium Molybdenum Oxide
Starting Material Purity: 99.9943%
Starting Material Lot No: 21410

CERTIFIED CONCENTRATION: 997 ± 3 µg/mL

The Certified Value is based upon the most precise method used to analyze this CRM. The following equations are used in the calculation of the certified value and the uncertainty:

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

$$\text{Uncertainty } (\pm) = \frac{2[(\sum s_i)^2]^{1/2}}{(n)^{1/2}}$$

($\bar{x}$) = mean x_i = individual results n = number of measurements $\sum s_i$ = The summation of all significant estimated errors.

Calculated Value: 1000 µg/mL
Method: Calculated, based on starting material.

Instrument Analysis: 997 ± 3 µg/mL (Average of 2 runs)
Method: Inductively Coupled Plasma Spectroscopy (ICP) vs NIST SRM 3134.

TRACE METALLIC IMPURITIES DETERMINED BY ICP-MS AND ICP-OES IN µg/mL:

Custom-Grade solutions tested for trace metallic impurities by ICP-MS were analyzed in an ULPA-Filtered Clean Room.
An ULPA-Filter is 99.9985% efficient for the removal of particles down to 0.3 µm.

<u>O</u>	Al	<0.050	<u>M</u>	Dy	<0.012	<u>O</u>	Li	<0.010	<u>O</u>	Pr	<0.10	<u>i</u>	Te	
<u>M</u>	Sb	0.0094	<u>M</u>	Er	<0.010	<u>M</u>	Lu	<0.00080	<u>M</u>	Re	<0.0020	<u>M</u>	Tb	<0.00060
<u>M</u>	As	<0.020	<u>M</u>	Eu	<0.0060	<u>O</u>	Mg	<0.050	<u>M</u>	Rh	<0.0020	<u>M</u>	Tl	<0.0020
<u>M</u>	Ba	<0.020	<u>M</u>	Gd	<0.0020	<u>M</u>	Mn	<0.0080	<u>M</u>	Rb	0.025	<u>M</u>	Th	<0.0020
<u>M</u>	Be	<0.0010	<u>M</u>	Ga	<0.0020	<u>i</u>	Hg		<u>M</u>	Ru	<0.0040	<u>M</u>	Tm	<0.00080
<u>M</u>	Bi	<0.00080	<u>M</u>	Ge	<0.012	<u>s</u>	Mo		<u>M</u>	Sm	<0.0020	<u>M</u>	Sn	<0.010
<u>O</u>	B	<0.50	<u>M</u>	Au	<0.0060	<u>O</u>	Nd	<0.050	<u>O</u>	Sc	<0.050	<u>O</u>	Ti	<0.0050
<u>O</u>	Cd	<0.50	<u>M</u>	Hf	<0.0040	<u>M</u>	Ni	<0.016	<u>M</u>	Se	<0.016	<u>M</u>	W	0.056*
<u>O</u>	Ca	0.00026	<u>M</u>	Ho	<0.0010	<u>O</u>	Nb	<0.10	<u>O</u>	Si	<0.10	<u>M</u>	U	<0.0040
<u>O</u>	Ce	<0.050	<u>M</u>	In	0.0024	<u>n</u>	Os		<u>M</u>	Ag	<0.0040	<u>M</u>	V	<0.0040
<u>M</u>	Cs	<0.00060	<u>M</u>	Ir	<0.010	<u>M</u>	Pd	<0.010	<u>O</u>	Na	<0.10	<u>M</u>	Yb	<0.0020
<u>M</u>	Cr	<0.010	<u>O</u>	Fe	<0.50	<u>i</u>	P		<u>M</u>	Sr	<0.0010	<u>M</u>	Y	<0.080
<u>M</u>	Co	<0.0060	<u>M</u>	La	<0.0010	<u>M</u>	Pt	<0.0040	<u>i</u>	S		<u>M</u>	Zn	<0.040
<u>M</u>	Cu	<0.012	<u>M</u>	Pb	<0.0060	<u>O</u>	K	0.0098	<u>M</u>	Ta	<0.014	<u>M</u>	Zr	<0.010

M - checked by ICP-MS O - checked by ICP-OES i - spectral interference n - not checked for s - solution standard element

* W impurity may be due to Mo doubly charged ion and represents the maximum concentration possible.

ANALYZED DENSITY OF SOLUTION (measured at 22°C): 0.999 g/mL

(over)

QA:KL Rev.0731020N



Inorganic Ventures, Inc.

195 Lehigh Avenue • Suite 4 • Lakewood, NJ 08701

Orders: 800-669-6799 • FAX (732) 901-1903

Technical Support: 800-569-6799

Quality Assurance Manager

EXPIRES
1st 2004



Certificate of Analysis

CUSTOM-GRADE SOLUTION

Catalog Number: CGP1-1, CGP1-2 and CGP1-5

1000 µg/mL Phosphorus in H₂O

Lot Number: T-P01118

Starting Material:
Starting Material Purity:
Starting Material Lot No:

Phosphoric Acid
99.999%
J18804

010247
INORGANIC LABS/RADCHEM LABS
DATE RECEIVED: 05/01/03
DATE EXPIRED: 06/01/2004
DATE OPENED: 05/01/03
INORG: 4049 PO: F52180

CERTIFIED CONCENTRATION: 1004 ± 2 µg/mL

The Certified Value is based upon the most precise method used to analyze this CRM. The following equations are used in the calculation of the certified value and the uncertainty:

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

$$\text{Uncertainty } (\pm) = \frac{2(\sum s_i^2)^{1/2}}{(n)^{1/2}}$$

($\bar{x}$) = mean

x_i = individual results

n = number of measurements

$\sum s_i$ = The summation of all significant estimated errors.

Classical Wet Assay: 1004 ± 2 µg/mL

Method: Acidimetric Titration vs NIST SRM 84k KHP.

Instrument Analysis: 1000 ± 3 µg/mL

Method: Inductively Coupled Plasma Spectroscopy (ICP) vs NIST SRM 3139a.

The independent samples t-test was used to determine if there is agreement between the above assay methods at the 95% confidence interval. Both methods were compared and showed agreement within the stated uncertainties. This agreement is a confirmation of the accuracy of this CRM.

TRACE METALLIC IMPURITIES DETERMINED BY ICP-MS AND ICP-OES IN µg/mL:

Custom-Grade solutions tested for trace metallic impurities by ICP-MS were analyzed in an ULPA-Filtered Clean Room.
An ULPA-Filter is 99.9985% efficient for the removal of particles down to 0.3 µm.

<u>O</u> Al <0.040	<u>M</u> Dy <0.00060	<u>M</u> Li <0.0010	<u>M</u> Pr <0.000030	<u>M</u> Te <0.0030
<u>M</u> Sb 0.012	<u>M</u> Er <0.00050	<u>M</u> Lu <0.000040	<u>M</u> Re <0.00010	<u>M</u> Tb <0.000030
<u>M</u> As <0.0010	<u>M</u> Eu <0.00030	<u>M</u> Mg <0.0030	<u>M</u> Rh <0.00010	<u>M</u> Tl <0.00010
<u>M</u> Ba <0.0010	<u>M</u> Gd <0.00010	<u>M</u> Mn <0.00040	<u>M</u> Rb <0.00010	<u>M</u> Th <0.00010
<u>M</u> Be <0.000050	<u>M</u> Ga 0.00070	<u>O</u> Hg <0.020	<u>M</u> Ru <0.00020	<u>M</u> Tm <0.000040
<u>M</u> Bi <0.000040	<u>M</u> Ge <0.00060	<u>M</u> Mo <0.00020	<u>M</u> Sm <0.00010	<u>M</u> Sn <0.00050
<u>M</u> B <0.0070	<u>M</u> Au <0.00030	<u>M</u> Nd <0.00020	<u>n</u> Sc	<u>n</u> Ti
<u>M</u> Cd <0.00030	<u>M</u> Hf <0.00020	<u>O</u> Ni <0.050	<u>O</u> Se <0.40	<u>M</u> W <0.0010
<u>O</u> Ca <0.010	<u>M</u> Ho <0.000050	<u>M</u> Nb <0.000050	<u>O</u> Si <0.020	<u>M</u> U <0.00020
<u>M</u> Ce <0.00050	<u>M</u> In <0.030	<u>n</u> Os	<u>M</u> Ag <0.00020	<u>M</u> V <0.00020
<u>M</u> Cs <0.000030	<u>M</u> Ir <0.00050	<u>M</u> Pd <0.00050	<u>O</u> Na <0.090	<u>M</u> Yb <0.00010
<u>M</u> Cr <0.00050	<u>O</u> Fe <0.050	<u>s</u> P	<u>M</u> Sr <0.000050	<u>M</u> Y <0.0040
<u>M</u> Co <0.00030	<u>M</u> La <0.000050	<u>M</u> Pt <0.00020	<u>n</u> S	<u>M</u> Zn 0.0035
<u>M</u> Cu 0.080	<u>M</u> Pb <0.00030	<u>n</u> K	<u>M</u> Ta <0.00070	<u>M</u> Zr <0.00050

M - checked by ICP-MS

O - checked by ICP-OES

i - spectral interference

n - not checked for

s - solution standard element

ANALYZED DENSITY OF SOLUTION (measured at 22°C): 0.999 g/mL

(over)

QA:KL Rev.122002DN

Inorganic Ventures, Inc.

195 Lehigh Avenue • Suite 4 • Lakewood, NJ 08701

Orders: 800-669-6799 • FAX (732) 901-1903

Technical Support: 800-569-6799

Quality Assurance Manager

EXPIRES

01/2004



inorganic ventures / iv labs

195 lehigh avenue, suite 4, lakewood, nj 08701 usa
phone: 800-669-6799 • 732-901-1900 • fax: 732-901-1903
e-mail: ivsales@ivstandards.com • website: www.ivstandards.com

certificate of analysis 010248

1.0 Inorganic Ventures / IV Labs is an ISO Guide 34-2000 Certified Reference Material (CRM) Manufacturer: Certificate #883-02. The certificate is designed and the certified value(s) and uncertainty(ies) are determined in accordance with ISO Guide 31-2000 (Reference Materials - Contents of certificates and label(s)), ISO Guide 34-2000 "Quality System Guidelines for the Production of Reference Materials," and ISO Guide 35-1989 "Certification of Reference Materials - General and Statistical Principles."

2.0 DESCRIPTION OF CRM Custom-Grade 1000 µg/mL Silicon in H₂O tr. HNO₃ tr. HF

Catalog Number: CGSI1-1, CGSI1-2, and CGSI1-5

Lot Number: W-SI02082

Starting Material: SiO₂

Starting Material Purity (%): 99.996367

Starting Material Lot No C05310C

Matrix: H₂O tr. HNO₃ tr. HF

INORGANIC LABS/RADCHEM LABS Pg. 1 of 2
DATE RECEIVED: 06/20/03
DATE EXPIRED: 07/01/2004
DATE OPENED: 06/23/03
INORG: 4153 PO: F52370

3.0 CERTIFIED VALUES AND UNCERTAINTIES

Certified Concentration: 1000 ± 5 µg/mL

Certified Density: 1.002 g/mL (measured at 22° C)

The Certified Value is based upon the most precise method used to analyze this CRM. The following equations are used in the calculation of the certified value and the uncertainty:

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

($\bar{x}$) = mean

x_i = individual results

n = number of measurements

$$\text{Uncertainty } (\pm) = \frac{2(s)}{(n)^{1/2}}$$

$\sum$ = The summation of all significant estimated errors.

(Most common are the errors from instrumental measurement, weighing, dilution to volume, and the fixed error reported on the NIST SRM certificate of analysis.)

4.0 TRACEABILITY TO NIST AND VALUES OBTAINED BY INDEPENDENT METHODS

□ "Property of the result of a measurement or the value of a standard whereby it can be related to stated references, usually national or international standards, through an unbroken chain of comparisons all having stated uncertainties." (ISO VIM, 2nd ed., 1993, definition 6.10)

□ This IV product is Traceable to NIST via direct comparison to NIST SRMs. The uncertainties for each certified value are reported, taking into account the SRM uncertainty error and the measurement, weighing and volume dilution errors.

4.1 Assay Method #1 1000 ± 5 µg/mL (Avg 2 runs)

ICP Assay NIST SRM 3150 Lot Number: 991108

Assay Method #2 1001 µg/mL

Gravimetric NIST SRM Lot Number: See Sec. 4.2

4.2 BALANCE CALIBRATION - All balances are checked daily using in-house procedure number 6-IMM-001. The weights used for testing are annually compared to Gerhart Scale Corporation's master weights and are traceable to the National Institute of Standards and Technology (NIST). The NIST Traceability numbers are 692476 - Class 1 and 692476A - Class 2. The NIST test number is 822/260017-98. All analytical balances are calibrated every 4 months by Gerhart Scale Corp. of South Amboy. The balances are calibrated with a class 1 and/or class 2 analytical weight set. These weights are tested annually by a NIST / NVLAP accredited calibration lab. The NIST test number is 822/260017-98.

4.3 THERMOMETER CALIBRATION - The thermometers used in the determination of the final densities are calibrated vs standard thermometer No. 903-2680 which was certified in accordance with the procedures outlined by ASTM E77-87 and NIST Monograph 150 using NIST Test Nos. and Std Nos.: 769543, 217368/769543, 217368/P14452, 176240/P14452, 176240. The in-house procedure No. is 2-QC-001. Thermometers which are not calibrated vs standard thermometer No. 903-2680 are traceable to NIST Identification Nos. 92564, 119016, 471047 and NIST test report Nos. 811/258522, 811/2557078, and 236090.

4.4 GLASSWARE CALIBRATION - In-house procedure 3-QC-002 is used to calibrate all Class A Glassware used in the manufacture and quality control of Custom Grade Standards.

5.0 TRACE METALLIC IMPURITIES (TMI) DETERMINED BY ICP/MS AND ICP-OES IN µg/mL

Custom-Grade solutions are tested for trace metallic impurities by Axial ICP-OES and ICP-MS. The result from the most sensitive method for each element, is reported below. Solutions tested by ICP-MS were analyzed in an ULPA-Filtered Clean Room. An ULPA-Filter is 99.9985% efficient for the removal of particles down to 0.3 µm.

<u>O</u> Al 0.02730	<u>M</u> Dy < 0.01358	<u>O</u> Li < 0.00009	<u>M</u> Pr < 0.00068	<u>M</u> Te < 0.06791
<u>M</u> Sb < 0.00113	<u>M</u> Er < 0.01132	<u>M</u> Lu < 0.00091	<u>M</u> Re < 0.00226	<u>M</u> Tb < 0.00068
<u>M</u> As < 0.02264	<u>M</u> Eu < 0.00679	<u>O</u> Mg < 0.04991	<u>M</u> Rh < 0.00226	<u>M</u> Tl < 0.00226
<u>M</u> Ba < 0.02264	<u>M</u> Gd < 0.00226	<u>M</u> Mn < 0.00906	<u>M</u> Rb < 0.00226	<u>M</u> Th < 0.00226
<u>O</u> Be < 0.00091	<u>M</u> Ga < 0.00226	<u>O</u> Hg < 0.04991	<u>M</u> Ru < 0.00453	<u>M</u> Tm < 0.00091
<u>M</u> Bi < 0.00091	<u>M</u> Ge < 0.01358	<u>M</u> Mo < 0.00453	<u>M</u> Sm < 0.00226	<u>M</u> Sn < 0.01132
<u>O</u> B 0.02409	<u>M</u> Au < 0.00679	<u>M</u> Nd < 0.00453	<u>O</u> Sc < 0.00091	<u>O</u> Ti 0.01325
<u>M</u> Cd < 0.00679	<u>M</u> Hf < 0.00453	<u>O</u> Ni < 0.01044	<u>M</u> Se < 0.01811	<u>M</u> W < 0.02264
<u>O</u> Ca 0.00135	<u>M</u> Ho < 0.00113	<u>M</u> Nb < 0.00113	<u>s</u> Si	<u>M</u> U < 0.00453
<u>M</u> Ce < 0.01132	<u>M</u> In < 0.02264	<u>n</u> Os	<u>M</u> Ag < 0.00453	<u>O</u> V < 0.00408
<u>M</u> Cs < 0.00068	<u>M</u> Ir < 0.01132	<u>M</u> Pd < 0.01132	<u>O</u> Na 0.02008	<u>M</u> Yb < 0.00226
<u>O</u> Cr < 0.00681	<u>O</u> Fe < 0.00499	<u>O</u> P < 0.02269	<u>O</u> Sr < 0.00032	<u>M</u> Y < 0.09055
<u>M</u> Co < 0.00679	<u>M</u> La < 0.00113	<u>M</u> Pt < 0.00453	<u>O</u> S < 0.11342	<u>M</u> Zn < 0.04528
<u>O</u> Cu < 0.00454	<u>M</u> Pb < 0.00679	<u>O</u> K < 0.00771	<u>M</u> Ta 0.00200	<u>M</u> Zr < 0.01132

M - Checked by ICP-MS O - Checked by ICP-OES i - Spectral Interference n - Not Checked For s - Solution Standard Element

6.0 INTENDED USE

For the calibration of analytical instruments including but not limited to the following:

ICP-MS, ICP-OES, FAAS, GFAA, XRF, and DCP

For the validation of analytical methods

For the preparation of "working reference samples"

For interference studies and the determination of correction coefficients

For detection limit and linearity studies

For additional intended uses, contact IV Technical Staff

7.0 INSTRUCTIONS FOR THE CORRECT USE OF THIS REFERENCE MATERIAL

Storage & Handling - Keep tightly sealed when not in use. Store and use at $20 \pm 4^\circ\text{C}$. Do not pipet from container. Do not return portions removed for pipetting to container.

Atomic Weight; Valence; Coordination Number; Chemical Form in Solution - 28.0855; +4; 6; $\text{Si}(\text{OH})_4(\text{F})_2$

010250

Chemical Compatibility - Soluble in HCl, HF, H_3PO_4 , H_2SO_4 and HNO_3 as the $\text{Si}(\text{OH})_4(\text{F})_2$. Avoid neutral to basic media. Unstable at ppm levels with metals that would pull F⁻ away (i.e. Do not mix with Alkaline or Rare Earths, or high levels of transition elements unless they are fluorinated). Stable with most inorganic anions with a tendency to hydrolyze forming silicic acid (silicic acid is soluble up to ~100 ppm in water) in all dilute acids except HF.

Stability - 2-100 ppb levels - stability unknown - (alone or mixed with all other metals) as the $\text{Si}(\text{OH})_4(\text{F})_2$. 1-10,000 ppm single element solutions as the $\text{Si}(\text{OH})_4(\text{F})_2$ chemically stable for years in 2-5 % HNO_3 / trace HF in a LDPE container.

Si Containing Samples (Preparation and Solution) - Metal (Soluble in 1:1:1 H_2O / HF / HNO_3) Oxide - SiO_2 , amorphous (Dissolve by heating in 1:1:1 H_2O / HF / HNO_3) Oxide - quartz (Fuse in Pt^{a} with Na_2CO_3); Geological Samples (Fuse in Pt^{a} with Na_2CO_3 followed by HCl solution of the fuseate); Organic Matrices containing silicates and non volatile silicon compounds (Dry ash at 450°C in Pt^{a} and dissolve by gently warming with 1:1:1 H_2O / HF / H_2SO_4 or fuse / ash with Na_2CO_3 and dissolve fuseate with HCl / H_2O); Silicone Oils - dimethyl silicones depolymerize to form volatile monomer units when heated (Measure directly in alcoholic KOH / xylene mixture where sample is treated first with the KOH at $60-100^\circ\text{C}$ to "unzip" the Si-O-Si polymeric structure or digest with concentrated $\text{H}_2\text{SO}_4/\text{H}_2\text{O}_2$ followed by cooling and dissolution of the dehydrated silica with HF.) Note that the direct analysis of silicone oils in an organic solvent will result in false high results due to high vapor pressure of volatile monomer units like hexamethylcyclotrisiloxane. The KOH forms the $\text{K}_2\text{Si}(\text{CH}_3)_2\text{O}^-$ salt which is not volatile at room temperature.

Atomic Spectroscopic Information (ICP-OES D.L.s are given as radial/axial view):

Technique/Line	Estimated D.L.	Order	Type	Interferences (underlined indicates severe at 10ppb)
ICP-OES 251.611 nm	0.012 / 0.003 $\mu\text{g/mL}$	1	ion	Ta, U, Zn, Th
ICP-OES 212.412 nm	0.02 / 0.01 $\mu\text{g/mL}$	1	ion	Hf, Os, <u>Mo</u> , Ta
ICP-OES 288.158 nm	0.03 / 0.004 $\mu\text{g/mL}$	1	ion	<u>Ta</u> , Ce, Cr, Cd, Th
ICP-MS 28 amu	4000 - 8000 ppt	n/a	M ⁺	¹⁴ N ₂ , ¹² C ¹⁸ O

HF Note: This standard should not be prepared or stored in glass.

8.0 HAZARDOUS INFORMATION - Please refer to the enclosed Material Safety Data sheet for information regarding this CRM.

9.0 HOMOGENEITY - This solution was mixed according to procedure IV-MPM-004 and is guaranteed to be homogeneous.

10.0 QUALITY STANDARD DOCUMENTATION

10.1 ISO 9001:2000 Quality Management System Registration - QMI Certificate Number 010105

Recognized by:

Registrar Accreditation Board (ANSI-RAB)

Standards Council of Canada (SCC)

Dutch Council for Accreditation (RVA)

Entidad Mexicana de Acreditacion, a.c.(EMA)

Members of IQ Net International Certification Network:

Argentina (IRAM), Australia (QAS), Austria (ÖQS), Belgium (Avinter), Brazil (FCAV), Canada (QMI), Hong Kong (HKQAA), Columbia (ICONTEC), Czech Republic (CQS), Denmark (DS), Finland (SFS), France (AFAQ), Germany (DQS), Greece (ELOT), Hungary (MSZT), Ireland (NSAI), Israel (SII), Italy (CISQ), Japan (JQA), Korea (KSA-QA), Netherlands (KEMA), Norway (NCS), Poland (PCBC), Portugal (APCER), Singapore (PSB), Slovenia (SIQ), Spain (AENOR), Switzerland (SQS)

10.2 ISO/IEC 17025 - 1999 "General Requirements for the Competence of Testing and Calibration"

- Chemical Testing - Accredited A2LA Certificate Number 883.01



10.3 ISO/IEC Guide 34 - 2000 "General Requirements for the Competence of Reference Material Producers"

- Reference Materials Production - Accredited A2LA Certificate Number 883.02

A2LA Mutual Recognition Agreement Partners:

Australia (NATA), Austria (Bmwa), Belgium (BELTEST) (BKO-OBE), Canada (SCC), Chinese Taipei (CNLA), Czech Republic (NAO), Denmark (DANAK), Finland (FINAS), France (COFRAC), Germany (DAR), Hong Kong (HKAS), Ireland (NAB), Italy (SIT) (SINAL), Japan (JAB) (JNLA), Republic of Korea (KOLAS), The Netherlands (RvA), New Zealand (IANZ), Norway (NA), Portugal (IPQ), Singapore (SAC-SINGLAS), Spain (ENAC), Sweden (SWEDAC), Switzerland (SAS), United Kingdom (UKAS) and United States (NVLAP) (ICBO ES)

10.4 10CFR50 Appendix B - Nuclear Regulatory Commission

- Domestic Licensing of Production and Utilization Facilities

10.5 10CFR21 - Nuclear Regulatory Commission - Reporting Defects and Non-Compliance

10.6 MIL-STD-45662A (Obsolete/Observed)

INORGANIC LABS/RADCHEM LABS Pg. 2 of 2

DATE RECEIVED: 06/20/03


DATE EXPIRED: 07/01/2004

DATE OPENED: 06/23/03

INORG: 4153 PO: F52370

11.0 DATE OF CERTIFICATION AND PERIOD OF VALIDITY

010251



11.1 IV Shelf Life - The period of time during which the concentration of the analyte(s) in a properly packaged, unopened, and unused standard stored under environmentally controlled and monitored conditions will remain within the specified uncertainty range. Shelf life is limited primarily by transpiration (loss of water from the solution) and infrequently, by chemical instability. Transpiration studies (P-SP01020) of chemically-stable solutions performed at Inorganic Ventures / IV Labs indicate a CRM shelf-life of four years for solutions packaged in 500-mL low density polyethylene bottles. When stored under special conditions that minimize transpiration and instability, the shelf life can be extended past this limit.

11.2 Expiration Date - The date after which a CRM should not be used. Routine laboratory use of a CRM increases transpiration losses and the chance of contamination which affect the integrity of the CRM and limit its useful life. Inorganic Ventures / IV Labs concurs with state and federal regulatory agencies' recommendations that solution standards be assigned a one-year expiration date.

Certification Date: February 24, 2003

Expiration Date:

EXPIRES
01/22/04

12.0 NAMES AND SIGNATURES OF CERTIFYING OFFICERS

Certificate Prepared By: JoAnn Struthers, QA Administrative Assistant

Certificate Approved By: Katalin Le, QC Supervisor

Certifying Officer: Paul Gaines, Chemist, Senior Technical Director

JoAnn Struthers

Katalin Le

Paul Gaines



Certificate of Analysis

CUSTOM-GRADE SOLUTION

Catalog Number: CGTI1-1, CGTI1-2 and CGTI1-5

1000 µg/mL Titanium in 1.4% HNO₃ (abs) tr. H₂O **010252**

This standard should not be prepared or stored in glass.

Lot Number: **T-TI02039**

Starting Material:
Starting Material Purity:
Starting Material Lot No:

Titanium Metal
99.999%
F29I14

INORGANIC LABS/RADCHEM LABS

DATE RECEIVED: 11/12/02
DATE EXPIRED: 12/11/2003 *DR*
DATE OPENED: 11/13/02
INORG: 3962 PO: F52061

CERTIFIED CONCENTRATION: 1010 ± 3 µg/mL

The Certified Value is based upon the most precise method used to analyze this CRM. The following equations are used in the calculation of the certified value and the uncertainty:

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

$$\text{Uncertainty } (\pm) = \frac{2[(\sum s_i)^2]^{1/2}}{(n)^{1/2}}$$

($\bar{x}$) = mean

x_i = individual results

n = number of measurements

$\sum s_i$ = The summation of all significant estimated errors.

Calculated Value: **1002 µg/mL**

Method: Calculated, based on starting material.

Instrument Analysis: **1010 ± 3 µg/mL** (Average of 3 runs)

Method: Inductively Coupled Plasma Spectroscopy (ICP) vs NIST SRM 3162a.

TRACE METALLIC IMPURITIES DETERMINED BY ICP-MS AND ICP-OES IN µg/mL:

Custom-Grade solutions tested for trace metallic impurities by ICP-MS were analyzed in an ULPA-Filtered Clean Room.

An ULPA-Filter is 99.9985% efficient for the removal of particles down to 0.3 µm.

<u>Q</u> Al <0.010	<u>M</u> Dy <0.00060	<u>M</u> Li <0.0010	<u>M</u> Pr <0.000030	<u>M</u> Te <0.0030
<u>M</u> Sb <0.000050	<u>M</u> Er <0.00050	<u>M</u> Lu <0.000040	<u>M</u> Re <0.00010	<u>M</u> Tb <0.000030
<u>M</u> As <0.0010	<u>M</u> Eu <0.00030	<u>Q</u> Mg <0.020	<u>M</u> Rh <0.00010	<u>M</u> Tl <0.00010
<u>M</u> Ba <0.0010	<u>M</u> Gd <0.00010	<u>M</u> Mn 0.0020	<u>M</u> Rb <0.00010	<u>M</u> Th <0.00010
<u>M</u> Be <0.000050	<u>M</u> Ga <0.00010	<u>Q</u> Hg <0.050	<u>M</u> Ru <0.00020	<u>M</u> Tm <0.000040
<u>M</u> Bi <0.000040	<u>M</u> Ge <0.00060	<u>M</u> Mo <0.00020	<u>M</u> Sm <0.00010	<u>M</u> Sn <0.00050
<u>Q</u> B <0.050	<u>M</u> Au <0.00030	<u>M</u> Nd <0.00020	<u>Q</u> Sc <0.0020	<u>s</u> Ti
<u>Q</u> Cd <0.00030	<u>M</u> Hf <0.00020	<u>Q</u> Ni <0.050	<u>Q</u> Se <0.40	<u>M</u> W <0.0010
<u>Q</u> Ca <0.010	<u>M</u> Ho <0.000050	<u>Q</u> Nb <0.000050	<u>Q</u> Si <0.010	<u>M</u> U <0.00020
<u>M</u> Ce <0.00050	<u>Q</u> In <0.020	<u>n</u> Os	<u>M</u> Ag <0.00020	<u>M</u> V <0.00020
<u>M</u> Cs <0.000030	<u>M</u> Ir <0.00050	<u>M</u> Pd <0.00050	<u>Q</u> Na 0.12	<u>M</u> Yb <0.00010
<u>M</u> Cr <0.00050	<u>Q</u> Fe <0.010	<u>i</u> P	<u>M</u> Sr <0.000050	<u>M</u> Y <0.0040
<u>M</u> Co <0.00030	<u>M</u> La <0.000050	<u>M</u> Pt <0.00020	<u>i</u> S	<u>M</u> Zn 0.19
<u>Q</u> Cu <0.040	<u>M</u> Pb <0.00030	<u>n</u> K 0.23	<u>M</u> Ta <0.00070	<u>M</u> Zr <0.00050

M - checked by ICP-MS

O - checked by ICP-OES

i - spectral interference

n - not checked for

s - solution standard element

ANALYZED DENSITY OF SOLUTION (measured at 22°C): **1.011 g/mL**

QA:KL Rev.060502DM

(over)

Inorganic Ventures, Inc.

195 Lehigh Avenue • Suite 4 • Lakewood, NJ 08701

Orders: 800-669-6799 • FAX (732) 901-1903

Technical Support: 800-569-6799

Quality Assurance Manager

EXPIRES
12/2003



inorganic ventures / iv labs

195 lehigh avenue, suite 4, lakewood, nj 08701 usa
phone: 800-669-6799 • 732-901-1900 • fax: 732-901-1903
e-mail: ivsales@ivstandards.com • website: www.ivstandards.com

010253

certificate of analysis

1.0 Inorganic Ventures / IV Labs is an ISO Guide 34-2000 Certified Reference Material (CRM) Manufacturer: Certificate #883-02. The certificate is designed and the certified value(s) and uncertainty(ies) are determined in accordance with ISO Guide 31-2000 (Reference Materials - Contents of certificates and label(s)), ISO Guide 34-2000 "Quality System Guidelines for the Production of Reference Materials," and ISO Guide 35-1989 "Certification of Reference Materials - General and Statistical Principles."

2.0 DESCRIPTION OF CRM Custom-Grade 1000 µg/mL Strontium in 0.1% (abs) HNO₃

Catalog Number: CGSR1-1, CGSR1-2, and CGSR1-5

Lot Number: T-SR01123

Starting Material: SrCO₃

Starting Material Purity (%): 99.9951

Starting Material Lot No 22593

Matrix: 0.1% (abs) HNO₃

INORGANIC LABS/RADCHEM LABS Pg. 1 of 2
DATE RECEIVED: 06/20/03
DATE EXPIRED: 07/01/2004
DATE OPENED: 06/23/03
INORG: 4154 PO: F52370

3.0 CERTIFIED VALUES AND UNCERTAINTIES

Certified Concentration: 998 ± 2 µg/mL

Certified Density: 1.000 g/mL (measured at 22° C)

The Certified Value is based upon the most precise method used to analyze this CRM. The following equations are used in the calculation of the certified value and the uncertainty:

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

($\bar{x}$) = mean

x_i = individual results

n = number of measurements

$$\text{Uncertainty } (\pm) = \frac{2[(\sum s_i^2)^{1/2}]}{(n)^{1/2}}$$

$\sum S$ = The summation of all significant estimated errors.

(Most common are the errors from instrumental measurement, weighing, dilution to volume, and the fixed error reported on the NIST SRM certificate of analysis.)

The independent samples t-test was used to determine if there is agreement between the above assay methods at the 95% confidence interval. Both methods were compared and showed agreement within the stated uncertainties. This agreement is a confirmation of the accuracy of this CRM.

4.0 TRACEABILITY TO NIST AND VALUES OBTAINED BY INDEPENDENT METHODS

• "Property of the result of a measurement or the value of a standard whereby it can be related to stated references, usually national or international standards, through an unbroken chain of comparisons all having stated uncertainties." (ISO VIM, 2nd ed., 1993, definition 6.10)

• This IV product is Traceable to NIST via direct comparison to NIST SRMs. The uncertainties for each certified value are reported, taking into account the SRM uncertainty error and the measurement, weighing and volume dilution errors.

4.1 Assay Method #1 998 ± 2 µg/mL

EDTA NIST SRM 928 Lot Number: 880710

Assay Method #2 1002 ± 8 µg/mL

ICP Assay NIST SRM 3153a Lot Number: 990906

4.2 BALANCE CALIBRATION - All balances are checked daily using in-house procedure number 6-IMM-001. The weights used for testing are annually compared to Gerhart Scale Corporation's master weights and are traceable to the National Institute of Standards and Technology (NIST). The NIST Traceability numbers are 692476 - Class 1 and 692476A - Class 2. The NIST test number is 822/260017-98. All analytical balances are calibrated every 4 months by Gerhart Scale Corp. of South Amboy. The balances are calibrated with a class 1 and/or class 2 analytical weight set. These weights are tested annually by a NIST / NVLAP accredited calibration lab. The NIST test number is 822/260017-98.

010254

4.3 THERMOMETER CALIBRATION - The thermometers used in the determination of the final densities are calibrated vs standard thermometer No. 903-2680 which was certified in accordance with the procedures outlined by ASTM E77-87 and NIST Monograph 150 using NIST Test Nos. and Std Nos.: 769543, 217368/769543, 217368/P14452, 176240/P14452, 176240. The in-house procedure No. is 2-QC-001. Thermometers which are not calibrated vs standard thermometer No. 903-2680 are traceable to NIST Identification Nos. 92564, 119016, 471047 and NIST test report Nos. 811/258522, 811/2557078, and 236090.

4.4 GLASSWARE CALIBRATION - In-house procedure 3-QC-002 is used to calibrate all Class A Glassware used in the manufacture and quality control of Custom Grade Standards.

5.0 TRACE METALLIC IMPURITIES (TMI) DETERMINED BY ICP/MS AND ICP-OES IN µg/mL

Custom-Grade solutions are tested for trace metallic impurities by Axial ICP-OES and ICP-MS. The result from the most sensitive method for each element, is reported below. Solutions tested by ICP-MS were analyzed in an ULPA-Filtered Clean Room. An ULPA-Filter is 99.9985% efficient for the removal of particles down to 0.3 µm.

<u>Q</u> Al < 0.00090	<u>M</u> Dy < 0.00600	<u>Q</u> Li < 0.00003	<u>M</u> Pr < 0.00030	<u>Q</u> Te < 0.10000
<u>M</u> Sb < 0.00050	<u>M</u> Er < 0.00500	<u>M</u> Lu < 0.00040	<u>M</u> Re < 0.00100	<u>M</u> Tb < 0.00030
<u>Q</u> As < 0.00500	<u>M</u> Eu < 0.00300	<u>Q</u> Mg 0.00037	<u>Q</u> Rh < 0.00600	<u>M</u> Tl < 0.00100
<u>M</u> Ba 0.04001	<u>M</u> Gd < 0.00100	<u>Q</u> Mn 0.00018	<u>i</u> Rb	<u>M</u> Th < 0.00100
<u>Q</u> Be < 0.00009	<u>M</u> Ga < 0.00100	<u>Q</u> Hg < 0.01500	<u>Q</u> Ru < 0.00300	<u>M</u> Tm < 0.00040
<u>M</u> Bi < 0.00040	<u>M</u> Ge < 0.00600	<u>M</u> Mo < 0.00200	<u>M</u> Sm < 0.00100	<u>M</u> Sn < 0.00500
<u>Q</u> B < 0.00060	<u>M</u> Au < 0.00300	<u>M</u> Nd < 0.00200	<u>M</u> Sc < 0.01000	<u>M</u> Ti < 0.05001
<u>M</u> Cd < 0.00300	<u>M</u> Hf < 0.00200	<u>Q</u> Ni < 0.00300	<u>Q</u> Se < 0.05000	<u>M</u> W < 0.01000
<u>Q</u> Ca 0.03600	<u>M</u> Ho < 0.00050	<u>M</u> Nb < 0.00050	<u>Q</u> Si 0.00056	<u>M</u> U < 0.00200
<u>M</u> Ce < 0.00500	<u>Q</u> In < 0.00200	<u>n</u> Os	<u>M</u> Ag < 0.00200	<u>M</u> V < 0.00200
<u>M</u> Cs < 0.00030	<u>M</u> Ir < 0.00500	<u>Q</u> Pd < 0.00400	<u>Q</u> Na 0.00520	<u>M</u> Yb < 0.00100
<u>Q</u> Cr < 0.00080	<u>Q</u> Fe 0.00080	<u>Q</u> P < 0.00480	<u>s</u> Sr	<u>Q</u> Y < 0.00004
<u>M</u> Co < 0.00300	<u>M</u> La < 0.00050	<u>M</u> Pt < 0.00200	<u>n</u> S	<u>M</u> Zn < 0.02000
<u>Q</u> Cu < 0.00140	<u>M</u> Pb < 0.00300	<u>Q</u> K < 0.00170	<u>M</u> Ta < 0.00700	<u>M</u> Zr < 0.00500

M - Checked by ICP-MS O - Checked by ICP-OES i - Spectral Interference n - Not Checked For s - Solution Standard Element

6.0 INTENDED USE

For the calibration of analytical instruments including but not limited to the following:
 ICP-MS, ICP-OES, FAAS, GFAA, XRF, and DCP
 For the validation of analytical methods
 For the preparation of "working reference samples"
 For interference studies and the determination of correction coefficients
 For detection limit and linearity studies
 For additional intended uses, contact IV Technical Staff

7.0 INSTRUCTIONS FOR THE CORRECT USE OF THIS REFERENCE MATERIAL

Storage & Handling - Keep tightly sealed when not in use. Store and use at 20 ± 4°C. Do not pipet from container. Do not return portions removed for pipetting to container.

Atomic Weight; Valence; Coordination Number; Chemical Form in Solution - 87.62; +2; 6; Sr(H₂O)₈·2

Chemical Compatibility - Soluble in HCl, and HNO₃. Avoid H₂SO₄, HF and neutral to basic media. Stable with most metals and inorganic anions forming insoluble silicate, carbonate, hydroxide, oxide, fluoride, sulfate, oxalate, chromate, arsenate and tungstate in neutral aqueous media.

Stability - 2-100 ppb levels stable for months in 1% HNO₃ / LDPE container. 1-10,000 ppm solutions chemically stable for years in 1 - 3.5% HNO₃ / LDPE container.

Sr Containing Samples (Preparation and Solution) - Metal (Best dissolved in diluted HNO₃); Ores (Carbonate fusion in Pt⁰ followed by HCl dissolution); Organic Matrices (Dry ash and dissolution in dilute HCl).

Atomic Spectroscopic Information (ICP-OES D.L.s are given as radial/axial view):

Technique/Line	Estimated D.L.	Order	Type	Interferences (underlined indicates severe at ≈ concs.)
ICP-OES 407.771 nm	0.0004 / 0.00006 µg/mL	1	ion	U, Ce
ICP-OES 421.552 nm	0.0008 / 0.00004 µg/mL	1	ion	Rb
ICP-OES 460.733 nm	0.07 / 0.003 µg/mL	1	atom	Ce
ICP-MS 88 amu	1200 ppt	n/a	M ⁺	⁷² Ge ¹⁸ O, ¹⁷⁶ Yb ² , ¹⁷⁶ Lu ² , ¹⁷⁶ Hf ²

8.0 **HAZARDOUS INFORMATION** - Please refer to the enclosed Material Safety Data sheet for information regarding this CRM.

9.0 **HOMOGENEITY** - This solution was mixed according to procedure IV-MPM-004 and is guaranteed to be homogeneous.

10.0 QUALITY STANDARD DOCUMENTATION

010255

10.1 ISO 9001:2000 Quality Management System Registration - QMI Certificate Number 010105

Recognized by:

Registrar Accreditation Board (ANSI-RAB)
Standards Council of Canada (SCC)
Dutch Council for Accreditation (RVA)
Entidad Mexicana de Acreditacion, a.c.(EMA)



Members of IQ Net International Certification Network:

Argentina (IRAM), Australia (QAS), Austria (ÖQS), Belgium (Avinter), Brazil (FCAV), Canada (QMI), Hong Kong (HKQAA), Columbia (ICONTEC), Czech Republic (CQS), Denmark (DS), Finland (SFS), France (AFAQ), Germany (DQS), Greece (ELOT), Hungary (MSZT), Ireland (NSAI), Israel (SII), Italy (CISQ), Japan (JQA), Korea (KSA-QA), Netherlands (KEMA), Norway (NCS), Poland(PCBC), Portugal (APCER), Singapore (PSB), Slovenia (SIQ), Spain (AENOR), Switzerland (SQS)

10.2 ISO/IEC 17025 - 1999 "General Requirements for the Competence of Testing and Calibration" - Chemical Testing - Accredited A2LA Certificate Number 883.01



10.3 ISO/IEC Guide 34 - 2000 "General Requirements for the Competence of Reference Material Producers" - Reference Materials Production - Accredited A2LA Certificate Number 883.02

A2LA Mutual Recognition Agreement Partners:

Australia (NATA), Austria (Bmwa), Belgium (BELTEST) (BKO-OBE), Canada (SCC), Chinese Taipei (CNLA), Czech Republic (NAO), Denmark (DANAK), Finland (FINAS), France (COFRAC), Germany (DAR), Hong Kong (HKAS), Ireland (NAB), Italy (SIT) (SINAL), Japan (JAB) (JNLA), Republic of Korea (KOLAS), The Netherlands (RvA), New Zealand (IANZ), Norway (NA), Portugal (IPQ), Singapore (SAC-SINGLAS), Spain (ENAC), Sweden (SWEDAC), Switzerland (SAS), United Kingdom (UKAS) and United States (NVLAP) (ICBO ES)

10.4 10CFR50 Appendix B - Nuclear Regulatory Commission - Domestic Licensing of Production and Utilization Facilities

10.5 10CFR21 - Nuclear Regulatory Commission - Reporting Defects and Non-Compliance

10.6 MIL-STD-45662A (Obsolete/Observed)

11.0 DATE OF CERTIFICATION AND PERIOD OF VALIDITY



11.1 **IV Shelf Life** - The period of time during which the concentration of the analyte(s) in a properly packaged, unopened, and unused standard stored under environmentally controlled and monitored conditions will remain within the specified uncertainty range. Shelf life is limited primarily by transpiration (loss of water from the solution) and infrequently, by chemical instability. Transpiration studies (P-SP01020) of chemically-stable solutions performed at Inorganic Ventures / IV Labs indicate a CRM shelf-life of four years for solutions packaged in 500-mL low density polyethylene bottles. When stored under special conditions that minimize transpiration and instability, the shelf life can be extended past this limit.

11.2 **Expiration Date** - The date after which a CRM should not be used. Routine laboratory use of a CRM increases transpiration losses and the chance of contamination which affect the integrity of the CRM and limit its useful life. Inorganic Ventures / IV Labs concurs with state and federal regulatory agencies' recommendations that solution standards be assigned a one-year expiration date.

Certification Date: October 31, 2002

Expiration Date: **EXPIRES**

01/22/04

INORGANIC LABS/RADCHEM LABS 79-2 of 2
DATE RECEIVED: 06/20/03
DATE EXPIRED: 07/01/2004
DATE OPENED: 06/23/03
INORG: 4154 PO: 752370

12.0 NAMES AND SIGNATURES OF CERTIFYING OFFICERS

Certificate Prepared By: Debbie Newman, QA Administrator

Certificate Approved By: Katalin Le, QC Supervisor

Certifying Officer: Paul Gaines, Chemist, Senior Technical Director

010256

Debbie Newman
Katalin Le

Paul Gaines



inorganic ventures / iv labs

195 lehigh avenue, suite 4, lakewood, nj 08701 usa
phone: 800-669-6799 • 732-901-1900 • fax: 732-901-1903
e-mail: ivsales@ivstandards.com • website: www.ivstandards.com

010257 certificate of analysis

1.0 Inorganic Ventures / IV Labs is an ISO Guide 34-2000 Certified Reference Material (CRM) Manufacturer: Certificate #883-02. The certificate is designed and the certified value(s) and uncertainty(ies) are determined in accordance with ISO Guide 31-2000 (Reference Materials - Contents of certificates and label(s)), ISO Guide 34-2000 "Quality System Guidelines for the Production of Reference Materials," and ISO Guide 35-1989 "Certification of Reference Materials - General and Statistical Principles."

2.0 DESCRIPTION OF CRM Custom-Grade 1000 µg/mL Tin In H₂O tr HNO₃ tr HF

Catalog Number: CGSN1-1, CGSN1-2, and CGSN1-5
Lot Number: T-SN01110
Starting Material: #1, #2 SN SHOT
Starting Material Purity (%): 99.9969
Starting Material Lot No AAJM H12H17, A23F09
Matrix: H₂O tr HNO₃ tr HF

INORGANIC LABS/RADCHEM LABS Pg 1 of 2
DATE RECEIVED: 3/26/03
DATE EXPIRED: 4/01/2004
DATE OPENED: 3/28/03
INORG: 4026 PD: F52167

3.0 CERTIFIED VALUES AND UNCERTAINTIES

Certified Concentration: 1003 ± 5 µg/mL

Certified Density: 1.000 g/mL (measured at 22° C)

The Certified Value is based upon the most precise method used to analyze this CRM. The following equations are used in the calculation of the certified value and the uncertainty:

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

$\bar{x}$ = mean

x_i = individual results

n = number of measurements

$$\text{Uncertainty } (\pm) = \frac{2[(\sum s_i)^2]^{1/2}}{(n)^{1/2}}$$

$\sum S$ = The summation of all significant estimated errors

(Most common are the errors from instrumental measurement, weighing, dilution to volume, and the fixed error reported on the NIST SRM certificate of analysis.)

4.0 TRACEABILITY TO NIST AND VALUES OBTAINED BY INDEPENDENT METHODS

☐ "Property of the result of a measurement or the value of a standard whereby it can be related to stated references, usually national or international standards, through an unbroken chain of comparisons all having stated uncertainties." (ISO VIM, 2nd ed., 1993, definition 6.10)

☐ This IV product is Traceable to NIST via direct comparison to NIST SRMs. The uncertainties for each certified value are reported, taking into account the SRM uncertainty error and the measurement, weighing and volume dilution errors.

4.1 Assay Method #1 1003 ± 5 µg/mL (Avg 2 runs)

ICP Assay NIST SRM 3161a Lot Number: 993107

Assay Method #2 1000 µg/mL

Gravimetric NIST SRM Lot Number: See Sec. 4.2

4.2 BALANCE CALIBRATION - All balances are checked daily using in-house procedure number 6-IMM-001. The weights used for testing are annually compared to Gerhart Scale Corporation's master weights and are traceable to the National Institute of Standards and Technology (NIST). The NIST Traceability numbers are 692476 - Class 1 and 692476A - Class 2. The NIST test number is 822/260017-98. All analytical balances are calibrated every 4 months by Gerhart Scale Corp. of South Amboy. The balances are calibrated with a class 1 and/or class 2 analytical weight set. These weights are tested annually by a NIST / NVLAP accredited calibration lab. The NIST test number is 822/260017-98.

010258

4.3 THERMOMETER CALIBRATION - The thermometers used in the determination of the final densities are calibrated vs standard thermometer No. 903-2680 which was certified in accordance with the procedures outlined by ASTM E77-87 and NIST Monograph 150 using NIST Test Nos. and Std Nos.: 769543, 217368/769543, 217368/P14452, 176240/P14452, 176240. The in-house procedure No. is 2-QC-001. Thermometers which are not calibrated vs standard thermometer No. 903-2680 are traceable to NIST Identification Nos. 92564, 119016, 471047 and NIST test report Nos. 811/258522, 811/2557078, and 236090.

4.4 GLASSWARE CALIBRATION - In-house procedure 3-QC-002 is used to calibrate all Class A Glassware used in the manufacture and quality control of Custom Grade Standards.

5.0 TRACE METALLIC IMPURITIES (TMI) DETERMINED BY ICP/MS AND ICP-OES IN µg/mL

Custom-Grade solutions are tested for trace metallic impurities by Axial ICP-OES and ICP-MS. The result from the most sensitive method for each element, is reported below. Solutions tested by ICP-MS were analyzed in an ULPA-Filtered Clean Room. An ULPA-Filter is 99.9985% efficient for the removal of particles down to 0.3 µm.

<u>Q</u> Al 0.00100	<u>M</u> Dy < 0.01217	<u>Q</u> Li < 0.00002	<u>M</u> Pr < 0.00061	<u>M</u> Te < 0.06087
<u>Q</u> Sb < 0.01000	<u>M</u> Er < 0.01014	<u>M</u> Lu < 0.00081	<u>M</u> Re < 0.00203	<u>M</u> Tb < 0.00061
<u>M</u> As < 0.02029	<u>M</u> Eu < 0.00609	<u>Q</u> Mg < 0.00003	<u>M</u> Rh < 0.00203	<u>M</u> Tl < 0.00203
<u>Q</u> Ba < 0.00070	<u>M</u> Gd < 0.00203	<u>M</u> Mn < 0.00812	<u>M</u> Rb < 0.00203	<u>M</u> Th < 0.00203
<u>M</u> Be < 0.00101	<u>M</u> Ga < 0.00203	<u>Q</u> Hg < 0.01500	<u>M</u> Ru < 0.00406	<u>M</u> Tm < 0.00081
<u>M</u> Bi < 0.00081	<u>M</u> Ge < 0.01217	<u>M</u> Mo < 0.00406	<u>M</u> Sm < 0.00203	<u>S</u> Sn
<u>Q</u> B < 0.01200	<u>M</u> Au < 0.00609	<u>M</u> Nd < 0.00406	<u>M</u> Sc < 0.02029	<u>M</u> Ti < 0.10144
<u>Q</u> Cd < 0.00500	<u>M</u> Hf < 0.00406	<u>Q</u> Ni < 0.01000	<u>M</u> Se < 0.01623	<u>M</u> W < 0.02029
<u>Q</u> Ca < 0.00150	<u>M</u> Ho < 0.00101	<u>M</u> Nb < 0.00101	<u>Q</u> Si 0.02700	<u>M</u> U < 0.00406
<u>M</u> Ce < 0.01014	<u>M</u> In < 0.02029	<u>n</u> Os	<u>M</u> Ag < 0.00406	<u>M</u> V < 0.00406
<u>M</u> Cs < 0.00061	<u>M</u> Ir < 0.01014	<u>M</u> Pd < 0.01014	<u>Q</u> Na 0.00325	<u>M</u> Yb < 0.00203
<u>M</u> Cr < 0.01014	<u>Q</u> Fe < 0.00110	<u>Q</u> P < 0.00500	<u>M</u> Sr < 0.00101	<u>M</u> Y < 0.08115
<u>Q</u> Co < 0.00200	<u>M</u> La < 0.00101	<u>M</u> Pt < 0.00406	<u>N</u> S	<u>M</u> Zn < 0.04058
<u>M</u> Cu < 0.01217	<u>M</u> Pb < 0.00609	<u>Q</u> K < 0.00200	<u>M</u> Ta < 0.01420	<u>M</u> Zr < 0.01014

M - Checked by ICP-MS O - Checked by ICP-OES I - Spectral Interference n - Not Checked For s - Solution Standard Element

6.0 INTENDED USE

For the calibration of analytical instruments including but not limited to the following:
 ICP-MS, ICP-OES, FAAS, GFAA, XRF, and DCP
 For the validation of analytical methods
 For the preparation of "working reference samples"
 For interference studies and the determination of correction coefficients
 For detection limit and linearity studies
 For additional intended uses, contact IV Technical Staff

7.0 INSTRUCTIONS FOR THE CORRECT USE OF THIS REFERENCE MATERIAL

010259

Storage & Handling - Keep tightly sealed when not in use. Store and use at $20 \pm 4^\circ\text{C}$. Do not pipet from container. Do not return portions removed for pipetting to container.

Atomic Weight; Valence; Coordination Number; Chemical Form in Solution - 118.710; +4; 4,5, 6,7,8 $\text{Sn}(\text{OH})_4\text{F}_2$

Chemical Compatibility - Soluble in HCl and dilute HF / HNO_3 . Avoid neutral to basic media. Unstable at ppm levels with metals that would pull F⁻ away. (i.e. Do not mix with Alkaline or Rare Earths or high levels of transition elements unless they are fluorinated.) Stable with most inorganic anions provided it is in the chemical form shown above.

Stability - 2-100 ppb levels stable (alone or mixed with all other metals that are at comparable levels) as the $\text{Sn}(\text{OH})_4\text{F}_2$ for 1 year in 1% HNO_3 / LDPE container. 1-10,000 ppm single element solutions as the $\text{Sn}(\text{OH})_4\text{F}_2$ chemically stable for years in 2-5% HNO_3 / trace HF in a LDPE container.

Sn Containing Samples (Preparation and Solution) - Metal (Soluble in HF / HNO_3 or HCl); Oxides - SnO (soluble in HCl), SnO_2 - very resistant to all acids including HF (Fusion with equal parts of Na_2CO_3 and S. It is then soluble in water or dilute acids as the stannate.); Alloys (Treat first 0.1 g with 10 mL conc. H_2SO_4 to boiling until the alloy disintegrates and nearly all of the sulfuric acid is expelled. Then add 100 mL H_2O , free water and 50 mL of conc HCl or transfer to a plastic container and add 1 mL HF in either case warming gently to bring about solution.) Organic Matrices (Volatility and precipitation of the insoluble stannic oxide are problems. Consultation of the literature should be made for individual matrices / Sn compounds.)

Atomic Spectroscopic Information (ICP-OES D.L.s are given as radial/axial view):

Technique/Line	Estimated D.L.	Order	Type	Interferences (underlined indicates severe at 100ppb)
ICP-OES 189.989 nm	0.03 / 0.003 $\mu\text{g/mL}$	1	ion	
ICP-OES 242.949 nm	0.1 / 0.01 $\mu\text{g/mL}$	1	atom	W, Mo, Rh, Ta, Co
ICP-MS 120 amu	5 ppt	n/a	M ⁺	¹¹⁸ Te, ¹¹⁸ Ru, ¹¹⁸ O, ¹¹⁸ Pd, ¹¹⁸ Os

HF Note: This standard should not be prepared or stored in glass.

8.0 HAZARDOUS INFORMATION - Please refer to the enclosed Material Safety Data sheet for information regarding this CRM.

9.0 HOMOGENEITY - This solution was mixed according to procedure IV-MPM-004 and is guaranteed to be homogeneous.

10.0 QUALITY STANDARD DOCUMENTATION

10.1 ISO 9001:2000 Quality Management System Registration - QMI Certificate Number 010105

Recognized by:

Registrar Accreditation Board (ANSI-RAB)

Standards Council of Canada (SCC)

Dutch Council for Accreditation (RVA)

Entidad Mexicana de Acreditacion, a.c.(EMA)

Members of IQ Net International Certification Network:

Argentina (IRAM), Australia (QAS), Austria (ÖQS), Belgium (Avinter), Brazil (FCAV), Canada (QMI), Hong Kong (HKQAA), Columbia (ICONTEC), Czech Republic (CQS), Denmark (DS), Finland (SFS), France (AFAQ), Germany (DQS), Greece (ELOT), Hungary (MSZT), Ireland (NSAI), Israel (SII), Italy (CISQ), Japan (JQA), Korea (KSA-QA), Netherlands (KEMA), Norway (NCS), Poland (PCBC), Portugal (APCER), Singapore (PSB), Slovenia (SIQ), Spain (AENOR), Switzerland (SQS)

10.2 ISO/IEC 17025 - 1999 "General Requirements for the Competence of Testing and Calibration"

- Chemical Testing - Accredited A2LA Certificate Number 883.01



10.3 ISO/IEC Guide 34 - 2000 "General Requirements for the Competence of Reference Material Producers"

- Reference Materials Production - Accredited A2LA Certificate Number 883.02

A2LA Mutual Recognition Agreement Partners:

Australia (NATA), Austria (BmWA), Belgium (BELTEST) (BKO-OBE), Canada (SCC), Chinese Taipei (CNLA), Czech Republic (NAO), Denmark (DANAK), Finland (FINAS), France (COFRAC), Germany (DAR), Hong Kong (HKAS), Ireland (NAB), Italy (SIT) (SINAL), Japan (JAB) (JNLA), Republic of Korea (KOLAS), The Netherlands (RvA), New Zealand (IANZ), Norway (NA), Portugal (IPQ), Singapore (SAC-SINGLAS), Spain (ENAC), Sweden (SWEDAC), Switzerland (SAS), United Kingdom (UKAS) and United States (NVLAP) (ICBO ES)

10.4 10CFR50 Appendix B - Nuclear Regulatory Commission

- Domestic Licensing of Production and Utilization Facilities

10.5 10CFR21 - Nuclear Regulatory Commission - Reporting Defects and Non-Compliance

10.6 MIL-STD-45662A (Obsolete/Observed)

INORGANIC LABS/RADCHEM LABS Pg. 2 of 2

DATE RECEIVED: 3/26/03

DATE EXPIRED: 4/01/2004

DATE OPENED: 3/28/03

INORG: 4026 PO: F52167

11.0 DATE OF CERTIFICATION AND PERIOD OF VALIDITY

010260



11.1 IV Shelf Life - The period of time during which the concentration of the analyte(s) in a properly packaged, unopened, and unused standard stored under environmentally controlled and monitored conditions will remain within the specified uncertainty range. Shelf life is limited primarily by transpiration (loss of water from the solution) and infrequently, by chemical instability. Transpiration studies (P-SP01020) of chemically-stable solutions performed at Inorganic Ventures / IV Labs indicate a CRM shelf-life of four years for solutions packaged in 500-mL low density polyethylene bottles. When stored under special conditions that minimize transpiration and instability, the shelf life can be extended past this limit.

11.2 Expiration Date - The date after which a CRM should not be used. Routine laboratory use of a CRM increases transpiration losses and the chance of contamination which affect the integrity of the CRM and limit its useful life. Inorganic Ventures / IV Labs concurs with state and federal regulatory agencies' recommendations that solution standards be assigned a one-year expiration date.

Certification Date: December 17, 2002

Expiration Date:

EXPIRES

12/2004


12.0 NAMES AND SIGNATURES OF CERTIFYING OFFICERS

Certificate Prepared By: Debbie Newman, QA Administrator

Certificate Approved By: Katalin Le, QC Supervisor

Certifying Officer: Paul Gaines, Chemist, Senior Technical Director

Debbie Newman
Katalin Le
Paul Gaines



inorganic ventures / iv labs

195 lehigh avenue, suite 4, lakewood, nj 08701 usa
phone: 800-669-6799 • 732-901-1900 • fax: 732-901-1903
e-mail: ivsales@ivstandards.com • website: www.ivstandards.com

certificate of analysis

010261

- 1.0 Inorganic Ventures / IV Labs is an ISO Guide 34-2000 Certified Reference Material (CRM) Manufacturer: Certificate #883-02.** The certificate is designed and the certified value(s) and uncertainty(ies) are determined in accordance with ISO Guide 31-2000 (Reference Materials - Contents of certificates and label(s)), ISO Guide 34-2000 "Quality System Guidelines for the Production of Reference Materials," and ISO Guide 35-1989 "Certification of Reference Materials - General and Statistical Principles."

2.0 DESCRIPTION OF CRM Custom-Grade 1000 µg/mL Bismuth in 3.5% (abs) HNO₃

Catalog Number: CGBI1-1 and CGBI1-5
Lot Number: **W-BI01089**
Starting Material: Bi needles
Starting Material Purity (%): 99.999090
Starting Material Lot No G25L16
Matrix: 3.5% (abs) HNO₃

INORGANIC LABS/RADCHEM LABS Pg. 1 of 2
DATE RECEIVED: 07/31/03
DATE EXPIRED: 08/01/2004
DATE OPENED: 08/01/03
INORG: 4200 PO: F52383

3.0 CERTIFIED VALUES AND UNCERTAINTIES

Certified Concentration: 1002 ± 4 µg/mL

Certified Density: 1.026 g/mL (measured at 22° C)

The Certified Value is based upon the most precise method used to analyze this CRM. The following equations are used in the calculation of the certified value and the uncertainty:

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

($\bar{x}$) = mean

x_i = individual results

n = number of measurements

$$\text{Uncertainty } (\pm) = \frac{2[(\sum s_i^2)^{1/2}]}{(n)^{1/2}}$$

$\sum s_i^2$ = The summation of all significant estimated errors

(Most common are the errors from instrumental measurement, weighing, dilution to volume, and the fixed error reported on the NIST SRM certificate of analysis.)

4.0 TRACEABILITY TO NIST AND VALUES OBTAINED BY INDEPENDENT METHODS

□ "Property of the result of a measurement or the value of a standard whereby it can be related to stated references, usually national or international standards, through an unbroken chain of comparisons all having stated uncertainties." (ISO VIM, 2nd ed., 1993, definition 6.10)

□ This IV product is Traceable to NIST via direct comparison to NIST SRMs. The uncertainties for each certified value are reported, taking into account the SRM uncertainty error and the measurement, weighing and volume dilution errors.

4.1 Assay Method #1 1002 ± 4 µg/mL (Avg 2 runs)

ICP Assay NIST SRM 3106 Lot Number: 991212

Assay Method #2 1002 µg/mL

Gravimetric NIST SRM Lot Number: See Sec. 4.2

4.2 BALANCE CALIBRATION - All balances are checked daily using in-house procedure number 6-IMM-001. The weights used for testing are annually compared to Gerhart Scale Corporation's master weights and are traceable to the National Institute of Standards and Technology (NIST). The NIST Traceability numbers are 692476 - Class 1 and 692476A - Class 2. The NIST test number is 822/260017-98. All analytical balances are calibrated every 4 months by Gerhart Scale Corp. of South Amboy. The balances are calibrated with a class 1 and/or class 2 analytical weight set. These weights are tested annually by a NIST / NVLAP accredited calibration lab. The NIST test number is 822/260017-98.

010262

4.3 THERMOMETER CALIBRATION - The thermometers used in the determination of the final densities are calibrated vs standard thermometer No. 903-2680 which was certified in accordance with the procedures outlined by ASTM E77-87 and NIST Monograph 150 using NIST Test Nos. and Std Nos.: 769543, 217368/769543, 217368/P14452, 176240/P14452, 176240. The in-house procedure No. is 2-QC-001. Thermometers which are not calibrated vs standard thermometer No. 903-2680 are traceable to NIST Identification Nos. 92564, 119016, 471047 and NIST test report Nos. 811/258522, 811/2557078, and 236090.

4.4 GLASSWARE CALIBRATION - In-house procedure 3-QC-002 is used to calibrate all Class A Glassware used in the manufacture and quality control of Custom Grade Standards.

5.0 TRACE METALLIC IMPURITIES (TMI) DETERMINED BY ICP/MS AND ICP-OES IN µg/mL

Custom-Grade solutions are tested for trace metallic impurities by Axial ICP-OES and ICP-MS. The result from the most sensitive method for each element, is reported below. Solutions tested by ICP-MS were analyzed in an ULPA-Filtered Clean Room. An ULPA-Filter is 99.9985% efficient for the removal of particles down to 0.3 µm.

<u>Q</u> Al 0.00012	<u>M</u> Dy < 0.01202	<u>Q</u> Li < 0.00002	<u>M</u> Pr < 0.00060	<u>M</u> Te < 0.06008
<u>M</u> Sb < 0.00100	<u>M</u> Er < 0.01001	<u>M</u> Lu < 0.00080	<u>M</u> Re < 0.00200	<u>M</u> Tb < 0.00060
<u>M</u> As < 0.02003	<u>M</u> Eu < 0.00601	<u>Q</u> Mg 0.00070	<u>M</u> Rh < 0.00200	<u>M</u> Tl < 0.00200
<u>M</u> Ba < 0.02003	<u>M</u> Gd < 0.00200	<u>Q</u> Mn < 0.00020	<u>M</u> Rb < 0.00200	<u>M</u> Th < 0.00200
<u>M</u> Be < 0.00100	<u>M</u> Ga < 0.00200	<u>Q</u> Hg < 0.01500	<u>M</u> Ru < 0.00401	<u>M</u> Tm < 0.00080
<u>s</u> Bi	<u>M</u> Ge < 0.01202	<u>M</u> Mo < 0.00401	<u>M</u> Sm < 0.00200	<u>M</u> Sn < 0.01001
<u>M</u> B < 0.14018	<u>M</u> Au < 0.00601	<u>M</u> Nd < 0.00401	<u>M</u> Sc < 0.02003	<u>M</u> Ti < 0.10013
<u>Q</u> Cd 0.00017	<u>M</u> Hf < 0.00401	<u>M</u> Ni < 0.01602	<u>M</u> Se < 0.01602	<u>M</u> W < 0.02003
<u>Q</u> Ca 0.00245	<u>M</u> Ho < 0.00100	<u>M</u> Nb < 0.00100	<u>Q</u> Si 0.00105	<u>M</u> U < 0.00401
<u>M</u> Ce < 0.01001	<u>Q</u> In 0.00105	<u>n</u> Os	<u>M</u> Ag < 0.00401	<u>M</u> V < 0.00401
<u>M</u> Cs < 0.00060	<u>M</u> Ir < 0.01001	<u>Q</u> Pd < 0.00400	<u>Q</u> Na 0.00240	<u>M</u> Yb < 0.00200
<u>Q</u> Cr 0.00020	<u>Q</u> Fe 0.00014	<u>Q</u> P < 0.01000	<u>M</u> Sr < 0.00100	<u>M</u> Y < 0.08011
<u>M</u> Co < 0.00601	<u>M</u> La < 0.00100	<u>M</u> Pt < 0.00401	<u>Q</u> S < 0.03000	<u>Q</u> Zn 0.00008
<u>Q</u> Cu 0.00014	<u>Q</u> Pb 0.00135	<u>Q</u> K 0.00039	<u>M</u> Ta < 0.01402	<u>M</u> Zr < 0.01001

M - Checked by ICP-MS

O - Checked by ICP-OES

i - Spectral Interference

n - Not Checked For

s - Solution Standard Element

6.0 INTENDED USE

For the calibration of analytical instruments including but not limited to the following:
ICP-MS, ICP-OES, FAAS, GFAA, XRF, and DCP

For the validation of analytical methods

For the preparation of "working reference samples"

For interference studies and the determination of correction coefficients

For detection limit and linearity studies

For additional intended uses, contact IV Technical Staff

7.0 INSTRUCTIONS FOR THE CORRECT USE OF THIS REFERENCE MATERIAL

010263

Storage & Handling - Keep tightly sealed when not in use. Store and use at $20 \pm 4^\circ\text{C}$. Do not pipet from container. Do not return portions removed for pipetting to container.

Atomic Weight; Valence; Coordination Number; Chemical Form in Solution - 208.9804; +3 6; $\text{Bi}(\text{O})(\text{H}_2\text{O})_5^{3+}$

Chemical Compatibility - Stable in HCl , HNO_3 , H_2SO_4 and HF . Avoid basic media forming insoluble hydroxide. Stable with most metals and inorganic anions in acidic media. Many salts that are insoluble in water are soluble in HCl , HNO_3 and HF . The major problem with Bi^{3+} is its tendency to hydrolyze at higher concentrations or in dilute acid. Nitric acid solutions should be 5% to hold the Bi in solution in the 100 to 10000 $\mu\text{g/mL}$ concentration range.

Stability - 2-100 ppb levels stable for months in 1% HNO_3 /LDPE container. 1-10,000 ppm solutions chemically stable for years in 5-7% HNO_3 /LDPE container.

Bi Containing Samples (Preparation and Solution) - Metal (soluble in HNO_3); Oxides (Soluble in HNO_3); Alloys (Dissolve in conc. 4:1 HCl/HNO_3 . Heating may be required.) Organic based (dry ash at 450°C and dissolve ash in HNO_3 or acid digestion with conc. hot sulfuric acid adding hydrogen peroxide dropwise and carefully until clear.)

Atomic Spectroscopic Information (ICP-OES D.L.s are given as radial/axial view):

Technique/Line	Estimated D.L.	Order	Type	Interferences (underlined indicates severe at all concs.)
ICP-OES 223.061 nm	0.04 / 0.005 $\mu\text{g/mL}$	1	atom	Th, Ir, Ti, Cu
ICP-OES 308.772 nm	0.08 / 0.01 $\mu\text{g/mL}$	1	atom	Th, U, Zr, Hf, Fe
ICP-OES 222.825 nm	0.1 / 0.02 $\mu\text{g/mL}$	1	atom	Cr, Hf, Ce, Os
ICP-MS 209 amu	2 ppt	n/a	M	<u>Bi</u> , <u>Pb</u> , <u>Mo</u>

8.0 HAZARDOUS INFORMATION - Please refer to the enclosed Material Safety Data sheet for information regarding this CRM.

9.0 HOMOGENEITY - This solution was mixed according to procedure IV-MPM-004 and is guaranteed to be homogeneous.

10.0 QUALITY STANDARD DOCUMENTATION

10.1 ISO 9001:2000 Quality Management System Registration - QMI Certificate Number 010105

Recognized by:

Registrar Accreditation Board (ANSI-RAB)

Standards Council of Canada (SCC)

Dutch Council for Accreditation (RVA)

Entidad Mexicana de Acreditacion, a.c.(EMA)

Members of IQ Net International Certification Network:

Argentina (IRAM), Australia (QAS), Austria (ÖQS), Belgium (Avinter), Brazil (FCAV), Canada (QMI), Hong Kong (HKQAA), Columbia (ICONTEC), Czech Republic (CQS), Denmark (DS), Finland (SFS), France (AFAQ), Germany (DQS), Greece (ELOT), Hungary (MSZT), Ireland (NSAI), Israel (SII), Italy (CISQ), Japan (JQA), Korea (KSA-QA), Netherlands (KEMA), Norway (NCS), Poland (PCBC), Portugal (APCER), Singapore (PSB), Slovenia (SIQ), Spain (AENOR), Switzerland (SQS)

10.2 ISO/IEC 17025 - 1999 "General Requirements for the Competence of Testing and Calibration"

- Chemical Testing - Accredited A2LA Certificate Number 883.01



10.3 ISO/IEC Guide 34 - 2000 "General Requirements for the Competence of Reference Material Producers"

- Reference Materials Production - Accredited A2LA Certificate Number 883.02

A2LA Mutual Recognition Agreement Partners:

Australia (NATA), Austria (BmWA), Belgium (BELTEST) (BKO-OBE), Canada (SCC), Chinese Taipei (CNLA), Czech Republic (NAO), Denmark (DANAK), Finland (FINAS), France (COFRAC), Germany (DAR), Hong Kong (HKAS), Ireland (NAB), Italy (SIT) (SINAL), Japan (JAB) (JNLA), Republic of Korea (KOLAS), The Netherlands (RvA), New Zealand (IANZ), Norway (NA), Portugal (IPQ), Singapore (SAC-SINGLAS), Spain (ENAC), Sweden (SWEDAC), Switzerland (SAS), United Kingdom (UKAS) and United States (NVLAP) (ICBO ES)

10.4 10CFR50 Appendix B - Nuclear Regulatory Commission

- Domestic Licensing of Production and Utilization Facilities

10.5 10CFR21 - Nuclear Regulatory Commission - Reporting Defects and Non-Compliance

10.6 MIL-STD-45662A (Obsolete/Observed)

INORGANIC LABS/RADCHEM LABS Pg. 2 of 2
 DATE RECEIVED: 07/31/03
 DATE EXPIRED: 08/01/2004
 DATE OPENED: 08/01/03
 INORG: 4200 PO: F52383

11.0 DATE OF CERTIFICATION AND PERIOD OF VALIDITY

010264



11.1 IV Shelf Life - The period of time during which the concentration of the analyte(s) in a properly packaged, unopened, and unused standard stored under environmentally controlled and monitored conditions will remain within the specified uncertainty range. Shelf life is limited primarily by transpiration (loss of water from the solution) and infrequently, by chemical instability. Transpiration studies (P-SP01020) of chemically-stable solutions performed at Inorganic Ventures / IV Labs indicate a CRM shelf-life of four years for solutions packaged in 500-mL low density polyethylene bottles. When stored under special conditions that minimize transpiration and instability, the shelf life can be extended past this limit.

11.2 Expiration Date - The date after which a CRM should not be used. Routine laboratory use of a CRM increases transpiration losses and the chance of contamination which affect the integrity of the CRM and limit its useful life. Inorganic Ventures / IV Labs concurs with state and federal regulatory agencies' recommendations that solution standards be assigned a one-year expiration date.

Certification Date: March 28, 2003

Expiration Date:

EXPIRES
1/2004

12.0 NAMES AND SIGNATURES OF CERTIFYING OFFICERS

Certificate Prepared By: JoAnn Struthers, QA Administrative Assistant

Certificate Approved By: Katalin Le, QC Supervisor

Certifying Officer: Paul Gaines, Chemist, Senior Technical Director

JoAnn Struthers
Katalin Le
Paul Gaines

**SOUTHWEST RESEARCH INSTITUTE
NUCLEAR PROJECT
CLIENT: Div. 20
TASK ORDER: 031106-12
SRR: 25200
SDG: 237578
CASE: CNWRA
VTSR: November 06, 2003
PROJECT#: 06002.01.141**

Pipette Calibrations

SwRI – Div. 01, Inorganic Labs' Fixed Volume Pipette Verification Log

(Space provide for Inorganic Laboratories' Fixed Volume Pipette Verification Spreadsheet)

010263*L. Wright
11/4/03***SwRI - Div. 01, Inorganic Labs' Fixed Volume Pipette Spreadsheet**

Eppendorf #	True Value (uL)	1st Reading (g)	2nd Reading (g)	3rd Reading (g)	Avg Wt (g)	% of True Value
Lab30	1000	1.0184	1.0102	1.0165	1.015	101.50
TMA1	1000	1.0187	1.0064	1.0019	1.009	100.90
TMA2	1000	1.0175	1.0126	1.0052	1.012	101.18
TMA3	1000	Out of Service				
TMA6	1000	1.0013	1.0015	1.0033	1.002	100.20
TMB1	900	0.9002	0.9001	0.9005	0.900	100.03
TMC1	800	0.803	0.7978	0.7995	0.800	100.01
TMDD1	750	0.7484	0.7491	0.7479	0.748	99.80
TMD1	700	0.7009	0.6906	0.6998	0.697	99.59
TMD2	700	0.7065	0.7082	0.7049	0.707	100.93
TME1	600	0.607	0.605	0.6051	0.606	100.95
TMF2	500	0.5077	0.5062	0.5026	0.506	101.10
TMF5	500	0.509	0.5061	0.501	0.505	101.07
ICF1	500	0.5012	0.5014	0.5003	0.501	100.19
L30-500	500	0.4997	0.4998	0.4989	0.499	99.89
TMG3	400	0.3979	0.3958	0.3995	0.398	99.43
TMH1	300	Out of Service				
TMH2	300	0.2993	0.2978	0.3009	0.299	99.78
TMJ1	250	0.25	0.2503	0.2499	0.250	100.03
TMJ2	250	0.2488	0.2492	0.2495	0.249	99.67
TMJ3	250	0.2504	0.2503	0.2519	0.251	100.35
TMK2	200	0.2	0.2001	0.2006	0.200	100.12
TML1	150	0.1486	0.1493	0.149	0.149	99.31
TMM1	120	0.1207	0.1202	0.1205	0.120	100.39
TMN3	100	0.1003	0.1001	0.1002	0.100	100.20
ICN1	100	0.1004	0.0999	0.0993	0.100	99.87
TMQ1	80	0.0799	0.08	0.0799	0.080	99.92
TMR1	70	Out of Service				
TMS1	60	Out of Service				
LAB-30A	50	0.0498	0.0497	0.0496	0.050	99.40
TMU1	40	0.0402	0.0403	0.0402	0.040	100.58
TMU2	40	0.04	0.0399	0.04	0.040	99.92
TMV1	30	0.0296	0.0298	0.0297	0.030	99.00
L30-20	20	0.0199	0.0198	0.0199	0.020	99.33
TMW1	25	0.025	0.025	0.0249	0.025	99.87
TMY1	15	Out of Service				

FRM-246 (Rev 1/Mar 03)

SwRI - Div. 01, Inorganic Labs' Fixed Volume Pipette Verification Log

Balance #: 16Thermometer #: G011diH2O Temperature (°C): 20.10266

Eppendorf #	True Value (µL)	1st Reading (g)	2nd Reading (g)	3rd Reading (g)
Lab30	1000	1.0184	1.0102	1.0165
TMA1	1000	1.0187	1.0064	1.0019
TMA2	1000	1.0175	1.0126	1.0052
~~TMA3~~	~~1000~~	~~out~~		
<u>TMA6</u> ~~out~~	1000 ^{10/13/03}	1.0414 1.0013	1.0220 1.0015	1.0105 1.0033
TMB1	900	0.9002	0.9001	0.9005
TMC1	800	0.8030	0.7978	0.7995
TMDD1	750	0.7484	0.7491	0.7479
TMD1	700	0.7009	0.6906	0.6998
TMD2	700	0.7065	0.7082	0.7049
TME1	600	0.6070	0.6050	0.6051
TMF2	500	0.5077	0.5062	0.5026
TMF5	500	0.5090	0.5061	0.5010
ICF1	500	0.5012	0.5014	0.5003
L30-500	500	0.4997	0.4998	0.4989
TMG3	400	0.3979	0.3958	0.3995
<u>TMH1</u> ~~out~~	300 ^{of service}	0.2192	0.2100	0.2182
TMH2	300	0.2993	0.2978	0.3009
TMJ1	250	0.2500	0.2503	0.2499
TMJ2	250	0.2488	0.2492	0.2495
TMJ3	250	0.2504	0.2503	0.2519
TMK2	200	0.2000	0.2001	0.2006
TML1	150	0.1486	0.1493	0.1490
TMM1	120	0.1207	0.1202	0.1205
TMN3	100	0.1003	0.1001	0.1002
ICN1	100	0.1004	0.0999	0.0993
TMQ1	80	0.0799	0.0800	0.0799
~~TMR1~~	~~70~~	~~out~~		
~~TMS1~~	~~60~~	~~out~~		
LAB-30A	50	0.0498	0.0497	0.0496
TMU1	40	0.0402	0.0403	0.0402
TMU2	40	0.0400	0.0399	0.0400
TMV1	30	0.0296	0.0298	0.0297
L30-20	20	0.0199	0.0198	0.0199
TMW1	25	0.0250	0.0250	0.0249
~~TMY1~~	~~15~~	~~out~~		

Analyst: Lana WrightDate: 10/21/03Reviewed by: Val W. AllenDate: 12/02/03

SwRI – Div. 01, Inorganic Labs' Adjustable Volume Pipette Verification Log

(Space provided for Inorganic Laboratories' Adjustable Volume Pipette Verification Spreadsheet)

010267

SwRI
11-14-03

SwRI – Div. 01, Inorganic Laboratory Adjustable Pipette Verification Spreadsheet

Eppendorf #	True Value (µL)	1st Reading (g)	2nd Reading (g)	3rd Reading (g)	Avg Wt (g)	% of True Value
	20	0.0198	0.0198	0.0199	0.020	99.17
ADJ200-A	100	0.0988	0.0980	0.0984	0.098	98.40
	200	0.1987	0.1981	0.1984	0.198	99.20
	20	0.0199	0.0198	0.0198	0.020	99.17
ADJ200-C	100	0.0984	0.0983	0.0983	0.098	98.33
	200	0.1973	0.1975	0.1980	0.198	98.80
	20	0.0202	0.0199	0.0199	0.020	100.00
ADJ200-D	100	0.0980	0.0981	0.0988	0.098	98.30
	200	0.1977	0.1989	0.1986	0.198	99.20
	20				0.000	0.00
ADJ200-G	100				0.000	0.00
	200				0.000	0.00
	20				0.000	0.00
ADJ200-H	100				0.000	0.00
	200				0.000	0.00
	20				0.000	0.00
ADJ200-J	100				0.000	0.00
	200				0.000	0.00
	20				0.000	0.00
ADJ200	100				0.000	0.00
	200				0.000	0.00
	20				0.000	0.00
ADJ200	100				0.000	0.00
	200				0.000	0.00
	20				0.000	0.00
ADJ200	100				0.000	0.00
	200				0.000	0.00

FRM-247a (Rev 3/Oct 03)

SwRI Div. 01 – Inorganic Laboratory Adjustable Pipette Verification Log **010268**

Balance #: 16

Thermometer #: G 011

diH2O Temperature (° C) 22

Eppendorf #	True Value (µL)	1 st Reading (g)	2 nd Reading (g)	3 rd Reading (g)
	20	0.0198	0.0198	0.0199
ADJ200-A	100	0.0988	0.0980	0.0984
34	200	0.1987	0.1981	0.1984
	20	0.0199	0.0198	0.0198
ADJ200-C	100	0.0984	0.0983	0.0983
39	200	0.1973	0.1975	0.1980
	20	0.0202	0.0199	0.0199
ADJ200-D	100	0.0980	0.0981	0.0988
27	200	0.1977	0.1989	0.1986
	20			
ADJ200-G	100			
	200			
	20			
ADJ200-H	100			
	200			
	20			
ADJ200-J	100			
	200			
	20			
ADJ200	100			
	200			
	20			
ADJ200	100			
	200			

pipette verification performed 11-7 by L.W., computer spreadsheet posted 11-14 by J. Wills

Analyst: Laura Wright
Reviewed by: John Wills

Date: 11/7/03
Date: 11-14-03

SwRI – Div. 01, Inorganic Labs' Adjustable Volume Pipette Verification Log

(Space provided for Inorganic Laboratories' Adjustable Volume Pipette Verification Spreadsheet)

010269

SwRI – Div. 01, Inorganic Laboratory Adjustable Pipette Verification Spreadsheet

Eppendorf #	True Value (μL)	1st Reading (g)	2nd Reading (g)	3rd Reading (g)	Avg Wt (g)	% of True Value
ADJ200-A	20				0.000	0.00
	100				0.000	0.00
	200				0.000	0.00
ADJ200-C	20				0.000	0.00
	100				0.000	0.00
	200				0.000	0.00
ADJ200-D	20				0.000	0.00
	100				0.000	0.00
	200				0.000	0.00
ADJ200-G	20	0.0200	0.0199	0.0199	0.020	99.67
	100	0.1010	0.1001	0.1007	0.101	100.60
	200	0.2001	0.1997	0.1990	0.200	99.80
ADJ200-H	20	0.0201	0.0202	0.0201	0.020	100.67
	100	0.0990	0.0989	0.0991	0.099	99.00
	200	0.1989	0.1995	0.1979	0.199	99.38
ADJ200-J	20	0.0201	0.0202	0.0202	0.020	100.83
	100	0.0985	0.0990	0.0991	0.099	98.87
	200	0.1989	0.1991	0.1983	0.199	99.38
ADJ200	20				0.000	0.00
	100				0.000	0.00
	200				0.000	0.00
ADJ200	20				0.000	0.00
	100				0.000	0.00
	200				0.000	0.00
ADJ200	20				0.000	0.00
	100				0.000	0.00
	200				0.000	0.00

FRM-247a (Rev 3/Oct 03)

SwRI Div. 01 – Inorganic Laboratory Adjustable Pipette Verification Log **010270**

Balance #: 34

Thermometer #: 6011

diH2O Temperature (°C) 21

20 μ L – 200 μ L	Eppendorf #	True Value (μ L)	1 st Reading (g)	2 nd Reading (g)	3 rd Reading (g)
		20			
	ADJ200-A	100			
		200			
		20			
	ADJ200-C	100			
		200			
		20			
	ADJ200-D	100			
		200			
		20	0.0200	0.0199	0.0199
	ADJ200-G	100	0.1010	0.1001	0.1007
		200	0.2001	0.1997	0.1990
		20	0.0201	0.0202	0.0201
	ADJ200-H	100	0.0996	0.0989	0.0991
		200	0.1989	0.1995	0.1979
		20	0.0201	0.0202	0.0202
	ADJ200-J	100	0.0985	0.0990	0.0991
		200	0.1989	0.1991	0.1983
		20			
	ADJ200	100			
		200			
		20			
	ADJ200	100			
		200			

Analyst: J. Van
 Reviewed by: AMS

Date: 11/11/03
 Date: 11/17/03

SwRI – Div. 01, Inorganic Labs' Adjustable Volume Pipette Verification Log

(Space provided for Inorganic Laboratories' Adjustable Volume Pipette Verification Spreadsheet)

010271

SwRI – Div. 01, Inorganic Laboratory Adjustable Pipette Verification Spreadsheet

Eppendorf #	True Value (μL)	1st Reading (g)	2nd Reading (g)	3rd Reading (g)	Avg Wt (g)	% of True Value
ADJ200-A	20	0.0204	0.0202	0.0203	0.020	101.50
	100	0.0999	0.0998	0.0993	0.100	99.67
	200	0.1996	0.1990	0.1996	0.199	99.70
ADJ200-C	20	0.0200	0.0203	0.0203	0.020	101.00
	100	0.0995	0.0993	0.0991	0.099	99.30
	200	0.1992	0.1993	0.1996	0.199	99.68
ADJ200-D	20	0.0202	0.0204	0.0200	0.020	101.00
	100	0.0999	0.0995	0.0997	0.100	99.70
	200	0.1993	0.1995	0.1997	0.200	99.75
ADJ200-G	20					
	100					
	200					
ADJ200-H	20					
	100					
	200					
ADJ200-J	20					
	100					
	200					
ADJ200	20					
	100					
	200					
ADJ200	20					
	100					
	200					
ADJ200	20					
	100					
	200					

FRM-247a (Rev 3/Oct 03)

FRM-244 (Rev 2/Sept 02)

SwRI Div. 01 – Inorganic Laboratory Adjustable Pipette Verification Log 010272

Balance #: 16 Thermometer #: 6011 diH2O Temperature (° C) 21

20 µL – 200 µL	Eppendorf #	True Value (µL)	1 st Reading (g)	2 nd Reading (g)	3 rd Reading (g)
		20	.0204	.0202	.0203
	ADJ200-A	100	.0999	.0998	.0993
		200	.1996	.1990	.1996
		20	.0200	.0203	.0203
	ADJ200-C	100	.0995	.0993	.0991
		200	.1992	.1993	.1996
		20	.0202	.0204	.0200
	ADJ200-D	100	.0999	.0995	.0997
		200	.1993	.1995	.1997
		20			
	ADJ200-G	100			
		200			
		20			
	ADJ200-H	100			
		200			
		20			
	ADJ200-J	100			
		200			
		20			
	ADJ200	100			
		200			
		20			
	ADJ200	100			
		200			

Analyst: John Wilks
Reviewed by: Val W. Allen

Date: 11-21-03
Date: 12/22/03

SwRI – Div. 01, Inorganic Labs' Adjustable Volume Pipette Verification Log

(Space provided for Inorganic Laboratories' Adjustable Volume Pipette Verification Spreadsheet)

010273
SwRI – Div. 01, Inorganic Laboratory Adjustable Pipette Verification Spreadsheet

Eppendorf #	True Value (μL)	1st Reading (g)	2nd Reading (g)	3rd Reading (g)	Avg Wt (g)	% of True Value
	20				0.000	0.00
ADJ200-A	100				0.000	0.00
	200				0.000	0.00
	20				0.000	0.00
ADJ200-C	100				0.000	0.00
	200				0.000	0.00
	20				0.000	0.00
ADJ200-D	100				0.000	0.00
	200				0.000	0.00
	20	0.0201	0.0201	0.0200	0.020	100.33
ADJ200-G	100	0.0989	0.0984	0.0989	0.099	98.73
	200	0.2001	0.2003	0.2050	0.202	100.90
	20	0.0202	0.0203	0.0204	0.020	101.50
ADJ200-H	100	0.1009	0.1010	0.1007	0.101	100.87
	200	0.1986	0.1984	0.1989	0.199	99.32
	20	0.0197	0.0198	0.0200	0.020	99.17
ADJ200-J	100	0.0991	0.0993	0.0991	0.099	99.17
	200	0.1989	0.1991	0.1992	0.199	99.53
	20				0.000	0.00
ADJ200	100				0.000	0.00
	200				0.000	0.00
	20				0.000	0.00
ADJ200	100				0.000	0.00
	200				0.000	0.00
	20				0.000	0.00
ADJ200	100				0.000	0.00
	200				0.000	0.00

FRM-247a (Rev 3/Oct 03)

SwRI Div. 01 – Inorganic Laboratory Adjustable Pipette Verification Log 010274

Balance #: 34 Thermometer #: 3011 diH2O Temperature (°C) 21

Eppendorf #	True Value (μL)	20 μL – 200 μL		
		1 st Reading (g)	2 nd Reading (g)	3 rd Reading (g)
	20	J.V. 11/21/03		
ADJ200-A	100			
	200			
	20			
ADJ200-C	100			
	200			
	20			
ADJ200-D	100			
	200			
	20	0.0201	0.0201	0.0200
ADJ200-G	100	0.0989	0.0984	0.0989
	200	0.2001	0.2003	0.2050
	20	0.0202	0.0203	0.0204
ADJ200-H	100	0.1009	0.1010	0.1007
	200	0.1986	0.1984	0.1989
	20	0.0197	0.0198	0.0200
ADJ200-J	100	0.0991	0.0993	0.0991
	200	0.1989	0.1991	0.1992
	20	J.V. 11/21/03		
ADJ200	100			
	200			
	20			
ADJ200	100			
	200			

Analyst: J. Vaw
Reviewed by: J. Vaw

Date: 11/21/03
Date: 12/22/03

SwRI – Div. 01, Inorganic Labs' Adjustable Volume Pipette Verification Log

(Space provided for Inorganic Laboratories' Adjustable Volume Pipette Verification Spreadsheet)

010275

SwRI – Div. 01, Inorganic Laboratory Adjustable Pipette Verification Spreadsheet

Eppendorf #	True Value (µL)	1st Reading (g)	2nd Reading (g)	3rd Reading (g)	Avg Wt (g)	% of True Value
ADJ1000-C	100	0.0989	0.0991	0.0990	0.099	99.00
	500	0.4910	0.4961	0.4982	0.495	99.02
	1000	0.9834	0.9869	0.9843	0.985	98.49
ADJ1000-D	100	0.0990	0.0993	0.0996	0.099	99.30
	500	0.4924	0.4995	0.4983	0.497	99.35
	1000	0.9866	0.9860	0.9870	0.987	98.65
ADJ1000-E	100	0.1010	0.1015	0.0993	0.101	100.60
	500	0.4985	0.4977	0.4986	0.498	99.65
	1000	0.9977	0.9985	0.9961	0.997	99.74
ADJ1000-F	100	0.0999	0.1001	0.0992	0.100	99.73
	500	0.4989	0.4963	0.4960	0.497	99.41
	1000	1.0067	0.9976	0.9996	1.001	100.13
ADJ1000-G	100				0.000	0.00
	500				0.000	0.00
	1000				0.000	0.00
ADJ1000-H	100				0.000	0.00
	500				0.000	0.00
	1000				0.000	0.00
ADJ1000-J	100				0.000	0.00
	500				0.000	0.00
	1000				0.000	0.00
ADJ1000	100				0.000	0.00
	500				0.000	0.00
	1000				0.000	0.00
ADJ1000	100				0.000	0.00
	500				0.000	0.00
	1000				0.000	0.00

FRM-247b (Rev 2/Oct 03)

FRM-244 (Rev 2/Sept 02)

SwRI Div. 01 – Inorganic Laboratory Adjustable Pipette Verification Log

Balance #: 16

Thermometer #: G011

diH2O Temperature (° C) 22

010276

Eppendorf #	True Value (µL)	1 st Reading (g)	2 nd Reading (g)	3 rd Reading (g)
	100	0.0989	0.0991	0.0990
ADJ1000-C	500	0.4910	0.4961	0.4982
27	1000	0.9834	0.9869	0.9843
hard turn	100	0.0990	0.0993	0.0996
ADJ1000-D	500	0.4924	0.4995	0.4983
27	1000	0.9866	0.9860	0.9870
	100	0.1010	0.1015	0.0993
ADJ1000-E	500	0.4985	0.4977	0.4986
34	1000	0.9977 ^W 1.0065 ^{11/7/03}	0.9985	0.9961
	100	0.0999	0.1001	0.0992
ADJ1000-F	500	0.4989	0.4963	0.4960
39	1000	1.0067	0.9976	0.9996
	100			
ADJ1000-G	500			
	1000			
	100			
ADJ1000-H	500			
	1000			
	100			
ADJ1000-J	500			
	1000			
	100			
ADJ1000	500			
	1000			
	100			
ADJ1000	500			
	1000			

pipette verification performed 11-7 by L.W., spreadsheet posted 11-14 by J. Willis

Analyst: Auna Wright
 Reviewed by: John Willis

Date: 11/7/03
 Date: 11-14-03

SwRI – Div. 01, Inorganic Labs' Adjustable Volume Pipette Verification Log

010277

(Space provided for Inorganic Laboratories' Adjustable Volume Pipette Verification Spreadsheet)

SwRI – Div. 01, Inorganic Laboratory Adjustable Pipette Verification Spreadsheet

Eppendorf #	True Value (μL)	1st Reading (g)	2nd Reading (g)	3rd Reading (g)	Avg Wt (g)	% of True Value
	100				0.000	0.00
ADJ1000-C	500				0.000	0.00
	1000				0.000	0.00
	100				0.000	0.00
ADJ1000-D	500				0.000	0.00
	1000				0.000	0.00
	100				0.000	0.00
ADJ1000-E	500				0.000	0.00
	1000				0.000	0.00
	100				0.000	0.00
ADJ1000-F	500				0.000	0.00
	1000				0.000	0.00
	100	0.1025	0.1015	0.1013	0.102	101.77
ADJ1000-G	500	0.4996	0.4995	0.4992	0.499	99.89
	1000	0.9911	0.9915	0.9920	0.992	99.15
	100	0.1011	0.1014	0.1011	0.101	101.20
ADJ1000-H	500	0.4975	0.4965	0.4935	0.496	99.17
	1000	0.9978	0.9984	0.9981	0.998	99.81
	100	0.1009	0.1011	0.1015	0.101	101.17
ADJ1000-J	500	0.4985	0.4981	0.4973	0.498	99.59
	1000	0.9982	0.9971	0.9965	0.997	99.73
	100				0.000	0.00
ADJ1000	500				0.000	0.00
	1000				0.000	0.00
	100				0.000	0.00
ADJ1000	500				0.000	0.00
	1000				0.000	0.00

FRM-247b (Rev 2/Oct 03)

SwRI Div. 01 – Inorganic Laboratory Adjustable Pipette Verification Log **010278**

Balance #: 34 Thermometer #: G 011 diH₂O Temperature (° C) 21

Eppendorf #	True Value (μL)	1 st Reading (g)	2 nd Reading (g)	3 rd Reading (g)
100 μL – 1000 μL	100	~~16V 11/11/03~~	~~16V 11/11/03~~	~~16V 11/11/03~~
	ADJ1000-C			
	500			
	1000			
	100			
	ADJ1000-D			
	500			
	1000			
	100			
	ADJ1000-E			
	500			
	1000			
100 μL – 1000 μL	100	0.1025	0.1015	0.1013
	ADJ1000-G	0.4996	0.4993	0.4992
	500	0.9911	0.9915	0.9920
	1000	0.1011	0.1014	0.1011
	ADJ1000-H	0.4975	0.4965	0.4935
	500	0.9978	0.9984	0.9981
	1000	0.1009	0.1011	0.1015
	ADJ1000-J	0.4985	0.4981	0.4973
	500	0.9982	0.9971	0.9965
	1000	~~16V 11/11/03~~	~~16V 11/11/03~~	~~16V 11/11/03~~
	ADJ1000			
	500			
	1000			
100 μL – 1000 μL	ADJ1000			
	500			
	1000			
	ADJ1000			

Analyst: [Signature]
Reviewed by: [Signature]

Date: 11/11/03
Date: 11/17/03

SwRI – Div. 01, Inorganic Labs' Adjustable Volume Pipette Verification Log

(Space provided for Inorganic Laboratories' Adjustable Volume Pipette Verification Spreadsheet)

010279

SwRI – Div. 01, Inorganic Laboratory Adjustable Pipette Verification Spreadsheet

Eppendorf #	True Value (µL)	1st Reading (g)	2nd Reading (g)	3rd Reading (g)	Avg Wt (g)	% of True Value
	100	0.1000	0.0993	0.1000	0.100	99.77
ADJ1000-C	500	0.4959	0.4983	0.4964	0.497	99.37
	1000	0.9941	0.9959	0.9933	0.994	99.44
	100	0.1003	0.1016	0.1007	0.101	100.87
ADJ1000-D	500	0.4953	0.4973	0.4966	0.496	99.28
	1000	0.9929	0.9916	0.9957	0.993	99.34
	100	0.1007	0.1000	0.1009	0.101	100.53
ADJ1000-E	500	0.4945	0.4978	0.4939	0.495	99.08
	1000	0.9944	0.9964	0.9986	0.996	99.65
	100	0.1020	0.1018	0.1008	0.102	101.53
ADJ1000-F	500	0.4958	0.4991	0.4973	0.497	99.48
	1000	0.9995	1.0077	0.9987	1.002	100.20
	100					
ADJ1000-G	500					
	1000					
	100					
ADJ1000-H	500					
	1000					
	100					
ADJ1000-J	500					
	1000					
	100					
ADJ1000	500					
	1000					
	100					
ADJ1000	500					
	1000					

FRM-247b (Rev 2/Oct 03)

SwRI Div. 01 – Inorganic Laboratory Adjustable Pipette Verification Log

Balance #: 116

Thermometer #: 6011

diH2O Temperature (° C) 21 **010280**

Eppendorf #	True Value (µL)	1 st Reading (g)	2 nd Reading (g)	3 rd Reading (g)
100 µL – 1000 µL	100	.1000	.0993	.1000
	ADJ1000-C	500	.4959	.4983
		1000	.9941	.9959
ADJ1000-D	100	.1003	.1016	.1007
	500	.4953	.4973	.4966
	1000	.9929	.9916	.9957
ADJ1000-E	100	.1007	.1000	.1009
	500	.4945	.4978	.4939
	1000	.9944	.9964	.9986
ADJ1000-F	100	.1020	.1018	.1008
	500	.4958	.4991	.4973
	1000	.9995	1.0077	.9987
ADJ1000-G	100			
	500			
	1000			
ADJ1000-H	100			
	500			
	1000			
ADJ1000-J	100			
	500			
	1000			
ADJ1000	100			
	500			
	1000			
ADJ1000	100			
	500			
	1000			

Analyst: John Wilks
 Reviewed by: [Signature]

Date: 11-21-03
 Date: 12/22/03

SwRI – Div. 01, Inorganic Labs' Adjustable Volume Pipette Verification Log

010281

(Space provided for Inorganic Laboratories' Adjustable Volume Pipette Verification Spreadsheet)

11/21/03

SwRI – Div. 01, Inorganic Laboratory Adjustable Pipette Verification Spreadsheet

Eppendorf #	True Value (µL)	1st Reading (g)	2nd Reading (g)	3rd Reading (g)	Avg Wt (g)	% of True Value
	100				0.000	0.00
ADJ1000-C	500				0.000	0.00
	1000				0.000	0.00
	100				0.000	0.00
ADJ1000-D	500				0.000	0.00
	1000				0.000	0.00
	100				0.000	0.00
ADJ1000-E	500				0.000	0.00
	1000				0.000	0.00
	100				0.000	0.00
ADJ1000-F	500				0.000	0.00
	1000				0.000	0.00
	100	0.0991	0.0990	0.0989	0.099	99.00
ADJ1000-G	500	0.4985	0.4989	0.4990	0.499	99.76
	1000	0.9963	0.9968	0.9972	0.997	99.68
	100	0.1021	0.1010	0.1015	0.102	101.53
ADJ1000-H	500	0.5051	0.5042	0.5012	0.504	100.70
	1000	0.9975	0.9985	0.9979	0.998	99.80
	100	0.1009	0.1004	0.1007	0.101	100.67
ADJ1000-J	500	0.5011	0.4997	0.4992	0.500	100.00
	1000	1.0013	1.0025	1.0015	1.002	100.18
	100				0.000	0.00
ADJ1000	500				0.000	0.00
	1000				0.000	0.00
	100				0.000	0.00
ADJ1000	500				0.000	0.00
	1000				0.000	0.00

FRM-247b (Rev 2/Oct 03)

FRM-244 (Rev 2/Sept 02)

Book/page: 04 080

SwRI Div. 01 – Inorganic Laboratory Adjustable Pipette Verification Log **010282**

Balance #: 34

Thermometer #: G 011

diH2O Temperature (°C) 21

Eppendorf #	True Value (µL)	1 st Reading (g)	2 nd Reading (g)	3 rd Reading (g)
ADJ1000-C	100			
	500			
	1000			
ADJ1000-D	100			
	500			
	1000			
ADJ1000-E	100			
	500			
	1000			
ADJ1000-F	100			
	500			
	1000			
ADJ1000-G	100	~~0.9990~~	0.0990	0.0989
	500	0.4985	0.4989	0.4990
	1000	0.9963	0.9968	0.9972
ADJ1000-H	100	0.1021	0.1010	0.1015
	500	0.5051	0.5042	0.5012
	1000	0.9975	0.9985	0.9979
ADJ1000-J	100	1.009	1.004	1.007
	500	0.5011	0.4997	0.4992
	1000	1.0013	1.0025	1.0015
ADJ1000	100			
	500			
	1000			
ADJ1000	100			
	500			
	1000			

Analyst: J. Vor

Reviewed by: [Signature]

Date: 11/21/03

Date: 12/22/03

SwRI – Div. 01, Inorganic Labs' Adjustable Volume Pipette Verification Log

(Space provided for Inorganic Laboratories' Adjustable Volume Pipette Verification Spreadsheet)

010283

SwRI – Div. 01, Inorganic Laboratory Adjustable Pipette Verification Spreadsheet

Eppendorf #	True Value (μL)	1st Reading (g)	2nd Reading (g)	3rd Reading (g)	Avg Wt (g)	% of True Value
ADJ5000-C	500	0.4949	0.4932	0.4938	0.494	98.79
	2500	2.4899	2.4903	2.4962	2.492	99.69
	5000	5.0189	5.0071	5.0011	5.009	100.18
ADJ5000-G	500	0.5091	0.5068	0.5080	0.508	101.59
	2500	2.5106	2.5156	2.5040	2.510	100.40
	5000	5.0016	5.0095	5.0121	5.008	100.15
ADJ5000-H	500	0.5033	0.5016	0.5088	0.505	100.91
	2500	2.4950	2.4843	2.4910	2.490	99.60
	5000	5.0141	4.9930	4.9982	5.002	100.04
ADJ5000-I	500	0.5003	0.4975	0.4956	0.498	99.56
	2500	2.4983	2.4824	2.4996	2.493	99.74
	5000	4.9885	4.9874	4.9773	4.984	99.69
ADJ5000-J	500				0.000	0.00
	2500				0.000	0.00
	5000				0.000	0.00
ADJ5000-K	500				0.000	0.00
	2500				0.000	0.00
	5000				0.000	0.00
ADJ5000-L	500				0.000	0.00
	2500				0.000	0.00
	5000				0.000	0.00
ADJ5000	500				0.000	0.00
	2500				0.000	0.00
	5000				0.000	0.00
ADJ5000	500				0.000	0.00
	2500				0.000	0.00
	5000				0.000	0.00
ADJ5000	500				0.000	0.00
	2500				0.000	0.00
	5000				0.000	0.00

FRM-247c (Rev 2/Mar 03)

FRM-244 (Rev 2/Sept 02)

SwRI Div. 01 – Inorganic Laboratory Adjustable Pipette Verification Log **010284**Balance #: 16Thermometer #: G0011diH₂O Temperature (° C) 22

Eppendorf #	True Value (μL)	1 st Reading (g)	2 nd Reading (g)	3 rd Reading (g)
	500	0.4949	0.4932	0.4938
ADJ5000-C	2500	2.4899	2.4903	2.4962
39	5000	5.0189	5.0071	5.0011
	500	0.5091	0.5068	0.5080
ADJ5000-G	2500	2.5106	2.5156	2.5040
27	5000	5.0016	5.0095	5.0121
	500	0.5033	0.5016	0.5088
ADJ5000-H	2500	2.4950	2.4843	2.4910
34	5000	5.0141	4.9930	4.9982
	500	0.5003	0.4975	0.4956
ADJ5000-I	2500	2.4983	2.4824	2.4996
27	5000	4.9885	4.9874	4.9773
	500			
ADJ5000-J	2500			
	5000			
	500			
ADJ5000-K	2500			
	5000			
	500			
ADJ5000-L	2500			
	5000			
	500			
ADJ5000	2500			
	5000			
	500			
ADJ5000	2500			
	5000			
	500			
ADJ5000	2500			
	5000			
	500			
ADJ5000	2500			
	5000			

500 μL – 5000 μL

pipet verification performed 11-7 by L.W., spreadsheet
 posted 11-14 by J.W.

Analyst: Laura Wright
 Reviewed by: John Wells

Date: 11/7/03
 Date: 11-14-03

SwRI – Div. 01, Inorganic Labs' Adjustable Volume Pipette Verification Log

(Space provided for Inorganic Laboratories' Adjustable Volume Pipette Verification Spreadsheet)

010285

SwRI – Div. 01, Inorganic Laboratory Adjustable Pipette Verification Spreadsheet

Eppendorf #	True Value (µL)	1st Reading (g)	2nd Reading (g)	3rd Reading (g)	Avg Wt (g)	% of True Value
	500				0.000	0.00
ADJ5000-C	2500				0.000	0.00
	5000				0.000	0.00
	500				0.000	0.00
ADJ5000-G	2500				0.000	0.00
	5000				0.000	0.00
	500				0.000	0.00
ADJ5000-H	2500				0.000	0.00
	5000				0.000	0.00
	500				0.000	0.00
ADJ5000-I	2500				0.000	0.00
	5000				0.000	0.00
	500	0.4975	0.4979	0.4971	0.498	99.50
ADJ5000-J	2500	2.4988	2.4983	2.4981	2.498	99.94
	5000	5.0075	5.0072	5.0061	5.007	100.14
	500	0.5011	0.5054	0.5073	0.505	100.92
ADJ5000-K	2500	2.5011	2.5081	2.5031	2.504	100.16
	5000	5.0625	5.0355	5.0335	5.044	100.88
	500	0.4994	0.4993	0.4981	0.499	99.79
ADJ5000-L	2500	2.4979	2.4963	2.4981	2.497	99.90
	5000	5.0521	5.0241	5.0341	5.037	100.74
	500				0.000	0.00
ADJ5000	2500				0.000	0.00
	5000				0.000	0.00
	500				0.000	0.00
ADJ5000	2500				0.000	0.00
	5000				0.000	0.00
	500				0.000	0.00
ADJ5000	2500				0.000	0.00
	5000				0.000	0.00
	500				0.000	0.00
ADJ5000	2500				0.000	0.00
	5000				0.000	0.00
	500				0.000	0.00
ADJ5000	2500				0.000	0.00
	5000				0.000	0.00

FRM-247c (Rev 2/Mar 03)

SwRI Div. 01 – Inorganic Laboratory Adjustable Pipette Verification Log

Balance #: 34Thermometer #: G 011diH2O Temperature (°C) 21

Eppendorf #	True Value (μL)	1 st Reading (g)	2 nd Reading (g)	3 rd Reading (g)
ADJ5000-C	500	J.V. 11/11/03		
	2500			
	5000			
ADJ5000-G	500	J.V. 11/11/03		
	2500			
	5000			
ADJ5000-H	500	J.V. 11/11/03		
	2500			
	5000			
ADJ5000-I	500	J.V. 11/11/03		
	2500			
	5000			
ADJ5000-J	500	0.4975	0.4979	0.4971
	2500	2.4988	2.4983	2.4981
	5000	5.0075	5.0072	5.0061
ADJ5000-K	500	0.5011	0.5054	0.5073
	2500	2.5011	2.5081	2.5031
	5000	5.0625	5.0355	5.0335
ADJ5000-L	500	0.4994	0.4993	0.4981
	2500	2.4979	2.4963	2.4981
	5000	5.0521	5.0241	5.0341
ADJ5000	500	J.V. 11/11/03		
	2500			
	5000			
ADJ5000	500	J.V. 11/11/03		
	2500			
	5000			
ADJ5000	500	J.V. 11/11/03		
	2500			
	5000			
ADJ5000	500	J.V. 11/11/03		
	2500			
	5000			

Analyst: [Signature]Reviewed by: [Signature]Date: 11/11/03Date: 11/17/03

SwRI – Div. 01, Inorganic Labs' Adjustable Volume Pipette Verification Log

(Space provided for Inorganic Laboratories' Adjustable Volume Pipette Verification Spreadsheet)

010287

SwRI – Div. 01, Inorganic Laboratory Adjustable Pipette Verification Spreadsheet

Eppendorf #	True Value (µL)	1st Reading (g)	2nd Reading (g)	3rd Reading (g)	Avg Wt (g)	% of True Value
	500	0.4950	0.4927	0.4943	0.494	98.80
ADJ5000-C	2500	2.4796	2.4815	2.4826	2.481	99.25
	5000	5.0028	4.9968	4.9972	4.999	99.98
	500					
ADJ5000-G	2500	OUT	OF	SERVICE		
	5000					
	500	0.5086	0.5059	0.5065	0.507	101.40
ADJ5000-H	2500	2.5045	2.5054	2.5020	2.504	100.16
	5000	5.0054	4.9968	5.0064	5.003	100.06
	500	0.5002	0.5043	0.5010	0.502	100.37
ADJ5000-I	2500	2.4903	2.4916	2.4963	2.493	99.71
	5000	4.9918	4.9871	4.9931	4.991	99.81
	500					
ADJ5000-J	2500					
	5000					
	500					
ADJ5000-K	2500					
	5000					
	500					
ADJ5000-L	2500					
	5000					
	500					
ADJ5000	2500					
	5000					
	500					
ADJ5000	2500					
	5000					
	500					
ADJ5000	2500					
	5000					
	500					
ADJ5000	2500					
	5000					

FRM-247c (Rev 2/Mar 03)

SwRI Div. 01 – Inorganic Laboratory Adjustable Pipette Verification Log

010288

Balance #: 16Thermometer #: G-011diH2O Temperature (°C) 21

Eppendorf #	True Value (µL)	1 st Reading (g)	2 nd Reading (g)	3 rd Reading (g)
	500	.4950	.4927	.4943
ADJ5000-C	2500	2.4796	2.4815	2.4826
	5000	5.0028	4.9968	4.9972
	500			
ADJ5000-G	2500	out	OF	SERVICE
	5000			
	500	.5086	.5059	.5065
ADJ5000-H	2500	2.5045	2.5054	2.5020
	5000	5.0054	4.9968	5.0064
	500	.5002	.5043	.5010
ADJ5000-I	2500	2.4903	2.4916	2.4963
	5000	4.9918	4.9871	4.9931
	500			
ADJ5000-J	2500			
	5000			
	500			
ADJ5000-K	2500			
	5000			
	500			
ADJ5000-L	2500			
	5000			
	500			
ADJ5000	2500			
	5000			
	500			
ADJ5000	2500			
	5000			
	500			
ADJ5000	2500			
	5000			
	500			
ADJ5000	2500			
	5000			

500 µL – 5000 µL

ADJ 5000 G was taken out of service 11-17-03, adjust set screw is dispensing 4.5 mLs @ all settings. JW

Analyst: John WellsDate: 11-21-03Reviewed by: Valerie OlyanDate: 12/22/03

SwRI - Div. 01, Inorganic Labs' Adjustable Volume Pipette Verification Log

(Space provided for Inorganic Laboratories' Adjustable Volume Pipette Verification Spreadsheet)

010289

JV
11/21/03

SwRI - Div. 01, Inorganic Laboratory Adjustable Pipette Verification Spreadsheet

Eppendorf #	True Value (µL)	1st Reading (g)	2nd Reading (g)	3rd Reading (g)	Avg Wt (g)	% of True Value
	500				0.000	0.00
ADJ5000-C	2500				0.000	0.00
	5000				0.000	0.00
	500				0.000	0.00
ADJ5000-G	2500				0.000	0.00
	5000				0.000	0.00
	500				0.000	0.00
ADJ5000-H	2500				0.000	0.00
	5000				0.000	0.00
	500				0.000	0.00
ADJ5000-I	2500				0.000	0.00
	5000				0.000	0.00
	500				0.000	0.00
ADJ5000-J	2500	0.5003	0.5005	0.5007	0.501	100.10
	5000	2.4911	2.4895	2.4925	2.491	99.64
	500	4.9910	4.9900	4.9899	4.990	99.81
ADJ5000-K	2500	0.5025	0.5012	0.5009	0.502	100.31
	5000	2.4975	2.4985	2.4990	2.498	99.93
	500	5.0255	5.0122	5.0890	5.042	100.84
ADJ5000-L	2500	0.4992	0.4995	0.4990	0.499	99.85
	5000	2.4979	2.4985	2.4992	2.499	99.94
	500	5.0131	5.0078	5.0100	5.010	100.21
ADJ5000	2500				0.000	0.00
	5000				0.000	0.00
	500				0.000	0.00
ADJ5000	2500				0.000	0.00
	5000				0.000	0.00
	500				0.000	0.00
ADJ5000	2500				0.000	0.00
	5000				0.000	0.00
	500				0.000	0.00
ADJ5000	2500				0.000	0.00
	5000				0.000	0.00
	500				0.000	0.00

FRM-247c (Rev 2/Mar 03)

FRM-244 (Rev 2/Sept 02)

SwRI Div. 01 – Inorganic Laboratory Adjustable Pipette Verification Log **010290**

Balance #: 34

Thermometer #: G 011

diH2O Temperature (° C) 21

Eppendorf #	True Value (µL)	1 st Reading (g)	2 nd Reading (g)	3 rd Reading (g)
500 µL – 5000 µL	500	~~_____~~	~~_____~~	~~_____~~
	ADJ5000-C			
	2500			
	5000	~~_____~~	~~_____~~	~~_____~~
	500			
	ADJ5000-G			
	2500	~~_____~~	~~_____~~	~~_____~~
	5000			
	500			
	ADJ5000-H	~~_____~~	~~_____~~	~~_____~~
	2500			
	5000			
	500	~~_____~~	~~_____~~	~~_____~~
	ADJ5000-I			
	2500			
	5000	~~_____~~	~~_____~~	~~_____~~
	500			
	ADJ5000-J			
	2500	0.5003	0.5005	0.5007
	5000	2.4911	2.4895	2.4925
	5000	4.9910	4.9900	4.9899
	500	0.5025	0.5012	0.5009
	ADJ5000-K	2.4975	2.4985	2.4990
	5000	5.0255	5.0122	5.0890
	500	0.4992	0.4995	0.4990
	ADJ5000-L	2.4979	2.4985	2.4992
	5000	5.0131	5.0078	5.0100
	500	~~_____~~	~~_____~~	~~_____~~
	ADJ5000			
	2500			
	5000	~~_____~~	~~_____~~	~~_____~~
	500			
	ADJ5000			
	2500	~~_____~~	~~_____~~	~~_____~~
	5000			
	500			
	ADJ5000	~~_____~~	~~_____~~	~~_____~~
	2500			
	5000			
	500	~~_____~~	~~_____~~	~~_____~~
	ADJ5000			
	2500			
	5000	~~_____~~	~~_____~~	~~_____~~
	500			
	ADJ5000			
	2500	~~_____~~	~~_____~~	~~_____~~
	5000			
	500			

Analyst: J. Voss

Date: 11/21/03

Reviewed by: [Signature]

Date: 12/22/03

SOUTHWEST RESEARCH INSTITUTE

NUCLEAR PROJECT

CLIENT: Div. 20

TASK ORDER: 031106-12

SRR: 25200

SDG: 237578

CASE: CNWRA

VTSR: November 06, 2003

PROJECT#: 06002.01.141

010291
JG
1.2.04

Balance Calibrations

SOUTHWEST RESEARCH INSTITUTE
Division 01
BALANCE VERIFICATION LOG

BALANCE #	LAB #:	SERIAL #:	TOLERANCE:	COMMENTS:
16	28	P37987	±0.0005	
Date	Std Wt (g)	Recorded Wt (g)	Operator	
11-11-03	2.0000	2.0000	JW	SW- JS0526-15
11-12-03	2.0000	2.0000	JW	"
11-13-03	2.0000	2.0000	JW	"
11-14-03	2.0000	2.0000	JW	"
11-17-03	2.0000	2.0000	JW	"

If balance is out of limits, clean the balance and re-calibrate using Class "S" weights.
If balance is still out of limits, place a "DO NOT USE" sign on it and call (x5896) for service.

Page # 10

FRM-112 (Rev 1/Dec 99)

✓ *ama*
11/17/03

010291

Southwest Research Institute
Division 01
BALANCE VERIFICATION LOG

BALANCE #	LAB #:	SERIAL #:	TOLERANCE:	COMMENTS:
16	28	P37987	±0.0005	
Date	Std Wt (g)	Recorded Wt (g)	Operator	
11-24-03	2.0000	2.0000	Jew	SN: J50526-15
11-25-03	2.0000	2.0000	Jew	"
11-26-03	2.0000	2.0001	Jew	"
12-1-03	2.0000	2.0000	Jew	"
12-2-03	2.0000	2.0000	Jew	"
12-3-03	2.0000	2.0000	Jew	"
12-4-03	2.0000	2.0000	Jew	"
12-6-03	2.0000	2.0001	adm	SN: 99-J50526-15
12-5-03	2.0000	2.0000	Jew	"
12-8-03	2.0000	2.0000	Jew	"

If balance is out of limits, clean the balance and re-calibrate using Class "S" weights.

If balance is still out of limits, place a "DO NOT USE" sign on it and call (x5896) for service.

Page # 12

FRM-112 (Rev 1/Dec 99)

010292

Southwest Research Institute
Division 01
BALANCE VERIFICATION LOG

BALANCE #	LAB #:	SERIAL #:	TOLERANCE:	COMMENTS:
19	27	0068597	±0.05	
Date	Std Wt (g)	Recorded Wt (g)	Operator	
11-11-03	10.00	10.00	Jew	SN: 99-J50624-S
11-12-03	10.00	10.00	Jew	"
11-13-03	10.00	10.01	Jew	"
11-14-03	10.00	10.00	Jew	"
11-17-03	10.00	10.00	Jew	"
11-18-03	10.00	10.00	Jew	"
11-19-03	10.00	10.01	Jew	"
11-20-03	10.00	10.00	Jew	"
11-21-03	10.00	10.00	Jew	"
11-24-03	10.00	10.00	Jew	"

If balance is out of limits, clean the balance and re-calibrate using Class "S" weights.

If balance is still out of limits, place a "DO NOT USE" sign on it and call (x5896) for service.

Page # 6

FRM-112 (Rev 1/Dec 99)

010293

Southwest Research Institute
Division 01
BALANCE VERIFICATION LOG

BALANCE #	LAB #:	SERIAL #:	TOLERANCE:	COMMENTS:
34	28	1116031935	±0.0005	
Date	Std Wt (g)	Recorded Wt (g)	Operator	
11-11-03	2.0000	1.9999	Jew	9-2 11-11-03 SN: 99-J450526 15
11-12-03	2.0000	2.0000	Jew	"
11-13-03	2.0000	2.0000	Jew	"
11-14-03	2.0000	2.0000	Jew	"
11-17-03	2.0000	1.9999	Jew	"
11-18-03	2.0000	2.0001	Jew	"
11-19-03	2.0000	2.0000	Jew	"
11-20-03	2.0000	2.0001	Jew	"
11-21-03	2.0000	2.0000	Jew	"
11-24-03	2.0000	2.0001	Jew	"

If balance is out of limits, clean the balance and re-calibrate using Class "S" weights.

If balance is still out of limits, place a "DO NOT USE" sign on it and call (x5896) for service.

Page # 12

FRM-112 (Rev 1/Dec 99)

010294

**SOUTHWEST RESEARCH INSTITUTE
NUCLEAR PROJECT
CLIENT: Div. 20
TASK ORDER: 031106-12
SRR: 25200
SDG: 237578
CASE: CNWRA
VTSR: November 06, 2003
PROJECT#: 06002.01.141**

DI Water Verification

010295

VDS 12/22/03