RECEIVED REGION 1

'n4 222 15 80 :10

Carolyn S. Bentivegna, Ph.D. Seton Hall University 400 South Orange Ave. South Orange, NJ 07079 December 14, 2004

Sattar Lodhi, Ph. D. Nuclear Materials Safety Branch 2 Division of Nuclear Materials Safety Region I King of Prussia, Pa 19406

03000886

Dear Dr. Lodhi:

This letter is to provide documentation for the decommissioning of Dr. Friedman's laboratory in room 81 of Duffy Hall at Seton Hall University. Dr. Friedman is leaving the university and her laboratory space will be used for other purposes not involving radioisotopes.

Dr. Friedman has used ³H, ⁴⁵Ca and ³⁵S at Seton Hall. She last used ³H on 12-10-04 just before I performed the laboratory swab tests. She last used ³⁵S on 4-01-02. She last used ⁴⁵Ca on 1-20-04. All sources and waste were removed from her laboratory and put into our radioactive materials storage area under my authority.

A comprehensive swab test was run on her lab. Samples were counted for ³H, ⁴⁵Ca and ³⁵S. ⁴⁵Ca (250 keV) has a similar emission spectrum to ³⁵S (168 keV), so both isotopes were counted together with the emission window set from 0 – 250 keV. A reference for ³H and background were run with the ³H samples to verify counter efficiency. A reference for ¹⁴C and background were run for ³⁵S/⁴⁵Ca as there are no current stocks of these isotopes at Seton Hall and ¹⁴C has a similar emission, 156 keV. The scintillation counter (TRI-CARB 19000CA, serial # 102312) was professionally set-up and calibrated on 4-6-04 by Tritech Field Engineering, Inc. Efficiencies of ³H and ¹⁴C were found to be within acceptable limits for sensitivity (see calculations below). All swab tests were negative.

Enclosed are the following:

1. Map of Dr. Friedman's lab showing areas where swab tests were performed.

2. Beta counter printout of swab test for ³H including standards for ³H and Background at positions 34 and 35, respectively.

Standard H-3 UNQ S/N HMQ0212 DPM 101700 REF 31Oct01 to 31Oct06

. 15 *j*n

REGION 1

136137

Standard BKG S/N BMQ4604

REF 31Oct01 to 31Oct06

3. Beta counter printout of swab test ³⁵S/⁴⁵Ca including ¹⁴C standard and background as above in positions 34 and 35, respectively.

Standard C-14 UNQ S/N CMQ3008 DPM 98900 REF 31Oct01 to 31Oct06

- 4. Swap test for the trap in the sink where equipment coming in contact with radioisotopes was washed.
- 5. Calculations for counting efficiency.
- 6. Service record for TRI-CARB scintillation counter

In addition, there is a Beckman scintillation counter in Dr. Friedman's laboratory, Model LS-5000TD, ser # 7050993. The source was removed on Tuesday, 12-14-04 by a local representative of Beckman Coulter, Inc., Mr. Henry McManus. The Beckman Coulter, Inc. license # is 0441-30. The serial # for the source is 598860 lot #1401. A letter of receipt is pending.

Address:

Beckman Coulter

4300 N. Harbour Blvd.

Fullerton, CA 92834-3100

Telephone:

McMannus 732-432-0056

The counter without the source will be moved to:

Muhlenberg Regional Medical Center Park Avenue and Randolph Road Plainfield, NJ 07060

The counter will be under the care of:

Martin Gizzi, MD, PhD New Jersey Neuroscience Institute 69 James Street Edison, NJ 08818 (732) 321-7950


Please contact me if you have any questions about these documents: e-mail address is bentivea@shu.edu and telephone number is 973-275-2113.

Thank you for your attention,

Carolyn S. Bentivegna, Ph.D.

Radiation Safety Officer

Seton Hall University



```
Protocol #:12
                                                       10-Dec-2004 00:39
                      Name: Tritium
Region A: LL-UL= 0.0-12.0
                                     \circ
                                       Bkq= 0.00
                                                    %2 Sigma=0.00
                            Lcr=
Region B: LL-UL=12.0-156.
                             LCr=
                                     O Bkg= 0.00
                                                    %2 Sigma=0.00
Region C: LL-UL= 0.0- 0.0 Lcr=
                                       Bkg= 0.00
                                     ()
                                                    %2 Sigma=0.00
                  QIP = tSIE
Time =
        1.00
                                       ES Terminator = Count
Conventional DPM
                         Nuclide 2 =
Nuclide 1 =
                  0
                                            0
  S#
       TIME
                CPMA
                                DFM1 FLAG
                         SIS
   1
        1.00
               14.00 69.274
       1.00
               11.00 99.161
RS232 Output Format
  No Spectrum Data
  Protocol Heading
  Format: P#, PID, S#, TIME, CPMA, DPM1, CPMB, DPM2
  No IPA Data Transmission
   ...;
               14.00 54.072
       1..00
       1.00
               13.00 82.791
   4
   :::;r
       1,,00
               21.00 77.380
       1.00
               14.00 91.714
   7
       1.00
               10.00 75.182
       1.00
   8
               12.00 72.319
   9
       1.,00
               14.00 68.664
       1.00
               13.00 71.308
  10
  1.1.
       1.00
               12.00 50,214
  12
       1.00
               36.00 37.097
                                   See Swab
pusitions on map
  1.3
       1.00
               13.00 103.59
  1.4
               17.00 53.864
       i.OO
  15
       1..00
               15.00 61.525
  1.6
       1.OO
               17.00 65.755
  1.7
       1.,00
               11.00 69.907
  1.8
       1.00
               14.00 68.896
  19
       1.00
               13.00 61.621
  20
       1.00
               19.00 60.119
  21
       1,00
               14,00 56.251
  27.2
       1.00
               12.00 106.31
  23
       i.00
               17.00 64.898
       1.00
  24
               19,00 53,607
  25
               14.00 60.990
       1.00
  26
       j.,()()
               14.00 62.595
  27
       1,,00
               14.00 76.383
  28
       1.00
               12.00 70.411
  29
       1..00
               16.00 93.304
  30
       1.00
               17.00 78.702
  51
       1.,()()
               29.00 41.556
  22
       1.00
               22.00 46.094
  33
               12.00 88.007
                             - Reference Background
       1.,00
               12.00 82.287
  34
       1.00
       1.00 51137.0 20.214 - 3H Standard
  FRONT PIN JAM FWD
```

```
Protocol #: 2
                      Name:355
                                                        12-Dec-2004 22:03
Region A: LL-UL= 0.0-250.
                              Lcr=
                                      0 Bkg= 0.00
                                                      %2 Sigma=0.00
Region B: LL-UL= 4.0-250.
                              L. (C) }" ==
                                      O Bkg= 0.00
                                                      %2 Sigma=0.00
Region C: LL-UL= 0.0- 0.0 Lcr=
                                      0 Bkg= 0.00
                                                      %2 Sigma=0.00
Time = 1.00
                   GIP = SIS
  S#
       TIME
                CPMA
                         SIS FLAG
   1.
       1..00
               29.00 108.20
   \mathbb{Z}
       1.00
               27.00 35.013
   3
       1..00
               41.00 82.207
   4
               32.00 72.977
       1.00
   ET;
       1.00
               36.00 112.97
   6
       1.00
               32.00 54.068
   7
       1.00
               39,00 34,257
  8
       1.00
               28.00 102.32
  9
       1.00
               28.00 80.422
 10
               42.00 74.136
       1.00
 11
       1.00
               35,00 87,771
 1.2
       1.00
               43.00 60.579
 1.3
       1.00
               33,00 83,655
 1,4
       1.00
               34.00 78.450
 1.5
       1.00
               32.00 90.883
                                           See Swals
positions on map
 16
       1.00
              30.00 60.027
 1.7
       1.00
              21.00 103.48
 18
       1.00
              32.00 31.398
 19
       1..00
              35.00 43.656
 20
      1.00
              27.00 97.548
 21
      1.00
             706.00
                     9.337
 22
      1.00
              38.00 45.827
 23
      1.00
              35.00 70.207
 24
      1.00
              32.00 101.37
 25
      1.00
              37.00 75.392
 26
      1..00
              34.00 58.771
 27
      1..00
              37.00 51.707
 28
      1.00
              39.00 75.435
29
      1.00
              49.00 50.591
30
      1.00
              35.00 103.36
31
      1.00
              32.00 58.131
32
      1.00
             297.00 27.933
1.00
              42.00 31.794
      1.00 95507.0 159.70 — Backgrand — Backgrand — Mc Pafererse
34
35
                        15618
              98900
                95%
```

1) Ocument 4

```
Protocol #:12
                 Name:Tritium
                                               13-Dec-2004 03:00
Region A: LL-UL= 0.0-12.0 Lcr=
                               O Bkg= 0.00 %2 Sigma=0.00
Region B: LL-UL=12.0-156. Lcr=
                               0 Bkg= 0.00
                                            %2 Sigma=0.00
Region C: LL-UL= 0.0- 0.0 Lcr=
                               O Bkg= 0.00 %2 Sigma=0.00
Time = 1.00 QIP = tSIE
                                 ES Terminator = Count
Conventional DPM
Nuclide 1 =
               0
                     Nuclide 2 =
 S#
      TIME
             CPMA
                     918
                            DPM1 FLAG
  1
      1..00
            15.00 60.434 1st trap sample 3 sink #1 postini
```

2

1.00

```
Protocol #: 2
                    Name:358
                                                    i3-Dec-2004 03:04
Region A: LL-UL= 0.0-250. Lcr=
                                  0 Bkg= 0.00
                                                  %2 Sigma=0.00
Region B: LL-UL= 4.0-250. Lcr= 0 Bkg= 0.00
                                                  %2 Sigma=0.00
Region C: LL-UL= 0.0- 0.0 Lcr= 0 Bkg= 0.00 %2 Sigma=0.00
                QIP = SIS
Time = 1.00
  S#
       TIME
               CPMA
                        SIS FLAG
              36.00 95.828 1st trop 3 5.h2 #1 position 34.00 80.590 ond trap 3 5.h2 #1 position
   1
       1.00
   2
       1.00
```

TRI-CARB 1900 CA Sensitivity

Minimum acceptable efficiency for 3/+ equals € 58%.

Corrected DPM: 101700 X . 8381 (de cay)
= 85235

minimum acceptusle efficiency for "C equals 95%

$$EM = \frac{9}{0} \frac{CPM}{DPM \delta} \frac{CPM}{VC standard}$$

$$= \frac{9}{98900} \frac{96.6 2 V}{98900}$$

Acceptable SIS equals 156 ±8.
Results: 159.7

Documen #6

O ENTERED

TRITECH Field Engineering, Inc.

600 Central Avenue East
Edgewater, Maryland 21037
800-886-7004 www.tritechinc.com
Field Service Report

Call Number		Manufacțure	f	Model	*	Serial N	lumber	Engineer
040406	-001	101 Parland		1900	WC17 /023/2		12_	Low Jemetro
Date Arrived	Time	Date Completed Time		Service Type PMC Contract Warranty Reimburs.		P.O. Number Tim Brackie		
4/6/64	11:30	4/6/15	13:33	PMC				1 Daying
Institution			Department		Roc	m No.	User	
Setar	Hall		Mike / F	en H.		12		Edwal
Address				State	ZIP	Pho	one	FAX
	(Dan	1	11		┙.		
Comments:								
Rosenital								

Resolution Description		Contract	Reimb	Assist	Unit Rate	Amou	nt
Recorda U unct	Travel Hours		1.0		NU	140	a
check out Ray	Labor Hours		20		280	280	۵
	Par Descri	11	Part lumber	Qty	Unit Price		
performance charle							
			n F C	E 1 10			
C1711-95.94			(U) -	1:9 200			
16 = 63.77 Bly normal.			Ву				
Bly normal.			<u> </u>		ļ	<u> </u>	
	-			<u></u>	ļ	ļ	-
				Tota	Net	<u> </u>	-
				T	ax		
				겍	oice ount	42	000

\bigcirc . \bigcirc	-	11/1	elea
Customer Signature: - 18 300	Kist were	Date: 40	G/ 07
Customer Signature:	7	- 7	

1	This is to acknowledge the rec 12 ((4 (2004)) includes an administrative rev	eipt of your letter/application dated, and to inform you that the initial processing which iew has been performed.
	There were no administrati technical reviewer. Please omissions or require addition	-of83(-of ve omissions. Your application was assigned to a note that the technical review may identify additional onal information.
	Please provide to this office	e within 30 days of your receipt of this card
	A copy of your action has been Branch, who will contact you s	n forwarded to our License Fee & Accounts Receivable separately if there is a fee issue involved.
	Your action has been assigned When calling to inquire about You may call us on (610) 337-	this action, please refer to this control number.
	NRC FORM 532 (RI) (6-96)	Sincerely, Licensing Assistance Team Leader

	: (FOR LFMS USE) : INFORMATION FROM LTS
BETWEEN:	:
License Fee Management Branch, ARM and Regional Licensing Sections	Program Code: 03620 Status Code: 0 Fee Category: EX 3M Exp. Date: 20060228 Fee Comments: 170.11(A)(4) Decom Fin Assur Reqd: N
LICENSE FEE TRANSMITTAL	
A. REGION	
1. APPLICATION ATTACHED Applicant/Licensee: SETON HALL UNIT Received Date: 20041215 Docket No: 3000886 Control No.: 136137 License No.: 29-09831-01 Action Type: Amendment	VERSITY
2. FEE ATTACHED Amount: Check No.:	
3. COMMENTS Signed _ Date _	M.a. lerbins
B. LICENSE FEE MANAGEMENT BRANCH (Check	when milestone 03 is entered //)
1. Fee Category and Amount:	
2. Correct Fee Paid. Application may Amendment Renewal License	be processed for:
3. OTHER	
Signed	

Date