

308

Q200506060008

Scientific Notebook No 697: Alloy 22
Repassivation Test (03/28/2005 through
05/31/2005)

LABORATORY NOTEBOOK

CNWRA/SwRI

CNWRA
CONTROLLED

COPY 697

NOTEBOOK NO. _____

ISSUED TO _____

ON _____ 20____


DEPARTMENT _____

RETURNED _____ 20____

Brian K. Deedy - *B.K.D.* - BKD

Steven B. Young - *Steve Young* - SBY

DARRELL S. DUNN *Darrell Dunn DD*



SCIENTIFIC NOTEBOOK COMPANY
2831 LAWRENCE AVENUE
STEVENSVILLE, MICHIGAN 49127
(800) 537-3028 - <http://www.snco.com>

TITLE _____

From Page No. _____

Initial Scientific notebook entry for crevice corrosion repassivation potential measurements

Title: Alloy 22 Repassivation Tests

Tests Performed by: Darrell S. Dunn, Div 20; Brian Derby, Div. 18. Steve Young. Div 18

Objectives: Measure Crevice Corrosion Repassivation Potentials for Alloy 22 as a function of solution chemistry, temperature and metallurgical condition

Equipment: Keithley 614/617. Potentiostat (Solartron 1287/1480, EG&G 273, PAR 263 or equivalent) and CorrView Software or equivalent, Electrochemical test cell.

Laboratory oven for exposure of test specimens at 600 to 900 °C, Thermocouple and thermocouple meter for thermally aged specimens

Laboratory oven for exposure of test specimens at 1100 to 1350 °C, Thermocouple and thermocouple meter for solution annealed specimens

Materials: Alloy C-22, heat to be identified in each test. Other materials and heats to be added and identified prior to testing.

Specimen specifications: Specimens will be equivalent to 20.01402.571.027 or a similar crevice corrosion specimen unless otherwise specified.

Measurement Parameters: Potential and Current of specimen during test. Temperature and time of exposure for thermally aged or heat treated specimens.

Required level of accuracy: Temperature ± 5 °C, Time of exposure ± 1 minute, Potentials ± 1 mV, Current ± 0.01 microamp.

Uncertainty and Sources of Error: Current measurement error can occur for localized corrosion processes because the actively corroding area is not the same as the surface area of the test specimen.

To Page No. _____


Witnessed & Understood by me,

Date

Invented by

Date

Recorded by



3/28/05



From Page No. _____			
			To Page No. _____

Witnessed & Understood by me,

Date

Invented by

Date

4/12/05

Recorded by



From Page No. _____			
			To Page No. _____

Witnessed & Understood by me,

Date

Invented by

Date

4/12/05

Recorded by

From Page No. _____

POTENTIAL SCAN AND HOLD

Objective: Same as page 1.

Alloy / Heat No. : C-22 Heat# 2277-3-3266

Specimen Preparation: Specimen machined to dimensions specified in CNWRA Drawing. Specimen surfaces polished to 600 Grit finish using SiC paper. Specimen cleaned in acetone and rinsed in DI water. PTFE crevice forming washers attached to specimen using insulated C-276 hardware. Hardware Torque to 50 in-oz.

Torque Screwdriver: Snap-on USA SN: 1001200319
Cal: 9/3/04 Due: 3/3/05

Initial Weight: 39.88041g Model: Sartorius Genius SN: 12809099
Final Weight: 39.79926g Cal: 11/10/04 Due: 5/10/05

Solution: 4M NaCl
467.51g NaCl Lot: 042966
+DI to 2L

Reagents measured with Model: OHAUS SN: 2883
Cal: 7/16/04 1/14/05 Due: 11/16/05 7/14/05
SBY 3/28/05 SBY 3/28/05
Initial pH: 8.751 Model: Orion EA 940 SN: 2330
Final pH: 6.847 Cal: 7/21/04 Due: 7/21/05
pH Probe: #13-620-296 SN: 4079126

Test Temperature: 60°C Measured with Hg Thermometer SN: E98-191
Cal: 01/07/05 Due: 07/07/05

Counter Electrode: Platinum Flag
Reference Electrode: Fisher 13-620-51 SN: 7282317

Gas: 99.999% Nitrogen Gas
Ecorr: -0.431 V_{SCE} Model: Keithley 614 SN: 467374
Ept: +0.066 V_{SCE} Cal: 12/2/04 Due: 12/2/2005
Potentiostat: Solartron 1480 SN: 00240551

DATA FILE: C-22 galvD1

Number of Crevice Corrosion Sites: 0/24 (24 max.)

No crevice corrosion.
Slight gold staining.
*Repolish for further testing.

To Page No. _____

Witnessed & Understood by me,

Date

Invented by


Date

Recorded by

Steve Young

2/28/05

From Page No. _____



Potential, V_{SCE}

To Page No. _____

Witnessed & Understood by me,

Date

Invented by

Date

Recorded by

Steve Young

2/28/05

From Page No. _____

POTENTIAL SCAN AND HOLD

Objective: Same as page 1.

Alloy / Heat No.: C-22 Heat# 2277-3-3266

Specimen Preparation: Specimen machined to dimensions specified in CNWRA Drawing. Specimen surfaces polished to 600 Grit finish using SiC paper. Specimen cleaned in acetone and rinsed in DI water. PTFE crevice forming washers attached to specimen using insulated C-276 hardware. Hardware Torque to 50 in-oz.

Torque Screwdriver: Snap-on USA SN: 1001200319
Cal: 9/3/04 Due: 3/3/05

Initial Weight: 40.82229g Model: Sartorius Genius SN: 12809099
Final Weight: 40.80486g Cal: 11/10/04 Due: 5/10/05

Solution: 4M NaCl
467.55g NaCl Lot: 042966
+ DI to 2L

Reagents measured with Model: OHAUS SN: 2883
Cal: 7/15/04 1/14/05 Due: 1/15/05 7/14/05
SBY 3/28/05 SBY 3/28/05
Initial pH: 8.474 Model: Orion EA 940 SN: 2330
Final pH: 7.265 Cal: 7/21/04 Due: 7/21/05
pH Probe: #13-620-296 SN: 4079126

Test Temperature: 60°C Measured with Hg Thermometer SN: E98-191
Cal: 01/07/05 Due: 07/07/05

Counter Electrode: Platinum Flag

Reference Electrode: Fisher 13-620-51 SN: 7282317

Gas: 99.999% Nitrogen Gas

Ecorr: -0.313V_{SCE} Model: Keithley 614 SN: 467374
Ept: -0.504V_{SCE} Cal: 12/2/04 Due: 12/2/2005

Potentiostat: Solartron 1480 SN: 00240551

DATA FILE: C-22 galvD2

Number of Crevice Corrosion Sites: 0/24 (24 max.)

No crevice corrosion.
Slight gold tint staining.
* Repolish for further testing.

To Page No. _____

Witnessed & Understood by me,

Date

Invented by


Date

Recorded by

Steve Young

2/28/05

From Page No. _____



To Page No. _____

Witnessed & Understood by me,

Date

Invented by

Date

Recorded by

Steve Young

2/28/05

From Page No. _____

POTENTIAL SCAN AND HOLD

Objective: Same as page 1.

Alloy / Heat No. : C-22 Heat# 2277-3-3266

Specimen Preparation: Specimen machined to dimensions specified in CNWRA Drawing. Specimen surfaces polished to 600 Grit finish using SiC paper. Specimen cleaned in acetone and rinsed in DI water. PTFE crevice forming washers attached to specimen using insulated C-276 hardware. Hardware Torque to 50 in-oz.

Torque Screwdriver: Snap-on USA SN: 1001200319
Cal: 9/3/04 Due: 3/3/05

Initial Weight: 40.75490g Model: Sartorius Genius SN: 12809099
Final Weight: 40.72257g Cal: 11/10/04 Due: 5/10/05

Solution: 4M NaCl
467.52g NaCl Lot: ^{SBY 2/28/05} 467: 042966
+ DI to 2L

Reagents measured with Model: OHAUS SN: 2883
Cal: ~~7/15/04~~ 1/14/05 Due: ~~4/15/05~~ 7/14/05
SBY 3/28/05
Initial pH: 8.609 Model: Orion EA 940 SN: 2330
Final pH: 6.511 Cal: 7/21/04 Due: 7/21/05
pH Probe: #13-620-296 SN: 4079126

Test Temperature: 60°C Measured with Hg Thermometer SN: 323007
Cal: 10/14/04 Due: 10/14/2005

Counter Electrode: Platinum Flag

Reference Electrode: Fisher 13-620-52 SN: 4028036

Gas: 99.999% Nitrogen Gas

Ecorr: -0.412V_{SCE} Model: Keithley 614 SN: 467374
Ept: +0.040V_{SCE} Cal: 12/2/04 Due: 12/2/2005

Potentiostat: Solartron 1480 SN: 00240551

DATA FILE: C-22 galvD3

Number of Crevice Corrosion Sites: 0 / 24 (24 max.)

No crevice corrosion.
Slight gold staining.
* Repolish for further testing.

To Page No. _____

Witnessed & Understood by me,

Date

Invented by


Date

Recorded by

Steve Young

2/28/05

From Page No. _____



Potential, V_{SCE}

To Page No. _____

Witnessed & Understood by me,

Date

Invented by

Date

Recorded by

Steve Young

2/28/05

From Page No. _____

POTENTIAL SCAN AND HOLD

Objective: Same as page 1.

Alloy / Heat No. : C-22 Heat# 2277-3-3266

Specimen Preparation: Specimen machined to dimensions specified in CNWRA Drawing. Specimen surfaces polished to 600 Grit finish using SIC paper. Specimen cleaned in acetone and rinsed in DI water. PTFE crevice forming washers attached to specimen using insulated C-276 hardware. Hardware Torque to 50 in-oz.

Torque Screwdriver:	Snap-on USA Cal: 9/3/04	SN: 1001200319 Due: 3/3/05
Initial Weight: 40.83315g	Model: Sartorius Genius	SN: 12809099
Final Weight: 40.83355g	Cal: 11/10/04	Due: 5/10/05

Solution: 4M NaCl
467.52g NaCl Lot: 042966
+ DI to 2L

Reagents measured with	Model: OHAUS Cal: 7/15/04 1/14/05 SBY 3/28/05	SN: 2883 Due: 7/14/05 7/14/05 SBY 3/28/05
Initial pH: 8.638	Model: Orion EA 940	SN: 2330
Final pH: 7.890	Cal: 7/21/04	Due: 7/21/05
	pH Probe: #13-620-296	SN: 4079126

Test Temperature: 60°C	Measured with Hg Thermometer	SN: 323007
	Cal: 10/14/04	Due: 10/14/2005

Counter Electrode: Platinum Flag

Reference Electrode: Fisher 13-620-52 SN: 4028036

Gas: 99.999% Nitrogen Gas

Ecorr: -0.371V _{SCE}	Model: Keithley 614	SN: 467374
Ept: -0.145V _{SCE}	Cal: 12/2/04	Due: 12/2/2005

Potentiostat: Solartron 1480 SN: 00240551


DATA FILE: C-22 galvD4

Number of Crevice Corrosion Sites: /24 (24 max.)

No crevice corrosion.
Slight gold tint staining.
* Repolish for further testing.

To Page No. _____

Witnessed & Understood by me,	Date	Invented by	Date	
		Recorded by <i>Steve Young</i>	2/28/05	



Witnessed & Understood by me,	Date	Invented by	Date	
		Recorded by <i>Steve Young</i>	2/28/05	

From Page No. _____

POTENTIAL SCAN AND HOLD

Objective: Same as page 1.

Alloy / Heat No. : C-22 Heat# 2277-3-3266

Specimen Preparation: Specimen machined to dimensions specified in CNWRA Drawing. Specimen surfaces polished to 600 Grit finish using SiC paper. Specimen cleaned in acetone and rinsed in DI water. PTFE crevice forming washers attached to specimen using insulated C-276 hardware. Hardware Torque to 50 in-oz.

Torque Screwdriver: Snap-on USA SN: 1001200319 Cal: 9/3/04 Due: 3/3/05

Initial Weight: 40.21623g Model: Sartorius Genius SN: 12809099 Final Weight: 40.21541g Cal: 11/10/04 Due: 5/10/05

Solution: 4M NaCl
467.52g NaCl Lot: 042966
+ DI to 2L

Reagents measured with Model: OHAUS SN: 2883 Cal: 7/15/04 1/14/05 Due: 1/15/05 7/14/05 SBY 3/28/05 Model: Orion EA 940 SN: 2330 Due: 7/21/05 pH Probe: #13-620-296 SN: 4079126

Test Temperature: 60°C Measured with Hg Thermometer SN: E98-191 Cal: 01/07/05 Due: 07/07/05

Counter Electrode: Platinum Flag Reference Electrode: Fisher 13-620-51 SN: 7282317

Gas: 99.999% Nitrogen Gas Ecorr: -0.320V_{SCE} Model: Keithley 614 SN: 467374 Ept: -0.072V_{SCE} Cal: 12/2/04 Due: 12/2/2005 Potentiostat: Solartron 1480 SN: 00240551

DATA FILE: C-22 galvD5

Number of Crevice Corrosion Sites: 0/24 (24 max.)

No crevice corrosion.
Slight gold tint staining.
*Repolish for further testing.

To Page No. _____

Witnessed & Understood by me,

Date

Invented by

Date


Recorded by

Steve Young

2/28/05

TITLE _____

From Page No. _____



To Page No. _____

Witnessed & Understood by me,

Date

Invented by

Date

Recorded by

Steve Young

2/28/05

From Page No. _____

POTENTIAL SCAN AND HOLD

Objective: Same as page 1.

Alloy / Heat No. : C-22 Heat# 2277-3-3266

Specimen Preparation: Specimen machined to dimensions specified in CNWRA Drawing. Specimen surfaces polished to 600 Grit finish using SiC paper. Specimen cleaned in acetone and rinsed in DI water. PTFE crevice forming washers attached to specimen using insulated C-276 hardware. Hardware Torque to 50 in-oz.

Torque Screwdriver: Snap-on USA SN: 1001200319
Cal: 9/3/04 Due: 3/3/05

Initial Weight: 40.58426g Model: Sartorius Genius SN: 12809099
Final Weight: 40.58407g Cal: 11/10/04 Due: 5/10/05

Solution: 4M NaCl
467.52g NaCl Lot: 042966
+ DI to 2L

Reagents measured with Model: OHAUS SN: 2883
Cal: 7/15/04 1/14/05 Due: 7/14/05
SBY 3/28/05
Initial pH: 8.544 Model: Orion EA 940 SN: 2330
Final pH: 7.793 Cal: 7/21/04 Due: 7/21/05
pH Probe: #13-620-296 SN: 4079126

Test Temperature: 60°C Measured with Hg Thermometer SN: E98-191
Cal: 01/07/05 Due: 07/07/05

Counter Electrode: Platinum Flag

Reference Electrode: Fisher 13-620-51 SN: 7282317

Gas: 99.999% Nitrogen Gas

Ecorr: -0.273V_{SCE} Model: Keithley 614 SN: 467374
Ept: +0.170V_{SCE} Cal: 12/2/04 Due: 12/2/2005

Potentiostat: Solartron 1480 SN: 00240551

DATA FILE: C-22 galvD6

Number of Crevice Corrosion Sites: 0/24 (24 max.)

No crevice corrosion.
Light staining.
* Repolish for further testing.

To Page No. _____

Witnessed & Understood by me,

Date

Invented by


Date

Recorded by

Steve Young

3/2/05

From Page No. _____



To Page No. _____

Witnessed & Understood by me,

Date

Invented by

Date

Recorded by

Steve Young

3/2/05

From Page No. _____

POTENTIAL SCAN AND HOLD

Objective: Same as page 1.

Alloy / Heat No. : C-22 Heat# 2277-3-3266

Specimen Preparation: Specimen machined to dimensions specified in CNWRA Drawing. Specimen surfaces polished to 600 Grit finish using SiC paper. Specimen cleaned in acetone and rinsed in DI water. PTFE crevice forming washers attached to specimen using insulated C-276 hardware. Hardware Torque to 50 in-oz.

Torque Screwdriver:	Snap-on USA Cal: 9/3/04	SN: 1001200319 Due: 3/3/05
Initial Weight:	40.70585g Model: Sartorius Genius	SN: 12809099
Final Weight:	40.70547g Cal: 11/10/04	Due: 5/10/05

Solution: 0.0001 M NaCl
0.013g NaCl Lot: 042966
+ DI to 2L

Reagents measured with	Model: OHAUS Cal: 7/15/04 1/14/05 SBY 3/28/05	SN: 2883 Due: 4/15/05 7/14/05 SBY 3/28/05
Initial pH:	5.744 Model: Orion EA 940	SN: 2330
Final pH:	7.103 Cal: 7/21/04	Due: 7/21/05
	pH Probe: #13-620-296	SN: 4079126

Test Temperature:	95°C Measured with Hg Thermometer	SN: 323007
	Cal: 10/14/04	Due: 10/14/2005

Counter Electrode: Platinum Flag

Reference Electrode: Fisher 13-620-52 SN: 4028036

Gas: 99.999% Nitrogen Gas

Ecorr:	- .554V _{SCE} Model: Keithley 614	SN: 467374
Ept:	+ .102V _{SCE} Cal: 12/2/04	Due: 12/2/2005

Potentiostat: Solartron 1480 SN: 00240551


DATA FILE: C-22 Trap 21a

Number of Crevice Corrosion Sites: 0/24 (24 max.)

No crevice corrosion.
Slight gold tint.
* Repolish for further testing.

To Page No. _____

Witnessed & Understood by me,	Date	Invented by	Date
		Recorded by <i>Steve Young</i>	3/2/05



Witnessed & Understood by me,	Date	Invented by	Date
		Recorded by <i>Steve Young</i>	3/2/05

From Page No. _____

POTENTIAL SCAN AND HOLD

Objective: Same as page 1.

Alloy / Heat No. : C-22 Heat# 2277-3-3266

Specimen Preparation: Specimen machined to dimensions specified in CNWRA Drawing. Specimen surfaces polished to 600 Grit finish using SiC paper. Specimen cleaned in acetone and rinsed in DI water. PTFE crevice forming washers attached to specimen using insulated C-276 hardware. Hardware Torque to 50 in-oz.

Torque Screwdriver: Proto 6104 SN: 139072 Cal: 9/29/04 Due: 3/29/05

Initial Weight: 39.61342g Model: Sartorius Genius SN: 12809099 Cal: 11/10/04 Due: 5/10/05
Final Weight: 39.61329g

Solution: 0.0001M NaCl
0.014g NaCl lot: 042966
+ DI to 2L

Reagents measured with Model: OHAUS SN: 2883 Cal: 1/14/05 Due: 7/14/05

Initial pH: 5.92 Model: Orion EA 940 SN: 2330 Cal: 7/21/04 Due: 7/21/05
Final pH: 7.61 pH Probe: #13-620-296 SN: 4079126

Test Temperature: 95°C Measured with Hg Thermometer SN: E98-191 Cal: 01/07/05 Due: 07/07/05

Counter Electrode: Platinum Flag

Reference Electrode: Fisher 13-620-51 SN: 7282317

Gas: 99.999% Nitrogen Gas

Ecorr: -576 mV_{SCE} Model: Keithley 614 SN: 467374 Cal: 12/2/04 Due: 12/2/2005
Ept: -358 mV_{SCE}

Potentiostat: Solartron 1480 SN: 00240551

DATA FILE: C-22 Trap 21b

Number of Crevice Corrosion Sites: 0/24 (24 max.)

No crevice corrosion.
Mild staining present.
* Repolish for further testing.

To Page No. _____

Witnessed & Understood by me, _____

Date _____

Invented by _____


Date _____

Recorded by _____

3/28/05

Steve Young

From Page No. _____



To Page No. _____

Witnessed & Understood by me, _____

Date _____

Invented by _____

Date _____

Recorded by _____

3/28/05

Steve Young

From Page No. _____

POTENTIAL SCAN AND HOLD

Objective: Same as page 1.

Alloy / Heat No. : C-22 Heat# 2277-3-3266

Specimen Preparation: Specimen machined to dimensions specified in CNWRA Drawing. Specimen surfaces polished to 600 Grit finish using SiC paper. Specimen cleaned in acetone and rinsed in DI water. PTFE crevice forming washers attached to specimen using insulated C-276 hardware. Hardware Torque to 50 in-oz.

Torque Screwdriver: Snap-on USA SN: 1001200319
Cal: 9/3/04 Due: 3/3/05

Initial Weight: 40.76356g Model: Sartorius Genius SN: 12809099
Final Weight: 40.76247g Cal: 11/10/04 Due: 5/10/05

Solution: 0.0005M NaCl
0.058g NaCl Lot: 042966
+ DI to 2L

Reagents measured with Model: OHAUS SN: 2883
Cal: ~~7/15/04~~ 1/14/05 Due: ~~1/15/05~~ 7/14/05
SBY 3/28/05
Initial pH: 5.821 Model: Orion EA 940 SN: 2330
Final pH: 6.823 Cal: 7/21/04 Due: 7/21/05
pH Probe: #13-620-296 SN: 4079126

Test Temperature: 95°C Measured with Hg Thermometer SN: 323007
Cal: 10/14/04 Due: 10/14/2005

Counter Electrode: Platinum Flag
Reference Electrode: Fisher 13-620-52 SN: 4028036

Gas: 99.999% Nitrogen Gas
Ecorr: -0.555V_{SCE} Model: Keithley 614 SN: 467374
Ept: -0.273V_{SCE} Cal: 12/2/04 Due: 12/2/2005
Potentiostat: Solartron 1480 SN: 00240551

DATA FILE: C-22 Trap 22a

Number of Crevice Corrosion Sites: 0/24 (24 max.)

No crevice corrosion.
Mild surface staining.
* Repolish for further testing.

To Page No. _____

Witnessed & Understood by me,

Date

Invented by

Date


Recorded by

Steve Young

3/3/05

TITLE _____

From Page No. _____



To Page No. _____

Witnessed & Understood by me,

Date

Invented by

Date

Recorded by

Steve Young

3/3/05

From Page No. _____

POTENTIAL SCAN AND HOLD

Objective: Same as page 1.

Alloy / Heat No. : C-22 Heat# 2277-3-3266

Specimen Preparation: Specimen machined to dimensions specified in CNWRA Drawing. Specimen surfaces polished to 600 Grit finish using SiC paper. Specimen cleaned in acetone and rinsed in DI water. PTFE crevice forming washers attached to specimen using insulated C-276 hardware. Hardware Torque to 50 in-oz.

Torque Screwdriver: Snap-on USA SN: 1001200319
Cal: 9/3/04 Due: 3/3/05

Initial Weight: 41.00339g Model: Sartorius Genius SN: 12809099
Final Weight: 41.00246g Cal: 11/10/04 Due: 5/10/05

Solution: 0.0005M NaCl
0.057g NaCl Lot: 042966
+ DI to 2L

Reagents measured with Model: OHAUS SN: 2883
Cal: ~~7/16/04~~ 1/14/05 Due: ~~4/16/05~~ 7/14/05
SBY 3/28/05
Model: Orion EA 940 SN: 2330
Cal: 7/21/04 Due: 7/21/05
pH Probe: #13-620-296 SN: 4079126

Initial pH: 5.855
Final pH: 6.921

Test Temperature: 95°C Measured with Hg Thermometer SN: E98-191
Cal: 01/07/05 Due: 07/07/05

Counter Electrode: Platinum Flag
Reference Electrode: Fisher 13-620-51 SN: 7282317

Gas: 99.999% Nitrogen Gas

Ecorr: -0.535V_{SCE} Model: Keithley 614 SN: 467374
Ept: +0.217V_{SCE} Cal: 12/2/04 Due: 12/2/2005

Potentiostat: Solartron 1480 SN: 00240551


DATA FILE: C-22 Trap 22b

Number of Crevice Corrosion Sites: 0/24 (24 max.)

No crevice corrosion.
Lightly blue stained.
* Repolish for further testing.

To Page No. _____

Witnessed & Understood by me,	Date	Invented by	Date	
		Recorded by <i>Steve Young</i>	3/3/05	



Witnessed & Understood by me,	Date	Invented by	Date	
		Recorded by <i>Steve Young</i>	3/3/05	

From Page No. _____

POTENTIAL SCAN AND HOLD

Objective: Same as page 1.

Alloy / Heat No. : C-22 Heat# 2277-3-3266

Specimen Preparation: Specimen machined to dimensions specified in CNWRA Drawing. Specimen surfaces polished to 600 Grit finish using SiC paper. Specimen cleaned in acetone and rinsed in DI water. PTFE crevice forming washers attached to specimen using insulated C-276 hardware. Hardware Torque to 50 in-oz

Torque Screwdriver: Snap-on USA SN: 1001200319
Cal: 9/3/04 Due: 3/3/05

Initial Weight: 40.89905g Model: Sartorius Genius SN: 12809099
Final Weight: 40.89808g Cal: 11/10/04 Due: 5/10/05

Solution: 0.001M NaCl
0.118g NaCl Lot: 042966
+ DI to 2L

Reagents measured with Model: OHAUS SN: 2883
Cal: ~~7/15/04~~ 1/14/05 Due: ~~4/15/05~~ 7/14/05
SBY 3/28/05
Initial pH: 5.684 Model: Orion EA 940 SN: 2330
Final pH: 6.202 Cal: 7/21/04 Due: 7/21/05
pH Probe: #13-620-296 SN: 4079126

Test Temperature: 95°C Measured with Hg Thermometer SN: 323007
Cal: 10/14/04 Due: 10/14/2005

Counter Electrode: Platinum Flag

Reference Electrode: Fisher 13-620-52 SN: 4028036

Gas: 99.999% Nitrogen Gas

Ecorr: -0.605V_{SCE} Model: Keithley 614 SN: 467374
Ept: -0.082V_{SCE} Cal: 12/2/04 Due: 12/2/2005

Potentiostat: Solartron 1480 SN: 00240551

DATA FILE: C-22 Trap 23a

Number of Crevice Corrosion Sites: 0/24 (24 max.)

No crevice corrosion.
Gold and blue staining.
* Repolish for further testing.

To Page No. _____

Witnessed & Understood by me,

Date

Invented by


Date

Recorded by

Steve Young

3/4/05

From Page No. _____



Potential, V_{SCE}

Witnessed & Understood by me,

Date

Invented by

Date

Recorded by

Steve Young

3/4/05

To Page No. _____

From Page No. _____

POTENTIAL SCAN AND HOLD

Objective: Same as page 1.

Alloy / Heat No.: C-22 Heat# 2277-3-3266

Specimen Preparation: Specimen machined to dimensions specified in CNWRA Drawing. Specimen surfaces polished to 600 Grit finish using SiC paper. Specimen cleaned in acetone and rinsed in DI water. PTFE crevice forming washers attached to specimen using insulated C-276 hardware. Hardware Torque to 50 in-oz.

Torque Screwdriver: Snap-on USA SN: 1001200319 Cal: 9/3/04 Due: 3/3/05

Initial Weight: 39.71435g **Model:** Sartorius Genius SN: 12809099
Final Weight: 39.71331g Cal: 11/10/04 Due: 5/10/05

Solution: 0.001M NaCl
 0.116g NaCl Lot: 042966
 + DI to 2L

Reagents measured with Model: OHAUS SN: 2883 Cal: 7/15/04 1/14/05 Due: 4/15/05 7/14/05
 SBY 3/22/05 SBY 3/22/05
 Model: Orion EA 940 SN: 2330 Cal: 7/21/04 Due: 7/21/05
 pH Probe: #13-620-296 SN: 4079126

Initial pH: 5.872
Final pH: 6.714

Test Temperature: 95°C Measured with Hg Thermometer SN: E98-191 Cal: 01/07/05 Due: 07/07/05

Counter Electrode: Platinum Flag
Reference Electrode: Fisher 13-620-51 SN: 7282317

Gas: 99.999% Nitrogen Gas

Ecorr: -582V_{SCE} Model: Keithley 614 SN: 467374 Cal: 12/2/04 Due: 12/2/2005
Ept: +212V_{SCE}

Potentiostat: Solartron 1480 SN: 00240551


DATA FILE: C-22 Trap 23b

Number of Crevice Corrosion Sites: 0/24 (24 max.)

No crevice corrosion.
 Mild blue and gold staining.
 * Repolish for further testing.

To Page No. _____

Witnessed & Understood by me,	Date	Invented by	Date
		Recorded by <i>Steve Young</i>	3/4/05



Witnessed & Understood by me,	Date	Invented by	Date
		Recorded by <i>Steve Young</i>	3/4/05

From Page No. _____

POTENTIAL SCAN AND HOLD

Objective: Same as page 1.

Alloy / Heat No. : C-22 Heat# 2277-3-3266

Specimen Preparation: Specimen machined to dimensions specified in CNWRA Drawing. Specimen surfaces polished to 600 Grit finish using SiC paper. Specimen cleaned in acetone and rinsed in DI water. PTFE crevice forming washers attached to specimen using insulated C-276 hardware. Hardware Torque to 50 in-oz.

Torque Screwdriver: Snap-on USA SN: 1001200319
Cal: 9/3/04 Due: 3/3/05

Initial Weight: 40.69642g Model: Sartorius Genius SN: 12809099
Final Weight: 40.69443g Cal: 11/10/04 Due: 5/10/05

Solution: 0.005M NaCl
0.584g NaCl Lot: 042771
+ DI to 2L

Reagents measured with Model: OHAUS SN: 2883
Cal: ~~7/15/04~~ 1/14/05 Due: ~~4/15/05~~ 7/14/05
SBY 3/28/05 SBY 3/28/05
Initial pH: 5.474 Model: Orion EA 940 SN: 2330
Final pH: 7.001 Cal: 7/21/04 Due: 7/21/05
pH Probe: #13-620-296 SN: 4079126

Test Temperature: 95°C Measured with Hg Thermometer SN: 323007
Cal: 10/14/04 Due: 10/14/2005

Counter Electrode: Platinum Flag
Reference Electrode: Fisher 13-620-52 SN: 4028036

Gas: 99.999% Nitrogen Gas
Ecorr: -0.406V_{SCE} Model: Keithley 614 SN: 467374
Ept: +0.188V_{SCE} Cal: 12/2/04 Due: 12/2/2005
Potentiostat: Solartron 1480 SN: 00240551

DATA FILE: C-22 Trap 24a

Number of Crevice Corrosion Sites: 0/24 (24 max.)

No crevice corrosion.
Blue and gold staining.
* Repolish for further testing.

To Page No. _____

Witnessed & Understood by me,

Date

Invented by


Date

Recorded by

Steve Young

3/8/05

From Page No. _____



To Page No. _____

Witnessed & Understood by me,

Date

Invented by

Date

Recorded by

Steve Young

3/8/05

From Page No. _____

POTENTIAL SCAN AND HOLD

Objective: Same as page 1.

Alloy / Heat No. : C-22 Heat# 2277-3-3266

Specimen Preparation: Specimen machined to dimensions specified in CNWRA Drawing. Specimen surfaces polished to 600 Grit finish using SiC paper. Specimen cleaned in acetone and rinsed in DI water. PTFE crevice forming washers attached to specimen using insulated C-276 hardware. Hardware Torque to 50 in-oz.

Torque Screwdriver: Snap-on USA SN: 1001200319
Cal: 9/3/04 Due: 3/3/05

Initial Weight: 40.49937g Model: Sartorius Genius SN: 12809099
Final Weight: 40.49781g Cal: 11/10/04 Due: 5/10/05

Solution: 0.005M NaCl
0.585g NaCl Lot: 042771
+ DI to 2L

Reagents measured with Model: OHAUS SN: 2883
Cal: 7/15/04 1/14/05 Due: 1/15/06 7/14/05
SBY 3/28/05 SBY 3/28/05
Model: Orion EA 940 SN: 2330
Cal: 7/21/04 Due: 7/21/05
pH Probe: #13-620-296 SN: 4079126

Initial pH: 5.621
Final pH: 6.745

Test Temperature: 95 °C Measured with Hg Thermometer SN: E98-191
Cal: 01/07/05 Due: 07/07/05

Counter Electrode: Platinum Flag
Reference Electrode: Fisher 13-620-51 SN: 7282317

Gas: 99.999% Nitrogen Gas
Ecorr: -0.414V_{SCE} Model: Keithley 614 SN: 467374
Ept: +0.229V_{SCE} Cal: 12/2/04 Due: 12/2/2005
Potentiostat: Solartron 1480 SN: 00240551


DATA FILE: C-22 Trap 24b

Number of Crevice Corrosion Sites: 0/24 (24 max.)

No crevice corrosion.
Blue and gold staining.
* Repolish for further testing.

To Page No. _____

Witnessed & Understood by me,	Date	Invented by	Date
		Recorded by <i>Steve Young</i>	3/8/05



Witnessed & Understood by me,	Date	Invented by	Date
		Recorded by <i>Steve Young</i>	3/8/05

From Page No. _____

POTENTIAL SCAN AND HOLD

Objective: Same as page 1.

Alloy / Heat No. : C-22 Heat# 2277-3-3266

Specimen Preparation: Specimen machined to dimensions specified in CNWRA Drawing. Specimen surfaces polished to 600 Grit finish using SiC paper. Specimen cleaned in acetone and rinsed in DI water. PTFE crevice forming washers attached to specimen using insulated C-276 hardware. Hardware Torque to 50 in-oz.

Torque Screwdriver: Snap-on USA SN: 1001200319
Cal: 9/3/04 Due: 3/3/05

Initial Weight: 38.87466g Model: Sartorius Genius SN: 12809099
Final Weight: 38.87236g Cal: 11/10/04 Due: 5/10/05

Solution: 0.01M NaCl
1.169g NaCl Lot: 042771
+ DI to 2L

Reagents measured with Model: OHAUS SN: 2883
Cal: 7/15/04 1/14/05 Due: 7/15/05 7/14/05
SBY 3/23/05 SBY 3/23/05
Model: Orion EA 940 SN: 2330
Cal: 7/21/04 Due: 7/21/05
pH Probe: #13-620-296 SN: 4079126

Initial pH: 5.775
Final pH: 6.573

Test Temperature: 95°C Measured with Hg Thermometer SN: 323007
Cal: 10/14/04 Due: 10/14/2005

Counter Electrode: Platinum Flag
Reference Electrode: Fisher 13-620-52 SN: 4028036

Gas: 99.999% Nitrogen Gas

Ecorr: -0.521V_{SCE} Model: Keithley 614 SN: 467374
Ept: +0.237V_{SCE} Cal: 12/2/04 Due: 12/2/2005

Potentiostat: Solartron 1480 SN: 00240551


DATA FILE: C-22 Trap 25a

Number of Crevice Corrosion Sites: 0/24 (24 max.)

No crevice corrosion.
Multi color staining.
* Repolish for further testing.

To Page No. _____

Witnessed & Understood by me,	Date	Invented by	Date
		Recorded by <i>Steve Young</i>	3/9/05



Witnessed & Understood by me,	Date	Invented by	Date
		Recorded by <i>Steve Young</i>	3/9/05

From Page No. _____

POTENTIAL SCAN AND HOLD

Objective: Same as page 1.

Alloy / Heat No. : C-22 Heat# 2277-3-3266

Specimen Preparation: Specimen machined to dimensions specified in CNWRA Drawing. Specimen surfaces polished to 600 Grit finish using SiC paper. Specimen cleaned in acetone and rinsed in DI water. PTFE crevice forming washers attached to specimen using insulated C-276 hardware. Hardware Torque to 50 in-oz.

Torque Screwdriver:	Snap-on USA Cal: 9/3/04	SN: 1001200319 Due: 3/3/05
Initial Weight:	39.69481g	Model: Sartorius Genius SN: 12809099
Final Weight:	39.69317g	Cal: 11/10/04 Due: 5/10/05
Solution:	0.01M NaCl 1.168g NaCl Lot: 042771 + DI to 2L	

Reagents measured with	Model: OHAUS Cal: 7/45/04 1/14/05 SBY 3/28/05	SN: 2883 Due: 4/16/05 7/14/05 SBY 3/28/05
Initial pH:	5.870	Model: Orion EA 940 SN: 2330
Final pH:	6.813	Cal: 7/21/04 Due: 7/21/05
	pH Probe: #13-620-296	SN: 4079126

Test Temperature:	95°C	Measured with Hg Thermometer Cal: 01/07/05	SN: E98-191 Due: 07/07/05
-------------------	------	---	------------------------------

Counter Electrode: Platinum Flag

Reference Electrode: Fisher 13-620-51 SN: 7282317

Gas: 99.999% Nitrogen Gas

Ecorr:	-0.610V _{SCE}	Model: Keithley 614 Cal: 12/2/04	SN: 467374 Due: 12/2/2005
Ept:	+0.249V _{SCE}		

Potentiostat: Solartron 1480 SN: 00240551


DATA FILE: C-22 Trap 25b

Number of Crevice Corrosion Sites: 0/24 (24 max.)

No crevice corrosion.
Multicolor staining.
* Repolish for further testing.

To Page No. _____

Witnessed & Understood by me,	Date	Invented by	Date
		Recorded by <i>Steve Young</i>	3/9/05



Witnessed & Understood by me,	Date	Invented by	Date
		Recorded by <i>Steve Young</i>	3/9/05

From Page No. _____

POTENTIAL SCAN AND HOLD

Objective: Same as page 1.

Alloy / Heat No. : C-22 Heat# 2277-3-3266

Specimen Preparation: Specimen machined to dimensions specified in CNWRA Drawing. Specimen surfaces polished to 600 Grit finish using SiC paper. Specimen cleaned in acetone and rinsed in DI water. PTFE crevice forming washers attached to specimen using insulated C-276 hardware. Hardware Torque to 50 in-oz.

Torque Screwdriver: Snap-on USA SN: 1001200319
Cal: 9/3/04 Due: 3/3/05

Initial Weight: 40.10622g Model: Sartorius Genius SN: 12809099
Final Weight: 40.10377g Cal: 11/10/04 Due: 5/10/05

Solution: 0.05 M NaCl
5.843g NaCl Lot: 042771
+ DI to 2L

Reagents measured with Model: OHAUS SN: 2883
Cal: ~~7/15/04~~ 1/14/05 Due: ~~7/15/05~~ 7/14/05
SBY 3/28/05 SBY 3/28/05
Initial pH: 5.550 Model: Orion EA 940 SN: 2330
Final pH: 7.321 Cal: 7/21/04 Due: 7/21/05
pH Probe: #13-620-296 SN: 4079126

Test Temperature: 95°C Measured with Hg Thermometer SN: 323007
Cal: 10/14/04 Due: 10/14/2005

Counter Electrode: Platinum Flag
Reference Electrode: Fisher 13-620-52 SN: 4028036

Gas: 99.999% Nitrogen Gas
Ecorr: -0.646V_{SCE} Model: Keithley 614 SN: 467374
Ept: +0.248V_{SCE} Cal: 12/2/04 Due: 12/2/2005
Potentiostat: Solartron 1480 SN: 00240551

DATA FILE: C-22 Trap 26a

Number of Crevice Corrosion Sites: 0 / 24 (24 max.)

No crevice corrosion.
Mild gold staining.
* Repolish for further testing.

To Page No. _____

Witnessed & Understood by me,

Date

Invented by


Date

Recorded by

Steve Young

3/10/05

From Page No. _____



To Page No. _____

Witnessed & Understood by me,

Date

Invented by

Date

Recorded by

Steve Young

3/10/05

From Page No. _____

POTENTIAL SCAN AND HOLD

Objective: Same as page 1.

Alloy / Heat No. : C-22 Heat# 2277-3-3266

Specimen Preparation: Specimen machined to dimensions specified in CNWRA Drawing. Specimen surfaces polished to 600 Grit finish using SiC paper. Specimen cleaned in acetone and rinsed in DI water. PTFE crevice forming washers attached to specimen using insulated C-276 hardware. Hardware Torque to 50 in-oz.

Torque Screwdriver: Snap-on USA SN: 1001200319
Cal: 9/3/04 Due: 3/3/05

Initial Weight: 38.72336g Model: Sartorius Genius SN: 12809099
Final Weight: 38.72000g Cal: 11/10/04 Due: 5/10/05

Solution: 0.05M NaCl
5.843g NaCl Lot: 042771
+ DI to 2L

Reagents measured with Model: OHAUS SN: 2883
Cal: 7/15/04 1/14/05 Due: 4/15/05 7/14/05
SBY 3/23/05
Initial pH: 5.766 Model: Orion EA 940 SN: 2330
Final pH: 7.277 Cal: 7/21/04 Due: 7/21/05
pH Probe: #13-620-296 SN: 4079126

Test Temperature: 95° C Measured with Hg Thermometer SN: E98-191
Cal: 01/07/05 Due: 07/07/05

Counter Electrode: Platinum Flag

Reference Electrode: Fisher 13-620-51 SN: 7282317

Gas: 99.999% Nitrogen Gas

Ecorr: -0.539V_{SCE} Model: Keithley 614 SN: 467374
Ept: +0.211V_{SCE} Cal: 12/2/04 Due: 12/2/2005

Potentiostat: Solartron 1480 SN: 00240551

DATA FILE: C-22 Trap 266

Number of Crevice Corrosion Sites: 1 /24 (24 max.)

Crevice corrosion present.
Mild gold staining.

To Page No. _____

Witnessed & Understood by me,

Date

Invented by


Date

Recorded by

Steve Young

3/10/05

From Page No. _____



To Page No. _____

Witnessed & Understood by me,

Date

Invented by

Date

Recorded by

Steve Young

3/10/05

From Page No. _____

POTENTIAL SCAN AND HOLD

Objective: Same as page 1.

Alloy / Heat No. : C-22 Heat# 2277-3-3266

Specimen Preparation: Specimen machined to dimensions specified in CNWRA Drawing. Specimen surfaces polished to 600 Grit finish using SiC paper. Specimen cleaned in acetone and rinsed in DI water. PTFE crevice forming washers attached to specimen using insulated C-276 hardware. Hardware Torque to 50 in-oz.

Torque Screwdriver: Proto 6104 SN: 139072
Cal: 9/29/04 Due: 3/29/05

Initial Weight: 40.61896g Model: Sartorius Genius SN: 12809099
Final Weight: 40.61854g Cal: 11/10/04 Due: 5/10/05

Solution: 0.1M NaCl
11.688g NaCl Lot: 042771
+ DI to 2L

Reagents measured with Model: OHAUS SN: 2883
Cal: 7/15/04 1/14/05 Due: 1/15/05 7/14/05
SBY 3/23/05 SBY 3/23/05
Initial pH: 5.96 Model: Orion EA 940 SN: 2330
Final pH: 8.06 Cal: 7/21/04 Due: 7/21/05
pH Probe: #13-620-296 SN: 4079126

Test Temperature: 95°C Measured with Hg Thermometer SN: 323007
Cal: 10/14/04 Due: 10/14/2005

Counter Electrode: Platinum Flag

Reference Electrode: Fisher 13-620-52 SN: 4028036

Gas: 99.999% Nitrogen Gas

Ecorr: -0.336V_{SCE} Model: Keithley 614 SN: 467374
Ept: +0.160V_{SCE} Cal: 12/2/04 Due: 12/2/2005

Potentiostat: Solartron 1480 SN: 00240551

DATA FILE: C-22 Trap 27a

Number of Crevice Corrosion Sites: 1/24 (24 max.)

Crevice corrosion present.
Mild gold staining.

To Page No. _____

Witnessed & Understood by me,

Date

Invented by

Date


Recorded by

Steve Young

3/23/05

TITLE _____

From Page No. _____



To Page No. _____

Witnessed & Understood by me,

Date

Invented by

Date

Recorded by

Steve Young

3/23/05

From Page No. _____

POTENTIAL SCAN AND HOLD

Objective: Same as page 1.

Alloy / Heat No. : C-22 Heat# 2277-3-3266

Specimen Preparation: Specimen machined to dimensions specified in CNWRA Drawing. Specimen surfaces polished to 600 Grit finish using SiC paper. Specimen cleaned in acetone and rinsed in DI water. PTFE crevice forming washers attached to specimen using insulated C-276 hardware. Hardware Torque to 50 in-oz.

Torque Screwdriver:	Proto 6104	SN: 139072
	Cal: 9/29/04	Due: 3/29/05
Initial Weight: 40.91883g	Model: Sartorius Genius	SN: 12809099
Final Weight: 40.91780g	Cal: 11/10/04	Due: 5/10/05

Solution: 0.1M NaCl
11.688g NaCl Lot: 042771
+ DI to 2L

Reagents measured with	Model: OHAUS	SN: 2883
	Cal: 7/15/04 1/14/05	Due: 4/16/05 7/14/05
	SBY 3/28/05	SBY 3/28/05
Initial pH: 6.15	Model: Orion EA 940	SN: 2330
Final pH: 8.20	Cal: 7/21/04	Due: 7/21/05
	pH Probe: #13-620-296	SN: 4079126

Test Temperature: 95°C	Measured with Hg Thermometer	SN: E98-191
	Cal: 01/07/05	Due: 07/07/05

Counter Electrode: Platinum Flag

Reference Electrode: Fisher 13-620-51 SN: 7282317

Gas: 99.999% Nitrogen Gas

Ecorr: -0.230V _{SCE}	Model: Keithley 614	SN: 467374
Ept: -0.034V _{SCE}	Cal: 12/2/04	Due: 12/2/2005

Potentiostat: Solartron 1480 SN: 00240551


DATA FILE: C-22 Trap 27b

Number of Crevice Corrosion Sites: 1 / 24 (24 max.)

Corrosion present at 1 site.
Mild gold staining.

To Page No. _____

Witnessed & Understood by me,	Date	Invented by	Date
		Recorded by <i>Steve Young</i>	3/23/05



Witnessed & Understood by me,	Date	Invented by	Date
		Recorded by <i>Steve Young</i>	3/23/05

From Page No. _____

POTENTIAL SCAN AND HOLD

Objective: Same as page 1.

Alloy / Heat No. : C-22 Heat# 2277-3-3266

Specimen Preparation: Specimen machined to dimensions specified in CNWRA Drawing. Specimen surfaces polished to 600 Grit finish using SIC paper. Specimen cleaned in acetone and rinsed in DI water. PTFE crevice forming washers attached to specimen using insulated C-276 hardware. Hardware Torque to 50 in-oz.

Torque Screwdriver:	Snap-on USA Cal: 9/3/04	SN: 1001200319 Due: 3/3/05
Initial Weight: 40.81802	Model: Sartorius Genius	SN: 12809099
Final Weight: 40.81376	Cal: 11/10/04	Due: 5/10/05

Solution: 0.5M NaCl
58.42g NaCl Lot: 042771
+ DI to 2L

Reagents measured with	Model: OHAUS Cal: 7/15/04 1/14/05 SBY 3/28/05	SN: 2883 Due: 1/15/05 7/14/05 SBY 3/28/05
Initial pH: 6.015	Model: Orion EA 940	SN: 2330
Final pH: 8.071	Cal: 7/21/04	Due: 7/21/05
	pH Probe: #13-620-296	SN: 4079126

Test Temperature: 95°C	Measured with Hg Thermometer	SN: 323007
	Cal: 10/14/04	Due: 10/14/2005

Counter Electrode: Platinum Flag

Reference Electrode: Fisher 13-620-52 SN: 4028036

Gas: 99.999% Nitrogen Gas

Ecorr: -207V _{SCE}	Model: Keithley 614	SN: 467374
Ept: -378V _{SCE}	Cal: 12/2/04	Due: 12/2/2005

Potentiostat: Solartron 1480 SN: 00240551


DATA FILE: C-22 Trap 28a

Number of Crevice Corrosion Sites: 2 / 24 (24 max.)

Crevice corrosion present.
Mild gold staining.

To Page No. _____

Witnessed & Understood by me,	Date	Invented by	Date
		Recorded by <i>Steve Young</i>	3/21/05



To Page No. _____

Witnessed & Understood by me,	Date	Invented by	Date
		Recorded by <i>Steve Young</i>	3/21/05

From Page No. _____

POTENTIAL SCAN AND HOLD

Objective: Same as page 1.

Alloy / Heat No. : C-22 Heat# 2277-3-3266

Specimen Preparation: Specimen machined to dimensions specified in CNWRA Drawing. Specimen surfaces polished to 600 Grit finish using SiC paper. Specimen cleaned in acetone and rinsed in DI water. PTFE crevice forming washers attached to specimen using insulated C-276 hardware. Hardware Torque to 50 in-oz.

Torque Screwdriver: Proto 6104 SN: 139072 Cal: 9/29/04 Due: 3/29/05

Initial Weight: 40.63081g Model: Sartorius Genius SN: 12809099
Final Weight: 40.61654g Cal: 11/10/04 Due: 5/10/05

Solution: 0.5M NaCl
58.45g NaCl Lot: 042771
+ DI to 2L

Reagents measured with Model: OHAUS SN: 2883 Cal: 1/14/05 Due: 7/14/05

Initial pH: 5.80 Model: Orion EA 940 SN: 2330
Final pH: 8.21 Cal: 7/21/04 Due: 7/21/05
pH Probe: #13-620-296 SN: 4079126

Test Temperature: 95°C Measured with Hg Thermometer SN: E98-191 Cal: 01/07/05 Due: 07/07/05

Counter Electrode: Platinum Flag

Reference Electrode: Fisher 13-620-51 SN: 7282317

Gas: 99.999% Nitrogen Gas

Ecorr: -414 mV_{SCE} Model: Keithley 614 SN: 467374
Ept: -172 mV_{SCE} Cal: 12/2/04 Due: 12/2/2005

Potentiostat: Solartron 1480 SN: 00240551

DATA FILE: C-22 Trap 28b

Number of Crevice Corrosion Sites: 2 / 24 (24 max.)

Two corrosion sites.
Mild gold staining.

To Page No. _____

Witnessed & Understood by me,

Date

Invented by


Date

Recorded by

Steve Young

3/28/05

From Page No. _____



To Page No. _____

Witnessed & Understood by me,

Date

Invented by

Date

Recorded by

Steve Young

3/28/05

From Page No. _____

POTENTIAL SCAN AND HOLD

Objective: Same as page 1.

Alloy / Heat No. : C-22 Heat# 2277-3-3266

Specimen Preparation: Specimen machined to dimensions specified in CNWRA Drawing. Specimen surfaces polished to 600 Grit finish using SiC paper. Specimen cleaned in acetone and rinsed in DI water. PTFE crevice forming washers attached to specimen using insulated C-276 hardware. Hardware Torque to 50 in-oz.

Torque Screwdriver: Proto 6104 SN: 139072 Cal: 9/29/04 Due: 3/29/05

Initial Weight: 40.49546g Model: Sartorius Genius SN: 12809099
Final Weight: 40.49171g Cal: 11/10/04 Due: 5/10/05

Solution: 1M NaCl
116.88g NaCl Lot: 037047A
+ DI to 2L

Reagents measured with Model: OHAUS SN: 2883 Cal: 1/14/05 Due: 7/14/05

Initial pH: 7.03 Model: Orion EA 940 SN: 2330
Final pH: 7.99 Cal: 7/21/04 Due: 7/21/05
pH Probe: #13-620-296 SN: 4079126

Test Temperature: 95°C Measured with Hg Thermometer SN: 323007 Cal: 10/14/04 Due: 10/14/2005

Counter Electrode: Platinum Flag

Reference Electrode: Fisher 13-620-52 SN: 4028036

Gas: 99.999% Nitrogen Gas

Ecorr: -317 mV_{SCE} Model: Keithley 614 SN: 467374
Ept: +44.5 mV_{SCE} Cal: 12/2/04 Due: 12/2/2005

Potentiostat: Solartron 1480 SN: 00240551

DATA FILE: C-22 Trap 29a

Number of Crevice Corrosion Sites: 1 / 24 (24 max.)

One crevice corrosion site.
Mild gold staining.

To Page No. _____

Witnessed & Understood by me,

Date

Invented by

Date


Recorded by

Steve Young

3/29/05

TITLE _____

From Page No. _____



To Page No. _____

Witnessed & Understood by me,

Date

Invented by

Date

Recorded by

Steve Young

3/29/05

From Page No. _____

POTENTIAL SCAN AND HOLD

Objective: Same as page 1.

Alloy / Heat No. : C-22 Heat# 2277-3-3266

Specimen Preparation: Specimen machined to dimensions specified in CNWRA Drawing. Specimen surfaces polished to 600 Grit finish using SiC paper. Specimen cleaned in acetone and rinsed in DI water. PTFE crevice forming washers attached to specimen using insulated C-276 hardware. Hardware Torque to 50 in-oz.

Torque Screwdriver: Proto 6104 SN: 139072 Cal: 9/29/04 Due: 3/29/05

Initial Weight: 39.96277g Model: Sartorius Genius SN: 12809099
Final Weight: 39.95558g Cal: 11/10/04 Due: 5/10/05

Solution: 1M NaCl
116.89g NaCl Lot: 037047A
+ DI to 2L

Reagents measured with Model: OHAUS SN: 2883 Cal: 1/14/05 Due: 7/14/05

Initial pH: 7.00 Model: Orion EA 940 SN: 2330
Final pH: 7.97 Cal: 7/21/04 Due: 7/21/05
pH Probe: #13-620-296 SN: 4079126

Test Temperature: 95°C Measured with Hg Thermometer SN: 323007
Cal: 10/14/04 Due: 10/14/2005

Counter Electrode: Platinum Flag

Reference Electrode: Fisher 13-620-52 SN: 4028036

Gas: 99.999% Nitrogen Gas

Ecorr: -431 mV_{SCE} Model: Keithley 614 SN: 467374
Ept: -49.3 mV_{SCE} Cal: 12/2/04 Due: 12/2/2005

Potentiostat: Solartron 1480 SN: 00240551

DATA FILE: C-22 Trap 29b


Number of Crevice Corrosion Sites: 2 / 24 (24 max.)

Two crevice corrosion sites present.
Mild gold staining.

To Page No. _____

Witnessed & Understood by me,	Date	Invented by	Date
		Recorded by <i>Steve Young</i>	3/29/05

From Page No. _____



To Page No. _____

Witnessed & Understood by me,	Date	Invented by	Date
		Recorded by <i>Steve Young</i>	3/29/05

From Page No. _____

POTENTIAL SCAN AND HOLD

Objective: Same as page 1.

Alloy / Heat No. : C-22 Heat# 2277-3-3266

Specimen Preparation: Specimen machined to dimensions specified in CNWRA Drawing. Specimen surfaces polished to 600 Grit finish using SiC paper. Specimen cleaned in acetone and rinsed in DI water. PTFE crevice forming washers attached to specimen using insulated C-276 hardware. Hardware Torque to 50 in-oz.

Torque Screwdriver: Proto 6104 SN: 139072 Cal: 9/29/04 Due: 3/29/05

Initial Weight: 39.62373g Model: Sartorius Genius SN: 12809099
Final Weight: 39.45674g Cal: 11/10/04 Due: 5/10/05

Solution: 4M NaCl
467.52g NaCl lot: 037047A
+ DI to 2L

Reagents measured with Model: OHAUS SN: 2883 Cal: 1/14/05 Due: 7/14/05

Initial pH: 8.91 Model: Orion EA 940 SN: 2330
Final pH: 7.63 Cal: 7/21/04 Due: 7/21/05
pH Probe: #13-620-296 SN: 4079126

Test Temperature: 95°C Measured with Hg Thermometer SN: 323007 Cal: 10/14/04 Due: 10/14/2005

Counter Electrode: Platinum Flag

Reference Electrode: Fisher 13-620-52 SN: 4028036

Gas: 99.999% Nitrogen Gas

Ecorr: -429mV_{SCE} Model: Keithley 614 SN: 467374
Ept: -120mV_{SCE} Cal: 12/2/04 Due: 12/2/2005

Potentiostat: Solartron 1480 SN: 00240551

DATA FILE: C-22 Trap 30a

Number of Crevice Corrosion Sites: 7 / 24 (24 max.)

Seven corrosion sites.
Mild staining.

To Page No. _____

Witnessed & Understood by me,

Date

Invented by


Date

Recorded by

Steve Young

3/30/05

From Page No. _____



To Page No. _____

Witnessed & Understood by me,

Date

Invented by

Date

Recorded by

Steve Young

3/30/05

From Page No. _____

POTENTIAL SCAN AND HOLD

Objective: Same as page 1.

Alloy / Heat No. : C-22 Heat# 2277-3-3266

Specimen Preparation: Specimen machined to dimensions specified in CNWRA Drawing. Specimen surfaces polished to 600 Grit finish using SiC paper. Specimen cleaned in acetone and rinsed in DI water. PTFE crevice forming washers attached to specimen using insulated C-276 hardware. Hardware Torque to 50 in-oz.

Torque Screwdriver: Proto 6104 SN: 139072 Cal: 9/29/04 Due: 3/29/05

Initial Weight: 39.81386g Model: Sartorius Genius SN: 12809099
Final Weight: 39.65486g Cal: 11/10/04 Due: 5/10/05

Solution: 4M NaCl
467.51g NaCl Lot: 037047A
+ DI to 2L

Reagents measured with Model: OHAUS SN: 2883 Cal: 1/14/05 Due: 7/14/05

Initial pH: 8.92 Model: Orion EA 940 SN: 2330
Final pH: 7.57 Cal: 7/21/04 Due: 7/21/05
pH Probe: #13-620-296 SN: 4079126

Test Temperature: 95°C Measured with Hg Thermometer SN: E98-191 Cal: 01/07/05 Due: 07/07/05

Counter Electrode: Platinum Flag

Reference Electrode: Fisher 13-620-51 SN: 7282317

Gas: 99.999% Nitrogen Gas

Ecorr: -461mV_{SCE} Model: Keithley 614 SN: 467374
Ept: +9.09mV_{SCE} Cal: 12/2/04 Due: 12/2/2005

Potentiostat: Solartron 1480 SN: 00240551

DATA FILE: C-22 Trap 30b

Number of Crevice Corrosion Sites: 20 /24 (24 max.)


Twenty corrosion sites.
Mild staining.

To Page No. _____

Witnessed & Understood by me,	Date	Invented by	Date
		Recorded by <i>Steve Young</i>	3/30/05

TITLE _____

From Page No. _____



To Page No. _____

Witnessed & Understood by me,	Date	Invented by	Date
		Recorded by <i>Steve Young</i>	3/30/05

From Page No. _____

POTENTIAL SCAN AND HOLD

Objective: Same as page 1.

Alloy / Heat No. : C-22 Heat# 2277-3-3266

Specimen Preparation: Specimen machined to dimensions specified in CNWRA Drawing. Specimen surfaces polished to 600 Grit finish using SiC paper. Specimen cleaned in acetone and rinsed in DI water. PTFE crevice forming washers attached to specimen using insulated C-276 hardware. Hardware Torque to 50 in-oz.

Torque Screwdriver: Proto 6104 SN: 139072 Due: 3/29/05
Cal: 9/29/04

Initial Weight: 39.57235 g Model: Sartorius Genius SN: 12809099
Final Weight: 39.56772 g Cal: 11/10/04 Due: 5/10/05

Solution: 5.5M NaCl
642.84g NaCl Lot: 042771
+ DI to 2L

Reagents measured with Model: OHAUS SN: 2883 Due: 7/14/05
Cal: 1/14/05

Initial pH: 8.60 Model: Orion EA 940 SN: 2330
Final pH: 7.87 Cal: 7/21/04 Due: 7/21/05
pH Probe: #13-620-296 SN: 4079126

Test Temperature: 60°C Measured with Hg Thermometer SN: 323007 Due: 10/14/2005
Cal: 10/14/04

Counter Electrode: Platinum Flag

Reference Electrode: Fisher 13-620-52 SN: 4028036

Gas: 99.999% Nitrogen Gas

Ecorr: -458 mV_{SCE} Model: Keithley 614 SN: 467374
Ept: -3.42 mV_{SCE} Cal: 12/2/04 Due: 12/2/2005

Potentiostat: Solartron 1480 SN: 00240551


DATA FILE: C-22 Trap 31a

Number of Crevice Corrosion Sites: 0/24 (24 max.)

No crevice corrosion.
Mild gold staining.
*Repolish for further testing.

To Page No. _____

Witnessed & Understood by me,	Date	Invented by	Date
		Recorded by <i>Steve Young</i>	3/31/05



Witnessed & Understood by me,	Date	Invented by	Date
		Recorded by <i>Steve Young</i>	3/31/05

From Page No. _____

POTENTIAL SCAN AND HOLD

Objective: Same as page 1.

Alloy / Heat No. : C-22 Heat# 2277-3-3266

Specimen Preparation: Specimen machined to dimensions specified in CNWRA Drawing. Specimen surfaces polished to 600 Grit finish using SiC paper. Specimen cleaned in acetone and rinsed in DI water. PTFE crevice forming washers attached to specimen using insulated C-276 hardware. Hardware Torque to 50 in-oz.

Torque Screwdriver: Proto 6104 SN: 139072
Cal: 9/29/04 Due: 3/29/05

Initial Weight: 38.90124g Model: Sartorius Genius SN: 12809099
Final Weight: 38.89820g Cal: 11/10/04 Due: 5/10/05

Solution: 5.5M NaCl
642.85g NaCl Lot: 042771
+ DI to 2L

Reagents measured with Model: OHAUS SN: 2883
Cal: 1/14/05 Due: 7/14/05

Initial pH: 8.56 Model: Orion EA 940 SN: 2330
Final pH: 7.84 Cal: 7/21/04 Due: 7/21/05
pH Probe: #13-620-296 SN: 4079126

Test Temperature: 60°C Measured with Hg Thermometer SN: E98-191
Cal: 01/07/05 Due: 07/07/05

Counter Electrode: Platinum Flag

Reference Electrode: Fisher 13-620-51 SN: 7282317

Gas: 99.999% Nitrogen Gas

Ecorr: -527mV_{SCE} Model: Keithley 614 SN: 467374
Ept: +25mV_{SCE} Cal: 12/2/04 Due: 12/2/2005

Potentiostat: Solartron 1480 SN: 00240551

DATA FILE: C-22 Trap 31b

Number of Crevice Corrosion Sites: 3 /24 (24 max.)

Three corrosion sites.
Mild staining.

To Page No. _____

Witnessed & Understood by me,

Date

Invented by


Date

Recorded by

Steve Young

3/31/05

From Page No. _____



To Page No. _____

Witnessed & Understood by me,

Date

Invented by

Date

Recorded by

Steve Young

3/31/05

From Page No. _____

POTENTIAL SCAN AND HOLD

Objective: Same as page 1.

Alloy / Heat No. : C-22 Heat# 2277-3-3266

Specimen Preparation: Specimen machined to dimensions specified in CNWRA Drawing. Specimen surfaces polished to 600 Grit finish using SiC paper. Specimen cleaned in acetone and rinsed in DI water. PTFE crevice forming washers attached to specimen using insulated C-276 hardware. Hardware Torque to 50 in-oz.

Torque Screwdriver: Proto 6104 SN: 139072 Due: 3/29/05
Cal: 9/29/04

Initial Weight: 38.73905g Model: Sartorius Genius SN: 12809099
Final Weight: 38.73351g Cal: 11/10/04 Due: 5/10/05

Solution: 5.5M NaCl
642.83g NaCl Lot: 042771
+ DI to 2L

Reagents measured with Model: OHAUS SN: 2883 Due: 7/14/05
Cal: 1/14/05

Initial pH: 8.34 Model: Orion EA 940 SN: 2330 Due: 7/21/05
Final pH: 7.61 Cal: 7/21/04 SN: 4079126
pH Probe: #13-620-296

Test Temperature: 80°C Measured with Hg Thermometer SN: 323007 Due: 10/14/2005
Cal: 10/14/04

Counter Electrode: Platinum Flag

Reference Electrode: Fisher 13-620-52 SN: 4028036

Gas: 99.999% Nitrogen Gas

Ecorr: -360mV_{SCE} Model: Keithley 614 SN: 467374 Due: 12/2/2005
Ept: +150mV_{SCE} Cal: 12/2/04

Potentiostat: Solartron 1480 SN: 00240551

DATA FILE: C-22 Trap 32a

Number of Crevice Corrosion Sites: 5 /24 (24 max.)

5 corrosion sites.
Mild staining.

To Page No. _____

Witnessed & Understood by me,

Date

Invented by


Date

Recorded by

Steve Young

4/8/05

From Page No. _____



To Page No. _____

Witnessed & Understood by me,

Date

Invented by

Date

Recorded by

Steve Young

4/8/05

From Page No. _____

POTENTIAL SCAN AND HOLD

Objective: Same as page 1.

Alloy / Heat No. : C-22 Heat# 2277-3-3266

Specimen Preparation: Specimen machined to dimensions specified in CNWRA Drawing. Specimen surfaces polished to 600 Grit finish using SiC paper. Specimen cleaned in acetone and rinsed in DI water. PTFE crevice forming washers attached to specimen using insulated C-276 hardware. Hardware Torque to 50 in-oz.

Torque Screwdriver: Proto 6104 SN: 139072 Cal: 9/29/04 Due: 3/29/05

Initial Weight: 39.69039g Model: Sartorius Genius SN: 12809099
Final Weight: 39.67656g Cal: 11/10/04 Due: 5/10/05

Solution: 5.5M NaCl
 642.83g NaCl Lot: 042771
 + DI to 2L

Reagents measured with Model: OHAUS SN: 2883 Cal: 1/14/05 Due: 7/14/05

Initial pH: 8.53 Model: Orion EA 940 SN: 2330
Final pH: 7.68 Cal: 7/21/04 Due: 7/21/05
 pH Probe: #13-620-296 SN: 4079126

Test Temperature: 80° C Measured with Hg Thermometer SN: E98-191 Cal: 01/07/05 Due: 07/07/05

Counter Electrode: Platinum Flag

Reference Electrode: Fisher 13-620-51 SN: 7282317

Gas: 99.999% Nitrogen Gas

Ecorr: -280mV_{SCE} Model: Keithley 614 SN: 467374
Ept: +110mV_{SCE} Cal: 12/2/04 Due: 12/2/2005

Potentiostat: Solartron 1480 SN: 00240551

DATA FILE: C-22 Trap 32b

Number of Crevice Corrosion Sites: 10 /24 (24 max.)

Crevice corrosion present.
Mild staining.

To Page No. _____

Witnessed & Understood by me,	Date	Invented by	Date	
		Recorded by <i>Steve Young</i>	4/8/05	

From Page No. _____

To Page No. _____

Witnessed & Understood by me,	Date	Invented by	Date	
		Recorded by <i>Steve Young</i>	4/8/05	

From Page No. _____

POTENTIAL SCAN AND HOLD

Objective: Same as page 1.

Alloy / Heat No. : C-22 Heat# 2277-3-3266

Specimen Preparation: Specimen machined to dimensions specified in CNWRA Drawing. Specimen surfaces polished to 600 Grit finish using SiC paper. Specimen cleaned in acetone and rinsed in DI water. PTFE crevice forming washers attached to specimen using insulated C-276 hardware. Hardware Torque to 50 in-oz.

Torque Screwdriver: Proto 6104 SN: 139072
Cal: ~~9/29/04~~ 4/11/05 Due: ~~3/20/05~~ 10/11/05
SBY 4/21/05
Model: Sartorius Genius SN: 12809099
Initial Weight: 39.49352g
Final Weight: 39.38727g Cal: ~~11/10/04~~ 11/10/04 Due: 5/10/05
SBY 4/21/05
Solution: 5.5M NaCl
642.85g NaCl Lot: 042771
+ DI to 2L

Reagents measured with Model: OHAUS SN: 2883
Cal: 1/14/05 Due: 7/14/05
Initial pH: 8.17 Model: Orion EA 940 SN: 2330
Final pH: 7.32 Cal: 7/21/04 Due: 7/21/05
pH Probe: #13-620-296 SN: 4079126

Test Temperature: 95°C Measured with Hg Thermometer SN: 323007
Cal: 10/14/04 Due: 10/14/2005

Counter Electrode: Platinum Flag
Reference Electrode: Fisher 13-620-52 SN: 4028036

Gas: 99.999% Nitrogen Gas
Ecorr: -222 mV_{SCE} Model: Keithley 614 SN: 467374
Ept: +268 mV_{SCE} Cal: 12/2/04 Due: 12/2/2005
Potentiostat: Solartron 1480 SN: 00240551

DATA FILE: C-22 Trap 33a

Number of Crevice Corrosion Sites: 15 /24 (24 max.)


Crevice corrosion present.
Mild staining.

To Page No. _____

Witnessed & Understood by me,	Date	Invented by	Date
	4/12/05	Recorded by Steve Young	4/12/05

TITLE _____

From Page No. _____



To Page No. _____

Witnessed & Understood by me,	Date	Invented by	Date
		Recorded by Steve Young	4/12/05

From Page No. _____

POTENTIAL SCAN AND HOLD

Objective: Same as page 1.

Alloy / Heat No. : C-22 Heat# 2277-3-3266

Specimen Preparation: Specimen machined to dimensions specified in CNWRA Drawing. Specimen surfaces polished to 600 Grit finish using SiC paper. Specimen cleaned in acetone and rinsed in DI water. PTFE crevice forming washers attached to specimen using insulated C-276 hardware. Hardware Torque to 50 in-oz.

Torque Screwdriver: Proto 6104 SN: 139072
Cal: 9/29/04 4/11/05 Due: 3/29/05 10/11/05
SBY 4/21/05
Model: Sartorius Genius SN: 12809099
Initial Weight: 40.41751g Due: 5/10/05
Final Weight: 40.38564g Cal: 11/10/04
Solution: 5.5M NaCl
642.83g NaCl Lot: 042771
+ DI to 2L

Reagents measured with Model: OHAUS SN: 2883
Cal: 1/14/05 Due: 7/14/05
Initial pH: 8.32 Model: Orion EA 940 SN: 2330
Final pH: 7.41 Cal: 7/21/04 Due: 7/21/05
pH Probe: #13-620-296 SN: 4079126

Test Temperature: 95 °C Measured with Hg Thermometer SN: E98-191
Cal: 01/07/05 Due: 07/07/05

Counter Electrode: Platinum Flag
Reference Electrode: Fisher 13-620-51 SN: 7282317

Gas: 99.999% Nitrogen Gas
Ecorr: -376 mV_{SCE} Model: Keithley 614 SN: 467374
Ept: +109 mV_{SCE} Cal: 12/2/04 Due: 12/2/2005
Potentiostat: Solartron 1480 SN: 00240551

DATA FILE: C-22 Trap 33b

Number of Crevice Corrosion Sites: 10 /24 (24 max.)

Crevice corrosion present.
Mild staining.

To Page No. _____

Witnessed & Understood by me,

Date

Invented by


Date

Recorded by

4/12/05

Steve Young

From Page No. _____



To Page No. _____

Witnessed & Understood by me,

Date

Invented by

Date

Recorded by

4/12/05

Steve Young

From Page No. _____

POTENTIAL SCAN AND HOLD

Objective: Same as page 1.

Alloy / Heat No.: C-22 Heat# 2277-3-3266

Specimen Preparation: Specimen machined to dimensions specified in CNWRA Drawing. Specimen surfaces polished to 600 Grit finish using SiC paper. Specimen cleaned in acetone and rinsed in DI water. PTFE crevice forming washers attached to specimen using insulated C-276 hardware. Hardware Torque to 50 in-oz.

Torque Screwdriver: Proto 6104 SN: 139072
Cal: 9/29/04 4/11/05 Due: 2/29/05 10/11/05
SBY 4/21/05
Model: Sartorius Genius SN: 12809099
Initial Weight: 39.46514g
Final Weight: 39.46365g Cal: 11/10/04 Due: 5/10/05

Solution: 3.5M MgCl₂ (7M Cl)
1423.20g MgCl₂ Lot: 041703
+ DI to 2L

Reagents measured with Model: OHAUS SN: 2883
Cal: 1/14/05 Due: 7/14/05
Initial pH: 3.30 Model: Orion EA 940 SN: 2330
Final pH: 5.14 Cal: 7/21/04 Due: 7/21/05
pH Probe: #13-620-296 SN: 4079126

Test Temperature: 60°C Measured with Hg Thermometer SN: 323007
Cal: 10/14/04 Due: 10/14/2005

Counter Electrode: Platinum Flag
Reference Electrode: Fisher 13-620-52 SN: 4028036

Gas: 99.999% Nitrogen Gas
Ecorr: -183 mV_{SCE} Model: Keithley 614 SN: 467374
Ept: +450 mV_{SCE} Cal: 12/2/04 Due: 12/2/2005
Potentiostat: Solartron 1480 SN: 00240551

DATA FILE: C-22 Trap 34a

Number of Crevice Corrosion Sites: 5 /24 (24 max.)

Five corrosion sites.
No staining.

To Page No. _____

Witnessed & Understood by me,

Date

Invented by


Date

Recorded by

Steve Young

4/14/05

From Page No. _____



To Page No. _____

Witnessed & Understood by me,

Date

Invented by

Date

Recorded by

Steve Young

4/14/05

From Page No. _____

POTENTIAL SCAN AND HOLD

Objective: Same as page 1.

Alloy / Heat No. : C-22 Heat# 2277-3-3266

Specimen Preparation: Specimen machined to dimensions specified in CNWRA Drawing. Specimen surfaces polished to 600 Grit finish using SiC paper. Specimen cleaned in acetone and rinsed in DI water. PTFE crevice forming washers attached to specimen using insulated C-276 hardware. Hardware Torque to 50 in-oz.

Torque Screwdriver: Proto 6104 SN: 139072
Cal: 9/29/04 4/11/05 Due: 3/29/05 10/11/05
SBY 4/21/05
Initial Weight: 39.50691 g Model: Sartorius Genius SN: 12809099
Final Weight: 39.50549 g Cal: 11/10/04 Due: 5/10/05

Solution: 3.5M MgCl₂ (7M Cl)
1423.17g MgCl₂ Lot: 041703
+ DI to 2L

Reagents measured with Model: OHAUS SN: 2883
Cal: 1/14/05 Due: 7/14/05

Initial pH: 3.71 Model: Orion EA 940 SN: 2330
Final pH: 5.14 Cal: 7/21/04 Due: 7/21/05
pH Probe: #13-620-296 SN: 4079126

Test Temperature: 60 °C Measured with Hg Thermometer SN: E98-191
Cal: 01/07/05 Due: 07/07/05

Counter Electrode: Platinum Flag

Reference Electrode: Fisher 13-620-51 SN: 7282317

Gas: 99.999% Nitrogen Gas

Ecorr: -211 mV_{SCE} Model: Keithley 614 SN: 467374
Ept: +97 mV_{SCE} Cal: 12/2/04 Due: 12/2/2005

Potentiostat: Solartron 1480 SN: 00240551

DATA FILE: C-22 Trap 34b

Number of Crevice Corrosion Sites: 4 /24 (24 max.)

Four corrosion sites.
No staining.

To Page No. _____

Witnessed & Understood by me,

Date

Invented by

Date


Recorded by

Steve Young

4/14/05

TITLE _____

From Page No. _____



To Page No. _____

Witnessed & Understood by me,

Date

Invented by

Date

Recorded by

Steve Young

4/14/05

From Page No. _____

POTENTIAL SCAN AND HOLD

Objective: Same as page 1.

Alloy / Heat No. : C-22 Heat# 2277-3-3266

Specimen Preparation: Specimen machined to dimensions specified in CNWRA Drawing. Specimen surfaces polished to 600 Grit finish using SiC paper. Specimen cleaned in acetone and rinsed in DI water. PTFE crevice forming washers attached to specimen using insulated C-276 hardware. Hardware Torque to 50 in-oz.

Torque Screwdriver: Proto 6104 SN: 139072
Cal: 9/29/04 4/11/05 SBY 4/21/05 Due: 3/29/05 10/11/05
Initial Weight: 40.25084 g Model: Sartorius Genius SN: 12809099
Final Weight: 40.21571 g Cal: 11/10/04 Due: 5/10/05

Solution: 5M MgCl₂ (10M Cl)
2033.05g MgCl₂ · 6H₂O Lot: 042821
+ DI to 2L

Reagents measured with Model: OHAUS SN: 2883
Cal: 1/14/05 Due: 7/14/05
Initial pH: 2.28 Model: Orion EA 940 SN: 2330
Final pH: 4.05 Cal: 7/21/04 Due: 7/21/05
pH Probe: #13-620-296 SN: 4079126

Test Temperature: 60°C Measured with Hg Thermometer SN: 323007
Cal: 10/14/04 Due: 10/14/2005

Counter Electrode: Platinum Flag
Reference Electrode: Fisher 13-620-52 SN: 4028036

Gas: 99.999% Nitrogen Gas
Ecorr: -149 mV_{SCE} Model: Keithley 614 SN: 467374
Ept: +475 mV_{SCE} Cal: 12/2/04 Due: 12/2/2005
Potentiostat: Solartron 1480 SN: 00240551

DATA FILE: C-22 Trap 35a


Number of Crevice Corrosion Sites: 20/24 (24 max.)

No corrosion under crevice washers,
only outside edges.
No staining.

To Page No. _____

Witnessed & Understood by me,	Date	Invented by	Date
		Recorded by <i>Steve Young</i>	4/19/05

From Page No. _____



To Page No. _____

Witnessed & Understood by me,	Date	Invented by	Date
		Recorded by <i>Steve Young</i>	4/19/05

From Page No. _____

POTENTIAL SCAN AND HOLD

Objective: Same as page 1.

Alloy / Heat No. : C-22 Heat# 2277-3-3266

Specimen Preparation: Specimen machined to dimensions specified in CNWRA Drawing. Specimen surfaces polished to 600 Grit finish using SiC paper. Specimen cleaned in acetone and rinsed in DI water. PTFE crevice forming washers attached to specimen using insulated C-276 hardware. Hardware Torque to 50 in-oz.

Torque Screwdriver: Proto 6104 SN: 139072
Cal: 9/29/04 4/11/05 Due: 3/20/06 10/11/05
SBY 4/21/05
Model: Sartorius Genius SN: 12809099
Cal: 11/10/04 Due: 5/10/05

Initial Weight: 40.5974g
Final Weight: 40.5515g

Solution: 5M MgCl₂ (10M Cl)
2033.09g MgCl₂ · 6H₂O Lot: 042821
+ DI to 2L

Reagents measured with Model: OHAUS SN: 2883
Cal: 1/14/05 Due: 7/14/05

Initial pH: 2.34 Model: Orion EA 940 SN: 2330
Final pH: 4.09 Cal: 7/21/04 Due: 7/21/05
pH Probe: #13-620-296 SN: 4079126

Test Temperature: 60°C Measured with Hg Thermometer SN: E98-191
Cal: 01/07/05 Due: 07/07/05

Counter Electrode: Platinum Flag
Reference Electrode: Fisher 13-620-51 SN: 7282317

Gas: 99.999% Nitrogen Gas

Ecorr: -345 mV_{SCE} Model: Keithley 614 SN: 467374
Ept: +211 mV_{SCE} Cal: 12/2/04 Due: 12/2/2005
Potentiostat: Solartron 1480 SN: 00240551


DATA FILE: C-22 Trap 35b

Number of Crevice Corrosion Sites: 22/24 (24 max.)

Minor corrosion under crevice washers.
No staining.

To Page No. _____

Witnessed & Understood by me,	Date	Invented by	Date	
		Recorded by <i>Steve Young</i>	4/19/05	



Witnessed & Understood by me,	Date	Invented by	Date	
		Recorded by <i>Steve Young</i>	4/19/05	

From Page No. _____

POTENTIAL SCAN AND HOLD

Objective: Same as page 1.

Alloy / Heat No. : C-22 Heat# 2277-3-3266

Specimen Preparation: Specimen machined to dimensions specified in CNWRA Drawing. Specimen surfaces polished to 600 Grit finish using SiC paper. Specimen cleaned in acetone and rinsed in DI water. PTFE crevice forming washers attached to specimen using insulated C-276 hardware. Hardware Torque to 50 in-oz.

Torque Screwdriver:	Proto 6104	SN: 139072
	Cal: 6/29/04 4/11/05	Due: 3/29/05 10/11/05
	SBY 4/21/05	SBY 4/21/05
Initial Weight: 39.54123g	Model: Sartorius Genius	SN: 12809099
Final Weight: 39.46246g	Cal: 11/10/04	Due: 5/10/05

Solution: 3.5M MgCl₂ (7M Cl)
1423.24g MgCl₂·6H₂O Lot: 04282
+ DI to 2L

Reagents measured with	Model: OHAUS	SN: 2883
	Cal: 1/14/05	Due: 7/14/05
Initial pH: 4.20	Model: Orion EA 940	SN: 2330
Final pH: 3.16	Cal: 7/21/04	Due: 7/21/05
	pH Probe: #13-620-296	SN: 4079126

Test Temperature: 80° C Measured with Hg Thermometer SN: 323007
Cal: 10/14/04 Due: 10/14/2005

Counter Electrode: Platinum Flag

Reference Electrode: Fisher 13-620-52 SN: 4028036

Gas: 99.999% Nitrogen Gas

Ecorr: -150 mV _{SCE}	Model: Keithley 614	SN: 467374
Ept: +480 mV _{SCE}	Cal: 12/2/04	Due: 12/2/2005

Potentiostat: Solartron 1480 SN: 00240551


DATA FILE: C-22 Trap 36a

Number of Crevice Corrosion Sites: 11 / 24 (24 max.)

Mild staining.

To Page No. _____

Witnessed & Understood by me, _____ Date _____
Invented by _____ Date _____
Recorded by Steve Young 4/20/05



Witnessed & Understood by me, _____ Date _____
Invented by _____ Date _____
Recorded by Steve Young 4/20/05

From Page No. _____

POTENTIAL SCAN AND HOLD

Objective: Same as page 1.

Alloy / Heat No. : C-22 Heat# 2277-3-3266

Specimen Preparation: Specimen machined to dimensions specified in CNWRA Drawing. Specimen surfaces polished to 600 Grit finish using SiC paper. Specimen cleaned in acetone and rinsed in DI water. PTFE crevice forming washers attached to specimen using insulated C-276 hardware. Hardware Torque to 50 in-oz.

Torque Screwdriver:	Proto 6104	SN: 139072
	Cal: 9/29/04 4/11/05	Due: 3/29/05 10/11/05
	SBY 4/21/05	SBY 4/21/05
Initial Weight: 40.37935 g	Model: Sartorius Genius	SN: 12809099
Final Weight: 40.35763 g	Cal: 11/10/04	Due: 5/10/05

Solution: 3.5M MgCl₂ (7M Cl)
1423.14g MgCl₂ · 6H₂O Lot: 043664
+ DI to 2L

Reagents measured with	Model: OHAUS	SN: 2883
	Cal: 1/14/05	Due: 7/14/05
Initial pH: 4.69	Model: Orion EA 940	SN: 2330
Final pH: 5.44	Cal: 7/21/04	Due: 7/21/05
	pH Probe: #13-620-296	SN: 4079126

Test Temperature: 80° C	Measured with Hg Thermometer	SN: E98-191
	Cal: 01/07/05	Due: 07/07/05

Counter Electrode: Platinum Flag		
Reference Electrode: Fisher 13-620-51		SN: 7282317

Gas: 99.999% Nitrogen Gas		
Ecorr: -380 mV _{SCE}	Model: Keithley 614	SN: 467374
Ept: +510 mV _{SCE}	Cal: 12/2/04	Due: 12/2/2005
Potentiostat: Solartron 1480		SN: 00240551

DATA FILE: C-22 Trap 36b


Number of Crevice Corrosion Sites: 10/24 (24 max.)

Mild staining.

To Page No. _____

Witnessed & Understood by me,	Date	Invented by	Date
		Recorded by <i>Steve Young</i>	4/20/05

From Page No. _____



To Page No. _____

Witnessed & Understood by me,	Date	Invented by	Date
		Recorded by <i>Steve Young</i>	4/20/05

From Page No. _____

POTENTIAL SCAN AND HOLD

Objective: Same as page 1.

Alloy / Heat No. : C-22 Heat# 2277-3-3266

Specimen Preparation: Specimen machined to dimensions specified in CNWRA Drawing. Specimen surfaces polished to 600 Grit finish using SiC paper. Specimen cleaned in acetone and rinsed in DI water. PTFE crevice forming washers attached to specimen using insulated C-276 hardware. Hardware Torque to 50 in-oz.

Torque Screwdriver: Proto 6104 SN: 139072 Cal: 4/11/05 Due: 10/11/05

Initial Weight: 39.04885 g Model: Sartorius Genius SN: 12809099
Final Weight: 39.02677 g Cal: 11/10/04 Due: 5/10/05

Solution: 5M MgCl₂ (10M Cl)
2033.11g MgCl₂ · 6H₂O Lot: 043664
+ DI to 2L

Reagents measured with Model: OHAUS SN: 2883 Cal: 1/14/05 Due: 7/14/05

Initial pH: 2.23 Model: Orion EA 940 SN: 2330
Final pH: 4.96 Cal: 7/21/04 Due: 7/21/05
pH Probe: #13-620-296 SN: 4079126

Test Temperature: 80°C Measured with Hg Thermometer SN: 323007 Cal: 10/14/04 Due: 10/14/2005

Counter Electrode: Platinum Flag

Reference Electrode: Fisher 13-620-52 SN: 4028036

Gas: 99.999% Nitrogen Gas

Ecorr: -370 mV_{SCE} Model: Keithley 614 SN: 467374
Ept: +460 mV_{SCE} Cal: 12/2/04 Due: 12/2/2005

Potentiostat: Solartron 1480 SN: 00240551

DATA FILE: C-22 Trap 37a


Number of Crevice Corrosion Sites: 11 / 24 (24 max.)

Corrosion initiating around outside edges of crevice washers and moving outward. Grain boundaries clearly defined. Dark staining.

To Page No. _____

Witnessed & Understood by me,	Date	Invented by	Date
		Recorded by <i>Steve Young</i>	

From Page No. _____



To Page No. _____

Witnessed & Understood by me,	Date	Invented by	Date
		Recorded by <i>Steve Young</i>	

From Page No. _____

POTENTIAL SCAN AND HOLD

Objective: Same as page 1.

Alloy / Heat No. : C-22 Heat# 2277-3-3266

Specimen Preparation: Specimen machined to dimensions specified in CNWRA Drawing. Specimen surfaces polished to 600 Grit finish using SiC paper. Specimen cleaned in acetone and rinsed in DI water. PTFE crevice forming washers attached to specimen using insulated C-276 hardware. Hardware Torque to 50 in-oz.

Torque Screwdriver: Proto 6104 SN: 139072 Cal: 4/11/05 Due: 10/11/05

Initial Weight: 39.27164g Model: Sartorius Genius SN: 12809099
Final Weight: 39.24825g Cal: 11/10/04 Due: 5/10/05

Solution: 5M MgCl₂ (10M Cl)
2033.10g MgCl₂·6H₂O Lot: 043664
+ DI to 2L

Reagents measured with Model: OHAUS SN: 2883 Cal: 1/14/05 Due: 7/14/05

Initial pH: 2.41 Model: Orion EA 940 SN: 2330
Final pH: 4.80 Cal: 7/21/04 Due: 7/21/05
pH Probe: #13-620-296 SN: 4079126

Test Temperature: 80° C Measured with Hg Thermometer SN: E98-191 Cal: 01/07/05 Due: 07/07/05

Counter Electrode: Platinum Flag

Reference Electrode: Fisher 13-620-51 SN: 7282317

Gas: 99.999% Nitrogen Gas

Ecorr: -410 mV_{SCE} Model: Keithley 614 SN: 467374
Ept: +470 mV_{SCE} Cal: 12/2/04 Due: 12/2/2005

Potentiostat: Solartron 1480 SN: 00240551


DATA FILE: C-22 Trap 37b

Number of Crevice Corrosion Sites: 13 / 24 (24 max.)

Corrosion initiating around outside edges of crevice washers and moving outward.
Grain boundaries clearly defined.
Dark staining.

To Page No. _____

Witnessed & Understood by me,	Date	Invented by	Date
		Recorded by <i>Steve Young</i>	4/22/05



Witnessed & Understood by me,	Date	Invented by	Date
		Recorded by <i>Steve Young</i>	4/22/05

From Page No. _____

POTENTIAL SCAN AND HOLD

Objective: Same as page 1.

Alloy / Heat No. : C-22 Heat# 2277-3-3266

Specimen Preparation: Specimen machined to dimensions specified in CNWRA Drawing. Specimen surfaces polished to 600 Grit finish using SiC paper. Specimen cleaned in acetone and rinsed in DI water. PTFE crevice forming washers attached to specimen using insulated C-276 hardware. Hardware Torque to 50 in-oz.

Torque Screwdriver: Proto 6104 SN: 139072
Cal: 4/11/05 Due: 10/11/05

Initial Weight: 39.42491g Model: Sartorius Genius SN: 12809099
Final Weight: 39.40042g Cal: 11/10/04 Due: 5/10/05

Solution: 3.5M MgCl₂ (7M Cl)
1423.15g MgCl₂ · 6H₂O Lot: 043664
+ DI to 2L

Reagents measured with Model: OHAUS SN: 2883
Cal: 1/14/05 Due: 7/14/05

Initial pH: 3.54 Model: Orion EA 940 SN: 2330
Final pH: 4.50 Cal: 7/21/04 Due: 7/21/05
pH Probe: #13-620-296 SN: 4079126

Test Temperature: 95°C Measured with Hg Thermometer SN: 323007
Cal: 10/14/04 Due: 10/14/2005

Counter Electrode: Platinum Flag
Reference Electrode: Fisher 13-620-52 SN: 4028036

Gas: 99.999% Nitrogen Gas
Ecorr: -365 mV_{SCE} Model: Keithley 614 SN: 467374
Ept: +456 mV_{SCE} Cal: 12/2/04 Due: 12/2/2005
Potentiostat: Solartron 1480 SN: 00240551

DATA FILE: C-22 Trap 38a

Number of Crevice Corrosion Sites: 24 /24 (24 max.)

Mild surface staining on all surfaces.
Crevice corrosion develops from outer edge
of crevice washers.

To Page No. _____

Witnessed & Understood by me,

Date

Invented by


Date

Recorded by

Steve Young

4/26/05

From Page No. _____



To Page No. _____

Witnessed & Understood by me,

Date

Invented by

Date

Recorded by

Steve Young

4/26/05

From Page No. _____

POTENTIAL SCAN AND HOLD

Objective: Same as page 1.

Alloy / Heat No. : C-22 Heat# 2277-3-3266

Specimen Preparation: Specimen machined to dimensions specified in CNWRA Drawing. Specimen surfaces polished to 600 Grit finish using SiC paper. Specimen cleaned in acetone and rinsed in DI water. PTFE crevice forming washers attached to specimen using insulated C-276 hardware. Hardware Torque to 50 in-oz.

Torque Screwdriver: Proto 6104 SN: 139072 Cal: 4/11/05 Due: 10/11/05

Initial Weight: 39.06284 g Model: Sartorius Genius SN: 12809099
Final Weight: 39.05949 g Cal: 11/10/04 Due: 5/10/05

Solution: 3.5M MgCl₂
1423.17g MgCl₂ · 6H₂O Lot: 043664
+ DI to 2L

Reagents measured with Model: OHAUS SN: 2883 Cal: 1/14/05 Due: 7/14/05

Initial pH: 3.59 Model: Orion EA 940 SN: 2330
Final pH: 5.23 Cal: 7/21/04 Due: 7/21/05
pH Probe: #13-620-296 SN: 4079126

Test Temperature: 95°C Measured with Hg Thermometer SN: E98-191 Cal: 01/07/05 Due: 07/07/05

Counter Electrode: Platinum Flag

Reference Electrode: Fisher 13-620-51 SN: 7282317

Gas: 99.999% Nitrogen Gas

Ecorr: -405 mV_{SCE} Model: Keithley 614 SN: 467374
Ept: +448 mV_{SCE} Cal: 12/2/04 Due: 12/2/2005

Potentiostat: Solartron 1480 SN: 00240551

DATA FILE: C-22 Trap 38b


Number of Crevice Corrosion Sites: 6/24 (24 max.)

Mild surface staining on all sides.
Corrosion from crevice washer outer edges only.

To Page No. _____

Witnessed & Understood by me,	Date	Invented by	Date
		Recorded by <i>Steve Young</i>	4/26/05

From Page No. _____



To Page No. _____

Witnessed & Understood by me,	Date	Invented by	Date
		Recorded by <i>Steve Young</i>	4/26/05

From Page No. _____

POTENTIAL SCAN AND HOLD

Objective: Same as page 1.

Alloy / Heat No. : C-22 Heat# 2277-3-3266

Specimen Preparation: Specimen machined to dimensions specified in CNWRA Drawing. Specimen surfaces polished to 600 Grit finish using SiC paper. Specimen cleaned in acetone and rinsed in DI water. PTFE crevice forming washers attached to specimen using insulated C-276 hardware. Hardware Torque to 50 in-oz.

Torque Screwdriver: Proto 6104 SN: 139072 Cal: 4/11/05 Due: 10/11/05

Initial Weight: 39.63893g Model: Sartorius Genius SN: 12809099
Final Weight: 39.52530g Cal: 11/10/04 Due: 5/10/05

Solution: 5M MgCl₂ (10M Cl)
2033.08g MgCl₂ · 6H₂O Lot: 043664
+ DI to 2L

Reagents measured with Model: OHAUS SN: 2883 Cal: 1/14/05 Due: 7/14/05

Initial pH: 3.10 Model: Orion EA 940 SN: 2330
Final pH: 5.28 Cal: 7/21/04 Due: 7/21/05
pH Probe: #13-620-296 SN: 4079126

Test Temperature: 95° C Measured with Hg Thermometer SN: 323007 Cal: 10/14/04 Due: 10/14/2005

Counter Electrode: Platinum Flag

Reference Electrode: Fisher 13-620-52 SN: 4028036

Gas: 99.999% Nitrogen Gas

Ecorr: -365 mV_{SCE} Model: Keithley 614 SN: 467374
Ept: +447 mV_{SCE} Cal: 12/2/04 Due: 12/2/2005

Potentiostat: Solartron 1480 SN: 00240551

DATA FILE: C-22 Trap 39a

Number of Crevice Corrosion Sites: 24/24 (24 max.)

Dark black staining on all surfaces.
Corrosion develops from outer edges
of crevice washers.

To Page No. _____

Witnessed & Understood by me, .

Date

Invented by


Date

Recorded by

Steve Young

4/27/05

From Page No. _____



To Page No. _____

Witnessed & Understood by me,

Date

Invented by

Date

Recorded by

Steve Young

4/27/05

From Page No. _____

POTENTIAL SCAN AND HOLD

Objective: Same as page 1.

Alloy / Heat No. : C-22 Heat# 2277-3-3266

Specimen Preparation: Specimen machined to dimensions specified in CNWRA Drawing. Specimen surfaces polished to 600 Grit finish using SiC paper. Specimen cleaned in acetone and rinsed in DI water. PTFE crevice forming washers attached to specimen using insulated C-276 hardware. Hardware Torque to 50 in-oz.

Torque Screwdriver: Proto 6104 SN: 139072
Cal: 4/11/05 Due: 10/11/05

Initial Weight: 39.72659g Model: Sartorius Genius SN: 12809099
Final Weight: 39.64435g Cal: 11/10/04 Due: 5/10/05

Solution: 5M MgCl₂ (10M Cl)
2033.13g MgCl₂ · 6H₂O Lot: 043664
+ DI to 2L

Reagents measured with Model: OHAUS SN: 2883
Cal: 1/14/05 Due: 7/14/05

Initial pH: 2.57 Model: Orion EA 940 SN: 2330
Final pH: 5.31 Cal: 7/21/04 Due: 7/21/05
pH Probe: #13-620-296 SN: 4079126

Test Temperature: 95 °C Measured with Hg Thermometer SN: E98-191
Cal: 01/07/05 Due: 07/07/05

Counter Electrode: Platinum Flag

Reference Electrode: Fisher 13-620-51 SN: 7282317

Gas: 99.999% Nitrogen Gas

Ecorr: -385 mV_{SCE} Model: Keithley 614 SN: 467374
Ept: +454 mV_{SCE} Cal: 12/2/04 Due: 12/2/2005

Potentiostat: Solartron 1480 SN: 00240551

DATA FILE: C-22 Trap 39b

Number of Crevice Corrosion Sites: 24 / 24 (24 max.)

Dark black staining on all surfaces.
Corrosion develops from outer edges of crevice washer.

To Page No. _____

Witnessed & Understood by me,

Date

Invented by


Date

Recorded by

Steve Young

4/27/05

From Page No. _____



To Page No. _____

Witnessed & Understood by me,

Date

Invented by

Date

Recorded by

Steve Young

4/27/05

From Page No. _____

POTENTIAL SCAN AND HOLD

Objective: Same as page 1.

Alloy / Heat No. : C-22 Heat# 2277-3-3266

Specimen Preparation: Specimen machined to dimensions specified in CNWRA Drawing. Specimen surfaces polished to 600 Grit finish using SiC paper. Specimen cleaned in acetone and rinsed in DI water. PTFE crevice forming washers attached to specimen using insulated C-276 hardware. Hardware Torque to 50 in-oz.

Torque Screwdriver: Proto 6104 SN: 139072 Cal: 4/11/05 Due: 10/11/05

Initial Weight: 39.73573g Model: Sartorius Genius SN: 12809099
Final Weight: 39.69930g Cal: 11/10/04 Due: 5/10/05

Solution: 3.5M MgCl₂ (7M Cl)
1423.15g MgCl₂ · 6H₂O Lot: 043664
+ DI to 2L

Reagents measured with Model: OHAUS SN: 2883 Cal: 1/14/05 Due: 7/14/05

Initial pH: 3.83 Model: Orion EA 940 SN: 2330
Final pH: 5.61 Cal: 7/21/04 Due: 7/21/05
pH Probe: #13-620-296 SN: 4079126

Test Temperature: 105 °C Measured with Hg Thermometer SN: 323007 Cal: 10/14/04 Due: 10/14/2005

Counter Electrode: Platinum Flag

Reference Electrode: Fisher 13-620-52 SN: 4028036

Gas: 99.999% Nitrogen Gas

Ecorr: -375 mV_{SCE} Model: Keithley 614 SN: 467374
Ept: +457 mV_{SCE} Cal: 12/2/04 Due: 12/2/2005

Potentiostat: Solartron 1480 SN: 00240551

DATA FILE: C-22 Trap 40a


Number of Crevice Corrosion Sites: 24/24 (24 max.)

Dark blue staining on all surfaces.
Corrosion develops from outer edge of crevice washers.

To Page No. _____

Witnessed & Understood by me,	Date	Invented by	Date
		Recorded by <i>Steve Young</i>	5/3/05

From Page No. _____



To Page No. _____

Witnessed & Understood by me,	Date	Invented by	Date
		Recorded by <i>Steve Young</i>	5/3/05

From Page No. _____

POTENTIAL SCAN AND HOLD

Objective: Same as page 1.

Alloy / Heat No. : C-22 Heat# 2277-3-3266

Specimen Preparation: Specimen machined to dimensions specified in CNWRA Drawing. Specimen surfaces polished to 600 Grit finish using SiC paper. Specimen cleaned in acetone and rinsed in DI water. PTFE crevice forming washers attached to specimen using insulated C-276 hardware. Hardware Torque to 50 in-oz.

Torque Screwdriver: Proto 6104 SN: 139072 Cal: 4/11/05 Due: 10/11/05

Initial Weight: 39.89529g Model: Sartorius Genius SN: 12809099
Final Weight: 39.85040g Cal: 11/10/04 Due: 5/10/05

Solution: 3.5M MgCl₂ (7M Cl)
1423.16g MgCl₂ · 6H₂O Lot: 043664
+ DI to 2L

Reagents measured with Model: OHAUS SN: 2883 Cal: 1/14/05 Due: 7/14/05

Initial pH: 3.45 Model: Orion EA 940 SN: 2330
Final pH: 5.35 Cal: 7/21/04 Due: 7/21/05
pH Probe: #13-620-296 SN: 4079126

Test Temperature: 105 °C Measured with Hg Thermometer SN: E98-191 Cal: 01/07/05 Due: 07/07/05

Counter Electrode: Platinum Flag

Reference Electrode: Fisher 13-620-51 SN: 7282317

Gas: 99.999% Nitrogen Gas

Ecorr: -381 mV_{SCE} Model: Keithley 614 SN: 467374
Ept: +457 mV_{SCE} Cal: 12/2/04 Due: 12/2/2005

Potentiostat: Solartron 1480 SN: 00240551

DATA FILE: C-22 Trap 40b

Number of Crevice Corrosion Sites: 24/24 (24 max.)

Dark blue tint staining on all surfaces.
Corrosion develops from outer edges of crevice washers.

To Page No. _____

Witnessed & Understood by me,

Date

Invented by


Date

Recorded by

Steve Young

5/3/05

From Page No. _____



To Page No. _____

Witnessed & Understood by me,

Date

Invented by

Date

Recorded by

Steve Young

5/3/05



From Page No. _____

Continued testing

In notebook

712

I have reviewed this scientific notebook and find it in compliance with QAP-001. There is sufficient information regarding procedures used for conducting tests, acquiring and analyzing data so that another qualified individual could repeat the activity.

[Signature] 6/6/05

To Page No. _____

Witnessed & Understood by me,

Date

Invented by

Date

Recorded by

[Signature]

5/3/05

ADDITIONAL INFORMATION FOR SCIENTIFIC NOTEBOOK NO. 697

Document Date:	03/28/2005
Availability:	Southwest Research Institute® Center for Nuclear Waste Regulatory Analyses 6220 Culebra Road San Antonio, Texas 78228
Contact:	Southwest Research Institute® Center for Nuclear Waste Regulatory Analyses 6220 Culebra Road San Antonio, TX 78228-5166 Attn.: Director of Administration 210.522.5054
Data Sensitivity:	<input checked="" type="checkbox"/> "Non-Sensitive" <input type="checkbox"/> Sensitive <input type="checkbox"/> "Non-Sensitive - Copyright" <input type="checkbox"/> Sensitive - Copyright
Date Generated:	06/04/2005
Operating System: (including version number)	Windows NT
Application Used: (including version number)	Excel and others unknown
Media Type: (CDs, 3 1/2, 5 1/4 disks, etc.)	1 - CD
File Types: (.exe, .bat, .zip, etc.)	xls, dat
Remarks: (computer runs, etc.)	Media contains two folders: galvanodynamic and potentiodynamic containing data files.