Appendix D

ESEM and SEM/EDS Data for Test #3 Day-30 Corrosion Products

List of Figures

Figure D-1: SEM image magnified 200 times for a Test #3 Day-30 powder on the
submerged rack. (T3~RackPowder013)
Figure D-2: EDS counting spectrum for the powder on the submerged rack shown in Figure
D-1. (T3RackPowder09)
Figure D-3: SEM image magnified 500 times for a Test #3 Day-30 powder on a submerged
galvanized steel coupon. (T3~Gal~Steel~Powder014)
Figure D-4: EDS counting spectrum for the white powder on galvanized steel shown in
Figure D-3. (CorPdct~GalStel10)
Figure D-5: EDS counting spectrum for the dark powder on galvanized steel shown in
Figure D-3. (CorPrdct~GalSteel11)
Figure D-6: SEM image magnified 500 times for a Test #3 Day-30 powder on a submerged
copper coupon. (T3~Copper~Powder015)
Figure D-7: EDS counting spectrum for the powder on copper shown in Figure D-6.
(CorPrdct~cu12)
Figure D-8: SEM image magnified 100 times for Test #3 Day-30 corrosion products on a
Figure D-8: SEM image magnified 100 times for Test #3 Day-30 corrosion products on a submerged aluminum coupon. (T3D30CorrPrdctSubmAl039)
submerged aluminum coupon. (T3D30CorrPrdctSubmAl039)

Figure D-14	4: EDS counting spectrum for the corrosion products (particles) shown in Figure	
	D-13. (t3cont30)	D-16
	List of Tables	
Table D-1.	Chemical Compositions for T3RackPowder09, Figure D-2.	. D-6
Table D-2.	Chemical Compositions for CorPrdct~GalSteel11, Figure D-5	. D-9
Table D-3.	Chemical Compositions for CorPrdct~cu12, Figure D-7	D-11
Table D-4.	Chemical Compositions for T3D30CorrPrdtAl22, Figure D-10	D-14


For ICET tests, one process of interest is the corrosion effect of metal and concrete coupons. One means of understanding the corrosion process is through direct examination of the corrosion products after the test is completed. For this purpose, corrosion products were collected when Test #3 was shut down (May 5, 2005). These corrosion products included (1) fine powders on a vertical piece of the submerged CPVC rack; (2) corrosion products on a submerged galvanized steel coupon; (3) corrosion products on a submerged copper coupon; (4) corrosion products on a submerged aluminum coupon; and (5) corrosion products on a submerged concrete coupon.

Corrosion products were collected by directly adhering the sample onto double-sided carbon tape suitable for SEM/EDS examination. After the samples were dried in air, a Au/Pd coating was applied to enhance the surface conductivity of the samples and to prevent possible charging problems during the SEM examination process. Based on EDS results, a semi-quantitative elemental analysis was performed after calibration of the x-ray signal using an internal standard of the microscopy. This appendix presents the SEM/EDS data that were generated on May 17, 2005. Available logbook entries for this laboratory session are included in this appendix as transcribed notes.

Transcribed Laboratory Log

Laboratory session from May 9, 2005.

Test #3 Day-30 Powder Deposits



Powder on Rack

Powder on Galvanized Steel

Powder on Copper


Image: T3~Copper~Powder015 $500 \times$ Figure D-6

EDS: CorPrdct_cu12 Powder on image 015 Figure D-7

Transcribed Laboratory Log

Laboratory session from May 17, 2005.

Test #3 Day-30 Metal Coupons



**Coat with Gold


Corrosion Product on Submerged Aluminum

Image:	T3D30CorrPrdctSubmAl039	100 ×		Figure D-8
	T3D30CorrPrdctSubmAl040	100 ×	Backscattering	Figure D-9
EDS:	T3D30CorrPrdtAl22		Particles on image 040	Figure D-10
Image:	T3D30CorrPrdctSubmAl041	1000 ×		Figure D-11

Transcribed Laboratory Log

Laboratory session from May 10, 2005.

Test #3 Day-30 ESEM.



Deposits on the Top of Concrete

Image:	T3Cont28	1000 ×		Figure D-12
	t3cont29	100 ×		Figure D-13
EDS:	t3cont30		Particles on image 29	Figure D-14




Figure D-1: SEM image magnified 200 times for a Test #3 Day-30 powder on the submerged rack. (T3~RackPowder013)




Figure D-2: EDS counting spectrum for the powder on the submerged rack shown in Figure D-1. (T3RackPowder09)

The results from the chemical composition analysis for T3RackPowder09 are given in Table D-1.

Table D-1. Chemical Compositions for T3RackPowder09, Figure D-2.

May 9 2005

```
Group
          : NRC
Sample
          : T3D30
                   ID# : 9
Comment
          : Powder on submerged Rack
Condition: Full Scale: 20KeV(10eV/ch, 2Kch)
            Live Time
                           60.000 sec
                                          Aperture #
            Acc. Volt
                        : 15.0 KV
                                          Probe Current: 1.608E-09 A
            Stage Point : X=23.811 Y=58.398 Z= 9.938
                        : Mon May 9 15:02:46 2005
            Acq. Date
Element
            Mode
                     ROI (KeV)
                                K-ratio(%)
                                            +/-
                                                    Net/Background
  O K
           Normal
                    0.25- 0.77
                                 16.4417
                                            0.0037
                                                       6188 /
                                                                    66
 Na K
                                                        448 /
           Normal
                    0.81- 1.27
                                  0.4147
                                            0.0011
                                                                   64
                                                       2708 /
 Si K
           Normal
                    1.50- 2.05
                                  1.8593
                                            0.0006
                                                                  288
  PK
           Normal
                    1.75- 2.38
                                  10.0274
                                            0.0058
                                                       9012 /
                                                                  160
                                                      15798 /
 Ca K
           Normal
                    3.39- 4.30
                                  22.3191
                                            0.0043
                                                                   22
  CK
           Normal
                    0.09- 0.46
                                                         19 /
                                  0.0345
                                            0.0005
                                                                  166
                             Chi_square = 55.4587
Element Mass%
                 Atomic%
                           ZAF
     0
         48.147
                 67.3988 1.8951 0.9734 1.9469 1.0000
    Na
          0.977
                  0.9516 1.5243 1.0269 1.4843 1.0000
    Si
          3.302
                  2.6332 1.1493 0.9767 1.1837 0.9941
     P
         13.197
                  9.5423 0.8517 1.1722 0.7283 0.9977
    Ca
         34.174
                 19.0962 0.9909 0.9956 0.9953 1.0000
          0.203
                  0.3778 3.7963 1.0208 3.7192 0.9999
Total
        100.000 100.0000
Normalization factor = 1.5452
        100.000 100.0000
Total
Normalization factor = 2.1120
```




Figure D-3: SEM image magnified 500 times for a Test #3 Day-30 powder on a submerged galvanized steel coupon. (T3~Gal~Steel~Powder014)




Figure D-4: EDS counting spectrum for the white powder on galvanized steel shown in Figure D-3. (CorPdct~GalStel10)




Figure D-5: EDS counting spectrum for the dark powder on galvanized steel shown in Figure D-3. (CorPrdct \sim GalSteel11)

The results from the chemical composition analysis for CorPrdct~GalSteel11 are given in Table D-2.

Table D-2. Chemical Compositions for CorPrdct~GalSteel11, Figure D-5.

May 9 2005

: NRC Group

Sample : T3D30 ID# : 11

Comment : Corrosion product of Gal-Steel Condition: Full Scale: 20KeV(10eV/ch,2Kch)
Live Time: 60.000 sec Aperture #: 1
Acc. Volt: 15.0 KV Probe Current: 1.

Probe Current: 1.608E-09 A

Stage Point : X=12.508 Y=58.337 Z= 9.938 Acq. Date : Mon May 9 15:23:47 2005

Elen	nent	Mode	ROI (KeV)	K-ratio(%)	+/-	Net/Back	gro	und
0	K	Normal	0.25- 0.77	24.6017	0.0044	9259	7	71
Mg	K	Normal	0.97- 1.57	0.1342	0.0003	211	1	424
Si	K	Normal	1.50- 2.05	4.3383	0.0008	6318	1	94
Ca	K	Normal	3.39- 4.30	1.4823	0.0016	1049	1	28
Zn	K	Normal	8.22-10.03	13.8053	0.0096	1501	1	6
Al	K	Normal	1.19- 1.83	0.1257	0.0003	196	1	77
C	K	Normal	0.09- 0.46	0.0000	0.0000	0	1	254

Chi square = 39.1864

Element	Mass*	Atomic%	ZAF	Z	A	F	
0	47.383	72.4077	0.9119	0.9354	0.9750	1.0000	
Mg	0.655	0.6591	2.3126	0.9284	2.4928	0.9993	
Si	13.143	11.4404	1.4344	0.9377	1.5298	0.9999	
Ca	3.059	1.8657	0.9770	0.9542	1.0247	0.9991	
Zn	35.285	13.1966	1.2102	1.2129	0.9978	1.0000	
Al	0.475	0.4304	1.7897	0.9594	1.8692	0.9979	
C	0.000	0.0000	5.2279	0.9811	5.3287	1.0000	

Total 100.000 100.0000

Normalization factor = 2.1120




Figure D-6: SEM image magnified 500 times for a Test #3 Day-30 powder on a submerged copper coupon. (T3~Copper~Powder015)




Figure D-7: EDS counting spectrum for the powder on copper shown in Figure D-6. (CorPrdct~cu12)

The results from the chemical composition analysis for CorPrdct~cu12 are given in Table D-3.

Table D-3. Chemical Compositions for CorPrdct~cu12, Figure D-7.

May 9 2005

```
Group
          : NRC
          : T3D30 ID# : 112
Sample
          : Corrosion product of Copper
Comment
Condition: Full Scale: 20KeV(10eV/ch,2Kch)
            Live Time
                       : 60.000 sec
                                         Aperture #
                                                        : 1
                        : 15.0 KV
                                         Probe Current: 1.609E-09 A
            Acc. Volt
            Stage Point : X=17.909 Y=70.722 Z= 9.938
                      : Mon May 9 15:34:51 2005
            Acq. Date
                                            +/-
Element
           Mode
                     ROI (KeV)
                                K-ratio(%)
                                                   Net/Background
                                                                   54
 OK
           Normal
                    0.25- 0.77
                                 16.9583
                                           0.0036
                                                       6386 /
Na K
           Normal
                    0.81- 1.27
                                  0.3712
                                           0.0011
                                                        402 /
                                                                   46
 Si K
           Normal
                    1.50- 2.05
                                  1.7009
                                           0.0006
                                                       2479 /
                                                                  254
                                                       7744 /
                                                                  164
 PK
           Normal
                    1.75- 2.38
                                  8.6109
                                           0.0054
                                           0.0040
                                 19.7450
                                                      13984 /
                                                                   28
 Ca K
           Normal
                    3.39- 4.30
                                                         0 /
                                                                  156
                    0.09- 0.46
                                 0.0000
                                           0.0000
 CK
           Normal
                             Chi_square = 59.8005
Element Mass%
                           ZAF
                 Atomic*
                 70.0693 1.8119 0.9751 1.8582 1.0000
    0
         51.054
          0.947
                  0.9045 1.5353 1.0287 1.4923 1.0000
   Na
                  2.5484 1.1533 0.9785 1.1852 0.9945
    Si
          3.260
                 8.6603 0.8538 1.1744 0.7286 0.9978
     P
         12.216
                 17.8176 0.9913 0.9976 0.9936 1.0000
    Ca
         32.523
          0.000
                 0.0000 3.6900 1.0225 3.6090 0.9999
        100.000 100.0000
```

Total 100.000 100.0000 Normalization factor = 1.6616



Figure D-8: SEM image magnified 100 times for Test #3 Day-30 corrosion products on a submerged aluminum coupon. (T3D30CorrPrdctSubmAl039)




Figure D-9: Backscattered SEM image magnified 100 times for Test #3 Day-30 corrosion products on a submerged aluminum coupon. (T3D30CorrPrdctSubmAl040)




Figure D-10: EDS counting spectrum for the corrosion products (particles) shown in Figure D-9. (T3D30CorrPrdtAl22)

The results from the chemical composition analysis for T3D30CorrPrdtAl22 are given in Table D-4.

Table D-4. Chemical Compositions for T3D30CorrPrdtAl22, Figure D-10.

May 17 2005

: NRC Group Sample : T3D30 ID# : 25 Comment : Corr Prodct on Submerged Al Condition: Full Scale: 20KeV(10eV/ch,2Kch) Live Time : 60.000 sec Aperture # : 2 Acc. Volt : 15.0 KV Probe Current: 7.686E-09 A Stage Point: X=17.100 Y=58.819 Z=11.000 Acq. Date : Tue May 17 16:07:30 2005 Element Mode ROI (KeV) K-ratio(%) +/-Net/Background Normal O K 0.25- 0.77 43.5604 0.0044 10399 / 29 Normal 0.81- 1.27 497 / Na K 0.7246 0.0010 45 15016 / Al K Normal 1.26- 1.78 15.1592 0.0021 184 1.50- 2.05 6.6963 0.0008 3.40- 4.30 4.2145 0.0082 0.00- 0.36 0.6943 0.0002 Si K Normal 6185 / 674 Ca K Normal 1892 / 14 вк Normal 10 69 / -----Chi square = 107.5732 Element Mass% Atomic% ZAF

Si 10.480 6.6618 1.3149 0.9816 1.3397 0.9999 Ca 4.923 2.1928 0.9814 0.9916 0.9896 1.0001 B 11.207 18.5068 13.5610 1.1342 11.9568 1.0000

52.728 58.8383 1.0170 0.9767 1.0413 1.0000

1.120 0.8699 1.2989 1.0311 1.2612 0.9988 19.542 12.9304 1.0831 0.9932 1.0934 0.9974

Total 100.000 100.0000

0 Na

Al

Normalization factor = 1.1902

Fe 2.399 1.5375 0.9896 1.0503 1.0060 0.9366

Total 100.000 100.0000

Normalization factor = 3.2486




Figure D-11: SEM image magnified 1000 times for Test #3 Day-30 corrosion products on a submerged aluminum coupon. (T3D30CorrPrdctSubmAl041)




Figure D-12: ESEM image magnified 1000 times for Test #3 Day-30 corrosion products on a submerged concrete coupon. (T3cont28)




Figure D-13: ESEM image magnified 100 times for Test #3 Day-30 corrosion products on a submerged concrete coupon. (t3cont29)




Figure D-14: EDS counting spectrum for the corrosion products (particles) shown in Figure D-13. (t3cont30)