



Q5

06-31343-01
03037855

MS-16

FAX

To: BETSY

From: DENISE

Fax #: 620-337-5269

Fax #: 203-288-3283

Date: 12-18-08

of Pages: 2

MESSAGE

Betsy,

Please let me know if you need any additional information. Happy Holidays.

Denise



142956

Device and Source Information

Sealed Sources for Quality Assurance and Calibration

ADAC Laboratories products require regular quality assurance (QA) and calibration to ensure their continued high quality performance. As noted in the previous subsection, some products require a specific sealed source for QA and calibration (e.g., Precedence requires the Gd-153 sources described on page 15; MCD and MCD-AC require the Ge-68 source described on page 16). For other products, depending what you have purchased, you may need to amend your RAM license to include one or more of the following sealed sources.

Co-57

Co-57 sealed sheet sources are used with nuclear medicine cameras as quality control and calibration sources to determine the intrinsic or extrinsic resolution of the camera detector.

The following approved supplier manufactures Co-57 sealed sheet sources according to ADAC's specifications.

SEALED SOURCE		
Manufacturer and Distributor	Isotope Products Laboratories (IPL) An Eckert & Ziegler Company 24937 Avenue Tibbitts Valencia, CA 91355	661-309-1010 661-257-8303 (fax)
Quantity	1	
Type	Co-57 flood source, rectangular, 24" x 36", encapsulated in plastic	
Model Number	FL Series Flood Sources, FL24R-057-10M	
Registration Number	CA0406S180S	
Radionuclide	Cobalt-57	
Single Source Strength	not to exceed 10 mCi	