

Nguyen, Janice

From: [REDACTED]
Sent: Friday, August 26, 2011 5:42 AM
To: Nguyen, Janice
Subject: Fwd: NRC
Attachments: LMP_Model14C_Survey_Meter_with_44-9_Pancake_Probe.pdf; Ludlum_14-C_Detection_Limits.pdf

Here is the survey meter efficiency information. I hope you have a great weekend!

I will be working in one of the doctors office today so if you have any questions, please feel free to use my cell phone
[REDACTED]

Regards,
Curtis

-----Original Message-----

From: Norweck, Jim <jnorweck@radiology-inc.com>
To: [REDACTED]
Sent: Fri, Aug 26, 2011 1:57 am
Subject: NRC and ACR

James Norweck, MS, DABR

Radiology Inc.
5221 US 60 East
Huntington, WV 25705
Office: 304-522-1550 x234
[REDACTED]

FAX: 304-522-0704
jnorweck@radiology-inc.com


This email communication and any attachments may contain confidential and privileged, or otherwise protected information for the use of the designated recipients named above only. If you are not the intended recipient, you are hereby notified that you have received this information in error and that any review, disclosure, dissemination, distribution or copying of it or its contents is prohibited and unauthorized. If you have received this communication in error, please notify me immediately by replying to this message and deleting it from your computer as well as destroying any copies that may have been created. Thank you

09/13/11
Rec'd in DRUMS

RADIATION SAFETY


Survey Meters (GM / Scintillation)

Model 14C



Part Number: 48-1611

Model 44-9



Part Number: 47-1539

Model 14C Survey Ratemeter with Pancake Probe

Introduction

This general purpose, handheld analog ratemeter supports operating two separate radiation detectors. A switch on the front panel allows the user to select between the internally mounted GM detector for detecting gamma exposure over a range of 0 - 2000 mR/hr or the external Model 44-9 GM pancake detector. The pancake detector is sensitive to alpha, beta and gamma and is the industry standard for detecting contamination. This survey meter additionally supports externally connected scintillation detectors in lieu of GM's. The Model 14C can be used in a wide range of Medical and Health Physics applications. With the addition of the Model 180-2 sample holder, it can be utilized to make a quick evaluation of wipe test surveys in the Nuclear Medicine department.

Specifications

Model 14C, Survey Meter

MULTIPLIERS: X0.1, X1, X10, X100, X1000

LINEARITY: Reading within $\pm 10\%$ of true value with detector connected

ENERGY RESPONSE: Within $\pm 15\%$ of true value between 60 keV - 3 MeV (internal detector only)

CONNECTOR: Series "C" (others available)

AUDIO: Built in unimorph speaker with ON/OFF switch (greater than 60 dB at 2 feet)

HIGH VOLTAGE: 900 V (setting can be checked on meter); THRESHOLD: 30 mV \pm 10 mV

RESPONSE: Toggle switch for FAST (4 seconds) or SLOW (22 seconds) from 10% to 90% of final reading

POWER: 2 each "D" cell batteries (housed in sealed compartment that is externally accessible)

BATTERY LIFE: Typically greater than 2000 hours with alkaline batteries

TEMPERATURE RANGE: -4°F to 122°F (-20° to 50°C)

SIZE: 6.5" H x 3.5" W x 8.5" L (16.5 x 8.9 x 21.6 cm)

WEIGHT: 3.5 lbs (1.6 kg) including batteries

Model 44-9, GM Pancake Detector

WINDOW: 1.7 ± 0.3 mg/cm² mica

WINDOW AREA: Active - 15 cm²; Open - 12 cm²

EFFICIENCY(2pi): Typically 5% ¹³⁷Cs; 22% ⁹⁰Sr/⁹⁰Y; 19% ¹³²Ic; 32% ³²P; 15% ²¹⁰Pb

SENSITIVITY: Typically 3300 cpm/mR hr (¹³⁷Cs gamma)

ENERGY RESPONSE: Energy dependent


DEAD TIME: Typically 80 μ s

Options:

Model 44-2, 1" x 1" NaI Gamma Detector: PN: 47-1532

1 μ Ci ¹³⁷Cs Check Source: PN: 01-5196

Check Source Holder: PN: 4062-166



Model 180-2
Sample Holder

RadCalc Ludlum Edition

File Edit Help

DLs - Surface Contamination

Ludlum Instrument | Ratemeter | Pre-Set Time

Background Count Rate (Rb) = 80 cpm

Inst. Model 14C Response Time: Fast Slow

Detector 44-9 Detector Area: Active Open Edit

Detector Area (A) = 15 cm² Detector Efficiency (Ed) = 5 %

Sample Count (Rs) = 0 cpm

Use MARSSIM method

Results

Lc (Critical Level) = 36.04 cpm above bkgd.

Ld (Detection Limit) = 81.08 cpm above bkgd.

MDA = 1622 dpm/detector

MDA = 1.081E4 dpm/100 cm²

MDA = 0.0007305 µCi/detector

MDA = 4.87E-5 µCi/100 cm²

All values calculated to 95% CL via MARSSIM methods:

Sample Activity = 0 dpm/detector

Sample Activity = 0 dpm/100 cm²

DETECTION LIMITS--SURFACE CONTAMINATION

INPUT DATA:

Background Count = 80 cpm
 Detector used: Ludlum Model 44-9
 Instrument used: Ludlum Model 14C on slow response
 Detector Area = 15 cm²
 Detector Efficiency = 5 % (worst case assumed)

RESULTS:

Critical Level (Lc) = 36.04 cpm above bkgd.
 Detection Limit (Ld) = 81.08 cpm above bkgd.
 Minimum Detectable Activity (MDA) = 1622 dpm/detector
 Minimum Detectable Activity (MDA) = 10,810 dpm/100 cm²
 Minimum Detectable Activity (MDA) = 0.0007305 µCi/detector
 Minimum Detectable Activity (MDA) = 4.87E-5 µCi/100 cm²

All values calculated to 95% CL via MARSSIM methods

Calculated by RadCalc Ludlum Edition version 1.0