

QTRLY HDR GYN Cylinder Size Check

Please place cylinders with 4 different sizes into the circles of proper size and verify that the diameter matches. Next, take x-ray of each cylinder with dummy source all the way in and measure the cylinder tip to center of dummy source distance. This distance should be within 1 mm of the nominal value.

QTR 1Date:

Width: L3.5cm _____ M 3.0cm _____ S 2.5cm _____ XS 2.0cm _____
Length(tip to dummy): 3.5 cm _____ 3.0cm _____ 2.5 cm _____ 2.0cm _____
Signature _____ Date _____

QTR 2Date:

Width: L3.5cm _____ M 3.0cm _____ S 2.5cm _____ XS 2.0cm _____
Length(tip to dummy): 3.5 cm _____ 3.0cm _____ 2.5 cm _____ 2.0cm _____
Signature _____ Date _____

QTR 3Date:

Width: L3.5cm _____ M 3.0cm _____ S 2.5cm _____ XS 2.0cm _____
Length(tip to dummy): 3.5 cm _____ 3.0cm _____ 2.5 cm _____ 2.0cm _____
Signature _____ Date _____

QTR 4Date:

Width: L3.5cm _____ M 3.0cm _____ S 2.5cm _____ XS 2.0cm _____
Length(tip to dummy): 3.5 cm _____ 3.0cm _____ 2.5 cm _____ 2.0cm _____
Signature _____ Date _____

