NUREG-1189 Vol. 2

Assessment of the Public Health Impact From the Accidental Release of UF6 at the Sequoyah Fuels Corporation Facility at Gore, Oklahoma

Docket No. 40-8027 License No. SUB-1010

Appendices

U.S. Nuclear Regulatory Commission

Ad Hoc Interagency Public Health Assessment Task Force



8603270335 860331 PDR ADDCK 04008027 B

NOTICE

Availability of Reference Materials Cited in NRC Publications

Most documents cited in NRC publications will be available from one of the following sources:

- The NRC Public Document Room, 1717 H Street, N.W. Washington, DC 20555
- The Superintendent of Documents, U.S. Government Printing Office, Post Office Box 37082, Washington, DC 20013-7082
- 3. The National Technical Information Service, Springfield, VA 22161

Although the listing that follows represents the majority of documents cited in NRC publications, it is not intended to be exhaustive.

Referenced documents available for inspection and copying for a fee from the NRC Public Document Room include NRC correspondence and internal NRC memoranda; NRC Office of Inspection and Enforcement bulletins, circulars, information notices, inspection and investigation notices; Licensee Event Reports; vendor reports and correspondence; Commission papers; and applicant and licensee documents and correspondence.

The following documents in the NUREG series are available for purchase from the GPO Sales Program: formal NRC staff and contractor reports, NRC-sponsored conference proceedings, and NRC booklets and brochures. Also available are Regulatory Guides, NRC regulations in the *Code of Federal Regulations*, and *Nuclear Regulatory Commission Issuances*.

Documents available from the National Technical Information Service include NUREG series reports and technical reports prepared by other federal agencies and reports prepared by the Atomic Energy Commission, forerunner agency to the Nuclear Regulatory Commission.

Documents available from public and special technical libraries include all open literature items, such as books, journal and periodical articles, and transactions. *Federal Register* notices, federal and state legislation, and congressional reports can usually be obtained from these libraries.

Documents such as theses, dissertations, foreign reports and translations, and non-NRC conference proceedings are available for purchase from the organization sponsoring the publication cited.

Single copies of NRC draft reports are available free, to the extent of supply, upon written request to the Division of Technical Information and Document Control, U.S. Nuclear Regulatory Commission, Washington, DC 20555.

Copies of industry codes and standards used in a substantive manner in the NRC regulatory process are maintained at the NRC Library, 7920 Norfolk Avenue, Bethesda, Maryland, and are available there for reference use by the public. Codes and standards are usually copyrighted and may be purchased from the originating organization or, if they are American National Standards, from the American National Standards Institute, 1430 Broadway, New York, NY 10018.

NUREG-1189 Vol. 2

Assessment of the Public Health Impact From the Accidental Release of UF₆ at the Sequoyah Fuels Corporation Facility at Gore, Oklahoma

Docket No. 40-8027 License No. SUB-1010

Appendices

Manuscript Completed: March 1986 Date Published: March 1986

Ad Hoc Interagency Public Health Assessment Task Force

U.S. Nuclear Regulatory Commission Washington, D.C. 20555



ABSTRACT

2

ļ

and

Following the accidental release of UF₆ from the Sequoyah Fuels Facility on January 4, 1986, an Ad Hoc Interagency Public Health Assessment Task Force was established. The Task Force consists of technical staff members from various agencies who have prepared this assessment of the public health impact associated with the accidental release.

The assessment consists of two volumes and is based on data from the accident available as of February 14, 1986.

Volume 1 of the report describes the effects from the intake of uranium and fluoride and summarizes the findings and recommendations of the Task Force.

Volume 2 of the report contains Appendices which provide more detailed information used in the assessment and support the discussion in Volume 1.

TABLE OF CONTENTS

Appendix		Page
2.2	Description of SFC's Facility Operation	1
3.1.1	Photos Showing the Location of the 14-Ton Cyclinder and the Rupture	11
3.3.1	NRC's Preliminary Notifications of Event or Unusual Occurrence Following the Accident Occurring on January 4, 1986	19
4.1.2	Meteorology and Atmospheric Dispersion	37
4.1.4.2	SFC's Investigation Report of the Effects on Cattle and Animals in the Potentially Affected Areas After the Accident	77
5.1.1	Air Sample Analysis by SFC and ORNL	131
5.2.1	Interorganizational Comparison Results	143
5.2.2	K-M Technical Center Laboratory Analysis Procedures	147
5.2.3	ORNL's Analytical Procedure	221
5.2.4	OSDH Sampling Procedure	247
5.2.5	SFC Sampling Procedure (2 Documents)	251
5.2.6	OSDH Laboratory Analysis Procedures	259
5.2.7	Environmental Sampling: Vegetation and Soil	265
5.3.1	Raw Surface Water Data	373
5.4.1	An Aerial Radiological Survey of the Kerr-McGee Facility and Surrounding Area	383
5.5.1	Results of HPGE In Situ Measurements	397
5.6.1	Medical Examination and Diagnosis for Onsite Workers	403
5.6.2	Uranium Bioassay Results for SFC Workers and Offsite Individuals	407
5.6.3	Intake Calculation Models and Assumptions	423
5.6.4	Chemical Bioassay Results for SFC Workers	431
5.6.5	Creatinine Correction Methodology and Factors	435
6.1.2.1	Summary of Dose Conversion Factors and Intake Parameters Used for Radiological Dose Calculation	443

ACKNOWLEDGEMENTS

The Ad Hoc Task Force acknowledges the assistance of the following individuals in providing data and analysis used in this assessment:

U. S. Nuclear Regulatory Commission	-	C.	Everett Jierree Spitzberg
Department of Energy & EG&G	-		Blauenstein Burson
Oak Ridge National Laboratory		R.	Eckerman Leggett Stewart
State of Oklahoma Department of Health	-	D.	Craig Morris Walker
Battelle Pacific Northwest Laboratory	-	D.	Fisher
Brookhaven National Laboratory	-	E.	Lessard

l

We also acknowledge the contributions of the following individuals:

U.S.	Nuclear Regulatory Commission	-		Congel Sturz
U.S.	Environmental Protection Agency	-	С.	Farland Gazda Makris

e

APPENDIX 2.2

DESCRIPTION OF SFC'S FACILITY OPERATION

APPENDIX 2.2 Description of SFC's Facility Operation

2.2.1 Receiving and sampling

Uranium feed material is delivered to the SFC facility in three forms: dry yellowcake (ore concentrate), wet yellowcake slurry, and wet UF₄ slurry. Processing of wet yellowcake slurry and UF₄ slurry will be discussed in Section 2.2.2.

The dry concentrate (yellowcake) is received from uranium mills in 55-gal steel drums. Each drum has an identifying number so that an accurate record may be maintained of the contents and analysis for purposes of uranium accountability and for billing or payment to customers. The drums are stored four to a pallet and up to three pallets high on an outside storage pad. The drums and pallets are strapped down for storm protection.

After weighing, each drum is emptied through a falling stream sampling unit. This unit consists of two samplers in series, each taking a small cut. The first cut provides an initial sample weight of about 1 percent or less of the contents of a drum. The sample is split down again by a factor of about 50 and is collected in trays. The material collected in the trays is processed to a final sample pulp by the operations of drying, pulverizing, riffling, sieving, and blending as needed.

The remaining ore concentrate from each drum can be fed directly into the two digester feed hoppers or can be redrummed for later use. The hoppers have a combined capacity of about 50 tons, which is one day's output of the sampling system.

All operations in the dry yellowcake sampling area of the main building involving open drums or uncontained yellowcake powders are performed under dust collection hoods.

2.2.2 Dry Yellowcake Digestion

Dry yellowcake digestion is the primary digestion process in the SFC facility. This is a batch process where the concentrate is reacted with preheated nitric acid in one of three 4000-gal digestion tanks to convert the uranium, present in the form of oxides or diuranates, to uranyl nitrate solution.

Initially, about 1400 gal of nitric acid of up to 60 percent concentration is pumped from the tank farm to a digester tank. Steam coils in the digester tank are used to raise the acid temperature to between 190 and 220°F. When the proper temperature is attained, dry concentrate is fed from the digester feed hoppers into the digester tank using screw-type feeders. About 12,000 lb of concentrate is added to the nitric acid. The normal batch volume is about 3100 gal. An agitator is operated continuously to promote both thorough mixing of the concentrate in the acid and thermal equilibrium. The reaction of the concentrate with nitric acid is exothermic, and water-cooled coils are used to remove excess heat to control the solution temperature. The reaction is also accompanied by evolution of nitrogen oxides. When digestion is complete, the slurry is transferred to an adjustment tank and the digester tank and transfer lines are flushed with water. The adjustment tank provides a holdup period for cooling of the slurry, for adjustment of the chemistry of the slurry to prevent uranium loss as precipitates, and to improve uranium recovery in the subsequent solvent extraction process. Composition adjustment may involve modifying the acidity of the uranyl nitrate slurry or the addition of special feed solutions from other digestion processes discussed below.

Wet yellowcake slurry containing approximately 38 wt percent of water is delivered from some uranium mills in stainless steel cargo tanks. The cargo tank arrives at the plant about 60 percent full. Pumps and piping enable transfer of the slurry to plant receiving tanks and recycle to the cargo tank.

In batch quantities, 60 percent HNO_3 is fed into the cargo tank until the solids in the yellowcake slurry are completely dissolved. This uranyl nitrate solution is then pumped to a receiving tank where the contents are weighed, sampled, and transferred to a 10,000-gal storage tank. Subsequently, this solution is pumped to the uranyl nitrate feed preparation tanks for final adjustment of solution chemistry before entering the solvent extraction operation (Fig. 2.2.1).

Tanks, pumps, and piping for processing the wet yellowcake slurry and the uranyl nitrate solutions are mostly enclosed in a separate building surrounded by curbs that can contain the total volume of the tanks in case of tank rupture. The parking area for the cargo tank is also curbed to contain the entire tank volume.

The UF₄ slurry with about 50 percent water is received in 55-gal drums with polyethylene liners. These drums are stored in a building with a curbed foundation that is adjacent to the west wall of the yellowcake slurry receiving building. The drums are emptied into a 4500-gal digester tank containing an $HNO_3-A1(NO_3)_3$ solution. The resulting uranyl nitrate solution can be stored or pumped to the feed preparation tanks for final adjustment before entering the solvent extraction operation (Fig. 2.2.1). Nitrogen oxide and hydrogen fluoride gases evolve during this digestion process.

A miscellaneous batch digester 1500 gal is used to recover uranium from fluorination ash and dusts collected in dust control systems by dissolution in a HNO_3 -Al(NO_3)₃ solution. Fluorination ash is derived from the exit filters of the fluorination reactors which are blown down periodically into ash receivers below the filters. Batch weight in the miscellaneous digester is approximately 2000 lb of ash [approximately 750 lb of uranium], and the resulting solution is blended with the primary digestion product from yellowcake digestion and fed to the solvent extraction system. Again, nitrogen oxides and hydrogen fluoride are produced.

The liquids of all digester tanks discussed above are eventually pumped into the solvent extraction system for purification of the uranium. There are no solid wastes from the digestion operations. Off-gases from the digester tanks are cooled with dilute HNO_3 and treated to remove most of the acid fumes, water vapor, and noxious nitrogen oxides. Condensed vapors are returned to the




Figure 2.2.1 Process diagram for uranium hexafluoride (UF₆) production from yellowcake (U_3O_8), including process inputs and waste product treatment

S

ES 6077

digestion tanks, caustic chemical-waste-containing fluoride is collected for disposal, and the remaining off-gases are piped to the HNO_3 recovery plant for recycle of remaining HNO_3 vapors.

The primary and miscellaneous digester tanks have spill curbs to contain any accidental leakage. Drainage from these curbed areas can be pumped into other digester tanks for processing.

2.2.3 Solvent Extraction

The uranyl nitrate slurry is fed from the adjustment tanks to the solvent extraction system. The uranyl nitrate is processed by countercurrent extraction in a series of pumper-decanters using 30 percent vol tributyl phosphate (TBP) in hexane solvent. Two control criteria have been established for the solvent extraction process. First, the digester slurry feed rate is adjusted at constant solvent flow to assure an acceptably low uranium loss in the raffinate waste stream. Second, a high uranium saturation level is maintained in the organic extraction product to obtain satisfactory separation of uranium from the impurities.

There are two liquid streams flowing from the series of pumper-decanters: the raffinate waste stream with very low uranium content and the organic extraction product that includes most of the uranium. The raffinate stream passes from the solvent extraction system to a raffinate decanter vessel where it is washed with hexane. The raffinate is then discharged to one of the four settling basins of clarifier for treatment and disposal.

The organic extraction product containing the uranium is then washed in a twostage mixer-settler scrub system using 0.01 M HNO_3 . The aqueous phase with impurities from the first-stage wash is recycled to the solvent extraction tanks, and the washed organic product is again contacted with 0.01 M HNO_3 in the re-extraction column where the uranium is stripped from the organic phase and transferred to an aqueous solution. The aqueous uranyl nitrate containing 80 to 100 g of uranium per liter is washed with hexane to remove TBP and transferred to the evaporation system feed tank.

The hexane from the raffinate wash, re-extraction column, and the aqueous uranium product wash is processed in the solvent rework system for recycle to the solvent extraction system.

2.2.4 Solvent Rework

The impure hexane solvent waste stream from the solvent extraction system is processed in the solvent rework (purification) system. Most of the residual uranium is recovered by washing the hexane with 25 percent ammonium sulfate. The TBP hydrolysis products, dibutyl phosphate and monobutyl phosphate, and remaining uranium are then removed as sodium salts by scrubbing the hexane with 5 percent sodium hydroxide (NaOH). The purified hexane is then adjusted for recycle to the solvent extraction system by addition of new hexane and/or TBP.

2.2.5 Evaporation and Boildown

From the evaporation system feed tank, the aqueous uranyl nitrate liquor flows through a feed condensate heat exchanger to a combination evaporator-condenser. In this equipment, the uranium concentration is increased to about 500 g/L as water and HNO_3 are evaporated. These vapors are compressed, recycled through the evaporator condenser to assist in the evaporation of the liquor, and condensed. The condensate collects in the recovered acidic condensate tank and is then pumped through the feed condensate heat exchanger to heat the incoming uranyl nitrate liquor. The acidic condensate tank is vented to dispose of small amounts of saturated steam, traces of hexane, and noncondensable gases. The acidic condensate is then stored for reprocessing and recovery of the HNO_3 .

The concentrated uranyl nitrate liquor ($\sim 500 \text{ g/L}$) is collected in a 2000-gal surge tank and held at 200°F until transferred to the boildown system.

The boildown system consists of three tanks in which batches of product from the evaporation system are heated using steam coils to boil off additional water and nitric acid, which increases the uranium concentration in the aqueous solution to about 1200 g/L. Vapor from the boildown tanks is condensed and the acidic condensate (0.1 M nitric acid) is transferred to a storage tank for subsequent recovery and recycle of the acid. The concentrated uranium product is pumped from the bottom of the boildown tanks into a recirculation piping system that permits feed to the denitration system.

2.2.6 Denitration

The concentrated product, uranyl nitrate hexahydrate (UNN), is bled from the boildown recirculation loop to the denitration system at the rate of 1.1 gpm. The UNH is heated in enclosed denitrator trays to a temperature of about 550°F with a slight negative pressure. The UNH is decomposed according to the following approximate reaction:

 $UO_2 (NO_3)_2 (n H_2 0) \rightarrow UO_3 + NO + NO_2 + O_2 + n H_2 0.$

The uranium trioxide (UO_3) produced in the denitrator trays consists of a free flowing, minus 8 mesh material that is discharged over an internal well and through seal legs to an enclosed screw-type system conveyor. The granular material is then ground in a pulverizer to a powder finer than 60 mesh (250 μ m) and discharged into the UO₃ storage bin.

The off-gases from the denitrator trays are scrubbed with 40 percent HNO_3 for removal of particulates and a portion of the nitrogen oxides and water vapor. The scrubbing fluid is cooled to $100^{\circ}F$ and recirculated to the off-gas scrubbers. Excess nitric and scrubbing fluids are returned to the HNO_3 storage tank.

The conveyor system, pulverizer, and storage bins are vented through a bin vent baghouse. Fine particulates of UO_3 trapped in the baghouse are periodically blown back into the storage bin by reverse air flow.

2.2.7 Reduction

The UO₃ is transferred by screw conveyor to the UO₃ feed bin which is slightly pressurized with nitrogen. The UO₃ is then moved by screw conveyor and charged into a two-stage fluid-bed reactor at about 2000 lb/h. Cracked ammonia (hydrogen and nitrogen) from an ammonia dissociator is introduced at the bottom of both stages. The reduction of UO₃ with hydrogen proceeds according to the following reaction:

 $UO_3 + H_2 \rightarrow UO_2 + H_2O$ (vapor) + heat.

Sufficient cracked ammonia is introduced to assure complete reduction. The thermal energy (heat) released in the reaction is used to help maintain the reaction temperature of 950-1100°F. Electrical furnaces are used to attain initial operation of the reactors and air-cooled coils are used to maintain temperature control.

Excess gas and UO_2 product from the reduction reactors are discharged to a UO_2 filter bin where the UO_2 is filtered from the gas by porous stainless steel filter elements that can remove 98 percent of all particles larger than 0.7 μ m. The filters are periodically cleansed by a reverse flow of nitrogen.

Solids collected in the filter bin are discharged to a screw conveyor and subsequently transported to the HF recovery reactor or directly to the $\rm UO_2$ seal bin for storage.

Exit gases from the filter bin are discharged through another stainless steel filter element to remove final traces of UO_2 product. The filtered off-gases, consisting of ammonia, nitrogen, water vapor, hydrogen, and hydrogen sulfide, are piped to a reduction hydrogen burner before disposal to the atmosphere.

The HF recovery reactor is a screw-type conveyor where contact between UO_2 powder and unreacted HF from the filter bin of the subsequent hydrofluorination stage can reduce loss of fluoride to the environment. The off-gas from the HF recovery reactor is filtered to remove particulates and cooled to condense the HF. The HF condensate and noncondensable gases are routed to the HF scrubber in the gaseous waste disposal system. The UO₂ powder partially reacted with HF is discharged into the UO₂ seal bin.

2.2.8 Hydrofluorination

Uranium dioxide is screw-fed from the UO_2 seal bin to the first stage of the hydrofluorination system where the powder is contacted with a countercurrent flow of HF in a stirred fluid-bed reactor. The reactor wall temperature is maintained at about 400-650°F using electrical heaters and air blowers. A mixture of UO_2 , UF_4 product, excess HF, and water vapor are discharged from the first-stage reactor at a temperature of about 600-700°F into the UO_2 - UF_4 filter bin where the solids are separated from the gases using porous carbon filter elements. These filters remove 98 percent of all particles larger than 0.7 µm and are periodically cleaned by back-flushing with hot nitrogen gas. Gases discharged from the UO_2 - UF_4 filter bin are routed to the HF recovery reactor discussed in the previous section. The solids in the UO_2 - UF_4 filter bins are transferred by screw conveyor to the second-stage hydrofluorinator for contact with anhydrous HF vapor heated to about 1000°F. The second-stage reactor is

maintained at about 900° F with an electrical furnace and a cooling air blower. The UF₄ product, excess HF vapor, water vapor, and other gases are discharged from the second-stage reactor into the UF₄ seal bin where most of the UF₄ solids are collected and most of the gases and entrained fine particles are piped to the first-stage fluorinator.

The UF₄ powder is transferred from the UF₄ seal bin by screw conveyor to a UF₄ nitrogen lift chamber where gaseous nitrogen is used to purge entrained HF vapor from the UF₄ powder. The UF₄ and all the gases then enter the UF₄ filter seal bin where the gases are separated and piped to the HF scrubber in the gaseous waste disposal system. The UF₄ powders are transferred by screw conveyor to UF₄ storage bins.

2.2.9 Fluorination

The UF₄ is transferred to a two-stage fluorination system by enclosed conveyors. The UF₄, charged into the first stage at a temperature of 100° F, is converted to UF₆ by reaction with elemental fluorine at a temperature of 850° F according to the following reaction:

 $UF_4 + F_2 \rightarrow UF_6 + thermal energy.$

Excess reaction heat is removed from the reaction vessels using steam-cooled coils. About 15 percent excess fluorine over the stoichiometric conversion requirement is introduced into the first-stage reactor to maximize conversion of UF₄ to UF₆. Approximately 1 percent of the feed material is collected as unreacted UF₄ or intermediate fluorides in ash containers at the bottom of the first-stage reactors. Each first-stage reactor is periodically shut down, and the ash is ground and returned to the UF₄ storage bin.

The outlet gases from the first-stage fluorination reactors are cooled to about 300°F, and entrained solids are removed in a cyclone separator and by porous Monel filter elements in the F₂ reactor filters. The filters remove 98 percent of all particles greater than 0.7 μ m and are periodically cleaned using a backflow of air. The particulates are returned to the miscellaneous digester for reprocessing.

The filtered gas, containing UF₆, F₂, HF, O₂, and N₂, is passed through the primary cold traps at 30°F where approximately 90 percent of the UF₆ is condensed as a solid. The noncondensed gases are heated to 800°F and forced into the second-stage fluorinators (cleanup reactors) for reaction with an excess of the stoichiometric requirement of UF₄ at a temperature above 750°F. About 10 percent of the UF₄ will be unreacted and is removed from the bottom of the cleanup reactor by screw conveyor where the product is cooled to 400°F. The excess UF₄ is subsequently returned to the UF₄ storage bin for feed to the fluorinators. The system is designed, monitored, and controlled so that the excess F₂ from the primary fluorinator will be completely consumed by the reactions with UF₄ and any UO₂ or UO₂F₂ in the feed material.

The gases and entrained solids from the cleanup reactor enter the cleanup reactor filter where 98 percent of the particles greater than 0.7 μ m are removed by porous Monel filters (primary and backup). These filters are also periodically cleaned using a backflow of air.

The filtered gas is passed through the secondary cold traps for condensation of UF_6 as a solid at -65 to -75°F. The noncondensable gases are pumped to the HF scrubber in the gaseous waste disposal system. Fluorine is monitored in the noncondensable gas discharge as an element in the process control.

2.2.10 UF6 handling and shipping

When any primary or secondary cold trap is full, the UF₆ is melted and drained by gravity through porous Monel filters into evacuated 10-ton or 14-ton shipping cylinders where it slowly solidifies at ambient temperature. The piping from the cold traps to the cylinders consists of a 2-in.-diam sloped, steam-traced transfer header. The transfer header and filter are connected to the cylinder valve through a 3/4-in.-diam copper tubing pigtail.

Two transfer headers and filling stations are in use. The headers are manifolded so that each filling station can handle UF₆ from either the primary or the secondary cold traps. The primary cold traps have a nominal capacity of 10 tons. After filling, the cylinders are weighed and transferred to the cylinder storage area in the yard. The cylinders are cooled for a minimum of 5 days before shipment.

APPENDIX 3.1.1

PHOTOS SHOWING THE LOCATION OF THE 14-TON CYLINDER AND THE RUPTURE

PHOTOS FOR APPENDIX 3.1.1

Figure	3.1.1.1A	The <u>Sequoyah Facility's Process Buildings Viewed from the</u> <u>South</u> . The sod has been stripped from the lawn in front of the building. Route 10 is at the upper right.
Figure	3.1.1.2A	The Sequoyah Facility, Viewed from the North-Northwest. The self-propelled crane is preparing to lift away the damaged steam chest from over the ruptured cylinder.
Figure	3.1.1.3A	Removing the Steam Chest from the Ruptured Cylinder. Both the cylinder and the steam chest are in the positions in which they were found after the cylinder ruptured. A temporary cover is in place over the rupture.
Figure	3.1.1.4A	The 52-inch Long Rupture. At its midpoint, the opening is about 8 inches wide. The cylinder wall is 5/8-inch thick steel. Water is draining out after rinsing out the cylinder.




Figure 3.1.1.1A The Sequoyah Facility's Process Buildings Viewed from the South




Figure 3.1.1.2A The Sequoyah Facility, Viewed from the North-Northwest.



Figure 3.1.1.3A Removing the Steam Chest from the Ruptured Cylinder

17




Figure 3.1.1.4A The 52-inch Long Rupture

APPENDIX 3.3.1

8

Ø

10.00

NRC'S PRELIMINARY NOTIFICATIONS OF EVENT OR UNUSUAL OCCURRENCE FOLLOWING THE ACCIDENT OCCURRING ON JANUARY 4, 1986

4

January 5, 1986

PRELIMINARY NOTIFICATION OF EVENT OR UNUSUAL OCCURRENCE -- PNO-IV-86-02

This preliminary notification constitutes EARLY notice of events of POSSIBLE safety or public interest significance. The information is as initially received without verification or evaluation, and is basically all that is known by the Region IV staff on this date.

FACILITY:	Licensee Emergency Classification:
Sequoyah Facility	Notification of Unusual Event
Sequoyah Fuels Corporation	Alert
(formerly Kerr-McGee Nuclear Corp.)	Site Area Emergency
Gore, Oklahoma	General Emergency
License No. SUB-1010	X Not Applicable

SUBJECT: RUPTURED URANIUM HEXAFLUORIDE CYLINDER

At approximately 11:30 a.m. on January 4, 1986, a cylinder containing approximately 29,500 pounds of uranium hexafluoride (UF6) ruptured, releasing its contents into the atmosphere and resulting in the death of one company employee from the inhalation of hydrogen fluoride (HF). As a result of this incident, the NRC has established an emergency response organization at the site to examine the resultant health and environmental impacts and an Augmented Investigation Team (AIT) to investigate the circumstances of the cylinder rupture.

According to statements from company officials, plant employees had overfilled the cylinder (maximum limit is 27,560 pounds) and were in the process of heating the cylinder to facilitate removing the excess uranium hexafluoride when the rupture occurred. An approximate three foot-long split occurred in the wall of the cylinder, allowing the contents of the cylinder to be released to the atmosphere in gaseous form. Company officials estimate that the contents of the cylinder emptied in about 30 to 40 minutes.

Contact with the air resulted in the formation of hydrogen fluoride (HF) and uranyl fluoride. Although it is believed, based on radiation surveys and sample results to date, that most of the uranyl fluoride was deposited in the immediate vicinity of the plant, the hydrogen fluoride, a toxic and corrosive chemical, formed a plume that was visible at more than a mile from the plant site before it was dispersed. Interstate 40, which runs east-west about a mile south of the plant, was closed for a short period of time after the incident by the Oklahoma Highway Patrol.

About 25 people, mostly company employees, were hospitalized as a result of the incident, but only a few were diagnosed as suffering from respiratory problems related to inhalation of hydrogen fluoride. Those exhibiting actual respiratory problems were transferred to a medical facility in Fort Smith, Arkansas, from a hospital in Sallisaw, Oklahoma. One person, a company employee, died enroute to Fort Smith. No additional deaths have been reported. Numerous company employees and local residents have reported to the Sallisaw

PNO-IV-86-02

hospital at the company's urging for initial screening. Many have been asked for urine samples which the company will analyze at its Technical Center in Oklahoma City.

The Company continues to perform radiation surveys at the plant and to perform surveys and take samples from the environment in the areas that the plume is believed to have traversed.

Although the rupture took place in a steam-heating facility outside the main process buildings, contamination was spread through most plant buildings by the building ventilation systems. The company continues to decontaminate those areas needed to manage its recovery activities.

Region IV's emergency response activities are being conducted in parallel with those of the Augmented Investigation Team. Both are headed by Region IV's Director of the Division of Radiation Safety and Safeguards. The complement of NRC personnel presently at the site include six technical staff Region IV (two from the Uranium Recovery Field Office), one each from NMSS and IE, and Region IV's Public Affairs Officer. The team is focusing on the events that led to the cylinder rupture, the resultant health and environmental impacts and the recovery activities that are currently taking place. Region IV's Incident Response Center and the Headquarters Operations Center were activated and staffed promptly after receiving notification of this incident. The company has agreed to leave the equipment and facilities related to this incident untouched and not to resume facility operations without NRC concurrence. The team is making arrangements to interview all plant employees who were at the facility at the time of the incident and to have those interviews transcribed.

Arrangements have been made through the Department of Energy for an aerial radiation survey. Both the state of Oklahoma and the state of Arkansas have been contacted regarding this event. The Oklahoma Department of Health dispatched two individuals to the plant site on Saturday. The NRC's response team has been coordinating its activities with representatives of the state of Oklahoma and the licensee.

National and local news media have run stories on this incident. The company and NRC have been responding to media inquiries since shortly after the incident. The company issued a news announcement at about 2:00 p.m. Saturday. A press conference was arranged by Congressman Mike Synar of Oklahoma and took place at about 2:00 p.m. Sunday in Muskogee. The NRC's response team leader and Oklahoma state officials took part in the press conference. The response team leader also briefed Congressman Synar separately on Sunday morning.

A hearing on a license amendment for this facility that was to have taken place on January 7 and 8 in Fort Smith, Arkansas, has been indefinitely postponed by the presiding Atomic Safety and Licensing Board Judge. PN0-IV-86-02

Region IV was notified of this incident at about noon January 4 by the Headquarters Operations Officer. This information is current as of 7:00 a.m. on Monday, January 6. This PN will be updated periodically.

-3-

CONTACT: E. F. Bates, FTS: 728-8100

LIMINARY NOTIFICATION OF EVENT OR UNUSUAL OCCURRENCE -- PNO-IV-86-02A

This preliminary notification constitutes EARLY notice of events of POSSIBLE safety or public interest significance. The information is as initially received without verification or evaluation, and is basically all that is known by the Region IV staff on this date.

FACILITY:	Sequoyah Facility Sequoyah Fuels Corporation (formerly Kerr-McGee Nuclear Corp.) Gore, Oklahoma License No. SUB-1010	Licensee Emergency Classification: Notification of Unusual Frent Alert Site Area Emergency General Emergency X Not Applicable
-----------	--	--

SUBJECT: RUPTURED URANIUM HEXAFLUORIDE CYLINDER (UPDATE)

The NRC's Augmented Investigation Team (AIT) and other response personnel continue to investigate the circumstances of the January 4, 1986, rupture of a cylinder containing approximately 15 tons of uranium hexafluoride (UF6) at the Sequoyah Fuels Corporation UF6 production facility near Gore, Oklahoma. The NRC team is evaluating both the events leading to the rupture and the radiological impact of the incident. The company has agreed not to resume processing operations without NRC concurrence. Nonradiological effluent processing systems are continuing to operate. This is an activity unrelated to the incident.

The investigation team began interviewing plant personnel on Monday, January 6, and has begun compiling a chronology of events. The team believes it is close to identifying the fundamental cause of the cylinder rupture. At a news briefing Tuesday, the company acknowledged that the plant's operating procedures prohibit heating a UF6 cylinder that has been overfilled.

An additional wisp of hydrogen fluoride gas was observed evolving from the residual UF6 in the ruptured cylinder but control measures by plant personnel have been successful in stopping it. No health or safety importance is attached to this release.

The company has formulated plans for cleaning up the processing portions of the facility and has already begun cleanup of the plant property itself, including the removal of sod from the lawn at the south end of the plant. The NRC team is reviewing the company's cleanup plans to ensure adequate protection of workers and ensure that no equipment essential to the investigation is disturbed. The team has asked the company to consider alternative methods of collecting fluids that will result from the washdown operations at the plant, rather than sending them through the plant's normal processing system and to the Illinois River.

Aerial radiation monitoring by EG&G, arranged through the Department of Energy, is expected to begin Wednesday. The survey, which will be conducted using a red, white and blue helicopter flying at an elevation of 150 feet, will cover an approximate

A/RA (PSCheck 1/7/86

12-square-mile area surrounding the plant. The company is taking steps to inform the public about the nature and purpose of these flights. The results of the aerial survey will be used to confirm ground survey data that appears to indicate that no significant amounts of uranium were deposited off the plant site. Initial review of environmental data available to the team at this time confirms early prediction that no significant offsite radiological impacts are expected.

The company reported today that all individuals hospitalized as a result of this incident have been released with the exception of one who remains hospitalized for reasons unrelated to the event.

The NRC team, now totaling 15 people at the site, is requesting additional personnel with expertise in instrument maintenance and surveillance, QA/QC and laboratory analytical procedures. These technicians are expected to arrive at the site tomorrow (January 8).

News briefings are being conducted daily in the vicinity of the facility. NRC, the company and the state of Oklahoma are represented at these briefings. The NRC team continues to work in cooperation with the company and the state of Oklahoma and to keep the state of Arkansas informed of the status of the recovery operations.

This information is current as of 1:00 p.m. on Tuesday, January 7, 1986. The team plans to update the information in this PN daily.

CONTACT: J. B. Baird (FTS: 728-8100)

PRELIMINARY NOTIFICATION OF EVENT OR UNUSUAL OCCURRENCE -- PNO-IV-86-028

This preliminary notification constitutes EARLY notice of events of POSSIBLE safety or public interest significance. The information is as initially received without verification or evaluation, and is basically all that is known by the Region IV staff on this date.

FACILITY:	Sequoyah Facility Sequoyah Fuels Corporation (formerly Kerr-McGee Nuclear Corp.)	Licensee Emergency Classification: Notification of Unusual Event Alert
	Gore, Oklahoma License No. SU8-1010	Site Area Emergency General Emergency

SUBJECT: RUPTURED URANIUM MEXAFLUORIDE CYLINDER (UPDATE)

The NRC's Augmented Investigation Team (AIT) and response team personnel continued to analyze information and data collected in the aftermath of the January 4, 1986 rupture of a cylinder containing approximately 15 tons of uranium hexafluoride (UF6). The rupture resulted in the death of one Sequoyah Fuels Company (SFC) employee. The AIT has completed interviews with 12 plant personnel but has temporarily discontinued interviewing to permit team members to analyze information gathered to date.

Based on the information obtained to date and visual observation of the ruptured cylinder, the best judgement of the team is that the cylinder failed by hydraulic rupture. The team has found no information which would indicate the cylinder was defective. On the contrary, the extensive deformation experienced by the cylinder before it ruptured indicates that it was subjected to internal pressures far beyond the pressures for which it was designed.

Such large pressures would be produced hydraulically when the liquid contents expanded to a volume greater than the volume of the vessel. The current state of the analysis at the site indicates that this condition occurred.

The company's cleanup of plant buildings and property continues. Decontamination of three areas of the plant process building -- the access areas, the change room and the maintenance shop -- has begun. Removal of contaminated sod from the lawn south of the main plant buildings continues and is expected to be completed soon. The NRC is reviewing the company's plans for disposition of the waste resulting from the cleanup. SFC is analyzing alternatives to discharging cleanup liquids (although expected to be within uranium concentration limits for discharge) to the Illincis River. Regarding contaminated sod, SFC and NRC are examining the possibility of transporting it to a company uranium mill tailings site in New Mexico.

The NRC team has completed a preliminary description of the probable trajectory of the plume. Using data from surveys and sample analysis, efforts are underway to define the extent of fluoride and uranium deposition. The company

Romartin:jc 1/8/88

continues to analyze the fluoride impact on animals in the path of the plume, but has seen no observable effects. Some fluoride damage to cedars has been observed by a specialist employed by the company. The NRC plans to produce a detailed report examining the radiological consequences of this incident. There is no data at this time to indicate any significant offsite radiological impacts.

The aerial radiation survey of the area surrounding the plant did not begin today as expected because of an FAA airworthiness directive that kept the helicopter grounded. This issue is expected to be resolved promptly.

The company reports today that all employees were were hospitalized as a result of this incident have been released from the hospital. Medical consultants to the company continue to monitor affected employees for evidence of kidney damage from uranium ingestion. Funeral services for the employee who died as a result of inhaling toxic hydrogen fluoride vapors took place today in Vian. Oklahoma.

News media interest in the incident has remained steady. Another briefing was held at 11:00 a.m. today in Gore. Oklahoma. The NRC team continues to work with the state of Oklahoma representatives at the site and to keep the Arkansas Department of Health informed of the status of the recovery operations. The team continues to brief Congressman Synar's staff daily. NRC, Oklahoma and company representatives have been requested to attend this evening's special meeting of the Gore Public Schools Board of Education to address concerns about location of Carlile School about one mile east of the plant site.

This information is current as of 1:00 p.m. on Wednesday, Jenuary 8, 1986.

CONTACT: J. 8. Baird (FTS: 728-8100)

January 9, 1986

PRELIMINARY NOTIFICATION OF EVENT OR UNUSUAL OCCURRENCE -- PNO-IV-86-02C

This preliminary notification constitutes EARLY notice of events of POSSIBLE safety or public interest significance. The information is as initially received without varification or evaluation, and is basically all that is known by the Region IV staff on this date.

FACILITY:	Sequoyah Facility	Licensee Emergency Classification:
	Sequoyah Fuels Corporation	Notification of Unusual Event
	(formerly Kerr-McGee Nuclear Corp.) Gore, Oklahoma	Aiert Site Area Emergency
	License No. SUB-1010	X Not Applicable

SUBJECT: RUPTURED URANIUM HEXAFLUORIDE CYLINDER (UPDATE)

The NRC's Augmented Investigation Team (AIT) is nearing the completion of its fact-finding activities and continues to evaluate the information that has been gathered to date. The site team leader announced to the press today the AIT's preliminary conclusion on the cause of the UF6 cylinder rupture. It appears that early next week the team can present estimates of the size and concentrations associated with the release and some conclusions at to the health and environmental impact of the incident.

The aerial surveys to monitor the vicinity of the plant for the presence of uranium and for any indication of fluoride damage to vegetation (using infrared photography) began today. The company announced that damage to offsite cedars previously reported to have been caused by fluoride was actually caused by a power company spraying a defoliant beneath its power lines.

The company continues to make progress in cleaning up the process areas and the plant property. Plans have been made to decontaminate the ruptured cylinder prior to its removal to an offsite location for a metallurgical examination. Arrangements are being made to obtain a complete inspection and test of the scale that was to weigh the cylinder.

Plant buildings were evacuated this morning when a fire alarm sounded. The company determined that the alarm was caused by low nitrogen pressure in the automatic fire suppression system in the solvent extraction building. The NRC is assisting the company in making arrangements to obtain necessary nitrogen supplies because the nitrogen must be delivered to an area of the facility that remains contaminated.

Senator Don Nickles (R-Oklahoma) visited the plant site this morning and was briefed by the NRC team leader and company officials. Senator Nickles, the NRC team leader, company and Oklahoma representatives then briefed the news media at an 11:00 a.m. news conference in Gore, Oklahoma. Last night in Gore, the NRC team leader appeared before a special meeting of the Gore Public Schools Board of Education that was called to assess the long-term safety of students in the Carlile elementary school, which is

RIV:RA RDMartin:jc 1/9/86 located about a mile east of the plant. The school was not in the direction of the plume from this incident and was surveyed in detail to confirm that no radioactivity from the incident had been deposited there. The school board delayed until Monday, January 13, a decision on whether to shut down the Carlile school and move students to other facilities.

This information is current as of 1:00 p.m. on Thursday, January 9, 1986.

CONTACT: J. B. Baird (FTS: 728-8100)

PRELIMINARY NOTIFICATION OF EVENT OR UNUSUAL OCCURRENCE -- PNO-IV-86-02D

This preliminary notification constitutes EARLY notice of events of POSSIBLE safety or public interest significance. The information is as initially received without verification or evaluation, and is basically all that is known by the Region IV staff on this date.

FACILITY:	Sequoyah Facility	Licensee Emergency Classification:
	Sequoyah Fuels Corporation (formerly Kerr-McGee Nuclear Corp.)	Notification of Unusual Event Alert
	Gore, Oklahoma License No. SUB-1010	General Emergency
		X Not Applicable

SUBJECT: RUPTURED URANIUM HEXAFLUORIDE CYLINDER (UPDATE)

The NRC's response team met with company officials this morning to discuss detailed plans for decontamination of the ruptured cylinder and the handling of any residual UF6 as the company prepares the cylinder for eventual transport to the Lawrence Livermore Laboratory for metallurgical examination. The principal concern of the team has been that adequate steps are taken to protect the health and safety of workers and others involved in this activity. These concerns have been satisfactorily addressed. The decontamination is underway.

Sequoyah Fuels Corporation (SFC) employees continue to decontaminate service areas of the plant process building. Removal of contaminated sod from the south lawn of the plant is complete. SFC officials assured the NRC this morning that decontamination activities over the weekend will involve only those portions of the facility for which NRC has already reviewed and approved plans. In addition, SFC plans to contract for additional health physics technician support to more carefully implement contamination control procedures as the cleanup of the plant and site continue. The NRC will maintain a presence at the facility over the weekend to monitor those activities.

The NRC's Augmented Investigation Team (AIT) has completed interviews with plant employees and has developed a preliminary sequence of events. The sequence is attached to this PN. The AIT's activities at the facility are virtually complete. Health physic coverage and environmental assessment activities will continue through the weekend. NRC presence at the site beyond early next week is expected to be limited to radiation protection specialists who will continue to monitor the decontamination of the plant.

Aerial surveys to monitor the vicinity of the plant for the presence of uranium and for any indication of fluoride damage to vegetation continued today. DOE expects the aerial radiation surveys to be completed sometime Saturday.

The Regional Administrator and other NRC staff members will brief the Commissioners on this incident at a public meeting beginning at 1:00 p.m. (EST) today. No news briefings at the plant site are scheduled today or over the weekend.

RIV:RA RDMartin:jc 1/10/86 The NRC response team continues to work in close cooperation with the state of Oklahoma Department of Health representatives at the plant site and to keep Arkansas Department of Health officials and Congressman Synar's staff informed of the status of recovery activities.

This information is current as of 1:00 p.m. on Friday, January 10, 1986. Additional PNs will be issued as necessary.

ATTACHMENT: Sequence of Events

CONTACT: J. B. Baird (FTS: 728-8100)

Sequoyah Facility, Sequoyah Fuels Corporation (formerly Kerr-McGee Nuclear Corp.), Gore, Oklahoma, License No. SUB-1010

SEQUENCE OF EVENTS

9/20/85

Cylinder No. E-2047 received at Sequoyah facility.

9/27/85 Cylinder passes 20-point quality inspection conducted by site engineer.

1/3/86 (10:00 am) Cylinder passes second 20-point quality inspection conducted by day shift operator. Cylinder is moved to south filling bay within the plant process building and has begun to be filled with the UF, product. The operator adjusts the counterweight of the scales on which the cylinder is positioned so that the scale dial is indicating no weight. That is, the counterweight is positioned to counterbalance the tare weight of the empty cylinder so that the scales will directly measure the net weight of the product.

> (3:30 pm) At the end of his shift, the day shift operator records that the cylinder has been loaded with 11,230 pounds of product and is still being filled.

(11:30 pm) At the end of his shift, the evening shift operator records that the cylinder has been loaded with an additional 12,200 pounds of product. The total weight of product in the cylinder is correctly noted as being 23,430 pounds.

1/4/86 (2:15 am) Midnight shift operator continues filling of cylinder. Although the targeted net weight of the cylinder is 27,500 pounds of product, the operator is unable to fill it beyond 26,400 pounds. Upon closer investigation, the operator observes that one wheel of the cart holding the cylinder is not fully on the scales which are set to indicate the net weight of the product in the cylinder. After restoring the cart to a position fully on the scales, the scale indicator pegs out at 29,500 pounds.

The operator informs the shift supervisor that the cylinder has been overfilled. The supervisor orders the operator to

begin removing product from the cylinder using a vapor evacuation technique.

- 1/4/86 (6:00 am) The operator moves the counterweight on the scales so that the scale dial indicator will be freed from its pegged position and will indicate the rate at which product is removed from the cylinder.
 - (6:15 am) Evacuation of product from the cylinder begins. The scale indicates that 150 pounds of product is removed in the succeeding 10 minutes. The operator can hear the flow of material through the cylinder valve.
 - (7:45 am) At the end of the midnight shift, evacuation of UF₆ from the cylinder is still in progress.

The operator informs the oncoming day shift operator of the problem.

(8:45 am) The day shift operator is unable to draw off any more product from the cylinder, presumably because the UF₆ has begun to solidify. The operator confers with the day shift supervisor who instructs the operator to transfer the cylinder to the southwest steam chest outside and north of the process building. The UF₆ will again be liquified in the steam chest before the cylinder is returned to the filling bay for further product extraction. The supervisor instructs the worker to leave the cylinder in the chest for 6 hours. (Heating an overfilled cylinder is prohibited by company procedures.)

Before transferring the cylinder from the scales, the counterweight is moved to its original position which had been demarked on the slide bar with a pen by the operator on the previous shift. The scale now again indicates 29,500 pounds.

- (9:15 am) The operator uses a fork lift truck to transfer the cylinder to the steam chest. Steam heating of the cylinder is begun.
- (11:30 am) The cylinder ruptures while in the steam chest. The UF quickly vaporizes forming UO_2F_2 and hydrofluoric acid. The operator in a scrubber building 50 feet from the steam

chest sustains lethal injuries caused by the acid fumes. Workers further away sustain less serious injuries. The plant is evacuated of approximately 40 workers.

The vapor release which continues for approximately 40 minutes is carried south by south-east by a gusting 25 mph wind.

January 13, 1986

PRELIMINARY NOTIFICATION OF EVENT OR UNUSUAL OCCURRENCE -- PNO-IV-86-02E

This preliminary notification constitutes EARLY notice of events of POSSIBLE safety or public interest significance. The information is as initially received without verification or evaluation, and is basically all that is known by NRC staff on this date.

FACILITY:	Sequoyah Facility Sequoyah Fuels Corporation	Licensee Emergency Classification: Notification of Unusual Event
	(formerly Kerr-McGee Nuclear Gore, Oklahoma License No. SUB-1010	

SUBJECT: RUPTURED URANIUM HEXAFLUORIDE CYLINDER (FINAL UPDATE)

The NRC's response team has begun to receive results of soil and vegetation samples taken from more than 100 locations in a broad area extending as far as 10 miles south of the facility (in the direction of the plume). Samples are being analyzed by the Kerr-McGee Technical Center, the state of Oklahoma Department of Health and the Oak Ridge National Laboratory (on behalf of NRC). Sample results to date indicate no deposition of uranium in detectable quantities beyond property owned by the licensee.

Preliminary results of aerial radiation surveys in a 25-square-mile area confirm the absence of uranium in concentrations above background except in the immediate areas of the facility. NRC and state of Oklahoma officials continue to analyze soil and vegetation samples for uranium and fluorine content. Company medical consultants also continue to analyze the results of uranium analyses done on the urine of both plant employees and members of the public.

The initial decontamination of the ruptured cylinder was completed without incident on Friday, January 10. As expected, a small amount of hydrogen fluoride vapor was released to the air during the operation of flushing residual materials from the UF6 cylinder. As a precaution, the state of Oklahoma Department of Health notified local officials prior to this activity. Other decontamination activities performed over the weekend focused on returning plant office space and the main gate area to normal conditions. The company plans to develop a program to periodically resurvey areas that have been cleaned to ensure that contamination is not spread during cleanup operations.



PNO-IV-86-02E

.

With the exception of radiation specialists who will maintain coverage during the decontamination of the facility, the remaining members of the response team plan to complete their onsite activities on Tuesday, January 14, and depart the site that afternoon.

The team plans a final news briefing at 3 p.m. Monday in Gore and a final briefing of U.S. Representative Mike Synar's staff on Tuesday. The team prepared a letter to the president of the Gore school board describing the results of radiation surveys of the Carlile School in preparation for a Monday evening meeting of the board to decide whether to move Carlile students to other facilities. All surveys at the school confirmed that no radiological contamination was present.

This information is current as of 1 p.m. on Monday, January 13, 1986. This is the final PN on this subject.

CONTACT: J. B. Baird (FTS: 728-8100)

APPENDIX 4.1.2

۱

METEOROLOGY AND ATMOSPHERIC DISPERSION

APPENDIX 4.1.2

ATMOSPHERIC DISPERSION AND METEOROLOGY

4.1.2.1 Atmospheric Dispersion Modeling

The atmospheric transport, diffusion, and deposition models used for assessing the impacts of the atmospheric releases are part of the Atmospheric Release Advisory Capability (ARAC) operated by the Lawrence Livermore National Laboratory under the auspices of the Department of Energy. The models^{7,8} used were the Mass-Adjusted Three Dimensional Wind Field (MATHEW) model and the Atmospheric Diffusion Particle-In-Cell (ADPIC) model. These models are fully described in References 7 and 8. The regional wind field was developed from information compiled from National Weather Service stations within 200 km, although data from the site dominated the wind field and consequently plume movement close to the faiclity. Data from two other stations affected the movement of the plume well-away from the facility. These stations are located in Fort Smith Arkansas (approximately 70 km east-southeast), and Page, Oklahoma (approximately 100 km south of the facility).

The influence of terrain and airflow trajectory is considered in determining the 3-dimensional wind field, with horizontal resolution of terrain of onetwentieth of the overall grid size and vertical resolution of terrain of 50 meters independent of grid size. Calculations of atmospheric dispersion were performed for three grid sizes: 10 km, 40 km, and 80 km. The total vertical depth modeled was 700 meters. Because of different terrain cell sizes and varying grid cell volumes close to the facility, contours of air concentration and deposition differ slightly in magnitude and location close to the facility for each computational grid used.

The ADPIC Model is based on a particle-in-cell technique with concentration represented by Lagrangian-marker particles inside a fixed Eulerian grid. From observations at the facility during the accident, the initial cloud was assumed to be an ellipsoidal volume 21 meters high (in the vertical direction), 76 meters wide (in the cross-wind direction), and 4 meters long (in the along wind direction). Uranyl fluoride was modeled as particles with a log-normal size cistribution, a median diameter of one micron, and a standard geometric deviation of 1.5 microns. A deposition velocity of 1 cm/sec was assumed. Hydrogen fluoride was modeled as a vapor with a deposition velocity of 1 cm/sec. The assumptions regarding the character and behavior of uranyl fluoride and hydrogen fluoride were developed to simplify for modeling purposes the complex chemical transformations and atmospheric processes related to condensation, precipitation, and sublimation. The assumptions appear reasonable given the nature of the release and meteorological conditions at the time of the accident; however, the relatively simple characterization of such complex processes is a major source of uncertainty in the atmospheric dispersion models used for this evaluation.

For example, when the cylinder ruptured in the steam chest, it released UF₆ in three physical states, i.e., some was a gas, some liquid, and some solid. It was reported that small pieces of solid UF₆ were scattered widely around the steam chest area. Some of the UF₆ reacted immediately with the steam and other water in the steam chest and in the surrounding air, but there was not

enough water present for an immediately reaction with all 13,400 kg UF6. The total water requirement would be about 1,370 kg water. Thus, the initial cloud (before hosing down by SFC workers) contained some unreacted UF6 as well as UO_2F_2 and HF. As the airborne UF₆ was carried downwind, it reacted with the moisture in the air. How far the UFs may have been transported before the reaction was essentially complete is not known, but even though the air had a relatively low moisture content, the UF6 would not have been transported more than a few km downwind before encountering sufficient atmospheric moisture to complete the reaction. If inhaled, the UF6 would quickly react with moisture in the lungs, bronchial passages, etc., to form UO2F2 and HF; the reaction also liberates heat. If deposited on surfaces (e.g., vegetation) the UF_6 would also react there with available moisture within a relatively short time period. The release has been modeled as the reaction products $\mathrm{UO}_2\mathrm{F}_2$ and HF; quantities of uranium and fluoride in UF6 are accurately presented as the sum of the quantities of UO2F2 and HF.

Furthermore, the cloud as formed within the first few seconds after cylinder rupture was exceptionally dense. This great density would promote, for example, rapid particle agglomeration and deposition in the vicinity of the plant buildings. No attempt was made to model such circumstances; instead, the dispersion modeling was designed to represent the cloud well when it had traveled a kilometer or more downwind.

However, even though the meteorological conditions were favorable to mathematical modeling of the cloud of released chemicals, the models provide a description based on experience and do not describe in second-by-second detail what actually occurred in the cloud of chemicals. Because the greater part of the chemicals was released in the first few minutes after cylinder rupture, it is possible that especially during this brief time period some small portions of the cloud carried substantially higher concentrations downwind a few kilometers before atmospheric mixing made the cloud more uniform and more dilute.

Based on initial observations of the accident and consideration of the process involved, type of cylinder, and the physical form of the material in the cylinder, 75 percent of the material was assumed to be released in the first 5 minutes following the rupture of the cylinder, with the remaining 25 percent of the material released during the next 40 minutes. Model calculations were performed for a unit release (1 gram) with a uniform release rate of 0.75 gram over 300 seconds (5 minutes) and of the remaining 0.25 gram over 2400 seconds (40 minutes). Use of a release normalized to one gram was chosen by the Ad Hoc Task Force for the dispersion modeling so that modeling results could contribute to the selection of the quantities to be used in the modeled exposure calculations to best represent their judgment of the quantities transported downwind to distances of a kilometer or more.

The dispersion simulation results produced by LLNL-ARAC are provided here in three forms: (1) contours of integrated air concentrations, (2) contours of cumulative deposition concentrations, and (3) instantaneous air concentrations. All three are presented for UO_2F_2 and HF separately; all three are presented normalized to a release of one gram of UO_2F_2 or one gram of HF. Thus to estimate meaningful exposures, these normalized values must be multiplied by the estimated source term, e.g., 5,900 kg UO_2F_2 . The contour plots for UO_2F_2 serve to provide an estimate for uranium exposures; to estimate fluoride exposures,

some combination of the results for both HF and UO_2F_2 must be used; the Task Force chose to use the sum of 1,500 kg HF and 5,900 kg UO_2F_2 . Results in the form of normalized quantities deposited per square meter are also provided; the same considerations apply. Figures 4.1.2.1A through 4.1.2.3A represent (in units of seconds per cubic meter) the ARAC-simulated normalized air concentrations for uranyl fluoride integrated over the total time of passage of the cloud. Figures 4.1.2.4A through 4.1.2.6A present the simulated normalized cumulative deposition (in units of inverse square meter). Figures 4.1.2.7A through 4.1.2.12A represent the same simulations for HF.

Figures 4.1.2.13A through 4.1.2.20A present contours of the simulated instantaneous normalized concentration of hydrogen fluoride at 10 minute intervals; Figure 4.1.2.13A begins the sequence at 11:40 a.m. local time, 10 minutes after the release began. The sequence illustrates the progression of the cloud and its transport and diffusion downwind. Figures 4.1.2.21A through 4.1.2.27A are the corresponding series for uranyl fluoride until the wind has blown the airborne material off the grid.


In addition, LLNL-ARAC has provided numerical results at selected points. Points selected are residences and sample locations near residences in the cloud path, towns downwind, and an array of other points downwind to 72 km. Figures 4.1.2.28A and 4.1.2.29A show the location of these points. The numerical values of the integrated normalized air concentrations and normalized cumulative depositions are presented in Tables 4.1.2.1A and 4.1.2.2A. Table 4.1.2.3A presents normalized air concentrations the peak 10-minute exposure period at selected points in the affected area north of the Robert S. Kerr Reservoir.


4.1.2.2 Meteorology


Meteorological conditions in the vicinity of the site at the time of the accident, with strong, gusty winds generally from the north-northwest and clear skies, reflected the passage of a cold front through the area approximately 6 hours earlier. At the time of the accident, onsite meteorological measurements indicated winds flowing from the north-northwest (resulting in transport to the southsoutheast) at about 18 mph with gusts up to about 30 mph. Because of the relatively high wind speeds and somewhat unsteady wind direction, atmospheric stability was considered to be neutral (Pasquill type "D"), which is typical of post-cold-front situations. This is also supported by analysis of the temperature structure in the planetary boundary layer from the nearest atmospheric sounding station at Monett, Missouri. Temperatures in the area were in the upper 40's (°F) and relative humidity was approximately 40 percent.

There is no official National Weather Service station in the immediate plant vicinity. The nearest station having similar topographic and climatological characteristics as the plant site is at Fort Smith, Arkansas, approximately 70 km (40 miles) east-southeast. The Sequoyah Fuels Corporation (SFC) maintains continuously recording wind speed and direction instrumentation at the plant. Although the onsite meteorological measurement system is located atop a sign on the roof of the process building and would likely be affected by airflow over the building and sign, the measurements at the time of the accident are consistent with the synoptic weather pattern and other measurements and observations in the region and, therefore, have been used for the assessments of atmospheric transport, diffusion, and deposition. The general transport direction away from the site was south-southeast, which is consistent with site and regional meteorological information. Near the facility, plume transport probably was influenced by the presence of buildings and initial release characteristics such as buoyancy and other release energetics.


Figure 4.1.2.30A is the surface weather map for Oklahoma and adjacent regions for 12:00 noon (18Z), Saturday, January 4, 1986, (from the U.S. Department of Commerce, NOAA/NWS/NMC Washington). It shows, among other things, the wind direction in east-central Oklahoma, and the cold front which by then had moved out of Oklahoma. Figure 4.1.2.31A is the 850 millibar weather map of the area at 6:00 a.m. (12Z), Saturday, January 4, 1986.







.






Figure 4.1.2.3A Integrated normalized air concentration of UO_2F_2




Figure 4.1.2.4A Normalized cumulative deposition of UO_2F_2








Figure 4.1.2.5A Normalized cumulative deposition of UO2F2










1



-31

Figure 4.1.2.8A Integrated normalized air concentration of MF




Figure 4.1.2.9A Integrated normalized air concentration of HF

1




Figure 4.1.2.10A Normalized cumulative deposition of HF




Figure 4.1.2.11A Normalized cumulative deposition of HF





Figure 4.1.2.12A Normalized cumulative deposition of HF



1.8

SEQUOYAH FUELS CORP. RELEASE

1

Figure 4.1.2.13A Instantaneous normalized air concentration of HF 10 minutes after cylinder rupture




Figure 4.1.2.14A Instantaneous normalized air concentration of HF 20 minutes after cylinder rupture




Figure 4.1.2.15A Instantaneous normalized air concentration of HF 30 minutes after cylinder rupture





Figure 4.1.2.16A

Instantaneous normalized air concentration of HF 40 minutes after cylinder rupture



-

SEQUOYAH FUELS CORP. RELEASE

-

l,

Figure 4.1.2.17A Instantaneous normalized air concentration of HF 50 minutes after cylinder rupture




Figure 4.1.2.18A Instantaneous normalized air concentration of HF 60 minutes after cylinder rupture




Figure 4.1.2.19A Instantaneous normalized air concentration of HF 70 minutes after cylinder rupture




Figure 4.1.2.20A Instantaneous normalized air concentration of HF 80 minutes after cylinder rupture




Figure 4.1.2.21A Instantaneous normalized air concentration of UO_2F_2 10 minutes after cylinder rupture




Figure 4.1.2.22A Instantaneous normalized air concentration of UO_2F_2 20 minutes after cylinder rupture




Figure 4.1.2.23A Instantaneous normalized air concentration of UO_2F_2 30 minutes after cylinder rupture




Figure 4.1.2.24A Instantaneous normalized air concentration of UO_2F_2 40 minutes after cylinder rupture




Figure 4.1.2.25A Instantaneous normalized air concentration of UO_2F_2 50 minutes after cylinder rupture






Figure 4.1.2.26A Instantaneous normalized air concentration of $\rm UO_2F_2$ 60 minutes after cylinder rupture




Figure 4.1.2.27A Instantaneous normalized air concentration of UO_2F_2 70 minutes after cylinder rupture




Figure 4.1.2.28A Locations, to 8 km downwind, of points for which numerical concentrations were calculated




Figure 4.1.2.29A Locations, to 72 km downwind, of points for which numerical concentrations were calculated

Table 4.1.2.1A Integrated normalized air concentration and deposition

The numerical values below represent normalized cumulative ground depositions and integrated normalized air concentration values. A 10 km by 10 km calculational grid was used. Seventy-five percent of the material was assumed to be released in the first 5 minutes with the rest being released in the following 40 minutes. Measures of wind values from the Sequoyah Fuels Plant were used which showed a 335 degree initial wind direction. The release point location is UTME - 311.25, UTMN - 3930.54.

Location					HF		2F2	
Lab	el		(Km)	UTME (km)	Int. Norm. (s/m ³)	Norm. Cum. (1/m ²)	Int. Norm. (s/m ³)	Norm. Cum (1/m ²)
1*	16		3928.4	311.0	0.0	0.0	0.0	0.0
2	17		3928.4	311.2	3.486E-9	0.0	2.410E-9	0.0
3	18		3928.4	311.4	9.386E-8	2.61E-10	1.017E-7	1.43E-10
4	19		3928.6	310.8	0.0	0.0	0.0	0.0
5	20		3928.6	311.7	2.725E-6	2.584E-8	2.704E-6	7.493E-9
6	21		3928.9	310.0	0.0	0.0	0.0	0.0
7	22		3929.1	311.1	0.0	0.0	0.0	0.0
8	23		3929.1	311.3	1.310E-7		1.388E-7	1.600E-9
9	24		3929.1	310.0	0.0	0.0	0.0	0.0
10	25	1	3929.3	311.1	0.0	0.0	0.0	0.0
11	26		3929.3	311.3	2.225E-7	3.258E-9	2.236E-7	1.483E-9
12	27		3929.4	310.6	0.0	0.0	0.0	0.0
13	28		3929.5	310.7	0.0	0.0	0.0	0.0
14	29		3929.5	311.0	0.0	0.0	0.0	0.0
15	30	-)	3929.6	310.8	0.0	0.0	0.0	0.0
16	P1		3928.6	311.3	4.247E-8	0.0	4.341E-8	0.0
17	P6		3927.8	311.9	8.534E-7	3.284E-9	8.597E-7	2.982E-10
18	P7		3928.4	312.5	1.146E-7	1.854E-10	1.144E-7	2.345E-10
19	P8		3927.9	312.6	2.937E-7	2.117E-9	9.112E-7	2.258E-9
20	P9		3928.2	312.9	5.492E-8	1.115E-12	7.418E-8	1.705E-10
21	P11		3928.7	313.5	0.0	0.0	0.0	0.0
22	P12		3928.9	313.8	0.0	0.0	0.0	0.0
23	P13		3928.3	313.5	2.893E-12	0.0	1.289E-11	0.0
24	P14		3928.1	313.4	2.894E-10	0.0	3.047E-10	0.0
25	P19		3926.7	315.6	0.0	0.0	0.0	0.0

Table 4.1.2.2A Integrated normalized air concentration and deposition

The numerical values represent normalized cumulative ground depositions and integrated normalized air concentration values. An 80 by 80 km grid was used. Seventy-five percent of the material was assumed to be released in the first 5 minutes with the rest being released in the following 40 minutes. Measures of wind values from the Sequoyah Fuels Plant were used which showed a 335 degree initial wind direction.

		Location		Н	F	UU	2F2
Labe	el	UTMN (km)	UTME (km)	Int. Norm. (s/m³)	Norm. Cum. (1/m²)	Int. Norm. (s/m³)	Norm. Cum (1/m ²)
26*	T.	3918.6	320.0	2.432E-7	1.297E-9	2.614E-7	2.093E-9
27	RH	3913.0	323.8	1.548E-7	4.795E-10	1.739E-7	5.165E-10
28	GA	3910.3	313.5	0.0	0.0	U.0	0.0
29	KE	3903.0	325.0	1.124E-7	1.032E-9	1.021E-7	1.114E-9
30	MC	3891.3	320.8	0.0	0.0	8.125E-13	0.0
31	IR	3901.0	322.1	4.485E-9	3.100E-11	8.220E-9	1.024E-10
32	BO	3895.3	337.3	3.760E-9	1.293E-10	3.785E-9	3.219E-11
33	PO	3880.0	352.3	2.380E-10	3.219E-13	1.230E-10	0.0
34	CO	3908.3	337.5	4.772E-13	0.0	3.218E-12	0.0
35	HE	3863.0	354.0	7.440E-10	0.0	1.277E-9	1.130E-11
36	WI	3870.0	344.0	1.376E-9	1.907E-12	3.665E-9	2.708E-12
37	FA	3868.0	332.0	9.963E-10	1.470E-13	3.680E-9	5.321E-11
38	SU	3864.4	335.0	4.343E-9	3.087E-12	1.226E-9	2.609E-11
39	\$7	3923.9	314.1	1.671E-07	9.173E-11	2.012E-07	1.195E-10
40	S12	3919.3	315.9	6.112E-08	5.281E-10	7.530E-08	4.661E-10
41	S24	3908.2	320.5	2.070E-08	7.619E-11	2.203E-08	2.894E-10
42	S40	3893.4	326.7	2.019E-08	1.098E-10	1.734E-08	1.389E-10
43	S56	3878.5	332.9	4.990E-09	0.0	6.488E-09	8.836E-12
44	S72	3863.6	339.0	1.418E-09	1.022E-11	7.378E-10	2.807E-12
45	E7	3925.4	316.4	1.128E-09	7.025E-11	5.076E-10	1.648E-12
46	E12	3921.9	319.9	1.353E-08	1.067E-10	1.248E-08	9.529E-11
47	E24	3913.5	328.3	7.421E-12	0.0	7.129E-12	0.0
48	E40	3902.1	339.7	2.172E-11	2.378E-13	5.823E-11	9.133E-13
49	E56	3890.7	351.1	8.950E-11	0.0	3.723E-13	0.0
50	E72	3879.3	362.5	1.746E-11	0.0	8.602E-12	0.0

đ

Location	Normalized Air Concentration Approximating Peak 10-Minute Period (1/m ³)			
	U02F2	HF		
2	1.9 × 10-12	5.3 × 10-12		
2 3 5 8	4.2 × 10-11	5.4 × 10-11		
5	1.5 × 10-9	1.6 x 10-9		
8	1.0×10^{-10}	1.3 × 10-10		
11	3.0 × 10-10	1.4 × 10-10		
16	1.8 × 10-11	4.3 × 10-11		
17	4.0 x 10-10	4.9 x 10-10		
18	1.0×10^{-10}	1.1 × 10-10		
19	5.0 x 10-10	7.0 x 10-10		
20	2.2×10^{-11}	2.6 x 10-11		
24		5.2 × 10-13		

Table 4.1.2.3A Normalized peak air concentrations

0

Locations 2 through 11 are sampling points. Locations 16, 17, 18, 19, 20 and 24 correspond respectively to residence locations 1, 6, 7, 8, 9 and 14 on Figure 4.1.3.1.

a b





Figure 4.1.2.30A The Surface Weather Map at 12:00 noon January 4, 1986



APPENDIX 4.1.4.2

SFC'S INVESTIGATION REPORT OF THE EFFECTS ON CATTLE AND ANIMALS IN THE POTENTIALLY AFFECTED AREAS AFTER THE ACCIDENT

SEQUOYAH FUELS CORPORATION

INTERNAL CORRESPONDENCE

TO	Joi	nn	Stauter
FROM	c.	L.	Couch

DATE January 31, 1986

SUBJECT Task Force (Interagency) Investigation: Data Items Inquiry for Health Impact Analysis

As you requested in your letter dated 1-27-86, the following responses to the Data Items Inquiry for Health Impact Analysis, Section E <u>Animal Population</u> were prepared. These questions were answered by Carol Couch, Environmental Engineer and Dr. John Miller DVM, Warner, Oklahoma.

Section E. Animal Population

- 1. Characterization of animal population at risk:
 - a. Description of sentinel herds
 - Species, number, age, diet, water source location on premises of fuel plant.
- RESPONSE (1) There were two herds of cattle (Bos taurus) that were exposed to the release on 1-4-86. Both herds belong to Kerr-McGee Corporation and were on facility owned lands. One herd consisted of 99 crossbred steers that were located southeast of the facility (Township 12N, Range 21E, parts of Sections 22 and 27). The other herd was 87 crossbred heifers located southsouthwest of the facility (Township 12N, Range 12E, parts of sections 21 and 28).

0

Both herds were 10 months to 1 year of age. They are allowed to graze on open pasture lands consisting of predominantly bermuda and fescue grasses. Their diets are supplemented with cattle cubes, hay, salt and mineral. They water from stock ponds located on the pasture areas.

Sequoyah Fuels produces an ammonium-nitrate fertilizer solution as a by-product of the process. During the summer growing season this product is applied to company owned lands under the direction of Oklahoma State University Consulting Agronomist, Dr. Billy Tucker. It is very closely regulated and monitored in accordance with NRC guidelines.

Based upon years of research and data collection the NRC has permitted Kerr-McGee unrestricted use of the lands for cattle grazing.

The hay the animals are fed is grown on the pastures during the summer. These pastures are part of the facility's annual fertilizer program. The hay was analyzed for As, B, Co, Cu, Fe, Mn, Mo, Ni, Pb, V, Zn, U, Th-230, Ra-226, and nitrates prior to its release for animal consumption. The hay must meet guidelines set by the NRC for these element concentrations or it is not released for use.

> Stock ponds in these pasture areas are monitored for N, As, Cu, Mo, Pb, Se, gross alpha, Ra-226s, Ra-226i, and U-238 as has been directed by the NRC to meet compliance of the annual fertilizer program.

- QUESTION (2) Is there a regular health monitoring program for these animals that is supervised by a veterinarian? If so, what information is available regarding baseline health status?
- RESPONSE (2) No. The ranch manager takes care of any immunization and regular herd health maintenance. If something should develop that he does not feel he can adequately handle a veterinarian is called.

There is not any baseline data on the herd. However they were received at the facility in good health and are growing and gaining weight very well.

- QUESTION (3) Are blood, urine, or tissue specimens collected and analyzed on a regular basis? If so, what are baseline fluoride levels in the animals?
- RESPONSE (3) No. The facility has very little if any effect on the surrounding environment. All facility discharges are well below standards set by the State, EPA, and NRC. There has not been any need prior, to the

Î

QUESTION (4) Current status of the herds.

RESPONSE (4) Dr. John Miller DVM, has monitored the herds since 1-6-86. He feels that they are in good general physical health.

> Urine samples were taken from 56 steers on 1/6-7/86. These samples were analyzed for uranium and fluoride. Analysis results are given in Item 1. Dr. Miller checked temperatures, eyes, nose, and mouth areas looking for any evidence of HF acid irritation, or toxin symptoms. Temperatures were normal. Their eyes did not show anything that would be acid related. There was not any excess mucus or lesions of the nose or mouth which would indicate a toxin. Their general health looked good.

The steers were held in a holding pen out of the plume area to prohibit any further contamination. However, two steers were lost due to the close confinement. Dr. Miller conducted necropsys on the two animals and found them to be stressed to death. The other steers "yead" them to death. This is common when strange cattle are introduced to a herd. (See Item 7)

The heifers were moved to an area well to the west of the plume path. Urine samples were taken

2

> from 12 of these animals on the 1-7-86. Dr. Miller examined the herd and found them to have normal temperature, no eye irritation or excessive mucus of the nose or mouth.

Heifer No. 320 was found dead on 1-15-86. Dr. Miller conducted a necropsy on the animal and bloat determined the cause of death to be bloat. The heifer's diet had been changed cogising her to bloat bloat

Analysis of the heifers are attached in Item 2.

There was one horse in the same area as the steers. This animal was checked on 1-7-86. His physical condition was fine. His lungs were clear, heart rate and temperature was normal, there was not any signs of irritation of the eyes or exposed skin. A urine sample was taken. (Results are given in Item 3.) Further testing of this animal was not warranted.

The ll heifers are being resampled weekly. Dr. Miller is checking the animals on a weekly basis. Twenty of the steers were selected for weekly testing. Dr. Miller is also keeping a watch on these animals. Testing of the urine for uranium and fluoride will continue until the level. have

ľ

dropped to those of area control cattle. The remaining animals have been returned to their respective pastures. The test animals are being held in separate pastures to allow ease of testing.

The ranch manager does a daily check on all the animals and does a count to assure all the animals are accounted for and are in good health.

Dr. John Miller feels the herd are in good general health. They have some very common herd problems such as ringworms, one case of rot foot, lice rub, watering of the eyes from hay and dust. However, none of these problems are abnormal to any herd in the area. Dr. Miller has not found any illnesses that have been a result of the release. His comments are given in Item 5. Dr. Billy Tucker, OSU consultant, evaluated the urine data, his comments are attached. (See Item 6)

1

TOPIC (b) Description of domestic animal population in area surrounding the fuel plant:

QUESTION 1) Species, number, location of livestock, particularly range or pasture animals.

RESPONSE (1) The response to the inquiry is only an estimation. To be able to be specific a survey would need to be conducted.

1

Typical domestic animals in the area would be as follows: dogs, Canis familearis; cats, Felis catus; horses, Equus przewalski; cattle, Bos taccus. There could be some swine (Sus) and possibly some goats (Capra aegagrus) in the area.

After the release on 1-4-86 herds of animals south of the facility, that could have had exposure to the release were surveyed. Also as persons, even those well out of the plume area, requested visits were made to their farms to check their animals. Item 4 is a summary of these off-site visits and the results of this investigation.

TOPIC (c) Description of wild animal population in area:

- 1) Species, number, location.
- Fish ponds in area or other sport or commercial fishing.
- RESPONSE This is not a question that personnel at the facility or Dr. John Miller are qualified to answer. The State Fish and Wildlife Department should have this information.

TOPIC (2) Estimation of dose to animals:

QUESTION a) Have the sentinel herds been monitored since the incident for signs of illness? If so, describe illnesses.

- RESPONSE (a) Yes. The ranch manager checks the herds daily and Dr. Miller checks the animals weekly. There have been no signs of illnesses that could or are related to the release.
- QUESTION (b) Have any blood, urine, or tissue samples been collected from the sentinel herds since the incident? If so, what tissues were collected, what are the results of any analysis that have been done, and are any of these tissues still available? Slaughter and necropsy of sentinel herd animals would not be appropriate.
- RESPONSE (b) Urine samples on 20 steers and 10 heifers are being collected weekly and this will continue until results fall within the range of control animals (See Item 4). The urine is being analyzed for uranium and fluoride.

No tissue was or is being collected.

QUESTION (c) What are the results of the survey of local veterinarians to determine if there have been any illnesses in livestock or pot animals that could be due to plume exposure? Provide a description of species affected illnesses seen, and location of animals in relation to fuel plant. Provide the results of any milk samples collected. Provide the results of any analysis of any milk samples collected.

RESPONSE (c)

Dr. John Miller has talked with Dr. Cindy Carter DVM, Gore, and she has not had any reported illnesses. Dr. Miller also visited with Dr. Gary Cox, DVM, Sallisaw. Dr. Cox had looked at some of resident 5R's swine and told him it was his opinion that the animals had swine dysentery and parasites. Dr. Cox did not feel the swine were ill as a result of the release. Dr. Cox did not report any other cases or illnesses.

Dr. Miller has responded to all the suspect cases that have been brought to the attention of the company. Dr. Miller has not found any illnesses that were a result of plume exposure.

The following types of domestic animals have been observed or examined and none have shown any illnesses as a result of the release: Cattle (steers, heifers, cows, bulls), horses (gelding and mares) goats, swine, ducks, geese, rabbits, dogs, and cats. The location of these animals ranged from the facility sight to 15 miles southeast of the facility.

0

No milk samples have been collected.

QUESTION (d) Has the State Fish and Wildlife Service surveyed wild animal populations near the plant to determine if there have been increases in illnesses or deaths? Have there been reports of fish die-off in the reservoir or the rivers near the plant?

RESPONSE (d) This question should be addressed to the appropriate State agencies. But, to our knowledge the answers to the above questions are no.

QUESTION (e) What steps have been enacted to control the movement of livestock and poultry in and out of or within (quarantine) the contaminated area.

RESPONSE (e) There is not nor is there any need for a quarantine of the area. The orly animals that have shown any contamination are the animals on the plant site that belong to Kerr-McGee and they are being monitored. They are showing no signs of illness. Off-site investigation of animals has not indicated any contamination or any plume related illnesses.

CLC:dd

88

Attachments

Kerr-McGee

Steer Urine Analyses 1-6-86

Sample No.	F mg/l	U pg/1
4	12	200
39	8.5	110
49	8.5 5.9 5.8 8.2 16	87
87	5.8	1100
35	8.2	2400
43	16	150
18	10	1400
95	13	410
13	3.6	83
5	11	6400
90	4.2	350
82	9.7	1000
68	11	350
76	8.2	830
38	14	740
75	13 3.6 11 4.2 9.7 11 8.2 14 13 6.9	160 750
11	6.9	750
46	6.9 7.5 11 4.1 5.6 10 5.9 7.5 8.1 6.4 2.8 4.8 9.8 2.7 11 21 7.5	590
71	7.5	130
36	11	290
100	4.1	150
84	5.6	1100
17	10	1400
30	5.9	220
56	7.5	220
61	8.1	390
34	6.4	210
42	2.8	75
44	4.8	39
44 15	9.8	330
77	2.7	120
77 73	11	730
85	21	580
58	7.5	160
86	2.5	710
48	9.1	190
96	2.5 9.1 3.1	37
55	5.7	110
20	6 5	220
20	6.5	220

*63	composites 3 samples	of	7.2	640
*24	composites 4 samples	of	7.3	130
*60	composites 6 samples	of	4.7	510
*88	composites 6 samples	of	3.9	540
КM	horse Ker	r-McGee se	3.2	22

ľ

Kerr-McGee

Heifer Urine Analyses

Sam	ple No.	F mg/l	U ug/1
369	yellow	4.2	770
371	yellow	1.4	1800
356	white	5.5	320
449	white	6.2	63
574		14	590
404		3.7	450
473		19	2200
74	white	16	7200
369	white	15	840
320	(deceased)	2.6	570
293		20	2700

Kerr-McGee Horse Urine Analyses

Sample No.	F mg/1	U ug/1
K-M horse	3.2	22

١

1

I

OFF SITE INVESTIGATION

INDEX:

1

1

Section	1:	Investigation of Animals South of
		Facility in Plume Path
Section	2:	Requested Investigation of Animals
Section	3:	Control Head Investigation

۱

1

SECTION 1

.

RESIDENCE: IR

DATE: January 7, 1986

LOCATION: Township 12N, Range 21E, Section 26

HERD TYPE: Cross bred and hereford cattle

HERD SIZE: Approximately 43 head

Dr. John Miller and I visited. ^{1R} and took urine samples from 4 of 1R's calves. 1R's herd size was approximately 43 head of herfords and crossbred cattle. Dr. Miller also examined two of 1R's young dogs. The dogs were in good health with no signs of illness or evidence of HF acid burn or toxin. The cattle were in good general health with signs common herd parasites. There was not any signs of acid or uranium related illnesses. Urine analysis results:

Sample No.	F mg/1	U mg/l
1 calf	1.3	Insufficient sample
2 calf	1.2	2
3 calf	0.9	8
4 calf	1.0	48

Urine analyses were normal and below those of the control group. Further testing was not warrented.

RESIDENCE:	2R (lease land, no house)
DATE:	January 8, 1986
LOCATION:	Township 12N, Range 21E, Section 27
HERD TYPE:	Horses (mares)
HERD SIZE:	3

Sample No.	F mg/l	U ug/1
Pk mare	0.3	2

CONCLUSION: General physical health was fine. Heart rate was normal, lungs were clear, temperature and urine pH was normal, no eye irratations or excessive mucus from the nose or mouth. Urine did not contain significant amounts of uranium or fluorides. Further testing was not necessary.

RESIDENCE: 3R

DATE:	January 7, 1986
LOCATION:	Township 12N, Range 21E, Section 27
HERD TYPE:	Cross bred cattle, horses
HERD SIZE:	Approximately 21 head cattle, 2 horses

Sample No.	F mg/1	U µg/1
17 cow	1.2	3
2 cow	1.2	4
10 cow	1.1	3
3 cow	0.6	2
80C cow	0.8	2
84C cow	9.7	2
Black Welch gelding	0.2	2

CONCLUSION: Dr. Miller examined the animals and found the temperatures were normal, there was not any signs of acid related skin or eye irritation, no excessive mucus which would indicate toxin. Dr. Miller's exam found the herd in good general physical condition. The horse had a normal pulse rate, clear lungs, no indication of acid burn to the eyes or exposed skin, no signs of excessive mucus. The horse was also in good general physical condition. No significant uranium or fluorine in the urine. Follow-up sampling was not warranted.

RESIDENCE: 4R

DATE: January 7, 1986


LOCATION: Township 12N, Range 21E, Section 26

HERD TYPE: Hereford cattle

HERD SIZE: Approximately 10

Sample No.	F mg/1	U µg/1
1 cow	2.8	2
2 cow	1.3	2
3 calf (6 wks. old)	0.6	4

CONCLUSION: Dr. Miller examined the animals and did not find any HF acid or uranium suspect illnesses. Eye, mouth, nose areas appeared normal. Temperatures were normal. Dr. Miller did not find anything except normal herd parasites. The urine contained no significant uranium or fluoride, further testing was not warranted.



RESIDENCE: 5R

DATE: January 15, 1986

LOCATION: Township 11N, Range 23E, Section 29 HERD TYPE: Swine, cattle, ducks, geese, dog, cat, rabbits HERD SIZE: (All estimates) - 25 swine, 5 cattle, 6 ducks, 4 geese, 1 dog, 1 cat, 6 rabbits

5R had lost 7 head of swine, he felt his animals were all ill due to the plant release. He also felt his illnesses were a result of the acid release.

No urine samples were taken at the time of this visit. Observation of the animals were as follows.

SWINE: The swine were of varying breeds and ages. The younger pigs were sick with diarhea and defecating at random. They were let run out in open dusty areas with little or no ground cover. The only feed available to the younger pigs was watered bread. The swine and cattle watered from a pond inhabited with water fowl (duck and geese). The fowl exhibited no signs of illness.

CATTLE: The cattle were also on the same barren pasture. They had been given loose hay. The cattle had some eye tearing which Dr. Miller felt was caused by the dust and hay. There was not excessive reddening of the eye area or skin irritations. The cattle did not show signs of excess mucus.

CAT:

DOG:

K,

The cat appeared in good general health. No symptoms of any illnesses.

The dog had a case of diarrhea. However in talking to 5R he explained that he had ran over the dog with a backhoe. The dog had lost control of his tail and the tip had been severed. 5R also stated that the animal had had severe swelling in the backend. Other than these problems and the fact that the dog appeared underweight he looked fine. There was not any eye irritation or excessive mucus from the nose or mouth.

RABBITS: The rabbits were located next to the pasture area the swine and cattle were located in. The rabbits did not have any eye irritation, or any signs of acid burn to the ears. They appeared in good general health. They looked under weight. They did not have food or water available to them. The rabbits were observed eating fur off of one another from adjoining cages. bR fed and watered the animals while we were there and the animals were very hungry and thirsty. 5R claims the animals were loosing weight since the release. Observation of the rabbits lends the impression that the rabbits were not being fed properly.

Dr. Miller suggested we purchase 3 swine and send two to OSU, Oklahoma Animal Diesease Diagnostic Laboratory. Three pigs were purchased and 5R gave the fourth pig, which was dying, for diagnoses. A full report from OSU is attached. See Item 9. Dr. Miller examined the other 2 pigs and found no evidence of HF or uranium related problems. Dr. Miller's findings were the same as those of OSU.

CONCLUSION: 5R's property was not in the plume path. Dr. Miller did not find any illnesses that were related to the facility release.

> 5R claimed his eyes were burned the day of the release and that he had holes burnt in his jean legs. However, his eyeglasses were not etched. The glass in the cap of his backhoe and vehicles showed no evidence of etching. Etching of glass should be quite evident and is characteristic of HF. None was observed.

Also, the rabbits should have shown some evidence of "burn" since they have extremely sensitive eyes and ears, but they did not.

In summary ^{5R} had several noted items indicating poor animal husbandry, infections, and parasites in his herds.

RESIDENCE: 6R

l

DATE: January 27, 1986

LOCATION: Township 11N, Range 22E, Section 3 HERD TYPE: Goat (capra aegagrus) HERD SIZE: 2 females

^{6R} called the facility 1-27-86. She had had a pregnant goat that had died the night before. She wanted someone to check the animal and see if its death was related to the plant release.

Dr. Miller conducted a necropsy on the animal. The body cavity of the animal was packed with fat. The liver was discolored to a yellowish brown and there was not any blood present in the liver. The heart showed signs of dead muscle. There was a layer of fat around the outside of the atrium areas. The goat was pregnant with 3 fetus and was approximately 1 - 1 1/2 months from delivery.

The stomach and intestine area was filled with parasitic "barber pole" worms.

Dr. Miller determined the cause of death to be pregnancy toxemia. This is the reason for the "cooked" liver. The goat had so much fat that she could not get proper circulation to allow her vital organs to function properly.

^{6R} also had another pregnant femal. Dr. Miller advised her to worm the animal and see that it was properly exercised each day.

6R witness the necropsy and was completely satisfied that the goat died of a physiological problem and not as a result of anything from the facility.

RESIDENCE: 7R

DATE: January 8, 1986

LOCATION: Township 12N, Range 21E, Section 36

HERD TYPE: Crossbred cattle, horse, cat

HERD SIZE: Approximately 15 cattle, 1 horse, 1 cat

Sample No.	F mg/1	<u>U ug/1</u>
609 cow	1.3	2
610 cow	1.8	3
611 cow	1.5	2
612 calf	0.7	2
P mare	0.8	2

CONCLUSION: Dr. Miller examined the 7R's cat and didn't find any release related symptoms. The cat did have a case of ear mites. The palomino mare was in good general physical health without any illness.

The sattle had normal temperatures. Eyes, nose, mouth and skin areas were clear and without irritation.

Dr. Miller examined one cow that had a conjested right lung. Dr. Miller did not think that this was related to the release or that it was a problem to the cow's health.

Sample results did not warrant any further testing.

The herds' general condition was good. However, 7R requested some additional urine testing and this was conducted on 1-16-86. Dr. Miller examined some additional cattle in 7R's herd and found them to be in good health. Some of the cows had a common uterin virus but there was not any illnesses that were related to the facility release. Analyses on five head of cattle sampled on 1-16-86 are still pending.

SECTION 3

CONTROL ANIMALS

RESIDENCE: 8R

۱

۱

į.

DATE: January 9, 1986

LOCATION: Rabbit Hill property Township 12N, Range 19E, Section 25

HERD TYPE: Cross bred steers and Holstein Heifer HERD SIZE: 12

Sample No.	F mg/1	U_ug/l
614	2.0	98
615	2.0	6
616	6.2	30
617	1.4	5
Holst. Heifer	2.8	22

CONTROL ANIMALS

RESIDENCE: ^{8R} DATE: January 9, 1986 LOCATION: Webbers Falls, Township 11N, Range 21E, Section 9 HERD TYPE: Cross bred cattle HERD SIZE: 12

Sample No.	F mg/l	<u>U µg/1</u>
Red cow	1.4	16
Bobtail	1.7	4
883	1.6	7
Blk. Pntr.	2.9	22
BWF COW	0.5	13

CONTROL ANIMALS

0

RESIDENC	E: K-N	M Corporation
DAT	E: Jan	nuary 9, 1986
LOCATIO	N: To	ownship 12N, Range 12E, Section 22
HERD TYP	E: C	rossbred steers
HERD SIZ	:E:	

Ē

Sample No.	F mg/1	U ng/1
14CG	1.6	47
27CG	1.2	13
34CG	0.8	12
160CG	2.0	9
154CG	1.5	23

CONTROL ANIMALS

RESIDENCE:	9R
DATE :	January 9, 1986
LOCATION:	Township 13N, Range 22E, Section 15.
HERD TYPE:	Limousin cattle, horses
HERD SIZE:	34 head of cattle, 3 horses

Sample No.	F mg/1	U ug/1
7961CG	0.6	< 2
Jenny	1.3	< 2
WCT8751	0.6	< 2
73BH27956	1.0	8
BH27961	0.5	10
CCCG	1.9	4
Poncho (horse)	1	< 2

ITEM 5

JOHN W. MILLER, D. V. M. TOWN & GOUNTRY VETERMARY OLINE BOR 447 WARNER, OKLA. 74460 TELEPHONE 018 - 403-2007

January 10, 1986

Ms. Carol Couch Kerr-McGee, Inc. Gore, Oklahoma 74435

Dear Ms. Couch:

As per your request, here are my observations on the conditions of the domestic animals I examined earlier this week.

Steers belonging to Kerr-McGee examined on 1-6-86 were in good general physical condition. Body temperatures and lung sounds of these examined were in normal ranges. Mucus membrane's of the eves, lips, gums, and nasal passages appeared normal on 30 steers. We did notice several cases of ring worm and skinned places on rear legs.

On Tuesday, 1-7-86, I examined about 190 head. They all exhibited symptoms seen regularly in a large animal practice. The black gelding on Kerr-McGee property revealed normal vital signs and urine pE. I observed no lesions to the mucus membranes.

On 1-8-86 the 7 year old Appaloosa mare examined appeared normal. No lesions of the eyes or mouth were noticed. Urine was collected and revealed a pH of 7. Two other mares in the pasture appeared in good general physical condition. The cattle we collected urine samples from exhibited normal vital signs with the exception of one black cow who exhibited congestion in the right lung. Her mucus membranes were normal as was her temperature. The palomino mare we collected urine from also appeared normal.

On 1-9-86 we examined and collected urine samples at four locations. All animals seen that day appeared normal.

Of the approximately 300 animals I have observed, examined, and/ or collected urine from, I have not observed a syndrome that I would consider out of the ordinary for a large animal practice in this area.

Sincerely,

John w. miller Jum.

John W. Miller. D.V.M.

ITEM 6

1212 N. Jardot Stillwater, Oklahoma 74074 January 27, 1986

Ms. Carol Couch Sequoyah Fuels Incorporated P. O. Box 610 Gore, OK 74435

Dear Carol:

I have reviewed the analyses from urine samples taken from the livestock owned by Sequoyah Fuels.

Even though concentrations of uranium and fluorine in the urine were higher in some of the animals than those generally found from other animals tested, these concentrations are still below suggested critical levels. Levels found in these animals pose no health problems to the animals and there is no risk of elevated levels being transmitted to the human chain from these animals or their products.


Sincerely yours,


Billy B. Tocker

Billy B. Tucker Agronomist and Soil Scientist

BBT jt

Enclosure





KEEP THIS SLIP FOR REFERENCE 5H 527 REDIPORM

ITEM 9

90

Oklahoma State University Oklahoma Animal Disease Diagnostic Laboratory

Case No: 112782

۱

1

5R's pig



Oklahoma State University

STILLWATER, OKLAHOMA 74078 (405) 624-6623

OKLAHOMA ANIMAL DISEASE DIAGNOSTIC LABORATORY

January 22, 1986

Mr. Butch Garver P.O. Box 393 Gore, Oklahoma 74435

Case No.: 112782 Owner: 5R

Dear Mr. Gar. T:

HISTOPATHOLOGIC EXAMINATION: Fig "A"- Colonic lesions consist of moderate to severe diffuse areas of superficial mucosal necrosis with associated accumulation of necrotic cellular debris and the presence of numerous bacterial colonies.

Mesenteric lymph nodes reveal moderate to heavy infiltrates of neutrophils, macrophages and occasional eosinophils within the paracortical areas with associated lymphoid hyperplasia of the cortical germinal centers.

There are no significant lesions in heart, lung, liver, spleen, kidney, stomach and intestines.

Pig "B" - Pulmonary findings include focally extensive infiltrates of neutrophils within significant portions of alveoli and air passages with the presence of several peribronchial and perivascular lymphoid infiltrates.

Colonic lesions reveal the presence of numerous trichura (whipworms) embedded within the superficial mucopa and extension into crypts.

There are no significant lesions in heart, liver, spleen, kidney, stomach and intestines.

Pig "C" - Colonic lesions reveal heavy numbers of <u>Trichura</u> <u>sp</u> (whipworus) embedded within the mucosa.

HISTOPATHOLOGIC CONCLUSION: Pig "A" - Severe necrotizing colitis; moderate suppurative lymphadenitis. Pig "B" - Moderate suppurative bronchopneumonia, severe trichuriasis. Pig "C" - severe trichuraisis.

Sincerely,

DVM. Patho logist





DIAGNOSTIC LABUMATORY α Cane No. ed cape (e.g. petred service) has b I a read seri is no apo 100 2000 OKLAHOMA STATE UNEVERSITY Schwater, Oldehome 74074 JAN 1 6 1986 Telephone: 405-624-8623 Condition of Specimen: C Excellent OKLAHOMA ANIMAL DISEASE Autonzeci C Other (Specify) THORAS NO SHO SABORN DEN PLEASE FILL OUT COMPLETELY WITH BLACK IN OX HILLER 5R VETERINARIAN, OWNER ADORERS ADDRESS_ 10:1 Der 2 DK 7:100 ZIP 773 8704 PHONE PHONE _ 10 - 86 Location of Case (County & States Ale DATE SUBLATTED Other (Specity) Lorries By Kige / Owner ... Mail Delivered by: Veterinarian. HISTORY sitting AGE 2 2 00 00 - SEX_ BREED SPECIES ____ WEIGHT 25 2 NUMBER DEAD NUMBER IN HERD FLOCK NUMBER SICK (excluding dead) ____ V RAISED ON OWNER'S PREMISES? YES NO. IF PURCHASED WHEN? ANY RECENT INTRODUCTIONS TO HERD (FLOCK)? YES . DATE INTRODUCED_ NÖ FIRST NOTICED SICK 1-14-86 CLINICAL SIGNS (INCLUDE TEMPERATURES) & POSTMORTEM FINDINGS en Actee regnozio tun trans 474429 418-6095527 VACCINATION HISTORY (DATES apress a perme 1M FIELD DIAGNOSIS MATERIALS SUBMITTED HOW MANY - REFRIGERATED FROZEN - RECEPTAN BY **EXAMINATIONS REQUESTED** Necropsy D Live Animal -Hetopethology" Dead Animal Fetus Cytology ġ Blood (Cloned) Mycoplasma (solation Blood (Anticoag) D-Bactenokogy Serum Anticiotic Sensitivity Kigney Virus Isolation & Identification** S UVOI 12 Parasitology floatation Thestine 1 P Mycology Lung n 0 Ventral Hernestology CBC Turnor Show boaton on model) -0 Puorescent Antibody Spleen 0 5 0 (Specify test ____ Brain 0 Serology (Complete form, C Lymph Node ____ 0 0 0 (ongin) 000000 (ada) Eyebal 0 Taxicology** MR 0 0 0 (Specily teel .. Urine 0 0 Cerebrospinel Fuid Win Path- Ca. Phos 0 0 0 Stomach Contents 5 1 -D Sweb __ 0 0 0 ma (ongen) Dorsal C Peed 0 0 0 Weter 0 0 0 0 Other n n 5 O Other . 0 0 Prozen Seaues not suitable. (Please record additional information on reverse.) ** Formatin Baad teause not suitable. Roused April, 1879

118

112782

SEROLOGY SAME

HA Liza to .

SAMPLE I

...

	Animal Identification					Lab Use	Leotospira (includes pomone, ictero, prippo, hardio,
ube	Name/No.	Species	Breed	Sex	Age	Only	Cantoola) (\$4.00) Anaplasmoete (\$4.00) Neonatal Isoerythrotysia (Complement-Lysis lest
-							requires neonatal blood in anticosquiar.) (\$4.00) Brucette abontus (Not official regulatory test except (evens) (\$4.00) Gen a4.005
-					+		500 (84.00) 500 (84.00) 500 (84.00)
-					-		Pis (\$4.00) Respiratory Syncytel Virus (\$4.00) Bovine Laukamie Virus (\$4.00)
					-		Chiemydia (54.00) Cattle Abortion Profile (Lapto, ISR, BVD, Bruceria) (59.00)
					-		Cattle Respiratory Disease Profile (IBR, Pts, BVD) (\$7 Swine Abortion Profile (Lepto, Brucella, Parvovirus, Pseudorables) (\$9.00)
							 Horse Abortion Profile (Lepto, Rhinopneumonitis, EIA-send EIA form) (\$9 00) Sheep Abortion Profile (Lepto, Brucella, Bluelongue, Chlamydia) (\$9 00) Swine Influenza Virus (\$4 00)
							Swine Parvovirus (\$4.00) Pseudorables (\$4.00) TGE (\$4.00)
	- C						Bruceila canis (\$4.00) Canine Distemper (Serum or CSF) (\$4.00) Feiline Leukemia (need unfixed blood smears) (\$4.00 El'A (\$5.00)
nt B-C	es for serology should b tubes are satisfactory contain toxic residues						EiA (\$5.00) Equine Ahinopheumonitis (\$4.00) Equine Influenza 1 (\$4.00) Equine Influenza 2 (\$4.00)

Laboratory charges are doubled for cases not owned by Cklahoma citizens.

	ADDITIONAL REMARKS		
4-5 fend ing	A		
		1	1
- Lance Estan	C.K. Arean	Jeraflor	le_
Joor doing	yuy?		

	and the second sec		
	and the second s	antenaariet ja 15 milie – 1	STAR 13
and the second se		1	



Oklahoma State University

STILLWATER, OKLAHOMA 74078 (405) 624-6623

OKLAHOMA ANIMAL DISEASE DIAGNOSTIC LABORATORY

January 22, 1986

Butch Garver P.O. Box 393 Gore, Oklahoma 74435 Case No.: 112782 Cwner: 5R Referral Lab.: Boren Veterinary Medical Teaching Hospital Oklahoma State University Stillwater, Oklahoma 74076

Dear Mr. Garver:

Attached hereto are the results of tests performed by an outside laboratory to which materials were referred from the above referenced case.

Sincerely,

Les Ray W. Ety, DVM, PhD Fathelegist

RWE/cas



CHUNH! NOTA NOTA C REMARKS Astron Dr. Miller V 1989 17 PALT TOUT UPwe -Bread/Bo JApa/Bez 122 M secting Hospits SwineB ANTUASE NA -K 10 PHOS 0 NOTES -T.P. BUN 10.7 1 CA -BOH CHOL BOOT a Boren Vetertnary Mer Oklahoma State Univ CPK BOPT CX348 BOOT CREAT DIGOXUN 00. LOH GLUC LLDH Tech: BE Total Chg. \$ 13.00 Date Done |-17-84 Reviewed By CHEMISTRY I 112782 TY 7 -RUE/OF STAT CHART COPY DATE PRICRITY STUDENT PLEASE CALL RESULTS 67 miller AT 171 REMARKS MADD Accress AC DEP RESULT TEST V RESULT erees St Ace Sex SIMAOL TEST V LIPASE ALB 4 MG Boren Veterinary Medical Tesching Hospilal Oklahoma State University A/G 127 NA AMYLASE ĸ BILL T 7.5 NOTES. PHOS BILL D. T.P BUN SAP CA SOH CHOL SGGT CL SGPT CPK SGOT CKMB DIGOXIN CREAT LOH co. LLOH Total Chg. 5 /3 CC GLUC CHEMISTRY I TAPA RE 0--- --- 1-17-Ri-Reviewed By MCac. C DIFF D BOUTTHE -Owec -T EMERGE CREC -IR 13 warden al OAL C REPEAT 0.900 5 RUEA anne 34 35 0 . -1 2 3 10 0 0 0 0 0 121

OKLAHOMA STATE UNIVERSITY	OKLAHOMA ANIMAL DISEASE DIAGNOSTIC LABORATORY STILLWATER, OKLAHOMA 74078 • (406) 624-6680
BACTERIOL	OGY REPORT
Preliminary Report	E Final Report
Veterinerian Dr. J. W. Miller	Owner_ 5R
	Received -16-16 Case No. 112712
Specimen(s) Pin A Isolate(s)	Specimen(s) Isolates(s
(Lung	B Ra B - (tive pig relemented)
Z Liver)	B Liver 2
Spleen } Now	E Silven (None
<pre>Spleen j Skidney</pre>	5 Small internet former Fare
S Lymph node - Not	
B Intestine (mail)	
D Brott Clan Normal Form	
Uterus	E Pac (+55-00 submitted)
Stomach contents	B Liver anteminents
Urine	E Intertine - Normal Flora
C Feces	
Muscle	0
Skin	
D'Turbinote such . Contaminents	
The bacteria isolated were The bacteria isolated were for antibiogram(s). Isolate referred to another An addendum report will for	
Intestinial samples men cultured	for Silmmella sop. and Trepresen
I intrine with negative near	4 * ·
hyperpenteries with negative real	
Date _ 1- 23-86	Bmonto



OKLAHOMA ANIMAL DISEASE DIAGNOSTIC LABORATORY

STILL WATER, OKLAHOMA 74078 . (408) 824-0623

Accession No. 86010563 Final x Preliminary Addendum		TOXICOLOGY REPORT	Case No. 112782 Related Cases Date Received 1-16-86		
Veterinarian	J.W. Miller, DVM	Owner	5R / Kerr-McGee		
SPEC	IMEN	ANALYTICAL METHODS	CONDITION		
x Liver x Kidney Brain Lung Fat Stomach contents Feces Eyeball	Serum CSF Milk Urine Feed	 Atomic Absorption Gas Chromatography E.C. Gas Chromatography Flame Gas Chromatography N-P Scanning UV Hp Liquid Chromatography X Ion Specific Electrode TLC Bench Chemistry Infrared Spectrophotometry Other 	 Sample Satisfactory Sample Unsatisfactory Sample Depleted Sample Referred Sample Disposed of 3 ample Stored* (one month) *Samples will be stored one month and will be discarded unless otherwise instructed. If conditions warrant, the submittor should collect and keep duplicate samples. 		

SPECIES: Swine

RESULTS AND COMMENTS:

Liver tissues were analyzed by ion specific electrode for fluoride.

liver A <1 ppm fluoride liver B <1 ppm fluoride

liver and kidney were analy:ad by atomic absorption spectrophotometry. Results are in ppm wet weight.

	Liver	Kidney
Chromiun	0.632	0.432
Cadmium	0.041	0.089
Lead	0.208	0.325
Arsenic	0.113	0.135

These levels are not clinically significant.

*Toxicology Charges ONLY \$ 60.00

Date Completed 1-22-86

*Please wait for statement before paying.

WCE/cas

William C. Edwards, DVM



Oklahoma State University

OKLAHOMA ANIMAL DISEASE DIAGNOSTIC LABORATORY

STILLWATER, OKLAHOMA 74078 (405) 624-6623

January 22, 1986

Butch Garver P.O. Box 393 Gore, Oklahoma 74435 Case No : 112782 Owner: 5R Referral Lab.: Boren Veterinary Medical Teaching Hospital Oklahoma State University Stillwater, Oklahoma 74178

Dear Mr. Garver:


Attached hereto are the results of tests performed by an outside laboratory to which materials were referred from the above referenced case.

Sincerely,

Ray W. Ely, DVM, PhD Pathologist

RWE/cas





TA Banker Stotesme to. statut ft.



Oklahoma State University

OKLAHOMA ANIMAL DISEASE DIAGNOSTIC LABORATORY

571LLWATER, CKLAHOMA 74078 (405) 624-6623

January 20, 1986

Mr. Butch Garver P. O. Box 393 Gore, Oklahoma 74435

Case No.: 112782 Owner: 5R

Dear Mr. Garver:

Two pigs are presented to necropsy with clinical history of several in the herr net doing well and in increased mertality rate.

POSTMORTEM EXAMINATION: The small white pig (A) appeared extremely emaciated with approximately 2-3 inches of the rectum prolapsed through the anus. Internal gross lesions revealed extensive patchy foci of necrosis involving the rectal mucosa with extension to involve the mucosa of the spiral colon. There were no other significant gross lesions noted. Fig (B) (2% month old Durce) revealed extremely heavy numbers of <u>Trichura</u> sp. (whipworms embedded within the mucosa of the rectum, spiral colon and cecum with the presence of a relatively fluid fecal material. Examination of the nasal cavity revealed a heavy catarrhal exudate present with associated complete loss.

POSTMORTEM CONCLUSION: Pig A - Emaciation, prolapsed rectum, moderate to severe necrotizing colitis; Pig B - severe whipworm infection, atrophic rhinitis of both the dorsal and ventral turbinates. There are no other significant gross lesions noted.

COMMENTS: Several pathologic processess are evident in both pigs of either infectious and parasitic origin and probably account for the increased mortality and "poor doers" in this herd. Other tests are being performed and further reports are pending.

Sincerely,

Ely, DVM. Ph Pathelogist

RWE:mf



DIAGNOSTIC LABORATORY ST LES 11 11. 10 bear 'AN 1 6 1986 Condition of Specimen: Telephone: 405-824-8623 C Excelent OKLAHOMA ANIMAL DISEASE D Autohzed C Other (Specify) THORNOSSE SHEDH BARDIN PLEASE FILL OUT COMPLETELY WITH BLAC 5R LILLER VETERINARIA I OWNER . ADDRES 14/2 547 9140 ADDRESS 710 8704 918 * PHONE . - 86 10 Location of Case (County & Stater al DATE SUSMITTED. Other (Spearly) Louis Owner Mail Delivered by: Veterinarian HISTORY AGE ZZ BARED X SPECIES _ · · SEX . WEIGHT da. 12 MBER IN HERD FLOCK NUMBER SICK (excluding dead) _ NUMBER DEAD RAISED ON OWNER'S PREMISES? YES 4 IF PURCHASED. WHEN? ANY RECENT INTRODUCTIONS TO HERD (FLOCK)? YES DATE INTRODUCED_ 14-86 1 -FIRST NOTICED SICK CLINICAL SIGNS (INCLUDE TEMPERATURES) & POSTMORTEM FINDINGS Lee rogracio 895527 VACCINATION HISTORY DATE & Dams 0 FIELD DIAGNOSIS . EXAMINATIONS REQUESTED -MATERIALS SUBMITTED - H POPMAN -Cry your RIGERATED #FROZEN Necropav E Live Animal E Matopathology" Dead Anmai Fetus Cytology Blood (Clotted) Mycopiasma isolation Biood (Anticoag) L Bacteriology Serum Antibiotic Senisitivity C Kidney Virus Isolation & Identification** -Over Farasnoogy floatation Billestine Mycology U Lung 5 0 Hernetotopy CBC 0 TUSTICS' (Show boston on model) 5 Puoresoent Antibody 5 -0 Scisen n Boacily test n 0 D 0 Brain 0 Lymph Node 5 õ ٥ 0 5 Eve 0 3 MR. 0 0 0 5. 0 UNING d Carabrospinal Puid D On P on Comencia . n 1 e ... r \mathbf{m} 0 1 C 0 0 · Froman Seau as not sub 5 C Other . n ** Formain fixed Beause not suited longi inform Laeseven ne note

Rovered April, 1979

127

	MOTAMACINE INFORMATION				NOLO	
					-	Lab Use
1	Animal Identification Name/No.	Scecies	Breed	Bex	Age	One
1					1	1
-						
			1	1	4	
				1		
-						
				1		1
-						1
			1		1	4
iami	ples for serology should b	e collected in ster	nie vacutair	er tubes	without	anticoa
t 8	-D tubes are satisfactor	y Bruceilosis tub	es and sor	me other	vacuun	n lubes
enti	y contain toxic residues					
Dor	atory charges are double	ed for cases not o	owned by (Okianoma	citizen	15
				ADDI	TIONA	L REMA
	·····	,				
	4-5	1 .110	12-			
	- per			1		+
	- per			- k	· · · ·	in
	Laine	eren	~ -	6.		
	- fine	Estan	~ -	6.6		
	- sine	Estan				
	Veor d	estar doing	n -	42 42		
	For a	estar loing	y ze	42		
	L'ane Voor 6	estar doing	n -	42 42		
	Trine Tool 6	estar doing	y'e	42 42		
	Veor a	ioing	y'	42		
	For a	estar loing	y'	42 42		
	Laine Voor a	estar doing	yer.	42		
	Laine Voor 4	estar doing	yre	42 42		
	Veor a	estar doing	yre	42 42		
	Laine Voor o	Loing C	yee	42 42		
	Laine Foor a	Loing t	yes	q2		
	Laine Foor a	estar doing	yee	q2		
	Laine Joor d	estar doing	y 2 ce	42 42		
	Laine Joor of	estar doing	222	42 42		
	Joor of	estar Loing	2 2 2 2 2 2 2 2 2 2 2 2 2 2 2 2 2 2 2	42 42		
	Joor of	Comp doing		42 42		
	Laine Joor	estar doing		42 42		
	Laine Joor d	estar doing		42		



Oklahoma State University

OKLAHOMA ANIMAL DISEASE DIACNOSTIC LABORATORY

STILLWATER. OKLAHOMA 74078 (405) 624-6623

January 21, 1986

Mr. Butch Garver P.O. Box 393 Gore, Oklahoma 74435 Case No.: 112782 Owner: 5R Referral Lab.: Boren Veterinary Medical Teaching Hospital Parasitology Laboratory Oklahoma State University Stillwater, Oklahoma 74078

Dear Mr. Garver:


Attached hereto are the results of tests performed by an outside laboratory to which materials were referred from the above referenced case.

Sincerely,

Ear 1 - m Ray W. Ely, DVM, Oth

Pathologist

RWE/cas



APPENDIX 5.1.1

AIR SAMPLE ANALYSES BY SFC AND ORNL

SEQUOYAH FUELS CORPORATION

POST OFFICE BOX 25851 + OKLAHOMA CITY OKLAHOMA 73125

January 27, 1986

FEDERAL EXPRESS

Mr. Edward Shum Uranium Fuel Licensing Branch Division of Fuel Cycle & Material Safety Office of Nuclear Material Safety & Safeguards United States Nuclear Regulatory Commission Willste Building 7915 Eastern Avenue Silver Spring, Maryland 20912

> Re: License SUB-1010 Docket 40-8027

Dear Mr. Shum:

During our telephone conversation on Monday, January 20, 1986, concerning NRC's evaluation of the UF release at the Sequoyah Facility on January 4, you asked for background information on several items. Sequoyah Fuels continues to assemble and evaluate the data obtained from sampling following the accident. As those data are finalized, they will be submitted to NRC for its use. We provide the information below:

 Describe the plant and residential potable water supply sources for the area south of the plant which may have been in the path of the plume.

The plant potable water supply is provided by a pipeline from the Lake Tenkiller reservoir. The subject area residences are served by Rural Water District #5 which buys water from the Sequoyah County Rural Water Supply Association. The water comes from the Tenkiller reservoir. Ron Barnett, Sequoyah County Health Department, indicates all residences south of the facility and I-40 are on the rural water system for potable water. They are in the process of confirming this.

2. Is the Robert S. Kerr Reservoir used as a potable water source?

To SFC's and the Sequoyah County Health Department's knowledge, it is not. The water is for recreational purposes.

 Please provide a map showing the location of residences south of I-40 in the probable downwind plume direction.

Attached is a map showing the residences, identified by family name. Additionally, known livestock assemblages in the area on the day of the release are noted. A key to the map is provided for reference noting family size and where sampling of livestock was done. Livestock sampling results will be submitted when completed. Mr. Edward Shum January 27, 1986 Page Two

> A map of the subject area, assembled from USGS quadrangle maps is also provided for your use. The USNRC on-site response team was provided with mylars of the subject study area, this map was one of them.

 Please provide facility weather data for January 4, 1986. Information concerning temperature, relative humidity and general description (visibility) besides the wind speed and direction would be helpful.

The plant wind speed and direction chart for the accident period are attached with reference to 11:30 AM, January 4, 1986 noted (1400 on the graph). Data from the Ft. Smith station concerning wind speed, direction, temperature and humidity are attached. The Oklahoma City station data for 11:00 AM, January 4, 1986 is also provided.

5. Please provide the results of uranium and fluoride analysis for the environmental air monitoring stations around the facility and off-site locations for the periods before, during and after the UF₆ release.

The data for the fenceline and off-site monitor stations are attached. Sampling was done according to those procedures prescribed in license SUB-1010, Chapter 5.

Please provide the results of urinalysis taken of employees and off-site residents.

This data was submitted on January 22, 1986, under separate cover to NRC Region IV. A copy was submitted to William Crow, NRC staff. Any updates to this submittal will be forwarded.

7. Have measurements been made to determine the amount of UF₆ that was potentially released to the environment?

At the time of the release emergency actions were taken which included water spray, with fog nozzles, to minimize the release. The resulting wash down water was routed to the facility emergency basin for positive containment. Analysis of the basin contents resulting from the response action, shows that approximately 9,250 pounds uranium were recovered. This is 13,700 pounds of UF₆ or about 46 percent of the material believed to have been in the cylinder. Further analysis of sod removed from in front of the building and facility decontamination solutions continues and these amounts will be added as they are reported. Edward Shum January 27, 1986 Page Three

Please contact me if you need further information at this time.

Sincerely,

J.C. Stanter

J.C. Stauter, Director Nuclear Licensing & Regulation

JCS/br

Attachments As Stated

cc: Robert L. Craig, Oklahoma State Department of Health

SEQUOYAH FACILITY FENCELINE AIR SAMPLING

	E-1	EAST	E-	2 WEST	Ε-	3 SOUTH	E-	4 NORTH
Sampling Time	MPC	Fluoride	MPC	Fluorlde	MPC	Fluorlde	MPC	Fluoride
01/01/86 (0001) to 01/01/86 (2400)	0.16		0.18		0.18		0.17	
01/02/86 (0001) to 01/02/86 (2400)	0.32		0.38		0.42		0.34	
01/03/86 (0001) to 01/03/86 (2400)	0.17		0.22		0.28		0.28	
01/04/86 (0001) to 01/04/86 (1800)	154.0	0.1*	0.15	0.0015*	0.28	0.018*	0.22	0.0028*
01/04/86 (1800) to 01/05/86 (1030)	0.90		0.40		0.46		0.50	
01/05/86 (1030) to 01/05/86 (2400)	0.19		0.19		0.28		0.19	
01/06/86 (0001) to 01/06/86 (1600)	0.19		0.20		0.20		0.15	
01/06/86 (1600) to 01/07/86 (1600)	<0.05		0.08		0.06		0.08	
01/07/86 (1600) to 01/08/86 (1600)	0.06		0.07		0.06		<0.05	
01/08/86 (1600) to 01/09/86 (1600)	<0.05		<0.05		<0.05		<0.05	
01/09/86 (1600) to 01/10/86 (1700)	2.00		<0.05		<0.05		<0.05	
01/10/86 (1700) to 01/11/86 (1800)	0.07		0.09		0.10		0.10	
01/11/86 (1800) to 01/12/86 (1500)	0.40		0.52		0.46		0.46	
01/12/86 (1500) to 01/13/86 (2400)	0.20		0.20		0.18		0.18	
01/14/86 (0001) to 01/14/86 (2400)	0.40	0.0021**	0.64	0.0020**	0.64	0.0005**	0.64	0.0022**
01/15/86 (0001) to 01/15/86 (2400)	0.18		0.40		0.48		0.30	

Fluoride results reported in ugm/1 Fluoride detection limit is 0.0005 ugm/1 * On 12/31/85 (1200) Off 01/04/86 (1800) ** On 01/04/86 (1800) Off 01/14/86 (1200) Sampling rate for fluoride is 0.2 cfm. Fluoride collected on calcium-impregnated filter paper.

Radionuclide results reported as fraction of MPC MPC = 5×10^{-12} uCl/ml Detection limit is 2.5×10^{-13} uCl/ml (daily) Sampling rate for radionuclides is 1 cfm.

SEQUOYAH FACILITY OFF-SITE ENVIRONMENTAL AIR SAMPLING

	2103		2105		2106		2107		2108	
Sampling Time	MPC	Fluoride								
12/10/85 (1200) to 12/17/85 (1200)	<0.007	<0.0005	<0.007	<0.0005	<0.007	<0.0005	<0.007	<0.0005	<0.007	<0.0005
12/17/85 (1200) to 12/24/85 (1200)	<0.007	<0.0005	<0.007	<0.0005	<0.007	<0.0005	<0.007	<0.0005	<0.007	<0.0005
12/24/85 (1200) to 12/31/85 (1200)	<0.007	<0.0005	<0.007	<0.0005	<0.007	<0.0005	<0.007	<0.0005	<0.007	0.00066
*12/31/85 (1200) to 01/04/86 (1800)	<0.007	0.0038	<0.007	0.00351	<0.007	0.00307	<0.007	0.00402	2.6	0.0038
*01/04/86 (1800) to 01/05/86 (1100)	<0.05		<0.05		<0.05		<0.05		<0.05	
*01/05/86 (1100) to 01/06/86 (1630)	<0.05		<0.05		<0.05		<0.05		<0.05	
01/06/86 (1630) to 01/07/86 (1530)	<0.05		<0.05		<0.05		<0.05		<0.05	
*01/07/86 (1530) to 01/08/86 (1530)	<0.05		<0.05		<0.05		<0.05		<0.05	
*01/08/86 (1530) to 01/09/86 (1530)	<0.05		<0.05		<0.05		<0.05		<0.05	
*01/09/86 (1530) to 01/10/86 (1700)	<0.05		<0.05		<0.05		NS		<0.05	
*01/10/86 (1700) to 01/11/86 (1700)	<0.05		<0.05		<0.05		<0.05		<0.05	
*01/11/86 (1700) to 01/12/85 (1500)	<0.05		<0.05		<0.05		<0.05		<0.05	
*01/12/86 (1500) to 01/14/86 (1200)	<0.05	0.00061	<0.05	0.00068	<0.05	<0.0005	<0.05	0.00051	<0.05	<0.0005

Fluoride results reported in ug/1 Sampling rate for fluoride is 0.2 cfm. Fluoride collected on calcium-impregnated filter paper.

Radionuclide results reported as fraction of MPC. Sampling rate for radionuclides is 1 cfm. MPC = 5 x 10⁻¹² uCi/mi Detection limit is 3.6 x 10⁻¹⁴ uCi/mi (weekly); 2.5 x 10⁻¹³ uCi/mi (daily) *Final count done on 01/23/86

.

Locations:

2103 -- East, 1,000 feet from plant 2105 -- $\frac{1}{2}$ Mile southwest of plant 2106 -- Carlisie School 2107 -- Highway 64 north of plant 2108 -- 1-40 south of plant

OAK RIDGE NATIONAL LABORATORY

OPERATED BY MARTIN MARIETTA ENERGY SYSTEMS, INC.

POST OFFICE BOX X OAK RIDGE, TENNESSEE 37831

February 19, 1986

Mr. Barry Zalcman US Nuclear Regulatory Commission 7915 Eastern Avenue Silver Springs, MD 20910

Dear Barry:

Enclosed are the particle size distribution measurements taken on the Sequoyah Nuclear Fuel air filters E-1 (ORNL #32129) and 2108 (Dec 31-Jan 4 sampling period). The tabular data shows the number of particles in each size range; the histogram indicates the same data in graphical form. The 2108 specimen was denoted "filter blank" by the microscopist, and two imaging modes were utilized to show the particles. The photographs are in pairs; one set is the conventional scanning electron image which shows all particles on the filter, while the greater intensity (brightness) of the backscattered SEM images indicates higher atomic numbered particles. From the backscattered images, we note that many of the particles were of low atomic numbered elements (aluminum and silicon) typical of airborne dirt. The 2108 filter contained relatively few uranium particles, compared to filter E-1.

We performed characteristic uranium x-ray "maps" of typical areas on the filter surfaces to verify the exact location and size of the uranium particles. Please note that an overlay of the SEM and the U maps identify the locations of the uranium particles.

Please notify me if we can be of further assistance in data interpretation.

Sincerely,

Joeph H. Stewart, Je

Joseph H. Stewart, Jr. CAPA Group

JHS:1p

cc: W. R. Laing

Sequoyah Air Sample - E-1 (ORNL #32129) Taken from 0001 Hours 1/4/86 - 1800 Hrs. 1/4/86

Class				Diam	eter,	u m		No.	of	Part	ticle	25
0 1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16 17 18 19 20 21				.32 .48 .64 .80 .96 1.11 1.27 1.43 1.59 1.75 1.91 2.07 2.23 2.39 2.55 2.71 2.87 3.02	32 48 64 80 96 -1.11 -1.27 -1.43 -1.59 -1.75 -1.91 -2.07 -2.23 -2.39 -2.55 -2.71 -2.87 -3.02 -3.18 -3.34				0 0 0 7 11 27 32 32 16 21 8 9 1 1 3 1 0 0 2			
	Class	No. Particles	Percent in Class									
	0123456789011231456789011231456718921	0000711752236128891131002	00035326680444001000	****		****	****	**** **** ***	***	***	***	

10.1

١

١

Sample - Filter Blank

Clas	s	Diameter, µm	No. of Particles
0 1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16 17 18 19 20 21		<.53 .5361 .6169 .6977 .7785 .8593 .93-1.01 1.01-1.09 1.09-1.17 1.17-1.25 1.25-1.33 1.33-1.41 1.41-1.49 1.49-1.57 1.57-1.64 1.64-1.72 1.72-1.80 1.80-1.88 1.88-1.96 1.96-2.04 2.04-2.12 >2.12	0 3 2 3 1 1 2 2 3 1 2 3 1 2 3 1 2 3 1 2 3 1 2 3 1 2 3 1 2 3 1 2 2 3 1 2 2 3 1 2 2 3 1 2 2 3 1 2 2 3 1 2 2 3 1 1 2 3 1 2 2 3 1 1 2 2 3 1 1 2 2 3 1 1 2 2 3 1 1 2 2 3 1 2 3 1 1 2 2 3 1 1 2 2 3 1 1 2 2 3 1 1 2 2 3 1 1 2 2 3 1 1 2 2 3 1 1 2 2 3 1 1 2 2 1 1 2 3 1 1 2 2 1 1 2 2 3 1 1 2 1 1 2 1 1 2 1 1 2 1 1 2 1 1 2 1 2 1 1 2 1 1 2 1 1 2 1 1 2 1 1 2 1 1 1 2 1 1 1 2 1 1 2 1 1 2 1 1 2 1 1 2 1 1 1 2 1 1 1 1 2 1 1 1 2 1 1 1 1 2 1 1 1 1 2 1 1 1 1 1 2 1
Class No. Particles	Percent in Class		
0 0 1 3 2 2 3 3 4 1 5 1 6 2 7 2 8 3 9 1 10 2 11 1 12 0 13 0 14 0 15 0 14 0 15 0 16 0 17 0 18 2 0 0 21 0	13 ***** 9 ***** 13 ***** 4 ***** 9 ***** 9 ***** 13 ***** 13 ***** 9 ***** 9 *****	********	****

Filter Blank Filter Blank (Backscatter Image) Sample Sample M&C 2000X M&C Mag: 2000X S --S Mag 2/18/86 DATE Zalcman 00 Barry SAMPLE NO. 7006 25KU X2000 3001 19 91 M&C 2 Sample M&C Filter Blank Filter Blank (Backscatter Image) Sample 14 S - M&C 2000X 2000X Mag Mag -4216 8-427. **DAK RIDGE NATIONAL LABORATORY** ANALYTICAL CHEMISTRY DIVISION SNE 5 RM Bldg Willste ACD REQ. NO. 51195 COMMENTS BUILDING 3005 X2000 10.00 1#& C 10,00 25K.U 25KU X2000 3004 M&C

PHOTOGRAPHIC REPORT SHEET

Other photos in this set were omitted because they were difficult to reproduce.

APPENDIX 5.2.1

ľ

INTERORGANIZATIONAL COMPARISON RESULTS

	OSDH		KMTC			ORNL			
	F	F	U	wt Loss	F	U	wt Loss		
Water	mg/l	mg/l	mg/l		µg/ml	ppm			
COMP-3	. 45	. 4	.018		. 38	<.03			
COMP-4	. 27	. 1	<.002		.06	₹.03			
COMP-5	. 24	. 2	<.002		. 18	₹.03			
Soil	mg/kg	ppm	ppm	%	µg/g	µg/g	%		
COMP-1	5.44	90	72	24.3	-	-	· · · ·		
COMP-2	2.76	70	23	26.0	2.9	8.4	26.2		
COMP-3	1.96	90	15	25.6	4.6	6.1	26.1		
COMP-4	1.48	140	<10	26.6	2.2	4.0	26.8		
COMP-5	5.6	80	15	25.4	2.6	5.2	26.2		
Vegetation	mg/kg	ppm	ppm	*	µg/g	µg∕g	%		
COMP-1	488	250	900	17.8	390	4970	16.8		
COMP-2	164	90	80	22.5	128	44	20.7		
COMP-3	238	120	160	19.0	144	22	18.3		
COMP-4	7.36	17	. 95	16.1	5		2 14.8		
COMP-5	10.8	18	3	11.8	18		9 10.6		

Interlaboratory Comparison Results

Note: Fluoride in soil from KMTC includes soluble and insoluble fluoride; fluoride in soil from OSDH and ORNL only includes soluble fluoride.

APPENDIX 5.2.2

K-M TECHNICAL CENTER LABORATORY ANALYSIS PROCEDURES

SEQUOYAH FUELS CORPORATION

POST OFFICE BOX 25861 + OKLAHOMA CITY OKLAHOMA 13125

January 27, 1986

FEDERAL EXPRESS

Donald A. Cool, Ph.D. Uranium Fuel Licensing Branch Division of Fuel Cycle & Material Safety C vice of Nuclear Material Safety & Safeguards United States Nuclear Regulatory Commission Willste Building 7915 Eastern Avenue Silver Spring, Maryland 20912

> Re: License SUB-1010 Docket 40-8027

Dear Dr. Cool:

Per your request of Friday, January 24, 1986, attached are the Kerr-McGee Technology Diviison, Technical Center analytical procedures, equipment list and calibration procedures applicable for the analysis of samples from the Sequoyah Facility.

Copies of these procedures were provided to Dr. Blair Nickolas, USNRC Region IV, during a visit to the Technical Center on January 7, 8, 9 and 10, 1986.

If you have any further questions concerning the analytical procedures or quality control/quality assurance protocols used, please contact Dr. Garet Van de Steeg at (405) 341-8551.

Sincerely,

J. C. Stauter

J.C. Stauter, Director Nuclear Licensing & Regulation

JCS/br

Attachments

APPENDIX I

EQUIPMENT LIST AND CALIBRATION PROCEDURES

	rumentation for the Detection of Nuclear Events
1.	Gamma Radiation
	a. Packard Auto-Gamma Multi-Channel Analyzer
	 Canberra Gemini Series 90 Gamma Pulse Height Analysis System
2.	Alpha and Beta Radiation
	a. Canberra Model 2400 Low-Background Gas Proportional Counter
	 b. Packard Tricarb Liquid Scintillation Spectrometer
	c. Ludlum Model 148 (GEV-1) Scintillation Counter
	d. Nuclear Measurements PC-4 Proportional Counter
	e. Random Scintillation Counter
	f. Princeton Gamma Tech Silicon Surface Barrier System
	g. Ludlum Model 248 (GEV-2) Scintillation Counter
	h. Capability
	i. Tracers
Indu	ctively Coupled Argon Plasma Emission Spectroscopy (ICP)
Aton	ic Absorption Spectroscopy (AAS)
Ator	
Ator	
	a. Perkin-Elmer (PE) Model 460
	a. Perkin-Elmer (PE) Model 460 b. Perkin-Elmer (PE) Model 280 y Radiation
	a. Perkin-Elmer (PE) Model 460 b. Perkin-Elmer (PE) Model 280
<u>X-R</u> 2	 a. Perkin-Elmer (PE) Model 460 b. Perkin-Elmer (PE) Model 280 by Radiation a. Siemens D-500 X-ray Diffraction (XRD) System
<u>X-R</u> 2	 a. Perkin-Elmer (PE) Model 460 b. Perkin-Elmer (PE) Model 280 by Radiation a. Siemens D-500 X-ray Diffraction (XRD) System b. Siemens SRS-300 X-ray Fluorescence (XRF) System
<u>X-R</u> 2	 a. Perkin-Elmer (PE) Model 460 b. Perkin-Elmer (PE) Model 280 a. Siemens D-500 X-ray Diffraction (XRD) System b. Siemens SRS-300 X-ray Fluorescence (XRF) System a. Perkin-Elmer (PE) Model 1750 (FT-IR) b. Varian Superscan I-BE, Beckman Model DU, Beckman
<u>X-R</u> 2	 a. Perkin-Elmer (PE) Model 460 b. Perkin-Elmer (PE) Model 280 by Radiation a. Siemens D-500 X-ray Diffraction (XRD) System b. Siemens SRS-300 X-ray Fluorescence (XRF) System a. Perkin-Elmer (PE) Model 1750 (FT-IR)

Appendix I

			Pa
Micro	800	<u>ру</u>	
1.	Opt	ical Microscopy	
	a.	Petrographic (Zeiss)	
	b.	Metallographic (Unitron)	
2.	Ele	ctron Microscopy	
Instru	umer	ntation for Analysis by Electrochemical Means	
	a.	Sargent Electrodeposition Unit and Mercury Cathode	
1	ь.	pH Meters and Specific Ion Meters	
	c.	Recording Titrators, Mettler and Metrohm/Brinkman	
Instru	umen	station for the Analysis of Gases and Liquids by	
Gas Cland Ic	hron	atography, High Pressure Liquid Chromatography Chromatography	
	R	Gas Chromatographs	
1	b.	High Pressure Liquid Chromatograph	
0		Ion Chromatograph	
d	1.	Perkin-Elmer Model Integrator	
٠		Hewlett-Packard Model 3390A Integrator	
Organi	LC A	nalysis	
		Coulometrics Carbon Analyzer	
b		Ramsbottom Carbon Residue Analyzer	
fuffle	Fu	rnaces	
Surfac	e A	rea Analyzer	
liscel	land	eous Equipment	
a		Frantz Isodynamic Separator, Andreason Pipettes, Sieves	
b		Sample Preparation	
c		Cleveland Open Cup and Pensky-Martins Closed Cup Tester	
d		Cameras	
e		Mechrolab Vapor Pressure Osmometer	
f		Brookfield Model HAT Viscometer	
g		pH Meters	
h		Specific Conductivity Meter	
1	2	Balances	

Appendix I

м.	Quality	Assurance - Quality	Control	18
	а.	Definition		19
	b.	Quality Assurance	Plan	19

Page

FUNCTIONS AND USE

A. Instrumentation for the Detection of Nuclear Events

- 1. Gamma Radiation
 - a. A Packard Auto-Gamma Multi-Channel Analyzer is a sodium iodide based gamma detection system. This analyzer is used for routine gamma spectral analysis where only limited spectral resolution is required. This instrument has a 300-sample capacity automatic sample changer, a 1024 channel multi-channel analyzer for storing the data, and a teletype for recording the spectral information. The entire system is situated in a controlled environment room where temperature is maintained at 21°C ± 1°C. The system and its essential components are to be maintained as follows:
 - Maintenance and repair of electronic components are to be performed by the manufacturer as required.
 - The counting system should be effficiency calibrated during each 300-sample run and energy calibrated weekly.
 - Outage of electricity sufficient to cause loss of bias to the detectors, or repair of detector or electronics is cause for recalibration prior to further analyses being performed.
 - b. A <u>Canberra Gemini Series 90 Gamma Pulse Height Analysis System</u> is used for qualitative and quantitative analysis. Four gamma detectors, three wide range Ge(Li)s and an intrinsic Ge are situated in low background gamma shields fabricated from 6-inch thick pre WW-II steel plate. This system is used for absolute measurement of gamma emitting radionuclides. The entire system is situated in a controlled environment room where the temperature is maintained at 21°C ± 1°C. The system and its essential components are to be maintained as follows:
 - Maintenance and repair for electronic components are to be performed by the manufacturer as required, but at least twice yearly.
 - Ge(Li) and Ge detector dewars are to be filled with liquid nitrogen weekly. The date of the liquid N₂ fill is to be logged into a file maintained at the detectors.
 - 3) The Ge(Li) and Ge detectors are calibrated at least quarterly, and checked at least weekly, for efficiency and energy calibrations using standards directly traceable to the U. S. National Bureau of Standards.

2. Alpha and Beta Radiation

A Canberra Model 2400 Low-Background Gas Proportional Counter a. is used for the detection of alpha and/or beta radiation. The thin-window (gold-coated mylar) detector is used to detect alpha and beta particles from the sample. The sample detector is one part of the detector system, the other part being an aluminum cased guard detector. The alpha and beta particles from the sample pass easily through the mylar window, but are stopped before entering the aluminum cased guard detector. The guard detector is used to detect background radiation from cosmic radiation and gamma radiation, which might register as a detected particle. The guard detector is used in an anticoincidence mode to reject unwanted background. The sample counted is mounted in a two inch planchet, allowing for the counting of precipitates and electrodeposits. The system has an automatic sample changer and microprocessor control for unattended operation. Up to five distinct programs may be entered, with variables such as sample number, preset time, bias voltage, discriminator window, efficiencies, and backgrounds entered for each program. The entire system is situated in a controlled environment room where temperature is maintained at 21°C ± 1°C. The system and its essential components are to be maintained as follows:

2

- Maintenance and repair of electronic components are to be performed by manufacturers as required.
- Efficiency calibration and plateau curves should be verified every three months, and following each repair.
- b. A Packard Tricarb Liquid Scintillation Spectrometer is used for the detection of alpha and/or beta radiation. This is a two-channel instrument which is capable of detecting two different radiation energies, or isotopes, simultaneously. The samples to be analyzed are first put in liquid form (in water, toluene or dioxane) and then mixed with special fluorescent materials. The nuclear radiation excites the fluorescent materials which, subsequently, fluoresce (emit light) and this light is detected by photomultiplier tubes. This is an automatic instrument capable of handling up to 150 samples at a time and can be operated overnight unattended. The entire system is situated in a controlled environment room where temperature is maintained at 21°C ± 1°C. The system and its essential components are to be maintained as follows:
 - Maintenance and repair of electronic components are to be performed by the manufacturer as required.
 - The counting system should be efficiency and energy calibrated at least once a year.

- The system should have its "normalization" checked every three months.
- Repair of detector or electronics is cause for recalibration prior to further analyses being performed.
- c. Two Ludlum Model 148 (GEV-1) Scintilation Counters are used for the detection of alpha radiation. These are five-station scintillation counters using photomultiplier tubes to count the light pulses emanating from the zinc sulfide detectors. The zinc sulfide detectors are disposable mylar disks coated with ZnS(Ag) which are placed in intimate contact with, and atop, the sample to be counted. The sample may be either a precipitate mounted on a filter paper or a deposit mounted on a stainless steel disk. The system is manually operated. The background is typically less than one count per 1000 minutes. The system is to be maintained as follows:
 - Maintenance and repair of electronic components are to be performed by the manufacturer as required.
 - 2) The counter is to have its efficiency calibrated every three months and a "plateau" curve checked yearly; each of these checks are to be run following any repair.
- d. A Nuclear Measurements PC-4 Proportional Counter is used for alpha or beta counting. This is a manually operated, windowless, P-10 gas, proportional counter which is used to calibrate reference alpha or beta standards. Additionally, this counter is used for routine, but infrequent, low-level alpha or beta counting of selected samples. The entire counter is situated in a controlled environment room where temperature is maintained at 21C° ± 1°C. The counter is to be maintained as follows:
 - Maintenance and repair of electronic components are to be performed by the manufacturer as required.
 - The counting chamber is to be removed, checked and cleaned (if necessary) every three months.
 - The counting efficiency should be recalibrated every three months following the maintenance described in 2) above.
- e. A <u>Random Scintillation Counter</u> is used for radium determinations by radon emanation. This photomultiplier tube based unit counts light pulses originating from within Lucas cells (containing radon) placed atop the PMT's. The interiors of the Lucas cells are coated with zinc sulfide which transforms the alphas emitted by radon and radon daughters to light pulses. The counter is to be maintained as follows:

- Maintenance and repair of electronic components are to be performed by the manufacturer as required.
- The counter is to have its "plateau" curve checked yearly or following any repair.
- The Lucas cells are to have their efficiency checked every three months or as required, whichever occurs first.
- f. A Princeton Gamma Tech Silicon Surface Barrier System, in conjunction with a Canberra Model 8100 Multi-Channel Analyzer, Canberra Bias Supplies and Amplifiers, and Houston Omnigraphic X-Y Recorder, is used for alpha pulse height analysis. Actinide elements are chemically separated into individual elements which are subsequently eectrodeposited onto stainless steel disks. The disks are then placed inside a vacuum chamber and beneath a silicon surface barrier detector (PGT-300-25-100). The instruments analyze the alpha emitted by the sample for energy, the energy being indicative of the isotope undergoing decay. The entire system is situated in a controlled environment room where temperature is maintained at 21°C ± 1°C. The system and its essential components are to be maintained as follows:
 - Maintenance and repair for electronic components are to be performed by the manufacturer as required.
 - 2) Detectors and associated electronics are to be calibrated for energy and efficiency at least quarterly using KMTC standard reference sources prepared from NBS or NBS-traceable standard reference materials.
 - Repair of detectors or electronics is cause for recalibration of affected components prior to further analyses being performed
- g. A Ludlum Model 248 (GEV-2) Scintillation Counter is used for the detection of alpha and beta radiation. This is a five-station scintillation counter using photomultiplier tubes to count the light pulses emanating from the zinc sulfide and beta phosphor detectors. The zinc sulfide detectors are disposable mylar disks coated with ZnS(Ag) which are placed in intimate contact with, and atop, the sample to be counted. The beta phosphors are 1/64 inch thick, Pilot-B beta phosphor disks. The sample is precipitate mounted on a filter paper. The system is manually operated. The background is typically less than one count per 1000 minutes for alpha and 10 counts per minute for beta. The counter is situated in a controlled environment room where the temperature is maintained at 21°C ± 1°C. The counter is to be maintained as follows:

- Maintenance and repair of electronic components are to be performed by the manufacturer as required.
- 2) The counter is to have its efficiency calibrated every three months and a "plateau" curve checked yearly; each of these checks are to be run following any repair.

h. Capability

With the present intsrumentation, we are capable of performing most types of radiometric analyses such as alpha, beta or gamma spectroscopy, Ra-226, Th-230 and total U, Pu, Th, etc. Additionally, virtually any radioactive isotope can be analyzed radiometrically should the need arise. Extremely low-detection levels are maintained such that we routinely analyze at levels of 1 x 10^{-15} curies.

i. Tracers

A program to use radioactive tracers in the study of problems relating to chemical processing, analysis and engineering is available. Radioactive tracers (isotopes) are used to simplify the often difficult problem of monitoring the disposition of materials in dilute systems or to simplify otherwise difficult or impossible analytical procedures. We are licensed to use any isotopes of atomic numbers 1 through 96 in any chemical form.

B. Inductively Coupled Argon Plasma Emission Spectroscopy (ICP)

An argon plasma is maintained by the interaction of an RF field with ionized argon gas and is reported to reach temperatures as high as 12,000°K near the top load coil. Analytical measurements are made at plasma observation heights that correspond to temperatues in the range of 5500-8000°K. These high temperatures ensure complete dissociation and effective atomization of the various sample matrix elements, thus minimizing chemical interferences.

In order to form an argon plasma, RF power is applied through induction coils which surround the quartz plasma torch assembly. This creates an oscillating magnetic field. Three separate argon streams ultimately determine the symmetry of the plasma. A tangential stream of argon (plasma gas) flows between the outer two quartz tubes of the ICF torch assembly. This stream produces a vortex stabilization of the plasma discharge. A secondary stream of argon (coolant gas) passes between the sample introduction jet and the adjacent concentric tube. The coolant gas flow lifts the discharge slightly above the tip of the sample introduction jet to protect the quartz jet from the intense heat. The sample aerosol is carried through the center quartz tube of the assembly via a third argon stream (carrier gas). The plasma is initiated when

the argon is exposed to the discharge of a Tesla coil, thus creating seed electrons and ions making the gas conductive. Inside the induced magnetic field, the charged particles are forced to flow in a closed annular path. As the particles meet resistance to their flow, ohmic heating takes place and additional ionization occurs. The process occurs almost instantaneously, and the plasma expands symmetrically.

The basis of all emission spectroscopy is that atoms or ions in an energized state, either directly or indirectly as a result of heating, spontaneously revert to a lower energy state and so doing emit a photon of energy. The ICP method of heating offers a much better chance for avoiding the self-absorption and self-reversal effects that are characteristic of both arcs and flames.

The light emitted by analyte atoms or ions must be converted to electrical signals that can be measured quantitatively. This is done by resolving the light into its component radiation by means of a diffraction grating and then measuring selected light intensities with photomultiplier tubes at specific analyte wavelengths.

After appropriate analyte emission wavelengths are chosen, the spectrometer can be calibrated using standard solutions. Then, analyte concentrations can be determined for all sample solutions.

ICP Instrumentation

Applied Research Laboratories (ARL) Model 3520

This ICP is a sequential instrument with a scanning 1.0 meter vacuum monochromator. It is equipped with a 1080 lines/mm concave reflecting grating (produced using laser interferometry) with a horizontal Paschen-Runge mounting. Light intensities are measured by two photomultiplier tubes (PMTs), one for UV and the other for visible up to IR. The two PMTs are mounted on a carriage. Metal-dielectric-metal, narrow bandpass filters allow the use of higher orders to provide higher resolution. The positioning of the primary slit, the selection of filters, the setting of PMT millivoltages, the selection PMTs, and the selection of PMT positions are microprocessor controlled.

The radio frequency (RF) generator is a crytal controlled oscillator operating at 27.12 MHz.

The central processor unit (CPU) is comprised of a PDP 11/23 computer with 16 bit parallel logic. It is equipped with a real time clock, 256K bytes MOS memory and bootstrap.

The computer is also equipped with an RLOI disk controller, a VT 100 display terminal, and an LA-120 printer.

The SAS/DPS software is a program package consisting of an RT-11 system and an SAP-11 Spectrometer Automation Program.

This ICP system is to be maintained as follows:

 Maintenance and repair of major components are to be performed by the instrument and computer manufacturer as required.

7

- There will be yearly preventive maintenance conducted by an ARL service engineer.
- A system log will be generated daily and compared to the specifications set by the manufacturer.
- 4) The nebulizer will be cleaned daily.
- Sensitivity will be evaluated monthly as prescribed by the manufacturer.

C. Atomic Absorption Specroscopy (AAS)

Atomic absorption is an analytical technique that is used for the quantitative determination of metal concentrations in a wide variety of materials.

Atomic absorption is the process that occurs when a ground state atom absorbs energy in the form of light of a specific wavelength and is thus elevated to an excited state. The amount of light energy absorbed at this wavelength will increase proportionately as the number of atoms of the selected element in the light path increases. Therefore, the linear relationship between the various analyte concentrations present in calibration standards and the light absorbed by those standards can be used to determine unknown analyte concentrations present in sample solutions.

The basic instrumentation for atomic absorption requires a primary light source, an atom source, a monochromator to isolate the specific wavelength of light to be used, a detector to measure the amount of light accurately, electronics to treat the signal and a data logging device to display the results. The light source used is either a hollow cathode lamp (HCL) or an electrodeless discharge lamp (EDL). The atom source used must produce free analyte atoms from the sample. The source of energy for free atom production is heat, most commonly in form of an air-acetylene or nitrous oxide-acetylene flame or a graphite furnace. The sample is introduced as an aerosol into the flame and as a liquid into the graphite furnace. The flame burner head or furnace cell is aligned so that the light beam passes through the flame or graphite tube where the light is absorbed.

Atomic Absorption Instruments

a. Two Perkin-Elmer (PE) Model 460s

These are microcomputerized double-beam atomic absorption spectrometers. They are equipped with a PE safety interlock burner

. . .

system for use with nitrous oxide-acetylene and air-acetylene flames. One of these is equipped with a PE Model HGA-500 graphite furnace for ultra-trace (low ppb) metals analyses. Hollow cathode and electrodeless discharge light sources are used. These systems are to be maintained as follows:

- Maintenance and repair are to be performed by the manufacturer as required.
- 2) The burner system is to be cleaned with a Contrad 70 cleaning solution and rinsed with 1+1 hydrochloric acid and deionized water. This is to be done on a daily basis or more frequently if required.
- 3) The nebulizer system is to be cleaned on a daily basis with a Contrad 70 cleaning solution and rinsed with 1+1 hydrochloric acid and deionized water. Replace the needle assembly and venturi every six months or sooner if required. Aspiration rate is to be adjusted to ~6 ml/min.
- Visual inspection of the optics for "clouding" should be performed every 3 months, or sooner if needed.
- The D₂ arc output and alignment should be checked monthly, or sooner if required.
- 6) The HGA-500 graphite furnace is to be checked as follows:
 - Internal gas flows check with bubble flow meter as described in the HGA-500 manual every 3 months, or sooner if required.
 - a server and the stated server all
 - Replace furnace cones every 3 months, or sooner if required.
 - Replace furnace quartz windows every 3 months, or sooner if required.
 - Furnace cooling system check every 3 months, or sooner if required.
 - Inspect and clean optical heat sensor before each analysis.
- Sensitivity checks must be conducted every 3 months, or sooner if required.
- A system log will be prepared on a daily basis and compared to the specifications set by the manufacturer.

b. Perkin Elmer (PE) Model 280

This single beam instrument is equipped with a PE safety interlock burner system for use with nitrous oxide-acetylene and air-acetylene flames. It is to be maintained using the same procedures as are used for the (PE) Model 460s.

9

The instrument was recently equipped with a 60 Hz mechanical chopper. Now EDLs can be used with the Model 280 specifically for hydride generation and cold vapor mercury analyses. The MHS-20 Hydride Generation Apparatus has been installed on the Model 280.

D. X-Ray Radiation

X-ray diffraction (XRD) and x-ray fluorescence (XRF) both utilize high energy electromagnetic radiation. This radiation is produced when electrons are "boiled off" the surface of a heated tungsten filament, and accelerated by a high voltage potential toward a suitable target. Primary x-rays are generated by the stepwise deceleration of high velocity electrons as they approach the atoms that make up the target of the x-ray tube. These x-rays, traveling in straight lines at the speed of light, exhibit properties of both waves and discreet particles. This duality allows them to interact with matter in several ways. Examples of wave properties are reflection, refraction, polarization, coherent scatter, velocity and diffraction. Examples of corpuscular properties are incoherent scatter, gas ionization, photoelectric absorption and scintillation production.

X-ray diffraction occurs when x-rays are scattered by a crystalline material. The scattered x-rays will interact by constructive interference to yield diffracted rays. If x-rays of known wavelength λ impinge on a crystal whose lattice planes are separated by the distance d, the radiation will be strongly reflected at specific angles θ . Since the wavelength λ is known and the 2θ angles of the reflections measured, the corresponding values of d (lattice plane spacings) can be calculated. This phenomenon was investigated by W. L. Bragg and his father W. H. Bragg, and is defined as Bragg's equation: (n λ = 2d sin θ).

X-ray fluorescence occurs when primary x-rays create a vacancy by expelling an electron from the inner orbitals of the atoms of a sample. Electrons from orbitals that are farther away from the nucleus can then fall into the vacancy created. The difference in energy between the orbitals is emitted as a secondary x-ray of a precise quanta. This radiation is dispersed into the individual spectral lines by reflection at an analyzer crystal. The amount of the element present is related to the intensity of the secondary fluorescent radiation detected from the sample.

a. X-ray Diffraction (XRD) (Siemens D-500)

X-ray diffraction permits the rapid qualtitative identification of crystalline solids by yielding an XRD pattern that is characteristic for that crystalline material. This pattern is compared to the patterns of known compounds in the JCPDS literature

Heres

(Joint Committee for Powder Diffraction Standards) for matching. This search of the literature can be performed either manually or by computer. The entire JCPDS data base (40,000 compounds, inorganic and organic) can be searched by computer while the Siemac microprocessor is independently running a scan.

Quantitative x-ray diffraction analysis is also possible by analyzing a series of carefully chosen peaks present in both the samples and a prepared standard for the phases of interest.

With the automated sample changer, up to 40 samples may be run under the control of the PDP 11/23+ computer. This enables overnight, unattended operation for both qualtitative and quantitative analysis.

Unit cell refinement is also possible using the program, APPLE, and determination of crystallite size by the program CRYSIZ.

These x-ray analyses are nondestructive in nature and the sample is available afterwards for analysis by other techniques.

The Siemens D-500 system consists of a KRISTALLOFLEX 805 generator, thallium doped NaI detector with graphite monochromator, copper targeted x-ray tube (molybdenum optional), DIFFRAC V software package with JCPDS data base, Tektronix 4027A graphics terminal, DEC LA-120 printer, Hewlett-Packard 7220C graphics plotter, and DEC PDP 11/23+ computer for control of the 40 position sample changer and the goniometer.

Routine and emergency service is provided by Siemens under a yearly service contract. The following procedure will check on the alignment of the system optics. For the "five fingers of quartz" XRD scan, analyze the novaculite standard under the following conditions and compare to previous scans.

Slit Number	Degrees	KV	MA	Det. Volt
I	1	40	30	1050
II III	1			= 1 second
IV V	.018			ed = 1/5 per min 2 cm/minute

Scale = 1×10^3

Start the scan at 67.5 degrees and end at 68.7 degrees.

An additional check on the overall alignment may be performed by scanning the NBS silicon standard at .01 degree two-theta steps and counting for one second at each step. The program ADR is then run to locate the peaks and determine their intensities. The results of the scans are compared to previous scans.

Slit Number	Degrees	KV	MA		
I	1	50	30		
II	1				
III	1	Time	constant	1	second
IV	.05				
V	.15				

Analyze the novaculite standard from 20.0 through 30.0 degrees two-theta and run the program ADR to locate the peaks and determine their intensities. Compare to previous scans.

Slit Number	Degrees	KV	MA		
I	1	50	30		
II	1				
III	1	Time	constant	= 1	second
IV	.05				
V	.15				

b. X-Ray Fluorescence (XRF) (Siemens SRS-300)

X-ray fluorescence is a rapid qualitative and quantitative analytical method for the determination of elemental concentrations. It is useful over a wide range of concentrations and for a wide variety of samples. The Siemens SRS-300 automated sequential x-ray fluorescence unit is capable of detecting elements from boron (Z = 5) through uranium (Z = 92), but is presently configured to analyze from sodium (Z = 11) on up.

Under the control of the DEC PDP 11/23+ computer, up to 10 samples can be scanned under a wide range of operating parameters. The atmosphere of the sample chamber, the voltages, analyzing crystals, detectors, counting times, pulse height windows, soller slits, deadtime correction, filters, mask sizes, and sample rotation are under control of the computer. These can be varied to provide optimum analyzing conditions for each element and each matrix. Sample preparation is usually minimal, with grinding and pelletizing handling the majority of the samples submitted. The method is nondestructive; however, the sample usually has a binder added to it as an aid in pelletizing.

Routine and emergency service is provided by Siemens under a yearly service contract.

Presently, a synthetic ilmenite standard is being used as a check on the instrument's performance. A more durable standard would be beneficial for the long term. A polished metal alloy of known composition is something we hope to obtain shortly for use in verifing the instrument's operating condition.

1

.

E. Absorption of Ultraviolet, Visible and Infrared Light

a. Perkin-Elmer (PE) Model 1750 FT-IR

Instrumentation has been described which utilizes the emission of light for analysis. The converse, the absorption of light, is more adaptable for the analysis of ionic and molecular species in solutions, gases, and some solids, and it can also be used for the analysis of metal atoms in the atomic or ground state.

Fourier-Transform Infrared Spectroscopy (FT-IR) is a light absorption technique that employs the lower energy region of that part of the electromagnetic spectrum referred to as "light". Samples in solid, liquid or gaseous form are subjected to infrared radiation in the 2.5-25 micrometer wavelength range. When the energy of the radiation coincides with the proper energy to raise a molecule to an excited molecular state, the radiation is absorbed. The absorption energy is characteristic of the vibration (bending), stretching or rotation mode of the molecule. Infrared is principally used for the investigation of organic molecules.

The Varian Superscan I-BE is a recording ultraviolet (UV) - visible b. (Vis) spectrophotometer. Its range is from 190 to 900 nanometers (nm). Ultraviolet-Vis absorption spectroscopy employs the middle and higher energy region of that part of the electromagnetic spectrum referred to as "light". Samples in either solid, liquid, or gaseous form are subjected to UV or Vis radiation. When the energy of the radiation coincides with the energy of an electronic transition, an electron is promoted to a higher energy level and radiation is absorbed. The electronic may be either molecular (i.e., promoting an electron normally thought of as being associated with a chemical bond such as the UV absorption of aromatic molecules) or atomic (i.e., promoting an electron normally thought of as being associated with an atom such as d-d transition for "transition-series" elements in solution). Ultraviolet-Vis spectroscopy is principally used for the investigation of inorganic ions and inorganic complexes of "transition-series metals".

Other instruments available of this type, but less flexible, are the Beckman Models DU and DB.

c. A Jarrell-Ash Fluorimeter Model 26-000 is available for the determination of uranium in solid samples. Uranium in a mixed sodium-lithium fluoride melt will fluoresce under the influence of near ultraviolet light. The amount of fluorescence is dependent upon the quantity of uranium-VI present. Typical sensitivities are on the order of 1 to 5 parts per billion uranium.

d. Calibration Procedure

The spectrophotometers are situated in several areas of the KMTC and, at times, are moved to other locations. They are to be maintained as follows:

- Maintenance and repair are to be performed by the manufacturer as required or at their direction.
- Efficiency and/or wavelength calibrations, as required, will be performed prior to analysis of each sample set.
- Wavelength calibrations will be checked at least annually.

F. Microscopy

1. Optical Microscopy

a. Petrographic Microscope

This microscope (Zeiss) is used for a variety of analyses. It is especially useful for studying rocks and minerals. Minerals can be identified in transmitted light by measuring optical properties such as indices of refraction. Thin sections of rocks may be studied to determine the relationship of ore minerals to gangue, types and degree of rock alteration, and other information useful in projects related to minerals exploration, ore recovery, and ore processing. This microscope also is used for measuring the particle size of powders, determining shapes of particles, degree of agglomeration, etc. Attachments are available for studying opaque materials in reflected light. Photomicroscopy accessories are also available.

b. Metallographic (Unitron)

The metallograph is used for studying metals, certain ceramics, and other opaque materials in reflected light. It can be used for identifying opaque ore minerals, measuring the grain size of metals and ceramics, studying inclusions in metals, determining the size and distribution of pores in ceramics, aiding in corrosion studies, etc. The metallograph is equipped with photomicroscopy accessories. Polishing equipment is available for preparing mounts.

2. Electron Microscopy

The Phillips EM-300 Electron Microscope is used to view details of interest that are too small to be seen with an optical microscope. The optical microscope theoretically resolves details down to 0.2 micrometers, but in actual practice very little information can be obtained in details less than 0.5 micrometers. The electron microscope is capable of resolving details down to about 0.5 nanometers. It is used to study fine-grained materials such as TiO₂ pigments, UO₂ powders, and MnO₂ powders. Small surface details on large size materials also can be examined. Particles less than 0.5 micrometers can be measured precise y by using a Zeiss particle-size analyzer in conjunction with photomicrographs of the materials of interest.

The electron microscope is to be calibrated annually as follows:

Take a series of photographs of a standard line grid (54,864 lines/inch) at taps 1, 4, 7, 10, 13, 16 and 19. Compare photographs with photos taken from previous calibration checks under identical conditions.

G. Instrumentation for Analysis by Electrochemical Means

Instruments available in our laboratory which employ electrochemical principles are:

a. Sargent Electrodeposition Unit and Mercury Cathode

This instrument allows electrodeposition of metals such as copper, lead and zinc on platinum electrodes. The deposited metal is weighed and calculated as percent of sample. The mercury cathode is used to separate unwanted metals such as Fe, Cu, Ni, Cr, Mn, Mo and V. It has been instrumental in our development of a referee method for the determination of vanadium in ferrophosphorus material.

b. pH Meters and Specific Ion Meters

Ion concentrations (activities) are measured in solutions by employing a potentiometer to detect potential differences between a glass or specific ion electrode and reference electrode. Several meters are available.

c. Recording Titrators, Mettler and Metrohm/Brinkman

These instruments were designed to record and control all practical electrometric titrations using either glass or metal indicating electrodes.

Because their primary function is to record the voltage between the electrodes, the electrodes and the chemical system determine the limitations of their use. The second function is to control the rate of addition of titrant to the cell, thereby avoiding reaction rate errors on lengthy titrations.

H. Instrumentation for the Analysis of Gases and Liquids by Gas Chromatography, High Pressure Liquid Chromatography and Ion Chromatography

Gas chromatography (GC) is useful for the separation and analysis of gases and liquids with boiling points up to 450°C. Organics with boiling points above 450°C and thermally degradable species are analyzed by high pressure liquid chromatography (HPLC). Both techniques employ carrier fluids that are used to elute organic sample components through

the columns at different rates, depending on how they partition between the stationary phase and carrier. Components are detected by various means as they emerge from the column. The detected components are displayed as peaks on a potentiometric strip chart recorder.

In GC, the carrier is a gas such as helium and samples are eluted as gases (liquids may be converted to gases in a heated compartment prior to entering the column). In HPLC, the carrier is a liquid such as methanol and components are eluted as liquids.

In GC, temperature programming is used to permit the analysis of both heavy and light components in a single sample. In HPLC, the temperature is usually low and constant and changes in column material and column polarity are used to effect the separation of heavy and light components in liquid samples.

a. Gas Chromatographs available are as follows:

Manufacturer	Model #	Thermal Conduct. Detector (TCD)	Flame Ionization Detector (FID)
Carle	SX-AGC	Yes	No
Hewlett-Packard	5830	No	Yes
Hewlett-Packard	5840	No	Yes

They are to be maintained and calibrated as follows:

- Maintenance and repair are to be performed by the manufacturer or according to his specifications as required.
- Standards are to be run prior to each sample set being analyzed to calibrate the instrument.
- b. High Pressue Liquid Chromatograph

One HPLC instrument is available; namely, the <u>Waters Model</u> <u>ALC-202/401</u> equipped with solvent programmer, refractive index and ultraviolet detectors and using the Fisher Recordall 2-pen recorder. It is to be maintained and calibrated as follows:

- Maintenance and repair are to be performed by the manufacturer or according to his specifications as required.
- Standards are to be run prior to each sample set being analyzed to calibrate the instrument.

c. Ion Chromatograph

A Dionex Model 16 ion chromatograph, capable of simulataneous anion and cation analysis, is available.

This technology, discovered by Branscomb at the National Bureau of Standards (1968), developed and reduced to practice by Dow Chemical Company (1975), is analogous to gas and high pressure liquid chromatography and is applicable to aqueous and select nonaqueous systems. Ions are separated, chromatographically, using proprietary ion-exchange resins; the eluted ions are then quantized by a specific conductance meter. Using this technique, most cations of moderate to strong bases or anions of moderate to strong acids may be determined, individually or in mixtures.

The ion chromatograph is to be maintained and calibrated as follows:

- Maintenance and repair are to be performed by the manufacturer or according to his specifications as required.
- Standards are to be run prior to each sample set being analyzed but at least once every eight hours.
- d. A <u>Perkin-Elmer Model 1 Integrator</u> and a Beckman 10" strip chart recorder are used to integrate peaks and trace chromatograms respectively.
- e. Two <u>Hewlett-Packard Model 3390A Integrators</u> are used to record and integrate peaks from gas, liquid, and ion chromatographs and provide trace chromatograms.

I. Organic Analysis

a. Coulometrics Carbon Analyzer

This method allows the determination of total inorganic carbon and/or total carbon and/or total organic carbon in aqueous and inorganic or organic solid samples. The carbon is oxidized to CO_2 or inorganic carbon liberated as CO_2 which is then measured coulometrically.

b. Ramsbottom Carbon Residue Analyzer

A <u>Ramsbottom Carbon Residue Analyzer</u> is available for determination of the amount of carbon residue remaining after evaporation and pyrolysis of an oil.

J. Muffle Furnaces

An assortment of muffle and tube furnaces is available.

K. Surface Area Analyzer

The Quantachrome Monosorb Model MS-8 is used to measure the total surface area of fine powders such as TiO₂, MnO₂, UO₂, etc. It performs this function by measuring the amount of nitrogen that is adsorbed on the surface. It is useful for characterizing materials and provides a means of establishing quality control on a variety of processed materials.

The instrument is to be calibrated quarterly as follows:

Analyze commercially-prepared standards and compare data to accepted value. Inspect sample tubes for chips and/or cracks, and check heating mantles for proper operation.

L. Miscellaneous Equipment

- a. A variety of small instruments, such as a Frantz Isodynamic Separator, Andreason Pipettes, Sieves, etc. are used to aid in phase separations, measure particle size, measure bulk densities and true density, etc. Many of the techniques used for characterizing materials do not require elaborate equipment.
- b. The equipment available includes a 4" Jaw Crusher, 6" Jaw Crusher. Cone Crusher, BiCO Disk Pulverizer, Rotap Screen Shaker, a complete complement of U. S. Standard Screens, assorted riffles and other equipment necessary for the preparation of most types of mineral ores and coal for laboratory analysis.

c. Cleveland Open Cup and Pensky-Martins Closed Cup Tester

This equipment is used to determine the flash and fire points of petroleum.

d. Cameras (Polaroid MP-3, Crown-Graphic, Ecc.)

Various cameras and related equipment are available for macrophotography and photomicroscopy. Photographs can be made for permanent records. Anything from coarse cores and large pieces of equipment down to the fine microscopic features on a variety of materials can be photographed.

e. Mechrolab Vapor Pressure Osmometer

The Mechrolab Vapor Pressure Osmometer is used for the determination of the molecular weight of organic molecules.

f. Brookfield Model HAT Viscometer

The <u>Brookfield Model HAT Viscometer</u> is used for measuring viscosities between 1-10 centipoise. Also available are <u>Cannon-Fenske Viscometers</u> for measuring liquids between 1-20,000 centipoise. The temperature range for viscosity measurements is 0-200°C.

The instrument is to be calibrated quarterly as follows:

Analyze DC-200 Silicone Fluid at several temperatures up to 500°F and compare data with that obtained from previous analysis under identical conditions. Also, check instrument at room temperature by analyzing 900 cp and 11,000 cp standards.

g. pH Meters

These instruments are used for pH and millivolt determination. Standards or buffers are to be run prior to each sample set being analyzed to calibrate the instrument. In addition, the millivolt scale will be tested annually using a known potential source.

h. Specific Conductivity Meter - YSI-31

This instrument is used for specific conductivity determination. Standards are to be run prior to each sample set being analyzed.

i. Balańces

Analytical balances are to be tested and certified biannually by an outside consultant and given internal checks annually.

M. Quality Assurance - Quality Control

A quality assurance - quality control program is maintained using the principles outlined in Quality Assurance Principles for Analytical Laboratories, Frederick M. Garfield, Association of Official Analytical Chemists, 1984.

a. Definition

Quality Control is defined as a planned system of activities whose purpose is to provide a quality product. Quality Assurance is defined as a planned system of activites whose purpose is to provide assurance that the quality control program is actually effective. For a quality assurance program to be effective, there must be a plan.

b. Quality Assurance Plan

In reality, a quality assurance plan is composed of three essential components, and costs are associated with each component:

- Prevention
- Assessment or appraisal
- Correction

<u>Prevention</u> requires an orderly program of planning and positive actions before or during analyses to ensure that analytical systems are functioning properly. Examples are quality control planning, training, calibration of instruments, instrument maintenance, and frequent standardization of "standard solutions", etc.

Assessment is a form of control that includes periodic checks on performance to determine precision and accuracy. Examples include analysis of duplicate and check samples, peer check of chart readings and calculations, and validation of methodology.

Correction is action taken to determine causes of quality defects and to restore proper functioning of the analytical system. This can involve trouble shooting to correct malfunctioning equipment, examination of more check samples, re-evaluations of methodology, retraining, etc.

KERR-MCGEE CORPORATION	KMTC - 121-F-2		
TECHNICAL CENTER	ISSUED7-1-74	REVISED76	
ANALYTICAL METHODS		PAGEOF3	
DETERMINATION OF	TTRE OF SAMPLE		
Fluoride	Ilmenite		
ANALYTICAL NETHOD	WRITTEN BY	APPROVED	
Specific Ion Electrode	JRJ/GEV	CHL	

Method:

Fusion of ilmenite in Na₂CO₃, separation of fluoride from impurities by distillation and, then, fluoride determination using specific-ion electrode.

Apparatus:

- 1. Orion 9409A Fluoride Electrode
- 2. Calomel or Silver/Silver Chloride Reference Electrode
- 3. Metrohm pH Meter Model 102 (Brinkmann) or Equivalent
- 4. Metrohm Single Probe pH Electrode or Equivalent
- 5. Platinum 50 ml Size Crucible with Cover
- 6. Plastic Beakers, 150 ml, 400 ml
- 7. 250 ml Polyethylene Volumetric Flasks
- 8. Fisher Burners
- 9. Distillation Apparatus Consisting of a 500 ml Round-Bottom Flask, Graham Coil-Type Glass Condenser, Distillation Arm, Powerstat, heating Mantle.

Reagents:

- 1. Sodium Carbonate (Na2CO3) Reagent Grade
- 2. Sulfuric Acid Solution (1:1)
- 3. Acetic Acid Reagent Grade
- 4. Sodium Chloride Reagent Grade
- 5. 10N Sodium Hydroxide
- 6. Buffer Solution Dissolve 58 g NaCl, 2 g CDTA (1,2-cyclohexylenedinitrilotetraacetic acid) in 500 ml fluoride-free deionized water. Add 57 ml conc. acetic acid. Adjust pH of solution to between 5.0 and 5.5 with 10N NaOH. Dilute to one liter.
- Fluoride Standard Stock Solution: Dissolve 2.211 gm of NaF (dried @ 105°C) in water and dilute to 1 liter. 1 ml = 1000 µgF.

DETERMINATION OF		
Fluoride		KMTC - 121-F-2
ANALYTICAL METHOD	TYPE OF SAMPLE	
Specific Ion Electrode	Ilmenite	Poge 2 3

- 8. Dilute 100 ml of stock fluoride solution to 1 liter. 1 ml = 100 µgF.
- 9. Dilute 10.0 ml of stock fluoride solution to 1 liter. 1 ml = 10 µgF.
- Prepare 0, 1, 5, 10, 100 and 1000 µgF fluoride standards by using 50 ml of buffer and diluting to 100 ml with distilled water.

Procedure:

- 1. Weigh 0.5 g of 100% pass 100-mesh sample into a 50 ml platinum crucible. Add 10 g Na_CO3 and mix with metal stirring rod. $^{\rm l}$
- Fuse at red heat for at least 30 minutes (1 hour preferred) over Fisher burner with cover and, then, cool.
- Leach the Na₂CO₃ fused cake with water until dissolved, transferring leach to a 100 ml plastic beaker, dilute to 100 ml volume.
- 4. Transfer to a 500 ml distillation flask.
- Neutralize sample (slowly) with 1 + 1 sulfuric acid, add a total of 150 ml of 1 + 1 sulfuric acid. Add another 50 ml of distilled water.
- Connect the distillation flask to the distilling column and distill until SO₃ fumes. (Distillate volume should be about 225 ml). Rinse condenser with several milliliters of distilled water and add rinses to the distillate.
- 7. Transfer distillate to a 250 ml volumetric flask, dilute to volume. Mix well.
- Transfer 50 ml of distillate to a 150 ml plastic beaker, add 50 ml of buffer (reagent No. 6).
- 9. Insert fluoride electrode into solution.
- 10. Wait 5 minutes or until equilibration is reached (electrode potential stabilizes).
- 11. Read mv, call this (S).
- 12. Add 50 ug fluoride (0.5 ml of 100 ug/mlF, reagent No. 8).
- 13. Allow to equilibrate.
- 14. Read my, call this (A).
- 15. Add 50 µg fluoride (0.5 ml of 100 µg/mlF, reagent No. 8).
- 16. Allow to equilibrate.
- 17. Read mv, call this (B).
- 18. Plot µgF vs. mv on semilogarithmic graph paper²,³.

DETERMINATION OF		
Fluoride		KMTC - 121-F-2
ANALYTICAL METHOD	TYPE OF SAMPLE	
Specific Ion Electrode	Ilmenite	Page of 3

19. From graph, record ug of S, A, and B.

20. Calculations:

- F, $\mu g/g = \frac{\mu g \text{ of } S \times \mu g \text{ of } F \text{ added } \times 250 \text{ ml}}{[\mu g \text{ of } (A + B) \mu g \text{ of } S] \times 0.5 \text{ g x } 50 \text{ ml}}$
- Note: ¹A blank of all reagents should be carried through the procedure for a blank dedetermination.

²The concentration of fluoride is plotted on the log axis and the electrode potential developed with the standard on the linear axis.

³The dilution-error from the standard addition results in a potential error of 1% in the resulting potential obtained. This error may be considered negligible.

Precision studies on four different ilmenites, for a total of fifteen determinations, indicated by averaging the fifteen relative percent standard deviations, the precision in the 30-35 ppm range is \pm 7.8% of the value. Accordingly, it is recommended that at 50 ppm F, two values from two different determinations should be considered in agreement if the average of the two lies within plus or minus 4 ppm and that the average be considered the assay of that sample if the sample is prepared and split in an approved manner. For a 50 ppm F standard, the value of plus or minus 4 ppm represents one standard deviation; which is to say, if a 50 ppm F absolute standard was assayed, a large number of times, 67% of the assays would fall within the range of one standard deviation or from 46 to 54 ppm (50 \pm 4 ppm), 95% of the assays would fall within the range of two standard deviations or from 42 to 58 ppm (50 \pm 8 ppm), and 99 \pm % of the assays would fall within the range of three standard deviations or from 38 to 62 ppm (50 \pm 12 ppm).

Manpower Requirements:

One sample of ilmenite, run in duplicate, from a split of the original, requires eight clock hours and five manhours of labor. This assessment assumes a dual distillation set-up which may be run simultaneously. Each additional sample, to be run in duplicate, will require an additional three-and-one-half manhours.

KERR-MCGEE CORPORATION	KMTC - 76-U-14		
KERR-MCGEE CORPORATION	ISSUED 1-14-71	REVISED	
ANALYTICAL METHODS		PAGEOF3	
DETERMINATION OF	TYPE OF SAMPLE		
Uranium	Urine		
ANALYTICAL METHOD	WRITTEN BY	APPROVED	
Fluorimetric	DMK/GMFIII	CHL	

Method

An aliquot of urine is deposited on a fluorimetric flux mixture (2% lithium in fluoride - 98% Sodium Fluoride) in platinum dishes. The pellets are dried, fused, and allowed to cool, forming buttons which slide easily from the dish. The fluoresence of the fluoride button is measured in a Jarrel-Ash fluorimeter, compared with a standard uranium curve, and reported as micrograms of uranium per liter of urine.

Limits of the Method

The fluorimeter will measure from one to 10,000 micrograms of uranium per liter of Urine. Background fluoresence reduces the accuracy of the analysis at very low concentrations (1-10µg per liter). Above 10µg per liter, the accuracy is within 10%. High concentrations (over 1000µg) are usually diluted prior to analysis to prevent saturation of the fluoresence in the pellet.

Cost

The labor for the determination of a set of 20 pellets would be 40 minutes. The number of samples, depend upon the number of standards and blank analyzed per set of 20 pellets.

Equipment and Reagents

- 1. Jarrel-Ash Galvenck Morrison Fluorimeter, Newtonville, Massachusetts.
- Fusion Burner Modified Fletcher Burner w/Nichrome gauze to accomodate 20 pellets.
- Platinum dishes Stamped from discs 0.015 inch thick and 0.75 inch in diameter purchased from J. T. Baker and Company, Newark, New Jersey.
- 4. Pellet Maker Custom fabricated from Stainless steel.
- 5. Timer for fluorimetric burner (TDAF)
- 6. Propane Regulator with fluorinator and needle valve.
- 7. Air regulator Air blower assemble w/booster pump.
- 8. Heat lamp.
- 9. Die for forming platinum fusion dishes.

DETERMINATION OF		
Uranium		
ANALYTICAL METHOD	TYPE OF SAMPLE	KMTC -76-U-14
Fluorimetric	Urine	Page 2 of 3

- 10. Sodium Fluoride, reagent grade. Meets fluorimetric specifications, lithium fluoride flux. A reagent jar containing 98% sodium fluoride and 2% lithium fluoride powders is agitated on a shaker and or dilutable rolling mill for a period of 36 hours. Aliquots of the mixture is sampled until a suitable blank is determined or the entire procedure repeated.
- 11. Uranium Standard Solution (chemically pure, NBS certified, Natural uranium).

Stock Solution "A" 500ug/ml

Dissolve 58.9 mg of NBS certified U_3O_8 in 2 ml conc HNO₃ and evaporate to dryness. Take up with 2 ml water containing 10 drops of conc HNO₃, transfer to a 100 ml Pyrex volumetric flask and dilute to volume.

Standard Solution "B" 0.5µg/ml

Transfer 0.1 ml of solution "A" to a 50 ml pyrex volumetric flask, add 5 drops of conc HNO₃ and dilute to volume. The standard solutions should be stored in polyethylene bottles.

Procedure

Sample Preparation

- 1. A clear sample may be used as received. The addition of 1% by volume of HCl at the time of sampling will prevent deterioration.
- If the sample is cloudy, add 1% by volume of conc HC1 and let stand overnight.
- 3. If the sample is clear, proceed with the analysis.
- 4. If the sample is still not clear, add 15 ml of conc HNO₃ and 2 ml of 30% H₂O₂. Carefully evaporate to dryness, cool, and make up to the original volume with distilled water.

Determination

- 1. Adjust the pellet maker to accommodate 400 mg of flux and add this amount onto each platinum disc.
- 2. Pipette 0.2 ml of urine onto each pellet. Using duplicate aliquots per sample.
- 3. Place samples under heat lamp and dry for 10 minutes.
- 4. Ignite the residue over the fusion burner using the following method.
 - a. Set the timer switch for 3 min., and ignite the burner. (switch on timer)b. Adjust the volumeter setting to 18, this will approximate a pressure
 - allowing flame temperature to fuse the flux at 900°C.
 - c. After two minutes, adjust the volumeter to 21, decreasing the temperature 850°C.
 - d. Allow the fusion to continue through the time period. The timer will cut off the air and the gas simultaneously.

DETERMINATION OF		
Uranium		KUTO TO THE
ANALYTICAL METHOD	TYPE OF SAMPLE	KMTC - 76-U-14
Fluorimetric	Urine	Page 3 of 3

5. When the fluoride melts have sole fied, remove the gauze and dishes from the burner and allow to cool for at least 15 minutes.

Measurement of Fluorescence

- 1. Push the sample slide to the forward stop, so that the empty aluminum holder is in the read position.
- 2. Adjust the 0.1 key to zero on the meter scale, with the background control.
- 3. Adjust the 0.1 key to zero with the fine adjustment. (Use a screwdriver).
- Pull the slide to the Front stop, depress the 0.1 key and check for zero, using the background control.
- 5. Invert each platinum dish and allow the fused button to slide out of the pan.
- 6. Using forceps, carefully arrange blanks, standards, and samples in preparation for measurement.
- Place each sample in the aluminum holder, insert the slide to the forward or read position, depress the scale 0.1 or 1, and read the scale deflection, recording for calculation.

Calculations

From the following equation determine micrograms of U or $U_3 O_8$ in the unknown urine sample.

Micrograms of U in unknown: µg of U308 in standard x (ave unknown - blank) ave standard - blank

Cleaning of Platinum Ware

- 1. Rinse in hot distilled H20.
- 2. Let stand overnight in 9N HNO3, and fuse with fluoride flux.
- 3. Fuse with potassium bisulfate, wash with tap water and soak in $9N \ HNO_3$ as above.
- 4. Rinse with distilled water and dry over flame.
- 5. Fuse with fluroide flux on each dish. Usual fusing time.
- Select 3 dishes and observe the reading on the meter scale. If any show a reading above a normal blank level the entire batch is re-cleaned, and three random samples again tested.

References

1. Fluorimetric Determinations of Uranium

Frederick A. Centanni, Arthur M. Ross, and Michael A. DeSesa, Raw Materials Development Laboratory, National Lead Co., Inc. Winchester, Mass. Analytical Chemistry, Vol 28, No. 11, November 1956.

KERR-MCGEE CORPORATION	KMTC - 149-EC-10	
KERR-MCGEE CORPORATION	ISSUED 12/1/77	None None
ANALYTICAL METHODS		AGEOF8
Thorium, Uranium and Plutonium	TYPE OF SAMPLE Water, So Products	oil, Vegetation and or By-Products
Alpha Pulse Height Analysis	GEVS	CHL

Method:

Except for natural uranium or natural thorium, environmentally important concentrations of actinides cannot be readily determined by ordinary chemical means. Additionally, chemical methods do not offer any information as to the isotopic composition of the element in question. Therefore, for samples where either an isotopic composition and/ or ultra low level of detection is required, alpha pulse height analysis is the method of choice.

This method utilizes the following technologies:

- In-situ formation of calcium phosphate is used to concentrate the actinides from a large aqueous volume.
- Anion exchange is used to permit separation of the actinides into relatively pure form.
 - a) For thorium, good separation is achieved from all elements except uranium. A separation factor in the vicinity of 10⁶ is achieved for most elements, but only about 10⁴ separation is achieved between uranium and thorium. As a result, some uranium will report to the electrodeposit when analyzing for thorium in a uranium rich material such as yellow cake.
 - b) For uranium, good separation is achieved from all elements except thorium and those which form stable anion complexes with chloride (i.e., iron, aluminum, etc.). For this reason, whenever a material rich in these metals is to be analyzed for uranium, such as soils or mill tailings, the uranium must first be solvent-extracted from a nitrate or sulfate system in order to eliminate these interferences.
 - c) For plutonium, good separation is achieved from all metals except for iron from iron-rich systems (i.e., liver tissue, iron or steel and their alloys, etc.). When analyzing an iron-rich system, the ion-exchange step should be repeated once. (Note: Three evaporations to dryness with concentrated nitric acid area required to completely remove traces of chloride and fluoride from the first ion-exchange strip liquor prior to the second ionexchange. Failure to remove the chloride or fluoride will result in poor recovery of the plutonium.).
- 3. Electrodeposition is utilized to prepare an "infinitesimally thin" deposit necessary for good alpha pulse height energy resolution. A thick deposit (> ten micrograms) causes alphas, from the bottom of the deposit, to degrade (lose energy) as they pass up through the deposit. As a result, the silicon surface barrier detector will see alphas with a lower energy which will produce

DETERMINATION OF		T
Thorium, Uranium and Plutonium		WHITE HIS SS
ANALYTICAL METHOD	Type of sample Water, Soil, Vegetation	KMTC - 149-EC-10
Alpha Pulse Height Analysis	and Products or By-Products	Page of 8

a spectral smear. Poor resolution of alpha energy will destroy the usefulness of the procedure. Alpha resolution should be < 50 KeV FWHM (full width at half maximum) and any spectra with an alpha resolution of > 100 KeV FWHM should result in rejection of the analysis.

Range and Precision:

Amounts of Th, U, or Pu as low as one femtocurie $(1 \times 10^{-15} \text{ curie or } 3.7 \times 10^{-5} \text{ disintegrations per second})$ can be detected over a counting time of 1000 minutes using silicon surface barrier alpha detectors in conjunction with a multichannel analyzer. The precision attained is a function of the activity in the sample and the recovery of the internal standard; however, typical precision, at one standard deviation, ranges from 3 to 50%.

Cost:

Cost (manhours) = k + nt

Where:

k = 2 manhours; t = 2 manhours, and n = number of samples.

Apparatus:

 Silicon surface barrier detector(s) (Princeton Gamma Tech, 300-25-100) is used in conjunction with a multichannel analyzer.

2. Millipore filter equipment including 0.45µ Millipore Filter.

Magnetic stirring hotplate(s).

- Ion-exchange column(s), Bio-Rad or equivalent, 1.5 cm I.D. X 15 cm in length, and accessories.
- 5. Ripple-free D. C. power supply capable of at least 1.5 amp @ 10 volts.
- Electrodeposition cells (Talvite, Anal. Chem. 1972) and polished stainless steel discs.
- 7. Assorted glass and plastic ware common to analytical laboratories.

Reagents:

- 1. a) Concentrated sulfuric acid (H₂SO₄).
 - b) Concentrated hydrochloric acid (HCl).
 - c) Concentrated nitric acid (HNO3).
 - d) Concentrated hydrofluoric acid (HF).
 - e) Concentrated phosphoric acid (H_PO_4).

Thorium, Uranium and Plutonium KMTC - 149-EC-10			
TICA	a Pulse Height Analysis	TYPE OF SAMPLE Water, Soil, Vegetation	Page 3 of 8
api	a ruise neight maxyors	and Products or By-Products	
	f. Concentrated ammonium	hydroxide (NH,OH).	
	g. Concentrated perchlori	c acid (HClO4).	
2.	50% calcium nitrate soluti	on [Ca(NO ₃) ₂].	
3.	7.2F nitric acid		
4.	Bio-Rad (or Dowex) AG 1 X	4, 50-100 mesh, anion exchange resid	n.
5.	1.5F hydrochloric acid.		
6.	0.1F silver nitrate.		
7.	8F hydrochloric acid.		
8.	0.5F hydrochloric acid.		
9.	0.36F hydrochloric acid,	0.008F in hydrofluoric acid.	
10.	Thymol blue indicator.		
11.	0.15F ammonium hydroxide.		
12.	95% or 100% ethanol (not	denatured).	
13.	0.18F sulfuric acid (H_2SO)	u).	
14.	Sodium nitrite.		
15.	Standard Pu-242, U-232 or Th-228 solution. The Pu-242 (99.9999 + %), prepared by cyclotron and purified by mass spectrometry (Alexander and Irene Dupzik, Lawrence Livermore Laboratory) is available from John Harley, New York Dept. of Health.		
	brated "in-nouse."		
San	mple Preparation:		
Α.	Aqueous Solutions		
	1. Filter sample through	1 0.45µ Millipore filter.	
	2. For samples < 100 ml,	, use 100-250-ml beaker.	
	a) Add 25 ml conc. H	HNO3, heat to 75-80°C while stirring	. Helling the second

- b) Add the spike (Pu-242, U-232 or Th-228). Make sure the exact amount of spike added is known and recorded. Add only that spike necessary for the analysis. For environmental samples, spikes should be on the order of 2-12 dpm.
- c) Continue the heating for at least 30 minutes (60 min for Pu).
- d) Evaporate the solution to dryness. CAUTION, DO NOT BAKE THE RESIDUE.

OLTERNINATION OF		
Thorium, Uranium and Plutonium		KMTC - 149-EC-10
ANALYTICAL METHOD	TYPE OF SAMPLE Water, Soil, Vegitation	
Alpha Pulse Height Analysis	and Products or By-Products	Page of 8

- e) If insols begin to form when the volume is ≤ 10 ml, then, to the dry residue, add 25 ml 7.2F HNO₃, heat to 100°C, cool, filter through 0.45µ Millipore, wash residue twice with 10 ml 7.2F HNO₃, discard filter and any remaining residue, and, evaporate filtrate to dryness. CAUTION, DO NOT BAKE THE THE RESIDUE.
- f) Save the residue for actinide determination.
- 3. For samples > 100 ml (or, at the option of the analyst samples \leq 100 ml), use 1000-2000 ml beaker.
 - a) Where necessary, dilute to ~500 ml with double deionized, filtered water.
 - b) Add 40 ml conc. HNO3 make sure pH < 2. Stir sample.
 - c) Add 0.5 ml 50% Ca(NO3) 2 solution with stirring.
 - d) Add 1 ml conc. H3PO4 with stirring.
 - e) Continue the stirring and add the spike (Pu-242, U-232 or Th-228). Make sure the exact amount of spike added is known and recorded. Add only that spike necessary for the analysis.
 - f) Heat the sample to 75 to 80°C on a stirring hot-plate; continue the stirring.
 - g) Continue the heating for at least 30 minutes (60 min for Pu). Do steps h) and i) while sample is at 75-80°C.
 - h) <u>Slowly</u>, and with vigorous stirring, add 200 ml conc. NH₄OH, at a rate of ~ 10 ml/minute. Check pH; if pH < 9.5 continue adding NH₄OH until pH > 9.5. Turn off heat.
 - i) Stir for 30 minutes.
 - Remove magnetic stirrer; allow sample to stand for at least 30 minutes to further cool and settle.
 - k) When sample is cool enough to handle, filter through 0.45µ Millipore. Filter by initially decanting supernatant liquor from top of beaker until ~ 100 ml remains. Slurry the precipitate in the remaining 100 ml and pour into the filter.
 - If precipitate remains in the beaker, or on the walls of the filter, wash with 1.5 NH₄OH and rinse into the filter. Discard filtrate.
 - m) Transfer the precipitate to a 50-100-ml beaker and ash the sample at 350°C in an O₂ atmosphere for 35 ± 5 minutes. Remove sample from oven and cool to room temperature.
 - n) Save the precipitate for actinide determination.

DETERMINATION OF		
Thorium, Uranium and Plutonium		KMTC - 149-EC-10
ANALYTICAL METHOD	TYPE OF SAMPLE Water, Soil, Vegetation	
Alpha Pulse Height Analysis	and Products or By-Products	Pageof8

- B. Non-Aqueous Samples: (0.1-10 gms dry weight)
 - 1. Dry ash the sample in the high temperature oven.
 - a) Start sample at 300-350°C, no 0_2 , and maintain for ~ 2 hours.
 - b) Turn on O₂ while maintaining sample at 350°C for ∿ 2 hours.
 - c) Increase temperature to 500°C, 02 on, and heat for for at least 4 hours or overnight (~ 16 hours).
 - d) Use platinum dish, Pyrex, or Vycor beaker.
 - Transfer the cool, dry ashed, sample to a Pt dish, cover with conc. HNO₃ (but at least 10 ml added very slowly) and mix well. Caution, vigorous gas evolution (CO₂) may occur. Add 5 to 10 ml of conc. HF and mix well.
 - Add spike (Pu-242, U-232 or Th-228). Make sure exact amount of spike added is known and recorded. Add only that spike necessary.
 - 4. Heat to ~ 100 °C and slowly evaporate slurry to a moist residue (should take ~ 1 hour). Remove from heat, repeat the HNO₃ and HF addition and reheat to evaporate slurry to a moist residue. Repeat until the residue becomes a white material.
 - Add 5 to 10 ml conc. HClO₄, fume to wet residue (1-3 ml volume), cool. (If Fu is to be determined, do not add HClO₄.)
 - Add 25 ml 7.2F HNO₃ and filter through 0.45µ Millipore filter. Wash residue twice with 5 ml each of 7.2F HNO₃. Discard Residue.
 - 7. Save the filtrate for actinide determination.

C. Thorium Analysis:

- To the residue (step 2.f) or 3.n) for aqueous samples or filtrate from step 7 for non-aqueous samples):
 - a) Dissolve or dilute to 80 ml with 7.2F HNO3.
 - b) Pass the sample through a charged resin column. (Dowex 1 X 4 or BioRad 1 X 4, NO₃ form, C1⁻ free by AgNO₃ test). Flow rate should be 1 to 3 ml per minute. Discard eluate.
 - c) Wash the column with 150 to 250 ml 7.2F HNO3 at 1-3 ml per minute. Discard washings.
 - d) Elute the thorium from the column with 80 ml 1.5F HCl at 1-3 ml per minute.
 - e) Save the HCl eluate for electrodeposition.

DETERMINATION OF		
Thorium, Uranium and Plutonium		KMTC - 149-EC-10
ANALYTICAL NETHOD	TYPE OF BAMPLE Water, Soil, Vegetation	
Alpha Pulse Height Analysis	and Products or By-Products	Page 0 of 8

- D. Uranium Analysis:
 - To the residue (step 2.f) or 3.n) for aqueous samples or step 7 for non-aqueous samples evaporated to moist residue), if Fe (III) content is > 1 ppm, proceed as follows:
 - a) Dissolve residue into 10 ml 2F HNO3/2F Al(NO3)3.
 - b) Extract uranium twice into 10 ml ethylacetate.
 - c) Transfer ethylacetate (20 ml) to clear, dry, beaker and evaporate to dryness. CAUTION, DO NOT BAKE RESIDUE.
 - 2. If < 1 ppm, proceed to next step.
 - a) Dissolve the residue in 80 ml 8F HCl.
 - b) Pass the sample through a charged resin column. (Dowex 1 X 4 or BioRad 1 X 4, Cl⁻form). The flow rate should be 1-3 ml per minute. Discard eluate.
 - c) Wash the column with 150 to 250 ml 8F HC1. Discard washings. Use flow rate of 1-3 ml per minute.
 - d) Elute uranium from column with 80 ml 0.5F HCl at a flow rate of 1-3 ml per minute.
 - e) Save the HCl eluate for electrodeposition.
- E. Plutonium Analysis:
 - To the residue (step 2.f) or 3.n) for aqueous samples or filtrate from step 7 for non-aqueous samples):
 - a) Dissolve residue in 80 ml. 7.2F HNO3 or dilute to 80 ml with 7.2F HNO3, and .5 g boric acid (H3BO3).
 - b) Add several milligrams NaNO₂, heat to 90°C, cool and let stand overnight (or heat to, and maintain at, 100°C for at least 1 hour).
 - c) Pass the sample through a charged resin column. (Dowex 1 X 4 or BioRad 1 X 4, NO, form, Cl free by AgNO, test). Flow rate should be 1-3 ml per minute. Discard eluate.
 - d) Wash column with 150-250 ml 7.2F HNO₃ at a flow rate of 1-3 ml per minute. Discard washings.
 - e) Elute plutonium from column with 80 ml of 0.36F HC1/0.008F HF at a flow rate of 1-3 ml per minute.
 - f) Save the HC1/HF eluate for electrodeposition.

Thorium, Uranium and Plutonium		KMTC - 149-EC-10
	TYPE OF SAMPLE Water, Soil, Vegetation	149-20-10
Alpha Pulse Height Analysis	and Products or By-Products	Page 7 of 8

F. Electrodeposition:

- To the eluates from Th, U, or Pu procedures (step C.1.e) for Th, step D.1.e) for U, and step E.1.f) for Pu):
 - a) Add 0.30 ml conc. H2SO4 and 5 ml conc. HNO3.
 - b) Evaporate to dense SO, fumes and cool.
 - c) Add 5.0 ml H₂O.
 - d) Add thymol blue indicator.
 - e) Adjust pH to 2 (first color change) with gaseous NH₃. Do not overshoot pH of 2. Obtain NH₃ by cutting inside stem off of polyethylene wash bottle flush with cap. Fill wash bottle 1/2 full with conc. NH₃OH. "Puff" NH₃ vapors from bottle. Direct NH₃ stream at, but above, surface of liquid sample.
 - f) Transfer to electrodeposition cell. Have clean, degreased stainless steel planchet in place in the cell.
 - g) Add 5 ml 0.18F H2SO4 use for beaker rinse.
 - h) Repeat step g) at once.
 - Adjust to pH of 2 (first color change) with gaseous NH₃. Do not overshoot pH of 2.
 - j) Place cell in electrodeposition apparatus and plate actinide onto the cathode (stainless steel planchet) for 1.2 amp/hour (0.6 amps for 2 hours or 1.2 amps for 1 hour).
 - k) With current flowing, add 5 ml conc. NH40H and continue current flow for 1-2 minutes.
 - With current still flowing, remove platinum anode from solution. Turn off current and disconnect cell.
 - m) Discard aqueous contents of cell.
 - wash cell three times with 0.15F NH4OH. Do not direct stream onto plated surface.
 - o) Wash cell twice with 95% ethanol.
 - p) Disassemble cell and carefully remove planchet. Do not touch plated surface.
 - q) Place planchet on hotplate at 250-300°C and heat for 5-10 minutes.
 - r) Remove planchet from hotplate and cool to room temperature.

DETERMINATION OF		1
Thorium, Uranium and Plutonium		KMTC - 149-EC-10
ANALYTICAL NETHOD	TYPE OF SAMPLE Water, Soil, Vegetation	147-20-10
Alpha Pulse Height Analysis	and Products or By-Products	Poge 8 of 8

G. Alpha Pulse Height Analysis:

- 1. Active area of surface barrier detector should be at least 2X's that of the plated surface of the planchet.
- Place planchet in close proximity to, but not touching, surface barrier detector. The surface barrier detector and sample must be inside a vacuum chamber with a pressure < lmm Hg.
- 3. Count the sample, over the energy range 3.8 to 5.8 MeV, for 1000 minutes.
- 4. Integrate the alpha peaks and identify each isotope by its alpha energy.
- 5. For the internal standard, the recovery is determined by the equation

counts per 1000 min found = Recovery Factor

- Divide each alpha integral by the recovery factor and divide the quotient by 1000. The result is the activity of that alpha (isotope) in the sample in units of disintegrations per minute (dpm).
- 7. To convert to picocuries per sample: dpm = 2.22 = pCi To convert to femtocuries per sample: dpm = .00222 = fCi
- 8. In alpha pulse height analysis, using silicon surface barrier detectors, the sample size should be chosen such that the alpha activity, on the planchet, exclusive of the internal standard, lies in the range of from 1.5 to 15 pCi. However, for low level samples it usually is impractical to exceed 1 liter or 10 grams of sample and, as a result, one often finds the alpha activity, exclusive of the internal standard, lies in the range of from < 0.001 to < 1.5 pCi.</p>
- For samples where the recovery of the internal standard is < 50%, based upon a measured counting geometry, the assay should be rejected.

KMITC - 200-XRE-1	
139000	PAGE OF
Solids and Lig	uids
WRITTEN BY	GEV
	Solids and Lig

METHOD:

Primary x-rays are generated by the bombardment of a suitable target material by electrons from a heated tungsten filament. These electrons are accelerated to a high velocity by a high voltage potential between the tungsten filament and the target material.

As the electrons deccelerate, they give up their energy as x-rays in a series of stepwise transitions.

Y

The primary x-rays from the rhodium-targeted tube are used to excite the inner electrons of a sample. These electrons (from the K, L and M shells primarily) are expelled from the inner shells of the atom. Electrons from orbitals farther away from the nucleus then fall into the vacancies created. The difference in energy between the orbitals is emitted as a secondary x-ray photon of an exact energy. Because the x-rays come from the inner electrons of the atoms, they are not related to any of the chemical properties of the elements nor to the compounds in which they may be present. These emitted wavelengths are detected by using various analyzing crystals with precisely known d-spacings (distances between the planes of the crystal). Coupled with the crystals is a flow proportional counter and thallium-doped NaI detector. These detectors are used both separately and in tandem.

RANGE :

The present configuration of the Siemens SRS-300 x-ray fluorescence unit provides analytical capabilities from Sodium (element 11) through Uranium (element 92) for a wide range of sample types. The detection limits range from a few ppm to 100 percent. This is dependent upon the element of interest, the matrix effects, the analyzing crystal used, the excitation conditions, the detector, sample preparation techniques, the use of air, vacuum or helium for the analysis conditions, etc. These instrumental parameters are best determined by optimizing the analysis conditions experimentally prior to analysis.

EQUIPMENT AND REAGENTS :

- Siemens SRS-300 sequential wavelength x-ray spectrometer and 10-position sample changer under computer control (DEC family of PDP computers).
- DEC VT-240 monochrome graphics display terminal and DEC Letterwriter 100 for printout of results plus graphic plots.
- 3. Retsch Micro-Mill for grinding samples prior to analysis.
- Sample binder and various dies for pressing samples into pellets using a hydraulic press.

DETERMINATION OF		
Quantitative Elemental Analysis		KMTC - 200-XRF-1
ANALYTICAL METHOD	TYPE OF SAMPLE	
General X-Ray Fluorescence	Solids and Liquids	Page of

- Sample cups with polypropylene covers plus sample holders for the instrument.
- Compressed air to operate the system pneumatics, helium for the sample chamber atmosphere, P-10 gas for the flow proportional counter, plus a Haskris heat exchanger for the x-ray tube and generator.

SAMPLE PREPARATION:

To obtain a sample that is as homogeneous as possible, materials should be ground to as fine a particle size as practical. Particle size is an important consideration, especially when dealing with the lighter elements and their longer and less energetic wavelengths.

The ground solids are mixed with a binder and pressed into a pellet. A smooth surface on the finished pellet is desirable for both accurary and reproducibility of the analysis.

Liquids are generally homogeneous and may be pipetted into the sample cup'as received. The cup and liquid should be set aside for approximately 30 minutes and then checked for leaks. As the sample is inverted directly over the x-ray tube, a leak may deposit sample directly on the beryllium window, causing erroneous results or possibly failure of the tube itself. Semisolid samples (i.e. ROSE samples) may be gently heated and poured into the sample cup. After cooling, they may be analyzed under a helium atmosphere like the liquids.

INSTRUMENTAL PARAMETERS :

These should be determined experimentally for each analyte and it's corresponding matrix.

The most effective wavelength for exciting an element is the wavelength approximately 0.2 Angstrom shorter than the absorption edge wavelength for that element. For example, the Cu K-alpha line at 1.54 Angstrom is very effective in exciting the iron K-spectrum with it's K absorption edge at approximately 1.7 Angstrom.

The kV and mA at which the x-ray tube is operated determine the wavelength and intesity of the primary radiation. These (kV and mA) can be optimized for the elements of interest. The analyzing crystals and detectors can also be optimized for each element of interest.

The use of air, helium, or vacuum in the spectrometer is dependent upon the characteristics of the sample itself.

ANALYSIS:

A qualitative scan can be quickly performed to yield information about the elements that make up the sample (from Z=11 to 92). A quantitative analysis may then be performed in a number of ways.

DETERMINATION OF		
Quantitative Elemental Analysis		KMTC - 200-XRF-1
ANALYTICAL METHOD	TYPE OF SAMPLE	
General X-Ray Fluorescence	Solids and Liquids	Page 3 3

- 1. By standard addition of known quantities of the elements of interest.
- By preparation of a set of standards that have a matrix that is very nearly the same as that of the sample.
- 3. By preparing a standard with known quantities of the elements of interest and using the fundamental parameters approach (Criss Software) to calculate the mass absorption coefficients based on theoretical considerations.

The fundamental parameters approach provides for the analysis of samples with widely divergent matrices, requiring the preparation of only one or two standards instead of an entire suite of standards.

INTERFERENCES :

There are a number of areas where problems may rise in the XRF analysis. Peak overlap can be a problem when analyzing elements that lie very close to a extremely strong peak. For example, the target of the tube is the element rhodium, Z=45; the peak for ruthenium, Z=44, is so close to the rhodium peak that it can not be resolved because of the massive cahracteristic radiation peak coming from the target itself.

Another problem is that the matrix may give rise to absorption-enhancement effects that may affect the analyte line intensity. Another is that samples and standards must be "infinitely thick", or at least the same thickness.

KERR-MCGEE CORPORATION	KMTC - 204-U-16
TECHNICAL CENTER	ISSUED 4/15/76 REVISED 1/22/85
ANALYTICAL METHODS	PAGEOF9
DETERMINATION OF	TYPE OF SAMPLE Rocks, Pulps, Vegetations,
Uranium U ₃ O ₈	Cuttings and Solutions
ANALYTICAL METHOD	WRITTEN BY APPROVED
Fluorimetric	JRJ GEV

Method:

It has been proven experimentally and in actual application that, with the proper care, quick, dependable assays can be turned out with a minimum of manpower. One person can turn out one hundred assays per day without difficulty, including blanks and standards. A crew of three people can handle 80-120 solid samples per 8-hour day on a regular basis. The procedure adopted by the Kerr-McGee Technical Center laboratory was developed by National Lead Co., Winchester Laboratory, and reported in <u>"Analytical Chemistry</u>," Vol. 28, Page 1651, November, 1956.

A new rotary fusion burner system has been installed, on which 22 fluoride fusions are performed simultaneously under reproducible fusion conditions. A fluorometric flux mixture of 2% lithium fluoride and 98% sodium fluoride is used and has been found to be superior to other fluxes in sensitivity, precision and convenience. The coefficient of variation of the fusion method (only) employing this new flux has been determined to be 0.7%. This constitutes at least a sixfold improvement in precision over that of the previous fusion system. This method has been proven to be rapid, accurate and economical.

A combination of dilution and solvent extraction is used to separate the uranium from possible interferences and quenchers. An aliquot of a liquid sample or of a solution of a solid sample is diluted and acidified so that the diluted sample contains approximately .005 mg of U_3O_8 per ml in a 5% nitric acid solution. One milliliter of the diluted sample is salted with saturated aluminum nitrate solution and the uranium is extracted into 10 ml of ethyl acetate. Aliquots of 0.2 ml are removed and transferred onto pellets of 2% lithium fluoride/98% sodium fluoride flux in platinum dishes. The pellets are dried, fused and allowed to cool, forming buttons which slide easily from the dish. The fluorescence of the fluoride button is measured with a Jarrell-Ash Model 26000 fluorometer.

Apparatus:

40-ml Vials - 4-1/4" in length and 1" in diameter.

Aluminum Foil-Lined Caps - purchased in quantities and used only once.

Machlett Autopipet - 10-ml capacity.

Platinum Dishes - These dishes are formed from discs 0.015" thick and 0.75" in diameter. Discs are available from Englehard Instruments, Inc., 113 Astor Street, Newark, New Jersey.

DETERMINATION OF		
Uranium U ₃ O ₈		KMTC 204-U-16
ANALYTICAL METHOD	Rocks, Pulps, Vegetations,	KMI C 204-0-10
Fluorimetric	Cuttings and Solutions	Page _2of9

Dish Forming Tool - Available from Jarrell-Ash Co., Waltham, Mass., Catalog No. 26100. Before forming, discs should be annealed at 850°C. After forming, they should be cleaned with 1:1 nitric acid and rinsed thoroughly with water.

Pellet Maker - Jarrell-Ash Co., Catalog No. 26110.

- Rotary Burner The burner is a rotary fusion burner. The burner is available from EDA Intruments Inc., 5151 Ward Road, Wheat Ridge, Colorado 80033. The system consists of fifteen Humboldt burners arranged in a circle, two flowmeters, temperature readout meter, gas igniter automatch, low heat timer, high heat timer.
- Drying Units Fabricated in house using four 500 watt infrared heat lamps mounted in two separate boxes of two lamps each 19" long, 12" deep and 10" high. The boxes are constructed of 1/2 inch thick transite and is closed all around except for a 19" x 10" opening in front.

Propane Regulator - Must be able to maintain 30 lbs output pressure.

- Shaking Machine Eberback Model 6010. Available from Scientific Products.
- Fluorometer Jarrell-Ash Co., Model 26000 fitted with reflectance attachment, Model 26050.
- Microliter Pipet Eppendorf 200-µl pipet with a supply of disposable pipet tips.
- Reagents The use of reagents that contains negligible uranium or other fluorescence impurities is essential. With the exception of sodium fluoride, it has bee found that Mallinckrodt reagent grade chemicals give the bes all-around performance.
- 98% Sodium FLuoride 2% lithium fluoride Fluorometric grade, EDA Instruments, Inc., 5151 Ward Road, Wheat Ridge, Denver, Colorado 80033.

Ethyl Acetate - Mallinckrodt reagent grade.

Aluminum Nitrate - Mallinckrodt reagent grade.

A solution of aluminum nitrate is made by dissolving 5 lb in 1100 ml deionized H_2O .

A reagent blank prepared from the above reagent should read no more than 15 units on the most sensitive meter range. Full scale meter reading - 10 µamps.

DETERMINATION OF		
Uranium U308		KMTC - 204-U-16
ANALYTICAL METHOD	Rocks, Pulps, Vegetations,	2
Fluorimetric	Cuttings and Solutions	Page of

<u>Acid</u> - Reagent grade nitric, perchloric, hydrochloric and hydrofluoric acids.

Mallinckrodt reagent grade.

Perchloric acid should only be used in approved hoods with washdown facilities.

Standard U30g Solution - 1.000 grams of 99.99% U30g is weighed into a 250-ml beaker. Add 10 ml HClO4 and 10 ml HNO3 and fume to 3-5 ml. Cool, rinse sides of beaker with DI H20, and 50 ml concentrated HNO3 and dilute to 1.00 liter.

1.00 ml of this stock solution diluted to 200 ml will give a standard solution of 5 μ g/ml U₃O₈ which should be made fresh once each month.

A. Sample Preparation

1. Aqueous Samples

Dilute all aqueous samples according to the estimated concentration so that 1 ml of the diluted sample contains approximately .005 mg of U₃Og. Since the diluted sample should contain 5 percent (v/v) free nitric acid, samples which are alkaline musc be neutralized before the proper amount of nitric acid is added.

Samples which are known to contain reduced uranium, such as Wet Phosphoric Acid (WPA), must be oxidized to uranium (VI). The preferred method is to pipet an aliquot into a 250-ml beaker. Add 5 ml HClO4 and 5 ml HNO3 and take to strong HClO4 acid fumes, cool, add HNO3 to make dilution 3-5% v/v HNO3 and dilute in mixing cylinder.

Organic samples, such as those from solvent extraction studies, are treated by stripping the uranium into aqueous carbonate solution before proceeding with the dilution and extraction. Pipet a suitable volume of the sample into an extraction vial and add 2 to 3 ml of carbon tetrachloride to make the organic layer heavier than water. In special cases, wet oxidation is the preferred treatment.

Add 10 ml of 5 percent (w/v) sodium carbonate solution from an automatic pipet, cap the vial, agitate on the shaker for 1 to 2 minutes, and cent ifuge to separate the phases completely. Treat the resulting carbonate solution according to the procedure described above for an aqueous sample.

2. Solids

Pulped samples are weighed, normally 2 grams, into a 250-ml beaker and treated with 10 ml each of perchloric acid and nitric acid. Place on hot plate and take to fumes, remove and add 5 ml perchloric, 10 ml of

DETERMINATION OF		in the second
Uranium U308		KMTC 204-U-16
ANALYTICAL METHOD	Rocks, Pulps, Vegetations,	
Fluorimetric	Cuttings and Solutions	Poge 4 of 7

hydrochloric and mix. When hydrochloric acid ceases to effervesce, add 1-2 ml hydrofluoric acid. Cover and heat on low burner for a few minutes and then fume to 3-5 ml. Cool sample, add 5 ml HNO_3 and 25-30 ml H_2O and transfer to a 100-ml volumetric flask. Dilute to mark. No filtration is necessary. Shake well before aliquot is made for analysis.

3. Vegetations

Dry sample at 105°C over 16 hours. Pulverize dried sample using a Wiley Mill down to approximately 8-10 mesh. Weigh 30 gms of the dry pulverized sample into a 250-ml pyrex beaker. Ash in an oxygen oven overnight by slowly increasing the heat 50°C per hour to 550°C with an adequate oxygen flow on.

Cool and transfer the ash to a 250-ml PTFE beaker with a 10% nitric acid (HNO3) rinse, using a rubber policeman to dislodge solids from sides of beaker. Add 10 ml nitric acid, 5 ml perchloric acid and 2 ml hydrofluoric acid to the sample. Digest the sample approximately 230°C on a hot-plate and evaporate to incipient dryness. Add 10 ml of nitric acid to the beaker and transfer contents to a 200-ml volumetric flask using distilled water. Dilute to mark with deionized water and shake well. Transfer a 2 or 3 ml aliquot into an extraction vial and proceed to the extraction section of this procedure. If the fluorescence is not greater than 1 unit above background, transfer a sufficient amount of sample to a 125-ml separatory funnel and proceed with the extraction step described in KMTC-125-U-17.

B. Extraction

Pipet the desired aliquot of the diluted sample into an extraction vial. If no dilution is required and the sample is alkaline, acidify by adding concentrated nitric acid dropwise to the aliquot taken for extraction. Add approximately 15 ml of saturated aluminum nitrate and exactly 10 ml of ethyl acetate from the automatic pipet. Cap the vial and agitate on the shaker for 3 to 4 minutes. To avoid dilution of the saturated aluminum nitrate in the extraction vial, which may lead to low uranium values, aliquots for extraction should not exceed 2 ml. Extract a reagent blank and standard with each set of samples.

C. Preparation of Pellets

Place 22 platinum dishes on wire rings and then ensure that the rings and the sample are level (see pattern Figure). Place the rings and dishes under the drying unit until all the dishes are dry. The dishes must be dry since placing flux pellets onto wet dishes will eventually cause a severe attack on the platinum.

Using the pellet maker, prepare and place a pellet of fluoride flux in each dry dish according to the following procedure: Fill a 250-ml beaker (or other suitable container) with flux to a depth of 2-1/2 to 3 inches. While rotating the beaker, pack the barrel of the pellet maker by plunging

DETERMINATION OF		an har bland a start a
Uranium U308		KMTC - 204-U-16
ANALYTICAL METHOD	Rocks, Pulps, Vegetations,	5 . 9
Fluorimetric	Cuttings and Solutions	Page of

it vertically through the bed of flux. While the pellet maker is still in a vertical position, draw it across the bottom of the beaker. Holding the filled pellet maker about three-quarters of an inch above the platinum dish, push down on the plunger and discharge the pellet. Rotation of the beaker while filling the pellet maker will assure a uniform pellet weight, and drawing the pellet maker across the bottom of the beaker will assure a flat bottom to each pellet. Uniformity of pellet size is necessary to maintain reproductivity. Using the above procedure, adjust the pellet maker to form pellets that weigh 0.5 grams \pm .02 grams.

Using a 200-µl Eppendorf automatic pipet with disposable pipet tips, pipet 200 µl of ethyl acetate containing 0.5 micrograms U_3O_8 per ml. The ethyl acetate solution is prepared by extracting 1 ml of a standard 5 µg/ml U_3O_8 aqueous solution using extraction method described above. Two pellets are reserved for the reagent blank. Volatilize the ethyl acetate from the pellets by placing the wire rings holding the dishes under the drying units for approximately 10 minutes.

D. Fusion

The wire ring holder with the twenty-two platinum dishes and pellets is placed on the burner assembly. Be sure that the rings and the sample dishes are level. Turn on the exhaust fan in the heat exhaust hood system. Place the power switch in the ON position. Place the disc rotation switch in the ON position. Be sure the low heat flowmeter is set at 30 ml/min before igniting the burners. Push the Start/Ignite button. The igniter tip can be manually adjusted to give proper spark gap. The igniter will continue to spark until the system has ignited. Adjust the low heat flowmeter to give the desired temperature. The Low Heat cycle is pre-adjusted to be approximately 30 seconds. When the High Heat Timer light goes on, adjust the High Heat Flowmeter to give the desired temperature of approximately 950°-1000°C (80 ml/min). This cycle is usually set at 2 minutes. The Pre-Cooling/Annealing cycle that follows the High Heat cycle will be the same temperature as the Low Heat cycle. After the 3 minute fusion cycle, the burner control unit will automatically shut off the gas supply. As soon as the melted flux has solidified, remove the wire screen and place on a heat-resistant pad to cool.

The fluorescence intonsity increases for the first 15 minutes and then remains constant for about an hour. The cooling rate affects the ultimate crystalline structure of the fluoride buttons, which, in turn, affects the fluorescence of the buttons. Therefore, it is necessary to follow the same procedure in cooling each set of dishes; i.e., place the ring assembly on a heat-resistant material away from drafts. Protect the fluoride buttons from exposure to ultra-violet light during the cool-down period. Exposure to UV light decreases the intensity of the fluorescence. Unshielded fluorescent light tubes give off enough UV light to effect the uranium fluorescence.

If, after cooling, it is noticed that the fused flux has a brown or pinkish discoloration, it is probably due to too high a fusion temperature during the first two minutes. Reduce the temperature by increasing the

DETERMINATION OF		
Uranium U308		KMTC 204-U-16
ANALYTICAL METHOD	Pocke Pulpe Vegetations	6 0
Fluorimetric	Rocks, Pulps, Vegetations, Cuttings and Solutions	Page of

gas flow 0.2 units or so on the flowmeter. Once gas and air settings are established for a burner, these settings will remain reasonably constant over long periods of time.

E. Measurement of Fluorescence

1. Instrument Adjustment

- a. Adjust the phototube voltage to obtain a standard meter reading with a fluoride disk containing 0.1 microgram of U_3O_8 . A suggested standard reading is 50 microamperes. If possible, all samples should be diluted to read approximately the same.
- b. Adjust the meter zero by placing the slide in load position and manipulating the zero control.
- c. Push the sample slide in, depress the lowest sensitivity scale key, and adjust the background current to zero.

2. Fluorescence Measurements

The fluorescence intensities of the first two standards are measured first. Then measure the samples and finish by measuring the last two standards.

- a. Pick up a platinum dish with a pair of forceps and tip it so that the fluoride button slides into the polished depression in the slide.
- b. Push the slide to the backstop. For buttons of unknown intensity, depress the keys in increasing order of sensitivity so as not to drive the meter needle off scale.
- c. Pull the slide forward and remove the button with forceps, cleaning any fluoride particles from the slide by aspiration.

F. Calculations

For each set of twenty buttons, average the readings of all duplicate samples and the readings of the four control standards and subtract the reagent blank. Calculate uranium concentration by the following formula:

mg/l U₃O₈ = µg/ml in Std x (Avg reading of unknown) x (vol. of dilution in ml) Avg Std Reading aliquots in milliliters

ppm U₃O₈ = <u>µg/ml in Std x (Avg reading of unknown) x (vol. of dilution in ml)</u> Avg Std Reading (weight of sample in gms) x (aliquot in milliliters)

e.g., ppm $U_3O_8 = \frac{5 \ \mu g}{50 \ \text{Units}} \times \frac{(50 \ \text{units}) \times (100 \ \text{ml dilution volume})}{(2 \ \text{ml aliquot}) \times (2g \ \text{sample weight})} = 125 \text{ppm}$

DETERMINATION OF		
Uranium U308		KMTC 204-U-16
ANALYTICAL METHOD	Rocks, Pulps, Vegetations,	
Fluorimetric	Cuttings and Solutions	Page of 9

The above calculation assumes that the same volume of ethyl acetate is used for the standards and the unknowns; e.g., 10 ml for extraction and 200 μ l aliquot on the fluoride pellet or 1 ml aq. x $\frac{5 \ \mu g}{ml} \times \frac{0.2 \ ml}{10 \ ml} \frac{\text{ETAC}}{10 \ ml}$

= 0.1 µg U30g actually on pellet.

G. Cleanliness

As with any other method employed for trace analysis, contamination must be carefully guarded against at all times. The fluorometric method of analysis is generally used in a laboratory or plant where uranium ore dust or even uranium concentrate dust is present in the atmosphere. Therefore, a room supplied with filtered air and used only for fluorometric analysis is a requirement. Each week all the floors and benches in the room should be washed and about once a month the walls and overhead fixtures should be vacuum cleaned.

Glassware and other equipment used for this analysis should not be used for other purposes and should be scrupulously cleaned after each use. People working in the fluorometer room should not use protective hand creams since these creams usually have a strong ultraviolet fluorescence, and anything touched by the hands will become contaminated.

For best quality ultra-trace analysis, only plastic and teflon ware should be used. All glassware and quartzware contains measureable levels of uranium which can bias high the final result when ultra-trace levels of uranium are being determined.

H. Costs

In groups of at least 15 rock or soil samples, the cost is 0.23 man-hours per sample, excluding sample preparation (crushing and pulping) cost.

I. Sensitivity, Precision and Accuracy

The lower limit is 1 ppm; upper recommended limit is 5000 ppm.

Accuracy is very dependent upon the nature of the sample and the quality of the standards. In no circumstances should accuracy be assumed to be equivalent to the precision. In a study on a Grand Junction analyzed sandstone certified at 0.107Z U30g, the following results were obtained over about a one-month period:

0.107% U308	-
0.101	$\bar{x} = 0.109$ U ₃ O ₈
0.112	S = 0.0052% U308
0.107	RSD = 4.72%
0.109	
0.117	
0.113	

DETERMINATION OF		the second s
Uranium U308		KMTC 204-U-16
ANALYTICAL METHOD	Rocks, Pulps, Vegetations,	NMIC 204-0-10
Fluorimetric	Cuttings and Solutions	Poge 8 of 9

J. Other Precision Data


	Sample 1 Sand Tail, ppm U ₃ O ₈	Sample 2 Slime Tail ppm U ₃ O ₈	Sample 3 Sandstone, ppm U ₃ O ₈	Sample 4 Sandstone, ppm U ₃ O ₈
	33	64	140	230
	36	66	150	240
	31	62	150	240
	35	70	140	230
	36		150	230
	36		140	230
*	34.5	65.5	145	233
S	2.07	3.4	5.5	5.2
RSD	6.0	5.22	3.8	2.2
RSD	6.0	5.22	3.8	2.2

From these five levels of U_3O_8 , precision(s) can be described by a linear equation:

± s = -1.48 + .0478 *.

KMTC-204-U-16

Page 9 of 9



KERR-MCGEE CORPORATION		KMTC - 208-EC-22
TECHNICAL CENTER	ISSUED 10/13/82	
ANALYTICAL METHODS		PAGE 1 OF 23
Classification of Air-Borne Particulates in	Air-Borne Particu	lates
In-Vitro	AKD/AKW	GEV

INTRODUCTION:

The manufacture of nuclear fuel elements, as at the Sequoyah Facility of Sequoyah Fuels Corporation, where uranium hexafluoride is produced from the uranium concentrate, involves the processing of several uranium compounds which at certain stages are present in finely divided form. The continuous handling of large quantities of these materials, and some escape of volatile fluorides of uranium, inevitably results in some air contamination.(1) By determining the concentration and size-spectrum of particulates within the working environment as well as at the nearest residence, the likely intake and retention of material in the human body (a measure of associated hazard due to inhalation) may be calculated using the Standard Man data of the International Commission on Radiological Protection, "ICRP (1959)" in association with the lung model proposed by the ICRP Task Group on Lung Dynamics (1966). This special task group was created by the ICRP Committee (2) for the purpose of reviewing the so-called lung model; a scheme for computing dust deposition as a function of particle size and clearance from the human respiratory tract thereby providing a basis for lung dosimetry from radionuclide deposited and the setting of exposure limits. The clearance of deposited material in the human body is carried out by numerous biological processes. The rates of these processes depend upon the solubility of the materials concerned in biological fluids. It has been known for some years that some materials insoluble in water have nevertheless a significant solubility in biological fluids. The ICRP commented on it in lung dynamics. Some initial solubility studies were performed with real plasma, but due to bacterial deterioration of the plasma, enormous difficulties were encountered. It was subsequently decided to use a simulated lung fluid as the solvent for solubility studies of air-borne uranium particles.

According to the ICRP lung clearance model, the deposited material in the lung is classified into three groups. The basis of this classification is the rate at which it leaves the lung. These three classes are D, W, and Y, corresponding to half-times in the lung of 0 to 10 days, 11 to 100 days, and >100 days, respectively. It was originally assumed that the rate at which the mass leaves the lung is proportional to the lung burden and a simple relationship is obtained which predicts that the lung burden will diminish exponentially with time. But, in many cases, the clearance of material from the lung has been found to be progressively slower than predicted by a single exponential process, and it becomes necessary to include additional exponential terms. The deposited material is then classified according to the weight fractions of D, W and Y components. Although biological processes (3) like endocytosis and ciliary-mucus transport are known to contribute to the lung clearance, the complete and thorough biological data are not available. As a result, the dissolution half-times for materials

DETERMINATION OF The Solubility and Di	ssolution Rate Classification of	
Air-Borne Particulates in Simula	ted Lung Fluid	KMTC - 208-EC-22
ANALTTICAL METHOD	TYPE OF SAMPLE	2
In-Vitro	Air-Borne Particulates	Page 2 of 23

(air-borne particles) in simulated 'ung fluid have been assumed to be approximately equal to their lung clearance half-times.(1,2) When the lung-clearance classification for a material is determined by the solubility test in simulated lung fluid, its transport rates between other anatomical compartments are automatically evaluated. From these values, the residence time of the material and the associated radiation dose in each compartment can be calculated.(4)

The Nuclear Regulatory Commission requires Kerr-McGee Corporation to carry out the dissolution-rate classification test of uranium on each quarterly composited dust sample collected near the uranium conversion plant. Samples from other locations or facilities may be required in the future. The method followed here for the dissolution test is that of D. R. Kalkwarf, (5)as required also by the NRC. The dissolution test is carried out at 37°C. Although the U.S. Environmental Protection Agency has recommended that all particulates with <10 µm be considered inhalable for hazard evaluations, the solubility test is performed on the whole sample which contains both respirable and non-respirable dust. Since the lung is believed to be a site for efficient dissolution, the solubility test is performed in a well-stirred suspension to achieve maximum clearance rates. As a result, these values should also approximate the lung clearance rates that include contributions from endocytosis and ciliary-mucus transport.

SAMPLE COLLECTION AND PREPARATION

For each quarter, air particulate samples are collected on filter papers. A high-volume air sampler is used for collection of air particulates. The filter paper used is replaced by a new one (8 X 10" Nucleopore or Millipore) every week during the quarter. The dust sample collected on the filter paper contains both respirable and non-respirable particulates.

Each filter paper containing air particulates is split. One-half of each filter paper for the quarter is either retained or sent to the U.S. Nuclear Regulatory Commission for performing the verification of uranium dissolution test in simulated lung fluid. The other half is again split; one-half of this is kept for gross alpha and uranium analysis, and the other half is utilized for the uranium solubility test.

In order to get a representative sample, filter paper strips of proper size are cut from each original filter paper containing the sample. The sizes (areas) of the strips are determined in such a way that the volume of air passing through each filter paper strip is approximately the same. These dust coated filter paper strips are dried in a desiccator over anhydrous calcium sulfate (Drierite) for 5 days. At the end of the drying period, the dust particles are 'vacuumed off' the surface of each filter paper strip with a vacuum line fitted with a plastic filter-holder (Millipore Corporation, Swinex 25) containing a 25-mm diameter membrane filter (Millipore, Type HA, pore size 0.45 µm). The dust collected on the membrane filter is sterilized by running ethylene oxide gas, from a small lecture bottle, through the membrane filter at a very slow rate for 4 hours. Sterilization is important so that the possible effects of bacterial growth in the suspension can be prevented. The dust collected on the membrane

Determination of The Solubility and Dissolution Rate Classification of

Alr-Borne Particulates in Simulated Lung Fluid		KMTC - 208-EC-22
ANALYTICAL METHOD TYPE OF SAMPLE		
In-Vitro	Air-Borne Particulates	Page 3_ot_23

filter is then removed by a camel-hair brush into a weighed glass reacti-vial. The weight of the dust particles in the vial is determined. The dust particles in the glass vial are homogenized by shaking. This sample is generally found to contain some long fibers of the original filter paper.

If a sufficient mass of the dust sample is available, then a portion (known amount) of it is set aside for checking uranium mass balance, before transferring it to the glass vial.

PREPARATION OF SIMULATED LUNG FLUID:

The simulated lung fluid is prepared according to the method described by D. R. Kalkwarf.(5) The composition of the simulated lung fluid is given in Table I.

Table I

Composition of the Simulated Lung Fluid

Compound	Concentration, g/l of Water	
Magnesium chloride, hexahydrate, MgCl ₂ *6H ₂ O	0.203	g/1
Sodium chloride, NaCl Potassium chloride, KCl	6.019	
Disodium hydrogen phosphate, heptahydrate.	0.290	
Na 2HPO 4 *7H 20	0.268	**
Sodium sulfate, Na ₂ SO ₄	0.071	**
Anhydrous calcium chloride, CaCl ₂	0.278	**
Sodium acetate, trihydrate, NaH3C2D2*3H2O	0.953	**
Sodium bicarbonate, NaHCO3	2.604	**
Sodium citrate, dihydrate, Na ₃ H ₃ C ₆ O ₇ •2H ₂ O	0.097	**

The pH of this aqueous solution is adjusted to 7.3-7.4 by adding IN HCl dropwise, with the help of a pH-meter. It has been found that the electrolytic composition of average human interstitial lung fluid is identical to that of the simulated lung fluid. For comparison, the electrolytic compositions of both lung fluids are given in Table II.

	issolution Rate Classification of	1
Air-Borne Particulates in Simula	KMTC - 208-EC-22	
ANALYTICAL METHOD	TYPE OF SAMPLE	4 23
In-Vitro	Air-Borne Particulates	Pogeof

Table II

Electrolytic Compositions of Average Human Lung Fluid and Simulated Lung Fluid

Ion	Fuman Lung Fluid	Simulated Lung Fluid
Ca ⁺² Mg ⁺² K	5.0 meg/1	5.0 meq/1
Mg ⁺ 2	2.0 "	2.0 "
ĸ	4.0 "	4.0 "
Na ⁺	145.0 "	145.0 "
Total Cations	156.0 meg/1	156.0 meq/1
HC0 3-	31.0 meq/1	31.0 meg/1
Cl	114.0 "	114.0 "
H ₅ C ₆ O ₇ (citrate)	_	1.0 "
H ₃ C ₂ O ₂ (acetate)	-	7.0 "
Total Organic Anion	7.0 "	-
Protein	1.0 "	-
	2.0 "	2.0 "
HPO 4 2	1.0 "	1.0 "
Total Anions	156.0 meg/1	156.0 meq/1
рH	7.3 - 7.4	7.3 - 7.4

In this simulated lung fluid, acetate is substituted for all organic acid anions present in the actual lung fluid. Another difference between these two kinds of lung fluids is the absence of protein component in the simulated lung fluid. The lung fluid proteins are poorly characterized. Some investigations indicated that even the presence of very small amounts of protein substitutes (0.4% W/V)(1) which have very high molecular weight, renders the fluid almost unfilterable. It also causes bacterial growth in the solution. So, the protein is omitted in the main body of this work. The protein is present in very small amounts in the human lung fluid and is also known to be a poor ligand for complexing metal ions like uranium. Hence, the absence or presence of a protein component in the simulated lung fluid will not cause a significant change in the solubility of uranium compound. The protein component of the human lung fluid has been replaced by an ionically equivalent amount of citrate in the simulated lung fluid, as suggested by Moss(7). It is justified by the fact that the citrate ligand and the protein have comparable complexing capabilities.

The actual lung fluid also contains traces of phospholipids. These are not added to the simulated lung fluid because of the filtering problem. One recent study (8) showed that the dissolution rate of yellow cake (U $_{3}O_{8}$) is not affected by the presence of dipalmitoyl lecithin (a phospholipid) in the simulated lung fluid.

The pH of the simulated lung fluid prepared in the laboratory is checked periodically during the dissolution experiment. It is found that the solution becomes slightly cloudy and shows a higher pH on standing for several days. Under these circumstances, the pH of the solution is adjusted to 7.3 to obtain a clear solution or a fresh solution is made.

offermination of The Solubility and D	issolution Rate Classification of	
Air-Borne Particulates in Simula	KMTC - 208-EC-22	
ANALYTICAL METHOD	TYPE OF BAMPLE	-
In-Vitro	Air-Borne Particulates	Page of23

DISSOLUTION TECHNIQUE:

Equipment and Apparatus Needed:

- Reacti-Therm Module with a low density heater which has a built-in stirring capability (Pierce Chemical).
- 5-ml Reacti-vials (Pierce Chemical) with teflon-lined screw caps as shown in Figure 1.
- Teflon-coated magnetic stirrers (Pierce Chemical).
- 4. Stainless steel filter holder (Millipore, Swinex).
- 5. Disposable plastic syringes.
- Membrane filters (Millipore, 13-mm diameter, GC, 0.22 µm pores).
- 7. Fisher Centrific Centrifuge (Model 225).
- 8. Thermometer

Dissolution of the dust sample, which is vacuumed off the surface of the filter paper, is carried out in stirred, 5-ml volumes of simulated lung fluid at 37°C. Before the start of the actual dissolution test, the power of the Reacti-Therm Module is turned on and the temperature of the heating block/stirrer assembly is set at 39°C. It takes about 24 hours for the desired temperature to get stabilized within $\pm 1^{\circ}$ C. Preliminary examination indicated that the actual temperature of the fluid itself in the Reacti-vial shows 37°C $\pm 1^{\circ}$ C, when the temperature of the heating/stirring block reads 39 $\pm 1^{\circ}$ C. The stirring rate is also kept constant during the experiment.

The Reacti-vial, containing a known amount of the dust sample, and another vial containing 5-ml of simulated lung fluid alone are both placed into the heating block for 5-10 minutes before they are mixed. When they reach the desired temperature (37°C), the lung fluid (5-ml) in the vial is withdrawn by a syringe and is added to the vial containing the sample. The time is noted and considered to be the starting (zero) time. The suspension in the vial is stirred with a Teflon-coated magnetic stirrer. After selected intervals of time, the vial is removed from the heating/stirring block, and centrifuged in the Fisher Centrific Centrifuge for 5 minutes. The supernatant liquid is drawn into a plastic syringe. A stainless steel filter holder (Millipore, Swinnex) containing a membrane filter (Millipore, 13 mm diameter, GC, 0.22 µm pores) is fitted on the end of the syringe. The solution is filtered into a plastic vial. Two drops of concentrated HNO3 are added to this filtrate (5-ml) to prevent any adsorption of the uranium compound on the walls of the vial. This acidified solution is stored for uranium analysis. The membrane filter is then removed from the stainless steel filter holder and 5-ml of the simulated lung fluid is taken into the plastic syringe. The filter holder with a syringe needle attached is fitted onto the

DETERMINATION OF The Solubility and Di	issolution Rate Classification of	
Air-Borne Particulates in Simulat	KMTC - 208-EC-22	
ANALYTICAL METHOD	TYPE OF SAMPLE	 The solution is status for an
In-Vitro	Air-Borne Particulates	Page 6 of 23

syringe. The small amount of solid left on the filter is washed down into the Reacti-vial with the jet of lung fluid from the syringe. The filter paper used is also stored for uranium analysis. The Reacti-vial is then capped, shakened to resuspend all the solid particles and replaced in the heating block. This process is repeated until all the samples are collected for selected dissolution periods. At the end of the dissolution trial, the residual sample is suspended in H₂O and this whole suspension is analyzed for uranium.

A blank is also run with the simulated lung fluid alone under identical conditions (37°C and the same stirring rate). The 5-ml samples are taken approximately at the same intervals of time throughout the experiment and analyzed for uranium to find out if there is any uranium leached out from the glass material of the Reacti-vial or present in the lung fluid itself.

Uranium Analysis:

The filtrates from the dissolution trial are analyzed for uranium by direct fluorometry. (9) In some cases where the uranium content is very low, preconcentration is carried out by extraction from aluminum nitrate with ethylacetate, followed by the fluorometric analysis on the extract. In the case of the solid residue and the membrane filters, they are first digested with concentrated HNO₃, diluted and then the fluorometric analysis is performed on these diluted samples.

Discussion:

A. Reduction of Data and Evaluation of Dissolution Half-Times

If the sample contains one pure component and the particles are present in a narrow size range (10), then the fraction of the component (sample) remaining undissolved will decrease exponentially with time. But in the case of a real world sample, many uranium components are usually present with different dissolution rates. As a result, there will be additional exponential terms and the dissolution data will be best fit by an equation of the following form:

 $F = f_1 \exp (-0.693\frac{t}{T_1}) + f_2 \exp (-0.693\frac{t}{T_2}) + \dots + f_n \exp (-0.693\frac{t}{T_n})$

where F is the fraction of total uranium remaining undissolved after time t, and f, is the initial weight fraction of uranium components in the sample with dissolution half-time T_i . The values of F are calculated by substracting the amount of uranium dissolved during any sampling period from the amount undissolved at the beginning of that period and dividing this quantity by the total amount of uranium in the sample. The weight of uranium in the sample portion used in the dissolution trial is calculated by summing the amounts of uranium dissolved during the dissolution time increments and the amount of uranium in

Page 7 of 23 KMTC-208-EC-22




FIGURE 1. Cross Section of Dissolution Chamber for Dust Samples

204

.

DETERMINATION OF The Solubility and	Dissolution Rate Classification of		
		KMTC -	208-EC-22
ANALYTICAL METHOD	TYPE OF SAMPLE		
In-Vitro	Air-Borne Particulates	Page	8 of 23

the undissolved residue. As an example, the values of F calculated during the dissolution test of the third quarter, 1981 sample are given in Table III.

Table III

Amounts of Uranium Dissolved and Fraction of Total Uranium Remaining Undissolved During the Dissolution Test of the Air Particulate Sample (3rd Quarter, 1981)


Time (Days)	Uranium Dissolved (µg)	Fraction Undissolved (F)
0.00	0.00	1.00
0.021	14	0.613
0.063	6.5	0.433
0.104	2.6	0.360
0.167	1.1	0.330
0.250	0.55	0.314
0.976	0.75	0.294
1.26	0.32	0.285
1.92	0.22	0.279
2.25	0.14	0.275
5.21	0.24	0.268
6.96	0.14	0.264
12.2	0.023	0.263
20.9	0.036	0.262
35.1	0.070	0.261
76.1	0.10	0.258

Solid residue and composite of filter papers used for = 9.32 µg of U filtration of different fractions

Total uranium (µg) in the sample used for the = 36.2 µg of U experiment

The values of F are plotted against t (days) on semi-log paper. The shape of this plot indicates that the dissolution involves mainly two rate processes; the first one is very fast, and the other is a slow dissolution process. The plot of F-values versus t(days) for the third quarter 1981 sample, is also shown in Figure 2. From this plot, the values of T_i and f_i can be determined manually from the slopes of the tangents drawn on the plot. The first portion of the curve is so steep that it is very difficult to draw a tangent accurately. As a result, the values of T_i and f_i determined manually can be associated with 50 to 100% error. Instead, the best values of f1, T1, f2, T2, f3 and T3 are obtained by fitting the experimental data into a 3-term equation of the form (see page 7). A program (11) for non-linear regression analysis was developed for five

Page 9- of 23



Fraction of Uranium Undissolved (F) vs. Time (days)

FIGURE 2

DETERMINATION OF The Solubility and D:	Issolution Rate Classification of	
Air-Borne Particulates in Simulat	ted Lung Fluid	KMTC - 208-EC-22
ANALYTICAL METHOD	TYPE OF SAMPLE	
In-Vitro	Air-Borne Particulates	Poge 10 of 23

independent variables. This program is called 'NEWBRIT' and its Fortran Coding is given in the Appendix. For determining the sixth variable, one condition has been introduced into the program, that the summation of the weight fractions of three components must be equal to 1, that is $f_1 + f_2 + f_3 = 1$. This particular program needs a "seed" error value associated with the feeding values of both F and t in order to start. The arbitrary values of 0.0001 and 0.01 are introduced as errors associated with t and F, respectively. In the iterative program, these values of error are varied in such a way that the best values of f_1 , T_1 , T_2 and T_3 are obtained with minimum errors. Then the value of f_3 is determined by difference.

The data is entered into the data file "ASOK" as in the following example:

.6689 .0092 .1383 1.2454 320 .1 .00001 .001 .01 50 1. 1. 1. 100 1E5 0 .0001 1. .01 0.0208 .0001 .724 .01 .0625 .0001 .647 .01 .1042 .0001 .613 .01 .1874 .0001 .583 .01 .2708 .0001 .563 .01 .9422 .0001 .526 .01 1.276 .0001 .508 .01 1.928 .0001 .498 .01 3.02 .0001 .489 .01 6.21 .0001 .48 .01 9.24 .0001 .473 .01 13.95 .0001 .468 .01 27.12 .0001 .459 .01 78.12 .0001 .349 .01

Lines 1-3 remain unchanged from sample to sample as these are the "seed" values used by the program to calculate the variables. Lines 4-18 are the actual data obtained experimentally. Column 1 is the time in days from the start of the test.

Columns 2 and 4 are the arbitrary error values discussed above. Column 3 is the fraction of uranium undissolved at the time listed in Column 1. After the data is entered, the program "NEWBRIT" is run manually. At each iteration the operator enters values for IT (number of iterations) = usually 1, IDEM (maximum number of deming iterations) = usually 1, and GAMACT (maximum feasible step from binding bounds) = any number from 10^{-5} to 10^{-5} . After the first iteration, the value for GAMACT is varied to slowly approach the SumDEL (sum of values of the data adjustments) value calculated for the previous step. This method, when taken through 30-50 iterations, produces reproducible values for f₁, f₂, T₁, T₂ and T₃ with errors on the order of < 1% of the value.

DETERMINATION OF The Solubi	lity and Dissolution Rate Classification of)f
Air-Borne Particulates in Simulated Lung Fluid		KMTC - 208-EC-22
ANALYTICAL METHOD	THE OF SAMPLE	
In-Vitro	Air-Borne Particulates	Poge 11 of 23

B. Dissolution Rate Classification

Solubility classification of the sample is based on the weight fractions and dissolution half-times of the uranium components. In order to make a solubility classification, let us consider the dissolution parameters (weight fraction, f, and dissolution half-times, T) for the third quarter, 1981 air particulate sample. These dissolution parameters are given in Table IV.

Table IV

Dissolution Parameters for Third Quarter, 1981 Air Particulate Sample

f1	T1 (days)	<u> f</u> 2	T ₂ (days)	<u> f</u> 3	T ₃ (days)
0.50±0.045	0.012±0.0014	0.223±0.044	0.083±0.018	0.277±0.063	547±232

From these values, it is inferred that 50% of the uranium in the third quarter, 1981 sample had a dissolution half-time of 0.012 days, 22% of the uranium with half-time of 0.083 days, and 28% of the uranium dissolved with half-time of 547 days. So the dissolution rate of the uranium in the sample is classified as 0.50D, 0.22 D, 0.28 Y. Since the sample contains two components in the D category of the ICRP task group lung model, these can be combined (added) to show the total weight fraction of uranium in this category.

Since the dissolution trial is carried out over a period of around 80 days, it is sufficient to categorize the uranium components into the solubility classes (1) used by the ICRP Task Group on Lung Dynamics. The classification scheme adopted is therefore:

- Soluble compounds with a solubility half-life <1-10 days, inclusive, equivalent to clearance class D.
- (2) Moderately Soluble compounds with an estimated solubility half life >10-100 days, inclusive, equivalent to clearance class W.
- (3) Relatively Insoluble compounds with an estimated solubility half-life >100 days, equivalent to clearance class Y.

On this basis, the third quarter 1981 sample contains 72% soluble uranium compound and 28% insoluble uranium compound.

N. Cooke and F. B. Holt (1) studied the solubility of some known uranium compounds in simulated lung fluid and classified them into solubility classes similar to clearance classes used by ICRP.

DETERMINATION OF The Solutility and	Dissolution Rate Classification of	1	
Air-Borne Particulates in Simul			208-EC-22
In-Vitro	Air-Borne Particulates	Poge	12 of 23

Class D

Class W

Class Y

Very soluble	Moderately soluble	Relatively insoluble
0-10 days	>10-100 days	>100 days
Uranium hexafluoride	Ammonium diuranate	Uranium dioxide
Uranyl fluoride	Uranium trioxide	Uranium octoxide
Jranyl nitrate	Uranium tetrafluoride	

This is a very general classification, as Steckel and West $(\underline{12})$ showed that the solubility in simulated lung fluid of a number of uranium process materials and laboratory-prepared oxides was related to the effect of thermal history on chemical composition and particle size. So, the solubility of a uranium compound is very much dependent on the condition (physical conditions) under which it has gone through.

Seveninivish of The Solubility and Di	ssolution Rate Classification of	1
Air-Borne Particulates in Simulated Lung Pluid		KMTC - 208-EC-22
In-Vitro	Air-Borne Particulates	Poge 13 of 23

REFERENCES :

- N. Cooke and F. B. Holt, "The Solubility of Some Uranium Compounds in Simulated Lung Fluid," Health Physics, 27: 69-77 (1974).
- Task Group on Lung Dynamics for Committee II of the International Commission on Radiological Protection, "Deposition and Retention Models for Internal Dosimetry of the Human Respiratory Tract," Health Physics, 12: 173-206 (1966).
- P. E. Morrow, F. R. Gibb, and L. Johnson, "Clearance of Insoluble Dust from the Lower Respiratory Tract," Health Physics, 10: 543-555 (1964).
- J. R. Houston, D. L. Strenge, and E. C. Watson, DACRIN A Computer Program for Calculating Organ Dose from Acute or Chronic Radionuclide Inhalation. BNWL-B-389, Pacific Northwest Laboratory, Richland, WA (1975).
- D. L. Kalkwarf, Solubility Classification of Airborne Uranium Products from LWR-Fuel Plants, NUREG/CR-1428, PNL-3411, Pacific Northwest Laboratory, Richland, WA 99352.
- F. J. Miller, D. E. Gardner, J. A. Graham, R. E. Lee, Jr., W. E. Wilson, and J. D. Bachmann, "Size Considerations for Establishing a Standard for Inhalable Particles," J. Air Poll. Control Assoc., 29: 610-615 (1979).
- O. R. Moss, "Simulants of Lung Interstitial Fluid," Health Physics, 36: 447-448 (1979).
- N. A. Dennis, "Dissolution Rates of Yellow Cake in Simulated Lung Fluid," M. S. Thesis, University of Pittsburgh, Pittsburgh, PA (1979).
- KM-TC-125-U-17, Determination of Uranium in Water, Technical Center Report, TCR-77003, Dec. 31, 1977, Kerr-McGee Corporation.
- T. T. Mercer, "On the Role of Particle Size in the Dissolution of Lung Burdens," Health Physics, 13: 1211-1221 (1967).
- Y. B. Bard, "Nonlinear Paramete: Estimation," Academic Press, New York, 1974, Sec. A-3, Programmer: R. F. Williams, Oct. 76.
- L. M. Steenel and C. M. West 1968 Union Carbide Corporation, Nuclear Division, Report Y-1544-A.

```
KMTC-208-EC-22
Page 14 of 23
```

```
1010**RUNH*= (ULIB) RLIB;LIBRARY/APPLIB,R
 1020 PARAMETER KP=14, MOVF=2, MDLP=1, NP=5, NXNF=(NF*(NF+1))/2, MXMP=(MDLP*(MDLP+1))/2
 1022 CALL SR4J ; CALL CALLSS('REMO CLEARFILES\')
 1024 CALL ACCESS('AF, /ASOK*33*#',$99)
 1030 IMPLICIT REAL*8(A-H+0-Z)
 1040 REAL*8 P(NF), DELP(NF), PFD(NF), LB(NF), UB(NF)
 1050 REAL*8 F(MDLP), FF(MDLP), FZDZ(MDLP)
 1060 REAL *8 Z(MQVF), ZM(MQVF), R(MQVF), DELZ(MQVF), ZWORK(MQVF)
1070 REAL *8 FFTWFF(NXNF), FT(NXNF), FFTWFI(NXNF)
 1080 REAL*8 FP(MDLP,NF),FZ(MDLP,MQVP),FPTW(NF,MDLP)
1090 REAL #8 W(MXMF)
1100 REAL*8 ZMK(MQVP,KP),RK(MQVP,KP),DELZK(MQVP,KP)
1110 REAL*8 FK(MDLF,KF),FZDZK(MDLF,KF)
1120 REAL*8 FPK(MDLP, NP, KP), FZK(MDLP, MQVP, KP), WK(MXMP, KP)
1130# definitions
1140* NF = N = # of constants to evaluated
1150*
                     P(N), DELP(N), PFD(N)
1160* MDLP = MDL = # of model equations.
1170*
                     F(MDL), FF(MDL), FZDZ(MDL)
1180* MQVP = MQV = # of measured values
1190*
                     Z(MQV), ZM(MQV), ZWORK(MQV), DELZ(MQV), R(MQV)
1200* NXN<sup>5</sup> = NXN = correlation matrix = N*(N+1)/2
1210*
                     FFTWFF(NXN), FT(NXN), FPTWFI(NXN)
1220* MXMP = MXM = weighting matrix = MDL*(MDL+1)/2
1230*
                      W(MXM)
1240* limits sum MQV 2 sum MDL > N
1250* Jacobian differentials FZ(MDL,MQV), FF(MDL,N)
1260 CALL FPARAM(1,132)
1270 CALL EXPDAT(KF, MOVF, MDLF, NF, MXMF, ZM, R, DELZ, F, FZDZ, FF, FZ, W,
12808
                       ZMK, RK, DELZK, FK, FZDZK, FFK, FZK, WK)
1290 KXF=KF ; N=NF ; MDL=MDLP ; MQV=MQVP ; NXN=NXNF ; MXM=MXMP
1300 SSND = 1D-4
1310* start data preparation
1320 READ(33,10)(P(I),I=1,NP)
1330 READ(33,10)(LE(I),I=1,NF)
1340 READ(33,10)(UR(I),I=1,NP)
1350 10 FORMAT(V)
1360 DO 20 JK=1+KF
1370 READ(33,10)ZM(1),R(2),ZM(1),R(1)
1380 CALL WEXFD(JK, ZM, R, DELZ, F, FZDZ, FF, FZ, W)
1390 20 CONTINUE
1400 CALL DATIM(AD8, HRS)
1410 WRITE(06,5000)HRS, ADS
1420 5000 FORMAT(1H1,/////,* Correlation run*, S0X, F10.2, 3X, A8)
1430 CALL LUEBRI(ITER: KXP: MQVP: MDLP: NF: NXNP: MXMP: F: FF: Z: ZM. ZWORK: F: R.W.
1440%FZ,FF,FPTW,DELF,DELZ,FZDZ,FFTWFF,FT,FFTWFI,SSND+FFD,LB,UB)
1450 97 STOP ; END
1460*<<<
1470 SUBROUTINE EXPLAT(KP, MOVP, MDLP, NF, MXMP, ZM, R, DELZ, F, FZDZ, FF, FZ, W,
14802
                       ZMK+RK+DELZK+FK+FZDZK+FFK+FZK+WK)
1490 REAL#8 ZM(MOVF), R(MOVF), DELZ(MOVF)
1500 REAL#8 F(MDLF),FZDZ(MDLF)
1510 REAL*8 FP(MDLP,NP),FZ(MDLP,MQVP),W(MXMP)
1520 REAL*8 ZMK(MQUP,KP),RK(MQUP,KP),DELZK(MQUP,KP)
```

KMTC-208-EC-22 Page 15 of 23

```
1540 REAL *8 FFK(MDLP,NF,KP),FZK(MDLP,MQVF,KF),WK(MXMF,KF)
1550 RETURN
1560 ENTRY REXFD(K, ZM, R, DELZ, F, FZDZ, FF, FZ, W)
1570 DO 10 I=1, MQVF
1580 ZM(I) = ZMK(I,K)
1590 R(I) = RK(I,K)
1600 \text{ DELZ(I)} = \text{DELZK(I,K)}
1610 DO 10 J=1, MDLF
1620 \ 10 \ FZ(J,I) = FZK(J,I,K)
1630 DO 12 J=1, MDLF
1640 F(J) = FK(J,K)
1650 FZDZ(J) = FZDZK(J,K)
1660 DO 12 I=1,NF
1670 \ 12 \ FF(J,I) = FFK(J,I,K)
1680 DO 14 I=1,MXMF
1690 \ 14 \ W(I) = WK(I,K)
1700 RETURN
1710 ENTRY WEXPD(K, ZM, R, DELZ, F, FZDZ, FF, FZ, W)
1720 DO 20 I=1, MOVF
1730 \text{ ZMK}(I,K) = \text{ZM}(I)
1740 \ RK(I,K) = R(I)
1750 \text{ DELZK}(I \cdot K) = \text{DELZ}(I)
1760 DO 20 J=1,MALF
1770 \ 20 \ FZK(J,I,K) = FZ(J,I)
1780 DO 22J=1, MDLF
1790 FK(J,K) = F(J)
1800 FZDZK(J,K) = FZDZ(J)
1810 DO 22 I=1,NP
1820 22 FPK(J,I,K) = FP(J,I)
1830 DO 24 I=1, MXMF
1840 24 WK(I,N) = W(I)
1850 RETURN
1860 END
1870* 344
1880 SUBROUTINE MODEL (MQV, Z, MDL, F, N, F)
1890 REAL*8 2(1)+F(1),F(1),X,XZ
1900 XZ = -.693147181D0*Z(2)
1910 X = ODO ; IF(XZ/F(2).GT.-S7DO) X = DEXF(XZ/F(2))
1920 F(1) = -Z(1) + P(1)*X
1930 X = 0D0 ; IF(XZ/P(4),GT,-87D0) X = DEXF(XZ/P(4))
1940 F(1) = F(1) + F(3)*X + (1D0 - P(1) - P(3))*DEXF(XZ/P(5))
1950 RETURN ; END
1960*<<<
1970 SUBROUTINE LUEBRI(ITER, KF, MOVF, MDLF, NF, NXNF, MXMF, F, FF, Z, ZM, ZWORK, F, R, W,
1980%FZ,FF,FPTW,DELF,DELZ,FZDZ,FFTWFF,FT,FFTWFI,SSND,PFD,LB,UE)
1990 IMPLICIT REAL*8(A-H,O-Z)
 2000 REAL#8 F(NP), DELP(NP), PFD(NP), LB(NP), UP(NP)
 2010 REAL*8 F(MDLP),FF(MD)P),FZDZ(MDLP)
 2020 REAL*8 Z(MOVP), ZM(MOVP), R(MOVP), DELZ(MOVP), ZWORK(MOVP)
 2030 REAL*8 FPTWFP(NXNP),FT(NXNP),FPTWFI(NXNP)
 2040 REAL*8 SSND, LAM, E(24)
 2050 REAL#8 FP(MDLP,NP),FZ(MULF,MOVP),FPTW(NP,MDLP)
```

1530 REAL*8 FK(MDLP,KP),FZDZK(MDLP,KP)

```
2060 REAL *8 W(MXMF)
 2070 DIMENSION JJJ(24)
 2080 STEP1 = 100+SSND
 2090 STEP2 = 1D0-SSND
 2100 5002 FORMAT(V)
 2110 IDEM = 1 ; EAMACT = .100
 2120 KXF=KF ; N=NF ; MDL=MDLP ; MQV=MQVF ; NXN=NXNF ; MXM=MXMP
 2130*
 2140*
         start iteration loop
 2150*
 2170 1000 WRITE(06,900)
 2180 900 FORMAT(" Input - IT, IDEM, GAMACT")
 2190 READ(05,5002) IT, IDEM, EAMACT
2200 DO 90 JN=1+N
2210 \text{ FFD(JN)} = 000
2220 DO 90 IN=1, JN
2230 KK = JN*(JN-1)/2 + IN
2240 90 FPTWFP(KK) = 000
2250*
2260*
         start experimental data loop 1
2270*
2280 DO 400 JK=1,KXF
2290 CALL REXFD(JK, ZM, R, DELZ, F, FZDZ, FF, FZ, W)
2300 DO 100 JQ=1, MQV
2310 IF(IDEM.EQ.1)DELZ(JQ) = ODO
2320 \ 100 \ Z(JQ) = ZM(JQ) - DELZ(JQ)
2330 140 CALL MODEL(MQV,Z,MDL,F,N,F)
2340* COMMENT
2350 STEP = STEP1
2360 \text{ DENOM} = 100
2370 UO 142 JM=1,MDL
2380 DO 142 JO=1, MOV
2390 142 FZ(JM, JQ) = 000
2400 144 DO 146 JO=1.MOV
2410 146 ZWORK(JQ) = Z(JQ)
2420 DO 143 JA=1, MAY
2430 IF(J0.GT.1)ZWORK(J0-1) = Z(J0-1)
2440 \text{ ZWORK}(J0) = \text{STEF}(J0)
2450 IF(Z(JQ).E0.0D0) ZWORK(JQ) = STEP - 1D0
2460 CALL MODEL (MRV, ZWORK, MDL, FF, N.F)
2470 A = Z(JQ)
2480 IF(A.EQ.000) A=100
2490 DO 148 JM=1, MDL
2500 148 FZ(JM, JQ) = ((FF(JM)-F(JM))/SSND/A - FZ(JM, JQ))/DENOM
2510 IF(STEP.NE.STEP1) GOTO 150
2520 STEP = STEP2
2530 \text{ DENOM} = -210
2540 GOTO 144
2550 150 CONTINUE
2560*
         start FF calc.
2570 STEP = STEP1
2580 \text{ DENDM} = 100
2590 DO 152 JM=1,MDL
```

KMTC-208-EC-22 Page 17 of 23

1

```
2600 DO 152 JN=1.N
2610 152 FP(JM, JN) = 0D0
2620 154 DO 156 JN=1,N
2630 156 DELP(JN) = P(JN)
2640 DO 158 JN=1,N
2650 IF(JN.GT.1) DELP(JN-1) = P(JN-1)
2660 DELP(JN) = STEP*P(JN)
2670 IF(P(JN).EQ.ODO) DELP(JN) = STEP - 100
2680 CALL MODEL (MQV, Z, MDL, FF, N, DELP)
2690 A = P(JN)
2700 IF(A.EQ.000)A=100
2710 DO 158 JM=1,MDL
2720 158 FP(JM, JN) = ((FF(JM)-F(JM))/SSND/A - FP(JM, JN))/DENOM
2730 IF(STEP.NE.STEP1) GOTO 160
2740 STEP = STEP2
2750 \text{ DENOM} = -200
2760 GOTO 154
2770 160 CONTINUE
2780*
          start W calc.
2790 DO 200 JM=1, MDL
2800 DO 200 IM=1, JM
2810 KK = JM*(JM-1)/2 + IM
2820 W(KK) = 0D0
2830 DO 200 JQ=1,MQV
2840 200 W(KK) = W(KK) + FZ(IM, JQ)*FZ(JM, JQ)*R(JQ)**2
2850 IF(MDL.EG.1)W(1) = 1D0/W(1)
2860 IF(MDL.GT.1)CALL DSINV(W, MDL, 1.E-12, IER)
2870 DO 210 JM=1, MDL
2880 DO 210 JN=1,N
2890 \text{ FPTW}(JN \cdot JM) = 000
2900 DO 210 IM=1.MUL
2910 KK = MAX(JM*(JM-1), IM*(IM-1))/2 + MIN(JM, IM)
2920 210 FFTW(JN, JM) = FFTW(JN, JM) + FF(IM, JN)*W(KK)
2730 DO 220 IN=1+N
2940 DO 220 JIN=1.IN
2950 KK = IN*(IN-1)/2 + JIN
2960 DO 220 IM=1. HDL
2970 220 FPTWFP(KK) = FPTWFP(KK) + FPTW(JIN, IM)*FP(IM, IN)
2980 DO 230 JM=1,MDL
2990 FZDZ(JM) = 0D0
3000 DO 240 JO=1. MOV
3010 240 FZDZ(JM) = FZDZ(JM) + FZ(JM+JQ)*DELZ(JQ)
 3020 230 FZDZ(JH) = FZDZ(JH) + F(JM)
 3030 DD 260 JN=1 . N
 3040 DO 260 JM=1, MDL
 3050 260 FFD(JN) = FFD(JN) + FFTW(JN,JM)*FZDZ(JM)
 3060 CALL WEXPD(JK, ZM, R, DELZ, F, FZDZ, FF, FZ, W)
 3070 400 CONTINUE
 3080*
3090* end experimental data loop 1
 3100*
 3110 LAM = 000
 3120 DO 300 JIN=1,10
```

.

1

1

```
KMTC-208-EC-22
Page 18 of 23
```

```
3130 \ J0 = 0
3140 SUH = 000
3150 DO 310 JN=1+N
3160 DO 310 IN=1, JN
3170 \ J0 = J0+1
3180 FT(J0) = FPTWFP(J0)
3190 IF(JN-IN)310,312,310
3200 312 IF(JIN.GT.1) FT(J0) = FT(J0) + LAM*DABS(FT(J0))
3210 IF(JIN.GT.1.AND.FT(JO).EQ.ODO) FT(JO) = (1D0+LAM)*AVG
3220 \text{ SUM} = \text{SUM} + \text{DABS}(FT(JO))
3230 310 CONTINUE
3240 AVG = SUM/N
3250 CALL DSINV(FT,N,1.E-12,IER)
3260 IF(IER.EQ.0) GOTO 320
3270 IF(JIN, EQ.1) LAM = 5.D-3
3280 300 LAM = LAM#10D0
3290 RETURN
3300 \ 320 \ J0 = 0
3310 DO 330 JN=1,N
3320 DELP(JN) = 000
3330 DD 330 IN=1, JN
3340 JO = JO+1
3350 FPTWFI(JO) = FT(JO)
3360 DELP(JN) = DELP(JN) + FPTWFI(J0)*PFD(IN)
3370 IF(JN.NE.IN) DELP(IN) = DELP(IN) + FPTWFI(J0) * PFD(JN)
3380 330 CONTINUE
3390 GAMMAX = 1038
3400 CALL GRADPJ(LB, UB, P, FPTWEI, DELP, E, JJJ(N+1), JJJ(2*N+1), JJJ(3*N+1), N, GAMMAX, IF
3410 GAMACT = MIN(EAMACT, GAMMAX)
3420 SUMDEL = 0D0
3430 DO 340 JN=1,N
3440 BASE = P(JN)
3450 IF(BASE.EQ.ODO) BASE = DELP(JN)
3460 IF(BASE, EQ. 000) BASE = 100
3470 SUMDEL = SUMDEL + (DELP(JN)/BASE) ##2
3480 IF(IT.E0.1) WRITE(06,5008) JN, P(JN), DELP(JN)
3490 \ 340 \ P(JN) = P(JN) - GAMACT*DELP(JN)
3500 SUMDEL = DSORT(SUMDEL/N)
3510 WRITE(06,5004) IT, SUMDEL, GAMACT
3520 5004 FORMAT(I4, * SUMDEL = *, IPE15.8, *, GAMACT = *, E15.8)
3530 WSUM = ODO ; SUMSQ = ODO ; RUMSQ = ODO
3540*
3550*
        start experimental data loop 2
3560*
3570 DO 410 JK=1,KXF
3580 CALL REXPD(JK, ZM, R, DELZ, F, FZDZ, FP, FZ, W)
3590 DO 350 JN=1 , N
3600 DO 350 JM=1+MDL
3610 350 FZDZ(JM) = FZDZ(JM) - GAMACT*FP(JM,JN)*DELP(JN)
3620 DO 360 JO=1, MOV
3630 \text{ ZWORK}(JQ) = 000 \text{ ; RUMSQ} = \text{RUMSQ} + 100/(R(JQ)) * * 2
3640 DO 370 JM=1,MDL
3650 DO 370 IM=1, MDL
```

KMTC-208-EC-22 Page 19 of 23

```
3660 KK = MAX(JM*(JM-1), IM*(IM-1))/2 + MIN(JM, IM)
3670 370 ZWORK(JQ) = ZWORK(JQ) + FZ(JM,JQ)*W(KK)*FZDZ(IM)
3680 \text{ DELZ(JQ)} = ZWORK(JQ) * R(JQ) * 2
3690 IF(IT.LE.O.DR.ABS(SUMDEL).LT.1.D-4) WRITE(06,5008) JQ, ZM(JQ), DELZ(JQ)
3700 360 SUMSQ = SUMSQ + (DELZ(JQ)/R(JQ))**2
3710 KK = 0 ; DO 380 IM=1, MDL ; KK = KK+IM
3720 380 WSUM = WSUM + W(KK)
3730 CALL WEXPD(JK, ZM, R, DELZ, F, FZDZ, FP, FZ, W)
7740 410 CONTINUE
37508
3760*
       end experimental data loop 2
3770*
3780 SIGMA = DSQRT((SUMSQ*KXF)/(WSUM*(KXP-N)))
3790 WRITE(06,5010)SUMSQ,WSUM,SIGMA
3800 5010 FORMAT(* SUMSG = *,1FE15.8,*, WSUM = *,E15.8,*, SIGMA = *,E15.8)
3810 WRITE(06,5018)
3820 5018 FORMAT(1H ,4X,***,3X,*Parameter*,8X,*Sigma*,12X,*Delta*)
3830 KK = 0 ; DO 5006 I=1:N ; KK = KK+I
3640 FT(KK) = DSORT(SUMSO*FPTWFI(KK)/(KXP-N))
3850 WRITE(06,5008)I, P(I), FT(KK), DELP(I)
3860 5008 FORMAT(16,1PE17.8,2E17.8)
3870 5006 CONTINUE
3880 5014 FORMAT(3X, 1FE12.4, 9E12.4)
3890 IF(ABS(SUMDEL).LT.1.D-8.AND.IDEM.EQ.O.OR.IDEM.LT.O)GOTO 1002
3900 IF(IT.GT.0) GOTO 1000
3910*
3920*
      end iteration loop
3930*
3940 1002 WRITE(06,5016)
3950 5016 FORMAT(* Correlation matrix *)
3960 WSUM = WSUM/100 ; KK = 0 ; DO 1004 I=1,N
3970 KK = KK+I ; KKK = KK-I+1
3980 1004 WRITE(06,5014)(FFTWFI(J), J=KKK, KK)
3990 IF(IPERR.NE.O) RETURN
4000 1012 CONTINUE
4010 SUMSOF = 000
4020 DO 420 JK=1,KXP
4030 CALL REXPD(JK, ZM, R, DELZ, F, FZDZ, FF, FZ, W)
4040 CALL MODEL (MOV, ZM, MUL, F, N, F)
4050 DO 430 JM=1, MDL
4060 430 SUMSOF = SUMSOF + F(JM)**2
4070 420 CONTINUE
4080 PRINT, ERR', SUMSOF
4090 RETURN
4100 END
4110*<<<
 4120*SERID BRITT CLASS A NAME DSINV FORTRAN 06/01/79 17:08:06
 4130*READ DSINV FORTRAN B1 TESTLB 5/04/79 12:37
 4140%$ $1 BY: BRITT DATE: 01/04/79 FIRST ABSORB
 4150
           SUBROUTINE DSINV(A, N, EPS, IER)
 4160
          DIMENSION A(1)
 4170
         REAL #8 A.DIN,WORK
 4180
         CALL DMFSD(A, N, EFS, IER)
```

KMTC-208-EC-22 Page 20 of 23

4190		IF(IER) 9,1,1
4200	1	IFIV=N*(N+1)/2
4210		IND=IFIV
4220		DO 6 I=1 + N
4230		DIN=1.DO/A(IFIV)
4240		A(IFIV)=DIN
4250		MIN=N
4260		KEND=I-1
4270		LANF=N-KEND
4280		IF (: END) 5,5,2
4290	2	J=IND
4300		DO 4 K=1,KEND
4310		WORK=0.DO
4320		MIN=MIN-1
4330		LHOR=IFIV
4340		LVER=J
4350		DO 3 L=LANF, MIN
4360		LVER=LVER+1
4370		LHOR=LHOR+L
4380	3	WORK=WORK+A(LVER)*A(LHOR)
4390		A(J)=-WORK*DIN
4400	4	J=J-MIN
4410	5	IFIV=IFIV-MIN
4420	6	IND=IND-1
4430		DO 8 I=1+N
4440		IPIV=IPIV+I
4450		J=IPIV
4460		DO B K=I+N
4470		WORK = 0.00
4480		LHOR=J
4490		DO 7 L=K+N
4500		LVER=LHOR+K-I
4510		WORK=WORK+A(LHOR) *A(LVER)
	7	LHOR=LHOR+L
4530		A(J)=WORK
	8	J= J+K
		RETURN
4560		END
4570#<<<		
4580#SER		BRITT CLASS A NAME DMFSD FORTRAN 06/01/79 17:09:18
4590*REA		
		BY: BRITT DATE: 01/04/79 FIRST ABSORB
4610		SUBROUTINE DMFSD(A, N, EFS, IER)
4620		DIMENSION A(1)
4630		REAL #8 DPIV.DSUM.A
4640		IF(N-1) 10,1,1
		IER=0
4660		KPIV=0
4670		
4680		DO 11 K=1,N
4680		KFIV=KFIV+K
		INDEKFIU
4700		
4710		TOL=ABS(EFS*SNGL(A(KFIV)))

KMTC-208-EC-22 Page 21 of 23

4720		DO 11 I=K+N
4730		DSUM=0.DO
4740		IF(LEND) 2,4,2
4750	2	DO 3 L=1,LEND
4760		LANF=KPIV-L
4770		LIND=IND-L
4780	3	DSUM=DSUM+A(LANF)*A(LIND)
4790	4	DSUM=A(IND)-DSUM
4800		IF(I-K) 10,5,10
4810	5	IF(SNGL(DSUM)-TOL) 6,6,9
4820	6	IF(DSUM) 12,12,7
4830	7	IF(IER) 8,8,9
4840	8	IER=K-1
4850	9	DPIV=DSQRT(DSUM)
4860		A(KFIV)=DFIV
4870		DPIV=1.DO/DPIV
4880		GO TO 11
4890	10	A(IND)=DSUM*DPIV
4900	11	IND=IND+I
4910		RETURN
4920	12	IER=-1
4930		RETURN
4940		END
4950*>	>>>	
4960		SUBROUTINE GRADPJ(LB, UB, P, R, S, E, M, K, IDX, N, SSMAX, IER)
4970		IMPLICIT REAL*8(A-H, 0-Z)
47.80		REAL #8 LB(N), UB(N), F(N), R(1), S(N), E(N, 1)
4990		INTEGER M(N),K(N),IDX(N)
5000*		
5010*		LOCATE SYMMETRIC STORAGE MODE ARRAY ELEMENT FOR IS=J
5020*		
5030		LOC(I,J) = I*(I-1)/2+J
5040*		
5050*		CHECK FOR ACTIVE BOUNDS AND SET UP M AND K ARRAYS
5060*		
5070		IER=0
5080		NAC=0
5090		DO 10 I=1,N
5100		IF(DABS(F(I)-LB(I)) .LT. 1D-6) GO TO 4
5110		IF(P(I) .LT. LB(I)\$1.00001D0) GD TO 4
5120		GO TO 5
5130*		LOWER BOUND ACTIVE
5140	1	NAC=NAC+1
5150		P(I) = LB(I)
5160		M(NAC)=1
5170		K(NAC)=1
5180		IDX(NAC)=I
5190		GO TO 10
5200		5 IF(DABS(P(I)-UB(I)) .LT. 10-6) GO TO 6
5210	1.1	IF(P(I) .GT. 0.7999900*UB(I)) GO TO 6
5220		GO TO 10
5230*		UPPER BOUND ACTIVE
5240		6 NAC=NAC+1

5250 P(I) = UB(I)5260 M(NAC) = -15270 K(NAC) = 15280 IDX(NAC)=I 5290 10 CONTINUE 5300* IF IN INTERIOR, NO PROJECTION 5310 IF(NAC.EQ.0) GO TO 60 5320* 5330* CONSTRUCT TABLEAU AND FIND BINDING BOUNDS 5340* DO 20 I=1, NAC 5350 5360 E(I, NAC+1) = -S(IDX(I))5370 DO 20 J=1,NAC 5380 IF(I,GE,J) = R(I,J) = R(LOC(IDX(I),IDX(J)))5390 IF(I.LT.J) E(I,J) = R(LOC(IDX(J), IDX(I)))5400 20 CONTINUE 5410 ICOUNT=0 5420 25 ICOUNT=ICOUNT+1 5430 IF(ICOUNT.LE.10*NAC) GO TO 27 5440 IER=1 5450 RETURN 5460 27 A=1D38 DO 30 I=1,NAC 5470 5480 AI=M(I)*K(I)*E(I,NAC+1) 5490 IF(AI.GE.A) GO TO 30 5500 A=AI 5510 IR=1 5520 30 CONTINUE 5530 IF(A.GT.-10-6) GO TO 40 5540 K(IR) = -15550 CALL GJPVT(NAC; NAC+1, N, IR, IR, E, IER) 5560 IF(IER.EQ.0) GO TO 25 5570 RETURN 5580# 5590* COMPUTE PROJECTION STEP 5600* 40 DO 50 J=1,NAC 5610 5620 IF(K(J).NE.-1) GO TO 50 5630 DO 45 I=1+N 5640 IF(I.GE.IDX(J)) S(I)=S(I)+R(LOC(I.IDX(J)))*E(J.NAC+1) 5650 IF(I.LT.IDX(J)) S(I)=S(I)+R(LOC(IDX(J),I))*E(J,NAC+1) 5660 45 CONTINUE 5670 50 CONTINUE 5680* FOR BINDING BOUNDS SET STEP TO IDENTICALLY 0.0 AND LOAN 5690* LAGRANGE MULTIPLIERS(SENSITIVITY DERIVATIVES) INTO THE FIFST 5700* COLUMN OF E FOR RETURN. SET ALL OTHER MULTIFLIERS TO 0.0 5710* 5720* 5730 DO 54 I=1+N 5740 54 E(I+1)=0.0 5750 DO 55 I=1,NAC 5760 IF(K(I).NE.-1) GO TO 55 5770 S(IDX(I))=0.0

KMTC-208-EC-22 Page 23 of 23

```
E(IDX(I),1)=-E(I,NAC+1)
5780
       55 CONTINUE
5790
5800*
             MAXIMUM FEASIBLE STEP SIZE AND EXIT
5810*
5820*
       60 SSMAX=1038
5830
          DO 65 I=1 / N
5840
          IF(S(I),LT.0.0) SSMAX=DMIN1(SSMAX,-(UB(I)-P(I))/S(I))
5850
          IF(S(I).GT.0.0) SSMAX=DMIN1(SSMAX,-(LB(I)-P(I))/S(I))
5860
5870
       65 CONTINUE
          M(1)=NAC
5880
           IF (NAC.NE.O) RETURN
5890
           DO 70 I=1,N
5900
       70 E(I,1)=0.0
5910
5920
          RETURN
          END
5930
           SUBROUTINE GJPVT(N, M, ND, I, J, B, IER)
5940
           REAL #8 B(1)
5950
5960*
           PERFORM GAUSS-JORDAN FIVOT ON MATRIX B(N,M)
5970*
               STORED IN AN ARRAY WITH ROW DIMENSION ND
5980*
               USING B(I, J) AS THE PIVOT ELEMENT
5990*
               ON RETURN, IER = 0 IF NO ERROR
6000*
                           IER = -1 IF B(I,J) = 0
6010*
6020*
           REFERENCE: BARD, Y., B., "NONLINEAR PARAMETER ESTIMATION, "ACADEMIC
6030*
                         PRESS, NEW YORK, 1974, SEC. A-3
6040*
6050*
           FROGRAMMER: R.F. WILLIAMS, DCT 76
6060*
60707
              L IS A FUNCTION TO LOCATE MATRIX ELEMENT B(IROW, JCOL) IN ARRAY
$080*
6090*
           L(IROW, JCOL) = ND*(JCOL-1) + IROW
6100
           IER = 0
6110
           IF(B(L(I,J)),NE.0.D0) GO TO 10
6120
           IER=-1
6130
           GO TO 200
6140
        10 DO 100 IF=1,N
6150
           100 100 JQ=1,M
 6160
           IF(IF.EQ.I.OF.JQ.EQ.J) GO TO 100
 6170
           B(L(IP,JQ))=B(L(IP,JQ))-B(L(I,JQ))*E(L(IP,J))/B(L(I,J))
 6180
      100 CONTINUE
 6190
           DO 120 IP=1,N
 6200
           IF(IF.E0.I) GO TO 120
 6210
           B(L(IF,J)) = -B(L(IF,J))/B(L(I,J))
 6220
 6230 120 CONTINUE
 6240
           10 140 JQ=1,M
           IF(JQ.EQ.J) GO TO 140
 6250
           B(L(I, JQ)) = B(L(I, JQ))/B(L(I, J))
 6260
 6270 140 CONTINUE
           B(L(I,J))=1.D0/B(L(I,J))
 5280
 6290
       200 RETURN
 6300
           END
```

.

APPENDIX 5.2.3

*

١

ORNL'S ANALYTICAL PROCEDURE

OAK PEORE NATIONAL LABORATORY OPENATED BY MARTIN MARIETTA ENERGY EVENTSAS, NO.

POST OFFICE BOX X CAA AROLE, TENNEBBEE \$7901

January 30, 1986

Mr. Barry Zaloman U. S. Muclaar Regulatory Commission Millste Bldg., Room 517 7915 Eastern Avenue Silver Springs, Maryland 20910

Dear Barry:

I enclose the copies of the specific analytical procedures used to measure the uranium and fluoride concentrations in the specimens received after the January θ_0 1986 Sequoyah Nuclear Fuels incident.

As we discussed by telephone, we utilized our neutron activation facility to determine the uranium concentrations in the vegetation, soils, and water specimens. The limits of detection for uranium are 0.030 micrograms per cilliliter in water, and 0.020 micrograms per gram of vegetation or soil. The confidence level is 10.0 ± 0.2 micrograms of uranium.

The uranium concentrations in the two shipments of urine specimens were detarmined by ultraviolet excitation and measurement of the characteristic yellowgreen uranium fluorescence. The limit of detection for uranium is 0.005 micrograms of uranium per milliliter of urine. The relative standard deviation is 7.6 percent at the 0.05 microgram per milliliter concentration level. The detailed procedure is included, for your information, together with a letter which details the quality assurance standards measured with your specimens. Please note that we measured more standards than specimens, and that these quality assurance concentrations covered the entire range experienced in the specimens.

The procedures for measuring fluoride concentration in vegetation and water are attached. The fluoride in the soil specimens was extracted with deionized fluoride-free water, the extract was filtered to remove soil particles, and the final determination of fluoride was measured with a fluoride specific ion electrode as detailed in the surface water procedure. We analyzed a 4.0 microgram per gram National Bureau of Standards orchard grass standard as a quality assurance test and recovered 3.4 micrograms of fluoride. We also prepared duplicate ORNL standards containing 20.0 micrograms of fluoride and subjected them to the entire analytical procedure, including distillation. We recovered 20.3 and 20.9 micrograms of fluoride to demonstrate excellent recover-

As requested by Mr. Jerry Swift, I also enclose the data sheet from the NRC shipment of the COMP water, soil, and vegetation specimens received at ORNL on January 22, 1986. I have taken the information written on the plastic specimen bags and included it on the data sheet that was in the shipping carton.

KIN KERR-MCGEE CORPORATION		NTC - 149-EC-10
ANALYTICAL METHODS	133060	VISED None
Thorium, Uranium and Plutonium	TYPE OF SAMPLE Water, Soi Products of	1, Vegetation and r By-Products
Alpha Pulse Height Analysis	GEVS	CHL

Method:

Except for natural uranium or natural thorium, environmentally important concentrations of actinides cannot be readily determined by ordinary chemical means. Additionally, chemical methods do not offer any information as to the isotopic composition of the element in question. Therefore, for samples where either an isotopic composition and/ or ultra low level of detection is required, alpha pulse height analysis is the method of choice.

This method utilizes the following technologies:

- In-situ formation of calcium phosphate is used to concentrate the actinides from a large aqueous volume.
- Anion exchange is used to permit separation of the actinides into relatively pure form.
 - a) For thorium, good separation is achieved from all elements except uranium. A separation factor in the vicinity of 10⁶ is achieved for most elements, but only about 10⁴ separation is achieved between uranium and thorium. As a result, some uranium will report to the electrodeposit when analyzing for thorium in a uranium rich material such as yellow cake.
 - b) For uranium, good separation is achieved from all elements except thorium and those which form stable anion complexes with chloride (i.e., iron, aluminum, etc.). For this reason, whenever a material rich in these metals is to be analyzed for uranium, such as soils or mill tailings, the uranium must first be solvent-extracted from a nitrate or sulfate system in order to eliminate these interferences.
 - c) For plutonium, good separation is achieved from all metals except for iron from iron-rich systems (i.e., liver tissue, iron or steel and their alloys, etc.). When analyzing an iron-rich system, the ion-exchange step should be repeated once. (Note: Three evalorations to dryness with concentrated nitric acid area required to completely remove traces of chloride and fluoride from the first ion-exchange strip liquor prior to the second ionexchange. Failure to remove the chloride or fluoride will result in poor recovery of the plutonium.).
- 3. Electrodeposition is utilized to prepare an "infinitesimally thin" deposit necessary for good alpha pulse height energy resolution. A thick deposit (> ten micrograms) causes alphas, from the bottom of the deposit, to degrade (lose energy) as they pass up through the deposit. As a result, the silicon surface barrier detector will see alphas with a lower energy which will produce

DETERMINATION OF		
Therium, Uranium and Plutonium		KMTC - 149-EC-10
ANALYTICAL METHOD	TYPE OF SAMPLE Water, Soil, Vegetation	149-00-10
Alpha Pulse Height Analysis	and Products or By-Products	Page of 8

a spectral smear. Poor resolution of alpha energy will destroy the usefulness of the procedure. Alpha resolution should be < 50 KeV FWHM (full width at half maximum) and any spectra with an alpha resolution of > 100 KeV FWHM should result in rejection of the analysis.

Range and Precision:

Amounts of Th, U, or Pu as low as one femtocurie $(1 \times 10^{-15} \text{ curie or } 3.7 \times 10^{-5} \text{ disintegrations per second})$ can be detected over a counting time of 1000 minutes using silicon surface barrier alpha detectors in conjunction with a multichannel analyzer. The precision attained is a function of the activity in the sample and the recovery of the internal standard; however, typical precision, at one standard deviation, ranges from 3 to 50%.

Cost:

Cost (manhours) = k + nt

Where:

k = 2 manhours; t = 2 manhours, and n = number of samples.

Apparatus:

- Silicon surface barrier detector(s) (Princeton Gamma Tech, 300-25-100) is used in conjunction with a multichannel analyzer.
- 2. Millipore filter equipment including 0.45µ Millipore Filter.

3. Magnetic stirring hotplate(s).

- Ion-exchange column(s), Bio-Rad or equivalent, 1.5 cm I.D. X 15 cm in length, and accessories.
- 5. Ripple-free D. C. power supply capable of at least 1.5 amp @ 10 volts.
- Electrodeposition cells (Talvite, Aral. Chem. 1972) and polished stainless steel discs.
- 7. Assorted glass and plastic ware common to analytical laboratories.

Reagents:

- 1. a) Concentrated sulfuric acid (H₂SO₄).
 - b) Concentrated hydrochloric acid (HCl).
 - c) Concentrated nitric acid (HNO3).
 - d) Concentrated hydrofluoric acid (HF).
 - e) Concentrated phosphoric acid (H₃PO₄).

January 30, 1985

Mr. Barry Zalcman

We are pleased that the Martin Marietta Energy Systems staff and facilities were able to respond to your needs during the incident.

-2-

Please notify me if we can be of further assistance.

Sincerely,

seph N. Stewart

J. H. Stewart, Jr. CAFA Group Leader Analytical Chemistry Division Oak Ridge National Laboratory

JHS:sdc

Enclosures

cc: W. R. Laing W. D. Shults J. R. Stokely, Jr.

NVIE	UN	IMENTAL ANALYSIS PROCEDURE	EC 410
	RIDE	(IN VEGETATION), ION-SELECTIVE ELECTRODE METHOD	PAGE 1 OF 3
1.0	5000	E AND APPLICATION	
	1.1	This method is applicable to vegetation such as g and forest floor litter.	grasses, leaves,
	1.2	The range of the method may be varied by sample a dilution. For a 3 gram sample, the lowest fluori	the and/or

2.0 SUMMARY OF METHOD

2.1 The sample is air dried, puiverized in a Wiley mill and digested in potassium hydroxide. Phosphoric acid is added and the fluoride determined by the ion-selective electrode method using the standard addition technique.

3.0 INTERFERENCES

3.1 The strong phosphoric acid solution in which the fluoride is measured eliminates the interference of complexing agents such as aluminum, iron, thorium, and silicon compounds.

4.0 SAMPLE HANDLING AND PRESERVATION

4.1 No special regulrements.

5.0 APPARATUS

- 5.1 Orion Model 801 digital pH/mV meter or similar instrument for use with ion-selective electrodes.
- 5 2 Fluoride ion-selective electrode
- 5.3 Reference electroda
- 5.4 pll Meter
- 5.5 Wiley Mili, Intermediate Model
- 5.6 Magnetic stirrers.

ANTROVED BY CO		DP T. Kuranochi	Monow	R.W.	4 -12
PCOP A.E. Simmerca ORNL FC 410	 FC ATO	teo at	Simmon		

EC	410			
DATE		1977		151
BUPERSE				
PAGE	-			
	2	QF	3	-

6.0 REAGENTS

- 6.1 Fluoride Standard. 2000 mg/liter: dissolve 2.2102 grams of dried reagent grade sodium fluoride in distilled water and dilute to 500 ml. One milliter = 2.0 milligrams fluoride.
- 6.2 Hydrochloric acid, concentrated.
- 6.3 Phosphoric acid. 18 N: dilute 407 ml reagent grade phosphoric acid. 85%, to one liter with distilled water.
- 6.4 Potassium hydroxide, 2 N: dissolve 112 g reagent grade potassium hydroxide in distilled water and dilute to one liter.

7.0 PROCEDURE

- 7.1 Grind the dried sample to minus 40 mesh in the Wiley mill. Weigh a 3-g sample, or less if more than 400 micrograms of fluoride is expected, and transfer to a 150-ml beaker.
- 7.2 Add 35 ml of 2 N potassium hydroxide, cover beaker with watch glass and place on hot plate adjusted to a low setting. Also, place a beaker containing 35 ml of 2 N potassium hydroxide on hot. plate and carry through the procedure as a reagent blank.
- 7.3 Digest sample for a total of 90 minutes, rinsing the beaker walls and cover with distilled water after 30 and 60 minutes of digestion.
- 7.4 After 90 minutes, remove backers from hot plate, rinse covers and walls and cool to room temperature.
- 7.5 Add 25 ml of 18 N phosphoric acid and adjust solutions to pH 1.0 by dropwise addition of concentrated hydrochloric acid.
- 7.6 Place beaker on a magnet stirrer. Insert the fluoride ionselective electrode and the reference electrode. Allow stirred semple to equilibrate and record the millivolt reading.
- 7.7 Add 0.20 ml of the 2.0 mg/ml fluoride standard. The volume of standard added must be small so as not to significantly affect the volume of the sample solution.
- 7.8 Allow stirred sample to equilibrate again and record the second millivolt reading. Determine the fluoride content of the sample using the calculations given below.

EC 4				
July	21,	1977		
SUPER	EDES			
PAGE	3	0*	3	

- 7.9 If the sample is found to contain more than 400 micrograms of fluoride, repeat the analysis using a smaller sample.
- 8.0 CALCULATION

8.1
$$C_0 = \frac{C_4}{(antilog - 4 mV)} - 1$$

where:

Co = total micrograms fluoride in sample solution
Ca = micrograms standard fluoride added
A my = difference in potential of sample solution and fluoride
 spiked sample solution
S = 59.2, the theoretical electrode slope value for S at 25°C

8.2 Fluoride concentration, $ug/g = \frac{C_0 - B}{W}$

where:

- Co = the micrograms F in sample solution, from (8.1).
- B = the micrograms F in reagent blank, and
- H = the sample weight in grams.

9.0 REFERENCES

9.1 Galloway, H. L., Shoaf, R. E., Skaggs, C. H., "A Rapid Method for the Vetermination of Fluoride in Vegetation"; <u>American Industrial</u> Hygiene Association Journal, <u>36</u> p. 721 (1975).

1 1 1

	UNI	MENTAL ANALYS	ION		EC 133 DATE March 24. 1977
	FLUOR	RIDE. ION-SELECTIVE ELEC	TRODE METHOD		P 446
1.0	SCOPE	AND APPLICATION	:	a	
	1.1	This method is applicab surface, and saline wat	le to the mea ers, domestic	surement of and industr	fluoride in drinking, fal wastes.
-44	1.2	Concentrations of fluor measured.	ide from 0.1	up to 1000 m	ng/liter may be
	1.3	For total or total diss must be performed on sa distillation must also to be free from axtreme	mples prior t be performed	when samples	analysis. The
	1.				
2.0	SUMMA	RY OF METHOD			
	2.1	The fluoride is determine alectrode in conjunction reference electrode and scale or an ion-select for fluoride.	on with a star a pH meter h	ndard singla having an ex	junction sleeve-type panded millivolt
	2.2	The fluoride electrode across which a potentia may be represented by / SCE/.	al is develop	ed by fluori	de ions. The cell
in i Reck	2.3	Using a 300-ml sample,	the lowest co	oncentration	reported is 0.1 mg/1
3.0	SAMPL	E HANDLING AND PRESERVA	TION		
			·		
9080	3.1	No special requirement	5.		and the second
1	-	SEEDENCES			
4.0	INIER	RFERENCES			
	(4.1	Extrames of pH and the and Al ⁺³ which form co of interferences depen cations, the concentra Sample pH should be in adjustment is accompli- below) containing a st complexes aluminum (th iron.	mplexes with ds upon the c tion of fluor itially adjus shed by addit rong, chelati	fluoride int oncentration ide, and the ted to betwe ion of pH 5. Ing agent wh	terfere. The degree n of the complexing a pH of the sample. een 5 and 9. Final .0 buffer (described ich preferentially
PPROVER	-	LAPPROV			-
Y-12	Rog	WManaw oro	OF T. Kwy	wohi	
ROYED	"Di	C ATTRO	LU.R.	La	EC 133

A. . .

0

1

N L: MB 8	EC 1	33	
DATE	March		1977
SUPERI	EPES		
PAGE			a street of the second
		91	

5.0 APPARATUS

- S.! Electrometer (pH meter), with expanded mV scale, or an ion-selective meter such as the Orion Model 801.
- 5.2 Fluoride Ion Activity Electrode, such as Orion No. 94-09
- 5.3 Reference clectrode, single junction, sleeve-type, such as Orion No. 90-01, Beckman No. 40454, or Corning No. 476010.

....

80.0

- 5.4 Magnetic Mixer, Teflon-coated stirring bar.
- 5.5 Direct distillation apparatus for fluoride, as shown in Figure 1.

6.0 REAGENTS

- 6.1 Buffer solution, pH 5.0-5.5: To approximately 500 ml of distilled water in a 1-liter beaker add 57 ml of glacial acetic acid, 59 g of sodium chloride and 2 g of CDTA (1, 2 cyclohexylene dinitrilo tetracetic acid or cyclohexane diamine tetracidic acid). Stir to dissolve and cool to room temperature. Adjust pH of solution to between 5.0 and 5.5 with 5 N sodium hydroxide (about 150 ml will be required). Transfer solution to a 1-liter volumetric flask and dilute to the mark with distilled water. For work with brines, additional NaCl should be added to raise the chloride level to twice the highest expected level of chloride in the sample.
- 6.2 Sodium fluoride, stock solution: 1.0 ml = 0.1 mg F. Dissolve 0.2210 g of sodium fluoride in distilled water and dilute to 1 liter in a volumetric flask. Store in chemical-resistant glass or polysthylene. Use anhydrous NAF in preparing this solution.
- 6.3 Sodium fluoride, standard solution: 1.0 ml = 0.01 mg F. Dilute 100.0 ml of sodium fluoride stock solution (6.2) to 1000 ml with distilled water.
- 6.4 Sulfuric acid, conc.

7.0 CALIBRATION

7.1 Prepare a series of appropriate standards using the fluoride standard solution (6.3). The range of the concentration of the standards should enclose the concentration of the samples.

EC 133

HUMATA

:h 24	, 19	77		
		Colorado da	-	-
			*	
				-
	0	0#	OF	OF


7.2 Calibration of Electrometer: Proceed as described in (8.2). Using semilogarithmic graph paper, plot the concentration of fluoride in mg/liter on the log axis vs the electrode potential developed by the standard on the linear axis, starting with the lowest concentration at the bottom of the scale. Calibration of an ion-selective meter: Follow the directions of the manufacturer for the operation of the instrument.

8.0 PROCEDURE

- 8.1 Distillation
 - 8.1.1 Place 400 ml distilled water in the distilling flask and carefully add 200 ml conc. H₂SO₄. Swirl until the flask contents are homogeneous. Add 25-35 glass beads and connect the apparatus as shown in Figure 1 making sure all joints are tight. Begin heating slowly at first, then as rapidly as the efficiency of the condenser will permit (the distillate must be cool) until the temperature of the flask tontants reaches 180°C. Discard the distillate. This process serves to remove fluoride contamination and to adjust the acid-water ratio for subsequent distillations.
 - 8.1.2 After cooling the acid mixture remaining from the steps outlined in 8.1.1 or previous distillations, to 120°C or below, add 250 ml of sample, mix thoroughly, and distill as before until the temperature reaches 180°C. To prevent sulfate carry-over, do not permit the temperature to exceed 180°C.
 - 8.1.3 Add silver sulfate to the distilling flask at the rate of 5 mg per milligram of chloride when high-chloride samples are distilled.
 - 8.1.4 Continue to use the sulfuric acid solution in the flask until the contaminants from the water samples accumulate to an extent that recovery is affected or interferences appear in the distillate. Check suitability of the acid periodically by distilling standard fluoride samples. After the distillation of high-fluoride samples, flush the still with 300 ml distilled water and combine the two fluoride distillates. If necessary, repeat the flushing operation until the fluoride content of the distillates is at a minimum. Include the additional fluoride recovered with that of the first distillation. After periods of inactivity, similarly flush the still and discard the distillate.

232

EC 133



Internal Correspondence

MARTIN MARIETTA ENERGY SYSTEMS, INC.

MARTIN MARIETTA

January 27, 1986

J. H. Stewart, 4500-5, MS-140, X-10 (4-4895)

Analysis of NRC Urine Samples from Kerr-McGee Plant

Two sets of urine samples from Kerr-McGee were analyzed for total uranium by the Y-12 Plant Laboratory. The first set received by the laboratory on January 8, 1986, totaled 218 samples, and the second set received on January 10, 1980, totaled 1/5 samples. For both sets, blanks and appropriate concentrations of calibration standards, internal controls and external controls were used. In addition, approximately every fourth sample was spiked with suitable aliquots of prepared standard solutions to ensure that no interfering substances were encountered in the urine. All samples and spikes were run in duplicate, with reruns made as necessary to confirm results.

Three sets of calibration standards were used to cover a range of 0 to 10 μ g/mL: 0 to 0.05 μ g/mL; 0 to 0.1 μ g/mL; and 0.05 to 10.0 μ g/mL. All samples were read with the appropriate calibration standards, along with internal and external controls. Calibration standards, internal controls, and spikes were prepared and diluted using 0.05 μ g/mL and 10.0 μ g/mL stock solutions. Listed below are the calibration sets and the controls that were used with each set. All values and readings are given in μ g/mL.

Α.	Calibration Standards	Internal Controls	External Controls
	0 and 0.050 0 and 0.1 0.050 and 10.0	0.025 0.025, 0.050 0.1	Blank, 0.020, 0.033 Blank, 0.020, 0.033 0.631, 1.26, 5.05
в.	Internal Controls	Total No. Read	Average Readings (95% Confidence Level)
	0.025 0.050 0.1	159 57 48	0.024 ± 0.006 0.052 ± 0.008 0.104 ± 0.020
	External Controls*		
	0.020 0.033 0.631 1.26 5.05	149 144 48 38 48	$\begin{array}{c} 0.020 \pm 0.006 \\ 0.033 \pm 0.008 \\ 0.660 \pm 0.052 \\ 1.30 \pm 0.19 \\ 4.98 \pm 0.71 \end{array}$

J. H. Stewart January 27, 1986 Page 2

*External control samples were prepared for the Plant Laboratory by the Y-12 Quality Division. The standard values for the external controls were calculated by the Quality Division.

Concentrations and amount of spikes used depended upon the readings of the sample. The samples were spiked so that the spike of sample readings were approximately twice that of the sample. For this purpose, various dilutions were prepared using the stock uranium solutions. The average recovery of all the spiked samples was $99 \pm 24\%$.

K. a. Puchells

K. D. Pickell, 9995, MS-1, Y-12 (4-2958) - NoRC

KDP:dh

cc: J. G. Dorsey, 9995, MS-2
 R. J. McElhaney, 9995, MS-2
 L. E. White, 9995, MS-1

INDUSTRIAL HYGIENE ANALYSIS PROCEDURE

day ist		UNION CARBIDE CORPORA		HUMBER IHA-490 DATE June 15, 1983
URANI	UM IN	URINE, FLUOROMETRIC METHO	D	Nov. 27, 1979
Anal	yte:	Uranium	Method No.:	IHA-490
Matr	ix:	Urfne	Range:	0.005 to 10.0 yg/mL
Proc	edure:	Fluorometric after fluoride pellet fusion	Precision:	15% RSD at 0.01 ug/mL 7.5% RSD at 0.05 ug/mL
1.0	Prin	tiple of the Method	*	
	1.1		is placed in a san infrared lam	small flat platinum dish and
	1.2		in the dish and	fused by controlled heat.
	1.3	The yellow-green uranium ation is measured using -fluorescence of known ur		fuced by ultraviolet fllumin- er and compared to the
2.0	Rance	and Sensitivity		
	2.1		ge of the method 10 ug/mL by usir	is 0.005 to 1.0 ug/mL. The higher concentration
	2.2	The limit of detection is	s 0.005 ug/mL.	
3.0	Inter	ference		
	3.1	No interferring substance concentrations normally e	s are encountere experienced.	d în urin at
4.0	Precis	ston and Accuracy		말 하는 것은 가슴 같아?
	4.1	The precision (limit of e ±30% (15% relative standa (7.5% RSD) at 0.05 µg/mL.	FO OPVIATION PS	confidence level is D) at 0.01 µg/mL and ±15%
POVED	LAD	AMS Ellaney AND NOVES	The state of the s	ALLA APPROVED BY UCCND-COORD .
	I'M x	2 Beam er oreon	IM N.J. Du	man culler
PEOP	in for	E. C. Serrer and APPROVED	Luon.m. Fu	IHA-490
	TIM	M. Jurk OV CANL	I'm Jones A. En	1HA-490

236

RUNSEA	IHA-4	190		
BATE	June	15,	1983	
SUPEPS		27.	1979	
PAGE	2	0=	5	

4.2 The accuracy of the method depends on how well the matrix of the calibration standards matches that of the sample being analyzed. No significant bias has been identified.

5.0 Advantages and Disadvantages of the Method

14 de 14 de

- 5.1 Advantages over previous methods include simplicity, speed, and sensitivity.
- 5.2 One disadvantage is that, without chemical separation of the uranium, the accuracy of the measurement depends on the amount of dissolved solids in the sample, and on how well this matches the calibration standard.

6.0 Apparatus

- 6.1 Fluorophotometer, modified Oak Ridge National Laboratory Model Q-1165, or equivalent.
- 6.2 Furnace, induction; or burner, modified Fletcher.
- 5.3_ Dishes, platinum fusion, 772-inch diameter, 1/16-inch lip, 1/8-inch deep.
- 6.4 Fusion dish holder.
- 6.5 Pelletizer, delivering 0.3 or 0.6-g pellet of flux.
- 6.6 Infrared lamp and stand, 375-w.
- 6.7 Pipet, 100-or 200-uL, semiautomatic.

7.0 Reagents

7.1 Sodium fluoride - 2% lithium fluoride flux: Put 1359 g (3 pounds) of fluorometric grade sodium fluoride and 27 grams of lithium fluoride into a clean 5-pound glass reagent jar. Place the jar on the mechanical roller and roll for 48 hours.

NUNSER	IHA-4	90	1
DATE	June	15,	1983
SUPERSI	NOV.	27,	1979
PAGE	3	0#	5

- 7.1.1 Alternate flux, sodium fluoride: the lithium fluoride additive is unnecessary when induction heating is used for the fusion.
- 7.2 Stock uranium solution: Dissolve 117.9 mg of natural U₃O₈ in 8 M HNO₃ in a 100-mL beaker. Transfer the solution to a 100-mL volumetric flask using 2 M HNO₃ to transfer. Dilute to volume with 2 M HNO₃ and mix thoroughly. (1 mL = 1 mg U.)
- 7.3 Calibration Standards:
 - 7:3.1 Dilute suitable aliquots of the stock uranium solution (7.2) with acidified, blanked, and blended urine to give at least two appropriate standards within the operating range of 0.005 to 10 µg/mL. The standards--usually a high and a low-- should bracket the working range needed.

8.0 Procedure

- 8.1 Cleaning Fusion Dishes
 - 8.1.1 Adequate cleaning of the platinum fusion dishes is critical to the prevention of cross-contamination. Thoroughly mix 2 tablespoons of boric acid with 150 ml of water in a Pyrex dish. Carefully add 100 mL of conc sulfuric acid and boil the platinum dishes gently in this solution for 2 hours. Drain and rinse the dishes with distilled water.
 - 8.1.2 Cover the dishes with 1:1 HNO₃ and slowly add 30% hydrogen peroxide until mild effervescence occurs. Soak overnight.
 - 8.1.3 Dishes can be tested for cleanliness in the fluorophotometer. Dishes that resist cleaning by the method described in 8.1.1 and 8.1.2 can be cleaned by fusing potassium pyrosulfate pellets in them and then soaking them in 1:1 HNO₃ overnight.
- 8.2 Pelletizing and Fusion
 - 8.2.1 Place the required number of platinum fusion dishes in the dish holder or on a nichrome screen.

IHA-490

	IHA-4	190			
DATE	June	15,	1983		
SUPER	NOV.	27,	1979		
PAGT				5	

- 8.2.2 Swirl each sample bottle to mix the sample and to suspend all solids. Process a blank and the calibration standards with each set of samples.
- 8.2.3 Quickly insert the semiautomatic pipet into the sample bottle and withdraw an appropriate aliquot of solution (usually 100 uL for a 0.3-g pellet or 200 uL for a 0.6-g pellet); discharge the contents of the pipet into a platinum dish.
- 8.2.4 Place the dish holder or screen and the platinum dishes under the infrared lamp and svaporate all samples to dryness. Then remove the holder or screen from under the lamp.
- 8.2.5 With the pelletizer, place a 0.3- or 0.6-g pellet of the fluoride flux into each platinum fusion dish.
- 8.2.6 Fusion using a modified Fletcher Burner.
 - 8.2.6.1 Place the screen over the fusion burner; fuse for 2 minutes at approximately 1100°C, then reduce to 950°C and fuse for 1 minute. (An optical pyrometer will provide temperature estimates.)
 - 8.2.6.2 Remove the screen containing the dishes from the burner and allow to cool.
- 8.2.7 Alternate fusion, using an induction Furnace:
 - 8.2.7.1 Place four dishes at a time in the induction furnace coil.
 - 8.2.7.2 Apply 600 ma current for 60 seconds.
 - 8.2.7.3 Transfer the dishes back to the holder to cool. (CAUTION: Use forceps.)
- 8.3 Fluorescence Measurement
 - 8.3.1 Place the platinum dishes containing the blank, uranium standards, and sample melts on the turn-table of the fluorophotometer.

HA-490

NUMBE	IHA-	490	
DATE	June	i5,	1983
UPERI	NOV.	27,	1979
PAGE	5		

- 8.3.2 Move the blank under the ultraviolet light and adjust the zero control knob to make the instrument read zero.
- 8.3.3 Move the high uranium standard melt under the ultraviolet light and adjust the calibration control until the instrument reads the known value of the standard.
- 8.3.4 Move the low uranium standard melt under the ultraviolet light. It should read the known value of the standard if the instrument is properly calibrated. The fluorophotometer is now calibrated for direct readout of uranium in the working range, which could extend from 0.005 to 10 ug/mL. Readout may be in appropriate units other than ug/mL.
- 8.3.5 Move the sample pellets under the ultraviolet light of the fluorophotometer and read the uranium content of the aliquots.

9.0 Calculations

- 9.1 If the readings are in ug/mL, no further calculations are necessary.
- 9.2 If the readings are in ug or ng, use the following formulas:

U.
$$ug/mL = \frac{vg}{A}$$

U. $ug/mL = \frac{ng}{A-1000}$

where: A = volume of sample (see 8.2.3) pipetted onto the pellet. mL.

9.3 If the ug U/day for an employee is to be reported, use the following formula:

where: V = volume excreted, mL T = time interval, hr 24 = hr/day

10:0 Reference

10.1 Centanni, F. A., Moss, A. M., and DeSesse, M. A.; "Fluorometric Determination of Uranium": Anal. Chem., 28. p. 1651 (1956).

240

THA-490

5 11924 5-22-80

Determination of Uranium by Neutron Activation and Delayed Neutron Counting

1.0 SCOPE

This method is applicable to the determination of ²³⁵U in liquid or solid samples. Total uranium analyses can be obtained by the method if the isotopic composition of uranium in the sample is known.

2.0 PRINCIPLE

This procedure assumes familiarity with the delayed neutron counting facility and pneumatic tube control panel at the Oak Ridge Research Reactor (ORR). After irradiation time and aliquot size are established, the sample is irradiated and the neutrons emitted following decay of certain short-lived fission products of 235 U are counted. A 235 U standard that corresponds to the configuration of the unknown sample is irradiated and counted exactly in the same manner as the sample. An empty rabbit is also irradiated and the neutron count rate is subtracted from the count rates of the sample and standard. The resulting count rates and weights of the standard and sample are compared to obtain the uranium content of the sample.

3.0 APPARATUS

Analytical balance Micro-pipettes Heat lamp Rabbits and inserts Soldering iron

4.0 SAFETY

The hood that contains the tabbit loading station is considered a contamination zone and general health physics rules apply. Care should be

5 11924-2 5-22-80

taken when loading a rabbit so as not to touch any contaminated equipment.

5.0 PROCEDURE

1

5.1 Sample Preparation

- 5.1.1. If a solid specimen is large and crystalline in structure, it should be crushed or ground to a fine consistency and blended well before weighing.
- 5.1.2. Transfer the sample to a previously prepared and tared insert and weigh the sample. Record the weight in the log book for delayed neutron measurements along with the sample identification. Cap the insert and heat seal the cap to the body of the insert. Load the insert into a rabbit for irradiation.
- 5.1.3. If a small insert is used, the finger of a plastic glove or other similar material should be put in the rabbit to keep the insert in place.
- 5.1.4. If the specimen is a liquid, pipet an aliquot (usually 1 ml) directly into an insert and evaporate the liquid under a heat lamp. After the solution has evaporated, allow it to cool and heat seal the insert. Transfer the insert to a rabbit for irradiation.
- 5.1.5. Some organic compounds that will not dry under ordinary heat lamp temperatures must be treated differently. Pipet these kinds of samples into an insert that contains enough accelerator (Fisher #9-906) to absorb the aliquot (usually 100 µl). Reat seal the insert and load it in a rabbit for irradiation.
- 5.2. Standardization of the Delayed Nautron Counting System
- 5.2.1. Take a 400 second background count of the empty system and record it in the log book. Place the ²⁴¹Am-boron neutron source into

242

5 11924-3 5-22-80

the access tube of the system and count it for 400 seconds. These two counts are not used in the actual analysis, but they are a good check as to the status of the system.

- 5.2.2. Set the control panel for a 20 second irradiation with the rabbit return selector switched for return to the delayed neutron counter. Set the decay timer for 12 seconds and the timer for the counter on 40 seconds.
- 5.2.3. Irradiate, count, and record counts for three blanks using this set up.
- 5.2.4. Irradiate, count, and record counts for two ²³⁵U standards and on:N.B.S. coal and fly ash (Note A).
- 5:3: Analysis of the Samples
- 5.3.1. Irradiate, count, and record the counts for each of the samples in the same manner as with the standards and quality control camples. The samples must correspond in configuration with the standards (standard in large insert - sample in large insert, etc.).
- 5.3.2. The uranium contents can be calculated by hand or by using the MONSTR program on the PDP-15 computer. Hand calculations are performed with the following equation:

235U, ppm = A8/CD Natural U, ppm = U/P

where:

for A = sample aliquot counts corrected for background,

NOTE A: The uranium standards are prepared beforehand exactly as the unknowns are prepared. They may be reused if the insert remains in good condition.

5 11924-4

H = weight of 235U in the standard, Mg.

C = aliquot weight, g,

D = standard counts corrected for background,

u = 235u, ppm.

P = abundance of 235U in sample U.

5.3.3. Store the data on DEC tape and save the line printer output for reporting the data.

6.0 PRECISION

1

6.1. The precision for uranium by this method is excellent if the count rate of the sample is significantly above the blank count rate. Precision can be improved by careful preparation of the unknowns and standards and careful attention to maintaining the same irradiation geometry for standards and samples. If the criteria are met, the accuracy is probably +2.52

7.0 INTERFERENCES

A specimen that is either depleted or enriched in ²³⁵U will of course give erroneous natural-U results. Any element present in high concentration with a high capture cross section could cause a negative bias in the results. Under a high gamma field, the BF, neutron detectors lose efficiency, and the resulting count rate may be somewhat low giving erroneous results.

8.0 REFERENCES

 Dyer, F. F., Emery, J. F., and Leddicotte, G. W., A Comprehensive Study of the Neutron Activation Analysis of Uranium by Delayed Neutron Counting. ORNL-3342 (1962).

5 11924-5 5-22-80

 "Reactor Operations and Experiments," ORNL Radiation Safety and Control Manual, Section 8-9.

Written by:

J. W. Wade

Approved by:

- 1. J. R. Stokely, Section Head
- 2. W. D. Shults, Director, ORNL Analytical Chemistry Division

APPENDIX 5.2.4

OSDH SAMPLING PROCEDURE

Joan K Leavier M.D. Commissioner

Bears of Health Edward H. Fiel, J. M.D. Aresiden Wollace I vel, M.D. Vice Aresi Jene Robert D. McCullough & D.O. Scottany

John B. Carmichael D.D.S. James A. Cox, Jr. M.D. Linde M. Johnson, M.D. Ernest D. Marc: Water Scott Misson, III WA, Toal: "Dator OKLAHOMA STATE

RO. BOX 53551 1000 N.E. TENTH OKLAHOMA CITY, OK 73152



AN EQUAL OPPORTUNITY EMPLOYER

January 27, 1986

MEHORANDUM

To: Sequoyah Fuese Incident File

From: Robert L. Craig, Director Radiation Protection Division

Subject: Sample Collection Technique

Samples collected after January 5 and before January 10 taken as follows:

Soil - collected by sticking a shovel into the ground about 2 to 4 inches and placing it in a plastic bag.

Vegetation - a quantity of leaves (either dead or living) were placed in a plastic bag. The vegetation was collected in the general area of the soil sample.

These samples are suitable for tracing the path of the plume but are not suitable for determining depositions of material per unit area.

Soil and vegetation samples collected after January 9 were taken by removing the grass and leaves from an area measuring one foot on each side. The vegetation was placed in a plastic bag. After removal of the vegetable matter, the soil from that square foot area was removed and placed in another plastic bag. An attempt was made to remove only the 1/2 to 1 inch of soil but in some cases as much as 2 inches of soil was taken.

All water samples were collected from standing surface water and placed in half-gallon plastic jugs.

Two samples of each type were taken from each sampling location. One was sent the NRC Laboratory and the other to the OSDH Laboratory for analysis.

APPENDIX 5.2.5

1

.

.

SFC SAMPLING PROCEDURE (2 DOCUMENTS)

*

4

SEQUOYAH FUELS CORPORATION

INTERNAL CORRESPONDENCE

TO W. L. Utnage

DATE January 21, 1986

FROM A. W. Norwood

SUBJECT Sequoyah Recovery Operation

Information concerning the survey and sampling techniques used, and the chronological sequence of actions taken by the Sequoyah Fuels' Cimarron Facility Response team has been requested by E.G.&G. personnel. With your concurrence this correspondence will be provided to them to assist in the documentation of the sample and survey results.

Cimarron survey and sampling team members at various points in time were:

A. W. Norwood - Cimarron Facility Manager
Ron Fine - Health Physics Supervisor
Will Rogers - Decontamination Supervisor
Joe Kegin - Maintenance Supervisor
Claude Thompson - Health Physics Specialist
Jack Andrews - Senior Health Physics Technician
Frank Murch - Senior Health Physics Technician
Don Rall - Decontamination Supervisor

All of the direct alpha survey readings provided on the survey map documents were taken with portable alpha survey instruments with 60 cm² probes. The sensivitiy of the portable alpha survey instruments allows good confidence in readings as low as 200 dpm/ 60 cm²; therefore, negative readings are shown as <200 dpm/60 cm² on the survey maps.

Off-site smears at the residences south of highway I-40 were counted with portable smear counters with a sensitivity that allows good confidence in readings as low as 6 dpm/100 cm². Most of the smears taken on site within the exclusion area fence were counted with portable alpha survey instruments with negative readings shown as <200 dpm/100 cm². All smear readings shown on the survey maps are per 100 cm².

The sampling criteria for the vegetation and soil grab samples taken around the residences by the Cimarron team was only that at least 100 grams of sample be taken to assure an ample quantity for analysis. This same criteria was used for the 180' grid grab samples taken on-site from the lawn between the process building and the south exclusion area fence as shown on the grid map. Five of these grab samples were taken immediately outside of the exclusion area fence as shown on the grid map. The five grab samples taken just outside of the exclusion area fence were split with NRC. Both vegetation and soil grab samples were taken on the 180' grid pattern and soil grab sampling was repeated in the areas where the grass was removed on the same grid pattern. W. L. Utnage January 21, 1986 Page 2

The soil, vegetation, and 12" core sampling of the area outside of the exclusion area south to highway I-40 was started on Wednesday, January 8. These samples consisted of all of the vegetation on one square foot and one square foot of soil (1/4" deep). The soil was taken from the same square foot from which the vegetation sample had been removed. The core samples were taken in the same general area.

Survey maps and sampling grid maps were submitted to document control covering the following chronological sequence of actions taken by the Cimarron response team.

Saturday, January 4, 1986

Supplemented on-site survey instrumentation with equipment from Cimarron, spot surveyed administration area and paint plant entrance area, checked plant air support systems and began a survey of eight residences located south of I-40. Direct. and smear surveys of out buildings, fences, vehicles, etc., were taken using alpha survey meters. Both vegetation and soil samples were also collected and submitted to the research center.

One two man team accompanied Scott Munson to the Sallisaw and Ft. Smith Hospitals to survey, clean and release the locations used to treat personnel. Medical staffs were also surveyed and cleared. Upon completion of the residence survey, the remaining team returned to the plant and began a survey of the administration area. Contamination was found, offices were cleaned, plastic was placed over the carpet to contain the in-place contamination, shoe covers were issued, and two contamination control points were established. One H-P technician remained at the guard station to survey personnel and vehicles exiting and entering the plant. Decontaminated H-P office for air sample service.

Sunday, January 5, 1986

One team of H-P technicians went to all the area hospitals to survey, clean and release. One H-P went to the Carlisle School with Bob Craig, OHD, to survey. Highway 10 was surveyed and closed. The south guard station was also surveyed and closed. The North guard station was established as the facility entry gate and surveyed of personnel and vehicles exiting the facility was initiated. One H-P stayed at the guard station. Additional equipment and supplies were brought from Cimarron. Surveyed residence.

Monday, January 6, 1986

Definition of the radiological dropout from the plume was conducted by surveying south of the facility and water samples were taken in indicated plume area. Cleaning of Highway 10 began. Samples of the defined plume area were January 21, 1986 Page 3

taken on 180' centers for vegetation and soil inside the exclusion area fence. The Tyler home south of I-40 was surveyed and samples were taken and submitted for analysis. A survey of the process area was made by Ron Fine. Don Majors was established as the focal point to document environmental surveys and sample results onto maps.

Tuesday, January 7, 1986

Orientations of contractor personnel for health and safety was given and stripping of the contaminated areas in front of the plant began to avoid further contamination spread. H-P monitors were used to control traffic through the high airborne area and to monitor operations. Definition of the plume area was completed between exclusionary fence and I-40. The ambulances used to transport the injured employees were tracked down, surveyed, and released. Washing of the contaminated grass outside the plant fence was initiated with water/soda ash solution. Cleaning of Highway 10 was completed. It was surveyed at 10° intervals and released.

Wednesday, January 8, 1996

Stripping operations continued by contractor personnel wearing respirators. The drivers were also wearing lapel samples. Began taking samples of the plume area (i.e., one on the center line, one on each boundary, and one at 1000 feet outside the plume on each side). These samples were taken at 1000 foot intervals with baseline samples. One square foot of grass, one square foot of soil (1/4") and a core sample at each location was taken for the NRC. Urine samples of the stripping operators were collected each day.

Thursday, January 9, 1986

Completed grass stripping. Moved the excavated soil to a diked storage area northwest of plant. Completed the soil and vegetation sampling of the plume grid and began taking NRC core samples. Began washdown operations of roadways and vehicles. Pulled "pinch" samples of soil where the grass had been removed for counting in the H-P proportional counter. Resampled the soil after excavation in accordance with the previous grid.

Friday, January 10, 1986

Completed the roadway and sidewalk survey at 25 foot intervals (NRC verified). Surveyed the O.G.&E. substation and recorded no cleaning was done. Received data from proportional counter on pinch sampling taken on Thursday indicating soil remaining W. L. Utnage January 21, 1986 Page 4

in front yard was 8 dpm/2" planchet. Resurveyed office area. Steam cleaned vehicles. Washed down front of building. Continued NRC core sampling.

Saturday, January 11, 1986

Contacted Quadrex for H-P assistance to monitor personnel exit of the restricted area in response to NRC's comments. Continued NRC core sampling. Lined sod storage with PVC liner for contaminated control. Cleaned and released south guard station, Highway 10 entry, and service road and reestablished south guard station operations. Watered Highway 10. Steam cleaned poles and signs at front of plant and released.

Sunday, January 12, 1986

Surveyed UF₄ construction site and released. Surveyed all roads inside of perimeter fence and released. Resurveyed plant administration area. Continued NRC core sampling.

Monday, January 13, 1986

Resodding began after orientation by H-P. Surveyed administration areas for final release of shoe cover requirements. Continued NRC core sampling. Wayne Norwood assumed additional responsibilities for process clean-up.

Tuesday, January 14, 1986

Met with NRC and released administration area of shoe cover requirements. Completed NRC core sampling (for the third time). Surveyed, cleaned and released lunchroom vendor's car. Follow-up to ensure fencing of sod and storage area was completed. Continued resodding. Resampled areas in plume area as follow-up to water samples taken on 1/6/86. Urine samples from the re-sodders will be collected upon completion of the work activity.

a. n. noused

AWN:dd

xc: C. A. Grosclaude Document Control

1/17/86 To Don Major's From J. E Charland Subject - Sompling Accorde Environmental Jomples To kan Scatt of I-400 for evaluation of impost as a result of the accident at sayangh Fules on 1/4/06. All songlos ware sont to the KM-Tech Carter for analysis for U and F. All Songlas were taken doing The wook of 1/4/86. Songle method were as follows. 1. Veg, tation - Where vegstation was present an area 1 42 in area was mortal off. Clipars were used to Cut De vegetion off al transfard to plastic Jacks. These ware then mortal has I don't becation and tragate to the plant. Junjar, Carton, bar and other songles other Then groond Cover were Songled near The lastion 03 a glib some. 2. Soil- The are directly bolow the Vegitation Somple was Collacked to a dept of appox 1/2" inch. These Songles were trons for to uplo ste Bog of B: Morked Low Idontification In one where no grow come was

ausholle a 125 area was collecte to a dest of approx 12". 3. Rediction Measoran ents. - direct Repha readings were made at each sompting site using a alpha Similforion probe al 200 PAC -1_ Scular. & A background measor nut Was made covering The morphe with The plastic Guard. Direct alphe measorenet was othen mode by placing the probe directly over the area to de messarel. Both Back your and direct measure Coont. were note I minute courts were Usal. 4 Water Joydes - Grob Jonglos ware to han at all draininges song I to on The county Road at locations inducto on map (CN-247) No preservities Of L'Itration was done. Jonglos wer tronsported to the Tech couler her andys Ar US F al Gross X. 258

APPENDIX 5.2.6

OSDH LABORATORY ANALYSIS PROCEDURES

January 13, 1986

MEMORANDUM

- TO: Dale McHard, Chief Radiation and Special Hazards
- FROM: Judith A. Duncan, Chief (141) State Environmental Laboratory Service
- SUBJECT: Summary of analytical methods used for analyses of soil, vegetation and water samples from the area of the recent Sequoyah Fuels Incident.

Bob Craig requested that we provide a written summary of our analytical methods and quality control procedures for the recent fluoride analyses which we performed on samples submitted in relation to the Sequoyah Fuels Incident earlier this month. It is my understanding that NRC requested this information.

Attached you will find an outline of our methods. Should you or anyone from NRC have any questions, please contact me.

SEDIMENTS:

5 grams of sediment were weighed into a 100 milliliter beaker. Equal amounts of Deionized water and TISAB (Total lonic Strength Adjustment Buffer) were added. Sample was then refluxed for 1 hour. Sample was brought back up to original volume with D.I. water if evaporation occurred. After sample cooled it was analyzed on an Orion Fluoride Meter. Fluoride value was determined by following calculation:

$$\frac{V_{W} \times C_{f}}{Wt} = C_{s} mg/kg$$

where

V_w = Volume of D.I. water added in ml.

- Cf = Concentration of fluoride read from meter
- wt = Weight of sample used in grams

Cs = Concentration of fluoride in sample

- FOLIAGE: The procedure for foliage differs from that of sediment only in the weight of sample used and the need to remove the foliage from the DI/TISAB solution to enable stir bar to rotate.
- QUALITY CONTROL: Duplicates and spikes were done every 10 samples. % Recovery ranged from 60% to 127%. The data measurements were made in hopes of approximating the passage of the plume by looking for gross differences in the fluoride concentration.

GROSS ALPHA AND GROSS BETA ANALYTICAL PROCEDURES

SAMPLE PREPARATION:

SOIL: The sample was dried on a hot plate and ground. 2.5 grams of sample was placed in a 50 ml capped centrifuge rube, 25 ml of deionized water was added, shaken for about 10 minutes and centrifuged. 10 ml of the supernatant was transferred to a 2 in. x 1/4 in. ringed SS planchette and dried under a heat lamp. The planchette was counted on a low background, thin window proportional counter for 10 minutes.

VEGETATION: 2 grams of the sample (as received) were placed in a 50 ml capped centrifuge tube, 25 ml of deionizied water was added and shaken for about 10 minutes. The liquid was decanted and centrifuged. 10 ml of the liquid was transferred to a 2 in. x 1/4 in. tared planchette and sufficient CaSO, solution was added to ensure that the solids on the planchette would be in excess of 10 milligrams. The planchette was dried under a heat lamp, weighed, and counted for 10 minutes.

CALCULATION: Gross alpha and gross beta concentrations were calculated as follows:

Activity (pCi/gm) = Net X 2.5 WT X EFF X 2.22

Where:

Net =	Net cpm of alpha or beta activity
WT =	Sample weight in grans
EFF =	Counting efficiency (cpm/dpm)
2.5 =	Ratio of volume of water added to sample
1.1.1	volume of sample analyzed.
2.22=	dpm/pCi

QUALITY CONTROL:

Counter background and counting efficiency were checked prior to counting each batch of samples. APPENDIX 5.2.7

)

Ą

0

-

ENVIRONMENTAL SAMPLING: VEGETATION AND SOIL

Preface

Shortly after the UF_6 release on January 4, 1986, an extensive program of environmental measurements and sampling was begun, with the objective of obtaining measurements sufficient to establish the nature and extent of contamination from the release. The effort involved staff from SFC, from the Oklahoma State Department of Health, from Headquarters and Region IV of NRC, and, under the auspices of the Department of Energy, from EG&G, Inc.

Acknowledgements

For the production of this data base through many long days of environmental sampling, an exemplary achievement under the circumstances, we must extend our thanks to the following:

- To Robert Craig, Dale McHard and other staff of the Oklahoma State Department of Health and local government agencies for their cooperation and generous efforts in extending and completing the sample collection of soil, vegetation and water,
- To John Stauter, Jim Cleveland, Don Majors, WAyne Norwood and the other staff of the Sequoyah Facility and Kerr-McGee for their cooperation in the extensive sampling effort,
- To John Doyle, Zolin Burson, Allen Fritzsche, Keith Roesner, Ronald Reiman, Marc Rivera, Hollis Berry and the 15 other people of EG&G, Inc., who performed the aerial photography, the Aerial Measurements System survey, the in situ surface level radioactive contaminated survey, the computerized sampling data management system, and the magnificent photographic record of the proceedings,
- To Marvin Dockter of the Department of Energy for his on-the-scene coordination of the EG&G efforts and to Joseph Deal and Fritz Wolff of DOE Headquarters for their support to this effort,
- To Jerry Everett, Blair Spitzberg, Harry Pettengill, and Candice Jierree of NRC's Region IV for their aid and guidance in establishing a standardized sampling system and an organized data base,
- To George Greenly, Kevin Foster, John Nostrom and Thomas Sullivan of LLNL-ARAC for their assistance on atmospheric dispersion,
- To Joseph Stewart, Susan Allay, and their staff at ORNL for their skillful and cooperative efforts on sample analyses,
- To Barry Zalcman of NRC Headquarters who provided essential on-the-scene coordination.

Sample Collection

Environmental sampling began on January 4, 1986, soon after the release. Initially, "grab" samples of soil, vegetation, and water were taken. These served to indicate roughly the concentrations and extent of contamination. However, comparisons of results from different locations and sampling organizations were difficult due to use of nonuniform sampling procedures. Therefore, procedures for more uniform sampling of the ground surface were designed, and put into practice by the SFC sampling teams on January 8, 1986, and by the OSDH sampling teams on January 10, 1986. Thus samples collected on those dates or later are of a more uniform quality.

Samples were collected by two organizations, the Oklahoma State Department of Health and by Sequoyah Fuels Corporation. Most samples were split, with the second portion being sent to Oak Ridge National Laboratory for confirmatory analysis. Thus analyses were performed by SFC or the Kerr-McGee Technical Center (designated K-M), by OSDH and by ORNL.

Maps of Sampling Locations and Results

Sampling locations were plotted at the SFC directly on overlays on two base maps, one being a 1" = 400' map of the site and its immediate surroundings, and the other a composite of 1:24,000 scale U.S. Geological Survey maps covering a larger area. These two base map areas have been designated the near-field and the far-field areas. They are also characterized by having different sampling location coordinate systems, the near field using the SFC plant grid and the far field using the land survey township-range-section grid found on the USGS maps. (Futher details on the coordinate systems are in a later section.)

The sampling locations and results on (in most cases) the original overlays have been superimposed on base maps and are included here (reduced in size) as Figures 5.2.7.1A through 5.2.7.24A. The significance of use of the originals is that they reflect the location plotting done at the time, at the SFC, and thus they include the contribution to accuracy made by input of the individuals who personally collected the samples. Figures 5.2.7.1A through 5.2.7.4A show the sampling locations for vegetation, soil, and water. Figures 5.2.7.5A through 5.2.7.24A show the analysis results for vegetation and soil (results for water sampling are in Appendix 5.3.1). Results for grab samples are italicized. Some grab sample locations were resampled under the systematic sampling program and thus one location may have two results, for example, fluoride on vegetation. For reproduction here, the overlays have been superimposed on basemaps redrawn to eliminate topography and other extraneous detail which tends to obscure the results.

Tables of Sampling Locations and Results

The plots of analysis results are followed by tables of the results. The tables provide the sample label, the coordinates (plant grid or township-range-section) of the sample location, a designation of the organization which performed the sample analysis, the date the sample was collected, the type of analysis (i.e., whether for uranium or for fluoride), the numerical result and its units, and comments which call out unusual characteristics of the sample. Some samples were ana'yzed for alpha and beta radioactivity (rather than for uranium per se) and are listed accordingly.

At the time the samples were collected, some, taken near residences, were labeled with the residents' names. To allow the residents some privacy, these labels have been changed.

Tables 5.2.7.1A through 5.2.7.17A are listings by type of analysis (e.g., U or F) and by the analysis-performing organization. Tables 5.2.7.18A through 5.2.7.25A are the same information sorted numerically and alphabetically by the sample label (to simplify looking up results by a location label); these tables also include water sample results as found in Appendix 5.3.1.

Coordinate Systems

The near-field figures, e.g., Figure 5.2.7.1A, show the plant grid; the grid is in feet, with the lines shown being 1,000 feet apart. The plant buildings are in Section 21 of T12N, R21E; Highway 10 runs north-south on the east boundary of Section 21. The southeast corner of Section 21 (which is in Highway 10) is the point on the plant grid which has coordinates N-10,000, E-10,000. The N-10,000 line of the plant grid lies along the section line (in Highway 10). The availability of maps with the plant grid superimposed made the plant grid a natural choice for the locations of near-field sampling.

Similarly, the availability of U.S. Geological Survey maps showing the sections and section lines made the township-range-section grid a natural choice for the locations of far-field sampling. In Oklahoma, as in much of the United States, the land was surveyed in sections which are a square, one mile on each side. Where the topography permits, many roads and fences run on the section line, or at right angles to it at 1/2 or 1/4 of the section. Thus it is relatively simple to determine one's location on a map which also shows the section lines. The far-field f gures, e.g., Figure 5.2.7.3A, show the section lines, the township and range lines, and the section numbers of the sections adjoining the intersections of the township and range lines. Figure 5.2.7.25A shows the section numbering in each 6-mile by 6-mile township, and the relationship of the numbers of the sections with the numbers of adjacent sections. For purposes of designating a sampling location more precisely than the nearest mile (section), an ad-hoc division of sections into 1/14-mile subdivisions was made, with 14 strips designated A through N from north to south, and 14 strips numbered 1 through 14 from west to east. This arrangement is shown in Figure 5.2.7.26A. This arrangement permits designating a sample location to the nearest 1/14 mile, roughly 100 meters.

For its atmospheric dispersion simulation calculations, LLNL-ARAC chose to use the Universal Transverse Mercator (UTM) coordinate system, which is also indicated on the U.S. Geological Survey maps, in kilometers. The advantage of the UTM coordinate system is that it is shown on the maps for long distances without significant corrections. (The plant grid is only shown for a few miles. The township-range-section grid contains offsets at the T8N, T9N line and at the T12N, T13N line, and is significantly offset at the Oklahoma-Arkansas border.) A difficulty with the UTM coordinate system is that it is set about 1.2 degrees of angle counterclockwise from the section lines at the plant site, and the angular difference changes slightly over long distances. For reference, the southeast corner of Section 21, T12N, R21E (N-10,000; E-10,000 on the plant grid) is approximately at UTM 3929.49N, 311.41E (coordinates in km). Figures 4.1.2.7, 4.1.3.1, 4.1.2.28A and 4.1.2.29A, show both the UTM coordinate grid and the township-range-section grid.




Figure 5.2.7.1A Near field locations sampled by OSDH




Figure 5.2.7.2A Near field locations sampled by SFC





Figure 5.2.7.3A Far field locations sampled by OSDH

Figure 5.2.7.4A Far field locations sampled by SFC



273




Figure 5.2.7.5A ORNL results, near field, fluoride on vegetation




Figure 5.2.7.6A OSDH results, near field, fluoride on vegetation





Figure 5.2.7.7A K-M results, near field, fluoride on vegetation






Figure 5.2.7.9A OSDH results, near field, fluoride in soil




Figure 5.2.7.10A K-M results, near field, fluoride in soil

Ģ




Figure 5.2.7.11A ORNL results, near field, uranium on vegetation




Figure 5.2.7.12A K-M results, near field, uranium on vegetation




Figure 5.2.7.13A ORNL results, near field, uranium in soil





Figure 5.2.7.14A K-M results, near field, uranium in soil





Figure 5.2.7.15A ORNL results, far field, fluoride on vegetation



.

Figure 5.2.7.16A OSDH results, far field, fluoride on vegetation









Figure 5.2.7.18A ORNL results, far field, fluoride in soil








Figure 5.2.7.20A K-M results, far field, fluoride in soil




Figure 5.2.7.21A ORNL results, far field, uranium on vegetation




Figure 5.2.7.22A K-M results, far field, uranium on vegetation




Figure 5.2.7.23A ORNL results, far field, uranium on soil




Figure 5.2.7.24A K-M results, far ield, uraniam in soil

	R20E	R21E					R21E	R22E
T12N	36	31	32	33	34	35	36	31
TIIN	1	6	5	4	3	2	1	6
	12	7	8	9	10	11	12	7
	13	18	17	16	15	14	13	18
	24	19	20	21	22	23	24	19
	25	30	29	28	27	26	25	30
TIIN	36	31	32	33	34	35	36	31
TION	1	6	5	4	3	2	1	6

Township and Range

Figure 5.2.7.25A Township-Range-Section grid




Figure 5.2.7.26A Grid for subsection designations

RAW SOIL AND VEGETATION DATA

NEAR FIELD ENVIRONMENTAL SAMPLING RESULTS BY MEDIA AND TYPE ANALYSIS

Table 5.2.7.1A

LABEL	NORTH	EAST I	LAB	DATE	TYPE	RESULT	UNITS	COMMENTS
105	11700	8450	KMTC	1-09	F	20.0	ppm	
108	10600	8750		1-09	F	30.0	ppm	
111	14250	8450 1		1-09	F	40.0	ppm	
116	10650	12300		1-09	F	10.0	ppm	
118	9750	12450		1-09	F	14.0	ppm	
120	9700	8850		1-09	F	20.0	ppm	
120	9600	9900 1		1-09	F	210.0	ppm	
122	9700	10700		1-09	F	100.0	ppm	
124	9750	11850		1-09	F	20.0	ppm	
128	10650	11300		1-09	F	35.0	ppm	
130	10650	10350		1-09	F	70.0	ppm	
130	11600	10650		1-09	F	80.0	ppm	
132	9250	13150		1-09	F	28.0	ppm	
134	9150	12050		1-09	F	16.0	ppm	
	9200		KMTC	1-09	F	90.0	ppm	
138	9150	9000		1-09	F	15.0	ppm	
142	10500	9700		1-09	F	50.0	ppm	
144	11600	9550		1-09	F	73.0	ppm	
146			KMTC	1-09	F	640.0	ppm	
150	11450	10100		1-05	F	1500.0	ppm	CEDAR
50	10750		KMTC	1-06	F	2100.0	ppm	CLOAR
A-01	13130			1-06	F	2600.0	ppm	
A-07	13130		KMTC	1-06	F	3000.0		
A-13	13130		KMTC		F	230.0	ppm	
A-19	13130	9760		1-06		35.0	ppm	
A-25	13130	9950		1-06	F		ppm	1/4 N.GATE
BG 1/4	14750	9500		1-06	F	40.0	ppm	CRES
CRES08	6150		KMTC	1-03	F	42.0	ppm	DRY LEAVES
DRES12	6100		KMTC	1-04	F	160.0	ppm	UNI LLAVES
E-01	13000		KMTC	1-06	F	700.0	ppm	
E-07	13000		KMTC	1-06	F	3100.0	ppm	
E-13	13000		KMTC	1-06	F	3700.0	ppm	
E-19	13000	9770		1-06	F	760.0	ppm	
E-25	13000	9950		1-06	F	80.0	ppm	
HRES02	6950	10650		1-04	F	24.0	ppm	
HRES03	6950	10650		1-04	F	80.0	ppm	
J-07	12840	9410		1-06	F	870.0	ppm	
J-13	12840	9600		1-06	F	1700.0	ppm	
J-19	12840	9780		1-06	F	1300.0	ppm	
J-25	12840	9950		1-06	F	140.0	ppm	TREE DATA DRCT
LRES04	6950	9900		1-04	F	3.0	ppm	TREE DATA DRUT
P-01	12670	9230		1-06	F	270.0	ppm	
P-07	12670		KMTC	1-06	F	2600.0	ppm	
P-12	12670		KMTC	1-06	F	450.0	ppm	AVEDAGE VALUE
P-13	12670	9600	KMTC	1-06	F	2900.0	ppm	AVERAGE VALUE
P-25	12670	9950	KMTC	1-06	F	320.0	ppm	
R-15	6050	13750	KMTC	1-04	F	32.0	ppm	ENTER
SM-21	7450	17250	KMTC	1-04	F	5.0	ppm	LEAVES
A-25	13130	9950	ORNL	1-06	F	33.0	ug/g	
BG 1/4	14750	9500	ORNL	1-06	F	7.3	ug/g	

Table 5.2.7.1A (Continued)

LABEL	NORTH	EAST LAB	DATE	TYPE	RESULT	UNITS	COMMENTS
CARLILE	14250	15100 ORNL	1-05	F	8.3	ug/g	
CRES	6150	10800 ORNL		F		ug/g	
E-25		9950 ORNL	1-06	F	16.2 58.0	ug/g	
J-25		9950 ORNL	1-06	F	150.0	ug/g	
P-19	12670	9780 ORNL	1-06 1-06	F	2900.0		
P-25	12670	9950 ORNL	1-06	F			
T1-01	12500	7200 ORNL	1-05	F	280.0		
T1-01	12500	7200 ORNL	1-05	F	1.7	ug/g	
T1-02	12500	8000 ORNL	1-05	r F	1.7	ug/g	
T1-03	12500	8550 ORNL	1-05	F	1.2	ug/g	
T1-04	12500	8900 ORNL	1-05	F	1.7 1.2 4.7 19.2	ug/g	
T1-05	12500	0300 ORNL	1-05	F	19.2	ug/g	
T1-06	12500	9200 ORNL 9500 ORNL	1-05	F	1800.0		
T1-07	12500	9500 URNL	1-05	F	930.0	~ ~	
	12500	9850 ORNL	1-05		750.0		
T1-08		10200 ORNL	1-05	F	44.6	ug/g	
T1-09		10600 ORNL 11300 ORNL	1-05	F	17.5	ug/g	
T1-11	Contraction of the contraction of the	11300 ORNL	1-05	F	8.8	ug/g	
T1-12		11900 ORNL	1-05	F	10.2	ug/g	
T2-01		6200 ORNL	1-05 1-05 1-05 1-05	F	4.7	ug/g	
T2-03	10650	7550 ORNL	1-05	F	5.2	ug/g	
T2-05	10800	9300 ORNL	1-05	F	250.0	ug/g	
T2-06		9950 ORNL	1-05	F	290.0	ug/g	
T2-08		11650 ORNL	1-06	F	2.3	ug/g	
T2-09	11350	12700 ORNL	1-06	F	2.3	ug/g	
T2-10	11400	13700 ORNL	1-06		2.7	ug/g	
T3-01	8300	6250 ORNL	1-06	F	4.7	ug/g	
T3-02	8600	6950 ORNL	1-06	F	3.0	ug/g	
T3-03	8650	9000 ORNL	1-06		7.2	ug/g	
T3-04	8350	10300 ORNL	1-06		420.0	ug/g	
T3-05	8350	11200 ORNL	1-06		67.0	ug/g	
T3-06	8400	12050 ORNL	1-06	F	7.0	ug/g	
T3-07	8550	13050 ORNL	1-06	F	3.6		
CARLILE	14250	15100 OSDH	1-05	F	4.3	mg/kg	GRASS
CARLILE	14100	15300 OSDH	1-14		2.56	mg/kg	URA33
COMP-1		9500 OSDH			488.0		LAB COMPARISON
COMP-2	8800	10100 OSDH	1-16	F	164.0	mg/kg	
COMP-3	10050	10700 OSDH	1-16	F	128.0	mg/kg	LAB COMPARISON
COMP-4	14450	14800 OSDH	1-16	F		mg/kg	LAB COMPARISON
COMP-5	12700	11200 OSDH	1-16	F	7.36	mg/kg	LAB COMPARISON
CRES	6150	10750 OSDH		F	1080	mg/kg	LAB COMPARISON
HRES	6950	10650 OSDH	1-05		32.0	mg/kg	GRASS
LR-01	7300		1-08	F	24.0	mg/kg	
LR-01		10450 OSDH	1-11	F	56.8	mg/kg	
	7300	10450 OSDH	1-11	F	56.8	mg/kg	
LR-02	6000	10000 OSDH	1-11	F	2.96	mg/kg	
LR-02	6000	10000 OSDH	1-11	F	2.96	mg/kg	
LR-03	6000	10800 OSDH	1-11	F	56.00	mg/kg	
LR-03	6000	10800 OSDH	1-11	F	56.00	mg/kg	
LR-04	5800	12100 OSDH	1-11	F	5.60	mg/kg	
LR-04	5800	12100 OSDH	1-11	F	5.60	mg/kg	
LR-05	5975	13950 OSDH	1-11	F	7.12	mg/kg	

Table 5.2.7.1A (Continued)

LABEL	NORTH	EAST	LAB	DATE	TYPE	RESULT	UNITS	COMMENTS
LR-05	5975	13950	OSDH	1-11	F	7.12	mg/kg	
LR-05A	6100	14450		1-14	F	13.20	mg/kg	
LR-05A	6100	14450		1-14	F	13.20	mg/kg	
LR-06	6000	15600		1-14	F	4.48	mg/kg	
R-01	11750	10700		1-11	F	33.60	mg/kg	
R-02	10800	10700		1-11	F	84.00	mg/kg	
R-03	10750	11550		1-11	F	6.80	mg/kg	
RD-01	8150	10550		1-16	F	40.00	mg/kg	
RD-02	7400	10750		1-16	F	46.40	mg/kg	
RD-03	6700	10750		1-16	F	53.60	mg/kg	
RD-04	6300	10750		1-16	F	64.00	mg/kg	
T1-01	12500		OSDH	1-05	F	1.2	mg/kg	EVERGREEN
T1-01	12500		OSDH	1-14	F	3.68	mg/kg	
T1-02	12500	8000	OSDH	1-05	F	1.8	mg/kg	GRASS
T1-03	12500	8550	OSDH	1-05	F	2.6	mg/kg	GRASS
T1-04	12500	8900	OSDH	1-05	F	10.4	mg/kg	EVERGREEN
T1-04	12500	8900	OSDH	1-14	F	4.96	mg/kg	
T1-05	12500	9200	OSDH	1-05	F	408.0	mg/kg	EVERGREEN
T1-05	12500	9200	OSDH	1-14	F	28.8	mg/kg	
T1-06	12500	9500	OSDH	1-05	F	500.0	mg/kg	GRASS
T1-06	12500	9500	OSDH	1-14	F	1320.0	mg/kg	
T1-07	12500	9850	OSDH	1-05	F	216.0	mg/kg	GRASS
T1-07	12500	9850	OSDH	1-14	F	240.0	mg/kg	
T1-08	12600	10200	OSDH	1-05	F	55.0	mg/kg	GRASS
T1-08	12600	10200	OSDH	1-14	F	20.0	mg/kg	
T1-09	12600	10600	OSDH	1-05	F	20.0	mg/kg	GRASS
T1-09	12600	10600	OSDH	1-14	F	17.6	mg/kg	
T1-11	12600	11300	OSDH	1-05	F	11.3	mg/kg	GRASS
T1-12	12550	11900	OSDH	1-05	F	9.9	mg/kg	EVERGREEN
T1-12	12550	11900	OSDH	1-14	F	7.36	mg/kg	
T1-13	12650	13200	OSDH	1-14	F	2.56	mg/kg	
T1-14	13200	14100	OSDH	1-14	F	5.92	mg/kg	
T2-01	10750		OSDH	1-05	F	4.0	mg/kg	
T2-02	10550		OSDH	1-14	F	2.16	mg/kg	
T2-03	10650	7550	OSDH	1-05	F	1.6	mg/kg	
T2-03	10650	7550	OSDH	1-14	F	3.12	mg/kg	
T2-04	10700		OSDH	1-14	F	4.72	mg/kg	
T2-05	10800		OSDH	1-14	F	21.6	mg/kg	
T2-05	10800	the case of the second	OSDH	1-05	F	91.3	mg/kg	
T2-06	11500		OSDH	1-14	F	584.0	mg/kg	
T2-06	11500		OSDH	1-06	F	122.2	mg/kg	
T2-07	11500	10650		1-14	F	22.4	mg/kg	
T2-08	11500	11650		1-06	F	3.1	mg/kg	
T2-08	11500	11650		1-14	F	18.8	mg/kg	
T2-09	11350	12700		1-06	F	<0.8	mg/kg	
T2-10	11400	13700		1-06	F	1.6	mg/kg	
T2-10	11400	13700	OSDH	1-14	F	5.52	mg/kg	

Table 5.2.7.1A (Continued)

LABEL	NORTH	EAST	LAB	DATE	TYPE	RESULT	UNITS	COMMENTS	
T3-01	8300		OSDH	1-06	F	4	mg/kg	EVERGREEN	
T3-01	8300	6250	OSDH	1-14	F	2.48	mg/k		
T3-02	8600		OSDH	1-06	F	2.4	mg/ky	EVERGREEN	
T3-02A	8250		OSDH	1-16	F	1.76	mg/kg		
T3-03	8650		OSDH	1-06	F	8.8	mg/kg	EVERGREEN	
T3-03	8650	900r	USDH	1-14	F	9.28	mg/kg		
T3-03A	8500	9550	OSDH	1-16	F	5.92	mg/kg		
T3-03B	8300	9950	OSDH	1-16	F	41.60	mg/kg		
T3-03C	8800	10050	OSDH	1-16	F	13.68	mg/kg		
T3-04	8350	10300	OSDH	1-06	F	102.0	mg/kg	EVRGRN/LEAVES	
T3-04	8350	10300	OSDH	1-14	F	4.00	mg/kg		
T3-04	8350	10700	OSDH	1-16	F	73.60	mg/kg		
T3-04B	8350	11000	OSDH	1-16	F	72.00	mg/kg		
T3-05	8350	11200	OSDH	1-06	F	93.0	mg/kg	EVERGREEN	
T3-05	8350	11200		1-14	F	30.4	mg/kg		
T3-06	8400	12050	OSDH	1-06	F	5.8	mg/kg	EVERGREEN	
T3-06	8400	12050	OSDH	1-14	F	5.84	mg/kg		
T3-07	8550	13050	OSDH	1-06	F	4.2	mg/kg	GRASS	
T3-07	8550	13050	OSDH	1-14	F	8.48	mg/kg		
T5-01	17700	8700	OSDH	1-15	F	5.92	mg/kg		
T5-01	17700		OSDH	1-08	F	2.48	mg/kg		
T5-02	16500	10000	OSDH	1-08	F	4.56	mg/kg		
T5-02	16500	10000	OSDH	1-08	F	1.76	mg/kg		
T5-03	16300	10350	OSDH	1-08	F	3.6	mg/kg		
T5-04	15500	8500	OSDH	1-08	F	2.56	mg/kg		
T6-01	13950	7700	OSDH	1-08	F	28.8	mg/kg		
T6-02	13900	8400	OSDH	1-08	F	48.0	mg/kg		
T6-03	14650	8300	OSDH	1-08	F	10.48	mg/kg		
T6-05	14450	9550	OSDH	1-08	F	15.84	mg/kg		

Table 5.2.7.2A

LABEL	NORTH	EAST	LAB	DATE	TYPE	RESULT	UNITS	COMMENTS
103	13800	10200	КМТС	1-09	F	120.0	ppm	
106	11700		KMTC	1-09	F	150.0	ppm	
109	10600	8750		1-09	F	160.0	ppm	
112	14250	8450		1-09	F	110.0	ppm	
115	11500	11600		1-09	F	100.0	ppm	
117	10650	12300		1-09	F	80.0	ppm	
119	9750	12450		1-09	F	110.0	ppm	
121	9700		KMTC	1-09	F	70.0	ppm	
123	9600		KMTC	1-09	F	130.0	ppm	
125	9700	10700		1-09	F	60.0	ppm	
127	9750	11850		1-09	F	90.0	ppm	
129	10650	11300		1-09	F	80.0	ppm	
131	10650	10350		1-09	F	70.0	ppm	
133	11600	10650		1-09	F	110.0	ppm	
135	9250	13150		1-09	F	90.0	ppm	
137	9150	12050		1-09	F	150.0	ppm	
139	9200	11000		1-09	F	120.0	ppm	
141	9300	10100		1-09	F	160.0	ppm	
143	9150		KMTC	1-09	F	80.0	ppm	
145	10500		KMTC	1-09	F	100.0	ppm	
147	11600	9550	KMTC	1-09	F	30.0	ppm	
149	11650	9900	KMTC	1-09	F	160.0	ppm	
187	12600	7200	KMTC	1-14	F	60.0	ppm	
190	14750	9500	KMTC	1-14	F	70.0	ppm	
193	12550	10200		1-14	F	80.0	ppm	
196	12600	10750		1-14	F	70.0	ppm	
199	12550	11850		1-14	F	80.0	ppm	
202	13150	10200		1-14	F	100.0	ppm	
205	13150	11250	KMTC	1-14	F	150.0	ppm	
209	12500	9900	KMTC	1-14	F	80.0	ppm	
218	12500	9350	KMTC	1-14	F	120.0	ppm	
221	12500	9100	KMTC	1-14	F	90.0	ppm	
A-01	13130	9220	KMTC	1-06	F	250.0	ppm	
A-07	13130	9400	KMTC	1-06	F	410.0	ppm	
A-07	13130	9400	KMTC	1-09	F	410.0	ppm	
A-13	13130	9580	KMTC	1-06	F	320.0	ppm	
A-13	13130	9580	KMTC	1-09	F	140.0	ppm	
A-25	13130	9950	KMTC	1-06	F	110.0	ppm	
BG 1/4	14750		KMTC	1-06	F	220.0	ppm	
CRES10	6150	10750	KMTC	1-04	F	190.0	ppm	CRES
DRES11	6100	11100	KMTC	1-04	F	40.0	ppm	
DRES13	6100	11100	KMTC	1-04	F	60.0	ppm	
E-01	13000	9230	KMTC	1-06	F	100.0	ppm	
E-01	13000	9230	KMTC	1-09	F	180.0	ppm	
E-07	13000	9410	KMTC	1-06	F	170.0	ppm	
E-07	13000	9410	KMTC	1-09	F	170.0	ppm	
E-13	13000	9600	KMTC	1-06	F	230.0	ppm	
E-13	13000	9600	KMTC	1-09	F	150.0	ppm	
E-19	13000	9770	KMTC	1-06	F	40.0	ppm	
E-25	13000		KMTC	1-06	F	170.0	ppm	
HRES01	6950	10650	KMIC	1-04	F	110.0	ppm	

LABEL	NORTH	EAST	LAB	DATE	TYPE	RESULT	UNITS	COMMENTS
J-07	12840	9410	KMTC	1-06	F	150.0	ppm	
J-07	12840		KMTC	1-09	F	180.0	ppm	
J-13	12840		KMTC	1-06	F	220.0	ppm	
J-13	12840		KMTC	1-09	F	200.0	ppm	
J-19	12840		KMTC	1-06	F	90.0	ppm	
J-25	12840		KMTC	1-06	F	110.0	ppm	
LRES05	6950		KMTC	1-04	F	120.0	ppm	
LRES06	6950		KMTC	1-04	F	60.0	ppm	
P-07	12670		KMTC	1-06	F	300.0	ppm	
P-07	12670		KMTC	1-09	F	170.0	ppm	
P-13	12670		KMTC	1-06	F	50.0	ppm	
P-13	12670		KMTC	1-09	F	290.0	ppm	
P-19	12670		KMTC	1-06	F	210.0	ppm	
P-25	12670		KMTC	1-06	F	90.0	ppm	
R-16	6050	13750		1-04	F	80.0	ppm	
SM-20	7450	17250		1-04	F	40.0	ppm	
A-25	13130		ORNL	1-06	F	9.0	ug/g	WET
A-25	13130		ORNL	1-06	F	11.4	ug/g	DRY
BG 1/4	14750		ORNL	1-06	F	3.7	ug/g	WET
BG 1/4	14750		ORNL	1-06	F	4.7	ug/g	DRY
CARLILE	14250	15100		1-05	F	1.2		DRY
CARLILE	14250	15100		1-05	F	0.8	ug/g	WET
CRES	6150	10750		1-05	F	1.9	ug/g	
CRES	6150	10750		1-05	F		ug/g	WET
E-25	13000		ORNL	1-06	F	2.8	ug/g	DRY
E-25	13000		ORNL	1-06	F	7.4	ug/g	WET
J-25	12840		ORNL	1-06	F	8.6	ug/g	DRY
J-25	12840		ORNL	1-06	F	8.3	ug/g	WET
P-19	12670	9600		1-06	F	11.2	ug/g	DRY
P-19	12670	9600		1-06	F	28.7	ug/g	WET
P-25	12670	9780		1-06	F	32.8	ug/g	DRY
P-25	12670	9780		1-06	F	14.2	ug/g	WET
T1-01	12500	7200		1-05	F	21.0	ug/g	DRY
T1-01	12500	7200		1-05	F F	3.1	ug/g	DRY
T1-02	12500	8000			F F	2.3	ug/g	WET
T1-02	12500	8000		1-05	F	1.4	ug/g	DRY
T1-03	12500	8550		1-05	F	1.2	ug/g	WET
T1-03	12500	8550		1-05	F	4.2	ug/g	DRY
11-04	12500	8900		1-05	F	3.4	ug/g	WET
T1-04	12500	8900		1-05		7.5	ug/g	DRY
T1-05		9200			F	5.7	ug/g	WET
T1-05	12500	9200		1-05	F	23.5	ug/g	DRY
T1-06	12500	9200		1-05	F	17.7	ug/g	WET
T1-06	12500			1-05	F	22.5	ug/g	WET
T1-07	12500	9500 9850		1-05	F	26.3	ug/g	DRY
T1-07				1-05	F	12.2	ug/g	WET
T1-07	12500	9850		1-05	F	22.8	ug/g	DRY
T1-08	12600	10200		1-05	F	6.6	ug/g	DRY
T1-08	12600	10200		1-05	F.	5.1	ug/g	WET
T1-09	12600	10600		1-05	F	2.1	ug/g	DRY
	12600	10600		1-05	F	1.7	ug/g	WET
T1-11	12600	11300	UKNL	1-05	F	7.6	ug/g	DRY

LABEL	NORTH	EAST	LAB	DATE	TYPE	RESULT	UNITS	COMM	IENTS
T1-11	12600	11300	ORNL	1-05	۲	4.6	ug/g	WET	
T1-12	12550	11900		1-05	F	2.4	ug/g	DRY	
T1-12	12550	11900		1-05	F	1.4	ug/g	WET	
T2-01	10750		ORNL	1-05	F	1.1	ug/g	DRY	
T2-01	10750		ORNL	1-05	F	0.9	ug/g	WET	
T2-02	10550		ORNL	1-05	F	2.9	ug/g	DRY	
T2-02	10550		ORNL	1-05	F	2.5	ug/g	WET	
T2-02	10700		ORNL	1-05	F	15.1	ug/g	DRY	
T2-04	10700		ORNL	1-05	F	13.9	ug/g	WET	
				1-05	F	1.4		DRY	
T2-05	10800		ORNL		F		ug/g	WET	
T2-05	10800		ORNL	1-05		1.2	ug/g	DRY	
T2-06	11450		ORNL	1-05	F	21.4	ug/g		
T2-06	11450		ORNL	1-05	F	20.2	ug/g	WET	
T2-07	11500	10650		1-06	F	14.9	ug/g	DRY	
T2-07	11500	10650		1-06	F	14.3	ug/g	WET	
T2-08	11500	11650		1-06	F	13.9	ug/g	DRY	
T2-08	11500	11650		1-06	F	10.2	ug/g	WET	
T2-09	11350	12700		1-06	F	0.6	ug/g	DRY	
T2-09	11350	12700		1-06	F	0.5	ug/g	WET	
T2-10	11400	13700	ORNL	1-06	F	4.1	ug/g	WET	
T2-10	11400	13700	ORNL	1-06	F	6.5	ug/g	DRY	
T3-01	8300	6250	ORNL	1-06	F	1.1	ug/g	DRY	
T3-01	8300	6250	ORNL	1-06	F	0.9	ug/g	WET	
T3-02	8600		ORNL	1-06	F	3.4	ug/g	DRY	
T3-02	8600		ORNL	1-06	F	2.9	ug/g	WET	
T3-03	8650		ORNL	1-06	F	3.2	ug/g	DRY	
T3-03	8650		ORNL	1-06	F	2.8	ug/g	WET	
T3-04	8350	10300		1-06	F	1.4	ug/g	DRY	
T3-04	8350	10300		1-06	F	1.3	ug/g	WET	
T3-05	8350	11200		1-06	F	2.2	ug/g	DRY	
T3-05	8350	11200		1-06	F	1.6	ug/g	WET	
	8400	12050		1-06	F	2.1	ug/g	DRY	
T3-06				1-06	F	1.7	ug/g	WET	
T3-06	8400	12050						DRY	
T3-07	8550	13050		1-06	F	2.0	ug/g		
T3-07		13050		1-06		1.4	ug/g	WET	
CARLILE	14100	15300		1-14	F	1.24	mg/kg		
CARLILE	14250	15100		1-05	F	<0.4	mg/kg		
COMP-1	11800	9500		1-16	F	5.44	mg/kg		COMPARISON
COMP-2	8800	10100		1-16	F	2.76	mg/kg		COMPARISON
COMP-3	10050	10700		1-16	F	1.96	mg/kg		COMPARISON
COMP-4	14450	14800		1-16	F	1.48	mg/kg		COMPARISON
COMP-5	12700	11200		1-16	F	5.6	mg/kg	LAB	COMPARISON
CRES	6150	10750	OSDH	1-05	F	0.6	mg/kg		
LR-01	7300	10450	OSDH	1-11	F	0.88	mg/kg		
LR-01	7300	10450		1-11	F	0.88	mg/kg		
LR-02	6000	10000		1-11	F	1.36	mg/kg		
LR-02	6000	10000		1-11	F	1.36	mg/kg		
LR-03	6000	10800		1-11	F	1.68	mg/kg		
LR-03	6000	10800		1-11	F	1.68	mg/kg		
LR-04	5800	12100		1-11	F	1.12	mg/kg		
LR-04		12100		1-11	F	1.12	mg/kg		
LR-04	5800	15100	0304	1-11	F	1.12	ing/ kg		

LABEL	NORTH	EAST LAB	DAIL	TYPE	RESULT	UNITS	COMMENTS
LR-05	5975	13950 OSDH	1-11	F	1.92	mg/kg	
LR-05	5975	13950 OSDH	1-11	F	1.92	mg/kg	
LR-05A	6100	14450 OSDH	1-14	F	1.36	mg/kg	
LR-06	6000	15600 OSDH	1-14	F	1.12	mg/kg	
R-01	11750	10700 OSDH	1-11	F	3.32	mg/kg	
R-02	10800	10700 OSDH	1-11	F	2.64	mg/kg	
R-03	10750	11550 OSDH	1-11	F	4.00	mg/kg	
RD-01	8150	10550 OSDH	1-16	F	5.20	mg/kg	
RD-02	7400	10750 OSDH	1-16	F	2.96	mg/kg	
RD-03	6700	10750 OSDH	1-16	F	0.48	mg/kg	
RD-04	6300	10750 OSDH	1-16	F	8.80	mg/kg	
T1-01	12500	7200 OSDH	1-05	F	0.8	mg/kg	
T1-02	12500	8000 OSDH	1-05	F	<0.4	mg/kg	
T1-02	12500	8000 OSDH	1-14	F	0.92	mg/kg	
T1-03	12500	8550 OSDH	1-05	F	1.3	mg/kg	
T1-04	12500	8900 OSDH	1-05	F	0.9	mg/kg	
T1-04	12500	8900 OSDH	1-14	F	1.84	mg/kg	
T1-05	12500	9200 OSDH	1-05	F	4.8	mg/kg	
T1-05	12500	9200 OSDH	1-14	F	4.08	mg/kg	
T1-06	12500	9500 OSDH	1-05	F	3.6	mg/kg	
T1-06	12500	9500 OSDH	1-14	F	3.64		
T1-07	12500	9850 OSDH	1-05	F	2.0	mg/kg	
T1-07	12500	9850 OSDH	1-14	F	3.00	mg/kg	
T1-08	12600	10200 OSDH	1-05	F	0.8	mg/kg	
T1-08	12600	10200 OSDH	1-14	F		mg/kg	
T1-09	12600	10600 OSDH	1-05	F	2.20	mg/kg	
T1-09	12600	10600 OSDH	1-14	F	0.6	mg/kg	
T1-11	12600	11300 OSDH	1-05	F	1.96	mg/kg	
T1-12	12550	11900 OSDH	1-05	F	1.0	mg/kg	
T1-12	12550	11900 OSDH	1-14	F	0.8	mg/kg	
T1-13	12650	13200 OSDH	1-14	F	3.76	mg/kg	
T1-14	13200	14100 OSDH	1-14	F	0.48	mg/kg	
T2-01	10750	6200 OSDH	1-05	F	0.96	mg/kg	
T2-02	10550	6800 OSDH	1-14	F	<0.4	mg/kg	
T2-02	10550	6800 OSDH	1-05	F		mg/kg	
T2-03	10650	7550 OSDH	1-14	F	0.48	mg/kg	
T2-03	10650	7550 OSDH	1-05	F	1.2	mg/kg	
T2-04	10700	8400 OSDH	1-14	F	1.6	mg/kg	
T2-04	10700	8400 OSDH	1-05	F	0.76	mg/kg	
T2-05	10800	9300 OSDH	1-14	F	1.96	mg/kg	
T2-05	10800	9300 OSDH	1-05	F	0.8	mg/kg	
T2-06	11500	9850 OSDH	1-14	F	0.4	mg/kg	
T2-06	11500	9850 OSDH		F	1.4	mg/kg	
T2-07	11500	10650 OSDH	1-06 1-14	F	1.12	mg/kg	
T2-07	11500	10650 OSDH	1-14	F	0.92	mg/kg	
T2-08	11500	11650 OSDH	1-14	F	6.72	mg/kg	
T2-08	11500	11650 OSDH		F	11.2	mg/kg	
T2-09	11350	12700 OSDH	1-06		3.12	mg/kg	
T2-10	11400	13700 OSDH	1-06 1-14	F	<0.8	mg/kg	
12 10	11400	13/00 0304	1-14	r	1.32	mg/kg	

LABEL	NORTH	EAST L	AB	DATE	TYPE	RESULT	UNITS	COMMENTS
T2-10	11400	13700 0	SDH	1-06	F	1.6	mg/kg	
T3-01	8300	6250 0	SDH	1-06	F	0.8	mg/kg	
T3-01	8300	6250 0	SDH	1-14	F	1.32	mg/kg	
T3-02	8600	6950 0	SDH	1-06	F	1.2	mg/kg	
T3-02A	8250	8750 0	SDH	1-16	F	3.00	mg/kg	
T3-03	8650	9000 0	SDH	1-06	F	2.5	mg/kg	
T3-03	8650	9000 0	SDH	1-14	F	2.36	mg/kg	
T3-03A	8500	9550 0	SDH	1-16	F	21.60	mg/kg	
T3-03B	8300	9950 0	SDH	1-16	F	5.56	mg/kg	
T3-03C	8800	10050 0	SDH	1-16	F	<0.40	mg/kg	
T3-04	8350	10300 0	SDH	1-06	F	<0.8	mg/kg	
T3-04	8350	10300 0	SDH	1-14	F	0.52	mg/kg	
T3-04A	8350	10700 0	SDH	1-16	F	8.80	mg/kg	
T3-04B	8350	11000 0	SDH	1-16	F	1.68	mg/kg	
T3-05	8350	11200 0	SDH	1-06	F	1.2	mg/kg	
T3-05	8350	11200 0	SDH	1-14	F	0.84	mg/kg	
T3-06	8400	12050 0	SDH	1-06	F	1.4	mg/kg	
T3-06	8400	12050 0	SDH	1-14	F	0.76	mg/kg	
T3-07	8550	13050 0	SDH	1-06	F	<0.8	mg/kg	
T3-07	8550	13050 0	SDH	1-14	F	0.64	mg/kg	
T5-01	17700	8700 0	SDH	1-15	F	1.84	mg/kg	
T5-01	17700	8700 0	SDH	1-08	F	0.68	mg/kg	
T5-03	16300	10350 0	SDH	1-08	F	<0.4	mg/kg	
T6-01	13950	7700 0	SDH	1-08	F	3.56	mg/kg	
T6-02	13900	8400 0	SDH	1-08	F	9.0	mg/kg	
T6-03	14650	8300 0	SDH	1-08	F	2.28	mg/kg	
T6-04	14650	8700 03	SDH	1-08	F	1.8	mg/kg	
							5 5	

Table 5.2.7.3A

NEAR FIELD ENVIRONMENTAL SAMPLING RESULTS URANIUM ON VEGETATION SAMPLES

LABEL	NORTH	EAST	LAB	DATE	TYPE	RESULT	UNITS	COMMENTS
105	11700	8450	KMTC	1-09	U	5.3	ppm	
108	10600	8750	KMTC	1-09	U	5.2	ppm	
111	14250	8450	KMTC	1-09	Ŭ	120.0	ppm	
116	10650	12300	KMTC	1-09	Ŭ	9.2	ppm	
118	9750	12450	KMTC	1-09	U	2.7		
							ppm	
120	9700	8850	KMTC	1-09	U	7.7	ppm	
122	9600	9900	KMTC	1-09	U	480.0	ppm	
124	9700	10700	KMTC	1-09	U	110.0	ppm	
126	9750	11850	KMTC	1-09	U	4.6	ppm	
128	10650	11300	KMTC	1-09	U	22.0	ppm	
130	1(650	10350	KMTC	1-09	U	71.0	ppm	
132	11600	10650	KMTC	1-09	U	25.0	ppm	
134	9250	13150	KMTC	1-09	U	5.4	ppm	
136	9150	12050	KMTC	1-09	U	7.2	ppm	
138	9200	11000	KMTC	1-09	U	69.0	ppm	
142	9150	9000	KMTC	1-09	U	3.2	ppm	
144	10500	9700	KMTC	1-09	U	92.0	ppm	
146	11600	9550	KMTC	1-09	U	130.0	ppm	
150	11450	9700	KMTC	1-09	U	1400.0	ppm	
198	12550	11850	KMTC	1-14	U	2.5	ppm	
201	13150	10200		1-14	U	9.1	ppm	
204	13150	11250		1-14	U	1.6	ppm	
208	12500		KMTC	1-14	Ŭ	540.0	ppm	
211	12500		KMTC	1-14	Ŭ	360.0	ppm	
214	12500		KMTC	1-14	U	1530.0		
50	10750	10100		1-06	Ŭ	17000.0	ppm	
A-01	13130		KMTC	1-06		4200.0	ppm	
A-07	13130		KMTC	1-06	U	8000.0	ppm	
A-13	13130		KMTC	1-06	U	6200.0	ppm	
A-19	13130		KMTC	1-06	U	7150.0	ppm	
A-25	13130		KMTC	1-06	U	49.0	ppm	
BG 1/4	14750		KMTC	1-06	U	64.0	ppm	
CRES08	6150	10750		1-04	U	12.0	ppm	
DRES12	6100	11100		1-04	U	220.0	ppm	
E-01	13000	9230		1-06	U	1340.0		
E-07	13000	9410	KMTC	1-06	U	14200.0	ppm	
E-13	13000	9600	KMTC	1-06	U	7420.0	ppm	
E-19	13000	9770	KMTC	1-06	U	21600.0	ppm	
E-25	13000	9950	KMTC	1-06	U	40.0	ppm	
HRES02	6950	10650		1-04	U	32.0	ppm	
HRES03	6950	10650		1-04	U	0.95	ppm	
J-07	12840		KMTC	1-06	U	1520.0	ppm	
J-13	12840		KMTC	1-06	Ŭ	6640.0	ppm	
J-19	12840		KMTC	1-06	Ŭ	3220.0	ppm	
J-25	12840		KMTC	1-06	Ŭ	210.0		
LRES04	6950		KMTC	1-08	U	<0.4	ppm	
							ppm	
P-01	12670		KMTC	1-06	U	5300.0	ppm	
P-07	12670		KMTC	1-06	U	4480.0	ppm	
P-12	12670		KMTC	1-06	U	5940.0	ppm	
P-13	12670	9600	KMTC	1-06	U	3800.0	ppm	

NEAR FIELD ENVIRONMENTAL SAMPLING RESULTS URANIUM ON VEGETATION SAMPLES

LABEL	NORTH	EAST	LAB	DATE	TYPE	RESULT	UNITS	COMMENTS
P-25	12670	9950	KMTC	1-06	U	360.0	ppm	
R-15	6050	13750	KMTC	1-04	U	0.6	opm	
SM-21	7450	17250	KMTC	1-04	U	0.6	ppm	
A-25	13130	9950	ORNL	1-06	U	11.8	ppm	
BG 1/4	14750	9500	ORNL	1-06	U	3.91	ppm	
CARLILE	14250	15100	ORNL	1-05	U	0.188	ppm	
CRES	6150	10800	ORNL	1-05	U	13.8	ppm	
E-25	13000	9950	ORNL	1-06	U	24.2	ppm	
J-25	12840	9950	ORNL	1-06	U	148.0	ppm	
P-19	12670	9780	ORNL	1-06	U	3600.0	ppm	
P-25	12670	9950	ORNL	1-06	U	186.0	ppm	
T2-05	10800	9300	ORNL	1-05	U	393.0	ppm	
T2-06	11450	9950	ORNL	1-05	U	285.0	ppm	
T2-08	11500	11650	ORNL	1-06	U	0.633	ppm	
T2-09	11350	12700	ORNL	1-06	U	1.53	ppm	
T2-10	11400	13700	ORNL	1-06	U	0.596	ppm	
T3-01	8300	6250	ORNL	1-06	U	5.03	ppm	
T3-02	8600	6950	ORNL	1-06	U	3.21	ppm	
T3-03	8650	9000	ORNL	1-06	U	3.39	ppm	
T3-04	8350	10300	ORNL	1-06	U	624.0	ppm	
T3-05	8350	11200	ORNL	1-06	U	177.0	ppm	
T3-06	8400	12050	ORNL	1-06	U	6.0	ppm	
T3-07	8550	13050	ORNL	1-06	U	35.6	ppm	

NEAR FIELD ENVIRONMENTAL SAMPLING RESULTS URANIUM IN SOIL SAMPLES

LABEL	NORTH	EAST LA	AB	DATE	TYPE	RESULT	UNITS	COMMENTS
103	13800	10200 KM	ITC	1-09	U	20.0	ppm	
106	11700	100 C 100 C 100 C 100 C	MTC	1-09	U	<10	ppm	
109	10600	and the second second	MTC	1-09	U	<10	ppm	
112	14250		MTC	1-09	U	290.0	ppm	
115	11500		ATC	1-09	U	<10	ppm	
	10650	12300 K		1-09	Ŭ	20.0	ppm	
117	9750	12450 K		1-09	Ŭ	20.0	ppm	
119	9700	8850 K		1-09	U	20.0	ppm	
121	9600	9900 KI		1-09	Ŭ	25.0	ppm	
123	9700	10700 K		1-09	Ŭ	20.0	ppm	
125		11850 K		1-09	U	<10	ppm	
127	9750	11300 K		1-09	Ŭ	<10	ppm	
129	10650		MTC	1-09	U	30.0	ppm	
131	10650			1-09	U	15.0	ppm	
133	11600	10650 K		1-09	U	<10	ppm	
135	9250	13150 K			U	<10	ppm	
137	9150	12050 K		1-09	U	<10		
139	9200	11000 K		1-09		15.0	ppm	
141	9300	10100 K		1-09	U		ppm	
143	9150	9000 K		1-09	U	<10	ppm	
145	10500	9700 K		1-09	U	30.0	ppm	
147	11600	9550 K		1-09	U	35.0	ppm	
149	11650	9900 K		1-09	U	95.0	ppm	
187	12600	7200 K		1-14	U	20.0	ppm	
190	14750	9500 K		1-14	U	50.0	ppm	
193	12550	10200 K		1-14	U	30.0	ppm	
196	12600	10750 K		1-14	U	20.0	ppm	
199	12550		(MTC	1-14	U	10.0	ppm	
202	13150		MTC	1-14	U	30.0	ppm	
205	13150		(MTC	1-14	U	<10	ppm	
209	12500		KMTC	1-14	U	140.0	ppm	
218	12500		KMTC	1-14	U	220.0	ppm	
221	12500	9100 H	KMTC	1-14	U	40.0	ppm	
A-01	13130	9220	KMTC	1-06	U	2180.0	ppm	
A-07	13130	9400	KMTC	1-09	U	590.0	ppm	
A-07	13130	9400	KMTC	1-06	U	730.0	ppm	
A-13	13130	9580	KMTC	1-06	U	655.0	ppm	
A-13	13130	9580 1	KMTC	1-09	U	200.0	ppm	
A-25	13130	9950	KMTC	1-06	U	25.0	ppm	
BG 1/4	14750	9500 1	KMTC	1-06	U	<10	ppm	
CRES10	6150	10750	KMTC	1-04	U	<10	ppm	
DRES11	6100	11100		1-04	U	15.0	ppm	
DRES13	6100	11100		1-04	U	13.0	ppm	
E-01	1:000	9230		1-06	U	210.0	ppm	
E-01	13000	9230		1-09	U	300.0	ppm	
E-07	13000	9410		1-06	U	465.0	ppm	
E-07	13000	9410		1-09	U	380.0	ppm	
E-13	13000	9600		1-06	U	130.0	ppm	
E-13	13000	9600		1-09	U	130.0	ppm	
E-19	13000	9770		1-06	U	15.0	ppm	
E-25	13000	9950		1-06	U	20.0	ppm	
HRESOI	6950	10650		1-04	U	<10	ppm	
11112-001	0000	20000			-			

LABEL	NORTH	EAST	LAB	DATE	TYPE	RESULT	UNITS	COMMENTS
J-07	12840	9410	КМТС	1-06	U	125.0	ppm	
J-07	12840		KMTC	1-09	U	250.0	ppm	
J-13	12840		KMTC	1-06	U	200.0	ppm	
J-13	12840		KMTC	1-09	U	125.0	ppm	
J-19	12840		KMTC	1-06	U	165.0	ppm	
J-25	12840		KMTC	1-06	U	20.0	ppm	
LRES05	6950		KMTC	1-04	Ŭ	<10	ppm	
LRES06	6950		KMTC	1-04	U	<10	ppm	
P-07	12670		KMTC	1-06	ŭ	235.0	ppm	
P-07	12670		KMTC	1-09	U	90.0	ppm	
P-13	12670		KMTC	1-05	U	55.0	ppm	
P-13	12670		KMTC	1-09	U	175.0	ppm	
P-19	12670	9780		1-06	U	90.0	ppm	
P-25	12670	9950		1-06	Ŭ	40.0	ppm	
R-16	6050	13750		1-04	Ŭ	<10	ppm	
SM-20	7450	17250		1-04	U	11.0	ppm	
A-25	13130	9760		1-06	Ŭ	10.5	ppm	WET
A-25	13130	9760		1-06	Ŭ	13.2	ppm	DRY
BG 1/4	14750	9500		1-06	Ŭ	6.1	ppm	WET
BG 1/4	14750	9500		1-06	Ŭ	7.88	ppm	DRY
CARLILE	14250	15100		1-05	Ŭ	2.7	ppm	WET
CARLILE	14250	15100		1-05	Ŭ	3.7	ppm	DRY
CRES	6150	10750		1-05	Ŭ	2.9	ppm	DRY
CRES	6150	10750		1-05	Ŭ	2.0	ppm	WET
E-25	13000	9950		1-06	Ŭ	8.9	ppm	WET
E-25	13000	9950		1-06	Ŭ	10.3	ppm	DRY
J-25	12840	9950		1-06	Ŭ	9.8	ppm	WET
J-25	12840	9950		1-06	Ŭ	13.1	ppm	DRY
P-19	12670	9600		1-06	Ŭ	135.0	ppm	WET
P-19	12670	9600		1-06	U	154.0	ppm	DRY
P-25	12670	9780		1-06	Ŭ	19.2	ppm	WET
P-25	12670	9780		1-06	Ŭ	28.4	ppm	DRY
T1-01	12500	7200		1-05	Ŭ	4.6	ppm	WET
T1-01	12500	7200		1-05	Ŭ	6.2	ppm	DRY
T1-02	12500	8000		1-05	Ŭ	5.7	ppm	DRY
T1-02	12500	8000		1-05	Ŭ	4.5	ppm	WET
T1-03	12500	8550		1-05	Ŭ	5.7	ppm	WET
T1-03	12500	8550		1-05	Ŭ	7.0	ppm	DRY
T1-04	12500	8900		1-05	Ŭ	13.0	ppm	WET
T1-04	12500	8900		1-05		17.3	ppm	DRY
T1-05	12500	9200		1-05	U	34.0	ppm	WET
T1-05	12500	9200		1-05	U	45.1	ppm	DRY
T1-06	12500	9500		1-05	U	90.7	ppm	DRY
T1-06	12500	9500		1-05	Ŭ	77.6	ppm	WET
T1-07	12500	9850		1-05	U	33.7	ppm	WET
T1-07	12500	9850		1-05	U	62.7	ppm	DRY
T1-08	12500	10200		1-05	U	7.3		DRY
T1-08	12600	10200		1-05	U	5.7	ppm	WET
T1-08	12600	10200		1-05	U	4.5	ppm	DRY
T1-09					U		ppm	WET
	12600	10600		1-05		3.6	ppm	DRY
T1-11	12600	11300	UKNL	1-05	U	6.0	ppm	DRT

NEAR FIELD ENVIRONMENTAL SAMPLING RESULTS URANIUM IN SOIL SAMPLES

LABEL	NORTH	EAST LAB	DATE	TYPE	RESULT	UNITS	COMMENTS
T1-11	12600	11300 ORNL	1-05	U	3.7	ppm	WET
T1-12	12550	11900 ORNL	1-05	U	6.0	ppm	DRY
T1-12	12550	11900 ORNL	1-05	U	3.5	ppm	WET
T2-01	10750	6200 ORNL	1-05	U	2.9	ppm	WET
T2-01	10750	6200 ORNL	1-05	U	3.5	ppm	DRY
T2-02	10550	6800 ORNL	1-05	U	2.6	ppm	DRY
T2-02	10550	6800 ORNL	1-05	U	2.2	ppm	WET
T2-03	10650	7550 ORNL	1-05	U	3.2	ppm	WET
T2-03	10650	7550 ORNL	1-05	U	3.9	ppm	DRY
T2-04	10700	8400 ORNL	1-05	U	4.1	ppm	DRY
T2-04	10700	8400 ORNL	1-05	U	3.2	ppm	WET
T2-05	10800	9300 ORNL	1-05	U	4.5	ppm	DRY
T2-05	10800	9300 ORNL	1-05	U	3.9	ppm	WET
T2-06	11450	9950 ORNL	1-05	U	25.9	ppm	DRY
T2-06	11450	9950 ORNL	1-05	U	24.4	ppm	WET
T2-07	11500	10650 ORNL	1-06	U	5.23	ppm	DRY
T2-07	11500	10650 ORNL	1-06	U	5.0	ppm	WET
T2-08	11500	11650 ORNL	1-06	U	5.21	ppm	DRY
T2-08	11500	11650 ORNL	1-06	U	3.8	ppm	WET
T2-09	11350	12700 ORNL	1-06	U	3.95	ppm	DRY
T2-09	11350	12700 ORNL	1-06	U	3.0	ppm	WET
T2-10	11400	13700 ORNL	1-06	U	2.4	ppm	WET
T2-10	11400	13700 ORNL	1-06	U	3.86	ppsi	DRY
T3-01	8300	6250 ORNL	1-06	U	4.34	ppm	DRY
T3-01	8300	6250 ORNL	1-06	U	3.3	ppm	WET
T3-02	8600	6950 ORNL	1-06	U	3.86	ppm	DRY
T3-02	8600	6950 ORNL	1-06	U	3.3	ppm	WET
T3-03	8650	9000 ORNL	1-06	U	5.34	ppm	DRY
T3-03	8650	9000 ORNL	1-06	U	4.7	ppm	WET
T3-04	8350	10300 ORNL	1-06	U	4.28	ppm	DRY
T3-04	8350	10300 ORNL	1-06	U	3.9	ppm	WET
T3-05	8350	11200 ORNL	1-06	U	4.0	ppm	DRY
T3-05	8350	11200 ORNL	1-06	U	3.0	ppm	WET
T3-06	8400	12050 ORNL	1-06	U	4.88	ppm	DRY
T3-06	8400	12050 ORNL	1-06	U	3.9	ppm	WET
T3-07	8550	13050 ORNL	1-06	U	3.9	ppm	DRY
T3-07	8550	13050 ORNL	1-06	U	2.8	ppm	WET

310

NEAR FIELD ENVIRONMENTAL SAMPLING RESULTS SOIL AND VEGETATION SAMPLES

LABEL	NORTH	EAST LAB	DATE	TYPE	RESULT	UNITS	COMMENTS
CRES09 CRE309	6150 6150	10750 KMT 10750 KMT		F U	90.0 <10	ppm ppm	MIX SOIL/VEG

NEAR FIELD ENVIRONMENTAL SAMPLING RESULTS ALPHA ON VEGETATION SAMPLES

-

LABEL	NORTH	EAST	LAB	DATE	TYPE	RESULT	UNITS	COMMENTS
CARLILE	14250	15100	OSDH	1-05	ALPHA	0.8	pCi/g	GRASS
CRES	6150	10750		1-05	ALPHA	1.1	pCi/g	GRASS
HRES	6950	10650		1-08	ALPHA	5.0	pCi/g	
T1-01	12500	7200		1-05	ALPHA	1.4	pCi/g	EVERGREEN
T1-01 T1-02	12500	8000		1-05	ALPHA	0.1	pCi/g	GRASS
	12500		OSDH	1-05	ALPHA	0.5	pCi/g	GRASS
T1-03	12500	8900	OSDH	1-05	ALPHA	1.3	pCi/g	EVERGREEN
T1-04	12500	9200	OSDH		ALPHA	4020.0	pCi/g	EVERGREEN
T1-05 T1-06	12500	9500	OSDH	1-05	ALPHA	1220.0	pCi/g	GRASS
	12500	9950	OSDH	1-05	ALPHA	124.0	pCi/g	GRASS
T1-07 T1-08	12600	10200	OSDH	1-05	ALPHA	1.0	pCi/g	GRASS
T1-08	12600		0.0011	3 05	ALPHA	0.0	pCi/g	
	12600	11300	OSDH	1-05	ALPHA		pCi/q	GRASS
T1-11 T1-12	12550	11900	OSDH	1-05	ALPHA		pCi/g	
T2-01	10750	6200	OSDH	1-05	ALPHA	-0.2	pCi/g	
T2-01	10650	7550	OSDH	1-05	ALPHA	-0.4	pCi/g	
T2-05	10800	9300	OSDH	1-05	ALPHA	78.0	pCi/g	
	11500	9850	OSDH	1-06	ALPHA	41.0	pCi/g	
T2-06		11650	OSDH	1-06	ALPHA	-0.5	pCi/g	
T2-08 T2-09			OSDH	1-06	ALPHA		pCi/g	
T2-09	11400		OSDH	1-06	ALPHA		pCi/g	
T3-01	5000		OSDH	1-06	ALPHA	1.6	pCi/g	EVERGREEN
T3-01	8600	6950	OSDH	1-06	ALPHA	1.0	pCi/g	
T3-02	8650	6950 9000	OSDH	1-06		0.6	pCi/g	EVERGREEM
T3-03	8350	10300	OSDH	1-06		1070.0		
T3-04	8350		OSDH	1-06		27.0	pCi/g	
T3-05	8400		OSDH	1-06		0.7	pCi/g	
T3-07	8550		OSDH	1-06	ALPHA	0.6	pCi/g	
T5-01	17700	8700		1-08	ALPHA	0.2	pCi/g	
T5-02	16500		OSDH	1-08		-0.2	pCi/g	
T5-02	16500		OSDH	1-08	ALPHA	0.0	pCi/g	
T5-02	16300		OSDH	1-08	ALPHA	0.0	pCi/g	
	15500	8500		1-08	ALPHA	0.2	pCi/g	
T5-04 T6-01	13950	7700	OSDH	1-08	ALPHA	1.2	pCi/g	
	13950	8400	O OSDH	1-08	ALPHA	2.5	pCi/g	
T6-02 T6-03	14650		0 OSDH	1-08	ALPHA	1.0	pCi/g	
	14650		0 OSDH	1-08	ALPHA	0.9	pCi/g	
T6-05	14450	333	0 03011	1 00	The THIT			

.

00

20

NEAR FIELD ENVIRONMENTAL SAMPLING RESULTS BETA ON VEGETATION SAMPLES

LABEL	NORTH	EAST	LAB	DATE	TYPE	RESULT	UNITS	COMMENTS
CARLILE	14250	15100	OSDH	1-05	BETA	0.7	pCi/g	GRASS
CRES	6150	10750	OSDH	1-05	BETA	0.2	pCi/g	GRASS
HRES	6950	10650	OSDH	1-08	BETA	1.9	pCi/g	
T1-01	12500	7200	OSDH	1-05	BETA	1.0		EVERGREEN
T1-02	12500	8000	OSDH	1-05	BETA	1.2	pCi/g	GRASS
T1-03	12500	8550 8900 9200	OSDH	1-05	BETA	1.1	pCi/g	GRASS
T1-04	12500	8900	OSDH	1-05	BETA	0.1	pCi/g	EVERGREEN
T1-05	12500	9200	OSDH	1-05 1-05	BETA	110.0	pCi/g	EVERGREEN
T1-06	12500	9500	OSDH	1-05	BETA	95.0	pCi/q	GRASS
T1-07	12500	9850	OSDH	1-05	BE ⁻ A	23.0	pCi/g	GRASS
T1-08	12600	10200	OSDH	1-05	BETA	5.0	pCi/g	GRASS
T1-09			OSDH				pCi/g	CRASS
T1-11			OSDH	1-05			pCi/g	GRASS
T1-12	12550	11900	OSDH	1-05	BETA		pCi/g	EVERGREEN
T2-01	10750	6200	OSDH	1-05	BETA	0.7	pCi/g	
T2-03	10650	7550	OSDH	1-05	BETA	3.1	pCi/g	
12-05	10800	6200 7550 9300 9850	OSDH	1-05 1-05 1-05	BETA	7.0	pCi/g	
T2-06	11500	9850	OSDH	1-06	BETA	6.0 6.0	pCi/g	
T2-08	11500	11650	OSDH	1-06	BETA	6.0	pCi/g	
T2-09	11350	12700	OSDH	1-06	BETA	2.3	pCi/g	
T2-10		13700						
T3-01		6250		1-06				
T3-02	8600		OSDH	1-06	BETA	0.4	pCi/g	
T3-03	8650			1-06 1-06 1-06	BETA	0.0	pCi/g	
T3-04		10300		1-06	BETA	77.0	pCi/g	EVRGRN/LEAVES
T3-05		11200		1-06	BETA	8 0	pCi/g	EVERGREEN
T3-06	8400		OSDH	1-06	BETA	-0.1	pCi/g	
T3-07	8550			1-06	BETA	-0.1	pCi/g	GRASS
T5-01		8700		1-08	BETA	-0.0	pCi/g	
T5-02		10000		1-08			pCi/g	
T5-02		10000		1-08			pCi/g	
T5-03	16300	10350	OSDH	1-08	BETA		pCi/g	
T5-04	15500	8500	OSDH	1-08	BETA	-0.1	pCi/g	
T5-04 T6-01 T6-02 T6-03	13950	7700	OSDH	1-08	BETA	0.7 1.1 0.3	pCi/g	
T6-02	13900	8400	OSDH	1-08	BETA	1.1	pCi/g	
T6-03	14650	8300	OSDH	1-08	BETA	0.3	pCi/g	
T6-05	14450	9550	OSDH	1-08 1-08 1-08 1-08	BETA	-0.5	pCi/g	

Table 5.2.7.8A

LABEL	NORTH	EAST	LAB	DATE	TYPE	RESULT	UNITS	COMMENTS
CARLILE	14250	15100	OSDH	1-05	ALPHA	0.1	pC1/g	
CRES	6150	10750	OSDH	1-05	ALPHA		pCi/g	
T1-01	12500	7200	OSDH	1-05	ALPHA	0.1	pCi/g	
T1-02		8000				0.1	pCi/g	
T1-03	12500	8550	OSDH	1-05	ALPHA	0.1 0.4	pCi/g	
T1-04		8900				0.9	pCi/g	
T1-05			OSDH					
T1-06	12500		OSDH					
T1-07		9850		1-05	ALPHA			
T1-08	12600	10200		1-05	ALPHA	-0.3		
T1-09	12600	10600	OSDH	1-05	ALPHA	1.6		
T1-11	12600	11300	OSDH	1-05			pCi/g	
T1-12	12550	11900	OSDH	1-05	ALPHA	0.4	pCi/g	
T2-01	10750	6200	OSDH	1-05	ALPHA	0.1	pCi/g	
T2-02	10550	6800	OSDH	1-05	ALPHA	0.1	pCi/g	
T2-03	10650	7550	OSDH	1-05		132.0	pCi/g	
T2-04		8400					pCi/g	
T2-05		9300				0.1	pCi/g	
T2-06		9850					pCi/g	
T2-07	11500		OSDH			3.0	pCi/g	
T2-08	11500		OSDH				pCi/g	
T2-09	11350		OSDH					
T2-10	11400	13700	OSDH	1-06	ALPHA	-0.3	pCi/g	
T3-01	8300	6250	OSDH	1-06	ALPHA		pCi/g	
T3-02	8600	6250 6950 9000	OSDH	1-06 1-06	ALPHA		pCi/g	
T3-03	8650	9000	OSDH	1-06	ALPHA ALPHA	-0.3	pCi/g	
T3-04	26 4 75 1 1	10.400	USUM	= 110	ALPHA	-0.3	pCi/g	
T3-05		11200	OSDH	1-06	ALPHA	0.9	pCi/g	
T3-06	8400	12050	OSDH	1-06	ALPHA	0.4	pCi/g	
T3-07	8550	13050	OSDH	1-06	ALPHA	0.1	pCi/g	
T5-01		8700		1-08	ALPHA	-0.3	pC1/g	
T5-03	16300		OSDH			-0.3	pCi/g	
T6-01		7700				0.1	pCi/g	
T6-02		8400				-0.3	pCi/g	
T6-03		8300					pCi/g	
T6-04	14650	8700	OSDH	1-08	ALPHA	0.4	pCi/g	

Table 5.2.7.9A

LABEL	NORTH	EAST LAB	DATE	TYPE	RESULT	UNITS	COMMENTS
CARLILE	14250	15100 OSDH	1-05	BETA	-0.3	pCi/g	
CRES	6150	10750 OSDH		BETA	-0.0	pCi/g	
T1-01	12500	7200 OSDH	1-05	BETA		pCi/g	
T1-02	12500	8000 OSDH	1-05	BETA		pCi/g	
T1-03		8550 OSDH	1-05			pCi/g	
T1-04		8900 OSDH				pCi/g	
T1-05	12500	9200 OSDH			0.2	pCi/g	
T1-06	12500	9500 OSDH		BETA	-0.2	pCi/g	
T1-07		9850 OSDH		BETA	-0.0	pCi/g	
T1-08	12600	10200 OSDH	1-05	BETA	-0.2	pCi/g	
T1-09	12600	10600 OSDH	1-05	BETA BETA BETA BETA	-0.5	pCi/g	
T1-11	12600	11300 OSDH	1-05	BETA	-0.4	pCi/g	
T1-12	12550	11900 OSDH	1-05	BETA	0.2	pCi/g	
T2-01	10750	6200 OSDH	1-05	BETA	-0.5	pCi/g	
T2-02	16550	6800 OSDH	1-05	BETA	0.1	pCi/g	
T2-03		7550 OSDH				pCi/g	
T2-04		8400 OSDH				pCi/g	
T2-05		9300 OSDH				pCi/g	
T2-06		9850 OSDH				pCi/g	
T2-07		10650 OSDH		BETA		pCi/g	
T2-08	11500	11650 OSDH	1-06	BETA		pCi/g	
T2-09	11350	12700 OSDH	1-06	BETA	0.3	pCi/g	
T2-10	11400	13700 OSDH	1-06	BETA	0.6	pCi/g	
T3-01	8300	6250 OSDH 6950 OSDH 9000 OSDH 10300 OSDH 11200 OSDH	1-06 1-06 1-06 1-06 1-06	BETA	-0.2	pCi/g	
T3-02	8600	6950 OSDH	1-06	BETA	1 0	pCi/g	
T3-03	8650	9000 OSDH	1-06	BETA	1.0	pCi/g	
T3-04	8350	10300 OSDH	1-06	BETA	0.1	pCi/g	
T3-05	8350	11200 OSDH	1-06	BETA	0.4	pCi/g	
T3-06	8400	12050 OSDH	1-06	BEIA	0.2	pCi/g	
T3-07	8550	13050 OSDH	1-06	BETA	-0.1	pCi/g	
T5-01		8700 OSDH	1-08	BETA	0.6	pCi/g	
T5-03	15300	10350 OSDH	1-08	BETA	-0.2	pCi/g	
T6-01		7700 OSDH			0.7	pCi/g	
T6-02		8400 OSDH				pCi/g	
T6-03		8300 OSDH				pCi/g	
T6-04	14650	8700 OSDH	1-08	BETA	-0.6	pCi/g	

FAR FIELD ENVIRONMENTAL SAMPLING RESULTS BY MEDIA AND TYPE ANALYSIS

7C 6 T11N 7 C6 1-08 KMTC 10.0 PPM J.GRASS 7F 7 T11N 7 F7 108 KMTC 20.0 PPM WHEAT 11F 7 T11N 11 F7 1-08 KMTC 10.0 PPM GAT 11F 7 T11N 11 F7 1-08 KMTC 10.0 PPM GAT LEAVES 12N 8 T11N 12 N8 1-08 KMTC 2.0 PPM CAT LEAVES 12N 8 T11N 12 N8 1-08 KMTC 2.0 PPM CAT TIPS 13B1 T11N 13 G11 109 KMTC 10.0 PPM CAT TIPS 13B1 T11N 14 B4 1-08 KMTC 5.0 PPM FESCUE TIPS 14B12 T11N 14 B12 1-08 KMTC 5.0 PPM FESCUE <td< th=""><th>LABEL</th><th>TOWNSHIP</th><th>SECTION</th><th>SUBSECT</th><th>DATE</th><th>LAB</th><th>RESULTS</th><th>UNITS</th><th>COMMENTS</th></td<>	LABEL	TOWNSHIP	SECTION	SUBSECT	DATE	LAB	RESULTS	UNITS	COMMENTS
7F 7 T11N 7 F7 1-08 KMTC 20.0 PPM WHEAT 7F 7 T11N 7 F7 1-08 KMTC 30.0 PPM CAT LEAVES 11F 7 T11N 11 F7 1-08 KMTC 43.0 PPM CAT LEAVES 11F 7 T11N 11 F7 1-08 KMTC 43.0 PPM WILLOW 12N 8 T11N 12 N8 1-08 KMTC 7.0 PPM WILLOW 12N 8 T11N 12 N8 1-08 KMTC 2.0 PPM CEDAR 13B1 T11N 13 B1 1-08 KMTC 10.0 PPM CEDAR 14B4 T11N 14 B4 1-08 KMTC 10.0 PPM FESCUE TIPS 14B12 T11N 14 B12 1-08 KMTC 0.0 PPM FESCUE TIPS 14B12 T11N 16 L12 1-10 KMTC 0.0 PPM FESCUE TIPS <tr< td=""><td>70 6</td><td>T11N</td><td>7</td><td>C6</td><td>1-08</td><td>KMTC</td><td>10.0</td><td>PPM</td><td>J. GRASS</td></tr<>	70 6	T11N	7	C6	1-08	KMTC	10.0	PPM	J. GRASS
7F 7 T11N 7 F7 1-08 KMTC 30.0 PPM 11F 7 T11N 11 F7 1-08 KMTC 10.0 PPM BEANS 11F 7 T11N 11 F7 1-08 KMTC 12.0 PPM BEANS 12N 8 T11N 12 N8 1-08 KMTC 5.0 PPM CAT LEAVES 12N 8 T11N 12 N8 1-08 KMTC 7.0 PPM CAT TIPS 13B 1 T11N 13 B1 1-08 KMTC 10.0 PPM CAT TIPS 13B 1 T11N 13 B1 1-08 KMTC 10.0 PPM CAT TIPS 13B 1 T11N 14 B1 1-08 KMTC 10.0 PPM CEDAR 13G11 T11N 14 B12 1-08 KMTC 5.0 PPM FESCUE TIP 14B12 T11N 16 B10 1-09 KMTC 9.0 PPM FESCUE TIP 16L12 T11N								PPM	WHEAT
11F 7 T11N 11 F7 1-08 KMTC 43.0 PPM CAT LEAVES 11F 7 T11N 11 F7 1-08 KMTC 43.0 PPM WILLOW 12N 8 T11N 12 N8 1-08 KMTC 5.0 PPM WILLOW 12N 8 T11N 12 N8 1-08 KMTC 7.0 PPM CAT TIPS 12N 8 T11N 12 N8 1-08 KMTC 2.0 PPM CEDAR 13B1 T11N 13 G11 1-09 KMTC 10.0 PPM CEDAR 14B8 4 T11N 14 B12 1-08 KMTC 4.0 PPM FESCUE TIP 14B12 T11N 14 B12 1-08 KMTC 0.0 PPM FESCUE TIP 14B12 T11N 16 B12 1-08 KMTC 0.0 PPM FESCUE TIP 16B10 T11N 18 H10 1-09 KMTC 0.0 PPM TATN TATN TA				F7				PPM	
11F 7 1108 KMTC 43.0 PPM BEANS 11F 7 111N 11 F7 1-08 KMTC 12.0 PPM WILLOW 12N 8 T11N 12 N8 1-08 KMTC 5.0 PPM CAT TIPS 12N 8 T11N 12 N8 1-08 KMTC 2.0 PPM CAT TIPS 13B 1 T11N 13 G11 1-08 KMTC 4.0 PPM CEDAR 13B 1 T11N 13 G11 1-08 KMTC 10.0 PPM CEDAR 148 4 T11N 14 B12 1-08 KMTC 10.0 PPM FESCUE TIP 14812 T11N 16 L12 1-08 KMTC 9.0 PPM FESCUE FESCUE TIP 16112 T11N 16 L12 1-08 KMTC 20.0 PPM FESCUE TSECUE TIP				F7				PPM	CAT LEAVES
11F 7 T11N 11 F7 1-08 KMTC 12.0 PPM WILLOW 12N 8 T11N 12 N8 1-08 KMTC 7.0 PPM CAT TIPS 12N 8 T11N 12 N8 1-08 KMTC 2.0 PPM CAT TIPS 13B1 T11N 13 B1 1-08 KMTC 2.0 PPM CEDAR 13G11 T1N 13 B1 1-08 KMTC 10.0 PPM 148 T11N 14 B4 1-08 KMTC 5.0 PPM 14812 T11N 14 B12 1-08 KMTC 9.0 PPM 16112 T11N 16 E12 1-0 KMTC 9.0 PPM 1884 T11N 18 H10 1-09 KMTC 10.0 PPM 1411 T11N 1 A11 1-09 KMTC 10.0 PPM 12614 T									
12N 8 T11N 12 N8 1-08 KMTC 5.0 PPM 12N 8 T11N 12 N8 1-08 KMTC 7.0 PPM CAT TIPS 12N 8 T11N 12 N8 1-08 KMTC 2.0 PPM CEDAR 13B 1 T11N 13 G11 1-08 KMTC 4.0 PPM 14B 4 T11N 13 G11 1-08 KMTC 10.0 PPM 14B 4 T11N 14 B12 1-08 KMTC 5.0 PPM 14B12 T11N 14 B12 1-08 KMTC 5.0 PPM FESCUE TIP 14B12 T1N 14 B12 1-08 KMTC 4.0 PPM FESCUE TIP 14B12 T1N 16 L12 1-08 KMTC 9.0 PPM FESCUE TIP 16L12 T1N 18 H10 1-09 KMTC 10.0 PPM TAT TIN 2.5 TAT TAT TAT TAT <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td>									
12N 8 T11N 12 N8 1-08 KMTC 7.0 PPM CAT TIPS 12N 8 T11N 12 N8 1-08 KMTC 2.0 PPM CEDAR 13B 1 T11N 13 G11 1-09 KMTC 10.0 PPM 14B 4 T11N 14 B4 1-08 KMTC 10.0 PPM 14B 5 T1N 14 B8 1-08 KMTC 6.0 PPM 14B12 T11N 14 B12 1-08 KMTC 5.0 PPM FESCUE 16B10 T1N 16 B10 1-09 KMTC 50.0 PPM FESCUE 18B4 T11N 18 B4 1-09 KMTC 30.0 PPI 188.0 1-08 KMTC 20.0 PPM 1414 T11N 1 111 111 14 1-08 KMTC 10.0 PPM 22614 111N 2 G14 1-10									
12N 8 T11N 12 N8 1-08 KMTC 2.0 PPM CEDAR 13B1 T11N 13 B1 1-08 KMTC 2.0 PPM 13G11 T1N 13 B1 1-09 KMTC 10.0 PPM 14B 4 T11N 14 B4 1-08 KMTC 12.0 PPM 14B12 T11N 14 B12 1-08 KMTC 4.0 PPM FESCUE TIP 14B12 T11N 14 B12 1-08 KMTC 5.0 PPM FESCUE TIP 16L12 T11N 16 L12 1-10 KMTC 30.0 PPM FESCUE 18H10 T11N 18 H10 1-09 KMTC 20.0 PPM 23A10 T1N 23 A10 1-08 KMTC 10.0 PPM 25H 3 T12N 25 H3 1-10 KMTC 10.0 PPM									CAT TIPS
138 1 111N 13 B1 1-08 KMTC 4.0 PPM 13611 T11N 13 G11 1-09 KMTC 12.0 PPM 148 4 T11N 14 B4 1-08 KMTC 6.0 PPM 148.2 T11N 14 B12 1-08 KMTC 6.0 PPM 148.2 T11N 14 B12 1-08 KMTC 4.0 PPM FESCUE TIP 148.2 T11N 14 B12 1-08 KMTC 9.0 PPM FESCUE FESCUE TIP 161.2 T11N 16 L12 1-10 KMTC 30.0 PPM FESCUE TIP 188.4 T11N 18 B4 1-09 KMTC 20.0 PPM T1N 14 T1N 14 10.0 NMTC 20.0 PPM T1N 14 11.0 14.0 14.0 NMTC 10.0 PPM T1N 14.1 10.0 KMTC 10.0 PPM T1N 14.1									
13G11 T11N 13 G11 1-09 KMTC 10.0 PPM 148 T11N 14 B4 1-08 KMTC 12.0 PPM 148 T11N 14 B4 1-08 KMTC 5.0 PPM 14812 T11N 14 B12 1-08 KMTC 5.0 PPM FESCUE TIP 14812 T11N 14 B12 1-08 KMTC 5.0 PPM FESCUE TIP 16112 T11N 16 B10 1-09 KMTC 50.0 PPM FESCUE TIP 1884 T11N 18 B4 1-09 KMTC 40.0 PPM T1N 14 111 109 KMTC 10.0 PPM 24614 111N 22 G14 1-10 KMTC 10.0 PPM 23A10 10.0 KMTC 10.0 PPM 25H 3 T12N 25 H3 1-10 KMTC 10.0									
148 4 T11N 14 B4 1-08 KMTC 12.0 PPM 148 8 T11N 14 B8 1-08 KMTC 6.0 PPM 14812 T11N 14 B12 1-08 KMTC 5.0 PPM FESCUE FIP 14812 T11N 14 B12 1-08 KMTC 4.0 PPM FESCUE FIP 16810 T11N 16 B10 1-09 KMTC 9.0 PPM FESCUE FESCUE 1884 T11N 16 L12 1-10 KMTC 9.0 PPM FESCUE FESCUE FESCUE T1N 14 140 9.0 PPM FESCUE									
148 8 T11N 14 B8 1-08 KMTC 6.0 PPM 14812 T11N 14 B12 1-08 KMTC 5.0 PPM FESCUE T1P 14812 T11N 14 B12 1-08 KMTC 4.0 PPM FESCUE FESCUE 16810 T11N 16 B10 1-09 KMTC 9.0 PPM 16112 T11N 16 L12 1-10 KMTC 30.0 PPM 18810 T1N 18 B4 1-09 KMTC 9.0 PPM 1A11 T1N 18 H10 1-09 KMTC 1.0 PPM 22614 T1N 22 G14 1-0 KMTC 1.0 PPM 23A10 T1N 23 A10 1-08 KMTC 10.0 PPM 25H 3 T12N 25 H3 1-10 KMTC 10.0 PPM 266 T12N 26 G8 1-10 KMTC 10.0 PPM									
14812 T11N 14 B12 1-08 KMTC 5.0 PPM FESCUE T1P 14812 T11N 14 B12 1-08 KMTC 4.0 PPM FESCUE 16810 T11N 16 B10 1-09 KMTC 9.0 PPM 18810 T11N 16 L12 1-10 KMTC 50.0 PPM 18810 T11N 16 L12 1-10 KMTC 9.0 PPM 1811 T11N 18 H10 1-09 KMTC 20.0 PPM 1411 T11N 1 H14 1-08 KMTC 1.0 PPM 23A10 T1N 23 A10 1-08 KMTC 10.0 PPM 25H 3 T12N 25 H3 1-10 KMTC 10.0 PPM 266 1 T12N 27 A13 1-08 KMTC 4.0 PPM 27A 13 T12N									
14812 T11N 14 B12 1-08 KMTC 4.0 PPM FESCUE 18810 T11N 16 B10 1-09 KMTC 9.0 PPM 16L12 T11N 16 L12 1-10 KMTC 50.0 PPM 1884 T11N 18 H10 1-09 KMTC 9.0 PPM 188110 T11N 18 H10 1-09 KMTC 9.0 PPM 1A11 T11N 1 A11 1-09 KMTC 20.0 PPM 22614 T11N 22 G14 1-10 KMTC 8.0 PPM 23A10 T11N 23 A10 1-08 KMTC 8.0 PPM 25H 3 T12N 25 H3 1-10 KMTC 10.0 PPM 26G 8 T12N 26 G1 1-10 KMTC 20.0 PPM 27A13 T12N 27 A13 1-08 KMTC 4.0 PPM SMILAX 27D 1 T12N 27 <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td>FESCUE TIP</td>									FESCUE TIP
16810 T11N 16 B10 1-09 KMTC 9.0 PPM 16L12 T11N 16 L12 1-10 KMTC 50.0 PPM 18B4 T11N 18 B4 1-09 KMTC 30.0 PPi 18H10 T11N 18 B4 1-09 KMTC 40.0 PPM 1A11 T11N 1 A11 1-09 KMTC 20.0 PPM 22G14 T11N 2 G14 1-10 KMTC 10.0 PPM 23A10 T1N 23 A10 1-08 KMTC 9.0 PPM 25H 3 T12N 25 H3 1-10 KMTC 10.0 PPM 266 1 T12N 25 H3 1-10 KMTC 0.0 PPM 27A13 T12N 27 A13 1-08 KMTC 4.0 PPM SMILAX 27H 9 T12N 27 A13 1-08 </td <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td>									
I6L12 T11N 16 L12 1-10 KMTC 50.0 PPM 18B 4 T11N 18 B4 1-09 KMTC 30.0 PP' 18H10 T11N 18 H10 1-09 KMTC 9.0 PPM 1A11 T11N 1 A11 1-09 KMTC 20.0 PPM 22G14 T11N 2 G14 1-10 KMTC 1.0 PPM 23A10 T11N 23 A10 1-08 KMTC 1.0 PPM 25E 8 T12N 25 H3 1-10 KMTC 10.0 PPM 25H 3 T12N 25 H3 1-10 KMTC 10.0 PPM 26G 8 T12N 26 G1 1-10 KMTC 10.0 PPM 27A13 T12N 27 A13 1-08 KMTC 4.0 PPM SMILAX 27D 1 T12N 27 D1 1-10<									LOODE
188 4 T11N 18 B4 1-09 KMTC 30.0 PP;' 18H10 T11N 18 H10 1-09 KMTC 9.0 PPM 1A11 T11N 1 A11 1-09 KMTC 40.0 PPM 1H14 T11N 1 H14 1-08 KMTC 20.0 PPM 22G14 T11N 22 G14 1-10 KMTC 1.0 PPM 23A10 T11N 23 A10 1-08 KMTC 1.0 PPM 25H 3 T12N 25 H3 1-10 KMTC 10.0 PPM 25H 3 T12N 25 H3 1-10 KMTC 10.0 PPM 26G 1 T-10 KMTC 10.0 PPM 266 1 1-10 KMTC 20.0 PPM 27A13 T12N 27 A13 1-08 KMTC 10.0 PM									
18H10 T11N 18 H10 1-09 KMTC 9.0 PPM 1A11 T11N 1 A11 1-09 KMTC 40.0 PPM 1H14 T11N 1 H14 1-08 KMTC 20.0 PPM 22G14 T11N 22 G14 1-10 KMTC 1.0 PPM 23A10 T11N 23 A10 1-08 KMTC 8.0 PPM 25H 3 T12N 25 H3 1-10 KMTC 10.0 PPM 26G 1 T12N 25 H3 1-10 KMTC 10.0 PPM 26G 3 T12N 25 H3 1-10 KMTC 0.0 PPM 26G 8 T12N 27 A13 1-08 KMTC 4.0 PPM SMILAX 27D 1 T12N 27 A13 1-08 KMTC 10.0 PPM 27H 9 T12N 27 H9 1-10 </td <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td>									
1A11 T11N 1 A11 1-09 KMTC 40.0 PPM 1H14 T11N 1 H14 1-08 KMTC 20.0 PPM 23A10 T11N 22 G14 1-10 KMTC 1.0 PPM 23A10 T11N 22 G14 1-10 KMTC 1.0 PPM 23A10 T11N 23 A10 1-08 KMTC 1.0 PPM 25E 8 T12N 25 H3 1-10 KMTC 10.0 PPM 25H 3 T12N 25 H3 1-10 KMTC 10.0 PPM 266 8 T12N 26 G1 1-10 KMTC 20.0 PPM 27A13 T12N 27 A13 1-08 KMTC 4.0 PPM SMILAX 27D 1 T12N 27 A13 1-08 KMTC 10.0 PPM 27H 9 T12N 27 H9 1-10 KMTC 20.0 PPM 27K 4 T12N 27 K4									
1H14 T11N 1 H14 1-08 KMTC 20.0 PPM 22G14 T11N 22 G14 1-10 KMTC 1.0 PPM 23A10 T11N 23 A10 1-08 KMTC 8.0 PPM 25E T12N 25 E8 1-10 KMTC 9.0 PPM 25H 3 T12N 25 H3 1-10 KMTC 10.0 PPM 26G 1 T12N 26 G1 1-10 KMTC 10.0 PPM 26G 3 T12N 26 G3 1-10 KMTC 20.0 PPM 27A13 T12N 27 A13 1-08 KMTC 4.0 PPM SMILAX 27D 1 T12N 27 D1 1-10 KMTC 20.0 PPM 27H 9 T12N 27 H9 1-10 KMTC 20.0 PPM 27H 9 T12N 27 H9 1-10 KMTC 20.0 PPM 286 4 T12N 28 G4									
22G14 T11N 22 G14 1-10 KMTC 1.0 PPM 23A10 T11N 23 A10 1-08 KMTC 8.0 PPM T11N 25E T12N 25 E8 1-10 KMTC 10.0 PPM 25H 3 T12N 25 H3 1-10 KMTC 10.0 PPM 26G 1 T12N 25 H3 1-10 KMTC 10.0 PPM 26G 1 T12N 26 G1 1-10 KMTC 10.0 PPM 26G 8 T12N 26 G8 1-10 KMTC 20.0 PPM 27A13 T12N 27 A13 1-08 KMTC 4.0 PPM SMILAX 27D T12N 27 D1 1-10 KMTC 10.0 PPM 27H 9 T12N 27 H9 1-10 KMTC 20.0 PPM 27K 4 T12N 27 K4 1-10 KMTC 20.0 PPM 2801 T12N 28									
23A10 T11N 23 A10 1-08 KMTC 8.0 PPM T11N 25E 8 T12N 25 E8 1-10 KMTC 10.0 PPM 25H 3 T12N 25 H3 1-10 KMTC 9.0 PPM 25H 3 T12N 25 H3 1-10 KMTC 9.0 PPM 266 1 T12N 25 H3 1-10 KMTC 10.0 PPM 266 1 T12N 26 G8 1-10 KMTC 20.0 PPM 27A13 T12N 27 A13 1-08 KMTC 4.0 PPM SMILAX 27D 1 T12N 27 D1 1-10 KMTC 80.0 PPM 27H 9 T12N 27 H9 1-10 KMTC 20.0 PPM 27H 9 T12N 27 K4 1-10 KMTC 20.0 PPM 286 4 T12N 28 G4 1-10 KMTC 20.0 <t< td=""><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td></t<>									
25E 8 T12N 25 E8 1-10 KMTC 10.0 PPM 25H 3 T12N 25 H3 1-10 KMTC 9.0 PPM 26G 1 T12N 25 H3 1-10 KMTC 10.0 PPM 26G 1 T12N 26 G1 1-10 KMTC 10.0 PPM 26G 8 T12N 26 G8 1-10 KMTC 20.0 PPM 27A13 T12N 27 A13 1-08 KMTC 4.0 PPM SMILAX 27D 1 T12N 27 D1 1-10 KMTC 10.0 PPM SMILAX 27D 1 T12N 27 D1 1-10 KMTC 10.0 PPM SMILAX 27H 9 T12N 27 H9 1-10 KMTC 20.0 PPM 286 4 T12N 27 K4 1-10 KMTC 20.0 PPM 280 2 T12N 28 G4 1-10 KMTC 20.0 PPM 280 3 T12N									T11N
25H 3 T12N 25 H3 1-10 KMTC 9.0 PPM 25H 3 T12N 25 H3 1-10 KMTC 10.0 PPM 26G 1 T12N 26 G1 1-10 KMTC 10.0 PPM 26G 8 T12N 26 G8 1-10 KMTC 20.0 PPM 27A13 T12N 27 A13 1-08 KMTC 6.0 PPM CEDAR 27A13 T12N 27 A13 1-08 KMTC 4.0 PPM SMILAX 27D 1 T12N 27 D1 1-10 KMTC 80.0 PPM 27H 9 T12N 27 H9 1-10 KMTC 20.0 PPM 27K 4 T12N 27 K4 1-10 KMTC 20.0 PPM 28G 4 T12N 28 G4 1-10 KMTC 20.0 PPM 28J11 T12N 28 J11 1-10 KMTC 20.0 PPM 28N 2 T12N 28 <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td>ITTH</td>									ITTH
25H 3 T12N 25 H3 1-10 KMTC 10.0 PPM 26G 1 T12N 26 G1 1-10 KMTC 10.0 PPM 26G 8 T12N 26 G8 1-10 KMTC 20.0 PPM 27A13 T12N 27 A13 1-08 KMTC 6.0 PPM CEDAR 27A13 T12N 27 A13 1-08 KMTC 4.0 PPM SMILAX 27D 1 T12N 27 D1 1-10 KMTC 80.0 PPM 27H 9 T12N 27 H9 1-10 KMTC 10.0 PPM 27H 9 T12N 27 H9 1-10 KMTC 20.0 PPM 27K 4 T12N 27 K4 1-10 KMTC 20.0 PPM 286 4 T12N 28 G4 1-10 KMTC 20.0 PPM 28N 2 T12N 28 N2 1-10 KMTC 20.0 PPM 28N 2 T12N 28 <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td>									
266 1 T12N 26 G1 1-10 KMTC 10.0 PPM 26G 8 T12N 26 G8 1-10 KMTC 20.0 PPM 27A13 T12N 27 A13 1-08 KMTC 6.0 PPM CEDAR 27A13 T12N 27 A13 1-08 KMTC 4.0 PPM SMILAX 27D 1 T12N 27 D1 1-10 KMTC 80.0 PPM 27H 9 T12N 27 H9 1-10 KMTC 10.0 PPM 27H 9 T12N 27 H9 1-10 KMTC 20.0 PPM 27K 4 T12N 27 H9 1-10 KMTC 20.0 PPM 286 4 T12N 28 G4 1-10 KMTC 10.0 PPM 28N 2 T12N 28 G4 1-10 KMTC 20.0 PPM 28N 2 T12N 28 N2 1-10 KMTC 20.0 PPM 28N 2 T12N 28 <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td>									
26G 8 T12N 26 G8 1-10 KMTC 20.0 PPM 27A13 T12N 27 A13 1-08 KMTC 6.0 PPM CEDAR 27A13 T12N 27 A13 1-08 KMTC 4.0 PPM SMILAX 27D 1 T12N 27 D1 1-10 KMTC 80.0 PPM 27H 9 T12N 27 H9 1-10 KMTC 20.0 PPM 27H 9 T12N 27 H9 1-10 KMTC 20.0 PPM 27K 4 T12N 27 K4 1-10 KMTC 20.0 PPM 286 4 T12N 28 G4 1-10 KMTC 20.0 PPM 28N 2 T12N 28 J11 1-10 KMTC 20.0 PPM 28N 2 T12N 28 N2 1-10 KMTC 20.0 PPM 29D 3 T12N 28 N2 1-10 KMTC 20.0 PPM 213 T1N 2									
27A13 T12N 27 A13 1-08 KMTC 6.0 PPM CEDAR 27A13 T12N 27 A13 1-08 KMTC 4.0 PPM SMILAX 27D 1 T12N 27 D1 1-10 KMTC 80.0 PPM SMILAX 27D 1 T12N 27 D1 1-10 KMTC 80.0 PPM 27H 9 T12N 27 H9 1-10 KMTC 10.0 PPM 27K 4 T12N 27 K4 1-10 KMTC 20.0 PPM 286 4 T12N 28 G4 1-10 KMTC 20.0 PPM 280 2 T12N 28 G4 1-10 KMTC 20.0 PPM 280 2 T12N 28 J11 1-10 KMTC 20.0 PPM 28N 2 T12N 28 N2 1-10 KMTC 20.0 PPM 28N 2 T12N 28 N2 1-10 KMTC 2.0 PPM 29D 3 T12N									
27A13 T12N 27 A13 1-08 KMTC 4.0 PPM SMILAX 27D 1 T12N 27 D1 1-10 KMTC 80.0 PPM 27H 9 T12N 27 H9 1-10 KMTC 10.0 PPM 27H 9 T12N 27 H9 1-10 KMTC 20.0 PPM 27K 4 T12N 27 K4 1-10 KMTC 20.0 PPM 286 4 T12N 28 G4 1-10 KMTC 20.0 PPM 280 4 T12N 28 G4 1-10 KMTC 20.0 PPM 2811 T12N 28 G4 1-10 KMTC 20.0 PPM 28N 2 T12N 28 N2 1-10 KMTC 20.0 PPM 28N 2 T12N 28 N2 1-10 KMTC 20.0 PPM 29D 3 T12N 29 D3 1-10 KMTC 20.0 PPM <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td>CEDAD</td>									CEDAD
27D 1 T12N 27 D1 1-10 KMTC 80.0 PPM 27H 9 T12N 27 H9 1-10 KMTC 10.0 PPM 27H 9 T12N 27 H9 1-10 KMTC 20.0 PPM 27K 4 T12N 27 K4 1-10 KMTC 20.0 PPM 28G 4 T12N 28 G4 1-10 KMTC 1.0 PPM 28G 4 T12N 28 G4 1-10 KMTC 20.0 PPM 28G 4 T12N 28 G4 1-10 KMTC 20.0 PPM 28J11 T12N 28 J11 1-10 KMTC 20.0 PPM 28N 2 T12N 28 N2 1-10 KMTC 20.0 PPM 29D 3 T12N 28 N2 1-10 KMTC 2.0 PPM 213 T11N 2 E14 1-09 KMTC 20.0 PPM 2M 8 T11N 2 M8 1-09 <t< td=""><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td></t<>									
27H 9 T12N 27 H9 1-10 KMTC 10.0 PPM 27H 9 T12N 27 H9 1-10 KMTC 20.0 PPM 27K 4 T12N 27 K4 1-10 KMTC 20.0 PPM 28G 4 T12N 28 G4 1-10 KMTC 1.0 PPM 28G 4 T12N 28 G4 1-10 KMTC 20.0 PPM 28G 4 T12N 28 G4 1-10 KMTC 20.0 PPM 28J11 T12N 28 J11 1-10 KMTC 20.0 PPM 28N 2 T12N 28 J11 1-10 KMTC 20.0 PPM 28N 2 T12N 28 N2 1-10 KMTC 10.0 PPM 29D 3 T12N 28 N2 1-10 KMTC 20.0 PPM 213 T11N 2 E14 1-09 KMTC 60.0 PPM 2M 8 T11N 2 M8 1-09									SMILAX
27H 9T12N27H91-10KMTC20.0PPM27K 4T12N27K41-10KMTC20.0PPM28G 4T12N28G41-10KMTC1.0PPM28G 4T12N28G41-10KMTC20.0PPM28J11T12N28J111-10KMTC20.0PPM28N 2T12N28N21-10KMTC20.0PPM28N 2T12N28N21-10KMTC2.0PPM28N 2T12N28N21-10KMTC2.0PPM29D 3T12N29D31-10KMTC20.0PPM2E14T11N2E141-09KMTC20.0PPM2L13T11N2L131-09KMTC60.0PPM30F 7T12N30F71-10KMTC6.0PPM30F 7T12N30F71-10KMTC6.0PPM31B10T12N31B101-10KMTC10.0PPM33B14T12N33B141-10KMTC8.0PPM33I 3T12N33131-10KMTC20.0PPM									
27K 4 T12N 27 K4 1-10 KMTC 20.0 PPM 28G 4 T12N 28 G4 1-10 KMTC 1.0 PPM 28G 4 T12N 28 G4 1-10 KMTC 20.0 PPM 28J11 T12N 28 G4 1-10 KMTC 20.0 PPM 28J11 T12N 28 J11 1-10 KMTC 20.0 PPM 28J11 T12N 28 J11 1-10 KMTC 20.0 PPM 28N 2 T12N 28 N2 1-10 KMTC 2.0 PPM 28N 2 T12N 28 N2 1-10 KMTC 2.0 PPM 29D 3 T12N 29 D3 1-10 KMTC 2.0 PPM 2E14 T1N 2 E14 1-09 KMTC 20.0 PPM 2M 8 T1N 2 M8 1-09 KMTC 40.0 PPM 30F 7 T12N 30 F7 1-10									
28G 4 T12N 28 G4 1-10 KMTC 1.0 PPM 28G 4 T12N 28 G4 1-10 KMTC 20.0 PPM 28J11 T12N 28 J11 1-10 KMTC 20.0 PPM 28J11 T12N 28 J11 1-10 KMTC 20.0 PPM 28N 2 T12N 28 N2 1-10 KMTC 2.0 PPM 28N 2 T12N 28 N2 1-10 KMTC 10.0 PPM 29D 3 T12N 29 D3 1-10 KMTC 2.0 PPM 2E14 T1N 2 E14 1-09 KMTC 20.0 PPM 2L13 T1N 2 L13 1-09 KMTC 60.0 PPM 2M 8 T11N 2 M8 1-09 KMTC 40.0 PPM 30F 7 T12N 30 F7 1-10 KMTC 6.0 PPM 31B10 T12N 31 B10 1-10									
28G 4 T12N 28 G4 1-10 KMTC 20.0 PPM 28J11 T12N 28 J11 1-10 KMTC 20.0 PPM 28N 2 T12N 28 N2 1-10 KMTC 2.0 PPM 28N 2 T12N 28 N2 1-10 KMTC 2.0 PPM 28N 2 T12N 28 N2 1-10 KMTC 10.0 PPM 29D 3 T12N 29 D3 1-10 KMTC 2.0 PPM 29L13 T11N 2 E14 1-09 KMTC 60.0 PPM 2L13 T11N 2 L13 1-09 KMTC 40.0 PPM 2M 8 T11N 2 M8 1-09 KMTC 6.0 PPM 30F 7 T12N 30 F7 1-10 KMTC 6.0 PPM 31B10 T12N 31 B10 1-10 KMTC 10.0 PPM 33B14 T12N 33 B14 1-10 <t< td=""><td></td><td></td><td></td><td></td><td></td><td>and the second second</td><td></td><td></td><td></td></t<>						and the second			
28J11 T12N 28 J11 1-10 KMTC 20.0 PPM 28N 2 T12N 28 N2 1-10 KMTC 2.0 PPM 28N 2 T12N 28 N2 1-10 KMTC 2.0 PPM 28N 2 T12N 28 N2 1-10 KMTC 10.0 PPM 29D 3 T12N 29 D3 1-10 KMTC 20.0 PPM 2E14 T11N 2 E14 1-09 KMTC 20.0 PPM 2L13 T11N 2 L13 1-09 KMTC 60.0 PPM 2M 8 T11N 2 M8 1-09 KMTC 40.0 PPM 30F 7 T12N 30 F7 1-10 KMTC 6.0 PPM 30F 7 T12N 30 F7 1-10 KMTC 10.0 PPM 31B10 T12N 31 B10 1-10 KMTC 10.0 PPM 33B14 T12N 33 B14 1-10 <									
28N 2 T12N 28 N2 1-10 KMTC 2.0 PPM 28N 2 T12N 28 N2 1-10 KMTC 10.0 PPM 29D 3 T12N 29 D3 1-10 KMTC 2.0 PPM 29D 3 T12N 29 D3 1-10 KMTC 2.0 PPM 2E14 T11N 2 E14 1-09 KMTC 20.0 PPM 2L13 T11N 2 L13 1-09 KMTC 60.0 PPM 2M 8 T11N 2 M8 1-09 KMTC 40.0 PPM 30F 7 T12N 30 F7 1-10 KMTC 6.0 PPM 30F 7 T12N 30 F7 1-10 KMTC 6.0 PPM 31B10 T12N 31 B10 1-10 KMTC 10.0 PPM 33814 T12N 33 B14 1-10 KMTC 20.0 PPM 33I 3 T12N 33 13 1-10									
28N 2 T12N 28 N2 1-10 KMTC 10.0 PPM 29D 3 T12N 29 D3 1-10 KMTC 2.0 PPM 2E14 T11N 2 E14 1-09 KMTC 20.0 PPM 2L13 T11N 2 L13 1-09 KMTC 60.0 PPM 2M 8 T11N 2 M8 1-09 KMTC 60.0 PPM 30F 7 T12N 30 F7 1-10 KMTC 6.0 PPM 30F 7 T12N 30 F7 1-10 KMTC 6.0 PPM 31B10 T12N 31 B10 1-10 KMTC 10.0 PPM 33B14 T12N 33 B14 1-10 KMTC 8.0 PPM 33I 3 T12N 33 I3 1-10 KMTC 20.0 PPM									
29D 3 T12N 29 D3 1-10 KMTC 2.0 PPM 2E14 T11N 2 E14 1-09 KMTC 20.0 PPM 2L13 T11N 2 L13 1-09 KMTC 60.0 PPM 2M 8 T11N 2 M8 1-09 KMTC 40.0 PPM 30F 7 T12N 30 F7 1-10 KMTC 6.0 PPM 30F 7 T12N 30 F7 1-10 KMTC 6.0 PPM 30F 7 T12N 30 F7 1-10 KMTC 6.0 PPM 30F 7 T12N 30 F7 1-10 KMTC 6.0 PPM 31B10 T12N 31 B10 1-10 KMTC 10.0 PPM 33B14 T12N 33 B14 1-10 KMTC 20.0 PPM 33I 3 T12N 33 13 1-10 KMTC 20.0 PPM									
2E14 T11N 2 E14 1-09 KMTC 20.0 PPM 2L13 T11N 2 L13 1-09 KMTC 60.0 PPM 2M 8 T11N 2 M8 1-09 KMTC 40.0 PPM 30F 7 T12N 30 F7 1-10 KMTC 6.0 PPM 30F 7 T12N 30 F7 1-10 KMTC 6.0 PPM 30F 7 T12N 30 F7 1-10 KMTC 6.0 PPM 30F 7 T12N 30 F7 1-10 KMTC 6.0 PPM 31B10 T12N 31 B10 1-10 KMTC 10.0 PPM 33B14 T12N 33 B14 1-10 KMTC 20.0 PPM 33I 3 T12N 33 13 1-10 KMTC 20.0 PPM									
2L13 T11N 2 L13 1-09 KMTC 60.0 PPM 2M 8 T11N 2 M8 1-09 KMTC 40.0 PPM 30F 7 T12N 30 F7 1-10 KMTC 6.0 PPM 30F 7 T12N 30 F7 1-10 KMTC 6.0 PPM 30F 7 T12N 30 F7 1-10 KMTC 6.0 PPM 30F 7 T12N 30 F7 1-10 KMTC 10.0 PPM 31B10 T12N 31 B10 1-10 KMTC 10.0 PPM 33B14 T12N 33 B14 1-10 KMTC 8.0 PPM 33I 3 T12N 33 13 1-10 KMTC 20.0 PPM									
2M 8 T11N 2 M8 1-09 KMTC 40.0 PPM 30F 7 T12N 30 F7 1-10 KMTC 6.0 PPM 30F 7 T12N 30 F7 1-10 KMTC 6.0 PPM 31B10 T12N 31 B10 1-10 KMTC 10.0 PPM 33B14 T12N 33 B14 1-10 KMTC 8.0 PPM 33I 3 T12N 33 I3 1-10 KMTC 20.0 PPM		T11N							
30F 7 T12N 30 F7 1-10 KMTC 6.0 PPM 30F 7 T12N 30 F7 1-10 KMTC 6.0 PPM 31B10 T12N 31 B10 1-10 KMTC 10.0 PPM 33B14 T12N 33 B14 1-10 KMTC 8.0 PPM 33I 3 T12N 33 I3 1-10 KMTC 20.0 PPM	2L13		2						
30F 7 T12N 30 F7 1-10 KMTC 6.0 PPM 31B10 T12N 31 B10 1-10 KMTC 10.0 PPM 33B14 T12N 33 B14 1-10 KMTC 8.0 PPM 33I 3 T12N 33 I3 1-10 KMTC 20.0 PPM									
31B10 T12N 31 B10 1-10 KMTC 10.0 PPM 33B14 T12N 33 B14 1-10 KMTC 8.0 PPM 33I 3 T12N 33 I3 1-10 KMTC 20.0 PPM	30F 7	T12N							
33B14 T12N 33 B14 1-10 KMTC 8.0 PPM 33I 3 T12N 33 I3 1-10 KMTC 20.0 PPM		T12N	30						
33B14 T12N 33 B14 1-10 KMTC 8.0 PPM 33I 3 T12N 33 I3 1-10 KMTC 20.0 PPM	31B10	T12N	31	B10	1-10	KMTC	10.0		
	33B14	T12N	33				8.0	PPM	
	33I 3	T12N	33	13	1-10	KMTC	20.0		
	33M14	T12N	33	M14	1-10	KMTC	6.0	PPM	

LABEL	TOWNSHIP	SECTION	SUBSECT	DATE	LAB	RESULTS	UNITS	COMMENTS
33M14	T12N	33	M14	1-10	KMTC	5.0	PPM	
34D 7	T12N	34	D7	1-10	KMTC	30.0	PPM	
35A 9	T12N	35	A9	1-10	KMTC	6.0	PPM	
35A 9	T12N	35	A9	1-10	KMTC	11.0	PPM	
35B 1	T12N	35	B1	1-10	KMTC	9.0	PPM	
35B14	T12N	35	814	1-10	KMTC	3.0	PPM	
35F10	T12N	35	F10	1-10	KMTC	20.0	PPM	
35L 5	T12N	35	L5	1-09	KMTC	5.0	PPM	
35L 5	T12N	35	L5	1-09	KMTC	8.0	PPM	
35L11	T12N	35	L11	1-10	KMTC	30.0	PPM	
35L11	T12N	35	L11	1-10	KMTC	<1	PPM	
36A13	T12N	36	A13	1-10	KMTC	20.0	PPM	
36G 5	T12N	36	G5	1-09	KMTC	20.0	PPM	
36H14	T12N	36	H14	1-10	KMTC	<1	PPM	
3F 7	TIIN	3	F7	1-09	KMTC	8.0	PPM	
4A 1	T11N	4	Al	1-10	KMTC	6.0	PPM	
4A 3	TIIN	4	A3	1-10	KMTC	20.0	PPM	
4A 6	TIIN	4	AG	1-10	KMTC	60.0	PPM	
4J14	T11N	4	J14	1-09	KMTC	15.0	PPM	
4J14	TIIN	4	J14	1-09	KMTC	6.0	PPM	
4N 7	TIIN	4	N7	1-09	KMTC	30.0	PPM	
4N 8	TIIN	4	N8	1-09	KMTC	6.0	PPM	
50/HWY	T12N	22	G1	1-06	KMTC	1500.0	PPM	CEDAR
5N 6	TIIN	5	NG	1-09	KMTC	30.0	PPM	CEDAR
6A 7	TIIN	6	A7	1-09	KMTC	30.0	PPM	
6C14	TIIN	6	C14	1-09	KMTC	10.0	PPM	
8A 1	TIIN	8	Al	1-09	KMTC	8.0	PPM	
8J 7	TIIN	8	J7	1-09	KMTC	10.0	PPM	
8M14	TIIN	8	M14	1-09	KMTC	9.0	PPM	
9B 2	TIIN	9	B2	1-09	KMTC	5.5	PPM	
BG 1/4	T12N	21	B13	1-06	KMTC	40.0	PPM	1/4 N.GATE
CRES08	T12N	27	K3	1-04	KMTC	42.0	PPM	CRES
DRES12	T12N	27	K4	1-04	KMTC	160.0	PPM	DRES/LEAVES
H-02	T12N	27	I2	1-04	KMTC	24.0	PPM	CEDAR
H-03	T12N	27		1-04	KMTC	80.0		HRES
LRES04	T12N	28	I14	1-04	KMTC	3.0	PPM	TREE
RKRES15		34	A13	1-04	KMTC	32.0	PPM	RKRES/CEDAR
RRES18	TIIN	34	B6	1-04	KMTC	7.0	PPM	RRES/CEDAR
SMRES21		26	K6	1-04	KMTC	5.0	PPM	SMRES/LEAVES
TRES53	TIIN	23	M10	1-06	KMTC	23.0	PPM	TRES
YRES51	TIIN	26	N3	1-05	KMTC	2600.0	PPM	YRES
BG 1/4	T12N	21	B13	1-06	ORNL	7.3	UG/G	TRES
CARLILE		22	D14	1-05	ORNL	8.3	UG/G	
T4-01E	TIIN	2	C11	1-06	ORNL	3.1	UG/G	
T4-02E	TIIN	2	B8	1-06	ORNL	2.1	UG/G	
T4-03E	T12N	35	N5	1-06	ORNL	2.1		
T4-04E	T12N	35	L3	1-06	ORNL	2.9	UG/G	
T4-04W	T12N	34	L13	1-06	ORNL	6.1	UG/G UG/G	
T4-05E	T12N	35	K1	1-06	ORNL			
14 052	1 2 2 11	55	N.L.	1-00	ORNE	9.3	UG/G	

AB-01 TIN 17 A14 1-15 OSDH 2.72 MG/KG AB-02 TIN 16 A11 1-15 OSDH 1.68 MG/KG AB-04 TIN 9 B1 1-15 OSDH 1.52 MG/KG AB-04 TIN 4 C8 1-15 OSDH 2.64 MG/KG AB-05 TIN 4 C8 1-15 OSDH 2.64 MG/KG AB-06 TI2N 31 A8 1-15 OSDH 2.64 MG/KG AB-07 T12N 13 G14 1-15 OSDH 2.96 MG/KG AB-08 T12N 13 G14 1-06 OSDH 4.3 MG/KG CARLLE T12N 12 D14 1-06 OSDH 2.96 MG/KG CARLIE T12N 28 E10 1-06 OSDH 2.9 MG/KG CNTRLO3 <t11n< td=""> 28 E10 1-06 OS</t11n<>	LABEL	TOWNSHIP	SECTION	SUBSECT	DATE	LAB	RESULTS	UNITS	COMMENTS
AB-03 T12N 32 N14 1-15 OSDH 3.28 MG/KG AB-04 T11N 9 B1 1-15 OSDH 1.52 MG/KG AB-05 T12N 29 N14 1-15 OSDH 2.64 MG/KG AB-07 T12N 31 A8 1-15 OSDH 1.60 MG/KG AB-07 T12N 18 M6 1-15 OSDH 1.60 MG/KG AB-07 T12N 18 M6 1-15 OSDH 1.60 MG/KG AB-09 T12N 13 G14 1-15 OSDH 2.96 MG/KG CARLILE T12N 22 D14 1-06 OSDH 2.96 MG/KG CARLILE T12N 22 D14 1-06 OSDH 2.9 MG/KG CARLILE T12N 22 D14 1-06 OSDH 2.9 MG/KG CNTRLO2 T10N 35 C14	AB-01	TIIN	17	A14	1-15	OSDH	2.72	MG/KG	
AB-04 T11N 9 B1 1-15 OSDH 1.52 MG/KG AB-05 T11N 4 C8 1-15 OSDH 2.64 MG/KG AB-06 T12N 31 A8 1-15 OSDH 2.64 MG/KG AB-07 T12N 31 A8 1-15 OSDH 2.64 MG/KG AB-07 T12N 13 G14 1-15 OSDH 2.60 MG/KG AB-09 T12N 13 G14 1-05 OSDH 4.3 MG/KG CARLILE T12N 22 D14 1-14 OSDH 2.56 MG/KG CARLILE T1N 22 D14 1-06 OSDH 2.96 MG/KG CNTRLO3 T1N 28 E10 1-06 OSDH 2.4 MG/KG CNTRLO3 T1N 28 E10 1-06 OSDH 2.5 MG/KG CNTRLO5 T1N 13 P F1	AB-02	TIIN	16	A11	1-15	OSDH	1.68	MG/KG	
AB-05 T11N 4 C8 1-15 OSDH 2.64 MG/KG AB-06 T12N 29 N14 1-15 OSDH 1.68 MG/KG AB-07 T12N 19 N8 1-15 OSDH 1.68 MG/KG AB-09 T12N 18 M6 1-15 OSDH 2.60 MG/KG AB-10 T12N 12 014 1-05 OSDH 4.60 MG/KG CARLILE T12N 22 D14 1-05 OSDH 2.96 MG/KG CARLILE T12N 22 D14 1-06 OSDH 2.96 MG/KG CARLICE T1N 22 D14 1-06 OSDH 2.96 MG/KG CNTRLO2 T1N 28 E10 1-06 OSDH 2.4 MG/KG CNTRLO5 T11N 19 F13 1-06 OSDH 2.4 MG/KG CNTRLO5 T11N 12 C12	AB-03	T12N	32	N14	1-15	OSDH	3.28	MG/KG	
AB-05 T11N 4 C8 1-15 OSDH 2.64 MG/KG AB-06 T12N 19 N14 1-15 OSDH 1.68 MG/KG AB-07 T12N 19 N8 1-15 OSDH 2.40 MG/KG AB-09 T12N 18 M6 1-15 OSDH 2.60 MG/KG AB-01 T12N 13 G14 1-15 OSDH 2.96 MG/KG CARLILE T12N 22 D14 1-06 OSDH 2.96 MG/KG CARLILE T12N 22 D14 1-06 OSDH 2.96 MG/KG CARLOT T1N 28 E10 1-06 OSDH 2.9 MG/KG CNTRLO3 T1N 28 E10 1-06 OSDH 2.4 MG/KG CNTRLO5 T1N 19 F13 1-06 OSDH 2.4 MG/KG CNTRLO5 T1N 19 F13 1-06 OSDH 2.7 MG/KG CNTRLO5 T1N 21 <t< td=""><td>AB-04</td><td>TIIN</td><td>9</td><td>B1</td><td>1-15</td><td>OSDH</td><td>1.52</td><td>MG/KG</td><td></td></t<>	AB-04	TIIN	9	B1	1-15	OSDH	1.52	MG/KG	
AB-07 T12N 31 A8 1-15 OSDH 2.40 MG/KG AB-08 T12N 19 N8 1-15 OSDH 3.60 MG/KG AB-09 T12N 13 G14 1-15 OSDH 2.96 MG/KG CARLILE T12N 13 G14 1-15 OSDH 2.96 MG/KG CARLILE T12N 22 D14 1-14 OSDH 2.96 MG/KG CARLILE T12N 22 D14 1-14 OSDH 2.96 MG/KG CNTRLO2 T10N 35 C14 1-06 OSDH 2.56 MG/KG CNTRLO3 T11N 28 E10 1-06 OSDH 2.5 MG/KG CNTRLO5 T11N 19 F13 1-06 OSDH 2.2 MG/KG CNTRLO5 T11N 13 106 OSDH 2.7 MG/KG CNTRL07 T11N 13 106 OSDH <td>AB-05</td> <td>T11N</td> <td></td> <td>C8</td> <td>1-15</td> <td>OSDH</td> <td>2.64</td> <td>MG/KG</td> <td></td>	AB-05	T11N		C8	1-15	OSDH	2.64	MG/KG	
AB-07 T12N 31 A8 1-15 OSDH 2.40 MG/KG AB-08 T12N 19 N8 1-15 OSDH 3.60 MG/KG AB-09 T12N 13 G14 1-15 OSDH 3.60 MG/KG CARLILE T12N 13 G14 1-15 OSDH 4.3 MG/KG CARLILE T12N 22 D14 1-06 OSDH 0.9 MG/KG CNTRL01 T09N 16 - 1-06 OSDH 1.9 MG/KG CNTRL02 T10N 35 C14 1-06 OSDH 2.5 MG/KG CNTRL05 T11N 19 F13 1-06 OSDH 2.2 MG/KG CNTRL05 T11N 13 B7 1-06 OSDH 2.7 MG/KG CNTRL05 T11N 21 C12 1-06 OSDH 1.9 MG/KG CNTRL07 T11N 21 C12			29				1.68	MG/KG	
AB-08 T12N 19 N8 1-15 OSDH 1.60 MG/KG AB-09 T12N 18 M6 1-15 OSDH 2.96 MG/KG CARLILE T12N 22 D14 1-05 OSDH 4.3 MG/KG CARLILE T12N 22 D14 1-05 OSDH 4.3 MG/KG CARLILE T12N 22 D14 1-06 OSDH 2.96 MG/KG CNTRLO2 T10N 35 C14 1-06 OSDH 2.56 MG/KG CNTRLO3 T11N 28 E10 1-06 OSDH 2.4 MG/KG CNTRLO5 T11N 19 F13 1-06 OSDH 2.2 MG/KG CNTRLO5 T11N 13 B7 1-06 OSDH 2.7 MG/KG CNTRLO8 T11N 21 C12 1-06 OSDH 2.7 MG/KG CNTRLO9 T11N 21 C12				A8				MG/KG	
AB-09 T12N 18 M6 1-15 OSDH 3.60 MG/KG AB-10 T12N 13 G14 1-15 OSDH 2.96 MG/KG CARLILE T12N 22 D14 1-05 OSDH 2.96 MG/KG CARLILE T12N 22 D14 1-14 OSDH 2.56 MG/KG CNTRLOZ TION 35 C14 1-06 OSDH 0.9 MG/KG CNTRLOZ TION 35 C14 1-06 OSDH 2.4 MG/KG CNTRLOS TIIN 28 E10 1-06 OSDH 2.4 MG/KG CNTRLOS TIIN 19 F13 1-06 OSDH 2.2 MG/KG CNTRLOS TIIN 13 B7 1-06 OSDH 2.7 MG/KG CNTRLOS TIIN 21 C12 1-06 OSDH 2.2 MG/KG CNTRLOT TIIN 21 C12 1-06 OSDH 1.9 MG/KG CNTRLOT TIIN 21<	A CONTRACT AND A DOMESTICS						1.60	MG/KG	
AB-10 T12N 13 G14 1-15 OSDH 2.96 MG/KG CARLILE T12N 22 D14 1-05 OSDH 4.3 MG/KG CARLILE T12N 22 D14 1-14 OSDH 2.56 MG/KG CNTRLO1 T09N 16 - 1-06 OSDH 0.9 MG/KG CNTRLO2 T10N 35 C14 1-06 OSDH 2.4 MG/KG CNTRLO3 T11N 28 E10 1-06 OSDH 2.4 MG/KG CNTRLO5 T11N 19 F13 1-06 OSDH 2.4 MG/KG CNTRL05 T11N 13 B7 1-06 OSDH 2.2 MG/KG CNTRL08 T10N 8 K8 1-06 OSDH 2.7 MG/KG CNTRL09 T11N 21 C12 1-06 OSDH 1.9 MG/KG CNTRL07 T11N 21 C12								MG/KG	
CARLILE T12N 22 D14 1-05 OSDH 4.3 MG/KG GRASS CARLILE T12N 22 D14 1-14 OSDH 2.56 MG/KG CNTRLO1 TOSN 16 - 1-06 OSDH 1.9 MG/KG CNTRLO2 TION 35 C14 1-06 OSDH 2.4 MG/KG CNTRLO3 TIIN 28 E10 1-06 OSDH 2.4 MG/KG CNTRLO5 TIIN 19 F13 1-06 OSDH 2.4 MG/KG CNTRLO5 TIIN 13 B7 1-06 OSDH 2.4 MG/KG CNTRLO7 TIN 13 B7 1-06 OSDH 2.7 MG/KG CNTRL07 TIN 21 C12 1-06 OSDH 2.7 MG/KG CNTRL07 TIN 21 C13 1-06 OSDH 1.2 MG/KG CNTRL13 TIN 2							2.96	MG/KG	
CARLILE T12N 22 D14 1-14 OSDH 2.56 MG/KG CNTRLO1 TOSN 16 - 1-06 OSDH 0.9 MG/KG CNTRLO3 T1N 28 E10 1-06 OSDH 2.4 MG/KG CNTRLO4 T1N 20 N10 1-06 OSDH 2.4 MG/KG CNTRLO5 T1N 19 F13 1-06 OSDH 2.4 MG/KG CNTRLO6 T09N 4 113 1-06 OSDH 2.7 MG/KG CNTRLO7 T1N 13 B7 1-06 OSDH 2.7 MG/KG CNTRLO9 T1N 21 C12 1-06 OSDH 2.7 MG/KG CNTRL09 T1N 21 C12 1-06 OSDH 1.9 MG/KG CNTRL10 T1N 21 I13 1-06 OSDH 1.4 MG/KG CNTRL10 T12N 21 I13								MG/KG	GRASS
CNTRL01 TO9N 16 - 1-06 OSDH 0.9 MG/KG CNTRL02 T10N 35 C14 1-06 OSDH 1.9 MG/KG CNTRL03 T11N 28 E10 1-06 OSDH 2.4 MG/KG CNTRL05 T11N 19 F13 1-06 OSDH 2.2 MG/KG CNTRL05 T09N 4 T13 1-06 OSDH 2.2 MG/KG CNTRL07 T11N 13 B7 1-06 OSDH 2.7 MG/KG CNTRL09 T11N 21 C12 1-06 OSDH 2.7 MG/KG CNTRL09 T11N 21 C12 1-06 OSDH 2.2 MG/KG CNTRL10 T11N 21 C12 1-06 OSDH 2.7 MG/KG CNTRL13 T11N 21 C12 1-06 OSDH 1.2 MG/KG CNTRL13 T11N 21 C12									
CNTRL02 TION 35 C14 1-06 OSDH 1.9 MG/KG CNTRL03 T11N 28 E10 1-06 OSDH 2.4 MG/KG CNTRL05 T11N 20 N10 1-06 OSDH 2.5 MG/KG CNTRL05 T11N 19 F13 1-06 OSDH 2.2 MG/KG CNTRL05 T11N 13 B7 1-06 OSDH 0.9 MG/KG CNTRL07 T11N 13 B7 1-06 OSDH 2.2 MG/KG CNTRL09 T11N 21 C12 1-06 OSDH 2.7 MG/KG CNTRL09 T11N 21 C12 1-06 OSDH 1.9 MG/KG CNTRL10 T11N 21 C12 1-06 OSDH 1.2 MG/KG CNTRL13 T1N 21 C12 1-06 OSDH 1.4 MG/KG CNTRL14 T1N Z1 C13									
CNTRL03 T11N 28 E10 1-06 OSDH 2.4 MG/KG CNTRL04 T11N 20 N10 1-06 OSDH 2.5 MG/KG CNTRL05 T11N 19 F13 1-06 OSDH 2.2 MG/KG CNTRL06 T09N 4 113 1-06 OSDH 0.9 MG/KG CNTRL07 T11N 13 B7 1-06 OSDH 2.7 MG/KG CNTRL09 T11N 21 C12 1-06 OSDH 2.2 MG/KG CNTRL09 T11N 21 C12 1-06 OSDH 1.9 MG/KG CNTRL13 T1N 21 C12 1-06 OSDH 1.4 MG/KG CNTRL13 T1N 21 C12 1-06 OSDH 1.4 MG/KG CNTRL13 T1N 21 C13 1-05 OSDH 2.0 MG/KG CNTRL14 T12N Z N1				C14					
CNTRL04 T11N 20 N10 1-06 OSDH 2.5 MG/KG CNTRL05 T11N 19 F13 1-06 OSDH 2.2 MG/KG CNTRL06 T09N 4 113 1-06 OSDH 0.9 MG/KG CNTRL07 T11N 13 B7 1-06 OSDH 2.7 MG/KG CNTRL09 T11N 21 C12 1-06 OSDH 2.2 MG/KG CNTRL09 T11N 21 C12 1-06 OSDH 2.2 MG/KG CNTRL10 T11N 21 C12 1-06 OSDH 3.0 MG/KG CNTRL13 T11N 21 C12 1-06 OSDH 1.4 MG/KG CNTRL13 T11N 3 C2 1-06 OSDH 1.4 MG/KG CRES T12N 26 E3 1-05 OSDH 3.0 MG/KG G-01 T12N 8 L5 <t< td=""><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td></t<>									
CNTRLO5 T11N 19 F13 1-06 OSDH 2.2 MG/KG CNTRL06 T09N 4 I13 1-06 OSDH 1.4 MG/KG CNTRL07 T11N 13 B7 1-06 OSDH 0.9 MG/KG CNTRL08 T10N 8 K8 1-06 OSDH 2.7 MG/KG CNTRL09 T11N 21 C12 1-06 OSDH 2.2 MG/KG CNTRL09 T11N 21 C12 1-06 OSDH 3.0 MG/KG CNTRL10 T11N 21 I13 1-06 OSDH 1.2 MG/KG CNTRL15 T09N 16 G13 1-06 OSDH 1.4 MG/KG G-01 T12N 8 L5 1-15 OSDH 3.0 MG/KG G-02 T12N 7 B13 1-15 OSDH 3.04 MG/KG G-01 T12N 7 K9 1-11 OSDH 2.48 MG/KG G-02 T12N <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td>									
CNTRLO6 TO9N 4 II3 1-06 OSDH 1.4 MG/KG CNTRLO7 TI1N 13 B7 1-06 OSDH 0.9 MG/KG CNTRLO8 TION 8 KB 1-06 OSDH 2.7 MG/KG CNTRLO9 TIIN 21 C12 1-06 OSDH 2.2 MG/KG CNTRLO9 TIIN 21 C12 1-06 OSDH 1.9 MG/KG CNTRLI3 TIIN 21 C12 1-06 OSDH 1.4 MG/KG CNTRLI3 TIIN 21 C12 1-06 OSDH 1.4 MG/KG CNTRLI3 TIIN 3 C2 I-06 OSDH 1.4 MG/KG CNTRLI5 TO9N 16 G13 1-06 OSDH 1.4 MG/KG GCN2 T12N 8 B3 1-15 OSDH 8.00 MG/KG G-01 T12N 7 K9									
CNTRLO7 T11N 13 B7 1-06 OSDH 0.9 MG/KG CNTRLO8 T10N 8 K8 1-05 OSDH 2.7 MG/KG CNTRLO9 T11N 21 C12 1-06 OSDH 2.2 MG/KG CNTRLO9 T11N 21 C12 1-06 OSDH 1.9 MG/KG CNTRL10 T11N 21 I13 1-06 OSDH 1.9 MG/KG CNTRL11 T12N 26 E3 1-06 OSDH 1.2 MG/KG CNTRL15 T09N 16 G13 1-06 OSDH 1.4 MG/KG GRES T12N 8 L5 1-15 OSDH 32.0 MG/KG G-02 T12N 8 B3 1-15 OSDH 2.08 MG/KG G-03 T12N 7 B13 1-15 OSDH 2.04 MG/KG G-01 T12N 8 N1 1-11 OSDH 2.06 MG/KG G-03									
CNTRLOB TION 8 K8 1-06 OSDH 2.7 MG/KG CNTRLO9 T11N 21 C12 1-06 OSDH 2.2 MG/KG CNTRLO9 T11N 21 C12 1-06 OSDH 1.9 MG/KG CNTRL10 T11N 21 C12 1-06 OSDH 3.0 MG/KG CNTRL11 T12N 26 E3 1-06 OSDH 1.2 MG/KG CNTRL13 T11N 3 C2 1-06 OSDH 1.4 MG/KG CNTRL13 T12N 26 E3 1-05 OSDH 1.4 MG/KG GCNTRL13 T12N 27 K3 1-05 OSDH 4.00 MG/KG G-01 T12N 8 B3 1-15 OSDH 2.08 MG/KG G-04 T12N 7 B13 1-15 OSDH 2.06 MG/KG HC-02 T11N 23 N1									
CNTRL09 T11N 21 C12 1-06 OSDH 2.2 MG/KG CNTRL09 T11N 21 C12 1-06 OSDH 1.9 MG/KG CNTRL10 T11N 21 T13 1-06 OSDH 1.9 MG/KG CNTRL13 T1N 3 C2 1-06 OSDH 1.2 MG/KG CNTRL15 T09N 16 G13 1-06 OSDH 1.4 MG/KG CRES T12N 27 K3 1-05 OSDH 32.0 MG/KG GRASS G-01 T12N 8 L5 1-15 OSDH 32.0 MG/KG GRASS G-01 T12N 8 B3 1-15 OSDH 3.04 MG/KG G-01 T12N 7 K9 1-15 OSDH 2.08 MG/KG G-01 T12N 16 F1 1-15 OSDH 2.04 MG/KG GL-01 T12N 7 </td <td>and the second second second second</td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td>	and the second second second second								
CNTRL09 T11N 21 C12 1-06 OSDH 1.9 MG/KG CNTRL10 T11N 21 I13 1-06 OSDH 3.0 MG/KG CNTRL11 T12N 26 E3 1-06 OSDH 1.2 MG/KG CNTRL13 T11N 3 C2 1-06 OSDH 1.4 MG/KG CNTRL15 T09N 16 G13 1-05 OSDH 1.4 MG/KG G-01 T12N 8 L5 1-15 OSDH 32.0 MG/KG G-02 T12N 8 B3 1-15 OSDH 8.00 MG/KG G-03 T12N 7 B13 1-15 OSDH 2.08 MG/KG G-04 T12N 7 K9 1-15 OSDH 2.04 MG/KG HC-01 T11N 23 N1 1-11 OSDH 2.48 MG/KG HC-03 T11N 20 N14 1-11 OSDH 2.00 MG/KG HC-07 T11N 21 A1									
CNTRL10 T11N 21 T13 1-06 OSDH 3.0 MG/KG CNTRL11 T12N 26 E3 1-06 OSDH 1.2 MG/KG CNTRL13 T11N 3 C2 1-06 OSDH 1.4 MG/KG CNTRL15 T09N 16 G13 1-06 OSDH 1.4 MG/KG CRES T12N Z7 K3 1-05 OSDH 32.0 MG/KG G-01 T12N 8 L5 1-15 OSDH 4.00 MG/KG G-02 T12N 8 B3 1-15 OSDH 2.08 MG/KG G-03 T12N 7 B13 1-15 OSDH 2.08 MG/KG GL-01 T12N 16 F1 1-15 OSDH 2.08 MG/KG HC-02 T11N 22 N1 1-11 OSDH 2.48 MG/KG HC-03 T11N 20 N14 1-11 <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td>a starte fille and</td> <td></td>								a starte fille and	
CNTRL11 T12N 26 E3 1-06 OSDH 1.2 MG/KG CNTRL13 T11N 3 C2 1-06 OSDH 1.4 MG/KG CNTRL15 T09N 16 G13 1-06 OSDH 1.4 MG/KG CRES T12N 27 K3 1-05 OSDH 32.0 MG/KG GRASS G-01 T12N 8 L5 1-15 OSDH 4.00 MG/KG G-03 T12N 7 B13 1-15 OSDH 2.08 MG/KG G-04 T12N 7 K9 1-15 OSDH 2.08 MG/KG GL-01 T12N 7 K9 1-15 OSDH 2.08 MG/KG GL-01 T12N 20 N1 1-11 OSDH 2.16 MG/KG HC-02 T11N 22 N1 1-11 OSDH 2.48 MG/KG HC-03 T11N 20 N14 1-11 OSDH 2.24 MG/KG HC-04 T11N 26									
CNTRL13 T11N 3 C2 1-06 OSDH 1.4 MG/KG CNTRL15 T09N 16 G13 1-06 OSDH 1.6 MG/KG GRASS GO1 T12N 8 L5 1-15 OSDH 32.0 MG/KG GRASS G-01 T12N 8 L5 1-15 OSDH 4.00 MG/KG GRASS G-02 T12N 8 B3 1-15 OSDH 8.00 MG/KG G-03 T12N 7 B13 1-15 OSDH 2.08 MG/KG G-04 T12N 7 K9 1-15 OSDH 2.04 MG/KG GL-01 T12N 16 F1 1-15 OSDH 2.16 MG/KG HC-02 T11N 23 N1 1-11 OSDH 2.00 MG/KG HC-03 T11N 20 N14 1-11 OSDH 2.00 MG/KG HC-04 T11N 21 A1 1-11 OSDH 2.24 MG/KG HC-07 T11N 2									
CNTRL15 TO9N 16 G13 1-06 OSDH 1.6 MG/KG CRES T12N 27 K3 1-05 OSDH 32.0 MG/KG GRASS G-01 T12N 8 L5 1-15 OSDH 4.00 MG/KG GRASS G-02 T12N 8 B3 1-15 OSDH 2.08 MG/KG G-03 T12N 7 B13 1-15 OSDH 2.08 MG/KG G-04 T12N 7 K9 1-15 OSDH 2.04 MG/KG GL-01 T12N 16 F1 1-15 OSDH 2.16 MG/KG HC-01 T11N 23 N1 1-11 OSDH 2.48 MG/KG HC-03 T11N 20 N14 1-11 OSDH 2.24 MG/KG HC-03 T11N 21 A1 1-11 OSDH 2.24 MG/KG HC-03 T11N 26									
CRES T12N 27 K3 1-05 OSDH 32.0 MG/KG GRASS G-01 T12N 8 L5 1-15 OSDH 4.00 MG/KG GRASS G-02 T12N 8 B3 1-15 OSDH 8.00 MG/KG G-03 T12N 7 B13 1-15 OSDH 2.08 MG/KG G-04 T12N 7 K9 1-15 OSDH 3.04 MG/KG GL-01 T12N 16 F1 1-15 OSDH 2.16 MG/KG HC-01 T11N 23 N1 1-11 OSDH 2.04 MG/KG HC-02 T11N 22 N1 1-11 OSDH 2.04 MG/KG HC-03 T11N 20 N14 1-11 OSDH 2.24 MG/KG HC-05 T11N 15 N5 1-11 OSDH 2.32 MG/KG HC-07 T1N 26 N1 1-11 OSDH 2.32 MG/KG HC-08 T1N									
G-01 T12N 8 L5 1-15 OSDH 4.00 MG/KG G-02 T12N 8 B3 1-15 OSDH 8.00 MG/KG G-03 T12N 7 B13 1-15 OSDH 2.08 MG/KG G-04 T12N 7 K9 1-15 OSDH 3.04 MG/KG GL-01 T12N 16 F1 1-15 OSDH 2.16 MG/KG HC-01 T11N 23 N1 1-11 OSDH 2.48 MG/KG HC-02 T11N 22 N1 1-11 OSDH 2.00 MG/KG HC-03 T11N 20 N14 1-11 OSDH 2.00 MG/KG HC-04 T1N 21 A1 1-11 OSDH 2.24 MG/KG HC-05 T11N 26 N1 1-11 OSDH 2.32 MG/KG HC-07 T1N 28 M11 1-11 OSDH 2.24 MG/KG HC-10 T10N 4 N1 <th< td=""><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td>GRASS</td></th<>									GRASS
G-02 T12N 8 B3 1-15 OSDH 8.00 MG/KG G-03 T12N 7 B13 1-15 OSDH 2.08 MG/KG G-04 T12N 7 K9 1-15 OSDH 3.04 MG/KG GL-01 T12N 16 F1 1-15 OSDH 2.16 MG/KG HC-01 T11N 23 N1 1-11 OSDH 2.48 MG/KG HC-02 T11N 22 N1 1-11 OSDH 2.48 MG/KG HC-03 T11N 20 N14 1-11 OSDH 2.00 MG/KG HC-04 T11N 21 A1 1-11 OSDH 2.00 MG/KG HC-05 T11N 15 N5 1-11 OSDH 2.24 MG/KG HC-07 T11N 26 N1 1-11 OSDH 2.32 MG/KG HC-09 T11N 32 A12 1-11 OSDH 2.32 MG/KG HC-10 T10N 4 N1									annoo
G-03 T12N 7 B13 1-15 OSDH 2.08 MG/KG G-04 T12N 7 K9 1-15 OSDH 3.04 MG/KG GL-01 T12N 16 F1 1-15 OSDH 2.16 MG/KG HC-01 T11N 23 N1 1-11 OSDH 2.48 MG/KG HC-02 T11N 22 N1 1-11 OSDH 2.48 MG/KG HC-03 T11N 20 N14 1-11 OSDH 2.00 MG/KG HC-04 T11N 21 A1 1-11 OSDH 2.24 MG/KG HC-05 T11N 15 N5 1-11 OSDH 2.200 MG/KG HC-07 T11N 26 N1 1-11 OSDH 2.24 MG/KG HC-08 T11N 28 M11 1-11 OSDH 2.32 MG/KG HC-10 T10N 4 N1 1-11 OSDH 2.24 MG/KG HC-11 T10N 4 N14			0						
G-04 T12N 7 K9 1-15 OSDH 3.04 MG/KG GL-01 T12N 16 F1 1-15 OSDH 2.16 MG/KG HC-01 T11N 23 N1 1-11 OSDH 2.48 MG/KG HC-02 T11N 22 N1 1-11 OSDH 1.04 MG/KG HC-03 T11N 20 N14 1-11 OSDH 2.00 MG/KG HC-04 T11N 21 A1 1-11 OSDH 2.24 MG/KG HC-05 T11N 15 N5 1-11 OSDH 2.24 MG/KG HC-07 T11N 26 N1 1-11 OSDH 2.32 MG/KG HC-08 T11N 28 M11 1-11 OSDH 2.32 MG/KG HC-10 T10N 4 N1 1-11 OSDH 2.24 MG/KG HC-11 T10N 4 N1 1-11 OSDH 2.52 MG/KG HC-11 T10N 4 N14			0						
GL-01 T12N 16 F1 1-15 OSDH 2.16 MG/KG HC-01 T11N 23 N1 1-11 OSDH 2.48 MG/KG HC-02 T11N 22 N1 1-11 OSDH 1.04 MG/KG HC-03 T11N 20 N14 1-11 OSDH 2.00 MG/KG HC-04 T11N 21 A1 1-11 OSDH 2.24 MG/KG HC-05 T11N 15 N5 1-11 OSDH 2.24 MG/KG HC-07 T11N 26 N1 1-11 OSDH 2.00 MG/KG HC-08 T11N 28 M11 1-11 OSDH 2.32 MG/KG HC-10 T10N 4 N1 1-11 OSDH 2.24 MG/KG HC-11 T10N 4 N1 1-11 OSDH 2.32 MG/KG HC-11 T10N 4 N14 1-11 OSDH 1.52 MG/KG HC-11 T10N 4 N14									
HC-01 T11N 23 N1 1-11 OSDH 2.48 MG/KG HC-02 T11N 22 N1 1-11 OSDH 1.04 MG/KG HC-03 T11N 20 N14 1-11 OSDH 2.00 MG/KG HC-04 T11N 21 A1 1-11 OSDH 1.28 MG/KG HC-05 T11N 15 N5 1-11 OSDH 2.24 MG/KG HC-07 T11N 26 N1 1-11 OSDH 2.00 MG/KG HC-07 T1N 26 N1 1-11 OSDH 2.00 MG/KG HC-07 T1N 28 M11 1-11 OSDH 2.00 MG/KG HC-09 T1N 32 A12 1-11 OSDH 2.32 MG/KG HC-10 T10N 4 N1 1-11 OSDH 1.52 MG/KG HC-11 T10N 4 N14 1-11 OSDH 1.60 MG/KG HC-12 T11N 22 A13									
HC-02 T11N 22 N1 1-11 OSDH 1.04 MG/KG HC-03 T11N 20 N14 1-11 OSDH 2.00 MG/KG HC-04 T11N 21 A1 1-11 OSDH 1.28 MG/KG HC-05 T11N 15 N5 1-11 OSDH 2.24 MG/KG HC-07 T11N 26 N1 1-11 OSDH 2.32 MG/KG HC-08 T11N 28 M11 1-11 OSDH 2.32 MG/KG HC-09 T1N 32 A12 1-11 OSDH 1.52 MG/KG HC-10 T10N 4 N1 1-11 OSDH 2.24 MG/KG HC-11 T10N 4 N1 1-11 OSDH 1.52 MG/KG HC-11 T10N 4 N14 1-11 OSDH 1.60 MG/KG HC-11 T10N 4 N14 1-11 OSDH 1.52 MG/KG HC-12 T11N 22 A13									
HC-03 T11N 20 N14 1-11 OSDH 2.00 MG/KG HC-04 T11N 21 A1 1-11 OSDH 1.28 MG/KG HC-05 T11N 15 N5 1-11 OSDH 2.24 MG/KG HC-07 T11N 26 N1 1-11 OSDH 2.00 MG/KG HC-08 T11N 28 M11 1-11 OSDH 2.32 MG/KG HC-09 T1N 32 A12 1-11 OSDH 1.52 MG/KG HC-10 T10N 4 N1 1-11 OSDH 2.24 MG/KG HC-11 T10N 4 N1 1-11 OSDH 2.32 MG/KG HC-11 T10N 4 N1 1-11 OSDH 1.52 MG/KG HC-11 T10N 4 N14 1-11 OSDH 1.60 MG/KG HC-12 T11N 22 A13 1-11 OSDH 0.90 MG/KG HC-13 T1N 26 A14									
HC-04 T11N 21 A1 1-11 OSDH 1.28 MG/KG HC-05 T11N 15 N5 1-11 OSDH 2.24 MG/KG HC-07 T11N 26 N1 1-11 OSDH 2.00 MG/KG HC-08 T11N 28 M11 1-11 OSDH 2.32 MG/KG HC-09 T11N 32 A12 1-11 OSDH 1.52 MG/KG HC-10 T10N 4 N1 1-11 OSDH 2.24 MG/KG HC-11 T10N 4 N1 1-11 OSDH 1.52 MG/KG HC-11 T10N 4 N14 1-11 OSDH 1.60 MG/KG HC-12 T11N 22 A13 1-11 OSDH 1.52 MG/KG HC-13 T11N 26 A14 1-11 OSDH 0.90 MG/KG HC-13 T11N 26 A14 1-11 OSDH 1.28 MG/KG HC-15 T10N 3 A1									
HC-05 T11N 15 N5 1-11 OSDH 2.24 MG/KG HC-07 T11N 26 N1 1-11 OSDH 2.00 MG/KG HC-08 T11N 28 M11 1-11 OSDH 2.32 MG/KG HC-09 T11N 32 A12 1-11 OSDH 1.52 MG/KG HC-10 T10N 4 N1 1-11 OSDH 2.24 MG/KG HC-11 T10N 4 N1 1-11 OSDH 1.52 MG/KG HC-11 T10N 4 N14 1-11 OSDH 1.60 MG/KG HC-12 T11N 22 A13 1-11 OSDH 1.52 MG/KG HC-12 T11N 22 A13 1-11 OSDH 0.90 MG/KG HC-13 T11N 26 A14 1-11 OSDH 3.60 MG/KG HC-14 T12N 25 A14 1-11 OSDH 1.28 MG/KG HC-15 T10N 3 A1 <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td>									
HC-07 T11N 26 N1 1-11 OSDH 2.00 MG/KG HC-08 T11N 28 M11 1-11 OSDH 2.32 MG/KG HC-09 T11N 32 A12 1-11 OSDH 1.52 MG/KG HC-10 T10N 4 N1 1-11 OSDH 2.24 MG/KG HC-11 T10N 4 N14 1-11 OSDH 1.60 MG/KG HC-11 T10N 4 N14 1-11 OSDH 1.52 MG/KG HC-12 T11N 22 A13 1-11 OSDH 1.52 MG/KG HC-13 T11N 26 A14 1-11 OSDH 0.90 MG/KG HC-13 T11N 26 A14 1-11 OSDH 3.60 MG/KG HC-14 T12N 25 A14 1-11 OSDH 1.28 MG/KG HC-15 T10N 3 A1 1-11 OSDH 4.08 MG/KG HC-16 T11N 33 N2 <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td>									
HC-08 T11N 28 M11 1-11 OSDH 2.32 MG/KG HC-09 T11N 32 A12 1-11 OSDH 1.52 MG/KG HC-10 T10N 4 N1 1-11 OSDH 2.24 MG/KG HC-11 T10N 4 N14 1-11 OSDH 1.60 MG/KG HC-11 T10N 4 N14 1-11 OSDH 1.52 MG/KG HC-11 T10N 4 N14 1-11 OSDH 1.52 MG/KG HC-12 T11N 22 A13 1-11 OSDH 0.90 MG/KG HC-13 T11N 26 A14 1-11 OSDH 3.60 MG/KG HC-14 T12N 25 A14 1-11 OSDH 1.28 MG/KG HC-15 T10N 3 A1 1-11 OSDH 4.08 MG/KG HC-16 T11N 33 N2 1-11 OSDH 5.76 MG/KG									
HC-09 T11N 32 A12 1-11 OSDH 1.52 MG/KG HC-10 T10N 4 N1 1-11 OSDH 2.24 MG/KG HC-11 T10N 4 N14 1-11 OSDH 1.60 MG/KG HC-11 T10N 4 N14 1-11 OSDH 1.60 MG/KG HC-11 T10N 4 N14 1-11 OSDH 1.52 MG/KG HC-12 T11N 22 A13 1-11 OSDH 0.90 MG/KG HC-13 T11N 26 A14 1-11 OSDH 3.60 MG/KG HC-14 T12N 25 A14 1-11 OSDH 1.28 MG/KG HC-15 T10N 3 A1 1-11 OSDH 4.08 MG/KG HC-16 T11N 33 N2 1-11 OSDH 5.76 MG/KG	and the second second								
HC-10 T10N 4 N1 1-11 OSDH 2.24 MG/KG HC-11 T10N 4 N14 1-11 OSDH 1.60 MG/KG HC-11 T10N 4 N14 1-11 OSDH 1.52 MG/KG HC-11 T10N 4 N14 1-11 OSDH 1.52 MG/KG HC-12 T11N 22 A13 1-11 OSDH 0.90 MG/KG HC-13 T11N 26 A14 1-11 OSDH 3.60 MG/KG HC-14 T12N 25 A14 1-11 OSDH 1.28 MG/KG HC-15 T10N 3 A1 1-11 OSDH 4.08 MG/KG HC-16 T11N 33 N2 1-11 OSDH 5.76 MG/KG									
HC-11 TION 4 N14 1-11 OSDH 1.60 MG/KG HC-11 TION 4 N14 1-11 OSDH 1.52 MG/KG HC-12 TIIN 22 A13 1-11 OSDH 0.90 MG/KG HC-13 TIIN 26 A14 1-11 OSDH 3.60 MG/KG HC-14 T12N 25 A14 1-11 OSDH 1.28 MG/KG HC-15 TION 3 A1 1-11 OSDH 4.08 MG/KG HC-16 TIIN 33 N2 1-11 OSDH 5.76 MG/KG									
HC-11 TION 4 N14 1-11 OSDH 1.52 MG/KG HC-12 T11N 22 A13 1-11 OSDH 0.90 MG/KG HC-13 T11N 26 A14 1-11 OSDH 3.60 MG/KG HC-14 T12N 25 A14 1-11 OSDH 1.28 MG/KG HC-15 T10N 3 A1 1-11 OSDH 4.08 MG/KG HC-16 T11N 33 N2 1-11 OSDH 5.76 MG/KG									
HC-12 T11N 22 A13 1-11 OSDH 0.90 MG/KG HC-13 T11N 26 A14 1-11 OSDH 3.60 MG/KG HC-14 T12N 25 A14 1-11 OSDH 1.28 MG/KG HC-15 T10N 3 A1 1-11 OSDH 4.08 MG/KG HC-16 T11N 33 N2 1-11 OSDH 5.76 MG/KG									
HC-13 T11N 26 A14 1-11 OSDH 3.60 MG/KG HC-14 T12N 25 A14 1-11 OSDH 1.28 MG/KG HC-15 T10N 3 A1 1-11 OSDH 4.08 MG/KG HC-16 T11N 33 N2 1-11 OSDH 5.76 MG/KG	CALIFORNIA CARDONNA								
HC-14 T12N 25 A14 1-11 OSDH 1.28 MG/KG HC-15 T10N 3 A1 1-11 OSDH 4.08 MG/KG HC-16 T11N 33 N2 1-11 OSDH 5.76 MG/KG									
HC-15 TION 3 A1 1-11 OSDH 4.08 MG/KG HC-16 TIIN 33 N2 1-11 OSDH 5.76 MG/KG									
HC-16 T11N 33 N2 1-11 OSDH 5.76 MG/KG	HC-14								
	HC-15	TION	3						
	HC-16	T11N	33	N2	1-11				
	HC-17	TION	9	A14	1-11	OSDH	2.80	MG/KG	

LABEL	TOWNSHIP	SECTION	SUBSECT	DATE	LAB	RESULTS	UNITS	COMMENTS	
HC-18	TION	10	A13	1-11	OSDH	3.60	MG/KG		
HC-19	T11N	29	A14	1-11		6.40	MG/KG		
HC-20	TIIN	30		1-11	OSDH	2.88	MG/KG		
P-01	T12N	26	A11	1-13	OSDH	2.24	MG/KG		
P-02	T12N	25	H4	1-13	OSDH	<0.1	MG/KG		
P-03	T12N	36	A13	1-13	OSDH	1.36	MG/KG		
P-04	T12N	36	D10	1-13	OSDH	<0.1			
P-05	TIIN	1	A11	1-13	OSDH		MG/KG		
P-06	TIIN	i	A7			1.76	MG/KG		
P-07		35	M10	1-13	OSDH	2.24	MG/KG		
P-08	T12N T12N	35		1-13	OSDH	9.68	MG/KG		
P-09			D10	1-13	OSDH	6.80	MG/KG		
	T12N	26	H6	1-13	OSDH	4.64	MG/KG		
S-01	TIIN	34	Al	1-11	OSDH	2.8	MG/KG		
S-02	TIIN	28	N1	1-11	OSDH		MG/KG		
S-03	TIIN	30	N14	1-11	OSDH		MG/KG		
S-04	T11N	25	N14	1-11	OSDH		MG/KG		
S-05	T11N	36	Al	1-11			MG/KG		
S-06	T11N	34	N14	1-11		3.12	MG/KG		
S-07	TION	2	Al	1-11			MG/KG		
S-08	TIIN	31	NI	1-11		5.44	MG/KG		
S-09	TION	6	A14	1-11		2.16	MG/KG		
S-10	TION	5	NI	1-11	OSDH	5.68	MG/KG		
S-11	TION	9	Al	1-11	OSDH	5.04	MG/KG		
S-12	TION	4	Al	1-11	OSDH	2.16	MG/KG		
S-13	T11N	33	N14	1-11	OSDH	4.96	MG/KG		
S-14	T10N	6	NI	1-11	OSDH	3.59	MG/KG		
S-15	TION	2	N14	1-11	OSDH	2.88	MG/KG		
S-16	TION	2	NI	1-11	OSDH		MG/KG		
T3-08	T12N	26	C12	1-14	OSDH		MG/KG		
T4-01E	T11N	2	C11	1-06	OSDH	2.6	MG/KG	EVERGREEN	
T4-02E	T11N	2	B8	1-06			MG/KG	EVERGREEN	
T4-03E	T12N	35	N5	1-06			MG/KG	EVERGREEN	
T4-04E	T12N	35	L3	1-06			MG/KG	EVERGREEN	
T4-04W	T12N	34	L13	1-06			MG/KG	LEAVES	
T4-05E	T12N	35	K1		OSDH			EVERGREEN	
WR-01	T11N	15	Al	1-15	OSDH	3.36	MG/KG	CYCHUNEEN	
WR-02	TIIN	17	A14	1-15	OSDH	5.36	MG/KG		
WR-03	TIIN	17	I14	1-15	OSDH	3.60	MG/KG		
WR-05	T11N	9	Al	1-15	OSDH	3.28	MG/KG		
WR-07	TIIN	6	N14	1-15	OSDH	1.68	MG/KG		
WR-08	TIIN	6	N7	1-15	OSDH	2.08			
WR-09	TIIN	6	Gl	1-15	OSDH	1.84	MG/KG		
WR-10	TIIN	17	Al	1-15	OSDH		MG/KG		
111 10	1 7 7 14	11	UT	1 10	USUN	2.32	MG/KG		

Table 5.2.7.11A

LABEL	TOWNSHIP	SECTION	SUBSECT	DATE	LAB	RESULTS	UNITS	COMMENTS
1H14	T11N	1	H14	1-10	КМТС	60.0	PPM	
2E14	TIIN	2	E14	1-10	KMTC	100.0	PPM	
4A 6	TIIN	4	A6	1-10	KMTC	170.0	PPM	
4J14	TIIN	4	J14	1-09	KMTC		PPM	
5N 6	TIIN	5	N6	1-09	KMTC	160.0	PPM	
6A 7	T11N		A7	1-09	KMTC	180.0	PPM	
6C14	TIIN	6	C14	1-09	KMTC	190.0	PPM	
6I11	TIIN		C14 I11 C6	1-09	KMTC	130 0		
7C 6	TIIN	7	C6	1-08	KMTC	130.0	PPM	
7F 7		7	F7	1-08	KMTC	120.0	PPM	
8A 1	TIIN	8	A1	1-09	KMTC	140.0		
8J 7	TIIN		J7	1-09	KMTC	50.0	PPM	
8M14	TIIN	8	A1 J7 M14	1-09	KMTC	30.0	PPM	
9B 2	T11N	9	B2	1-09	KMTC	250.0	PPM	
9I 8		9	18	1-09	KMTC	90.0	PPM	
	TIIN		B8	1-09	KMTC	30.0	PPM	
	TIIN			1-08		30.0	PPM	
	TIIN		M10	1-08		30.0	PPM	
	TIIN		F7	1-08		50.0	PPM	
	TIIN			1-08		50.0	PPM	
	TIIN			1-09		50.0	PPM	
	TIIN			1-08			PPM	
	TIIN		B1	1-08			PPM	
13G11	TIIN	13					PPM	
	TIIN		H2	1-08		60.0	PPM	
14B 4	TIIN	14					PPM	
14B 8	TIIN	14	B8				PPM	
	TIIN		B12				PPM	
	TIIN	16	B10				PPM	
	TIIN	16	L2				PPM	
	TIIN	16					PPM	
	TIIN	18	B4				PPM	
18H10	TIIN	18	H10					
1A11	TIIN	1	A11					
22G14			G14			70.0	PPM	
23A10	TIIN	23	A10	1-08	KMTC	50.0	PPM	T11N
25E 8	T12N	25	E8	1-10	KMTC	60.0	PPM	
25H 3	T12N	25	H3	1-10	KMTC	60.0	PPM	
25H14	T12N	25	H14	1-10	KMTC	80.0	PPM	
26G 1	T12N	26	G1	1-10	KMTC	60.0	PPM	
26G 8	T12N	26	G8	1-10	KMTC	70.0	PPM	
27A13	T12N	27	A13	1-08	KMTC	90.0	PPM	
27D 1	T12N	27	D1	1-10	KMTC	90.0	PPM	
27H 9	T12N	27	H9	1-10	KMTC	180.0	PPM	
27K 4	T12N	27	K4	1-10	KMTC	50.0	PPM	
28G 4	T12N	28	G4	1-10	KMTC	100.0	PPM	
28G 4 28J11	T12N	28	J11	1-10	KMTC	140.0	PPM	
28511 28N 2	T12N	28	N2	1-10	KMTC	170.0	PPM	
29D 3	T12N	29	D3	1-10	KMTC	50.0	PPM	
250 5	1 7714	6.5	00					

LABEL	TOWNSHIP	SECTION	SUBSECT	DATE	LAB	RESULTS	UNITS	COMMENTS
2L13	TIIN	2	L13	1-09	КМТС	50.	PPM	
2M 8	TIIN		M8	1-09		50.0	PPM	
30F 7	T12N		F7	1-10		60.0	PPM	
31B10	T12N		B10	1-10	KMTC	80.0	PPM	
	T12N			1-10	KMTC	170.0	PPM	
	T12N		13	1-10	KMTC			
33M14		33	M14	1-10		90.0	PPM	
	T12N			1-10	KMTC	90.0	PPM	
	T12N			1-10	KMTC	120.0	PPM	
	TI2N			1-10	KMIC	60.0	PPM	
	T12N		DI	1-10	KMTC	70.0	PPM	
	T12N		B14	1-10	KMTC	70.0	PPM	
35FI0	T12N		F10	1-10	KMTC			
351 3	T12N	35	L5	1-09	KMTC	70.0	PPM	
35111	TION	35	L11	1-10	KMTC	210.0	PPM	
	T12N		A13	1-10 1-09	KMTC	100.0	PPM	
	T12N		G5	1-09	KMTC	70 0	PPM	
36H14	T12N	36	H14	1-10	KMTC	160.0	PPM	
35 /	T11N T11N	3	F7	1-10 1-09 1-10 1-09 1-06 1-04 1-04	KMTC	90.0	PPM	
4A 1	TIIN	4	Al	1-10	KMTC	140.0	PPM	
4N /	T11N	4	N7	1-09	KMTC	90.0	PPM	
BG 1/4	T12N	21	B13	1-06	KMTC	220.0	PPM	
CRES10	T12N	27	K3	1-04	KMTC	190.0	PPM	CRES
DRES11	T12N	27	K4	1-04	KMTC	40.0	PPM	DRES
DRES13	T12N	27	14.7			60.0	PPM	DRES
HRES01	T12N	27	12	1-04	KMTC	110.0	PPM	HRES
LRES05	T12N	28	I14	1-04	KMTC	120.0	PPM	LRES
LRES06	T12N	28	I14	1-04	KMTC	60.0	PPM	LRES
RKRES16	T12N	34	A13	1-04	KMTC	80.0		RKRES
RRES17	TIIN	34	B6	1-04	KMTC	90.0		RRES
WALCIG	111N	34	B6	1-04	KMTC		PPM	
SHRES14	T12N T12N	27	N6	1-04	KMTC			SHRES
SMRES20	T12N	26	K6	1-04	KMTC			SMRES
TRES54	T11N	23	M10	1-06	KMTC			TRES
YRES52	T11N	26	N3	1-05	KMTC			
BG 1/4	T12N	21			ORNL	4.7	UG/G	DRY
BG 1/4	T12N	21	B13	1-06	ORNL	3.7	UG/G	WET
CARLILE		22	D14	1-05	ORNL	1.2	UG/G	DRY
CARLILE		22	D14	1-05	ORNL	0.8	UG/G	WET
CRES	T12N	27	K3	1-05	ORNL	2.8	UG/G	DRY
CRES	T12N	27	K3	1-05	ORNL	1.9	UG/G	WET
T4-01E	TIIN	2	C11	1-06	ORNL	1.2	UG/G	DRY
	T11N	2	C11	1-06	ORNL	1.0	UG/G	WET
T4-01W	T12N	27	N1	1-06	ORNL	0.7		
T4-01W	T12N	27	NI	1-06	ORNL	0.6	UG/G	DRY
T4-02E	TIIN	2	B8	1-06	ORNL		UG/G	WET
T4-02E	TIIN	2	88	1-06	ORNL	0.6	UG/G	DRY
T4-02W	T12N	33	L14	1-06		0.5	UG/G	WET
T4-02W	T12N	33	L14	1-06	ORNL	0.7	UG/G	DRY
	1 4.6.11	55	C 1.4	1.00	ORNL	0.6	UG/G	WET

LABEL	TOWNSHIP	SECTION	SUBSECT	DATE	LAB	RESULTS	UNITS	COMME	NTS
T4-03E	T12N	35	N5	1-06	ORNL	1.3	UG/G	DRY	
T4-03W	T12N	34	86	1-06	ORNL	0.8	UG/G	DRY	
T4-03W	T12N	34	B6	1-06	ORNL	0.7	UG/G	WET	
T4-04E	T12N	35	L3	1-06	ORNL		UG/G	DRY	
T4-04E	T12N	35	L3	1-06	ORNL		UG/G	WET	
T4-04W	T12N	34	L13	1-06	ORNL	1.0	UG/G	DRY	
T4-04W	T12N	34	L13	1-06	ORNL	0.8	UG/G	WET	
T4-05E	T12N	35	K1	1-06	ORNL		UG/G	DRY	
T4-05E	T12N	35	K1	1-06	ORNL		UG/G	WET	
AB-01	TIIN	17	A14	1-15	OSDH		MG/KG	HEI	
AB-02	TIIN	16	A11	1-15	OSDH		MG/KG		
AB-03	T12N	32	N14	1-15	OSDH	0.44	MG/KG		
AB-03	TIIN	9	B1	1-15	OSDH	0.56	MG/KG		
AB-04 AB-05	TIIN	4	C8	1-15	OSDH		MG/KG		
							MG/KG		
AB-06	T12N	29	N14	1-15	OSDH				
AB-07	T12N	31	A8	1-15	OSDH		MG/KG		
AB-08	T12N	19	N8	1-15	OSDH		MG/KG		
	T12N	18	M6	1-15	OSDH		MG/KG		
	T12N	13	G14	1-15			MG/KG		
CARLILE		22	D14	1-14	OSDH		MG/KG		
CARLILE		22	D14	1-05	OSDH		MG/KG		
CNTRL01		16	5. U	1-06	OSDH	0.4	MG/KG		
CNTRL02		35	C14	1-06	OSDH		MG/KG		
CNTRL03		28	E10	1-06	OSDH		MG/KG		
CNTRL04		20	N10	1-06	OSDH	0.7	MG/KG		
CNTRL05	T11N	19	F13	1-06	OSDH	0.9	MG/KG		
CNTRL07	T11N	13	B7	1-06	OSDH	<0.4	MG/KG		
CNTRL08	TION	8	K8	1-06	OSDH	0.6	MG/KG		
CNTRL09	T11N	21	C12	1-06	OSDH	1.3	MG/KG		
CNTRL10	T11N	21	I13	1-06	OSDH	1.3	MG/KG		
CNTRL11	T12N	26	E3	1-06	OSDH	1.5	MG/KG		
CNTRL12		14	N6	1-06	OSDH	0.7	MG/KG		
CNTRL13		3	C2	1-06	OSDH	0.4	MG/KG		
CNTEL14		-	1	1-06	OSDH	0.5	MG/KG	S OF	SALLISAW
CNTRL15		16	G13	1-06	OSDH		MG/KG		
CRES	T12N	27	KS	1-05	OSDH	0.6	MG/KG		
G-01	T12N	8	L5	1-15	OSDH	0.72	MG/KG		
G-02	T12N	8	B3	1-15	OSDH	1.84	MG/KG		
G-03	T12N	7	B13	1-15	OSDH	0.64	MG/KG		
G-04	T12N	7	K9	1-15	OSDH	1.52	MG/KG		
GL-01	T12N	16	F1	1-15	OSDH	0.52	MG/KG		
			A11		OSDH	0.88	MG/KG		
P-01	T12N	26	H4	1-13	OSDH	0.50	MG/KG		
P-02	T12N	25		1-13		0.68	MG/KG		
P-03	T12N	36	A13	1-13	OSDH				
P-04	T12N	36	D10	1-13	OSDH	0.96	MG/KG		
P-05	T11N	1	A11	1-13	OSDH	5.56	MG/KG		
P-06	TIIN	1	A7	1-13	OSDH	0.52	MG/KG		
P-07	T12N	35	M10	1-13	OSDH	<0.1	MG/KG		
P-08	T12N	35	D10	1-13	OSDH	0.48	MG/KG		

LABEL	TOWNSHIP	SECTION	SUBSECT	DATE	LAB	RESULTS	UNITS	COMMENTS
P-09	T12N	26	H6	1-13	OSDH	1.16	MG/KG	
T3-08	T12N	26	C12	1-14	OSDH	1.24	MG/KG	
T4-01E	T11N	2	C11	1-06	OSDH	1.1	MG/KG	
T4-01W	T12N	27	N1	1-06	OSDH	1.1	MG/KG	
T4-02E	T11N	2	B8	1-06	OSDH	<0.8	MG/KG	
T4-02W	T12N	33	L14	1-06	OSDH	2.3	MG/KG	
T4-03E	T12N	35	N5	1-06	OSDH	<0.8	MG/KG	
T4-03W	T12N	34	B6	1-06	OSDH	1.6	MG/KG	
T4-04E	T12N	35	L3	1-06	OSDH	<0.8	MG/KG	
T4-04W	T12N	34	L13	1-06	OSDH	1.4	MG/KG	
T4-05E	T12N	35	K1	1-06	OSDH	1.2	MG/KG	
W8-01	T11N	15	A1	1-15	OSDH	1.40	MG/KG	
WR-02	T11N	17	A14	1-15	OSDH	1.64	MG/KG	
WR-03	T11N	17	I14	1-15	OSDH	1.12	MG/KG	
WR-05	T11N	9	A1	1-15	OSDH	1.28	MG/KG	
WR-07	T11N	6	N14	1-15	OSDH	<0.40	MG/KG	
WR-08	T11N	6	N7	1-15	OSDH	1.08	MG/KG	
WR-09	T11N	6	G1	1-15	OSDH	0.52	MG/KG	
WR-10	TIIN	17	A1	1-15	OSDH	<0.40	MG/KG	

Table 5.2.7.12A

LABEL	TOWNSHIP	SECTION	SUBSECT	DATE	LAB	RESULTS	UNITS	COMMENTS
7C 6	T11N	7	C6	1-08	КМТС	1.2	PPM	J. GRASS
7F 7	T11N	7	F7	1-08	KMTC	3.9	PPM	WHEAT
7F 7	T11N	7	F7	1-08	KMTC	3.5	PPM	
11F 7	T11N	11	F7	1-08	KMTC	9.3	PPM	CAT LEAVES
11F 7	TIIN	11	F7	1-08	KMTC	31.0	PPM	BEANS
11F 7	TIIN	11	F7	1-08	KMTC	18.0	PPM	WILLOW
12N 8	TIIN	12	N8	1-08	KMTC	5.9	PPM	CEDAR
12N 8	TIIN	12	N8	1-08	KMTC	3.7	PPM	CAT TIPS
12N 8	TIIN	12	N8	1-08	KMTC	2.6	PPM	onti Trito
13B 1	TIIN	13	B1	1-08	KMTC	2.2	PPM	
13G11	TIIN	13	G11	1-09	KMTC	1.9	PPM	
14B 4	TIIN	14	B4	1-08	KMTC	1.5	PPM	
14B 8	TIIN	14	B8	1-08	KMTC	0.6	PPM	
14B12	TIIN	14	B12	1-08	KMTC	8.1	PPM	FESCUE
14B12	TIIN	14	B12	1-08	KMTC	7.8	PPM	FESCUE TIP
16B10	TIIN	16	B10	1-09	KMTC	0.6	PPM	ILSCOL III
16L12	TIIN	16	L12	1-10	KMTC	1.1	PPM	
	TIIN			1-09	KMTC	0.7	PPM	
18B 4		18	B4		KMTC	<0.5	PPM	
18H10	T11N	18	H10	1 09		0.7	PPM	
1A11	TIIN	1	A11	1-09	KMTC		PPM	
1H14	TIIN	1	H14	1-08	KMTC	1.4		
22G14	TIIN	22	G14	1-10	KMTC	<0.5	PPM	7111
23A10	T11N	23	A10	1-08	KMTC	28.0	PPM	T11N
25E 8	T12N	25	E8	1-10	KMTC	1.2	PPM	
25H 3	T12N	25	H3	1-10	KMTC	<0.5	PPM	
25H 3	T12N	25	H3	1-10	KMTC	<0.5	PPM	
26G 1	T12N	26	G1	1-10	KMTC	<0.5	PPM	
26G 8	T12N	26	G8	1-10	KMTC	1.1	PPM	
27A13	T12N	27	A13	1-08	KMTC	4.5	PPM	SMILAX
27A13	T12N	27	A13	1-08	KMTC	5.7	PPM	CEDAR
27D 1	T12N	27	D1	1-10	KMTC	52.0	PPM	
27H 9	T12N	27	Н9	1-10	KMTC	2.3	PPM	
27H 9	T12N	27	H9	1-10	KMTC	1.1	PPM	
27K 4	T12N	27	K4	1-10	KMTC	38.0	PPM	
28G 4	T12N	28	G4	1-10	KMTC	1.2	PPM	
28G 4	T12N	28	G4	1-10	KMTC	<0.5	PPM	
28J11	T12N	28	J11	1-10	KMTC	1.7	PPM	
28N 2	T12N	28	N2	1-10	KMTC	<0.5	PPM	
28N 2	T12N	28	N2	1-10	KMTC	6.6	PPM	
29D 3	T12N	29	D3	1-10	KMTC	<0.5	PPM	
2E14	TIIN	2	E14	1-09	KMTC	<0.5	PPM	
2L13	TIIN	2	L13	1-09	KMTC	5.5	PPM	
2M 8	TIIN	2	M8	1-09	KMTC	55.0	PPM	
30F 7	T12N	30	F7	1-10	KMTC	<0.5	PPM	
30F 7	T12N	30	F7	1-10	KMTC	0.5	PPM	
31B10	T12N	31	B10	1-10	KMTC	<0.5	PPM	
33B14	T12N	33	B14	1-10	KMTC	<0.5	PPI'	
33I 3	T12N	33	13	1-10	KMTC	<0.5	F	
	T12N	33	M14	1-10	KMTC	<0.5	F.	
33M14	17.514	22	117.4	1 10	MITC	10.0		

LABEL	TOWNSHIP	SECTION	SUBSECT	DATE	LAB	RESULTS	UNITS	COMMENTS
33M14	TIEN	33	M14	1-10	КМТС	<0.5	PPM	
34D 7	T12N	34	D7	1-10		16.0	PPM	
35A 9	T12N	35	A9	1-10		<0.5	PPM	
35A 9	T12N	35	A9	1-10		0.8	PPM	
	T12N	35	B1	1-10		<0.5	PPM	
35B14	T12N	35	B14	1-10			PPM	
	T12N	35	F10	1-10			PPM	
	T12N	35	L5	1-09			PPM	
	T12N	35	L5	1-09		and the second	PPM	
	T12N	35	L11	1-10			PPM	
	T12N	35	L11	1-10		1.2	PPM	
	T12N	36	A13	1-10			PPM	
	T12N	36	G5	1-09		0.9	PPM	
	T12N	36		1-10		0.5	PPM	
	TIIN	3		1-09		0.6	PPM	
	TIIN	4		1-10			PPM	
4A 3	TIIN	4		1-10			PPM	
4A 6	TIIN	4		1-10			PPM	
4J14	TIIN	4		1-09		1.0	PPM	
4J14	T11N	4		1-09			PPM	
4N 7	TIIN	4	N7	1-09			PPM	
	T11N	4	N8	1-09			PPM	
50/HWY	T12N	22	G1	1-06				
5N 6	T11N	5	N6	1-09		0.5	PPM	
	T11N	6	A7	1-09			PPM	
6C14	T11N	6	C14	1-09		1.8	PPM	
8A 1	TIIN	8	Al	1-09		1.7	PPM	
8J 7	TIIN	8	J7	1-09	KMTC	1.1	PPM	
8N14	T11N	8	N14	1-09		<0.5	PPM	
9B 2	T11N	8	82	1-09	KMTC	1.7	PPM	
BG 1/4	T12N	21	B13	1-06	KMTC	64.0	PPM	
	T12N	27	K3	1-04	KMTC	12.0	PPM	CRES
DRES12	T12N	27	K4	1-04	KMTC	220.0	PPM	DRES
H-02	T12N	27	12	1-04	KMTC	32.0	PPM	CEDAR
H-03	T12N					0.95	PPM	HRES
LRES04	T12N	28	I14	1-04	KMTC	<0.4	PPM	LRES
RKRES15		34	A13	1-04	KMTC	0.6	PPM	RKRES
	T11N	34	86	1-04	KMTC	2.7	PPM	RRES
SMRES21	T12N	26	K6	1-04	KMTC	0.6	PPM	SMRES
TRES53	T11N	23	M10	1-06	KMTC	16.0	PPM	STILLS
YRES51	TIIN	26	N3	1-05	KMTC	8.1	PPM	
BG 1/4	T12N	21	B13	1-06	ORNL	3.91	PPM	
T4-01E	T11N	2	C11	1-06	ORNL	5.29	PPM	
T4-02E	TIIN	2	B8	1-06	ORNL	1.99	PPM	
T4-03E	T12N	35	N5	1-06	ORNL	1.0	PPM	
T4-04E	T12N	35	L3	1-06	ORNL	0.807	PPM	
T4-04W	T12N	34	L13	1-06	ORNL	3.68	PPM	
T4-05E	T12N	35	K1	1-06	ORNL	19.6	PPM	

Table 5.2.7.13A

FAR FIELD ENVIRONMENTAL SAMPLING RESULTS FOR URANIUM IN SOIL SAMPLES

LABEL	TOWNSHIP	SECTION	SUBSECT	DATE	LAB	RESULTS	UNITS	COMMENTS
1H14	T11N	1	H14	1-10	KMTC	<10	PPM	
	TIIN	2	E14	1-10	KMTC	<10	PPM	
2E14				1-10	KMTC	<10	PPM	
4A 6	TIIN	4	A6	1-09	KMTC	<10	PPM	
4J14	TIIN	4	J14			<10	PPM	
5N 6	TIIN	5 6	N6	1-09	KMTC		PPM	
6A 7	TIIN	0	A7	1-09		<10	PPM	
6014	TIIN	6	C14	1-09		<10	PPM	
6111	TIIN	6 7	I11	1-09		<10	PPM	
70 6	T11N		C6	1-08		<10	PPM	
7F 7	T11N	7	F7	1-08		<10		
8A 1	TIIN	8	Al	1-09		<10	PPM	
8J 7	TIIN	8	J7	1-09		<10	PPM	
8N14	TIIN	8	M14	1-09	KMTC	20.0	PPM	
9B 2	T11N	9	B2	1-09	KMTC	<10	PPM	
9I 8	TIIN	9	18	1-09	KMTC	15.0	PPM	
10B 8	TIIN	10	B8	1-09	KMTC	<10	PPM	
10L 8	TIIN	10	L8	1-08	KMTC	<10	PPM	
10M10	TIIN	10	M10	1-08	KMTC	<10	PPM	
11F 7	TIIN	11	F7	1-08	KMTC	<10	PPM	
11N 6	TIIN	11	N6	1-08	KMTC	<10	PPM	
12F 7	TIIN	12	F7	1-09	KMTC	<10	PPM	
12N 8	TIIN	12	N8	1-08	KMTC	<10	PPM	
13B 1	TIIN	13	81	1-08	KMTC	15.0	PPM	
13G11	TIIN	13	G11	1-09	KMTC	15.0	PPM	
13H 2	TIIN	13	H2	1-08	KMTC	15.0	PPM	
14B 4	T11N	14	B4	1-08	KMTC	<10	PPM	
14B 8	TIIN	14	B8	1-08	KMTC	15.0	PPM	
14812	TIIN	14	B12	1-08	KMTC	15.0	PPM	
16810	TIIN	16	B10	1-09	KMTC	<10	PPM	
16L 2	T11N	16	L2	1-10	KMTC	<10	PPM	
16L12	TIIN	16	L12	1-10	KMTC	<10	PPM	
18B 4	TIIN	18	B4	1-09		<10	PPM	
18H10	T11N	18	H10	1-09		<10	PPM	
1A11	TIIN	1	A11	1-09	KMTC	<10	PPM	
22G14	T11N	22	G14	1-10	KMTC	<10	PPM	
23A10	TIIN	23	A10	1-08	KMTC	10.0	PPM	TIIN
25E 8	T12N	25	E8	1-10	KMTC	15.0	PPM	
25H 3	T12N	25	H3	1-10	KMTC	15.0	PPM	
25H14	T12N	25	H14	1-10	KMTC	15.0	PPM	
26G 1	T12N	26	G1	1-10	KMTC	15.0	PPM	
26G 8	T12N	26	G8	1-10	KMTC	10.0	PPM	
27A13	T12N	27	A13	1-08	KMTC	10.0	PPM	
27D 1	T12N	27	D1	1-10	KMTC	15.0	PPM	
27H 9	T12N	27	H9	1-10	KMTC	<10	PPM	
27K 4	T12N	27	K4	1-10	KMTC	15.0	PPM	
28G 4	T12N	28	G4	1-10	KMTC	<10	PPM	
28J11	T12N	28	J11	1-10	KMTC	<10	PPM	
28N 2	T12N	28	N2	1-10	KMTC	<10	PPM	
29D 3	T12N	29	D3	1-10	KMTC	15.0	PPM	

LABEL	TOWNSHIP	SECTION	SUBSECT	DATE	LAB	RESULTS	UNITS	COMMENTS
	TIIN	2	L13	1-09	KMTC	<10	PPM	
	TIIN	2	L13	1-09	KMTC	<10	PPM	
	T11N	2	M8	1-09	KMTC	<10	PPM	
	T12N	30	F7	1-10	KMTC	15.0	PPM	
	T12N	31	B10	1-10	KMTC	<10	PPM	
	T12N	33	B14	1-10	KMTC	<10	PPM	
	T12N		13	1-10	KMTC	<10	PPM	
	T12N	33		1-10	KMTC	<10	PPM	
	T12N	34	D7	1-10	KMTC	<10	PPM	
	T12N	35	A9	1-10		15.0	PPM	
	T12N	35	B1	1-10	KMTC	10.0	PPM	
	T12N	35	B14	1-10	KMTC	10.0	PPM	
	T12N	35	F10	1-10	KMTC	15.0	PPM	
	T12N	35	L5	1-09	KMTC	10.0	PPM	
	T12N	35	L11	1-10	KMTC	<10	PPM	
	T12N	36	A13	1-10	KMTC	<10	PPM	
	T12N	36	G5	1-09	KMTC	<10	PPM	
	T12N	36	H14	1-10	KMTC	<10	PPM	
	TIIN	3	F7	1-09	KMTC	<10	PPM	
4A 1	T11N	4	A1	1-10	KMTC	<10	PPM	
4N 7	T11N	4	N7	1-09	KMTC	<10	PPM	
BG 1/4	T12N	21	B13	1-06	KMTC	<10	PPM	
CRES09	T12N	27	K3	1-04		<10	PPM	CRES
CRES10	T12N	27	K3	1-04	KMTC	<10	PPM	CRES
DRES11	T12N	27	K4	1-04		15.0	PPM	DRES
DRES13	T12N	27	K4	1-04		13.0	PPM	DRES
HRES01	T12N	27	12	1-04		<10	PPM	HRES
LRES05	T12N	28	I14	1-04		<10	PPM	LRES
LRES06	T12N	28	I14	1-04		<10	PPM	LRES
RKRES16	T12N	34	A13	1-04		<10	PPM	RKRES
RRES17	TIIN	34	B6	1-04		<10	PPM	RRES
RRES19	T11N	34	B6	1.04	KMTC	<10	PPM	RRES
SHRES14	T12N	27		1-04			PPM	SHRES
SMRES20	T12N	26	К6	1-04	KMTC			SMRES
	TIIN	23	M10	1-06		<10	PPM	
YRES52	TIIN	26	N3	1-05		<10	PPM	
BG 1/4	T12N	21	B13	1-06	ORNL	7.88	PPM	DRY
BG 1/4	T12N	21	B13	1-06		6.1	PPM	WET
CARLILE	T12N	22	D14	1-05		3.7	PPM	DRY
CARLILE	T12N	22	D14	1-05		2.7	PPM	WET
CRES	T12N	27	K3	1-05	ORNL	2.9	PPM	DRY
CRES	T12N	27	K3	1-05		2.0	PPM	WET
T4-01E	T11N	2	C11	1-06	ORNL	3.72	PPM	DRY
T4-01E	TIIN	2	C11	1-06	ORNL	3.1	PPM	WET
T4-01W	T12N	27	N1	1-06	ORNL	1.52	PPM	DRY
	T12N	27	N1	1-06	ORNL	1.3	PPM	WET
	TIIN	2	B8	1-06	ORNL	3.59	PPM	DRY
	TIIN	2	88	1-06		3.0	PPM	WET
	T12N	33	L14	1-06	ORNL	1.5	PPM	WET
				1 00	UNITE	4.0		HL I

LABEL	TOWNSHIP	SECTION	SUBSECT	DATE	LAB	RESULTS	UNITS	COMMENTS
T4-02W	T12N	33	L14	1-06	ORNL	1.67	PPM	DRY
T4-03E	T12N	35	N5	1-06	ORNL	3.8	PPM	DRY
T4-03E	T12N	35	N5	1-06	ORNL	2.9	PPM	WET
T4-03W	T12N	34	B6	1-06	OKNL	2.76	PPM	DRY
T4-03W	T12N	34	B6	1-06	ORNL	2.2	PPM	WET
T4-04E	T12N	35	L3	1-06	ORNL	3.73	PPM	DRY
T4-04E	T12N	35	L3	1-06	ORNL	3.0	PPM	WET
T4-04W	T12N	34	L13	1-06	ORNL	2.45	PPM	DRY
T4-04W	T12N	34	L13	1-06	ORNL	2.0	PPM	WET
T4-05E	T12N	35	K1	1-06	ORNL	2.29	PPM	DRY
T4-05E	T12N	35	K1	1-06	ORNL	1.9	PPM	WET

Table 5.2.7.14A

LABEL	TOWNSHIP	SECTION	SUBSECT	DATE	LAB	RESULTS	UNITS	COMMENTS
CARLILE	T12N	22	D14	1-05	OSDH	0.8	PCI/G	GRASS
CNTRL01	TOON	16	-	1-06	OSDH	-0.4	PCI/G	
CNTRL02	TION	35	C14	1-06	OSDH	0.2	PCI/G	
CNTRL03	TIIN	28	E10	1-06	OSDH	0.1	PCI/G	
CNTRL04	T11N	20	N10	1-06	OSDH	0.5	PCI/G	
CNTRL05	TIIN	19	F13	1-06	OSDH	-0.3	PCI/G	
CNTRL06	T09N	4	I13	1-06	OSDH	-0.4	PCI/G	
CNTRL07	T11N	13	B7	1-06	OSDH	-0.3	PCI/G	
CNTRL08	TION	8	K8	1-06	OSDH	1.0	PCI/G	
CNTRL09	TIIN	21	C12	1-06	OSDH	1.0	PCI/G	
CNTRL09	TIIN	21	C12	1-06	OSDH	0.4	PCI/G	
CNTRL10	TIIN	21	I13	1-06	OSDH	0.1	PCI/G	
CNTRL11	T12N	26	E3	1-06	OSDH	0.2	PCI/G	
CNTRL13	TIIN	3	C2	1-06	OSDH	0.2	PCI/G	
CNTRL15	T09N	16	G13	1-06	OSDH	-0.3	PCI/G	
CRES	T12N	27	K3	1-05	OSDH	1.1	PCI/G	GRASS
HRES	T12N	27	12		OSDH	5.0	PCI/G	
T4-01E	T11N	2	C11	1-06	OSDH	0.3	PCI/G	EVERGREEN
T4-02E	TIIN	2	88	1-06	OSDH	-0.2	PCI/G	EVERGREEN
T4-03E	T12N	35	N5	1-06	OSDH	-0.1	PCI/G	EVERGREEN
T4-04E	T12N	35	L3	1-06	OSDH	0.0	PCI/G	EVERGREEN
T4-04W	T12N	34	L13	1-06	CSDH	0.1	PCI/G	LEAVES
T4-05E	T12N	35	K1	1-06	OSDH	4.5	PCI/G	EVERGREEN

Table 5.2.7.15A

1	LABEL	TOWNSHIP	SECTION	SUBSECT	DATE	LAB	RESULTS	UNITS	COMMENTS
1	CARLILE	T12N	22	D14	1-05	OSDH	0.7	PCI/G	GRASS
	CNTRL01		16	-	1-06	OSDH	-0.1	PCI/G	
	CNTRL02		35	C14	1-06	OSDH	2.0	PCI/G	
	CNTRL03		28	E10	1-06	OSDH	0.3	PCI/G	
	CNTRL04		20	N10	1-06	OSDH	0.3	PCI/G	
	CNTRL05	and the second sec	19	F13	1-06	OSDH	0.2	PCI/G	
	CNTRL06		4	I13	1-06	OSDH	0.6	PCI/G	
	CNTRL07		13	B7	1-06	OSDH	0.8	PCI/G	
	CNTRL08		8	K8	1-06	OSDH	0.7	PCI/G	
	CNTRL09		21	C12	1-06	OSDH	7.0	PCI/G	
1	CNTRL09	TIIN	21	C12	1-06	OSDH	-0.3	PCI/G	
(CNTRL10	T11N	21	I13	1-06	OSDH	-0.2	PCI/G	
(CNTRL11	T12N	26	E3	1-06	OSDH	0.1	PCI/G	
1	CNTRL13	TIIN	3	C2	1-06	OSDH	3.0	PCI/G	
1	CNTRL15	TOON	16	G13	1-06	OSDH	0.3	PCI/G	
1	CRES	T12N	27	K3	1-05	OSDH	0.2	PCI/G	GRASS
1	HRES	T12N	27	I2	1/5	OSDH	1.9	PCI/G	
	T4-01E	TIIN	2	C11	1-06	OSDH	-0.2	PCI/G	EVERGREEN
	T4-02E	T11N	2	B8	1-06	OSDH	-0.3	PCI/G	EVERGREEN
2	T4-03E	T12N	35	N5	1-06	OSDH	0.2	PCI/G	EVERGREEN
1	T4-04E	T12N	35	L3	1-06	OSDH	-0.2	PCI/G	EVERGREEN
	T4-04W	T12N	34	L13	1-06	OSDH	0.4	PCI/G	LEAVES
1	T4-05E	T12N	35	K1	1-06	OSDH	0.8	PCI/G	EVERGREEN

Table 5.2.7.16A

LABEL	TOWNSHIP	SECTION	SUBSECT	DATE	LAB	RESULTS	UNITS	COMMENTS
CARLTIE	TION	22	014	1.05	0.000			
CARLILE		22	D14	1-05	OSDH	0.1	PCI/G	
CNTRL01		16	-	1-06	OSDH	0.8	PCI/G	
CNTRL02		35	C14	1-06	OSDH	0.4	PCI/G	
CNTRL03		28	E10	1-06	OSDH	-0.2	PCI/G	
CNTRL04		20	N10	1-06	OSDH	1.1	PCI/G	
CNTRL05		19	F13	1-06	OSDH	-0.3	PCI/G	
CNTRL07		13	87	1-06	OSDH	-0.3	PCI/G	
CNTRL08		8	K8	1-06	OSDH	0.4	PCI/G	
CNTRL09		21	C12	1-06	OSDH	0.1	PCI/G	
CNTRL10	T11N	21	I13	1-06	OSDH	-0.3	PCI/G	
CNTRL11	T12N	26	E3	1-06	OSDH	-0.3	PCI/G	
CNTRL12		14	N6	1-06	OSDH	0.5	PCI/G	
CNTRL13		3	C2	1-06	OSDH	0.4	PCI/G	
CNTRL14	TIIN		-	1-06	OSDH	0.7	PCI/G	S OF SALLISAW
CNTRL15	T09N	16	G13	1-06	OSDH	-0.3	PCI/G	
CRES	T12N	27	K3	1-05	OSDH	-0.1	PCI/G	
T4-01E	T11N	2	C11	1-06	OSDH	-0.3	PCI/G	
T4-01W	T12N	27	N1	1-06	OSDH	-0.2	PCI/G	
T4-02E	T11N	2	B8	1-06	OSDH	0.4	PCI/G	
T4-02W	T12N	33	L14	1-06	OSDH	0.7	PCI/G	
T4-03E	T12N	35	N5	1-06	OSDH	0.1	PCI/G	
T4-03W	T12N	34	B6	1-06	OSDH	0.1	PCI/G	
T4-04E	T12N	35	L3	1-06	OSDH	0.1	PCI/G	
T4-04W	T12N	34	L13	1-06	OSDH	0.1	PCI/G	
T4-05E	T12N	35	K1	1-06	OSDH	-0.3	PCI/G	

Table 5.2.7.17A

LABEL	TOWNSHIP	SECTION	SUBSECT	DATE	LAB	RESULTS	UNITS	COMMENTS
CARLILE	T12N	22	D14	1-05	OSDH	-0.3	PCI/G	
CNTRL01		16	-	1-06	OSDH	-0.3	PCI/G	
CNTRL02		35	C14	1-06	OSDH	0.3	PCI/G	
CNTRL03		28	E10	1-06	OSDH	-0.0	PCI/G	
CNTRL04		20	N10	1-06	OSDH	-0.3	PCI/G	
CNTRL05		19	F13	1-06	OSDH	-0.0	PCI/G	
CNTRL07		13	B7	1-06	OSDH	0.5	PCI/G	
CNTRL08		8	K8	1-06	OSDH	-0.4	PCI/G	
CNTRL09		21	C12	1-06	OSDH	-0.3	PCI/G	
CNTRL10		21	I13	1-06	OSDH	0.3	PCI/G	
CNTRL11		26	E3	1-06	OSDH	-0.2	PCI/G	
CNTRL12		14	N6	1-06	OSDH	0.9	PCI/G	
CNTRL13		3	C2	1-06	OSDH	0.5	PCI/G	
CNTRL14		-	-	1-06	OSDH	0.6	PCI/G	S OF SALLISAW
CNTRL15	T09N	16	G13	1-06	OSDH	-0.4	PC1/G	
CRES	T12N	27	K3	1-05	OSDH	-0.0	PCI/G	
T4-01E	T11N	2	C11	1-06	OSDH	0.1	PCI/G	
T4-01W	T12N	27	N1	1-06	OSDH	-0.2	PCI/G	
T4-02E	T11N	2	B8	1-06	OSDH	-0.2	PCI/G	
T4-02W	T12N	33	L14	1-06	OSDH	0.1	PCI/G	
T4-03E	T12N	35	N5	1-06	OSDH	0.6	PCI/G	
T4-03W	T12N	34	B6	1-06	OSDH	-0.6	PCI/G	
T4-04E	T12N	35	L3	1-06	OSDH	0.1	PCI/G	
T4-04W	T12N	34	L13	1-06	OSDH	-0.4	PCI/G	
T4-05E	T12N	35	K1	1-06	OSDH	-0.4	PCI/G	

NEAR FIELD ENVIRONMENTAL SAMPLING RESULTS BY TYPE OF SAMPLE SORTED BY LABEL

Table 5.2.7.18A

NEAR FIELD ENVIRONMENTAL SAMPLING RESILTS FOR FLUORIDE SAMPLES

LABEL	NORTH	EAST	LAB	DATE	MEDIA	RESULT	UNITS	COMMENTS
1W	14500	14800	KMTC	1-06	WATER	0.1	mg/1	
2W	13250	13950		1-06	WATER	0.1	mg/1	
3W	12750	11200		1-06	WATER	0.2	mg/1	
4W	11700	11150		1-06	WATER	0.2	mg/1	
5W	10100	10950		1-06	WATER	0.4	mg/1	
6W	9850		KMTC	1-06		0.1	mg/1	
7W	10600		KMTC	1-06	WATER	0.2	mg/1	
8W	11800		KMTC	1-06		0.3	mg/1	
10W	11800		KMTC	1-06		0.1	mg/1	
11W	16350		KMTC	1-06		<0.1	mg/1	
12W	6750	11050	the state of the s	1-06		0.1	mg/1	
13W	7475	13350		1-06		<0.1	mg/1	
14W	18450		KMTC	1-06		<0.1	mg/1	
14W	18450		KMTC	1-08		<0.1	mg/1	
15W	12900		KMTC	1-06		<0.1	mg/1	
15W	12900		KMTC	1-08				
17W	7650		KMTC	1-06			mg/1	
18W	15100	10300		1-08		0.1	mg/1	
19W	6850	9050		1-08		0.2	mg/1	
20W	7750	17000		1-08		0.1		
103	13800	10200		1-09	SOIL		ppm	
105	11700		KMTC	1-09	VEG	20.0	ppm	
106	11700		KMTC	1-09	SOIL		ppm	
108	10600		KMTC	1-09	VEG	30.0	ppm	
109	10600		KMTC	1-09			ppm	
111	14250		KMTC	1-09	VEG		ppm	
112	14250		KMTC	1-09	SOIL	110.0	ppm	
115	11500	11600		1-09	SOIL	100.0	ppm	
116	10650	12300		1-09	VEG	10.0	ppm	
117	10650	12300	KMTC	1-09	SOIL	80.0	ppm	
118	9750	12450	KMTC	1-09	VEG	14.0	ppm	
119	9750	12450	KMTC	1-09	SOIL	110.0	ppm	
120	9700		KMTC	1-09	VEG	20.0	ppm	
121	9700	8850	KMTC	1-09	SOIL	70.0	ppm	
122	9600	9900	KMTC	1-09	VEG	210.0	ppm	
123	9600		KMTC	1-09	SCIL	130.0	ppm	
124	9700	10700	KMTC	1-09	VEG	100.0	ppm	
125	9700	10700	KMTC	1-09	SOIL	60.0	ppm	
126	9750	11850	KMTC	1-09	VEG	20.0	ppm	
127	9750	11850	KMTC	1-09	SOIL	90.0	ppm	
128	10650	11300		1-09	VEG	35.0	ppm	
129	10650	11300		1-09	SOIL	80.0	ppm	
130	10650	10350		1-09	VEG	70.0	ppm	
131	10650	10350		1-09	SOIL	70.0	ppm	
132	11600	10650		1-09	VEG	80.0	ppm	
133	11600	10650		1-09	SOIL	110.0	ppm	
134	9250	13150		1-09	VEG	28.0	ppm	
135	9250	13150		1-09	SOIL	90.0	ppm	

LABEL	NORTH	EAST	LAB	DATE	MEDIA	RESULT	UNITS	COMMENTS
136	9150	12050	КМТС	1-09	VEG	16.0	ppm	
137	9150	12050	KMTC	1-09	SOIL	150.0	ppm	
138	9200	11000	KMTC	1-09	VEG	90.0	ppm	
139	9200	11000	KMTC	1-09	SOIL	120.0	ppm	
141	9300	10100	KMTC	1-09	SOIL	160.0	ppm	
142	9150	9000	KMTC	1-09	VEG	15.0	ppm	
143	9150	9000	KMTC	1-09	SOIL	80.0	ppm	
144	10500	9700	KMTC	1-09	VEG	50.0	ppm	
145	10500	9700	KMTC	1-09	SOIL	100.0	ppm	
146	11600	9550	KMTC	1-09	VEG	73.0	ppm	
147	11600	9550	KMTC	1-09	SOIL	30.0	ppm	
149	11650		KMTC	1-09	SOIL		ppm	
150	11450		KMTC	1-09	VEG		ppm	
187	12600		KMTC	1-14	SOIL	60.0	ppm	
190	14750		KMTC	1-14	SOIL		ppm	
193	12550	10200		1-14	SOIL		ppm	
196	12600	10750		1-14	SOIL		ppm	
198	12550	11850		1-14	VEG	I	ppm	
199	12550	11850		1-14			ppm	
201	13150	10200		1-14		I	ppm	
202	13150	10200		1-14		100.0	ppm	
204	13150	11250		1-14			ppm	
205	13150	11250		1-14			ppm	
208		9900		1-14			ppm	
209		9900		1-14			ppm	
211	12500		KMTC	1-14			ppm	
214	12500		KMTC	1-14			ppm	
218	12500		KMTC	1-14			ppm	
221	12500		KMTC	1-14				
300	12450		KMTC	1-06				
301	12450		KMTC	1-14		0.2	mg/1	
302	12550			1-06			mg/1	
303	12550			1-14	WATER			
304	12550			1-06	WATER			
305	12550			1-14		4.8	mg/1	
306	11450		KMTC	1-06	WATER		mg/1	
307	11450		KMTC	1-14	WATER	0.5		
50	10750		KMTC	1-06	VEG	1500.0	ppm	CEDAR
A-01	13130		KMTC	1-06	VEG	2100.0	ppm	
A-01	13130		KMTC	1-06	SOIL	250.0	ppm	
A-07	13130		KMTC	1-06	VEG	2600.0	ppm	
A-07	13130		KMTC	1-06	SOIL	410.0	ppm	
A-07	13130		KMTC	1-09	SOIL	410.0	ppm	
A-13	13130		KMTC	1-06	SOIL	320.0	ppm	
A-13	13130		KMTC	1-06	VEG	3000.0	ppm	
A-13	13130		KMTC	1-09	SOIL	140.0	ppm	
A-19	13130		KMTC	1-06	VEG	230.0	ppm	
A-25	13130		KMTC	1-06	VEG	35.0	ppm	
1 20	10100	5550	MILLE	1 00	VLU	55.0	ppm	

LABEL	NORTH	EAST	LAB	DATE	MEDIA	RESULT	UNITS	COMMENTS
A~25	13130	9950	KMTC	1-06	SOIL	110.0	ppm	
A-25	13130	9760		1-06	SOIL	9.0	ug/g	WET
A-25	13130	9760		1-06	SOIL	11.4	ug/g	DRY
A-25	13130	9950		1-06	VEG	33.0	ug/g	UNI
BG 1/4	14750	9500		1-06	VEG	7.3	ug/g	
BG 1/4	14750	9500		1-06	SOIL	3.7	ug/g	WET
BG 1/4	14750	9500		1-06	SOIL	4.7	ug/g	DRY
BG 1/4	14750	9500		1-06	VEG	40.0	ppm	1/4 N. GATE
BG 1/4	14750	9500		1-06	SOIL	220.0	ppm	4/4 H. UNIC
CARLILE	14250	15100		1-05	SOIL	1.2	ug/g	DRY
CARLILE	14250	15100		1-05	SOIL	0.8	ug/g	WET
CARLILE	14250	15100		1-05	VEG	8.3	ug/g	HE I
CARLILE	14250	15100		1-05	VEG	4.3	mg/kg	GRASS
CARLILE	14100	15300		1-14	VEG	2.56	mg/kg	
CARLILE	14100	15300		1-14	SOIL	1.24	mg/kg	
CARLILE	14250	15100		1-05	SOIL	<0.4	mg/kg	
COMP-1	11800	9500		1-16	VEG	488.0	mg/kg	LAB COMPARISON
COMP-1	11800	9500		1-16	SOIL	5.44	mg/kg	LAB COMPARISON
COMP-2	8800	10100		1-16	VEG	164.0	mg/kg	LAB COMPARISON
COMP-2	8800	10100		1-16	SOIL	2.76	mg/kg	LAB COMPARISON
COMP-3	10050	10700		1-16	VEG	128.0	mg/kg	LAB COMPARISON
COMP-3	10050	10700		1-16	SOIL	1.96	mg/kg	LAB COMPARISON
COMP-3	10050	10700		1-16	WATER	0.45	mg/1	LAB COMPARISON
COMP-4	14450	14800		1-16	VEG	7.36	mg/kg	LAB COMPARISON
COMP-4	14450	14800		1-16	SOIL	1.48	mg/kg	LAB COMPARISON
CCMP-4	14450	14800		1-16	WATER	0.27	mg/1	LAB COMPARISON
COMP-5	12700	11200		1-16	VEG	1080	mg/kg	LAB COMPARISON
COMP-5	12700	11200		1-16	SOIL	5.6	mg/kg	LAB COMPARISON
COMP-5	12700	11200		1-16	WATER	0.24	mg/1	LAB COMPARISON
CRES	6150	10750		1-05	VEG	32.0	mg/kg	GRASS
CRES	6150	10800		1-05	WATER	0.19	mg/1	
CRES	6150	10800		1-05	WATER	0.09	ug/ml	
CRES	6150	10750		1-05	SOIL	1.9	ug/g	WET
CRES	6150	10750		1-05	SOIL	2.8	ug/g	DRY
CRES	6150	10800		1-05	VEG	16.2	ug/g	
CRES	6150	10750		1-05	SOIL	0.6	mg/kg	
CRES08	6150	10750		1-03	VEG	42.0	ppm	CRES
CRES09	6150	10750		1-04	S-V	90.0	ppm	MIX SOIL/VEG
CRES10	6150	10750		1-04	SOIL	190.0	ppm	CRES
DRES11	6100	11100		1-04	SOIL	40.0	ppm	
DRES12	6100	11100		1-04	VEG	160.0	ppm	DRY LEAVES
DRES13	6100	11100		1-04	SOIL	60.0	ppm	
E-01	13000	9230		1-06	VEG	700.0	ppm	
E-01	13000	9230		1-06	SOIL	100.0	ppm	
E-01	13000	9230		1-09	SOIL	180.0	ppm	
E-07	13000	9410		1-06	VEG	3100.0	ppm	
E-07	13000	9410		1-06	SOIL	170.0	ppm	
	20000	0110	an u	2 00			p p at	

0

١

1

.

LABEL	NORTH	EAST LAB	DATE	MEDIA	RESULT	UNITS	COMMENTS
E-07	13000	9410 KMTC	1-09	SOIL	170.0	ppm	
E-13	13000	9600 KMTC	1-06	VEG	3700.0	ppm	
E-13	13000	9600 KMTC	1-06	SOIL	230.0	ppm	
E-13	13000	9600 KMTC	1-09	SOIL	150.0	ppm	
E-19	13000	9770 KMTC	1-06	SOIL	40.0	ppm	
E-19	13000	9770 KMTC	1-06	VEG	760.0	ppm	
E-25	13000	9950 KMTC	1-06	VEG	80.0	ppm	
E-25	13000	9950 KMTC	1-06	SOIL	170.0	ppm	
E-25	13000	9950 ORNL	1-06	SOIL	7.4	ug/g	WET
E-25	13000	9950 ORNL	1-06	SOIL	8.6	ug/g	DRY
E-25	13000	9950 ORNL	1-06	VEG	58.0	ug/g	
HRES	6950	10650 OSDH	1-08	VEG	24.0	mg/kg	
HRESOI	6950	10650 KMTC	1-04	SOIL	110.0	ppm	
HRES02	6950	10650 KMTC	1-04	VEG	24.0	ppm	
HRES02	6950	10650 KMTC	1-04	VEG	80.0	ppm	
J-07	12840	9410 KMTC	1-04	VEG	870.0	ppm	
	12840	9410 KMTC	1-06	SOIL	150.0	ppm	
J-07			1-09	SOIL	180.0	ppm	
J-07	12840			SOIL	220.0		
J-13	12840	9600 KMTC	1-06			ppm	
J-13	12840	9600 KMTC	1-06	VEG	1700.0 200.0	ppm	
J-13	12840	9600 KMTC	1-09	SOIL		ppm	
J-19	12840	9780 KMTC	1-06	VEG	1300.0	ppm	
J-19	12840		1-06	SOIL	90.0	ppm	
J-25	12840		1-06	VEG	140.0	ppm	
J-25	12840		1-06	SOIL	110.0	ppm	
J-25	12840	9950 ORNL	1-06	SOIL	8.3	ug/g	WET
J-25	12840	9950 ORNL	1-06	SOIL	11.2	ug/g	DRY
J-25	12840	9950 ORNL	1-06	VEG	150.0	ug/g	
KM H20	13250	9400 OSDH	1-17	WATER	0.24	mg/1	
LAGOON	12800	9200 OSDH	1-17	WATER	4.40	mg/1	
LAGOON	12800	9200 KMTC	1-17	WATER	3.9	mg/1	
LR-01	7300	10450 OSDH	1-11	SOIL	0.88	mg/kg	
LR-01	7300	10450 OSDH	1-11	VEG	56.8	mg/kg	
LR-01	7300	10450 OSDH	1-11	VEG	56.8	mg/kg	
LR-01	7300	10450 OSDH	1-11	SOIL	0.88	mg/kg	
LR-02	6000	10000 OSDH	1-11	SOIL	1.36	mg/kg	
LR-02	6000	10000 OSDH	1-11	VEG	2.96	mg/kg	
LR-02	6000	10000 OSDH	1-11	SOIL	1.36	mg/kg	
LR-02	6000	10000 OSDH	1-11	VEG	2.96	mg/kg	
LR-03	6000	10800 OSDH	1-11	SOIL	1.68	mg/kg	
LR-03	6000	10800 OSDH	1-11	VEG	56.00	mg/kg	
LR-03	6000	10800 OSDH	1-11	SOIL	1.68	mg/kg	
LR-03	6000	10800 OSDH	1-11	VEG	56.00	mg/kg	
LR-04	5800	12100 OSDH	1-11	SOIL	1.12	mg/kg	
LR-04	5800	12100 OSDH	1-11	VEG	5.60	mg/kg	
	5800	12100 OSDH	1-11	SOIL	1.12	mg/kg	
LR-04		12100 OSDH	1-11	VEG	5.60	mg/kg	
LR-04	5800			SOIL	1.92	mg/kg	
LR-05	5975	13950 OSDH	1-11	SUIL	1.96	mg, kg	

NEAR FIELD ENVIRONMENTAL SAMPLING RESULTS FOR FLUORIDE SAMPLES

LABEL	NORTH	EAST	LAB	DATE	MEDIA	RESULT	UNITS	COMMENTS
LR-05	5975	13950	OSDH	1-11	VEG	7.12	mg/kg	
LR-05	5975	13950		1-11	SOIL	1.92	mg/kg	
LR-05	5975	13950		1-11	VEG	7.12	mg/kg	
LR-05A	6100	14450		1-14	VEG	13.20	mg/kg	
LR-05A	6100	14450		1-14	VEG	13.20	mg/kg	
LR-05A	6100	14450		1-14	SOIL	1.36	mg/kg	
LR-06	6000	15600		1-14	VEG	4.48	mg/kg	
LR-06	6000	15600		1-14	SOIL	1.12	mg/kg	
LRES04	6950		KMTC	1-04	VEG	3.0	ppm	TREE ?DATA DRCT
LRES05	6950		KMTC	1-04	SOIL	120.0	ppm	THE FORTH ONOT
LRES06	6950		KMTC	1-04	SOIL	60.0	ppm	
LRES07	6950		KMTC	1-04	WATER	<0.1	mg/1	RAIN BARREL
P-01	12670		KMTC	1-06	VEG	270.0	ppm	RATH DARACE
P-07	12670		KMTC	1-06	VEG	2600.0	ppm	
P-07	12670		KMTC	1-06	SOIL		ppm	
P-07	12670		KMTC	1-09	SOIL	170.0	ppm	
P-12	12670		KMTC	1-06	VEG	450.0	ppm	
P-13	12670		KMTC	1-06	SOIL	50.0	ppm	
P-13	12670		KMTC	1-06	VEG	2900.0	ppm	AVERAGE VALUE
P-13	12670		KMTC	1-09	SOIL	290.0	19 J. 19	AVERAUE VALUE
P-19	12670		KMTC	1-06	SOIL	210.0	ppm ppm	
P-19	12670		ORNL	1-06	SOIL	28.7	ug/g	WET
P-19	12670		ORNL	1-06	SOIL	32.8	ug/g	DRY
P-19	12670		ORNL	1-06	VEG	2900.0	ug/g	UKI
P-19A	12670		KMTC	1-06	VEG	450.0		
P-25	12670		KMTC	1-06	VEG	320.0	ppm ppm	
P-25	12670		KMTC	1-06	SOIL	90.0	ppm	
P-25	12670		ORNL	1-06	SOIL	14.2	ug/g	WET
P-25	12670		ORNL	1-06	SOIL	21.0	ug/g	DRY
P-25	12670		ORNL	1-06	VEG	280.0		UKI
PW 1632	13250		KMTC	1-15	WATER	<0.1	ug/g	
R-01	11750	10700		1-11	VEG	33.60	mg/l	
R-01	11750	10700		1-11	SOIL	3.32	mg/kg	
R-02	10800	10700		1-11	VEG	84.00	mg/kg	
R-02	10800	10700		1-11	SOIL	2.64	mg/kg	
R-03	10750	11550		1-11	VEG	6.80	mg/kg	
R-03	10750	11550		1-11	SOIL	4.00	mg/kg	
R-15	6050	13750		1-04	VEG	32.0	mg/kg	
R-16	6050	13750		1-04	SOIL		ppm	
RD-01	8150	10550			VEG	80.0	ppm	
RD-01	8150	10550		1-16 1-16	SOIL	40.00	mg/kg	
RD-01						5.20	mg/kg	
RD-02	7400	10750		1-16	VEG	46.40	mg/kg	
RD-02 RD-03	7400		OSDH	1-16	SOIL	2.96	mg/kg	
	6700		OSDH	1-16	VEG	53.60	mg/kg	
RD-03	6700		OSDH	1-16	SOIL	0.48	mg/kg	
RD-03	6600		OSDH	1-16	WATER	0.22	mg/1	
RD-04	6300	10/50	OSDH	1-16	VEG	64.00	mg/kg	

2 1.00

					1.1.1.1.1.1	and the second second	distant 1	
LABEL	NORTH	EAST	LAB	DATE	MEDIA	RESULT	UNITS	COMMENTS
RD-04	6300	10750	OSDH	1-16	SOIL	8.80	mg/kg	
RD-04	6350	10750		1-16	WATER	0.30	mg/l	
RD-05	6200	10850		1-16	WATER	0.17	mg/1	
RD-06	6900	10200		1-16	WATER	0.19	mg/1	
RD-07	10100	10700		1-16	WATER	0.32	mg/1	
RD-08	9900		OSDH	1-16	WATER	0.20	mg/1	
SM-20	7450	17250		1-04	SOIL	40.0	ppm	
SM-21	7450	17250		1-04	VEG	5.0		LEAVES
STRM GH	12650		OSDH	1-16	WATER	7.00	ppm	LEAVES
T1-01	12500		OSDH	1-15			mg/1	
					SOIL	0.8	mg/kg	THEREPERT
T1-01	12500		OSDH	1-05	VEG	1.2	mg/kg	EVERGREEN
T1-01	12500		OSDH	1-14	VEG	3.68	mg/kg	0.014
T1-01	12500		ORNL	1-05	SOIL	3.1	ug/g	DRY
T1-01	12500		ORNL	1-05	SOIL	2.3	ug/g	WET
T1-01	12500		ORNL	1-05	VEG	1.7	ug/g	
T1-01	12500		ORNL	1-05	VEG	1.7	ug/g	
T1-02	12500		OSDH	1-05	SOIL	<0.4	mg/kg	
T1-02	12500		OSDH	1-05	VEG	1.8	mg/kg	GRASS
T1-02	12500		OSDH	1-14	SOIL	0.92	mg/kg	
T1-02	12500		ORNL	1-05	SOIL	1.4	ug/g	DRY
T1-02	12500		ORNL	1-05	SOIL	1.2	ug/g	WET
T1-02	12500	8000	ORNL	1-05	VEG	1.2	ug/g	
T1-03	12500	8550	OSDH	1-05	SOIL	1.3	mg/kg	
T1-03	12500	8550	OSDH	1-05	VEG	2.6	mg/kg	GRASS
T1-03	12500	8550	ORNL	1-05	SOIL	4.2	ug/g	DRY
T1-03	12500	8550	ORNL	1-05	SOIL	3.4	ug/g	WET
T1-03	12500	8550	ORNL	1-05	VEG	4.7	ug/g	
T1-03&4	12500	8750	OSDH	1-14	WATER	0.23	mg/1	
T1-04	12500		OSDH	1-05	SOIL	0.9	mg/kg	
T1-04	12500		OSDH	1-05	VEG	10.4	mg/kg	EVERGREEN
T1-04	12500		OSDH	1-14	VEG	4.96	mg/kg	
T1-04	12500		OSDH	1-14	SOIL	1.84	mg/kg	
T1-04	12500		ORNL	1-05	SOIL	7.5	ug/g	DRY
T1-04	12500		ORNL	1-05	SOIL	5.7	ug/g	WET
T1-04	12500		ORNL	1-05	VEG		ug/g	HL.I
T1-05	12500		OSDH	1-05	SOIL	4.8	mg/kg	
T1-05	12500		OSDH	1-05	VEG	408.0		EVERGREEN
T1-05	12500		OSDH	1-14	VEG		mg/kg	EVENUNEEN
						28.8	mg/kg	
T1-05	12500		OSDH	1-14	SOIL	4.08	mg/kg	0.0V
T1-05	12500		ORNL	1-05	SOIL	23.5	ug/g	DRY
T1-05	12500		ORNL	1-05	SOIL	17.7	ug/g	WET
T1-05	12500		ORNL	1-05	VEG	1800.0	ug/g	
T1-06	12500		OSDH	1-05	SOIL	3.6	mg/kg	
T1-06	12500	9500	OSDH	1-05	VEG	500.0	mg/kg	GRASS

LABEL	NORTH	EAST LAB	DATE	MEDIA	RESULT	UNITS	COMMENTS
T1-06	12500	9500 OSDH	1-14	VEG	1320.0	mg/kg	
T1-06	12500	9500 OSDH	1-14	SOIL	3.64	mg/kg	
T1-06	12500	9500 OSDH	1-14	WATER	1.80	mg/1	
T1-06	12500	9500 ORNL	1-05	SOIL	22.5	ug/g	WET
T1-06	12500	9500 ORNL	1-05	SOIL	26.3	ug/g	DRY
T1-06	12500	9500 ORNL	1-05	VEG	930.0	ug/g	
T1-07	12500	9850 OSDH	1-05	SOIL	2.0	mg/kg	
T1-07	12500	9850 OSDH	1-05	VEG	216.0	mg/kg	GRASS
T1-07	12500	9850 OSDH	1-14	VEG	240.0	mg/kg	
T1-07	12500	9850 OSDH	1-14	SOIL	3.00	mg/kg	
T1-07	12500	9850 ORNL	1-05	SOIL	12.2	ug/g	WET
T1-07	12500	9850 ORNL		SOIL	22.8	ug/g	DRY
T1-07	12500	9850 ORNL	1-05	VEG	750.0	ug/g	
T1-08	12600	10200 OSDH	1-05	SOIL	0.8	mg/kg	
T1-08	12600	10200 OSDH	1-05	VEG	55.0	mg/kg	GRASS
T1-08	12600	10200 OSDH	1-14	VEG	20.0	mg/kg	
T1-08	12600	10200 OSDH	1-14	SOIL	2.20	mg/kg	
T1-08	12600	10200 ORNL	1-05	SOIL	6.6	ug/g	DRY
T1-08	12600	10200 ORNL	1-05	SOIL	5.1	ug/g	WET
T1-08	12600	10200 ORNL	1-05	VEG	44.6	ug/g	
T1-09	12600	10600 OSDH	1-05	SOIL	0.6	mg/kg	and do
T1-09	12600	10600 OSDH	1-05	VEG	20.0	mg/kg	GRASS
T1-09	12600	10600 OSDH	1-14	VEG	17.6	mg/kg	
T1-09	12600	10600 OSDH	1-14	SOIL	1.96	mg/kg	0.014
T1-09	12600	10600 ORNL	1-05	SOIL	2.1	ug/g	DRY
T1-09	12600	10600 ORNL	1-05	SOIL	1.7	ug/g	WET
T1-09	12600	10600 ORNL	1-05	VEG	17.5	ug/g	
T1-10	12750	11100 OSDH	1-05	WATER	0.2	mg/1	
T1-10	12750	11100 ORNL	1-05	WATER	0.21	ug/ml	
T1-10	12750	11100 ORNL	1-05	WATER	0.19	ug/ml	
T1-11	12600	11300 OSDH	1-05	SOIL	1.0	mg/kg	COACC
T1-11	12600	11300 OSDH	1-05	VEG	11.3	mg/kg	GRASS
T1-11	12600	11300 ORNL	1-05	SOIL	7.6	ug/g	DRY WET
T1-11	12600	11300 ORNL	1-05	SOIL	4.6	ug/g	WEI
T1-11 T1-12	12600			VEG			
T1-12		11900 OSDH				mg/kg	EVERGREEN
T1-12				VEG		mg/kg	EVERGREEN
T1-12				VEG		mg/kg	
T1-12 T1-12			1-14			mg/kg	DRY
			1-05		2.4	ug/g ug/g	WET
T1-12 T1-12	12550		1-05	SOIL VEG		ug/g ug/g	ALL I
T1-12 T1-13	12550 12650	13200 OSDH				mg/kg	
T1-13	12650		1-14			mg/kg	
		13200 050H	1-14	WATED	0.48	mg/kg mg/l	
T1-13 T1-14	12200 13200	14100 OSDH	1-14	VEC	5.92	mg/kg	
T1-14	13200	14100 OSDH	1-14	SOIL	0.96	mg/kg	
T1-14	13250	14100 OSDH	1-14	WATER	0.29	mg/1	
11-14	19520	14000 0304	1 14	MATCH	0.23	ing/ i	

LABEL	NORTH	EAST	LAB	DATE	MEDIA	RESULT	UNITS	COMMENTS
T2-01	10750	6200	ORNL	1-05	SOIL	1.1	ug/g	DRY
T2-01	10750	6200	ORNL	1-05	SOIL		ug/g	WET
T2-01	10750	6200	ORNL	1-05	VEG		ug/g	
T2-01	10750	6200	OSDH	1-05			mg/kg	
T2-01	10750	6200	OSDH	1-05			mg/kg	
T2-02	10550	6800	ORNL	1-05			ug/g	DRY
T2-02	10550	6800	ORNL	1-05				WET
T2-02	10550	6800	OSDH	1-14			mg/kg	
T2-02	10550		OSDH	1-05		0.48		
T2-02	10550	6800	OSDH	1-14			mg/kg	
T2-03	10650		OSDH	1-05			mg/kg	
T2-03	10650		ORNL	1-05			ug/g	
T2-03	10650		ORNL				ug/g	WET
T2-03	10650		ORNL				ug/g	DRY
T2-03	10650		OSDH					
T2-03	10650		OSDH		SOIL		mg/kg	
T2-03	10650		OSDH	1-14		3.12	mg/kg	
T2-04	10350		ORNL	1-05	WATER		ug/ml	
T2-04	10700		ORNL	1-05				DRY
T2-04	10700		ORNL					WET
T2-04	10700		OSDH					
T2-04	10700		OSDH				mg/kg	
T2-04	10700		OSDH		VEG	4.72	mg/kg	
T2-04	10350		OSDH	1-14	WATER		mg/l	
T2-04	10350		OSDH	1-05			mg/1	
T2-05	10800		ORNL	1-05			ug/g	DRY
T2-05	10800		ORNL	1-05	SOIL	1.2	ug/g	WET
T2-05	10800		ORNL	1-05	VEG	250.0	ug/g	
T2-05	10800		OSDH	1-14	SOIL	0.8	mg/kg	
T2-05	10800		OSDH	1-05		0.4	mg/kg	
T2-05	10800		OSDH	1-14				
T2-05	10800		OSDH	1-05	VEG	91.3	mg/kg	
T2-06	11450		ORNL	1-05	WATER			
T2-06	11450		ORNL	1-05	SOIL			DRY
T2-06		9950	ORNL		SOIL		ug/g	
T2-06	11450		ORNL	1-05	VEG	290.0	ug/g	
T2-06	11500		OSDH	1-14	SOIL	1.4	mg/kg	
T2-06	11500	9850	OSDH	1-06	SOIL	1.12	mg/kg	
T2-06	11500		OSDH	1-14	VEG	584.0	mg/kg	
T2-06	11500	9850	OSDH	1-06	VEG	122.2	mg/kg	
T2-06	11500	9950	OSDH	1-06	WATER	0.37	mg/1	
T2-07	11700	11100	ORNL	1-06	WATER	0.09	ug/ml	
T2-07	11500	10650		1-06	SOIL	14.9	ug/g	DRY
T2-07	11500	10650		1-06	SOIL	14.3	ug/g	WET
T2-07	11500	10650	OSDH	1-14	SOIL	0.92	mg/kg	
T2-07	11500	10650	OSDH	1-06	SOIL	6.72	mg/kg	

LABEL	NORTH	EAST	LAB	DATE	MEDIA	RESULT	UNITS	COMMENTS
T2-07	11500	10650	0cnu	1-14	VEG	22.4	na/ka	
T2-07	11500	10650		1-05	WATER	0.23	mg/kg	
T2-08	11500	11650		1-05	VEG	3.1	mg/l	
T2-08	11500	11650		1-16	WATER	0.17	mg/kg mg/l	
T2-08	11500	11650		1-06	SOIL	13.9	ug/g	DRY
T2-08	11500	11650		1-06	SOIL	10.2	ug/g	WET
T2-08	11500	11650		1-06	VEG	2.3	ug/g	MEI
T2-08	11500	11650		1-14	SOIL	11.2	mg/kg	
T2-08	11500	11650		1-06	SOIL	3.12	mg/kg	
T2-08	11500	11650		1-14	VEG	18.8	mg/kg	
T2-09	11350	12700		1-06	VEG	<0.8	mg/kg	
T2-09	11350	12700		1-06	SOIL	0.6	ug/g	DRY
T2-09	11350	12700		1-06	SOIL	0.5	ug/g	WET
T2-09	11350	12700		1-06	VEG	2.3	ug/g	
T2-09	11350	12700	OSDH	1-06	SOIL	<0.8	mg/kg	
T2-10	11400	13700		1-06	VEG	1.6	mg/kg	
T2-10	11400	13700	OSDH	1-16	WATER	0.16	mg/1	
T2-10	11400	13700	ORNL	1-06	SOIL	4.1	ug/g	WET
T2-10	11400	13700	ORNL	1-06	SOIL	6.5	ug/g	DRY
T2-10	11400	13700		1-06	VEG	2.7	ug/g	
T2-10	11400	13700	OSDH	1-14	SOIL	1.32	mg/kg	
T2-10	11400	13700		1-06	SOIL	1.6	mg/kg	
T2-10	11400	13700		1-14	VEG	5.52	mg/kg	
T3-01	8300		OSDH	1-06	SOIL	0.8	mg/kg	
T3-01	8300		OSDH	1-06	VEG	4	mg/kg	EVERGREEN
T3-01	8300		OSDH	1-14	VEG	2.48	mg/kg	
T3-01	8300		OSDH	1-14	SOIL	1.32	mg/kg	-
T3-01	8300		ORNL	1-06	SOIL	1.1	ug/g	DRY
T3-01	8300		ORNL	1-06	SOIL	0.9	ug/g	WET
T3-01	8300	6250	ORNL	1-06	VEG	4.7	ug/g	
T3-02	8600		OSDH	1-06	SOIL	1.2	mg/kg	CUEDODEEN
T3-02	8600		OSDH	1-06	VEG	2.4e	mg/kg	EVERGREEN
T3-02	8600	6950		1-06	SOIL	3.4	ug/g	DRY
T3-02	8600		ORNL	1-06	SOIL	2.9	ug/g	WET
T3-02	8600		ORNL	1-06	VEG .	3.0	ug/g	
T3-02A	8250		OSDH	1-16	VEG	1.76	mg/kg	
T3-02A	8250		OSDH	1-16	SOIL	3.00	mg/kg	
T3-03	8650		OSDH	1-06	SOIL	2.5	mg/kg	CUEDODEEN.
T3-03	8650		OSDH	1-06	VEG	8.8	mg/kg	EVERGREEN
T3-03	8650		OSDH	1-14	VEG	9.28	mg/kg	
T3-03	8650		OSDH	1-14	SOIL	2.36	mg/kg	DRV
T3-03	8650	9000		1-06	SOIL	3.2	ug/g	DRY
T3-03	8650		ORNL	1-06	SOIL	2.8	ug/g	WET
T3-03	8650		ORNL	1-06	VEG VEG	7.2 5.92	ug/g	
T3-03A	8500		OSDH OSDH	1-16		21.60	mg/kg	
T3-03A T3-03B	8500			1-16	SOIL VEG	41.60	mg/kg	
13-038	8300	3320	OSDH	1-16	VCO	41.00	mg/kg	

NEAR FIELD ENVIRONMENTAL SAMPLING RESULTS FOR FLUORIDE SAMPLES

LABEL	NORTH	EAST LAB	DATE	MEDIA	RESULT	UNITS	COMMENTS
T3-03B	8300	9950 OSDH	1-16	SOIL	5.56	mg/kg	
T3-03C	8800	10050 OSDH	1-16	VEG	13.68	mg/kg	
T3-03C	8800	10050 OSDH	1-16	SOIL	<0.40	mg/kg	
T3-04	8350	10300 OSDH	1-06	SOIL	<0.8	mg/kg	
T3-04	8350	10300 OSDH	1-06	VEG	102.0	mg/kg	EVRGRN/LEAVES
T3-04	8350	10300 OSDH	1-14	VEG	4.00	mg/kg	
T3-04	8350	10300 OSDH	1-14	SOIL	0.52	mg/kg	
T3-04	8350	10700 OSDH	1-16	VEG	73.60	mg/kg	
T3-04	8350	10300 ORNL	1-06	SOIL	1.4	ug/g	DRY
T3-04	8350	10300 ORNL	1-06	SOIL	1.3	ug/g	WET
T3-04	8350	10300 ORNL	1-06	VEG	420.0	ug/g	
T3-04A	8350	10700 OSDH	1-16	SOIL	8.80	mg/kg	
T3-04A	8150	10800 OSDH	1-16	WATER	0.12	mg/1	
T3-04B	8350	11000 OSDH	1-16	VEG	72.00	mg/kg	
T3-04B	8350	11000 OSDH	1-16	SOIL	1.68	mg/kg	
T3-05	8350	11200 OSDH	1-06	SOIL	1.2	mg/kg	
T3-05	8350	11200 OSDH	1-06	VEG	93.0	mg/kg	EVERGREEN
T3-05	8350	11200 OSDH	1-14	VEG	30.4	mg/kg	
T3-05	8350	11200 OSDH	1-14	SOIL	0.84	mg/kg	
T3-05	8350	11200 ORNL	1-06	SOIL	2.2	ug/g	DRY
T3-05	8350	11200 ORNL	1-06	SOIL	1.6	ug/g	WET
T3-05	8350	11200 ORNL	1-06	VEG	67.0	ug/g	
T3-06	8400	12050 OSDH	1-06	SOIL	1.4	mg/kg	
T3-06	8400	12050 OSDH	1-06	VEG	5.8	mg/kg	EVERGREEN
T3-06	8400	12050 OSDH	1-14	VEG	5.84	mg/kg	
T3-06	8400	12050 OSDH	1-14	SOIL	0.76	mg/kg	
T3-06	8400	12050 ORNL	1-06	SOIL	2.1	ug/g	DRY
T3-06	8400	12050 ORNL	1-06	SOIL	1.7	ug/g	WET
T3-06	8400	12050 ORNL	1-06	VEG	7.0	ug/g	
T3-07	8550	13050 OSDH	1-06	SOIL	<0.8	mg/kg	
T3-07	8550	13050 OSDH	1-06	VEG	4.2	mg/kg	GRASS
T3-07	8550	13050 OSDH	1-14	VEG	8.48	mg/kg	
T3-07	8550	13050 OSDH	1-14	SOIL	0.64	mg/kg	
T3-07	8550	13050 ORNL	1-06	SOIL	2.0	ug/g	DRY
T3-07	8550	13050 ORNL	1-06	SOIL	1.4	ug/g	WET
T3-07	8550	13050 ORNL	1-06	VEG	3.6	ug/g	
T5-01	17700	8700 OSDH	1-15	VEG	5.92	mg/kg	
T5-01	17700	8700 OSDH	1-15	SOIL	1.84	mg/kg	
T5-01	17700	8700 OSDH	1-08	SOIL	0.68	mg/kg	
T5-01	17700	8700 OSDH	1-08	VEG	2.48	mg/kg	
T5-02	16500	10000 OSDH	1-08	VEG	4.56	mg/kg	
T5-02	16500	10000 OSDH	1-08	VEG	1.76	mg/kg	
T5-03	16300	10350 OSDH	1-08	SOIL	<0.4	mg/kg	
T5-03	16300	10350 OSDH	1-08	VEG	3.6	mg/kg	
T5-04	15500	8500 OSDH	1-08	VEG	2.56	mg/kg	
T5-04	1550	8500 OSDH	1-08	WATER	0.32	mg/l	
T6-01	13950	7700 OSDH	1-08	SOIL	3.56	mg/kg	
T6-01	13950	7700 OSDH	1-08	VEG	28.8	mg/kg	

LABEL	NORTH	EAST	LAB	DATE	MEDIA	RESULT	UNITS	COMMENTS
T6-02	13900	8400	OSDH	1-08	SOIL	9.0	mg/kg	
T6-02	13900	8400	OSDH	1-08	VEG	48.0	mg/kg	
T6-03	14650	8300	OSDH	1-08	SOIL	2.28	mg/kg	
T6-03	14650	8300	OSDH	1-08	VEG	10.48	mg/kg	
T6-04	14650	8700	OSDH	1-08	SOIL	1.8	mg/kg	
T6-05	14450	9550	OSDH	1-08	VEG	15.84	mg/kg	

Table 5.2.7.19A

LABEL	NORTH	EAST	LAB	DATE	MEDIA	RESULT	UNITS	COMMENTS
1W	14500	14800	KMTC	1-06	WATER	0.11	mg/1	
2W	13250	13950		1-06	WATER	0.23	mg/1	
3W	12750	11200		1-06	WATER	0.19		
4W	11700	11150		1-06	WATER	<0.002	mg/1	
5W	10100	10950		1-06			mg/1	
6W	9850					0.011	mg/1	
7W			KMTC	1-06		<0.002	mg/1	
	10600		KMTC	1-06		0.32	mg/1	
8W	11800		KMTC	1-06	California de la calegaria de la	0.6	mg/1	
10W	11800		KMTC	1-06			mg/1	
11W	16350		KMTC	1-06		<0.002	mg/1	
12W	6750	11050		1-06		0.003	mg/1	
13W	7475	13350		1-06		<0.002	mg/1	
14W	18450		KMTC	1-06		<0.002	mg/1	
14W	18450		KMTC	1-08		0.003	mg/1	
15W	12900		KMTC	1-06		0.018	mg/1	
15W	12900		KMTC	1-08	WATER	0.013	mg/1	
17W	7650	4350	KMTC	1-06	WATER	0.006	mg/1	
18W	15100	10300	KMTC	1-08	WATER	0.002	mg/1	
19W	6850	9050	KMTC	1-08	WATER	0.004	mg/1	
20W	7750	17000	KMTC	1-08			mg/1	
103	13800	10200		1-09			ppm	
105		8450		1-09			ppm	
106		8450		1-09			ppm	
108		8750		1-09			ppm	
109	10600	8750	KMTC	1-09			ppm	
111		8450		1-09	VEG		ppm	
112	14250		KMTC	1-09	SOIL		ppm	
115	11500	11600		1-09				
116	10650	12300		1-09			ppm	
117	10650	12300		1-09			ppm	
118	9750	12450		1-09		20.0	ppm	
119				1-09	COTI	2. /	ppm	
120		12450		1-09	SOIL	20.0	ppm	
120		8850		1-09		7.7	ppm	
	2100	0000		1-09	SOIL	20.0	ppm	
122		9900		1-09		480.0		
123	9600		KMTC	1-09	SOIL		ppm	
124	9700	10700		1-09	VEG	110.0	ppm	
125	9700	10700		1-09	SOIL	20.0	ppm	
126	9750	11850		1-09	VEG	4.6	ppm	
127	9750	11850		1-09	SOIL	<10	ppm	
128	10650	11300	KMTC	1-09	VEG	22.0	ppm	
129	10650	11300		1-09	SOIL	<10	ppm	
130	10650	10350	KMTC	1-09	VEG	71.0	ppm	
131	10650	10350	KMTC	1-09	SOIL	30.0	ppm	
132	11600	10650	KMTC	1-09	VEG	25.0	ppm	
133	11600	10650		1-09	SOIL	15.0	ppm	
134	9250	13150		1-09	VEG	5.4	ppm	
135	9250	13150		1-09	SOIL	<10	ppm	
	100 million (100 m						be berry	

LABEL	NORTH	EAST	LAB	DATE	MEDIA	RESULT	UNITS	COMMENTS
136	9150	12050	KMTC	1-09	VEG	7.2	ppm	
137	9150		KMTC				ppm	
138	9200	11000		1-09			ppm	
139	9200		KMTC				ppm	
135	9300		KMTC				ppm	
141	9150		KMTC				ppm	
142	9150		KMTC				ppm	
143	10500		KMTC				ppm	
145	10500		KMTC				ppm	
145	11600			1-09			ppm	
140	11600		KMTC				ppm	
149	11650		KMTC				ppm	
149	11450		KMTC				ppm	
187	12600		KMTC				ppm	
190		9500		1-14			ppm	
193	12550		KMTC				ppm	
195	12600		KMTC				ppm	
198	12550	11850		1-14			ppm	
199	12550	11850		1-14			ppm	
201	13150	10200		1-14			ppm	
202	13150		KMTC				ppm	
202	13150		KMTC				ppm	
204	13150	11250		1-14			ppm	
208	12500		KMTC	1-14	VEG		ppm	
209	12500		KMTC	1-14	SOIL		ppm	
211	12500		KMTC	1-14			ppm	
214	12500		KMTC	1-14		1530.0	ppm	
218	12500		KMTC	1-14	SOIL		ppm	
221	12500		KMTC	1-14	SOIL	40.0	ppm	
300	12450		KMTC	1-06	WATER		mg/1	
301	12450		KMTC	1-14	WATER		mg/1	
302	12550		KMTC	1-06	WATER		mg/1	
303	12550		KMTC	1-14	WATER		mg/1	
303	12550		KMTC	1-06	WATER	1.2	mg/1	
305	12550		KMTC				mg/1	
306	11450		KMTC	1-14	WATER	0.25	mg/1	
306	11450		KMTC	1-06	WATER	0.25	mg/1	
307	11450		KMTC	1-14	WATER	0.21	mg/1	
50	10750	10100		1-06	VEG	17000.0	ppm	
			KMTC	1-06	VEG	4200.0	ppm	
A-01	13130		KMTC	1-06	SOIL	2180.0		
A-01 A-07	13130		KMTC	1-06	VEG	8000.0	ppm ppm	
A-07	13130 13130		KMTC	1-08	SOIL	590.0	ppm	
			KMTC	1-09	SOIL	730.0	ppm	
A-07	13130		KMTC	1-06	SOIL	655.0		
A-13	13130		KMTC	1-06	VEG	6200.0	ppm	
A-13	13130			1-08	SOIL	200.0	ppm	
A-13	13130		KMTC				ppm	
A-19	13130	3100	KMTC	1-06	VEG	7150.0	ppm	

LABEL	NORTH	EAST	LAB	DATE	MEDIA	RESULT	UNITS	COMMENTS
		6,10,1		LITTL.	116.0111	HESOEI	01110	CONFICTING S
A-25	13130	9950	KMTC	1-06	VEG	49.0	ppm	
A-25	13130	9950	KMTC	1-06	SOIL	25.0	ppm	
A-25	13130	9760	ORNL	1-06		10.5	ppm	WET
A-25	13130		ORNL	1-06	SOIL	13.2	ppm	DRY
A-25	13130		ORNL	1-06		11.8	ppm	
BG 1/4	14750		ORNL	1-06	VEG	3.91	ppm	
BG 1/4	14750		ORNL	1-06	SOIL	6.1	ppm	WET
BG 1/4	14750		ORNL	1-06		7.88	ppm	DRY
BG 1/4			KMTC	1-06		64.0	ppm	DRT
BG 1/4			KMTC	1-06		<10	ppm	
CARLILE	14250	15100		1-05	the set of the	2.7	ppm	WET
CARLILE	14250	15100		1-05		3.7		DRY
CARLILE	14250	15100		1-05		0.188	ppm	DRT
COMB ST	12600		KMTC	1-07		0.29	ppm	
COMB ST	12600		KMTC	1-06			mg/1	
COMB ST	12600		KMTC	1-05	WATER	0.95	mg/1	
COMB ST	12600		KMTC			0.36	mg/1	
COMB ST	12600		KMTC	1-04	WATER	0.36	mg/1	
COMB ST				1-03	WATER	0.61	mg/1	
COMB ST	12600		KMTC	1-02	WATER	0.57	mg/1	
		8650		1-01	WATER	0.43	mg/1	
CRES	6150	10800		1-05	WATER	0.05	ug/ml	
CRES	6150	10750		1-05	SOIL	2.9	ppm	DRY
CRES	6150	10750		1-05	SOIL	2.0	ppm	WET
CRES	6150	10800		1-05	VEG	13.8	ppm	
CRES08	6150	10750		1-04	VEG	12.0	ppm	
CRES09	6150	10750		1-04	S-V	<10	ppm	
CRES10	6150	10750		1-04	SOIL	<10	ppm	
DRES11	6100	11100		1-04	SOIL	15.0	ppm	
DRES12	6100	11100		1-04	VEG	220.0	ppm	
DRES13	6100	11100		1-04	SOIL	13.0	ppm	
E-01	13000		KMTC	1-06	VEG	1340.0	ppm	
E-01	13000		KMTC	1-06	SOIL	210.0	ppm	
E-01	13000		KMTC	1-09	SOIL	300.0	ppm	
E-07	13000		KMTC	1-06	VEG	14200.0	ppm	
E-07	13000	9410		1-06	SOIL	465.0	ppm	
E-07	13000		KMTC	1-09	SOIL	380.0	ppm	
E-13	13000		KMTC	1-06	VEG	7420.0	ppm	
E-13	13000		KMTC	1-06	SOIL	130.0	ppm	
E-13	13000		KMTC	1-09	SOIL	130.0	ppm	
E-19	13000		KMTC	1-06	SOIL	15.0	ppm	
E-19	13000		KMTC	1-06	VEG	21600.0	ppm	
E-25	13000		KMTC	1-06	VEG	40.0	ppm	
E-25	13000		KMTC	1-06	SOIL	20.0	ppm	
E-25	13000	9950	ORNL	1-06	SOIL	8.9	ppm	WET
E-25	13000	9950	ORNL	1-06	SOIL	10.3	ppm	DRY
E-25	13000		ORNL	1-06	VEG	24.2	ppm	
HRES01	6950		KMTC	1-04	SOIL	<10	ppm	
HRES02	6950		KMTC	1-04	VEG	32.0	ppm	

LABEL	NORTH	EAST	LAB	DATE	MEDIA	RESULT	UNITS	COMMENTS
HRES03	6950	10650	KMTC	1-04	VEG	0.95	ppm	
J-07	12840		KMTC	1-06	VEG	1520.0	ppm	
J-07	12840		KMTC	1-06	SOIL	125.0	ppm	
J-07	12840		KMTC	1-09		250.0	ppm	
J-13	12840		KMTC	1-06	SOIL	200.0	ppm	
J-13	12840		KMTC	1-06	VEG	6640.0	ppm	
J-13	12840		KMTC	1-09	SOIL	125.0	ppm	
J-19	12840		KMTC	1-06	VEG	3220.0	ppm	
J-19	12840		KMTC	1-06	SOIL	165.0	ppm	
J-25	12840		KMTC	1-06	VEG	210.0	ppm	
J-25	12840		KMTC	1-06	SOIL	20.0	ppm	
J-25	12840		ORNL	1-06	SOIL	9.8	ppm	WET
J-25	12840		ORNL	1-06	SOIL	13.1	ppm	DRY
J-25	12840		ORNL	1-06	VEG	148.0	ppm	
LAGOON	12800		KMTC	1-17	WATER	I	mg/1	
LRES04	6950		KMTC	1-04	VEG	<0.4	ppm	
LRES05	6950		KMTC	1-04	SOIL	<10	ppm	
LRES06	6950		KMTC	1-04	SOIL	<10	ppm	
LRES07	6950		KMTC	1-04	WATER	<0.002	mg/1	FRONT PORCH
P-01	12670		KMTC	1-06	VEG	5300.0	ppm	
P-07	12670		KMTC	1-06	VEG	4480.0	ppm	
P-07	12670		KMTC	1-06	SOIL	235.0	ppm	
P-07	12670		KMTC	1-09	SOIL	90.0	ppm	
P-12	12670		KMTC	1-06	VEG	5940.0	ppm	
P-13	12670		KMTC	1-06	VEG	3800.0	ppm	
P-13	12670		KMTC	1-06	SOIL	55.0	ppm	
P-13	12670		KMTC	1-09	SOIL	175.0	ppm	
P-19	12670		KMTC	1-06	SOIL	90.0	ppm	
P-19	12670		ORNL	1-06	SOIL	135.0	ppm	WET
P-19	12670		ORNL	1-06	SOIL	154.0	ppm	DRY
P-19	12670	9780		1-06	VEG	3600.0	ppm	
P-25	12670		KMTC	1-06	VEG	360.0	ppm	
P-25	12670		KMTC	1-06	SOIL	40.0	ppm	
P-25	12670	9780		1-06	SOIL	19.2	ppm	WET
P-25	12670	9780		1-06	SOIL	28.4	ppm	DRY
P-25	12670		ORNL	1-06	VEG	186.0	ppm	
PW 1632	13250		KMTC	1-15	WATER	<0.002	mg/1	
R-15	6050	13750		1-04	VEG	0.6	ppm	
R-16	6050	13750		1-04	SOIL	<10	ppm	
SM-20	7450	17250		1-04	SOIL	11.0	ppm	
SM-21	7450	17250		1-04	VEG	0.6	ppm	
T1-01	12500		ORNL	1-05	SOIL	4.6	ppm	WET
T1-01	12500		ORNL	1-05	SOIL	6.2	ppm	DRY
T1-02	12500		ORNL	1-05	SOIL	5.7	ppm	DRY
T1-02	12500		ORNL	1-05	SOIL	4.5	ppm	WET
T1-03	12500		ORNL	1-05	SOIL	5.7	ppm	WET
12 00	20000	0000	Sector Sector	A. 16.10.				

LABEL	NORTH	EAST	LAB	DATE	MEDIA	RESULT	UNITS	COMMENTS
T1-03	12500	8550	ORNL	1-05	SOIL	7.0	ppm	DRY
T1-04	12500		ORNL	1-05	SOIL	13.0	ppm	WET
T1-04	12500		ORNL	1-05	SOIL	17.3	ppm	DRY
T1-05	12500		ORNL	1-05	SOIL		ppm	WET
T1-05	12500		ORNL	1-05	SOIL	45.1	ppm	DRY
T1-06	12500		ORNL	1-05		90.7	ppm	DRY
T1-06	12500		ORNL	1-05			ppm	WET
T1-07	12500		ORNL				ppm	WET
T1-07	12500		ORNL		SOIL		ppm	DRY
T1-08	12600	10200		1-05	SOIL	7.3	ppm	DRY
T1-08	12600		ORNL		SOIL	5.7	ppm	WET
T1-09	12600	10600		1-05		4.5	ppm	DRY
T1-09	12600	10600		1-05	SOIL	3.6	ppm	WET
T1-10	12750	11100		1-05	WATER		ug/ml	WLI
T1-10	12750	11100		1-05	WATER		ug/ml	
T1-11	12600	11300		1-05	SOIL	6.0	ppm	DRY
T1-11	12600	11300		1-05	SOIL	3.7	ppm	WET
T1-12	12550	11900		1-05	SOIL	6.0	ppm	DRY
T1-12	12550	11900		1-05	SOIL	3.5	ppm	WET
T2-01	10750		ORNL	1-05	SOIL	2.9		WET
T2-01	10750		ORNL	1-05	SOIL	3.5	ppm	DRY
T2-02	10550		ORNL	1-05	SOIL	2.6	ppm	DRY
T2-02	10550		ORNL	1-05	SOIL	2.2	ppm ppm	WET.
T2-03	10650		ORNL	1-05	SOIL	3.2	ppm	WET
T2-03	10650		ORNL	1-05	SOIL	3.9	ppm	DRY
T2-04	10350		ORNL	1-05		0.003	ug/ml	UNI
T2-04	10700		ORNL	1-05	SOIL	4.1	ppm	DRY
T2-04	10700		OKNL	1-05	SOIL	3.2	ppm	WET
T2-05	10800		ORNL	1-05	SOIL	4.5	ppm	DRY
T2-05	10800		ORNL	1-05	SOIL	3.9	ppm	WET
T2-05	10800		ORNL	1-05	VEG	393.0	ppm	WL I
T2-06	11450		ORNL	1-05	WATER	0.27	ug/ml	
T2-06	11450		ORNL	1-05	SOIL	25.9	ppm	DRY
T2-06	11450		ORNL	1-05	SOIL	24.4	ppm	WET
T2-06	11450		ORNL	1-05	VEG	285.0	ppin	WL I
T2-07	11700	11100		1-06	WATER	0.03	ug/m1	
T2-07	11500	10650		1-06	SOIL	5.23	ppm .	DRY
T2-07	11500	10650		1-06	SOIL	5.0	ppm	WET
T2-08	11500	11650		1-06	SOIL	5.21	ppm	DRY
T2-08	11500	11650		1-06	SOIL	3.8	ppm	WET
T2-08	11500	11650		1-06	VEG	0.633	ppm	HL I
T2-09	11350	12700		1-06	SOIL	3.95	ppm	DRY
T2-09	11350	12700		1-06	SOIL	3.0	ppm	WET
T2-09	11350		ORNL	1-06	VEG	1.53	ppm	
T2-10	11400		ORNL	1-06	SOIL	2.4	ppm	WET
T2-10	11400		ORNL	1-06	SOIL	3.86	ppm	DRY
T2-10	11400		ORNL	1-06	VEG	0.596	ppm	with the second s
T3-01	8300		ORNL	1-06	SOIL	4.34	ppm	DRY

LABEL	NORTH	EAST I	LAB	DATE	MEDIA	RESULT	UNITS	COMMENTS
T3-01	8300	6250 0	ORNL	1-06	SOIL	3.3	ppm	WET
T3-01	8300	6250 0	ORNL	1-06	VEG	5.03	ppm	
T3-02	8600	6950 (ORNL	1-06	SOIL	3.86	ppm	DRY
T3-02	8600		ORNL	1-06	SOIL	3.3	ppm	WET
T3-02	8600		ORNL	1-06	VEG	3.21	ppm	
T3-03	8650		ORNL	1-06	SOIL	5.34	ppm	DRY
T3-03	8650		ORNL	1-06	SOIL	4.7	ppm	WET
T3-03	8650		ORNL	1-06	VEG	3.39	ppm	
T3-04	8350		ORNL	1-06	SOIL	4.28	ppm	DRY
T3-04	8350		ORNL	1-06	SOIL	3.9	ppm	WET
T3-04	8350		ORNL	1-06	VEG	624.0	ppm	
T3-05	8350	11200 (ORNL	1-06	SOIL	4.0	ppm	DRY
T3-05	8350		ORNL	1-06	SOIL	3.0	ppm	WET
T3-05	8350		ORNL	1-06	VEG	177.0	ppm	
T3-06	8400		ORNL	1-06	SOIL	4.88	ppm	DRY
T3-06	8400		ORNL	1-06	SOIL	3.9	ppm	WET
T3-06	8400		ORNL	1-06	VEG	6.0	ppm	
T3-07	8550	the second s	ORNL	1-06	SOIL	3.9	ppm	DRY
T3-07	8550		ORNL	1-06	SOIL	2.8	ppm	WET
T3-07	8550		ORNL	1-06	VEG	35.6	ppm	

Table 5.2.7.20A

LABEL	NORTH	EAST LAB	DATE	MEDIA	RESULT	UNITS	COMMENTS
CARLILE	14250	15100 OSDH	1-05	VEG	0.8	pCi/g	GRASS
CARLILE	14250	15100 OSDH	1-05	SOIL			
CRES	6150	10750 OSDH	1-05	SOIL	-0.1	pCi/g	
CRES		10750 OSDH	1-05	VEG	1.1	pCi/g	GRASS
CRES	6150	10800 OSDH	1-05	WATER	0.0	pCi/1	
HRES	6950	10650 OSDH	1-08	VEG	5.0	pCi/q	
T1-01	12500	7200 OSDH	1-05	SOIL	0.1	pCi/q	
T1-01	12500	7200 OSDH	1-05	VEG	1.4	pCi/g	EVERGREEN
T1-02	12500	8000 OSDH	1-05	SOIL	0.1	pCi/g	
T1-02	12500	8000 OSDH	1-05	VEG	0.1	pCi/a	GRASS
T1-03	12500	8550 OSDH	1-05	SOIL	0.4	pCi/a	
T1-03	12500	8550 OSDH	1-05	VEG	0.5	pCi/a	GRASS
T1-04	12500	8900 OSDH	1-05	SOIL	0.9	pCi/a	
T1-04	12500	8900 OSDH	1-05	VEG	1.3	pCi/a	EVERGREEN
T1-05	12500	9200 OSDH	1-05	SOIL	0.4	pCi/a	
T1-05	12500	9200 OSDH	1-05	VEG	4020.0	pCi/g	EVERGREEN
T1-06	12500	9500 OSDH	1-05	SOIL	1.8	pCi/o	
T1-06	12500	9500 OSDH	1-05	VEG	1220.0	pCi/g	GRASS
T1-07	12500	9850 OSDH	1-05	SOIL	-0.3	pCi/g	
T1-07	12500	9850 OSDH	1-05	VEG	124.0	pCi/g	GRASS
T1-08	12600	10200 OSDH	1-05	SOIL	-0.3	pCi/g	
T1-08	12600	10200 OSDH	1-05	VEG	1.0	nCi/a	GRASS
T1-09	12600	10600 OSDH	1-05	SOIL	1.6	pCi/g	
T1-09	12600	10600 OSDH	1-05	VEG	1.6	pCi/g	GRASS
T1-10	12750	11100 OSDH	1-05	WATER	100.0	pCi/1	
T1-11	12600	11300 OSDH	1-05	SOIL	-0.3	pCi/g	
T1-11	12600	11300 OSDH	1-05	VEG	1.3	pCi/g	GRASS
T1-12	12550	11300 OSDH 11300 OSDH 11900 OSDH 11900 OSDH 6200 OSDH 6200 OSDH 6800 OSDH 7550 OSDH 7550 OSDH 8400 OSDH	1-05	SOIL	0.4	pCi/g	
T1-12	12550	11900 OSDH	1-05	VEG	0.0	pCi/g	EVERGREEN
T2-01	10750	6200 OSDH	1-05	SOIL	0.1	pCi/g	
T2-01	10750	6200 OSDH	1-05	VEG	-0.2	pCi/g	
T2-02	10550	6800 OSDH	1-05	SOIL	0.1	pCi/g	
T2-03	10650	7550 OSDH	1-05	SOIL	132.0	pCi/g	
T2-03	10650	7550 OSDH	1-05	VEG	-0.4	pCi/g	
T2-04	****	0100 00011	* ***	of the sk too	57 - A.	hourd	
T2-04	10350	8400 OSDH	1-05	WATER	150.0	pCi/1	
T2-05	10800	9300 OSDH	1-05	SOIL	0.1	pCi/g	
T2-05	10800	9300 OSDH	1-05	VEG	78.0	pCi/g	
T2-06	11500	9850 OSDH	1-06	SOIL	1.9	pCi/g	
T2-06	11500	9850 OSDH	1-06	VEG	41.0	pCi/g	
T2-06	11500	9950 OSDH	1-06	WATER	124.0	pCi/1	
T2-07	11500	10650 OSDH	1-06	SOIL	3.0	pCi/g	
T2-08	11500	11650 OSDH	1-06	SOIL	0.1	pCi/g	
T2-08	11500	11650 OSDH	1-06	VEG	-0.5	pCi/g	
T2-09	11350	12700 OSDH	1-06	SOIL	1.2	pCi/g	
T2-09	11350	12700 OSDH	1-06	VEG	0.1	pCi/g	
T2-10	11400	13700 OSDH	1-06	SOIL	-0.3	pCi/g	
T2-10	11400	13700 OSDH	1-06	VEG	0.3	pCi/g	

LABEL	NORTH	EAST LAB	DATE	MEDIA	RESULT	UNITS	COMMENTS
LADEL	NORTH	CHOI CHO	DATE	TIL WATT			
T3-01	8300	6250 OSDH	1-06	SOIL	-0.3	pCi/g	
T3-01	8300	6250 OSDH	1-06	VEG	1.6	pCi/g	EVERGREEN
T3-02	8600	6950 OSDH	1-06	SOIL	3.0	pCi/g	
T3-02	8600	6950 OSDH	1-06	VEG	1.0	pCi/g	EVERGREEN
T3-03	8650	9000 OSDH	1-06	SOIL	-0.3	pCi/g	
T3-03	8650	9000 OSDH	1-06	VEG	0.6	pCi/g	EVERGREEN
T3-04	8350	10300 OSDH	1-06	SOIL	-0.3	pCi/g	
T3-04	8350	10300 OSDH	1-06	VEG	1070.0	pCi/g	EVRGRN/LEAVES
T3-05	8350	11200 OSDH	1-06	SOIL	0.9	pCi/g	
T3-05	8350	11200 OSDH	1-06	VEG	27.0	pCi/g	EVERGREEN
T3-06	8400	12050 OSDH	1-06	SOIL	0.4	pCi/g	
T3-06	8400	12050 OSDH	1-06	VEG	0.7	pCi/g	EVERGREEN
T3-07	8550	13050 OSDH	1-06	SOIL	0.1	pCi/g	
T3-07	8550	13050 OSDH	1-06	VEG	0.6	pCi/g	GRASS
T5-01	17700	8700 OSDH	1-08	SOIL	-0.3	pCi/g	
T5-01	17700	8700 OSDH	1-08	VEG	0.2	pCi/g	
T5-02	16500	10000 OSDH	1-08	VEG	-0.2	pCi/g	
T5-02	16500	10000 OSDH	1-08	VEG	0.0	pCi/g	
T5-03	16300	10350 OSDH	1-08	SOIL	-0.3	pCi/g	
T5-03	16300	10350 OSDH	1-08	VEG	0.0	pCi/g	
T5-04	15500	8500 OSDH	1-08	VEG	0.2	pCi/g	
T5-04	1550	8500 OSDH	1-08	WATER	7.0	pCi/1	
T6-01	13950	7700 OSDH	1-08	SOIL	0.1	pCi/g	
T6-01	13950	7700 OSDH	1-08	VEG	1.2	pCi/g	
T6-02	13900	8400 OSDH	1-08	SOIL	-0.3	pCi/g	
T6-02	13900	8400 OSDH	1-08	VEG	2.5	pCi/g	
T6-03	14650	8300 OSDH	1-08	SOIL	0.4	pCi/g	
T6-03	14650	8300 OSDH	1-08	VEG	1.0	pCi/g	
T6-04	14650	8700 OSDH	1-08	SOIL	0.4	pCi/g	
T6-05	14450	9550 OSDH	1-08	VEG	0.9	pCi/g	

Table 5.2.7.21A

NEAR FIELD ENVIRONMENTAL SAMPLING RESULTS FOR BETA SAMPLES

LABEL	NORTH	EAST	LAB	DATE	MEDIA	RESULT	UNITS	COMMENTS
CARLILE		15100	OSDH	1-05		0.7		GRASS
CARLILE			OSDH	1-05	SOIL	-0.3	pCi/g	
	6150	10750	OSDH	1-05	SOIL	-0.0	pCi/g	
CRES			OSDH	1-05	VEG	0.2	pCi/g	GRASS
CRES	6150	10800	OSDH	1-05	WATER	-1.0	pCi/1	
HRES	6950	10650	OSDH	1-08	VEG	1.9	pC1/g	
T1-01	12500	7200	OSDH	1-05	SOIL	0.6	pCi/g	
T1-01	12500	7200	OSDH	1-05	VEG	1.0	pCi/g	EVERGREEN
T1-02	12500	8000	OSDH	1-05	SOIL	0.4	pCi/g	
T1-02	12500	8000	OSDH	1-05	VEG	1.2	pCi/g	GRASS
T1-03		8550	OSDH	1-05	SOIL	-0.0	pCi/g	
T1-03	12500	8550	OSDH	1-05	VEG	1.1	pCi/g	GRASS
T1-04					SOIL	-0.1	pCi/g	
T1-04		8900	OSDH	1-05	VEG	0.1	pCi/q	EVERGREEN
T1-05		9200	OSDH	1-05	SOIL	0.2	pCi/g	
T1-05	12500	9200	OSDH	1-05	VEG	110.0	pCi/g	EVERGREEN
T1-06	12500	9500	OSDH	1-05	SOIL	-0.2	pCi/q	
T1-06				1-05	VEG	95.0	pCi/g	GRASS
T1-07				1-05	SOIL	-0.0	pCi/g	
T1-07				1-05	VEG	23.0	pCi/g	GRASS
T1-08		10200	OSDH	1-05	SOIL	-0.2	pCi/a	
T1-08	12600	10200	OSDH	1-05	VEG	5.0	pCi/a	GRASS
T1-09	12600	10600	OSDH	1-05	SOIL	-0.5	pCi/a	
T1-09	12600	10600	OSDH	1-05	VEG	0.5	pCi/q	GRASS
T1-10	12750	11100	OSDH	1-05	WATER	12.0	pCi/l	
T1-11			OSDH	1-05	SOIL	-0.4	pCi/g	
T1-11			OSDH	1-05	VEG	1.7	pCi/g	GRASS
T1-12			OSDH	1-05	VEG SOIL VEG	0.2	pCi/g	
T1-12				1-05	VEG	0.3	pCi/g	EVERGREEN
T2-01	10750			1-05	SOTI	-0.5	nCila	
T2-01				1-05	VEG SOIL SOIL	0.7	pCi/g	
T2-02	10550	6800	OSDH	1-05	SOIL	0.1	pCi/g	
T2-03				1-05	SOIL	23.0	pCi/g	
T2-03	10650			1-05	VEG	3.1	pCi/g	
T2-04	10700	8400	OSDH	1-05	SOIL	0.6	pCi/g	
T2-04	10350		OSDH	1-05	WATER	9.0	pCi/1	
T2-05	10800		OSDH	1-05	SOIL	0.7	pCi/g	
T2-05	10800		OSDH	1-05	VEG	7.0	pCi/g	
T2-06	11500		OSDH	1-06	SOIL	-0.7	pCi/g	
T2-06	11500		OSDH	1-06	VEG	6.0	pCi/g	
T2-06	11500		OSDH	1-06	WATER	26.0	pCi/l	
T2-07	11500	10650		1-06	SOIL	0.7	pCi/g	
T2-08	11500		OSDH	1-06	SOIL	0.4	pCi/g	
T2-08	11500		OSDH	1-06	VEG	6.0	pCi/g	
T2-09	11350		OSDH	1-06	SOIL	0.3	pCi/g	
T2-09	11350		OSDH	1-06	VEG	2.3	pCi/g	
T2-10	11400	13700		1-06	SOIL	0.6	pCi/g	
T2-10	11400	13700	OSDH	1-06	VEG	0.3	pCi/g	

NEAR FIELD ENVIRONMENTAL SAMPLING RESULTS FOR BETA SAMPLES

1

LABEL	NORTH	EAST LAB	DATE	MEDIA	RESULT	UNITS	COMMENTS
T3-01	8300	6250 OSDH	1-06	SOIL	-0.2	pCi/g	
T3-01	8300	6250 OSDH	1-06	VEG	0.6	pCi/g	EVERGREEN
T3-02	8600	6950 OSDH	1-06	SOIL	1.0	pCi/g	
T3-02	8600	6950 OSDH	1-06		0.4	pCi/g	EVERGREEN
T3-03	8650	9000 OSDH	1-06		1.0	pCi/g	
T3-03	8650	9000 OSDH	1-06	VEG	0.0	pCi/g	EVERGREEN
T3-04	8350	10300 OSDH	1-06	SOIL	0.1	pCi/g	
T3-04	8350	10300 OSDH	1-06		77.0	pCi/g	EVRGRN/LEAVES
T3-05	8350	11200 OSDH	1-06		0.4	pCi/g	
T3-05	8350	11200 OSDH	1-06		8.0	pCi/g	EVERGREEN
T3-06	8400	12050 OSDH	1-06		0.2	pCi/g	
T3-06	8400	12050 OSDH	1-06		-0.1	pCi/g	EVERGREEN
T3-07	8550	13050 OSDH	1-06		-0.1	pCi/g	Salara (Teller)
T3-07	8550	13050 OSDH	1-06	VEG	-0.1	pCi/g	GRASS
T5-01	17700	8700 OSDH	1-08	SOIL	0.6	pCi/g	
T5-01	17700	8700 OSDH	1-08		-0.0	pCi/g	
T5-02	16500	10000 OSDH	1-08	VEG	-0.2	pCi/g	
T5-02	16500	10000 OSDH	1-08	VEG	-0.1	pCi/g	
T5-03	16300	10350 OSDH	1-08	SOIL	-0.2	pCi/g	
T5-03	16300	10350 OSDH	1-08	VEG	-0.3	pCi/g	
T5-04	15500	8500 OSDH	1-08	VEG	-0.1	pCi/g	
T5-04	1550	8500 OSDH	1-08	WATER	5.0	pCi/1	
T6-01	13950	7700 OSDH	1-08	SOIL	0.7	pCi/g	
T6-01	13950	7700 OSDH	1-08	VEG	0.7	pCi/g	
T6-02	13900	8400 OSDH	1-08	SOIL	-0.2	pCi/g	
T6-02	13900	8400 OSDH	1-08	VEG	1.1	pCi/g	
T6-03	14650	8300 OSDH	1-08	SOIL	-0.2	pCi/g	
T6-03	14650	8300 OSDH	1-08	VEG	0.3	pCi/g	
T6-04	14650	8700 OSDH	1-08	SOIL	-0.6	pCi/g	
T6-05	14450	9550 OSDH	1-08	VEG	-0.5	pCi/g	
a second second							

FAR FIELD ENVIRONMENTAL SAMPLING RESULTS BY SAMPLE ANALYSIS SORTED BY LABEL

Table 5.2.7.22A

FAR FIELD ENVIRONMENTAL SAMPLING RESULTS FOR FLUORIDE SAMPLES

LABEL	TOWNSHP	RANGE	SECTI	ON SUBSE	ECT LAB	DATE	MEDIA	RESUL	TS UNITS	COMMENTS
16W	T11N	R21E	34	B6	KMTC	1-06	WATER	0.2	MG/L	
1H14	TIIN	R21E	1	H14	KMTC	1-10	SOIL	60.0	PPM	
2E14	TIIN	R21E	2	E14	KMTC	1-10	SOIL	100.0	PPM	
4A 6	T11N	R21E	4	A6	KMTC	1-10	SOIL	170.0	PPM	
4J14	T11N	R21E	4	J14	KMTC	1-09	SOIL	220.0	PPM	
5N 6	T11N	R22E	5	N6	KMTC	1-09	SOIL	160.0	PPM	
6A 7	T11N	R22E	6	A7	KMTC	1-09	SOIL	130.0	PPM	
6C14	TIIN	R22E	6	C14	KMTC	1-09	SOIL	190.0	PPM	
6I11	T11N	R22E	6	I11	KMTC	1-09	SOIL	130.0	PPM	
7C 6	TIIN	R22E	7	C6	KMTC	1-08	SOIL	130.0	PPM	
7C 6	T11N	R22E	7	C6	KMTC	1-08	VEG	10.0	PPM	J. GRASS
7F 7	TIIN	R22E	7	F7	KMTC	1-08	SOIL	120.0	PPM	
7F 7	TIIN	R22E	7	F7	KMTC	1-08	VEG	20.0	PPM	WHEAT
7F 7	T11N	R22E	7	F7	KMTC	1-08	VEG	30.0	PPM	
8A 1	TIIN	R22E	8	Al	KMTC	1-09	SOIL	140.0	PPM	
8J 7	TIIN	R21E	8	J7	KMTC	1-09	SOIL	50.0	PPM	
8M14	T11N	R22E	8	M14	KMTC	1-09	SOIL	30.0	PPM	
9B 2	T11N	R21E	9	B2	KMTC	1-09	SOIL	250.0	PPM	
9I 8	TIIN	R21E	9	18	KMTC	1-09	SOIL	90.0	PPM	
108 8	T11N	R21E	10	88	KMTC	1-09	SOIL	30.0	PPM	
10L 8	TIIN	R21E	10	L8	KMTC	1-08	SOIL	30.0	PPM	
10M10	TIIN	R21E	10	M10	KMTC	1-08	SOIL	30.0	PPM	
11F 7	TIIN	R21E	11	F7	KMTC	1-08	SOIL	50.0	PPM	
11F 7	TIIN	R21E	11	F 7	KMTC	1-08	VEG	10.0	PPM	CAT LEAVES
11F 7	T11N	R21E	11	F7	KMTC	1-08	VEG	43.0	PPM	BEANS
11F 7	TIIN	R21E	11	F7	KMTC	1-08	VEG	12.0	PPM	WILLOW
11N 6	TIIN	R21E	11	N6	KMTC	1-08	SOIL	50.0	PPM	
12F 7	TIIN	R21E	12	F7	KMTC	1-09	SOIL	50.0	PPM	
12N 8	TIIN	R21E	12	N8	KMTC	1-08	SOIL	130.0	PPM	
12N 8	TIIN	R21E	12	N8	KMTC	1-08	VEG	5.0	PPM	
12N 8	TIIN	R21E	12	N8	KMTC	1-08	VEG	7.0	PPM	CAT TIPS
12N 8	TIIN	R21E	12	N8	KMTC	1-08	VEG	2.0	PPM	CEDAR
138 1	T11N	R21E	13	81	KMTC	1-08	SOLL	50.0	PPM	
13B 1	TIIN	R21E	13	81	KMTC	1-08	VEG	4.0	PPM	
13G11	T11N		13	G11	KMTC	1-09		80.0	PPM	
13G11	TIIN	R21E	13	G11	KMTC	1-09	VEG	10.0	PPM	
13H 2	TIIN	R21E	13	H2	KMTC	1-08	SOIL	60.0	PPM	
14B 4	TIIN	R21E	14	84	KMTC	1-08	SOIL	100.0	PPM	
14B 4	TIIN	R21E	14	84	KMTC	1-08	VEG	12.0	PPM	
14B 8	TIIN	R21E	14	88	KMTC	1-08	SOIL	70.0	PPM	
14B 8	TIIN	R21E	14	88	KMTC	1-08	VEG	6.0	PPM	
14B12	TIIN	R21E	14	B12	KMTC	1-08	SOIL	60.0	PPM	
14812	TIIN	R21E	14	812	KMTC	1-08	VEG	5.0	PPM	FESCUE TIP
14B12	TIIN	R21E	14	812	KMTC	1-08	VEG	4.0	PPM	FESCUE
16810	TIIN	R21E	16	810	KMTC	1-09	SOIL	40.0	PPM	
16810	TIIN	R21E	16	810	KMTC	1-09	VEG	9.0	PPM	
16L 2	TIIN	R21E	16	L2	KMTC	1-10	SOIL	280.0	PPM	
16L12	TIIN	R21E	16	L12	KMTC	1-10	SOIL	120.0	PPM	
16L12	TIIN	R21E R21E	16	L12	KMTC	1-10	VEG	50.0	PPM	
18B 4	TIIN	R22E	18	84	KMTC	1-09	SOIL	60.0	PPM	

FAR FIELD ENVIRONMENTAL SAMPLING RESULTS FOR FLUORIDE SAMPLES

LABEL	TOWNSHP	RANGE	SECTION	SUBSECT	LAB	DATE	MEDIA	RESULTS	UNITS	COMMENT	S
18B 4	T11N	R22E	18	84	KMTC	1-09	VEG	30.0	PPM		
18H10	T11N	R22E	18	H10	KMTC	1-09	SOIL	90.0	PPM		
18H10	T11N	R22E	18	H10	KMTC	1-09	VEG	9.0	PPM		
1A11	T11N	R21E	1	A11	KMTC	1-09	SOIL	90.0	PPM		
1A11	T11N	R21E	1	A11	KMTC	1-09	VEG	40.0	PPM		
1H14	TIIN	R21E	1	H14	KMTC	1-08	VEG	20.0	PPM		
22G14	TIIN	R21E	22	G14	KMTC	1-10	SOIL	70.0	PPM		
22G14	T11N	R21E	22	G14	KMTC	1-10	VEG	1.0	PPM		
23A10	TIIN	R21E	23	A10	KMTC	1-08	SOIL	50.0	PPM	TIIN	
23A10	T11N	R21E	23	A10	KMTC	1-08	VEG	8.0	PPM	TIIN	
25E 8	T12N	R21E	25	E8	KMTC	1-10	SOIL	60.0	PPM	TTH	
25E 8	T12N	R21E	25	E8	KMTC	1-10	VEG	10.0	PPM		
25H 3	TI2N	R21E	25	H3	KMTC	1-10	SOIL	60.0	PPM		
25H 3	TI2N	R21E	25	H3	KMTC	1-10					
25H 3	T12N	R21E	25	H3	KMTC		VEG	9.0	PPM		
25H14	TI2N	R21E	25			1-10	VEG	10.0	PPM		
25114	T12N			H14	KMTC	1-10	SOIL	80.0	PPM		
26E 7		R21E	25	L14	KMTC	1-14	WATER	<0.1	MG/L		
26F 6	T12N	R21E	26	E7	KMTC	1-14	WATER	<0.1	MG/L		
	T12N	R21E	26	F6	KMTC	1-14	WATER	<0.1	MG/L	SURFACE	WATER
26G 1	T12N	R21E	26	G1	KMTC	1-10	SOIL	60.0	PPM		
26G 1	T12N	R21E	26	G1	KMTC	1-10	VEG	10.0	PPM		
26G 8	T12N	R21E	26	G8	KMTC	1-10	SOIL	70.0	PPM		
26G 8	T12N	R21E	26	G8	KMTC	1-10	VEG	20.0	PPM		
26L 1	T12N	R21E	26	Ll	KMTC	1-14	WATER	<0.1	MG/L	SURFACE	(POND)
27A13	T12N	R21E	27	A13	KMTC	1-08	SOIL	90.0	PPM		
27A13	T12N	R21E	27	A13	KMTC	1-08	VEG	6.0	PPM	CEDAR	
27A13	T12N	R21E	27	A13	KMTC	1-08	VEG	4.0	PPM	SMILAX	
270 1	T12N	R21E	27	01	KMTC	1-10	SOIL	90.0	PPM		
27D 1	T12N	R21E	27	01	KMTC	1-10	VEG	80.0	PPM		
27H 9	T12N	R21E	27	H9	KMTC	1-10	SOIL	180.0	PPM		
27H 9	T12N	R21E	27	H9	KMTC	1-10	VEG	10.0	PPM		
27H 9	T12N	R21E	27	H9	KMTC	1-10	VEG	20.0	PPM		
27I 0	T12N	R21E	27	10	KMTC	1-14	WATER	<0.1	MG/L	SURFACE	WATER
27K 4	T12N	R21E	27	K4	KMTC	1-10	SOIL	50.0	PPM		
27K 4	T12N	R21E		K4	KMTC	1-10		20.0	PPM		
27L 1	T12N	R21E	27	L1	KMTC	1-14	WATER	0.1	MG/L		
27L 4	T12N	R21E	27	14	KMTC	1-14	WATER	0.1	MG/L	SURFACE	WATER
27L 8	T12N	R21E	27	L8	KMTC	1-14	WATER	<0.1	MG/L	SURFACE	
28G 4	T12N	R21E	28	G4 .	KMTC	1-10	SOIL	100.0	PPM	JUNIAGE	MAILA
28G 4	T12N	R21E	28	G4	KMTC	1-10	VEG	1.0	PPM		
28G 4	T12N	R21E	28	G4	KMTC	1-10	VEG	20.0	PPM		
28H14	T12N	R21E	28	H14	KMTC	1-14	WATER		MG/L	SURFACE	WATER
28J11	T12N	R21E	28	J11	KMTC	1-10		0.1	PPM	SURFACE	MATCH
28J11	T12N	R21E	28	J11	KMTC		SOIL	140.0			
28N 2	T12N	R21E R21E				1-10	VEG	20.0	PPM		
28N 2			28	N2	KMTC	1-10	SOIL	170.0	PPM		
28N 2	T12N	R21E	28	N2	KMTC	1-10	VEG	2.0	PPM		
	T12N	R21E	28	N2	KMTC	1-10	VEG	10.0	PPM		
290 3	T12N	R22E	29	03	KMTC	1-10	SOIL	50.0	PPM		
29D 3	T12N	R22E	29	D3	KMTC	1-10	VEG	2.0	PPM		

ø

ø

FAR FIELD ENVIRONMENTAL SAMPLING RESULTS FOR FLUORIDE SAMPLES

LABEL	TOWNSHP	RANGE	SECTION	SUBSECT	LAB	DATE	MEDIA	RESULTS	UNITS	COMMENTS
2E14	TIIN	R21E	2	E14	KMTC	1-09	VEG	20.0	PPM	
	TIIN	R21E	2	L13	KMTC	1-09	SOIL	50.	PPM	
2L13 2L13	TIIN	R21E	2	L13	KMTC	1-09	VEG	60.0	PPM	a shekara a s
		R21E	2	M8	KMTC	1-09	SOIL	50.0	PPM	
2M 8	TIIN	R21E R21E	2	MB	KMTC	1-09	VEG	40.0	PPM	
2M 8	TIIN	R22E	30	F7	KMTC	1-10	SOIL	60.0	PPM	
30F 7	T12N		30	F7	KMTC	1-10	VEG	6.0	PPM	
30F 7	T12N	R22E		F7	KMTC	1-10	VEG	6.0	PPM	
30F 7	T12N	R22E	30	B10	KMTC	1-10	SOIL	80.0	PPM	
31B10	T12N	R22E	31	B10	KMTC	1-10	VEG	10.0	PPM	
31810	T12N	R22E	31	B14	KMTC	1-10	SOIL	170.0	PPM	1.1.1
33814	T12N	R21E	33		KMTC	1-10	VEG	8.0	PPM	
33B14	T12N	R21E	33	814	KMTC	1-10	SOIL	130.0	PPM	1.6.7
331 3	T12N	R21E	33	13	KMTC	1-10	VEG	20.0	PPM	
33I 3	T12N	R21E	33	13	KMTC	1-10	SOIL	90.0	PPM	이 가지 않는 것이 없다.
33M14	T12N	R21E	33	M14		1-10	VEG	6.0	PPM	
33M14	T12N	R21E	33	M14	KMTC		VEG	5.0	PPM	이 이 집에 같은 것이 같이 없다.
33M14	T12N	R21E	33	M14	KMTC	1-10	SOIL	120.0	PPM	
34D 7	T12N	R21E	34	D7	KMTC	1-10		30.0	PPM	
34D 7	T12N	R21E	34	D7	KMTC	1-10	VEG	60.0	PPM	2. 1. Active 1
35A 9	T12N	R21E	35	A9	KMTC	1-10	SOIL		PPM	1. 1. 1. 1. 1. 1. 1. 1. 1. 1. 1. 1. 1. 1
35A 9	T12N	R21E	35	A9	KMTC	1-10	VEG	6.0	PPM	1 N. 1984
35A 9	T12N	R21E	35	A9	KMTC	1-10	VEG	11.0	PPM	
35B 1	T12N	R21E	35	81	KMTC	1-10	SOIL	70.0	PPM	
35B 1	T12N	R21E	35	81	KMTC	1-10	VEG	9.0		
35814	T12N	R21E	35	814	KMTC	1-10	SOIL	70.0	PPM	1943
35814	T12N	R21E	35	814	KMTC	1-10	VEG	3.0	PPM	
35F10	T12N	R21E	35	F10	KMTC	1-10	SOIL	60.0	PPM	
35F10	T12N	R21E	35	F10	KMTC	1-10	VEG	20.0	PPM	
35L 5	T12N	R21E	35	L5	KMTC	1-09	SOIL	70.0		
35L 5	T12N	R21E	35	L5	KMTC	1-09	VEG	5.0	PPM	
35L 5	T12N	R21E	35	L5	KMTC	1-09	VEG	8.0	PPM	
35L11	T12N	R21E	35	L11	KMTC	1-10	SOIL	210.0		1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1
35L11	T12N	R21E	35	L11	KMTC	1-10	VEG	30.0	PPM	1 A A A A A A A A A A A A A A A A A A A
35L11	T12N	R21E	35	L11	KMTC	1-10	VEG	<1	PPM	
36A13	T12N	R21E	36	A13	KMTC	1-10	SOIL	100.0		
36A13	T12N	R21E	36	A13	KMTC	1-10	VEG	20.0	PPM	
36011	T12N	R21E	36	011	KMTC	1-14	WATER	<0.1	MG/L	
36G 5	T12N	R21E	36	G5	KMTC	1-09	SOIL	70.0	PPM	
36G 5	T12N	R21E	36	G5	KMTC	1-09	VEG	20.0	PPM	
36H14	T12N	R21E	36	H14	KMTC	1-10	SOIL	160.0	PPM	
36H14	T12N	R21E	36	H14	KMTC	1-10	VEG	<1	PPM	
36L 5	T12N	R21E	36	L5	KMTC	1-14	WATER	<0.1	MG/L	CUDEACE WATER
36M10	T12N	R21E	36	M10	KMTC	1-14	WATER	<0.1	MG/L	SURFACE WATER
3F 7	T11N	R21E	3	F7	KMTC	1-09	SOIL	90.0	PPM	
3F 7	TIIN	R21E	3	F7	KMTC	1-09	VEG	8.0	PPM	
4A 1	TIIN	R21E	4	A1	KMTC	1-10	SOIL	140.0	PPM	
4A 1	TIIN	R21E	4	A1	KMTC	1-10	VEG	6.0	PPM	
4A 3	TIIN	R21E	4	A3	KMTC	1-10	VEG	20.0	PPM	
4A 6	TIIN	R21E	4	A6	KMTC	1-10	VEG	60.0	PPM	

FAR FIELD ENVIRONMENTAL SAMPLING RESULTS FOR FLUORIDE SAMPLES

LABEL	TOWNSHP	RANGE	SECTION	SUBSECT	LAB	DATE	MEDIA	RESULTS	UNITS	COMMENTS
4J14	T11N	R21E	4	J14	KMTC	1-09	VEG	15.0	PPM	
4J14	T11N	R21E	4	J14	KMTC	1-09	VEG	6.0	PPM	
4N 7	T11N	R21E	4	N7	KMTC	1-09	SOIL	90.0	PPM	
4N 7	T11N	R21E	4	N7	KMTC	1-09	VEG	30.0	PPM	
4N 8	T11N	R21E	4	N8	KMTC	1-09	VEG	6.0	PPM	
50/HWY	T12N	R21E	22	G1	KMTC	1-06	VEG	1500.0	PPM	CEDAR
5N 6	TIIN	R22E	5	N6	KMTC	1-09	VEG	30.0	PPM	020111
6A 7	TIIN	R22E	6	A7	KMTC	1-09	VEG	30.0	PPM	
6C14	TIIN	R22E	6	C14	KMTC	1-09	VEG	10.0	PPM	
8A 1	TIIN	R22E	8	Al	KMTC	1-09	VEG	8.0	PPM	
8J 7	TIIN	R21E	8	J7	KMTC	1-09	VEG	10.0	PPM	
8M14	TIIN	R22E	8	M14	KMTC	1-09	VEG	9.0	PPM	
9B 2	TIIN	R21E	9	B2	KMTC	1-09	VEG	5.5	PPM	
AB-01	TIIN	R21E	17	A14	OSDH	1-15	SOIL	0.64	MG/KG	
			17	A14		1-15	VEG	2.72	MG/KG	
AB-01	T11N	R21E			OSDH				MG/KG	
AB-02	TIIN	R21E	16	A11	OSDH	1-15	SOIL	0.68		
AB-02	TIIN	R21E	16	A11	OSDH	1-15	VEG	1.68	MG/KG	
AB-03	T12N	R21E	32	N14	OSDH	1-15	SOIL	0.44	MG/KG	
AB-03	T12N	R21E	32	N14	OSDH	1-15	VEG	3.28	MG/KG	
AB-04	TIIN	R21E	9	B1	OSDH	1-15	SOIL	0.56	MG/KG	
AB-04	TIIN	R21E	9	81	OSDH	1-15	VEG	1.52	MG/KG	
AB-05	T11N	R21E	4	C8	OSDH	1-15	SOIL	0.52	MG/KG	
AB-05	T11N	R21E	4	C8	OSDH	1-15	VEG		MG/KG	
AB-06	T12N	R21E	29	N14	OSDH	1-15	SOIL		MG/KG	
AB-06	T12N	R21E	29	N14	OSDH	1-15	VEG		MG/KG	
AB-07	T12N	R21E	31	A8	OSDH	1-15	SOIL		MG/KG	
AB-07	T12N	R21E	31	A8	OSDH	1-15	VEG		MG/KG	
AB-08	T12N	R21E	19	N8	OSDH	1-15	SOIL		MG/KG	
AB-08	T12N	R21E	19	N8	OSDH	1-15	VEG	1.60	MG/KG	
AB-09	T12N	R21E	18	M6	OSDH	1-15	SOIL	1.04	MG/KG	
AB-09	T12N	R21E	18	M6	OSDH	1-15	VEG	3.60	MG/KG	
AB-10	T12N	R20E	13	G14	OSDH	1-15	SOIL	0.52	MG/KG	
AB-10	T12N	R20E	13	G14	USDH	1-15	VEG	2.96	MG/KG	
BG 1/4	T12N	R21E	21	B13	ORNL	1-06	SOIL	4.7	UG/G	DRY
BG 1/4	T12N	R21E	21	813	ORNL	1-06	SOIL	3.7	UG/G	WET
BG 1/4	T12N	R21E	21	813	KMTC	1-06	SOIL	220.0	PPM	
BG 1/4	T12N	R21E	21	813	ORNL	1-06	VEG	7.3	UG/G	
BG 1/4	T12N	R21E	21	B13	KMTC	1-06	VEG	40.0	PPM	1/4 N.GATE
CARLILE		R21E	22	014	ORNL	1-05	SOIL	1.2	UG/G	DRY
CARLILE		R21E	22	014	ORNL	1-05	SOIL	0.8	UG/G	WET
CARLILE		R21E	22	014	OSDH	1-14	SOIL	1.24	MG/KG	
CARLILE		R21E	22	014	OSDH	1-05	SOIL	<0.4	MG/KG	
CARLILE		R21E	22	014	OSDH	1-05	VEG	4.3	GRASS	
CARLILE		R21E	22	D14	ORNL	1-05	VEG	8.3	UG/G	
CARLILE		R21E	22	D14	OSDH	1-14	VEG	2.56	MG/KG	
CNTRL01		R24E	16	-	OSDH	1-06	SOIL	0.4	MG/KG	
CNTRL01		R24E	16	*	OSDH	1-06	VEG	0.9	MG/KG	
CNTRL02		R22E	35	C14	OSDH	1-06	SOIL	0.5	MG/KG	
CNTRL02		R22E	35	C14	OSDH	1-06	VEG	1.9	MG/KG	
UNITAL UZ	1.7.014	N. C. C. L.	33	014	0301	1.00	ALC	1.3	mu/ nu	

FAR FIELD ENVIRONMENTAL SAMPLING RESULTS FOR FLUORIDE SAMPLES

LADEL	TOWNCUD	DANCE	SECTION	SUBSECT	LAR	DATE	MEDIA	RESULTS	UNITS	COMMENTS
LABEL	TOWNSHP	RANGE	SECTION	SUBSECT	LAD	UNIC	ALD IN	HEODEIO	Siters	COLUMN TELEVIS
CNTRL03	T11N	R22E	28	E10	OSDH	1-06	SOIL	<0.4	MG/KG	
CNTRL03		R22E	28	E10	OSDH	1-06	VEG	2.4	MG/KG	
CNTRL04		R22E	20	N10	OSDH	1-06	SOIL	0.7	MG/KG	
CNTRL04		R22E	20	N10	OSDH	1-06	VEG	2.5	MG/KG	
CNTRL05		R22E	19	F13	OSDH	1-06	SOIL	0.9	MG/KG	
CNTRL05		R22E	19	F13	0 SDH	1-06	VEG	2.2	MG/KG	
CNTRL06		R22E	4	I13	OSDH	1-06	VEG	1.4	MG/KG	
CNTRL07		R21E	13	87	OSDH	1-06	SOIL	<0.4	MG/KG	
CNTRL07		R21E	13	87	OSDH	1-06	VEG	0.9	MG/KG	
CNTRL08		R22E	8	K8	OSDH	1-06	SOIL	0.6	MG/KG	
CNTRL08		R22E	8	K8	OSDH	1-06	VEG	2.7	MG/KG	
CNTRL09	TIIN	R20E	21	C12	OSDH	1-06	SOIL	1.3	MG/KG	
CNTRL09	TIIN	R20E	21	C12	OSDH	1-06	VEG	2.2	MG/KG	
CNTRL09		R20E	21	C12	OSDH	1-06	VEG	1.9	MG/KG	
CNTRL10		R20E	21	113	OSDH	1-06	SOIL	1.3	MG/KG	
CNTRL10		R20E	21	I13	OSDH	1-06	VEG	3.0	MG/KG	
CNTRL11	T12N	R19E	26	E3	OSDH	1-06	SOIL	1.5	MG/KG	
CNTRL11	T12N	R19E	26	E3	OSDH	1-06	VEG	1.2	MG/KG	
CNTRL12	TIIN	R20E	14	N6	OSDH	1-06	SOIL	0.7	MG/KG	
CNTRL13	TIIN	R23E	3	C2	OSDH	1-06	SOIL	0.4	MG/KG	
CNTRL13		R23E	3	C2	OSDH	1-06	VEG	1.4	MG/KG	
CNTRL14	TIIN	R23E	-	-	OSDH	1-06	SOIL	0.5	MG/KG	S OF SALLISAW
CNTRL15		R20E	16	G13	OSDH	1-06	SOIL	0.5	MG/KG	
CNTRL15	T09N	R20E	16	G13	OSDH	1-06	VEG	1.6	MG/KG	
CRES	T12N	R21E	27	K3	OSDH	1-05	SOIL	0.6	MG/KG	
CRES	T12N	R21E	27	K3	ORNL	1-05	SOIL	2.8	UG/G	DRY
CRES	T12N	R21E	27	K3	ORNL	1-05	SOIL	1.9	UG/G	WET
CRES	T12N	R21E	27	K3	OSDH	1-05	VEG	32.0	MG/KG	GRASS
CRES	T12N	R21E	27	K3	OSDH	1-05	WATER	0.19	MG/L	
CRES	T12N	R21E	27	K3	ORNL	1-05	WATER	0.09	UM/ML	ODEC
CRES08	T12N	R21E	27	К3	KMTC	1-04	VEG	42.0	PPM	CRES
CRES09	T12N	R21E	27	K3	KMTC	1-04	S/V	90.0	PPM	MIX LEAVES/DIF
CRES10	T12N	R21E	27	K3	KMTC	1-04	SOIL	190.0	PPM	CRES
DIRTYCR		R21E	6	F9	OSDH	1-14	WATER	0.32	MG/L	DREC
DRES11		R21E	27	K4	KMTC	1-04	SOIL		PPM	DRES DRES/LEAVES
DRES12	T12N	R21E	27	K4	KMTC	1-04	VEG	160.0	PPM	DRES
DRES13	T12N	R21E	27	K4	KMTC	1-04	SOIL	60.0	MG/KG	ORES
G-01	T12N	R21E	8	L5	OSDH	1-15	SOIL	0.72	MG/KG	
G-01	T12N	R21E	8	15	OSDH	1-15	VEG	4.00	MG/KG	
G-02	T12N	R21E	8	83	OSDH	1-15	SOIL	1.84	MG/KG	
G-02	T12N	R21E	8	83	OSDH	1-15	VEG	8.00	MG/KG	
G-03	T12N	R21E	7	813	OSDH	1-15	SOIL	0.64	MG/KG	
G-03	T12N	R21E	/	813	OSDH	1-15	VEG	2.08	MG/KG	
G-04	T12N	R21E	/	K9	OSDH	1-15	SOIL	1.52	MG/KG	
G-04	T12N	R21E	7	K9	OSDH	1-15	VEG	3.04	MG/KG	
GL-01	T12N	R21E	16	F1	OSDH	1-15	SOIL	0.52	MG/KG	
GL-01	T12N	R21E	16	F1	OSDH	1-15	VEG	2.16	PPM	CEDAR
H-02	T12N	R21E	27	12	KMTC	1-04	VEG	24.0	PPM	HRES
H-03	T12N	R21E	27	12	KMTC	1-04	VEG	80.0	r.m	THE S

à.

	LABEL	TOWNSHP	RANGE	CECTION	CURCECT	1.40	DATE	NEDIL			
	LADEL	TUWNSHP	MANGE	SECTION	SUBSECT	LAB	DATE	MEDIA	RESULTS	UNITS	COMMENTS
	HC-01	TIIN	R21E	23	N1	OSDH	1-11	VEG	2.48	MG/KG	
	HC-02	TIIN	R21E	2.2	NI	OSDH		VEG	1.04	MG/KG	
	HC-03	TIIN	R21E	20	N14	OSDH		VEG	2.00	MG/KG	
1.	HC-04	T11N	R21E	21	Al	OSDH		VEG	1.28	MG/KG	
	HC-05	T11N	R21E	15	N5	OSDH		VEG	2.24	MG/KG	
11	HC-07	TIIN	R21E	26	N1		1-11	VEG	2.00	MG/KG	
	HC-08	TIIN	R21E		M11	OSDH	1-11	VEG	2.32	MG/KG	
		TIIN	R21E		A12	OSDH	1-11	VEG	1.52	MG/KG	
	HC-10	TION	R21E			OSDH	1-11	VEG			
	HC-11	TION	R21E	4		OSDH	1-11	VEG	2.24	MG/KG	
	HC-11	TION	R21E	4			1-11	VEG	1.60	MG/KG	
	HC-12	TIIN		22	A13	OSDH	1-11	VEG.	1.52	MG/KG	
	HC-13	TIIN	R21E				1-11	VEG	0.90	MG/KG	
	HC-14	T12N			A14	OSDH	1-11	VEG	3.60	MG/KG	
	HC-14		R21E		A14	OSDH	1-11	VEG	1.28	MG/KG	
	HC-16	TION	R21E		Al	OSDH	1-11	VEG	4.08	MG/KG	
		TIIN	R21E	33	N2	OSDH	1-11	VEG	5.76	MG/KG	
	HC-17	TION	R22E		A14	OSDH	1-11	VEG	2.80	MG/KG	
	HC-18	TION	R22E		A13	OSDH	1-11	VEG	3.60	MG/KG	
		TIIN	R22E	29	A14	OSDH	1-11	VEG	6.40	MG/KG	
		TIIN	R22E		A14			VEG	2.88	MG/KG	
		T12N	R21E		12		1-04	SOIL	110.0	PPM	HRES
		T12N	R21E	28	I14	KMTC	1-04	VEG	3.0	PPM	TREE
		T12N	R21E		I14	KMTC	1-04	SOIL	120.0	PPM	LRES
		T12N	R21E	28	I14	KMTC	1-04	SOIL	60.0	PPM	LRES
		T12N	R21E	28	I14	KMTC	1-04	WATER	<0.1	MG/L	RAIN BARREL
		T12N	R21E	26	A11	OSDH	1-13	SOIL	0.88	MG/KG	
		T12N	R21E	26	A11		1-13	VEG	2.24	MG/KG	
1	P-01	T12N	R21E	26	A11		1-13	WATER	0.23	MG/L	
1	P-02	T12N	R21E	25	H4		1-13	SOIL	0.50	MG/KG	
1	P-02	T12N	R21E	25	H4		1-13	VEG	<0.1	MG/KG	
1	P-03	T12N	R21E		A13		1-13	SOIL	0.68	MG/KG	
1	P-03	T12N	R21E	36	A13	OSDH.			1.36	MG/KG	
	P-04		R21E	36	D10	OSDH		SOIL	0.96	MG/KG	
F	P-04	T12N	R21E	36	D10 D10	OSDH		VEG	<0.1	MG/KG	
				1	A11	OSDH	1-13	SATI	5 56	MC/KC	
	P-05	T11N	R21E	ĩ	A11	OSDH	1-13	VEG	1.76		
	-06	T11N	R21E	ĩ	A7	OSDH	1-13	SOIL		MG/KG	
	-06	TIIN	R21E	î	A7	OSDH	1-13		0.52	MG/KG	
	-07	T12N	R21E		M10	OSDH		VEG	2.24	MG/KG	
	-07	T12N	R21E		M10		1-13	SOIL	<0.1	MG/KG	
	-08	T12N	R21E			OSDH	1-13	VEG	9.68	MG/KG	
	-08	T12N	R21E	35	D10	OSDH	1-13	SOIL	0.48	MG/KG	
	-09	T12N			010	OSDH	1-13	VEG	6.80	MG/KG	
	-09		R21E		HG	OSDH	1-13	SOIL	1.16	MG/KG	
	-10	T12N	R21E		H6	OSDH	1-13	VEG	4.64	MG/KG	
		T12N	R21E	28	F14	OSDH	1-13	WATER	0.49	MG/L	
	REDHILL		R21E		H4	OSDH	1-14	WATER	0.15	MG/L	
	RKRES15		R21E	34	A13	KMTC	1-04	VEG	32.0	PPM	RKRES/CEDAR
	RKRES16		R21E		A13	KMTC	1-04	SOIL	80.0	PPM	RKRES
	RES17	TIIN	R21E	34	86	K'ATC	1-04	SOIL	90.0	PPM	RRES

FAR FIELD ENVIRONMENTAL SAMPLING RESULTS FOR FLUORIDE SAMPLES

4

)

LABEL	TOWNSHP	RANGE	SECTION	SUBSECT	LAB	DATE	MEDIA	RESULTS	UNITS	COMMENTS
RRES18	TIIN	R21E	34	86	KMTC	1-04	VEG	7.0	PPM	RRES/CEDAR
RRES19	TIIN	R21E	34	B6	KMTC	1-04	SOIL	100.0	PPM	RRES
S-01	TIIN	R22E	34	Al	OSDH	1-11	VEG	2.8	MG/KG	
S-02	TIIN	R22E	28	NI	OSDH	1-11	VEG	6.64	MG/KG	
5-02	TIIN	R22E	30	N14	OSDH	1-11	VEG	3.92	MG/KG	
5-04	TIIN	R21E	25	N14	OSDH	1-11	VEG	2.56	MG/KG	
	TIIN	R21E	36	Al	OSDH	1-11	VEG	3.52	MG/KG	
S-05			34	N14	OSDH	1-11	VEG	3.12	MG/KG	
S-06	TIIN	R21E			OSDH	1-11	VEG	3.04	MG/KG	
S-07	TION	R21E	2	A1 .	OSDH	1-11	VEG	5.44	MG/KG	
S-08	TIIN	R22E	31	NI			VEG	2.16	MG/KG	
S-09	TION	R22E	6	A14	OSDH	1-11			MG/KG	
S-10	TION	R22E	5	NI	OSDH	1-11	VEG	5.68		
S-11	TION	R22E	9	A1	OSDH	1-11	VEG	5.04	MG/KG	
S-12	TION	R22E	4	A1	OSDH	1-11	VEG	2.16	MG/KG	
S-13	TIIN	R22E	33	N14	OSDH	1-11	VEG	4.96	MG/KG	
S-14	TION	R22E	6	NI	OSDH	1-11	VEG	3.68	MG/KG	
S-15	TION	R21E	2	N14	OSDH	1-11	VEG	2.88	MG/KG	
S-16	TION	R21E	2	NI	OSDH	1-11	VEG	3.12	MG/KG	
SHRES14		R21E	27	N6	KMTC	1-04	SOIL	80.0	PPM	SHRES
SMRES20	T12N	R21E	26	K6	KMTC	1-04	SOIL	40.0	PPM	SMRES
SMRES21	T12N	R21E	26	K6	KMTC	1-04	VEG	5.0	PPM	SMRES/LEAVES
STRM GH	TIIN	R21E	4	KO	OSDH	1-16	WATER	7.00	MG/L	
SWR PND	T11N	R21E	4	KO	OSDH	1-14	WATER	0.18	MG/L	
T3-08	T12N	R21E	26	C12	OSDH	1-14	SOIL	1.24	MG/KG	
T3-08	T12N	R21E	26	C12	OSDH	1-14	VEG	2.40	MG/KG	
T4-01E	T11N	R21E	2	C11	OSDH	1-06	SOIL	1.1	MG/KG	
T4-01E	T11N	R21E	2	C11	ORNL	1-06	SOIL	1.2	UG/G	DRY
T4-01E	T11N	R21E	2	C11	ORNL	1-06	SOIL	1.0	UG/G	WET
T4-01E	T11N	R21E	2	C11	OSDH	1-06	VEG	2.6	MG/KG	EVERGREEN
T4-01E	T11N	R21E	2	C11	ORNL	1-06	VEG	3.1	UG/G	
T4-01W	T12N	R21E	27	N1	OSDH	1-06	SOIL	1.1	MG/KG	
T4-01W	T12N	R21E	27	N1	ORNL	1-06	SOIL	0.7	UG/G	DRY
T4-01W	T12N	R21E	27	N1	ORNL	1-06	SOIL	0.6	UG/G	WET
T4-02E	TIIN	R21E	2	B8	OSDH	1-06	SOIL	<0.8	MG/KG	
T4-02E		R21E	2	B8	ORNL	1-06	SOIL	0.6	UG/G	DRY
		R21E		88	ORNL	1-06	SOIL	0.5	UG/G	WET
T4-02E	TIIN	R21E	2	88	OSDH	1-06	VEG	2.7	MG/KG	EVERGREEN
T4-02E	TIIN		2	B8	ORNL	1-06	VEG	2.1	UG/G	E.F.E.C.B.C.B.E.F.
T4-02E	TIIN	R21E			OSDH	1-06	SOIL	2.3	MG/KG	
T4-02W	T12N	R21E	33	L14		1-06	SOIL	0.7	UG/G	DRY
T4-02W	T12N	R21E	33	L14	ORNL			0.6	UG/G	WET
T4-02W	T12N	R21E	33	L14	ORNL	1-06	SOIL		MG/KG	WL 1
T4-03E	T12N	R21E	35	N5	OSDH	1-06	SOIL	<0.8		WET
T4-03E	T12N	R21E	35	N5	ORNL	1-06	SOIL	?	UG/G	DRY
T4-03E	T12N	R21E	35	N5	ORNL	1-06	SOIL	1.3	UG/G	
T4-03E	T12N	R21E	35	N5	OSDH	1-06	VEG	3.5	MG/KG	EVERGREEN
T4-03E	T12N	R21E	35	N5	ORNL	1-06	VEG	2.9	UG/G	
T4-03W	T12N	R21E	34	86	OSDH	1-06	SOIL	1.6	MG/KG	D D L
T4-03W	T12N	R21E	34	B6	ORNL	1-06	SOIL	0.8	UG/G	DRY
T4-03W	T12N	R21E	34	B6	ORNL	1-06	SOIL	0.7	UG/G	WET

FAR FIELD ENVIRONMENTAL SAMPLING RESULTS FOR FLUORIDE SAMPLES

T4-03W T12N R21E 34 C7 OSDH 1-06 WATER 0.17 MG/KG T4-04E T12N R21E 35 L3 OSDH 1-06 SOIL 0.9 UG/G DRY T4-04E T12N R21E 35 L3 ORNL 1-06 SOIL 0.7 UG/G WET T4-04E T12N R21E 35 L3 OSDH 1-06 VEG 1.8 MG/KG EVERGREEN T4-04E T12N R21E 34 L13 OSDH 1-06 SOIL 1.4 MG/KG EVERGREEN T4-04W T12N R21E 34 L13 OSDH 1-06 SOIL 0.8 UG/G DEY T4-04W T12N R21E 34 L13 OSDH 1-06 SOIL 0.8 UG/G DEY T4-04W T12N R21E 34 L13 OSDH 1-06 SOIL 0.2 UG/G WET T4-04D T12N R21E 35 K1 OSDH 1-06	LABEL	TOWNSHP	RANGE	SECTION	SUBSECT	LAB	DATE	MEDIA	RESULTS	UNITS A	COMMENTS
T4-04E T12N R21E 35 L3 ORNL 1-06 SOIL 0.9 UG/G DRY T4-04E T12N R21E 35 L3 ORNL 1-06 VEG 1.8 MG/KG EVERGREEN T4-04E T12N R21E 35 L3 ORNL 1-06 VEG 2.2 UG/G WET T4-04W T12N R21E 34 L13 OSDH 1-06 SOIL 1.4 MG/KG EVERGREEN T4-04W T12N R21E 34 L13 OSDH 1-06 SOIL 0.8 UG/G WET T4-04W T12N R21E 34 L13 ORNL 1-06 SOIL 0.8 UG/G WET T4-04W T12N R21E 34 L13 ORNL 1-06 SOIL 0.8 UG/G WET T4-05E T12N R21E 35 K1 OSNL 1-06 SOIL 1.2 MG/KG UE/KG T4-05E T12N R21E 35 K1 ORNL <t< td=""><td>T4-03W</td><td>T12N</td><td>R21E</td><td>34</td><td>C7</td><td>OSDH</td><td>1-06</td><td>WATER</td><td>0.17</td><td>MG/KG</td><td></td></t<>	T4-03W	T12N	R21E	34	C7	OSDH	1-06	WATER	0.17	MG/KG	
T4-04E T12N R21E 35 L3 ORNL 1-06 SOIL 0.9 UG/G DRY T4-04E T12N R21E 35 L3 ORNL 1-06 VEG 1.8 MG/KG EVERGREEN T4-04E T12N R21E 35 L3 ORNL 1-06 VEG 2.2 UG/G WET T4-04W T12N R21E 34 L13 OSDH 1-06 SOIL 1.4 MG/KG EVERGREEN T4-04W T12N R21E 34 L13 OSDH 1-06 SOIL 0.8 UG/G WET T4-04W T12N R21E 34 L13 ORNL 1-06 SOIL 0.8 UG/G WET T4-04W T12N R21E 34 L13 ORNL 1-06 SOIL 0.8 UG/G WET T4-05E T12N R21E 35 K1 OSNL 1-06 SOIL 1.2 MG/KG UE/KG T4-05E T12N R21E 35 K1 ORNL <t< td=""><td>T4-04E</td><td>T12N</td><td>R21E</td><td>35</td><td>L3</td><td>OSDH</td><td>1-06</td><td>SOIL</td><td><0.8</td><td>MG/KG</td><td></td></t<>	T4-04E	T12N	R21E	35	L3	OSDH	1-06	SOIL	<0.8	MG/KG	
T4-04E T12N R21E 35 L3 ORNL 1-06 SOIL 0.7 UG/G WET T4-04E T12N R21E 35 L3 OSDH 1-06 VEG 1.8 MG/KG EVERGREEN T4-04W T12N R21E 34 L13 OSDH 1-06 SOIL 1.4 MG/KG T4-04W T12N R21E 34 L13 OSNL 1-06 SOIL 1.4 MG/KG T4-04W T12N R21E 34 L13 OSNL 1-06 SOIL 0.8 UG/G WET T4-04W T12N R21E 34 L13 OSNL 1-06 VEG 54.0 MG/KG LEAVES T4-04W T12N R21E 35 K1 OSNL 1-06 SOIL 0.7 UG/G DRY T4-05E T12N R21E 35 K1 ORNL 1-06 SOIL 0.7 UG/G DRY T4-05E T12N R21E 35 K1 ORNL 1-06 VEG	T4-04E	T12N	R21F	35	13	ORNI	1-06	SOTI	0 9	UG/G	DRY
T4-04W T12N R21E 34 L13 ORNL 1-06 SOIL 1.0 UG/G DRY T4-04W T12N R21E 34 L13 ORNL 1-06 SOIL 0.8 UG/G WET T4-04W T12N R21E 34 L13 ORNL 1-06 VEG 54.0 MG/KG LEAVES T4-04W T12N R21E 35 K1 OSDH 1-06 SOIL 0.7 UG/G DRY T4-05E T12N R21E 35 K1 ORNL 1-06 SOIL 0.7 UG/G WET T4-05E T12N R21E 35 K1 ORNL 1-06 VEG 11.2 MG/KG EVERGREEN T4-05E T12N R21E 35 K1 ORNL 1-06 VEG 23.0 PPM TRES TRES53 T11N R21E 23 M10 KMTC 1-06 SOIL 1.40 MG/KG WR-01 T11N R22E 15 A1 OSDH 1-15	T4-04E	T12N	R21E	35	L3	ORNL	1-06	SOIL	0.7	UG/G	WET
T4-04W T12N R21E 34 L13 ORNL 1-06 SOIL 1.0 UG/G DRY T4-04W T12N R21E 34 L13 ORNL 1-06 SOIL 0.8 UG/G WET T4-04W T12N R21E 34 L13 ORNL 1-06 VEG 54.0 MG/KG LEAVES T4-04W T12N R21E 35 K1 OSDH 1-06 SOIL 0.7 UG/G DRY T4-05E T12N R21E 35 K1 ORNL 1-06 SOIL 0.7 UG/G WET T4-05E T12N R21E 35 K1 ORNL 1-06 VEG 11.2 MG/KG EVERGREEN T4-05E T12N R21E 35 K1 ORNL 1-06 VEG 23.0 PPM TRES TRES53 T11N R21E 23 M10 KMTC 1-06 SOIL 1.40 MG/KG WR-01 T11N R22E 15 A1 OSDH 1-15	T4-04E	T12N	R21E	35	L3	OSDH	1-06	VEG	1.8	MG/KG	EVERGREEN
T4-04W T12N R21E 34 L13 ORNL 1-06 SOIL 1.0 UG/G DRY T4-04W T12N R21E 34 L13 ORNL 1-06 SOIL 0.8 UG/G WET T4-04W T12N R21E 34 L13 ORNL 1-06 VEG 54.0 MG/KG LEAVES T4-04W T12N R21E 35 K1 OSDH 1-06 SOIL 0.7 UG/G DRY T4-05E T12N R21E 35 K1 ORNL 1-06 SOIL 0.7 UG/G WET T4-05E T12N R21E 35 K1 ORNL 1-06 VEG 11.2 MG/KG EVERGREEN T4-05E T12N R21E 35 K1 ORNL 1-06 VEG 23.0 PPM TRES TRES53 T11N R21E 23 M10 KMTC 1-06 SOIL 1.40 MG/KG WR-01 T11N R22E 15 A1 OSDH 1-15	T4-04E	T12N	R21E	35	13	ORNL	1-06	VEG	2.2	UG/G	
T4-04W T12N R21E 34 L13 ORNL 1-06 SOIL 1.0 UG/G DRY T4-04W T12N R21E 34 L13 ORNL 1-06 SOIL 0.8 UG/G WET T4-04W T12N R21E 34 L13 ORNL 1-06 VEG 54.0 MG/KG LEAVES T4-04W T12N R21E 35 K1 OSDH 1-06 SOIL 0.7 UG/G DRY T4-05E T12N R21E 35 K1 ORNL 1-06 SOIL 0.7 UG/G WET T4-05E T12N R21E 35 K1 ORNL 1-06 VEG 11.2 MG/KG EVERGREEN T4-05E T12N R21E 35 K1 ORNL 1-06 VEG 23.0 PPM TRES TRES53 T11N R21E 23 M10 KMTC 1-06 SOIL 1.40 MG/KG WR-01 T11N R22E 15 A1 OSDH 1-15	T4-04W	T12N	R21E	34	L13	OSDH	1-06	SOIL	1.4	MG/KG	
T4-04W T12N R21E 34 L13 ORNL 1-06 VEG 6.1 UG/G T4-05E T12N R21E 35 K1 OSDH 1-06 SOIL 1.2 MG/KG T4-05E T12N R21E 35 K1 ORNL 1-06 SOIL 0.7 UG/G DRY T4-05E T12N R21E 35 K1 ORNL 1-06 SOIL 0.6 UG/G WET T4-05E T12N R21E 35 K1 ORNL 1-06 VEG 11.2 MG/KG EVERGREEN T4-05E T12N R21E 23 M10 KMTC 1-06 VEG 23.0 PPM TRES TRES53 T11N R21E 23 M10 KMTC 1-06 SOIL 90.0 PPM TRES WR-01 T11N R22E 15 A1 OSDH 1-15 SOIL 1.40 MG/KG WR-01 T11N R22E 17 A14 OSDH 1-15 SOIL 1.64 MG	T4-04W	T12N	R21E	34	L13	ORNL	1-06	SOIL	1.0	UG/G	DRY
T4-04W T12N R21E 34 L13 ORNL 1-06 VEG 6.1 UG/G T4-05E T12N R21E 35 K1 OSDH 1-06 SOIL 1.2 MG/KG T4-05E T12N R21E 35 K1 ORNL 1-06 SOIL 0.7 UG/G DRY T4-05E T12N R21E 35 K1 ORNL 1-06 SOIL 0.6 UG/G WET T4-05E T12N R21E 35 K1 ORNL 1-06 VEG 11.2 MG/KG EVERGREEN T4-05E T12N R21E 23 M10 KMTC 1-06 VEG 23.0 PPM TRES TRES53 T11N R21E 23 M10 KMTC 1-06 SOIL 90.0 PPM TRES WR-01 T11N R22E 15 A1 OSDH 1-15 SOIL 1.40 MG/KG WR-01 T11N R22E 17 A14 OSDH 1-15 SOIL 1.64 MG		T12N	R21E	34	L13	ORNL	1-06	SOIL	0.8	UG/G	WET
T4-04W T12N R21E 34 L13 ORNL 1-06 VEG 6.1 UG/G T4-05E T12N R21E 35 K1 OSDH 1-06 SOIL 1.2 MG/KG T4-05E T12N R21E 35 K1 ORNL 1-06 SOIL 0.7 UG/G DRY T4-05E T12N R21E 35 K1 ORNL 1-06 SOIL 0.6 UG/G WET T4-05E T12N R21E 35 K1 ORNL 1-06 VEG 11.2 MG/KG EVERGREEN T4-05E T12N R21E 23 M10 KMTC 1-06 VEG 23.0 PPM TRES TRES53 T11N R21E 23 M10 KMTC 1-06 SOIL 90.0 PPM TRES WR-01 T11N R22E 15 A1 OSDH 1-15 SOIL 1.40 MG/KG WR-01 T11N R22E 17 A14 OSDH 1-15 SOIL 1.64 MG	T4-04W	T12N	R21E	34	L13	OSDH	1-06	VEG	54.0	MG/KG	LEAVES
WR-01 T11N R22E 15 A1 03DH 1-15 S01L 1.40 MG/KG WR-01 T11N R22E 15 A1 03DH 1-15 VEG 3.36 MG/KG WR-01 T11N R22E 9 N13 0SDH 1-15 WATER 0.18 MG/KG WR-02 T11N R22E 17 A14 0SDH 1-15 S01L 1.64 MG/KG WR-02 T11N R22E 17 A14 0SDH 1-15 WATER 0.26 MG/KG WR-02 T11N R22E 9 N2 0SDH 1-15 WATER 0.26 MG/KG WR-03 T11N R22E 17 114 0SDH 1-15 SOIL 1.12 MG/KG WR-03 T11N R22E 17 114 0SDH 1-15 VEG 3.60 MG/KG WR-03 T11N R22E 8 L14 0SDH 1-15 WATER 0.1 MG/L WR-05 T11N R22E <t< td=""><td>T4-04W</td><td>T12N</td><td>R21E</td><td>34</td><td>L13</td><td>ORNL</td><td>1-06</td><td>VEG</td><td>6.1</td><td>UG/G</td><td></td></t<>	T4-04W	T12N	R21E	34	L13	ORNL	1-06	VEG	6.1	UG/G	
WR-01 T11N R22E 15 A1 03DH 1-15 S01L 1.40 MG/KG WR-01 T11N R22E 15 A1 03DH 1-15 VEG 3.36 MG/KG WR-01 T11N R22E 9 N13 0SDH 1-15 WATER 0.18 MG/KG WR-02 T11N R22E 17 A14 0SDH 1-15 S01L 1.64 MG/KG WR-02 T11N R22E 17 A14 0SDH 1-15 WATER 0.26 MG/KG WR-02 T11N R22E 9 N2 0SDH 1-15 WATER 0.26 MG/KG WR-03 T11N R22E 17 114 0SDH 1-15 SOIL 1.12 MG/KG WR-03 T11N R22E 17 114 0SDH 1-15 VEG 3.60 MG/KG WR-03 T11N R22E 8 L14 0SDH 1-15 WATER 0.1 MG/L WR-05 T11N R22E <t< td=""><td>T4-05E</td><td>T12N</td><td>R21E</td><td>35</td><td>K1</td><td>OSDH</td><td>1-06</td><td>SOIL</td><td>1.2</td><td>MG/KG</td><td></td></t<>	T4-05E	T12N	R21E	35	K1	OSDH	1-06	SOIL	1.2	MG/KG	
WR-01 T11N R22E 15 A1 03DH 1-15 S01L 1.40 MG/KG WR-01 T11N R22E 15 A1 03DH 1-15 VEG 3.36 MG/KG WR-01 T11N R22E 9 N13 0SDH 1-15 WATER 0.18 MG/KG WR-02 T11N R22E 17 A14 0SDH 1-15 S01L 1.64 MG/KG WR-02 T11N R22E 17 A14 0SDH 1-15 WATER 0.26 MG/KG WR-02 T11N R22E 9 N2 0SDH 1-15 WATER 0.26 MG/KG WR-03 T11N R22E 17 114 0SDH 1-15 SOIL 1.12 MG/KG WR-03 T11N R22E 17 114 0SDH 1-15 VEG 3.60 MG/KG WR-03 T11N R22E 8 L14 0SDH 1-15 WATER 0.1 MG/L WR-05 T11N R22E <t< td=""><td>T4-05E</td><td>T12N</td><td>R21E</td><td>35</td><td>K1</td><td>ORNL</td><td>1-05</td><td>SOIL</td><td>0.7</td><td>UG/G</td><td>DRY</td></t<>	T4-05E	T12N	R21E	35	K1	ORNL	1-05	SOIL	0.7	UG/G	DRY
WR-01 T11N R22E 15 A1 03DH 1-15 S01L 1.40 MG/KG WR-01 T11N R22E 15 A1 03DH 1-15 VEG 3.36 MG/KG WR-01 T11N R22E 9 N13 0SDH 1-15 WATER 0.18 MG/KG WR-02 T11N R22E 17 A14 0SDH 1-15 S01L 1.64 MG/KG WR-02 T11N R22E 17 A14 0SDH 1-15 WATER 0.26 MG/KG WR-02 T11N R22E 9 N2 0SDH 1-15 WATER 0.26 MG/KG WR-03 T11N R22E 17 114 0SDH 1-15 SOIL 1.12 MG/KG WR-03 T11N R22E 17 114 0SDH 1-15 VEG 3.60 MG/KG WR-03 T11N R22E 8 L14 0SDH 1-15 WATER 0.1 MG/L WR-05 T11N R22E <t< td=""><td>T4-05E</td><td>T12N</td><td>R21E</td><td>35</td><td>K1</td><td>ORNL</td><td>1-06</td><td>SOIL</td><td>0.6</td><td>UG/G</td><td>WET</td></t<>	T4-05E	T12N	R21E	35	K1	ORNL	1-06	SOIL	0.6	UG/G	WET
WR-01 T11N R22E 15 A1 03DH 1-15 S01L 1.40 MG/KG WR-01 T11N R22E 15 A1 03DH 1-15 VEG 3.36 MG/KG WR-01 T11N R22E 9 N13 0SDH 1-15 WATER 0.18 MG/KG WR-02 T11N R22E 17 A14 0SDH 1-15 S01L 1.64 MG/KG WR-02 T11N R22E 17 A14 0SDH 1-15 WATER 0.26 MG/KG WR-02 T11N R22E 9 N2 0SDH 1-15 WATER 0.26 MG/KG WR-03 T11N R22E 17 114 0SDH 1-15 SOIL 1.12 MG/KG WR-03 T11N R22E 17 114 0SDH 1-15 VEG 3.60 MG/KG WR-03 T11N R22E 8 L14 0SDH 1-15 WATER 0.1 MG/L WR-05 T11N R22E <t< td=""><td>T4-05E</td><td>T12N</td><td>R21E</td><td>35</td><td>K1</td><td>OSDH</td><td>1-06</td><td>VEG</td><td>11.2</td><td>MG/KG</td><td>EVERGREEN</td></t<>	T4-05E	T12N	R21E	35	K1	OSDH	1-06	VEG	11.2	MG/KG	EVERGREEN
WR-01 T11N R22E 15 A1 03DH 1-15 S01L 1.40 MG/KG WR-01 T11N R22E 15 A1 03DH 1-15 VEG 3.36 MG/KG WR-01 T11N R22E 9 N13 0SDH 1-15 WATER 0.18 MG/KG WR-02 T11N R22E 17 A14 0SDH 1-15 S01L 1.64 MG/KG WR-02 T11N R22E 17 A14 0SDH 1-15 WATER 0.26 MG/KG WR-02 T11N R22E 9 N2 0SDH 1-15 WATER 0.26 MG/KG WR-03 T11N R22E 17 114 0SDH 1-15 SOIL 1.12 MG/KG WR-03 T11N R22E 17 114 0SDH 1-15 VEG 3.60 MG/KG WR-03 T11N R22E 8 L14 0SDH 1-15 WATER 0.1 MG/L WR-05 T11N R22E <t< td=""><td>T4-05E</td><td>T12N</td><td>R21E</td><td>35</td><td>K1</td><td>ORNL</td><td>1-06</td><td>VEG</td><td>9.3</td><td>UG/G</td><td></td></t<>	T4-05E	T12N	R21E	35	K1	ORNL	1-06	VEG	9.3	UG/G	
WR-01 T11N R22E 15 A1 03DH 1-15 S01L 1.40 MG/KG WR-01 T11N R22E 15 A1 03DH 1-15 VEG 3.36 MG/KG WR-01 T11N R22E 9 N13 0SDH 1-15 WATER 0.18 MG/KG WR-02 T11N R22E 17 A14 0SDH 1-15 S01L 1.64 MG/KG WR-02 T11N R22E 17 A14 0SDH 1-15 WATER 0.26 MG/KG WR-02 T11N R22E 9 N2 0SDH 1-15 WATER 0.26 MG/KG WR-03 T11N R22E 17 114 0SDH 1-15 SOIL 1.12 MG/KG WR-03 T11N R22E 17 114 0SDH 1-15 VEG 3.60 MG/KG WR-03 T11N R22E 8 L14 0SDH 1-15 WATER 0.1 MG/L WR-05 T11N R22E <t< td=""><td>TRES53</td><td>T11N</td><td>R21E</td><td>23</td><td>M10</td><td>KMTC</td><td>1-06</td><td>VEG</td><td>23.0</td><td>PPM</td><td>TRES</td></t<>	TRES53	T11N	R21E	23	M10	KMTC	1-06	VEG	23.0	PPM	TRES
WR-01 T11N R22E 15 A1 03DH 1-15 S01L 1.40 MG/KG WR-01 T11N R22E 15 A1 03DH 1-15 VEG 3.36 MG/KG WR-01 T11N R22E 9 N13 0SDH 1-15 WATER 0.18 MG/KG WR-02 T11N R22E 17 A14 0SDH 1-15 S01L 1.64 MG/KG WR-02 T11N R22E 17 A14 0SDH 1-15 WATER 0.26 MG/KG WR-02 T11N R22E 9 N2 0SDH 1-15 WATER 0.26 MG/KG WR-03 T11N R22E 17 114 0SDH 1-15 SOIL 1.12 MG/KG WR-03 T11N R22E 17 114 0SDH 1-15 VEG 3.60 MG/KG WR-03 T11N R22E 8 L14 0SDH 1-15 WATER 0.1 MG/L WR-05 T11N R22E <t< td=""><td>TRES54</td><td>TIIN</td><td>R21E</td><td>23</td><td>M10</td><td>KMTC</td><td>1-06</td><td>SOIL</td><td>90.0</td><td>PPM</td><td>TRES</td></t<>	TRES54	TIIN	R21E	23	M10	KMTC	1-06	SOIL	90.0	PPM	TRES
WR-02 111N R22E 17 A14 OSDH 1-15 VEG 5.36 MG/KG WR-02 T11N R22E 9 N2 OSDH 1-15 WATER 0.26 MG/L WR-03 T11N R22E 17 I14 OSDH 1-15 SOIL 1.12 MG/KG WR-03 T11N R22E 17 I14 OSDH 1-15 VEG 3.60 MG/KG WR-03 T11N R22E 17 I14 OSDH 1-15 VEG 3.60 MG/KG WR-04 T11N R22E 8 L14 OSDH 1-15 WATER 0.1 MG/L WR-05 T11N R22E 9 A1 OSDH 1-15 SOIL 1.28 MG/KG WR-05 T11N R22E 9 A1 OSDH 1-15 VEG 3.28 MG/KG WR-06 T11N R22E 5 N3 OSDH 1-15 SOIL <0.40	WR-01	TIIN	R22E	15	A1	USUN	1 13	2011	1.40	PRU/ NU	
WR-02 111N R22E 17 A14 OSDH 1-15 VEG 5.36 MG/KG WR-02 T11N R22E 9 N2 OSDH 1-15 WATER 0.26 MG/L WR-03 T11N R22E 17 I14 OSDH 1-15 SOIL 1.12 MG/KG WR-03 T11N R22E 17 I14 OSDH 1-15 VEG 3.60 MG/KG WR-03 T11N R22E 17 I14 OSDH 1-15 VEG 3.60 MG/KG WR-04 T11N R22E 8 L14 OSDH 1-15 WATER 0.1 MG/L WR-05 T11N R22E 9 A1 OSDH 1-15 SOIL 1.28 MG/KG WR-05 T11N R22E 9 A1 OSDH 1-15 VEG 3.28 MG/KG WR-06 T11N R22E 5 N3 OSDH 1-15 SOIL <0.40	WR-01	T11N	R22E	15	A1	OSDH	1-15	VEG	3.36	MG/KG	
WR-02 111N R22E 17 A14 OSDH 1-15 VEG 5.36 MG/KG WR-02 T11N R22E 9 N2 OSDH 1-15 WATER 0.26 MG/L WR-03 T11N R22E 17 I14 OSDH 1-15 SOIL 1.12 MG/KG WR-03 T11N R22E 17 I14 OSDH 1-15 VEG 3.60 MG/KG WR-03 T11N R22E 17 I14 OSDH 1-15 VEG 3.60 MG/KG WR-04 T11N R22E 8 L14 OSDH 1-15 WATER 0.1 MG/L WR-05 T11N R22E 9 A1 OSDH 1-15 SOIL 1.28 MG/KG WR-05 T11N R22E 9 A1 OSDH 1-15 VEG 3.28 MG/KG WR-06 T11N R22E 5 N3 OSDH 1-15 SOIL <0.40	WR-01	TIIN	R22E	9	N13	OSDH	1-15	WATER	0.18	MG/L	
WR-02 111N R22E 17 A14 OSDH 1-15 VEG 5.36 MG/KG WR-02 T11N R22E 9 N2 OSDH 1-15 WATER 0.26 MG/L WR-03 T11N R22E 17 I14 OSDH 1-15 SOIL 1.12 MG/KG WR-03 T11N R22E 17 I14 OSDH 1-15 VEG 3.60 MG/KG WR-03 T11N R22E 17 I14 OSDH 1-15 VEG 3.60 MG/KG WR-04 T11N R22E 8 L14 OSDH 1-15 WATER 0.1 MG/L WR-05 T11N R22E 9 A1 OSDH 1-15 SOIL 1.28 MG/KG WR-05 T11N R22E 9 A1 OSDH 1-15 VEG 3.28 MG/KG WR-06 T11N R22E 5 N3 OSDH 1-15 SOIL <0.40	WR-02	TIIN	R22E	17	A14	OSDH	1-15	SOIL	1.64	MG/KG	
WR-08 T11N R22E 6 N7 0SDH 1-15 SOIL 1.08 MG/KG WR-08 T11N R22E 6 N7 0SDH 1-15 SOIL 1.08 MG/KG WR-08 T11N R22E 6 N7 0SDH 1-15 VEG 2.08 MG/KG WR-09 T11N R22E 6 G1 0SDH 1-15 SOIL 0.52 MG/KG	WR-02	T11N	R22E	17	A14	OSDH	1-15	VEG	5.36	MG/KG	
WR-08 T11N R22E 6 N7 0SDH 1-15 SOIL 1.08 MG/KG WR-08 T11N R22E 6 N7 0SDH 1-15 SOIL 1.08 MG/KG WR-08 T11N R22E 6 N7 0SDH 1-15 VEG 2.08 MG/KG WR-09 T11N R22E 6 G1 0SDH 1-15 SOIL 0.52 MG/KG		TIIN	R22E	9	N2	OSDH	1-15	WATER	0.26	MG/L	
WR-08 T11N R22E 6 N7 0SDH 1-15 SOIL 1.08 MG/KG WR-08 T11N R22E 6 N7 0SDH 1-15 SOIL 1.08 MG/KG WR-08 T11N R22E 6 N7 0SDH 1-15 VEG 2.08 MG/KG WR-09 T11N R22E 6 G1 0SDH 1-15 SOIL 0.52 MG/KG		TIIN	R22E	17	114	OSDH	1-15	SOIL	1.12	MG/KG	
WR-08 T11N R22E 6 N7 0SDH 1-15 SOIL 1.08 MG/KG WR-08 T11N R22E 6 N7 0SDH 1-15 SOIL 1.08 MG/KG WR-08 T11N R22E 6 N7 0SDH 1-15 VEG 2.08 MG/KG WR-09 T11N R22E 6 G1 0SDH 1-15 SOIL 0.52 MG/KG		TIIN	R22E	17	I14	OSDH	1-15	VEG	3.60	MG/KG	
WR-08 T11N R22E 6 N7 0SDH 1-15 SOIL 1.08 MG/KG WR-08 T11N R22E 6 N7 0SDH 1-15 SOIL 1.08 MG/KG WR-08 T11N R22E 6 N7 0SDH 1-15 VEG 2.08 MG/KG WR-09 T11N R22E 6 G1 0SDH 1-15 SOIL 0.52 MG/KG		TIIN	R22E	8	L14	OSDH	1-15	WATER	0.1	MG/L	
WR-08 T11N R22E 6 N7 0SDH 1-15 SOIL 1.08 MG/KG WR-08 T11N R22E 6 N7 0SDH 1-15 SOIL 1.08 MG/KG WR-08 T11N R22E 6 N7 0SDH 1-15 VEG 2.08 MG/KG WR-09 T11N R22E 6 G1 0SDH 1-15 SOIL 0.52 MG/KG		TIIN	R22E	9	A1	OSDH	1-15	SOIL	1.28	MG/KG	
WR-08 T11N R22E 6 N7 0SDH 1-15 SOIL 1.08 MG/KG WR-08 T11N R22E 6 N7 0SDH 1-15 SOIL 1.08 MG/KG WR-08 T11N R22E 6 N7 0SDH 1-15 VEG 2.08 MG/KG WR-09 T11N R22E 6 G1 0SDH 1-15 SOIL 0.52 MG/KG		TIIN	R22E	9	Al	OSDH	1-15	VEG	3.28	MG/KG	
WR-08 T11N R22E 6 N7 0SDH 1-15 SOIL 1.08 MG/KG WR-08 T11N R22E 6 N7 0SDH 1-15 SOIL 1.08 MG/KG WR-08 T11N R22E 6 N7 0SDH 1-15 VEG 2.08 MG/KG WR-09 T11N R22E 6 G1 0SDH 1-15 SOIL 0.52 MG/KG		TIIN	R22E	5	N3	OSDH	1-15	WATER	0.26	MG/L	
WR-08 T11N R22E 6 N7 0SDH 1-15 SOIL 1.08 MG/KG WR-08 T11N R22E 6 N7 0SDH 1-15 SOIL 1.08 MG/KG WR-08 T11N R22E 6 N7 0SDH 1-15 VEG 2.08 MG/KG WR-09 T11N R22E 6 G1 0SDH 1-15 SOIL 0.52 MG/KG		TIIN	R22E	6	N14	OSDH	1-15	SOIL	<0.40	MG/KG	
WR-09 T11N R22E 6 G1 OSDH 1-15 SOTL 0.52 MG/KG		TIIN	R22E		N14	OSDH	1-15	VEG	1.68	rna/ nu	
WR-09 T11N R22E 6 G1 OSDH 1-15 SOTL 0.52 MG/KG		TIIN	R22E	6	N7	OSDH	1-15	SOIL	1.08		
WR-09 T11N R22E 6 G1 OSDH 1-15 SOIL 0.52 MG/KG WR-09 T11N R22E 6 G1 OSDH 1-15 VEG 1.84 MG/KG WR-10 T11N R22E 17 A1 OSDH 1-15 SOIL <0.40 MG/KG WR-10 T11N R22E 17 A1 OSDH 1-15 VEG 2.32 MG/KG YRES51 T11N R21E 26 N3 KMTC 1-05 VEG 2600.0 PPM YRES YRES52 T11N R21E 26 N3 KMTC 1-05 SOIL 100.0 PPM YRES		TIIN	R22E	6	N7	OSDH	1-15	VEG	2.08		
WR-09 T11N R22E 6 G1 OSDH 1-15 VEG 1.84 MG/KG WR-10 T11N R22E 17 A1 OSDH 1-15 SOIL <0.40		TIIN	R22E	6	Gl	OSDH	1-15	SOIL	0.52	MG/KG	
WR-10 T11N R22E 17 A1 OSDH 1-15 SOIL <0.40 MG/KG WR-10 T11N R22E 17 A1 OSDH 1-15 VEG 2.32 MG/KG YRES51 T11N R21E 26 N3 KMTC 1-05 VEG 2600.0 PPM YRES YRES52 T11N R21E 26 N3 KMTC 1-05 SOIL 100.0 PPM YRES		TIIN	RZZE	6	Gl	OSDH	1-15	VEG	1.84	MG/KG	
WR-10 T11N R22E 17 A1 OSDH 1-15 VEG 2.32 MG/KG YRES51 T11N R21E 26 N3 KMTC 1-05 VEG 2600.0 PPM YRES YRES52 T11N R21E 26 N3 KMTC 1-05 S01L 100.0 PPM YRES	WR-10	TIIN	R22E	17	Al	OSDH	1-15	SOIL	<0.40	MG/KG	
YRES51 TIIN R21E 26 N3 KMTC 1-05 VEG 2600.0 PPM YRES YRES52 TIIN R21E 26 N3 KMTC 1-05 S01L 100.0 PPM YRES	WR-10	TIIN	RZZE	17	Al	OSDH	1-15	VEG	2.32	MG/KG	
YRESSZ TIIN RZIE 26 N3 KMTC 1-05 SOIL 100.0 PPM YRES	YRES51	TIIN	R21E	26	N3	KMTC	1-05	VEG	2600.0	PPM	YRES
	YRES52	TIIN	RZIE	26	N3	KMTC	1-05	SOIL	100.0	MAA	YRES

]

· .

۱

FAR FIELD ENVIRONMENTAL SAMPLING RESULTS FOR URANIUM SAMPLES

LABEL	TOWNSHP	RANGE	SECT	ION SUBSE	CT LAB	DATE	MEDIA	RESUL	TS UNITS	COMMENTS
16W	TIIN	R21E	34	86	KMTC	1-06	WATER	0.006	MG/L	
1H14	TIIN	R21E	1	H14	KMTC	1-10	SOIL	<10	PPM	
2E14	TIIN	R21E	2	E14	KMTC	1-10	SOIL	<10	PPM	
4A 6	TIIN	R21E	4	A6	KMTC	1-10	SOIL	<10	PPM	
4J14	TIIN	R21E	4	J14	KMTC	1-09	SOIL	<10	PPM	
5N 6	TIIN	R22E	5	NG	KMTC	1-09	SOIL	<10	PPM	
6A 7	TIIN	R22E	6	A7	KMTC	1-09	SOIL	< 10	PPM	
6C14	T11N	R22E	6	C14	KMTC	1-09	SOIL	<10	PPM	
6111	TIIN	R22E	6	I11	KMTC	1-09	SOIL	<10	PPM	
70 6	TIIN	R22E	7	C6	KMTC	1-08	SOIL	<10	PPM	
70 6	TIIN	R22E	7	C6	KMTC	1-08	VEG	1.2	PPM	J. GRASS
7F 7	TIIN	R22E	7	F7	KMTC	1-08	SOIL	<10	PPM	
7F 7	TIIN	R22E	7	F7	KMTC	1-08	VEG	3.9	PPM	WHEAT
7F 7	TIIN	R22E	7	F7	KMTC	1-08	VEG	3.5	PPM	WITTER TT
8A 1	TIIN	R22E	8	Al	KMTC	1-09	SOIL	< 10	PPM	
8J 7	TIIN	R21E	8	J7	KMTC	1-09	SOIL	<10	PPM	
8N14	TIIN	R22E	8	M14	KMTC	1-09	SOIL	20.0	PPM	
9B 2	TIIN	R21E	9	B2	KMTC	1-09	SOIL	<10	PPM	
91 8	TIIN	R21E	9	18	KMTC	1-09	SOIL	15.0	PPM	
	TIIN				KMTC	1-09	SOIL	<10	PPM	
108 8		R21E	10	88	KMTC	1-09	SOIL	<10	PPM	
10L 8	TIIN	R21E	10	LS		1-08	SOIL	<10	PPM	
10M10	TIIN	R21E	10	M10	KMTC			<10	PPM	
11F 7	TIIN	R21E	11	F7	KMTC	1-08	SOIL		PPM	CAT LEAVES
11F 7	TIIN	R21E	11	F7	KMTC	1-08	VEG	9.3		BEANS
11F 7	TIIN	R21E	11	F7	KMTC	1-08	VEG	31.0	PPM	
11F 7	TIIN	R21E	11	F7	KMTC	1-08	VEG	18.0	PPM	WILLOW
11N 6	T11N	R21E	11	N6	KMTC	1-08	SOIL	<10	PPM	
12F 7	T11N	R21E	12	F7	KMTC	1-09	SOIL	<10	PPM	
12N 8	TIIN	R21E	12	N8	KMTC	1-08	SOIL	< 10	PPM	00040
12N 8	TIIN	R21E	12	N8	KMTC	1-08	VEG	5.9	PPM	CEDAR
12N 8	T11N	R21E	12	N8	KMTC	1-08	VEG	3.7	PPM	CAT TIPS
12N 8	TIIN	R21E	12	N8	KMTC	1-08	VEG	2.6	PPM	
13B 1	T11N	R21E	13	81	KMTC	1-08	SOIL	15.0	PPM	
13B 1	T11N	R21E	13	81	KMTC	1-08	VEG	2.2	PPM	
13G11		R21E	13	G11	KMTC	1-09		15.0	PPM	
13G11	TIIN	R21E	13	G11	KMTC	1-09	VEG	1.9	PPM	
13H 2	TIIN	R21E	13	H2	KMTC	1-08	SOIL	15.0	PPM	
148 ;	TIIN	R21E	14	B4	KMTC	1-08	SOIL	< 10	PPM	
148 4	TIIN	R21E	14	B4	KMTC	1-08	VEG	1.5	PPM	
148 8	TIIN	R21E	14	88	KMTC	1-08	SOIL	15.0	PPM	
14B 8	T11N	R21E	14	88	KMTC	1-08	VEG	0.6	PPM	
14B12	T11N	R21E	14	B12	KMTC	1-08	SOIL	15.0	PPM	
14812	TIIN	R21E	14	B12	KMTC	1-08	VEG	8.1	PPM	FESCUE
14812	TIIN	R21E	14	B12	KMTC	1-08	VEG	7.8	PPM	FESCUE TIP
16810	TIIN	R21E	16	B10	KMTC	1-09	SOIL	<10	PPM	
16B10	T11N	R21E	16	810	KMTC	1-09	VEG	0.6	PPM	
16L 2	T11N	R21E	16	L2	KMTC	1-10	SOIL	<10	PPM	
16L12	TIIN	R21E	16	L12	KMTC	1-10	SOIL	< 10	PPM	
16L12	TIIN	R21E	16	L12	KMTC	1-10	VEG	1.1	PPM	

.

FAR FIELD ENVIRONMENTAL SAMPLING RESULTS FOR URANIUM SAMPLES

۱

١

LABEL	TOWNSHP	RANGE	SECTION	SUBSECT	LAB	DATE	MEDIA	RESULTS	UNITS	COMMENTS	
18B 4	TIIN	R22E	18	84	KMTC	1-09	SOIL	< 10	PPM		
18B 4	TIIN	R22E	18	84	KMTC	1-09	VEG	0.7	PPM		
18H10	TIIN	R22E	18	H10	KMTC	1-09	SOIL	<10	PPM		
18H10	TIIN	R22E	18	H10	KMTC	1-09	VEG	<0.5	PPM		
1A11	TIIN	R21E	1	A11		1-09	SOIL	<10	PPM		
1A11	TIIN	R21E	1	A11	KMTC	1-09	VEG	0.7	PPM		
1H14	TIIN	R21E	î	H14		1-08	VEG	1.4	PPM		
22G14	TIIN	R21E	22	G14		1-10	SOIL	<10	PPM		
22G14	TIIN	R21E	22	G14	KMTC	1-10	VEG	<0.5	PPM		
23A10	TIIN	R21E	23	A10		1-08	SOIL	10.0	PPM	T11N	
23A10	TIIN	R21E	23	A10		1-08	VEG	28.0	PPM	TIIN	
25E 8	T12N	R21E	25	E8		1-10	SOIL	15.0	PPM	1 2 2 1	
25E 8	T12N	R21E	25	E8		1-10	VEG	1.2	PPM		
25H 3	T12N	R21E	25	H3		1-10	SOIL	15.0	PPM		
25H 3	T12N	R21E	25	H3		1-10	VEG	<0.5	PPM		
25H 3	T12N	R21E	25	H3		1-10	VEG -	<0.5	PPM		
25H14	T12N	R21E	25	H14		1-10	SOIL	15.0	PPM		
25L14	T12N	R21E	25	L14		1-14	WATER	<0.002	MG/L		
26E 7	T12N	R21E	26	E7	KMTC	1-14	WATER	0.006	MG/L		
26F 6	T12N	R21E	26	F6		1-14	WATER	0.003	MG/L	SURFACE	WATER
26G 1	T12N	R21E	26	G1		1-10	SOIL	15.0	PPM	JUNIAUL	MATER
26G 1	T12N	R21E	26	GI	KMTC	1-10	VEG	<0.5	PPM		
26G 8		R21E	26	G8		1-10	SOIL	10.0	PPM		
	T12N			G8	KMTC	1-10	VEG	1.1	PPM		
26G 8	T12N	R21E	26					0.003	MG/L		
26L 1	T12N	R21E	26	L1	KMTC	1-14	WATER	10.0	PPM		
27A13	T12N	R21E	27	A13	KMTC	1-08	SOIL	4.5	PPM	SMILAX	
27A13	T12N	R21E	27	A13	KMTC	1-08	VEG	5.7	PPM	CEDAR	
27A13	T12N	R21E	27	A13	KMTC	1-08	VEG		PPM	CEDAR	
270 1	T12N	R21E	27	01	KMTC	1-10	SOIL	15.0	PPM		
270 1	T12N	R21E	27	D1	KMTC	1-10	VEG	52.0	PPM		
27H 9	T12N	R21E	27	H9		1-10	SOIL	< 10			
27H 9	T12N	R21E	27	H9		1-10	VEG	2.3	PPM		
27H 9	T12N	R21E	27	H9		1-10	VEG	1.1	PPM	CUDEACE	WATER
271 0	T12N	R21E	27	10	KMTC	1-14	WATER	<0.002	MG/L	SURFACE	WATER
		R21E				1-10					
27K 4	T12N	R21E	27	K4	KMTC	1-10	VEG	38.0	PPM		
27L 1	T12N	R21E	27	Ll	KMTC	1-14	WATER	<0.002	MG/L	CUDEACE	WATER
27L 4	T12N	R21E	27	L4	KMTC	1-14	WATER	<0.002	MG/L	SURFACE	
27L 8	T12N	R21E	27	L8	KMTC	1-14	WATER	0.004	MG/L	SURFACE	WATER
28G 4	T12N	R21E	28	G4	KMTC	1-10	SOIL	<10	PPM		
28G 4	T12N	R21E	28	G4	KMTC	1-10	VEG	1.2	PPM		
28G 4	T12N	R21E	28	G4	KMTC	1-10	VEG	<0.5	PPM	CUDENOS	UNTER
28H14	T12N	R21E	28	H14	KMTC	1-14	WATER	<0.002	MG/L	SURFACE	WATER
28J11	T12N	R21E	28	J11	KMTC	1-10	SOIL	< 10	PPM		
28J11	T12N	R21E	28	J11	KMTC	1-10	VEG	1.7	PPM		
28N 2	T12N	R21E	28	N2	KMTC	1-10	SOIL	< 10	PPM		
28N 2	T12N	R21E	28	N2	KMTC	1-10	VEG	<0.5	PPM		
28N 2	T12N	R21E	28	N2	KMTC	1-10	VEG	6.6	PPM		
290 3	T12N	R22E	29	D3	KMTC	1-10	SOIL	15.0	PPM		

FAR FIELD ENVIRONMENTAL SAMPLING RESULTS FOR URANIUM SAMPLES

LABEL	TOWNSHP	RANGE	SECTION	SUBSECT	LAB	DATE	MEDIA	RESULTS	UNITS	COMMENTS
29D 3	T12N	R22E	29	D3	KMTC	1-10	VEG	<0.5	PPM	
2E14	TIIN	R21E	2	E14	KMTC	1-09	VEG	<0.5	PPM	
2L13	TIIN	R21E	2	L13	KMTC	1-09	SOIL	<10	PPM	
2L13	TIIN	R21E	2	L13	KMTC	1-09	SOIL	<10	PPM	
2L13	TIIN	R21E	2	L13	KMTC	1-09	VEG	5.5	PPM	
2M 8	TIIN	R21E	2	M8	KMTC	1-09	SOIL	<10	PPM	
2M 8	TIIN	R21E	2	M8	KMTC	1-09	VEG	55.0	PPM	
30F 7	TI2N	R22E	30	F7	KMTC	1-10	SOIL	15.0	PPM	
30F 7	T12N	R22E	30	F7	KMTC	1-10	VEG	<0.5	PPM	
		R22E	30	F7	KMTC	1-10	VEG	0.5	PPM	
30F 7	T12N				KMTC	1-10	SOIL	<10	PPM	
31B10	T12N	R22E	31	B10			VEG	<0.5	PPM	
31B10	T12N	R22E	31	B10	KMTC	1-10			PPM	
33814	T12N	R21E	33	B14	KMTC	1-10	SOIL	<10		
33B14	T12N	R21E	33	B14	KMTC	1-10	VEG	<0.5	PPM	
33I 3	T12N	R21E	33	13	KMTC	1-10	SOIL	<10	PPM	
33I 3	T12N	R21E	33	13	KMTC	1-10	VEG	<0.5	PPM	
33M14	T12N	R21E	33	M14	KMTC	1-10	SOIL	<10	PPM	
33M14	T12N	R21E	33	M14	KMTC	1-10	VEG	<0.5	PPM	
33M14	T12N	R21E	33	M14	KMTC	1-10	VEG	<0.5	PPM	
34D 7	T12N	R21E	34	D7	KMTC	1-10	SOIL	<10	PPM	
34D 7	T12N	R21E	34	D7	KMTC	1-10	VEG	16.0	PPM	
35A 9	T12N	R21E	35	A9	KMTC	1-10	SOIL	15.0	PPM	
35A 9	T12N	R21E	35	A9	KMTC	1-10	VEG	<0.5	PPM	
35A 9	T12N	R21E	35	A9	KMTC	1-10	VEG	0.8	PPM	
35B 1	T12N	R21E	35	B1	KMTC	1-10	SOIL	10.0	PPM	
35B 1	T12N	R21E	35	B1	KMTC	1-10	VEG	<0.5	PPM	
35B14	T12N	R21E	35	B14	KMTC	1-10	SOIL	10.0	PPM	
35B14	T12N	R21E	35	B14	KMTC	1-10	VEG	<0.5	PPM	
35F10	T12N	R21E	35	F10	KMTC	1-10	SOIL	15.0	PPM	
35F10	T12N	R21E	35	F10	KMTC	1-10	VEG	0.6	PPM	
35L 5	T12N	R21E	35	L5	KMTC	1-09	SOIL	10.0	PPM	
	T12N	R21E	35	L5	KMTC	1-09	VEG	<0.5	PPM	
35L 5			35	L5	KMTC	1-09	VEG	1.0	PPM	
35L 5	T12N	R21E		L11	KMTC	1-10	SOIL	<10	PPM	
35L11	T12N	R21E	35					<0.5	PPM	
35L11		P.21E	35	L11	KMTC			1.2	PPM	
35L11	T12N	R21E	35	L11	KMTC	1-10	VEG		PPM	
36A13	T12N	R21E	36	A13	KMTC	1-10	SOIL	<10	PPM	
36A13	T12N	R21E	36	A13	KMTC	1-10	VEG	<0.5		
36D11	T12N	R21E	36	D11	KMTC	1-14	WATER	0.004	MG/L	
36G 5	T12N	R21E	36	G5	KMTC	1-09	SOIL	<10	PPM	
36G 5	T12N	R21E	36	G5	KMTC	1-09	VEG	0.9	PPM	
36H14	T12N	R21E	36	H14	KMTC	1-10	SOIL	<10	PPM	
36H14	T12N	R21E	36	H14	KMTC	1-10	VEG	0.5	PPM	
36L 5	T12N	R21E	36	L5	KMTC	1-14	WATER	0.003	MG/L	
36M10	T12N	R21E	36	M10	KMTC	1-14	WATER	<0.002	MG/L	
3F 7	TIIN	R21E	3	F7	KMTC	1-09	SOIL	<10	PPM	
3F 7	TIIN	R21E	3	F7	KMTC	1-09	VEG	0.6	PPM	
4A 1	TIIN	R21E	4	Al	KMTC	1-10	SOIL	<10	PPM	
4A 1	TIIN	R21E	4	Al	KMTC	1-10	VEG	<0.5	PPM	
4A 3	TIIN	R21E	4	A3	KMTC	1-10	VEG	<0.5	PPM	
40.0	TTH	NE de la	-	110	1					

¥

FAR FIELD ENVIRONMENTAL SAMPLING RESULTS FOR URANIUM SAMPLES

											1.0
LABEL	TOWNSHP	RANGE	SECTION	SUBSECT	LAB	DATE	MEDIA	RESULTS	UNITS	COMMEN	TS
4A 6	TIIN	R21E	4	A6	KMTC	1-10	VEG	<0.5	PPM		
4J14	T11N	R21E	4	J14	KMTC	1-09	VEG	1.0	PPM		1. A
4J14	TIIN	R21E	4	J14	KMTC	1-09	VEG	<0.5	PPM		- F. 13
4N 7	T11N	R21E	4	N7	KMTC	1-09	SOIL	<10	PPM		
4N 7	TIIN	R21E	4	N7	KMTC	1-09	VEG	<0.5	PPM		
4N 8	TIIN	R21E	4	N8	KMTC	1-09	VEG	<0.5	PPM		
50/HWY	T12N	R21E	22	G1	KMTC	1-06	VEG	17000.0	PPM		1.1
5N 6	T11N	R22E	5	N6	KMTC	1-09	VEG	0.5	PPM		
6A 7	T11N	R22E	6	A7	KMTC	1-09	VEG	<0.5	PPM		
6C14	TIIN	R22E	6	C14	KMTC	1-09	VEG	1.8	PPM		
8A 1	T11N	R22E	8	A1	KMTC	1-09	VEG	1.7	PPM		
8J 7	T11N	R21E	8	J7	KMTC	1-09	VEG	1.1	PPM		
8N14	T11N	R22E	8	N14	KMTC	1-09	VEG	<0.5	PPM		1.1
9B 2	T11N	R21E	8	82	KMTC	1-09	VEG	1.7	PPM		1.10
	T12N	R21E	21	B13	ORNL	1-06	SOIL	7.88	PPM	DRY	1.1
	T12N	R21E	21	B13	ORNL	1-06	SOIL	6.1	PPM	WET	1.1
	T12N	R21E	21	B13	KMTC	1-06	SOIL	<10	PPM		1111
	T12N	R21E	21	B13	ORNL	1-06	VEG	3.91	PPM		1.11
	T12N	R21E	21	B13	KMTC	1-06	VEG	64.0	PPM		
CARLILE		R21E	22	D14	ORNL	1-05	SOIL	3.7	PPM	DRY	1.1
CARLILE		R21E	22	D14	ORNL	1-05	SOIL	2.7	PPM	WET	
CRES	T12N	R21E	27	K3	ORNL	1-05	SOIL	2.9	PPM	DRY	
CRES	T12N	R21E	27	K3	ORNL	1-05	SOIL	2.0	PPM	WET	1.1
CRES	T12N	R21E	27	K3	ORNL	1-05	WATER	0.05	UM/ML	HL I	- 112
CRES08	T12N	R21E	27	K3	KMTC	1-04	VEG	12.0	PPM	CRES	120
	T12N	R21E	27	K3	KMTC	1-04	SOIL	<10	PPM	CRES	1.5
CRES10	T12N	R21E	27	K3	KMTC	1-04	SOIL	<10	PPM	CRES	
	TI2N	R21E	27	K4	KMTC	1-04	SOIL	15.0	PPM	DRES	1.5.5
DRES12	T12N	R21E	27	K4	KMTC	1-04	VEG	220.0	PPM	DRES	1 C 1
DRES12	T12N	R21E	27	K4	KMTC	1-04	SOIL	13.0	PPM	DRES	
H-02	T12N	R21E	27	12	KMTC	1-04	VEG	32.0	PPM	CEDAR	1.11
	T12N	R21E	27	12	KMTC	1-04	VEG	0.95	PPM	HRES	
	T12N	R21E	27	12	KMTC	1-04	SOIL	<10	PPM	HRES	1.12
	T12N	R21E	28	I14	KMTC	1-04	VEG	<0.4	PPM	LRES	19.24
LRES05		R21E	28	I14 I14	KMTC	1-04		<10	PPM	LRES	
	T12N	R21E	28	I14 I14	KMTC	1-04		<10	PPM	LRES	
	T12N	R21E	28	114	KMTC	1-04	SOIL	<0.002	MG/L	FRONT	DODCH
RKRES15											PURCH
RKRES16		R21E	34	A13	KMTC	1-04	VEG	0.6	PPM	RKRES	
		R21E	34	A13	KMTC	1-04	SOIL	<10	PPM	RKRES	1.1.1
RRES17		R21E	34	B6	KMTC	1-04	SOIL	<10	PPM	RRES	1.1
	T11N	R21E	34	B6	KMTC	1-04	VEG	2.7	PPM	RRES	1.1.1.1
	TIIN	R21E	34	B6	KMTC	1-04	SOIL	<10	PPM	RRES	14.73
SHRES14		R21E	27	N6	KMTC	1-04	SOIL	<10	PPM	RRES	1.1
SMRES20		R21E	26	K6	KMTC	1-04	SOIL	11.0	PPM	SMRES	
SMRES21		R21E	26	K6	KMTC	1-04	VEG	0.6	PPM	SMRES	
T4-01E	T11N	R21E	2	C11	ORNL	1-06	SOIL	3.72	PPM	DRY	
T4-01E	TIIN	R21E	2	C11	ORNL	1-06	SOIL	3.1	PPM	WET	
T4-01E	TIIN	R21E	2	C11	ORNL	1-06	VEG	5.29	PPM		
T4-01W	T12N	R21E	27	N1	ORNL	1-06	SOIL	1.52	PPM	DRY	
T4-01W	T12N	R21E	27	NI	ORNL	1-06	SOIL	1.3	PPM	WET	
T4-02E	T11N	R21E	2	B8	CRNL	1-06	SOIL	3.59	PPM	DRY	

LABEL	TOWNSHP	RANGE	SECTION	SUBSECT	LAB	DATE	MEDIA	RESULTS	UNITS	COMMENTS
T4-02E	TIIN	R21E	2	B8	ORNL	1-06	SOIL	3.0	PPM	WET
T4-02E	T11N	R21E	2	B8	ORNL	1-06	VEG	1.99	PPM	
T4-02W	T12N	R21E	33	L14	ORNL	1-06	SOIL	1.5	PPM	WET
T4-02W	T12N	R21E	33	L14	ORNL	1-06	SOIL	1.67	PPM	DRY
T4-03E	T12N	R21E	35	N5	ORNL	1-06	SOIL	3.8	PPM	DRY
T4-03E	T12N	R21E	35	N5	ORNL	1-06	SOIL	2.9	PPM	WET
T4-03E	T12N	R21E	35	N5	ORNL	1-06	VEG	1.0	PPM	
T4-03W	T12N	R21E	34	B6	ORNL	1-06	SOIL	2.76	PPM	DRY
T4-03W	T12N	R21E	34	B6	ORNL	1-06	SOIL	2.2	PPM	WET
T4-04E	T12N	R21E	35	L3	ORNL	1-06	SOIL	3.73	PPM	DRY
T4-04E	T12N	R21E	35	L3	ORNL	1-06	SOIL	3.0	PPM	WET
T4-04E	T12N	R21E	35	L3	ORNL	1-06	VEG	0.807	PPM	
T4-04W	T12N	R21E	34	L13	ORNL	1-06	SOIL	2.45	PPM	DRY
T4-04W	T12N	R21E	34	L13	ORNL	1-06	SOIL	2.0	PPM	WET
T4-04W	T12N	R21E	34	L13	ORNL	1-06	VEG	3.68	PPM	
T4-05E	T12N	R21E	35	K1	ORNL	1-06	SOIL	2.29	PPM	DRY
T4-05E	T12N	R21E	35	K1	ORNL	1-06	SOIL	1.9	PPM	WET
T4-05E	T12N	R21E	35	K1	OP.NL	1-06	VEG	19.6	PPM	
TRES53	TIIN	R21E	23	M10	KMTC	1-06	VEG	16.0	PPM	
TRES54	T11N	R21E	23	M10	KMTC	1-06	SOIL	<10	PPM	
YRES51	T11N	R21E	26	N3	KMTC	1-05	VEG	8.1	PPM	
YRES52	T11N	R21E	26	N3	KMTC	1-05	SOIL	<10	PPM	

Table 5.2.7.24A

FAR FIELD ENVIRONMENTAL SAMPLING RESULTS FOR ALPHA PARTICLES

LABEL	TOWNSHP	RANGE	SECTION	SUBSECT	LAB	DATE	MEDIA	RESULTS	UNITS	COMMENTS
CARLILE	T12N	R21E	22	D14	OSDH	1-05	SOIL	0.1	PCI/G	
	T12N		22		OSDH	1-05	VEG	0.8	PCI/G	GRASS
CNTDLOI	TOON	DOAE				1-05	COTI	0.8		URASS
CNIRLUI	TOON	R24E	16	-	OSDH	1-06	SOIL	0.8	PCI/G	
CNIRLOI	TOON	R24E	16		OSDH	1-06	VEG	-0.4	PCI/G	
CNTRL02	TION	R22E	35	C14	OSDH	1-06	SOIL	0.4	PCI/G	
CNTRL02	TION	R22E	35	C14	OSDH	1-06	VEG	0.2	PCI/G	
ALL THE R. L. A. A.	100.00 00 00 00	A A A A	28	E10	OSDH	1-06	SOIL	-0.2	PCI/G	
CNTRL03	T11N	R22E	28	E10	OSDH	1-06	VEG	0.1	PCI/G	
CNTRI 04	TIIN	R22E	20	N10	OSDH	1-06	SOTI	1 1	PCI/G	
CNTRI 04	T11N T11N T11N T11N T11N T11N T09N T11N T11N	R22E	20	N10	0SDH	1-06	VEG	-0.4 0.2 -0.2 0.1 1.1 0.5	PCI/G	
CNTDLOS	TIIN	R22E	10	513	ACOU	1.00	COTI	0.0		
CNTRLOS	TITN	RECE	19	F13	OSDH	1-06	SUIL	-0.3	PCI/G	
CNIKLUS	TILN	R22E	19	F13	USUH	1-06	VEG	-0.3	PCI/G	
CNIRLOB	TOON	R22E	4	113	OSDH	1-06	VEG	-0.4	PCI/G	
CNTRL07	TIIN	R21E	19 19 4 13	87	OSDH	1-06	SOIL	-0.3 -0.3 -0.4 -0.3 -0.3	PCI/G	
CNTRL07	T11N	R21E	13	87	OSDH	1-06	VEG	-0.3	PCI/G	
CNTRL08	T11N T10N	R22E	13 8 8 21	K8	OSDH	1-06	SOIL	0.4	PCI/G	
CNIRL08	TION	RZZE	8	K8	OSDH	1-06 1-06 1-06	VEG	1.0	PCI/G	
CNTRL09	T11N T11N	R20E	21	C12	OSDH	1-06	SOIL	0.1	PCI/G	
CNTRIOS	TIIN	R20E	21	C12	OSDH	1-06	VEG	1.0	PCI/G	
CNTRLOG	TIIN	R20E	21	C12 C12	OCDH	1-06	VEG	0.4		
CNTRL 10	TIIN	REDE	21	112	OCDU	1-06 1-06	VEG	0.4		
CNIRLIU	TIIN	R2OE	<u> </u>	172	USUN	1-06	SOIL	-0.3	PCI/G	
CNIRLIO	TIIN	R20E	21	113		1-06	VEG	0.1	PCI/G	
CNTRLII	T12N	R19E	26	E3	OSDH		SOIL VEG	-0.3	PCI/G	
CNTRL11	T12N	R19E	26	E3		1-06	VEG	0.2	PCI/G	
CNTRL12	TIIN	R20E	14	N6	OSDH	1-06	SOIL	0.5	PCI/G	
CNTRL13	T11N	R23E	3	C2	OSDH	1-06	SOIL	0.4	PCI/G	
CNTRL13	TIIN	R23E	3	C2	OSDH	1-06	VEG	0.2	PCI/G	
CNTRL14	T11N	R23E	-		OSDH	1-06	SOIL		PCI/G	S OF SALLISAW
CNTRL15		R20E	16		OSDH	1-06	SOIL	-0.3	PCI/G	o or origination
CNTRL15		R20E	16		OSDH	1-06	VEG	-0.3	PCI/G	
CRES		R21E								
CRES	TION	RZIE	27	K3		1-05	SOIL		PCI/G	00100
CRES	TION	R21E	27	K3	USDH	1-05	VEG			GRASS
CRES	112N	R21E	27	K.3	OSDH	1-05	WATER		PCI/L	
HRES	T12N	R21E	27	12	OSDH		VEG	5.0	PCI/G	
T4-01E	TIIN	R21E	2	C11	OSDH	1-06	SOIL	-0.3	PCI/G	
T4-01E	T11N	R21E	2	C11	OSDH	1-06	VEG	0.3	PCI/G	EVERGREEN
T4-01W	T12N	R21E	27	NI	OSDH	1-06	SOIL	-0.2	PCI/G	
T4-02E	TIIN	R21E	2	88	OSDH	1-06	SOIL	0.4	PCI/G	
T4-02E	T11N	R21E	2	B8	OSDH	1-06	VEG	-0.2	PCI/G	EVERGREEN
T4-02W	T12N	R21E	33	L14	OSDH	1-06	SOIL	0.7	PCI/G	
T4-03E	T12N	R21E	35	N5	OSDH	1-06	SOIL	0.1	PCI/G	
T4-03E	TI2N	R21E	35	N5						EVERGREEN
					OSDH	1-06	VEG	-0.1	PCI/G	EVERGREEN
T4-03W	T12N	R21E	34	B6	OSDH	1-06	SOIL	0.1	PCI/G	
T4-04E	T12N	R21E	35	L3	OSDH	1-06	SOIL	0.1	PCI/G	
T4-04E	T12N	R21E	35	L3	OSDH	1-06	VEG	0.0	PCI/G	EVERGREEN
T4-04W	T12N	R21E	34	L13	OSDH	1-06	SOIL	0.1	PCI/G	
T4-04W	T12N	R21E	34	L13	OSDH	1-06	VEG	0.1	PCI/G	LEAVES
T4-05E	T12N	R21E	35	K1	OSDH	1-06	SOIL	-0.3	PCI/G	
T4-05E	T12N	R21E	35	K1	OSDH	1-06	VEG	4.5	PCI/G	EVERGREEN
		The second second			0.0011	1 00	1.0	1.5		to The I the I the be I t

Table 5.2.7.25A

FAR FIELD ENVIRONMENTAL SAMPLING RESULTS FOR BETA PARTICLES

LABEL	TOWNSHP	RANGE	SECTION	SUBSECT	LAB	DATE	MEDIA	RESULTS	UNITS	COMMENTS
CARLILE	T12N	R21E	22	D14	OSDH	1-05	SOIL	-0.3	PCI/G	
CARLILE		R21E	22	D14	OSDH	1-05	VEG	0.7	PCI/G	GRASS
CNTRL01		R24E	16	-	OSDH	1-06	SOIL	-0.3	PCI/G	GIUIOS
CNTRL01		R24E	16	-	OSDH	1-06	VEG	-0.1	PCI/G	
CNTRL01		R22E	35	C14		1-06	SOIL	0.3	PCI/G	
					OSDH				PC1/G	
CNTRL02		R22E	35	C14	OSDH	1-06	VEG	2.0		
CNTRL03		R22E	28	E10	OSDH	1-06	SOIL	-0.0	PCI/G	
CNTRL03	TIIN	R22E	28	E10	OSDH	1-06	VEG	0.3	PCI/G	
			20	N10	OSDH	1-06	SOIL	-0.3	PCI/G	
CNTRL04		R22E	20	N10	OSDH	1-06	VEG	0.3	PCI/G	
CNTRL05		R22E	19	F13	OSDH	1-06	SOIL	-0.0	PCI/G	
CNTRL05		R22E	19	F13	OSDH	1-06	VEG	0.2	PCI/G	
CNTRL06		R22E	4	I13	OSDH	1-06	VEG	0.6	PCI/G	
CNTRL07	TIIN	R21E	13	B7	OSDH	1-06	SOIL	0.5	PCI/G	
CNTRL07	TIIN	R21E	13	B7	OSDH	1-06	VEG	0.8	PCI/G	
CNTRL08	TION	R22E	8	K8	OSDH	1-05	SOIL	-0.4	PCI/G	
CNTRL08	T10N	R22E	8	K8	OSDH	1-06	VEG	0.7	PCI/G	
CNTRL09		R20E	21	C12	OSDH	1-06	SOIL	-0.3	PCI/G	
CNTRL09		R20E	21	C12	OSDH	1-06	VEG	7.0	PCI/G	
CNTRL09		R20E	21	C12	OSDH	1-06	VEG	-0.3	PCI/G	
CNTRL10		R20E	21	I13	OSDH	1-06	SOIL	0.3	PCI/G	
CNTRL10		R20E	21	I13	OSDH	1-06	VEG	-0.2	PCI/G	
CNTRL11		R19E	26	E3	OSDH	1-06	SOIL	-0.2	PCI/G	
			26	E3		1-06	VEG	0.1	PCI/G	
CNTRL11		R19E			OSDH					
CNTRL12		R20E	14	N6	OSDH	1-06	SOIL	0.9	PCI/G	
CNTRL13		R23E	3	C2	OSDH	1-06	SOIL	0.5	PCI/G	
CNTRL13		R23E	3	C2	OSDH	1-06	VEG	3.0	PCI/G	
CNTRL14		R23E	*	-	OSDH	1-06	SOIL	0.6	PCI/G	S OF SALLISAW
CNTRL15		R20E	16	G13	OSDH	1-06	SOIL	-0.4	PCI/G	
CNTRL15		R20E	16	G13	OSDH	1-06	VEG	0.3	PCI/G	
CRES	T12N	R21E	27	K3	OSDH	1-05	SOIL	-0.0	PCI/G	
CRES	T12N	R21E	27	K3	OSDH	1-05	VEG	0.2	PCI/G	GRASS
CRES	T12N	R21E	27	K3	OSDH	1-05	WATER	-1.0	PCI/L	
HRES	T12N	R21E	27	12	OSDH	1/5	VEG	1.9	PCI/G	
T4-01E	T11N	R21E	2	C11	OSDH	1-06	SOIL	0.1	PCI/G	
T4-01E	T11N	R21E	2	C11	OSDH	1-06	VEG	-0.2	PCI/G	EVERGREEN
	T12N	R21E	27	N1	OSDH	1-05	SOIL	-0.2	PCI/G	
	T11N	R21E		B8	OSDH	1-06	SOIL	-0.2	PCI/G	
	TIIN	R21E		B8	OSDH	1-06	VEG	-0.3	PCI/G	EVERGREEN
	T12N	R21E		L14	OSDH	1-06	SOIL	0.1	PCI/G	
	T12N	R21E		N5	OSDH	1-06	SOIL	0.6	PCI/G	
	T12N	R21E		N5	OSDH	1-06	VEG	0.2	PCI/G	EVERGREEN
		R21E		86	OSDH	1-06	SOIL	-0.6	PCI/G	LILMUNLLI
	T12N							0.1	PCI/G	
	T12N	R21E		L3	OSDH	1-06	SOIL			EVEDODEEN
	T12N	R21E	35	L3	OSDH	1-06	VEG	-0.2	PCI/G	EVERGREEN
	T12N	R21E		L13	OSDH	1-06	SOIL	-0.4	PCI/G	LEAVES
	T12N	R21E		L13	OSDH	1-06	VEG	0.4	PCI/G	LEAVES
	T12N	R21E	35	K1	OSDH	1-06	SOIL	-0.4	PCI/G	EUE DODE EU
T4-05E	T12N	R21E	35	K1	OSDH	1-06	VEG	0.8	PCI/G	EVERGREEN

APPENDIX 5.3.1

RAW SURFACE WATER DATA

APPENDIX 5.3.1

Surface waters were sampled by both SFC and OSDH. Sampling locations are shown in Appendix 5.2.1 and in Section 5.3.1. Sampling locations and results are listed in the tables which follow.

NEAR FIELD ENVIRONMENTAL SAMPLING RESULTS FLUORIDE IN WATER SAMPLES

LABEL	NORTH	EAST	LAB	DATE	TYPE	RESULT	UNITS	COM	MENTS
1W	14500	14800	КМТС	1-06	F	0.1	mg/1		
2W	13250	13950	KMTC	1-06	F	0.1	mg/l		
3W	12750	11200	KMTC	1-06	F	0.2	mg/1		
44	11700	11150	KMTC	1-06	F	0.2	mg/1		
5W	10100	10950	KMTC	1-06	F	0.4	mg/1		
6W	9850	9450	KMTC	1-06	F	0.1	mg/1		
7W	10600	7000	KMTC	1-06	F	0.2	mg/1		
8W	11800	7000	KMTC	1-06	F	0.3	mg/1		
10W	11800	6200	KMTC	1-06	F	0.1	mg/1		
11W	16350	9450	KMTC	1-06	F	<0.1	mg/1		
12W	6750	11050	KMTC	1-06	F	0.1	mg/1		
13W	7475	13350	KMTC	1-06	F	<0.1	mg/l		
14W	18450	6950	KMTC	1-06	F	<0.1	mg/1		
14W	18450	6950	KMTC	1-08	F	<0.1	mg/1		
15W	12900	5950	KMTC	1-06	F	<0.1	mg/1		
15W	12900	5950	KMTC	1-08	F	<0.1	mg/1		
17W	7650	4350	KMTC	1-06	F	0.2	mg/1		
18W	15100	10300	KMTC	1-08	F	0.1	mg/1		
19W	6850	9050	KMTC	1-08	F	0.2	mg/1		
20W	7750	17000	KMTC	1-08	F	0.1	mg/1		
300	12450	9900	KMTC	1-06	F	0.3	mg/l		
301	12450	9900	KMTC	1-14	F	0.2	mg/l		
302	12550	9550	KMTC	1-06	F	1.6	mg/l		
303	12550	9550	KMTC	1-14	F	1.5	mg/1		
304	12550	9400	KMTC	1-06	F	2.0	mg/l		
305	12550	9400	KMTC	1-14	F	4.8	mg/l		
306	11450	9550	KMTC	1-06	F	0.5	mg/1		
307	11450	9550	KMTC	1-14	F	0.5	mg/1		
LAGOON	12800	9200	KMTC	1-17	F	3.9	mg/l		
LRES07	6950	9900	KMTC	1-04	F	<0.1	mg/1	RAT	N BARREL
PW 1632	13250	9400	KMTC	1-15	F	<0.1	mg/1		or o
CRES	6150	10800	ORNL	1-05	F	0.09	ug/ml		
T1-10	12750	11100	ORNL	1-05	F	0.21	ug/ml		
T1-10	12750	11100	ORNL	1-05	F	0.19	ug/ml		
T2-04	10350	8400	ORNL	1-05	F	0.15	ug/ml		
T2-06	11450	9950	ORNL	1-05	F	0.46	ug/ml		
T2-07	11700	11100	ORNL	1-06	F	0.09	ug/ml		
COMP-3	10050	10700	OSDH	1-16	F	0.45	mg/l	LAR	COMPARISON
COMP-4	14450	14800	OSDH	1-16	F	0.27	mg/1		COMPARISON
COMP-5	12700	11200	OSDH	1-16	F	0.24	mg/1		COMPARISON
CRES	6150	10800	OSDH	1-05	F	0.19	mg/1	LHO	COM ANT DOM
KM H20	13250	9400	OSDH	1-17	F	0.24	mg/1		
LAGOON	12800	9200	OSDH	1-17	F	4.40	mg/1		
RD-03	6600	10800	OSDH	1-16	F	0.22	mg/1		
RD-04	6350	10750	OSDH	1-16	F	0.22	1997 I		
RD-05	6200	10850	OSDH	1-16	F	0.30	mg/1 mg/1		
RD-06	6900	10200	OSDH	1-16	F	0.17			
10 00	0300	10200	USUN	1-10	r	0.19	mg/1		

NEAR FIELD ENVIRONMENTAL SAMPLING RESULTS FLUORIDE IN WATER SAMPLES

LABEL	NORTH	EAST	LAB	DATE	TYPE	RESULT	UNITS	COMMENTS
RD-07	10100	10700	OSDH	1-16	F	0.32	mg/1	
RD-08	9900	9450	OSDH	1-16	F	0.20	mg/1	
STRM GH	12650	9550	OSDH	1-16	F	7.00	mg/	
T1-03&4	12500	8750	OSDH	1-14	F	0.23	mg/l	
T1-06	12500	9500	OSDH	1-14	F	1.80	mg/l	
T1-10	12750	11100	OSDH	1-05	F	0.2	mg/1	
T1-13	12200	13150	OSDH	1-14	F	0.41	mg/1	
T1-14	13250	14000	OSDH	1-14	F	0.29	mg/1	
T2-04	10350	8650	OSDH	1-14	F	0.36	mg/1	
T2-04	10350	8400	OSDH	1-05	F	0.26	mg/1	
T2-06	11500	9950	OSDH	1-06	F	0.37	mg/1	
T2-07	11500	10650	OSDH	1-05	F	0.23	mg/1	
T2-08	11500	11650	OSDH	1-16	F	0.17	mg/1	
T2-10	11400	13700	OSDH	1-16	F	0.16	mg/1	
T3-04A	8150	10800	OSDH	1-16	F	0.12	mg/1	
T5-04	1550	8500	OSDH	1-08	F	0.32	mg/1	

1

1

NEAR FIELD ENVIRONMENTAL SAMPLING RESULTS URANIUM IN WATER SAMPLES

LABEL	NORTH	EAST	LAB	DATE	TYPE	RESULT	UNITS	COMMENTS
1W	14500	14800	KMTC	1-06	U	0.11	mg/1	
2W	13250	13950	KMTC	1-06	U	0.23	mg/1	
3W	12750	11200	KMTC	1-06	U	0.19	mg/1	
4W	11700	11150	KMTC	1-06	U	<0.002	mg/1	
5W	10100	10950	KMTC	1-06	U	0.011	mg/l	
6W	9850	9450	KMTC	1-06	U	<0.002	mg/1	
7W	10600	7000	KMTC	1-06	U	0.32	mg/1	
8W	11800	7000	KMTC	1-06	U	0.6	mg/1	
10W	11800	6200	KMTC	1-06	U	<0.002	mg/1	
11W	16350	9450	KMTC	1-06	U	<0.002	mg/1	
12W	6750	11050	KMTC	1-06	U	0.003	mg/1	
13W	7475	13350	KMTC	1-06	U	<0.002	mg/1	
14W	18450	6950	KMTC	1-06	U	<0.002	mg/1	
14W	18450	6950	KMTC	1-08	U	0.003	mg/1	
15W	12900	5950	KMTC	1-06	U	0.018	mg/1	
15W	12900	5950	KMTC	1-08	U	0.013	mg/1	
17W	7650	4350	KMTC	1-06	U	0.006	mg/1	
18W	15100	10300	KMTC	1-08	U	0.002	mg/1	
19W	6850	9050	KMTC	1-08	U	0.004	mg/1	
20W	7750	17000	KMTC	1-08	U	<0.002	mg/1	
300	12450	9900	KMTC	1-06	U	0.042	mg/1	
301	12450	9900	KMTC	1-14	U	0.017	mg/1	
302	12550	9550	KMTC	1-06	U	0.94	mg/1	
303	12550	9550	KMTC	1-14	U	1.1	mg/1	
304	12550	9400	KMTC	1-06	U	1.2	mg/1	
305	12550	9400	KMTC	1-14	U	0.25	mg/1	
306	11450	9550	KMTC KMTC	1-14 1-06	UU	0.25	mg/1	
306	11450	9550	KMTC	1-14	U	0.21	mg/1	
307	11450	9550	KMTC	1-07	U	0.29	mg/1 mg/1	
COMB ST	12600	8650	KMTC	1-06	U	0.95	mg/1	
COMB ST	12600	8650	KMTC	1-05	U	0.36	mg/1	
COMB ST COMB ST	12600 12600	8650 8650	KMTC	1-03	U	0.36	mg/1	
		8650	KMTC	1-04		0.61	mg/1	
COMB ST COMB ST	12600		KMTC	1-02	U	0.57	mg/1	
COMB ST	12600 12600	8650 8650	KMTC	1-02	U	0.43	mg/1	
	12800	9200	KMTC	1-17	U	I.45	mg/1	
LAGOON LRES07	6950	9900	KMTC	1-04	U	<0.002	mg/1	FRONT PORCH
		9400	KMTC	1-15	U	<0.002	mg/1	PRONT PORCH
PW 1632 CRES	13250 6150	10800	ORNL	1-15	U	0.05	ug/ml	
		11100	ORNL	1-05	U	0.15	ug/ml	
T1-10 T1-10	12750 12750	11100	ORNL	1-05	U	0.026	ug/ml	
T2-04	10350	8400	ORNL	1-05	U	0.003	ug/ml	
T2-04	11450	9950	ORNL	1-05	U	0.003	ug/ml	
T2-06	11450	11100	ORNL	1-05	U	0.03	ug/ml	
12-07	11/00	11100	UNNE	1-00	0	0.05	ug/mi	

NEAR FIELD ENVIRONMENTAL SAMPLING RESULTS ALPHA IN WATER SAMPLES

LABEL	NORTH	EAST	LAB	DATE	TYPE	RESULT	UNITS	COMMENTS
CRES	6150	10800	OSDH	1-05	ALPHA	0.0	pCi/1	
T1-10	12750	11100	OSDH	1-05	ALPHA	100.0	pCi/1	
T2-04	10350	8400	OSDH	1-05	ALPHA	150.0	pCi/1	
T2-06	11500	9950	OSDH	1-06	ALPHA	124.0	pCi/1	
T5-04	1550	8500	OSDH	1-08	ALPHA	7.0	pCi/1	

NEAR FIELD ENVIRONMENTAL SAMPLING RESULTS BETA IN WATER SAMPLES

LABEL	NORTH	EAST	LAB	DATE	TYPE	RESULT	UNITS	COMMENTS
CRES	6150	10800	OSDH	1-05	BETA	-1.0	pCi/1	
T1-10	12750	11100	OSDH	1-05	BETA	12.0	pCi/1	
T2-04	10350	8400	OSDH	1-05	BETA	9.0	pCi/1	
T2-06	11500	9950	OSDH	1-06	BETA	26.0	pCi/1	
T5-04	1550	8500	OSDH	1-08	BETA	5.0	pCi/1	

ľ

FAR FIELD ENVIRONMENTAL SAMPLING RESULTS FOR FLUORIDE IN WATER SAMPLES

LABEL	TOWNSHIP	SECTION	SUBSECT	DATE	LAB	RESULTS	UNITS	COMMENTS
16W	TIIN	34	B6	1-06	кмтс	0.2	MG/L	
25L14	T12N	25	L14	1-14	KMTC	<0.1	MG/L	
26E 7	T12N	26	E7	1-14	KMTC	<0.1	MG/L	CUDEACE MATER
26F 6	T12N	26	F6	1-14	KMTC	<0.1	MG/L	SURFACE WATER
26L 1	T12N	26	L1	1-14	KMTC	<0.1	MG/L	SURFACE (POND)
271 0	T12N	27	10	1-14	KMTC	<0.1	MG/L	SURFACE WATER
27L 1	T12N	27	L1	1-14	KMTC	0.1	MG/L	
27L 4	T12N	27	L4	1-14	KMTC	0.1	MG/L	SURFACE WATER
27L 8	T12N	27	L8	1-14	KMTC	<0.1	MG/L	SURFACE WATER
28H14	T12N	28	H14	1-14	KMTC	0.1	MG/L	SURFACE WATER
36D11	T12N	36	D11	1-14	KMTC	<0.1	MG/L	
36L 5	T12N	36	L5	1-14	KMTC	<0.1	MG/L	
36M10	T12N	36	M10	1-14	KMTC	<0.1	MG/L	SURFACE WATER
LRES07	T12N	28	I14	1-04	KMTC	<0.1	MG/L	RAIN BARREL
CRES	T12N	27	K3	1-05	ORNL	0.09	UM/ML	
CRES	T12N	27	K3	1-05	OSDH	0.19	MG/L	
DIRTYCR	TIIN	6	F9	1-14	OSDH	0.32	MG/L	
P-01	T12N	26	A11	1-13	OSDH	0.23	MG/L	
P-10	T12N	28	F14	1-13	OSDH	0.49	MG/L	
REDHILL	TIIN	9	H4	1-14	OSDH	0.15	MG/L	
STRM GH	TIIN	4	KO	1-16	OSDH	7.00	MG/L	
SWR PND	TIIN	4	KO	1-14	OSDH	0.18	MG/L	
T4-03W	T12N	34	C7	1-06	OSDH	0.17	MG/KG	
WR-01	TIIN	9	N13	1-15	OSDH	0.18	MG/L	
WR-02	TIIN	9	N2	1-15	OSDH	0.26	MG/L	
WR-04	TIIN	8	L14	1-15	OSDH	0.1	MG/L	
WR-06	TIIN	5	N3	1-15	OSDH	0.26	MG/L	
111 00	1 7 7 14	5	1152	1 10	00011			

1.

FAR FIELD ENVIRONMENTAL SAMPLING RESULTS FOR URANIUM IN WATER SAMPLES

LABEL	TOWNSHIP	SECTION	SUBSECT	DATE	LAB	RESULTS	UNITS	COMMENTS	S	
16W	T11N	34	B6	1-06	KMTC	0.006	MG/L			
25L14	T12N	25	L14	1-14	KMTC	<0.002	MG/L			
26E 7	T12N	26	E7	1-14	KMTC	0.006	MG/L			
26F 6	T12N	26	F6	1-14	KMTC	0.003	MG/L	SURFACE	WATER	
26L 1	T12N	26	L1	1-14	KMTC	0.003	MG/L			
271 0	T12N	27	IO	1-14	KMTC	<0.002	MG/L	SURFACE	WATER	
27L 1	T12N	27	L1	1-14	KMTC	<0.002	MG/L			
27L 4	T12N	27	L4	1-14	KMTC	<0.002	MG/L	SURFACE	WATER	
27L 8	T12N	27	L8	1-14	KMTC	0.004	MG/L	SURFACE	WATER	
28H14	T12N	28	H14	1-14	KMTC	<0.002	MG/L	SURFACE	WATER	
36D11	T12N	36	D11	1-14	KMTC	0.004	MG/L			
36L 5	T12N	36	L5	1-14	KMTC	0.003	MG/L			
36M10	T12N	36	M10	1-14	KMTC	<0.002	MG/L			
LRES07	T12N	28	114	1-04	KMTC	<0.002	MG/L			
CRES	T12N	27	K3	1-05	ORNL	0.05	UM/ML			


Page No.	1		D ENVIRON OR ALPHA		a secolar and the second		TS	
LABEL	TOWNSHIP	SECTION	SUBSECT	DATE	LAB	RESULTS	UNITS	COMMENTS
CRES	T12N	27	К3	1-05	OSDH	0.0	PCI/L	
			D ENVIRON FOR BETA				TS	
LABEL	TOWNSHIP	SECTION	SUBSECT	DATE	LAB	RESULTS	UNITS	COMMENTS
CRES	T12N	27	K3	1-05	OSDH	-1.0	PCI/L	

APPENDIX 5.4.1

AN AERIAL RADIOLOGICAL SURVEY OF THE KERR-MCGEE FACILITY AND SURROUNDING AREA ERERGY MEASUREMENTS GROUP

EG&G SURVEY REPORT NRC-8105 JULY 1981

ł.



N.S.

AN AERIAL RADIOLOGICAL SURVEY OF THE

3

KERR - McGEE FACILITY

AND SURROUNDING AREA

GORE, OKLAHOMA DATES OF SURVEY: 25-30 JULY 1980

DISCLAIMER

This report was prepared as an account of work sponsored by the United States Government. Neither the United States nor the United States Department of Energy, nor any of their employees, makes any warranty, express or implied, or assumes any legal liability or responsibility for the accuracy, completeness, or usefulness of any information, apparatus, product, or process disclosed, or represents that its use would not infringe privately owned rights. Reference herein to any specific commercial product, process, or service by trade name, mark, manufacturer, or otherwise, does not necessarily constitute or imply its endorsement, recommendation, or favoring by the United States Government or any agency thereof. The views and opinions of authors expressed herein do not necessarily state or reflect those of the United States Government or any agency thereof.

Printed in the United States of America.

Available from:

National Technical Information Service U.S. Department of Commerce 5285 Port Royal Road Springfield, Virginia 22161

Price:

Printed Copy \$5.00 Microfiche \$3.50 ENERGY MEASUREMENTS GROUP

EG&G Survey Report NRC-8105 JULY 1981

AN AERIAL RADIOLOGICAL SURVEY OF THE

KERR - MCGEE FACILITY

AND SURROUNDING AREA

GORE, OKLAHOMA

DATES OF SURVEY: 25-30 JULY 1980

E. L. Feimster Project Scientist

REVIEWED BY

W. J. Tipton, Head Radiation Sciences Section

This Document is UNCLASSIFIED

20 Stoke

G. P. Stobie Classification Officer

This work was performed by EG&G for the United States Nuclear Regulatory Commission through an EAO transfer of funds to Contract Number DE-AC08-76NV01183 with the United States Department of Energy.

ABSTRACT

An aerial radiological survey was conducted over a 17 km² area, centered on the Kerr-McGee facility at Gore, Oklahoma, during 25-30 July 1980. An isoradiation contour map was produced that shows the distribution of exposure rates (at the 1 m level) due to terrestrial gamma ray activity within the surveyed area. These results showed several areas directly over the Kerr-McGee facility with radiation levels significantly above natural background. Spectral analysis revealed the presence of protactinium-234, a daughter of separated uranium-238 (²²⁶Ra and daughters removed). Bismuth-214, a daughter of the natural uranium-238 decay chain, was also present. The estimated maximum exposure rate range over this facility was 220-320 μ R/h. The typical background ranges were 4 to 6.5 μ R/h over vegetation and 6.5 to 9.0 μ R/h over cultivated land areas. These exposure rate ranges are from terrestrial sources only. An additional 4 μ R/h due to cosmic ray contributions must be added to obtain the total external exposure rate.

CONTENTS

3

đ

2 Abstract

ľ

Sections

- 4 1.0 Introduction
- 4 2.0 Site Description
- 4 3.0 Survey Procedures and Equipment
- 5 4.0 Discussion and Results

Figures

- 6 1 Exposure Rate Isoradiation Contours Overlaid on an Aerial Photograph
- 7 2 Gamma Pulse Height Spectrum Typical of Natural Background
- 7 3 Gamma Pulse Height Spectrum Indicating the Presence of Anomalous Concentrations of Separated U-238 and Radium Daughters from the Normal Uranium-238 Decay Chain
- 8 Appendix—Survey Parameters
- 9 References
- 10 Distribution

The United States Department of Energy (DOE) maintains an aerial surveillance operation called the Aerial Measuring Systems (AMS). AMS is operated for the DOE by EG&G, Inc. This continuing nationwide program, started in 1958, involves surveys to monitor and document background radiation levels throughout the United States. At the request of the DOE or other federal and state agencies (such as the United States Nuclear Regulatory Commission), the AMS is deployed for various aerial survey operations. This report is the result of a survey requested by the NRC for the Kerr-McGee facility located at Gore, Oklahoma.

Aerial radiological detection systems average the radiation levels due to gamma-emitting radionuclides existing over an area of several acres. The systems are capable of detecting anomalous gamma count rates and determining the specific radionuclides causing the anomalies; however, because of averaging, they tend to underestimate the magnitude of localized sources as compared with ground-based readings. As such, the indicated radiation levels in the vicinity of anomalies are not definitive. Ground surveys are required for accurate definition of the extent and intensity of such anomalies.

The results of the survey are reported as radiation exposure rates in microroentgens per hour (μ R/h) at 1 meter above the ground surface. Approximate anñúal radiation dose levels expressed as millirem per year (mrem/y) are obtained by multiplying μ R/h by 8.76. This conversion number applies only to the external radiation dose component.

Natural background radiation originates from radioactive elements present in the earth and cosmic rays entering the earth's atmosphere from space. The terrestrial gamma rays originate primarily from the uranium decay chain, the thorium decay chain, and radioactive potassium. Local concentrations of these nuclides produce radiation levels at the surface of the earth ranging from 1 to 15μ R/h (9 to 130 mrem/y). Some areas with high uranium and thorium concentrations in surface minerals exhibit even higher radiation levels, especially in the western states. For example, in the Colorado Plateau area the average radition level is above 200 mrem/y. At some locations in Brazil and India, the natural radiation level is above 1000 mrem/y. One member of each of the uranium and thorium decay chains is a noble gas which can both diffuse through soil and be borne by air to other locations. Thus, the level of this airborne radiation depends on the meteorological conditions, the mineral content of the soil, the soil permeability, and other conditions existing at each location at any particular time. The airborne radiation contributes from 1 to 10% of the natural background radiation levels.

Cosmic rays, the space component, interact in a complicated manner with the elements of the earth's atmosphere and the soil. These interactions produce an additional natural source of gamma radiation. Radiation levels due to cosmic rays vary with altitude and geomagnetic latitude: they range from 3.7 to 23 μ R/h (up to 200 mrem/y).¹

2.0 SITE DESCRIPTION

The Kerr-McGee facility is located approximately 4.5 km southeast of the towns of Webber Falls and Gore, Oklahoma; its boundaries cover approximately 2000 acres. Most of the land area is covered with vegetation or cultivated for agricultural use. The site is bordered on the west by the Illinois River (Figure 1).

The primary operation at this facility is the conversion of the uranium compound U_3O_8 (referred to as "yellow cake") to UF₆, uranium hexaflouride. The raw material, U_3O_8 , is stored on site, however, the refined p oduct, UF₆, is shipped immediately to its various users.

3.0 SURVEY PROCEDURES AND EQUIPMENT

An aerial photograph was used to define the area to be surveyed. Parallel lines spaced at 76 m intervals were flown at 76 m altitude.

The survey vehicle, a Hughes H-500 helicopter, carried a pilot, an equipment operator and a lightweight version of a specialized data recording system. Two detector pods were mounted on the sides of the helicopter; each pod contained ten sodium iodide, Nal (TI), detectors. The crystal in each detector was 12.7 cm in diameter and 5.1 cm in height. Gamma ray signals from the twenty detectors were summed and routed through an analog-to-digital converter and pulse height analyzer. Gamma ray counting rates and energy spectral data were accumulated in 1-second intervals and recorded on magnetic tape.

The helicopter position was established with two systems: a Microwave Ranging System (MRS) and a radio altimeter. The MRS master station, mounted in the helicopter, interrogated two remote transceivers mounted at stationary positions overlooking the survey area. By measuring the round trip propagation time between the master and remote stations, the master computed the distance to each. These distances were also recorded on magnetic tape each second. In subsequent computer processing, they were converted to position coordinates.

The radio altimeter determined the aircraft altitude by connecting the time cf a round trip pulsed signal to distance between the aircraft and the ground. These data were also recorded on magnetic tape so that any variations in gamma signal strength caused by altitude fluctuation could be compensated accurately.

The detectors and electronic systems that accumulate and record the data are described in detail in previous reports.^{2,3}

Data processing was done primarily with a computer based analysis laboratory mounted in a mobile van. An extensive collection of software routines was available for data processing. The first data processing that was accomplished produced terrestrial exposure rate isoradiation contours. These contours were constructed from gross count rate numbers, which refer to integral count rates in that portion of the gamma ray energy spectrum between 0.05 and 3.0 MeV. These count rates were converted to exposure rates at 1 m above the ground using a factor of 940 counts per second per μ R/h, a number obtained from calibration data over a test range. The

exposure rate isoradiation contours shown in Figure 1 are from terrestrial sources only. A total exposure rate can be obtained by adding an average cosmic ray value of 4 μ R/h.


4.0 DISCUSSION AND RESULTS

Results of the aerial survey are shown in Figure 1 in the form of exposure rate contours superimposed on a photograph of the survey area.

Several areas of localized activity significantly higher than background activity were encountered directly over the Kerr-McGee facility. The rapidly changing radiation levels over the site preclude accurate extrapolation of the aerial data to ground level exposure rates in the near vicinity of the site. The maximum estimated exposure rate (220-320 µR/h) over the localized areas may be an underestimate of the actual exposure rate as measured on the ground. In areas where heavy vegetation, such as trees, covered the terrain, the background exposure rate was in the range of 4-6.5 µR/h. In bare ground areas, cultivated fields, and areas with minimum vegetation the background was 6.5-9.0µR/h. One small area of cultivated land on the western boundary of the survey showed an exposure rate in the range of 9-13 µR/h. A gamma ray spectrum typical of the area is shown in Figure 2.


Spectral analysis of the anomalous areas indicated the presence of large quantities of natural uranium and seperated uranium (226Ra and daughters removed). Spectra accumulated over all the elevated activity areas showed photopeaks of protactinium-234 and bismuth-214 (Figure 3); these are daughter products associated with seperated uranium-238 and normal uranium-238, respectively.

The exposure rates quoted here and shown in Figure 1 are from terrestrial sources only. Cosmic ray radiation contributes an additional exposure rate of approximately 4μ R/h to the total exposure rate.



392

ų



393

1

1.0

-*10*

<

-x102

-x10²

+.8

ILess TB+12

APPENDIX: SURVEY PARAMETERS

Location:	Gore, Oklahoma and surrounding area, centered on the Kerr-McGee facility
Survey Coverage:	17 km ²
Survey Date:	25-30 July 1980
Project Scientist:	E. L. Feimster
Survey Altitude:	76 m (250 feet)
Line Spacing:	76 m (250 feet)
Detector Array:	20 Na. (TI)
Acquisition System:	REDAR III
Aircraft:	Hughes H-500 Helicopter
Data Dassessies	

Data Processing:

A. Gross Counts

Window: 0.05 - 3.0 MeV Conversion Factor: 940 cps/μR/h Cosmic ray contribution: 4 μR/h

REFERENCES

- Klement, A.W.; Miller, C.R.; Min, R.P.; Shleren, B. August 1972. Estimate of Ionizing Radiation Doses in the United States 1960-2000. R.S. EPA Report ORP/CD72-1. Washington D.C.: Environmental Protection Agency.
- Boyns, P.K. 1976. The Aerial Radiological Measuring System (ARMS): Systems, Procedures and Sensitivity. Report No. EGG-1183-1691. Las Vegas, NV: EG&G.
- 3. Jobst, J.E. 1979. "The Aerial Measuring Systems Program." Nuclear Safety 20:136-47.

DISTRIBUTION

NRC

-

L. K. Cohen (15)

DOE/OES

L. J. Deal (1)

EG&G

Z. G.	Burson, LVAO (1)	
J.F.	Doyle, LVAO (1)	
E. L.	Feimster, WAMD (1)	
H.A.	Lamonds, SBO (1)	
R. E.	Lounsbury, WAMD (4)	
L.G.	Sasso, LVAO (1)	
T. P.	Stuart, LVAO (1)	
W.J.	Tipton, LVAO (1)	

LIBRARIES

DOE/NV

J. K. Magruder (1)

AMO (8) Las Vegas (1) Santa Barbara (1)

KERR-MCGEE FACILITY NRC-8105 DATES OF SURVEY: 25-30 JULY 1980 DATE OF REPORT: JULY 1981

-


APPENDIX 5.5.1

1

RESULTS OF HPGE IN SITU MEASUREMENTS

APPENDIX 5.5.1

EG&G, Inc. performed a survey of locations using an intrinsic germanium detector, as described in Section 5.5.1. Results are shown in Figure 5.5.1. Figure 5.5.1A in this Appendix shows the measurement locations with their labels. Table 5.5.1A is a listing of the measurement results; it lists them as results for Pa-234m samples because Pa-234m peaks in the gamma spectrum were used for the measurement.





Page No. 1 2/27/86

NEAR FIELD ENVIRONMENTAL SAMPLING RESULTS FOR Pa-234m SAMPLES

L	ABEL	NORTH	EAST	LAB	DATE	MEDIA	RESULT	UNITS	COMMENTS
Z	01	8750	10100	EG&G	1-07	DEPO	< 90	nCi/m2	
	02		12600		1-07	DEPO	< 90	nCi/m2	
	03	9550	16050	EG&G	1-07	DEPO		nCi/m2	
	04	8850	7450	FG&G	1-07		< 90	nCi/m2	
	05	14350	14950	EG&G	1-08	DEPO	< 90	nCi/m2	
	06	11500	10200	FG&G	1-08	DEPO	95.3	nCi/m2	
	07	14200	10200 10100	EG&G	1-08	DEPO	< 90	nCi/m2	
		12900	10100	FG&G	1-08 1-08 1-09	DEPO	126.5	nCi/m2	
	09	8300	10100 10500	FG&G	1-09	DEPO	<90	nCi/m2	
	10	8350	11050 11050 11450	FG&G	1-09	DEPO		nCi/m2	
	11	8350	11450	EG&G	1-09	DEPO		nCi/m2	
	12		10100	EG&G			< 90	nCi/m2	
7	13	9100	10100 10100	EG&G	1-09			nCi/m2	
	14		10200	FG&G	1-09	DEPO	<90	nCi/m2	
	15	9950	10100	EG&G	1-09	DEPO	< 90	nCi/m2	
	16	10400	10100	EG&G	1-09	DEPO	94.9		
7	17	10900	10050	EGRG	1-09	DEPO	106.8		
7	18	12050					108.2		
7	19	12000	10100	EG&G	1-09 1-09	DEPO	96.8		
			10100		1-09	DEPO			
	20	11550	10600	EG&G	1-10		< 90	nCi/m2	
	21	11550	11200			DEPO		nCi/m2	
	22		11800			DEPO			
	23	11550		EG&G	1-10				
	24		9450	EG&G		DEPO			
	25		6750			DEPO			
	26		9700						
	27		10800		1-10				
	28	10750			1-10		96.4		
	29	10550			1-10				
		6300			1-10				
		6050			1-10				
	32		10800		1-11				
	33		11300		1-11				
	34	10800			1-11				
		10800	11850	EG&G	1-11	DEPO	<90	nCi/m2	
	36	12550	10650	EG&G	1-11	DEPO	<90	nCi/m2	
Z	37	12600	11550	EG&G	1-11	DEPO	109.5	nCi/m2	
Z	38	12550	12200	EG&G	1-11	DEPO	106.0	nCi/m2	
	39	12500	9850	EG&G	1-11	DEPO	210.5	nCi/m2	
Ζ	40	12550	9750	EG&G	1-11	DEPO	297.0	nCi/m2	
Z	41	12500	9700	EG&G	1-11	DEPO	383.0	nCi/m2	
Ζ	42	12500	9600	EG&G	1-11	DEPO	652.5	nCi/m2	
Z	43	12500	9450	EG&G	1-11	DEPO	465.0	nCi/m2	
	44	12500	9350	EG&G	1-11	DEPO	332.5	nCi/m2	
	45	12500	9250	EG&G	1-11	DEPO	395.0	nCi/m2	
Z	1.00								

Page No. 02/27/86

2

NEAR FIELD ENVIRONMENTAL SAMPLING RESULTS FOR Pa-234m SAMPLES

L	ABEL	NORTH	EAST	LAB	DATE	MEDIA	RESULT	UNITS	COMMENTS
ZZZ	47 48 49 50	15300 15400 14800 12500	8850 8100 7550 8650	EG&G EG&G EG&G EG&G	1-12 1-12 1-12 1-12	DEPO DEPO DEPO DEPO	111.5 95.8 <90 <90	nCi/m2 nCi/m2 nCi/m2 nCi/m2	
Z	51	12500	8950	EG&G	1-12	DEPO	<90	nCi/m2	WASHED DWN PRIOR
Z	52	12550	9100	EG&G	1-12	DEPO	<90	nCi/m2	WASHED DWN PRIOR

APPENDIX 5.6.1

MEDICAL EXAMINATION AND DIAGNOSIS FOR ONSITE WORKERS

SEQUOYAH FUELS CORPORATION

P.O. BOX 610 . GORE, OKLAHOMA 74435

February 13, 1986

Dr. Donald A. Cool Uranium Fuels Licensing Branch Division of Fuel Cycle & Material Safety U.S. Nuclear Regulatory Commission Washington, D.C. 20555

> RE: Interagency Meeting February 6, 1986

Dear Dr. Cool:

At the Interagency Meeting of February 6, 1986 your subcommittee on UF Assessment Task Group requested the following information on individuals by code designation:

- 1. Identification of employees who responded to the emergency to cool the ruptured UF cylinder and minimize the spread of UF $_6$ into the atmosphere.
- Identify employees sent to the hospital(s) and indicate the treatment received.
- Identify subcontractor and/or other outside individuals sent to the hospital(s) for examination and treatment.

In response to this request, we supplied the following information:

 Employees responding to the emergency who used water spray to cool and control UF₆ vapors. (Response was from upwind from the UF₆ vapors.)

5E, 7E, 8E, 11E, 16E, 20E, 25E, and 29E.

A total of 8 employees.

II. Employees examined at the Sequoyah Memorial Hospital in Sallisaw, Okla. and retained overnight for observation.

14E, 22E, 27E, and 31E.

A total of 4 employees.

A SUBSIDIARY OF KERR MCGEE CORPORATION

Dr. Donald A. Cool February 13, 1986 Page 2

> III. Employees examined at the Sequoyah Memorial Hospital in Sallisaw and retained more than 24 hours for observation and "Alka-Seltzer" treatment.

> > 1E, 2E, 6E, 8E, 9E, 12E, 16E, 17E, 29E, 23E, 25E, 26E, 28E, and 29E.

A total of 14 employees.

IV. Employees examined first at Sequoyah Memorial Hospital, Sallisaw, Okla. then transferred to Sparks Memorial Hospital, Ft. Smith, Arkansas for observation relative to possible lung irritation and/or damage from HF exposure.

5E, 15E, 18E.

A total of three employees.

V. Subcontractor and other outside personnel examined at Sequoyah Memorial Hospital, Sallisaw, Okla. and retained overnight for observation. (i.e. same category as employees in II above.)

5C, 6C, 9C, and 52P.

A total of 4 persons.

Only one person, 25E, was treated for skin irritation and burn from HF.

In addition to the above categories 10 employees were examined at Sequoyah Memorial Hospital and immediately released and 6 subcontractor personnel upwind of the UF₆ cylinder rupture, but on site, were examined at the same hospital and immediately released.

Your subcommittee also asked for the medical diagnosis for the employees in Category IV above and we received the following information from the hospital for all three individuals:

Diagnosis: Toxic inhalation exposure; medical code 506.3.

The employee fatality is not included in the above summary; the cause of death was also medical code 506.3.

Sincerely,

Wellion of Utroge

W. L. Utnage Sequoyah Facility Manager

WLU:dd

APPENDIX 5.6.2

URANIUM BIOASSAY RESULTS FOR SFC WORKERS AND OFFSITE INDIVIDUALS

U (ug/1) FECAL(ppm)	K ORNL	5240 5300 994 5300	3100 6800 3100 6800 .48				
(1/bn) N							
(1/bn) N		994 994	000	-			
(i/bn) N	¥		31	810	10		1500
n	PEMBROOK	950	6500 1500	825	13		1000 200
						3500	
				15		9400 3100 1400	
	SFC	6900 990 84 35 35 32 32 32	7400 2800 91 448 34 5	983 110 110 122 123 122 122 122 122	00288222	9200 2600 1300 210 110 82	1300 230 19
VOID	TIME	1548 0400 0400 0512 0612 0612 0612 0612 0612 0745 0745 0745 0745	1600 0400 0500 0500 0600 1500 1000	1915 0600 0600 0600 0530 0730 0600	2000 2200 2200 2200 2200 2200 0600 0600	1630 2100 0600 0700 0630 0605	1600 0400 0700 0630
	SAMPLE DATE	01/04/86 01/05/86 01/05/86 01/08/86 01/09/86 01/11/86 01/11/86 01/11/86	01/04/86 01/05/85 01/05/85 01/07/86 01/09/86 01/10/86 01/11/86	01/04/86 01/06/86 01/07/86 01/08/86 01/10/86 01/11/8 01/11/8	01/04/86 01/05/86 01/05/86 01/05/86 01/02/86 01/03/86 01/03/86	01/04/86 01/04/86 01/05/86 01/03/86 01/08/86 01/09/86	01/04/86 01/05/86 01/08/86 01/09/86
	10 NO:	¥	3£	36	ц а	56	6E

I

0)	AVERAGE	12.03	1.34	6.35		1.60	6.79
Intake (mg-U)	ORNL	9.6	1.2	5.6			5.5
Intal	- NI	15.1	1.60	5.1		2	8.3
	BNL	11.4	1.21	5.96		1.80	6.56
FECAL (ppm)		.28					
ļ			200	1200		091	
	ORNI	757	290	1210		472	2280
(1/bn) N	PEMBROOK	\$11 511	250	1650 600		200	625
				73			
	SFC	3800 820 15 54 36	5555 3366 86228 86288 86288 86288 86288 86288 8628 86288 866	1600 6400 180 130 34 34	640 130 22 10 10	270 270 113 115 115 115 115 115 115	1600 630 28 57
	VOID	1700 0800 22000 0600 0600 1200	1500 2015 2315 0300 1100 0630 0630 0630	1950 0500 0630 0630 0630 0630 0645	1700 1700 0830 0630 0600 0700	2230 0400 0600 0600 0600 0600 0600 0600 06	1835 0600 2130 2240
	SAMPLE DATE	01/04/86 01/05/86 01/05/86 01/07/86 01/07/86 01/12/86 01/12/86	01/02/86 01/04/86 01/04/86 01/05/86 01/05/86 01/05/86 01/07/86 01/07/86 01/12/86	01/04/86 01/05/86 01/05/86 01/07/86 01/08/86 01/09/86 01/10/66	01/04/86 01/05/86 01/07/86 01/08/86 01/10/86 01/10/86	01/04/86 01/05/86 01/05/86 01/06/86 01/06/86 01/06/86 01/06/86 01/09/86 01/12/86 01/12/86	01/04/86 01/05/86 01/06/86 01/07/86
	10 NO:	7.E	BE	36	10E	311	32.
					410		

410

-

				U	(ug/1)			FECAL(ppm)		Int	ake (mg-	-U)
10 NO;	SAMPLE DATE	VOID	SFC		PEMBROOK	ORNL			BNL	PNL	ORNL	AVERAGI
13E	01/04/86	1740	700		775	486			2.31	2.7	2.1	2.37
	01/05/86	0400	300			400			c. 31	2.1	2.1	2.31
	01/05/86	1330	95									
	01/06/86	1300	63									
	01/08/86	1000	27									
	01/09/86	0830	17									
	01/09/86	1230	10									
	01/10/86	0930	23									
	01/10/86	1230	5									
14E	01/04/86	2215	2500		2150	3030			8.18	7.6	5.8	7.19
	01/05/86	0400	390		400	495			0.10	1.0	2.0	1.15
	01/05/86	1330	52									
	01/06/86	1900	39									
	01/08/86	0630	28									
	01/09/86	0630	18									
	01/10/86	0600	15									
	01/11/86	0630	9									
	01/11/86	1000	25									
	01/13/86	0630	8									
15E	01/04/86	2300	190		200	153			1.04	1.5	0.99	1.18
	01/05/86		18									
	01/06/86	1.000	11									
	01/07/86	0600	9									
	01/07/86	0730	17									
	01/07/86	2315	6									
	01/08/86	0430	7									
	01/08/86		14									
	01/09/86	1500	6									
	01/10/86	1030	<2									
					2010				12.00	1.1		
16E	01/04/86	1845	3500		3250	2520			14.3	14.2	8.9	12.47
	01/05/86	0400	720		400	740						
	01/07/86	1800	42									
	01/10/86	1000	15									
	01/11/86	0900	10									
	01/13/86	0530	17					and the state of the second				
	01/12/86	0700						. 19				
17E	01/04/86	1600	11000	13000	13000	2250	2200		26.5	34.2	22.2	27.63
	01/05/86	0400	3200			3940						
	01/06/86	0800	110									
	01/07/86	0700	120									
	01/08/86	1500	93									
	01/09/86	0800	36									
	01/10/86	1039	50									
	01/11/86	0700	67									

					U (ug/1)	-		FECAL(ppm)		Inte	ke (mg	-U)
10 NO;	SAMPLE DATE	VOID	SFC		PEMBROOK	ORNL			BNL	PNL	ORNL	AVERAGE
18E	01/02/86 01/04/86 01/06/86 01/07/86 01/08/86 01/09/86 01/10/86 01/11/86 01/11/86 01/12/86 01/13/86	2300 1600 0800 0800 0800 0700 0800 0700 0800 0700	<2 3200 312 131 6 46 8 21 31 8	2700	3100	2540	2500		4.20	11,2	9.6	8.33
19E	01/04/86 01/05/86 01/07/86 01/09/86 01/10/86 01/13/86	1700 1300 0600 0600 0600 0600	620 64 21 3 29 <2						1.37	1.6	1,1	1.36
20E	01/04/86 01/04/86 01/05/86 01/08/86 01/08/86 01/09/86 01/10/86 01/11/86 01/11/86	1700 2245 0400 0730 0600 0600 0600 1600	250 140 170 27 18 12 13 23 10		180 110	124 166	120 170		1.03	1.6	2.0	1.54
21E	01/04/86 01/04/86 01/05/86 01/06/86 01/06/86 01/07/86 01/12/86	1700 2000 1300 0615 2100 2100 0630	62 47 3 5 17 9 3	67	50	84			0.138	0.32	0.32	0,26
22E	01/04/86 01/04/86 01/05/86 01/05/86 01/05/86 01/05/86 01/05/86 01/05/86 01/06/86 01/07/86 01/09/86	1700 1950 0400 1300 1400 0600 0700 0700	730 200 38 610 12 44 42 6		140 17	700 150			1.25	1.8	2.1	1.72
	01/10/86 01/11/86 01/12/86	0822 0600 0600	10 5									

412

	D NO: SAMPLE	23E 01/04/ 00/04/ 00/00	24E 01/04/86 01/05/86 01/06/86	25E 01/04/86 01/05/86 01/06/86 01/09/86 01/09/86 01/09/86	2/5 E 01/04/86 01/05/86 01/06/86 01/08/86 01/11/86 01/11/86	7E 01/04/86 01/04/86 01/04/86 01/05/86 01/05/86 01/02/86 01/02/86 01/12/86 01/12/86	28E 01/04/86 01/05/86 01/08/86
	DATE	888866 88866 88666666666 8666666666666	86 86 86	86 86 86 86 86	(86 (86 (86 (86 (86 (86	186 186 1866 1866 1866 1866 1866 1866 1	/86 /86 /86
	VOID	1945 2000 2330 0300 0700 0700 0700 0700 0805 0700 0805 0700 0805 0700 0805 0700 0805 2335	1900	1745 0800 1545 0800 0810 1647	2130 0400 1320 1315 1425	0600 1700 2000 0400 0600 0600 0600 0600 0600 06	1700 1300 1300
	SFC	4400 5500 3900 11100 147 71 71 71 71 71 71 71 71 71 71 71 71 71	90 8	760 370 52 12 12	7400 1500 110 73	2500 1100 780 39 33 33 33 33 33 33 33 33	1500 22 11 8
(1/bn) N	PEMBROOK	3500	ε	700 200	8000	1650	1400
1	ORNL	5340	68 28	742 286	5960	1360 526	1510
1					6000	1400	
FECAL(ppm)		4.1			14.		
	BNL	23.1	0.455	3.28	6.97	6.25	2.93
Inta	PNL	24.8	0.550	й.6	9.1	6.1	4,2
Intake (mq-U)	ORNL	16.1	0.33	3.0	15.0	<i>1</i> 7	2.4
(1)	AVERAGE	21.33	0.45	3.63	10.35	5.92	3.18

()	AVERAGE	5.9 19.20	0, 73	3.66
-bm) ax	ORML	5.9	0.68 0.99 0.73	9.4
Intake (mg-U)	PNL	24.0	0.68	<i>n</i> .1
	BNL	27.7	0.509	2.29
FECAL (ppm)				
		3600		1000
	ORNL	3620 3600 560	114	1010 1000
(1/bn) 0	PEMBROOK	6000	24	1200
		00th	151	1400
	SFC	13 4600 430 22 29 29 29	100 150 150 191 19	2300 1800 16 72 69
	VOID	1620 0400 1600	1700 1900 0400 1300 2130 2730	1700 1730 1600 0700 0700
	SAMPLE DATE	/03/ /05//05//05//05//05//05//00//00//00//0	01/11/86 01/04/86 01/05/86 01/05/86 01/05/86 01/05/86 01/13/86	01/04/86 01/04/86 01/04/86 01/05/86 01/05/86 01/07/86 01/08/86
	10 NO:	29£	30£	316

414

				n	(1/bn) N	1		FECAL (ppm)		intake (mg-U)	(n-bm)		
:ON DI	SAMPLE DATE	VOID	SFC		PEMBROOK	ORML			BNL	PNL	ORNI	AVERAGE	
10	01/04/86 01/05/86	1930 0400	48		49 65	<i>~~~</i>			0.0201	0.0244	0.0154	0.0200	
20	01/04/86 01/05/86		55			Ş			ł	1	0.0595	\$650.	
3C	RESULTS PENDING	NG											
40	01/04/86 01/05/86	1547	2600		1050	7560	1600						
50	01/04/86 01/04/86 01/05/86	0090	1230 1700 690										
60	01/04/86	1535		1400	1500	850			5.78	1.748	9.66	1.73	
	01/05/86 01/05/86 01/05/86		1300		500	1250							
70	01/04/86 01/05/86	1930	# (Y			41			0.0201	0.0245	0.0154	0.0200	
8C	01/04/86 01/04/86 01/05/86 01/05/86 01/05/86 01/07/86 01/07/86	1700 1800 0400 1300 1700 1500 1330	120 98 41 7 7 7		64 1	æ							
- 26	01/04/86 01/05/86	1555	230			230							
100	01/04/86 01/05/86	1515	22		Ş	ŝ			0.0196	0.0768	0.0768 0.0167	0.0377	

				U (u	(I/p	and the second	FECAL(ppm)		Intake	(mg-U)	-
10 NO:	SAMPLE DATE	VOID	SFC	P	EMBROOK	ORNL		BNL	PNL	ORNL	AVERAGE
1P	01/07/86	1535	3								
2P	01/07/86	1535	11								
3 P	01/05/86	2200	12								
4 P	01/05/86	0400	<2		<5	<5					
5P	01/05/86	1200	<2								
6P	01/05/86		12								
7 P	01/07/86	0345	<2								
8P	01/05/86		<2								
9P	01/06/86	1200	<2								
10P			<3								
11P			3								
12P	01/04/86 01/05/86	1530 0500 e	est. 93	69 26	110 30	92 39		0.429	0.463	0.291	0.39
13P	01/04/86 01/05/86	1530 0400es	240 st.110		240 60	170 83		0.862	1.183	0.729	0.92
14P	01/04/86 01/04/86	1930	<2 4		9	<5					
15P	01/04/86	1930	37 15 <2		<5	<5		0.189	0.227	0.142	0.19
	01/13/86	0630									
16P	01/07/86		<2								
17P	NO RESULTS										
18P	01/06/86		11								
19P	01/04/86		6 4		<5	<5					
20P	01/04/86 01/05/86 01/08/86	1540 1710	7 16 <2		15 11	10 <5					
21P	01/04/86	2400	<2								

I

				U (ug/1)	1		FECAL(ppm)		Intak	e (mq-U)	
ID NO:	SAMPLE DATE	VOID	SFC	PEMBR	ROOK	ORNL		BNL	PNL	ORNL	AVERAG
22P	01/07/86		<2								
23P	01/06/86	1400	<2								
24P	01/07/86	0935	<2								
25 P	01/05/86		53		32	20					
26P	01/05/86		47		19	20					
27P	01/07/86	1530	2								
28P	01/07/86	1530	2								
29P	01/07/86	1530	<2								
30P			10								
31P			9								
32P	01/05/86	1930	8								
33P	01/07/86		<2								
34P	01/06/86		<3								
35P	01/07/86		<2								
36P	01/05/86		<2								
37P	01/05/86		3								
38P	01/04/86	1620	<2								
39P	01/06/86	1715	<2								
40P			0								
41P	01/04/86	1900	<4	5		6					
42P	01/06/86		8								
43P	01/06/86		8		<5	<5					
44P			19								
45P	01/05/86 01/06/86	2315 1000	<3 18								

Liste d

				Anna and Andrews	U (ug/1)	in the second second	FECAL(ppm)		Intak	e (mg-U)	
	10 NO;	SAMPLE DATE	VOID	SFC	PEMBROOK	ORNL		BNL	PNL	ORNL	AVERAG
	46P	01/05/86 01/06/86 01/06/86 01/07/86	2330 0900 1830 0615	<2 <2 <2 3		<5 <5					
	479	01/05/86	1045	7	<5	<5					
	48P	01/08/86 01/10/86		<2 <2							
	49P	01/06/86	1015	7							
	50P	01/06/86	1015	<2							
	51P	01/05/86	1630	6							
	52P	01/04/86 01/04/86	1700 1900	<7 <7	<5	<5					
	53P	01/04/86 01/05/86	2245 1600	14 15	<5	17 <5					
	54P		0400 2040	<2 12							
110	55P	01/05/86	1845	10							
	56P	01/07/86	1115	<2							
	57P	01/07/86	1115	<2							
	58P	01/08/86		<2							
	59P	01/06/86	1955	<2							
	60P	NO RESULTS									
	61P	NO RESULTS									
	62P	01/06/86	1030	4							
	63P	01/06/86	0030	6							
	64P	01/07/86		<2							
	65P	01/06/86 01/06/86	1520 1520	<2 <2							
	66P			7							

				U (ug/1)		FECAL(ppm)		Intake	(mg-U)	
10 NO:	SAMPLE_DATE	VOID	SFC	PEMBROOK	ORNL		BNL	PNL	ORNL	AVERAGE
67P	01/07/86	0820	10							
68P	01/07/86	1300	6							
69P	01/06/86	1700	<2							
70P	01/07/86		<2							
71P	01/04/86	1930	<2 4	6	<5					
72P	01/04/86	1930	27 6	<5	-5					
73P	NO RESULTS									
74P	NO RESULTS									
75P	NO RESULTS									
76P	01/05/86		49							
77P	01/07/86	0600	14							
78P	01/04/86 01/05/86	0400	15 27	13	16 7					
79P	01/04/86 01/05/86	1745 0400	15 17	8 7	<5					
80P	01/06/86	1410	4							
81P	01/08/86		77							
82P	01/07/86		<2							
83P	01/04/86	1515	<7		<5		0.0326	0.0374	0.0184	0.03
84P	NO RESULTS									
85P	01/06/86	1030	8							
86P	01/07/86	1535	<2							
87P	01/09/86	1300	<2							
88 P	01/04/86 01/05/86	0400	17 <2		<5					
89P	01/04/36 01/05/86	1600 1600	13 5	<5	23 <5					

.

					U (ug/1)		FECAL(ppm)		Intak	e (mg~U)	
	ID NO:	SAMPLE DATE	VOID	SFC	PEMBROOK	ORNL		BNL	PNL	ORNL	AVERAGE
	90P	01/04/86 01/04/86	1650	<4 2	<5	<5					
	91P	01/08/86		.7							
	92P	01/04/86 01/05/86	2315 0420	<2 <2	<5 <5	<5 <12					
	93P	01/06/86	1215	7							
	94P	01/04/86	0400 1830	13 <4							
	95P	01/08/86		<2							
	96P	01/06/86	0950	5							
	978	NO RESULTS									
	98P	NO RESULTS									
	99P			<2							
		01/04/86		<2 <4	<5	<5					
420	100P	01/04/86 01/04/86		<2 4	<5 <5	<5 82					
	101P	01/04/86 01/04/86		<2 4	<5	<5					
	102P	01/04/86 01/05/86	1930 0400	4 <2	<5	14 <5					
	103P	01/07/86	1440	<2							
	104 P	01/07/86		<2							
	105P	01/07/86	1735	<2							
	106P	01/04/86	1900	<14	8 <5	11					
	107P	01/04/86	1930	<4	18 <5	<5					
	108P	01/04/86 01/04/86	1900	<2 8	<5 <5	<5 <5					
	109P	01/04/86 01/04/86	1900	<2 8	<5	<5					

-

	-			
	1	1	ť	Ľ.
	100.000	1	1	٢
	1	ii	ú	a
	1			
	3		2	٥
	1	ş		ļ
	100			7
	1			2
	3			2
	3			
	3			
	8			-
	1			
				-
	and a second sec		-	The second secon
	and a second sec	- TAN		
	and the second s	- TAN		
	and the second s	- CALL		The second secon
	and the second s	- TANK		a state of the sta
	and the second se	L CALL		The second secon
	Table - States	T L CALL		"HOLL
	Contraction of the local division of the loc	C T L CALL		THE PARTY OF THE P
	A TANK OF A DATE	A T LOAL	A DE	
	Contraction of the Contraction o	A T DALL	A DE	THE PARTY OF THE P
	And the second se	CAL VAL	A DECK DECK DECK DECK DECK DECK DECK DECK	
	The second se	CAL VAL	T ALL LUNC L	CONT DUN -
	A LOUDE DATE OF THE OWNER OWNER OF THE OWNER OWNE	LOAT LOAL		LUNI LUNI
	The second se	LOAT LOAL.		I LONI PURE
	A LOUDE DISTANCE OF THE REAL O	TANT VALL		I CONT FORM
	A DESCRIPTION OF THE PARTY OF T	CINKTINKI, C		I LUMITUM.
	A COLORADO CONTRACTOR OF A COLORADO	CLURTINE.		T LUM TUNE
	A DESCRIPTION OF THE PARTY OF T	CINATION.		I LUMITUM.
	き してきましてきのうちの 白いろしい		A DECK T DECK T DECK T	
	き しはあま うちまんひとう れたいと シー	THE PART OF THE PA	A REAL PROPERTY AND A REAL	
	Contraction of the Contraction of the Contraction of the	CINET DALL		
	Contraction of the Contraction of the Contraction of the	CINET DALL		
	Contraction of the Contraction of the Contraction of the	CINET DALL		
	Contraction of the Contraction of the Contraction of the	CINET DALL		
	Contraction of the Contraction of the Contraction of the	CINET DALL		
	ましい います ひろうろう ちちち ひろう しろう	CINET DALL		
	Contraction of the Contraction of the Contraction of the	CINET DALL		
	Contraction of the Contraction of the Contraction of the	CINET DALL		
	「日本書」」「「「「「」」」」」」」		- NIT I TITI	
	「日本書」」「「「「「」」」」」」」		- NIT I TITI	
	「日本書」」「「「「」」」」」」」」」」」」」」」」」」」」」」」」」」」」」」	CINET DALL	- NIT I TITI	

				n	(1/bn) N		FECAL(ppm)	m)	Intake	Intake (mg-U)	1
ID NO:	SAMPLE DATE	VOID	SFC		PEMBROOK	ORNL		BNL	PNL	ORNL	AVERAGE
110P	01/04/86 01/04/86	1930	6 A		\$	\$					
111P	01/04/86	1900	\$\$		<5	<5					
112P	01/07/86		9								
113P	01/04/86 01/05/86	1818 0400	52		99	99					
114P	01/04/86		6								
1156	01/04/86 01/05/86	2230	10 6		33	99					
116P	01/04/86	2230	5		15	\$					
4211	01/05/86	1000	19		<5>	\$					
118P			16								
119P	98/50/10	1600	<5		<5	\$					
120P	01/07/86		4								
121P	01/04/86	1600	43	35	20	111		0.204	0.229	0.116	0.18
122P	01/04/86	1800	20		31	20	16	0.0988	0.122	0.0663	3 0.10
123P	01/04/86	1645	20 40		31	16 20		0.192	0.213	0.116	0.17
124P			51		10	36					
125P	01/05/86		12								
126P	01/05/86		4								
1278	01/06/86	1310	<2								
128P			4								
129P	01/01/86		2								
130P	01/04/86	1920	6.6		9	\$					
131P	01/04/86	1910	10		<2	<2					
132P	01/06/86 01/07/86	1600	<5 6								

APPENDIX 5.6.3

INTAKE CALCULATION MODELS AND ASSUMPTIONS



BROOKHAVEN NATIONAL LABORATORY

ASSOCIATED UNIVERSITIES, INC.

Upton, Long Island, New York 11973

(516) 282 FTS 666 4250

Safety & Environmental Protection Division

January 17, 1986

Dr. Allen Brodsky Senior Health Physicist Division of Radiation Programs and Earth Sciences Office of Nuclear Regulatory Research Washington, D. C. 20555

SUBJECT: Estimate of Intake Post UF6 Release at Sequoyah Fuels Corporation

Dear Dr. Brodsky:

This is a letter to report results of estimates of intake, based on preliminary bioassay data sent to me on January 14, 1986. This letter briefly outlines the methods used to interpret the bioassay measurements, and additional details will be reported in the monthly report for FIN No. A-3289. The bioassay data are that reported to the Nuclear Regulatory Commission by Sequoyah Fuels Corporation, and the estimates of intake, which I computed from this bioassay data, are tabulated on the attachment.

In order to interpret the data, I assumed that urine was excreted at the rate of 1.4 liters per day, and from this assumption, I estimated the accumulated amount of uranium in urine. I based this on the reported times of sample, or the reported time of the single exposure to UF6, from which I estimated the volume of urine in each sample. Based on this volume and the reported uranium concentration, I estimated the micrograms of uranium excreted, from the point in time of a previous urine sample, up to the time of the sample under consideration. From this, I estimated the total micrograms of uranium excreted, from exposure up to the time of sample, and it was this estimate which led me to compute the estimated intake, which is reported in column 5 of the attachment.

The estimate of intake, in column 5, was based on the fraction of intake excreted to urine, up to the time of sample, and the estimate of total uranium excreted through the urine pathway, up to the time of sample. The fraction of intake was based on lung, GI tract, and systemic models, which were given in ICRP Publication 30. An inhalation intake of class D aerosols, with an AMAD of 1 micron, was assumed. An exposed person's estimate of intake showed some variation from sample to sample. Some of this variation was due to biological variables, and some was due to physical variables, such as those associated with sample collection and measurement, and those associated with computation and transcription of the results. The last column, on the attachment, represents the best estimate of intake, and it was based on the least squares method of fitting data.

2

If you have any questions, please do not hesitate to call.

Sincerely,

Edward T. Lesiard

Edward T. Lessard Program Manager, Interpretation of Bioassay Measurements

ETL/cjl

cc: J. Baum

Attachment*

*Note: The data in this attachment has been summarized in Appendix 5.6.2.

Colorizion of Intakes (estimated) and Organ Dosas for UF at Km.

/sourcestons:	
1. Time & intake = 1130 hr	5, Jan. 4, 1986.
2. 11% = Class D solubility	
3. Analised ICR P. 30 meter	slie vote constants, and assumed an
Massl with AMAD = 1.	
	(ICRP-23). Doily unnary oniput= 1.4 L
C. This withly sis weighted	equilly. For duplicate samples, the man
ring wes used as a si	
6. Contrasion Sectors:	
D. 5036 uci / a Unot.	2. 543 X104 Bg/g Unot.
7. Frind of Criculation: C	reputer code GENMOD, by J.R. Johnson
and D. W. Dun ford, Chell	River Nuchar Lolos, Onterio, Canada.
1. Julisots events 1.4 L unn	
Fronture:	
1. Tires (tenths Staday) pos	t-exposure were obtervised for each unnulyeit
2. A vielae of Uranium exerciti	m (mg/dry) was determined for each sample

2. A volue of Uronium exercition (mg/dry) was determined for each sample comming a drity wrincing exercition of 1.4 L.

3. The expected traction of uranium in urine for various times post experiments of the assay (10/100) to the organized dramma Generation Generation Generation, was determined. The avery a Generation for all bioresays of an individual was coloulated. This was its assumed intake, in my Unot.

4. The asserved intoke for each individual who converted to Becquerd white, and obser values were obtained from the GENMOD output fields.

Dowell P. Fiely 1-16-1

Makel by J. Toknism of Chark River Canada

Allen - A copy with the references added and some bust changes -

KEEckerman

5

ESTIMATES OF INTAKE OF URANIUM BY PERSONS NEAR ACCIDENTAL RELEASE AT KERR-McGEE SITE

> K. F. Eckerman and R. W. Leggett Health and Safety Research Division Oak Ridge National Laboratory

Estimates of uranium intake have been made for several persons on or near the site of the recent accident at the Kerr-McGee plant in Oklahoma. These estimates were based on uranium concentrations in urine of the subjects at times from a few hours to a few days after exposure.

Various excretion models have been applied to the general problem of estimating uranium uptake to blood based on the level of uranium in urine. Bernard and Struxness (1957) determined a best-fitting power function for urinary excretion data for patients injected with uranium nitrate:

$$Y_{-}(t) = 34.3 t^{-1.5}$$

where Y is the excretion rate in $\frac{\pi}{h}$ h and t is measured in hours. This function is not a good approximation to the data for the first day after injection and hence is not suitable for the present analysis.

The alternate function

$$Y_{u}(t) = 20\% t^{-1.5}$$
,

where t is in days, was later adopted by the ICRP (1968). This function is not suitable for use in the present situation because of the lack of detail offered for the first day after exposure. Moreover, this function may overestimate urinary excretion at early times.

A third approach is to combine the uranium retention functions and the generic clearance rate from the transfer compartment as given in ICRP 30 (1979). This approach yields seemingly reasonable excretion rates at early times when combined with the ICRP lung model for class D material (the classification for uranium hexafluoride given in ICRP 30). This is fortuitous, however, since the ICRP generic clearance rate from the transfer compartment appears to be a poor approximation for uranium, and the ICRP 30 retention models for uranium do not include consideration of the early times after exposure. (Errors involved in the blood clearance rate may be offset somewhat by errors in the lung model as applied to soluble uranium.) In general, the use of ICRP 30 retention models for estimation of excretion rates should be avoided because the ICRP 30 retention models were not designed or intended for analysis of short-term retention of radionuclides.

We appealed to the original injection data for the "Boston patients" considered by Bernard and Struxness (1957), combined with injection data on healthy persons and persons with osteoporosis (see Hursh and Spoor 1973) to derive an excretion function for uranium that has reached blood. The function

 $Y(t) = 7.3 e^{-0.2t} + 2.4 e^{-0.075t} + 0.08 e^{-0.011t} + 0.004 e^{-0.0014t}$

where Y_{th} is the %/h excreted in urine and t is in hours, appears to approximate these combined data closely for times from 30 minutes to

Received 2/4/86 a Grodsty

several days after injection.

The data provided to us on the Kerr-McGee subjects include only the concentration of uranium and the date/time of the urine sample. Some assumptions had to be made in order to relate our theoretical blood-tourine clearance curve to the concentrations of uranium in urine reported for various times of voiding. We assumed that the average urinary excretion rate for each subject is 1.4 liters of urine per day, or 0.058 liters per hour. For a given sample it was assumed that the measured concentration represented an integrated value for all urine formed since the time of the preceding sample, unless that was more than 5 hours. If the time since the preceding reported sample was more than 5 hours, we assumed that the measured concentration represented an integrated value for urine formed over the preceding 3 hours.

There remained the troublesome problem of determining the delay time of UF6 and resulting compounds in the respiratory tract. Analysis of data of Boback (1975), Chalabreysse (1970), and Moore and Kathren (1985) for persons accidentally exposed to uranium herafluoride indicate that any of three different patterns of excretion may be observed. The excretion rates for some persons appear to be the same as if there were almost no delay of material in the lungs, since the pattern of arinary excretion was very similar to our blood-to-urine model given above. Other persons exhibited a slight extention of urinary excretion of uranium, as if there were a brief delay (with $T_{1/2}$ of perhaps 0.5 day or greater) in the lungs. For a third small group of persons, there appeared a nearly bimodal distribution of the uranium excretion rates, with high initial rates falling quickly toward zero after a few days and later rising abruptly, sometimes to levels approaching the initial excretion rates. The bimodel distribution of excretion rates appears to be fairly unusual and was not considered in the present analysis.

In ICRP 30 it is recommended that inhaled uranium hexafluoride be considered as a class D compound. The ICRP lung model for class D compounds is intended only as a first approximation for a large class of soluble materials that actually have a fairly wide distribution of clearance times, varying up to several days. It is evident from the published data on uranium hexafluoride mentioned above that a substantial portion of inhaled uranium hexafluoride (much more than would be predicted by the ICRP lung model) must reach the blocd fairly rapidly.

We specialized the ICRP lung model for class D materials to uranium hexafluoride and to the present exposure situation by choosing two fractions, one clearing rapidly and the other more slowly, in such a way as to "stabilize" our estimates of intake. (In agreement with the ICRP model for class D material, "rapid" was taken to mean $T_{1/2}=0.01$ days and "more slowly" was taken to mean $T_{1/2}=0.5$ days.) Suppose, for example, that there are six urine samples for subject A, taken at various times. Given a lung retention model and our blood-to-urine model, the concentration of uranium in each sample can be used to estimate an intake at time 0. Uncertainties associated with our blood-to-urine model appear small compared with those associated with the delay in the lungs. Thus, a large variance in the six estimates would suggest to D: that we have incorrectly estimated delay in the respiratory tract, while stable estimates (a small variance) would suggest that the delay has been estimated reasonably well. There is no single lung model that would yield most stable estimates for all subjects. For most subjects, estimates were

stabilized reasonably well by assuming that 30% of the activity enters blood with half-time 0.01 days and 20% enters blood with half-time 0.5 days; the total amount assumed to enter blood (50%) agrees with the ICRP lung model for class D material. Our specialized lung model was used for all subjects. We note that estimates were stabilized nearly as well and fairly similar estimates were obtained by assuming no delay at all in the lungs and 50% clearance to blood. Estimates became very unstable if long-term delays in the lung were assumed.

Estimates for all 45 subjects are summarized in the attached table. The best choice of a lung model, and the subsequent estimates of intake, are expected to vary somewhat as more data on urinary excretion are collected.

REFERENCES

- Bernard SR and Struxness EG, 1957, <u>A Study of the Distribution and Excretion of Uranium in Man</u>, ORNL 2304.
- [2] Boback MW, 1975, A review of uranium excretion and clinical urinalysis data in accidental exposure cases, Report ERDA 93, <u>Confer-</u> ence on <u>Occupational Health Experience with Uranium</u>, 28-29 April 1975, pp. 226-231.
- [3] Hursh JB and Spoor NL, 1973, Data on man, in Hodge HC, Stannard JN, Hursh JB (eds.), <u>Uranium</u>, <u>Plutonium</u>, <u>Transplutonic Elements</u>, New York: Springer-Verlag, pp. 197-239.
- [4] International Commission on Radiological Protection, 1968, Publication 10, <u>Report on Evaluation of Radiation Doses to Body Tissues</u> from Internal Contamination due to Occupational Exposure. Oxford: Permagon Press.
- [5] International Commission on Radiological Protection, 1979, Publication 30, <u>Limits for Intakes of Radionuclides by Workers</u>, Oxford: Permagon Press.
- [6] Moore RH and Kathren RL, 1985, A World War II uranium hexafluoride inhalation event with pulmonary implications for today, J. Occup. Med. 27, 753-756.
- [7] Chalabreysse J, 1970, Etude et resultats d'examens effectues a la suite d'une inhalation de composes dits solubles d'uranium naturel, <u>Radioprotection</u> 5, 1-17.

APPENDIX 5.6.4

103

1

CHEMICAL BIOASSAY RESULTS FOR SFC WORKERS

2,3

Name (1D)	Sex	Weight	Date	Time	Osmolality (MOS/KG)	Creatinine (MG/DL)	Protein (MG/DL)	Glucose (MG/DL)	LDH (MU/24HR)	Albumin (MG/L)	(β -M) (ng/ml)	(NAG) (nmol/min/ml)
		11-5	01/23/86	1510	965	473.1	5	18	21	<68	244	9.05
1E	м	145	01/24/86	0715	513	96.6	ND	14	5	<68	247	3.86
			01/25/86	0800	467	96.1	1	2.9	3	<68	278	2.69
	м	208	01/23/86	0800	562	135.4	1	6	5 4	<68	160	7.78
SE		200	01/24/86	0600	299	83.6	ND	4	4	<68	144	7.44
	2.5		01/02/06	0000	998	182.9	ND	9	9	<68	244	10.23
3E	м	260	01/23/86 01/24/86	0900	714	146.3	1	6	5 8	<68	133	5.02
			01/25/86	0700	728	152.5	1	6	8	<68	173	6.38
4E	м		01/23/86	1115	258	34.5	ND	2	2	<68	40 57	1.24
41	14		01/24/86	1210	261	36.5	ND	2	1	<68	142	4.54
			01/25/85	0900	629	142.1	1	5	8	<68		
6.5	м	207	01/23/86	0800	1021	267.0	14	11	10	<68	189	5.72 2.98
5E	11	201	01/24/86	1120	858	161.1	ND	7	1	<68 <68	156 169	1.61
			01/25/86	0600	1029	234.0	2	1	10	~00		
		202	01/23/86	0800	237	46.6	4	9	2	<68	20	5.25
6E	F	202	01/24/86	0730	410	34.1	1	1462	2	<68	80	3.91 4.06
			01/25/86	0730	228	38.2	2	6	1	<68	20	
		1.70	01/23/86	?	496	137.0	5	0.2	7	<68	460	7.41
7E	м	170	01/24/86	2	605	144.3	9	6	5	<68	889	4.55 2.35
			01/25/86	1000	352	100.0	9	3.2	6	<68	256	
		100	01/23/86	1030	778	268.0	2	8	11	< 68	120	6.23
8E	M	190	01/24/86	1100	361	105.9	ND	14	5	<68	84	3.03 9.23
			01/25/86	0900	800	312.6	5	11	13	<68	187	
	~	1.0.2	01/22/86	0845	281	31.3	ND	1	1	<68	20	1.08
9E	F	123	01/23/86 01/24/86	0905	463	134.6	ND	8	6	<68	289	4.01
			01/25/86	0750	398	123.2	1	4.9	5	<68	511.11	3.06
			01/01/06	1145	875	226.0	ND	13	12	<68	100	3.06
10E	м	190	01/23/86 01/24/86	1357	743	253.6	3	14	11	<68	256	2.28
			01/25/86	1025	400	69.4	ND	13.4	.3	<68	51	0.63
		1.6.6	01/23/86	1000	750	148.0	3	5	8	<68	120	4.25
11E	м	145	01/24/86	1100	1.20	146.0	ND	7	7	<68	196	5.29 6.12
			01/25/86	1000	789	169.3	1	6	8	<68	187	0.12
			01/02/04	1020	1114	446.0	ND	13	19	<68	156	12.66
12E	м	210	01/23/86	1030	990	352.8	3 .	9	12	<68	550	6.78
			01/24/86 01/25/86	7	890	181.0	2	7	11	<68	164	2.88
		179	01/23/86	1530	246	35.5	ND -	3	2	<68	80	3.56
14E	м	178	01/24/86	and the same same	371	65.2	ND	5	5	<68	256	6.72
			01/25/86		153	15.6	ND	0.4		<68	36	
		100	01/22/86	?	215	35.5	5	2	3	<68	30	3.64
15E	F	155	01/23/86 01/24/86	?	320	43.0	1	12	5	<68	109	1.98
			01/25/86		433	60.0	1	0.1	14	<68	151	1.95
			01/22/06	0920	1014	212.0	4	20	11	<68	267	8.78
16E	м	155	01/23/86 01/24/86		982	244.4	4	18	14	<68	700	13.65
			01/24/00	Sars								

	Name (1D)	Sex	Weight	Date	Time	Osmolality (MOS/KG)	Creatinine (MG/DL)	Protein (MG/DL)	Glucose (MG/DL)	LDH (MU/24HR)	Albumin (MG/L)	(B-M) (ng/ml)	(NAG) (nmol/min/ml)
	17E	м	210	01/23/86 01/24/86 01/25/86	0830 1000 0600	218 260 481	29.0 27.3 69.4	ND ND 1	3 2 2.9	1 2 4	<68 <68 <68	43 43 311	1.29 1.09 2.35
	18E	м	170	01/23/86 01/24/86 01/25/86	1605 1530 0900	347 232 437	60.0 18.4 158.1	ND ND 1	3 2 5	5 0 7	<68 <68 <68	66 32 222	1.16 0.57 5.91
	19E	м	210	01/23/86 01/24/86 01/25/86	0923 0850 0935	459 870 538	69.9 319.2 50.5	2 9 2	4 15 3	16 4	<68 69 <68	89 236 58	1.64 21 4.64
	20E	м	155	01/23/86 01/24/86 01/25/86	0925 0855 1145	382 395 192	55.8 64.8 17.2	ND ND ND	4 3 1	3 4 1	<68 <68 <68	116 89 36	1.16 2.19 1.29
	21E	м	180	01/23/86 01/24/86 01/25/86	1118 0800 0730	756 316 858	140.7 42.7 151.9	1	7 3 5.7	6 3 6	<68 <68 <68	124 72 367	3.59 1.43 6.46
	22E	м	187	01/23/86 01/24/86 01/25/86	1100 0855 0875	494 462 579	74.5 80.6 88.1	ND ND	4 5 4.2	543	<68 <68 <68	22 33 98	1.64 0.64 1.89
	23E	м	190	01/23/86 01/24/86 01/25/86	1145 1235 1040	730 828 389	327.0 237.8 168.0	7 ND 1	15 67 15	14 10 9	<68 <68 <68	278 256 278	10.39 4.97 2.91
434	25E	м	274	01/23/86 01/24/86 01/25/86	1100 0905 0400	715 966 362	241.8 368.1 100.3	ND 5 1	5 9 3.5	9 16 7	<68 <68 <68	109 267 151	2.43 2.62 2.32
	26E	F	182	01/23/86 01/24/86 01/25/86	1120 1130 0430	937 360 497	223.0 101.7 116.0	17 ND 1	22 9 6.1	24 18 17	<68 <68 <68	100 200 322	10.63 7.52 6.75
	27E	м	172	01/23/86 01/24/86 01/25/86	1050 0905 0715	319 430 741	43.1 73.7 154.4	3 ND 1	1 7 12	3 3 8	<58 <68 <68	40 67 129	1.29 1.76 2.53
	28£	м	210	01/23/86 01/24/86 01/25/86	0900 1100 7	806 612	158.6 156.6	2	9 6	8 8	<68 <68 <68	105 189	2.53 3.61
	29E	м	190	01/23/86 01/24/86 01/25/86	? ? ?	718 769 903	128.3 161.4 186.8	3 5 5	7 8 7.4	579	<ó8 <68 <68	61 156 344	0.95 3.96 2.98
	30E	м	118	01/23/86 01/24/86 01/25/86	? 1025 0830	693 410 471	272.0 84.0	ND ND 1	12 2 3	13 3	<68 <68 <68	116 82 128	9.79 2.15 4.64
	31E	м	153	01/23/86 01/24/86 01/25/86	0915 0930 0915	946 1019 927	175.4 232.7 228.0	15 1 4	7 8 6.1	11 18 14	<68 <68 <68	10 124 324	0.74 4.35 3.61
	8C	м	220	01/23/86 01/24/86 01/25/86	? 1230 1015	679 716 763	158.3 209.8 189.9	ND 4 ND	15 15 18.2	7 12 12	<68 <68 <68	129 249 356	7.22 5.74 5.88

APPENDIX 5.6.5

CREATININE CORRECTION METHODOLOGY AND FACTORS

Creatinine Correction

The following calculations of 24-hour creatinine excretion values are based on a 1.25 l/day average urine volume, the measured creatinine values, and the body weight data provided. On the basis of the mg creatinine/day/kg body weight estimates, all urinary indicator values were normalized by determining the ratio of each estimate to a normal value of 23 mg creatinine/day/kg body <u>whenever</u> the estimate lay outside a normal range of creatinine excretion values, viz., 16 to 30 mg creatinine/day/kg. Each ratio is listed as the creatinine correction. Two examples of the application of the creatinine correction to the urinary-injury indicator measurements follow the tabulation.

Date	Subject	mg/1	mg/day	mg/day/kg body wt.	creatinine correction
1/23/86	1E	4731	5914	90	0.26
1/24/86		966	1208	18	(1)
1/25/86		961	1201	18	(1)
	2E	1354	1693	18	(1)
		836	1045	11	2.1
	3E	1829	2286	19	(1)
		1463	1829	16	(1)
		1525	1906	16	(1)
	4E	345	431		
		365	456		
		1421	1776		
	5E	2670	3338	35	0.66
		1611	2014	22	(1)
		2340	2925	32	0.72
	6E	466	583	6	3.8
		341	426	5	4.6
	··· N	382	478	5	4.6
	7E	1370	1713	22	(1)
		1443	1804	23	(1)
		1000	1250	16	(1)
	8E	2680	3350	39	0.60
		1059	1324	15	1.5
		3129	3911	45	2.0
	9E	313	391	7	3.3
		1346	1683	30	(1)
		1232	1540	28	(1)
	10E	2260	2825	33	0.70
		2536	3170	37	0.62
		694	868	10	2.3

ł

Creatinine Data

		Creatinin	ne (continued)	
11E	1480	1850	28	(1)
	1460	1825	28	(1)
	1693	2116	32	0.72
12E	4460	5575	59	0.39
	3528	4410	46	0.50
	1810	2263	24	(1)
13E				
		- -		
14E	355	444	6	3.8
	652	815	10	2.3
	156	195	2	11
15E	355	444	6	3.8
	430	538	8	2.9
	600	750	11	2.1
16E	2120	2650	39	0.60
	2440	3050	44	0.52
17E	290	363	4	5.8-
	273	341	4	5.8
	694	668	9	2.6
18E	600	750	10	2.3
	184	230	3	7.7
	1581	1976	26	(1)
19E	699	874	9	2.6
	3192	3990	42	0.55
	505	631	7	3.3
20E	558	698	10	2.3
	648	810	12	1.9
	172	215	3	7.7
21E	1407	1759	21	(1)
	427	534	7	3.3
	1519	1899	23	(1)

22E 745 931 11	
	2.1
806 1008 12	1.9
881 1101 13	1.8
23E 3270 4088 48	0.48
2378 2973 35	1.5
1680 2100 24	(1)
24E	**
25E 2418 3023 24	(1)
3681 4601 37	0.62
1003 1254 10	2.3
26E 2230 2788 34	0.68
1017 1271 15	(1)
1160 1450 17	(1)
27E 431 539 7	3.3
737 921 12	1.9
1544 1930 25	(1)
28E 1586 1983 21	(1)
1566 1958 21	(1)
29E 1283 1604 19	(1)
1614 2018 23	(1)
1868 2335 27	(1)
30E 2720 3400 63	0.37
840 1050 19	(1)
31E 1754 2193 31	0.74
2327 2909 42	0.55
2280 2850 41	0.56
8C 1583 1979 20	(1)
2098 2623 26	(1)
1899 2374 24	(1)

Examples of the Application of the Creatinine Correction to Urinary

Measurements from Subject 5E and 6E

-

.

Subject 5E

Measured Values	Re-expressed Values	Creatinine Corrected Values
NAG 1.24 nmoles/min/ml	1.24 nmoles/min/ml	0.82 nmoles/min/ml
1.71	1.71	1.71
4.54	4.54	3.27
BMG 189 ng/ml	189 µg/1	125 µg/1
156	156	156
169	169	122
Protein 4 mg/dl	0.53 mg/day/kg	0.35 mg/day/kg
2	0.27	0.19
Glucose 11 mg/dl	110 mg/1	73 mg/1
7	70	70
7	70	50
Subject 6E		
Measured Values	Re-expressed Values	Creatinine Corrected Values
NAG 5.25 nmoles/min/ml	5.25 nmoles/min/ml	20.0 nmoles/min/ml
3.91	3.91	18.0
4.06	4.06	18.7
BMG 20 ng/ml	20 µg/1	70 µg/1
80	80	368 µg/1
20	20	92 µg/1
Protein 4 mg/dl	0.55 mg/day/kg	2.1 mg/day/kg
1	0.13	0.6
2	0.27	1.2
Glucose 9 mg/dl	90 mg/1	342 mg/1
1462	14620	67000
6	60	276

APPENDIX 6.1.2.1

١

۱

.

SUMMARY OF DOSE CONVERSION FACTORS AND INTAKE PARAMETERS USED FOR RADIOLOGICAL DOSE CALCULATION

External Dose Conversion Factors

Releases of radioactive gases and particulates to the atmosphere may result in external doses by exposure to and/or immersion in the plume and by exposure to contaminated land surfaces. The dose conversion factors are summarized by Kocher in ORNL/MUREG-70, and those used in this report are shown in Table 6.1.2.1.1. Usage factors are given in Table 6.1.2.1.4.

		Or	gan	
Radionuclide	Total body	Bone	Kidney	Lungs
E	xposure to ground s	urfaces (millir	rem∕year per µCi	/cm ²)
234U	7.1×10^2	3.0×10^{2}	1.0×10^{2}	1.7×10^{2}
2350	1.5×10^{5}	2.1×10^{5}	1.3 x 10 ⁵	1.4 x 10 ⁵
238U	5.7 x 10 ²	2.1×10^{2}	5.9 x 10 ¹	1.2×10^{2}
	Immersion in ai	r (millirem/yea	ar per µCi∕cm³)	
234U	6.8×10^{5}	7.1×10^{5}	3.7×10^{5}	4.1×10^{5}
235U	6.8×10^{8}	9.4×10^{8}	5.9×10^{8}	1.2×10^{8}
2381	4.6×10^{5}	4.5×10^{5}	2.2×10^{5}	2.5×10^{5}

Table 6.1.2.1.1 Dose conversion factors for external exposure pathways

Source: D. C. Kocher, Dose-Rate Conversion Factors for External Exposure to Photons and Electrons, ORNL/NUREG-79, Oak Ridge National Laboratory, August 1981.

Con	nmitted dose eq	uivalent (rem/µ	iCi)
Total body	Bone	Kidneys	Lungs
	Class D		
2.5	7.7 × 10 ¹	1.6×10^{2}	9.4 × 10-1
2.2	6.9 x 10 ¹	1.5×10^{1}	8.5 × 10-1
2.1	6.8 × 10 ¹	1.5 x 10 ¹	8.3 × 10-1
	Total body 2.5 2.2	Total body Bone Class D 2.5 7.7 x 10 ¹ 2.2 6.9 x 10 ¹	Class D 2.5 7.7×10^{1} 1.6×10^{2} 2.2 6.9×10^{1} 1.5×10^{1}

Table 6.1.2.1.2 Dose conversion factors for inhalation exposure pathways--AMAD = 1 μm

Source: D. E. Dunning, Jr., G. G. Killough, S. R. Bernard, J. G. Pleasant, and P. J. Walsh, Estimates of Internal Dose Equivalent to 22 Target Organs for Radionuclide Occurring in Routine Releases From Nuclear Fuel Cycle Facilities. Vol. III, OTNL/NUREG/TM-190/V3, Oak Ridge National Laboratory, October 1981.

Table 6.1.2.1.1 Dose conversion factors for external exposure pathways

	Cor	nmitted dose ed	quivalent (rem/)	uCi)
Radionuclide	Total body	Bone	Kidneys	Lungs
		Class D		
234U	2.6×10^{-1}	7.8	1.7	1.7×10^{-2}
235U	2.2×10^{-1}	7.1	1.5	1.6×10^{-2}
238U	2.2×10^{-1}	7.0	1.5	1.5×10^{-2}

Source: D. E. Dunning, Jr., G. G. Killough, S. R. Bernard, J. G. Pleasant, and P. J. Walsh, Estimates of Internal Dose Equivalent to 22 Target Organs for Radionuclide Occurring in Routine Releases From Nuclear Fuel Cycle Facilities. Vol. III, OTNL/NUREG/TM-190/V3, Oak Ridge National Laboratory, October 1981.

Pathway	Exposed individual
Vegetables, kg/year	94
Milk, L/year	310
Meat, kg/year	110
Inhalation, m ³ /year	8000

Table 6.1.2.1.4 Intake parameters (adult) used in lieu of site-specific data

U.S. NUCLEAR REGULATORY COMMISSION	1. REPORT NUMBER (Assigned by DDC)
BIBLIOGRAPHIC DATA SHEET	NUREG-1189, Volume 2
Assessment of the Public Health Impact From The Accidental	and the second
Release of UF6 at the Sequoyah Fuels Corporation Facility at Gore, Oklahoma	3. RECIPIENT'S CCESSION NO.
AUTHOR(S)	5. DATE REPORT COMPLETED
Ad Hoc Interagency Public Health Assessment Task Force	March YEAR March 1986
PERFORMING ORGANIZATION NAME AND MAILING ADDRESS (Include Zip Code)	DATE REPORT ISSUED
U.S. Nuclear Regulatory Commission Washington, DC 20555	Mirch 1986 6. (Lawe blank)
	8. (Leave blank)
2. SPONSORING ORGANIZATION NAME AND MAILING ADDRESS (Include Zip Code)	
Same as above.	10. PROJECT/TASK/WORK UNIT NO.
	11. FIN NO.
3. TYPE OF REPORT PERIOD ODVER	RED (Inclusive dates)
Task Force Report	
	14. (Leave blank)
Reference Docket 40-8027 and License SUB-1010 16. ABSTRACT (200 words or less) Following the accidental release of UEs from the Sequeval	Fuels Facility on January 4,
Reference Docket 40-8027 and License SUB-1010	Fuels Facility on January 4, Force was established. The us agencies who have prepared th the accidental release. le as of February 14, 1986, and intake of uranium and fluoride. ake of uranium and fluoride sk Force.
16. ABSTRACT (200 words or less) Following the accidental release of UF6 from the Sequoyah 1986, an Ad Hoc Interagency Public Health Assessment Task Task Force consists of technical staff members from various this assessment of the public health impact associated with The assessment is based on data from the accident availabil describes the chemical and radiological effects from the Volume 1 of the report describes the effects from the inter and summarizes the findings and recommendations of the Task Volume 2 of the report consists of Appendices which provide	Fuels Facility on January 4, Force was established. The us agencies who have prepared th the accidental release. le as of February 14, 1986, and intake of uranium and fluoride. ake of uranium and fluoride sk Force.
Reference Docket 40-8027 and License SUB-1010 16. ABSTRACT (200 words or less) Following the accidental release of UE6 from the sequoyah 1986, an Ad Hoc Interagency Public Health Assessment Task Task Force consists of technical staff members from variou this assessment of the public health impact associated with The assessment is based on data from the accident available describes the chemical and radiological effects from the inter Volume 1 of the report describes the effects from the inter and summarizes the findings and recommendations of the Tail Volume 2 of the report consists of Appendices which providused in the assessment.	Fuels Facility on January 4, Force was established. The us agencies who have prepared th the accidental release. le as of February 14, 1986, and intake of uranium and fluoride. ake of uranium and fluoride sk Force.
Reference Docket 40-8027 and License SUB-1010 16 ABSTRACT (200 words or (ess) Following the accidental release of UK6 from the Sequoyah 1986, an Ad Hoc Interagency Public Health Assessment Task Task Force consists of technical staff members from variou this assessment of the public health impact associated with The assessment is based on data from the accident availabil describes the chemical and radiological effects from the intri- and summarizes the findings and recommendations of the Tail Volume 1 of the report consists of Appendices which providues in the assessment. 17. KEY WORDS AND DOCUMENT ANALYSIS assessment fluoride environmental impact exposure 174. DESCRIPTE	Fuels Facility on January 4, Force was established. The us agencies who have prepared th the accidental release. le as of February 14, 1986, and intake of uranium and fluoride. ake of uranium and fluoride sk Force.
Reference Docket 40-8027 and License SUB-1010 6. ABSTRACT (200 words or less) Following the accidental release of UE6 from the Sequoyah 1986, an Ad Hoc Interagency Public Health Assessment Task Task Force consists of technical staff members from variou this assessment of the public health impact associated with The assessment is based on data from the accident availab describes the chemical and radiological effects from the intra and summarizes the findings and recommendations of the Tail Volume 1 of the report consists of Appendices which providused in the assessment. 17. KEY WORDS AND DOCUMENT ANALYSIS assessment uranium fluoride environmental impact exposure release 17a. DESCRIPTE 17b. IDENTIFIERS/OPEN-ENDED TERMS	Fuels Facility on January 4, Force was established. The us agencies who have prepared th the accidental release. le as of February 14, 1986, and intake of uranium and fluoride. ake of uranium and fluoride sk Force.

Î

UNITED STATES NUCLEAR REGULATORY COMMISSION WASHINGTON, D.C. 20555

OFFICIAL BUSINESS PENALTY FOR PRIVATE USE, \$300

1

SPECIAL FOURTH CLASS RATE POSTAGE & FEES PAID USNRC WASH 0.C PERMIT No. G 67

4

.

.