## CONTAINMENT SYSTEMS

## 3/4.6.4 COMBUSTIBLE GAS CONTROL

## HYDROGEN ANALYZERS

#### LIMITING CONDITION FOR OPERATION

3.6.4.1 Two independent containment hydrogen analyzers shall be OPERABLE.

APPLICABILITY: MODES 1 and 2.

#### ACTION:

. .

- a. With one containment hydrogen analyzer inoperable, recore the inoperable analyzer to OPERABLE status within 30 days or the in at least HOT STANDBY within the next 6 hours.
- b. With both containment hydrogen analyzers inoperable, restore at least one analyzer to OPERABLE status within 72 hours and comply with the requirements of ACTION a, or be in at least HOT STANDBY within the next 6 hours.

#### SURVEILLANCE REQUIREMENTS

- 4.6.4.1 Each hydrogen analyzer shall be demonstrated OPERABLE by the performance of a CHANNEL FUNCTIONAL TEST at least once per 31 days, and at least once per 92 days on a STAGGERED TEST BASIS by performing a CHANNEL CALIBRATION using sample gases containing:
  - a. Zero volume percent hydrogen, balance nitrogen.
  - b. 9.5 volume percent hydrogen, balance nitrogen.

ATTACHMENT B

## CONTAINMENT SYSTEMS

## 3/4.6.4 COMBUSTIBLE GAS CONTROL

HYDROGEN ANALYZERS

#### LIMITING CONDITION FOR OPERATION

3.6.4.1 Two independent containment hydrogen analyzers shall be OPERABLE.

APPLICABILITY: MODES 1 and 2.

### ACTION:

- a. With one containment hydrogen analyzer inoperable, restore the inoperable analyzer to OPERABLE status within 30 days or be in at least HOT STANDBY within the next 6 hours.
- b. With both containment hydrogen analyzers inoperable, restore at least one analyzer to OPERABLE status within 72 hours and comply with the requirements of ACTION a, or be in at least HOT STANDBY within the next 6 hours.

# SURVEILLANCE REQUIREMENTS

- 4.6.4.1 Each hydrogen analyzer shall be demonstrated OPERABLE by the performance of a CHANNEL FUNCTIONAL TEST at least once per 31 days, and at least once per 92 days on a STAGGERED TEST BASIS by performing a CHANNEL CALIBRATION using sample gases containing:
  - a. Zero volume percent hydrogen, balance nitrogen.
  - b. 9.5 volume percent hydrogen, balance nitrogen.

Delete 1

4.6.4.1 AT LEAST ONCE PER 31 days each Hydroson ANALYZER Shall be DEMONSTRATED OPERABLE by PERFORMING A CHANNEL CALIBRATION USING A SANGLE GAS.

LDD-1