



Calibration Control Kit

Part No.: PCD4000

Store at room temperature away from light

Français: Kit de Contrôle de Calibration - No. Art: PCD4000
Tenir aux température ambiante sans lumière

Portugués: Kit de Contrôle Calibragem - No.Peça: PCD4000
Manter em temperatura ambiente sem luz

Español: Kit de control de Calibración - No.Pieza:PCD4000
Preservar en temperatura ambiental sin luz

Italiano: Kit di Controllo Calibrato - Cod.Comp.:PCD4000
Conservare in temperatura ambiente senza luce

Deutsch: Kalibrierungskontroll Kit - Art.Nr.:PCD4000
Lagern auf Raumtemperatur ohne Licht

Nederlands: Calibratiecontrole Kit - Art.Nr.:PCD4000
Bewaren bij kamertemperatuur zonder licht

中文: 校准质控试剂盒 - 编号:PCD4000
室温避光保存

Manufactured for Hygiena by Medical Packaging Corp.

ID: LBL0262
Rev. 0

3" x 2.125"

PLEASE READ YOUR PROOF CAREFULLY BEFORE APPROVAL

Once a proof is signed off as "OK" and is printed, We will **NOT** be held **RESPONSIBLE** for errors and no credit will be given.

* **NOTE:** Please do not be alarmed if some or all of the text look jagged. This is normal for outlined text in a pdf proof.

Label Color:



Black

**DIELINE
DOES NOT
PRINT**

← IMPORTANT!!!

Please Verify Color Swatches Are Correct.

Please indicate the correct unwind position

Unwind Positions

1 <input type="checkbox"/>	2 <input type="checkbox"/>	3 <input type="checkbox"/>	4 <input type="checkbox"/>	5 <input type="checkbox"/>
Top-Off	Bottom-Off	Right-Off	Left-Off	Hand Apply

FOR PMS COLOR, PLEASE REFER TO THE PANTONE COLOR FORMULA GUIDE.

FOR 4 COLOR PROCESS, PLEASE REQUEST FOR A DIGITAL PROOF.

Stock:	Adhesive:	It is the sole responsibility of the client to approve art and copy for accuracy prior to reproduction. Before signing below , please check for: <ul style="list-style-type: none"> • Typographical or graphic error or omissions • Barcode correctness (if applicable) • Compliance to any and all government rules and regulations • Compliance to originally supplied customer art or electronic disk graphics • Die cut to copy registration as well as accuracy of Graphics and Type 	
Client: Hygiena	Issue: 1		<input type="checkbox"/> OK TO PRINT AS IS
File: LBL0262 Calibration Control Kit			AUTHORIZED SIGNATURE _____ DATE _____
Date: 09-22-21			<input type="checkbox"/> PLEASE SUBMIT NEW PROOF
Designer: K.C.			

THE INK SPOT, INC. Phone: (626) 338-4500 Fax: (626) 338-4533